




Laboratory Reagents Handbook



Trusted solutions uncompromised analysis

Our online catalog features access to an expanded range of content to help guide you through any application, on any instrument:

- Find products online
- Access to method development support
- Information on application resources
- Links to chromatography videos



www.thermoscientific.com/catalog



2012-2013

Thermo Scientific Chromatography
Columns and Consumables
Catalog

**Expanded interactive
online version now
available!**



www.thermoscientific.com/catalog

Thermo
SCIENTIFIC
Part of Thermo Fisher Scientific

Your partner of choice for chemicals and reagents

Welcome to the new Laboratory Reagents catalogue which reflects the range of choice we offer you :

- Over 4400 Fisher Chemical products dedicated to many analytical applications
- More than 1000 Fisher BioReagents products for molecular biology research, biochemistry and cellular biology
- Including more than 250 new products as Optima® LC/MS Solvents, high purity Acids for Trace Elemental analysis, microbiology media and modifying enzymes

Fisher Chemical

- Acids & Bases for Trace Elemental analysis see page 14
- Buffers Standard Solutions see page 65
- Indicators see page 158
- Karl Fischer Titration Aqualine™ see page 170
- Pharmacopeia grade : Reagents & Solvents see page 234
- Solvents, Blends and Reagents for Chromatography see page 308
- Standard Elemental solution for Trace Metal analysis see page 314
- Standard Volumetric solution see page 323



Fisher BioReagents

- Molecular Biology see page 385
- Protein Chemistry see page 438
- Cell Biology see page 452
- Core BioReagents & complementary products see page 479
- Index by product name see page 499
- Index by manufacturer number see page 501



Whatever your field of activity: industry, production, quality control, research, analysis or development, our teams provide you with the best service and support available.

Your chemical specialist to answer your needs

Smart chemistry solutions

Chemistry experience and expertise

- Large field sales force with highly experienced specialists to advise and support you
- Huge warehouse space ensures high product availability
- An experienced Customer service team dedicated to making your contact with Fisher Scientific easy and efficient

Bottle labelling meets GHS requirements and includes:

- Certificate of analysis
- Grade and application
- Physical, chemical and transportation data

Specification and safety data sheets, online certificates of analysis

- Specification and safety data sheets can be downloaded based on the product reference
- Certificate of analysis can be downloaded based on the product reference and the lot number at www.acros.com

Packaging for quality, safety and convenience

- Solvents are available in plastic-coated glass bottles that confine the liquid in the event of bottle breakage
- Innovative, environmentally friendly packaging: pre-formed cardboard with integrated absorbent that are entirely recyclable and meet UN regulations for transportation.

Our Customer Service Team are on hand to make placing your order a quick and easy experience.

Request a catalogue
fisherchemical@thermofisher.com



Acros Organics Catalogue



Thermo Scientific Catalogue
Chromatography Columns and Consumables

2 classifications to simplify your searches

To easily navigate through this catalog, you can access product information by using 2 modes of reading:

Classification by applications

Acids & bases for trace elemental analysis	see page 14	Pharmacopeia grade: reagents & solvents	see page 234
Buffers standard solutions	see page 65	Solvents, blends and reagents for chromatography	see page 308
Indicators	see page 158	Standard elemental solutions for trace metal analysis	see page 314
Karl Fischer titration Aqualine™	see page 170	Standard volumetric solutions	see page 323

Alphabetical classification and quality in decreasing order

The adopted nomenclature is in compliance with the International Union of Pure and Applied Chemistry rules. We retained some routinely used denominations.

- 1 colour coded application :
- 2 pack size & packaging
- 3 physical and chemical data
- 4 transportation data
- 5 hazard symbol
- 6 specifications
- 7 Catalogue number and unit pack size code

	Code suffix	Pack size		Code suffix	Pack size
Liquids	07	= 250mL	Solids	48	= 100g
	08	= 500mL		50	= 250g
	15	= 1L		53	= 500g
	17	= 2.5L		60	= 1Kg
	21	= 5L		63	= 3Kg
	24	= 10L		65	= 5Kg
	25	= 25L		70	= 25Kg

Acetonitrile NEW PACK Certified, Eur.Ph		
HPLC for gradient analysis		
packaging	Cat. No	Mfr. No
1 L Glass bottle	1079-4741	A/0627/15
2.5 L Glass bottle	1066-0131	A/0627/17
2.5 L Coated Glass bottle	1063-0131	A/0627/PB17
5 L Aluminium Can	1050-0911	A/0627/21
10 L Returnable Drum	1265-6667	A/0627/10R55
25 L Returnable Drum	1266-6667	A/0627/25R55
30 L Returnable Drum	1134-1842	A/0627/30R55
C, H, N UN: 1648 ADR/RID: 3, II		
CAS: 75-05-8 MDL: 41-04 UN: 1648 ADR/RID: 3, II		
Density at 20°C: 0.781 Boiling point: 81 to 82°C Boiling pressure: 760 mmHg Melting point: -46°C Flash point: 2°C		
Assay (GC): >99.9%		
Guaranteed analysis		
Acidity / Alkalinity	<0.0005 meq/g	Residue after evaporation
Water	<0.01%	<2 ppm
Absorbance specification		
Max. Absorbance (1cm cell v HPLC water)		
195 nm	<0.1 A.U.	230 nm
200 nm	<0.05 A.U.	240 nm
210 nm	<0.02 A.U.	250 nm
220 nm	<0.01 A.U.	Gradient (max eluted peak)
230 nm	<0.01 A.U.	<0.01 A.U.
Gradient analysis:		
Absorbance of the largest eluted peak <0.01 A.U.		
Conditions		
Column: C18 Reverse Phase		
Flow rate: 2 mL / min		
Wavelength: 210 nm		
Absorbance range: 0-0.32 A.U.		
Solvent A: Water		
Solvent B: Acetonitrile		
Gradient: linear 95% / 5% B 100% B during 20 min.		
Lot analysis on the label		
Filtered to 0.2 µm.		

Solvents, reagents, standard solutions and titrants,
A product to suit every application



For specialist solvents, solutions, specialised packaging
for a customised specification,
contact us on: scs.emea@thermofisher.com

4 Fisher Chemical product grades

for a broader offer

Application

For Analysis

Specified

Titration

Optima® LC/MS

Ultra high-purity solvent for use with LC-MS systems. Low background noise in MS detection. Very low metal ion impurities content. Very low detection limit for the development of methods. Reduction in plastic contaminant peaks.

Certified Advanced HPLC Gradient

Advanced HPLC gradient grade. HPLC gradient guarantee with a very low level of gradient baseline drift. Includes lot analysis and absorbance curve on the label.

HPLC Electrochemical

HPLC solvents for Electrochemical and UV detectors. Guaranteed for low electrochemical activity and low UV absorbance/high transmission. Includes lot analysis and absorbance curve on label.

Certified HPLC-MS

HPLC solvents specified for use with LC/MS systems. Guaranteed for low level of trace metals and non-volatile residue. Low level of absorbance, performance under gradient conditions.

CertiFied™ HPLC Gradient

HPLC solvents suitable for gradient analysis. Guaranteed for low absorbance/high UV transmission and low concentration of non-volatile impurities. In some instance may be suitable for fluorescence detection. Includes lot analysis and absorbance curve on the label.

GPC

For gel permeation chromatography. Filtered to 0.2µm. Low water, residue and colour. Unique chemical range – Actual lot analysis on the pack label.

UHPLC Gradient Grade

Solvent certified for UHPLC analysis with high UV transmission. Low background noise at 210nm and 254nm. Filtered at 0.1µm for ultra low particulates.

HPLC Fluorescence

HPLC solvents for Fluorescence and UV detectors. Guaranteed for low fluorescence between 250nm and 750nm emission & excitation wavelengths.

Distol

For pesticide and petroleum residue analysis. Specified for ECD, NPD and FID detectors. GC range supplied with certificate of analysis.

Optima® Grade

ICP-MS Ultra TraceMetal® analysis. Ultrapure acids, bases & water. Tested for up to 65 parameters at 1-100 ppt level. The highest performance acid range.

For Analysis

Certified reagents for analytical application. Tested for up to 18 guaranteed parameters. Actual lot analysis on the pack label.

Buffers

Standard (colour coded) NIST certified solutions for pH. Ready to use, with an accuracy factor of ±0.02 pH at 20°C. Also available as concentrates, packaged in ampoules.

Trace Metal® Grade

TraceMetal® for ICP trace metal analysis. Acids tested for up to 65 parameters at ppb levels.

For Analysis Conform Eur.Ph.

Certified reagents for analytical application. Many of these products also comply with various pharmacopoeias, including BP, EP and USP. Tested for up to 18 guaranteed parameters. Actual lot analysis on the pack label.

Karl Fischer - Aqualine®

Karl Fischer reagents for the determination of moisture Volumetric and coulometric reagents and standards. Pyridine free, rapid titration and a stable end-point. Supplied in single packs or in ampoules.

Primar Plus® Grade

Acids for trace elemental analysis. Acids are tested for up to 40 parameters. Levels to ppb.

Volumetric solution

Standard solutions for volumetric analysis. Accuracy factor up to 0.999 - 1.001 NIST traceability. Ready to use.

Dried Distilled

For organic synthesis. Guaranteed low water content.

Specified Laboratory Reagents (SLR)

Specified Laboratory Reagents for general laboratory applications. Tested for up to 13 parameters.

Solutrate

Concentrated solutions for volumetric analysis. Supplied in singles or pack of six sealed ampoules.



Our analytical reagents and solvents
are produced at a manufacturing site registered ISO 9001:2000 and ISO 14001

Fisher Chemical packaging

for quality, safety and convenience

Plastic bottle HDPE

Volume, mL	Thread	Description	Dimensions, mm	Mass, g
500	DIN45	Natural HDPE	177 x Ø 78	52
1000	DIN45	Natural HDPE	212 x 97 x 97	112
1000	DIN45	Black HDPE	213 x 97 x 97	112
2500	DIN45	Natural HDPE	300 x Ø 134	190
2500	DIN45	Black HDPE	301 x Ø 134	190
5000	DIN45	Natural HDPE	339 x 168 x 168	321
5000	DIN45	Black HDPE	340 x 168 x 168	321

Glass bottle

Volume, mL	Thread	Description	Dimensions, mm	Mass, g
100	DIN31	Amber Glass	122 x 50 x 50	140
250	DIN31	Amber Glass	160 x 63 x 63	270
500	DIN54	Amber Glass	183 x Ø 83	430
1000	DIN45	Amber Glass	230 x Ø 102	685
2500	DIN45	Amber Glass	298 x Ø 134	1300
100	DIN31	Coated Amber Glass	122 x 51 x 51	157
250	DIN31	Coated Amber Glass	160 x 64 x 64	287
500	DIN54	Coated Amber Glass	183 x Ø 84	453
1000	DIN45	Coated Amber Glass	230 x Ø 103	720
2500	DIN45	Coated Amber Glass	298 x Ø 135	1360

HDPE Plastic bottle solids (wide opening)

Volume, mL	Thread	Description	Dimensions, mm	Mass, g
50	DIN28	HDPE silver	88.5 x 36 x 36	13,5
100	DIN32	HDPE silver	101 x 45 x 45	18
250	DIN45	HDPE silver	129.5 x 59 x 59	30
350	DIN45	HDPE silver	164 x 6 x 60	36.5
500	DIN54	HDPE silver	162 x 75 x 75	47
750	DIN54	HDPE silver	188 x 83 x83	58.5
1000	DIN54	HDPE silver	201 x 85 x 85	77
1500	DIN80	HDPE silver	209 x 108 x 108	100
2500	DIN80	HDPE silver	260 x 122 x 122	130
4000	DIN80	HDPE silver	294 x 142 x 142	250

Glass bottle solids

Volume, mL	Thread	Description	Dimensions, mm	Mass, g
30	DIN25	Amber Glass	66 x Ø 33	40
100	DIN32	Amber Glass	110 x 45 x 45	300
250	DIN45	Amber Glass	150 x 60 x 60	750
500	DIN54	Amber Glass	175 x 70 x 70	1500
1000	DIN54	Amber Glass	215 x 90 x 90	3000

Totally recyclable packaging

to respect the environment

Other packaging

Volume	Description	Dimensions, mm	Mass, g
150mL	Plastic Ampoule		15
5000mL	Aluminium can DIN45	300 x Ø 175	360
5000mL	florinated HDPE plastic bottle	339 x 168 x 168	321
5000mL	Florinated HDPE container (stackable)	255 x 193 x 146	255
5000mL	Metal can	200 x Ø 230	1200
10L	Cubitainer	220 h x 230 w	460
10L	Plastic Jerrican HDPE transparent and black	325 x 195 x 230	450
25L	Mild steel drum (2" & 3/4")	473 x Ø 285	3840
25L	Plastic drum	445 x 286 x 286	1350
15L	Plastic keg	297 x Ø 338	810
20L	Plastic keg	377 x Ø 338	1080
25L	Plastic keg	457 x Ø 338	1370
50L	Fibre keg	525 x Ø 335	1868
100L	Fibre keg	775 x Ø 395	3834



Our environmental approach is designed to deal with the main environmental issues

- Chemicals are delivered in totally recyclable pre-formed cardboard with integrated absorbent meeting UN regulations for transportation
- Solvents are available in plastic-coated bottles that confine the liquid in the event of breakage
- All our shipping partners adhere to a sustainable development charter
- Storage, packaging and processing solutions adapted to various types of waste
- Compaction of damaged packaging plastic and films for reprocessing and recovery
- Respect for the environment code and the regulations of Installations Classified for the Protection of the Environment
- Adherence to the EU Regulatory Directives for DEEE, RoHS, REACH



Various types of packaging
to meet your specifications



Respect our planet
this is a reality; it is part of our daily lives

Harmonised labelling system

meeting GHS requirements



New labels

Complying with the Globally Harmonized System (GHS) requirements for classification and labelling of chemicals. See key opposite

New hazard pictograms

Symbols are now more visible and easier to understand. Former symbols printed in black on orange squares are replaced with nine new danger pictograms including black symbols printed in a white lozenge bordered in red. See page 9.

New regulations

Regulation (EC) No. 1272/2008, designated by the GHS or CLP acronyms came into effect on January 20, 2009. This regulation governs the classification, labelling and packaging of substances and mixtures. The English CLP acronym stands for "Classification, Labelling and Packaging of substances and mixtures". As from 2015, this regulation will replace and supersede European Directive 67/548/CEE defining substances and Directive 1999/45/CE defining preparations.

Header-Key

- 1 New GHS symbols
- 2 Fisher Chemical Brand logo
- 3
 - Code for re-ordering
 - Batch(lot) number
 - Packsize
 - Formula & Molecular weight
- 4
 - Chemical name
 - Grade & application
 - Lot analysis showing actual QC data
- 5 General application information
- 6 The H+P Statements
 - Packing date
 - Expiration product information (including date opened)
 - UN & CAS number
 - Storage condition
- 7 Information about the manufacturing site
- 8 Physical product data
- 9 Data Matrix Code (2D Code) for the fastest possible logistic processes

Hazard warning symbols and precautionary phrases

complying with European standards

Physical Hazards

Hazard classes and hazard categories*	Label elements NEW**	Label elements OLD
Explosives • Unstable explosives • Explosives, divisions 1.1 to 1.3	H200 H201, H202, H203	(R2, R3)
Self-reactive substances, mixtures, types A, B Organic peroxides, types A, B	H240, H241, H242, H243	
Explosives, division 1.4	H204	No classification
Flammable gases, category 1 Flammable aerosols, category 1 Flammable liquids, category 1	H220, H222, H224	(R12), (R12), (R12)
Flammable liquids, category 2 Flammable solids, category 1 Flammable solids, category 2	H225, H228, H226	(R11), (R11), (R11)
Flammable aerosols, category 2 Flammable liquids, category 3	H223, H226	No symbol No classification Flashpoint 56-60°C
Pyrophoric liquids, category 1 Pyrophoric solids, category 1 Substances, mixtures which in contact with water emit flammable gases, categories 1, 2 and category 3	H250, H251, H252, H253	(R17), (R17), (R17)
Self-reactive substances, mixtures, type B Self-reactive substances, mixtures, types C, D and types E, F Self-heating substances, mixtures, category 1 and category 2	H241, H242, H243, H251, H252	(R12), (R12)
Organic peroxides, type B Organic peroxides, types C, D Organic peroxides, types E, F	H241, H242, H243	(R7), (R7)
Oxidising gases, category 1 Oxidising liquids, categories 1, 2 and category 3 Oxidising solids, categories 1, 2 and category 3	H270, H271, H272, H273	(R8, R8, R8), (R8, R8, R8)
Gases under pressure • Compressed gases • Liquefied gases • Refrigerated liquefied gases • Dissolved gases	H280, H281, H282	No classification
Corrosive to metals, category 1	H290	No classification

Health Hazards

Hazard classes and hazard categories*	Label elements NEW**	Label elements OLD
Acute toxicity, categories 1, 2 • Oral • Dermal • Inhalation	H300, H310, H330	(R23, R24, R25)
Acute toxicity, category 3 • Oral • Dermal • Inhalation	H302, H312, H332	(R22, R23, R24)
Skin corrosion, categories 1A, 1B Corrosion, categories 1A, 1B Reproductive toxicity, categories 1A, 1B STOT***, single exposure, category 1 STOT***, repeated exposure, category 1	H314, H335, H360D1, H360D2, H373	(R34, R35, R36, R37, R38)
Respiratory sensitisation, category 1 Respiratory sensitisation, category 2	H334, H335	(R34, R35)
Skin sensitisation, category 1 Skin sensitisation, category 2 Reproductive toxicity, category 2 STOT***, single exposure, category 2 STOT***, repeated exposure, category 2	H335, H360D1, H360D2, H373	(R34, R35, R36, R37, R38)
Acute toxicity, category 4 • Oral • Dermal • Inhalation	H302, H312, H332	(R22, R23, R24)
Skin corrosion, categories 1A, 1B, 1C	H314	(R34, R35)
Skin irritation, category 1	H315	(R36, R37)
Skin irritation, category 2 Eye irritation, category 2 Skin sensitisation, category 1 STOT***, single exposure, category 3 • Respiratory tract irritation • Mucous membranes	H315, H319, H335, H336	(R36, R37, R38)

Environmental Hazards

Hazardous to the aquatic environment, category 1	H400 H410	Warning		H400 H410
Hazardous to the aquatic environment, category 2	H411			H411

Detail of the H & P & EUH phrases ► see p. 496

New hazard pictograms, clearer and safer



www.acros.com
to access all GHS information

Your simplified approach to regulated products

Products subject to legal provisions

Ethanol: ethyl alcohol

For any purchase order bound for UK and UK Overseas Departments, please send us a copy of the "Advance Notification of Occupation" or the warehouse certification (UT No or excise duty No) certified by the Customs Department.
A copy of this document must be sent to our Customer Relations Department.

Chemical weapons precursors, drug precursors and regulated export products

They are subject to preliminary declaration. According to the classification of the products, their quality, quantity and the country of destination, the enforcement authorities issue a delivery authorisation.

Substances which weaken the ozone layer

The use of substances which weaken the ozone layer (ODS: Ozone-Depleting Substances) imposes new delivery requirements:

- Production of your identification number that you must request from the European Commission. This number can be obtained by visiting the website: <http://ec.europa.eu/clima/policies/ods> found on the Laboratory-ODS Database tab.
- The e-mail response from the European Commission showing the number must be forwarded to our Customer Services Department
- Subsequently fill out and sign the "Notification" document downloadable from the www.fisher.co.uk website. This notification must be renewed along with each order.

REACH

Isolated intermediaries

In accordance with the (EC) N 1907/2006 "REACH" regulation, some substances are listed as isolated intermediates.

These substances are only be delivered to our customers upon presentation of a certificate of use in compliance with the strictly controlled conditions as specified under Article 18 (4) of the present regulation.

Regulation

The implementation of REACH calls for a close collaboration among Thermo Fisher Scientific, its suppliers and its customers, which allows us to guarantee continuity in terms of availability of our chemicals on the basis of the REACH requirements. The chemicals that are affected by the REACH regulation are registered. We follow very closely all the substances which are or will be listed and will require autorisation at a later stage.

For further information

AUSTRALIA	auinfo@thermofisher.com
CHINA	cnchemical@thermofisher.com
INDIA	qfc.customer@thermofisher.com
JAPAN	sales.bid.jp@thermofisher.com
KOREA	sales.fsk@thermofisher.com
MALAYSIA	enquiry.my@thermofisher.com
NEW ZEALAND	nzinfo@thermofisher.com

Your specialised chemical services partner We support your chemistry

Specialised Chemical Services

When you need a product that is not readily available in a chemical catalogue, we can help.

We are in a position to offer you solutions adapted to your needs through our knowledge and experience in this area:

- a portfolio of over 100,000 products
- a global network of more than 3,500 suppliers and manufacturers

Our Specialised Chemical Services team is ready directly to discuss your specific chemical needs

- Semi-bulk and bulk chemical service
- Tailored solvents and solvent blends & solutions
- Testing services
- Customised packaging and labelling

For further information

AUSTRALIA	auinfo@thermofisher.com
CHINA	cnchemical@thermofisher.com
INDIA	qfc.customer@thermofisher.com
JAPAN	sales.bid.jp@thermofisher.com
KOREA	sales.fsk@thermofisher.com
MALAYSIA	enquiry.my@thermofisher.com
NEW ZEALAND	nzinfo@thermofisher.com



A personalized advice
to optimize the processing of your order



For custom solutions
our team of scientists is at your disposal

- Acids & bases for trace elemental analysis see page 14
- Buffers standard solutions see page 65
- Indicators see page 158
- Karl Fischer titration Aqualine™ see page 170
- Pharmacopeia grade : reagents & solvents see page 234
- Solvents, blends and reagents for chromatography see page 308
- Standard elemental solution for trace metal analysis see page 314
- Standard volumetric solution see page 322

Acetaldehyde

specified reagents for general laboratory work- extra pure

SLR

packaging	Mfr. No
500 mL Glass bottle	A/0080/08
2.5 L Glass bottle	A/0080/17
C ₂ H ₄ O CAS: 75-07-0 MW: 44.05 Density at 20°C: 0.785 Boiling point: 21°C Melting point: -123°C Flash point: -27°C ADR/RID: 3, I	EINECS: 200-836-8 UN: 1089 GHS: Flam. Liq. 1, Carc. 2, Eye Irrit. 2, STOT SE 3 H224, H319, H351, H335 P261, P280, P305+P351+P338, P210, P243, P403+P235

Assay (CG): >99%

Guaranteed analysis

Acidity <0.05 meq/g

Acetamide

specified reagents for general laboratory work- extra pure

SLR

packaging	Mfr. No
100 g Plastic jar	A/0160/48
C ₂ H ₅ NO CAS: 60-35-5 MW: 59.07 Density at 20°C: 1.16 Boiling point: 221°C	Melting point: 79 to 81°C EINECS: 200-473-5 GHS: Carc. 2 H351 P301+P310, P281

Assay: >98%

Acetanilide

specified reagents for general laboratory work- extra pure

SLR

packaging	Mfr. No
250 g Plastic jar	A/0200/50
500 g Plastic jar	A/0200/53
C ₆ H ₅ NO CAS: 103-84-4 MW: 135.17 Boiling point: 304°C Melting point: 111 to 115°C Flash point: 173°C	EINECS: 203-150-7 GHS: Acute Tox. 4, Eye Irrit. 2 H302, H319 P301+P312, P280, P305+P351+P338

Guaranteed analysis

Residue after ignition <0.05%

Acetic acid

for LC/MS analysis

NEW

Optima

packaging	Mfr. No
1 mL Ampoule	A113-1AMP
10x1 mL Ampoule	A113-10X1AMP
50 mL HDPE Plastic bottle	A113-50
C ₂ H ₄ O ₂ CAS: 64-19-7 MW: 60.05 Density at 20°C: 1.048 Boiling point: 117 to 118°C Melting point: 16 to 16.5°C Flash point: 40°C ADR/RID: 8, 3, II	EINECS: 200-580-7 UN: 2789 GHS: Flam. Liq. 3, Skin Corr. 1A H226, H314 P280, P305+P351+P338, P302+P352, P304+P340, P210, P309+P311

Assay: 99.7%

Guaranteed analysis

Appearance Colorless, clear liquid	Magnesium (Mg) <10 ppb
Colour 10 APHA	Manganese (Mn) <10 ppb
Water <0.1%	Nickel (Ni) <10 ppb
Residue after evaporation <5 ppm	Potassium (K) <20 ppb
Also meets ACS specifications : Pass Test	Propan--1-ol <10 ppb
Aluminium (Al) <20 ppb	Sodium (Na) <50 ppb
Calcium (Ca) <50 ppb	Sulfate (SO ₄) <0.00004 ppm
Copper (Cu) <10 ppb	n-Butyl formate <0.00004 ppm
Iron (Fe) <10 ppb	Zinc (Zn) <10 ppb
Lead (Pb) <10 ppb	

Absorbance specification

Max. Absorbance (1cm cell v HPLC water)

255 nm	: <1.0 A.U.	LC/MS Gradient Suitability (at 0.1% aqueous)	: Pass Test
280 nm	: <0.05 A.U.		
350 nm	: <0.02 A.U.		

This product has been tested on an LC-MS Thermo Scientific instrument to ensure suitability.

Acetic acid glacial

for HPLC

Certified

packaging	Mfr. No
500 mL Coated Glass bottle	A/0406/PB08
1 L Coated Glass bottle	A/0406/PB15
C ₂ H ₄ O ₂ CAS: 64-19-7 MW: 60.05 Density at 20°C: 1.048 Boiling point: 117 to 118°C Melting point: 16 to 16.5°C Flash point: 40°C ADR/RID: 8, 3, II	EINECS: 200-580-7 UN: 2789 GHS: Flam. Liq. 3, Skin Corr. 1A H226, H314 P280, P305+P351+P338, P302+P352, P304+P340, P210, P309+P311

Assay: >99.7%

Absorbance specification

Max. Absorbance (1cm cell v HPLC water)

260 nm (sol. aqu. 1 M)	: <0.1 A.U.	280 nm (sol. aqu. 1 M)	: <0.02 A.U.
270 nm (sol. aqu. 1 M)	: <0.02 A.U.		

Acetic acid glacial (d=1.08) is equivalent to a strength of 17.4M (17.4N).
Lot analysis on the label.
Filtered to 0.2 µm.

Acetic acid glacial

HPLC for electrochemical detection

Certified

packaging	Mfr. No
500 mL Coated Glass bottle	A/0407/PB08
C ₂ H ₄ O ₂ CAS: 64-19-7 MW: 60.05 Density at 20°C: 1.048 Boiling point: 117 to 118°C Melting point: 16 to 16.5°C Flash point: 40°C ADR/RID: 8, 3, II	EINECS: 200-580-7 UN: 2789 GHS: Flam. Liq. 3, Skin Corr. 1A H226, H314 P280, P305+P351+P338, P302+P352, P304+P340, P210, P309+P311

Assay: >99.7%

Guaranteed analysis

Electrochemical test	Pass	Lead (Pb)	<5 ppm
Iron (Fe)	<5 ppm	n-Butyl formate	<0.0005%

Absorbance specification

Max. Absorbance (1cm cell v HPLC water)

260 nm (sol. aqu. 1 M)	: <0.1 A.U.	280 nm (sol. aqu. 1 M)	: <0.02 A.U.
270 nm (sol. aqu. 1 M)	: <0.02 A.U.		

Tested by HPLC-EC to ensure low electrochemical response.

Detector - ESA Coulochem II

Mode-oxidative

Electrode- Graphite.

Acetic acid glacial (d=1.08) is equivalent to a strength of 17.4M (17.4N).

Lot analysis on the label.

Filtered to 0.2 µm.

We Support Your Chemistry For

- Semi-bulk and bulk requirements
 - Tailored solvents and solvent blends or customized solutions
 - Testing services or customized packaging and labeling
- Please email your contact details and a brief description of your requirements to begeel.sha.quotes@thermofisher.com



Acetic acid glacial min.99%
for ultra trace elemental analysis

Optima

packaging		Mfr. No
250 mℓ	Pre-cleaned Teflon FEP bottle	A465-250
500 mℓ	Pre-cleaned Teflon FEP bottle	A465-500
1 ℓ	Pre-cleaned Teflon FEP bottle	A465-1
C ₂ H ₄ O ₂ CAS: 64-19-7 MW: 60.05 Density at 20°C: 1.048 Boiling point: 117 to 118°C Melting point: 16 to 16.5°C Flash point: 40°C ADR/RID: 8, 3, II		EINECS: 200-580-7 UN: 2789 GHS: Flam. Liq. 3, Skin Corr. 1A H226, H314 P280, P305+P351+P338, P302+P352, P304+P340, P210, P309+P311



Assay: >=99%			
Guaranteed analysis			
Zirconium (Zr)	<10 ppt	Manganese (Mn)	<10 ppt
Gallium (Ga)	<10 ppt	Molybdenum (Mo)	<10 ppt
Terbium (Tb)	<1 ppt	Neodymium (Nd)	<1 ppt
Indium (In)	<1 ppt	Nickel (Ni)	<50 ppt
Dysprosium (Dy)	<1 ppt	Potassium (K)	<50 ppt
Aluminium (Al)	<50 ppt	Praseodymium (Pr)	<1 ppt
Antimony (Sb)	<50 ppt	Rhenium (Re)	<10 ppt
Arsenic (As)	<50 ppt	Rhodium (Rh)	<50 ppt
Barium (Ba)	<10 ppt	Ruthenium (Ru)	<50 ppt
Beryllium (Be)	<10 ppt	Samarium (Sm)	<1 ppt
Bismuth (Bi)	<10 ppt	Scandium (Sc)	<10 ppt
Cadmium (Cd)	<10 ppt	Propan--1-ol	<50 ppt
Calcium (Ca)	<50 ppt	Sodium (Na)	<100 ppt
Cerium (Ce)	<10 ppt	Strontium (Sr)	<10 ppt
Chromium (Cr)	<10 ppt	Tellurium (Te)	<1 ppt
Cobalt (Co)	<10 ppt	Thallium (Tl)	<10 ppt
Copper (Cu)	<10 ppt	Thorium (Th)	<1 ppt
Erbium (Er)	<1 ppt	Thulium (Tm)	<1 ppt
Europium (Eu)	<1 ppt	Tin (Sn)	<50 ppt
Gadolinium (Gd)	<1 ppt	Titanium (Ti)	<10 ppt
Germanium (Ge)	<10 ppt	Uranium (U)	<1 ppt
Hafnium (Hf)	<10 ppt	Vanadium (V)	<10 ppt
Holmium (Ho)	<1 ppt	Ytterbium (Yb)	<1 ppt
Iron (Fe)	<50 ppt	Yttrium (Y)	<1 ppt
Lanthanum (La)	<1 ppt	Zinc (Zn)	<50 ppt
Lead (Pb)	<10 ppt	Platinum (Pt)	<50 ppt
Lithium (Li)	<10 ppt	Rubidium (Rb)	<10 ppt
Lutetium (Lu)	<10 ppt	Tungsten (W)	<10 ppt
Magnesium (Mg)	<50 ppt		

Bottled in class 10 cleanroom conditions. Individually double-bagged in class 100 cleanroom conditions. Each bottle is individually packaged. Certificate of analysis delivered with each bottle.

Acetic acid glacial min.99%
for trace metal analysis

Trace Metal

packaging		Mfr. No
500 mℓ	HDPE plastic bottle	A507-P500
1 ℓ	HDPE plastic bottle	A507-P1
2.5 ℓ	HDPE plastic bottle	A507-P212
C ₂ H ₄ O ₂ CAS: 64-19-7 MW: 60.05 Density at 20°C: 1.048 Boiling point: 117 to 118°C Melting point: 16 to 16.5°C Flash point: 40°C ADR/RID: 8, 3, II		EINECS: 200-580-7 UN: 2789 GHS: Flam. Liq. 3, Skin Corr. 1A H226, H314 P280, P305+P351+P338, P302+P352, P304+P340, P210, P309+P311



Assay: >=99%			
Guaranteed analysis			
Colour	<10 APHA	Cesium (Cs)	<0.1 ppb
Aluminium (Al)	<1 ppb	Chromium (Cr)	<1 ppb
Antimony (Sb)	<0.5 ppb	Cobalt (Co)	<0.1 ppb
Arsenic (As)	<0.5 ppb	Copper (Cu)	<0.5 ppb
Barium (Ba)	<0.5 ppb	Dysprosium (Dy)	<0.1 ppb
Beryllium (Be)	<0.1 ppb	Erbium (Er)	<0.1 ppb
Bismuth (Bi)	<0.1 ppb	Europium (Eu)	<0.1 ppb
Cadmium (Cd)	<0.5 ppb	Gadolinium (Gd)	<0.1 ppb
Calcium (Ca)	<1 ppb	Gallium (Ga)	<0.1 ppb
Cerium (Ce)	<0.1 ppb		

Germanium (Ge)	<0.5 ppb	Scandium (Sc)	<0.1 ppb
Hafnium (Hf)	<0.1 ppb	Selenium (Se)	<1 ppb
Holmium (Ho)	<0.1 ppb	Silver (Ag)	<1 ppb
Indium (In)	<0.1 ppb	Sodium (Na)	<1 ppb
Iron (Fe)	<1 ppb	Strontium (Sr)	<0.5 ppb
Lanthanum (La)	<0.1 ppb	Sulfate (SO ₄)	<0.5 ppb
Lead (Pb)	<0.1 ppb	Tellurium (Te)	<0.5 ppb
Lithium (Li)	<0.1 ppb	Terbium (Tb)	<0.1 ppb
Lutetium (Lu)	<0.1 ppb	Thallium (Tl)	<0.1 ppb
Magnesium (Mg)	<0.5 ppb	Thorium (Th)	<0.1 ppb
Manganese (Mn)	<0.5 ppb	Thulium (Tm)	<0.1 ppb
Mercury (Hg)	<1 ppb	Tin (Sn)	<0.5 ppb
Molybdenum (Mo)	<0.5 ppb	Titanium (Ti)	<0.5 ppb
Neodymium (Nd)	<0.1 ppb	Total chloride (Cl)	<1 ppb
Nickel (Ni)	<0.5 ppb	Total phosphorus (P)	<1 ppb
Platinum (Pt)	<0.5 ppb	Tungsten (W)	<0.5 ppb
Potassium (K)	<1 ppb	Uranium (U)	<0.1 ppb
Praseodymium (Pr)	<0.1 ppb	Vanadium (V)	<0.5 ppb
Rhenium (Re)	<0.1 ppb	Ytterbium (Yb)	<0.1 ppb
Rhodium (Rh)	<0.5 ppb	Yttrium (Y)	<0.1 ppb
Rubidium (Rb)	<0.1 ppb	Zinc (Zn)	<1 ppb
Ruthenium (Ru)	<0.5 ppb	Zirconium (Zr)	<0.1 ppb
Samarium (Sm)	<0.1 ppb	Dichromate demand	Conform

Bottled in class 10 cleanroom conditions. Lot analysis on the label.

Acetic acid glacial
for trace metal analysis

Primar Plus

packaging		Mfr. No
250 mℓ	HDPE Plastic bottle	A/0411/PB07
500 mℓ	HDPE Plastic bottle	A/0411/PB08
1 ℓ	HDPE Plastic bottle	A/0411/PB15
2.5 ℓ	HDPE Plastic bottle	A/0411/PB17
C ₂ H ₄ O ₂ CAS: 64-19-7 MW: 60.04 Density at 20°C: 1.048 Boiling point: 117 to 118°C Melting point: 16 to 16.5°C Flash point: 40°C ADR/RID: 8, 3, II		EINECS: 200-580-7 UN: 2789 GHS: Flam. Liq. 3, Skin Corr. 1A H226, H314 P280, P305+P351+P338, P302+P352, P304+P340, P210, P309+P311



Assay: >99.8%			
Guaranteed analysis			
Aluminium (Al)	<0.005 ppm	Magnesium (Mg)	<0.005 ppm
Antimony (Sb)	<0.005 ppm	Manganese (Mn)	<0.002 ppm
Arsenic (As)	<0.005 ppm	Mercury (Hg)	<0.005 ppm
Barium (Ba)	<0.005 ppm	Molybdenum (Mo)	<0.002 ppm
Beryllium (Be)	<0.005 ppm	Nickel (Ni)	<0.002 ppm
Bismuth (Bi)	<0.005 ppm	Potassium (K)	<0.5 ppm
Boron (B)	<0.005 ppm	Selenium (Se)	<0.005 ppm
Cadmium (Cd)	<0.002 ppm	Strontium (Sr)	<0.002 ppm
Calcium (Ca)	<0.1 ppm	Thallium (Tl)	<0.005 ppm
Chromium (Cr)	<0.05 ppm	Tin (Sn)	<0.005 ppm
Cobalt (Co)	<0.002 ppm	Titanium (Ti)	<0.005 ppm
Copper (Cu)	<0.005 ppm	Total chloride (Cl)	<0.0001%
Gallium (Ga)	<0.005 ppm	Total phosphorus (P)	<0.5 ppm
Germanium (Ge)	<0.005 ppm	Total sulphur (S)	<0.1 ppm
Indium (In)	<0.005 ppm	Vanadium (V)	<0.002 ppm
Iron (Fe)	<0.05 ppm	Zinc (Zn)	<0.01 ppm
Lead (Pb)	<0.002 ppm	Zirconium (Zr)	<0.005 ppm
Lithium (Li)	<0.002 ppm	Residue after ignition	<0.0002%

Safe, environmentally friendly and effective packaging

Fisher Chemical outer packaging is UN compliant and 100% recyclable. It's an easy and safe way to reduce your waste management costs and contribute to a cleaner environment.

Acetic acid glacial
for analysis

Certified AR, Eur.Ph., BP ,USP

packaging		Mfr. No
500 mℓ	HDPE Plastic bottle	A/0400/PB08
1 ℓ	HDPE Plastic bottle	A/0400/PB15
1 ℓ	Coated Glass bottle	A/0400/PC15
2.5 ℓ	HDPE Plastic bottle	A/0400/PB17
2.5 ℓ	Coated Glass bottle	A/0400/PC17
25 ℓ	HDPE Plastic container	A/0400/25
45 ℓ	Plastic container	A/0400/26
C ₂ H ₄ O ₂ CAS: 64-19-7 MW: 60.05 Density at 20°C: 1.048 Boiling point: 117 to 118°C Melting point: 16 to 16.5°C Flash point: 40°C ADR/RID: 8, 3, II		EINECS: 200-580-7 UN: 2789 GHS: Flam. Liq. 3, Skin Corr. 1A H226, H314 P280, P305+P351+P338, P302+P352, P304+P340, P210, P309+P311



Assay: 99.7 to 100.5%			
Guaranteed analysis			
Colour	<10 APHA	Magnesium (Mg)	<0.05 ppm
Water	<0.2%	Manganese (Mn)	<0.01 ppm
Residual solvents test	: Pass	Mercury (Hg)	<0.02 ppm
Residue after evaporation	<5 ppm	Molybdenum (Mo)	<0.01 ppm
Substances Reducing K ₂ Cr ₂ O ₇	<0.005%	Nickel (Ni)	<0.02 ppm
Formic acid (H. COOH)	<0.05%	Potassium (K)	<0.5 ppm
Heavy metals	<5 ppm	Silver (Ag)	<0.01 ppm
Oxygen absorbed test	: Pass	Sodium (Na)	<1 ppm
Conforms to BP, EP & USP	: Pass	Strontium (Sr)	<0.01 ppm
Aluminium (Al)	<0.01 ppm	Tin (Sn)	<0.01 ppm
Arsenic (As)	<0.01 ppm	Titanium (Ti)	<0.01 ppm
Barium (Ba)	<0.01 ppm	Total chloride (Cl)	<0.0001%
Cadmium (Cd)	<0.01 ppm	Total nitrogen (N)	<2 ppm
Calcium (Ca)	<1 ppm	Total phosphorus (P)	<1 ppm
Chromium (Cr)	<0.2 ppm	Total silicon (Si)	>0 ppm
Cobalt (Co)	<0.01 ppm	Total sulphur (S)	<1 ppm
Copper (Cu)	<0.02 ppm	Vanadium (V)	<0.01 ppm
Iron (Fe)	<0.2 ppm	Zinc (Zn)	<0.05 ppm
Lead (Pb)	<0.01 ppm	Acetic anhydride	<0.2%
Lithium (Li)	<0.01 ppm		

Lot analysis on the label. Acetic acid glacial (d=1.08) is equivalent to a strength of 17.4M (17.4N).

Acetic acid glacial 99-100%
specified reagents for general laboratory work- extra pure

SLR, Eur.Ph.,BP, USP

packaging		Mfr. No
500 mℓ	HDPE Plastic bottle	A/0360/PB08
1 ℓ	HDPE Plastic bottle	A/0360/PB15
2.5 ℓ	HDPE Plastic bottle	A/0360/PB17
25 ℓ	HDPE Plastic container	A/0360/25
200 ℓ	HDPE Plastic container	A/0360/27
C ₂ H ₄ O ₂ CAS: 64-19-7 MW: 60.05 Density at 20°C: 1.048 Boiling point: 117 to 118°C Melting point: 16 to 16.5°C Flash point: 40°C ADR/RID: 8, 3, II		EINECS: 200-580-7 UN: 2789 GHS: Flam. Liq. 3, Skin Corr. 1A H226, H314 P280, P305+P351+P338, P302+P352, P304+P340, P210, P309+P311



Assay: >99%			
Guaranteed analysis			
Calcium (Ca)	<10 ppm	Potassium (K)	<10 ppm
Copper (Cu)	<1 ppm	Sodium (Na)	<20 ppm
Iron (Fe)	<2 ppm	n-Butyl formate	<0.0005%
Lead (Pb)	<2 ppm	Total sulphur (S)	<5 ppm
Magnesium (Mg)	<2 ppm	Zinc (Zn)	<1 ppm

Acetic acid glacial (d=1.08) is equivalent to a strength of 17.4M (17.4N).

Acetic acid glacial
pure, For non aqueous titrations

packaging		Mfr. No
2.5 ℓ	Coated Glass bottle	A/0420/PB17
C ₂ H ₄ O ₂ CAS: 64-19-7 MW: 60.05 Density at 20°C: 1.048 Boiling point: 117 to 118°C Melting point: 16 to 16.5°C Flash point: 40°C ADR/RID: 8, 3, II		EINECS: 200-580-7 UN: 2789 GHS: Flam. Liq. 3, Skin Corr. 1A H226, H314 P280, P305+P351+P338, P302+P352, P304+P340, P210, P309+P311



Guaranteed analysis			
Water	0.01 to 0.1%	Identity test	: Pass

Acetic acid glacial
pure, aldehyde free

packaging		Mfr. No
500 mℓ	HDPE Plastic bottle	A/0430/PB08
C ₂ H ₄ O ₂ CAS: 64-19-7 MW: 60.05 Density at 20°C: 1.048 Boiling point: 117 to 118°C Freezing point: >16.1°C Melting point: 16 to 16.5°C Flash point: 40°C		ADR/RID: 8, 3, II EINECS: 200-580-7 UN: 2789 GHS: Flam. Liq. 3, Skin Corr. 1A H226, H314 P280, P305+P351+P338, P302+P352, P304+P340, P210, P309+P311



Absorbance specification
Max. Absorbance (1cm cell v HPLC water)
570 nm <0.02 A.U.
Suitable for use in determination of total cholesterol in blood plasma. Ref. Macyntre & Ralston, Biochem J.1954.56,XLIII.
Acetic acid glacial (d=1.08) is equivalent to a strength of 17.4M (17.4N).

Acetic acid solution 0.1M (0.1N)
NIST Standard Concentrate, For Volumetric analysis

packaging		Mfr. No
1	Ampoule	J/0551C/05
6	Ampoule	J/0551C/90
EINECS: 200-580-7		

Each ampoule makes 1 ℓ of solution. See the full range of Standard Volumetric solution page 323.

Acetic acid solution 1M (1N approx.6%)
pure, ready to use

packaging		Mfr. No
2.5 ℓ	HDPE Plastic bottle	J/0550/17
EINECS: 200-580-7 GHS: Skin Corr. 1A H314 P280, P305+P351+P338, P301+P330+P331, P310		



Manufactured to a tolerance ± 0.002

Acetic acid solution 2M (2N approx.12%)
pure, ready to use

packaging		Mfr. No
2.5 ℓ	HDPE Plastic bottle	J/0549/PB17
CAS: 64-19-7 Density at 20°C: 1.0 Boiling point: 100°C ADR/RID: 8, III EINECS: 200-580-7		UN: 2790 GHS: Skin Irrit. 2, Eye Irrit. 2 H319, H315 P305+P351+P338, P280



Manufactured to a tolerance ± 0.002

Acetic anhydride for analysis		Certified AR	
packaging		Mfr. No	
500 mℓ	HDPE Plastic bottle	A/0480/PB08	
2.5 ℓ	HDPE Plastic bottle	A/0480/PB17	
C ₄ H ₆ O ₃	UN: 1715		
CAS: 108-24-7	GHS: Flam. Liq. 3, Acute		
MW: 102.09	Tox. 4, Skin Corr. 1B		
Density at 20°C: 1.087	H226, H330, H302, H314, EUH071		
Boiling point: 140°C	P301+P330+P331, P280, P305+P351+P338, P310, P304+P340, P261, P210		
Melting point: -73.1°C			
Flash point: 54°C			
ADR/RID: 8, 3, II			
EINECS: 203-564-8			
Assay: >98%			
Guaranteed analysis			
Colour	<10 APHA	Nickel (Ni)	<0.02 ppm
Substances		Potassium (K)	<0.5 ppm
reducing KMnO ₄	<0.01%	Sodium (Na)	<1 ppm
Calcium (Ca)	<0.5 ppm	Total chloride (Cl)	<0.0002%
Chromium (Cr)	<0.02 ppm	Total nitrogen (N)	<5 ppm
Copper (Cu)	<0.02 ppm	Total phosphorus (P)	<5 ppm
Iron (Fe)	<0.1 ppm	Total sulphur (S)	<2 ppm
Lead (Pb)	<0.02 ppm	Zinc (Zn)	<0.1 ppm
Magnesium (Mg)	<0.05 ppm		
Typical analysis			
Aluminium (Al)	<0.1 ppm	Manganese (Mn)	<0.01 ppm
Arsenic (As)	<0.01 ppm	Molybdenum (Mo)	<0.01 ppm
Barium (Ba)	<0.01 ppm	Strontium (Sr)	<0.01 ppm
Cadmium (Cd)	<0.01 ppm	Titanium (Ti)	<0.01 ppm
Cobalt (Co)	<0.01 ppm	Vanadium (V)	<0.01 ppm
Lithium (Li)	<0.01 ppm		
Lot analysis on the label.			

Acetic anhydride specified reagents for general laboratory work- extra pure		SLR
packaging		Mfr. No
500 mℓ	HDPE Plastic bottle	A/0440/PB08
2.5 ℓ	HDPE Plastic bottle	A/0440/PB17
C ₄ H ₆ O ₃	UN: 1715	
CAS: 108-24-7	GHS: Flam. Liq. 3, Acute	
MW: 102.09	Tox. 4, Skin Corr. 1B	
Density at 20°C: 1.087	H226, H330, H302, H314,	
Boiling point: 140°C	EUH071	
Melting point: -73.1°C	P301+P330+P331, P280,	
Flash point: 54°C	P305+P351+P338, P310,	
ADR/RID: 8, 3, II	P304+P340, P261, P210	
EINECS: 203-564-8		
Assay: >97%		
Guaranteed analysis		
Colour	<20 APHA	Sodium (Na) <20 ppm
Calcium (Ca)	<10 ppm	n-Butyl formate <0.005%
Copper (Cu)	<10 ppm	Total nitrogen (N) <50 ppm
Iron (Fe)	<10 ppm	Total phosphorus (P) <50 ppm
Lead (Pb)	<10 ppm	Total sulphur (S) <50 ppm
Magnesium (Mg)	<10 ppm	Zinc (Zn) <10 ppm
Potassium (K)	<10 ppm	

High Purity Acids for Elemental Analysis

Fisher Chemical offers a complete range of acids to meet the most challenging applications. Choose from our Optima®, TraceMetal® and PrimarPlus® grades to ensure that the sensitivity you require is always available.

Read more page 24



Acetone NEW PACK. for HPLC			Certified	
packaging		Mfr. No		
1 ℓ	Glass bottle	A/0606/15		
2.5 ℓ	Glass bottle	A/0606/17		
10 ℓ	Returnable Drum	A/0606/10R55		
25 ℓ	Returnable Drum	A/0606/25R55		
30 ℓ	Returnable Drum	A/0606/30R55		
C ₃ H ₆ O	EINECS: 200-662-2			
CAS: 67-64-1	UN: 1090			
MW: 58.08	GHS: Flam. Liq. 2, Eye Irrit.			
Density at 20°C: 0.790	2, STOT SE 3			
Boiling point: 56°C	H225, H336, H319, EUH066			
Melting point: -95°C	P261, P280,			
Flash point: -20°C	P305+P351+P338, P210, P240			
ADR/RID: 3, II				
Assay (CG): >99.8%				
Guaranteed analysis				
Acidity / Alkalinity		<0.0003 meq/g	Residue after evaporation	<2 ppm
Water		<0.1%		
Absorbance specification				
Max. Absorbance (1cm cell v HPLC water)				
330 nm		<1 A.U.	350 nm	<0.005 A.U.
335 nm		<0.3 A.U.	360 nm	<0.005 A.U.
340 nm		<0.09 A.U.	370 nm	<0.005 A.U.
Lot analysis on the label.				
Filtered to 0.2 µm.				

Acetone

for trace metal analysis

Primar

C ₃ H ₆ O	EINECS: 200-662-2		
CAS: 67-64-1	UN: 1090		
MW: 58.08	GHS: Flam. Liq. 2, Eye Irrit.		
Density at 20°C: 0.790	2, STOT SE 3		
Boiling point: 56°C	H225, H336, H319, EUH066		
Melting point: -95°C	P261, P280,		
Flash point: -20°C	P305+P351+P338, P210, P240		
ADR/RID: 3, II			
Assay (CG): >99.9%			
Guaranteed analysis			
Acidity / Alkalinity	<0.0005 meq/g	Manganese (Mn)	<0.005 ppm
Colour	<5 APHA	Molybdenum (Mo)	<0.005 ppm
Water	<0.2%	Nickel (Ni)	<0.005 ppm
Residue after evaporation	<2 ppm	Potassium (K)	<0.05 ppm
Substances		Propan--1-ol	<0.005 ppm
reducing KMnO ₄	<0.0005%	Sodium (Na)	<0.2 ppm
Aluminium (Al)	<0.02 ppm	Strontium (Sr)	<0.005 ppm
Barium (Ba)	<0.005 ppm	Titanium (Ti)	<0.005 ppm
Cadmium (Cd)	<0.005 ppm	Total phosphorus (P)	<0.02 ppm
Calcium (Ca)	<0.05 ppm	Total silicon (Si)	<0.02 ppm
Chromium (Cr)	<0.005 ppm	Total sulphur (S)	<0.1 ppm
Cobalt (Co)	<0.005 ppm	Vanadium (V)	<0.005 ppm
Copper (Cu)	<0.01 ppm	Zinc (Zn)	<0.01 ppm
Iron (Fe)	<0.05 ppm	Di-iso-propylether	<0.005%
Lead (Pb)	<0.005 ppm	Methanol	<0.05%
Lithium (Li)	<0.005 ppm	Methylpentadienes	<0.02%
Magnesium (Mg)	<0.02 ppm		

Acetone for residue and pesticide analysis			Distol
packaging		Mfr. No	
1 ℓ	Glass bottle	A/0603/15	
2.5 ℓ	Glass bottle	A/0603/17	
C ₃ H ₆ O	EINECS: 200-662-2		
CAS: 67-64-1	UN: 1090		
MW: 58.08	GHS: Flam. Liq. 2, Eye Irrit.		
Density at 20°C: 0.790	2, STOT SE 3		
Boiling point: 56°C	H225, H336, H319, EUH066		
Melting point: -95°C	P261, P280,		
Flash point: -20°C	P305+P351+P338, P210, P240		
ADR/RID: 3, II			

Guaranteed analysis			
Acidity / Alkalinity	<0.0005%	GC (ECD) Analysis	Pass
Water	<0.2%	GC (FID) Analysis	Pass
Residue after evaporation	<2 ppm	GC (NPD) analysis	Pass
GC Analysis			
Pesticides analysis: No impurity peak that has a retention time greater than Lindane will have a peak height greater than produced by 2500 pg / mℓ Heptachlorepoxide (ECD) and 10 pg / mℓ Parathion (NPD).			
Hydrocarbon analysis: no impurity chain lenght is greater than Dodecane, is present at a concentration greater than 0.01 ppm (by FID).			
Conditions : Column: J&W DB 5; 30m x 0.32mm x 0.1 μm dF. Temperature profile: (1) 50°C, 2 min (2) 40°C / min up to 320°C (3) 320°C, 30 min Injection: 5 μℓ cold on-column.			

Acetone		
extra dry, for synthesis		
packaging		Mfr. No
1 ℓ	Glass bottle	A/0602/15
2.5 ℓ	Glass bottle	A/0602/17
C ₃ H ₆ O	EINECS: 200-662-2	
CAS: 67-64-1	UN: 1090	
MW: 58.08	GHS: Flam. Liq. 2, Eye Irrit.	
Density at 20°C: 0.790	2, STOT SE 3	
Boiling point: 56°C	H225, H336, H319, EUH066	
Melting point: -95°C	P261, P280,	
Flash point: -20°C	P305+P351+P338, P210, P240	
ADR/RID: 3, II		
Guaranteed analysis		
Water	<0.01%	Residue after evaporation <2 ppm

Acetone			
for electronic use (MOS)			
packaging		Mfr. No	
2.5 ℓ	Glass bottle	A/0601/17	
5 ℓ	Aluminium can	A/0601/21	
25 ℓ	Metal drum	A/0601/25	
C ₃ H ₆ O		EINECS: 200-662-2	
CAS: 67-64-1		UN: 1090	
MW: 58.08		GHS: Flam. Liq. 2, Eye Irrit.	
Density at 20°C: 0.790		2, STOT SE 3	
Boiling point: 56°C		H225, H336, H319, EUH066	
Melting point: -95°C		P261, P280,	
Flash point: -20°C		P305+P351+P338, P210, P240	
ADR/RID: 3, II			
Assay (CG): >99.8%			
Guaranteed analysis			
Acidity / Alkalinity	<0.0005 meq/g	Magnesium (Mg)	<0.02 ppm
Colour	<5 APHA	Manganese (Mn)	<0.02 ppm
Water	<0.2%	Molybdenum (Mo)	<0.02 ppm
Residue after evaporation	<5 ppm	Nickel (Ni)	<0.02 ppm
Substances		Potassium (K)	<0.1 ppm
reducing KMnO ₄	<0.0005%	Sodium (Na)	<0.2 ppm
Aluminium (Al)	<0.02 ppm	Strontium (Sr)	<0.02 ppm
Barium (Ba)	<0.02 ppm	Titanium (Ti)	<0.02 ppm
Cadmium (Cd)	<0.02 ppm	Total phosphorus (P)	<0.02 ppm
Calcium (Ca)	<0.1 ppm	Total silicon (Si)	<0.02 ppm
Chromium (Cr)	<0.02 ppm	Vanadium (V)	<0.02 ppm
Cobalt (Co)	<0.02 ppm	Zinc (Zn)	<0.02 ppm
Copper (Cu)	<0.01 ppm	Di-iso-propylether	<0.005%
Iron (Fe)	<0.05 ppm	Methanol	<0.05%
Lead (Pb)	<0.02 ppm	Methylpentadienes	<0.02%
Lithium (Li)	<0.02 ppm		

Acetone NEW PACK. for analysis		Certified AR Eur.Ph	
packaging		Mfr. No	
500 ml	Coated Glass bottle	A/0600/PB08	
1 ℓ	Glass bottle	A/0600/15	
2.5 ℓ	Glass bottle	A/0600/17	
2.5 ℓ	Coated Glass bottle	A/0600/PB17	
2.5 ℓ	HDPE Plastic bottle	A/0600/PC17	
5 ℓ	Fluorinated HDPE Plastic bottle	A/0600/21	
5 ℓ	HDPE Plastic bottle	A/0600/PC21	
25 ℓ	Metal drum	A/0600/25	
200 ℓ	Metal drum	A/0600/27	
C ₃ H ₆ O	EINECS: 200-662-2		
CAS: 67-64-1	UN: 1090		
MW: 58.08	GHS: Flam. Liq. 2, Eye Irrit.		
Density at 20°C: 0.790	2, STOT SE 3		
Boiling point: 56°C	H225, H336, H319, EUH066		
Melting point: -95°C	P261, P280,		
Flash point: -20°C	P305+P351+P338, P210, P240		
ADR/RID: 3, II			
Assay (CG): >99.8%			
Guaranteed analysis			
Acidity / Alkalinity	<0.0005 meq/g	Iron (Fe)	<0.05 ppm
Colour	<5 APHA	Lead (Pb)	<0.02 ppm
Water	<0.2%	Magnesium (Mg)	<0.02 ppm
Residue after evaporation	<5 ppm	Potassium (K)	<0.1 ppm
Substances		Sodium (Na)	<0.2 ppm
reducing KMnO ₄	<0.0005%	Total phosphorus (P)	<0.02 ppm
Calcium (Ca)	<0.1 ppm	Zinc (Zn)	<0.02 ppm
Copper (Cu)	<0.01 ppm		
Typical analysis			
Aluminium (Al)	<0.02 ppm	Molybdenum (Mo)	<0.02 ppm
Barium (Ba)	<0.02 ppm	Nickel (Ni)	<0.02 ppm
Cadmium (Cd)	<0.02 ppm	Total silicon (Si)	<0.02 ppm
Chromium (Cr)	<0.02 ppm	Strontium (Sr)	<0.02 ppm
Cobalt (Co)	<0.02 ppm	Titanium (Ti)	<0.02 ppm
Lithium (Li)	<0.02 ppm	Vanadium (V)	<0.02 ppm
Manganese (Mn)	<0.02 ppm		
Lot analysis on the label.			

Acetone

NEW PACK.

specified reagents for general laboratory work- extra pure

SLR

packaging

Mfr. No

500 mℓ

Glass bottle

A/0560/08

1 ℓ

Glass bottle

A/0560/15

2.5 ℓ

Glass bottle

A/0560/17

2.5 ℓ

Coated Glass bottle

A/0560/PB17

2.5 ℓ

HDPE Plastic bottle

A/0560/PC17

5 ℓ

Fluorinated HDPE Plastic bottle

A/0560/21

5 ℓ

Fluorinated HDPE square Drum

A/0560/FP21

5 ℓ

HDPE Plastic bottle

A/0560/PC21

10 ℓ

HDPE Plastic container

A/0560/PC24

25 ℓ

Metal drum

A/0560/25

25 ℓ

HDPE Plastic container

A/0560/PC25

200 ℓ

Metal drum

A/0560/27

C₃H₆O

EINECS: 200-662-2

UN: 1090

GHS: Flam. Liq. 2, Eye Irrit.

2, STOT SE 3

H225, H336, H319, EUH066

P261, P280,

P305+P351+P338, P210, P240

CAS: 67-64-1

MW: 58.08

Density at 20°C: 0.790

Boiling point: 56°C

Melting point: -95°C

Flash point: -20°C

ADR/RID: 3, II

Assay (CG): >99%

Guaranteed analysis

Acidity / Alkalinity

<0.001 meq/g

Copper (Cu)

<2 ppm

Colour

<10 APHA

Iron (Fe)

<10 ppm

Water

<0.5%

Lead (Pb)

<5 ppm

Residue after evaporation

<50 ppm

Magnesium (Mg)

<1 ppm

Substances

Potassium (K)

<5 ppm

reducing KMnO₄

<0.005%

Sodium (Na)

<5 ppm

Calcium (Ca)

<1 ppm

Zinc (Zn)

<1 ppm

Acetone Technical		
packaging	Mfr. No	
2.5 ℓ Glass bottle	A/0520/17	
2.5 ℓ Coated Glass bottle	A/0520/PB17	
25 ℓ Metal drum	A/0520/25	
C ₃ H ₆ O	EINECS: 200-662-2	
CAS: 67-64-1	UN: 1090	
MW: 58.08	GHS: Flam. Liq. 2, Eye Irrit.	
Density at 20°C: 0.790	2, STOT SE 3	
Boiling point: 56°C	H225, H336, H319, EUH066	
Melting point: -95°C	P261, P280,	
Flash point: -20°C	P305+P351+P338, P210, P240	
ADR/RID: 3, II		
Assay (CG): > 98%		

Acetonitrile for LC/MS analysis	Optima
packaging	Mfr. No
500 ml Glass bottle	A955-500
1 ℓ Glass bottle	A955-1
2.5 ℓ Glass bottle	A955-212
4 ℓ Glass bottle	A955-4
C ₂ H ₃ N	EINECS: 200-835-2
CAS: 75-05-8	UN: 1648
MW: 41.04	GHS: Flam. Liq. 2, Acute
Density at 20°C: 0.785	Tox. 4, Eye Irrit. 2
Boiling point: 81 to 82°C	H225, H319, H332, H302, H312
Melting point: -46°C	P280, P302+P352,
Flash point: 2°C	P301+P312, P304+P340,
ADR/RID: 3, II	P305+P351+P338, P210, P240
Assay (CG): <99.9%	

Guaranteed analysis			
Acidity	0.008 meq/g	Iron (Fe)	<5 ppb
Alkalinity	<0.0006 meq/g	Lead (Pb)	<5 ppb
Water	<0.01%	Magnesium (Mg)	<10 ppb
Residue after evaporation	<0.8 ppm	Manganese (Mn)	<5 ppb
Aluminium (Al)	<25 ppb	Nickel (Ni)	<5 ppb
Barium (Ba)	<5 ppb	Potassium (K)	<10 ppb
Cadmium (Cd)	<5 ppb	Propan--1-ol	<5 ppb
Calcium (Ca)	<25 ppb	Sodium (Na)	<50 ppb
Chromium (Cr)	<5 ppb	Tin (Sn)	<5 ppb
Cobalt (Co)	<5 ppb	Zinc (Zn)	<10 ppb
Copper (Cu)	<5 ppb		

Absorbance specification			
Max. Absorbance (1cm cell v HPLC water)			
190 nm	<1 A.U.	280 nm	<0.005 A.U.
195 nm	<0.15 A.U.	LC/MS at Negative Mode as Aldicarb . . .	<
200 nm	<0.05 A.U.	50 ppb	
205 nm	<0.04 A.U.	LC/MS at Positive Mode as Reserpine . . .	<
210 nm	<0.03 A.U.	50 ppb	
215 nm	<0.025 A.U.	LC/UV Gradient (single max peak) at	
220 nm	<0.015 A.U.	210nm	<0.002 A.U.
225 nm	<0.015 A.U.	LC/UV Gradient (single max peak) at	
230 nm	<0.01 A.U.	254nm	<0.0005 A.U.
254 nm	<0.005 A.U.		

This product has been tested on an LC-MS Thermo Scientific instrument to ensure suitability.
Filtered to 0.2µm.

LC/MS Solvents Starter kit for LC/MS analysis	Optima
packaging	Mfr. No
Box Glass bottle	LCMSKIT
Containing acetonitrile, methanol, water and iso-propanol (1ℓ of each) for LC/MS analysis plus pre-treated empty bottle for blend	

Acetonitrile for HPLC-MS	Certified
packaging	Mfr. No
1 ℓ Glass bottle	A/0638/15
2.5 ℓ Glass bottle	A/0638/17
C ₂ H ₃ N	EINECS: 200-835-2
CAS: 75-05-8	UN: 1648
MW: 41.04	GHS: Flam. Liq. 2, Acute
Density at 20°C: 0.781	Tox. 4, Eye Irrit. 2
Boiling point: 81 to 82°C	H225, H319, H332, H302, H312
Melting point: -46°C	P280, P302+P352,
Flash point: 2°C	P301+P312, P304+P340,
ADR/RID: 3, II	P305+P351+P338, P210, P240

Assay (CG): >99.9%			
Guaranteed analysis			
Acidity	<0.0002 meq/g	Iron (Fe)	<0.01 ppm
Alkalinity	<0.0001 meq/g	Lead (Pb)	<0.01 ppm
Water	<0.01%	Magnesium (Mg)	<0.01 ppm
Residue after evaporation	<2 ppm	Manganese (Mn)	<0.01 ppm
Aluminium (Al)	<0.05 ppm	Nickel (Ni)	<0.01 ppm
Barium (Ba)	<0.01 ppm	Potassium (K)	<0.01 ppm
Cadmium (Cd)	<0.01 ppm	Propan--1-ol	<0.01 ppm
Calcium (Ca)	<0.01 ppm	Sodium (Na)	<0.05 ppm
Chromium (Cr)	<0.01 ppm	Tin (Sn)	<0.01 ppm
Cobalt (Co)	<0.01 ppm	Zinc (Zn)	<0.01 ppm
Copper (Cu)	<0.01 ppm		

Absorbance specification			
Max. Absorbance (1cm cell v HPLC water)			
200 nm	<0.03 A.U.	Gradient (max eluted peak). .	<0.2 mA.U.
230 nm	<0.005 A.U.	LC/MS Gradient Suitability	: Pass test
254 nm	<0.005 A.U.		

Gradient analysis:
Absorbance of the largest eluted peak: <0.2 mA.U.
Conditions :
Column: C18 Reverse Phase.
Flow rate: 2 ml / min
Wavelength 254 nm
Solvent A: Water
Solvent B: Acetonitrile
Gradient: linear 95% /5% B 100% B during 20 min.
Lot analysis on the label.
This product has been tested on an LC-MS instrument to ensure suitability.
Filtered to 0.2 µm.

Acetonitrile for UHPLC gradient grade analysis	Certified UHPLC
packaging	Mfr. No
1 ℓ Coated Glass bottle	A/0650/PB15
2.5 ℓ Coated Glass bottle	A/0650/PB17
C ₂ H ₃ N	EINECS: 200-835-2
CAS: 75-05-8	UN: 1648
MW: 41.04	GHS: Flam. Liq. 2, Acute
Density at 20°C: 0.781	Tox. 4, Eye Irrit. 2
Boiling point: 81 to 82°C	H225, H319, H332, H302, H312
Melting point: -46°C	P280, P302+P352,
Flash point: 2°C	P301+P312, P304+P340,
ADR/RID: 3, II	P305+P351+P338, P210, P240

Assay (CG): >99.9%			
Guaranteed analysis			
Acidity / Alkalinity	0.0002 meq/g	Residue after evaporation	<1 ppm
Water	<0.01%		

Absorbance specification			
Max. Absorbance (1cm cell v HPLC water)			
195 nm	<0.1 A.U.	230 nm	<0.005 A.U.
200 nm	<0.05 A.U.	240 nm	<0.005 A.U.
210 nm	<0.02 A.U.	250 nm	<0.005 A.U.
220 nm	<0.01 A.U.		

Gradient analysis:
Absorbance of the largest eluted peak at 210nm: <2mA.U.
Absorbance of the largest eluted peak at 220nm: <0.1mA.U.
Baseline Drift at 220nm: <8mA.U.
Conditions:
Column: C18 Reverse phase
Wavelength:210nm
Wavelength:220nm
Solvent A: Water

Solvent B: Acetonitrile
Gradient: Linear Gradient from 95% A/5%B to 100% B
Suitable for UHPLC instrument .
Lot analysis on the label.
Filtered to 0.1µm

UHPLC Solvents starter kit		NEW	Certified UHPLC
for UHPLC gradient grade analysis			
packaging			Mfr. No
Box	Coated Glass bottle		UHPLCKIT
Containing acetonitrile, methanol, water (1ℓ of each) for UHPLC gradient analysis			
also includes UHPLC solvents brochure			
Lot analysis on the label			

Acetonitrile HPLC for gradient analysis	Certified Advanced HPLC
packaging	Mfr. No
1 ℓ Coated Glass bottle	A/0632/PB15
2.5 ℓ Coated Glass bottle	A/0632/PB17
C ₂ H ₃ N	EINECS: 200-835-2
CAS: 75-05-8	UN: 1648
MW: 41.04	GHS: Flam. Liq. 2, Acute
Density at 20°C: 0.781	Tox. 4, Eye Irrit. 2
Boiling point: 81 to 82°C	H225, H319, H332, H302, H312
Melting point: -46°C	P280, P302+P352,
Flash point: 2°C	P301+P312, P304+P340,
ADR/RID: 3, II	P305+P351+P338, P210, P240

Assay (CG): >99.9%			
Guaranteed analysis			
Acidity / Alkalinity	<0.0008 meq/g	Residue after evaporation	<2 ppm
Water	<0.01%		

Absorbance specification			
Max. Absorbance (1cm cell v HPLC water)			
195 nm	<0.1 A.U.	230 nm	<0.005 A.U.
200 nm	<0.05 A.U.	240 nm	<0.005 A.U.
210 nm	<0.02 A.U.	250 nm	<0.005 A.U.
220 nm	<0.01 A.U.	Drift in gradient baseline	<8 mA.U.

Gradient analysis (30% to 95% of the solvent):
Absorbance of the largest eluted peak: <0.1 mA.U.
Conditions :
Column: C18 Reverse Phase.
Flow rate: 3 ml / min
Wavelength 220 nm
Solvent A: Water
Solvent B: Acetonitrile
Gradient: 10 min at 30% solvent B followed by a linear gradient to 95% solvent B over 15 min.
Lot analysis on the label.Filtered to 0.2µ.

Acetonitrile HPLC for gradient analysis	NEW PACK.	Certified, Eur.Ph
packaging	Mfr. No	
1 ℓ Glass bottle	A/0627/15	
2.5 ℓ Glass bottle	A/0627/17	
2.5 ℓ Coated Glass bottle	A/0627/PB17	
5 ℓ Aluminium can	A/0627/21	
10 ℓ Returnable Drum	A/0627/10RSS	
25 ℓ Returnable Drum	A/0627/25RSS	
30 ℓ Returnable Drum	A/0627/30RSS	
C ₂ H ₃ N	EINECS: 200-835-2	
CAS: 75-05-8	UN: 1648	
MW: 41.04	GHS: Flam. Liq. 2, Acute	
Density at 20°C: 0.781	Tox. 4, Eye Irrit. 2	
Boiling point: 81 to 82°C	H225, H319, H332, H302, H312	
Melting point: -46°C	P280, P302+P352,	
Flash point: 2°C	P301+P312, P304+P340,	
ADR/RID: 3, II	P305+P351+P338, P210, P240	

Assay (CG): >99.9%			
Guaranteed analysis			
Acidity / Alkalinity	<0.0008 meq/g	Residue after evaporation	<2 ppm
Water	<0.01%		

Absorbance specification			
Max. Absorbance (1cm cell v HPLC water)			
195 nm	<0.1 A.U.	230 nm	<0.005 A.U.
200 nm	<0.05 A.U.	240 nm	<0.005 A.U.
210 nm	<0.02 A.U.	250 nm	<0.005 A.U.
220 nm	<0.01 A.U.	Gradient (max eluted peak). .	<0.01 A.U.


Gradient analysis:
Absorbance of the largest eluted peak <0.01 AU.
Conditions :
Column: C18 Reverse Phase.
Flow rate: 2 ml / min
Wavelength 210 nm
Absorbance range: 0-0.32 AU.
Solvent A: Water
Solvent B: Acetonitrile
Gradient: linear 95% /5% B 100% B during 20 min.
Lot analysis on the label.
Filtered to 0.2 µm.

Acetonitrile for HPLC-RMN with low propionitrile level	Certified
packaging	Mfr. No
1 ℓ Glass bottle	A/0636/15
2.5 ℓ Glass bottle	A/0636/17
C ₂ H ₃ N	EINECS: 200-835-2
CAS: 75-05-8	UN: 1648
MW: 41.04	GHS: Flam. Liq. 2, Acute
Density at 20°C: 0.781	Tox. 4, Eye Irrit. 2
Boiling point: 81 to 82°C	H225, H319, H332, H302, H312
Melting point: -46°C	P280, P302+P352,
Flash point: 2°C	P301+P312, P304+P340,
ADR/RID: 3, II	P305+P351+P338, P210, P240

Assay (CG): >99.9%			
Guaranteed analysis			
Acidity / Alkalinity	<0.0008 meq/g	Residue after evaporation	<2 ppm
Water	<0.01%	Propionitrile	<0.0005%

Absorbance specification			
Max. Absorbance (1cm cell v HPLC water)			
195 nm	<0.1 A.U.	230 nm	<0.005 A.U.
200 nm	<0.05 A.U.	240 nm	<0.005 A.U.
210 nm	<0.02 A.U.	250 nm	<0.005 A.U.
220 nm	<0.01 A.U.		

Gradient analysis:
Absorbance of the largest eluted peak: <0.01A.U.
Conditions :
Column: C18 Reverse Phase.
Flow rate: 2 ml / min
Wavelength 210 nm
Solvent A: Water
Solvent B: Acetonitrile
Gradient: linear 95% /5% B 100% B during 20 min.
Lot analysis on the label.
Filtered to 0.2 µm.



Acetonitrile **NEW PACK.** Certified
for HPLC

packaging		Mfr. No
1 ℓ	Glass bottle	A/0626/15
2.5 ℓ	Glass bottle	A/0626/17
2.5 ℓ	Coated Glass bottle	A/0626/PB17
5 ℓ	Aluminium can	A/0626/21
10 ℓ	Returnable Drum	A/0626/10RSS
25 ℓ	Returnable Drum	A/0626/25RSS
30 ℓ	Returnable Drum	A/0626/30RSS
C ₂ H ₃ N CAS: 75-05-8 MW: 41.04 Density at 20°C: 0.781 Boiling point: 81 to 82°C Melting point: -46°C Flash point: 2°C ADR/RID: 3, II		EINECS: 200-835-2 UN: 1648 GHS: Flam. Liq. 2, Acute Tox. 4, Eye Irrit. 2 H225, H319, H332, H302, H312 P280, P302+P352, P301+P312, P304+P340, P305+P351+P338, P210, P240
Assay (CG): >99.8%		
Guaranteed analysis		
Acidity / Alkalinity	<0.0008 meq/g	Residue after evaporation <2 ppm
Water	<0.05%	
Absorbance specification		
Max. Absorbance (1cm cell v HPLC water)		
230 nm	<0.05 A.U.	270 nm <0.005 A.U.
240 nm	<0.015 A.U.	280 nm <0.005 A.U.
250 nm	<0.01 A.U.	290 nm <0.005 A.U.
260 nm	<0.005 A.U.	
Lot analysis on the label.		
Filtered to 0.2 µm.		

Assay (CG): >99.8%	
Guaranteed analysis	
Acidity / Alkalinity	<0.0008 meq/g
Water	<0.05%
Absorbance specification	
Max. Absorbance (1cm cell v HPLC water)	
230 nm	<0.05 A.U.
240 nm	<0.015 A.U.
250 nm	<0.01 A.U.
260 nm	<0.005 A.U.
Lot analysis on the label.	
Filtered to 0.2 µm.	

Acetonitrile Certified
HPLC for fluorescence detection

packaging		Mfr. No
1 ℓ	Coated Glass bottle	A/0630/PB15
2.5 ℓ	Coated Glass bottle	A/0630/PB17
C ₂ H ₃ N CAS: 75-05-8 MW: 41.04 Density at 20°C: 0.781 Boiling point: 81 to 82°C Melting point: -46°C Flash point: 2°C ADR/RID: 3, II		EINECS: 200-835-2 UN: 1648 GHS: Flam. Liq. 2, Acute Tox. 4, Eye Irrit. 2 H225, H319, H332, H302, H312 P280, P302+P352, P301+P312, P304+P340, P305+P351+P338, P210, P240
Assay (CG): >99.9%		
Guaranteed analysis		
Water	<0.02%	Residue after evaporation <2 ppm
Absorbance specification		
Max. Absorbance (1cm cell v HPLC water)		
200 nm	<0.1 A.U.	230 nm <0.01 A.U.
205 nm	<0.05 A.U.	240 nm <0.005 A.U.
210 nm	<0.035 A.U.	250 nm <0.005 A.U.
220 nm	<0.025 A.U.	Fluorescence Pass
Excitation in the range 250-750nm produces no fluorescence emission greater than the solvent natural Raman signal at 350nm excitation. This is equivalent to the emission maximum of approximately 1ppb quinine sulfate.		
Lot analysis on the label.		
Filtered to 0.2 µm (excluding fluorescence scan).		

Assay (CG): >99.9%	
Guaranteed analysis	
Water	<0.02%
Residue after evaporation	<2 ppm

Absorbance specification	
Max. Absorbance (1cm cell v HPLC water)	
200 nm	<0.1 A.U.
205 nm	<0.05 A.U.
210 nm	<0.035 A.U.
220 nm	<0.025 A.U.

Excitation in the range 250-750nm produces no fluorescence emission greater than the solvent natural Raman signal at 350nm excitation. This is equivalent to the emission maximum of approximately 1ppb quinine sulfate.

Lot analysis on the label.

Filtered to 0.2 µm (excluding fluorescence scan).

Acetonitrile Certified
HPLC for electrochemical detection

packaging		Mfr. No
1 ℓ	Coated Glass bottle	A/0635/PB15
2.5 ℓ	Coated Glass bottle	A/0635/PB17
C ₂ H ₃ N CAS: 75-05-8 MW: 41.04 Density at 20°C: 0.781 Boiling point: 81 to 82°C Melting point: -46°C Flash point: 2°C ADR/RID: 3, II		EINECS: 200-835-2 UN: 1648 GHS: Flam. Liq. 2, Acute Tox. 4, Eye Irrit. 2 H225, H319, H332, H302, H312 P280, P302+P352, P301+P312, P304+P340, P305+P351+P338, P210, P240
Assay (CG): >99.6%		

Guaranteed analysis	
Acidity / Alkalinity	<0.0008 meq/g
Water	<0.05%
Residue after evaporation	<2 ppm
Electrochemical test	Pass

Absorbance specification	
Max. Absorbance (1cm cell v HPLC water)	
230 nm	<0.1 A.U.
240 nm	<0.015 A.U.
250 nm	<0.01 A.U.

Purified then tested by HPLC-EC to show low electrochemical response.

Detector - ESA Coulochem II

Mode-oxidative

Graphite Electrode.

Lot analysis on the label.

Filtered to 0.2 µm.

Acetonitrile
for trace metal analysis

packaging		Mfr. No
2.5 ℓ	Glass bottle	A/0624/17
C ₂ H ₃ N CAS: 75-05-8 MW: 41.04 Density at 20°C: 0.781 Boiling point: 81 to 82°C Melting point: -46°C Flash point: 2°C ADR/RID: 3, II		EINECS: 200-835-2 UN: 1648 GHS: Flam. Liq. 2, Acute Tox. 4, Eye Irrit. 2 H225, H319, H332, H302, H312 P280, P302+P352, P301+P312, P304+P340, P305+P351+P338, P210, P240

Assay (CG): >99.9%	
Guaranteed analysis	
Acidity / Alkalinity	<0.0005 meq/g
Colour	<5 APHA
Water	<0.05%
Residue after evaporation	<2 ppm
Aluminium (Al)	<0.02 ppm
Barium (Ba)	<0.005 ppm
Cadmium (Cd)	<0.005 ppm
Calcium (Ca)	<0.05 ppm
Chromium (Cr)	<0.005 ppm
Cobalt (Co)	<0.005 ppm
Copper (Cu)	<0.01 ppm
Iron (Fe)	<0.05 ppm
Lead (Pb)	<0.005 ppm
Lithium (Li)	<0.005 ppm
Magnesium (Mg)	<0.02 ppm
Manganese (Mn)	<0.005 ppm
Molybdenum (Mo)	<0.005 ppm
Nickel (Ni)	<0.005 ppm
Potassium (K)	<0.05 ppm
Propan--1-ol	<0.005 ppm
Sodium (Na)	<0.2 ppm
Strontium (Sr)	<0.005 ppm
Titanium (Ti)	<0.005 ppm
Total phosphorus (P)	<0.02 ppm
Total silicon (Si)	<0.02 ppm
Total sulphur (S)	<0.1 ppm
Vanadium (V)	<0.005 ppm
Zinc (Zn)	<0.01 ppm
Allyl alcohol	<0.005%
Methanol	<0.005%
Propionitrile	<0.05%

Acetonitrile
for residue and pesticide analysis

packaging		Mfr. No
1 ℓ	Glass bottle	A/0623/15
2.5 ℓ	Glass bottle	A/0623/17
C ₂ H ₃ N CAS: 75-05-8 MW: 41.04 Density at 20°C: 0.781 Boiling point: 81 to 82°C Melting point: -46°C Flash point: 2°C ADR/RID: 3, II		EINECS: 200-835-2 UN: 1648 GHS: Flam. Liq. 2, Acute Tox. 4, Eye Irrit. 2 H225, H319, H332, H302, H312 P280, P302+P352, P301+P312, P304+P340, P305+P351+P338, P210, P240

Guaranteed analysis	
Water	<0.1%
Residue after evaporation	<2 ppm
GC (FID) Analysis	Pass
GC (ECD) Analysis	Pass

GC Analysis

Analysis of pesticides: from the retention time of Lindane, after concentration of 2500 times, no impurity peak higher than that obtained with 10 pg / mL heptachloro-oxide (ECD) and 5 pg / mL Parathion (NPD).

Analysis of hydrocarbons: no impurity chain longer than the dodecane is present at a concentration greater than 0.01 ppm (FID).

Conditions :

Column: J&W DB 5: 30m x 0.32mm x 0.1µm df.

Temperature profile:

(1) 50°C for 2 min

(2) 40°C / min until 320°C

(3) 320°C for 30 min

Injection: 5µl cold on the column.

Acetonitrile
for DNA synthesis and protein sequencing

packaging		Mfr. No
1 ℓ	Glass bottle	A/0628/15
2.5 ℓ	Glass bottle	A/0628/17
2.5 ℓ	Coated Glass bottle	A/0628/PB17
C ₂ H ₃ N CAS: 75-05-8 MW: 41.04 Density at 20°C: 0.781 Boiling point: 81 to 82°C Melting point: -46°C Flash point: 2°C ADR/RID: 3, II		EINECS: 200-835-2 UN: 1648 GHS: Flam. Liq. 2, Acute Tox. 4, Eye Irrit. 2 H225, H319, H332, H302, H312 P280, P302+P352, P301+P312, P304+P340, P305+P351+P338, P210, P240

Assay (CG): >99.9%	
Guaranteed analysis	
Alkalinity	<0.00006 meq/g
Water	<0.003%
Acidity	<0.00015 meq/g

Absorbance specification	
Max. Absorbance (1cm cell v HPLC water)	
190 nm	<1 A.U.
200 nm	<0.05 A.U.
230 nm	<0.01 A.U.
250 nm	<0.005 A.U.

Acetonitrile
extra dry, for synthesis

packaging		Mfr. No
1 ℓ	Glass bottle	A/0622/PB15
2.5 ℓ	Glass bottle	A/0622/PB17
C ₂ H ₃ N CAS: 75-05-8 MW: 41.04 Density at 20°C: 0.781 Boiling point: 81 to 82°C Melting point: -46°C Flash point: 2°C ADR/RID: 3, II		EINECS: 200-835-2 UN: 1648 GHS: Flam. Liq. 2, Acute Tox. 4, Eye Irrit. 2 H225, H319, H332, H302, H312 P280, P302+P352, P301+P312, P304+P340, P305+P351+P338, P210, P240

Guaranteed analysis	
Water	<0.01%
Residue after evaporation	<2 ppm

Acetonitrile **NEW PACK.** Certified AR
for analysis

packaging		Mfr. No
1 ℓ	Coated Glass bottle	A/0621/PB15
2.5 ℓ	Coated Glass bottle	A/0621/PB17
5 ℓ	Fluorinated HDPE Plastic bottle	A/0621/21
25 ℓ	Metal drum	A/0621/25
200 ℓ	Metal drum	A/0621/27
C ₂ H ₃ N CAS: 75-05-8 MW: 41.04 Density at 20°C: 0.781 Boiling point: 81 to 82°C Melting point: -46°C Flash point: 2°C ADR/RID: 3, II		EINECS: 200-835-2 UN: 1648 GHS: Flam. Liq. 2, Acute Tox. 4, Eye Irrit. 2 H225, H319, H332, H302, H312 P280, P302+P352, P301+P312, P304+P340, P305+P351+P338, P210, P240

Assay (CG): >99.5%	
Guaranteed analysis	
Acidity / Alkalinity	<0.0008 meq/g
Colour	<10 APHA
Water	<0.1%
Residue after evaporation	<10 ppm
Calcium (Ca)	<0.05 ppm
Copper (Cu)	<0.05 ppm
Iron (Fe)	<0.2 ppm
Lead (Pb)	<0.05 ppm
Magnesium (Mg)	<0.05 ppm
Potassium (K)	<0.05 ppm
Sodium (Na)	<1 ppm
Zinc (Zn)	<0.5 ppm

Typical analysis	
Aluminium (Al)	<0.02 ppm
Barium (Ba)	<0.02 ppm
Cadmium (Cd)	<0.02 ppm
Chromium (Cr)	<0.02 ppm
Cobalt (Co)	<0.02 ppm
Lithium (Li)	<0.02 ppm
Manganese (Mn)	<0.02 ppm
Molybdenum (Mo)	<0.02 ppm
Nickel (Ni)	<0.02 ppm
Total silicon (Si)	<0.02 ppm
Strontium (Sr)	<0.02 ppm
Titanium (Ti)	<0.02 ppm
Vanadium (V)	<0.02 ppm

Lot analysis on the label.

Acetonitrile **NEW PACK.** SLR
specified reagents for general laboratory work- extra pure

packaging		Mfr. No
500 mL	Coated Glass bottle	A/0620/PB08
2.5 ℓ	Coated Glass bottle	A/0620/PB17
5 ℓ	Fluorinated HDPE Plastic bottle	A/0620/21
25 ℓ	Metal drum	A/0620/25
200 ℓ	Metal drum	A/0620/27
C ₂ H ₃ N CAS: 75-05-8 MW: 41.04 Density at 20°C: 0.781 Boiling point: 81 to 82°C Melting point: -46°C Flash point: 2°C ADR/RID: 3, II		EINECS: 200-835-2 UN: 1648 GHS: Flam. Liq. 2, Acute Tox. 4, Eye Irrit. 2 H225, H319, H332, H302, H312 P280, P302+P352, P301+P312, P304+P340, P305+P351+P338, P210, P240

Assay (CG): >99%	
Guaranteed analysis	
Acidity / Alkalinity	<0.001 meq/g
Colour	<20 APHA
Water	<0.3%
Residue after evaporation	<50 ppm
Calcium (Ca)	<1 ppm
Copper (Cu)	<2 ppm
Iron (Fe)	<2 ppm
Lead (Pb)	<2 ppm
Magnesium (Mg)	<0.5 ppm
Potassium (K)	<1 ppm
Sodium (Na)	<5 ppm
Zinc (Zn)	<2 ppm

Acetonitrile with 0.1% formic acid (v/v) **NEW** Optima
for LC/MS analysis

packaging		Mfr. No
500 mL	Glass bottle	LS120-500
1 ℓ	Glass bottle	LS120-1
2.5 ℓ	Glass bottle	LS120-212
4 ℓ	Glass bottle	LS120-4
Density at 20°C: 0.781 Flash point: 6°C ADR/RID: 3, II EINECS: 200-835-2 UN: 1648 GHS: Flam. Liq. 2, Acute		Tox. 4, Eye Irrit. 2 H225, H302, H312, H332, H319 P280, P302+P352, P301+P312, P304+P340, P305+P351+P338, P210, P240

Guaranteed analysis	
Assay (v/v)	: 0.095 to 0.105%
Colour	<10 APHA
Water	<0.01%
Residue after evaporation	<1 ppm
Aluminium (Al)	<25 ppb
Calcium (Ca)	<50 ppb
Copper (Cu)	<10 ppb
Iron (Fe)	<10 ppb
Lead (Pb)	<10 ppb
Magnesium (Mg)	<10 ppb
Manganese (Mn)	<10 ppb
Nickel (Ni)	<10 ppb
Potassium (K)	<20 ppb
Propan--1-ol	<10 ppb
Sodium (Na)	<50 ppb
Zinc (Zn)	<20 ppb

Absorbance specification	
Max. Absorbance (1cm cell v HPLC water)	
210 nm	<1.3 A.U.
220 nm	<1.25 A.U.
230 nm	<0.75 A.U.
254 nm	<0.03 A.U.
LC/MS Gradient Suitability : Pass test	

This product has been tested on an LC-MS Thermo Scientific instrument to ensure suitability.

Filtered to 0.2µm.

We Support Your Chemistry For

- Semi-bulk and bulk requirements
- Tailored solvents and solvent blends or customized solutions
- Testing services or customized packaging and labeling

Please email your contact details and a brief description of your requirements to begei.shi.quotes@thermofisher.com

Acetonitrile with 0.1% formic acid (v/v)

NEW

Certified

HPLC for gradient analysis

packaging	Mfr. No
2.5 ℓ Glass bottle	FE/1100/17
Density at 20°C: 0.781 Flash point: 6°C ADR/RID: 3, II EINECS: 200-835-2 UN: 1648 GHS: Flam. Liq. 2, Acute	Tox. 4, Eye Irrit. 2 H225, H302, H312, H332, H319 P280, P302+P352, P301+P312, P304+P340, P305+P351+P338, P210, P240
Guaranteed analysis Residue after evaporation <=1 ppm Formic acid (% v/v) : 0.095 to 0.105%	
Absorbance specification Max. Absorbance (1cm cell v HPLC water) 210 nm <1.3 A.U. 230 nm <0.8 A.U. 240 nm <0.3 A.U.	
Lot analysis on the label. Filtered to 0.2µm.	

Acetonitrile with 0.1% trifluoroacetic acid (v/v)

NEW

Optima

for LC/MS analysis

packaging	Mfr. No
500 mℓ Glass bottle	LS121-500
1 ℓ Glass bottle	LS121-1
2.5 ℓ Glass bottle	LS121-212
4 ℓ Glass bottle	LS121-4
CAS: 76-05-1 Density at 20°C: 0.781 Boiling point: 82°C Melting point: -50°C Flash point: 6°C ADR/RID: 3, II EINECS: 200-929-3	UN: 1648 GHS: Flam. Liq. 2, Acute Tox. 4, Eye Irrit. 2 H225, H302, H312, H319, H332 P280, P301+P312, P304+P341, P302+P352, P305+P351+P338, P210, P240
Guaranteed analysis Assay (v/v) : 0.095 to0.105% Colour <10 APHA Water <0.01% Residue after evaporation <1 ppm Aluminium (Al) <25 ppb Calcium (Ca) <50 ppb Copper (Cu) <10 ppb Iron (Fe) <10 ppb	
Absorbance specification Max. Absorbance (1cm cell v HPLC water) 210 nm <0.6 A.U. 220 nm <0.55 A.U. 230 nm <0.4 A.U.	
This product has been tested on an LC-MS Thermo Scientific instrument to ensure suitability. Filtered to 0.2µm.	

Acetonitrile with 0.1% trifluoroacetic acid (v/v)

NEW

Certified

HPLC for gradient analysis

packaging	Mfr. No
2.5 ℓ Glass bottle	FE/0200/17
CAS: 76-05-1 Density at 20°C: 0.781 Boiling point: 82°C Melting point: -50°C Flash point: 6°C ADR/RID: 3, II EINECS: 200-929-3	UN: 1648 GHS: Flam. Liq. 2, Acute Tox. 4, Eye Irrit. 2 H225, H302, H312, H319, H332 P280, P301+P312, P304+P341, P302+P352, P305+P351+P338, P210, P240
Guaranteed analysis Residue after evaporation <=1 ppm Gradient analysis at 254 nm. : <0.01 A.U.	
Absorbance specification Max. Absorbance (1cm cell v HPLC water) 190 nm : <1 A.U. 200 nm <0.5 A.U. 210 nm <0.36 A.U. 230 nm <0.3 A.U.	
Lot analysis on the label. Filtered to 0.2µm.	

Acetophenone

specified reagents for general laboratory work- extra pure

SLR

packaging	Mfr. No
500 Coated Glass bottle	A/0640/PB08
C ₆ H ₅ COCH ₃ CAS: 98-86-2 MW: 120.15 Density at 20°C: 1.03 Boiling point: 202°C Melting point: 19 to 20°C	Flash point: 76°C EINECS: 202-708-7 GHS: Acute Tox. 4, Eye Irrit. 2 H302, H319 P280, P302+P352, P301+P312, P305+P351+P338
Assay (CG): >99%	

2-Acetoxybenzoic acid

See Acetylsalicylic acid, page 23

Acetylacetone

for analysis

Certified AR

packaging	Mfr. No
100 mℓ Coated Glass bottle	A/0802/PB05
C ₅ H ₈ O ₂ CAS: 123-54-6 MW: 100.12 Density at 20°C: 0.975 Boiling point: 133 to 135°C Melting point: -23°C Flash point: 35.5°C	ADR/RID: 3, 6.1, III EINECS: 204-634-0 UN: 2310 GHS: Flam. Liq. 3, Acute Tox. 4 H226, H302 P280, P301+P312, P210
Assay (CG): >99%	
Guaranteed analysis Colour <20 APHA Water <0.1% Residue after evaporation <500 ppm Calcium (Ca) <0.5 ppm Copper (Cu) <1 ppm Iron (Fe) <1 ppm	
Typical analysis Aluminium (Al) <0.05 ppm Barium (Ba) <0.05 ppm Cadmium (Cd) <0.05 ppm Chromium (Cr) <0.05 ppm Cobalt (Co) <0.05 ppm Lithium (Li) <0.05 ppm Manganese (Mn) <0.05 ppm	
Lot analysis on the label	

Acetylacetone

specified reagents for general laboratory work- extra pure

SLR

packaging	Mfr. No
500 mℓ Coated Glass bottle	A/0800/PB08
C ₅ H ₈ O ₂ CAS: 123-54-6 MW: 100.12 Density at 20°C: 0.975 Boiling point: 133 to 135°C Melting point: -23°C Flash point: 35.5°C	ADR/RID: 3, 6.1, III EINECS: 204-634-0 UN: 2310 GHS: Flam. Liq. 3, Acute Tox. 4 H226, H302 P280, P301+P312, P210
Assay (CG): >98%	
Guaranteed analysis Colour <30 APHA Water <0.5% Residue after evaporation <100 ppm Calcium (Ca) <2 ppm Copper (Cu) <5 ppm Iron (Fe) <5 ppm	
Lead (Pb) <5 ppm Magnesium (Mg) <0.5 ppm Potassium (K) <0.5 ppm Sodium (Na) <10 ppm Zinc (Zn) <0.5 ppm	

Acetyl chloride

for analysis

Certified AR

packaging	Mfr. No
100 mℓ Coated Glass bottle 500 mℓ Coated Glass bottle	A/0920/PB05 A/0920/PB08
C ₂ H ₃ ClO CAS: 75-36-5 MW: 78.50 Density at 20°C: 1.1 Boiling point: 51°C Melting point: -112°C Flash point: 4°C ADR/RID: 3, 8, II EINECS: 200-865-6	UN: 1717 GHS: Flam. Liq. 2, Skin Corr. 1B, Eye Dam. 1 H225, H314, EUH014 P301+P330+P331, P280, P305+P351+P338, P310, P402+P404, P301+P312, P210, P240
Assay: >99%	
Guaranteed analysis Colour <20 APHA Residue after evaporation <0.01% Calcium (Ca) <2 ppm Copper (Cu) <2 ppm Iron (Fe) <2 ppm Lead (Pb) <2 ppm	
Magnesium (Mg) <2 ppm Potassium (K) <2 ppm Sodium (Na) <2 ppm Total phosphorus (P) <10 ppm Total sulphur (S) <2 ppm Zinc (Zn) <2 ppm	
Lot analysis on the label.	

Acetyl chloride

specified reagents for general laboratory work- extra pure

SLR

packaging	Mfr. No
500 mℓ Coated Glass bottle	A/0880/PB08
C ₂ H ₃ ClO CAS: 75-36-5 MW: 78.50 Density at 20°C: 1.1 Boiling point: 51°C Melting point: -112°C Flash point: 4°C ADR/RID: 3, 8, II EINECS: 200-865-6	UN: 1717 GHS: Flam. Liq. 2, Skin Corr. 1B, Eye Dam. 1 H225, H314, EUH014 P301+P330+P331, P280, P305+P351+P338, P310, P402+P404, P301+P312, P210, P240
Assay: >98%	
Guaranteed analysis Residue after evaporation <0.05% Calcium (Ca) <10 ppm Copper (Cu) <10 ppm Iron (Fe) <10 ppm Lead (Pb) <10 ppm Magnesium (Mg) <10 ppm	
Potassium (K) <10 ppm Sodium (Na) <10 ppm Total phosphorus (P) <50 ppm Total sulphur (S) <20 ppm Zinc (Zn) <10 ppm	

Acetyl chloride

technical

packaging	Mfr. No
5 Ampoule	A/0881/99
C ₂ H ₃ ClO CAS: 75-36-5 MW: 78.50 Density at 20°C: 1.1 Boiling point: 51°C Melting point: -112°C Flash point: 4°C ADR/RID: 3, 8, II EINECS: 200-865-6	UN: 1717 GHS: Flam. Liq. 2, Skin Corr. 1B, Eye Dam. 1 H225, H314, EUH014 P301+P330+P331, P280, P305+P351+P338, P310, P402+P404, P301+P312, P210, P240

Standard Volumetric solution

For your convenience we offer a range of NIST traceable solutions. A full list of our products can be found on page 323. Can't see what you need? Please email your contact details and a brief description of your requirements to begei.sha.quotes@thermofisher.com

Acetylsalicylic acid

specified reagents for general laboratory work- extra pure

SLR

packaging	Mfr. No
500 g Plastic jar	A/0960/53
C ₉ H ₈ O ₄ CAS: 50-78-2 MW: 180.16 Density at 20°C: 1.35 Melting point: 136 to 140°C Flash point: 250°C ADR/RID: 6.1, III EINECS: 200-064-1	UN: 2811 GHS: Acute Tox. 3, Skin Irrit. 2, Eye Irrit. 2, STOT SE 3 H301, H335, H315, H319 P301+P310, P261, P302+P352, P280, P305+P351+P338
Assay: >99.5%	
Guaranteed analysis Residue after ignition <0.1% Arsenic (As) <2 ppm	
Lead (Pb) <10 ppm Salicylic acid <0.05%	



OPTIMA® grade for Ultra Trace Elemental analysis

as low as 1-100 ppt level detection, tested up to 65 parameters
For ICP-MS & ICP equipment
Bottled in class 10 cleanroom conditions. Individually double-bagged in class 100 cleanroom conditions.
Certificate of analysis delivered with each bottle.

designation	packsize	package type	Mfr. No	page
Acetic acid glacial	250 mL	Pre-cleaned Teflon FEP bottle	A465-250	14
	500 mL	Pre-cleaned Teflon FEP bottle	A465-500	14
	1 L	Pre-cleaned Teflon FEP bottle	A465-1	14
Ammonia solution 20%	250 mL	Pre-cleaned HDPE bottle	A470-250	33
	500 mL	Pre-cleaned HDPE bottle	A470-500	33
	1 L	Pre-cleaned HDPE bottle	A470-1	33
		Pre-cleaned HDPE bottle	A470-1	33
Hydrobromic acid 44-49%	500 mL	Pre-cleaned Teflon FEP bottle	A471-500	146
Hydrochloric acid 32-35%	250 mL	Teflon PFA bottle	A466-250	147
	500 mL	Pre-cleaned Teflon PFA bottle	A466-500	147
	1 L	Pre-cleaned Teflon PFA bottle	A466-1	147
Hydrofluoric acid 47-51%	250 mL	Pre-cleaned Teflon PFA bottle	A463-250	152
	500 mL	Pre-cleaned Teflon PFA bottle	A463-500	152
	1 L	Pre-cleaned Teflon PFA bottle	A463-1	152
Hydrogen peroxyde 30-32%	500 mL	Pre-cleaned Teflon FEP bottle	P170-500	154
Nitric acid 67-69%	250 mL	Pre-cleaned Teflon FEP bottle	A467-250	215
	500 mL	Pre-cleaned Teflon FEP bottle	A467-500	215
	1 L	Pre-cleaned Teflon FEP bottle	A467-1	215
Perchloric acid 65-71%	250 mL	Pre-cleaned Teflon FEP bottle	A469-250	228
	500 mL	Pre-cleaned Teflon FEP bottle	A469-500	228
	1	Pre-cleaned Teflon FEP bottle	A469-1	228
Sulfuric acid 93-98%	250 mL	Pre-cleaned Teflon FEP bottle	A468-250	336
	500 mL	Pre-cleaned Teflon FEP bottle	A468-500	336
	1 L	Pre-cleaned Teflon FEP bottle	A468-1	336
Water	500 mL	LDPE bottle	W9-500	368
	1 L	LDPE bottle	W9-1	368
	2 L	LDPE bottle	W9-2	368

PRIMAR PLUS® grade for Trace Metal analysis **NEW PACK.**

at 1-10 ppb level detection, tested up to 40 parameters
For AAS equipment
Lot analysis on the label.

designation	packsize	package type	Mfr. No	page
Acetic acid glacial	250 mL	HDPE Plastic bottle	A/0411/PB07	14
	500 mL	HDPE Plastic bottle	A/0411/PB08	14
	1 L	HDPE Plastic bottle	A/0411/PB15	14
	2.5 L	HDPE Plastic bottle	A/0411/PB17	14
Hydrochloric acid 37% d=1.18	250 mL	HDPE Plastic bottle	H/1196/PB07	148
	500 mL	HDPE Plastic bottle	H/1196/PB08	148
	1 L	HDPE Plastic bottle	H/1196/PB15	148
	2.5 L	HDPE Plastic bottle	H/1196/PB17	148
Nitric acid 68% d =1.42	500 mL	HDPE Plastic bottle	N/2272/PB08	216
	1 L	HDPE Plastic bottle	N/2272/PB15	216
	2.5 L	HDPE Plastic bottle	N/2272/PB17	216
Sulfuric acid 94-98% d=1.83	250 mL	HDPE Plastic bottle	S/9231/PB07	337
	500 mL	HDPE Plastic bottle	S/9231/PB08	337
	1 L	HDPE Plastic bottle	S/9231/PB15	337
	2.5 L	HDPE Plastic bottle	S/9231/PB17	337

TRACE METAL® grade for Trace Metal analysis **NEW**

at 1 ppb level detection, tested up to 65 parameters
For ICP-MS & ICP equipment
Bottled in class 10 cleanroom conditions.
Lot analysis on the label.

designation	packsize	package type	Mfr. No	page
Acetic acid glacial	500 mL	HDPE plastic bottle	A507-P500	14
	1 L	HDPE plastic bottle	A507-P1	14
	2.5 L	HDPE plastic bottle	A507-P212	14
Ammonia solution 20%	500 mL	HDPE plastic bottle	A512-P500	33
Hydrochloric acid 34-37%	500 mL	HDPE plastic bottle	A508-P500	147
	1 L	HDPE plastic bottle	A508-P1	147
	2.5 L	HDPE plastic bottle	A508-P212	147
Hydrofluoric acid 47-51%	500 mL	HDPE plastic bottle	A513-P500	152
Nitric acid 67-69%	500 mL	HDPE plastic bottle	A509-P500	216
	1 L	HDPE plastic bottle	A509-P1	216
	2.5 L	HDPE plastic bottle	A509-P212	216
Perchloric acid 65-71%	500 mL	HDPE plastic bottle	A511-P500	228
		Coated Glass bottle	P/1292/PB08	228
	1 L	HDPE plastic bottle	A511-P1	228
		Coated Glass bottle	P/1292/PB15	228
		HDPE plastic bottle	A511-P212	228
Sulfuric acid 93-98%	500 mL	HDPE plastic bottle	A510-P500	336
	1 L	HDPE plastic bottle	A510-P1	336
	2.5 L	HDPE plastic bottle	A510-P212	336

Acrylamide
for electrophoresis, unstabilized

packaging	Mfr. No
250 g Plastic jar	A/0977/50
1 kg Plastic jar	A/0977/60
25 kg Plastic keg	A/0977/70

C₃H₅NO
CAS: 79-06-1
MW: 71.08
Density at 30°C: 1.122
Boiling point: 125°C
Melting point: 82 to 86°C
Flash point: 138°C
ADR/RID: 6.1, III
EINECS: 201-173-7
UN: 2074
GHS: Carc. 1B, Muta. 1B,

Repr. 2, Acute Tox. 3, Skin Irrit. 2, Eye Irrit. 2, Skin Sens. 1, STOT RE 1
H315, H319, H301, H317, H332, H350, H340, H372, H361f, H312
P301+P310, P280, P302+P350, P201, P308+P313, P260, P304+P340, P305+P351+P338



Acrylamide
extra pure, unstabilized

SLR

packaging	Mfr. No
250 g Plastic jar	A/0976/50

C₃H₅NO
CAS: 79-06-1
MW: 71.08
Density at 30°C: 1.122
Boiling point: 125°C
Melting point: 82 to 86°C
Flash point: 138°C
ADR/RID: 6.1, III
EINECS: 201-173-7
UN: 2074
GHS: Carc. 1B, Muta. 1B,

Repr. 2, Acute Tox. 3, Skin Irrit. 2, Eye Irrit. 2, Skin Sens. 1, STOT RE 1
H315, H319, H301, H317, H332, H350, H340, H372, H361f, H312
P301+P310, P280, P302+P350, P201, P308+P313, P260, P304+P340, P305+P351+P338



Assay (CG): >98.5%

Adipic acid
pure

packaging	Mfr. No
500 g Plastic jar	A/1000/53

C₆H₁₀O₄
CAS: 124-04-9
MW: 146.15
Boiling point: 265°C
Melting point: 151 to 153°C

Flash point: 196°C
EINECS: 204-673-3
GHS: Eye Irrit. 2
H319
P280, P305+P351+P338



Assay: >99%

Adipoyl chloride
pure

packaging	Mfr. No
5 mL Glass bottle	A/1020/99

C₆H₈Cl₂O₂
CAS: 111-50-2
MW: 183.04
Density at 20°C: 1.25
Boiling point: 105 to 107°C
Flash point: 72°C
ADR/RID: 8, 6.1, II
EINECS: 203-876-4
UN: 2922

GHS: Acute Tox. 4, Skin Corr. 1B, Eye Dam. 1
H314, H312, H332, H302, EUH029, EUH014
P280, P301+P330+P331, P305+P351+P338, P310, P302+P352, P402+P404, P301+P312, P304+P340



UHPLC

Filtered to 0.1 µm, Fisher Chemical UHPLC gradient grade solvents (Acetonitrile, Methanol and Water) are perfect for those working at high pressure.
See page 309 for more information



Agar
technical, granular powder

packaging	Mfr. No
100 g Plastic jar	A/1080/48
500 g Plastic jar	A/1080/53
10 kg Plastic keg	A/1080/68

CAS: 9002-18-0
EINECS: 232-658-1

Guaranteed analysis
Identity test : Pass

L-Alanine
pure

packaging	Mfr. No
100 g Plastic jar	A/1162/48

C₃H₇NO₂
CAS: 56-41-7
MW: 89.09

Density at 20°C: 1.401
Melting point: 314 to 316°C
EINECS: 200-273-8

Assay: >99%

Guaranteed analysis
Loss on drying <0.1% Iron (Fe) <10 ppm
Specific rotation 14° to 16° Lead (Pb) <10 ppm
Residue after ignition <0.1%

DL-Alanine
pure

packaging	Mfr. No
100 g Plastic jar	A/1160/48

C₃H₇NO₂
CAS: 302-72-7
MW: 89.09

Melting point: 272 to 275°C
EINECS: 206-126-4

Guaranteed analysis
Complies with standard IR spectrum: Pass

Albumin, egg
specified reagents for general laboratory work- extra pure

SLR

packaging	Mfr. No
100 g Plastic jar	A/1280/48
250 g Plastic jar	A/1280/50

CAS: 9006-50-2
Density at 20°C: 1.0
Melting point: 61°C

Guaranteed analysis
Identity test : Pass

Albumin, bovine (fraction V)
extra pure, powder

SLR

packaging	Mfr. No
25 g Plastic jar	A/1278/46

CAS: 9048-46-8
EINECS: 232-936-2

Guaranteed analysis
Loss on drying <6% pH (aqueous solution 1%) 6 to 8
Protein content : >96%

pH Buffers

NIST Certified Standard Solution & Concentrate
A full range to discover on page 63



Alizarin C.I. 58000
pure, indicator grade

packaging	Mfr. No
10 g Glass jar	A/1460/45

C₁₄H₈O₄
CAS: 72-48-0
MW: 240.22
Boiling point: 430°C
Melting point: 287 to 289°C

EINECS: 200-782-5
GHS: Acute Tox. 4
H302
P301+P312



Guaranteed analysis
pH colour change test : Pass Complies with standard IR spectrum: Pass
pH 11.0-13.0 pink to violet.

Alizarin fluorine blue

See 3-Aminomethylalizarin-NN-diacetic acid, page 32

Alizarin red S, C.I.58005
pure

packaging	Mfr. No
25 g Glass jar	A/1520/46

C₁₄H₇NaO₇S
CAS: 130-22-3
MW: 342.26
EINECS: 204-981-8

Guaranteed analysis
pH colour change test : Pass Complies with standard IR spectrum: Pass

Alkyd resin
pure, solution 70 % wt

packaging	Mfr. No
1 l Glass bottle	A/1552/15

Density at 20°C: 0.96
Flash point: 41°C
ADR/RID: 3, III
UN: 1866
GHS: Flam. Liq. 3, Asp. Tox.

1, Aquatic Chronic 2
H226, H304, H411
P301+P310, P331, P280, P210, P240, P273



Guaranteed analysis
Complies with standard IR spectrum: Pass

Allylthiourea
pure, Determination of biological oxygen demand in water

packaging	Mfr. No
50 g Plastic jar	A/1566/47

C₃H₈N₂S
CAS: 109-57-9
MW: 116.19
Melting point: >70°C
ADR/RID: 6.1, III

EINECS: 203-683-5
UN: 2811
GHS: Acute Tox. 3
H301
P264, P270, P301+P310, P280



Guaranteed analysis
Melting Point tested >70°C

Almond oil

SLR

packaging	Mfr. No
100 mL Glass bottle	A/1565/05

Guaranteed analysis
Identity test : Pass

alpha-Amylase bacterial
technical, 50-100 units/mg

packaging	Mfr. No
5 g Glass jar	A/7190/44

CAS: 9000-85-5
EINECS: 232-560-9
GHS: Resp. Sens. 1

H334
P261, P304+P341, P342+P311



Guaranteed analysis
Complies with standard IR spectrum: Pass

Alumina

See "Aluminium oxide, page 30

Aluminium standard metal solution 1000 ppm
for AAS, solution HNO₃ 1M

packaging	Mfr. No
100 mL HDPE Plastic bottle	J/8000/05
500 mL HDPE Plastic bottle	J/8000/08

ADR/RID: 8, II
EINECS: 236-751-8
UN: 2031
GHS: Skin Corr. 1B, Eye Dam. 1

H314
P280, P305+P351+P338, P301+P330+P331, P310



Manufactured to a tolerance of ± 0.5%.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.

Aluminium standard metal solution 1000 ppm
for ICP, solution HNO₃ 1M

packaging	Mfr. No
500 mL HDPE Plastic bottle	J/8206/08

ADR/RID: 8, II
EINECS: 236-751-8
UN: 2031
GHS: Skin Corr. 1B, Eye Dam. 1

H314
P280, P305+P351+P338, P301+P330+P331, P310



Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.

Aluminium standard metal solution 10 000 ppm
for ICP, solution HNO₃ 1M

packaging	Mfr. No
100 mL HDPE Plastic bottle	J/8205/05

ADR/RID: 8, II
UN: 2031
GHS: Skin Corr. 1B, Eye Dam. 1

H314
P280, P305+P351+P338, P301+P330+P331, P310



Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.

Safe, environmentally friendly and effective packaging

Fisher Chemical outer packaging is UN compliant and 100% recyclable. It's an easy and safe way to reduce your waste management costs and contribute to a cleaner environment.



Aluminium standard metal solution 1000 ppm

for ICP solution in HNO₃ 2%, NIST certified

packaging	Mfr. No
100 mL Plastic bottle	SPPLAL22Y
250 mL Plastic bottle	SPPLAL22T
500 mL Plastic bottle	SPPLAL22X
ADR/RID: 8, II	H315, H319
EINECS: 236-751-8	P280, P305+P351+P338,
UN: 2031	P301+P330+P331, P310
GHS: Skin Corr. 1B, Eye Dam. 1	



Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Aluminium standard metal solution 1000 ppm

for ICP solution in HCL 2%, NIST certified

packaging	Mfr. No
500 mL Plastic bottle	SPPLAL12X

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Aluminium standard metal solution 10 000 ppm

for ICP solution in HNO₃ 5%, NIST certified

packaging	Mfr. No
100 mL Plastic bottle	SPPLAL23Y
500 mL Plastic bottle	SPPLAL23X
H314	



Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Aluminium standard metal solution 10 000 ppm

for ICP solution in HCL 5%, NIST certified

packaging	Mfr. No
500 mL Plastic bottle	SPPLAL13X

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Aluminium standard metal solution 1000 ppm

for ICP-MS solution in HNO₃ 2%, NIST certified

packaging	Mfr. No
100 mL Plastic bottle	SPCLAL22Y
H315, H319	



Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Aluminium foil

extra pure, approx. 0.1 mm thick

SLR

packaging	Mfr. No
250 g Plastic jar	A/1570/50
500 g Plastic jar	A/1570/53
Al	Boiling point: 2327°C
CAS: 7429-90-5	Melting point: 660°C
MW: 26.97	ADR/RID: 4.3/II
Density at 20°C: 2.7	EINECS: 231-072-3

Guaranteed analysis

Identity test : Pass

Aluminium thin sheet

extra pure, approx. 0.5 mm thick

SLR

packaging	Mfr. No
250 g Plastic jar	A/1690/50

Al	Boiling point: 2327°C
CAS: 7429-90-5	Melting point: 660°C
MW: 26.97	ADR/RID: 4.3/II
Density at 20°C: 2.7	EINECS: 231-072-3

Guaranteed analysis

Identity test : Pass

Aluminium powder

extra-pure, powder fine

SLR

packaging	Mfr. No
250 g Glass jar	A/1605/50
500 g Glass jar	A/1605/53

Al	EINECS: 231-072-3
CAS: 7429-90-5	UN: 1396
MW: 26.97	GHS: Flam. Sol. 1, Water-react. 2
Density at 20°C: 27.02	H228, H261
Boiling point: 2327°C	P402+P404, P280,
Melting point: 660°C	P231+P232, P210, P240
ADR/RID: 4.3, II	



Guaranteed analysis

Identity test : Pass

Aluminium powder

extra-pure, coarse (grease free)

SLR

packaging	Mfr. No
250 g Glass jar	A/1610/50
500 g Glass jar	A/1610/53

Al	ADR/RID: 4.3, II
CAS: 7429-90-5	EINECS: 231-072-3
MW: 26.97	UN: 1396
Density at 20°C: 27.02	GHS: Pyr. Sol. 1, Water-react. 2
Boiling point: 2327°C	H261, H250
Melting point: 660°C	P210, P370+P378, P402



Guaranteed analysis

Identity test : Pass

Aluminium metal flake

technical

SLR

packaging	Mfr. No
250 g Plastic jar	A/1600/50

Al	CAS: 7429-90-6
MW: 26.97	Boiling point: 2327°C

Guaranteed analysis

Identity test : Pass

Aluminium metal wire

for analysis, 1.0 mm diameter - 25 mtr

Certified AR

packaging	Mfr. No
COIL Plastic jar	A/1810/99
Al	Boiling point: 2327°C
CAS: 7429-90-5	Melting point: 660°C
MW: 26.97	EINECS: 231-072-3
Density at 20°C: 2.7	

Assay: >99.9%

Guaranteed analysis

Copper (Cu)	<50 ppm	Total nitrogen (N)	<10 ppm
Iron (Fe)	<40 ppm	Total silicon (Si)	<50 ppm
Manganese (Mn)	<20 ppm	Insoluble Acids	<0.005%

Lot analysis on the label

Aluminium ammonium sulfate dodecahydrate

for analysis

Certified AR

packaging	Mfr. No
500 g Plastic jar	A/1880/53
NH ₄ Al(SO ₄) ₂ ·12H ₂ O	EINECS: 232-055-3
CAS: 7784-26-1	GHS: Skin Irrit. 2, Eye Irrit. 2,
MW: 453.33	STOT SE 3
Density at 20°C: 1.64	H319, H315, H335
Boiling point: 200°C	P280, P305+P351+P338
Melting point: 93.5°C	



Assay: 99 to 100.5%

Guaranteed analysis

pH (aqueous solution 5%)	>3	Nitrate (NO ₃)	<0.001%
Calcium (Ca)	<20 ppm	Potassium (K)	<200 ppm
Copper (Cu)	<5 ppm	Sodium (Na)	<50 ppm
Iron (Fe)	<5 ppm	n-Butyl formate	<0.002%
Lead (Pb)	<5 ppm	Total silicon (Si)	<20 ppm
Magnesium (Mg)	<10 ppm	Zinc (Zn)	<10 ppm

Lot analysis on the label.

Aluminium ammonium sulfate dodecahydrate

specified reagents for general laboratory work- extra pure

SLR

packaging	Mfr. No
1 kg Plastic jar	A/1840/60
NH ₄ Al(SO ₄) ₂ ·12H ₂ O	EINECS: 232-055-3
CAS: 7784-26-1	GHS: Skin Irrit. 2, Eye Irrit. 2,
MW: 453.33	STOT SE 3
Density at 20°C: 1.64	H319, H315, H335
Boiling point: 200°C	P280, P305+P351+P338
Melting point: 93.5°C	



Assay: >99%

Guaranteed analysis

Calcium (Ca)	<50 ppm	Potassium (K)	<2000 ppm
Copper (Cu)	<20 ppm	Sodium (Na)	<100 ppm
Iron (Fe)	<20 ppm	n-Butyl formate	<0.005%
Lead (Pb)	<20 ppm	Zinc (Zn)	<50 ppm
Magnesium (Mg)	<50 ppm		

Aluminium chloride anhydrous

pure

packaging	Mfr. No
500 g Glass bottle	A/2000/53
AlCl ₃	EINECS: 231-208-1
CAS: 7446-70-0	UN: 1726
MW: 133.35	GHS: Skin Corr. 1B, Eye Dam. 1
Density at 20°C: 2.44	H314
Melting point: 194°C	P280, P305+P351+P338
ADR/RID: 8, II	



Assay: >98%

Guaranteed analysis

Iron (Fe)	<500 ppm	Sulfate (SO ₄)	<0.01%
-----------	----------	----------------------------	--------

Aluminium chloride hexahydrate

SLR, Eur.Ph., BP, USP

specified reagents for general laboratory work- extra pure

packaging	Mfr. No
500 g Plastic jar	A/1960/53
AlCl ₃ ·6H ₂ O	UN: 1726
CAS: 7784-13-6	GHS: Skin Irrit. 2, Eye Irrit. 2
MW: 241.44	H315, H319
Density at 20°C: 2.39	P302+P352, P280,
ADR/RID: 8/II	P305+P351+P338
EINECS: 231-208-1	



Assay: >99%

Guaranteed analysis

Calcium (Ca)	<100 ppm	Magnesium (Mg)	<100 ppm
Copper (Cu)	<50 ppm	Potassium (K)	<100 ppm
Iron (Fe)	<50 ppm	Sodium (Na)	<100 ppm
Lead (Pb)	<50 ppm	Zinc (Zn)	<100 ppm

Aluminium hydroxide dried

specified reagents for general laboratory work- extra pure

SLR

packaging	Mfr. No
1 kg Plastic jar	A/2080/60
Al(OH) ₃	EINECS: 244-492-7
CAS: 21645-51-2	GHS: Eye Irrit. 2
MW: 78.00	H319
Density at 20°C: 2.42	P280, P305+P351+P338
Melting point: 300°C	



Guaranteed analysis

Iron (Fe)	<100 ppm	n-Butyl formate	<0.01%
Sodium (Na)	<3000 ppm		

Aluminium nitrate nonahydrate

for analysis

Certified AR

packaging	Mfr. No
500 g Plastic jar	A/2202/53
1 kg Plastic jar	A/2202/60
Al(NO ₃) ₃ ·9H ₂ O	UN: 1438
CAS: 7784-27-2	GHS: Ox. Sol. 2, Skin Irrit. 2,
MW: 375.13	Eye Irrit. 2
Melting point: 73°C	H272, H315, H319
ADR/RID: 5.1, III	P302+P352, P280,
EINECS: 236-751-8	P305+P351+P338, P210



Assay: 99 to 101%

Guaranteed analysis

pH (aqueous solution 5%)	>2	Potassium (K)	<10 ppm
Ammonia (NH ₃)	<0.02%	Sodium (Na)	<50 ppm
Calcium (Ca)	<10 ppm	Total chloride (Cl)	<0.001%
Copper (Cu)	<5 ppm	Total phosphorus (P)	<20 ppm
Iron (Fe)	<20 ppm	Total silicon (Si)	<5 ppm
Lead (Pb)	<10 ppm	Total sulphur (S)	<50 ppm
Magnesium (Mg)	<5 ppm	Zinc (Zn)	<5 ppm

Typical analysis

Silver (Ag)	<2 ppm	Lithium (Li)	<2 ppm
Arsenic (As)	<1 ppm	Molybdenum (Mo)	<2 ppm
Barium (Ba)	<2 ppm	Nickel (Ni)	<2 ppm
Cadmium (Cd)	<2 ppm	Strontium (Sr)	<2 ppm
Chromium (Cr)	<2 ppm	Titanium (Ti)	<2 ppm
Cobalt (Co)	<2 ppm	Vanadium (V)	<2 ppm

Lot analysis on the label.

High Volume Solvent Delivery Systems

High-volume solvent delivery systems are available in 10L to 1000L, in Europe only. Please email your contact details and a brief description of your requirements to begei.sha.quotes@thermofisher.com



Aluminium nitrate nonahydrate		SLR
specified reagents for general laboratory work- extra pure		
packaging		Mfr. No
250 g	Plastic jar	A/2200/50
500 g	Plastic jar	A/2200/53
1 kg	Plastic jar	A/2200/60
5 kg	Plastic jar	A/2200/65
Al(NO ₃) ₃ ·9H ₂ O	UN: 1438	
CAS: 7784-27-2	GHS: Ox. Sol. 2, Skin Irrit. 2,	
MW: 375.13	Eye Irrit. 2	
Melting point: 73°C	H272, H315, H319	
ADR/RID: 5.1, III	P302+P352, P280,	
EINECS: 236-751-8	P305+P351+P338, P210	
Assay: >98%		

Guaranteed analysis		
Calcium (Ca)	<50 ppm	Sodium (Na)
Copper (Cu)	<20 ppm	n-Butyl formate
Iron (Fe)	<50 ppm	Total phosphorus (P)
Lead (Pb)	<20 ppm	Total silicon (Si)
Magnesium (Mg)	<50 ppm	Total sulphur (S)
Potassium (K)	<100 ppm	Zinc (Zn)

Aluminium oxide, acid Brockmann type I
For chromatography, "Camag" 100-250 mesh

packaging		Mfr. No
500 g	Plastic jar	A/2370/53
Al ₂ O ₃		Boiling point: 2980°C
CAS: 1344-28-1		Melting point: 2000°C
MW: 101.96		EINECS: 215-691-6
Guaranteed analysis		
pH (10% aq. suspension)		: 4 to 5
Acid, pH (4.0- 5.0)		

Aluminium oxide, basic Brockmann type I
For chromatography, "Camag" 100-250 mesh

packaging		Mfr. No
500 g	Plastic jar	A/2372/53
Al ₂ O ₃	Boiling point: 2980°C	
CAS: 1344-28-1	Melting point: 2000°C	
MW: 101.96	EINECS: 215-691-6	
Guaranteed analysis		
pH (10% aq. suspension)		: 9 to 10
Alkaline, pH (9.0- 10.0)		

Aluminium oxide, neutral Brockmann type I
For chromatography, "Camag" 100-250 mesh

packaging		Mfr. No
500 g	Plastic jar	A/2374/53
Al ₂ O ₃	Boiling point: 2980°C	
CAS: 1344-28-1	Melting point: 2000°C	
MW: 101.96	EINECS: 215-691-6	
Guaranteed analysis		
pH (10% aq. suspension)		: 6.5 to 7.5
Neutral, pH (6.5-7.5)		

Aluminium oxide, basic Brockmann type II
For chromatography, "Camag" 100-250 mesh

packaging		Mfr. No
500 g	Plastic jar	A/2380/53
Al ₂ O ₃	Boiling point: 2980°C	
CAS: 1344-28-1	Melting point: 2000°C	
MW: 101.96	EINECS: 215-691-6	
Guaranteed analysis		
pH (10% aq. suspension)		: 9 to 10
Alkaline, pH (9-10)		

Aluminium oxide			Certified AR
for analysis, calcined			
Packaging			Mfr. No
1 kg Plastic jar			A/2360/60
Al ₂ O ₃		Boiling point: 2980°C	
CAS: 1344-28-1		Melting point: 2000°C	
MW: 101.96		EINECS: 215-691-6	
Guaranteed analysis			
Loss on ignition	<0.5%	Sulfate (SO ₄)	<0.005%
Iron (Fe)	<50 ppm	n-Butyl formate	<0.005%
Lot analysis on the label.			

Aluminium oxide
pure, Desiccating agent

Packaging		Mfr. No
500 g	Plastic jar	A/2320/53
Al ₂ O ₃		Boiling point: 2980°C
CAS: 1344-28-1		Melting point: 2000°C
MW: 101.96		EINECS: 215-691-6
Guaranteed analysis		
Identity test		: Pass
Particle size 3-6 mm.		

Aluminium oxide calcined
specified reagents for general laboratory work- extra pure

Packaging		Mfr. No
500 g	Plastic jar	A/2280/53
1 kg	Plastic jar	A/2280/60
Al ₂ O ₃		Boiling point: 2980°C
CAS: 1344-28-1		Melting point: 2000°C
MW: 101.96		EINECS: 215-691-6
Guaranteed analysis		
Identity test		: Pass

Aluminium potassium sulfate dodecahydrate
specified reagents for general laboratory work- extra pure

Packaging		Mfr. No
500 g	Plastic jar	A/2400/53
1 kg	Plastic jar	A/2400/60
3 kg	Plastic jar	A/2400/63
AlK(SO ₄) ₂ ·12H ₂ O		
CAS: 7784-24-9		
MW: 474.39		
Density at 20°C: 1.757		
Melting point: 92°C		
Assay: >99%		
Guaranteed analysis		
Heavy metals	<20 ppm	Sodium (Na)
Ammonium (NH ₄)	<0.2%	n-Butyl formate
Calcium (Ca)	<100 ppm	Total nitrogen (N)
Copper (Cu)	<20 ppm	Total phosphorus (P)
Iron (Fe)	<50 ppm	Total silicon (Si)
Lead (Pb)	<20 ppm	Zinc (Zn)
Magnesium (Mg)	<100 ppm	

Aluminium silicate (hydrated)

See Kaolin - light, page 169

Aluminium sulfate hexadecahydrate		Certified AR	
for analysis			
packaging		Mfr. No	
500 g	Plastic jar	A/2600/53	
Al ₂ (SO ₄) ₃ ·16H ₂ O	1, STOT SE 3, Aquatic Chronic 3		
CAS: 16828-11-8	H412, H335, H318, H315		
MW: 630.39	P280, P273,		
EINECS: 233-135-0	P305+P351+P338,		
GHS: Skin Irrit. 2, Eye Dam.	P302+P352, P261		
Assay: 98 to 102%			
Guaranteed analysis			
pH (aqueous solution 5%)	>2.2	Sodium (Na)	<500 ppm
Calcium (Ca)	<50 ppm	Total chloride (Cl)	<0.001%
Copper (Cu)	<2 ppm	Total phosphorus (P)	<50 ppm
Iron (Fe)	<50 ppm	Total silicon (Si)	<20 ppm
Lead (Pb)	<10 ppm	Zinc (Zn)	<10 ppm
Magnesium (Mg)	<10 ppm	Total nitrogen (N)	<200 ppm
Potassium (K)	<20 ppm		
Typical analysis			
Silver (Ag)	<1 ppm	Manganese (Mn)	<1 ppm
Arsenic (As)	<5 ppm	Molybdenum (Mo)	<1 ppm
Barium (Ba)	<1 ppm	Nickel (Ni)	<1 ppm
Cadmium (Cd)	<1 ppm	Strontium (Sr)	<1 ppm
Chromium (Cr)	<1 ppm	Titanium (Ti)	<1 ppm
Cobalt (Co)	<1 ppm	Vanadium (V)	<1 ppm
Lithium (Li)	<1 ppm		
Lot analysis on the label.			

Aluminium sulfate hexadecahydrate
specified reagents for general laboratory work- extra pure

packaging		Mfr. No	
500 g	Plastic jar	A/2560/53	
1 kg	Plastic jar	A/2560/60	
5 kg	Plastic keg	A/2560/65	
Al ₂ (SO ₄) ₃ ·16H ₂ O			
CAS: 16828-11-8			
MW: 630.39			
EINECS: 233-135-0			
GHS: Skin Irrit. 2, Eye Dam.			
Assay: >94%			
Guaranteed analysis			
Calcium (Ca)	<500 ppm	Sodium (Na)	<1000 ppm
Copper (Cu)	<20 ppm	n-Butyl formate	<0.01%
Iron (Fe)	<100 ppm	Total nitrogen (N)	<1000 ppm
Lead (Pb)	<100 ppm	Total phosphorus (P)	<500 ppm
Magnesium (Mg)	<100 ppm	Total silicon (Si)	<200 ppm
Potassium (K)	<200 ppm	Zinc (Zn)	<100 ppm

Aluminium sulfate hexadecahydrate
Technical

packaging	Mfr. No	
2 kg Plastic jar	A/2520/61	
Al ₂ (SO ₄) ₃ ·16H ₂ O	1, STOT SE 3, Aquatic Chronic 3	
CAS: 16828-11-8	H412, H335, H318, H315	
MW: 630.39	P280, P273,	
EINECS: 233-135-0	P305+P351+P338,	
GHS: Skin Irrit. 2, Eye Dam.	P302+P352, P261	

Amberjet® 1000H cation
for Ion chromatography, Gel-strongly acidic, Hydrogen form (size 0.62-0.80mm),

packaging		Mfr. No
500 g	Plastic jar	IE/1012/53
CAS: 69011-20-7	H318	
Boiling point: 100°C	P280, P305+P351+P338, P310	
GHS: Eye Dam. 1		
Guaranteed analysis		
Identity test		: Pass

Amberlite® IR-120, Cation exchange resin
for Ion chromatography, Gel type (Hydrogen form , strongly acidic)

Packaging		Mfr. No
100 g	Plastic jar	IE/1015/48
500 g	Plastic jar	
Guaranteed analysis		
Identity test		: Pass
Styrene-DVB Matrix structure.		
<i>Premium quality conventional strongly acidic cation exchange resin with 8% DVB by weight, amber in colour and contains a high perfect bead content. The quality of this resin is such that it can be used in standard operations as well as where there is a need for a premium grade gel cation exchange resin.</i>		

Amberlite® IR-120, Cation exchange resin
for Ion chromatography, Gel type (Sodium form , strongly acidic)

Packaging	Mfr. No
500 g Plastic jar	IE/1020/53
Guaranteed analysis	
Identity test	: Pass
<i>Styrene-DVB Matrix structure.</i>	
<i>Premium quality conventional strongly acidic cation exchange resin with 8% DVB by weight, amber in colour and contains a high perfect bead content. The quality of this resin is such that it can be used in standard operations as well as where there is a need for a premium grade gel cation exchange resin.</i>	

Amberlite® 200, Cation exchange resin
for Ion chromatography, Macroreticular type (Sodium form, strongly acidic)

Packaging	Mfr. No
600 g Plastic jar	IE/1040/53
Guaranteed analysis	
Identity test	: Pass
<i>Styrene-DVB matrix structure.</i>	
<i>Excellent physical stability and resistance to oxidation. Contains approximately 20%DVB by weight, and can be expected to provide at least three times greater resistance to oxidation than conventional gel type cation exchange resin.</i>	

Amberlite® IRA-67 chloride
for Ion chromatography, Gel type, free base: weakly basic

packaging		Mfr. No
500 g	Plastic jar	IE/1080/53
Guaranteed analysis		
Identity test		: Pass

Amberlite® IRA-402, Anion exchange resin
for Ion chromatography, Gel type (Chloride form, strongly basic)

Packaging	Mfr. No
500 g Plastic jar	IE/1050/53
Guaranteed analysis	
Identity test	: Pass
<i>Styrene-DVB Matrix structure.</i>	
<i>Combined with a strong acid cation exchanger it reduces both strong and weak acid concentrations to extremely low levels.Its main use is water demineralisation.</i>	

Amberlite® MB-1 (Mixed resin)
for Ion chromatography

Packaging		Mfr. No
600 g	Plastic jar	IE/1505/53
Guaranteed analysis		
Identity test		: Pass
<i>Can be used to remove trace ionic impurities from acrylamide before polymerisation.</i>		

Amberlite® MB-6113 (Mixed resin)
for Ion chromatography

packaging	Mfr. No
500 g Plastic jar	IE/1502/53
GHS: Eye Dam. 1 H318 P280, P305+P351+P338, P310	
Guaranteed analysis	
Identity test : Pass	

Amberlite® FPC3500 cation
for Ion chromatography, Carboxylic functionality, Hydrogen form (size 0.45-0.65mm),

packaging	Mfr. No
100 g Plastic jar	IE/1046/48
CAS: 950182-46-4	
Guaranteed analysis	
Identity test : Pass	

Amberlyst® 15, Amberlyst exchange resin
for Ion chromatography, Macroreticular type (cation, strongly acidic), dry

packaging	Mfr. No
500 g Plastic jar	IE/2500/53
GHS: Skin Irrit. 2, Eye Irrit. 2, STOT SE 3 H319, H315, H335 P261, P280, P301+P312, P302+P352, P305+P351+P338	

Guaranteed analysis	
Identity test : Pass	
Amberlyst 15 Matrix structure	

Amidosulfonic acid See Sulfamic acid, page 334

4-Aminophenazone Certified AR
for analysis

packaging	Mfr. No
25 g Plastic jar	A/3042/46
C ₁₁ H ₁₃ N ₃ O CAS: 83-07-8 MW: 203.25 Density at 20°C: 0.800 Melting point: 105.5 to 110°C EINECS: 201-452-3	
GHS: Acute Tox. 4, Skin Irrit. 2, Eye Irrit. 2, STOT SE 3 H302, H315, H335, H319 P261, P301+P312, P302+P352, P280, P305+P351+P338	
Assay: >99%	
Guaranteed analysis	
Loss on drying <0.5% Solution : Pass	
Residue after ignition <0.1% Complies with standard IR spectrum: Pass	
Lot analysis on the label.	

4-Aminophenazone SLR
specified reagents for general laboratory work- extra pure

packaging	Mfr. No
25 g Plastic jar	A/3040/46
C ₁₁ H ₁₃ N ₃ O CAS: 83-07-8 MW: 203.25 Density at 20°C: 0.800 Melting point: 105.5 to 110°C EINECS: 201-452-3	
GHS: Acute Tox. 4, Skin Irrit. 2, Eye Irrit. 2, STOT SE 3 H302, H315, H335, H319 P261, P301+P312, P302+P352, P280, P305+P351+P338	
Assay: >98%	

Aminobenzene See Aniline, page 48

4-Aminobenzenesulfonamide

See Sulfanilamide, page 335

4-Aminobenzenesulfonic acid

See Sulfanilic acid 99+ % , page 335

4-Aminobenzoic acid SLR
specified reagents for general laboratory work- extra pure

packaging	Mfr. No
100 g Plastic jar	A/2840/48
C ₇ H ₇ NO ₂ CAS: 150-13-0 MW: 137.14 Melting point: 186 to 189°C EINECS: 205-753-0	
GHS: Skin Irrit. 2, Eye Irrit. 2, STOT SE 3 H315, H335, H319 P261, P302+P352, P280, P305+P351+P338, P264	
Assay: >99%	

3-Amino-7-dimethylamino-2-methylphenazine hydrochloride

See Neutral red C.I 50040, page 212

4-Amino-2.3-dimethyl-1-phenyl-3-pyrazolin-5-one

See 4-Aminophenazone, page 32

2-Aminoethanesulfonic acid

See Taurine 99+ % (2-aminoethane-sulfonic acid), page 342

2-Aminoethanol See Ethanolamine, page 123

3-Aminoethylbenzene

See 2-Phenylethylamine 98.5+ % , page 239

2-Amino-2-(hydroxymethyl)-1.3-propanediol

See Tris(hydroxymethyl) methylamine, page 362

L-2-Amino-3-(4-imidazolyl)propionic acid

See L-Histidine, page 144

3-Aminomethylalizarin-N,N-diacetic acid
pure, indicator grade

packaging	Mfr. No
1 g Glass jar	A/2930/43
C ₁₆ H ₁₅ NO ₈ ·2H ₂ O CAS: 3952-78-1 MW: 421.35	
Melting point: 190°C EINECS: 223-544-2	
Guaranteed analysis	
Complies with standard IR spectrum: Pass	
Conforming to the IR spectrum of reference.	



(2S,3S)-2-Amino-3-methylpentanoic acid

See L-iso-Leucine 99 % , page 182

2-Amino-2-methyl-propan-1-ol SLR
specified reagents for general laboratory work- extra pure

packaging	Mfr. No
500 mL Glass bottle	A/2960/08
C ₄ H ₁₁ NO CAS: 124-68-5 MW: 89.14 Density at 20°C: 0.930 Boiling point: 165°C Melting point: 31 to 32°C Flash point: 67°C ADR/RID: 8, III	
EINECS: 204-709-8 UN: 3259 GHS: Skin Irrit. 2, Eye Irrit. 2, Aquatic Chronic 3 H412, H315, H319 P273, P302+P352, P280, P305+P351+P338	
Assay (CG): >93%	
Guaranteed analysis	
Residue after ignition <0.005%	

3-Aminophthalhydrazide (Luminol) SLR
specified reagents for general laboratory work- extra pure

packaging	Mfr. No
1 g Vial	A/3150/43
5 g Glass jar	A/3150/44
C ₈ H ₇ N ₃ O ₂ CAS: 521-31-3 MW: 177.16 Melting point: 319 to 320°C EINECS: 208-309-4 GHS: Acute Tox. 4, Skin Irrit.	
2, Eye Irrit. 2, STOT SE 3 H335, H302, H315, H319 P261, P301+P312, P302+P352, P280, P305+P351+P338	
Guaranteed analysis	
Complies with standard IR spectrum: Pass	

L-2-Aminosuccinamic acid

See L-Asparagine, page 465

L-Aminosuccinic acid See L-Aspartic acid, page 51

Ammonia solution 20-22% NEW Optima
for ultra trace elemental analysis

packaging	Mfr. No
250 mL Pre-cleaned HDPE bottle	A470-250
500 mL Pre-cleaned HDPE bottle	A470-500
1 L Pre-cleaned HDPE bottle	A470-1
NH ₃ ·H ₂ O CAS: 1336-21-6 MW: 35.05 Density at 20°C: 0.920 ADR/RID: 8, III EINECS: 215-647-6	
UN: 2672 GHS: Skin Corr. 1B, STOT SE 3 H314, H335 P301+P330+P331, P305+P351+P338, P280, P261	
Assay: 20-22%	

Guaranteed analysis	
Zirconium (Zr) <10 ppt	Cadmium (Cd) <10 ppt
Gallium (Ga) <10 ppt	Calcium (Ca) <10 ppt
Terbium (Tb) <10 ppt	Cerium (Ce) <10 ppt
Indium (In) <10 ppt	Chromium (Cr) <10 ppt
Dysprosium (Dy) <10 ppt	Cobalt (Co) <10 ppt
Aluminium (Al) <20 ppt	Copper (Cu) <10 ppt
Antimony (Sb) <10 ppt	Erbium (Er) <10 ppt
Arsenic (As) <10 ppt	Europium (Eu) <10 ppt
Barium (Ba) <10 ppt	Gadolinium (Gd) <10 ppt
Beryllium (Be) <10 ppt	Germanium (Ge) <10 ppt
Bismuth (Bi) <10 ppt	

Gold (Au) <10 ppt	Samarium (Sm) <10 ppt
Holmium (Ho) <10 ppt	Scandium (Sc) <10 ppt
Iron (Fe) <10 ppt	Propan--1-ol <10 ppt
Lanthanum (La) <10 ppt	Sodium (Na) <20 ppt
Lead (Pb) <10 ppt	Strontium (Sr) <10 ppt
Lithium (Li) <10 ppt	Tellurium (Te) <10 ppt
Lutetium (Lu) <10 ppt	Thallium (Tl) <10 ppt
Magnesium (Mg) <10 ppt	Thorium (Th) <10 ppt
Manganese (Mn) <10 ppt	Thulium (Tm) <10 ppt
Mercury (Hg) <200 ppt	Tin (Sn) <10 ppt
Molybdenum (Mo) <10 ppt	Titanium (Ti) <10 ppt
Neodymium (Nd) <10 ppt	Uranium (U) <10 ppt
Nickel (Ni) <10 ppt	Vanadium (V) <10 ppt
Niobium (Nb) <10 ppt	Ytterbium (Yb) <10 ppt
Potassium (K) <10 ppt	Yttrium (Y) <10 ppt
Praseodymium (Pr) <10 ppt	Zinc (Zn) <10 ppt
Rhodium (Rh) <10 ppt	Rubidium (Rb) <10 ppt

Bottled in class 10 cleanroom conditions. Individually double-bagged in class 100 cleanroom conditions. Each bottle is individually packaged. Certificate of analysis delivered with each bottle.

Ammonia solution 20-22% NEW Trace Metal
for trace metal analysis

packaging	Mfr. No
500 mL HDPE plastic bottle	A512-P500
NH ₃ ·H ₂ O CAS: 1336-21-6 MW: 35.05 Density at 20°C: 0.921 ADR/RID: 8, III EINECS: 215-647-6	
UN: 2672 GHS: Skin Corr. 1B, STOT SE 3 H314, H335 P301+P330+P331, P305+P351+P338, P280, P261	

Assay: 20-22%	
Guaranteed analysis	
Colour <10 APHA	Nickel (Ni) <0.5 ppb
Aluminium (Al) <1 ppb	Niobium (Nb) <0.1 ppb
Antimony (Sb) <0.5 ppb	Palladium (Pd) <1 ppb
Arsenic (As) <1 ppb	Platinum (Pt) <1 ppb
Barium (Ba) <0.1 ppb	Potassium (K) <1 ppb
Beryllium (Be) <0.1 ppb	Praseodymium (Pr) <0.1 ppb
Bismuth (Bi) <0.1 ppb	Rhenium (Re) <1 ppb
Cadmium (Cd) <0.5 ppb	Rhodium (Rh) <0.5 ppb
Calcium (Ca) <1 ppb	Rubidium (Rb) <0.1 ppb
Cerium (Ce) <0.1 ppb	Ruthenium (Ru) <1 ppb
Cesium (Cs) <0.1 ppb	Samarium (Sm) <0.1 ppb
Chromium (Cr) <0.5 ppb	Scandium (Sc) <0.1 ppb
Cobalt (Co) <0.5 ppb	Selenium (Se) <1 ppb
Copper (Cu) <0.5 ppb	Silver (Ag) <0.5 ppb
Dysprosium (Dy) <0.1 ppb	Sodium (Na) <1 ppb
Erbium (Er) <0.1 ppb	Strontium (Sr) <0.1 ppb
Europium (Eu) <0.1 ppb	Sulfate (SO ₄) <1 ppm
Gadolinium (Gd) <0.1 ppb	Tellurium (Te) <0.1 ppb
Gallium (Ga) <0.1 ppb	Terbium (Tb) <0.1 ppb
Germanium (Ge) <0.1 ppb	Thallium (Tl) <0.1 ppb
Gold (Au) <0.5 ppb	Thorium (Th) <0.1 ppb
Hafnium (Hf) <1 ppb	Thulium (Tm) <0.1 ppb
Holmium (Ho) <0.1 ppb	Tin (Sn) <0.5 ppb
Indium (In) <0.1 ppb	Titanium (Ti) <0.5 ppb
Iron (Fe) <1 ppb	Total chloride (Cl) <0.5 ppm
Lanthanum (La) <0.1 ppb	Total phosphorus (P) <0.5 ppm
Lead (Pb) <0.1 ppb	Tungsten (W) <0.1 ppb
Lithium (Li) <0.1 ppb	Uranium (U) <0.1 ppb
Lutetium (Lu) <0.1 ppb	Vanadium (V) <0.5 ppb
Magnesium (Mg) <1 ppb	Ytterbium (Yb) <0.1 ppb
Manganese (Mn) <0.5 ppb	Yttrium (Y) <0.1 ppb
Mercury (Hg) <0.2 ppb	Zinc (Zn) <0.5 ppb
Molybdenum (Mo) <0.5 ppb	Zirconium (Zr) <0.1 ppb
Neodymium (Nd) <0.1 ppb	

Bottled in class 10 cleanroom conditions. Lot analysis on the label.



Ammonia solution 20% d=0.92		Certified AR
for analysis		
packaging	Mfr. No	
1 ℓ Glass bottle	A/3367/15	
2.5 ℓ Glass bottle	A/3367/17	
NH ₃ .H ₂ O	UN: 2672	
CAS: 1336-21-6	GHS: Skin Corr. 1B, STOT SE 3	
MW: 35.05	H314, H335	
Density at 20°C: 0.920	P301+P330+P331,	
ADR/RID: 8, III	P305+P351+P338, P280, P261	
EINECS: 215-647-6		

Assay: 20 to 22%	
Guaranteed analysis	
Colour	<10 APHA
Density at 20°C	: 0.917 to 0.923
Substances	
reducing KMnO ₄	<0.001%
Calcium (Ca)	<1 ppm
Carbonate (Na ₂ CO ₃)	<0.002%
Copper (Cu)	<0.05 ppm
Iron (Fe)	<0.2 ppm
Lead (Pb)	<0.05 ppm
Magnesium (Mg)	<0.1 ppm
Potassium (K)	<1 ppm
Sodium (Na)	<2 ppm
Total chloride (Cl)	<0.00005%
Total phosphorus (P)	<0.5 ppm
Total silicon (Si)	<1 ppm
Total sulphur (S)	<1 ppm
Zinc (Zn)	<0.2 ppm
Pyridine & homologues	<2 ppm

Lot analysis on the label.

Ammonia Solution 20% d=0.92

specified reagents for general laboratory work- extra pure

packaging	Mfr. No
1 ℓ Glass bottle	A/3366/15
5 ℓ HDPE Plastic bottle	A/3366/21
NH ₃ .H ₂ O	UN: 2672
CAS: 1336-21-6	GHS: Skin Corr. 1B, STOT SE 3
MW: 35.05	H314, H335
Density at 20°C: 0.920	P301+P330+P331,
ADR/RID: 8, III	P305+P351+P338, P280, P261
EINECS: 215-647-6	

Assay: 20 to 22%	
Guaranteed analysis	
Density at 20°C	: 0.91 to 0.93
Calcium (Ca)	<10 ppm
Carbonate (Na ₂ CO ₃)	<0.01%
Copper (Cu)	<1 ppm
Iron (Fe)	<1 ppm
Lead (Pb)	<2 ppm
Magnesium (Mg)	<1 ppm
Potassium (K)	<10 ppm
Sodium (Na)	<20 ppm
n-Butyl formate	<0.001%
Total phosphorus (P)	<20 ppm
Total sulphur (S)	<20 ppm
Zinc (Zn)	<2 ppm

Ammonia solution 25% d=0.91

for analysis

packaging	Mfr. No
1 ℓ Coated Glass bottle	A/3360/PB15
2.5 ℓ Coated Glass bottle	A/3360/PB17
NH ₃ .H ₂ O	UN: 2672
CAS: 1336-21-6	GHS: Skin Corr. 1B, Eye
MW: 35.05	Dam. 1, STOT SE 3, Aquatic
Density at 20°C: 0.910	Acute 1
Boiling point: 38°C	H335, H410, H314
Melting point: -57°C	P301+P330+P331, P280,
ADR/RID: 8, III	P305+P351+P338, P310,
EINECS: 215-647-6	P261, P273

Assay: 20 to 30%	
Guaranteed analysis	
Colour	<10 APHA
Substances	
reducing KMnO ₄	<0.001%
Calcium (Ca)	<1 ppm
Carbonate (Na ₂ CO ₃)	<0.002%
Copper (Cu)	<0.05 ppm
Iron (Fe)	<0.2 ppm
Lead (Pb)	<0.05 ppm
Magnesium (Mg)	<0.1 ppm
Potassium (K)	<1 ppm
Sodium (Na)	<5 ppm
Total chloride (Cl)	<0.00005%
Total phosphorus (P)	<0.5 ppm
Total silicon (Si)	<1 ppm
Total sulphur (S)	<1 ppm
Zinc (Zn)	<0.2 ppm
Pyridine & homologues	<2 ppm

Typical analysis	
Aluminium (Al)	<0.1 ppm
Barium (Ba)	<0.01 ppm
Cadmium (Cd)	<0.01 ppm
Chromium (Cr)	<0.02 ppm
Cobalt (Co)	<0.01 ppm
Lithium (Li)	<0.01 ppm
Manganese (Mn)	<0.01 ppm
Molybdenum (Mo)	<0.01 ppm
Nickel (Ni)	<0.02 ppm
Strontium (Sr)	<0.01 ppm
Titanium (Ti)	<0.01 ppm
Vanadium (V)	<0.01 ppm

Lot analysis on the label.
Ammonia is equivalent to a strength of 13.4M (13.4N).

Ammonia solution 25% d=0.91

specified reagents for general laboratory work- extra pure

packaging	Mfr. No
1 ℓ Coated Glass bottle	A/3320/PB15
2.5 ℓ Coated Glass bottle	A/3320/PB17
5 ℓ HDPE Plastic bottle	A/3320/PB21
NH ₃ .H ₂ O	UN: 2672
CAS: 1336-21-6	GHS: Skin Corr. 1B, Eye
MW: 35.05	Dam. 1, STOT SE 3, Aquatic
Density at 20°C: 0.910	Acute 1
Boiling point: 38°C	H335, H410, H314
Melting point: -57°C	P301+P330+P331, P280,
ADR/RID: 8, III	P305+P351+P338, P310,
EINECS: 215-647-6	P261, P273

Assay: 20 to 30%	
Guaranteed analysis	
Colour	<20 APHA
Calcium (Ca)	<10 ppm
Carbonate (Na ₂ CO ₃)	<0.01%
Copper (Cu)	<1 ppm
Iron (Fe)	<1 ppm
Lead (Pb)	<1 ppm
Magnesium (Mg)	<2 ppm
Potassium (K)	<10 ppm
Sodium (Na)	<20 ppm
n-Butyl formate	<0.001%
Total phosphorus (P)	<20 ppm
Total sulphur (S)	<20 ppm
Zinc (Zn)	<2 ppm

Ammonia is equivalent to a strength of 13.4M (13.4N).

Ammonia solution 35% d=0.88

for HPLC

packaging	Mfr. No
100 mℓ Coated Glass bottle	A/3295/PB05
NH ₃ .H ₂ O	GHS: Skin Corr. 1B, Eye
CAS: 1336-21-6	Dam. 1, STOT SE 3, Aquatic
MW: 35.05	Acute 1
Density at 20°C: 0.88	H335, H314, H400
ADR/RID: 8, III	P301+P330+P331, P280,
EINECS: 215-647-6	P305+P351+P338, P310,
UN: 2672	P261, P273

Assay: 32 to 36%	
Guaranteed analysis	
Iron (Fe)	<0.2 ppm
Lead (Pb)	<0.05 ppm

Absorbance specification	
Max. Absorbance (1cm cell v HPLC water)	
254 nm	<0.01 A.U.

Ammonia is equivalent to a strength of 18.1M (18.1N).
Lot analysis on the label.
Filtered to 0.2 μm.

Optima LC/MS

Fisher Chemical OPTIMA® LC/MS solvents, blends, additives and modifiers have set the standard of excellence for consistent, reproducible performance in the mobile phase of LC/MS.

Find out more on page 308

Ammonia solution 35% d=0.88

for trace metal analysis

packaging	Mfr. No
500 mℓ HDPE Plastic bottle	A/3290/PB08
1 ℓ HDPE Plastic bottle	A/3290/PB15
NH ₃ .H ₂ O	GHS: Skin Corr. 1B, Eye
CAS: 1336-21-6	Dam. 1, STOT SE 3, Aquatic
MW: 35.05	Acute 1
Density at 20°C: 0.88	H335, H314, H400
ADR/RID: 8, III	P301+P330+P331, P280,
EINECS: 215-647-6	P305+P351+P338, P310,
UN: 2672	P261, P273

Assay: 32 to 36%	
Guaranteed analysis	
Other organic impurities	: <1 ppm
Oxygen absorbed	: <0.0002%
Aluminium (Al)	<0.05 ppm
Arsenic (As)	<0.005 ppm
Barium (Ba)	<0.01 ppm
Cadmium (Cd)	<0.001 ppm
Calcium (Ca)	<0.2 ppm
Carbonate (Na ₂ CO ₃)	<0.001%
Chromium (Cr)	<0.005 ppm
Cobalt (Co)	<0.005 ppm
Copper (Cu)	<0.005 ppm
Iron (Fe)	<0.05 ppm
Lead (Pb)	<0.002 ppm
Lithium (Li)	<0.005 ppm
Magnesium (Mg)	<0.05 ppm
Manganese (Mn)	<0.005 ppm
Mercury (Hg)	<0.01 ppm
Molybdenum (Mo)	<0.005 ppm
Nickel (Ni)	<0.005 ppm
Potassium (K)	<0.05 ppm
Propan--1-ol	<0.05 ppm
Strontium (Sr)	<0.005 ppm
Tin (Sn)	<0.005 ppm
Titanium (Ti)	<0.005 ppm
n-Butyl formate	<0.00001%
Total phosphorus (P)	<0.05 ppm
Total sulphur (S)	<0.1 ppm
Vanadium (V)	<0.005 ppm
Zinc (Zn)	<0.01 ppm

Ammonia is equivalent to a strength of 18.1M (18.1N).

Ammonia solution 35% d=0.88

for AAS

packaging	Mfr. No
2.5 ℓ Coated Glass bottle	A/3285/PB17
NH ₃ .H ₂ O	GHS: Skin Corr. 1B, Eye
CAS: 1336-21-6	Dam. 1, STOT SE 3, Aquatic
MW: 35.05	Acute 1
Density at 20°C: 0.88	H335, H314, H400
ADR/RID: 8, III	P301+P330+P331, P280,
EINECS: 215-647-6	P305+P351+P338, P310,
UN: 2672	P261, P273

Assay: 32 to 36%	
Guaranteed analysis	
Cadmium (Cd)	<0.02 ppm
Calcium (Ca)	<0.5 ppm
Cobalt (Co)	<0.02 ppm
Copper (Cu)	<0.05 ppm
Iron (Fe)	<0.2 ppm
Lead (Pb)	<0.05 ppm
Magnesium (Mg)	<0.02 ppm
Manganese (Mn)	<0.01 ppm
Nickel (Ni)	<0.01 ppm
Potassium (K)	<0.2 ppm
Sodium (Na)	<0.5 ppm
Zinc (Zn)	<0.1 ppm

Ammonia solution 35% d=0.88

for electronic use (MOS)

packaging	Mfr. No
2.5 ℓ Coated Glass bottle	A/3222/PB17
5 ℓ HDPE Plastic bottle	A/3222/21
NH ₃ .H ₂ O	GHS: Skin Corr. 1B, Eye
CAS: 1336-21-6	Dam. 1, STOT SE 3, Aquatic
MW: 35.05	Acute 1
Density at 20°C: 0.88	H335, H314, H400
ADR/RID: 8, III	P301+P330+P331, P280,
EINECS: 215-647-6	P305+P351+P338, P310,
UN: 2672	P261, P273

Assay: 32 to 36%	
Iron (Fe)	<0.2 ppm
Lead (Pb)	<0.05 ppm

Guaranteed analysis	
Colour	<10 APHA
Density at 20°C	: 0.88 to 0.89
Substances	
reducing KMnO ₄	<0.001%
Aluminium (Al)	<0.1 ppm
Barium (Ba)	<0.01 ppm
Cadmium (Cd)	<0.01 ppm
Calcium (Ca)	<1 ppm
Carbonate (Na ₂ CO ₃)	<0.002%
Chromium (Cr)	<0.02 ppm
Cobalt (Co)	<0.01 ppm
Copper (Cu)	<0.05 ppm
Iron (Fe)	<0.2 ppm
Lead (Pb)	<0.05 ppm
Lithium (Li)	<0.01 ppm
Magnesium (Mg)	<0.1 ppm
Manganese (Mn)	<0.01 ppm
Molybdenum (Mo)	<0.01 ppm
Nickel (Ni)	<0.02 ppm
Potassium (K)	<1 ppm
Sodium (Na)	<5 ppm
Strontium (Sr)	<0.01 ppm
Titanium (Ti)	<0.01 ppm
n-Butyl formate	<0.00005%
Total phosphorus (P)	<0.5 ppm
Total sulphur (S)	<1 ppm
Vanadium (V)	<0.01 ppm
Zinc (Zn)	<0.2 ppm
Pyridine & homologues	<2 ppm

Ammonia solution 0.25M

for HPLC

packaging	Mfr. No
500 mℓ Coated Glass bottle	J/0100/PB08
NH ₃ .H ₂ O	GHS: Skin Corr. 1B, Eye
CAS: 1336-21-6	Dam. 1, STOT SE 3, Aquatic
MW: 35.05	Acute 1, Aquatic Chronic 1
Density at 20°C: 0.98	H314, H335, H410
Melting point: -2.9°C	P301+P330+P331, P280,
ADR/RID: 8, III	P305+P351+P338, P310,
EINECS: 215-647-6	P261, P273
UN: 2672	

Absorbance specification	
Max. Absorbance (1cm cell v HPLC water)	
Drift in gradient baseline	<40 mAU.
Absorbance of the largest eluted peak: <1.5 mAU.	

Ammonia solution 35% d=0.88

for analysis

packaging	Mfr. No
1 ℓ Coated Glass bottle	A/3280/PB15
2.5 ℓ Coated Glass bottle	A/3280/PB17
5 ℓ HDPE Plastic bottle	A/3280/21
25 ℓ HDPE Plastic container	A/3280/25
NH ₃ .H ₂ O	GHS: Skin Corr. 1B, Eye
CAS: 1336-21-6	Dam. 1, STOT SE 3, Aquatic
MW: 35.05	Acute 1
Density at 20°C: 0.88	H335, H314, H400
ADR/RID: 8, III	P301+P330+P331, P280,
EINECS: 215-647-6	P305+P351+P338, P310,
UN: 2672	P261, P273

Assay: 32 to 36%	
Guaranteed analysis	
Colour	<10 APHA
Density at 20°C	: 0.88 to 0.89
Substances	
reducing KMnO ₄	<0.001%
Calcium (Ca)	<1 ppm
Carbonate (Na ₂ CO ₃)	<0.002%
Copper (Cu)	<0.05 ppm
Iron (Fe)	<0.2 ppm
Lead (Pb)	<0.05 ppm
Magnesium (Mg)	<0.1 ppm
Potassium (K)	<1 ppm
Sodium (Na)	<5 ppm
Total chloride (Cl)	<0.00005%
Total phosphorus (P)	<0.5 ppm
Total silicon (Si)	<1 ppm
Total sulphur (S)	<1 ppm
Zinc (Zn)	<0.2 ppm
Pyridine & homologues	<2 ppm

Typical analysis	
Aluminium (Al)	<0.1 ppm
Barium (Ba)	<0.01 ppm
Cadmium (Cd)	<0.01 ppm
Chromium (Cr)	<0.02 ppm
Cobalt (Co)	<0.01 ppm
Lithium (Li)	<0.01 ppm
Manganese (Mn)	<0.01 ppm
Molybdenum (Mo)	<0.01 ppm
Nickel (Ni)	<0.02 ppm
Strontium (Sr)	<0.01 ppm
Titanium (Ti)	<0.01 ppm
Vanadium (V)	<0.01 ppm

Lot analysis on the label.
Ammonia is equivalent to a strength of 18.1M (18.1N).

Ammonia solution 35% d=0.88
specified reagents for general laboratory work- extra pure

SLR

packaging	Mfr. No
1 ℓ Coated Glass bottle	A/3240/PB15
2.5 ℓ Coated Glass bottle	A/3240/PB17

NH₃,H₂O
CAS: 1336-21-6
MW: 35.05
Density at 20°C: 0.88
ADR/RID: 8, III
EINECS: 215-647-6
UN: 2672

GHS: Skin Corr. 1B, Eye Dam. 1, STOT SE 3, Aquatic Acute 1
H335, H314, H400
P301+P330+P331, P280, P305+P351+P338, P310, P261, P273

!

Corrosion

Environment

Assay: 32 to 36%

Guaranteed analysis

Colour	<20 APHA	Magnesium (Mg)	<2 ppm
Density at 20°C	: 0.88 to 0.89	Potassium (K)	<10 ppm
Calcium (Ca)	<10 ppm	Sodium (Na)	<20 ppm
Carbonate (Na ₂ CO ₃)	<0.01%	n-Butyl formate	<0.001%
Copper (Cu)	<1 ppm	Total phosphorus (P)	<20 ppm
Iron (Fe)	<1 ppm	Total sulphur (S)	<20 ppm
Lead (Pb)	<1 ppm	Zinc (Zn)	<2 ppm

Ammonia is equivalent to a strength of 18.1M (18.1N).

Ammonium acetate
for LC/MS analysis

NEW

Optima

packaging	Mfr. No
50 g Glass bottle	A114-50

C₂H₇NO₂
CAS: 631-61-8
MW: 77.08

Melting point: 110 to 112°C
Flash point: 136°C
EINECS: 211-162-9

Assay: >99%

Guaranteed analysis

Appearance	White fine crystals w/ lumps	Nickel (Ni)	<10 ppb
Water	<1%	Nitrate (NO ₃)	<0.001%
Aluminium (Al)	<20 ppb	Potassium (K)	<20 ppb
Calcium (Ca)	<50 ppb	Propan--1-ol	<10 ppb
Copper (Cu)	<10 ppb	Sodium (Na)	<50 ppb
Iron (Fe)	<10 ppb	n-Butyl formate	<0.0005%
Lead (Pb)	<10 ppb	Zinc (Zn)	<20 ppb
Magnesium (Mg)	<10 ppb	Residue after ignition	<0.01%
Manganese (Mn)	<10 ppb		

Absorbance specification

Max. Absorbance (1cm cell v HPLC water)

LC/MS Gradient Suitability : Pass test

This product has been tested on an LC-MS Thermo Scientific instrument to ensure suitability.

Ammonia solution
0.88 SG, 1:2 dilution

packaging	Mfr. No
500 mℓ HDPE Plastic bottle	J/9015/08

ADR/RID: 8, III
UN: 2672
GHS: Skin Corr. 1B, Eye Dam. 1, STOT SE 3, Aquatic Acute 1

H335, H410, H314
P301+P330+P331, P280, P305+P351+P338, P310, P261, P273

!

Corrosion

Environment

Ammonia-Nitrogen standard 1000 ppm
For Ion chromatography in H₂O, NIST certified

packaging	Mfr. No
100 mℓ Plastic bottle	SPASNH3N92Y

CAS: 12125-02-9

Manufactured to a tolerance of ±0.3%.

Check with your local distributor for availability of Spex product

Ammonium standard
for Ion pair chromatography, Ready-to-use

packaging	Mfr. No
100 mℓ HDPE Plastic bottle	J/4542/05

Manufactured to a tolerance of ± 1%.

High Purity Acids
for Elemental Analysis

Fisher Chemical offers a complete range of acids to meet the most challenging applications. Choose from our Optima®, TraceMetal® and PrimarPlus® grades to ensure that the sensitivity you require is always available.

Read more page 24

Ammonium acetate
HPLC for electrochemical detection

Certified

packaging	Mfr. No
250 g Plastic jar	A/3450/50

C₂H₇NO₂
CAS: 631-61-8
MW: 77.08

Melting point: 110 to 112°C
Flash point: 136°C
EINECS: 211-162-9

Assay: >97%

Guaranteed analysis

pH (aqueous solution 5%)	6.5 to 7.5	Lead (Pb)	<1 ppm
Electrochemical test	Pass	n-Butyl formate	<0.0005%
Iron (Fe)	<1 ppm		

Absorbance specification

Max. Absorbance (1cm cell v HPLC water)

254 nm (sol. aqu. 1 M) : <0.02 A.U.

Tested by HPLC-EC to ensure low electrochemical response.
Detector - ESA Coulochem II
Mode-oxidative
Graphite Electrode.
Lot analysis on the label.
Filtered to 0.2 µm.

Ammonium acetate
For soil analysis

packaging	Mfr. No
500 g Plastic jar	A/3445/53
1 kg Plastic jar	A/3445/60

C₂H₇NO₂
CAS: 631-61-8
MW: 77.08

Melting point: 110 to 112°C
Flash point: 136°C
EINECS: 211-162-9

Assay: >97%

Guaranteed analysis

pH (aqueous solution 5%)	6.5 to 7.5	Magnesium (Mg)	<5 ppm
Residue after ignition	<0.01%	Nitrate (NO ₃)	<0.001%
Aluminium (Al)	<5 ppm	Phosphate (PO ₄)	<1 ppm
Calcium (Ca)	<10 ppm	Potassium (K)	<2 ppm
Copper (Cu)	<5 ppm	n-Butyl formate	<0.0005%
Iron (Fe)	<1 ppm	Total sulphur (S)	<5 ppm
Lead (Pb)	<1 ppm		

Ammonium acetate
for analysis

Certified AR

packaging	Mfr. No
250 g Plastic jar	A/3440/50
500 g Plastic jar	A/3440/53
1 kg Plastic jar	A/3440/60
5 kg Plastic keg	A/3440/65

C₂H₇NO₂
CAS: 631-61-8
MW: 77.08

Melting point: 110 to 112°C
Flash point: 136°C
EINECS: 211-162-9

Assay: >99%

Guaranteed analysis

pH (aq. Solution 5%) at 25°C	: 6.7 to 7.3	Sodium (Na)	<5 ppm
Insoluble matter	<0.005%	Total chloride (Cl)	<0.0005%
Calcium (Ca)	<2 ppm	Total phosphorus (P)	<2 ppm
Copper (Cu)	<1 ppm	Total silicon (Si)	<5 ppm
Iron (Fe)	<1 ppm	Total sulphur (S)	<10 ppm
Lead (Pb)	<2 ppm	Zinc (Zn)	<1 ppm
Magnesium (Mg)	<1 ppm	Sulfate (SO ₄)	<0.001%
Nitrate (NO ₃)	<0.001%	Residue after ignition	<0.01%
Potassium (K)	<10 ppm		

Typical analysis

Aluminium (Al)	<1 ppm	Lithium (Li)	<1 ppm
Silver (Ag)	<1 ppm	Manganese (Mn)	<1 ppm
Arsenic (As)	<2 ppm	Molybdenum (Mo)	<1 ppm
Barium (Ba)	<1 ppm	Nickel (Ni)	<1 ppm
Cadmium (Cd)	<1 ppm	Strontium (Sr)	<1 ppm
Chromium (Cr)	<1 ppm	Titanium (Ti)	<1 ppm
Cobalt (Co)	<1 ppm	Vanadium (V)	<1 ppm

Lot analysis on the label.

Ammonium acetate
extra pure, crystals

SLR

packaging	Mfr. No
500 g Plastic jar	A/3400/53
1 kg Plastic jar	A/3400/60
5 kg Plastic keg	A/3400/65

C₂H₇NO₂
CAS: 631-61-8
MW: 77.08

Melting point: 110 to 112°C
Flash point: 136°C
EINECS: 211-162-9

Assay: >97%

Guaranteed analysis

pH (aqueous solution 5%)	6.5 to 7.5	Sodium (Na)	<20 ppm
Calcium (Ca)	<20 ppm	n-Butyl formate	<0.001%
Copper (Cu)	<5 ppm	Total phosphorus (P)	<50 ppm
Iron (Fe)	<5 ppm	Total silicon (Si)	<50 ppm
Lead (Pb)	<5 ppm	Total sulphur (S)	<100 ppm
Magnesium (Mg)	<10 ppm	Zinc (Zn)	<10 ppm
Potassium (K)	<50 ppm		

Ammonium acetate solution 0.25M
HPLC for gradient analysis

Certified

packaging	Mfr. No
500 mℓ HDPE plastic bottle	J/0200/PB08

C₂H₇NO₃
CAS: 631-61-8
MW: 77.08

EINECS: 211-162-9
P262

Absorbance specification

Max. Absorbance (1cm cell v HPLC water)

Drift in gradient baseline : <40 mA.U.

Absorbance of the largest eluted peak: <1.5 mA.U.

Manufactured to a tolerance of ±0.02M

Conditions:

Column: C18 Reverse Phase

Flow Rate: 3ml/min

Wavelength: 254nm

Solvent A: Water + 10%Buffer

Solvent B: Acetonitrile + 10% Buffer

Gradient: 10 min. At 5% solvent B followed by a linear gradient to 95% solvent B over 15min.

Ammonium bromide
for analysis

Certified AR

packaging	Mfr. No
500 g Plastic jar	A/3560/53

NH₄Br
CAS: 12124-97-9
MW: 97.95
Density at 20°C: 2.42
Boiling point: 235°C

Melting point: 452°C
EINECS: 235-183-8
GHS: Eye Irrit. 2
H319
P280, P305+P351+P338

Assay: >99.5%

Guaranteed analysis

pH (aqueous solution 5%)	4.8 to 6	Magnesium (Mg)	<5 ppm
Loss on drying	<0.1%	Potassium (K)	<50 ppm
Bromate (BrO ₃)	<0.001%	Sodium (Na)	<20 ppm
Calcium (Ca)	<10 ppm	Total chloride (Cl)	<0.2%
Copper (Cu)	<1 ppm	Total phosphorus (P)	<10 ppm
Iodide (I)	<0.005%	Total silicon (Si)	<10 ppm
Iron (Fe)	<1 ppm	Total sulphur (S)	<50 ppm
Lead (Pb)	<2 ppm	Zinc (Zn)	<5 ppm

Typical analysis

Aluminium (Al)	<2 ppm	Lithium (Li)	<2 ppm
Silver (Ag)	<2 ppm	Manganese (Mn)	<2 ppm
Arsenic (As)	<2 ppm	Molybdenum (Mo)	<2 ppm
Barium (Ba)	<2 ppm	Nickel (Ni)	<2 ppm
Cadmium (Cd)	<2 ppm	Strontium (Sr)	<2 ppm
Chromium (Cr)	<2 ppm	Titanium (Ti)	<2 ppm
Cobalt (Co)	<2 ppm	Vanadium (V)	<2 ppm

Ammonium bromide
specified reagents for general laboratory work- extra pure

SLR

packaging	Mfr. No
500 g Plastic jar	A/3520/53

NH₄Br
CAS: 12124-97-9
MW: 97.95
Density at 20°C: 2.42
Boiling point: 235°C

Melting point: 452°C
EINECS: 235-183-8
GHS: Eye Irrit. 2
H319
P280, P305+P351+P338

Assay: >99%

Guaranteed analysis

Calcium (Ca)	<50 ppm	Sodium (Na)	<20 ppm
Copper (Cu)	<10 ppm	n-Butyl formate	<0.5%
Iron (Fe)	<10 ppm	Total phosphorus (P)	<100 ppm
Lead (Pb)	<10 ppm	Total silicon (Si)	<100 ppm
Magnesium (Mg)	<50 ppm	Total sulphur (S)	<500 ppm
Potassium (K)	<500 ppm	Zinc (Zn)	<50 ppm

36

Application For analysis Specified Titration

37

Ammonium carbonate		Certified
for HPLC		
packaging	Mfr. No	
250 g Plastic jar	A/3686/50	
NH ₄ HCO ₃ +NH ₂ .COO.NH ₄	EINECS: 208-058-0	
CAS: 506-87-6	GHS: Acute Tox. 4	
MW: 96.09	H302	
Melting point: 58°C	P301+P312	
Assay: >31%		
Guaranteed analysis		
Iron (Fe)	<1 ppm n-Butyl formate	<0.0001%
Lead (Pb)	<1 ppm	
Absorbance specification		
Max. Absorbance (1 cm cell v HPLC water)		
254 nm (sol. aqu. 1 M)		: <0.02 A.U.
Lot analysis on the label.		
Filtered to 0.2 µm.		

Ammonium carbonate		Certified AR
for analysis, powder		
packaging	Mfr. No	
500 g Plastic jar	A/3680/53	
1 kg Plastic jar	A/3680/50	
NH ₄ HCO ₃ .NH ₂ .COO.NH ₄	EINECS: 208-058-0	
CAS: 506-87-6	GHS: Acute Tox. 4	
MW: 96.09	H302	
Melting point: 58°C	P301+P312	
Assay: >31%		
Guaranteed analysis		
Calcium (Ca)	<10 ppm Sodium (Na)	<5 ppm
Copper (Cu)	<0.5 ppm Total chloride (Cl)	<0.0005%
Iron (Fe)	<10 ppm Total phosphorus (P)	<5 ppm
Lead (Pb)	<5 ppm Total silicon (Si)	<10 ppm
Magnesium (Mg)	<5 ppm Total sulphur (S)	<10 ppm
Potassium (K)	<5 ppm Zinc (Zn)	<5 ppm
Typical analysis		
Aluminium (Al)	<1 ppm Lithium (Li)	<1 ppm
Silver (Ag)	<1 ppm Manganese (Mn)	<1 ppm
Arsenic (As)	<1 ppm Molybdenum (Mo)	<1 ppm
Barium (Ba)	<1 ppm Nickel (Ni)	<1 ppm
Cadmium (Cd)	<1 ppm Strontium (Sr)	<1 ppm
Chromium (Cr)	<1 ppm Titanium (Ti)	<1 ppm
Cobalt (Co)	<1 ppm Vanadium (V)	<1 ppm
Lot analysis on the label.		

Ammonium carbonate		SLR
extra pure, powder		
packaging	Mfr. No	
500 g Plastic jar	A/3640/53	
1 kg Plastic jar	A/3640/60	
5 kg Plastic keg	A/3640/65	
NH ₄ HCO ₃ +NH ₂ .COO.NH ₄	EINECS: 208-058-0	
CAS: 506-87-6	GHS: Acute Tox. 4	
MW: 96.09	H302	
Melting point: 58°C	P301+P312	
Assay: >30%		
Guaranteed analysis		
Calcium (Ca)	<50 ppm Sodium (Na)	<50 ppm
Copper (Cu)	<5 ppm n-Butyl formate	<0.001%
Iron (Fe)	<20 ppm Total phosphorus (P)	<50 ppm
Lead (Pb)	<5 ppm Total silicon (Si)	<50 ppm
Magnesium (Mg)	<20 ppm Total sulphur (S)	<50 ppm
Potassium (K)	<20 ppm Zinc (Zn)	<50 ppm

Ammonium cerium(IV) nitrate		Certified AR
for analysis		
packaging	Mfr. No	  
100 g Plastic jar	A/3730/48	
250 g Plastic jar	A/3730/50	
(NH ₄) ₂ Ce(NO ₃) ₆	GHS: Ox. Sol. 2, Acute Tox. 4, Skin Irrit. 2, Eye Dam. 1, STOT SE 3	
CAS: 16774-21-3	H272, H335, H318, H302, H315	
MW: 548.23	P280, P305+P351+P338, P310, P210, P261, P301+P312	
Melting point: 107 to 108°C		
ADR/RID: 5.1, II		
EINECS: 240-827-6		
UN: 1477		
Assay: >99%		
Guaranteed analysis		
Loss on drying	<3% n-Butyl formate	<0.001%
Copper (Cu)	<10 ppm Total phosphorus (P)	<50 ppm
Iron (Fe)	<20 ppm Total silicon (Si)	<20 ppm
Lead (Pb)	<10 ppm Total sulphur (S)	<500 ppm
Magnesium (Mg)	<50 ppm Zinc (Zn)	<20 ppm
Potassium (K)	<50 ppm Insoluble Acids	<0.005%
Sodium (Na)	<50 ppm	
Lot analysis on the label.		

Ammonium cerium(IV) nitrate		SLR
specified reagents for general laboratory work- extra pure		
packaging	Mfr. No	  
100 g Plastic jar	A/3720/48	
250 g Plastic jar	A/3720/50	
(NH ₄) ₂ Ce(NO ₃) ₆	GHS: Ox. Sol. 2, Acute Tox. 4, Skin Irrit. 2, Eye Dam. 1, STOT SE 3	
CAS: 16774-21-3	H272, H335, H318, H302, H315	
MW: 548.23	P280, P305+P351+P338, P310, P210, P261, P301+P312	
Melting point: 107 to 108°C		
ADR/RID: 5.1, II		
EINECS: 240-827-6		
UN: 1477		
Assay: >98%		
Guaranteed analysis		
Loss on drying	<5% Sodium (Na)	<500 ppm
Copper (Cu)	<100 ppm n-Butyl formate	<0.005%
Iron (Fe)	<100 ppm Total phosphorus (P)	<500 ppm
Lead (Pb)	<100 ppm Total silicon (Si)	<200 ppm
Magnesium (Mg)	<200 ppm Total sulphur (S)	<500 ppm
Potassium (K)	<500 ppm Zinc (Zn)	<200 ppm

Ammonium cerium(IV) sulfate		Certified AR
for analysis, low in rare earths		
packaging	Mfr. No	
100 g Plastic jar	A/3771/48	
(NH ₄) ₄ Ce(SO ₄) ₄ .2H ₂ O	Melting point: 130°C	
CAS: 10378-47-9	EINECS: 231-567-4	
MW: 632.56		
Assay: >99%		
Guaranteed analysis		
Copper (Cu)	<10 ppm n-Butyl formate	<0.001%
Iron (Fe)	<50 ppm Total phosphorus (P)	<50 ppm
Lead (Pb)	<20 ppm Total silicon (Si)	<20 ppm
Magnesium (Mg)	<50 ppm Zinc (Zn)	<20 ppm
Potassium (K)	<50 ppm Other rare earths	<2%
Sodium (Na)	<50 ppm	
Lot analysis on the label.		



Ammonium cerium(IV) sulfate		SLR
specified reagents for general laboratory work- extra pure		
packaging	Mfr. No	
100 g Plastic jar	A/3760/48	
(NH ₄) ₄ Ce(SO ₄) ₄ .2H ₂ O	Melting point: 130°C	
CAS: 10378-47-9	EINECS: 231-567-4	
MW: 632.56		
Assay: 90 to 105%		
Guaranteed analysis		
Copper (Cu)	<500 ppm Sodium (Na)	<500 ppm
Iron (Fe)	<500 ppm n-Butyl formate	<0.05%
Lead (Pb)	<500 ppm Total phosphorus (P)	<500 ppm
Magnesium (Mg)	<500 ppm Total silicon (Si)	<500 ppm
Potassium (K)	<500 ppm Zinc (Zn)	<500 ppm

Ammonium chloride		Certified AR, Ph.Eur., BP
for analysis		
packaging	Mfr. No	
500 g Plastic jar	A/3920/53	
1 kg Plastic jar	A/3920/60	
2 kg Plastic jar	A/3920/61	
5 kg Plastic keg	A/3920/65	
25 kg Fibre Keg	A/3920/70	
NH ₄ Cl	GHS: Acute Tox. 4, Eye Irrit. 2	
CAS: 12125-02-9	H319, H302	
MW: 53.49	P270, P264, P280,	
Melting point: 340°C	P305+P351+P338, P301+P312	
EINECS: 235-186-4		
Assay: 99.5 to 100.5%		
Guaranteed analysis		
pH (aqueous solution 5%)	4.5 to 5.5 Nitrate (NO ₃)	<0.0005%
Heavy metals	<10 ppm Potassium (K)	<50 ppm
Sulfated ash	<0.1% Sodium (Na)	<50 ppm
Calcium (Ca)	<5 ppm Total phosphorus (P)	<5 ppm
Copper (Cu)	<1 ppm Total silicon (Si)	<2 ppm
Iodide (I)	<0.001% Total sulphur (S)	<10 ppm
Iron (Fe)	<2 ppm Zinc (Zn)	<2 ppm
Lead (Pb)	<2 ppm Pyridine & homologues	<10 ppm
Magnesium (Mg)	<1 ppm	
Typical analysis		
Aluminium (Al)	<1 ppm Manganese (Mn)	<1 ppm
Silver (Ag)	<1 ppm Molybdenum (Mo)	<1 ppm
Barium (Ba)	<1 ppm Nickel (Ni)	<1 ppm
Cadmium (Cd)	<1 ppm Strontium (Sr)	<1 ppm
Chromium (Cr)	<1 ppm Titanium (Ti)	<1 ppm
Cobalt (Co)	<1 ppm Vanadium (V)	<1 ppm
Lithium (Li)	<1 ppm	

Ammonium chloride		SLR
specified reagents for general laboratory work- extra pure		
packaging	Mfr. No	
500 g Plastic jar	A/3880/53	
1 kg Plastic jar	A/3880/60	
2 kg Plastic jar	A/3880/61	
5 kg Plastic keg	A/3880/65	
25 kg Plastic keg	A/3880/70	
NH ₄ Cl	GHS: Acute Tox. 4, Eye Irrit. 2	
CAS: 12125-02-9	H319, H302	
MW: 53.49	P270, P264, P280,	
Melting point: 340°C	P305+P351+P338, P301+P312	
EINECS: 235-186-4		
Assay: >99%		
Guaranteed analysis		
Calcium (Ca)	<50 ppm Sodium (Na)	<200 ppm
Copper (Cu)	<5 ppm Total phosphorus (P)	<20 ppm
Iron (Fe)	<10 ppm Total silicon (Si)	<20 ppm
Lead (Pb)	<10 ppm Total sulphur (S)	<100 ppm
Magnesium (Mg)	<10 ppm Zinc (Zn)	<50 ppm
Potassium (K)	<200 ppm	

Ammonium chloride		
Technical, powder		
packaging	Mfr. No	
2 kg Plastic jar	A/3800/61	
NH ₄ Cl	GHS: Acute Tox. 4, Eye Irrit. 2	
CAS: 12125-02-9	H319, H302	
MW: 53.49	P270, P264, P280,	
Melting point: 340°C	P305+P351+P338, P301+P312	
EINECS: 235-186-4		
Guaranteed analysis		
Identity test		: Pass

tri-Ammonium citrate		Certified AR
for analysis		
packaging	Mfr. No	
250 g Plastic jar	A/4082/50	
500 g Plastic jar	A/4082/53	
25 kg Fibre Keg	A/4082/70	
(NH ₄) ₃ C ₆ H ₅ O ₇	GHS: Skin Irrit. 2, Eye Irrit. 2,	
CAS: 3458-72-8	STOT SE 3	
MW: 243.22	H315, H319, H335	
Melting point: 185°C	P280, P305+P351+P338	
EINECS: 222-394-5		
Assay: 99 to 101%		
Guaranteed analysis		
pH (aqueous solution 5%)	6 to 7 Sodium (Na)	<50 ppm
Calcium (Ca)	<20 ppm n-Butyl formate	<0.001%
Copper (Cu)	<1 ppm Total phosphorus (P)	<20 ppm
Iron (Fe)	<5 ppm Total silicon (Si)	<20 ppm
Lead (Pb)	<5 ppm Total sulphur (S)	<20 ppm
Magnesium (Mg)	<5 ppm Zinc (Zn)	<2 ppm
Oxalate (C ₂ O ₄)	<0.01% Readily carbonisable substances test: Pass	
Potassium (K)	<20 ppm Reducing substance	Pass

tri-Ammonium citrate		SLR
specified reagents for general laboratory work- extra pure		
packaging	Mfr. No	
500 g Plastic jar	A/4080/53	
1 kg Plastic jar	A/4080/60	
25 kg Fibre Keg	A/4080/70	
(NH ₄) ₃ C ₆ H ₅ O ₇	GHS: Skin Irrit. 2, Eye Irrit. 2,	
CAS: 3458-72-8	STOT SE 3	
MW: 243.22	H315, H319, H335	
Melting point: 185°C	P280, P305+P351+P338	
EINECS: 222-394-5		
Assay: >98%		
Guaranteed analysis		
pH (aqueous solution 5%)	6 to 7 Sodium (Na)	<500 ppm
Calcium (Ca)	<50 ppm n-Butyl formate	<0.005%
Copper (Cu)	<10 ppm Total phosphorus (P)	<100 ppm
Iron (Fe)	<20 ppm Total silicon (Si)	<100 ppm
Lead (Pb)	<10 ppm Total sulphur (S)	<100 ppm
Magnesium (Mg)	<20 ppm Zinc (Zn)	<20 ppm
Potassium (K)	<100 ppm	

Ammonium copper(II) chloride dihydrate		SLR
specified reagents for general laboratory work- extra pure		
packaging	Mfr. No	
500 g Plastic jar	A/4240/53	
(NH ₄) ₂ CuCl ₄ .2H ₂ O		
CAS: 15610-76-1		
MW: 277.47		
EINECS: 239-690-5		
Assay: >98%		
Guaranteed analysis		
Calcium (Ca)	<50 ppm Potassium (K)	<500 ppm
Iron (Fe)	<100 ppm Sodium (Na)	<500 ppm
Lead (Pb)	<500 ppm Zinc (Zn)	<500 ppm
Magnesium (Mg)	<500 ppm	

Ammonium dichromate NEW		SLR
specified reagents for general laboratory work- extra pure		
packaging		Mfr. No
250 g	Plastic jar	A/4420/50
500 g	Plastic jar	A/4420/53
1 kg	Plastic jar	A/4420/60
(NH₄)CrO₇ CAS: 7789-09-5 MW: 252.70 Density at 20°C: 2.15 Melting point: 170°C Flash point: 218°C ADR/RID: 5.1, II EINECS: 232-143-1 UN: 1439 GHS: Carc. 1B, Muta. 1B, Repr. 1B, Acute Tox. 2, Skin		Corr. 1B, Resp. Sens. 1, Skin Sens. 1, STOT RE 1, Aquatic Acute 1, Aquatic Chronic 1 H272, H350, H340, H360FD, H330, H301, H372, H312, H314, H334, H317, H410 P280, P301+P330+P331, P302+P352, P305+P351+P338, P304+P340, P342+P311 

Assay: > 99.5%		
Guaranteed analysis		
Appearance	orange red crystals	Chloride (Cl) <0.005%
Loss on drying (105°C /2h)	: <3.0%	Iron (Fe) <0.002%
Insoluble matter	<0.005%	Sodium (Na) <0.005%
Calcium (Ca)	<0.002%	Sulfate (SO ₄) <0.01%

Ammonium dihydrogen orthophosphate			Certified AR
for analysis			
packaging			Mfr. No
500 g	Plastic jar		A/4560/53
1 kg	Plastic jar		A/4560/60
(NH ₄)H ₂ PO ₄		Melting point: 190°C	
CAS: 7722-76-1		EINECS: 231-764-5	
MW: 115.03			
Assay: >99%			
Guaranteed analysis			
pH (aqueous solution 5%)	4 to 4.5	Potassium (K)	<20 ppm
Calcium (Ca)	<10 ppm	Sodium (Na)	<50 ppm
Copper (Cu)	<5 ppm	Total chloride (Cl)	<0.0005%
Iron (Fe)	<5 ppm	Total silicon (Si)	<10 ppm
Lead (Pb)	<5 ppm	Total sulphur (S)	<20 ppm
Magnesium (Mg)	<10 ppm	Zinc (Zn)	<5 ppm
Nitrate (NO ₃)	<0.001%		
Typical analysis			
Aluminium (Al)	<2 ppm	Lithium (Li)	<2 ppm
Silver (Ag)	<2 ppm	Manganese (Mn)	<2 ppm
Arsenic (As)	<2 ppm	Molybdenum (Mo)	<2 ppm
Barium (Ba)	<2 ppm	Nickel (Ni)	<2 ppm
Cadmium (Cd)	<2 ppm	Strontium (Sr)	<2 ppm
Chromium (Cr)	<5 ppm	Titanium (Ti)	<2 ppm
Cobalt (Co)	<2 ppm	Vanadium (V)	<2 ppm
Lot analysis on the label.			

Ammonium dihydrogen orthophosphate		SLR
specified reagents for general laboratory work- extra pure		
packaging		Mfr. No
500 g	Plastic jar	A/4520/53
(NH ₄)H ₂ PO ₄ CAS: 7722-76-1 MW: 115.03	Melting point: 190°C EINECS: 231-764-5	
Guaranteed analysis		
Identity test		: Pass

Ammonium fluoride		Certified AR	
for analysis			
packaging		Mfr. No	
500 g	Plastic jar	A/4960/53	
NH ₄ F	UN: 2505		
CAS: 12125-01-8	GHS: Acute Tox. 3		
MW: 37.04	H311, H331, H301		
ADR/RID: 6.1, III	P260, P280, P302+P350,		
EINECS: 235-185-9	P304+P340, P301+P310		
Assay: >96%			
Guaranteed analysis			
Reaction (0.5% solution)	: 6 to 7	Sodium (Na)	<10 ppm
Calcium (Ca)	<5 ppm	n-Butyl formate	<0.0005%
Copper (Cu)	<2 ppm	Total phosphorus (P)	<10 ppm
Iron (Fe)	<5 ppm	Total silicon (Si)	<1000 ppm
Lead (Pb)	<5 ppm	Total sulphur (S)	<50 ppm
Magnesium (Mg)	<2 ppm	Zinc (Zn)	<2 ppm
Potassium (K)	<20 ppm		
Lot analysis on the label.			

Ammonium fluoride		SLR	
specified reagents for general laboratory work- extra pure			
packaging		Mfr. No	
500 g	Plastic jar	A/4920/53	
NH ₄ F	UN: 2505		
CAS: 12125-01-8	GHS: Acute Tox. 3		
MW: 37.04	H311, H331, H301		
ADR/RID: 6.1, III	P260, P280, P302+P350,		
EINECS: 235-185-9	P304+P340, P301+P310		
Assay: >95%			
Guaranteed analysis			
Reaction (0.5% solution)	: 6 to 7	Sodium (Na)	<100 ppm
Calcium (Ca)	<20 ppm	n-Butyl formate	<0.002%
Copper (Cu)	<10 ppm	Total phosphorus (P)	<50 ppm
Iron (Fe)	<20 ppm	Total silicon (Si)	<2000 ppm
Lead (Pb)	<20 ppm	Total sulphur (S)	<200 ppm
Magnesium (Mg)	<10 ppm	Zinc (Zn)	<10 ppm
Potassium (K)	<100 ppm		

Ammonium formate NEW		Optima
for LC/MS analysis		
packaging		Mfr. No
50 g	Glass bottle	A115-50
HCOONH ₄ CAS: 540-69-2 MW: 63.06 Melting point: 115 to 120°C EINECS: 208-753-9	GHS: Skin Irrit. 2, Eye Irrit. 2, STOT SE 3 H315, H319, H335 P261, P302+P352, P280, P305+P351+P338	
		
Assay: 99%		
Guaranteed analysis		
Appearance. White fine crystals w/ lumps	Nickel (Ni)	<10 ppb
Water	Potassium (K)	<20 ppb
Aluminium (Al)	Propan--1-ol	<10 ppb
Calcium (Ca)	Sodium (Na)	<50 ppb
Copper (Cu)	Sulfate (SO ₄)	<0.005%
Iron (Fe)	n-Butyl formate	<0.001%
Lead (Pb)	Zinc (Zn)	<20 ppb
Magnesium (Mg)	Residue after ignition	<0.005%
Manganese (Mn)		
Absorbance specification		
Max. Absorbance (1cm cell v HPLC water)		
LC/MS Gradient Suitability	: Pass test	
Protease (0.1% aqueous) : Not Detected		
This product has been tested on an LC-MS Thermo Scientific instrument to ensure suitability.		

Ammonium formate		Certified AR
for analysis		
packaging		Mfr. No
500 g	Plastic jar	A/5080/53
HCOONH ₄	GHS: Skin Irrit. 2, Eye Irrit. 2, STOT SE 3	
CAS: 540-69-2	H315, H319, H335	
MW: 63.06	P261, P302+P352, P280,	
Melting point: 115 to 120°C	P305+P351+P338	
EINECS: 208-753-9		
Assay: >99%		
Guaranteed analysis		
pH (aqueous solution 5%)	6 to 7	Sodium (Na) <5 ppm
Calcium (Ca)	<5 ppm	Total chloride (Cl) <0.001%
Copper (Cu)	<1 ppm	Total phosphorus (P) <10 ppm
Iron (Fe)	<2 ppm	Total silicon (Si) <10 ppm
Lead (Pb)	<2 ppm	Total sulphur (S) <50 ppm
Magnesium (Mg)	<2 ppm	Zinc (Zn) <2 ppm
Potassium (K)	<10 ppm	
Typical analysis		
Aluminium (Al)	<1 ppm	Lithium (Li) <1 ppm
Silver (Ag)	<1 ppm	Manganese (Mn) <1 ppm
Arsenic (As)	<2 ppm	Molybdenum (Mo) <1 ppm
Barium (Ba)	<1 ppm	Nickel (Ni) <1 ppm
Cadmium (Cd)	<1 ppm	Strontium (Sr) <1 ppm
Chromium (Cr)	<1 ppm	Titanium (Ti) <1 ppm
Cobalt (Co)	<1 ppm	Vanadium (V) <1 ppm
Lot analysis on the label.		

Ammonium formate		SLR	
specified reagents for general laboratory work- extra pure			
packaging		Mfr. No	
500 g Plastic jar		A/5040/53	
HCOONH ₄ CAS: 540-69-2 MW: 63.06 Melting point: 115 to 120°C EINECS: 208-753-9	GHS: Skin Irrit. 2, Eye Irrit. 2, STOT SE 3 H315, H319, H335 P261, P302+P352, P280, P305+P351+P338		
			
Assay: >98%			
Guaranteed analysis			
pH (aqueous solution 5%)	6 to 7	Sodium (Na)	<50 ppm
Calcium (Ca)	<20 ppm	n-Butyl formate	<0.005%
Copper (Cu)	<10 ppm	Total phosphorus (P)	<50 ppm
Iron (Fe)	<10 ppm	Total silicon (Si)	<50 ppm
Lead (Pb)	<10 ppm	Total sulphur (S)	<200 ppm
Magnesium (Mg)	<10 ppm	Zinc (Zn)	<10 ppm
Potassium (K)	<50 ppm		

Ammonium hydrogen carbonate		Certified AR	
for analysis			
packaging		Mfr. No	
250 g	Plastic jar	A/5160/50	
1 kg	Plastic jar	A/5160/60	
25 kg	Fibre Keg	A/5160/70	
NH ₄ HCO ₃	EINECS: 213-911-5		
CAS: 1066-33-7	GHS: Acute Tox. 4		
MW: 79.06	H302		
Melting point: 105°C	P264, P270, P301+P312, P330		
Assay: >99%			
Guaranteed analysis			
Calcium (Ca)	<5 ppm	Sodium (Na)	<10 ppm
Copper (Cu)	<2 ppm	Total chloride (Cl)	<0.0005%
Iron (Fe)	<1 ppm	Total phosphorus (P)	<10 ppm
Lead (Pb)	<2 ppm	Total silicon (Si)	<10 ppm
Magnesium (Mg)	<20 ppm	Total sulphur (S)	<10 ppm
Potassium (K)	<10 ppm	Zinc (Zn)	<5 ppm
Typical analysis			
Aluminium (Al)	<1 ppm	Manganese (Mn)	<1 ppm
Silver (Ag)	<1 ppm	Molybdenum (Mo)	<1 ppm
Barium (Ba)	<1 ppm	Nickel (Ni)	<1 ppm
Cadmium (Cd)	<1 ppm	Strontium (Sr)	<1 ppm
Chromium (Cr)	<1 ppm	Titanium (Ti)	<1 ppm
Cobalt (Co)	<1 ppm	Vanadium (V)	<1 ppm
Lithium (Li)	<1 ppm		
Lot analysis on the label.			

Ammonium hydrogen carbonate		SLR	
specified reagents for general laboratory work- extra pure			
packaging	Mfr. No		
1 kg Plastic jar	A/5120/60		
NH ₄ HCO ₃ CAS: 1066-33-7 MW: 79.06 Melting point: 105°C	EINECS: 213-911-5 GHS: Acute Tox. 4 H302 P264, P270, P301+P312, P330		
Assay: >98%			
Guaranteed analysis			
Calcium (Ca)	<20 ppm	Sodium (Na)	<50 ppm
Copper (Cu)	<10 ppm	n-Butyl formate	<0.002%
Iron (Fe)	<10 ppm	Total phosphorus (P)	<50 ppm
Lead (Pb)	<10 ppm	Total silicon (Si)	<50 ppm
Magnesium (Mg)	<50 ppm	Total sulphur (S)	<50 ppm
Potassium (K)	<50 ppm	Zinc (Zn)	<20 ppm

di-Ammonium hydrogen citrate		SLR
specified reagents for general laboratory work- extra pure		
packaging		Mfr. No
500 g Plastic jar		A/5200/53
C ₆ H ₁₄ N ₂ O ₇ CAS: 3012-65-5 MW: 226.19 Density at 20°C: 1.5	EINECS: 221-146-3 GHS: Eye Irrit. 2, STOT SE 3 H319, H335 P261, P280, P305+P351+P338	
Assay: >98%		
Guaranteed analysis		
Lead (Pb)	<50 ppm	n-Butyl formate
Sulfate (SO ₄)	<0.005%	<0.005%

Ammonium hydrogen difluoride			SLR	
specified reagents for general laboratory work- extra pure				
packaging			Mfr. No	
500 g	Plastic jar		A/5280/53	
2.5 kg	Plastic jar		A/5280/62	
25 kg	Fibre Keg		A/5280/70	
NH ₄ HF ₃ CAS: 1341-49-7 MW: 57.05 Density at 20°C: 1.50 Boiling point: 230°C Melting point: 125°C ADR/RID: 8, II		EINECS: 215-676-4 UN: 1727 GHS: Acute Tox. 3, Skin Corr. 1B, Eye Dam. 1 H314, H301 P301+P330+P331, P280, P305+P351+P338, P310		
				
				
Assay: >97%				
Guaranteed analysis				
Residue after ignition		<0.02%	Lead (Pb)	<100 ppm
Solution		: Pass	n-Butyl formate	<0.002%
Iron (Fe)		<100 ppm		

Ammonium hydrogen difluoride		
Technical		
packaging		Mfr. No
3 kg	Plastic jar	A/5240/63
NH ₄ HF ₂	EINECS: 215-676-4	
CAS: 1341-49-7	UN: 1727	
MW: 57.05	GHS: Acute Tox. 3, Skin	
Density at 20°C: 1.50	Corr. 1B, Eye Dam. 1	
Boiling point: 230°C	H314, H301	
Melting point: 125°C	P301+P330+P331, P280,	
ADR/RID: 8, II	P305+P351+P338, P310	
Guaranteed analysis		
Identity test : Pass		

di-Ammonium hydrogen orthophosphate Certified AR
for analysis

packaging		Mfr. No		
500 g	Plastic jar	A/5360/53		
1 kg	Plastic jar	A/5360/60		
(NH ₄) ₂ HPO ₄ CAS: 7783-28-0 MW: 132.06 Density at 20°C: 1.619 EINECS: 231-987-8		GHS: Skin Irrit. 2, Eye Irrit. 2, STOT SE 3 H335, H315, H319 P261, P302+P352, P280, P305+P351+P338		
Assay: >98%				
Guaranteed analysis				
pH (aqueous solution 5%)	7.8 to 8.2	Potassium (K)	<10 ppm	
Calcium (Ca)	<5 ppm	Sodium (Na)	<10 ppm	
Copper (Cu)	<5 ppm	Total chloride (Cl)	<0.0005%	
Iron (Fe)	<5 ppm	Total silicon (Si)	<50 ppm	
Lead (Pb)	<5 ppm	Total sulphur (S)	<50 ppm	
Magnesium (Mg)	<5 ppm	Zinc (Zn)	<5 ppm	
Nitrate (NO ₃)	<0.001%			
Typical analysis				
Aluminium (Al)	<1 ppm	Lithium (Li)	<1 ppm	
Silver (Ag)	<1 ppm	Manganese (Mn)	<1 ppm	
Arsenic (As)	<1 ppm	Molybdenum (Mo)	<1 ppm	
Barium (Ba)	<1 ppm	Nickel (Ni)	<1 ppm	
Cadmium (Cd)	<1 ppm	Strontium (Sr)	<1 ppm	
Chromium (Cr)	<1 ppm	Titanium (Ti)	<1 ppm	
Cobalt (Co)	<1 ppm	Vanadium (V)	<1 ppm	

di-Ammonium hydrogen orthophosphate SLR
specified reagents for general laboratory work- extra pure

Packaging		Mfr. No	
500 g	Plastic jar	A/5320/53	
1 kg	Plastic jar	A/5320/60	
(NH ₄) ₂ HPO ₄ CAS: 7783-28-0 MW: 132.06 Density at 20°C: 1.619 EINECS: 231-987-8		GHS: Skin Irrit. 2, Eye Irrit. 2, STOT SE 3 H335, H315, H319 P261, P302+P352, P280, P305+P351+P338	
			
Assay: >97%			
Guaranteed analysis			
Calcium (Ca)	<20 ppm	Sodium (Na)	<200 ppm
Copper (Cu)	<20 ppm	n-Butyl formate	<0.002%
Iron (Fe)	<20 ppm	Total silicon (Si)	<200 ppm
Lead (Pb)	<20 ppm	Total sulphur (S)	<200 ppm
Magnesium (Mg)	<20 ppm	Zinc (Zn)	<20 ppm
Potassium (K)	<50 ppm		

Ammonium hydrogen sulfate SLR
specified reagents for general laboratory work- extra pure

packaging		Mfr. No	
500 g	Plastic jar	A/5400/53	
NH ₄ SO ₄ CAS: 7803-63-6 MW: 115.11 Boiling point: 350°C Melting point: 147°C ADR/RID: 8, II	EINECS: 232-265-5 UN: 2506 GHS: Skin Corr. 1B, Eye Dam. 1 H314 P301+P330+P331, P280, P305+P351+P338, P310		
Assay: >98%			
Guaranteed analysis			
Calcium (Ca)	<50 ppm	Magnesium (Mg)	<100 ppm
Copper (Cu)	<50 ppm	Potassium (K)	<300 ppm
Iron (Fe)	<50 ppm	Sodium (Na)	<100 ppm
Lead (Pb)	<50 ppm	Zinc (Zn)	<50 ppm

Ammonium iron (III) citrate SLR
specified reagents for general laboratory work- extra pure

Packaging		Mfr. No	
500 g	Plastic jar	A/4640/53	
C ₆ H ₁₁ FeNO ₇ CAS: 1185-57-5 MW: 265.00 INECS: 214-686-6			
Assay: 20.5 to 22.5%			
Guaranteed analysis			
Lead (Pb)	<50 ppm	Zinc (Zn)	<1000 ppm

Ammonium Iron(II) sulfate hexahydrate Certified AR
for analysis

Packaging			Mfr. No
50 g	Plastic jar		A/4880/50
100 g	Plastic jar		A/4880/53
1 kg	Plastic jar		A/4880/60
3 kg	Plastic jar		A/4880/63
(NH ₄) ₂ Fe(SO ₄) ₂ ·6H ₂ O		GHS: Skin Irrit. 2, Eye Irrit. 2,	
CAS: 7783-85-9		STOT SE 3	
MW: 392.14		H315, H319, H335	
Melting point: 100°C		P261, P302+P352, P280,	
INECS: 233-151-8		P305+P351+P338	
Assay: >99%			
Guaranteed analysis			
pH (aqueous solution 5%)	3 to 5	Potassium (K)	<50 ppm
Calcium (Ca)	<10 ppm	Sodium (Na)	<50 ppm
Copper (Cu)	<10 ppm	Total chloride (Cl)	<0.001%
Iron (III)	<0.03%	Total phosphorus (P)	<30 ppm
Lead (Pb)	<10 ppm	Total silicon (Si)	<10 ppm
Magnesium (Mg)	<50 ppm	Zinc (Zn)	<50 ppm
Typical analysis			
Aluminium (Al)	<2 ppm	Manganese (Mn)	<0.02 %
Silver (Ag)	<2 ppm	Molybdenum (Mo)	<2 ppm
Arsenic (As)	<5 ppm	Nickel (Ni)	<0.001 %
Barium (Ba)	<2 ppm	Strontium (Sr)	<2 ppm
Chromium (Cr)	<2 ppm	Titanium (Ti)	<2 ppm
Cobalt (Co)	<0.001 %	Vanadium (V)	<2 ppm
Lithium (Li)	<2 ppm		
Lot analysis on the label.			

Ammonium iron(II) sulfate hexahydrate SLR
specified reagents for general laboratory work- extra pure

Packaging		Mfr. No	
500 g	Plastic jar	A/4840/53	
1 kg	Plastic jar	A/4840/60	
3 kg	Plastic jar	A/4840/63	
25 kg	Plastic keg	A/4840/70	
$(\text{NH}_4)_2\text{Fe}(\text{SO}_4)_2 \cdot 6\text{H}_2\text{O}$ CAS: 7783-85-9 MW: 392.14 Melting point: 100°C EINECS: 233-151-8		GHS: Skin Irrit. 2, Eye Irrit. 2, STOT SE 3 H315, H319, H335 P261, P302+P352, P280, P305+P351+P338	
			
Assay: >98%			
Guaranteed analysis			
Calcium (Ca)	<50 ppm	Sodium (Na)	<200 ppm
Copper (Cu)	<50 ppm	n-Butyl formate	<0.005%
Iron (III)	<0.05%	Total phosphorus (P)	<100 ppm
Lead (Pb)	<50 ppm	Total silicon (Si)	<50 ppm
Magnesium (Mg)	<200 ppm	Zinc (Zn)	<200 ppm
Potassium (K)	<200 ppm		



Ammonium iron(III) sulfate dodecahydrate Certified AR
for analysis

Packaging		Mfr. No
500 g	Plastic jar	A/4760/53
1 kg	Plastic jar	A/4760/60
NH ₄ Fe(SO ₄) ₂ ·12H ₂ O		Boiling point: 230°C
CAS: 7783-83-7		Melting point: 39 to 41°C
MW: 482.19		EINECS: 233-382-4
Density at 20°C: 1.71		
Assay: >99%		
Guaranteed analysis		
pH (aqueous solution 5%)	<3	Potassium (K) <20 ppm
Calcium (Ca)	<10 ppm	Sodium (Na) <20 ppm
Copper (Cu)	<5 ppm	Total chloride (Cl) <0.0005%
Iron (II)	<0.001%	Total phosphorus (P) <20 ppm
Lead (Pb)	<10 ppm	Total silicon (Si) <10 ppm
Magnesium (Mg)	<5 ppm	Zinc (Zn) <20 ppm
Typical analysis		
Aluminium (Al)	<2 ppm	Manganese (Mn) <0.001 %
Silver (Ag)	<2 ppm	Molybdenum (Mo) <2 ppm
Arsenic (As)	<5 ppm	Nickel (Ni) <2 ppm
Barium (Ba)	<2 ppm	Strontium (Sr) <2 ppm
Chromium (Cr)	<2 ppm	Titanium (Ti) <2 ppm
Cobalt (Co)	<2 ppm	Vanadium (V) <2 ppm
Lithium (Li)	<2 ppm	
Lot analysis on the label.		

Ammonium iron(III) sulfate dodecahydrate SLR
specified reagents for general laboratory work- extra pure

Packaging		Mfr. No	
500 g	Plastic jar	A/4720/53	
1 kg	Plastic jar	A/4720/60	
NH ₄ Fe(SO ₄) ₂ ·12H ₂ O		Boiling point: 230°C	
CAS: 7783-83-7		Melting point: 39 to 41°C	
MW: 482.19		EINECS: 233-382-4	
Density at 20°C: 1.71			
Assay: >98%			
Guaranteed analysis			
Calcium (Ca)	<50 ppm	Sodium (Na)	<100 ppm
Copper (Cu)	<20 ppm	n-Butyl formate	<0.002%
Iron (II)	<0.005%	Total phosphorus (P)	<100 ppm
Lead (Pb)	<50 ppm	Total silicon (Si)	<50 ppm
Magnesium (Mg)	<20 ppm	Zinc (Zn)	<100 ppm
Potassium (K)	<100 ppm		

Ammonium iron (II) sulfate solution 0.1M (0.1N)
NIST Standard solution ready to use, For Volumetric analysis

Packaging		Mfr. No
1 ℓ	Glass bottle	J/1101/15
2.5 ℓ	HDPE Plastic bottle	J/1101/PB17
CAS: 7783-85-9		
Density at 20°C: 1.04		
Manufactured to a tolerance of ± 0.5%.		
See the full range of Standard Volumetric solution page 323.		

Ammonium iron(II) sulfate solution 0.025M
Standard solution ready to use for volumetric analysis

Packaging		Mfr. No
2.5 ℓ	Glass bottle	J/1100/17
CAS: 7783-85-9		
Density at 20°C: 1.04		
Manufactured to a tolerance of ± 0.5%.		
See the full range of Standard Volumetric solution page 323.		

Ammonium metavanadate SLR
specified reagents for general laboratory work- extra pure

Packaging		Mfr. No	
100 g	Plastic jar	A/5680/48	
250 g	Plastic jar	A/5680/50	
25 kg	Fibre Keg	A/5680/70	
NH ₄ VO ₃ CAS: 7803-55-6 MW: 116.98 Melting point: 200°C ADR/RID: 6.1, II EINECS: 232-261-3		UN: 2859 GHS: Acute Tox. 3, Skin Irrit. 2, Eye Irrit. 2, STOT SE 3 H315, H335, H319, H301, H332 P301+P310, P302+P352, P280, P305+P351+P338, P260	
			
Assay: >99%			
Guaranteed analysis			
Calcium (Ca)	<50 ppm	Sodium (Na)	<100 ppm
Copper (Cu)	<10 ppm	n-Butyl formate	<0.002%
Iron (Fe)	<20 ppm	Total phosphorus (P)	<50 ppm
Lead (Pb)	<20 ppm	Total silicon (Si)	<50 ppm
Magnesium (Mg)	<10 ppm	Total sulphur (S)	<50 ppm
Potassium (K)	<100 ppm	Zinc (Zn)	<10 ppm

Ammonium molybdate tetrahydrate Certified AR
for analysis

Packaging		Mfr. No	
100 g	Plastic jar	A/5760/48	
250 g	Plastic jar	A/5760/50	
500 g	Plastic jar	A/5760/53	
1 kg	Plastic jar	A/5760/60	
25 kg	Fibre Keg	A/5760/70	
(NH ₄) ₆ Mo ₇ O ₂₄ ·4H ₂ O CAS: 12054-85-2 MW: 123.86 Density at 20°C: 2.49 Melting point: 190°C EINECS: 234-722-4		GHS: Acute Tox. 4, Skin Irrit. 2, Eye Irrit. 2, STOT SE 3 H315, H302, H335, H319 P261, P301+P312, P302+P352, P280, P305+P351+P338	
Assay: 99 to 101%			
Guaranteed analysis			
Calcium (Ca)	<5 ppm	Sodium (Na)	<1000 ppm
Copper (Cu)	<10 ppm	n-Butyl formate	<0.0005%
Iron (Fe)	<5 ppm	Total sulphur (S)	<50 ppm
Lead (Pb)	<5 ppm	Zinc (Zn)	<5 ppm
Magnesium (Mg)	<5 ppm	Arsenate, phosphate and silicate	<0.001%
Potassium (K)	<20 ppm		
Lot analysis on the label.			

Ammonium molybdate tetrahydrate SLR
specified reagents for general laboratory work- extra pure

Packaging		Mfr. No	
100 g	Plastic jar	A/5720/48	
250 g	Plastic jar	A/5720/50	
500 g	Plastic jar	A/5720/53	
(NH₄)₆Mo₇O₂₄·4H₂O CAS: 12054-85-2 MW: 123.86 Density at 20°C: 2.49 Melting point: 190°C EINECS: 234-722-4		GHS: Acute Tox. 4, Skin Irrit. 2, Eye Irrit. 2, STOT SE 3 H315, H302, H335, H319 P261, P301+P312, P302+P352, P280, P305+P351+P338 	
Assay: >98%			
Guaranteed analysis			
Calcium (Ca)	<20 ppm	Sodium (Na)	<5000 ppm
Copper (Cu)	<50 ppm	n-Butyl formate	<0.002%
Iron (Fe)	<20 ppm	Total silicon (Si)	<200 ppm
Lead (Pb)	<20 ppm	Total sulphur (S)	<200 ppm
Magnesium (Mg)	<20 ppm	Zinc (Zn)	<20 ppm
Potassium (K)	<100 ppm		

Ammonium nickel(II) sulfate hexahydrate		SLR
specified reagents for general laboratory work- extra pure		
packaging	Mfr. No	
250 g Plastic jar	A/5800/50	
500 g Plastic jar	A/5800/53	
(NH ₄) ₂ Ni(SO ₄) ₂ ·6H ₂ O CAS: 15699-18-0 MW: 394.97 Density at 20°C: 1.9 ADR/RID: 9, III EINECS: 239-793-5 UN: 3077 GHS: Carc. 1B, Muta. 2, Repr. 1B, Acute Tox. 4,	Resp. Sens. 1, Skin Sens. 1, STOT RE 1, Aquatic Acute 1, Aquatic Chronic 1 H302, H317, H332, H334, H341, H350i, H372, H410, H360D P273, P304+P340, P301+P312, P307+P311, P201, P302+P352, P280	
Guaranteed analysis		
Identity test : Pass		

Ammonium nitrate		Certified AR
for analysis		
packaging	Mfr. No	
500 g Plastic jar	A/5920/53	
1 kg Plastic jar	A/5920/60	
2.5 kg Plastic jar	A/5920/62	
NH ₄ NO ₃ CAS: 6484-52-2 MW: 80.04 Density at 20°C: 1.72 Boiling point: 210°C Melting point: 169°C ADR/RID: 5.1, III	EINECS: 229-347-8 UN: 1942 GHS: Ox. Sol. 2, Skin Irrit. 2, Eye Irrit. 2, STOT SE 3 H272, H315, H319, H335 P261, P302+P352, P280, P305+P351+P338, P210	
Assay: >99.5%		
Guaranteed analysis		
pH (aqueous solution 5%) : 4.5 to 6 Calcium (Ca) : <5 ppm Copper (Cu) : <1 ppm Iron (Fe) : <1 ppm Lead (Pb) : <2 ppm Magnesium (Mg) : <1 ppm Nitrite (NO ₂) : <0.0005%		Potassium (K) : <10 ppm Sodium (Na) : <5 ppm Total chloride (Cl) : <0.0005% Total phosphorus (P) : <5 ppm Total silicon (Si) : <5 ppm Total sulphur (S) : <10 ppm Zinc (Zn) : <1 ppm
Typical analysis		
Aluminium (Al) : <1 ppm Silver (Ag) : <1 ppm Arsenic (As) : <2 ppm Barium (Ba) : <1 ppm Cadmium (Cd) : <1 ppm Chromium (Cr) : <1 ppm Cobalt (Co) : <1 ppm		Lithium (Li) : <1 ppm Manganese (Mn) : <1 ppm Molybdenum (Mo) : <1 ppm Nickel (Ni) : <1 ppm Strontium (Sr) : <1 ppm Titanium (Ti) : <1 ppm Vanadium (V) : <1 ppm
Lot analysis on the label.		

Ammonium nitrate		SLR
specified reagents for general laboratory work- extra pure		
packaging	Mfr. No	
500 g Plastic jar	A/5880/53	
1 kg Plastic jar	A/5880/60	
5 kg Plastic keg	A/5880/65	
NH ₄ NO ₃ CAS: 6484-52-2 MW: 80.04 Density at 20°C: 1.72 Boiling point: 210°C Melting point: 169°C ADR/RID: 5.1, III	EINECS: 229-347-8 UN: 1942 GHS: Ox. Sol. 2, Skin Irrit. 2, Eye Irrit. 2, STOT SE 3 H272, H315, H319, H335 P261, P302+P352, P280, P305+P351+P338, P210	
Assay: >98%		
Guaranteed analysis		
Calcium (Ca) : <50 ppm Copper (Cu) : <20 ppm Iron (Fe) : <10 ppm Lead (Pb) : <10 ppm Magnesium (Mg) : <20 ppm Potassium (K) : <100 ppm		Sodium (Na) : <50 ppm n-Butyl formate : <0.002% Total phosphorus (P) : <10 ppm Total silicon (Si) : <50 ppm Total sulphur (S) : <100 ppm Zinc (Zn) : <20 ppm

Ammonium oxalate monohydrate		Certified AR
for analysis		
packaging	Mfr. No	
100 g Plastic jar	A/6040/48	
500 g Plastic jar	A/6040/53	
1 kg Plastic jar	A/6040/60	
(NH ₄) ₂ C ₂ O ₄ ·H ₂ O CAS: 6009-70-7 MW: 142.11 Density at 20°C: 1.5 Melting point: 70°C	ADR/RID: 6.1, III UN: 2811 GHS: Acute Tox. 4 H312, H302 P280, P302+P352, P301+P312	
Assay: 99 to 101%		
Guaranteed analysis		
Reaction (2.5% solution) : 6.3 to 6.7 Calcium (Ca) : <5 ppm Copper (Cu) : <5 ppm Iron (Fe) : <5 ppm Lead (Pb) : <5 ppm Magnesium (Mg) : <5 ppm Potassium (K) : <10 ppm		Sodium (Na) : <10 ppm Total chloride (Cl) : <0.0005% Total phosphorus (P) : <20 ppm Total silicon (Si) : <10 ppm Total sulphur (S) : <50 ppm Zinc (Zn) : <5 ppm
Typical analysis		
Aluminium (Al) : <2 ppm Silver (Ag) : <2 ppm Arsenic (As) : <2 ppm Barium (Ba) : <2 ppm Cadmium (Cd) : <2 ppm Chromium (Cr) : <2 ppm Cobalt (Co) : <2 ppm		Lithium (Li) : <2 ppm Manganese (Mn) : <2 ppm Molybdenum (Mo) : <2 ppm Nickel (Ni) : <2 ppm Strontium (Sr) : <2 ppm Titanium (Ti) : <2 ppm Vanadium (V) : <2 ppm
Lot analysis on the label.		

Ammonium oxalate monohydrate		SLR
specified reagents for general laboratory work- extra pure		
packaging	Mfr. No	
500 g Plastic jar	A/6000/53	
(NH ₄) ₂ C ₂ O ₄ ·H ₂ O CAS: 6009-70-7 MW: 142.11 Density at 20°C: 1.5 Melting point: 70°C	ADR/RID: 6.1, III UN: 2811 GHS: Acute Tox. 4 H312, H302 P280, P302+P352, P301+P312	
Assay: >99%		
Guaranteed analysis		
Calcium (Ca) : <50 ppm Copper (Cu) : <50 ppm Iron (Fe) : <20 ppm Lead (Pb) : <50 ppm Magnesium (Mg) : <50 ppm Potassium (K) : <100 ppm		Sodium (Na) : <100 ppm n-Butyl formate : <0.005% Total phosphorus (P) : <200 ppm Total silicon (Si) : <100 ppm Total sulphur (S) : <500 ppm Zinc (Zn) : <50 ppm

Ammonium persulfate		
for electrophoresis		
packaging	Mfr. No	
25 g Plastic jar	A/P470/46	
H ₈ N ₂ O ₈ S ₂ CAS: 7727-54-0 MW: 228.19 Density at 20°C: 1.98 Melting point: 120°C ADR/RID: 5.1, III EINECS: 231-786-5 UN: 1444 GHS: Ox. Sol. 3, Acute Tox.	4, Skin Irrit. 2, Eye Irrit. 2, Resp. Sens. 1, Skin Sens. 1, STOT SE 3 H272, H334, H335, H315, H317, H319, H302 P342+P311, P280, P301+P312, P302+P352, P305+P351+P338, P210	
Assay: >99%		
Guaranteed analysis		
Insoluble matter : <0.003% Copper (Cu) : <20 ppm Iron (Fe) : <10 ppm		Lead (Pb) : <20 ppm Manganese (Mn) : <1 ppm n-Butyl formate : <0.0005%

Ammonium persulfate		Certified AR
for analysis		
packaging	Mfr. No	
500 g Plastic jar	A/6160/53	
1 kg Plastic jar	A/6160/60	
(NH ₄) ₂ S ₂ O ₈ CAS: 7727-54-0 MW: 228.19 Density at 20°C: 1.98 Melting point: 120°C ADR/RID: 5.1, III EINECS: 231-786-5 UN: 1444 GHS: Ox. Sol. 3, Acute Tox.	4, Skin Irrit. 2, Eye Irrit. 2, Resp. Sens. 1, Skin Sens. 1, STOT SE 3 H272, H334, H335, H315, H317, H319, H302 P342+P311, P280, P301+P312, P302+P352, P305+P351+P338, P210	
Assay: >99%		
Guaranteed analysis		
Calcium (Ca) : <5 ppm Copper (Cu) : <5 ppm Iron (Fe) : <5 ppm Lead (Pb) : <5 ppm Magnesium (Mg) : <5 ppm Potassium (K) : <20 ppm		Sodium (Na) : <100 ppm Total chloride (Cl) : <0.0005% Total phosphorus (P) : <50 ppm Total silicon (Si) : <5 ppm Zinc (Zn) : <5 ppm
Typical analysis		
Aluminium (Al) : <1 ppm Silver (Ag) : <1 ppm Arsenic (As) : <5 ppm Barium (Ba) : <1 ppm Cadmium (Cd) : <1 ppm Chromium (Cr) : <1 ppm Cobalt (Co) : <1 ppm		Lithium (Li) : <1 ppm Manganese (Mn) : <1 ppm Molybdenum (Mo) : <1 ppm Nickel (Ni) : <1 ppm Strontium (Sr) : <1 ppm Titanium (Ti) : <1 ppm Vanadium (V) : <1 ppm
Lot analysis on the label.		

Ammonium persulfate		SLR
specified reagents for general laboratory work- extra pure		
packaging	Mfr. No	
1 kg Plastic jar	A/6120/60	
3 kg Plastic jar	A/6120/63	
(NH ₄) ₂ S ₂ O ₈ CAS: 7727-54-0 MW: 228.19 Density at 20°C: 1.98 Melting point: 120°C ADR/RID: 5.1, III EINECS: 231-786-5 UN: 1444 GHS: Ox. Sol. 3, Acute Tox.	4, Skin Irrit. 2, Eye Irrit. 2, Resp. Sens. 1, Skin Sens. 1, STOT SE 3 H272, H334, H335, H315, H317, H319, H302 P342+P311, P280, P301+P312, P302+P352, P305+P351+P338, P210	
Assay: >98%		
Guaranteed analysis		
Calcium (Ca) : <20 ppm Copper (Cu) : <50 ppm Iron (Fe) : <50 ppm Lead (Pb) : <50 ppm Magnesium (Mg) : <20 ppm Potassium (K) : <200 ppm		Sodium (Na) : <500 ppm n-Butyl formate : <0.002% Total phosphorus (P) : <500 ppm Total silicon (Si) : <50 ppm Zinc (Zn) : <50 ppm

Ammonium polysulfide 8% solution		Certified AR
for analysis		
packaging	Mfr. No	
500 mL Coated Glass bottle	A/6560/PB08	
2.5 L Coated Glass bottle	A/6560/PB17	
(NH ₄) ₂ S _x CAS: 12259-92-6 Density at 20°C: 1.07 ADR/RID: 8, 6.1, II EINECS: 235-512-5	UN: 2818 GHS: Skin Corr. 1B, Eye Dam. 1 H314, H400, EUH031 P280, P305+P351+P338, P273	
Assay: >6%		
Guaranteed analysis		
Residue after ignition : <0.01%		Carbonate (Na ₂ CO ₃) : <0.005%
Lot analysis on the label.		

Ammonium purpurate		
pure, Metal indicator		
packaging	Mfr. No	
5 g Glass jar	A/6200/44	
25 g Glass jar	A/6200/46	
C ₈ H ₈ N ₄ O ₆ CAS: 3051-09-0 MW: 284.19	Melting point: 300°C EINECS: 221-266-6	
Guaranteed analysis		
Complies with standard IR spectrum: Pass		

Ammonium pyrrolidine dithiocarbamate		
for AAS		
packaging	Mfr. No	
25 g Plastic jar	A/6220/46	
(C ₄ H ₈ N) ₂ CSSNH ₄ CAS: 5108-96-3 MW: 164.29 Melting point: 150 to 155°C ADR/RID: 6.1, III EINECS: 225-834-4 UN: 2811	GHS: Acute Tox. 4, Skin Irrit. 2, Eye Irrit. 2, STOT SE 3 H315, H319, H302, H335 P261, P301+P312, P302+P352, P280, P305+P351+P338	
Guaranteed analysis		
Residue after ignition : <0.05% Cobalt (Co) : <0.2 ppm Copper (Cu) : <0.2 ppm Iron (Fe) : <0.5 ppm		Lead (Pb) : <0.5 ppm Nickel (Ni) : <0.2 ppm Zinc (Zn) : <0.5 ppm
For the analysis by atomic absorption (AAS). Contains ammonium carbonate to delay decomposition.		

Ammonium sodium hydrogen orthophosphate tetrahydrate		SLR
specified reagents for general laboratory work- extra pure		
packaging	Mfr. No	
500 g Plastic jar	A/6280/53	
Na(NH ₄)HPO ₄ ·4H ₂ O CAS: 13011-54-6 MW: 209.07	Density at 20°C: 1.54 Melting point: 79°C EINECS: 235-860-8	
Assay: >95%		
Guaranteed analysis		
Calcium (Ca) : <50 ppm Copper (Cu) : <50 ppm Iron (Fe) : <50 ppm Lead (Pb) : <50 ppm		Magnesium (Mg) : <100 ppm Potassium (K) : <100 ppm Zinc (Zn) : <50 ppm

Ammonium sulfamate		
pure		
packaging	Mfr. No	
250 g Plastic jar	A/6400/50	
NH ₄ OSO ₂ NH ₂ CAS: 7773-06-0 MW: 114.12 Melting point: 132 to 135°C	EINECS: 231-871-7 GHS: Acute Tox. 4 H302 P280, P264, P270, P301+P312	
Assay: >99%		

Reagents, Solvents & Volumetric solutions comply with the European Pharmacopeia

Read more on page 234

Ammonium sulfate
for analysis

packaging	Mfr. No
500 g Plastic jar	A/6480/53
1 kg Plastic jar	A/6480/60
2.5 kg Plastic jar	A/6480/62
5 kg Plastic keg	A/6480/65
25 kg Fibre Keg	A/6480/70

(NH ₄) ₂ SO ₄ CAS: 7783-20-2 MW: 132.13	Melting point: 280°C EINECS: 231-984-1
---	---

Assay: >99.5%

Guaranteed analysis		
pH (aqueous solution 5%)	5 to 6	Nitrate (NO ₃)
pH(aq. Solution 5%) at 25°C	: 5 to 6	Potassium (K)
Loss on drying	<0.1%	Sodium (Na)
Heavy metals	<5 ppm	Total chloride (Cl)
Insoluble matter	<0.005%	Total phosphorus (P)
Calcium (Ca)	<5 ppm	Total silicon (Si)
Copper (Cu)	<1 ppm	Zinc (Zn)
Iron (Fe)	<1 ppm	Phosphate (PO ₄)
Lead (Pb)	<2 ppm	Residue after ignition
Magnesium (Mg)	<2 ppm	

Typical analysis			
Aluminium (Al)	<2 ppm	Lithium (Li)	<2 ppm
Silver (Ag)	<2 ppm	Manganese (Mn)	<2 ppm
Arsenic (As)	<2 ppm	Molybdenum (Mo)	<2 ppm
Barium (Ba)	<2 ppm	Nickel (Ni)	<2 ppm
Cadmium (Cd)	<2 ppm	Strontium (Sr)	<2 ppm
Chromium (Cr)	<2 ppm	Titanium (Ti)	<2 ppm
Cobalt (Co)	<2 ppm	Vanadium (V)	<2 ppm

Lot analysis on the label.

Ammonium sulfate
specified reagents for general laboratory work- extra pure

packaging	Mfr. No
500 g Plastic jar	A/6440/53
1 kg Plastic jar	A/6440/60
2.5 kg Plastic jar	A/6440/62
5 kg Plastic keg	A/6440/65

(NH ₄) ₂ SO ₄ CAS: 7783-20-2 MW: 132.13	Melting point: 280°C EINECS: 231-984-1
---	---

Assay: >99%

Guaranteed analysis			
pH (aqueous solution 5%)	5 to 6	Potassium (K)	<200 ppm
Calcium (Ca)	<50 ppm	Sodium (Na)	<200 ppm
Copper (Cu)	<5 ppm	n-Butyl formate	<0.002%
Iron (Fe)	<10 ppm	Total phosphorus (P)	<50 ppm
Lead (Pb)	<20 ppm	Total silicon (Si)	<100 ppm
Magnesium (Mg)	<20 ppm	Zinc (Zn)	<100 ppm
Nitrate (NO ₃)	<0.01%		

Ammonium thiocyanate
for analysis

packaging	Mfr. No
250 g Plastic jar	A/6680/50
500 g Plastic jar	A/6680/53
1 kg Plastic jar	A/6680/60

NH ₄ SCN CAS: 1762-95-4 MW: 76.12 Density at 20°C: 1.3 Melting point: 149°C Flash point: 190°C EINECS: 217-175-6	GHS: Acute Tox. 4, Aquatic Chronic 3 H302, H312, H412, H332, EUH032 P280, P302+P352, P273, P233, P301+P312, P304+P340
---	--

Assay: >99%

Guaranteed analysis			
pH (aqueous solution 5%)	4.8 to 5.8	Potassium (K)	<10 ppm
Iodine consuming substances	: <0.02%	Sodium (Na)	<10 ppm
Calcium (Ca)	<20 ppm	Sulfate (SO ₄)	<0.03%
Copper (Cu)	<2 ppm	Total chloride (Cl)	<0.005%
Iron (Fe)	<1 ppm	Total phosphorus (P)	<50 ppm
Lead (Pb)	<5 ppm	Total silicon (Si)	<10 ppm
Magnesium (Mg)	<5 ppm	Zinc (Zn)	<5 ppm

Typical analysis			
Sulfide (S)	<0.001%	Lithium (Li)	<1 ppm
Aluminium (Al)	<1 ppm	Manganese (Mn)	<1 ppm
Silver (Ag)	<1 ppm	Molybdenum (Mo)	<1 ppm
Arsenic (As)	<2 ppm	Nickel (Ni)	<1 ppm
Barium (Ba)	<1 ppm	Strontium (Sr)	<1 ppm
Cadmium (Cd)	<1 ppm	Titanium (Ti)	<1 ppm
Chromium (Cr)	<1 ppm	Vanadium (V)	<1 ppm
Cobalt (Co)	<1 ppm		

Lot analysis on the label.

Ammonium thiocyanate
specified reagents for general laboratory work- extra pure

packaging	Mfr. No
250 g Plastic jar	A/6640/50
500 g Plastic jar	A/6640/53
1 kg Plastic jar	A/6640/60

NH ₄ SCN CAS: 1762-95-4 MW: 76.12 Density at 20°C: 1.3 Melting point: 149°C Flash point: 190°C EINECS: 217-175-6	GHS: Acute Tox. 4, Aquatic Chronic 3 H302, H312, H412, H332, EUH032 P280, P302+P352, P273, P233, P301+P312, P304+P340
---	--

Assay: >98%

Guaranteed analysis			
Calcium (Ca)	<100 ppm	Sodium (Na)	<200 ppm
Copper (Cu)	<10 ppm	n-Butyl formate	<0.02%
Iron (Fe)	<5 ppm	Total phosphorus (P)	<200 ppm
Lead (Pb)	<20 ppm	Total silicon (Si)	<50 ppm
Magnesium (Mg)	<20 ppm	Zinc (Zn)	<20 ppm
Potassium (K)	<50 ppm		

Ammonium thiocyanate solution 0.1M (0.1N)
NIST Standard solution ready to use, For Volumetric analysis

packaging	Mfr. No
1 ℓ HDPE Plastic bottle	J/1500/PB15
2.5 ℓ HDPE Plastic bottle	J/1500/PB17

Manufactured to a tolerance of ± 0.2%.
See the full range of Standard Volumetric solution page 323.

Ammonium thiocyanate solution 0.1M (0.1N)
NIST Standard Concentrate, For Volumetric analysis

packaging	Mfr. No
1 Ampoule	J/1500C/05
6 Ampoule	J/1500C/90

CAS: 1762-95-4 EINECS: 217-175-6 EUH032 P233

Each ampoule makes 1 ℓ of solution.
See the full range of Standard Volumetric solution page 323.



Ammonium thiosulfate
specified reagents for general laboratory work- extra pure

packaging	Mfr. No
500 g Plastic jar	A/6720/53
25 kg Fibre Keg	A/6720/70

(NH ₄) ₂ S ₂ O ₃ CAS: 7783-18-8 MW: 148.20	Melting point: 150°C EINECS: 231-982-0
---	---

Guaranteed analysis
Identity test : Pass

n-Amyl acetate
extra pure, mix of 3 and 2-Methylbutyl acetate

packaging	Mfr. No
500 mℓ Glass bottle	A/6800/08

CH ₃ COO(CH ₂) ₄ CH ₃ CAS: 628-63-7 MW: 130.19 Density at 20°C: 0.870 Boiling point: 149°C Melting point: -70.8°C Flash point: 24°C	ADR/RID: 3, III EINECS: 211-047-3 UN: 1104 GHS: Flam. Liq. 3 H226, EUH066 P210
--	---

Assay (CG): >99%

iso-Amyl acetate
for analysis, Mix of 2- &3- Methylbutyl acetate

packaging	Mfr. No
500 mℓ Glass bottle	A/6920/08
2.5 ℓ Glass bottle	A/6920/17

CH ₃ COOCH ₂ CH ₂ CH(CH ₃) ₂ CAS: 123-92-2 MW: 130.19 Density at 20°C: 0.874 Boiling point: 142°C Melting point: -78°C Flash point: 25°C	ADR/RID: 3, III EINECS: 204-662-3 UN: 1104 GHS: Flam. Liq. 3 H226, EUH066 P210
--	---

Assay (CG): >99%

Guaranteed analysis		
Acidity	<0.002 meq/g	Lead (Pb)
Colour	<10 APHA	Magnesium (Mg)
Water	<0.1%	Potassium (K)
Residue after evaporation	<20 ppm	Sodium (Na)
Calcium (Ca)	<1 ppm	Total phosphorus (P)
Copper (Cu)	<0.05 ppm	Total sulphur (S)
Iron (Fe)	<0.2 ppm	Zinc (Zn)

Lot analysis on the label.

iso-Amyl acetate
extra pure, Mix of 2- &3- Methylbutyl acetate

packaging	Mfr. No
500 mℓ Glass bottle	A/6880/08
2.5 ℓ Glass bottle	A/6880/17

CH ₃ COOCH ₂ CH ₂ CH(CH ₃) ₃ CAS: 123-92-2 MW: 130.19 Density at 20°C: 0.874 Boiling point: 142°C Melting point: -78°C Flash point: 25°C	ADR/RID: 3, III EINECS: 204-662-3 UN: 1104 GHS: Flam. Liq. 3 H226, EUH066 P210
--	---

Assay (CG): >98%

Guaranteed analysis		
Acidity	<0.01 meq/g	Lead (Pb)
Colour	<20 APHA	Magnesium (Mg)
Water	<0.5%	Potassium (K)
Residue after evaporation	<200 ppm	Sodium (Na)
Calcium (Ca)	<5 ppm	Total sulphur (S)
Copper (Cu)	<1 ppm	Zinc (Zn)
Iron (Fe)	<2 ppm	

n-Amyl alcohol See Pentan-1-ol, page 226

iso-Amyl alcohol
specified reagents for general laboratory work- extra pure

packaging	Mfr. No
500 mℓ Glass bottle	A/6960/08
2.5 ℓ Coated Glass bottle	A/6960/PB17

(CH ₃) ₂ CHCH ₂ CH ₂ OH CAS: 123-51-3 MW: 88.15 Density at 20°C: 0.800 Boiling point: 130°C Melting point: -117°C Flash point: 43°C	ADR/RID: 3, III EINECS: 204-633-5 UN: 1105 GHS: Flam. Liq. 3, Acute Tox. 4, STOT SE 3 H226, H335, H332, EUH066 P261, P304+P340, P210
--	--

Assay (CG): >98%

Guaranteed analysis		
Acidity / Alkalinity	<0.01 meq/g	Lead (Pb)
Colour	<20 APHA	Magnesium (Mg)
Water	<0.5%	Potassium (K)
Residue after evaporation	<100 ppm	Sodium (Na)
Calcium (Ca)	<1 ppm	Total phosphorus (P)
Copper (Cu)	<1 ppm	Total sulphur (S)
Iron (Fe)	<2 ppm	Zinc (Zn)


iso-Amyl alcohol
extra pure, for determination of fat in milk

packaging	Mfr. No
500 mℓ Glass bottle	A/7000/08
2.5 ℓ Coated Glass bottle	A/7000/PB17

(CH ₃) ₂ CHCH ₂ CH ₂ OH CAS: 123-51-3 MW: 88.15 Density at 20°C: 0.800 Boiling point: 130°C Melting point: -117°C Flash point: 43°C	ADR/RID: 3, III EINECS: 204-633-5 UN: 1105 GHS: Flam. Liq. 3, Acute Tox. 4, STOT SE 3 H226, H335, H332, EUH066 P261, P304+P340, P210
--	--

Assay (CG): >98%

Guaranteed analysis		
Acidity / Alkalinity	<0.01 meq/g	Lead (Pb)
Colour	<20 APHA	Magnesium (Mg)
Water	<0.5%	Potassium (K)
Residue after evaporation	<100 ppm	Sodium (Na)
Calcium (Ca)	<1 ppm	Total phosphorus (P)
Copper (Cu)	<1 ppm	Total sulphur (S)
Iron (Fe)	<2 ppm	Zinc (Zn)



Aniline <i>for analysis</i>		Certified AR
packaging	Mfr. No	
500 mℓ Coated Glass bottle	A/7280/PB08	
2.5 ℓ Coated Glass bottle	A/7280/PB17	
C ₆ H ₇ N CAS: 62-53-3 MW: 93.13 Density at 20°C: 1.021 Boiling point: 181 to 185°C Melting point: -6.2°C Flash point: 76°C ADR/RID: 6.1, II EINECS: 200-539-3 UN: 1547	GHS: Carc. 2, Muta. 2, Acute Tox. 3, Eye Dam. 1, Skin Sens. 1, STOT RE 1, Aquatic Acute 1 H351, H372, H341, H311, H301, H318, H331, H317, H400 P301+P310, P280, P302+P350, P305+P351+P338, P304+P340, P260, P273	   

Assay (CG): >99.8%		
Guaranteed analysis		
Water	<0.1%	Potassium (K)
Residue after evaporation	<100 ppm	Sodium (Na)
Calcium (Ca)	<1 ppm	Total phosphorus (P)
Copper (Cu)	<0.5 ppm	Total sulphur (S)
Iron (Fe)	<0.5 ppm	Zinc (Zn)
Lead (Pb)	<0.1 ppm	Nitrobenzene (C ₆ H ₅ NO ₂)
Magnesium (Mg)	<0.5 ppm	

Aniline <i>specified reagents for general laboratory work- extra pure</i>		SLR
packaging	Mfr. No	
500 mℓ Coated Glass bottle	A/7240/PB08	
2.5 ℓ Coated Glass bottle	A/7240/PB17	
C ₆ H ₇ N CAS: 62-53-3 MW: 93.13 Density at 20°C: 1.021 Boiling point: 181 to 185°C Melting point: -6.2°C Flash point: 76°C ADR/RID: 6.1, II EINECS: 200-539-3 UN: 1547	GHS: Carc. 2, Muta. 2, Acute Tox. 3, Eye Dam. 1, Skin Sens. 1, STOT RE 1, Aquatic Acute 1 H351, H372, H341, H311, H301, H318, H331, H317, H400 P301+P310, P280, P302+P350, P305+P351+P338, P304+P340, P260, P273	   

Assay (CG): >99%		
Guaranteed analysis		
Water	<0.5%	Magnesium (Mg)
Residue after evaporation	<500 ppm	Potassium (K)
Calcium (Ca)	<5 ppm	Sodium (Na)
Copper (Cu)	<2 ppm	Total phosphorus (P)
Iron (Fe)	<2 ppm	Total sulphur (S)
Lead (Pb)	<1 ppm	Zinc (Zn)

Aniline blue C.I. 42755 <i>pure, water soluble</i>	
packaging	Mfr. No
25 g Glass jar	A/7300/46
C ₁₂ H ₂₅ N ₃ Na ₃ O ₉ S ₃ CAS: 28631-66-5 MW: 737.71 EINECS: 249-113-9 GHS: Skin Irrit. 2, Eye Irrit. 2,	
Guaranteed analysis	
Complies with standard IR spectrum: Pass	

Aniline hydrochloride <i>specified reagents for general laboratory work- extra pure</i>		SLR
packaging	Mfr. No	
100 g Glass jar	A/7360/48	
500 g Glass jar	A/7360/53	
C ₆ H ₇ N.HCl CAS: 142-04-1 MW: 129.59 Boiling point: 245°C Melting point: 196 to 199°C Flash point: 193°C ADR/RID: 6.1, III EINECS: 205-519-8 UN: 1548	GHS: Carc. 2, Muta. 2, Acute Tox. 3, Eye Dam. 1, Skin Sens. 1, STOT RE 1, Aquatic Acute 1 H351, H341, H331, H311, H301, H318, H317, H400 P280, P301+P310, P305+P351+P338, P302+P350, P304+P340, P310	   

Assay: >99%	
Guaranteed analysis	
Residue after ignition	<0.05%

8-Anilidonaphthalene-1-sulfonic acid ammonium salt <i>pure</i>	
packaging	Mfr. No
250 g Plastic jar	A/7486/50
C ₆ H ₄ NHC ₁₀ H ₆ SO ₃ NH ₄ CAS: 28836-03-5 MW: 316.38 EINECS: 249-265-6	 GHS: Skin Irrit. 2, Eye Irrit. 2 H319, H315 P280, P302+P352, P305+P351+P338
Guaranteed analysis	
Complies with standard IR spectrum: Pass	

Anion Eluent Concentrate <i>for Ion pair chromatography</i>	
packaging	Mfr. No
100 mℓ HDPE Plastic bottle	J/4534/05
0.0028M sodium bicarbonate / 0.0022M sodium carbonate solution. 100mℓ makes 5 ℓ of above solution.	

Anion standard mixed concentrate <i>for Ion pair chromatography (for column check)</i>	
packaging	Mfr. No
100 mℓ HDPE Plastic bottle	J/4552/05
Guaranteed analysis	
Fluoride (F)	29.7 to 30.3 ppm Sulfate (SO ₄)
Nitrate (NO ₃)	297 to 303 ppm
Manufactured to a tolerance of ± 1%. 100mℓ, to make 1000mℓ consisting of: 3ppm F, 6ppm Cl, 25ppm Br, 30ppm NO ₃ ⁻ , 60ppm PO ₄ ³⁻ , 30ppm SO ₄ ²⁻	

p-Anisaldehyde <i>specified reagents for general laboratory work- extra pure</i>		SLR
packaging	Mfr. No	
100 mℓ Coated Glass bottle	A/7490/PB05	
C ₈ H ₈ O ₂ CAS: 123-11-5 MW: 136.15 Density at 20°C: 1.121 Boiling point: 248 to 249°C Melting point: -1°C	Flash point: 116°C EINECS: 204-602-6 GHS: Acute Tox. 4 H302 P301+P330+P331	
Assay (CG): >98%		
Guaranteed analysis		
Acidity		

Anisole <i>pure</i>	
packaging	Mfr. No
500 mℓ Coated Glass bottle	A/7600/PB08
C ₇ H ₈ O CAS: 100-66-3 MW: 108.14 Density at 20°C: 0.990 Boiling point: 154°C Melting point: -37°C Flash point: 43°C	ADR/RID: 3, III EINECS: 202-876-1 UN: 2222 GHS: Flam. Liq. 3, Acute Tox. 4, Skin Irrit. 2 H226, H332, H315 P304+P340, P302+P352, P210
Assay (CG): >99%	
Guaranteed analysis	
Complies with standard IR spectrum: Pass	

Anthranilic acid <i>specified reagents for general laboratory work- extra pure</i>		SLR
packaging	Mfr. No	
100 g Glass jar	A/7680/48	
C ₇ H ₇ NO ₂ CAS: 118-92-3 MW: 137.14 Density at 20°C: 1.41 Melting point: 148°C Flash point: 150°C	EINECS: 204-287-5 GHS: Skin Irrit. 2, Eye Irrit. 2, STOT SE 3 H315, H335, H319 P261, P302+P352, P280, P305+P351+P338	
Assay: >98%		

Anthrone <i>specified reagents for general laboratory work- extra pure</i>		SLR
packaging	Mfr. No	
25 g Plastic jar	A/7760/46	
C ₁₄ H ₁₀ O CAS: 90-44-8 MW: 194.23 Boiling point: 721°C Melting point: 152 to 158°C	EINECS: 201-994-0 GHS: Skin Irrit. 2, Eye Irrit. 2, STOT SE 3 H315, H319, H335 P280, P305+P351+P338, P261	
Guaranteed analysis		
Residue after ignition	<0.2%	Sensitivity to sucrose
Pass		

Anti-bumping granules <i>extra pure, fused alumina</i>		SLR
packaging	Mfr. No	
250 g Plastic jar	A/7770/50	
500 g Plastic jar	A/7770/53	
1 kg Plastic jar	A/7770/60	
Al ₂ O ₃ CAS: 1344-28-1 MW: 101.96	Boiling point: 2980°C Melting point: 2000°C EINECS: 215-691-6	
Guaranteed analysis		
Identity test		

Antimony standard metal solution 1000 ppm <i>for AAS, solution HCl 1M</i>	
packaging	Mfr. No
100 mℓ HDPE Plastic bottle	J/8002/05
500 mℓ HDPE Plastic bottle	J/8002/08
ADR/RID: 8, III UN: 1789 GHS: Met. Corr. 1	H290 P234, P390
	
Manufactured to a tolerance of ± 0.5%. See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.	

Antimony standard metal solution 1000 ppm <i>for ICP, solution HCl 1M</i>	
packaging	Mfr. No
500 mℓ HDPE Plastic bottle	J/8215/08
ADR/RID: 8, III UN: 1789 GHS: Met. Corr. 1	H290 P234, P390
	

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
All ICP standards are manufactured from acids, metal & metal salts Ultra pure quality and contain less than <10 ppm total metal impurities.
See also the full range of "Standard solution" in the Standard solution section.

Antimony standard metal solution 10 000 ppm <i>for ICP, solution HCl 1M</i>	
packaging	Mfr. No
100 mℓ HDPE Plastic bottle	J/8210/05
ADR/RID: 8, III UN: 1789 GHS: Met. Corr. 1	H290 P234, P390
	

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.

Antimony standard metal solution 1000 ppm <i>for ICP-MS solution in H₂O/0.6% Tartaric Acid/HNO₃ trace, NIST certified</i>	
packaging	Mfr. No
100 mℓ Plastic bottle	SPCLS872Y
Manufactured to a tolerance of ± 0.5%. All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities. See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314. Check with your local distributor for availability of Spex product	

Antimony standard metal solution 1000 ppm <i>for ICP solution in H₂O/0.6% Tartaric Acid/HNO₃ trace, NIST certified</i>	
packaging	Mfr. No
100 mℓ Plastic bottle	SPPLSB72Y
250 mℓ Plastic bottle	SPPLSB72T
500 mℓ Plastic bottle	SPPLSB72X
Manufactured to a tolerance of ± 0.5%. All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities. See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314. Check with your local distributor for availability of Spex product	

Antimony standard metal solution 10 000 ppm <i>for ICP solution in H₂O/0.6% Tartaric Acid/HNO₃ trace, NIST certified</i>	
packaging	Mfr. No
100 mℓ Plastic bottle	SPPLSB73Y
500 mℓ Plastic bottle	SPPLSB73X
Manufactured to a tolerance of ± 0.5%. All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities. See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314. Check with your local distributor for availability of Spex product	

Antimony standard metal solution 10 000 ppm
for ICP solution in HCL 20%, NIST certified

packaging	Mfr. No
500 mℓ Plastic bottle	PLSB53X
H315, H319, H335, H351	 

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Antimony metal min.99+ %
pure

packaging	Mfr. No
500 g Plastic jar	A/7840/53
Sb	 
CAS: 7440-36-0	EINECS: 231-146-5
MW: 121.75	UN: 2871
Boiling point: 1635°C	GHS: Acute Tox. 4, Aquatic Chronic 2
Melting point: 630°C	H302, H332, H411
ADR/RID: 6.1, III	P273, P301+P312, P304+P340
Assay: >99%	

Antimony trichloride
for analysis

packaging	Mfr. No
250 g Plastic jar	A/8202/50
SbCl ₃	 
CAS: 10025-91-9	UN: 1733
MW: 228.11	GHS: Skin Corr. 1B, Aquatic Chronic 2
Boiling point: 223°C	H314, H411
Melting point: 73°C	P301+P330+P331, P280,
ADR/RID: 8, II	P305+P351+P338, P310,
EINECS: 233-047-2	P304+P340, P273
Assay: >99%	
Guaranteed analysis	
Insoluble in chloroform	: <0.05%
Calcium (Ca)	<100 ppm
Copper (Cu)	<5 ppm
Iron (Fe)	<10 ppm
Lead (Pb)	<20 ppm
Magnesium (Mg)	<10 ppm
Potassium (K)	<50 ppm
Sodium (Na)	<100 ppm
Total phosphorus (P)	<200 ppm
Total silicon (Si)	<20 ppm
Total sulphur (S)	<50 ppm
Zinc (Zn)	<50 ppm

Lot analysis on the label.

Antimony trichloride
specified reagents for general laboratory work- extra pure

packaging	Mfr. No
250 g Plastic jar	A/8200/50
SbCl ₃	 
CAS: 10025-91-9	UN: 1733
MW: 228.11	GHS: Skin Corr. 1B, Aquatic Chronic 2
Boiling point: 223°C	H314, H411
Melting point: 73°C	P301+P330+P331, P280,
ADR/RID: 8, II	P305+P351+P338, P310,
EINECS: 233-047-2	P304+P340, P273
Assay: >99%	

Guaranteed analysis	
Calcium (Ca)	<100 ppm
Copper (Cu)	<50 ppm
Iron (Fe)	<50 ppm
Lead (Pb)	<200 ppm
Magnesium (Mg)	<100 ppm
Potassium (K)	<200 ppm
Sodium (Na)	<500 ppm
Total silicon (Si)	<200 ppm
Total sulphur (S)	<200 ppm
Zinc (Zn)	<200 ppm

Antimony trioxide
specified reagents for general laboratory work- extra pure

packaging	Mfr. No
500 g Plastic jar	A/8240/53
Sb ₂ O ₃	
CAS: 1309-64-4	EINECS: 215-175-0
MW: 291.42	GHS: Carc. 2
Boiling point: 1550°C	H351
Melting point: 656°C	P281, P273, P260
Assay: >99%	

L(+)Arabinose
specified reagents for general laboratory work- extra pure

packaging	Mfr. No
25 g Plastic jar	A/8280/46
C ₅ H ₁₀ O ₅	Melting point: 155°C
CAS: 5328-37-0	EINECS: 226-214-6
MW: 150.13	
Guaranteed analysis	
Complies with standard IR spectrum: Pass	

L-Arginine
specified reagents for general laboratory work- extra pure

packaging	Mfr. No
25 g Plastic jar	A/8300/46
C ₆ H ₁₄ N ₄ O ₂	
CAS: 74-79-3	EINECS: 200-811-1
MW: 174.20	GHS: Eye Irrit. 2
Density at 20°C: 1.3	H319
Melting point: 223 to 224°C	P280, P305+P351+P338
Assay: >99%	
Guaranteed analysis	
Loss on drying	<0.3%
Specific rotation	26.5° to 27.5°
Residue after ignition	<0.1%
Ammonium (NH ₄)	<0.02%
Iron (Fe)	<50 ppm
Lead (Pb)	<20 ppm

Arsenic standard metal solution 1000 ppm
for AAS, solution HCl 1M

packaging	Mfr. No
100 mℓ HDPE Plastic bottle	J/8003/05
500 mℓ HDPE Plastic bottle	J/8003/08
ADR/RID: 8, III	H350
EINECS: 215-481-4	P281, P308+P313, P302+P352
UN: 1789	

Manufactured to a tolerance of ± 0.5%.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.

Arsenic standard metal solution 1000 ppm
for ICP solution in HNO₃ 2%, NIST certified

packaging	Mfr. No
100 mℓ Plastic bottle	SPPLAS22Y
250 mℓ Plastic bottle	SPPLAS22T
500 mℓ Plastic bottle	SPPLAS22X
H315, H319	

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Arsenic standard metal solution 1000 ppm
for ICP solution in HCL 2%, NIST certified

packaging	Mfr. No
500 mℓ Plastic bottle	SPPLAS12X
H350	

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Arsenic standard metal solution 10 000 ppm
for ICP solution in HNO₃ 5%, NIST certified

packaging	Mfr. No
100 mℓ Plastic bottle	SPPLAS23Y
500 mℓ Plastic bottle	SPPLAS23X
H314, H412	

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Arsenic standard metal solution 1000 ppm
for ICP-MS solution in HNO₃ 2%, NIST certified

packaging	Mfr. No
100 mℓ Plastic bottle	SPCLAS22Y
H315, H319	

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Arsenic trioxide
specified reagents for general laboratory work- extra pure

packaging	Mfr. No
500 g Plastic jar	A/8440/53
As ₂ O ₃	   
CAS: 1327-53-3	GHS: Acute Tox. 2, Skin Corr. 1B, Eye Dam. 1, Aquatic Acute 1, Aquatic Chronic 1
MW: 197.84	H314, H410, H350, H300
Melting point: 312.3°C	P280, P301+P330+P331,
ADR/RID: 6.1, II	P305+P351+P338, P201,
EINECS: 215-481-4	P308+P313, P273
UN: 1561	

Assay: >99%	
Guaranteed analysis	
Calcium (Ca)	<200 ppm
Copper (Cu)	<100 ppm
Iron (Fe)	<50 ppm
Lead (Pb)	<200 ppm
Magnesium (Mg)	<200 ppm
Potassium (K)	<500 ppm
Sodium (Na)	<500 ppm
n-Butyl formate	<0.05%
Total phosphorus (P)	<500 ppm
Total silicon (Si)	<200 ppm
Total sulphur (S)	<500 ppm
Zinc (Zn)	<200 ppm

L-Ascorbic acid
for analysis

packaging	Mfr. No
100 g Plastic jar	A/8882/48
500 g Plastic jar	A/8882/53
C ₆ H ₈ O ₆	Melting point: 190 to 192°C
CAS: 50-81-7	EINECS: 200-066-2
MW: 176.12	

Assay: >99.7%	
Guaranteed analysis	
Specific rotation	20.5° to 21.5°
Calcium (Ca)	<10 ppm
Copper (Cu)	<2 ppm
Iron (Fe)	<2 ppm
Lead (Pb)	<2 ppm
Magnesium (Mg)	<2 ppm
Potassium (K)	<5 ppm
Sodium (Na)	<5 ppm
Total chloride (Cl)	<0.005%
Total phosphorus (P)	<20 ppm
Total silicon (Si)	<20 ppm
Total sulphur (S)	<20 ppm
Zinc (Zn)	<2 ppm

Typical analysis	
Aluminium (Al)	<1 ppm
Silver (Ag)	<1 ppm
Arsenic (As)	<2 ppm
Barium (Ba)	<1 ppm
Cadmium (Cd)	<1 ppm
Chromium (Cr)	<1 ppm
Cobalt (Co)	<1 ppm
Lithium (Li)	<1 ppm
Manganese (Mn)	<1 ppm
Molybdenum (Mo)	<1 ppm
Nickel (Ni)	<1 ppm
Strontium (Sr)	<1 ppm
Titanium (Ti)	<1 ppm
Vanadium (V)	<1 ppm

Lot analysis on the label.

L-Ascorbic acid
extra pure, 99+%

packaging	Mfr. No
100 g Plastic jar	A/8880/48
500 g Plastic jar	A/8880/53
1 kg Plastic jar	A/8880/60
C ₆ H ₈ O ₆	Density at 20°C: 1.7
CAS: 50-81-7	Melting point: 190 to 192°C
MW: 176.13	EINECS: 200-066-2

Assay: >99%	
Guaranteed analysis	
Specific rotation	20° to 22°
Calcium (Ca)	<20 ppm
Copper (Cu)	<20 ppm
Iron (Fe)	<20 ppm
Lead (Pb)	<20 ppm
Magnesium (Mg)	<20 ppm
Potassium (K)	<50 ppm
Sodium (Na)	<50 ppm
n-Butyl formate	<0.02%
Zinc (Zn)	<20 ppm

Vitamin C

L-Aspartic acid
pure

packaging	Mfr. No
500 g Plastic jar	A/8961/53
C ₄ H ₇ NO ₄	Melting point: >230°C
CAS: 56-84-8	EINECS: 200-291-6
MW: 133.11	

Assay: >99%	
Guaranteed analysis	
Loss on drying	<0.5%
Specific rotation	24° to 26°
Residue after ignition	<0.1%
Iron (Fe)	<10 ppm
Lead (Pb)	<10 ppm

Azomethine H
pure, For colorimetric determination of boron

packaging	Mfr. No
10 g Plastic jar	A/9080/45
C ₁₇ H ₁₂ NO ₈ S ₂ Na	
CAS: 5941-07-1	
MW: 445.39	
EINECS: 227-698-1	
Guaranteed analysis	
Complies with standard IR spectrum: Pass	

Barbitone
extra pure, 99.5+ %

SLR

packaging	Mfr. No	
250 g Plastic jar	B/0052/50	
C ₈ H ₁₂ N ₂ O ₃ CAS: 57-44-3 MW: 184.20 Density at 20°C: 1.2 Melting point: 190°C	EINECS: 200-331-2 GHS: Acute Tox. 4 H302, H332 P261, P280, P301+P312, P305+P351+P338, P302+P352	
Assay: >99.5%		
Guaranteed analysis		
Copper (Cu)	<2 ppm	
Iron (Fe)	<5 ppm	
Lead (Pb)	<5 ppm	
n-Butyl formate		<0.05%
Residue after ignition		<0.05%

Product restriction only supplied to UK customers

Barbitone sodium
for analysis

Certified AR, BP

packaging	Mfr. No	
100 g Plastic jar	B/0102/48	
C ₈ H ₁₁ N ₂ O ₃ Na CAS: 144-02-5 MW: 206.18 EINECS: 205-613-9	GHS: Acute Tox. 4 H302, H332 P261, P280, P304+P340, P301+P312	
Assay: >99%		
Guaranteed analysis		
Loss on drying	<0.5%	
Free barbitone	: <0.1%	
Calcium (Ca)	<50 ppm	
Copper (Cu)	<5 ppm	
Iron (Fe)	<5 ppm	
Lead (Pb)	<2 ppm	
Magnesium (Mg)	<5 ppm	
Potassium (K)		<100 ppm
n-Butyl formate		<0.005%
Total phosphorus (P)		<20 ppm
Total silicon (Si)		<20 ppm
Total sulphur (S)		<20 ppm
Zinc (Zn)		<10 ppm

Typical analysis

Aluminium (Al)	<1 ppm
Arsenic (As)	<2 ppm
Barium (Ba)	<1 ppm
Cadmium (Cd)	<1 ppm
Chromium (Cr)	<1 ppm
Cobalt (Co)	<1 ppm
Lithium (Li)	<1 ppm
Manganese (Mn)	<1 ppm
Nickel (Ni)	<1 ppm
Titanium (Ti)	<1 ppm
Vanadium (V)	<1 ppm

Product restriction only supplied to UK customers. Lot analysis on the label.

Barbitone sodium
extra pure, 98.5+ %

SLR

packaging	Mfr. No	
100 g Plastic jar	B/0100/48	
500 g Plastic jar	B/0100/53	
C ₈ H ₁₁ N ₂ O ₃ Na CAS: 144-02-5 MW: 206.18 EINECS: 205-613-9	GHS: Acute Tox. 4 H302, H332 P261, P280, P304+P340, P301+P312	
Assay: >98.5%		
Guaranteed analysis		
Loss on drying	<1%	
Calcium (Ca)	<200 ppm	
Copper (Cu)	<20 ppm	
Iron (Fe)	<20 ppm	
Lead (Pb)	<10 ppm	
Magnesium (Mg)	<20 ppm	
Potassium (K)		<500 ppm
n-Butyl formate		<0.01%
Total phosphorus (P)		<100 ppm
Total silicon (Si)		<100 ppm
Total sulphur (S)		<100 ppm
Zinc (Zn)		<50 ppm

Product restriction only supplied to UK customers

Barbituric acid
specified reagents for general laboratory work- extra pure

SLR

packaging	Mfr. No
100 g Plastic jar	B/0120/48
C ₄ H ₄ N ₂ O ₃ CAS: 67-52-7 MW: 128.09	Melting point: 248 to 255°C EINECS: 200-658-0
Assay: >99%	

Barium standard metal solution 1000 ppm
for AAS, solution HNO₃ 1M

packaging	Mfr. No
100 mℓ HDPE Plastic bottle	J/8004/05
500 mℓ HDPE Plastic bottle	J/8004/08
Density at 20°C: 1.0	H314
ADR/RID: 8, II	P301+P330+P331, P280, P305+P351+P338, P310, P303+P361+P353
UN: 2031	
GHS: Skin Corr. 1B	

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.

Barium standard metal solution 1000 ppm
for ICP, solution HNO₃ 1M

packaging	Mfr. No
500 mℓ HDPE Plastic bottle	J/8221/08
ADR/RID: 8, II	P301+P330+P331, P280, P305+P351+P338, P310, P303+P361+P353
UN: 2031	
GHS: Skin Corr. 1B	
H314	

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
All ICP standards are manufactured from acids, metal & metal salts Ultra pure quality and contain less than <10 ppm total metal impurites.
See also the full range of "Standard solution" in the Standard solution section.

Barium standard metal solution 10 000 ppm
for ICP, solution HNO₃ 1M

packaging	Mfr. No
100 mℓ HDPE Plastic bottle	J/8220/05
ADR/RID: 8, II	P301+P330+P331, P280, P305+P351+P338, P310, P303+P361+P353
UN: 2031	
GHS: Skin Corr. 1B, Eye Dam. 1	
H314	

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.

Barium standard metal solution 1000 ppm
for ICP solution in HNO₃ 2%, NIST certified

packaging	Mfr. No
100 mℓ Plastic bottle	SPPLBA22Y
250 mℓ Plastic bottle	SPPLBA22T
500 mℓ Plastic bottle	SPPLBA22X
H315, H319	

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product



Barium standard metal solution 1000 ppm
for ICP solution in HCL 2%, NIST certified

packaging	Mfr. No
500 mℓ Plastic bottle	PLBA12X
H314	

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurites.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Barium standard metal solution 10 000 ppm
for ICP solution in HNO₃ 5%, NIST certified

packaging	Mfr. No
100 mℓ Plastic bottle	SPPLBA23Y
500 mℓ Plastic bottle	SPPLBA23X
H314	

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurites.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Barium standard metal solution 1000 ppm
for ICP-MS solution in HNO₃ 2%, NIST certified

packaging	Mfr. No
100 mℓ Plastic bottle	SPCLBA22Y
H315, H319	

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurites.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Barium
pure, metal in liquid paraffin

packaging	Mfr. No
25 g Aluminium jar	B/0150/46
Ba CAS: 7440-39-3 MW: 137.34 Density at 20°C: 3.51 Boiling point: 1640°C Melting point: 725°C ADR/RID: 4.3, II EINECS: 231-149-1 UN: 1400	GHS: Flam. Sol. 2, Water-react. 2, Skin Irrit. 2, Eye Irrit. 2, STOT SE 3 H261, H228, H335, H315, H319, EUH014 P210, P280, P302+P352, P305+P351+P338, P304+P340, P261

Standard Volumetric solution

For your convenience we offer a range of NIST traceable solutions. A full list of our products can be found on page 323. Can't see what you need? Please email your contact details and a brief description of your requirements to begei.sha.quotes@thermofisher.com

Barium carbonate
for analysis

Certified AR

packaging	Mfr. No	
100 g Plastic jar	B/0400/48	
500 g Plastic jar	B/0400/53	
BaCO ₃ CAS: 513-77-9 MW: 197.35 Density at 20°C: 4.43 Boiling point: 1450°C Melting point: 811°C	ADR/RID: 6.1, III EINECS: 208-167-3 UN: 1564 GHS: Acute Tox. 4 H302 P280, P301+P310	
Assay: >99%		
Guaranteed analysis		
Calcium (Ca)	<200 ppm	
Copper (Cu)	<10 ppm	
Iron (Fe)	<10 ppm	
Lead (Pb)	<10 ppm	
Magnesium (Mg)	<20 ppm	
Potassium (K)	<50 ppm	
Sodium (Na)		<100 ppm
n-Butyl formate		<0.002%
Total phosphorus (P)		<10 ppm
Total sulphur (S)		<10 ppm
Zinc (Zn)		<10 ppm

Lot analysis on the label.

Barium carbonate
specified reagents for general laboratory work- extra pure

SLR

packaging	Mfr. No	
500 g Plastic jar	B/0350/53	
BaCO ₃ CAS: 513-77-9 MW: 197.35 Density at 20°C: 4.43 Boiling point: 1450°C Melting point: 811°C	ADR/RID: 6.1, III EINECS: 208-167-3 UN: 1564 GHS: Acute Tox. 4 H302 P280, P301+P310	
Assay: >99%		
Guaranteed analysis		
Calcium (Ca)	<500 ppm	
Copper (Cu)	<50 ppm	
Iron (Fe)	<50 ppm	
Lead (Pb)	<50 ppm	
Magnesium (Mg)	<100 ppm	
Potassium (K)	<200 ppm	
Sodium (Na)		<500 ppm
n-Butyl formate		<0.01%
Total phosphorus (P)		<50 ppm
Total sulphur (S)		<50 ppm
Zinc (Zn)		<50 ppm

Barium chloride dihydrate
for analysis

Certified AR

packaging	Mfr. No	
500 g Plastic jar	B/0550/53	
1 kg Plastic jar	B/0550/60	
3 kg Plastic jar	B/0550/63	
5 kg Plastic jar	B/0550/65	
25 kg Plastic keg	B/0550/70	
BaCl ₂ ·2H ₂ O CAS: 10326-27-9 MW: 244.28 Boiling point: 1560°C Melting point: 962°C ADR/RID: 6.1, III	EINECS: 233-788-1 UN: 1564 GHS: Acute Tox. 3 H301, H332 P301+P310, P304+P340	
Assay: 99 to 101%		
Guaranteed analysis		
pH (aqueous solution 5%)	5 to 8	
Calcium (Ca)	<20 ppm	
Copper (Cu)	<2 ppm	
Iron (Fe)	<1 ppm	
Lead (Pb)	<2 ppm	
Magnesium (Mg)	<10 ppm	
Potassium (K)	<20 ppm	
Sodium (Na)		<50 ppm
Total nitrogen (N)		<20 ppm
Total phosphorus (P)		<10 ppm
Total silicon (Si)		<10 ppm
Total sulphur (S)		<20 ppm
Zinc (Zn)		<5 ppm

Typical analysis

Aluminium (Al)	<2 ppm
Silver (Ag)	<1 ppm
Arsenic (As)	<2 ppm
Cadmium (Cd)	<1 ppm
Chromium (Cr)	<2 ppm
Cobalt (Co)	<1 ppm
Lithium (Li)	<1 ppm
Manganese (Mn)	<1 ppm
Molybdenum (Mo)	<1 ppm
Nickel (Ni)	<1 ppm
Strontium (Sr)	<0.01 %
Titanium (Ti)	<1 ppm
Vanadium (V)	<1 ppm

Lot analysis on the label.

Barium chloride dihydrate			SLR
specified reagents for general laboratory work- extra pure			
packaging			Mfr. No
500 g	Plastic jar		B/0500/53
1 kg	Plastic jar		B/0500/60
5 kg	Plastic jar		B/0500/65
25 kg	Plastic keg		B/0500/70
BaCl ₂ ·2H ₂ O		EINECS: 233-788-1	
CAS: 10326-27-9		UN: 1564	
MW: 244.28		GHS: Acute Tox. 3	
Boiling point: 1560°C		H301, H332	
Melting point: 962°C		P301+P310, P304+P340	
ADR/RID: 6.1, III			
Assay: 99 to 101%			
Guaranteed analysis			
Calcium (Ca)	<200 ppm	Sodium (Na)	<200 ppm
Copper (Cu)	<20 ppm	Total phosphorus (P)	<100 ppm
Iron (Fe)	<10 ppm	Total silicon (Si)	<50 ppm
Lead (Pb)	<20 ppm	Total sulphur (S)	<200 ppm
Magnesium (Mg)	<10 ppm	Zinc (Zn)	<50 ppm
Potassium (K)	<100 ppm		

Barium chloride solution 1M (2N)
pure, ready to use

packaging	Mfr. No
1 ℓ Coated Glass bottle	J/1700/PB15
ADR/RID: 6.1, III UN: 3287	

Manufactured to a tolerance of ± 0.5%.

Barium diphenylamine-4-sulfonate
pure, Redox indicator

packaging	Mfr. No
5 g Plastic jar	B/0620/44
C ₂₄ H ₂₀ BaN ₂ O ₆ S ₂ CAS: 6211-24-1 MW: 633.89 Melting point: 300°C ADR/RID: 6.1, II	EINECS: 228-278-0 UN: 1564 GHS: Acute Tox. 4 H302, H332 P301+P312, P304+P340
Guaranteed analysis	
Identity test : Pass	

Barium hydroxide octahydrate
for analysis

packaging		Mfr. No	
500 g	Plastic jar	B/0700/53	
1 kg	Plastic jar	B/0700/60	
Ba(OH) ₂ ·8H ₂ O		UN: 3262	
CAS: 12230-71-6		GHS: Acute Tox. 4, Skin	
MW: 315.48		Corr. 1B, Eye Dam. 1	
Density at 20°C: 2.18		H332, H302, H314	
Melting point: 78°C		P280, P305+P351+P338,	
ADR/RID: 8, III		P308+P313	
EINECS: 241-234-5			
Assay: 98 to 102%			
Guaranteed analysis			
Calcium (Ca)	<50 ppm	Sodium (Na)	<100 ppm
Carbonate (Na ₂ CO ₃)	<2%	Total chloride (Cl)	<0.001%
Copper (Cu)	<5 ppm	Total phosphorus (P)	<10 ppm
Iron (Fe)	<5 ppm	Total silicon (Si)	<500 ppm
Lead (Pb)	<5 ppm	Total sulphur (S)	<10 ppm
Magnesium (Mg)	<5 ppm	Zinc (Zn)	<5 ppm
Potassium (K)	<50 ppm		

Typical analysis			
Aluminium (Al)	<2 ppm	Manganese (Mn)	<1 ppm
Silver (Ag)	<1 ppm	Molybdenum (Mo)	<1 ppm
Arsenic (As)	<2 ppm	Nickel (Ni)	<2 ppm
Cadmium (Cd)	<1 ppm	Strontium (Sr)	<0.5 %
Chromium (Cr)	<2 ppm	Titanium (Ti)	<1 ppm
Cobalt (Co)	<1 ppm	Vanadium (V)	<1 ppm
Lithium (Li)	<1 ppm		

Lot analysis on the label.

Barium hydroxide octahydrate
specified reagents for general laboratory work- extra pure

packaging			Mfr. No
500 g	Plastic jar		B/0650/53
1 kg	Plastic jar		B/0650/60
Ba(OH) ₂ ·8H ₂ O		UN: 3262	
CAS: 12230-71-6		GHS: Acute Tox. 4, Skin	
MW: 315.48		Corr. 1B, Eye Dam. 1	
Density at 20°C: 2.18		H332, H302, H314	
Melting point: 78°C		P280, P305+P351+P338,	
ADR/RID: 8, III		P308+P313	
EINECS: 241-234-5			
Assay: 97 to 103%			
Guaranteed analysis			
Calcium (Ca)	<200 ppm	Potassium (K)	<500 ppm
Carbonate (Na ₂ CO ₃)	<2%	Sodium (Na)	<1000 ppm
Copper (Cu)	<50 ppm	n-Butyl formate	<0.005%
Iron (Fe)	<50 ppm	Total phosphorus (P)	<100 ppm
Lead (Pb)	<50 ppm	Total sulphur (S)	<500 ppm
Magnesium (Mg)	<50 ppm	Zinc (Zn)	<50 ppm

Barium nitrate
for analysis

packaging		Mfr. No	
250 g	Plastic jar	B/0850/50	
500 g	Plastic jar	B/0850/53	
1 kg	Plastic jar	B/0850/60	
25 kg	Fibre Keg	B/0850/70	
Ba(NO ₃) ₂		EINECS: 233-020-5	
CAS: 10022-31-8		UN: 1446	
MW: 261.35		GHS: Ox. Sol. 2, Acute Tox. 4	
Density at 20°C: 3.23		H272, H302, H332	
Melting point: 592°C		P210, P301+P312, P304+P340	
ADR/RID: 5.1, 6.1, II			
Assay: >99%			
Guaranteed analysis			
pH (aqueous solution 5%)	5 to 8	Potassium (K)	<50 ppm
Calcium (Ca)	<10 ppm	Sodium (Na)	<50 ppm
Copper (Cu)	<5 ppm	Total phosphorus (P)	<20 ppm
Iron (Fe)	<2 ppm	Total silicon (Si)	<20 ppm
Lead (Pb)	<5 ppm	Total sulphur (S)	<20 ppm
Magnesium (Mg)	<10 ppm	Zinc (Zn)	<20 ppm

Typical analysis			
Aluminium (Al)	<1 ppm	Manganese (Mn)	<1 ppm
Silver (Ag)	<1 ppm	Molybdenum (Mo)	<1 ppm
Arsenic (As)	<2 ppm	Nickel (Ni)	<1 ppm
Cadmium (Cd)	<1 ppm	Strontium (Sr)	<0.05 %
Chromium (Cr)	<1 ppm	Titanium (Ti)	<1 ppm
Cobalt (Co)	<1 ppm	Vanadium (V)	<1 ppm
Lithium (Li)	<1 ppm		

Lot analysis on the label.

UHPLC

Filtered to 0.1 µm, Fisher Chemical UHPLC gradient grade solvents (Acetonitrile, Methanol and Water) are perfect for those working at high pressure.

See page 309 for more information



Barium nitrate			SLR
specified reagents for general laboratory work- extra pure			
packaging			Mfr. No
500 g	Plastic jar		B/0800/53
25 kg	Plastic keg		B/0800/70
Ba(NO ₃) ₂		EINECS: 233-020-5	
CAS: 10022-31-8		UN: 1446	
MW: 261.35		GHS: Ox. Sol. 2, Acute Tox. 4	
Density at 20°C: 3.23		H272, H302, H332	
Melting point: 592°C		P210, P301+P312, P304+P340	
ADR/RID: 5.1, 6.1, II			
Assay: >95%			

Barium peroxide			SLR
specified reagents for general laboratory work- extra pure			
packaging			Mfr. No
500 g	Plastic jar		B/1080/53
BaO ₂	EINECS: 215-128-4	 	
CAS: 1304-29-6	UN: 1449		
MW: 169.36	GHS: Ox. Sol. 2, Acute Tox. 4		
Density at 20°C: 4.96	H272, H332, H302		
Melting point: 450°C	P210, P304+P340, P301+P312		
ADR/RID: 5.1, 6.1, II			
Assay: >85%			
Guaranteed analysis			
Calcium (Ca)	<10000 ppm	n-Butyl formate	<0.04%
Iron (Fe)	<200 ppm	Insoluble Acids	<3%

Barium sulfate precipitated			SLR
<i>specified reagents for general laboratory work- extra pure</i>			
packaging			Mfr. No
500 g	Plastic jar		B/1120/53
2 kg	Plastic jar		B/1120/61
BaSO ₄		Density at 20°C: 4.5	
CAS: 7727-43-7		Melting point: 1580°C	
MW: 233.40		EINECS: 231-784-4	
Guaranteed analysis			
Identity test : Pass			

Benedict's reagent quantitative
pure

packaging	Mfr. No
1 ℓ Glass bottle	J/2150/15
GHS: Aquatic Chronic 3 H412, EUH032	

50 mL of solution is equivalent to 0.1 g glucose.
Blue to white colour change.

Benedict's reagent qualitative
pure

packaging	Mfr. No
1 ℓ Glass bottle	J/2100/15
2.5 ℓ Glass bottle	J/2100/17
Density at 20°C: 1.21 GHS: Eye Irrit. 2, Aquatic Chronic 3 H412, H319 P273, P280, P305+P351+P338	

Bentonite		SLR
<i>specified reagents for general laboratory work- extra pure</i>		
packaging		Mfr. No
1 kg	Plastic jar	B/1300/60
CAS: 1302-78-9		
EINECS: 215-108-5		
Guaranteed analysis		
Identity test : Pass		

Benzaldehyde for analysis		Certified AR	
packaging		Mfr. No	
250 ml	Coated Glass bottle	B/1450/PB07	
C ₇ H ₆ O	ADR/RID: 9, III		
CAS: 100-52-7	EINECS: 202-860-4		
MW: 106.13	UN: 1990		
Density at 20°C: 1.043	GHS: Acute Tox. 4		
Boiling point: 179°C	H302		
Melting point: -26°C	P280, P301+P310		
Flash point: 64°C			
Assay (CG): >99%			
Guaranteed analysis			
Colour	<40 APHA	Lead (Pb)	<1 ppm
Sulfated ash	<0.005%	Magnesium (Mg)	<1 ppm
Acidity	<0.02 meq/g	Potassium (K)	<20 ppm
Calcium (Ca)	<5 ppm	Sodium (Na)	<20 ppm
Ionized chloride	<0.005%	Total phosphorus (P)	<1 ppm
Copper (Cu)	<1 ppm	Total sulphur (S)	<1 ppm
Iron (Fe)	<1 ppm	Zinc (Zn)	<1 ppm

Lot analysis on the label

Benzaldehyde			SLR
specified reagents for general laboratory work- extra pure			
packaging			Mfr. No
500 ml Coated Glass bottle			B/1400/PB08
C ₇ H ₆ O	ADR/RID: 9, III		
CAS: 100-52-7	EINECS: 202-860-4		
MW: 106.13	UN: 1990		
Density at 20°C: 1.043	GHS: Acute Tox. 4		
Boiling point: 179°C	H302		
Melting point: -26°C	P280, P301+P310		
Flash point: 64°C			
Assay (CG): >99%			
Guaranteed analysis			
Colour	<50 APHA	Lead (Pb)	<5 ppm
Sulfated ash	<0.05%	Magnesium (Mg)	<5 ppm
Acidity	<0.1 meq/g	Potassium (K)	<100 ppm
Calcium (Ca)	<20 ppm	Sodium (Na)	<100 ppm
Ionized chloride	<0.05%	Total phosphorus (P)	<5 ppm
Copper (Cu)	<5 ppm	Total sulphur (S)	<5 ppm
Iron (Fe)	<5 ppm	Zinc (Zn)	<5 ppm

Benzamide			SLR
specified reagents for general laboratory work- extra pure			
packaging			Mfr. No
100 g	Plastic jar		B/1460/48
C ₇ H ₇ NO	Flash point: 180°C		
CAS: 55-21-0	EINECS: 200-227-7		
MW: 121.14	GHS: Muta. 2, Acute Tox. 4		
Boiling point: 228°C	H341, H302		
Melting point: 125 to 129°C	P280, P301+P312, P260		
Assay: >98%			

1.3-Benzenediol See Resorcinol, page 269

1.2-Benzenediol *See Catechol, page 80*

Benzenesulfonyl chloride SLR
specified reagents for general laboratory work- extra pure

packaging		Mfr. No
250 mℓ	Coated Glass bottle	B/1800/PB07
C ₆ H ₅ ClO ₂ S	EINECS: 202-636-6	
CAS: 98-09-9	UN: 2225	
MW: 176.63	GHS: Acute Tox. 4, Skin	
Density at 20°C: 1.384	Corr. 1B, Eye Dam. 1	
Boiling point: 251 to 252°C	H314, H302, EUH029	
Melting point: 14.5°C	P301+P330+P331, P280,	
Flash point: 128°C	P305+P351+P338, P310,	
ADR/RID: 8, III	P402+P404, P301+P312	

Benzoic acid Certified AR
for analysis

packaging		Mfr. No
100 g	Plastic jar	B/1950/48
500 g	Plastic jar	B/1950/53
1 kg	Plastic jar	B/1950/60
C ₇ H ₆ O ₂	2, Eye Irrit. 2, Resp. Sens. 1,	
CAS: 65-85-0	Skin Sens. 1, STOT SE 3	
MW: 122.12	H334, H317, H302, H319,	
Boiling point: 249°C	H315, H335	
Melting point: 121 to 123°C	P261, P304+P341,	
Flash point: 121°C	P342+P311, P280,	
EINECS: 200-618-2	P301+P312, P302+P352,	
GHS: Acute Tox. 4, Skin Irrit.	P305+P351+P338	
Assay: 99.9 to 100.5%		
Guaranteed analysis		
Water	<0.7%	Magnesium (Mg)
Substances		<2 ppm
reducing KMnO ₄	<0.01%	Potassium (K)
Freezing Point	: 121 to 123°C	<5 ppm
Heavy metals	<10 ppm	Sodium (Na)
Calcium (Ca)	<5 ppm	<5 ppm
Chlorine (Cl) compounds	<0.005%	n-Butyl formate
Copper (Cu)	<2 ppm	<0.0005%
Iron (Fe)	<2 ppm	Total phosphorus (P)
Lead (Pb)	<2 ppm	<20 ppm
		Calcium (Ca)
		<10 ppm
		Total silicon (Si)
		<10 ppm
		Total sulphur (S)
		<20 ppm
		Copper (Cu)
		<2 ppm
		Iron (Fe)
		<2 ppm
		Readily carbonisable substances test: Pass
		Residue after ignition
		<0.05%

Lot analysis on the label.

Benzoic acid SLR
specified reagents for general laboratory work- extra pure

packaging		Mfr. No
500 g	Plastic jar	B/1900/53
1 kg	Plastic jar	B/1900/60
25 kg	Fibre Keg	B/1900/70
C ₇ H ₆ O ₂	2, Eye Irrit. 2, Resp. Sens. 1,	
CAS: 65-85-0	Skin Sens. 1, STOT SE 3	
MW: 122.12	H334, H317, H302, H319,	
Boiling point: 249°C	H315, H335	
Melting point: 121 to 123°C	P261, P304+P341,	
Flash point: 121°C	P342+P311, P280,	
EINECS: 200-618-2	P301+P312, P302+P352,	
GHS: Acute Tox. 4, Skin Irrit.	P305+P351+P338	
Assay: >99.5%		
Guaranteed analysis		
Calcium (Ca)	<20 ppm	Sodium (Na)
Copper (Cu)	<10 ppm	<20 ppm
Iron (Fe)	<10 ppm	n-Butyl formate
Lead (Pb)	<5 ppm	<0.002%
Magnesium (Mg)	<10 ppm	Total phosphorus (P)
Potassium (K)	<20 ppm	<100 ppm
		Total silicon (Si)
		<20 ppm
		Total sulphur (S)
		<100 ppm
		Zinc (Zn)
		<10 ppm

Benzoin SLR
specified reagents for general laboratory work- extra pure

Packaging		Mfr. No
250 g	Plastic jar	B/2050/50
$C_{14}H_{12}O_2$		Boiling point: 194°C
CAS: 119-53-9		Melting point: 133 to 137°C
MW: 212.25		EINECS: 204-331-3
Guaranteed analysis		
Residue after ignition	<0.2%	

p-Benzoquinone SLR
specified reagents for general laboratory work- extra pure

packaging		Mfr. No
100 g	Glass jar	B/2150/48
C ₆ H ₄ O ₂	GHS: Acute Tox. 3, Skin Irrit.	
CAS: 106-51-4	2, Eye Irrit. 2, STOT SE 3,	
MW: 108.10	Aquatic Acute 1	
Density at 20°C: 1.31	H331, H315, H335, H319,	
Melting point: 112 to 116°C	H400, H301	
ADR/RID: 6.1, II	P304+P340, P261, P273,	
EINECS: 203-405-2	P302+P352, P280,	
UN: 2587	P305+P351+P338	
Assay: >98%		
Guaranteed analysis		
Water	<1%	Residue after ignition
		<0.5%

Benzoyl chloride Certified AR
for analysis

packaging		Mfr. No
500 mℓ	Coated Glass bottle	B/2350/PB08
C ₇ H ₅ ClO	UN: 1736	
CAS: 98-88-4	GHS: Acute Tox. 4, Skin	
MW: 140.57	Corr. 1B, Skin Sens. 1	
Density at 20°C: 1.21	H317, H302, H314, H312, H332	
Boiling point: 198°C	P280, P302+P350,	
Melting point: -1°C	P301+P330+P331,	
Flash point: 93°C	P305+P351+P338,	
ADR/RID: 8, II	P304+P340, P309+P311	
EINECS: 202-710-8		
Assay: >99.5%		
Guaranteed analysis		
Colour	<20 APHA	Magnesium (Mg)
Residue after ignition	<0.02%	<2 ppm
Calcium (Ca)	<5 ppm	Potassium (K)
Copper (Cu)	<0.1 ppm	<5 ppm
Iron (Fe)	<0.5 ppm	Sodium (Na)
Lead (Pb)	<0.1 ppm	<10 ppm
		Total phosphorus (P)
		<10 ppm
		Total sulphur (S)
		<10 ppm
		Zinc (Zn)
		<1 ppm

Benzoyl chloride SLR
specified reagents for general laboratory work- extra pure

packaging		Mfr. No
250 mℓ	Coated Glass bottle	B/2300/PB07
500 mℓ	Coated Glass bottle	B/2300/PB08
C ₆ H ₅ COCl	UN: 1736	
CAS: 98-88-4	GHS: Acute Tox. 4, Skin	
MW: 140.57	Corr. 1B, Skin Sens. 1	
Density at 20°C: 1.21	H317, H302, H314, H312, H332	
Boiling point: 198°C	P280, P302+P350,	
Melting point: -1°C	P301+P330+P331,	
Flash point: 93°C	P305+P351+P338,	
ADR/RID: 8, II	P304+P340, P309+P311	
EINECS: 202-710-8		
Assay: >98%		
Guaranteed analysis		
Colour	<30 APHA	Magnesium (Mg)
Residue after ignition	<0.05%	<10 ppm
Calcium (Ca)	<250 ppm	Potassium (K)
Copper (Cu)	<1 ppm	<250 ppm
Iron (Fe)	<1 ppm	Sodium (Na)
Lead (Pb)	<1 ppm	<50 ppm
		Total phosphorus (P)
		<100 ppm
		Total sulphur (S)
		<100 ppm
		Zinc (Zn)
		<5 ppm

Benzyl alcohol SLR, Eur.Ph., BP, USP
specified reagents for general laboratory work- extra pure

Packaging		Mfr. No
500 mℓ	Glass bottle	B/2550/08
2.5 ℓ	Coated Glass bottle	B/2550/PB17
<chem>c1ccccc1CO</chem>	Flash point: 96°C	
CAS: 100-51-6	EINECS: 202-859-9	
MW: 108.14	GHS: Acute Tox. 4	
Density at 20°C: 1.044	H332, H302	
Boiling point: 205°C	P280, P304+P340, P301+P312	
Melting point: -15°C		
Assay (CG): >98%		
Guaranteed analysis		
Residue after evaporation	<500 ppm	Benzaldehyde
Chlorine (Cl) compounds	<0.05%	<0.1%

Benzylamine SLR
extra pure, 98+%, (acidimetric)

packaging		Mfr. No
100 mℓ	Coated Glass bottle	B/2600/PB05
C ₇ H ₉ N	EINECS: 202-854-1	
CAS: 100-46-9	UN: 2735	
MW: 107.16	GHS: Acute Tox. 4, Skin Corr. 1B	
Density at 20°C: 0.980	H318, H312, H302	
Boiling point: 182 to 185°C	P280, P301+P330+P331,	
Melting point: -30°C	P305+P351+P338,	
Flash point: 72°C	P302+P352, P301+P312	
ADR/RID: 8, II		
Assay: >98%		

Benzyl chloroformate pure

packaging		Mfr. No
100 mℓ	Coated Glass bottle	B/2710/PB05
C ₇ H ₅ ClO ₂	Dam. 1, Aquatic Acute 1,	
CAS: 501-53-1	Aquatic Chronic 1	
MW: 170.60	H410, H314	
ADR/RID: 8, I	P301+P330+P331, P280,	
EINECS: 207-925-0	P305+P351+P338, P310,	
UN: 1739	P261, P273	
GHS: Skin Corr. 1B, Eye		

Beryllium standard metal solution 1000 ppm *for ICP solution in HNO₃ 2%, NIST certified*

packaging		Mfr. No
100 mℓ	Plastic bottle	SPPLBE22Y
250 mℓ	Plastic bottle	SPPLBE22T
500 mℓ	Plastic bottle	SPPLBE22X
H315, H319, H350		

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Beryllium standard metal solution 10 000 ppm *for ICP solution in HNO₃ 5%, NIST certified*

packaging		Mfr. No
100 mℓ	Plastic bottle	SPPLBE23Y
500 mℓ	Plastic bottle	SPPLBE23X
H302, H311, H314, H331, H350		

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Beryllium standard metal solution 1000 ppm *for ICP-MS solution in HNO₃ 2%, NIST certified*

packaging		Mfr. No
100 mℓ	Plastic bottle	SPCLBE22Y
H315, H320, H350		

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Bicarbonate indicator solution *pure, Concentrated stock reagent*

packaging		Mfr. No
250 mℓ	Glass bottle	J/2300/07
500 mℓ	Glass bottle	J/2300/08

Biphenyl pure

packaging		Mfr. No
500 g	Plastic jar	B/2975/53
C ₁₂ H ₁₀	EINECS: 202-163-5	
CAS: 92-52-4	UN: 3077	
GHS: Skin Irrit. 2, Eye Irrit. 2,	STOT SE 3, Aquatic Acute 1,	
Density at 20°C: 0.991	Aquatic Chronic 1	
Boiling point: 255°C	H315, H410, H319, H335	
Melting point: 68 to 72°C	P261, P273, P302+P352,	
Flash point: 113°C	P280, P305+P351+P338	
ADR/RID: 9, III		
Guaranteed analysis		
Complies with standard IR spectrum: Pass		

Bis(2-hydroxyethyl)-imino-tris-(hydroxymethyl)-methane

SLR

specified reagents for general laboratory work- extra pure

packaging	Mfr. No		
100 g Plastic jar	B/3010/48		
C ₈ H ₁₉ NO ₃ CAS: 6976-37-0 MW: 209.24 Melting point: 103 to 104°C ADR/RID: 8, III EINECS: 230-237-7 UN: 3259		GHS: Met. Corr. 1, Skin Irrit. 2, Eye Irrit. 2, STOT SE 3 H290, H315, H319, H335 P234, P390, P261, P302+P352, P280, P305+P351+P338	 
Guaranteed analysis Complies with standard IR spectrum: Pass			

Bismuth standard metal solution 1000 ppm

for ICP solution in HNO₃ 10%, NIST certified

packaging	Mfr. No	
100 mL Plastic bottle	SPPLBI42Y	
500 mL Plastic bottle	SPPLBI42X	
H314		

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Bismuth standard metal solution 10 ppm

for ICP-MS solution in HNO₃ 2%, NIST certified

packaging	Mfr. No	
100 mL Plastic bottle	SPCLBI21AY	
H315, H320		

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Bismuth(III) chloride

SLR

specified reagents for general laboratory work- extra pure

packaging	Mfr. No	
25 g Glass jar	B/3110/46	
BiCl ₃ CAS: 7787-60-2 MW: 315.33 Density at 20°C: 4.75 Boiling point: 447°C Melting point: 230 to 232°C		EINECS: 232-123-2 GHS: Skin Irrit. 2, Eye Irrit. 2, STOT SE 3 H315, H319, H335 P261, P302+P352, P280, P305+P351+P338
Guaranteed analysis Identity test : Pass		

High Volume Solvent Delivery Systems

High-volume solvent delivery systems are available in 10L to 1000L, in Europe only. Please email your contact details and a brief description of your requirements to begel.sha.quotes@thermofisher.com



Bismuth nitrate pentahydrate

Certified AR

for analysis

Packaging		Mfr. No	
100 g	Plastic jar	B/3200/48	
Bi(NO₃)₃·5H₂O CAS: 10035-06-0 MW: 485.07 Density at 20°C: 2.83 Boiling point: 80°C Melting point: 30°C ADR/RID: 5.1, II		EINECS: 233-791-8 UN: 1477 GHS: Ox. Sol. 2, Skin Irrit. 2, Eye Irrit. 2, STOT SE 3 H272, H315, H319, H335 P210, P261, P302+P352, P280, P305+P351+P338 	
Assay: >98%			
Guaranteed analysis			
Ammonia (NH ₃)	<0.005%	Sodium (Na)	<100 ppm
Calcium (Ca)	<20 ppm	n-Butyl formate	<0.001%
Copper (Cu)	<20 ppm	Total phosphorus (P)	<100 ppm
Iron (Fe)	<50 ppm	Total silicon (Si)	<50 ppm
Lead (Pb)	<10 ppm	Total sulphur (S)	<200 ppm
Magnesium (Mg)	<10 ppm	Zinc (Zn)	<10 ppm
Potassium (K)	<20 ppm		
Not analysis on the label.			

Lot analysis on the label.

3,3'-Bis[N,N-di(carboxymethyl)aminomethyl]thymolsulfonephthalein, sodium salt

See Thymolsulfonephthaleindi-(methyl-iminodiacetic acid) sodium salt, page 349

Bispyrazolone (3,3'-dimethyl-1,1'-diphenyl(4,4'-bi-2-pyrazoline)-5,5'-dione) pure

packaging	Mfr. No	
5 g Plastic jar	B/3620/44	
C ₂₀ H ₁₈ N ₄ O ₂ CAS: 7477-67-0 MW: 346.39		Melting point: 355°C EINECS: 231-282-5
Guaranteed analysis Complies with standard IR spectrum: Pass		

Borax See di-Sodium tetraborate decahydrated, page 305

Boric acid

for electrophoresis

packaging	Mfr. No	
500 g Plastic jar	B/P560/53	
1 kg Plastic jar	B/P560/60	
5 kg Plastic keg	B/P560/65	
H ₃ BO ₃ CAS: 10043-35-3 MW: 61.83 Melting point: 169°C EINECS: 233-139-2		GHS: Repr. 1B H360FD P310, P201, P308+P313, P304+P340
Assay: > 99.8%		

We Support Your Chemistry For

- Semi-bulk and bulk requirements
- Tailored solvents and solvent blends or customized solutions
- Testing services or customized packaging and labeling

Please email your contact details and a brief description of your requirements to begel.sha.quotes@thermofisher.com



Boric acid

for electronic use (MOS)

Packaging		Mfr. No	
500 g	Plastic jar	B/3820/53	
H ₃ BO ₃ CAS: 10043-35-3 MW: 61.83 Melting point: 169°C EINECS: 233-139-2	GHS: Repr. 1B H360FD P310, P201, P308+P313, P304+P340		
			
Assay: >99.5%			
Guaranteed analysis			
Water insoluble matter	<0.005%	Phosphate (PO ₄)	<10 ppm
Calcium (Ca)	<50 ppm	Sodium (Na)	<10 ppm
Copper (Cu)	<2 ppm	Sulfate (SO ₄)	<0.0002%
Iron (Fe)	<2 ppm	n-Butyl formate	<0.0002%
Lead (Pb)	<2 ppm		

Boric acid

for analysis

Certified AR

packaging		Mfr. No	
500 g	Plastic jar	B/3800/53	
1 kg	Plastic jar	B/3800/60	
2.5 kg	Plastic jar	B/3800/62	
25 kg	Fibre Keg	B/3800/70	
H ₃ BO ₃ CAS: 10043-35-3 MW: 61.83 Melting point: 169°C EINECS: 233-139-2		GHS: Repr. 1B H360FD P310, P201, P308+P313, P304+P340	
			
Assay: >99.8%			
Guaranteed analysis			
Water insoluble matter	<0.005%	Potassium (K)	<10 ppm
Insoluble in methanol	: <0.005%	Sodium (Na)	<20 ppm
Non volatile with Methanol	: <0.05%	Total chloride (Cl)	<0.0002%
Calcium (Ca)	<10 ppm	Total phosphorus (P)	<10 ppm
Copper (Cu)	<2 ppm	Total silicon (Si)	<10 ppm
Iron (Fe)	<2 ppm	Total sulphur (S)	<10 ppm
Lead (Pb)	<2 ppm	Zinc (Zn)	<5 ppm
Magnesium (Mg)	<2 ppm	Phosphate (PO ₄)	<10 ppm
Non-esterifiable impurities	<0.05%	Sulfate (SO ₄)	<0.01%
Typical analysis			
Aluminium (Al)	<2 ppm	Lithium (Li)	<1 ppm
Silver (Ag)	<1 ppm	Manganese (Mn)	<1 ppm
Arsenic (As)	<1 ppm	Molybdenum (Mo)	<1 ppm
Barium (Ba)	<1 ppm	Nickel (Ni)	<1 ppm
Cadmium (Cd)	<1 ppm	Strontium (Sr)	<1 ppm
Chromium (Cr)	<1 ppm	Titanium (Ti)	<1 ppm
Cobalt (Co)	<1 ppm	Vanadium (V)	<1 ppm
Lot analysis on the label.			

Lot analysis on the label.

Boric acid

extra pure, granular

SLR, Eur.Ph

packaging		Mfr. No	
500 g	Plastic jar	B/3750/53	
1 kg	Plastic jar	B/3750/60	
5 kg	Plastic keg	B/3750/65	
H ₃ BO ₃ CAS: 10043-35-3 MW: 61.83 Melting point: 169°C EINECS: 233-139-2		GHS: Repr. 1B H360FD P310, P201, P308+P313, P304+P340 	
Assay: >99.5%			
Guaranteed analysis			
Calcium (Ca)	<100 ppm	Sodium (Na)	<200 ppm
Copper (Cu)	<20 ppm	n-Butyl formate	<0.01%
Iron (Fe)	<20 ppm	Total phosphorus (P)	<100 ppm
Lead (Pb)	<20 ppm	Total silicon (Si)	<100 ppm
Magnesium (Mg)	<20 ppm	Total sulphur (S)	<200 ppm
Potassium (K)	<100 ppm	Zinc (Zn)	<50 ppm

Boric acid

extra pure, powder

SLR, Eur. Ph.

Packaging		Mfr. No	
600 g	Plastic jar	B/3760/53	
1 kg	Plastic jar	B/3760/60	
5 kg	Plastic keg	B/3760/65	
25 kg	Fibre Keg	B/3760/70	
H₃BO₃ CAS: 10043-35-3 MW: 61.83 Melting point: 169°C EINECS: 233-139-2		GHS: Repr. 1B H360FD P310, P201, P308+P313, P304+P340 	
Assay: >99.5%			
Guaranteed analysis			
Calcium (Ca)	<100 ppm	Sodium (Na)	<200 ppm
Copper (Cu)	<20 ppm	n-Butyl formate	<0.01%
Iron (Fe)	<20 ppm	Total phosphorus (P)	<100 ppm
Lead (Pb)	<20 ppm	Total silicon (Si)	<100 ppm
Magnesium (Mg)	<20 ppm	Total sulphur (S)	<200 ppm
Potassium (K)	<100 ppm	Zinc (Zn)	<50 ppm

Boron standard metal solution 1000 ppm

for AAS, solution aqueous

packaging	Mfr. No	
100 mL HDPE Plastic bottle	J/8006/05	
500 mL HDPE Plastic bottle	J/8006/08	
Density at 20°C: 1		

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.

Boron standard metal solution 1000 ppm

for ICP, solution aqueous

packaging	Mfr. No	
500 mL HDPE Plastic bottle	J/8231/08	
Density at 20°C: 1		

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.

Boron standard metal solution 10 000 ppm

for ICP, solution NH₃ 1%

packaging	Mfr. No	
100 mL HDPE Plastic bottle	J/8230/05	

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.

Boron standard metal solution 1000 ppm

for ICP solution in H₂O, NIST certified

packaging	Mfr. No	
100 mL Plastic bottle	SPPLB92Y	
250 mL Plastic bottle	SPPLB92T	
500 mL Plastic bottle	SPPLB92X	

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Boron standard metal solution 10 000 ppm
for ICP solution in H₂O, NIST certified

packaging	Mfr. No
100 mℓ Plastic bottle	SPPLB93Y
500 mℓ Plastic bottle	SPPLB93X

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Boric oxideSLR
specified reagents for general laboratory work- extra pure

packaging	Mfr. No
500 g Plastic jar	B/3850/53
B ₂ O ₃	ADR/RID: 8/III
CAS: 1303-86-2	EINECS: 215-125-8
MW: 69.61	UN: 3260
Density at 20°C: 2.46	GHS: Repr. 1B
Boiling point: 1860°C	H360FD
Melting point: 450°C	P201, P308+P313
Assay: >98%	

Bromate standard 1000 ppm
For Ion chromatography in H₂O, NIST certified

packaging	Mfr. No
100 mℓ Plastic bottle	SPASBRO392Y
Manufactured to a tolerance of ± 0.3%.	
Check with your local distributor for availability of Spex product	

Bromide standard solution 1000 ppm (1 mg/mℓ)
for Ion pair chromatography, Ready-to-use

packaging	Mfr. No
100 mℓ HDPE Plastic bottle	J/4544/05
Manufactured to a tolerance of ± 1%.	

BrominePrimar
for trace metal analysis

packaging	Mfr. No
100 mℓ Coated Glass bottle	B/4055/PB05
Br ₂	UN: 1744
CAS: 7726-95-6	GHS: Acute Tox. 2, Skin
MW: 159.82	Corr. 1A, Eye Dam. 1, Aquatic Acute 1
Density at 20°C: 3.111	H314, H330, H400
Boiling point: 58.7°C	P301+P330+P331, P280,
Melting point: -7.2°C	P305+P351+P338, P310,
ADR/RID: 8, 6.1, I	P304+P340, P273
EINECS: 231-778-1	

Assay: >99.5%	
Guaranteed analysis	
Other organic impurities test	: Pass
Aluminium (Al)	<0.02 ppm
Arsenic (As)	<1 ppm
Barium (Ba)	<0.005 ppm
Cadmium (Cd)	<0.005 ppm
Calcium (Ca)	<0.5 ppm
Chromium (Cr)	<0.005 ppm
Cobalt (Co)	<0.005 ppm
Copper (Cu)	<0.005 ppm
Iodine	<0.5 ppm
Iron (Fe)	<0.05 ppm
Lead (Pb)	<0.01 ppm
Lithium (Li)	<0.005 ppm
Magnesium (Mg)	<0.05 ppm
Manganese (Mn)	<0.005 ppm
Molybdenum (Mo)	<0.005 ppm
Nickel (Ni)	<0.005 ppm
Potassium (K)	<0.1 ppm
Propan--1-ol	<0.005 ppm
Sodium (Na)	<0.5 ppm
Strontium (Sr)	<0.005 ppm
Titanium (Ti)	<0.005 ppm
n-Butyl formate	<0.05%
Total phosphorus (P)	<0.05 ppm
Total silicon (Si)	<0.05 ppm
Total sulphur (S)	<0.5 ppm
Vanadium (V)	<0.005 ppm
Zinc (Zn)	<0.02 ppm

BromineCertified AR
for analysis, 99.5+ %, liquid

packaging	Mfr. No
100 mℓ Coated Glass bottle	B/4050/PB05
500 mℓ Coated Glass bottle	B/4050/PB08
Br ₂	UN: 1744
CAS: 7726-95-6	GHS: Acute Tox. 2, Skin
MW: 159.82	Corr. 1A, Eye Dam. 1, Aquatic Acute 1
Density at 20°C: 3.111	H314, H330, H400
Boiling point: 58.7°C	P301+P330+P331, P280,
Melting point: -7.2°C	P305+P351+P338, P310,
ADR/RID: 8, 6.1, I	P304+P340, P273
EINECS: 231-778-1	

Assay: >99.5%	
Guaranteed analysis	
Other organic impurities test	: Pass
Calcium (Ca)	<0.5 ppm
Copper (Cu)	<0.2 ppm
Iodine	<1 ppm
Iron (Fe)	<0.2 ppm
Lead (Pb)	<0.2 ppm
Magnesium (Mg)	<0.5 ppm
Potassium (K)	<0.5 ppm
Sodium (Na)	<1 ppm
n-Butyl formate	<0.05%
Total phosphorus (P)	<0.5 ppm
Total silicon (Si)	<0.5 ppm
Total sulphur (S)	<20 ppm
Zinc (Zn)	<0.2 ppm
Typical analysis	
Aluminium (Al)	<0.2 ppm
Arsenic (As)	<1 ppm
Barium (Ba)	<0.2 ppm
Cadmium (Cd)	<0.2 ppm
Chromium (Cr)	<0.2 ppm
Cobalt (Co)	<0.2 ppm
Lithium (Li)	<0.2 ppm
Manganese (Mn)	<0.2 ppm
Molybdenum (Mo)	<0.2 ppm
Nickel (Ni)	<0.2 ppm
Silver (Ag)	<0.2 ppm
Strontium (Sr)	<0.2 ppm
Titanium (Ti)	<0.2 ppm
Vanadium (V)	<0.2 ppm

Lot analysis on the label

BromineSLR
specified reagents for general laboratory work- extra pure

packaging	Mfr. No
100 mℓ Coated Glass bottle	B/4000/PB05
500 mℓ Coated Glass bottle	B/4000/PB08
Br ₂	UN: 1744
CAS: 7726-95-6	GHS: Acute Tox. 2, Skin
MW: 159.82	Corr. 1A, Eye Dam. 1, Aquatic Acute 1
Density at 20°C: 3.111	H314, H330, H400
Boiling point: 58.7°C	P301+P330+P331, P280,
Melting point: -7.2°C	P305+P351+P338, P310,
ADR/RID: 8, 6.1, I	P304+P340, P273
EINECS: 231-778-1	

Assay: >99%	
Guaranteed analysis	
Calcium (Ca)	<5 ppm
Copper (Cu)	<2 ppm
Iron (Fe)	<2 ppm
Lead (Pb)	<2 ppm
Magnesium (Mg)	<5 ppm
Potassium (K)	<5 ppm
Sodium (Na)	<10 ppm
n-Butyl formate	<0.5%
Total phosphorus (P)	<5 ppm
Total silicon (Si)	<5 ppm
Total sulphur (S)	<200 ppm
Zinc (Zn)	<2 ppm

Bromobenzene
pure

packaging	Mfr. No
250 mℓ Glass bottle	B/4300/07
C ₆ H ₅ Br	ADR/RID: 3, III
CAS: 108-86-1	EINECS: 203-623-8
MW: 157.02	UN: 2514
Density at 20°C: 1.49	GHS: Flam. Liq. 3, Skin Irrit.
Boiling point: 156°C	2, Aquatic Chronic 2
Melting point: -31°C	H226, H411, H315
Flash point: 51°C	P273, P302+P352, P210

Assay (CG): >99%



1-Bromobutane
pure

packaging	Mfr. No
100 mℓ Coated Glass bottle	B/4310/PB05
C ₄ H ₉ Br	EINECS: 203-691-9
CAS: 109-65-9	UN: 1126
MW: 137.03	GHS: Flam. Liq. 2, Skin Irrit.
Density at 20°C: 1.27	2, Eye Irrit. 2, STOT SE 3
Boiling point: 100 to 104°C	H225, H319, H335, H315
Melting point: -112°C	P261, P302+P352, P280,
Flash point: 10°C	P305+P351+P338, P210, P240
ADR/RID: 3, II	
Assay (CG): >98%	

Bromocresol green
pure, indicator grade

packaging	Mfr. No
1 g Glass jar	B/4320/43
5 g Glass jar	B/4320/44
25 g Glass jar	B/4320/46
C ₂₁ H ₁₄ Br ₄ O ₅ S	GHS: Skin Irrit. 2, Eye Irrit. 2,
CAS: 76-60-8	STOT SE 3
MW: 698.01	H335, H315, H319
Melting point: 225°C	P261, P302+P352, P280,
EINECS: 200-972-8	P305+P351+P338
Guaranteed analysis	
pH colour change test	: Pass
Complies with standard IR spectrum: Pass	

Bromocresol green
pure, solution 0.04%, indicator grade

packaging	Mfr. No
500 mℓ Glass bottle	B/4320L/08
ADR/RID: 3/II	GHS: Flam. Liq. 3
EINECS: 200-972-8	H226
UN: 1993	P210
Guaranteed analysis	
pH colour change test	: Pass

Bromocresol purple
pure, solid, indicator grade

packaging	Mfr. No
5 g Glass jar	B/4330/44
25 g Glass jar	B/4330/46
C ₂₁ H ₁₆ Br ₂ O ₅ S	GHS: Skin Irrit. 2, Eye Irrit. 2,
CAS: 115-40-2	STOT SE 3
MW: 540.21	H315, H319, H335
Melting point: 245°C	P261, P302+P352, P280,
EINECS: 204-087-8	P305+P351+P338
Guaranteed analysis	
pH colour change test	: Pass
Complies with standard IR spectrum: Pass	

Bromocresol purple
pure, solution 0.04%, indicator grade

packaging	Mfr. No
500 mℓ Glass bottle	B/4330L/08
Density at 20°C: 0.95	GHS: Flam. Liq. 3, Skin Corr. 1B
Flash point: 32°C	H226, H314
ADR/RID: 3, III	P301+P330+P331, P280,
EINECS: 204-087-8	P305+P351+P338, P310, P210
UN: 1993	

Guaranteed analysis	
pH colour change test	: Pass

Bromoethane
pure

packaging	Mfr. No
250 mℓ Coated Glass bottle	B/4350/PB07
500 mℓ Coated Glass bottle	B/4350/PB08
C ₂ H ₅ Br	EINECS: 200-825-8
CAS: 74-96-4	UN: 1891
MW: 108.97	GHS: Flam. Liq. 2, Carc. 2,
Density at 20°C: 1.46	Acute Tox. 4
Boiling point: 37 to 40°C	H225, H351, H302, H332
Melting point: -119°C	P281, P301+P312,
Flash point: -23°C	P304+P340, P210, P240
ADR/RID: 6.1, II	

Assay (CG): >98%

Bromoform
pure

packaging	Mfr. No
100 mℓ Coated Glass bottle	B/4400/PB05
CHBr ₃	UN: 2515
CAS: 75-25-2	GHS: Acute Tox. 3, Skin Irrit.
MW: 252.73	2, Eye Irrit. 2, Aquatic Chronic 2
Density at 20°C: 2.63	H315, H411, H331, H319, H302
Boiling point: 150 to 151°C	P304+P340, P273,
Melting point: 8°C	P301+P312, P302+P352,
ADR/RID: 6.1, III	P280, P305+P351+P338
EINECS: 200-854-6	

Guaranteed analysis	
Density at 20°C	: 2.58 to 2.68
Ethanol 90 % v/v	3 to 5%

Bromoform
technical, unstabilized

packaging	Mfr. No
250 mℓ Coated Glass bottle	B/4450/PB07
CHBr ₃	UN: 2515
CAS: 75-25-2	GHS: Acute Tox. 3, Skin Irrit.
MW: 252.73	2, Eye Irrit. 2, Aquatic Chronic 2
Density at 20°C: 2.63	H315, H411, H331, H319, H302
Boiling point: 150 to 151°C	P304+P340, P273,
Melting point: 8°C	P301+P312, P302+P352,
ADR/RID: 6.1, III	P280, P305+P351+P338
EINECS: 200-854-6	

Guaranteed analysis	
Density at 20°C	: 2.86 to 2.9

2-Bromo-2-methylpropane
pure

packaging	Mfr. No
100 mℓ Coated Glass bottle	B/4540/PB05
C ₄ H ₉ Br	EINECS: 208-065-9
CAS: 507-19-7	UN: 2342
MW: 137.03	GHS: Flam. Liq. 2, Skin Irrit.
Density at 20°C: 1.18	2, Eye Irrit. 2, STOT SE 3
Boiling point: 72 to 74°C	H225, H319, H335, H315
Melting point: -20°C	P261, P302+P352, P280,
Flash point: 18°C	P305+P351+P338, P210, P240
ADR/RID: 3, II	

Assay (CG): >97%



1-Bromonaphthalene

specified reagents for general laboratory work- extra pure

SLR

packaging

100 mℓ Coated Glass bottle

Mfr. No

B/4550/PB05

C₁₀H₇Br

Flash point: 110°C

CAS: 90-11-9

EINECS: 201-965-2

MW: 207.08

GHS: Acute Tox. 4

Density at 20°C: 1.489

H302

Boiling point: 279 to 281°C

P301+P312

Melting point: -1°C

Absorbance specification

Max. Absorbance (1cm cell v HPLC water)

Refractive index

1.657 to 1.66

!

Bromophenol red

pure, indicator grade

packaging

1 g Glass jar

Mfr. No

B/4660/43

C₁₉H₁₂Br₂O₃S

CAS: 2800-80-8

MW: 512.19

EINECS: 220-538-1

Guaranteed analysis

pH colour change test : Pass

Bromothymol blue

pure, indicator grade

packaging

5 g Glass jar

Mfr. No

B/4690/44

C₂₇H₂₈Br₂O₃S

CAS: 76-59-5

MW: 624.38

Melting point: 204°C

EINECS: 200-971-2

Guaranteed analysis

Complies with standard IR spectrum: Pass

Bromothymol blue

pure, solution 0.04%, indicator grade

packaging

500 mℓ Glass bottle

Mfr. No

B/4690L/08

Density at 20°C: 0.8

UN: 1170

Boiling point: 63.9°C

GHS: Flam. Liq. 2

Flash point: 11.1°C

H225

ADR/RID: 3, III

P210, P240

EINECS: 200-971-2

Guaranteed analysis

pH colour change test : Pass

🔥

Bromure standard 1000 ppm

For Ion chromatography in H₂O, NIST certified

packaging

100 mℓ Plastic bottle

Mfr. No


SPASBR92Y

500 mℓ Plastic bottle

SPASBR92X

Manufactured to a tolerance of ± 0.3%.

Check with your local distributor for availability of Spex product



Bromophenol blue

pure, indicator grade

packaging

5 g Glass jar

Mfr. No

B/4630/44

25 g Glass jar

B/4630/46

C₁₉H₁₀Br₄O₃S

Melting point: 279°C

CAS: 115-39-9

EINECS: 204-086-2

MW: 669.96

Guaranteed analysis

pH colour change test : Pass

Complies with standard IR spectrum: Pass

Bromophenol blue

pure, solution 0.2% in Methylated spirits

packaging

500 mℓ Glass bottle

Mfr. No

B/4631L/08

CAS: 115-39-9

EINECS: 204-086-2

Density at 20°C: 0.8

UN: 1170

Boiling point: 63.9°C

GHS: Flam. Liq. 2

Flash point: 11.1°C

H225

ADR/RID: 3, II

P210, P240

Guaranteed analysis

pH colour change test : Pass

🔥

Bromophenol blue

pure, solution 0.04%, indicator grade

packaging

100 mℓ Glass bottle

Mfr. No

B/4630L/05

500 mℓ Glass bottle

B/4630L/08

CAS: 115-39-9

GHS: Flam. Liq. 3

Flash point: 32°C

H226

ADR/RID: 3, III

P210, P280, P233, P240,

UN: 1170

P303+P361+P353

Guaranteed analysis

pH colour change test : Pass

🔥

NIST Buffer Standard solution ready to use for pH measurement					
Manufactured to a tolerance of ± 0.02 pH Units at 20°C					
pH 2.00 to 7.00 buffer solution are stabilized with 10 ppm of Mercury (II) chloride					
designation	packsize	package type	Mfr. No	page	
Buffer solution pH 2.00 (glycine)	500 mℓ	HDPE Plastic bottle	J/2800/08	65	
	1 L	HDPE Plastic bottle	J/2800/15	65	
	2.5 L	HDPE Plastic bottle	J/2800/17	65	
Buffer solution pH 3.00 (phthalate)	500 mℓ	HDPE Plastic bottle	J/2810/08	65	
	1 L	HDPE Plastic bottle	J/2810/15	65	
	2.5 L	HDPE Plastic bottle	J/2810/17	65	
Buffer solution pH 4.00 (phthalate)	500 mℓ	HDPE Plastic bottle	J/2820/08	65	
	1 L	HDPE Plastic bottle	J/2820/15	65	
	2.5 L	HDPE Plastic bottle	J/2820/17	65	
	5 L	HDPE Plastic bottle	J/2820/21	65	
	10 L	Bag in box	J/2820/24	65	
Buffer solution pH 4.00 BS formulated	1 L	HDPE Plastic bottle	J/2821/15	65	
Buffer solution pH 5.00 (phthalate)	500 mℓ	HDPE Plastic bottle	J/2830/08	65	
	1 L	HDPE Plastic bottle	J/2830/15	65	
	2.5 L	HDPE Plastic bottle	J/2830/17	65	
Buffer solution pH 6.00 (phosphate)	500 mℓ	HDPE Plastic bottle	J/2840/08	65	
	1 L	HDPE Plastic bottle	J/2840/15	65	
	2.5 L	HDPE Plastic bottle	J/2840/17	65	
Buffer solution pH 6.87 BS formulated	1 L	HDPE Plastic bottle	J/2845/15	65	
Buffer solution pH 7.00 (phosphate)	500 mℓ	HDPE Plastic bottle	J/2850/08	66	
	1 L	HDPE Plastic bottle	J/2850/15	66	
	2.5 L	HDPE Plastic bottle	J/2850/17	66	
	5 L	HDPE Plastic bottle	J/2850/21	66	
	10 L	Bag in box	J/2850/24	66	
Buffer solution pH 8.00 (borate)	500 mℓ	HDPE Plastic bottle	J/2860/08	66	
	1 L	HDPE Plastic bottle	J/2860/15	66	
	2.5 L	HDPE Plastic bottle	J/2860/17	66	
	10 L	Bag in box	J/2860/24	66	
Buffer solution pH 9.00 (borate)	2.5 L	HDPE Plastic bottle	J/2865/17	66	
Buffer solution pH 9.2 (borate)	500 mℓ	HDPE Plastic bottle	J/2870/08	66	
	1 L	HDPE Plastic bottle	J/2870/15	66	
	2.5 L	HDPE Plastic bottle	J/2870/17	66	
	5 L	HDPE Plastic bottle	J/2870/21	66	
	10 L	Bag in box	J/2870/24	66	
Buffer solution pH 9.225 BS formulated	1 L	HDPE Plastic bottle	J/2871/15	66	
Buffer solution pH 9.23 (borate)	1 L	HDPE Plastic bottle	J/2875/15	66	
Buffer solution pH 10.00 (borate)	500 mℓ	HDPE Plastic bottle	J/2880/08	66	
	1 L	HDPE Plastic bottle	J/2880/15	66	
	2.5 L	HDPE Plastic bottle	J/2880/17	66	
	5 L	HDPE Plastic bottle	J/2880/21	66	

NIST Buffer colour coded standard solution ready to use for pH measurement					
Manufactured to a tolerance of ± 0.02 pH Units at 20°C					
pH 4.00 & 7.00 buffer solution are stabilized with 10 ppm of Mercury (II) chloride					
designation	packsize	package type	Mfr. No	page	
Buffer colour coded solution pH 4.00 (phthalate) Red	500 mℓ	HDPE Plastic bottle	J/2825/08	65	
	1 L	HDPE Plastic bottle	J/2825/15	65	
	2.5 L	HDPE Plastic bottle	J/2825/17	65	
	5 L	HDPE Plastic bottle	J/2825/21	65	
Buffer colour coded solution pH 7.00 (phosphate) Yellow	10 L	Bag in box	J/2825/24	65	
	500 mℓ	HDPE Plastic bottle	J/2855/08	66	
	1 L	HDPE Plastic bottle	J/2855/15	66	
	2.5 L	HDPE Plastic bottle	J/2855/17	66	
	5 L	HDPE Plastic bottle	J/2855/21	66	
Buffer colour coded solution pH 10.00 (borate) Blue	10 L	Bag in box	J/2855/24	66	
	500 mℓ	HDPE Plastic bottle	J/2885/08	67	
	1 L	HDPE Plastic bottle	J/2885/15	67	
	2.5 L	HDPE Plastic bottle	J/2885/17	67	
	5 L	HDPE Plastic bottle	J/2885/21	67	
	10 L	Bag in box	J/2885/24	67	

NIST Buffer concentrated solution for pH measurement

Each ampoule makes up 500mℓ of solution
pH 2.00 to 7.00 buffer solution are stabilized with 10 ppm of Mercury (II) chloride

designation	packsize	package type		Mfr. No	page
Buffer concentrated solution pH 2.00 (glycine)	1	Ampoule		J/2800C/05	66
	6	Ampoule		J/2800C/90	66
Buffer concentrated solution pH 4.00 (phthalate)	6	Ampoule		J/2820C/90	66
Buffer concentrated solution pH 7.00 (phosphate)	1	Ampoule		J/2850C/05	67
	6	Ampoule		J/2850C/90	67
Buffer concentrated solution pH 8.00 (borate)	1	Ampoule		J/2860C/05	67
	6	Ampoule		J/2860C/90	67
Buffer concentrated solution pH 9.2 (borate)	1	Ampoule		J/2870C/05	67
	6	Ampoule		J/2870C/90	67
Buffer concentrated solution pH 10.00 (borate)	6	Ampoule		J/2880C/90	68

NIST Buffer colour coded concentrated solution for pH measurement

Each ampoule makes up 500mℓ of solution
Stabilized with 10 ppm of Mercury (II) chloride (except J/2885C)

designation	packsize	package type		Mfr. No	page
Buffer colour coded concentrated solution pH 4.00 (phthalate) Red	1 ampoule	Ampoule		J/2825C/05	65
	6 ampoule	Ampoule		J/2825C/90	65
Buffer colour coded concentrated solution pH 7.00 (phosphate) Yellow	1 ampoule	Ampoule		J/2855C/05	66
	6 ampoule	Ampoule		J/2855C/90	66
Buffer colour coded concentrated solution pH 10.00 (borate) Blue	1 ampoule	Ampoule		J/2885C/05	67
	6 ampoule	Ampoule		J/2885C/90	67

Buffer Tablets for pH measurement

Each ampoule makes up 100mℓ of solution

designation	packsize	package type		Mfr. No	page
Buffer tablets, pH 4 (phthalate)	50 PC	Glass jar		B/4760/77	65
Buffer tablets, pH 7 (phosphate)	50 PC	Glass jar		B/4765/77	66
Buffer tablets, pH 9.2 (borate)	50 PC	Glass jar		B/4770/77	66

Buffer solution pH 2.00 (glycine)
NIST Standard solution ready to use for pH measurement

packaging	Mfr. No
500 mℓ HDPE Plastic bottle	J/2800/08
1 ℓ HDPE Plastic bottle	J/2800/15
2.5 ℓ HDPE Plastic bottle	J/2800/17

Manufactured to a tolerance of ± 0.02 pH Units at 20°C
Stabilised with 10 ppm of Mercury (II) chloride

Buffer concentrated solution pH 2.00 (glycine)
NIST Standard for pH measurement

packaging	Mfr. No
1 Ampoule	J/2800C/05
6 Ampoule	J/2800C/90

GHS: Skin Corr. 1B, Eye Dam. 1

H314

P301+P330+P331, P280, P305+P351+P338, P310



Each ampoule makes up 500mℓ of solution.
Stabilised with 10 ppm of Mercury (II) chloride

Buffer solution pH 3.00 (phthalate)
NIST Standard solution ready to use for pH measurement

packaging	Mfr. No
500 mℓ HDPE Plastic bottle	J/2810/08
1 ℓ HDPE Plastic bottle	J/2810/15
2.5 ℓ HDPE Plastic bottle	J/2810/17

Manufactured to a tolerance of ± 0.02 pH Units at 20°C
Stabilised with 10 ppm of Mercury (II) chloride

Buffer solution pH 4.00 (phthalate)
NIST Standard solution ready to use for pH measurement

packaging	Mfr. No
500 mℓ HDPE Plastic bottle	J/2820/08
1 ℓ HDPE Plastic bottle	J/2820/15
2.5 ℓ HDPE Plastic bottle	J/2820/17
5 ℓ HDPE Plastic bottle	J/2820/21
10 ℓ Bag in box	J/2820/24

Manufactured to a tolerance of ± 0.02 pH Units at 20°C
Stabilised with 10 ppm of Mercury (II) chloride

Buffer concentrated solution pH 4 (phthalate)

packaging	Mfr. No
6 Ampoule	J/2820C/90

Buffer colour coded solution pH 4.00 (phthalate) Red
NIST Standard solution ready to use for pH measurement

packaging	Mfr. No
500 mℓ HDPE Plastic bottle	J/2825/08
1 ℓ HDPE Plastic bottle	J/2825/15
2.5 ℓ HDPE Plastic bottle	J/2825/17
5 ℓ HDPE Plastic bottle	J/2825/21
10 ℓ Bag in box	J/2825/24

Manufactured to a tolerance of ± 0.02 pH Units at 20°C
Stabilised with 10 ppm of Mercury (II) chloride

Buffer colour coded concentrated solution pH 4.00 (phthalate) Red
NIST Standard for pH measurement

packaging	Mfr. No
1 Ampoule	J/2825C/05
6 Ampoule	J/2825C/90

Each ampoule makes up 500mℓ of solution.
Stabilised with 10 ppm of Mercury (II) chloride.

Buffer solution pH 4.00 BS formulated
Standard solution ready to use for pH measurement

packaging	Mfr. No
1 ℓ HDPE Plastic bottle	J/2821/15

See the full range of Standard Volumetric solution page 323.

Buffer concentrated solution BS pH 4.00 (phtalate)
Standard for pH measurement

packaging	Mfr. No
1 Ampoule	J/2821C/05

Each ampoule makes up 500mℓ of solution
Stabilised with 10ppm mercury chloride

Buffer tablets, pH 4 (phthalate)
for pH measurement

packaging	Mfr. No
50 PC Glass jar	B/4760/77

Guaranteed analysis
pH of 10 tablets made up to 1000mℓ : pH of tablet made up to 100mℓ.: 3.95 to 4.05
3.95 to 4.05

Each tablet makes up 100mℓ of solution.

Buffer solution pH 5.00 (phthalate)
NIST Standard solution ready to use for pH measurement

packaging	Mfr. No
500 mℓ HDPE Plastic bottle	J/2830/08
1 ℓ HDPE Plastic bottle	J/2830/15
2.5 ℓ HDPE Plastic bottle	J/2830/17

EINECS: 212-889-4

Manufactured to a tolerance of ± 0.02 pH Units at 20°C
Stabilised with 10 ppm of Mercury (II) chloride

Buffer solution pH 6.00 (phosphate)
NIST Standard solution ready to use for pH measurement

packaging	Mfr. No
500 mℓ HDPE Plastic bottle	J/2840/08
1 ℓ HDPE Plastic bottle	J/2840/15
2.5 ℓ HDPE Plastic bottle	J/2840/17

EINECS: 231-913-4

Manufactured to a tolerance of ± 0.02 pH Units at 20°C
Stabilised with 10 ppm of Mercury (II) chloride

Buffer solution pH 6.87 BS formulated
Standard solution ready to use for pH measurement

packaging	Mfr. No
1 ℓ HDPE Plastic bottle	J/2845/15

Manufactured to a tolerance of ± 0.02 pH Units at 20°C
Stabilised with 10 ppm of Mercury (II) chloride

Buffer concentrated solution BS pH 6.87 (phosphate)
Standard for pH measurement

packaging	Mfr. No
6 Ampoule	J/2845C/90

Each ampoule makes up 500mℓ of solution
Stabilised with 10ppm mercury chloride

Buffer solution pH 7.00 (phosphate)

NIST Standard solution ready to use for pH measurement

packaging	Mfr. No
500 mℓ HDPE Plastic bottle	J/2850/08
1 ℓ HDPE Plastic bottle	J/2850/15
2.5 ℓ HDPE Plastic bottle	J/2850/17
5 ℓ HDPE Plastic bottle	J/2850/21
10 ℓ Bag in box	J/2850/24

Manufactured to a tolerance of ± 0.02 pH Units at 20°C
Stabilised with 10 ppm of Mercury (II) chloride

Buffer concentrated solution pH 7.00 (phosphate)

NIST Standard for pH measurement

packaging	Mfr. No
1 Ampoule	J/2850C/05
6 Ampoule	J/2850C/90

Each ampoule makes up 500 mℓ of solution.
Stabilised with 10 ppm of Mercury (II) Chloride.

Buffer colour coded solution pH 7.00 (phosphate) Yellow

NIST Standard solution ready to use for pH measurement

packaging	Mfr. No
500 mℓ HDPE Plastic bottle	J/2855/08
1 ℓ HDPE Plastic bottle	J/2855/15
2.5 ℓ HDPE Plastic bottle	J/2855/17
5 ℓ HDPE Plastic bottle	J/2855/21
10 ℓ Bag in box	J/2855/24

Manufactured to a tolerance of ± 0.02 pH Units at 20°C
Stabilised with 10 ppm of Mercury (II) chloride

Buffer colour coded concentrated solution pH 7.00 (phosphate) Yellow

NIST Standard for pH measurement

packaging	Mfr. No
1 Ampoule	J/2855C/05
6 Ampoule	J/2855C/90

Each ampoule makes up 500mℓ of solution.
Stabilised with 10 ppm of Mercury (II) chloride.

Buffer tablets, pH 7 (phosphate)

for pH measurement

packaging	Mfr. No
50 PC Glass jar	B/4765/77

Guaranteed analysis

pH of 10 tablets made up to 1000mℓ : 6.95 to 7.05

Each tablet makes up 100mℓ of solution.

Buffer solution pH 8.00 (borate)

NIST Standard solution ready to use for pH measurement

packaging	Mfr. No
500 mℓ HDPE Plastic bottle	J/2860/08
1 ℓ HDPE Plastic bottle	J/2860/15
2.5 ℓ HDPE Plastic bottle	J/2860/17
10 ℓ Bag in box	J/2860/24

Manufactured to a tolerance of ± 0.02 pH Units at 20°C

Buffer concentrated solution pH 8.00 (borate)

NIST Standard for pH measurement

packaging	Mfr. No
1 Ampoule	J/2860C/05
6 Ampoule	J/2860C/90

Each ampoule makes up 500mℓ of solution.

Buffer solution pH 9.00 (borate)

NIST Standard solution ready to use for pH measurement

packaging	Mfr. No
2.5 ℓ HDPE Plastic bottle	J/2865/17

Manufactured to a tolerance of ± 0.02 pH Units at 20°C

Buffer solution pH 9.2 (borate)

NIST Standard solution ready to use for pH measurement

packaging	Mfr. No
500 mℓ HDPE Plastic bottle	J/2870/08
1 ℓ HDPE Plastic bottle	J/2870/15
2.5 ℓ HDPE Plastic bottle	J/2870/17
5 ℓ HDPE Plastic bottle	J/2870/21
10 ℓ Bag in box	J/2870/24

Manufactured to a tolerance of ± 0.02 pH Units at 20°C

Buffer concentrated solution pH 9.2 (borate)

NIST Standard for pH measurement

packaging	Mfr. No
1 Ampoule	J/2870C/05
6 Ampoule	J/2870C/90

Each ampoule makes up 500mℓ of solution.

Buffer solution pH 9.225 BS formulated

Standard solution ready to use for pH measurement

packaging	Mfr. No
1 ℓ HDPE Plastic bottle	J/2871/15

Manufactured to a tolerance of ± 0.02pH Units at 20°C
Stabilised with 10ppm mercury chloride

Buffer tablets, pH 9.2 (borate)

for pH measurement

packaging	Mfr. No
50 PC Glass jar	B/4770/77

GHS: Repr. 1B

H360FD

P201, P281, P308+P313

Guaranteed analysis

pH of tablet made up to 100mℓ.: 9.15 to 9.25

Each tablet makes up 100mℓ of solution.

Buffer solution pH 9.23 (borate)

NIST Standard solution ready to use for pH measurement

packaging	Mfr. No
1 ℓ HDPE Plastic bottle	J/2875/15
H360FD	

Manufactured to a tolerance of ± 0.02 pH Units at 20°C

Buffer solution pH 10.00 (borate)

NIST Standard solution ready to use for pH measurement

packaging	Mfr. No
500 mℓ HDPE Plastic bottle	J/2880/08
1 ℓ HDPE Plastic bottle	J/2880/15
2.5 ℓ HDPE Plastic bottle	J/2880/17
5 ℓ HDPE Plastic bottle	J/2880/21

Manufactured to a tolerance of ± 0.02 pH Units at 20°C

Buffer concentrated solution pH 10 (borate)

packaging	Mfr. No
6 Ampoule	J/2880C/90

GHS: Acute Tox. 4, Eye Irrit. 2

H332, H319

P304+P340, P280, P305+P351+P338

Buffer colour coded solution pH 10.00 (borate) Blue

NIST Standard solution ready to use for pH measurement

packaging	Mfr. No
500 mℓ HDPE Plastic bottle	J/2885/08
1 ℓ HDPE Plastic bottle	J/2885/15
2.5 ℓ HDPE Plastic bottle	J/2885/17
5 ℓ HDPE Plastic bottle	J/2885/21
10 ℓ Bag in box	J/2885/24

Manufactured to a tolerance of ± 0.02 pH Units at 20°C

Buffer colour coded concentrated solution pH 10.00 (borate) Blue

NIST Standard for pH measurement

packaging	Mfr. No
1 Ampoule	J/2885C/05
6 Ampoule	J/2885C/90

GHS: Skin Irrit. 2, Eye Irrit. 2

H315, H319

P280, P305+P351+P338

Each ampoule makes up 500 mℓ of solution.

Full range indicator, pH range 1-13

for pH measurement

packaging	Mfr. No
100 mℓ Coated Glass bottle	U/0025/PB05
500 mℓ Coated Glass bottle	U/0025/PB08

Density at 20°C: 0.8

Flash point: 12°C

ADR/RID: 3, II

EINECS: 200-578-6

UN: 1170

GHS: Flam. Liq. 2, Acute

Tox. 4, STOT SE 2

H225, H371, H312, H332, H302

P280, P302+P352, P260,

P304+P340, P210, P240

Guaranteed analysis

pH colour change test : Pass

Universal indicator, pH range 4-10

for pH measurement

packaging	Mfr. No
100 mℓ Coated Glass bottle	U/0020/PB05
500 mℓ Coated Glass bottle	U/0020/PB08
2.5 ℓ Coated Glass bottle	U/0020/PB17

Flash point: 12°C

ADR/RID: 3, II

EINECS: 200-578-6

UN: 1170

GHS: Flam. Liq. 2, Acute

Tox. 4, STOT SE 2

H225, H371, H312, H332, H302

P280, P302+P352, P260,

P304+P340, P210, P240

Buffer solution ammonia-ammonium chloride

pure, For water hardness testing

packaging	Mfr. No
1 ℓ HDPE Plastic bottle	J/2500/PB15
2.5 ℓ HDPE Plastic bottle	J/2500/17

Density at 20°C: 0.93

ADR/RID: 8, III

EINECS: 235-186-4

UN: 2672

GHS: Skin Corr. 1B, STOT SE 3

H314, H335

P280, P301+P330+P331,

P310, P305+P351+P338,

P304+P340, P302+P352

1-Butane sulfonic acid sodium salt

for HPLC

Certified

packaging	Mfr. No
25 g Plastic jar	B/4791/46

CH₃(CH₂)₃SO₃Na

CAS: 2386-54-1

MW: 160.16

Melting point: >300°C

EINECS: 219-201-1

GHS: Skin Irrit. 2, Eye Irrit. 2,

STOT SE 3

H315, H335, H319

P261, P302+P352, P280,

P305+P351+P338

Assay: >98%

Lot analysis on the label.

Filtered to 0.2 µm.

iso-Butyric acid

pure

packaging	Mfr. No
500 mℓ Coated Glass bottle	B/5550/PB08

C₄H₈O₂

CAS: 79-31-2

MW: 88.11

Density at 20°C: 0.950

Boiling point: 153 to 154°C

Melting point: -47°C

Flash point: 59°C

ADR/RID: 3, 8, III

EINECS: 201-195-7

UN: 2529

GHS: Acute Tox. 4

H302, H312

P301+P310, P280,

P302+P350, P210

Assay: >99%

Butan-1-ol

For chromatography, for the determination of 17-Ketosteroids

packaging	Mfr. No
500 mℓ Coated Glass bottle	B/4860/PB08

C₄H₁₀O

CAS: 71-36-3

MW: 74.12

Density at 20°C: 0.810

Boiling point: 117.6°C

Melting point: -89°C

Flash point: 35°C

ADR/RID: 3, III

EINECS: 200-751-6

UN: 1120

GHS: Flam. Liq. 3, Acute

Tox. 4, Skin Irrit. 2, Eye

Dam. 1, STOT SE 3

H226, H315, H302, H336,

H318, H335

P280, P305+P351+P338,

P301+P312, P304+P340,

P210, P303+P361+P353

Assay (CG): >99%

Safe, environmentally
friendly and effective
packaging

Fisher Chemical outer packaging
is UN compliant and 100% recyclable.
It's an easy and safe way to reduce
your waste management costs and
contribute to a cleaner environment.



Butan-1-ol
for analysis

Certified AR

packaging		Mfr. No
500 mℓ	Glass bottle	B/4850/08
1 ℓ	Glass bottle	B/4850/15
2.5 ℓ	Glass bottle	B/4850/17
2.5 ℓ	Coated Glass bottle	B/4850/PB17
25 ℓ	Metal drum	B/4850/25

C ₄ H ₁₀ O	UN: 1120	
CAS: 71-36-3	GHS: Flam. Liq. 3, Acute	
MW: 74.12	Tox. 4, Skin Irrit. 2, Eye	
Density at 20°C: 0.810	Dam. 1, STOT SE 3	
Boiling point: 117.6°C	H226, H315, H302, H336,	
Melting point: -89°C	H318, H335	
Flash point: 35°C	P280, P305+P351+P338,	
ADR/RID: 3, III	P301+P312, P304+P340,	
EINECS: 200-751-6	P210, P303+P361+P353	

Assay (CG): >99.5%

Guaranteed analysis			
Acidity	<0.0005 meq/g	Magnesium (Mg)	<0.1 ppm
Colour	<10 APHA	Potassium (K)	<0.1 ppm
Water	<0.1%	Sodium (Na)	<1 ppm
Residue after evaporation	<10 ppm	Total phosphorus (P)	<0.2 ppm
Calcium (Ca)	<0.5 ppm	Total silicon (Si)	<0.05 ppm
Copper (Cu)	<0.02 ppm	Total sulphur (S)	<0.5 ppm
Iron (Fe)	<0.1 ppm	Zinc (Zn)	<0.2 ppm
Lead (Pb)	<0.05 ppm	Butan -2-ol	<0.05%

Typical analysis			
Aluminium (Al)	<0.2 ppm	Manganese (Mn)	<0.02 ppm
Barium (Ba)	<0.1 ppm	Molybdenum (Mo)	<0.02 ppm
Cadmium (Cd)	<0.05 ppm	Nickel (Ni)	<0.02 ppm
Chromium (Cr)	<0.02 ppm	Strontium (Sr)	<0.05 ppm
Cobalt (Co)	<0.02 ppm	Titanium (Ti)	<0.02 ppm
Lithium (Li)	<0.02 ppm	Vanadium (V)	<0.02 ppm

Lot analysis on the label.

Butan-1-ol
specified reagents for general laboratory work- extra pure

SLR

packaging		Mfr. No
500 mℓ	Glass bottle	B/4800/08
1 ℓ	Glass bottle	B/4800/15
1 ℓ	Metal can	B/4800/MC15
2.5 ℓ	Glass bottle	B/4800/17
2.5 ℓ	Coated Glass bottle	B/4800/PB17
5 ℓ	Fluorinated HDPE Plastic bottle	B/4800/21
25 ℓ	Metal drum	B/4800/25

C ₄ H ₁₀ O	UN: 1120	
CAS: 71-36-3	GHS: Flam. Liq. 3, Acute	
MW: 74.12	Tox. 4, Skin Irrit. 2, Eye	
Density at 20°C: 0.810	Dam. 1, STOT SE 3	
Boiling point: 117.6°C	H226, H315, H302, H336,	
Melting point: -89°C	H318, H335	
Flash point: 35°C	P280, P305+P351+P338,	
ADR/RID: 3, III	P301+P312, P304+P340,	
EINECS: 200-751-6	P210, P303+P361+P353	

Assay (CG): >99%

Guaranteed analysis			
Colour	<20 APHA	Magnesium (Mg)	<2 ppm
Water	<0.5%	Potassium (K)	<2 ppm
Residue after evaporation	<50 ppm	Sodium (Na)	<5 ppm
Acidity	<0.005 meq/g	Total phosphorus (P)	<2 ppm
Calcium (Ca)	<5 ppm	Total silicon (Si)	<2 ppm
Copper (Cu)	<2 ppm	Total sulphur (S)	<2 ppm
Iron (Fe)	<5 ppm	Zinc (Zn)	<2 ppm
Lead (Pb)	<5 ppm		

sec-Butanol (Butan-2-ol)
for analysis

Certified AR

packaging		Mfr. No
1 ℓ	Glass bottle	B/5210/15
2.5 ℓ	Coated Glass bottle	B/5210/PB17

C ₄ H ₁₀ O	EINECS: 201-158-5	
CAS: 78-92-2	UN: 1120	
MW: 74.12	GHS: Flam. Liq. 3, Eye Irrit.	
Density at 20°C: 0.800	2, STOT SE 3	
Boiling point: 99°C	H226, H319, H335, H336	
Melting point: -115°C	P261, P280,	
Flash point: 24°C	P305+P351+P338, P210	
ADR/RID: 3, III		

Assay (CG): >99.5%

Guaranteed analysis			
Acidity / Alkalinity	<0.0005 meq/g	Potassium (K)	<0.1 ppm
Colour	<10 APHA	Sodium (Na)	<0.5 ppm
Water	<0.2%	Total phosphorus (P)	<0.1 ppm
Residue after evaporation	<20 ppm	Total silicon (Si)	<0.02 ppm
2-Methylpropan-2-ol	: <0.05%	Total sulphur (S)	<1 ppm
Calcium (Ca)	<0.2 ppm	Zinc (Zn)	<0.1 ppm
Copper (Cu)	<0.02 ppm	Butan--2-one	<0.05%
Iron (Fe)	<0.1 ppm	Di-n-butyl ether	<0.05%
Lead (Pb)	<0.05 ppm	Propan--2-ol	<0.1%
Magnesium (Mg)	<0.05 ppm		

Typical analysis			
Manganese (Mn)	<0.01 ppm	Lithium (Li)	<0.01 ppm
Aluminium (Al)	<0.05 ppm	Molybdenum (Mo)	<0.01 ppm
Barium (Ba)	<0.02 ppm	Nickel (Ni)	<0.02 ppm
Cadmium (Cd)	<0.01 ppm	Strontium (Sr)	<0.01 ppm
Chromium (Cr)	<0.01 ppm	Titanium (Ti)	<0.01 ppm
Cobalt (Co)	<0.01 ppm	Vanadium (V)	<0.01 ppm

sec-Butanol (Butan-2-ol)
specified reagents for general laboratory work- extra pure

SLR

packaging		Mfr. No
1 ℓ	Glass bottle	B/5200/15
2.5 ℓ	Coated Glass bottle	B/5200/PB17
25 ℓ	Metal drum	B/5200/25

C ₄ H ₁₀ O	EINECS: 201-158-5	
CAS: 78-92-2	UN: 1120	
MW: 74.12	GHS: Flam. Liq. 3, Eye Irrit.	
Density at 20°C: 0.800	2, STOT SE 3	
Boiling point: 99°C	H226, H319, H335, H336	
Melting point: -115°C	P261, P280,	
Flash point: 24°C	P305+P351+P338, P210	
ADR/RID: 3, III		

Assay (CG): >99%

Guaranteed analysis			
Acidity / Alkalinity	<0.001 meq/g	Magnesium (Mg)	<2 ppm
Colour	<20 APHA	Potassium (K)	<2 ppm
Water	<0.5%	Sodium (Na)	<5 ppm
Residue after evaporation	<50 ppm	Total phosphorus (P)	<2 ppm
Calcium (Ca)	<5 ppm	Total silicon (Si)	<2 ppm
Copper (Cu)	<1 ppm	Total sulphur (S)	<5 ppm
Iron (Fe)	<2 ppm	Zinc (Zn)	<2 ppm
Lead (Pb)	<1 ppm		

iso-Butanol
for analysis

Certified AR

packaging		Mfr. No
500 mℓ	Glass bottle	B/5150/08
2.5 ℓ	Coated Glass bottle	B/5150/PB17
25 ℓ	Metal drum	B/5150/25

C ₄ H ₁₀ O	EINECS: 201-148-0	
CAS: 78-83-1	UN: 1212	
MW: 74.12	GHS: Flam. Liq. 3, Skin Irrit.	
Density at 20°C: 0.800	2, Eye Dam. 1, STOT SE 3	
Boiling point: 108°C	H226, H318, H336, H315, H335	
Melting point: -108°C	P280, P305+P351+P338,	
Flash point: 28°C	P310, P304+P340, P261,	
ADR/RID: 3, III	P302+P352, P210	



Assay (CG): >99.5%			
Guaranteed analysis			
Acidity / Alkalinity	<0.005 meq/g	Potassium (K)	<0.1 ppm
Colour	<10 APHA	Sodium (Na)	<1 ppm
Water	<0.2%	Total phosphorus (P)	<0.1 ppm
Residue after evaporation	<20 ppm	Total silicon (Si)	<0.02 ppm
2-Methylpropan-2-ol	: <0.5%	Total sulphur (S)	<1 ppm
Calcium (Ca)	<0.5 ppm	Zinc (Zn)	<0.2 ppm
Copper (Cu)	<0.1 ppm	Butan -1-ol	<0.5%
Iron (Fe)	<0.2 ppm	Butan -2-ol	<0.5%
Lead (Pb)	<0.2 ppm	Butan--2-one	<0.1%
Magnesium (Mg)	<0.1 ppm	Butyraldehyde	<0.2%

Typical analysis			
Aluminium (Al)	<0.05 ppm	Manganese (Mn)	<0.02 ppm
Barium (Ba)	<0.02 ppm	Molybdenum (Mo)	<0.02 ppm
Cadmium (Cd)	<0.02 ppm	Nickel (Ni)	<0.02 ppm
Chromium (Cr)	<0.02 ppm	Strontium (Sr)	<0.02 ppm
Cobalt (Co)	<0.02 ppm	Titanium (Ti)	<0.02 ppm
Lithium (Li)	<0.02 ppm	Vanadium (V)	<0.02 ppm

Lot analysis on the label.

iso-Butanol
specified reagents for general laboratory work- extra pure

SLR

packaging		Mfr. No
2.5 ℓ	Coated Glass bottle	B/5100/PB17
25 ℓ	Metal drum	B/5100/25
200 ℓ	Metal drum	B/5100/27

C ₄ H ₁₀ O	EINECS: 201-148-0	
CAS: 78-83-1	UN: 1212	
MW: 74.12	GHS: Flam. Liq. 3, Skin Irrit.	
Density at 20°C: 0.800	2, Eye Dam. 1, STOT SE 3	
Boiling point: 108°C	H226, H318, H336, H315, H335	
Melting point: -108°C	P280, P305+P351+P338,	
Flash point: 28°C	P310, P304+P340, P261,	
ADR/RID: 3, III	P302+P352, P210	

Assay (CG): >99%

Guaranteed analysis			
Acidity / Alkalinity	<0.05 meq/g	Magnesium (Mg)	<2 ppm
Colour	<20 APHA	Potassium (K)	<2 ppm
Water	<0.5%	Sodium (Na)	<5 ppm
Residue after evaporation	<50 ppm	Total phosphorus (P)	<2 ppm
Calcium (Ca)	<5 ppm	Total silicon (Si)	<2 ppm
Copper (Cu)	<1 ppm	Total sulphur (S)	<5 ppm
Iron (Fe)	<2 ppm	Zinc (Zn)	<2 ppm
Lead (Pb)	<1 ppm		

tert-Butanol
for analysis

Certified AR

packaging		Mfr. No
500 mℓ	HDPE Plastic bottle	B/5251/08
2.5 ℓ	HDPE Plastic bottle	B/5251/17

C ₄ H ₁₀ O	EINECS: 200-889-7	
CAS: 75-65-0	UN: 1120	
MW: 74.12	GHS: Flam. Liq. 2, Acute	
Density at 20°C: 0.780	Tox. 4, Eye Irrit. 2, STOT SE 3	
Boiling point: 83°C	H225, H319, H335, H332	
Melting point: 25 to 25.5°C	P261, P304+P340, P280,	
Flash point: 11°C	P305+P351+P338, P210, P240	
ADR/RID: 3, II		

Assay (CG): >99.5%

Guaranteed analysis			
Acidity / Alkalinity	<0.005 meq/g	Potassium (K)	<0.1 ppm
Colour	<10 APHA	Sodium (Na)	<2 ppm
Water	<0.1%	Total phosphorus (P)	<0.1 ppm
Residue after evaporation	<10 ppm	Total silicon (Si)	<0.2 ppm
Calcium (Ca)	<0.5 ppm	Total sulphur (S)	<1 ppm
Copper (Cu)	<0.2 ppm	Zinc (Zn)	<0.2 ppm
Iron (Fe)	<0.2 ppm	Butan -2-ol	<0.2%
Lead (Pb)	<0.1 ppm	Butan--2-one	<0.01%
Magnesium (Mg)	<0.1 ppm	Butyraldehyde	<0.01%

Typical analysis			
Aluminium (Al)	<0.1 ppm	Manganese (Mn)	<0.02 ppm
Barium (Ba)	<0.02 ppm	Molybdenum (Mo)	<0.02 ppm
Cadmium (Cd)	<0.05 ppm	Nickel (Ni)	<0.02 ppm
Chromium (Cr)	<0.02 ppm	Strontium (Sr)	<0.02 ppm
Cobalt (Co)	<0.02 ppm	Titanium (Ti)	<0.02 ppm
Lithium (Li)	<0.02 ppm	Vanadium (V)	<0.02 ppm

Lot analysis on the label.

tert-Butanol
specified reagents for general laboratory work- extra pure

SLR

packaging		Mfr. No
500 mℓ	HDPE Plastic bottle	B/5250/08
2.5 ℓ	HDPE Plastic bottle	B/5250/17
25 ℓ	HDPE Plastic container	B/5250/25

C ₄ H ₁₀ O	EINECS: 200-889-7	
CAS: 75-65-0	UN: 1120	
MW: 74.12	GHS: Flam. Liq. 2, Acute	
Density at 20°C: 0.780	Tox. 4, Eye Irrit. 2, STOT SE 3	
Boiling point: 83°C	H225, H319, H335, H332	
Melting point: 25 to 25.5°C	P261, P304+P340, P280,	
Flash point: 11°C	P305+P351+P338, P210, P240	
ADR/RID: 3, II		

Assay (CG): >99%

Guaranteed analysis			
Acidity / Alkalinity	<0.005 meq/g	Magnesium (Mg)	<2 ppm
Colour	<20 APHA	Potassium (K)	<2 ppm
Water	<0.3%	Sodium (Na)	<10 ppm
Residue after evaporation	<50 ppm	Total phosphorus (P)	<2 ppm
Calcium (Ca)	<5 ppm	Total silicon (Si)	<2 ppm
Copper (Cu)	<1 ppm	Total sulphur (S)	<5 ppm
Iron (Fe)	<2 ppm	Zinc (Zn)	<2 ppm
Lead (Pb)	<1 ppm		

2-n-Butoxyethanol
specified reagents for general laboratory work- extra pure

SLR

packaging		Mfr. No
500 mℓ	Coated Glass bottle	B/4900/PB08
2.5 ℓ	Coated Glass bottle	B/4900/PB17
25 ℓ	Metal drum	B/4900/25

C ₆ H ₁₄ O ₂	EINECS: 203-905-0	
CAS: 111-76-2	UN: 2810	
MW: 118.18	GHS: Acute Tox. 4, Skin Irrit.	
Density at 20°C: 0.901	2, Eye Irrit. 2	
Boiling point: 171°C	H319, H315, H302, H312, H332	
Melting point: -70°C	P280, P307+P311,	
Flash point: 62°C	P302+P350, P305+P351+P338	
ADR/RID: 6.1/III		

Assay (CG): >98%

iso-Butyl acetate
specified reagents for general laboratory work- extra pure

SLR

|--|

n-Butyl acetate
for electronic use (MOS)

packaging		Mfr. No
2.5 ℓ	Glass bottle	B/4960/17
5 ℓ	Aluminium can	B/4960/21
C ₆ H ₁₂ O ₂	ADR/RID: 3, III	
CAS: 123-86-4	EINECS: 204-658-1	
MW: 116.16	UN: 1123	
Density at 20°C: 0.880	GHS: Flam. Liq. 3, STOT SE 3	
Boiling point: 127 to 126°C	H226, H336, EUH066	
Melting point: -77.9°C	P304+P340, P261, P210, P240	
Flash point: 22°C		
Assay (CG): >99%		
Guaranteed analysis		
Acidity / Alkalinity	<0.001 meq/g	Nickel (Ni)
Colour	<10 APHA	Potassium (K)
Water	<0.1%	Sodium (Na)
pH colour change test	: Pass	Strontium (Sr)
Residue after evaporation	<10 ppm	Titanium (Ti)
Aluminium (Al)	<0.02 ppm	Total phosphorus (P)
Barium (Ba)	<0.01 ppm	Total silicon (Si)
Cadmium (Cd)	<0.01 ppm	Total sulphur (S)
Calcium (Ca)	<0.2 ppm	Vanadium (V)
Chromium (Cr)	<0.02 ppm	Zinc (Zn)
Cobalt (Co)	<0.02 ppm	Substances darkened
Copper (Cu)	<0.2 ppm	by H ₂ SO ₄
Iron (Fe)	<0.1 ppm	Di-n-butyl ether
Lead (Pb)	<0.05 ppm	iso-Butyl acetate
Lithium (Li)	<0.01 ppm	Butan -1-ol
Magnesium (Mg)	<0.02 ppm	n-Butyl formate
Manganese (Mn)	<0.01 ppm	Butyl propanoate
Molybdenum (Mo)	<0.01 ppm	

n-Butyl acetate
for analysis

packaging		Mfr. No
1 ℓ	Glass bottle	B/4951/15
2.5 ℓ	Glass bottle	B/4951/17
10 ℓ	Metal can	B/4951/24
200 ℓ	Metal drum	B/4951/27
C ₆ H ₁₂ O ₂	ADR/RID: 3, III	
CAS: 123-86-4	EINECS: 204-658-1	
MW: 116.16	UN: 1123	
Density at 20°C: 0.880	GHS: Flam. Liq. 3, STOT SE 3	
Boiling point: 127 to 126°C	H226, H336, EUH066	
Melting point: -77.9°C	P304+P340, P261, P210, P240	
Flash point: 22°C		
Assay (CG): >99%		
Guaranteed analysis		
Acidity / Alkalinity	<0.001 meq/g	Total phosphorus (P)
Colour	<10 APHA	Total silicon (Si)
Water	<0.1%	Total sulphur (S)
Residue after evaporation	<10 ppm	Zinc (Zn)
Calcium (Ca)	<0.2 ppm	Butan -1-ol
Copper (Cu)	<0.2 ppm	iso-Butyl acetate
Iron (Fe)	<0.1 ppm	n-Butyl formate
Lead (Pb)	<0.05 ppm	Butyl propanoate
Magnesium (Mg)	<0.02 ppm	Di-n-butyl ether
Potassium (K)	<0.2 ppm	Substances darkened
Sodium (Na)	<0.2 ppm	by H ₂ SO ₄
Typical analysis		
Aluminium (Al)	<0.02 ppm	Manganese (Mn)
Barium (Ba)	<0.01 ppm	Molybdenum (Mo)
Cadmium (Cd)	<0.01 ppm	Nickel (Ni)
Chromium (Cr)	<0.01 ppm	Strontium (Sr)
Cobalt (Co)	<0.01 ppm	Titanium (Ti)
Lithium (Li)	<0.01 ppm	Vanadium (V)

Lot analysis on the label.

n-Butyl acetate
specified reagents for general laboratory work- extra pure

packaging		Mfr. No
500 mℓ	Glass bottle	B/4950/08
2.5 ℓ	Glass bottle	B/4950/17
5 ℓ	Fluorinated HDPE Plastic bottle	B/4950/21
10 ℓ	Metal can	B/4950/24
25 ℓ	Metal drum	B/4950/25
200 ℓ	Metal drum	B/4950/27
C ₆ H ₁₂ O ₂	ADR/RID: 3, III	
CAS: 123-86-4	EINECS: 204-658-1	
MW: 116.16	UN: 1123	
Density at 20°C: 0.880	GHS: Flam. Liq. 3, STOT SE 3	
Boiling point: 127 to 126°C	H226, H336, EUH066	
Melting point: -77.9°C	P304+P340, P261, P210, P240	
Flash point: 22°C		
Assay (CG): >98%		
Guaranteed analysis		
Acidity / Alkalinity	<0.005 meq/g	Potassium (K)
Colour	<20 APHA	Sodium (Na)
Water	<0.5%	Total phosphorus (P)
Residue after evaporation	<50 ppm	Total silicon (Si)
Calcium (Ca)	<5 ppm	Total sulphur (S)
Copper (Cu)	<5 ppm	Zinc (Zn)
Iron (Fe)	<2 ppm	Substances darkened
Lead (Pb)	<2 ppm	by H ₂ SO ₄
Magnesium (Mg)	<2 ppm	

n-Butylamine
specified reagents for general laboratory work- extra pure

packaging		Mfr. No
250 mℓ	Coated Glass bottle	B/5260/PB07
C ₄ H ₁₁ N	UN: 1125	
CAS: 109-73-9	GHS: Flam. Liq. 2, Acute	
MW: 73.14	Tox. 3, Skin Corr. 1A, Eye	
Density at 20°C: 0.740	Dam. 1, STOT SE 3	
Boiling point: 78°C	H225, H314, H302, H311,	
Melting point: -49°C	H331, H335	
Flash point: -7°C	P280, P305+P351+P338,	
ADR/RID: 3, 8, II	P310, P210, P260, P264	
EINECS: 203-699-2		
Assay (CG): >98%		

n-Butyl bromide

See 1-Bromobutane, page 61

tert-Butyl bromide

See 2-Bromo-2-methylpropane, page 61

sec-Butyl chloride

See 2-Chlorobutane, page 85

n-Butyl chloride

See 1-Chlorobutane, page 85

Optima LC/MS

Fisher Chemical OPTIMA® LC/MS solvents, blends, additives and modifiers have set the standard of excellence for consistent, reproducible performance in the mobile phase of LC/MS.

Find out more on page 308



Butyl digol
specified reagents for general laboratory work- extra pure

packaging		Mfr. No
500 mℓ	Glass bottle	B/5300/08
2.5 ℓ	Glass bottle	B/5300/17
C ₈ H ₁₈ O ₃	Flash point: 100°C	
CAS: 112-34-5	EINECS: 203-961-6	
MW: 162.23	GHS: Eye Irrit. 2	
Density at 20°C: 0.955	H319	
Boiling point: 231°C	P280, P305+P351+P338	
Melting point: -68°C		
Assay (CG): >98%		

Butyl iodide

See 1-Iodobutane, page 161

tert-Butyl methyl ether
for HPLC

packaging		Mfr. No
2.5 ℓ	Glass bottle	M/4496/17
10 ℓ	Returnable Drum	M/4496/10RSS
25 ℓ	Returnable Drum	M/4496/25RSS
30 ℓ	Returnable Drum	M/4496/30RSS
C ₅ H ₁₂ O	ADR/RID: 3, II	
CAS: 1634-04-4	EINECS: 216-653-1	
MW: 88.15	UN: 2398	
Density at 20°C: 0.740	GHS: Flam. Liq. 2, Skin Irrit. 2	
Boiling point: 54 to 56°C	H225, H315	
Boiling point: -110°C	P302+P352, P210, P240	
Flash point: -28°C		
Assay (CG): >99.5%		
Guaranteed analysis		
Acidity / Alkalinity	<0.0002 meq/g	Residue after evaporation
Water	<0.02%	<2 ppm
Absorbance specification		
Max. Absorbance (1cm cell v HPLC water)		
220 nm	<0.25 A.U.	250 nm
230 nm	<0.15 A.U.	260 nm
240 nm	<0.1 A.U.	270 nm

Lot analysis on the label.

Filtered to 0.2 µm.

tert-Butyl methyl ether
specified reagents for general laboratory work- extra pure

packaging		Mfr. No
500 mℓ	Glass bottle	M/4494/08
2.5 ℓ	Glass bottle	M/4494/17
2.5 ℓ	Coated Glass bottle	M/4494/PB17
5 ℓ	Fluorinated HDPE Plastic bottle	M/4494/21
25 ℓ	Metal drum	M/4494/25
200 ℓ	Metal drum	M/4494/27
C ₅ H ₁₂ O	ADR/RID: 3, II	
CAS: 1634-04-4	EINECS: 216-653-1	
MW: 88.15	UN: 2398	
Density at 20°C: 0.740	GHS: Flam. Liq. 2, Skin Irrit. 2	
Boiling point: 54 to 56°C	H225, H315	
Melting point: -110°C	P302+P352, P210, P240	
Flash point: -28°C		
Assay (CG): >99.5%		
Guaranteed analysis		
Water	<0.05%	Residue after evaporation
		<20 ppm

n-Butyric acid
specified reagents for general laboratory work- extra pure

packaging		Mfr. No
500 mℓ	Coated Glass bottle	B/5500/PB08
C ₄ H ₈ O ₂	EINECS: 203-532-3	
CAS: 107-92-6	UN: 2820	
MW: 88.11	GHS: Skin Corr. 1B	
Density at 20°C: 0.960	H314	
Boiling point: 162 to 165°C	P280, P312, P302+P350,	
Melting point: -5--7°C	P301+P330+P331,	
Flash point: 69°C	P305+P351+P338	
ADR/RID: 8, III		

Cadmium standard metal solution 1000 ppm
for AAS, solution HNO₃ 1M

packaging		Mfr. No
100 mℓ	HDPE Plastic bottle	J/8005/05
500 mℓ	HDPE Plastic bottle	J/8005/08
Density at 20°C: 1.0	H350, H314, H312, H302	
ADR/RID: 6.1, 8, II	P280, P305+P351+P338,	
UN: 3289	P301+P330+P331, P201,	
GHS: Carc. 1B, Acute Tox. 4,	P308+P313, P302+P352	
Skin Corr. 1B		

Manufactured to a tolerance of ± 0.5%.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.

Cadmium standard metal solution 1000 ppm
for ICP, solution HNO₃ 1M

packaging		Mfr. No
500 mℓ	HDPE Plastic bottle	J/8236/08
Density at 20°C: 1.0	H350, H314, H312, H302	
ADR/RID: 6.1, 8, II	P280, P305+P351+P338,	
UN: 3289	P301+P330+P331, P201,	
GHS: Carc. 1B, Acute Tox. 4,	P308+P313, P302+P352	
Skin Corr. 1B		

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
All ICP standards are manufactured from acids, metal & metal salts Ultra pure quality and contain less than <10 ppm total metal impurities

See also the full range of "Standard solution" in the Standard solution section

High Purity Acids
for Elemental Analysis

Fisher Chemical offers a complete range of acids to meet the most challenging applications. Choose from our Optima®, TraceMetal® and PrimarPlus® grades to ensure that the sensitivity you require is always available.

Read more page 24



Cadmium standard metal solution 10 000 ppm
for ICP, solution HNO₃ 1M

packaging	Mfr. No
100 mℓ HDPE Plastic bottle	J/8235/05
ADR/RID: 6.1, 8, II UN: 3289 GHS: Carc. 1B, Acute Tox. 4, Skin Corr. 1B	H350, H314, H312, H302 P280, P305+P351+P338, P301+P330+P331, P201, P308+P313, P302+P352



Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.

Cadmium standard metal solution 1000 ppm
for ICP solution in HNO₃ 2%, NIST certified

packaging	Mfr. No
100 mℓ Plastic bottle	SPPLCD22Y
250 mℓ Plastic bottle	SPPLCD22T
500 mℓ Plastic bottle	SPPLCD22X
H315, H319, H350	



Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Cadmium standard metal solution 1000 ppm
for ICP solution in HCL 2%, NIST certified

packaging	Mfr. No
500 mℓ Plastic bottle	PLCD12X
H350	



Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Cadmium standard metal solution 10 000 ppm
for ICP solution in HNO₃ 5%, NIST certified

packaging	Mfr. No
100 mℓ Plastic bottle	SPPLCD23Y
500 mℓ Plastic bottle	SPPLCD23X
H314, H341, H350, H412, H373	



Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Cadmium standard metal solution 1000 ppm
for ICP-MS solution in HNO₃ 2%, NIST certified

packaging	Mfr. No
100 mℓ Plastic bottle	SPCLCD22Y
H315, H319, H350	



Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Cadmium
extra pure, metal granulated

SLR

packaging	Mfr. No
100 g Plastic jar	C/0080/48
Cd	GHS: Carc. 1B, Muta. 2, Repr. 2, Acute Tox. 2, STOT RE 1, Aquatic Acute 1, Aqua- tic Chronic 1
CAS: 7440-43-9	H330, H350, H341, H372, H410, H361fd
MW: 112.40	P304+P340, P307+P311, P260, P285, P201, P273, P280
Density at 25°C: 8.64	
Boiling point: 765°C	
Melting point: 321°C	
ADR/RID: 6.1, I	
EINECS: 231-152-8	
UN: 2570	



Assay: > 99%

Guaranteed analysis

Arsenic (As)	<2 ppm	Lead (Pb)	<200 ppm
Copper (Cu)	<30 ppm	Zinc (Zn)	<100 ppm
Iron (Fe)	<30 ppm		

Cadmium acetate
extra pure, 98+%

SLR

packaging	Mfr. No
500 g Plastic jar	C/0200/53
C ₄ H ₆ CdO ₄ ·2H ₂ O	GHS: Acute Tox. 4, Aquatic Acute 1, Aquatic Chronic 1
CAS: 543-90-8	H350, H302, H312, H332, H410
MW: 266.52	P280, P302+P352, P273, P308+P313
Density at 20°C: 2.01	
ADR/RID: 6.1, III	
UN: 2570	



Assay: >98%

Guaranteed analysis

Calcium (Ca)	<100 ppm	Potassium (K)	<100 ppm
Copper (Cu)	<100 ppm	Sodium (Na)	<200 ppm
Iron (Fe)	<100 ppm	n-Butyl formate	<0.01%
Lead (Pb)	<100 ppm	Total sulphur (S)	<200 ppm
Magnesium (Mg)	<100 ppm	Zinc (Zn)	<200 ppm

Reagents, Solvents & Volumetric
solutions comply with the
European Pharmacopeia

Read more on page 234

Cadmium chloride hydrate
specified reagents for general laboratory work- extra pure

SLR

packaging	Mfr. No
100 g Plastic jar	C/0280/48
CdCl ₂ ·H ₂ O	Repr. 1B, Acute Tox. 2, STOT RE 1, Aquatic Acute 1, Aquatic Chronic 1
CAS: 7790-78-5	H350, H340, H360FD, H330, H301, H372, H400, H410
MW: 183.32	P201, P280, P308+P313, P301+P312, P304+P340, P260, P273
Boiling point: 967°C	
Melting point: 34°C	
ADR/RID: 6.1, I	
EINECS: 233-296-7	
UN: 2570	
GHS: Carc. 1B, Muta. 1B,	



Assay: >79%

Guaranteed analysis

Calcium (Ca)	<100 ppm	Potassium (K)	<500 ppm
Copper (Cu)	<100 ppm	Sodium (Na)	<500 ppm
Iron (Fe)	<100 ppm	Total sulphur (S)	<500 ppm
Lead (Pb)	<100 ppm	Zinc (Zn)	<500 ppm

Cadmium sulfate octahydrate
for analysis

Certified AR

packaging	Mfr. No
100 g Plastic jar	C/0520/48
3CdSO ₄ ·8H ₂ O	Repr. 1B, Acute Tox. 2, STOT RE 1, Aquatic Acute 1, Aquatic Chronic 1
CAS: 7790-84-3	H350, H340, H360FD, H330, H301, H372, H410
MW: 769.52	P280, P301+P310, P310, P308+P313, P304+P340, P260, P273
Melting point: 41°C	
ADR/RID: 6.1, III	
EINECS: 233-331-4	
UN: 2570	
GHS: Carc. 1B, Muta. 1B,	



Assay: >99%

Guaranteed analysis

pH (aqueous solution 5%)	3 to 6	Sodium (Na)	<10 ppm
Calcium (Ca)	<20 ppm	n-Butyl formate	<0.001%
Copper (Cu)	<5 ppm	Total nitrogen (N)	<50 ppm
Iron (Fe)	<5 ppm	Total phosphorus (P)	<20 ppm
Lead (Pb)	<20 ppm	Total silicon (Si)	<20 ppm
Magnesium (Mg)	<5 ppm	Zinc (Zn)	<5 ppm
Potassium (K)	<2 ppm		

Lot analysis on the label.

Cadmium sulfate octahydrate
specified reagents for general laboratory work- extra pure

SLR

packaging	Mfr. No
100 g Plastic jar	C/0480/48
500 g Plastic jar	C/0480/53
3CdSO ₄ ·8H ₂ O	Repr. 1B, Acute Tox. 2, STOT RE 1, Aquatic Acute 1, Aquatic Chronic 1
CAS: 7790-84-3	H350, H340, H360FD, H330, H301, H372, H410
MW: 769.52	P280, P301+P310, P310, P308+P313, P304+P340, P260, P273
Melting point: 41°C	
ADR/RID: 6.1, III	
EINECS: 233-331-4	
UN: 2570	
GHS: Carc. 1B, Muta. 1B,	



Assay: >98%

Guaranteed analysis

Calcium (Ca)	<200 ppm	Sodium (Na)	<200 ppm
Copper (Cu)	<100 ppm	n-Butyl formate	<0.005%
Iron (Fe)	<100 ppm	Total phosphorus (P)	<100 ppm
Lead (Pb)	<200 ppm	Total silicon (Si)	<100 ppm
Magnesium (Mg)	<50 ppm	Zinc (Zn)	<100 ppm
Potassium (K)	<50 ppm		

Caffeine
specified reagents for general laboratory work- extra pure

SLR

packaging	Mfr. No
100 g Plastic jar	C/0800/48
C ₈ H ₁₀ N ₄ O ₂	EINECS: 200-362-1 UN: 1544
CAS: 58-08-2	GHS: Acute Tox. 4 H302
MW: 194.19	P301+P312, P280
Density at 20°C: 1.23	
Melting point: 234 to 239°C	
ADR/RID: 6.1, III	

Guaranteed analysis

Residue after ignition	<0.1%
------------------------------	-------

Calcium standard metal solution 1000 ppm
for AAS, solution HNO₃ 1M

packaging	Mfr. No
100 mℓ HDPE Plastic bottle	J/8010/05
500 mℓ HDPE Plastic bottle	J/8010/08
ADR/RID: 8, II	H314
UN: 2031	P280, P305+P351+P338, P301+P330+P331, P310
GHS: Skin Corr. 1B, Eye Dam. 1	



Manufactured to a tolerance of ± 0.5%.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.

Calcium standard metal solution 1000 ppm
for ICP, solution HNO₃ 1M

packaging	Mfr. No
500 mℓ HDPE Plastic bottle	J/8241/08
ADR/RID: 8, II	H314
UN: 2031	P280, P305+P351+P338, P301+P330+P331, P310
GHS: Skin Corr. 1B, Eye Dam. 1	



Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.

Calcium standard metal solution 10 000 ppm
for ICP, solution HNO₃ 1M

packaging	Mfr. No
100 mℓ HDPE Plastic bottle	J/8240/05
Density at 20°C: 1.0	H314
ADR/RID: 8, II	P280, P305+P351+P338, P301+P330+P331, P310
UN: 2031	
GHS: Skin Corr. 1B	



Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.

Calcium standard metal solution 1000 ppm
for ICP solution in HNO₃ 2%, NIST certified

packaging	Mfr. No
100 mℓ Plastic bottle	SPPLCA22Y
250 mℓ Plastic bottle	SPPLCA22T
H315, H319	



Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Calcium standard metal solution 1000 ppm
for ICP solution in HCL 2%, NIST certified

packaging	Mfr. No
500 mℓ Plastic bottle	SPPLCA12X

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurites.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Calcium standard metal solution 10 000 ppm
for ICP solution in HNO₃ 5%, NIST certified

packaging	Mfr. No
100 mℓ Plastic bottle	SPPLCA23Y
250 mℓ Plastic bottle	SPPLCA23T
500 mℓ Plastic bottle	SPPLCA23X
H314	



Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurites.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Calcium standard metal solution 10 000 ppm
for ICP solution in HCL 5%, NIST certified

packaging	Mfr. No
500 mℓ Plastic bottle	SPPLCA13X
H318	



Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurites.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Calcium standard metal solution 1000 ppm
for ICP-MS solution in HNO₃ 2%, NIST certified

packaging	Mfr. No
100 mℓ Plastic bottle	SPCLCA22Y
H315, H320	



Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurites.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Calcium
extra pure, metal granulated

SLR

packaging	Mfr. No
25 g Plastic jar	C/0840/46
100 g Plastic jar	C/0840/48
250 g Plastic jar	C/0840/50
Ca	UN: 1401
CAS: 7440-70-2	GHS: Water-react. 2
MW: 40.07	H261
Boiling point: 1484°C	P280, P231+P232, P223,
Melting point: 845°C	P335+P334, P370+P378,
ADR/RID: 4.3, II	P402+P404
EINECS: 231-179-5	



Guaranteed analysis
Identity test : Pass

Calcium hardness indicator tablets
pure, no murexide

packaging	Mfr. No
100 PC Glass jar	C/0885/79

Guaranteed analysis
Identity test : Pass

Calcium acetate
extra pure, dried

SLR

packaging	Mfr. No
250 g Plastic jar	C/0920/50
500 g Plastic jar	C/0920/53

(CH₃COO)₂Ca.xH₂O
CAS: 62-54-4
MW: 158.16
EINECS: 200-540-9

Assay: >95%

Loss on drying 200°C	: 5 to 7%	Sulfate (SO ₄)	<0.2%
Iron (Fe)	<20 ppm	n-Butyl formate	<0.005%
Lead (Pb)	<10 ppm		

Calcium bromide dihydrate
specified reagents for general laboratory work- extra pure

SLR

packaging	Mfr. No
250 g Plastic jar	C/0960/50

CaBr₂.2H₂O
CAS: 22208-73-7
MW: 235.93
Density at 20°C: 3.35
Melting point: 760°C

EINECS: 232-164-6
GHS: Eye Irrit. 2
H319
P280, P305+P351+P338



Guaranteed analysis
Identity test : Pass

Calcium carbide
Technical

packaging	Mfr. No
500 g Metal can	C/1000/53

Approx. CaC₂
CAS: 75-20-7
MW: 64.10
Boiling point: 2300°C
Melting point: >447°C
ADR/RID: 4.3, II

EINECS: 200-848-3
UN: 1402
GHS: Water-react. 1
H260
P223, P231+P232, P280



Guaranteed analysis
Identity test : Pass

Calcium carbonate
for analysis

Certified AR

packaging	Mfr. No
250 g Plastic jar	C/1120/50
500 g Plastic jar	C/1120/53
1 kg Plastic jar	C/1120/60
25 kg Fibre Keg	C/1120/70

CaCO₃
CAS: 471-34-1
MW: 100.09
Melting point: 825°C
EINECS: 207-439-9

GHS: Skin Irrit. 2, Eye Irrit. 2
H315, H319
P302+P352, P280,
P305+P351+P338



Assay: >99%

Guaranteed analysis			
Copper (Cu)	<5 ppm	Total chloride (Cl)	<0.01%
Iron (Fe)	<10 ppm	Total nitrogen (N)	<100 ppm
Lead (Pb)	<5 ppm	Total phosphorus (P)	<50 ppm
Magnesium (Mg)	<200 ppm	Total silicon (Si)	<20 ppm
Potassium (K)	<100 ppm	Total sulphur (S)	<100 ppm
Sodium (Na)	<200 ppm	Zinc (Zn)	<10 ppm

Typical analysis			
Aluminium (Al)	<0.001 %	Lithium (Li)	<1 ppm
Silver (Ag)	<1 ppm	Manganese (Mn)	<1 ppm
Arsenic (As)	<2 ppm	Molybdenum (Mo)	<1 ppm
Barium (Ba)	<0.005 %	Nickel (Ni)	<1 ppm
Cadmium (Cd)	<1 ppm	Strontium (Sr)	<0.05 %
Chromium (Cr)	<1 ppm	Titanium (Ti)	<1 ppm
Cobalt (Co)	<1 ppm	Vanadium (V)	<1 ppm

Lot analysis on the label.

Calcium carbonate
specified reagents for general laboratory work- extra pure

SLR, Eur.Ph., BP

packaging	Mfr. No
1 kg Plastic jar	C/1040/60
3 kg Plastic jar	C/1040/63
5 kg Plastic keg	C/1040/65
25 kg Fibre Keg	C/1040/70

CaCO₃
CAS: 471-34-1
MW: 100.09
Melting point: 825°C
EINECS: 207-439-9

GHS: Skin Irrit. 2, Eye Irrit. 2
H315, H319
P302+P352, P280,
P305+P351+P338



Assay: >98%

Guaranteed analysis			
Copper (Cu)	<100 ppm	n-Butyl formate	<0.05%
Iron (Fe)	<200 ppm	Total phosphorus (P)	<500 ppm
Lead (Pb)	<10 ppm	Total silicon (Si)	<1000 ppm
Magnesium (Mg)	<2000 ppm	Total sulphur (S)	<500 ppm
Potassium (K)	<500 ppm	Zinc (Zn)	<100 ppm
Sodium (Na)	<2000 ppm		

Calcium chloride
extra pure, fused, granular

SLR

packaging	Mfr. No
500 g Plastic jar	C/1400/53
1 kg Plastic jar	C/1400/60
2.5 kg Plastic jar	C/1400/62

CaCl₂
CAS: 10043-52-4
MW: 110.98
Boiling point: 1600°C
Melting point: 782°C

EINECS: 233-140-8
GHS: Eye Irrit. 2
H319
P301+P312, P280,
P305+P351+P338



Guaranteed analysis
Identity test : Pass

Particle size 1.0-2.0 mm.

Calcium chloride
Technical, fused, porous lump

packaging	Mfr. No
1 kg Plastic jar	C/1320/60
3 kg Plastic jar	C/1320/63
25 kg Fibre Keg	C/1320/70

CaCl₂
CAS: 10043-52-4
MW: 110.98
Boiling point: 1600°C
Melting point: 782°C

EINECS: 233-140-8
GHS: Eye Irrit. 2
H319
P301+P312, P280,
P305+P351+P338



Guaranteed analysis
Identity test : Pass

Calcium chloride dihydrate
for analysis

Certified AR

packaging	Mfr. No
250 g Plastic jar	C/1500/50
500 g Plastic jar	C/1500/53
1 kg Plastic jar	C/1500/60
5 kg Plastic keg	C/1500/65
25 kg Fibre Keg	C/1500/70

CaCl₂.2H₂O
CAS: 10035-04-8
MW: 147.02
Density at 20°C: 0.830
Melting point: 175°C

EINECS: 233-140-8
GHS: Eye Irrit. 2
H319
P280, P305+P351+P338



Assay: 99 to 102%

Guaranteed analysis

pH (aqueous solution 5%)	5 to 7	Sodium (Na)	<100 ppm
Copper (Cu)	<2 ppm	Total nitrogen (N)	<50 ppm
Iron (Fe)	<2 ppm	Total phosphorus (P)	<50 ppm
Lead (Pb)	<5 ppm	Total silicon (Si)	<50 ppm
Magnesium (Mg)	<50 ppm	Total sulphur (S)	<50 ppm
Potassium (K)	<50 ppm	Zinc (Zn)	<20 ppm

Typical analysis

Aluminium (Al)	<2 ppm	Lithium (Li)	<1 ppm
Silver (Ag)	<1 ppm	Manganese (Mn)	<1 ppm
Arsenic (As)	<2 ppm	Molybdenum (Mo)	<1 ppm
Barium (Ba)	<5 ppm	Nickel (Ni)	<1 ppm
Cadmium (Cd)	<1 ppm	Strontium (Sr)	<0.02 %
Chromium (Cr)	<1 ppm	Titanium (Ti)	<1 ppm
Cobalt (Co)	<1 ppm	Vanadium (V)	<1 ppm

Lot analysis on the label.

Calcium chloride hexahydrate
extra pure, 98.5%+

NEW

SLR

packaging	Mfr. No
500 g Plastic jar	C/1250/53
1 kg Plastic jar	C/1250/60
5 kg Plastic jar	C/1250/65

CaCl₂.6H₂O
CAS: 7774-34-7
MW: 219.08
EINECS: 233-140-8

GHS: Eye Irrit. 2
H319
P280, P305+P351+P338



Assay: >98.5%

Guaranteed analysis

Appearance. colourless, hygroscopic crystals	Iron (Fe)	<0.0002%
pH (aqueous solution 5%)	Magnesium (Mg)	<0.05%
Substances insoluble in Water	Potassium (K)	<0.01%
Heavy metals (as Pb)	Sodium (Na)	<0.2%
Barium (Ba)	Sulfate (SO ₄)	<0.005%

Calcium gluconate
specified reagents for general laboratory work- extra pure

SLR, Eur.Ph., USP

packaging	Mfr. No
500 g Plastic jar	C/1610/53

C₁₂H₂₂CaO₁₄.H₂O
CAS: 299-28-5
MW: 448.35

Melting point: 120°C
Flash point: >170°C
EINECS: 206-075-8

Assay: 99 to 104%

Guaranteed analysis

Arsenic (As)	<2 ppm	Sulfate (SO ₄)	<0.15%
Lead (Pb)	<5 ppm	n-Butyl formate	<0.08%



Calcium hydride
technical

packaging		Mfr. No
100 g	Plastic jar	C/1620/48
CaH ₂	UN: 1404	
CAS: 7789-78-8	GHS: Water-react. 1	
MW: 42.09	H260	
Density at 20°C: 1.9	P280, P223, P231+P232,	
Melting point: 675°C	P335+P334, P402+P404,	
ADR/RID: 4.3, I	P370+P378	
EINECS: 232-189-2		
Assay: >95%		

Calcium hydroxide
for analysis

Certified AR

packaging		Mfr. No
250 g	Plastic jar	C/1642/50
500 g	Plastic jar	C/1642/53
1 kg	Plastic jar	C/1642/60
Ca(OH) ₂	ADR/RID: 8/III	
CAS: 1305-62-0	EINECS: 215-137-3	
MW: 74.09	UN: 3262	
Density at 20°C: 2.24	GHS: Eye Dam. 1	
Boiling point: 2850°C	H318	
Melting point: 580°C	P280, P305+P351+P338, P310	
Assay: >97%		

Guaranteed analysis

Copper (Cu)	<10 ppm	Total nitrogen (N)	<50 ppm
Iron (Fe)	<50 ppm	Total phosphorus (P)	<50 ppm
Lead (Pb)	<10 ppm	Total silicon (Si)	<100 ppm
Magnesium (Mg)	<1000 ppm	Total sulphur (S)	<100 ppm
Potassium (K)	<500 ppm	Zinc (Zn)	<100 ppm
Sodium (Na)	<500 ppm	Insoluble Acids	<0.02%
Total chloride (Cl)	<0.1%		

Typical analysis

Aluminium (Al)	<0.001%	Lithium (Li)	<1 ppm
Silver (Ag)	<1 ppm	Manganese (Mn)	<5 ppm
Arsenic (As)	<2 ppm	Molybdenum (Mo)	<1 ppm
Barium (Ba)	<0.002%	Nickel (Ni)	<2 ppm
Cadmium (Cd)	<1 ppm	Strontium (Sr)	<0.02%
Chromium (Cr)	<1 ppm	Titanium (Ti)	<1 ppm
Cobalt (Co)	<1 ppm	Vanadium (V)	<1 ppm

Lot analysis on the label.

Calcium hydroxide

SLR

specified reagents for general laboratory work- extra pure

packaging		Mfr. No
1 kg	Plastic jar	C/1640/60
5 kg	Plastic keg	C/1640/65
Ca(OH) ₂	ADR/RID: 8/III	
CAS: 1305-62-0	EINECS: 215-137-3	
MW: 74.09	UN: 3262	
Density at 20°C: 2.24	GHS: Eye Dam. 1	
Boiling point: 2850°C	H318	
Melting point: 580°C	P280, P305+P351+P338, P310	
Assay: >95%		

Guaranteed analysis

Copper (Cu)	<100 ppm	n-Butyl formate	<0.1%
Iron (Fe)	<500 ppm	Total phosphorus (P)	<500 ppm
Lead (Pb)	<100 ppm	Total silicon (Si)	<5000 ppm
Magnesium (Mg)	<5000 ppm	Total sulphur (S)	<1000 ppm
Potassium (K)	<2000 ppm	Zinc (Zn)	<1000 ppm
Sodium (Na)	<500 ppm		

Calcium hypochlorite

Technical

packaging		Mfr. No
500 g	Plastic jar	C/1720/53
1 kg	Plastic jar	C/1720/60
Ca(OCl) ₂	GHS: Ox. Sol. 2, Acute Tox. 4, Skin Corr. 1B, Eye Dam. 1, Aquatic Acute 1	   
CAS: 7778-54-3	H272, H302, H314, H400, EUH031	
MW: 142.98	P280, P210,	
Density at 20°C: 2.35	P301+P330+P331, P310,	
Melting point: 100°C	P305+P351+P338, P273	
ADR/RID: 5.1, 8, II		
EINECS: 231-908-7		
UN: 3485		

Guaranteed analysis

Identity test : Pass

Calcium lactate pentahydrate

SLR, Eur.Ph., USP, BP

specified reagents for general laboratory work- extra pure

packaging		Mfr. No
500 g	Plastic jar	C/1820/53
5 kg	Plastic keg	C/1820/65
C ₆ H ₁₀ CaO ₆ ·5H ₂ O	Melting point: 240°C	
CAS: 5743-47-5	EINECS: 212-406-7	
MW: 308.29		

Assay: 98 to 102%

Guaranteed analysis

Acidity / Alkalinity	Pass	Barium (Ba)	<10 ppm
Loss on drying	22 to 27%	Iron (Fe)	<50 ppm
Heavy metals	<10 ppm	Lead (Pb)	<10 ppm
Magnesium & alkali salts	: <1%	n-Butyl formate	<0.02%
Arsenic (As)	<2 ppm		

Calcium nitrate tetrahydrate

Certified AR

for analysis

packaging		Mfr. No
500 g	Plastic jar	C/1882/53
1 kg	Plastic jar	C/1882/60
Ca(NO ₃) ₂ ·4H ₂ O	EINECS: 233-332-1	 
CAS: 13477-34-4	UN: 1454	
MW: 236.15	GHS: Ox. Sol. 2, Skin Irrit. 2,	
Density at 20°C: 1.86	Eye Irrit. 2, STOT SE 3	
Boiling point: 132°C	H272, H319, H315, H335	
Melting point: 45°C	P280, P305+P351+P338, P210	
ADR/RID: 5.1, III		

Assay: 99 to 101%

Guaranteed analysis

pH (aqueous solution 5%)	5 to 7	Sodium (Na)	<50 ppm
Copper (Cu)	<5 ppm	Total chloride (Cl)	<0.001%
Iron (Fe)	<5 ppm	Total phosphorus (P)	<10 ppm
Lead (Pb)	<5 ppm	Total silicon (Si)	<10 ppm
Magnesium (Mg)	<200 ppm	Total sulphur (S)	<50 ppm
Potassium (K)	<50 ppm	Zinc (Zn)	<5 ppm

Typical analysis

Aluminium (Al)	<2 ppm	Lithium (Li)	<1 ppm
Silver (Ag)	<1 ppm	Manganese (Mn)	<1 ppm
Arsenic (As)	<2 ppm	Molybdenum (Mo)	<1 ppm
Barium (Ba)	<0.002 %	Nickel (Ni)	<1 ppm
Cadmium (Cd)	<1 ppm	Strontium (Sr)	<0.01 %
Chromium (Cr)	<1 ppm	Titanium (Ti)	<1 ppm
Cobalt (Co)	<1 ppm	Vanadium (V)	<1 ppm
Tin (Sn)	<1 ppm		

Lot analysis on the label.

Calcium nitrate tetrahydrate

SLR

specified reagents for general laboratory work- extra pure

packaging		Mfr. No
500 g	Plastic jar	C/1880/53
1 kg	Plastic jar	C/1880/60
3 kg	Plastic jar	C/1880/63
Ca(NO ₃) ₂ ·4H ₂ O	EINECS: 233-332-1	 
CAS: 13477-34-4	UN: 1454	
MW: 236.15	GHS: Ox. Sol. 2, Skin Irrit. 2,	
Density at 20°C: 1.86	Eye Irrit. 2, STOT SE 3	
Boiling point: 132°C	H272, H319, H315, H335	
Melting point: 45°C	P280, P305+P351+P338, P210	
ADR/RID: 5.1, III		

Assay: >99%

Guaranteed analysis

pH (aqueous solution 5%)	4 to 7	Sodium (Na)	<500 ppm
Copper (Cu)	<20 ppm	n-Butyl formate	<0.005%
Iron (Fe)	<20 ppm	Total phosphorus (P)	<50 ppm
Lead (Pb)	<10 ppm	Total silicon (Si)	<50 ppm
Magnesium (Mg)	<500 ppm	Total sulphur (S)	<200 ppm
Potassium (K)	<200 ppm	Zinc (Zn)	<50 ppm

tri-Calcium di-orthophosphate

pure, precipitated

packaging		Mfr. No
500 g	Plastic jar	C/1960/53
Ca ₃ (PO) ₂	Melting point: 1670°C	
CAS: 7758-87-4	EINECS: 231-840-8	
MW: 310.17		

Guaranteed analysis

Identity test : Pass

Calcium oxide

extra pure, lump

SLR

packaging		Mfr. No
1 kg	Plastic jar	C/2100/60
CaO	UN: 1910	 
CAS: 1305-78-8	GHS: Acute Tox. 4, Skin	
MW: 56.08	Corr. 1A, Eye Dam. 1	
Boiling point: 2850°C	H314, H302, EUH014	
Melting point: 2570°C	P301+P330+P331, P280,	
ADR/RID: 8/III	P305+P351+P338,	
EINECS: 215-138-9	P301+P312, P402+P404	

Guaranteed analysis

Identity test : Pass

Calcium oxide

extra pure, powder

SLR

packaging		Mfr. No
1 kg	Plastic jar	C/2120/60
CaO	UN: 1910	 
CAS: 1305-78-8	GHS: Acute Tox. 4, Skin	
MW: 56.08	Corr. 1A, Eye Dam. 1	
Boiling point: 2850°C	H314, H302, EUH014	
Melting point: 2570°C	P301+P330+P331, P280,	
ADR/RID: 8/III	P305+P351+P338,	
EINECS: 215-138-9	P301+P312, P402+P404	

Guaranteed analysis

Identity test : Pass

pH Buffers
NIST Certified Standard
Solution & Concentrate
A full range to discover on page 63

Calcium oxide

pure, powder low in chloride and sulfide

packaging		Mfr. No
250 g	Plastic jar	C/2160/50
CaO	UN: 1910	 
CAS: 1305-78-8	GHS: Acute Tox. 4, Skin	
MW: 56.08	Corr. 1A, Eye Dam. 1	
Boiling point: 2850°C	H314, H302, EUH014	
Melting point: 2570°C	P301+P330+P331, P280,	
ADR/RID: 8/III	P305+P351+P338,	
EINECS: 215-138-9	P301+P312, P402+P404	

Guaranteed analysis

Sulfide (S) <0.005% n-Butyl formate <0.005%

Calcium stearate

Technical

packaging		Mfr. No
500 g	Plastic jar	C/2240/53
(C ₁₇ H ₃₅ -COO) ₂ Ca	EINECS: 216-472-8	
CAS: 1592-23-0	GHS: Eye Irrit. 2, STOT SE 3	
MW: 607.00	H319, H335	
Density at 20°C: 1.03	P261, P280,	
Melting point: 145 to 160°C	P305+P351+P338, P304+P341	

Guaranteed analysis

Complies with standard IR spectrum: Pass

Calcium sulfate anhydrous

extra pure, powder

SLR

packaging		Mfr. No
500 g	Plastic jar	C/2440/53
1 kg	Plastic jar	C/2440/60
CaSO ₄	Density at 20°C: 29.6	
CAS: 7778-18-9	Melting point: 1450°C	
MW: 136.15	EINECS: 231-900-3	

Guaranteed analysis

Identity test : Pass


Calcium sulfate anhydrous

pure, granular, Desiccating agent

packaging		Mfr. No
500 g	Plastic jar	C/2400/53
CaSO ₄	Density at 20°C: 29.6	
CAS: 7778-18-9	Melting point: 1450°C	
MW: 136.15	EINECS: 231-900-3	

Guaranteed analysis

Identity test : Pass



Calcium sulfate dihydrate
for analysis

Certified AR

Packaging		Mfr. No.	
250 g	Plastic jar	C/2360/50	
500 g	Plastic jar	C/2360/53	
1 kg	Plastic jar	C/2360/60	
CaSO ₄ ·2H ₂ O		Density at 20°C: 2.32	
CAS: 10101-41-4		EINECS: 231-900-3	
MW: 172.17			
Assay: 99 to 101%			
Guaranteed analysis			
Acidity	<0.002 meq/g	Sodium (Na)	<500 ppm
Copper (Cu)	<10 ppm	Total chloride (Cl)	<0.005%
Iron (Fe)	<10 ppm	Total nitrogen (N)	<10 ppm
Lead (Pb)	<10 ppm	Total phosphorus (P)	<50 ppm
Magnesium (Mg)	<100 ppm	Total silicon (Si)	<50 ppm
Potassium (K)	<100 ppm	Zinc (Zn)	<10 ppm
Typical analysis			
Aluminium (Al)	<0.005 %	Manganese (Mn)	<1 ppm
Silver (Ag)	<5 ppm	Molybdenum (Mo)	<2 ppm
Barium (Ba)	<0.005 %	Nickel (Ni)	<5 ppm
Cadmium (Cd)	<5 ppm	Strontium (Sr)	<0.02 %
Chromium (Cr)	<5 ppm	Titanium (Ti)	<5 ppm
Cobalt (Co)	<5 ppm	Vanadium (V)	<5 ppm
Lithium (Li)	<1 ppm		
Lot analysis on the label.			

Calcium sulfate dihydrate
extra pure, precipitated

SLR

packaging		Mfr. No	
500 g	Plastic jar	C/2320/53	
CaSO ₄ ·2H ₂ O		Density at 20°C: 2.32	
CAS: 10101-41-4		EINECS: 231-900-3	
MW: 172.17			
Assay: >98%			
Guaranteed analysis			
Copper (Cu)	<50 ppm	n-Butyl formate	<0.05%
Iron (Fe)	<50 ppm	Total phosphorus (P)	<200 ppm
Lead (Pb)	<50 ppm	Total silicon (Si)	<200 ppm
Magnesium (Mg)	<500 ppm	Zinc (Zn)	<50 ppm
Potassium (K)	<500 ppm		

Calcium sulfate hemihydrate
Technical

packaging		Mfr. No
3 kg	Plastic jar	C/2280/63
CaSO ₄ ·1/2H ₂ O		Melting point: 1450°C
CAS: 7778-18-9		EINECS: 231-900-3
MW: 145.15		
Guaranteed analysis		
Identity test		: Pass

Calcium tetrahydrogen di-orthophosphate monohydrate

SLR

specified reagents for general laboratory work- extra pure

packaging		Mfr. No	
500 g	Plastic jar	C/2520/53	
CaH ₄ O ₈ P ₂		GHS: Eye Dam. 1	
CAS: 10031-30-8		H318	
MW: 252.07		P280, P305+P351+P338, P310	
EINECS: 231-837-1			
Assay: 97 to 102%			
Guaranteed analysis			
Arsenic (As)	<2 ppm	n-Butyl formate	<0.005%
Iron (Fe)	<150 ppm		

DL-Camphor
extra pure, origin synthetic

SLR

Packaging		Mfr. No
250 g	Glass jar	C/2560/50
$C_{10}H_{16}O$		
CAS: 76-22-2		
MW: 152.24		
Boiling point: 204°C		
ADR/RID: 4.1, III		
EINECS: 200-945-0		
UN: 2717		
Guaranteed analysis		
Complies with standard IR spectrum: Pass		

Canada balsam
pure, dried in Xylene

packaging		Mfr. No	
100 mℓ	Glass bottle	J/2980/05	
Density at 20°C: 0.86		Tox. 4, Skin Irrit. 2	
Flash point: 24°C		H226, H312, H315, H332	
ADR/RID: 3, III		P210, P282, P302+P350,	
UN: 1993		P260, P304+P340	
GHS: Flam. Liq. 3, Acute			

Canada balsam
For microscopy, filtered

Packaging		Mfr. No
100 g	Glass bottle	C/2640/48
CAS: 8007-47-4		
Density at 20°C: 0.98		
Guaranteed analysis		
Complies with standard IR spectrum: Pass		

Caproic acid See n-Hexanoic acid, page 144

Capryl alcohol See Octan-1-ol, page 220

Carbon standard solution 1000 ppm
for ICP solution in H₂O, NIST certified

packaging		Mfr. No	
100 mℓ	Plastic bottle	SPPLC92Y	
500 mℓ	Plastic bottle	SPPLC92X	

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Standard Volumetric solution

For your convenience we offer a range of NIST traceable solutions. A full list of our products can be found on page 323.

Can't see what you need?
Please email your contact details and a brief description of your requirements to begeel.sha.quotes@thermofisher.com

Carbon disulfide
for analysis

Certified AR

packaging		Mfr. No	
500 mℓ	Coated Glass bottle	C/2920/PB08	
2.5 ℓ	Coated Glass bottle	C/2920/PB17	
CS ₂		UN: 1131	
CAS: 75-15-0		GHS: Flam. Liq. 2, Repr. 2,	
MW: 76.14		Skin Irrit. 2, Eye Irrit. 2,	
Density at 20°C: 1.262		STOT RE 1	
Boiling point: 46°C		H225, H315, H319, H372,	
Melting point: -111°C		H361fd	
Flash point: -30°C		P260, P302+P352, P280,	
ADR/RID: 3, 6.1, I		P305+P351+P338, P210,	
EINECS: 200-843-6		P240, P243	

Assay (CG): >99.8%			
Guaranteed analysis			
Acidity / Alkalinity	<0.0001 meq/g	Magnesium (Mg)	<0.02 ppm
Colour	<10 APHA	Potassium (K)	<0.05 ppm
Water	<0.05%	Sodium (Na)	<0.1 ppm
Residue after evaporation	<20 ppm	Total phosphorus (P)	<0.1 ppm
Calcium (Ca)	<0.2 ppm	Total silicon (Si)	<0.01 ppm
Copper (Cu)	<0.02 ppm	Zinc (Zn)	<0.1 ppm
Iron (Fe)	<0.05 ppm	Benzene	<0.01%
Lead (Pb)	<0.05 ppm		

Typical analysis			
Aluminium (Al)	<0.05 ppm	Manganese (Mn)	<0.01 ppm
Barium (Ba)	<0.02 ppm	Molybdenum (Mo)	<0.01 ppm
Cadmium (Cd)	<0.01 ppm	Nickel (Ni)	<0.01 ppm
Chromium (Cr)	<0.01 ppm	Strontium (Sr)	<0.01 ppm
Cobalt (Co)	<0.01 ppm	Titanium (Ti)	<0.01 ppm
Tin (Sn)	<0.01 ppm	Vanadium (V)	<0.01 ppm
Lithium (Li)	<0.01 ppm		

Carbon disulfide
extra pure, 99+%

SLR

packaging		Mfr. No	
500 mℓ	Coated Glass bottle	C/2880/PB08	
2.5 ℓ	Coated Glass bottle	C/2880/PB17	
CS ₂		UN: 1131	
CAS: 75-15-0		GHS: Flam. Liq. 2, Repr. 2,	
MW: 76.13		Skin Irrit. 2, Eye Irrit. 2,	
Density at 20°C: 1.262		STOT RE 1	
Boiling point: 46°C		H225, H315, H319, H372,	
Melting point: -111°C		H361fd	
Flash point: -30°C		P260, P302+P352, P280,	
ADR/RID: 3, 6.1, I		P305+P351+P338, P210,	
EINECS: 200-843-6		P240, P243	

Assay (CG): >99%			
Guaranteed analysis			
Acidity / Alkalinity	<0.0005 meq/g	Lead (Pb)	<5 ppm
Colour	<20 APHA	Magnesium (Mg)	<2 ppm
Water	<0.1%	Potassium (K)	<2 ppm
Residue after evaporation	<100 ppm	Sodium (Na)	<5 ppm
Calcium (Ca)	<5 ppm	Total phosphorus (P)	<2 ppm
Copper (Cu)	<2 ppm	Total silicon (Si)	<2 ppm
Iron (Fe)	<5 ppm	Zinc (Zn)	<2 ppm

Carborundum
extra pure, powder, fine about 180 grit

SLR

packaging		Mfr. No
500 g	Plastic jar	C/3092/53
CAS: 409-21-2		
Density at 20°C: 3.2		
Melting point: 2700°C		
EINECS: 206-991-8		
Guaranteed analysis		
Identity test		: Pass

Carborundum
extra pure, fine powder, extra fine about 300 grit

SLR

Packaging		Mfr. No
500 g	Plastic jar	C/3093/53
CAS: 409-21-2		
Density at 20°C: 3.2		
Melting point: 2700°C		
EINECS: 206-991-8		
Guaranteed analysis		
Identity test : Pass		

Carborundum
extra pure, powder, coarse, about 24 grit

SLR

Packaging		Mfr. No
500 g	Plastic jar	C/3090/53
CAS: 409-21-2		
Density at 20°C: 3.2		
Melting point: 2700°C		
EINECS: 206-991-8		
Guaranteed analysis		
Identity test		: Pass

Carborundum
pure

Packaging		Mfr. No
500 g	Plastic jar	C/3091/53
CAS: 409-21-2		
Density at 20°C: 3.2		
Melting point: 2700°C		
EINECS: 206-991-8		
Guaranteed analysis		
Identity test		: Pass

Carboxymethyl cellulose
extra pure, low viscosity

SLR

Packaging		Mfr. No
500 g	Plastic jar	C/3100/53
CAS: 9004-32-4		
Melting point: >300°C		
Guaranteed analysis		
pH at 20°C	: 6 to 8	Lead (Pb)
Viscosity	29 to 70 cP	<10 ppm

Carmine
pure

packaging		Mfr. No	
5 g	Glass jar	C/3160/44	
CAS: 1390-65-4			
EINECS: 215-724-4			

Guaranteed analysis			
Complies with standard IR spectrum: Pass			

Casein
extra pure, light white soluble

SLR

Packaging		Mfr. No
500 g	Plastic jar	C/3240/53
CAS: 9000-71-9		
Density at 20°C: 1.25		
Melting point: 280°C		
EINECS: 232-555-1		
Guaranteed analysis		
Complies with standard IR spectrum: Pass		

Castor oil			
pure			
packaging		Mfr. No	
500 ml	Glass bottle	C/3280/08	
2.5 l	Glass bottle	C/3280/17	
CAS: 8001-79-4		Melting point: -10°C	
Density at 20°C: 0.957		Flash point: 229°C	
Boiling point: 229°C		EINECS: 232-293-8	
Guaranteed analysis			
Complies with standard IR spectrum: Pass			

Catechol		SLR
specified reagents for general laboratory work- extra pure		
packaging		Mfr. No
100 g	Plastic jar	C/3320/48
500 g	Plastic jar	C/3320/53
C ₆ H ₆ O ₂		
CAS: 120-80-9		
UN: 2811		
GHS: Acute Tox. 4, Skin Irrit. 2, Eye Irrit. 2		
H319, H315, H302, H312		
P280, P312, P302+P350,		
P302+P352, P305+P351+P338		
P302+P352, P305+P351+P338		
P302+P352, P305+P351+P338		
Flash point: 131°C		
ADR/RID: 6.1, III		
Assay (CG): >99%		

Catecholsulfonphthalein	
<i>pure, Metal indicator</i>	
packaging	
5 g	Glass jar
C/3360/44	
C ₁₉ H ₁₄ O ₇ S	
CAS: 115-41-3	
MW: 386.36	
EINECS: 204-088-3	
Guaranteed analysis	
Complies with standard IR spectrum: Pass	

Cation standard mixed concentrate	
<i>for Ion pair chromatography (for check column)</i>	
packaging	
100 mℓ	HDPE Plastic bottle
J/4554/05	
Guaranteed analysis	
Lithium (Li)	28 to 32 ppm
Sodium (Na)	48 to 52 ppm
Potassium (K)	95 to 105 ppm
<i>Manufactured to a tolerance of ± 1%.</i>	
<i>100mℓ, to make 1000mℓ consisting of: 3ppm Li⁺, 5ppm Na⁺, 5ppm NH₄⁺, 10ppm K⁺</i>	

Cedarwood oil	
<i>For microscopy, thickened</i>	
packaging	
100 mℓ	Glass bottle
C/3440/05	
CAS: 8000-27-9	
Flash point: 110°C	
Guaranteed analysis	
Complies with standard IR spectrum: Pass	

Cedarwood oil	
<i>For microscopy, natural</i>	
packaging	
500 mℓ	Glass bottle
C/3400/08	
CAS: 8000-27-9	
Flash point: 110°C	
Guaranteed analysis	
Complies with standard IR spectrum: Pass	

Cellosolve *See 2-Ethoxyethanol, page 123*

Cellosolve acetate *See 2-Ethoxyethyl acetate, page 124*

Cerium standard metal solution 1000 ppm	
<i>for ICP solution in HNO₃ 2%, NIST certified</i>	
packaging	
100 mℓ	Plastic bottle
500 mℓ	Plastic bottle
SPPLCE22Y	
SPPLCE22X	
H315, H319	
	

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Cerium standard metal solution 10 000 ppm	
<i>for ICP solution in HNO₃ 5%, NIST certified</i>	
packaging	
100 mℓ	Plastic bottle
500 mℓ	Plastic bottle
SPPLCE23Y	
SPPLCE23X	
H314	
	

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Cerium(III) nitrate		SLR
extra pure, low in other rare earths		
packaging		Mfr. No
100 g	Plastic jar	C/3780/48
Ce(NO ₃) ₃ ·6H ₂ O		 
CAS: 10294-41-4		
MW: 434.23		
Density at 200°C: 1.67		
ADR/RID: 5.1, II		
EINECS: 233-297-2		
UN: 1477 GHS: Ox. Sol. 2, Skin Irrit. 2, Eye Irrit. 2, Aquatic Chronic 3 H272, H315, H412, H319 P210, P273, P302+P352, P280, P305+P351+P338		
Assay: >99%		
Guaranteed analysis		
Sulfate (SO ₄)	<0.05%	Other rare earths
n-Butyl formate	<0.02%	<1%

Cerium(IV) oxide		SLR
extra pure, low in other rare earths		
packaging		Mfr. No
100 g	Plastic jar	C/3680/48
CeO ₂	Melting point: 2600°C	
CAS: 1306-38-3	EINECS: 215-150-4	
MW: 172.12		
Guaranteed analysis		
Identity test : Pass		

Cerium(IV) sulfate tetrahydrate			SLR	
extra pure, low in other rare earths				
packaging			Mfr. No	
100 g	Plastic jar		C/3720/48	
Ce(SO ₄) ₂ ·4H ₂ O		EINECS: 237-029-5		
CAS: 10294-42-5		UN: 3085		
MW: 404.29		GHS: Ox. Sol. 2, Skin Corr. 1A		
Density at 20°C: 3.91		H272, H318		
Melting point: 350°C		P301+P330+P331, P280,		
ADR/RID: 5.1, 8, II		P305+P351+P338, P310, P210		
Assay: 98 to 102%				
Guaranteed analysis				
n-Butyl formate		<0.01%	Other rare earths	<1%

Cerium(IV) sulfate solution 0.1, (0.1N)		Eur.Ph. , BP
NIST Standard solution ready to use, For Volumetric analysis		
packaging		Mfr. No
1 ℓ	Glass bottle	J/3000/15
2.5 ℓ	Glass bottle	J/3000/17
10 ℓ	Bag in box	J/3000/24
CAS: 13590-82-4		
Density at 20°C: 1.1		
ADR/RID: 8, II		
EINECS: 237-029-5		
UN: 3264 GHS: Skin Irrit. 2, Eye Irrit. 2 H315, H319 P280, P305+P351+P338		
Manufactured to a tolerance of ± 0.2%.		
See the full range of Standard Volumetric solution page 323.		

Cerium(IV) sulfate, concentrated solution 0.05M, (0.05N)	
<i>Standard concentrate, For volumetric analysis</i>	
packaging	
6	Ampoule
J/3010C/90	
ADR/RID: 8, II	
UN: 2796	
GHS: Skin Corr. 1A, Eye Dam. 1	
H314	
	
<i>Each ampoule makes up 500mℓ of solution</i>	

Cesium standard metal solution 1000 ppm	
<i>for ICP solution in HNO₃ 2%, NIST certified</i>	
packaging	
100 mℓ	Plastic bottle
500 mℓ	Plastic bottle
SPPLCS22Y	
SPPLCS22X	
H315, H319	
	

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Cesium standard metal solution 10 000 ppm	
<i>for ICP solution in HNO₃ 5%, NIST certified</i>	
packaging	
100 mℓ	Plastic bottle
500 mℓ	Plastic bottle
SPPLCS23Y	
SPPLCS23X	
H314	
	

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Cesium chloride	
<i>for ultra centrifuge work , 99.9+%</i>	
packaging	
100 g	Plastic jar
C/0682/48	
CsCl ₂	
CAS: 7647-17-8	
MW: 168.36	
Density at 20°C: 3.983	
Assay: >99.9%	
Guaranteed analysis	
Chromium (Cr)	<0.5 ppm
Copper (Cu)	<0.5 ppm
Iron (Fe)	<0.5 ppm
Lead (Pb)	<1 ppm
Manganese (Mn)	<0.5 ppm
Zinc (Zn)	<0.5 ppm
Absorbance specification	
Max. Absorbance (1cm cell v HPLC water)	
260 nm (sol. aqu. 60%)	<0.05 A.U.
280 nm (sol. aqu. 60%)	<0.05 A.U.

Cesium chloride		Certified AR	
for analysis			
packaging		Mfr. No	
25 g	Plastic jar	C/0680/46	
100 g	Plastic jar	C/0680/48	
250 g	Plastic jar	C/0680/50	
1 kg	Plastic jar	C/0680/60	
CsCl ₂		Boiling point: 1290°C Melting point: 646°C EINECS: 231-600-2	
CAS: 7647-17-8			
MW: 168.36			
Density at 20°C: 3.983			
Assay: >99.5%			
Guaranteed analysis			
pH (aqueous solution 5%)	5 to 8	Sodium (Na)	<100 ppm
Calcium (Ca)	<10 ppm	Total nitrogen (N)	<20 ppm
Copper (Cu)	<2 ppm	Total phosphorus (P)	<10 ppm
Iron (Fe)	<2 ppm	Total silicon (Si)	<20 ppm
Lead (Pb)	<2 ppm	Total sulphur (S)	<50 ppm
Magnesium (Mg)	<5 ppm	Zinc (Zn)	<2 ppm
Potassium (K)	<20 ppm		
Typical analysis			
Aluminium (Al)	<2 ppm	Lithium (Li)	<1 ppm
Silver (Ag)	<1 ppm	Manganese (Mn)	<1 ppm
Arsenic (As)	<2 ppm	Molybdenum (Mo)	<1 ppm
Barium (Ba)	<0.002 %	Nickel (Ni)	<1 ppm
Cadmium (Cd)	<1 ppm	Titanium (Ti)	<1 ppm
Chromium (Cr)	<1 ppm	Vanadium (V)	<1 ppm
Cobalt (Co)	<1 ppm		
Lot analysis on the label.			



Cesium chloride		SLR
<i>specified reagents for general laboratory work- extra pure</i>		
packaging		Mfr. No
25 g	Plastic jar	C/0640/46
100 g	Plastic jar	C/0640/48
CsCl ₂	Boiling point: 1290°C	
CAS: 7647-17-8	Melting point: 646°C	
MW: 168.36	EINECS: 231-600-2	
Density at 20°C: 3.983		
Assay: >99%		

Guaranteed analysis			
pH (aqueous solution 5%)	5 to 8	Potassium (K)	<100 ppm
Calcium (Ca)	<50 ppm	Sodium (Na)	<500 ppm
Copper (Cu)	<10 ppm	Total phosphorus (P)	<50 ppm
Iron (Fe)	<10 ppm	Total silicon (Si)	<100 ppm
Lead (Pb)	<10 ppm	Total sulphur (S)	<200 ppm
Magnesium (Mg)	<20 ppm	Zinc (Zn)	<20 ppm

Cetyl trimethylammonium bromide

for Ion pair chromatography

Certified

packaging		Mfr. No
25 g	Plastic jar	C/3961/46
C ₁₉ H ₄₂ BrN	GHS: Acute Tox. 4, Skin Irrit.	
CAS: 57-09-0	2, Eye Irrit. 2, STOT SE 3,	
MW: 364.46	Aquatic Acute 1, Aquatic	
Melting point: 230°C	Chronic 1	
Flash point: 244°C	H319, H302, H410, H315, H335	
ADR/RID: 6.1, III	P261, P273, P301+P312,	
EINECS: 200-311-3	P302+P352, P280,	
UN: 2811	P305+P351+P338	
Assay: >98%		

Cetyl trimethylammonium bromide

for analysis

Certified AR

packaging		Mfr. No
100 g	Plastic jar	C/3962/48
500 g	Plastic jar	C/3962/53
C ₁₉ H ₄₂ BrN		GHS: Acute Tox. 4, Skin Irrit. 2, Eye Irrit. 2, STOT SE 3, Aquatic Acute 1, Aquatic Chronic 1
CAS: 57-09-0		
MW: 364.46		
Melting point: 230°C		
Flash point: 244°C		
ADR/RID: 6.1, III		
EINECS: 200-311-3		
UN: 2811		
Assay: >99%		

Guaranteed analysis			
Loss on drying	<1%	Potassium (K)	<10 ppm
Calcium (Ca)	<20 ppm	Sodium (Na)	<20 ppm
Copper (Cu)	<5 ppm	Total phosphorus (P)	<10 ppm
Iron (Fe)	<5 ppm	Total silicon (Si)	<20 ppm
Lead (Pb)	<5 ppm	Total sulphur (S)	<200 ppm
Magnesium (Mg)	<5 ppm	Zinc (Zn)	<5 ppm

Lot analysis on the label.

Cetyl trimethylammonium bromide

pure

packaging		Mfr. No
500 g	Plastic jar	C/3960/53
25 kg	Fibre Keg	C/3960/70
C ₁₉ H ₄₂ BrN		GHS: Acute Tox. 4, Skin Irrit.
CAS: 57-09-0		2, Eye Irrit. 2, STOT SE 3,
MW: 364.46		Aquatic Acute 1, Aquatic
Melting point: 230°C		Chronic 1
Flash point: 244°C		H319, H302, H410, H315, H335
ADR/RID: 6.1, III		P261, P273, P301+P312,
EINECS: 200-311-3		P302+P352, P280,
UN: 2811		P305+P351+P338
Guaranteed analysis		
Identity test : Pass		

Charcoal

extra pure, lump, wood

Packaging		Mfr. No
1 kg	Plastic jar	C/4020/60
CAS: 7440-44-0	EINECS: 231-153-3	
MW: 12.00	UN: 1361	
ADR/RID: 4.2/III		
Guaranteed analysis		
Identity test		: Pass

Charcoal

pure, powder wood

Packaging		Mfr. No
1 kg	Plastic jar	C/4010/60
CAS: 7440-44-0	EINECS: 231-153-3	
MW: 12.00	UN: 1361	
ADR/RID: 4.2/III		
Guaranteed analysis		
Identity test		: Pass

Charcoal

pure, powder, animal

Packaging		Mfr. No
1 kg	Plastic jar	C/3990/60
CAS: 7440-44-0		
MW: 12.00		
INECS: 231-153-3		
Guaranteed analysis		
Identity test		: Pass

Charcoal

pure, animal, granular for filters

Packaging		Mfr. No
1 kg	Plastic jar	C/4000/60
CAS: 7440-44-0		
MW: 12.00		
INECS: 231-153-3		
Guaranteed analysis		
Identity test		: Pass

Charcoal

pure, decolourising powder, activated

packaging	Mfr. No	
500 g Plastic jar	C/4040/53	
1 kg Plastic jar	C/4040/60	
5 kg Plastic keg	C/4040/65	
25 kg Fibre Keg	C/4040/70	
CAS: 7440-44-0	ADR/RID: 4.2/III	
MW: 12.00	EINECS: 231-153-3	
Density at 20°C: 1.8 to 2.1	UN: 1362	
Melting point: 3652°C		

Guaranteed analysis		
Identity test : Pass		

Charcoal

pure, decolourising powder,activated, acid washed

Packaging		Mfr. No
50 g	Plastic jar	C/4080/50
CAS: 7440-44-0	EINECS: 231-153-3	
MW: 12.00	UN: 1361	
ADR/RID: 4/2III		
Guaranteed analysis		
Identity test	: Pass	

Charcoal

pure, activated for gas adsorption, approx.18 mesh (1mm), granular

packaging		Mfr. No
500 g	Plastic jar	C/4120/53
1 kg	Plastic jar	C/4120/60
C		ADR/RID: 4.2/III
CAS: 7440-44-0		EINECS: 231-153-3
MW: 12.00		UN: 1361
Guaranteed analysis		
Identity test : Pass		

Chemical binder spillage treatment

pure

packaging		Mfr. No
500 g	Plastic jar	S/7435/53
Guaranteed analysis		
Identity test		: Pass

Chlorate standard 1000 ppm

For Ion chromatography in H₂O, NIST certified

packaging	Mfr. No	
100 mℓ Plastic bottle	SPASCLO392Y	
500 mℓ Plastic bottle	SPASCLO392X	

Manufactured to a tolerance of ± 0.3%.

Check with your local distributor for availability of Spex product

Chloral hydrate

SLR, Eur.Ph., BP

specified reagents for general laboratory work- extra pure

packaging		Mfr. No
500 g	Plastic jar	C/4280/53
10 kg	Plastic keg	C/4280/68
C ₂ H ₃ Cl ₃ O ₂		EINECS: 206-117-5
CAS: 302-17-0		UN: 2811
MW: 165.41		GHS: Acute Tox. 3, Skin Irrit.
Density at 20°C: 1.908		2, Eye Irrit. 2
Boiling point: 97°C		H301, H319, H315
Melting point: 57°C		P301+P310, P302+P352,
ADR/RID: 6.1, II		P280, P305+P351+P338
Assay: >99%		
Guaranteed analysis		
Residue after ignition	<0.1%	n-Butyl formate
		<0.01%

Chloramine-T Sodium salt

SLR, Eur.Ph., BP

specified reagents for general laboratory work- extra pure

packaging	Mfr. No	
250 g Plastic jar	C/4320/50	
500 g Plastic jar	C/4320/53	
C ₇ H ₇ ClNNaO ₂ .S.3H ₂ O	UN: 3263	
CAS: 7080-50-4	GHS: Acute Tox. 4, Skin	
MW: 281.69	Corr. 1B, Resp. Sens. 1	
Melting point: 170 to 177°C	H314, H302, H334, EUH031	
Flash point: 192°C	P301+P330+P331, P280,	
ADR/RID: 8, III	P305+P351+P338, P261,	
EINECS: 204-854-7	P304+P341, P301+P312	

Assay: >98%		
Guaranteed analysis		
sodium chloride	<1.5%	

Chloramphenicol

specified reagents for general laboratory work- extra pure

packaging		Mfr. No
50 g	Plastic jar	C/4322/47
C ₁₁ H ₁₂ Cl ₂ N ₂ O ₅	EINECS: 200-287-4	
CAS: 56-75-7	GHS: Carc. 1B	
MW: 323.14	H350	
Melting point: 149 to 153°C	P201, P308+P313	
Assay: >99%		

Chloride standard 100 ppm

For Ion chromatography in H₂O, NIST certified

packaging		Mfr. No
100 mℓ	Plastic bottle	SPASCL91Y
500 mℓ	Plastic bottle	SPASCL91X
Manufactured to a tolerance of ± 0.3%.		
Check with your local distributor for availability of Spex product		

Chloride standard solution 1000 ppm (1 mg/mℓ)

for Ion pair chromatography, Ready-to-use

packaging	Mfr. No	
100 mℓ HDPE Plastic bottle	J/4546/05	
CAS: 7447-40-7		
EINECS: 231-211-8		

Manufactured to a tolerance of ± 1%.

Chloride standard 1000 ppm

For Ion chromatography in H₂O, NIST certified

packaging		Mfr. No
100 mℓ	Plastic bottle	SPASCL92Y
500 mℓ	Plastic bottle	SPASCL92X
Manufactured to a tolerance of ± 0.3%.		
Check with your local distributor for availability of Spex product		

Chloride Colour Reagent

pure

packaging	Mfr. No	
4.5 ℓ HDPE Plastic bottle	C/4345/20	
Density at 20°C: 0.992 to 1.002	H226, H302, H312, H332, H370	
ADR/RID: 3, III	P280, P301+P312,	
UN: 1993	P302+P350, P304+P340,	
GHS: Flam. Liq. 3, Acute	P210, P311	
Tox. 4, STOT SE 1		

Made up to 1ℓ with Analytical Reagent water.

Chlorite standard 1000 ppm

For Ion chromatography in H₂O, NIST certified

packaging		Mfr. No
100 mL	Plastic bottle	SPASCL0292Y
500 mL	Plastic bottle	SPASCL0292X
Manufactured to a tolerance of $\pm 0.3\%$.		
Check with your local distributor for availability of Spex product		

Chloroacetic acid		Certified AR
for analysis		
packaging	Mfr. No	
500 g Glass jar	C/4400/53	
25 kg Plastic keg	C/4400/70	
CICH ₂ COOH	GHS: Acute Tox. 2, Skin	
CAS: 79-11-8	Corr. 1B, Eye Dam. 1, STOT	
MW: 94.50	SE 3, Aquatic Acute 1	
Density at 20°C: 1.58	H301, H335, H311, H314,	
Boiling point: 189°C	H400, H330	
Melting point: 61 to 63°C	P301+P310, P280,	
Flash point: 126°C	P302+P350,	
ADR/RID: 6.1, 8, II	P301+P330+P331,	
EINECS: 201-178-4	P305+P351+P338,	
UN: 1751	P304+P340, P261, P273	
Assay: >99.5%		

Guaranteed analysis			
Calcium (Ca)	<20 ppm	Potassium (K)	<10 ppm
Chromium (Cr)	<1 ppm	Sodium (Na)	<100 ppm
Copper (Cu)	<5 ppm	Total nitrogen (N)	<50 ppm
Iron (Fe)	<5 ppm	Total phosphorus (P)	<20 ppm
Lead (Pb)	<5 ppm	Total silicon (Si)	<10 ppm
Magnesium (Mg)	<5 ppm	Total sulphur (S)	<50 ppm
Nickel (Ni)	<5 ppm	Zinc (Zn)	<5 ppm
Typical analysis			
Aluminium (Al)	<1 ppm	Lithium (Li)	<1 ppm
Silver (Ag)	<1 ppm	Manganese (Mn)	<1 ppm
Arsenic (As)	<2 ppm	Molybdenum (Mo)	<1 ppm
Barium (Ba)	<1 ppm	Strontium (Sr)	<1 ppm
Cadmium (Cd)	<1 ppm	Titanium (Ti)	<1 ppm
Cobalt (Co)	<1 ppm	Vanadium (V)	<1 ppm

Chloroacetic acid SLR
specified reagents for general laboratory work- extra pure

packaging		Mfr. No	
250 g	Glass jar	C/4360/50	
500 g	Glass jar	C/4360/53	
ClCH ₂ COOH		GHS: Acute Tox. 2, Skin	
CAS: 79-11-8		Corr. 1B, Eye Dam. 1, STOT	
MW: 94.50		SE 3, Aquatic Acute 1	
Density at 20°C: 1.58		H301, H335, H311, H314,	
Boiling point: 189°C		H400, H330	
Melting point: 61 to 63°C		P301+P310, P280,	
Flash point: 126°C		P302+P350,	
ADR/RID: 6.1, 8, II		P301+P330+P331,	
EINECS: 201-178-4		P305+P351+P338,	
UN: 1751		P304+P340, P261, P273	
Assay: >99%			
Guaranteed analysis			
Calcium (Ca)	<100 ppm	Sodium (Na)	<500 ppm
Copper (Cu)	<20 ppm	Total phosphorus (P)	<100 ppm
Iron (Fe)	<20 ppm	Total silicon (Si)	<50 ppm
Lead (Pb)	<20 ppm	Total sulphur (S)	<200 ppm
Magnesium (Mg)	<20 ppm	Zinc (Zn)	<20 ppm
Potassium (K)	<50 ppm		

4-Chloroaniline SLR
specified reagents for general laboratory work- extra pure

packaging	Mfr. No	
100 g Plastic jar	C/4600/48	
C ₆ H ₄ ClN	UN: 2018	
CAS: 106-47-8	GHS: Carc. 1B, Acute Tox. 3,	
MW: 127.57	Skin Sens. 1, Aquatic Acute	
Density at 20°C: 1.169	1, Aquatic Chronic 1	
Boiling point: 232°C	H350, H331, H410, H317,	
Melting point: 68 to 72°C	H311, H301	
Flash point: >188°C	P280, P302+P350, P201,	
ADR/RID: 6.1, II	P308+P313, P304+P340, P273	
EINECS: 203-401-0		
Assay (CG): >99%		

Chloroauric (III) acid		SLR
specified reagents for general laboratory work- extra pure		
packaging	Mfr. No	
1 g Vial	C/4640/43	
HAuCl ₄ .xH ₂ O	GHS: Acute Tox. 4, Skin	
CAS: 16903-35-8	Corr. 1B, Eye Dam. 1	
MW: 339.78	H302, H314	
Melting point: 30°C	P301+P330+P331, P280,	
ADR/RID: 8, III	P305+P351+P338, P310,	
EINECS: 240-948-4	P301+P312	
UN: 3260		
Assay: >49%		
The price charged for this item is subject to alteration and will be that ruling at time of despatch.		

4-Chlorobenzenamine See 4-Chloroaniline, page 84

Chlorobenzene Certified AR
for analysis

packaging	Mfr. No	
1 ℓ Glass bottle	C/4681/15	
2.5 ℓ Glass bottle	C/4681/17	
2.5 ℓ Coated Glass bottle	C/4681/PB17	
C ₆ H ₅ Cl	ADR/RID: 3, III	
CAS: 108-90-7	EINECS: 203-628-5	
MW: 112.56	UN: 1134	
Density at 20°C: 1.108	GHS: Flam. Liq. 3, Acute	
Boiling point: 131°C	Tox. 4, Aquatic Chronic 2	
Melting point: -45°C	H226, H411, H332	
Flash point: 23°C	P273, P304+P340, P210	

Assay (CG): >99.8%			
Guaranteed analysis			
Acidity / Alkalinity	<0.0005 meq/g	Potassium (K)	<0.02 ppm
Colour	<10 APHA	Sodium (Na)	<0.5 ppm
Water	<0.05%	Total phosphorus (P)	<0.5 ppm
Residue after evaporation	<10 ppm	Total silicon (Si)	<0.02 ppm
Calcium (Ca)	<0.2 ppm	Total sulphur (S)	<0.5 ppm
Copper (Cu)	<0.02 ppm	Zinc (Zn)	<0.1 ppm
Iron (Fe)	<0.2 ppm	Benzene	<0.02%
Lead (Pb)	<0.05 ppm	o-Dichlorobenzene	<0.01%
Magnesium (Mg)	<0.2 ppm	p-Dichlorobenzene	<0.02%
Typical analysis			
Aluminium (Al)	<0.02 ppm	Manganese (Mn)	<0.01 ppm
Barium (Ba)	<0.01 ppm	Molybdenum (Mo)	<0.01 ppm
Cadmium (Cd)	<0.01 ppm	Nickel (Ni)	<0.01 ppm
Chromium (Cr)	<0.01 ppm	Strontium (Sr)	<0.01 ppm
Cobalt (Co)	<0.01 ppm	Titanium (Ti)	<0.01 ppm
Lithium (Li)	<0.01 ppm	Vanadium (V)	<0.01 ppm
Lot analysis on the label.			

Chlorobenzene SLR
specified reagents for general laboratory work- extra pure

packaging	Mfr. No	
1 ℓ Glass bottle	C/4680/15	
2.5 ℓ Glass bottle	C/4680/17	
2.5 ℓ Coated Glass bottle	C/4680/PB17	
25 ℓ Metal drum	C/4680/25	
C ₆ H ₅ Cl	ADR/RID: 3, III	
CAS: 108-90-7	EINECS: 203-628-5	
MW: 112.56	UN: 1134	
Density at 20°C: 1.108	GHS: Flam. Liq. 3, Acute	
Boiling point: 131°C	Tox. 4, Aquatic Chronic 2	
Melting point: -45°C	H226, H411, H332	
Flash point: 23°C	P273, P304+P340, P210	

Assay (CG): >99%

Guaranteed analysis			
Acidity / Alkalinity	<0.005 meq/g	Magnesium (Mg)	<2 ppm
Colour	<20 APHA	Potassium (K)	<2 ppm
Water	<0.5%	Sodium (Na)	<5 ppm
Residue after evaporation	<100 ppm	Total phosphorus (P)	<2 ppm
Calcium (Ca)	<5 ppm	Total silicon (Si)	<2 ppm
Copper (Cu)	<2 ppm	Total sulphur (S)	<2 ppm
Iron (Fe)	<2 ppm	Zinc (Zn)	<2 ppm
Lead (Pb)	<2 ppm		

1-Chlorobutane
for HPLC

packaging		Mfr. No
1 ℓ	Glass bottle	C/4756/15
2.5 ℓ	Glass bottle	C/4756/17
C ₄ H ₉ Cl		
CAS: 109-69-3		
EINECS: 203-696-6		
UN: 1127		
GHS: Flam. Liq. 2		
H225		
MW: 92.57		P210, P240, P304+P340
Density at 20°C: 0.880		
Boiling point: 77 to 78°C		
Melting point: -123°C		
Flash point: -12°C		
Assay (CG): >99.7%		
Guaranteed analysis		
Acidity / Alkalinity	<0.0003 meq/g	Residue after evaporation
Water	<0.02%	<2 ppm
Absorbance specification		
Max. Absorbance (1cm cell v HPLC water)		
220 nm	<1 A.U.	260 nm
230 nm	<0.2 A.U.	270 nm
240 nm	<0.05 A.U.	Wavelength at 1 A.U.
250 nm	<0.01 A.U.	: <220 nm
Lot analysis on the label.		
Filtered to 0.2 µm.		

1-Chlorobutane
for protein sequencing

packaging		Mfr. No
2.5 ℓ	Glass bottle	C/4758/17
C ₄ H ₉ Cl	ADR/RID: 3, II	
CAS: 109-69-3	EINECS: 203-696-6	
MW: 92.57	UN: 1127	
Density at 20°C: 0.880	GHS: Flam. Liq. 2	
Boiling point: 77 to 78°C	H225	
Melting point: -123°C	P210, P240, P304+P340	
Flash point: -12°C		
Assay (CG): >99.5%		
Guaranteed analysis		
Water	<0.02%	Residue after evaporation <2 ppm
Absorbance specification		
Max. Absorbance (1cm cell v HPLC water)		
220 nm	<1 A.U.	290 nm (sol. aqu. 1%) : <0.005 A.U.
250 nm	<0.01 A.U.	

1-Chlorobutane SLR
specified reagents for general laboratory work- extra pure

packaging		Mfr. No
100 mL	Glass bottle	C/4750/05
C ₄ H ₉ Cl	ADR/RID: 3, II	
CAS: 109-69-3	EINECS: 203-696-6	
MW: 92.57	UN: 1127	
Density at 20°C: 0.880	GHS: Flam. Liq. 2	
Boiling point: 77 to 78°C	H225	
Melting point: -123°C	P210, P240, P304+P340	
Flash point: -12°C		
Assay (CG): >98%		
Guaranteed analysis		
Residue after evaporation		<50 ppm

2-Chlorobutane		SLR
specified reagents for general laboratory work- extra pure		
packaging	Mfr. No	
100 mℓ Glass bottle	C/4759/05	
C ₄ H ₉ Cl	ADR/RID: 3, II	
CAS: 78-86-4	EINECS: 201-151-7	
MW: 92.57	UN: 1127	
Density at 20°C: 0.870	GHS: Flam. Liq. 2	
Boiling point: 68 to 70°C	H225	
Melting point: -140°C	P210, P240	
Flash point: -15°C		
Assay (CG): >95%		

2-Chloroethanol
pure

packaging	Mfr. No	
500 mℓ Coated Glass bottle	C/4880/PB08	
C ₂ H ₅ ClO	EINECS: 203-459-7	
CAS: 107-07-3	UN: 1135	
MW: 80.52	GHS: Flam. Liq. 3, Acute	
Density at 20°C: 1.2	Tox. 1, Eye Dam. 1	
Boiling point: 129 to 130°C	H226, H330, H310, H300, H318	
Melting point: -63°C	P280, P361, P302+P350,	
Flash point: 55°C	P310, P304+P340, P210	
ADR/RID: 6.1, 3, I		
Assay (CG): >99%		

Chloroform stabilized with amylene NEW PACK

for HPLC Certified

Packaging		Mfr. No
1 ℓ	Glass bottle	C/4966/15
2.5 ℓ	Glass bottle	C/4966/17
2.5 ℓ	Coated Glass bottle	C/4966/PB17
10 ℓ	Returnable Drum	C/4966/10RSS
25 ℓ	Returnable Drum	C/4966/25RSS
30 ℓ	Returnable Drum	C/4966/30RSS
CHCl ₃		EINECS: 200-663-8
CAS: 67-66-3		UN: 1888
MW: 119.38		GHS: Carc. 2, Acute Tox. 4,
Density at 20°C: 1.48		Skin Irrit. 2, STOT RE 2
Boiling point: 61°C		H315, H351, H302, H373
Melting point: -63°C		P281, P260,
ADR/RID: 6.1, III		P301+P312, P302+P352
		
		
Guaranteed analysis		
Acidity / Alkalinity	<0.0002 meq/g	Residue after evaporation
Water	<0.05%	<5 ppm
Absorbance specification		
Max. Absorbance (1cm cell v HPLC water)		
250 nm	<0.5 A.U.	270 nm
260 nm	<0.075 A.U.	280 nm
<0.01 A.U.		
<0.005 A.U.		
Lot analysis on the label.		
Filtered to 0.2 µm.		

Chloroform stabilized with amylene Certified
for Gel Permeation Chromatography (GPC)

Packaging		Mfr. No
2.5 ℓ	Coated Glass bottle	C/4970/PB17
CHCl ₃	EINECS: 200-663-8	
CAS: 67-66-3	UN: 1888	
MW: 119.38	GHS: Carc. 2, Acute Tox. 4,	
Density at 20°C: 1.48	Skin Irrit. 2, STOT RE 2	
Boiling point: 61°C	H315, H351, H302, H373	
Melting point: -63°C	P281, P260,	
ADR/RID: 6.1, III	P301+P312, P302+P352	
Guaranteed analysis		
Colour	<5 APHA	Residue after evaporation <5 ppm
Water	<0.02%	

Chloroform stabilized with amylene
for trace metal analysis

Primar

Packaging		Mfr. No.	
2.5 ℓ	Glass bottle	C/4961/17	
CHCl ₃	EINECS: 200-663-8		
CAS: 67-66-3	UN: 1888		
MW: 119.38	GHS: Carc. 2, Acute Tox. 4,		
Density at 20°C: 1.48	Skin Irrit. 2, STOT RE 2		
Boiling point: 61°C	H315, H351, H302, H373		
Melting point: -63°C	P281, P260,		
ADR/RID: 6.1, III	P301+P312, P302+P352		
Assay (CG): >99%			
Guaranteed analysis			
Acidity / Alkalinity	<0.0002 meq/g	Nickel (Ni)	<0.005 ppm
Colour	<5 APHA	Potassium (K)	<0.05 ppm
Water	<0.01%	Propan--1-ol	<0.005 ppm
Residue after evaporation	<2 ppm	Sodium (Na)	<0.2 ppm
Free Chlorine (Cl ₂)	<100 ppm	Strontium (Sr)	<0.005 ppm
Aluminium (Al)	<0.02 ppm	Titanium (Ti)	<0.005 ppm
Barium (Ba)	<0.005 ppm	Total phosphorus (P)	<0.02 ppm
Cadmium (Cd)	<0.005 ppm	Total silicon (Si)	<0.02 ppm
Calcium (Ca)	<0.05 ppm	Total sulphur (S)	<0.1 ppm
Chromium (Cr)	<0.005 ppm	Vanadium (V)	<0.005 ppm
Cobalt (Co)	<0.005 ppm	Zinc (Zn)	<0.01 ppm
Copper (Cu)	<0.01 ppm	Substances darkened	
Iron (Fe)	<0.05 ppm	by H ₂ SO ₄	<5 APHA
Lead (Pb)	<0.005 ppm	Carbon tetrachloride	<0.002%
Lithium (Li)	<0.005 ppm	Dichloromethane	<0.002%
Magnesium (Mg)	<0.02 ppm	Tetrachloroethylene	<0.002%
Manganese (Mn)	<0.005 ppm	Trichloroethylene	<0.002%
Molybdenum (Mo)	<0.005 ppm		

Chloroform stabilized with amylene
for residue and pesticide analysis

Distol

packaging			Mfr. No
1 ℓ	Glass bottle		C/4963/15
2.5 ℓ	Glass bottle		C/4963/17
CHCl ₃	EINECS: 200-663-8		
CAS: 67-66-3	UN: 1888		
MW: 119.38	GHS: Carc. 2, Acute Tox. 4,		
Density at 20°C: 1.48	Skin Irrit. 2, STOT RE 2		
Boiling point: 61°C	H315, H351, H302, H373		
Melting point: -63°C	P281, P260,		
ADR/RID: 6.1, III	P301+P312, P302+P352		
Guaranteed analysis			
Water	<0.05%	GC (FID) Analysis	Pass
Residue after evaporation	<2 ppm	GC (NPD) analysis	Pass
GC (ECD) Analysis	Pass		

GC Analysis
Analysis of pesticides: from the retention time of Lindane, after concentration of 2500 times, no impurity peak higher than that obtained with 10 pg / mL heptachloro-éroxide (ECD) and 5 pg / mL Parathion (NPD).
Analysis of hydrocarbons: no impurity chain longer than the dodecane is present at a concentration greater than 0.01 ppm (FID).
Conditions :
Column: J&W DB 5: 30m x 0.32mm x 0.1µm df.
Temperature profile:
(1) 50°C for 2 min
(2) 40°C / min until 320°C
(3) 320°C for 30 min
Injection: 5µl cold on the column.

Chloroform stabilized with amylene
extra dry, for synthesis

packaging		Mfr. No	
1 ℓ	Glass bottle	C/4962/15	
2.5 ℓ	Glass bottle	C/4962/17	
CHCl ₃	EINECS: 200-663-8		
CAS: 67-66-3	UN: 1888		
MW: 119.38	GHS: Carc. 2, Acute Tox. 4,		
Density at 20°C: 1.480	Skin Irrit. 2, STOT RE 2		
Boiling point: 61°C	H315, H351, H302, H373		
Melting point: -63°C	P281, P260,		
ADR/RID: 6.1, III	P301+P312, P302+P352		

Guaranteed analysis

Water <0.01% Residue after evaporation <2 ppm

Chloroform stabilised with amylene
for analysis

Certified AR

packaging		Mfr. No	
500 mℓ	Coated Glass bottle	C/4960/PB08	
1 ℓ	Glass bottle	C/4960/15	
1 ℓ	Coated Glass bottle	C/4960/PB15	
2.5 ℓ	Glass bottle	C/4960/17	
2.5 ℓ	Coated Glass bottle	C/4960/PB17	
25 ℓ	Metal drum	C/4960/25	
200 ℓ	Metal drum	C/4960/27	
CHCl ₃		EINECS: 200-663-8	
CAS: 67-66-3		UN: 1888	
MW: 119.38		GHS: Carc. 2, Acute Tox. 4,	
Density at 20°C: 1.48		Skin Irrit. 2, STOT RE 2	
Boiling point: 61°C		H315, H351, H302, H373	
Melting point: -63°C		P281, P260,	
ADR/RID: 6.1, III		P301+P312, P302+P352	
Assay (CG): >99.8%			
Guaranteed analysis			
Acidity / Alkalinity	<0.0005 meq/g	Lead (Pb)	<0.02 ppm
Colour	<5 APHA	Magnesium (Mg)	<0.05 ppm
Water	<0.02%	Potassium (K)	<0.1 ppm
Residue after evaporation	<5 ppm	Sodium (Na)	<0.2 ppm
Substances darkened		Total phosphorus (P)	<0.1 ppm
by (H ₂ SO ₄)	<5 APHA	Total silicon (Si)	<0.05 ppm
Acetone & Aldehyde	: Pass	Total sulphur (S)	<0.5 ppm
Acid & Chloride	: Pass	Zinc (Zn)	<0.05 ppm
Calcium (Ca)	<0.1 ppm	Carbon tetrachloride	<0.005%
Free Chlorine (Cl ₂)	<5 ppm	Dichloromethane	<0.005%
Copper (Cu)	<0.02 ppm	Tetrachloroethylene	<0.005%
Iron (Fe)	<0.1 ppm	Trichloroethylene	<0.005%

Typical analysis			
Aluminium (Al)	<0.05 ppm	Manganese (Mn)	<0.01 ppm
Barium (Ba)	<0.01 ppm	Molybdenum (Mo)	<0.01 ppm
Cadmium (Cd)	<0.01 ppm	Nickel (Ni)	<0.01 ppm
Chromium (Cr)	<0.01 ppm	Selenium (Se)	<0.01 ppm
Cobalt (Co)	<0.01 ppm	Strontium (Sr)	<0.01 ppm
Tin (Sn)	<0.01 ppm	Titanium (Ti)	<0.01 ppm
Lithium (Li)	<0.01 ppm	Vanadium (V)	<0.01 ppm

Suitable for use with dithizone.
Lot analysis on the label.
Stabilised with 100 ppm amylene.

Chloroform stabilised with amylene

NEW PACK.

SLR

specific reagents for general laboratory work- extra pure		
packaging		Mfr. No
500 ml	Glass bottle	C/4920/08
1 ℓ	Glass bottle	C/4920/15
2.5 ℓ	Glass bottle	C/4920/17
2.5 ℓ	Coated Glass bottle	C/4920/PB17
5 ℓ	HDPE Plastic bottle	C/4920/21
5 ℓ	HDPE Plastic bottle	C/4920/PC21
25 ℓ	Metal drum	C/4920/25
200 ℓ	Metal drum	C/4920/27
CHCl ₃	EINECS: 200-663-8	
CAS: 67-66-3	UN: 1888	
MW: 119.38	GHS: Carc. 2, Acute Tox. 4,	
Density at 20°C: 1.48	Skin Irrit. 2, STOT RE 2	
Boiling point: 61°C	H315, H351, H302, H373	
Melting point: -63°C	P281, P260,	
ADR/RID: 6.1, III	P301+P312, P302+P352	
Assay (CG): >99%		

Guaranteed analysis			
Acidity / Alkalinity	<0.001 meq/g	Lead (Pb)	<5 ppm
Colour	<10 APHA	Magnesium (Mg)	<2 ppm
Water	<0.1%	Potassium (K)	<2 ppm
Residue after evaporation	<20 ppm	Sodium (Na)	<5 ppm
Substances darkened		Total phosphorus (P)	<2 ppm
by (H ₂ SO ₄)	<20 APHA	Total silicon (Si)	<2 ppm
Calcium (Ca)	<5 ppm	Total sulphur (S)	<2 ppm
Copper (Cu)	<2 ppm	Zinc (Zn)	<2 ppm
Iron (Fe)	<5 ppm		

Stabilised with 100ppm amylene

2-Chloro-2-methyl propane
specified reagents for general laboratory work- extra pure

SLR

packaging		Mfr. No
500 mL	Glass bottle	C/4990/08
C ₄ H ₉ Cl	ADR/RID: 3, II	
CAS: 507-20-0	EINECS: 208-066-4	
MW: 92.57	UN: 1127	
Density at 20°C: 0.850	GHS: Flam. Liq. 2	
Boiling point: 51 to 52°C	H225	
Melting point: -25°C	P210, P240	
Flash point: -27°C		
Guaranteed analysis		
Complies with standard IR spectrum: Pass		

2-Chlorophenol
specified reagents for general laboratory work- extra pure

SLR

packaging		Mfr. No
500 ml	Coated Glass bottle	C/5040/PB08
C ₆ H ₅ ClO	UN: 2021	
CAS: 95-57-8	GHS: Acute Tox. 4, Aquatic	
MW: 128.56	Chronic 2	
Density at 20°C: 1.262	H411, H312, H302, H332	
Boiling point: 173 to 175°C	P280, P301+P330+P331,	
Melting point: 8°C	P305+P351+P338, P310,	
Flash point: 63°C	P304+P340, P302+P352,	
ADR/RID: 6.1, III	P273, P301+P312	
EINECS: 202-433-2		
Guaranteed analysis		
Complies with standard IR spectrum: Pass		

N-Chloro-p-toluenesulfonamide sodium salt

See Chloramine-T, page 83

Chlorotrimethylsilane
for gas chromatography

packaging		Mfr. No
250 g	Glass bottle	T/3520/50
C ₃ H ₉ ClSi	GHS: Flam. Liq. 2, Acute	  
CAS: 75-77-4	Tox. 3, Skin Corr. 1A, Eye Dam. 1	
MW: 108.65	H225, H314, H332, H301,	
Density at 20°C: 0.850	H312, EUH014, EUH071	
Boiling point: 57°C	P301+P310, P280,	
Melting point: -58°C	P301+P330+P331,	
Flash point: -28°C	P305+P351+P338, P310,	
ADR/RID: 3, 8, II	P302+P352, P402+P404,	
EINECS: 200-900-5	P304+P340, P210, P240	
UN: 1298		
Assay (CG): >98%		

Cholesterol
specified reagents for general laboratory work- extra pure

SLR

packaging		Mfr. No
100 g	Plastic jar	C/5360/48
C ₂₇ H ₄₆ O	Boiling point: 360°C	
CAS: 57-88-5	Melting point: 147 to 150°C	
MW: 386.67	Flash point: 250°C	
Density at 20°C: 1.06	EINECS: 200-353-2	
Guaranteed analysis		
Water	<0.5%	Residue after ignition
Specific rotation	-40.5° to -38.5°	<0.1%

Chromic Acid
pure, Mix detergent

packaging		Mfr. No	
1 ℓ	Coated Glass bottle	J/3500/PC15	
Density at 20°C: 1.84	Dam. 1, Resp. Sens. 1, Skin		
Boiling point: 290 to 338°C	Sens. 1, STOT RE 1		
Melting point: 10°C	H350, H340, H360FD,		
ADR/RID: 8, I	H372, H314, H334, H317		
EINECS: 231-639-5	P280, P301+P330+P331,		
UN: 2240	P310, P304+P341,		
GHS: Carc. 1B, Muta. 1B,	P303+P361+P353,		
Repr. 1B, Skin Corr. 1A, Eye	P305+P351+P338		

Chromium standard metal solution 1000 ppm
for AAS, solution HNO₃ 1M

packaging		Mfr. No	
100 mL	HDPE Plastic bottle	J/8015/05	
500 mL	HDPE Plastic bottle	J/8015/08	
ADR/RID: 8, II	P280, P201,		
UN: 3264	P301+P330+P331, P310,		
GHS: Muta. 1B, Skin Corr.	P305+P351+P338,		
1B, Eye Dam. 1	P303+P361+P353		
H314, H340, H350			

Manufactured to a tolerance of ± 0.5%.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.

Chromium standard metal solution 1000 ppm
for ICP, solution HNO₃ 1M

packaging		Mfr. No	
500 mL	HDPE Plastic bottle	J/8246/08	
ADR/RID: 8, II	P280, P201,		
UN: 3264	P301+P330+P331, P310,		
GHS: Muta. 1B, Skin Corr.	P305+P351+P338,		
1B, Eye Dam. 1	P303+P361+P353		
H314, H340, H350			

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.

High Volume Solvent Delivery Systems

High-volume solvent delivery systems are available in 10L to 1000L, in Europe only. Please email your contact details and a brief description of your requirements to begei.sha.quotes@thermofisher.com

Chromium standard metal solution 10 000 ppm
for ICP, solution HNO₃ 1M

packaging		Mfr. No
100 mℓ	HDPE Plastic bottle	J/8245/05
Density at 20°C: 1.0	H350, H340, H314, H334, H317, EUH044	 
ADR/RID: 8, II		
UN: 3264	P280, P301+P330+P331, P310, P303+P361+P353, P305+P351+P338, P308+P313	
GHS: Muta. 1B, Skin Corr. 1B, Resp. Sens. 1, Skin Sens. 1		

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.

Chromium standard metal solution 1000 ppm
for ICP solution in HNO₃ 2%, NIST certified

packaging		Mfr. No
100 mℓ	Plastic bottle	SPPLCR22Y
250 mℓ	Plastic bottle	SPPLCR22T
500 mℓ	Plastic bottle	SPPLCR22X
H315, H319		

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Chromium standard metal solution 1000 ppm
for ICP solution in HCL 2%, NIST certified

packaging		Mfr. No
500 mℓ	Plastic bottle	SPPLCR12X
Manufactured to a tolerance of ± 0.5%. All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities. See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314. Check with your local distributor for availability of Spex product		

Chromium standard metal solution 1000 ppm
for ICP solution in H₂O, NIST certified

packaging		Mfr. No
500 mℓ	Plastic bottle	SPPLCR92X
H317, H333, H334, H340, H350		

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Optima LC/MS

Fisher Chemical OPTIMA® LC/MS solvents, blends, additives and modifiers have set the standard of excellence for consistent, reproducible performance in the mobile phase of LC/MS.

Find out more on page 308



Chromium standard metal solution 10 000 ppm
for ICP solution in HNO₃ 5%, NIST certified

packaging		Mfr. No
100 mℓ	Plastic bottle	SPPLCR23Y
500 mℓ	Plastic bottle	SPPLCR23X
H314		

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Chromium standard metal solution 10 000 ppm
for ICP solution in HCL 5%, NIST certified

packaging		Mfr. No
500 mℓ	Plastic bottle	PLCR13X
Manufactured to a tolerance of ± 0.5%. All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities. See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314. Check with your local distributor for availability of Spex product		

Chromium standard metal solution 10 000 ppm
for ICP solution in H₂O, NIST certified

packaging		Mfr. No
500 mℓ	Plastic bottle	SPPLCR93X
H315, H317, H319, H334, H340, H350, H360, H373, H412		 

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Chromium standard metal solution 1000 ppm
for ICP-MS solution in HNO₃ 2%, NIST certified

packaging		Mfr. No
100 mℓ	Plastic bottle	SPCLCR22Y
H315, H320		

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Chromium(III) chloride hexahydrate

specified reagents for general laboratory work- extra pure

packaging		Mfr. No
100 g	Glass jar	C/5480/48
500 g	Glass jar	C/5480/53
CrCl ₃ ·6H ₂ O CAS: 10060-12-5 MW: 266.45 Melting point: 83°C ADR/RID: 8, 6.1, III EINECS: 233-038-3 UN: 2923		
GHS: Acute Tox. 4, Skin Corr. 1B, Eye Dam. 1 H312, H314, H332 P280, P301+P330+P331, P305+P351+P338, P310, P302+P352, P304+P340		 

Assay: >95%			
Guaranteed analysis			
Calcium (Ca)	<200 ppm	Sodium (Na)	<500 ppm
Copper (Cu)	<200 ppm	Total phosphorus (P)	<50 ppm
Iron (Fe)	<200 ppm	Total silicon (Si)	<100 ppm
Lead (Pb)	<100 ppm	Total sulphur (S)	<500 ppm
Magnesium (Mg)	<100 ppm	Zinc (Zn)	<100 ppm
Potassium (K)	<200 ppm		

Chromium(III) nitrate monohydrate

specified reagents for general laboratory work- extra pure

packaging			Mfr. No
500 g	Glass jar		C/5560/53
Cr(NO ₃) ₃ ·9H ₂ O	UN: 2720		
CAS: 7789-02-8	GHS: Ox. Sol. 2, Skin Irrit. 2,		
MW: 400.15	Eye Irrit. 2		
Melting point: 60°C	H272, H315, H319		
ADR/RID: 5.1, III	P280, P305+P351+P338,		
EINECS: 236-921-1	P302+P352, P210		
Assay: >95%			
Guaranteed analysis			
Iron (Fe)	<500 ppm	n-Butyl formate	<0.003%
Sulfate (SO ₄)	<0.01%		

Chromium(III) oxide, green

specified reagents for general laboratory work- extra pure

packaging		Mfr. No
500 g	Plastic jar	C/5680/53
Cr ₂ O ₃ CAS: 1308-38-9 MW: 151.99 Density at 20°C: 5.21	Boiling point: 4000°C Melting point: 2435°C EINECS: 215-160-9	
Guaranteed analysis		
Identity test : Pass		

Chromium(IV) oxide
for analysis

Certified AR

packaging		Mfr. No
250 g	Plastic jar	C/5920/50
500 g	Plastic jar	C/5920/53
CrO ₃ CAS: 1333-82-0 MW: 99.99 Density at 20°C: 2.7 Melting point: 196°C ADR/RID: 5.1, 6.1. 8, II EINECS: 215-607-8 UN: 1463 GHS: Ox. Sol. 1, Muta. 1B, Repr. 2, Acute Tox. 2, Skin Corr. 1A, Eye Dam. 1, Resp.	Sens. 1, Skin Sens. 1, STOT SE 3, STOT RE 1, Aquatic Acute 1, Aquatic Chronic 1 H271, H350, H340, H361f, H330, H301, H372, H314, H334, H317, H410, H335, H310 P280, P302+P350, P301+P330+P331, P305+P351+P338, P304+P340, P210, P260, P273	   

Assay: >99%			
Guaranteed analysis			
Calcium (Ca)	<50 ppm	n-Butyl formate	<0.005%
Copper (Cu)	<10 ppm	Total nitrogen (N)	<100 ppm
Iron (Fe)	<50 ppm	Total phosphorus (P)	<50 ppm
Lead (Pb)	<50 ppm	Total silicon (Si)	<100 ppm
Magnesium (Mg)	<20 ppm	Total sulphur (S)	<300 ppm
Potassium (K)	<50 ppm	Zinc (Zn)	<20 ppm
Sodium (Na)	<1000 ppm		

Typical analysis

Aluminium (Al)	<0.005%	Molybdenum (Mo)	<0.001%
Barium (Ba)	<0.001%	Nickel (Ni)	<0.005%
Cadmium (Cd)	<0.002%	Silver (Ag)	<0.001%
Cobalt (Co)	<0.002%	Strontium (Sr)	<0.001%
Lithium (Li)	<0.001%	Titanium (Ti)	<0.001%
Manganese (Mn)	<0.002%	Vanadium (V)	<0.001%

Lot analysis on the label

Chromium(VI) oxide

specified reagents for general laboratory work- extra pure

SLR

packaging		Mfr. No
250 g	Plastic jar	C/5880/50
500 g	Plastic jar	C/5880/53
CrO ₃ CAS: 1333-82-0 MW: 99.99 Density at 20°C: 2.7 Melting point: 196°C ADR/RID: 5.1, 6.1. 8, II EINECS: 215-607-8 UN: 1463 GHS: Ox. Sol. 1, Muta. 1B, Repr. 2, Acute Tox. 2, Skin Corr. 1A, Eye Dam. 1, Resp.	Sens. 1, Skin Sens. 1, STOT SE 3, STOT RE 1, Aquatic Acute 1, Aquatic Chronic 1 H271, H350, H340, H361f, H330, H301, H372, H314, H334, H317, H410, H335, H310 P280, P302+P350, P301+P330+P331, P305+P351+P338, P304+P340, P210, P260, P273	    

Assay: >99%

Guaranteed analysis

Calcium (Ca)	<100 ppm	Sodium (Na)	<2000 ppm
Copper (Cu)	<200 ppm	n-Butyl formate	<0.14%
Iron (Fe)	<200 ppm	Total phosphorus (P)	<200 ppm
Lead (Pb)	<100 ppm	Total silicon (Si)	<100 ppm
Magnesium (Mg)	<140 ppm	Total sulphur (S)	<500 ppm
Potassium (K)	<200 ppm	Zinc (Zn)	<140 ppm

Chromium(VI) oxide

Technical

packaging		Mfr. No
500 g	Plastic jar	C/5840/53
2.5 kg	Plastic jar	C/5840/62
CrO ₃ CAS: 1333-82-0 MW: 99.99 Density at 20°C: 2.7 Melting point: 196°C ADR/RID: 5.1, 6.1. 8, II EINECS: 215-607-8 UN: 1463 GHS: Ox. Sol. 1, Muta. 1B, Repr. 2, Acute Tox. 2, Skin Corr. 1A, Eye Dam. 1, Resp.	Sens. 1, Skin Sens. 1, STOT SE 3, STOT RE 1, Aquatic Acute 1, Aquatic Chronic 1 H271, H350, H340, H361f, H330, H301, H372, H314, H334, H317, H410, H335, H310 P280, P302+P350, P301+P330+P331, P305+P351+P338, P304+P340, P210, P260, P273	    

Reagents, Solvents & Volumetric solutions comply with the European Pharmacopeia

Read more on page 234



Chromium(III) potassium sulfate dodecahydrate SLR
specified reagents for general laboratory work- extra pure

Packaging		Mfr. No	
500 g	Plastic jar	C/5720/53	
KCr(SO ₄) ₂ ·12H ₂ O	GHS: Skin Irrit. 2, Eye Irrit. 2		
CAS: 7788-99-0	H315, H319		
MW: 499.39	P302+P352, P280,		
Density at 20°C: 1.83	P305+P351+P338		
Melting point: 89°C			
Assay: >98%			
Guaranteed analysis			
Calcium (Ca)	<1000 ppm	Sodium (Na)	<500 ppm
Copper (Cu)	<200 ppm	n-Butyl formate	<0.01%
Iron (Fe)	<200 ppm	Total phosphorus (P)	<200 ppm
Lead (Pb)	<100 ppm	Total silicon (Si)	<100 ppm
Magnesium (Mg)	<100 ppm	Zinc (Zn)	<100 ppm

Chromotropic acid sodium salt dihydrate SLR
extra pure, for formaldehyde determination

packaging		Mfr. No
25 g	Plastic jar	C/6002/46
C ₁₀ H ₆ O ₈ S ₂ Na ₂ ·2H ₂ O	GHS: Skin Irrit. 2, Eye Irrit. 2,	
CAS: 5808-22-0	STOT SE 3	
MW: 400.28	H335, H315, H319	
Melting point: 300°C	P261, P302+P352, P280,	
EINECS: 204-972-9	P305+P351+P338	
Guaranteed analysis		
Complies with standard IR spectrum: Pass		

Cinnamaldehyde SLR
specified reagents for general laboratory work- extra pure

packaging		Mfr. No
500 mL	Glass bottle	C/6098/08
C ₉ H ₈ O		
CAS: 104-55-2		
MW: 132.16		
EINECS: 203-213-9		
Assay (CG): > 98%		
Guaranteed analysis		
Complies with standard IR spectrum: Pass		

Cinnamic acid
technical, 99+%

packaging		Mfr. No
100 g	Plastic jar	C/6100/48
C ₉ H ₈ O ₂	EINECS: 210-708-3	
CAS: 621-82-9	GHS: Skin Irrit. 2, Eye Irrit. 2,	
MW: 148.05	STOT SE 3	
Boiling point: 300°C	H319, H315, H335	
Melting point: 133°C	P261, P280, P302+P352,	
Flash point: 100°C	P305+P351+P338, P304+P341	
Assay: >99%		

Citric acid, anhydrous SLR, Eur.Ph., BP, USP
specified reagents for general laboratory work- extra pure

packaging		Mfr. No
500 g	Plastic jar	C/6240/53
C ₆ H ₈ O ₇	EINECS: 201-069-1	
CAS: 77-92-9	GHS: Eye Irrit. 2	
MW: 192.13	H319	
Melting point: 153°C	P264, P305+P351+P338,	
Flash point: 100°C	P337+P313, P280	
Assay: >99%		

Citric acid monohydrate Certified
HPLC for electrochemical detection

Packaging		Mfr. No	
500 g	Plastic jar	C/6230/53	
C ₆ H ₈ O ₇ ·H ₂ O	GHS: Eye Irrit. 2		
CAS: 5949-29-1	H319		
MW: 210.14	P264, P280,		
Melting point: 135 to 152°C	P305+P351+P338, P337+P313		
Flash point: 174°C			
Assay: >99.8%			
Guaranteed analysis			
Electrochemical test	Pass	Lead (Pb)	<1 ppm
Iron (Fe)	<1 ppm	n-Butyl formate	<0.0005%
Absorbance specification			
Max. Absorbance (1cm cell v HPLC water)			
254 nm (sol. aqu. 1 M)	<0.5 A.U.	270 nm (sol. aqu. 1 M)	<0.1 A.U.

Tested by HPLC-EC to ensure low electrochemical response.
Detector - ESA Coulochem II
Mode-oxidative
Graphite Electrode.
Lot analysis on the label.
Filtered to 0.2 µm.

Citric acid monohydrate
For foodstuffs analysis, low in lead

Packaging		Mfr. No
500 g	Plastic jar	C/6180/53
C ₆ H ₈ O ₇ ·H ₂ O	GHS: Eye Irrit. 2	
CAS: 5949-29-1	H319	
MW: 210.14	P264, P280,	
Melting point: 135 to 152°C	P305+P351+P338, P337+P313	
Flash point: 174°C		
Guaranteed analysis		
Lead (Pb)	<0.5 ppm	

Citric acid monohydrate
for amino acid analysis

Packaging		Mfr. No
500 g	Plastic jar	C/P330/53
1 kg	Plastic jar	C/P330/60
<chem>C6H8O7.H2O</chem>		
CAS: 5949-29-1		
MW: 210.14		
Melting point: 135 to 152°C		
Flash point: 174°C		
Assay: >99.9%		
Guaranteed analysis		
Tartrate	: Pass	Lead (Pb)
Ninhydrin reaction test	: Pass	n-Butyl formate
Copper (Cu)	<0.5 ppm	Total phosphorus (P)
Iron (Fe)	<1 ppm	Total sulphur (S)
Component of electrophoresis buffers.		

Citric acid monohydrate Certified AR ,Eur.Ph., USP
for analysis

Packaging		Mfr. No
500 g	Plastic jar	C/6200/53
1 kg	Plastic jar	C/6200/60
2.5 kg	Plastic jar	C/6200/62
5 kg	Plastic keg	C/6200/65
25 kg	Fibre Keg	C/6200/70
C₆H₈O₇·H₂O CAS: 5949-29-1 MW: 210.14 Melting point: 135 to 152°C Flash point: 174°C		GHS: Eye Irrit. 2 H319 P264, P280, P305+P351+P338, P337+P313 
Assay: 99.8 to 102%		

Guaranteed analysis			
Insoluble matter	<0.005%	Total chloride (Cl)	<0.0005%
Solution	: Pass	Total phosphorus (P)	<3 ppm
Calcium (Ca)	<50 ppm	Total silicon (Si)	<5 ppm
Copper (Cu)	<0.5 ppm	Total sulphur (S)	<7 ppm
Iron (Fe)	<1 ppm	Zinc (Zn)	<1 ppm
Lead (Pb)	<1 ppm	Phosphate (PO ₄)	<0.001%
Magnesium (Mg)	<5 ppm	Sulfate (SO ₄)	<0.002%
Oxalate (C ₂ O ₄)	<0.005%	Readily carbonisable substances test:	Pass
Potassium (K)	<10 ppm	Residue after ignition	<0.02%
Sodium (Na)	<10 ppm		

Typical analysis			
Aluminium (Al)	<2 ppm	Lithium (Li)	<1 ppm
Silver (Ag)	<1 ppm	Manganese (Mn)	<1 ppm
Arsenic (As)	<2 ppm	Molybdenum (Mo)	<1 ppm
Barium (Ba)	<1 ppm	Nickel (Ni)	<1 ppm
Cadmium (Cd)	<1 ppm	Strontium (Sr)	<2 ppm
Chromium (Cr)	<1 ppm	Titanium (Ti)	<1 ppm
Cobalt (Co)	<1 ppm	Vanadium (V)	<1 ppm

Lot analysis on the label.

Citric acid monohydrate SLR
extra pure, granular

packaging		Mfr. No
500 g	Plastic jar	C/6160/53
1 kg	Plastic jar	C/6160/60
2.5 kg	Plastic jar	C/6160/62
5 kg	Plastic keg	C/6160/65
25 kg	Fibre Keg	C/6160/70
C ₆ H ₈ O ₇ ·H ₂ O	GHS: Eye Irrit. 2	
CAS: 5949-29-1	H319	
MW: 210.14	P264, P280,	
Melting point: 135 to 152°C	P305+P351+P338, P337+P313	
Flash point: 174°C		
Assay: >99.5%		

Cleve's acid See 1-Naphthylamine-7- sulfonic acid, page 211

Clove oil SLR
extra-pure

packaging		Mfr. No
100 mL	Glass bottle	C/6280/05
CAS: 8000-34-8	H315, H319	
Boiling point: 250°C	P280, P261, P302+P352,	
GHS: Skin Irrit. 2, Eye Irrit. 2	P305+P351+P338	
Guaranteed analysis		
Refractive index		1.535 to 1.545

Cobalt standard metal solution 1000 ppm
for AAS, solution HNO₃ 1M

packaging		Mfr. No
100 mL	HDPE Plastic bottle	J/8020/05
500 mL	HDPE Plastic bottle	J/8020/08
Density at 20°C: 1.0	H314	
ADR/RID: 8, II	P280, P305+P351+P338,	
UN: 2031	P301+P330+P331, P310	
GHS: Skin Corr. 1B, Eye Dam. 1		

See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.

Cobalt standard metal solution 1000 ppm
for ICP, solution HNO₃ 1M

packaging		Mfr. No
500 mℓ	HDPE Plastic bottle	J/8251/08
Density at 20°C: 1.0	H314	
ADR/RID: 8, II	P280, P305+P351+P338,	
UN: 2031	P301+P330+P331, P310	
GHS: Skin Corr. 1B, Eye Dam. 1		
Manufactured to a tolerance of ± 0.5%.		

All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
All ICP standards are manufactured from acids, metal & metal salts Ultra pure quality and contain less than <10 ppm total metal impurities
.

Cobalt standard metal solution 10 000 ppm
for ICP, solution HNO₃ 1M

packaging		Mfr. No
100 mL	HDPE Plastic bottle	J/8250/05
Density at 20°C: 1.0	H314	
ADR/RID: 8, II	P280, P301+P330+P331,	
UN: 2031	P310, P305+P351+P338,	
GHS: Skin Corr. 1B, Eye Dam. 1	P303+P361+P353	
Manufactured to a tolerance of ± 0.5%.		

All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.

Cobalt standard metal solution 1000 ppm
for ICP solution in HNO₃ 2%, NIST certified

packaging		Mfr. No
100 mL	Plastic bottle	SPPLCO22Y
250 mL	Plastic bottle	SPPLCO22T
500 mL	Plastic bottle	SPPLCO22X
H315, H319		
Manufactured to a tolerance of ± 0.5%.		
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.		
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.		
Check with your local distributor for availability of Spex product		

Cobalt standard metal solution 1000 ppm
for ICP solution in HCL 2%, NIST certified

packaging		Mfr. No
500 mL	Plastic bottle	SPPLCO12X
H350		
Manufactured to a tolerance of ± 0.5%.		
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.		
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.		
Check with your local distributor for availability of Spex product		



Cobalt standard metal solution 10 000 ppm
for ICP solution in HNO₃ 5%, NIST certified

packaging	Mfr. No
100 mL Plastic bottle	SPPLCO23Y
500 mL Plastic bottle	SPPLCO23X
H314, H317, H334	

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Cobalt standard metal solution 1000 ppm
for ICP-MS solution in HNO₃ 2%, NIST certified

packaging	Mfr. No
100 mL Plastic bottle	SPCLCO22Y
H315, H320	

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Cobalt
extra pure, metal powder

packaging	Mfr. No
100 g Plastic jar	C/6360/48
Co	GHS: Flam. Sol. 2, Resp. Sens. 1, Skin Sens. 1, Aquatic Chronic 4
CAS: 7440-48-4	
MW: 58.93	
Boiling point: 2870°C	H228, H317, H334, H413
ADR/RID: 4.1, III	P261, P304+P341, P342+P311, P280, P302+P352, P210
EINECS: 231-158-0	
UN: 3089	
Guaranteed analysis	
Identity test	: Pass

Cobalt(II) acetate tetrahydrate
pure

packaging	Mfr. No
250 g Plastic jar	C/6400/50
C ₄ H ₆ CoO ₄ ·4H ₂ O	Repr. 1B, Resp. Sens. 1, Skin Sens. 1, Aquatic Acute 1, Aquatic Chronic 1
CAS: 6147-53-1	
MW: 249.08	
Melting point: 140°C	H317, H341, H334, H410, H350i, H360F
ADR/RID: 9, III	P261, P280, P342+P311, P201, P308+P313, P273
UN: 3077	
GHS: Carc. 1B, Muta. 2,	
Guaranteed analysis	
Identity test	: Pass

Cobalt(II) carbonate
technical

packaging	Mfr. No
100 g Plastic jar	C/6459/48
CO ₃ Co	Repr. 1B, Resp. Sens. 1, Skin Sens. 1, Aquatic Acute 1, Aquatic Chronic 1
CAS: 513-79-1	
MW: 118.94	
ADR/RID: 9, III	H350i, H341, H360F, H317, H334, H410
EINECS: 208-169-4	
UN: 3077	
GHS: Carc. 1B, Muta. 2,	P261, P280, P308+P313, P304+P341, P302+P352
Guaranteed analysis	
Assay as Cobalt (Co)	: 43 to 47%

Cobalt(II) chloride anhydrous
specified reagents for general laboratory work- extra pure

packaging	Mfr. No
100 g Plastic jar	C/6600/48
CoCl ₂	Repr. 1B, Acute Tox. 4, Resp. Sens. 1, Skin Sens. 1, Aquatic Acute 1, Aquatic Chronic 1
CAS: 7646-79-9	
MW: 129.84	
Boiling point: 1049°C	H350i, H341, H360F, H302, H334, H317, H410
Melting point: 735°C	
ADR/RID: 6.1, III	P280, P301+P310, P304+P341, P308+P313, P261, P273
EINECS: 231-589-4	
UN: 3288	
GHS: Carc. 1B, Muta. 2,	
Guaranteed analysis	
Identity test	: Pass

Cobalt(II) chloride hexahydrate
for analysis

packaging	Mfr. No
100 g Plastic jar	C/6560/48
500 g Plastic jar	C/6560/53
1 kg Plastic jar	C/6560/60
5 kg Plastic jar	C/6560/65
CoCl ₂ ·6H ₂ O	STOT SE 2, Aquatic Acute 1, Aquatic Chronic 1
CAS: 7791-13-1	
MW: 237.93	
Melting point: 86°C	H410, H302, H341, H317, H334, H350i, H360F
ADR/RID: 9, III	P261, P304+P341, P342+P311, P280, P281, P201, P308+P313, P260, P273, P301+P312
UN: 3077	
GHS: Carc. 1B, Muta. 2, Repr. 1B, Acute Tox. 4, Resp. Sens. 1, Skin Sens. 1,	
Assay: 99 to 101%	
Guaranteed analysis	
pH (aqueous solution 5%)	4 to 6
Calcium (Ca)	<20 ppm
Copper (Cu)	<5 ppm
Iron (Fe)	<20 ppm
Lead (Pb)	<10 ppm
Magnesium (Mg)	<20 ppm
Potassium (K)	<20 ppm
Sodium (Na)	<100 ppm
Total nitrogen (N)	<50 ppm
Total phosphorus (P)	<20 ppm
Total silicon (Si)	<20 ppm
Total sulphur (S)	<50 ppm
Zinc (Zn)	<20 ppm

Lot analysis on the label.

Standard Volumetric solution

For your convenience we offer a range of NIST traceable solutions. A full list of our products can be found on page 323.

Can't see what you need? Please email your contact details and a brief description of your requirements to bege.shi.quotes@thermofisher.com

Cobalt(II) chloride hexahydrate
specified reagents for general laboratory work- extra pure

packaging	Mfr. No
100 g Plastic jar	C/6520/48
500 g Plastic jar	C/6520/53
1 kg Plastic jar	C/6520/60
CoCl ₂ ·6H ₂ O	STOT SE 2, Aquatic Acute 1, Aquatic Chronic 1
CAS: 7791-13-1	
MW: 237.93	
Melting point: 86°C	H410, H302, H341, H317, H334, H350i, H360F
ADR/RID: 9, III	P261, P304+P341, P342+P311, P280, P281, P201, P308+P313, P260, P273, P301+P312
UN: 3077	
GHS: Carc. 1B, Muta. 2, Repr. 1B, Acute Tox. 4, Resp. Sens. 1, Skin Sens. 1,	
Assay: >98%	
Guaranteed analysis	
Calcium (Ca)	<100 ppm
Copper (Cu)	<20 ppm
Iron (Fe)	<50 ppm
Lead (Pb)	<50 ppm
Magnesium (Mg)	<100 ppm
Potassium (K)	<500 ppm
Sodium (Na)	<500 ppm
Total phosphorus (P)	<100 ppm
Total silicon (Si)	<100 ppm
Total sulphur (S)	<200 ppm
Zinc (Zn)	<100 ppm

Cobalt(II) chloride hexahydrate
technical

packaging	Mfr. No
100 g Plastic jar	C/6500/48
500 g Plastic jar	C/6500/53
CoCl ₂ ·6H ₂ O	STOT SE 2, Aquatic Acute 1, Aquatic Chronic 1
CAS: 7791-13-1	
MW: 237.93	
Melting point: 86°C	H410, H302, H341, H317, H334, H350i, H360F
ADR/RID: 9, III	P261, P304+P341, P342+P311, P280, P281, P201, P308+P313, P260, P273, P301+P312
UN: 3077	
GHS: Carc. 1B, Muta. 2, Repr. 1B, Acute Tox. 4, Resp. Sens. 1, Skin Sens. 1,	
Assay: >96.5%	

Cobalt(II) nitrate hexahydrate
for analysis

packaging	Mfr. No
100 g Plastic jar	C/6680/48
250 g Plastic jar	C/6680/50
Co(NO ₃) ₂ ·6H ₂ O	GHS: Carc. 1B, Muta. 2, Repr. 1B, Resp. Sens. 1, Skin Sens. 1, Aquatic Acute 1, Aquatic Chronic 1
CAS: 10026-22-9	
MW: 291.03	
Density at 20°C: 1.88	
Melting point: 55 to 56°C	H317, H334, H350i, H341, H360F, H400, H410
ADR/RID: 9, III	P261, P281, P307+P311, P273
EINECS: 233-402-1	
UN: 3077	
Assay: 99 to 101%	
Guaranteed analysis	
Ammonium (NH ₄)	<0.05%
Calcium (Ca)	<20 ppm
Copper (Cu)	<10 ppm
Iron (Fe)	<20 ppm
Lead (Pb)	<10 ppm
Magnesium (Mg)	<100 ppm
Potassium (K)	<20 ppm
Sodium (Na)	<500 ppm
n-Butyl formate	<0.001%
Total phosphorus (P)	<20 ppm
Total silicon (Si)	<20 ppm
Total sulphur (S)	<50 ppm
Zinc (Zn)	<20 ppm

Lot analysis on the label.

Cobalt(II) nitrate hexahydrate
specified reagents for general laboratory work- extra pure

packaging	Mfr. No
100 g Plastic jar	C/6640/48
500 g Plastic jar	C/6640/53
Co(NO ₃) ₂ ·6H ₂ O	GHS: Carc. 1B, Muta. 2, Repr. 1B, Resp. Sens. 1, Skin Sens. 1, Aquatic Acute 1, Aquatic Chronic 1
CAS: 10026-22-9	
MW: 291.03	
Density at 20°C: 1.88	
Melting point: 55 to 56°C	H317, H334, H350i, H341, H360F, H400, H410
ADR/RID: 9, III	P261, P281, P307+P311, P273
EINECS: 233-402-1	
UN: 3077	
Assay: >98%	
Guaranteed analysis	
Calcium (Ca)	<100 ppm
Copper (Cu)	<50 ppm
Iron (Fe)	<50 ppm
Lead (Pb)	<50 ppm
Magnesium (Mg)	<200 ppm
Potassium (K)	<100 ppm
Sodium (Na)	<1000 ppm
n-Butyl formate	<0.005%
Total phosphorus (P)	<100 ppm
Total silicon (Si)	<100 ppm
Total sulphur (S)	<500 ppm
Zinc (Zn)	<100 ppm

Cobalt(II) sulfate heptahydrate
for analysis

packaging	Mfr. No
100 g Plastic jar	C/6880/48
CoSO ₄ ·7H ₂ O	Resp. Sens. 1, Skin Sens. 1, Aquatic Acute 1, Aquatic Chronic 1
CAS: 10026-24-1	
MW: 281.09	
Density at 20°C: 2.03	
Melting point: 98°C	H317, H410, H334, H350i, H341, H360F, H302
ADR/RID: 9, III	P261, P304+P341, P342+P311, P280, P201, P308+P313, P273
UN: 3077	
GHS: Carc. 1B, Muta. 2, Repr. 1B, Acute Tox. 4,	
Assay: 99 to 101%	
Guaranteed analysis	
pH (aqueous solution 5%)	4 to 6
Calcium (Ca)	<20 ppm
Copper (Cu)	<10 ppm
Iron (Fe)	<20 ppm
Lead (Pb)	<10 ppm
Magnesium (Mg)	<10 ppm
Potassium (K)	<20 ppm
Sodium (Na)	<50 ppm
n-Butyl formate	<0.001%
Total nitrogen (N)	<20 ppm
Total phosphorus (P)	<50 ppm
Total silicon (Si)	<20 ppm
Zinc (Zn)	<20 ppm
Typical analysis	
Aluminium (Al)	<2 ppm
Barium (Ba)	<1 ppm
Cadmium (Cd)	<1 ppm
Chromium (Cr)	<1 ppm
Lithium (Li)	<1 ppm
Manganese (Mn)	<0.001%
Molybdenum (Mo)	<1 ppm
Nickel (Ni)	<0.01%
Silver (Ag)	<1 ppm
Strontium (Sr)	<1 ppm
Titanium (Ti)	<1 ppm
Vanadium (V)	<1 ppm

Lot analysis on the label

Label

Fisher Chemical labels show GHS pictograms clearly. Together with actual lot analysis shown on the label, plus a wealth of additional information, we help you to work safely to achieve accurate and consistent results in your analytical applications.

Code: A/0627/17
Lot: 1298711
Net: 2.5 L
GHS N
MW: 41.04
CAS: 75-05-2
EC: 200-835-2
UN: 1150

Cobalt(II) sulfate heptahydrate
specified reagents for general laboratory work- extra pure

SLR

packaging	Mfr. No
250 g Plastic jar	C/6840/50
500 g Plastic jar	C/6840/53
CoSO ₄ ·7H ₂ O CAS: 10026-24-1 MW: 281.09 Density at 20°C: 2.03 Melting point: 98°C ADR/RID: 9, III UN: 3077 GHS: Carc. 1B, Muta. 2, Repr. 1B, Acute Tox. 4,	Resp. Sens. 1, Skin Sens. 1, Aquatic Acute 1, Aquatic Chronic 1 H317, H410, H334, H350i, H341, H360F, H302 P261, P304+P341, P342+P311, P280, P201, P308+P313, P273
  	
Assay: >98%	
Guaranteed analysis	
Calcium (Ca)<100 ppm	Sodium (Na)<500 ppm
Copper (Cu)<50 ppm	n-Butyl formate<0.005%
Iron (Fe)<50 ppm	Total phosphorus (P)<100 ppm
Lead (Pb)<50 ppm	Total silicon (Si)<100 ppm
Magnesium (Mg)<50 ppm	Zinc (Zn)<100 ppm
Potassium (K)<500 ppm	

2.4,6-Collidine
pure

packaging	Mfr. No
100 mℓ Coated Glass bottle	C/6940/PB05
C ₈ H ₁₁ N CAS: 108-75-8 MW: 121.18 Density at 20°C: 0.910 Boiling point: 171 to 172°C Melting point: -43°C Flash point: 57°C ADR/RID: 3, 6.1, III EINECS: 203-613-3	UN: 1992 GHS: Flam. Liq. 3, Acute Tox. 4, Skin Irrit. 2, Eye Irrit. 2, STOT SE 3 H226, H315, H312, H319, H332, H335 P280, P302+P352, P261, P304+P340, P305+P351+P338, P210
Assay (CG): >98%	

Collodion
pure

packaging	Mfr. No
500 mℓ Glass bottle	C/7000/08
CAS: 9004-70-0 H224, H302, H336, UN: 2059 GHS: Flam. Liq. 1, Acute	Tox. 4, STOT SE 3 EUH019, EUH066 P280, P210, P243, P301+P310

Congo red C.I. 22120
pure, pH and adsorption indicator

packaging	Mfr. No
25 g Glass jar	C/7020/46
C ₂₂ H ₂₂ N ₆ Na ₂ O ₆ S ₂ CAS: 573-58-0 MW: 696.65 Melting point: >360°C ADR/RID: 4.1, III	EINECS: 209-358-4 UN: 1325 GHS: Carc. 1B, Repr. 2 H350, H361d P201, P308+P313, P281, P233
Guaranteed analysis	
pH colour change test : Pass	

Coomassie Blue G250
for electrophoresis, (page blue G90), C.I. 42655

packaging	Mfr. No
25 g Glass jar	C/P541/46
C ₄₇ H ₄₈ N ₃ NaO ₇ S ₂ CAS: 6104-58-1 MW: 854.02	Density at 20°C: >> 1.0 Melting point: 100°C EINECS: 228-058-4

Coomassie Blue R250
for electrophoresis, (page blue 83) C.I. 42660

packaging	Mfr. No
25 g Glass jar	C/P540/46
C ₄₅ H ₄₄ N ₃ NaO ₇ S ₂ CAS: 6104-59-1 MW: 825.96	Melting point: 174 to 180°C EINECS: 228-060-5

Copper standard metal solution 1000 ppm
for ICP, solution HNO₃ 1M

packaging	Mfr. No
500 mℓ HDPE Plastic bottle	J/8256/08
ADR/RID: 8, II EINECS: 221-838-5 UN: 2031 GHS: Skin Corr. 1B, Eye Dam. 1	H314 P280, P305+P351+P338, P301+P330+P331, P310
Manufactured to a tolerance of ± 0.5%. All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities. See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.	

Copper standard metal solution 1000 ppm
for AAS, solution HNO₃ 1M

packaging	Mfr. No
100 mℓ HDPE Plastic bottle	J/8025/05
500 mℓ HDPE Plastic bottle	J/8025/08
ADR/RID: 8, II EINECS: 221-838-5 UN: 2031 GHS: Skin Corr. 1B, Eye Dam. 1	H314 P280, P305+P351+P338, P301+P330+P331, P310
Manufactured to a tolerance of ± 0.5%. See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.	

Copper standard metal solution 10 000 ppm
for ICP, solution HNO₃ 1M

packaging	Mfr. No
100 mℓ HDPE Plastic bottle	J/8255/05
ADR/RID: 8, II UN: 2031 GHS: Skin Corr. 1B	H314 P280, P305+P351+P338, P301+P330+P331, P310
Manufactured to a tolerance of ± 0.5%. All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities. See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.	

Copper standard metal solution 1000 ppm
for ICP solution in HNO₃ 2%, NIST certified

packaging	Mfr. No
100 mℓ Plastic bottle	SPPLCU22Y
250 mℓ Plastic bottle	SPPLCU22T
500 mℓ Plastic bottle	SPPLCU22X
H315, H319	

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurites.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Copper standard metal solution 1000 ppm
for ICP solution in HCL 2%, NIST certified

packaging	Mfr. No
500 mℓ Plastic bottle	SPPLCU12X
Manufactured to a tolerance of ± 0.5%. All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurites. See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314. Check with your local distributor for availability of Spex product	

Copper standard metal solution 10 000 ppm
for ICP solution in HNO₃ 5%, NIST certified

packaging	Mfr. No
100 mℓ Plastic bottle	SPPLCU23Y
500 mℓ Plastic bottle	SPPLCU23X
H314	
Manufactured to a tolerance of ± 0.5%. All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurites. See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314. Check with your local distributor for availability of Spex product	

Copper standard metal solution 10 000 ppm
for ICP solution in HCL 5%, NIST certified

packaging	Mfr. No
500 mℓ Plastic bottle	SPPLCU13X
Manufactured to a tolerance of ± 0.5%. All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurites. See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314. Check with your local distributor for availability of Spex product	

Copper standard metal solution 1000 ppm
for ICP-MS solution in HNO₃ 2%, NIST certified

packaging	Mfr. No
100 mℓ Plastic bottle	SPCLCU22Y
H315, H319	
Manufactured to a tolerance of ± 0.5%. All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurites. See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314. Check with your local distributor for availability of Spex product	

Copper
extra pure, metal powder

SLR

packaging	Mfr. No
100 g Glass jar	C/7120/48
250 g Glass jar	C/7120/50
500 g Glass jar	C/7120/53
Cu CAS: 7440-50-8 MW: 63.54 Density at 20°C: 89.2 Boiling point: 2580°C Melting point: 1083°C ADR/RID: 4.1, III	EINECS: 231-159-6 UN: 3089 GHS: Flam. Sol. 2, Skin Irrit. 2, Eye Irrit. 2, STOT SE 3 H228, H315, H319, H335 P261, P302+P352, P280, P305+P351+P338, P233, P210
Guaranteed analysis	
Identity test : Pass	

Copper
extra pure, metal granular

SLR

packaging	Mfr. No
250 g Plastic jar	C/7170/50
Cu CAS: 7440-50-8 MW: 63.54	Boiling point: 2595°C Melting point: 1083°C EINECS: 231-159-6
Guaranteed analysis	
Identity test : Pass	

Copper
extra pure, metal turning

SLR

packaging	Mfr. No
500 g Plastic jar	C/7160/53
1 kg Plastic jar	C/7160/60
Cu CAS: 7440-50-8 MW: 63.54	Boiling point: 2595°C Melting point: 1083°C EINECS: 231-159-6
Guaranteed analysis	
Identity test : Pass	

Copper
Technical, metal foil (about 0.1 mm)

packaging	Mfr. No
250 g Plastic jar	C/7050/50
500 g Plastic jar	C/7050/53
Cu CAS: 7440-50-8 MW: 63.54	Boiling point: 2595°C Melting point: 1083°C EINECS: 231-159-6
Guaranteed analysis	
Identity test : Pass	

Copper
pure, metal wire about 0.91 mm(20 swg)

packaging	Mfr. No
500 g Plastic jar	C/7210/53
Cu CAS: 7440-50-8 MW: 63.54	Boiling point: 2595°C Melting point: 1083°C EINECS: 231-159-6
Guaranteed analysis	
Identity test : Pass	

Copper(II) acetate monohydrate
specified reagents for general laboratory work- extra pure

SLR

packaging		Mfr. No
250 g	Plastic jar	C/7760/50
500 g	Plastic jar	C/7760/53
C ₄ H ₆ CuO ₄ CAS: 142-71-2 MW: 199.65 Melting point: 115°C ADR/RID: 9, III EINECS: 205-553-3 UN: 3077 GHS: Acute Tox. 4, Skin Irrit.		2, Eye Irrit. 2, STOT SE 3, Aquatic Acute 1, Aquatic Chronic 1 H335, H315, H319, H302, H410 P261, P273, P301+P312, P302+P352, P280, P305+P351+P338
Assay: >98%		
Guaranteed analysis		
Calcium (Ca)	<200 ppm	n-Butyl formate <0.005%
Iron (Fe)	<100 ppm	Total phosphorus (P) <100 ppm
Lead (Pb)	<200 ppm	Total silicon (Si) <100 ppm
Magnesium (Mg)	<50 ppm	Total sulphur (S) <200 ppm
Potassium (K)	<50 ppm	Zinc (Zn) <100 ppm
Sodium (Na)	<200 ppm	

Copper(II) bromide
specified reagents for general laboratory work- extra pure

SLR

packaging		Mfr. No
250 g	Glass jar	C/7840/50
CuBr ₂ CAS: 7789-45-9 MW: 233.36 Density at 20°C: 47.7 Boiling point: 900°C Melting point: 498°C ADR/RID: 8, III EINECS: 232-167-2		UN: 1759 GHS: Acute Tox. 4, Skin Corr. 1B, Eye Dam. 1 H314, H302 P280, P260, P301+P330+P331, P305+P351+P338, P310, P264
Guaranteed analysis		
Identity test		: Pass

Copper(II) carbonate basic
specified reagents for general laboratory work- extra pure

SLR

packaging		Mfr. No
250 g	Plastic jar	C/7880/50
500 g	Plastic jar	C/7880/53
25 kg	Fibre Keg	C/7880/70
CuCO ₃ ·Cu(OH) ₂ CAS: 12069-69-1 MW: 221.10 Density at 20°C: 4.0 Melting point: 200°C ADR/RID: 6.1, III EINECS: 235-113-6		UN: 3288 GHS: Acute Tox. 3, Skin Irrit. 2, Eye Irrit. 2, STOT SE 3 H335, H319, H315, H301 P264, P270, P261, P302+P352, P280, P305+P351+P338
Assay: 54 to 56%		
Guaranteed analysis		
Iron (Fe)	<2000 ppm	Sulfate (SO ₄) <0.2%

Copper(II) carbonate basic
technical

packaging		Mfr. No
500 g	Plastic jar	C/7879/53
3 kg	Plastic jar	C/7879/63
CuCO ₃ ·Cu(OH) ₂ CAS: 12069-69-1 MW: 221.10 Density at 20°C: 4.0 Boiling point: 240°C Melting point: 200°C ADR/RID: 6.1, III EINECS: 235-113-6		UN: 3288 GHS: Acute Tox. 3, Skin Irrit. 2, Eye Irrit. 2, STOT SE 3 H335, H319, H315, H301 P264, P270, P261, P302+P352, P280, P305+P351+P338
Guaranteed analysis		
Identity test		: Pass

Copper (I) chloride
extra pure, min. 96%

SLR

packaging		Mfr. No
100 g	Plastic jar	C/8680/48
CuCl CAS: 7758-89-6 MW: 99.00 Density at 20°C: 4.14 Boiling point: 1490°C Melting point: 430°C ADR/RID: 8, III		EINECS: 231-842-9 UN: 2802 GHS: Acute Tox. 4, Aquatic Acute 1, Aquatic Chronic 1 H302, H410, H400 P301+P310, P273
Assay: >96%		
Guaranteed analysis		
Calcium (Ca)	<200 ppm	Sodium (Na) <500 ppm
Iron (Fe)	<100 ppm	Total phosphorus (P) <100 ppm
Lead (Pb)	<200 ppm	Total silicon (Si) <100 ppm
Magnesium (Mg)	<50 ppm	Zinc (Zn) <100 ppm
Potassium (K)	<50 ppm	

Copper(II) chloride dihydrate
for analysis

Certified AR

packaging		Mfr. No
100 g	Plastic jar	C/7960/48
1 kg	Plastic jar	C/7960/60
CuCl ₂ ·2H ₂ O CAS: 10125-13-0 MW: 170.48 Density at 20°C: 2.54 Melting point: 100°C ADR/RID: 8, III UN: 2802 GHS: Met. Corr. 1, Acute		Tox. 3, Skin Irrit. 2, Eye Irrit. 2, STOT SE 3, Aquatic Acute 1, Aquatic Chronic 1 H301, H335, H319, H315, H410, H290 P261, P273, P302+P352, P280, P305+P351+P338, P301+P310
Assay: >99%		
Guaranteed analysis		
pH (aqueous solution 5%)	3 to 4	Sodium (Na) <20 ppm
Calcium (Ca)	<20 ppm	Total nitrogen (N) <50 ppm
Iron (Fe)	<20 ppm	Total phosphorus (P) <20 ppm
Lead (Pb)	<50 ppm	Total silicon (Si) <20 ppm
Magnesium (Mg)	<5 ppm	Total sulphur (S) <50 ppm
Potassium (K)	<5 ppm	Zinc (Zn) <20 ppm
Typical analysis		
Aluminium (Al)	<1 ppm	Lithium (Li) <2 ppm
Arsenic (As)	<2 ppm	Manganese (Mn) <2 ppm
Barium (Ba)	<2 ppm	Nickel (Ni) <5 ppm
Cadmium (Cd)	<2 ppm	Titanium (Ti) <2 ppm
Chromium (Cr)	<2 ppm	Vanadium (V) <2 ppm
Cobalt (Co)	<2 ppm	

Lot analysis on the label.

Copper(II) chloride dihydrate
specified reagents for general laboratory work- extra pure

SLR

packaging		Mfr. No
100 g	Plastic jar	C/7920/48
250 g	Plastic jar	C/7920/50
500 g	Plastic jar	C/7920/53
5 kg	Plastic keg	C/7920/65
CuCl ₂ ·2H ₂ O CAS: 10125-13-0 MW: 170.48 Density at 20°C: 2.54 Melting point: 100°C ADR/RID: 8, III UN: 2802 GHS: Met. Corr. 1, Acute		Tox. 3, Skin Irrit. 2, Eye Irrit. 2, STOT SE 3, Aquatic Acute 1, Aquatic Chronic 1 H301, H335, H319, H315, H410, H290 P261, P273, P302+P352, P280, P305+P351+P338, P301+P310
Assay: >98%		

Guaranteed analysis			
Calcium (Ca)	<200 ppm	Sodium (Na)	<200 ppm
Iron (Fe)	<100 ppm	Total phosphorus (P)	<100 ppm
Lead (Pb)	<200 ppm	Total silicon (Si)	<100 ppm
Magnesium (Mg)	<50 ppm	Total sulphur (S)	<200 ppm
Potassium (K)	<50 ppm	Zinc (Zn)	<100 ppm

Copper(I) iodide
specified reagents for general laboratory work- extra pure

SLR

packaging		Mfr. No
100 g	Glass jar	C/8800/48
CuI CAS: 7681-65-4 MW: 190.44 Boiling point: 1290°C Melting point: 605°C ADR/RID: 9, III EINECS: 231-674-6 UN: 3077		GHS: Acute Tox. 4, Skin Irrit. 2, Eye Irrit. 2, STOT SE 3, Aquatic Acute 1, Aquatic Chronic 1 H302, H335, H410, H319, H315 P261, P273, P301+P312, P302+P352, P280, P305+P351+P338
Guaranteed analysis		
Identity test		: Pass

Copper(II) nitrate trihydrate
specified reagents for general laboratory work- extra pure

SLR

packaging		Mfr. No
250 g	Glass jar	C/8080/50
500 g	Glass jar	C/8080/53
Cu(NO ₃) ₂ ·3H ₂ O CAS: 10031-43-3 MW: 241.60 Melting point: 114°C ADR/RID: 5.1, II EINECS: 221-838-5 UN: 1477 GHS: Ox. Sol. 2, Acute Tox.		4, Skin Irrit. 2, Eye Irrit. 2, Aquatic Acute 1, Aquatic Chronic 1 H272, H315, H319, H302, H400, H410 P210, P301+P312, P273, P302+P352, P280, P305+P351+P338
Assay: >95%		
Guaranteed analysis		
Calcium (Ca)	<200 ppm	n-Butyl formate <0.005%
Iron (Fe)	<100 ppm	Total phosphorus (P) <100 ppm
Lead (Pb)	<200 ppm	Total silicon (Si) <100 ppm
Magnesium (Mg)	<50 ppm	Total sulphur (S) <200 ppm
Potassium (K)	<50 ppm	Zinc (Zn) <100 ppm
Sodium (Na)	<200 ppm	

Copper(I) oxide
specified reagents for general laboratory work- extra pure

SLR

packaging		Mfr. No
250 g	Plastic jar	C/8840/50
Cu ₂ O CAS: 1317-39-1 MW: 143.08 Density at 20°C: 6.0 Boiling point: 1800°C Melting point: 1232°C ADR/RID: 9, III		EINECS: 215-270-7 UN: 3077 GHS: Acute Tox. 4, Aquatic Acute 1, Aquatic Chronic 1 H410, H302 P301+P312, P273
Guaranteed analysis		
Identity test		: Pass

UHPLC

Filtered to 0.1 µm, Fisher Chemical UHPLC gradient grade solvents (Acetonitrile, Methanol and Water) are perfect for those working at high pressure.

See page 309 for more information

Copper(II) oxide
for analysis

Certified AR

packaging		Mfr. No
250 g	Plastic jar	C/8360/50
CuO CAS: 1317-38-0 MW: 79.54 Melting point: 1326°C ADR/RID: 9, III EINECS: 215-269-1		UN: 3077 GHS: Acute Tox. 4, Aquatic Acute 1 H302, H400 P273, P301+P312
Assay: >98%		
Guaranteed analysis		
Calcium (Ca)	<200 ppm	n-Butyl formate <0.005%
Iron (Fe)	<500 ppm	Total nitrogen (N) <0.005 ppm
Lead (Pb)	<200 ppm	Total phosphorus (P) <100 ppm
Magnesium (Mg)	<50 ppm	Total silicon (Si) <200 ppm
Potassium (K)	<50 ppm	Total sulphur (S) <100 ppm
Sodium (Na)	<200 ppm	Zinc (Zn) <200 ppm
Typical analysis		
Aluminium (Al)	<2 ppm	Manganese (Mn) <2 ppm
Barium (Ba)	<5 ppm	Nickel (Ni) <0.001%
Cadmium (Cd)	<2 ppm	Tin (Sn) <5 ppm
Chromium (Cr)	<2 ppm	Titanium (Ti) <2 ppm
Cobalt (Co)	<2 ppm	Vanadium (V) <2 ppm
Lithium (Li)	<2 ppm	

Lot analysis on the label

Copper(II) oxide
Technical, powder

packaging		Mfr. No
250 g	Glass jar	C/8200/50
500 g	Glass jar	C/8200/53
CuO CAS: 1317-38-0 MW: 79.54 Melting point: 1326°C ADR/RID: 9, III EINECS: 215-269-1		UN: 3077 GHS: Acute Tox. 4, Aquatic Acute 1 H302, H400 P273, P301+P312
Assay: >98%		

Copper(II) oxide
pure, wire form

packaging		Mfr. No
100 g	Plastic jar	C/8320/48
500 g	Plastic jar	C/8320/53
CuO CAS: 1317-38-0 MW: 79.54 Melting point: 1326°C ADR/RID: 9, III EINECS: 215-269-1		UN: 3077 GHS: Acute Tox. 4, Aquatic Acute 1 H302, H400 P273, P301+P312
Guaranteed analysis		
Identity test		: Pass

Copper(II) sulfate anhydrous
specified reagents for general laboratory work- extra pure

SLR

packaging		Mfr. No
250 g	Plastic jar	C/8600/50
1 kg	Plastic jar	C/8600/60
CuSO ₄ CAS: 7758-98-7 MW: 159.61 Density at 20°C: 3.6 Melting point: 200°C ADR/RID: 6.1, III EINECS: 231-847-6		UN: 3288 GHS: Acute Tox. 4, Skin Irrit. 2, Eye Irrit. 2, Aquatic Acute 1, Aquatic Chronic 1 H410, H319, H302, H315 P280, P302+P352, P273, P301+P312, P305+P351+P338
Assay: >98%		
Guaranteed analysis		
Iron (Fe)	<1000 ppm	n-Butyl formate <0.01%
Lead (Pb)	<100 ppm	

Copper(II) sulfate pentahydrate
for analysis

Certified AR

packaging		Mfr. No
250 g	Plastic jar	C/8560/50
500 g	Plastic jar	C/8560/53
1 kg	Plastic jar	C/8560/60
3 kg	Plastic jar	C/8560/63
5 kg	Plastic jar	C/8560/65
25 kg	Plastic keg	C/8560/70
CuSO ₄ ·5H ₂ O		2, Eye Irrit. 2, Aquatic Acute
CAS: 7758-99-8		1, Aquatic Chronic 1
MW: 249.68		H315, H410, H302, H319
Melting point: 110°C		P273, P301+P312,
ADR/RID: 9, III		P302+P352, P280,
UN: 3077		P305+P351+P338
GHS: Acute Tox. 4, Skin Irrit.		
Assay: >99.5%		
Guaranteed analysis		
pH (aqueous solution 5%)	3.5 to 4	Sodium (Na)
Calcium (Ca)	<20 ppm	Total chloride (Cl)
Iron (Fe)	<20 ppm	Total nitrogen (N)
Lead (Pb)	<50 ppm	Total phosphorus (P)
Magnesium (Mg)	<5 ppm	Total silicon (Si)
Potassium (K)	<5 ppm	Zinc (Zn)
Typical analysis		
Aluminium (Al)	<2 ppm	Lithium (Li)
Arsenic (As)	<2 ppm	Manganese (Mn)
Barium (Ba)	<2 ppm	Nickel (Ni)
Cadmium (Cd)	<2 ppm	Titanium (Ti)
Chromium (Cr)	<2 ppm	Vanadium (V)
Cobalt (Co)	<2 ppm	

Lot analysis on the label.

Copper(II) sulfate pentahydrate
specified reagents for general laboratory work- extra pure

SLR

Packaging		Mfr. No	
250 g	Plastic jar	C/8520/50	
500 g	Plastic jar	C/8520/53	
1 kg	Plastic jar	C/8520/60	
3 kg	Plastic jar	C/8520/63	
5 kg	Plastic jar	C/8520/65	
CuSO ₄ ·5H ₂ O		2, Eye Irrit. 2, Aquatic Acute	
CAS: 7758-99-8		1, Aquatic Chronic 1	
MW: 249.68		H315, H410, H302, H319	
Melting point: 110°C		P273, P301+P312,	
ADR/RID: 9, III		P302+P352, P280,	
UN: 3077		P305+P351+P338	
GHS: Acute Tox. 4, Skin Irrit.		 	
Assay: >98%			
Guaranteed analysis			
Calcium (Ca)	<200 ppm	Sodium (Na)	<200 ppm
Iron (Fe)	<100 ppm	n-Butyl formate	<0.005%
Lead (Pb)	<200 ppm	Total phosphorus (P)	<100 ppm
Magnesium (Mg)	<50 ppm	Total silicon (Si)	<100 ppm
Potassium (K)	<50 ppm	Zinc (Zn)	<100 ppm

Copper(II) sulfate pentahydrate
Technical

Packaging		Mfr. No
500 g	Plastic jar	C/8480/53
1 kg	Plastic jar	C/8480/60
3 kg	Plastic jar	C/8480/63
5 kg	Plastic jar	C/8480/65
CuSO ₄ ·5H ₂ O		2, Eye Irrit. 2, Aquatic Acute
CAS: 7758-99-8		1, Aquatic Chronic 1
MW: 249.68		H315, H410, H302, H319
Melting point: 110°C		P273, P301+P312,
ADR/RID: 9, III		P302+P352, P280,
UN: 3077		P305+P351+P338
GHS: Acute Tox. 4, Skin Irrit.		
Guaranteed analysis		
Identity test		: Pass

Creatinine
pure

Packaging		Mfr. No
25 g	Plastic jar	C/7240/46
C ₄ H ₇ N ₃ O		Melting point: 255°C
CAS: 60-27-5		Flash point: +290°C
MW: 113.12		EINECS: 200-466-7
Assay: >98%		

o-Cresol
pure

Packaging		Mfr. No
500 mL	Coated Glass bottle	C/7320/PB08
C ₇ H ₈ O	UN: 3455	
CAS: 95-48-7	GHS: Acute Tox. 3, Skin	
MW: 108.14	Corr. 1B, Eye Dam. 1	
Density at 20°C: 1.04	H301, H314, H311	
Boiling point: 191°C	P301+P310, P280, P312,	
Melting point: 30 to 32°C	P302+P350,	
Flash point: 81°C	P301+P330+P331,	
ADR/RID: 6.1, 8, II	P305+P351+P338, P310	
EINECS: 202-423-8		
Guaranteed analysis		
Others impurities		<2%

m-Cresol
pure

Packaging		Mfr. No
600 mℓ	Coated Glass bottle	C/7360/PB08
2.5 ℓ	Coated Glass bottle	C/7360/PB17
25 ℓ	Mild steel drum	C/7360/25
C ₇ H ₈ O	EINECS: 203-577-9	
CAS: 108-39-4	UN: 2076	
MW: 108.14	GHS: Acute Tox. 3, Skin	
Density at 20°C: 1.03	Corr. 1B, Eye Dam. 1	
Boiling point: 203°C	H311, H314, H301	
Melting point: 8 to 10°C	P280, P312, P302+P350,	
Flash point: 86°C	P301+P330+P331,	
ADR/RID: 6.1, 8, II	P305+P351+P338, P310	
Assay (CG): >98%		
Guaranteed analysis		
Other isomers		<2%

Cresol purple
pure, solution 0.04 %, pH indicator

packaging		Mfr. No	
100 mL	Glass bottle	C/7450L/05	
CAS: 2303-01-7			
EINECS: 218-960-6			

Cresol red
pure, solid, indicator grade

Packaging		Mfr. No
5 g	Glass jar	C/7460/44
C ₂₁ H ₁₈ O ₅ S		Melting point: 290°C
CAS: 1733-12-6		EINECS: 217-064-2
MW: 382.42		
Guaranteed analysis		
pH colour change test		: Pass
		Complies with standard IR spectrum: Pass

Cresol red
pure, solution 0.02 % pH indicator

Packaging		Mfr. No
500 mL	Glass bottle	C/7460L/08
Flash point: 32°C		
ADR/RID: 3, III		
UN: 1170		
GHS: Flam. Liq. 3		
P210, P280, P303+P361+P353, P403+P235		
Guaranteed analysis		
pH colour change test : Pass		

o-Cresolphthaleindi-(methyliminodi-acetic acid)
pure, metal indicator

Packaging		Mfr. No
5 g	Glass jar	C/7440/44
$C_{32}H_{32}N_2O_{12}$		Melting point: 181°C
CAS: 2411-89-4		EINECS: 219-318-8
MW: 636.61		
Guaranteed analysis		
Complies with standard IR spectrum: Pass		

o-Cresolsulfonphthaleindi-(methyl-iminodiacetic acid) sodium salt
pure, Special reagent for metals (Xylenol orange)

Packaging		Mfr. No
1 g	Glass jar	C/7480/43
5 g	Glass jar	C/7480/44
C ₃₁ H ₃₁ N ₂ NaO ₁₃ S		Melting point: 222°C EINECS: 216-553-8
CAS: 1611-35-4		
MW: 694.65		
Guaranteed analysis		
Complies with standard IR spectrum: Pass		

18-Crown-6-ether (1.4.7.10.13.16-hexaoxacyclooctadecane) SLR
specified reagents for general laboratory work- extra pure

Packaging		Mfr. No
10 g	Glass jar	C/7551/45
$C_{12}H_{24}O_6$		
CAS: 17455-13-9		
MW: 264.32		
Melting point: 36 to 40°C		
Flash point: >110°C		
Assay (CG): > 98%		
Guaranteed analysis		
Complies with standard IR spectrum: Pass		

Cupferron
for analysis

Certified AR

Packaging		Mfr. No
100 g	Glass jar	C/7720/48
<chem>C6H9N3O2</chem>		
GHS: Skin Irrit. 2, Eye Irrit. 2, STOT SE 3		
CAS: 135-20-6		
MW: 155.16		
H335, H319, H350, H315		
P201, P308+P313, P261, P302+P352, P280, P305+P351+P338		
Melting point: 153 to 163°C		
ADR/RID: 6.1, III		
EINECS: 205-183-2		
UN: 2811		
Guaranteed analysis		
Extinction at 281 nm (E 1%/1cm) : >515		Residue after ignition <0.05%
Stabilized with ammonium carbonate. Suitable for precipitation of iron.		

Curcumin
pure, crystalline

Packaging		Mfr. No
1 g	Vial	C/8850/43
C₂₁H₂₀O₆ CAS: 458-37-7 MW: 368.39 Melting point: 183°C EINECS: 207-280-5		GHS: Skin Irrit. 2, Eye Irrit. 2, STOT SE 3 H319, H315, H335 P261, P302+P352, P280, P305+P351+P338 
Guaranteed analysis		
Identity test		: Pass

2.5-Cyclohexadiene-1.4-dione
See p-Benzoquinone (dried) , page 56

Cyclohexane NEW PACK.
for HPLC

Certified

Packaging		Mfr. No
1 ℓ	Glass bottle	C/8936/15
2.5 ℓ	Glass bottle	C/8936/17
2.5 ℓ	Coated Glass bottle	C/8936/PB17
10 ℓ	Returnable Drum	C/8936/10RSS
25 ℓ	Returnable Drum	C/8936/25RSS
30 ℓ	Returnable Drum	C/8936/30RSS
C ₆ H ₁₂		UN: 1145
CAS: 110-82-7		GHS: Flam. Liq. 2, Skin Irrit.
MW: 84.16		2, STOT SE 3, Asp. Tox. 1,
Density at 20°C: 0.770		Aquatic Acute 1, Aquatic
Boiling point: 81°C		Chronic 1
Melting point: 6.5°C		H225, H336, H304, H315, H410
Flash point: -18°C		P261, P301+P310, P273,
ADR/RID: 3, II		P304+P340, P302+P352, P210
INECS: 203-806-2		

Cyclohexane
for residue and pesticide analysis

Distol

packaging		Mfr. No
1 ℓ	Glass bottle	C/8933/15
2.5 ℓ	Glass bottle	C/8933/17
C ₆ H ₁₂ CAS: 110-82-7 MW: 84.16 Density at 20°C: 0.770 Boiling point: 81°C Melting point: 6.5°C Flash point: -18°C ADR/RID: 3, II EINECS: 203-806-2		UN: 1145 GHS: Flam. Liq. 2, Skin Irrit. 2, STOT SE 3, Asp. Tox. 1, Aquatic Acute 1, Aquatic Chronic 1 H225, H336, H304, H315, H410 P261, P301+P310, P273, P304+P340, P302+P352, P210
		
		
		
		

Guaranteed analysis			
Water	<0.02%	GC (FID) Analysis	Pass
Residue after evaporation	<2 ppm	GC (NPD) analysis	Pass
GC (ECD) Analysis	Pass		

GC Analysis
Analysis of pesticides: from the retention time of Lindane, after concentration of 2500 times, no impurity peak higher than that obtained with 10 pg / ml héptachlo-réoxide (ECD) and 5 pg / ml Parathion (NPD).
Analysis of hydrocarbons: no impurity chain longer than the dodecane is present at a concentration greater than 0.01 ppm (FID).
Conditions :
Column: J&W DB 5: 30m x 0.32mm x 0.1µm df.
Temperature profile: (1) 50°C, 2 min
(2) 40°C / min until 320°C
(3) 320°C for 30 min
Injection: 5µℓ cold on the column.

Cyclohexane
for analysis

Certified AR

packaging		Mfr. No
1 ℓ	Glass bottle	C/8921/15
2.5 ℓ	Glass bottle	C/8921/17
2.5 ℓ	Coated Glass bottle	C/8921/PB17
5 ℓ	Fluorinated HDPE Plastic bottle	C/8921/21
25 ℓ	Metal drum	C/8921/25
200 ℓ	Metal drum	C/8921/27
C ₆ H ₁₂ CAS: 110-82-7 MW: 84.16 Density at 20°C: 0.770 Boiling point: 81°C Melting point: 6.5°C Flash point: -18°C ADR/RID: 3, II EINECS: 203-806-2		UN: 1145 GHS: Flam. Liq. 2, Skin Irrit. 2, STOT SE 3, Asp. Tox. 1, Aquatic Acute 1, Aquatic Chronic 1 H225, H336, H304, H315, H410 P261, P301+P310, P273, P304+P340, P302+P352, P210
		
		
		
		

Assay (CG): >99.8%			
Guaranteed analysis			
Colour	<10 APHA	Sodium (Na)	<1 ppm
Water	<0.01%	Total phosphorus (P)	<0.1 ppm
Residue after evaporation	<10 ppm	Total silicon (Si)	<0.1 ppm
Calcium (Ca)	<0.5 ppm	Zinc (Zn)	<0.2 ppm
Copper (Cu)	<0.02 ppm	Benzene	<0.005%
Iron (Fe)	<0.2 ppm	Cyclohexene	<0.005%
Lead (Pb)	<0.05 ppm	Substances darkened by H ₂ SO ₄	<20 APHA
Magnesium (Mg)	<0.1 ppm		
Potassium (K)	<0.2 ppm		
Typical analysis			
Aluminium (Al)	<0.1 ppm	Manganese (Mn)	<0.02 ppm
Barium (Ba)	<0.02 ppm	Molybdenum (Mo)	<0.02 ppm
Cadmium (Cd)	<0.05 ppm	Nickel (Ni)	<0.05 ppm
Chromium (Cr)	<0.02 ppm	Strontium (Sr)	<0.02 ppm
Cobalt (Co)	<0.02 ppm	Titanium (Ti)	<0.02 ppm
Lithium (Li)	<0.02 ppm	Vanadium (V)	<0.02 ppm

Lot analysis on the label.

Cyclohexane
specified reagents for general laboratory work- extra pure

SLR

packaging		Mfr. No
500 mℓ	Coated Glass bottle	C/8920/PB08
2.5 ℓ	Glass bottle	C/8920/17
2.5 ℓ	Coated Glass bottle	C/8920/PB17
5 ℓ	Fluorinated HDPE Plastic bottle	C/8920/21
5 ℓ	Fluorinated HDPE square Drum	C/8920/FP21
25 ℓ	Metal drum	C/8920/25
C ₆ H ₁₂ CAS: 110-82-7 MW: 84.16 Density at 20°C: 0.770 Boiling point: 81°C Melting point: 6.5°C Flash point: -18°C ADR/RID: 3, II EINECS: 203-806-2		UN: 1145 GHS: Flam. Liq. 2, Skin Irrit. 2, STOT SE 3, Asp. Tox. 1, Aquatic Acute 1, Aquatic Chronic 1 H225, H336, H304, H315, H410 P261, P301+P310, P273, P304+P340, P302+P352, P210
		
		
		
		

Assay (CG): >99%			
Guaranteed analysis			
Colour	<20 APHA	Magnesium (Mg)	<1 ppm
Water	<0.05%	Potassium (K)	<2 ppm
Residue after evaporation	<50 ppm	Sodium (Na)	<5 ppm
Calcium (Ca)	<2 ppm	Total phosphorus (P)	<2 ppm
Copper (Cu)	<2 ppm	Total silicon (Si)	<2 ppm
Iron (Fe)	<5 ppm	Zinc (Zn)	<2 ppm
Lead (Pb)	<5 ppm		

Cyclohexanol
pure

packaging		Mfr. No
500 mℓ	Coated Glass bottle	C/9000/PB08
2.5 ℓ	Coated Glass bottle	C/9000/PB17
C ₆ H ₁₂ O CAS: 108-93-0 MW: 100.16 Density at 20°C: 0.960 Boiling point: 161°C Melting point: 23°C Flash point: 67°C		EINECS: 203-630-6 GHS: Acute Tox. 4, Skin Irrit. 2, STOT SE 3 H315, H302, H335, H332 P310, P304+P340, P261, P301+P312, P302+P352
		
Assay (CG): >98%		

Cyclohexanone
for analysis

Certified AR

Packaging		Mfr. No	
600 mℓ	Coated Glass bottle	C/9050/PB08	
2.5 ℓ	Coated Glass bottle	C/9050/PB17	
$C_6H_{10}O$ CAS: 108-94-1 MW: 98.15 Density at 20°C: 0.947 Boiling point: 155°C Melting point: -47°C Flash point: 43°C		ADR/RID: 3, III EINECS: 203-631-1 UN: 1915 GHS: Flam. Liq. 3, Acute Tox. 4 H226, H332 P280, P302+P350, P304+P340, P210	
			
			
Assay (CG): >99.5%			
Guaranteed analysis			
Acidity / Alkalinity	<0.003 meq/g	Lead (Pb)	<0.1 ppm
Colour	<10 APHA	Total phosphorus (P)	<2 ppm
Water	<0.2%	Total sulphur (S)	<5 ppm
Copper (Cu)	<0.05 ppm	Cyclohexanol	<0.2%
Iron (Fe)	<0.5 ppm	Cyclohexane	<0.01%

Lot analysis on the label.

Cyclohexanone
specified reagents for general laboratory work- extra pure

SLR

packaging		Mfr. No	
500 mℓ	Coated Glass bottle	C/9040/PB08	
2.5 ℓ	Coated Glass bottle	C/9040/PB17	
5 ℓ	Fluorinated HDPE Plastic bottle	C/9040/21	
25 ℓ	Metal drum	C/9040/25	
C ₆ H ₁₀ O		ADR/RID: 3, III	
CAS: 108-94-1		EINECS: 203-631-1	
MW: 98.15		UN: 1915	
Density at 20°C: 0.947		GHS: Flam. Liq. 3, Acute Tox. 4	
Boiling point: 155°C		H226, H332	
Melting point: -47°C		P280, P302+P350,	
Flash point: 43°C		P304+P340, P210	
Assay (CG): >99%			
Guaranteed analysis			
Acidity / Alkalinity	<0.03 meq/g	Potassium (K)	<2 ppm
Colour	<20 APHA	Sodium (Na)	<50 ppm
Calcium (Ca)	<5 ppm	Total phosphorus (P)	<10 ppm
Copper (Cu)	<1 ppm	Total silicon (Si)	<1 ppm
Iron (Fe)	<2 ppm	Total sulphur (S)	<20 ppm
Lead (Pb)	<1 ppm	Zinc (Zn)	<1 ppm
Magnesium (Mg)	<1 ppm		

Cyclohexene
specified reagents for general laboratory work- extra pure

SLR

packaging		Mfr. No
250 mℓ	Glass bottle	C/9080/07
500 mℓ	Glass bottle	C/9080/08
C ₆ H ₁₀ CAS: 110-83-8 MW: 82.15 Density at 20°C: 0.810 Boiling point: 83°C Melting point: -104°C Flash point: -20°C ADR/RID: 3, II		EINECS: 203-807-8 UN: 2256 GHS: Flam. Liq. 2, Acute Tox. 4, Asp. Tox. 1, Aquatic Chronic 2 H225, H304, H312, H302, H411 P280, P302+P352, P301+P310, P331, P210, P240
		
		
		
		
Assay (CG): >99%		

L-Cysteine
pure

packaging			Mfr. No
25 g	Plastic jar		C/9150/46
C ₃ H ₇ NO ₂ S		2, Eye Irrit. 2, STOT SE 3	
CAS: 52-90-4		H302, H319, H315, H335	
MW: 121.16		P261, P301+P312,	
Melting point: 220°C		P302+P352, P280,	
EINECS: 200-158-2		P305+P351+P338	
GHS: Acute Tox. 4, Skin Irrit.			
Assay: >98%			
Guaranteed analysis			
Specific rotation	6.5° to 8°	Iron (Fe)	<10 ppm
Residue after ignition	<0.1%	Lead (Pb)	<10 ppm

High Volume Solvent Delivery Systems

High-volume solvent delivery systems are available in 10L to 1000L, in Europe only. Please email your contact details and a brief description of your requirements to begei.shi.quotes@thermofisher.com



L-Cysteine hydrochloride monohydrate
pure

Packaging			Mfr. No
100 g	Plastic jar		C/9152/48
$C_3H_7NO_2S.HCl.H_2O$ CAS: 7048-04-6 MW: 175.63 Melting point: 175 to 177.8°C EINECS: 200-157-7		GHS: Skin Irrit. 2, Eye Irrit. 2, STOT SE 3 H315, H319, H335 P280, P305+P351+P338	
Assay: >98.5%			
Guaranteed analysis			
Specific rotation	5.5° to 7°	Iron (Fe)	<30 ppm
Residue after ignition	<0.1%	Lead (Pb)	<10 ppm
Ammonium (NH ₄)	<0.03%	Sulfate (SO ₄)	<0.02%
Arsenic (As)	<1 ppm		

L-Cystine
pure

Packaging		Mfr. No	
25 g	Plastic jar	C/9160/46	
C₆H₁₂N₂O₄S₂ CAS: 56-89-3 MW: 240.29 Melting point: 260 to 261°C EINECS: 200-296-3 GHS: Acute Tox. 4, Skin Irrit. 2, Eye Irrit. 2, STOT SE 3 H315, H319, H302, H335 P261, P301+P312, P302+P352, P280, P305+P351+P338			
Assay: >99%			
Guaranteed analysis			
Specific rotation	<-210°	Iron (Fe)	<20 ppm
Residue after ignition	<0.5%	Lead (Pb)	<10 ppm

Dansyl chloride
specified reagents for general laboratory work- extra pure

SLR

packaging		Mfr. No
5 g	Plastic jar	D/0008/44
C ₁₂ H ₁₂ ClNO ₂ S CAS: 605-65-2 MW: 269.75 Melting point: 71 to 74°C ADR/RID: 8, II EINECS: 210-092-6 UN: 3261		GHS: Acute Tox. 4, Skin Corr. 1B, Eye Dam. 1 H302, H314, EUH029 P301+P330+P331, P280, P305+P351+P338, P310, P402+P404, P301+P312
		
		
Assay: >99%		

Decahydronaphthalene
Technical

packaging		Mfr. No
2.5 ℓ	Glass bottle	D/0050/17
C ₁₀ H ₁₈ CAS: 91-17-8 MW: 138.25 Density at 20°C: 0.890 Boiling point: 187°C Melting point: -31°C Flash point: 58°C ADR/RID: 3, III EINECS: 202-046-9		UN: 1147 GHS: Flam. Liq. 3, Acute Tox. 3, Skin Irrit. 2, Eye Irrit. 2, STOT SE 3 H226, H319, H331, H335, H315 P304+P340, P261, P302+P352, P280, P305+P351+P338, P210
		
		

UHPLC

Filtered to 0.1 µm, Fisher Chemical UHPLC gradient grade solvents (Acetonitrile, Methanol and Water) are perfect for those working at high pressure. See page 309 for more information



n-Decane		
Technical		
packaging	Mfr. No	
100 mℓ Glass bottle	D/0020/05	
500 mℓ Glass bottle	D/0020/08	
C ₁₀ H ₂₂ CAS: 124-18-5 MW: 142.29 Density at 20°C: 0.735 Boiling point: 174°C Melting point: -30°C Flash point: 46°C	ADR/RID: 3, III EINECS: 204-686-4 UN: 2247 GHS: Flam. Liq. 3, Asp. Tox. 1 H226, H304, EUH066 P301+P310, P331, P210	 
Assay (CG): >98%		

1-Decanesulfonic acid sodium salt		Certified
for Ion pair chromatography with electrochemical detection		
packaging	Mfr. No	
5 g Glass jar	D/0035/44	
C ₁₀ H ₂₁ NaO ₃ S CAS: 13419-61-9 MW: 244.31	Melting point: >300°C EINECS: 236-525-9	
Assay: >98%		
Guaranteed analysis		
Electrochemical test Pass		
Absorbance specification		
Max. Absorbance (1cm cell v HPLC water)		
254 nm (sol. aqu. 0.005 M) : <0.05 A.U.		
Lot analysis on the label.		
EC Specification:		
Tested by HPLC-EC to ensure low electrochemical response.		
Detector - ESA Coulochem II		
Mode-Oxidative		
Electrode graphite		

1-Decanesulfonic acid sodium salt		Certified
for Ion pair chromatography		
packaging	Mfr. No	
25 g Glass jar	D/0030/46	
C ₁₀ H ₂₁ NaO ₃ S CAS: 13419-61-9 MW: 244.31	Melting point: >300°C EINECS: 236-525-9	
Assay: >98%		
Absorbance specification		
Max. Absorbance (1cm cell v HPLC water)		
254 nm (sol. aqu. 0.005 M) : <0.05 A.U.		
Lot analysis on the label.		

Dekalin		
pure		
packaging	Mfr. No	
2.5 ℓ Glass bottle	D/0052/17	
C ₁₀ H ₁₈ CAS: 91-17-8 MW: 138.25 Density at 20°C: 0.890 Boiling point: 187°C Melting point: -31°C Flash point: 58°C ADR/RID: 3, III EINECS: 202-046-9	UN: 1147 GHS: Flam. Liq. 3, Acute Tox. 3, Skin Irrit. 2, Eye Irrit. 2, STOT SE 3 H226, H319, H331, H335, H315 P304+P340, P261, P302+P352, P280, P305+P351+P338, P210	 
Assay (CG): >98%		

Devarda's alloy		
pure, Silylation reagent		
packaging	Mfr. No	
250 g Plastic jar	D/0100/50	
CAS: 8049-11-4 Density at 20°C: 5.8 Melting point: 550°C ADR/RID: 4.1, II	UN: 3089 GHS: Flam. Sol. 2 H228 P210, P240, P241, P280	
Guaranteed analysis		
Identity test : Pass		

Dextran 150		
pure, M.W. 100.000-200.000		

packaging	Mfr. No
100 g Plastic jar	D/0122/48
CAS: 9004-54-0	
EINECS: 232-677-5	
Guaranteed analysis	
Identity test	: Pass

Dextran 70		
pure, M.W. 60.000-80.000		

packaging	Mfr. No
100 g Plastic jar	D/0124/48
CAS: 9004-54-0	
EINECS: 232-677-5	
Guaranteed analysis	
Complies with standard IR spectrum: Pass	

Dextrin		
Technical		

packaging		Mfr. No
100 g	Plastic jar	D/0130/48
C ₆ H ₁₀ O ₅ CAS: 9004-53-9 MW: 162.07 EINECS: 232-675-4		
Guaranteed analysis		
Complies with standard IR spectrum: Pass		
Conforming to the IR spectrum of reference.		

Dextrose See D(+)-Glucose anhydrous, page 135

Diacetone alcohol		SLR
specified reagents for general laboratory work- extra pure		

packaging	Mfr. No	
500 mℓ Glass bottle	D/0150/08	
C ₆ H ₁₂ O ₂ CAS: 123-42-2 MW: 116.16 Density at 20°C: 0.939 Boiling point: 166°C Melting point: -42.8°C Flash point: 58°C	ADR/RID: 3, III EINECS: 204-626-7 UN: 1148 GHS: Eye Irrit. 2 H319 P264, P280, P305+P351+P338, P337+P313	
Guaranteed analysis		
Water <0.2% Residue after evaporation <200 ppm		
Refractive index at 20°C : <1.425		

Diacetyl		SLR
specified reagents for general laboratory work- extra pure		
packaging	Mfr. No	
100 mℓ Glass bottle	D/0250/05	
C ₄ H ₆ O ₂ CAS: 431-03-8 MW: 86.09 Density at 20°C: 0.985 Boiling point: 88°C Melting point: -4 to -2°C Flash point: 7°C ADR/RID: 3, II	EINECS: 207-069-8 UN: 2346 GHS: Flam. Liq. 2, Acute Tox. 3, Skin Irrit. 2, Eye Irrit. 2 H225, H302, H315, H319, H331 P280, P302+P352, P261, P310, P305+P351+P338, P210	 
Assay (CG): >98%		

1.2-Diaminoethane anhydrous		SLR
specified reagents for general laboratory work- extra pure		

packaging	Mfr. No	
500 mℓ Coated Glass bottle	D/0370/PB08	
2.5 ℓ Coated Glass bottle	D/0370/PB17	
25 ℓ Metal drum	D/0370/25	
C ₂ H ₈ N ₂ CAS: 107-15-3 MW: 60.10 Density at 20°C: 0.898 Boiling point: 117 to 118°C Melting point: 11°C Flash point: 38°C ADR/RID: 8, 3, II EINECS: 203-468-6 UN: 1604	GHS: Flam. Liq. 3, Acute Tox. 3, Skin Corr. 1B, Eye Dam. 1, Resp. Sens. 1, Skin Sens. 1 H226, H314, H302, H311, H334, H317 P280, P301+P330+P331, P310, P305+P351+P338, P303+P361+P353, P210	   

Assay: >98%

1.2-Di(2-aminoethoxy)ethane-N,N,N,N-tetraacetic acid		
pure		

packaging	Mfr. No	
25 g Plastic jar	D/0755/46	
C ₁₄ H ₂₄ N ₂ O ₁₀ CAS: 67-42-5 MW: 380.35	Melting point: 240 to 245°C Flash point: >200°C EINECS: 200-651-2	
Assay: >98%		

1.6-Diaminohexane		SLR
specified reagents for general laboratory work- extra pure		

packaging	Mfr. No	
25 g Glass jar	D/0760/46	
C ₆ H ₁₆ N ₂ CAS: 124-09-4 MW: 116.21 Boiling point: 204 to 205°C Melting point: 38 to 41°C Flash point: 81°C ADR/RID: 8, III EINECS: 204-679-6	UN: 2280 GHS: Acute Tox. 4, Skin Corr. 1B, Eye Dam. 1, STOT SE 3 H312, H314, H335, H302 P280, P301+P330+P331, P305+P351+P338, P310, P302+P352, P261	 
Assay: 99 to 101%		

1.6-Diaminohexane		
pure, solution 5% in water		

packaging	Mfr. No	
250 mℓ Coated Glass bottle	D/0765/PB07	
ADR/RID: 8, III UN: 3267 GHS: Skin Corr. 1B, Eye Dam. 1	H314 P301+P330+P331, P280, P305+P351+P338, P310	

L(+)-2.6-Diaminohexanoic acid hydrochloride	
See L-Lysine monohydrochloride, page 186	

3.7-Diamino-5-phenylphenazinium chloride	
See Phenosafranine, page 241	

Diastase	
pure, from malt, mixed alpha and beta amylase	

packaging	Mfr. No	
100 g Plastic jar	D/0800/48	
CAS: 9000-92-4 EINECS: 232-567-7 GHS: Resp. Sens. 1	H334 P261, P304+P341, P342+P311	
Guaranteed analysis		
Activity <4 EU/mg Reducing sugar Pass		

1.2-Dibromoethane		SLR
specified reagents for general laboratory work- extra pure		

packaging	Mfr. No	
500 mℓ Coated Glass bottle	D/1250/PB08	
C ₂ H ₄ Br ₂ CAS: 106-93-4 MW: 187.88 Density at 20°C: 2.173 Boiling point: 131 to 132°C Melting point: 9 to 10°C Flash point: >104°C ADR/RID: 6.1, I EINECS: 203-444-5	UN: 1605 GHS: Carc. 1B, Acute Tox. 3, Skin Irrit. 2, Eye Irrit. 2, STOT SE 3, Aquatic Chronic 2 H319, H315, H350, H411, H311, H331, H335, H301 P280, P302+P350, P308+P313, P304+P340, P261, P273, P305+P351+P338	  

Assay (CG): >98%

Guaranteed analysis	
Residue after evaporation <200 ppm	

5',5''-Dibromophenolsulfonphthalein	
See Bromophenol red, page 62	

3',3''-Dibromothymolsulfonephthalein	
See Bromothymol blue, page 62	

Di-n-butylamine		Certified
for HPLC		

packaging	Mfr. No	
250 mℓ Coated Glass bottle	D/1338/PB07	
C ₈ H ₁₉ N CAS: 111-92-2 MW: 129.25 Density at 20°C: 0.760 Boiling point: 159°C Melting point: -62°C Flash point: 39°C ADR/RID: 8, 3, II EINECS: 203-921-8	UN: 2248 GHS: Flam. Liq. 3, Acute Tox. 2, Skin Corr. 1B, Eye Dam. 1 H226, H311, H314, H330, H301 P280, P312, P302+P350, P301+P330+P331, P305+P351+P338, P310, P304+P340, P210	  

Assay (CG): >99%

Absorbance specification	
Max. Absorbance (1cm cell v HPLC water)	
254 nm (sol. aqu. 0.02 M) ... : <0.01 A.U.	
Lot analysis on the label.	
Filtered to 0.2 μm.	

Di-n-butylamine <i>pure</i>		
packaging	Mfr. No	
500 mℓ Coated Glass bottle	D/1340/PB08	
C ₈ H ₁₉ N CAS: 111-92-2 MW: 129.25 Density at 20°C: 0.760 Boiling point: 159°C Melting point: -62°C Flash point: 39°C ADR/RID: 8, 3, II EINECS: 203-921-8	UN: 2248 GHS: Flam. Liq. 3, Acute Tox. 2, Skin Corr. 1B, Eye Dam. 1 H226, H311, H314, H330, H301 P280, P312, P302+P350, P301+P330+P331, P305+P351+P338, P310, P304+P340, P210	  
Assay (CG): >98%		

Di-n-butylamine phosphate <i>for HPLC, solution 1.0 M</i>		Certified
packaging	Mfr. No	
100 mℓ Glass bottle	D/1341/05	
(CH ₃ (CH ₂) ₃) ₂ NH.H ₃ PO ₄ Density at 20°C: 1.06		
Absorbance specification Max. Absorbance (1cm cell v HPLC water) 254 nm <0.05 A.U. <i>Filtered to 0.2 µm.</i>		

2.6-Di-tert-butyl-p-cresol <i>pure</i>		
packaging	Mfr. No	
500 g Plastic jar	D/1345/53	
C ₁₅ H ₂₄ O CAS: 128-37-0 MW: 220.36 Boiling point: 265°C Melting point: 69 to 71°C Flash point: 127°C ADR/RID: 9, III EINECS: 204-881-4	UN: 3077 GHS: Acute Tox. 4, Skin Irrit. 2, Eye Irrit. 2, STOT SE 3, Aquatic Acute 1, Aquatic Chronic 1 H315, H335, H410, H319, H302 P261, P273, P302+P352, P280, P305+P351+P338	 
Guaranteed analysis Complies with standard IR spectrum: Pass		

Di-n-butyl ether <i>pure</i>		
packaging	Mfr. No	
2.5 ℓ Glass bottle 25 ℓ Metal drum 200 ℓ Metal drum	D/1350/17 D/1350/25 D/1350/27	
C ₈ H ₁₈ O CAS: 142-96-1 MW: 130.23 Density at 20°C: 0.768 Boiling point: 141°C Melting point: -95°C Flash point: 25°C ADR/RID: 3, III	EINECS: 205-575-3 UN: 1149 GHS: Flam. Liq. 3, Skin Irrit. 2, Eye Irrit. 2, STOT SE 3, Aquatic Chronic 3 H226, H319, H412, H335, H315 P261, P273, P302+P352, P280, P305+P351+P338, P210	 
Assay (CG): >99%		

UHPLC
Filtered to 0.1 µm, Fisher Chemical UHPLC gradient grade solvents (Acetonitrile, Methanol and Water) are perfect for those working at high pressure.
See page 309 for more information



Di-n-butyl phthalate <i>specified reagents for general laboratory work- extra pure</i>			SLR
packaging	Mfr. No		
500 mℓ Glass bottle 2.5 ℓ Glass bottle 200 ℓ Metal drum	D/1500/08 D/1500/17 D/1500/27		
C ₁₆ H ₂₂ O ₄ CAS: 84-74-2 MW: 278.35 Density at 20°C: 1.043 Boiling point: 340°C Melting point: -35°C Flash point: 188°C	ADR/RID: 9, III EINECS: 201-557-4 UN: 3082 GHS: Repr. 1B, Aquatic Acute 1 H360Df, H400 P201, P308+P313, P273, P304+P340	 	
Assay (CG): >99%			
Guaranteed analysis Acidity <0.003 meq/g Residue after ignition <0.01%			

Dichloroacetic acid <i>specified reagents for general laboratory work- extra pure</i>			SLR
packaging	Mfr. No		
100 mℓ Coated Glass bottle 1 ℓ Coated Glass bottle	D/1550/PB05 D/1550/PB15		
C ₂ H ₂ Cl ₂ O ₂ CAS: 79-43-6 MW: 128.94 Density at 20°C: 1.56 Boiling point: 194°C Freezing point: 11 to 13°C Melting point: 9 to 11°C Flash point: >112°C ADR/RID: 8, II	EINECS: 201-207-0 UN: 1764 GHS: Skin Corr. 1A, Aquatic Acute 1 H314, H400 P280, P312, P302+P350, P301+P330+P331, P305+P351+P338, P273	 	
Assay: >99%			
Guaranteed analysis Sulfate (SO ₄) <0.001% <i>Suitable for DNA synthesis.</i>			

1.2-Dichlorobenzene <i>for analysis</i>			Certified AR
packaging	Mfr. No		
500 mℓ Glass bottle	D/1610/08		
C ₆ H ₄ Cl ₂ CAS: 95-50-1 MW: 147.01 Boiling point: 179 to 180°C Melting point: -15°C Flash point: 67°C ADR/RID: 6.1, III EINECS: 202-425-9 UN: 1591	GHS: Acute Tox. 4, Skin Irrit. 2, Eye Irrit. 2, STOT SE 3, Aquatic Acute 1, Aquatic Chronic 1 H319, H302, H335, H410, H315 P261, P273, P301+P312, P302+P352, P280, P305+P351+P338	 	
Assay (CG): >99.5%			
Guaranteed analysis Acidity / Alkalinity <0.0005 meq/g Potassium (K) <0.2 ppm Colour <10 APHA Sodium (Na) <1 ppm Water <0.02% Total phosphorus (P) <0.2 ppm Residue after evaporation <20 ppm Total silicon (Si) <0.2 ppm Calcium (Ca) <1 ppm Total sulphur (S) <0.5 ppm Copper (Cu) <0.05 ppm Zinc (Zn) <0.5 ppm Iron (Fe) <0.5 ppm 1.4-Dichlorobenzene <0.1% Lead (Pb) <0.05 ppm 1.3-Dichlorobenzene <0.1% Magnesium (Mg) <0.2 ppm			
<i>Lot analysis on the label.</i>			

pH Buffers
NIST Certified Standard Solution & Concentrate
A full range to discover on page 63



1.2-Dichlorobenzene <i>specified reagents for general laboratory work- extra pure</i>			SLR
packaging	Mfr. No		
1 ℓ Glass bottle 2.5 ℓ Glass bottle 2.5 ℓ Coated Glass bottle	D/1600/15 D/1600/17 D/1600/PB17		
C ₆ H ₄ Cl ₂ CAS: 95-50-1 MW: 147.01 Density at 20°C: 1.30 Boiling point: 180°C Melting point: -15°C Flash point: 67°C ADR/RID: 6.1, III EINECS: 202-425-9	UN: 1591 GHS: Acute Tox. 4, Skin Irrit. 2, Eye Irrit. 2, STOT SE 3, Aquatic Acute 1, Aquatic Chronic 1 H319, H302, H335, H410, H315 P261, P273, P301+P312, P302+P352, P280, P305+P351+P338	 	
Assay (CG): >98%			
Guaranteed analysis Acidity / Alkalinity <0.005 meq/g Magnesium (Mg) <1 ppm Colour <20 APHA Potassium (K) <1 ppm Water <0.1% Sodium (Na) <5 ppm Residue after evaporation <200 ppm Total phosphorus (P) <1 ppm Calcium (Ca) <5 ppm Total silicon (Si) <1 ppm Copper (Cu) <1 ppm Total sulphur (S) <2 ppm Iron (Fe) <2 ppm Zinc (Zn) <2 ppm Lead (Pb) <1 ppm Other isomers <2%			

1.2-Dichloroethane <i>for HPLC</i>			Certified
packaging	Mfr. No		
1 ℓ Glass bottle 2.5 ℓ Glass bottle	D/1756/15 D/1756/17		
C ₂ H ₄ Cl ₂ CAS: 107-06-2 MW: 98.96 Density at 20°C: 1.25 Boiling point: 81 to 85°C Melting point: -35°C Flash point: 13°C ADR/RID: 3, 6.1, II EINECS: 203-458-1	UN: 1184 GHS: Flam. Liq. 2, Carc. 1B, Acute Tox. 4, Skin Irrit. 2, Eye Irrit. 2, STOT SE 3 H225, H350, H319, H335, H315, H302 P308+P313, P261, P302+P352, P280, P305+P351+P338, P210	  	
Assay (CG): >99.5%			
Guaranteed analysis Acidity / Alkalinity <0.0015 meq/g Residue after evaporation <2 ppm Water <0.02%			
Absorbance specification Max. Absorbance (1cm cell v HPLC water) 240 nm <0.15 A.U. 260 nm <0.005 A.U. 250 nm <0.05 A.U. 270 nm <0.005 A.U.			
<i>Lot analysis on the label.</i> <i>Filtered to 0.2 µm.</i>			

1.2-Dichloroethane <i>for analysis</i>			Certified AR
packaging	Mfr. No		
1 ℓ Glass bottle 2.5 ℓ Glass bottle 200 ℓ Metal drum	D/1751/15 D/1751/17 D/1751/27		
C ₂ H ₄ Cl ₂ CAS: 107-06-2 MW: 98.96 Density at 20°C: 1.25 Boiling point: 81 to 85°C Melting point: -35°C Flash point: 13°C ADR/RID: 3, 6.1, II EINECS: 203-458-1	UN: 1184 GHS: Flam. Liq. 2, Carc. 1B, Acute Tox. 4, Skin Irrit. 2, Eye Irrit. 2, STOT SE 3 H225, H350, H319, H335, H315, H302 P308+P313, P261, P302+P352, P280, P305+P351+P338, P210	  	
Assay (CG): >99.5%			

Guaranteed analysis			
Acidity / Alkalinity <0.0005 meq/g	Magnesium (Mg) <0.1 ppm		
Colour <10 APHA	Potassium (K) <0.2 ppm		
Water <0.03%	Sodium (Na) <0.2 ppm		
Residue after evaporation <20 ppm	Substances reducing KMnO ₄ <0.001%		
Calcium (Ca) <0.5 ppm	Total phosphorus (P) <0.2 ppm		
Free Chlorine (Cl ₂) <1 ppm	Total silicon (Si) <0.05 ppm		
Copper (Cu) <0.05 ppm	Total sulphur (S) <0.5 ppm		
Iron (Fe) <1 ppm	Zinc (Zn) <0.2 ppm		
Lead (Pb) <0.02 ppm			
Typical analysis			
Aluminium (Al) <0.02 ppm	Manganese (Mn) <0.02 ppm		
Barium (Ba) <0.02 ppm	Molybdenum (Mo) <0.02 ppm		
Cadmium (Cd) <0.02 ppm	Nickel (Ni) <0.02 ppm		
Chromium (Cr) <0.02 ppm	Strontium (Sr) <0.01 ppm		
Cobalt (Co) <0.02 ppm	Titanium (Ti) <0.01 ppm		
Lithium (Li) <0.01 ppm	Vanadium (V) <0.01 ppm		
<i>Lot analysis on the label.</i>			

1.2-Dichloroethane <i>specified reagents for general laboratory work- extra pure</i>			SLR
packaging	Mfr. No		
500 mℓ Glass bottle 2.5 ℓ Glass bottle 5 ℓ Fluorinated HDPE Plastic bottle 25 ℓ Metal drum	D/1750/08 D/1750/17 D/1750/21 D/1750/25		
C ₂ H ₄ Cl ₂ CAS: 107-06-2 MW: 98.96 Density at 20°C: 1.25 Boiling point: 81 to 85°C Melting point: -35°C Flash point: 13°C ADR/RID: 3, 6.1, II EINECS: 203-458-1	UN: 1184 GHS: Flam. Liq. 2, Carc. 1B, Acute Tox. 4, Skin Irrit. 2, Eye Irrit. 2, STOT SE 3 H225, H350, H319, H335, H315, H302 P308+P313, P261, P302+P352, P280, P305+P351+P338, P210	  	

Assay (CG): >99%			
Guaranteed analysis Acidity / Alkalinity <0.0005 meq/g Iron (Fe) <2 ppm Colour <20 APHA Lead (Pb) <1 ppm Water <0.1% Magnesium (Mg) <0.5 ppm Residue after evaporation <50 ppm Potassium (K) <1 ppm Free Chlorine (Cl ₂) <10 ppm Sodium (Na) <2 ppm Substances reducing KMnO ₄ <0.01% Total phosphorus (P) <1 ppm Calcium (Ca) <0.02 ppm Total silicon (Si) <0.2 ppm Copper (Cu) <1 ppm Total sulphur (S) <2 ppm Zinc (Zn) <1 ppm			

Dichloro (R) fluorescein <i>Technical</i>		
packaging	Mfr. No	
5 g Glass jar	D/1820/44	
C ₂₀ H ₁₀ Cl ₂ O ₃ CAS: 76-54-0 MW: 401.20 Density at 20°C: 0.790 Melting point: 280°C EINECS: 200-968-6	GHS: Skin Irrit. 2, Eye Irrit. 2, STOT SE 3 H315, H319, H335 P261, P302+P352, P280, P305+P351+P338	
Guaranteed analysis Complies with standard IR spectrum: Pass		

Optima LC/MS
Fisher Chemical OPTIMA® LC/MS solvents, blends, additives and modifiers have set the standard of excellence for consistent, reproducible performance in the mobile phase of LC/MS.
Find out more on page 308



Dichloromethane unstabilised **NEW PACK.** Certified for HPLC

packaging		Mfr. No
1 ℓ	Glass bottle	D/1857/15
2.5 ℓ	Glass bottle	D/1857/17
2.5 ℓ	Coated Glass bottle	D/1857/PB17
10 ℓ	Returnable Drum	D/1857/10RSS
25 ℓ	Returnable Drum	D/1857/25RSS
30 ℓ	Returnable Drum	D/1857/30RSS
CH ₂ Cl ₂ CAS: 75-09-2 MW: 84.93 Density at 20°C: 1.325 Boiling point: 39 to 40°C Melting point: -97°C	ADR/RID: 6.1, III EINECS: 200-838-9 UN: 1593 GHS: Carc. 2 H351 P281	
Assay (CG): >99.8%		
Guaranteed analysis		
Acidity / Alkalinity	<0.0001 meq/g	Residue after evaporation
Water	<0.02%	Free Chlorine (Cl ₂)
Absorbance specification		
Max. Absorbance (1cm cell v HPLC water)		
240 nm	<0.2 A.U.	270 nm
250 nm	<0.02 A.U.	280 nm
260 nm	<0.005 A.U.	Wavelength at 1 A.U.
Lot analysis on the label.		
Filtered to 0.2 µm.		

Dichloromethane stabilised with amylene **NEW PACK.** Certified for HPLC

packaging		Mfr. No
1 ℓ	Glass bottle	D/1856/15
2.5 ℓ	Glass bottle	D/1856/17
2.5 ℓ	Coated Glass bottle	D/1856/PB17
5 ℓ	Aluminium can	D/1856/21
10 ℓ	Returnable Drum	D/1856/10RSS
25 ℓ	Returnable Drum	D/1856/25RSS
30 ℓ	Returnable Drum	D/1856/30RSS
CH ₂ Cl ₂ CAS: 75-09-2 MW: 84.93 Density at 20°C: 1.325 Boiling point: 39 to 40°C Melting point: -97°C	ADR/RID: 6.1, III EINECS: 200-838-9 UN: 1593 GHS: Carc. 2 H351 P281	
Assay (CG): >99.8%		
Guaranteed analysis		
Acidity / Alkalinity	<0.0001 meq/g	Residue after evaporation
Water	<0.02%	Free Chlorine (Cl ₂)
Absorbance specification		
Max. Absorbance (1cm cell v HPLC water)		
233 nm	<1 A.U.	270 nm
240 nm	<0.2 A.U.	280 nm
250 nm	<0.02 A.U.	Wavelength at 1 A.U.
260 nm	<0.005 A.U.	
Lot analysis on the label.		
Filtered to 0.2 µm.		

Dichloromethane stabilised with methanol Certified for HPLC

packaging		Mfr. No
1 ℓ	Glass bottle	D/1859/15
2.5 ℓ	Glass bottle	D/1859/17
CH ₂ Cl ₂ CAS: 75-09-2 MW: 84.93 Density at 20°C: 1.325 Boiling point: 39 to 40°C Melting point: -97°C	ADR/RID: 6.1, III EINECS: 200-838-9 UN: 1593 GHS: Carc. 2 H351 P281	
Assay (CG): >99.8%		
Guaranteed analysis		
Acidity / Alkalinity	<0.0001 meq/g	Residue after evaporation
Water	<0.02%	Free Chlorine (Cl ₂)


Absorbance specification			
Max. Absorbance (1cm cell v HPLC water)			
240 nm	<0.2 A.U.	270 nm	<0.005 A.U.
250 nm	<0.02 A.U.	280 nm	<0.005 A.U.
260 nm	<0.005 A.U.		
Lot analysis on the label.			
Filtered to 0.2 µm.			

Dichloromethane stabilised with amylene Certified for Gel Permeation Chromatography (GPC)

packaging		Mfr. No
2.5 ℓ	Coated Glass bottle	D/1861/PB17
CH ₂ Cl ₂ CAS: 75-09-2 MW: 84.93 Density at 20°C: 1.325 Boiling point: 39 to 40°C Melting point: -97°C	ADR/RID: 6.1, III EINECS: 200-838-9 UN: 1593 GHS: Carc. 2 H351 P281	
Guaranteed analysis		
Colour	<5 APHA	Residue after evaporation
Water	<0.02%	

Dichloromethane stabilised with amylene Primar for trace metal analysis

packaging		Mfr. No
2.5 ℓ	Glass bottle	D/1855/17
CH ₂ Cl ₂ CAS: 75-09-2 MW: 84.93 Density at 20°C: 1.325 Boiling point: 39 to 40°C Melting point: -97°C	ADR/RID: 6.1, III EINECS: 200-838-9 UN: 1593 GHS: Carc. 2 H351 P281	
Assay (CG): >99.9%		
Guaranteed analysis		
Acidity / Alkalinity	<0.0002 meq/g	Manganese (Mn)
Colour	<5 APHA	Molybdenum (Mo)
Water	<0.01%	Nickel (Ni)
Residue after evaporation	<2 ppm	Potassium (K)
Chloroform stabilized with amylene	< 0.002%	Propan-1-ol
Free Chlorine (Cl ₂)	<100 ppm	Sodium (Na)
Aluminium (Al)	<0.02 ppm	Strontium (Sr)
Barium (Ba)	<0.005 ppm	Titanium (Ti)
Cadmium (Cd)	<0.005 ppm	Total phosphorus (P)
Calcium (Ca)	<0.05 ppm	Total silicon (Si)
Chromium (Cr)	<0.005 ppm	Total sulphur (S)
Cobalt (Co)	<0.005 ppm	Vanadium (V)
Copper (Cu)	<0.01 ppm	Zinc (Zn)
Iron (Fe)	<0.05 ppm	Substances darkened
Lead (Pb)	<0.005 ppm	by H ₂ SO ₄
Lithium (Li)	<0.005 ppm	Carbon tetrachloride
Magnesium (Mg)	<0.02 ppm	1.1-Dichloroethane



Dichloromethane stabilised with amylene Distol for residue and pesticide analysis

packaging		Mfr. No
1 ℓ	Glass bottle	D/1853/15
2.5 ℓ	Glass bottle	D/1853/17
CH ₂ Cl ₂ CAS: 75-09-2 MW: 84.93 Density at 20°C: 1.325 Boiling point: 39 to 40°C Melting point: -97°C	ADR/RID: 6.1, III EINECS: 200-838-9 UN: 1593 GHS: Carc. 2 H351 P281	
Guaranteed analysis		
Water	<0.02%	GC (FID) Analysis
Residue after evaporation	<2 ppm	GC (NPD) analysis
GC (ECD) Analysis	Pass	
GC Analysis		
Analysis of pesticides: from the retention time of Lindane, after concentration of 2500 times, no impurity peak higher than that obtained with 10 pg / ml héptachloréoxi (ECD) and 5 pg / ml Parathion (NPD).		
Analysis of hydrocarbons: no impurity chain longer than the dodecane is present at a concentration greater than 0.01 ppm (FID).		
Conditions :		
Column: J&W DB 5: 30m x 0.32mm x 0.1µm df.		
Temperature profile:		
(1) 50°C for 2 min		
(2) 40°C / min until 320°C		
(3) 320°C for 30 min		
Injection: 5µl cold on the column.		

Dichloromethane stabilised with amylene for nucleotide and peptide synthesis

packaging		Mfr. No
1 ℓ	Glass bottle	D/1858/15
2.5 ℓ	Glass bottle	D/1858/17
2.5 ℓ	Coated Glass bottle	D/1858/PB17
CH ₂ Cl ₂ CAS: 75-09-2 MW: 84.93 Density at 20°C: 1.325 Boiling point: 39 to 40°C Melting point: -97°C	ADR/RID: 6.1, III EINECS: 200-838-9 UN: 1593 GHS: Carc. 2 H351 P281	
Assay (CG): >99.9%		
Guaranteed analysis		
Acidity / Alkalinity	<0.0001 meq/g	Residue after evaporation
Water	<0.01%	Free Chlorine (Cl ₂)
Absorbance specification		
Max. Absorbance (1cm cell v HPLC water)		
240 nm	<0.1 A.U.	260 nm
250 nm	<0.01 A.U.	270 nm

Dichloromethane stabilised with amylene extra dry, for synthesis

packaging		Mfr. No
1 ℓ	Coated Glass bottle	D/1854/PB15
CH ₂ Cl ₂ CAS: 75-09-2 MW: 84.93 Density at 20°C: 1.325 Boiling point: 39 to 40°C Melting point: -97°C	ADR/RID: 6.1, III EINECS: 200-838-9 UN: 1593 GHS: Carc. 2 H351 P281	
Guaranteed analysis		
Water	<0.01%	Residue after evaporation

Dichloromethane stabilised with amylene **NEW PACK.** Certified AR for analysis

packaging		Mfr. No
1 ℓ	Glass bottle	D/1852/15
1 ℓ	Coated Glass bottle	D/1852/PB15
2.5 ℓ	Glass bottle	D/1852/17
2.5 ℓ	Coated Glass bottle	D/1852/PB17
5 ℓ	Fluorinated HDPE Plastic bottle	D/1852/21
10 ℓ	Metal can	D/1852/24
25 ℓ	Metal drum	D/1852/25
200 ℓ	Metal drum	D/1852/27
CH ₂ Cl ₂ CAS: 75-09-2 MW: 84.93 Density at 20°C: 1.325 Boiling point: 39 to 40 to 40°C Melting point: -97°C	ADR/RID: 6.1, III EINECS: 200-838-9 UN: 1593 GHS: Carc. 2 H351 P281	
Assay (CG): >99.8%		
Guaranteed analysis		
Acidity / Alkalinity	<0.0003 meq/g	Iron (Fe)
Colour	<5 APHA	Lead (Pb)
Water	<0.02%	Magnesium (Mg)
Residue after evaporation	<10 ppm	Potassium (K)
Substances darkened		Sodium (Na)
by (H ₂ SO ₄)	<5 APHA	Total phosphorus (P)
Chloroform stabilized with amylene	<0.01 %	Total silicon (Si)
Calcium (Ca)	<0.5 ppm	Total sulphur (S)
Free Chlorine (Cl ₂)	<0.3 ppm	Zinc (Zn)
Copper (Cu)	<0.02 ppm	Carbon tetrachloride
		1.1-Dichloroethane
Typical analysis		
Aluminium (Al)	<0.2 ppm	Manganese (Mn)
Barium (Ba)	<0.01 ppm	Molybdenum (Mo)
Cadmium (Cd)	<0.02 ppm	Nickel (Ni)
Chromium (Cr)	<0.02 ppm	Strontium (Sr)
Cobalt (Co)	<0.02 ppm	Titanium (Ti)
Lithium (Li)	<0.01 ppm	Vanadium (V)

Lot analysis on the label.

Dichloromethane stabilised with amylene **NEW PACK.** SLR specified reagents for general laboratory work- extra pure

packaging		Mfr. No
500 mℓ	Glass bottle	D/1850/08
1 ℓ	Glass bottle	D/1850/15
1 ℓ	Metal can	D/1850/MC15
2.5 ℓ	Glass bottle	D/1850/17
2.5 ℓ	Coated Glass bottle	D/1850/PB17
5 ℓ	HDPE Plastic bottle	D/1850/21
5 ℓ	HDPE Plastic bottle	D/1850/PC21
25 ℓ	Metal drum	D/1850/25
200 ℓ	Metal drum	D/1850/27
CH ₂ Cl ₂ CAS: 75-09-2 MW: 84.93 Density at 20°C: 1.325 Boiling point: 39 to 40°C Melting point: -97°C	ADR/RID: 6.1, III EINECS: 200-838-9 UN: 1593 GHS: Carc. 2 H351 P281	
Assay (CG): >99%		
Guaranteed analysis		
Acidity / Alkalinity	<0.001 meq/g	Lead (Pb)
Colour	<10 APHA	Magnesium (Mg)
Water	<0.05%	Potassium (K)
Residue after evaporation	<500 ppm	Sodium (Na)
Substances darkened		Total phosphorus (P)
by (H ₂ SO ₄)	<25 APHA	Total silicon (Si)
Calcium (Ca)	<2 ppm	Total sulphur (S)
Copper (Cu)	<1 ppm	Zinc (Zn)
Iron (Fe)	<2 ppm	

Dichloromethane stabilised with ethanol NEW SLR
specified reagents for general laboratory work- extra pure

packaging		Mfr. No
5 ℓ	Fluorinated HDPE Plastic bottle	D/1860/PC21
CH ₂ Cl ₂	ADR/RID: 6.1, III	
CAS: 75-09-2	EINECS: 200-838-9	
MW: 84.93	UN: 1593	
Density at 20°C: 1.325	H351	
Boiling point: 39 to 40°C	P281	
Melting point: -97°C		
Assay (CG): > 99%		
Guaranteed analysis		
Water	<0.05%	Residue after evaporation <500 ppm

Diethanolamine SLR
specified reagents for general laboratory work- extra pure

packaging		Mfr. No
500 mℓ	Glass bottle	D/2050/08
2.5 ℓ	Glass bottle	D/2050/17
C ₄ H ₁₁ NO ₂	EINECS: 203-868-0	
CAS: 111-42-2	UN: 3259	
MW: 105.14	GHS: Acute Tox. 4, Skin Irrit. 2, Eye Dam. 1, STOT RE 2	
Density at 20°C: 1.096	H318, H302, H373, H315	
Boiling point: 268°C	P280, P305+P351+P338, P310, P260,	
Melting point: 27 to 30°C	P301+P312, P302+P352	
Flash point: 179°C		
ADR/RID: 8, III		

Assay (CG): >99%		
Guaranteed analysis		
Colour	<20 APHA	Potassium (K) <20 ppm
Water	<0.5%	Sodium (Na) <50 ppm
Calcium (Ca)	<20 ppm	Total phosphorus (P) <2 ppm
Copper (Cu)	<5 ppm	Total silicon (Si) <2 ppm
Iron (Fe)	<5 ppm	Total sulphur (S) <5 ppm
Lead (Pb)	<2 ppm	Zinc (Zn) <20 ppm
Magnesium (Mg)	<5 ppm	

Diethylamine Certified AR
for analysis

packaging		Mfr. No
500 mℓ	Glass bottle	D/2110/08
2.5 ℓ	Glass bottle	D/2110/17
C ₄ H ₁₁ N	EINECS: 203-716-3	
CAS: 109-89-7	UN: 1154	
MW: 73.14	GHS: Flam. Liq. 2, Acute	
Density at 20°C: 0.710	Tox. 3, Skin Corr. 1A, Eye Dam. 1	
Boiling point: 55 to 58°C	H225, H314, H332, H302, H311	
Melting point: -50°C	P280, P305+P351+P338,	
Flash point: -23°C	P210, P260, P302+P352, P310	
ADR/RID: 3, 8, II		

Assay (CG): >99.5%		
Guaranteed analysis		
Colour	<10 APHA	Potassium (K) <0.2 ppm
Residue after evaporation	<100 ppm	Sodium (Na) <1 ppm
Calcium (Ca)	<1 ppm	Total phosphorus (P) <0.2 ppm
Copper (Cu)	<0.05 ppm	Total silicon (Si) <0.05 ppm
Iron (Fe)	<0.5 ppm	Total sulphur (S) <0.5 ppm
Lead (Pb)	<0.05 ppm	Zinc (Zn) <0.5 ppm
Magnesium (Mg)	<0.5 ppm	
Typical analysis		
Aluminium (Al)	<0.1 ppm	Manganese (Mn) <0.05 ppm
Barium (Ba)	<0.05 ppm	Molybdenum (Mo) <0.05 ppm
Cadmium (Cd)	<0.05 ppm	Nickel (Ni) <0.05 ppm
Chromium (Cr)	<0.05 ppm	Strontium (Sr) <0.05 ppm
Cobalt (Co)	<0.05 ppm	Titanium (Ti) <0.05 ppm
Lithium (Li)	<0.02 ppm	Vanadium (V) <0.05 ppm

Lot analysis on the label.

Diethylamine SLR
specified reagents for general laboratory work- extra pure

packaging		Mfr. No
500 mℓ	Glass bottle	D/2100/08
2.5 ℓ	Glass bottle	D/2100/17
C ₄ H ₁₁ N	EINECS: 203-716-3	
CAS: 109-89-7	UN: 1154	
MW: 73.14	GHS: Flam. Liq. 2, Acute	
Density at 20°C: 0.710	Tox. 3, Skin Corr. 1A, Eye Dam. 1	
Boiling point: 55 to 58°C	H225, H314, H332, H302, H311	
Melting point: -50°C	P280, P305+P351+P338,	
Flash point: -23°C	P210, P260, P302+P352, P310	
ADR/RID: 3, 8, II		

Assay (CG): >99%		
Guaranteed analysis		
Colour	<20 APHA	Potassium (K) <5 ppm
Residue after evaporation	<500 ppm	Sodium (Na) <10 ppm
Calcium (Ca)	<5 ppm	Total phosphorus (P) <2 ppm
Copper (Cu)	<1 ppm	Total silicon (Si) <2 ppm
Iron (Fe)	<2 ppm	Total sulphur (S) <5 ppm
Lead (Pb)	<1 ppm	Zinc (Zn) <5 ppm
Magnesium (Mg)	<2 ppm	

Diethylammonium diethyldithiocarbamate pure

packaging		Mfr. No
100 g	Plastic jar	D/2250/48
C ₅ H ₁₁ NS ₂ ·C ₄ H ₁₁ N	GHS: Skin Irrit. 2, Eye Irrit. 2, STOT SE 3	
CAS: 1518-58-7	H319, H315, H335	
MW: 222.42	P261, P280, P302+P352,	
Melting point: 8185°C	P305+P351+P338, P304+P341	
EINECS: 216-175-3		

Guaranteed analysis		
Complies with standard IR spectrum: Pass		

Diethylene glycol Certified AR
for analysis

packaging		Mfr. No
500 mℓ	Glass bottle	D/2951/08
2.5 ℓ	Glass bottle	D/2951/17
C ₄ H ₁₀ O ₃	Flash point: 143°C	
CAS: 111-46-6	EINECS: 203-872-2	
MW: 106.12	GHS: Acute Tox. 4, STOT RE 2	
Density at 20°C: 1.11	H302, H373	
Boiling point: 245°C	P280, P301+P310	
Melting point: -10°C		

Assay (CG): >99.5%		
Guaranteed analysis		
Colour	<10 APHA	Magnesium (Mg) <1 ppm
Water	<0.2%	Potassium (K) <1 ppm
Residue after ignition	<0.01%	Sodium (Na) <5 ppm
Calcium (Ca)	<2 ppm	Total phosphorus (P) <0.5 ppm
Copper (Cu)	<0.5 ppm	Total silicon (Si) <0.1 ppm
Iron (Fe)	<1 ppm	Total sulphur (S) <1 ppm
Lead (Pb)	<0.5 ppm	Zinc (Zn) <0.5 ppm

Lot analysis on the label.

Reagents, Solvents & Volumetric solutions comply with the European Pharmacopeia

Read more on page 234

Diethylene glycol SLR
specified reagents for general laboratory work- extra pure

packaging		Mfr. No
500 mℓ	Glass bottle	D/2950/08
2.5 ℓ	Glass bottle	D/2950/17
25 ℓ	HDPE Plastic container	D/2950/25
C ₄ H ₁₀ O ₃	Flash point: 143°C	
CAS: 111-46-6	EINECS: 203-872-2	
MW: 106.12	GHS: Acute Tox. 4, STOT RE 2	
Density at 20°C: 1.11	H302, H373	
Boiling point: 245°C	P280, P301+P310	
Melting point: -10°C		

Assay (CG): >99%		
Guaranteed analysis		
Colour	<30 APHA	Potassium (K) <5 ppm
Water	<0.5%	Sodium (Na) <20 ppm
Calcium (Ca)	<10 ppm	Total phosphorus (P) <2 ppm
Copper (Cu)	<2 ppm	Total silicon (Si) <1 ppm
Iron (Fe)	<5 ppm	Total sulphur (S) <5 ppm
Lead (Pb)	<2 ppm	Zinc (Zn) <5 ppm
Magnesium (Mg)	<5 ppm	Residue after ignition <0.05%

Diethylene glycol monoethyl ether pure

packaging		Mfr. No
2.5 ℓ	Glass bottle	E/1300/17
C ₆ H ₁₄ O ₃	Flash point: 94°C	
CAS: 111-90-0	EINECS: 203-919-7	
MW: 134.18	GHS: Acute Tox. 4, Eye Irrit. 2	
Density at 20°C: 0.990	H319, H302	
Boiling point: 197°C	P301+P312, P280,	
Melting point: -80°C	P305+P351+P338	

Guaranteed analysis		
Water	<0.05%	

Diethylenetriaminepenta-acetic acid sodium salt pure, solution about 40% w/w

packaging		Mfr. No
500 mℓ	Glass bottle	D/2350/08
GHS: Skin Irrit. 2, Eye Irrit. 2, STOT SE 3		
H315, H319, H335		
P280, P302+P352, P304+P341, P305+P351+P338, P261		

Diethyl ether stabilised with ethanol NEW PACK. Certified
for HPLC

packaging		Mfr. No
1 ℓ	Glass bottle	D/2506/15
2.5 ℓ	Glass bottle	D/2506/17
2.5 ℓ	Coated Glass bottle	D/2506/PB17
C ₄ H ₁₀ O	EINECS: 200-467-2	
CAS: 60-29-7	UN: 1155	
MW: 74.12	GHS: Flam. Liq. 1, Acute	
Density at 20°C: 0.714	Tox. 4, STOT SE 3	
Boiling point: 34.6°C	H224, H336, H302,	
Melting point: -116°C	EUH019, EUH066	
Flash point: -45°C	P261, P233, P301+P312,	
ADR/RID: 3, I	P210, P240, P243, P403+P235	

Assay (CG): >99%		
Guaranteed analysis		
Acidity / Alkalinity	<0.0002 meq/g	Residue after evaporation <2 ppm
Water	<0.02%	

Absorbance specification		
Max. Absorbance (1cm cell v HPLC water)		
205 nm	<1 A.U.	250 nm <0.05 A.U.
210 nm	<0.3 A.U.	260 nm <0.025 A.U.
220 nm	<0.2 A.U.	270 nm <0.01 A.U.
230 nm	<0.125 A.U.	280 nm <0.005 A.U.
240 nm	<0.075 A.U.	Wavelength at 1 A.U. : <205 nm

Lot analysis on the label.
Filtered to 0.2 µm.

Diethyl ether Primar
for trace metal analysis, stabilised with copper gauze

packaging		Mfr. No
2.5 ℓ	Glass bottle	D/2455/17
C ₄ H ₁₀ O	EINECS: 200-467-2	
CAS: 60-29-7	UN: 1155	
MW: 74.12	GHS: Flam. Liq. 1, Acute	
Density at 20°C: 0.714	Tox. 4, STOT SE 3	
Boiling point: 34.6°C	H224, H336, H302,	
Melting point: -116°C	EUH019, EUH066	
Flash point: -45°C	P261, P233, P301+P312,	
ADR/RID: 3, I	P210, P240, P243, P403+P235	

Assay (CG): >99.9%		
Guaranteed analysis		
Acidity / Alkalinity	<0.0002 meq/g	Molybdenum (Mo) <0.005 ppm
Colour	<5 APHA	Nickel (Ni) <0.005 ppm
Water	<0.01%	Potassium (K) <0.05 ppm
Residue after evaporation	<2 ppm	Propan--1-ol <0.005 ppm
Acetone	<0.002%	Sodium (Na) <0.2 ppm
Aluminium (Al)	<0.02 ppm	Strontium (Sr) <0.005 ppm
Barium (Ba)	<0.005 ppm	Titanium (Ti) <0.005 ppm
Cadmium (Cd)	<0.005 ppm	Total phosphorus (P) <0.02 ppm
Calcium (Ca)	<0.05 ppm	Total silicon (Si) <0.02 ppm
Chromium (Cr)	<0.005 ppm	Total sulphur (S) <0.1 ppm
Cobalt (Co)	<0.005 ppm	Vanadium (V) <0.005 ppm
Copper (Cu)	<0.01 ppm	Zinc (Zn) <0.01 ppm
Iron (Fe)	<0.05 ppm	Substances darkened
Lead (Pb)	<0.005 ppm	by H ₂ SO ₄ <5 APHA
Lithium (Li)	<0.005 ppm	Ethanol 90 % v/v <0.002%
Magnesium (Mg)	<0.02 ppm	Ether peroxide <1 ppm
Manganese (Mn)	<0.005 ppm	Methanol <0.005%

Diethyl ether Distol
for residue and pesticide analysis, stabilised with copper gauze

packaging		Mfr. No
1 ℓ	Glass bottle	D/2503/15
2.5 ℓ	Glass bottle	D/2503/17
C ₄ H ₁₀ O	EINECS: 200-467-2	
CAS: 60-29-7	UN: 1155	
MW: 74.12	GHS: Flam. Liq. 1, Acute	
Density at 20°C: 0.714	Tox. 4, STOT SE 3	
Boiling point: 34.6°C	H224, H336, H302,	
Melting point: -116°C	EUH019, EUH066	
Flash point: -45°C	P261, P233, P301+P312,	
ADR/RID: 3, I	P210, P240, P243, P403+P235	

Guaranteed analysis		
Water	<0.02%	GC (FID) Analysis Pass
Residue after evaporation	<2 ppm	GC (NPD) analysis Pass
GC (ECD) Analysis	Pass	

We Support Your Chemistry For

- Semi-bulk and bulk requirements
- Tailored solvents and solvent blends or customized solutions
- Testing services or customized packaging and labeling

Please email your contact details and a brief description of your requirements to begeel.sha.quotes@thermofisher.com

Diethyl ether stabilised with ethanol **NEW PACK.**
extra dry, for synthesis

packaging		Mfr. No
1 ℓ	Glass bottle	D/2502/15
2.5 ℓ	Glass bottle	D/2502/17
C ₄ H ₁₀ O		EINECS: 200-467-2
CAS: 60-29-7		UN: 1155
MW: 74.12		GHS: Flam. Liq. 1, Acute
Density at 20°C: 0.714		Tox. 4, STOT SE 3
Boiling point: 34.6°C		H224, H336, H302,
Melting point: -116°C		EUH019, EUH066
Flash point: -45°C		P261, P233, P301+P312,
ADR/RID: 3, I		P210, P240, P243, P403+P235
Guaranteed analysis		
Water	<0.01%	Peroxide <0.00005%
Residue after evaporation	<2 ppm	

Diethyl ether stabilised with BHT **NEW PACK.** Certified AR, Eur.Ph.
for analysis

packaging		Mfr. No
500 mℓ	Coated Glass bottle	D/2450/PB08
1 ℓ	Glass bottle	D/2450/15
1 ℓ	Coated Glass bottle	D/2450/PB15
2.5 ℓ	Glass bottle	D/2450/17
2.5 ℓ	Coated Glass bottle	D/2450/PB17
5 ℓ	Mild steel drum	D/2450/MS21
25 ℓ	Metal drum	D/2450/25
200 ℓ	Metal drum	D/2450/27
C ₄ H ₁₀ O		EINECS: 200-467-2
CAS: 60-29-7		UN: 1155
MW: 74.12		GHS: Flam. Liq. 1, Acute
Density at 20°C: 0.714		Tox. 4, STOT SE 3
Boiling point: 34.6°C		H224, H336, H302,
Melting point: -116°C		EUH019, EUH066
Flash point: -45°C		P261, P233, P301+P312,
ADR/RID: 3, I		P210, P240, P243, P403+P235
Assay (CG): >99.5%		
Guaranteed analysis		
Acidity / Alkalinity	<0.0002 meq/g	Substances
Colour	<5 APHA	reducing KMnO ₄
Water	<0.03%	Total phosphorus (P)
Residue after evaporation	<10 ppm	Total silicon (Si)
Carbonyl compounds	: <0.001%	Total sulphur (S)
Calcium (Ca)	<0.2 ppm	Zinc (Zn)
Copper (Cu)	<0.05 ppm	Acetone
Iron (Fe)	<0.1 ppm	Ethanol 90 % v/v
Lead (Pb)	<0.05 ppm	Ether peroxide
Magnesium (Mg)	<0.05 ppm	Methanol
Potassium (K)	<0.5 ppm	Substances darkened
Sodium (Na)	<0.5 ppm	by H ₂ SO ₄
Typical analysis		
Aluminium (Al)	<0.05 ppm	Manganese (Mn)
Barium (Ba)	<0.01 ppm	Molybdenum (Mo)
Cadmium (Cd)	<0.01 ppm	Nickel (Ni)
Chromium (Cr)	<0.01 ppm	Strontium (Sr)
Cobalt (Co)	<0.01 ppm	Titanium (Ti)
Lithium (Li)	<0.01 ppm	Vanadium (V)

Lot analysis on the label.

Safe, environmentally friendly and effective packaging
Fisher Chemical outer packaging is UN compliant and 100% recyclable. It's an easy and safe way to reduce your waste management costs and contribute to a cleaner environment.

Diethyl ether stabilised with BHT **NEW PACK.** SLR
specified reagents for general laboratory work- extra pure

packaging		Mfr. No
500 mℓ	Coated Glass bottle	D/2400/PB08
2.5 ℓ	Glass bottle	D/2400/17
2.5 ℓ	Coated Glass bottle	D/2400/PB17
5 ℓ	Mild steel drum	D/2400/MS21
25 ℓ	Metal drum	D/2400/25
200 ℓ	Metal drum	D/2400/27
C ₄ H ₁₀ O		EINECS: 200-467-2
CAS: 60-29-7		UN: 1155
MW: 74.12		GHS: Flam. Liq. 1, Acute
Density at 20°C: 0.714		Tox. 4, STOT SE 3
Boiling point: 34.6°C		H224, H336, H302,
Melting point: -116°C		EUH019, EUH066
Flash point: -45°C		P261, P233, P301+P312,
ADR/RID: 3, I		P210, P240, P243, P403+P235

Assay (CG): >99%

Acidity / Alkalinity	<0.001 meq/g	Lead (Pb)	<1 ppm
Colour	<10 APHA	Magnesium (Mg)	<2 ppm
Water	<0.2%	Potassium (K)	<2 ppm
Residue after evaporation	<20 ppm	Sodium (Na)	<5 ppm
Substances		Total phosphorus (P)	<1 ppm
reducing KMnO ₄	<0.005%	Total silicon (Si)	<1 ppm
Calcium (Ca)	<5 ppm	Total sulphur (S)	<1 ppm
Copper (Cu)	<1 ppm	Zinc (Zn)	<1 ppm
Iron (Fe)	<2 ppm	Ether peroxide	<10 ppm

Diethyl ether stabilised with pyrogallol
pure, for synthesis

Packaging		Mfr. No
600 mℓ	Glass bottle	D/2500/08
1 ℓ	Glass bottle	D/2500/15
2.5 ℓ	Glass bottle	D/2500/17
2.5 ℓ	Coated Glass bottle	D/2500/PB17
C ₄ H ₁₀ O		EINECS: 200-467-2
CAS: 60-29-7		UN: 1155
MW: 74.12		GHS: Flam. Liq. 1, Acute
Density at 20°C: 0.714		Tox. 4, STOT SE 3
Boiling point: 34.6°C		H224, H336, H302,
Melting point: -116°C		EUH019, EUH066
Flash point: -45°C		P261, P233, P301+P312,
ADR/RID: 3, I		P210, P240, P243, P403+P235
		
Guaranteed analysis		
Water	<0.02%	Peroxide <0.00005%
Residue after evaporation	<20 ppm	

1.2-Dihydroxyanthraquinone
See Alizarin, page 27

(3.4-Dihydroxy-2-anthraquinonyl)methylimino-diacetic acid dihydrate
See 3-Aminomethylalizarin-NN-diacetic acid, page 32

1.3-Dihydroxybenzene See Resorcinol, page 269

UHPLC
Filtered to 0.1 µm, Fisher Chemical UHPLC gradient grade solvents (Acetonitrile, Methanol and Water) are perfect for those working at high pressure.
See page 309 for more information

1.2-Dihydroxybenzene-3.5-disulfonic acid, disodium salt
pure, Metal indicator (Tiron) for Iron and Titanium.

packaging		Mfr. No
25 g	Plastic jar	D/3000/46
(OH) ₂ C ₆ H ₂ (SO ₃ Na) ₂ ·H ₂ O		GHS: Skin Irrit. 2, Eye Irrit. 2, STOT SE 3
CAS: 149-45-1		H319, H315, H335
MW: 332.22		P261, P302+P352, P280, P305+P351+P338
Melting point: >300°C		
EINECS: 205-741-5		

Guaranteed analysis
Complies with standard IR spectrum: Pass

1.8-Dihydroxynaphthalene-3.6-disulfonic acid, disodium salt dihydrate
See Chromotropic acid sodium salt dihydrate formaldehyde determination, page 90

di-(2-Hydroxyphenylimino)ethane
pure, for calcium determination

packaging		Mfr. No
10 g	Plastic jar	D/3100/45
C ₁₄ H ₁₂ N ₂ O ₂		GHS: Skin Irrit. 2, Eye Irrit. 2, STOT SE 3
CAS: 1149-16-2		H335, H319, H315
MW: 240.26		P261, P302+P352, P280, P305+P351+P338
Melting point: 201 to 205°C		
EINECS: 214-560-0		

Guaranteed analysis
Complies with standard IR spectrum: Pass

L(+)-Dihydroxysuccinic acid
See (+)-Tartaric acid, page 341

Di-iodomethane SLR
specified reagents for general laboratory work- extra pure

packaging		Mfr. No
25 mℓ	Glass bottle	D/3200/04
CH ₂ I ₂		UN: 2810
CAS: 75-11-6		GHS: Acute Tox. 4, Skin Irrit. 2, Eye Irrit. 2, STOT SE 3
MW: 267.84		H302, H312, H319, H332, H315, H335
Density at 20°C: 3.32		P302+P352, P280, P261, P301+P312, P304+P340, P305+P351+P338
Boiling point: 181°C		
Melting point: 6°C		
Flash point: >110°C		
ADR/RID: 6.1, III		
EINECS: 200-841-5		

Di-iso-propyl ether stabilised with BHT **NEW PACK.** Certified AR
for analysis

packaging		Mfr. No
1 ℓ	Glass bottle	D/5002/15
2.5 ℓ	Glass bottle	D/5002/17
5 ℓ	Fluorinated HDPE Plastic bottle	D/5002/21
25 ℓ	Metal drum	D/5002/25
C ₆ H ₁₄ O		ADR/RID: 3, II
CAS: 108-20-3		EINECS: 203-560-6
MW: 102.18		UN: 1159
Density at 20°C: 0.720		GHS: Flam. Liq. 2, STOT SE 3
Boiling point: 68°C		H225, H336, EUH019, EUH066
Melting point: -85.5°C		P261, P233, P210, P240
Flash point: -29°C		

Assay (CG): >99%

Guaranteed analysis		
Acidity / Alkalinity	<0.002 meq/g	Magnesium (Mg)
Colour	<10 APHA	Potassium (K)
Water	<0.1%	Sodium (Na)
Residue after evaporation	<20 ppm	Total phosphorus (P)
Calcium (Ca)	<0.5 ppm	Total silicon (Si)
Copper (Cu)	<0.05 ppm	Total sulphur (S)
Iron (Fe)	<0.1 ppm	Zinc (Zn)
Lead (Pb)	<0.02 ppm	Peroxide

Di-iso-propyl ether stabilised with BHT **NEW PACK.** SLR
specified reagents for general laboratory work- extra pure

packaging		Mfr. No
1 ℓ	Glass bottle	D/5000/15
2.5 ℓ	Glass bottle	D/5000/17
2.5 ℓ	Coated Glass bottle	D/5000/PB17
5 ℓ	Fluorinated HDPE Plastic bottle	D/5000/21
5 ℓ	Fluorinated HDPE square Drum	D/5000/FP21

C ₆ H ₁₄ O	ADR/RID: 3, II
CAS: 108-20-3	EINECS: 203-560-6
MW: 102.18	UN: 1159
Density at 20°C: 0.720	GHS: Flam. Liq. 2, STOT SE 3
Boiling point: 68°C	H225, H336, EUH019, EUH066
Melting point: -85.5°C	P261, P233, P210, P240
Flash point: -29°C	

Assay (CG): >99%

Guaranteed analysis		
Acidity / Alkalinity	<0.002 meq/g	Magnesium (Mg)
Colour	<20 APHA	Potassium (K)
Water	<0.2%	Sodium (Na)
Residue after evaporation	<100 ppm	Total phosphorus (P)
Calcium (Ca)	<5 ppm	Total silicon (Si)
Copper (Cu)	<1 ppm	Total sulphur (S)
Iron (Fe)	<2 ppm	Zinc (Zn)
Lead (Pb)	<1 ppm	

Di-iso-propylnaphthalene SLR
extra pure, mixed isomers

packaging		Mfr. No
2.5 ℓ	Glass bottle	D/5010/17

C ₁₆ H ₂₀	Boiling point: 290 to 299°C
CAS: 38640-62-9	Melting point: -47°C
MW: 212.33	EINECS: 254-052-6

Dimedone
pure

packaging		Mfr. No
25 g	Plastic jar	D/3220/46

C ₈ H ₁₂ O ₂	Melting point: 147 to 150°C
CAS: 126-81-8	EINECS: 204-804-4
MW: 140.18	

Guaranteed analysis
Identity test : Pass

1.2-Dimethoxyethane SLR
specified reagents for general laboratory work- extra pure

packaging		Mfr. No
500 mℓ	Coated Glass bottle	D/3280/PB08
C ₄ H ₁₀ O ₂		EINECS: 203-794-9
CAS: 110-71-4		UN: 2252
MW: 90.12		GHS: Flam. Liq. 2, Repr. 1B, Acute Tox. 4
Density at 20°C: 0.867		H225, H332, H360, EUH019
Boiling point: 84 to 86°C		P201, P308+P313, P233, P304+P340, P210
Melting point: -69°C		
Flash point: -6°C		
ADR/RID: 3, II		

Assay (CG): >99%

Dimethoxymethane

pure

packaging		Mfr. No
500 mℓ	Glass bottle	D/3400/08
C ₃ H ₈ O ₂	EINECS: 203-714-2	
CAS: 109-87-5	UN: 1234	
MW: 76.09	GHS: Flam. Liq. 2, Skin Irrit.	
Density at 20°C: 0.859 to 0.862	2, Eye Irrit. 2, STOT SE 3	
Boiling point: 41 to 42°C	H225, H319, H315, H335	
Melting point: -105°C	P261, P302+P352, P280,	
Flash point: -18°C	P305+P351+P338, P210, P240	
ADR/RID: 3, II		

N,N-Dimethylacetamide

for analysis

Certified AR

packaging		Mfr. No
500 mℓ	Glass bottle	D/3411/08
2.5 ℓ	Glass bottle	D/3411/17
2.5 ℓ	Coated Glass bottle	D/3411/PB17
C ₄ H ₉ NO	ADR/RID: 6.1/III	
CAS: 127-19-5	EINECS: 204-826-4	
MW: 87.12	UN: 2810	
Density at 20°C: 0.937	GHS: Repr. 1B, Acute Tox. 4	
Boiling point: 164 to 166°C	H312, H332, H360D	
Melting point: -20°C	P280, P201, P308+P313,	
Flash point: 70°C	P304+P340, P302+P352	

Assay (CG): >99.5%

Guaranteed analysis

Acidity / Alkalinity	<0.02 meq/g	Magnesium (Mg)	<0.5 ppm
Colour	<10 APHA	Potassium (K)	<1 ppm
Water	<0.1%	Sodium (Na)	<5 ppm
Residue after evaporation	<100 ppm	Total phosphorus (P)	<2 ppm
Calcium (Ca)	<2 ppm	Total silicon (Si)	<0.05 ppm
Copper (Cu)	<1 ppm	Total sulphur (S)	<5 ppm
Iron (Fe)	<2 ppm	Zinc (Zn)	<0.5 ppm
Lead (Pb)	<1 ppm		

Lot analysis on the label.

N,N-Dimethylacetamide

specified reagents for general laboratory work- extra pure

SLR

packaging		Mfr. No
500 mℓ	Glass bottle	D/3410/08
2.5 ℓ	Glass bottle	D/3410/17
C ₄ H ₉ NO	ADR/RID: 6.1/III	
CAS: 127-19-5	EINECS: 204-826-4	
MW: 87.12	UN: 2810	
Density at 20°C: 0.937	GHS: Repr. 1B, Acute Tox. 4	
Boiling point: 164 to 166°C	H312, H332, H360D	
Melting point: -20°C	P280, P201, P308+P313,	
Flash point: 70°C	P304+P340, P302+P352	

Assay (CG): >99%

Guaranteed analysis

Acidity / Alkalinity	<0.05 meq/g	Magnesium (Mg)	<2 ppm
Colour	<20 APHA	Potassium (K)	<5 ppm
Water	<0.5%	Sodium (Na)	<20 ppm
Residue after evaporation	<500 ppm	Total phosphorus (P)	<10 ppm
Calcium (Ca)	<10 ppm	Total silicon (Si)	<1 ppm
Copper (Cu)	<5 ppm	Total sulphur (S)	<20 ppm
Iron (Fe)	<10 ppm	Zinc (Zn)	<5 ppm
Lead (Pb)	<5 ppm		

High Volume Solvent Delivery Systems

High-volume solvent delivery systems are available in 10L to 1000L, in Europe only. Please email your contact details and a brief description of your requirements to begel.sha.quotes@thermofisher.com



Dimethylamine

pure, solution 40% w/v

packaging		Mfr. No
500 mℓ	Glass bottle	D/3450/08
C ₂ H ₇ N	EINECS: 204-697-4	
CAS: 124-40-3	UN: 1160	
MW: 45.09	GHS: Flam. Liq. 2, Acute	
Density at 20°C: 0.890	Tox. 4, Skin Corr. 1B	
Boiling point: 54°C	H225, H332, H315	
Melting point: -37°C	P210, P280,	
Flash point: -18°C	P305+P351+P338, P310	
ADR/RID: 3, 8, II		

4-Dimethylaminobenzaldehyde

for analysis

Certified AR

packaging		Mfr. No
25 g	Plastic jar	D/3600/46
100 g	Plastic jar	D/3600/48
C ₉ H ₁₁ NO	EINECS: 202-819-0	
CAS: 100-10-7	GHS: Acute Tox. 4, Aquatic	
MW: 149.19	Chronic 3	
Boiling point: 176 to 177°C	H412, H302	
Melting point: 70 to 75°C	P273, P301+P312	
Flash point: 164°C		

Assay: >99%

Guaranteed analysis

Insoluble in alcohol	: <0.01%	Total chloride (Cl)	<0.005%
Calcium (Ca)	<20 ppm	Total phosphorus (P)	<20 ppm
Copper (Cu)	<5 ppm	Total silicon (Si)	<20 ppm
Iron (Fe)	<5 ppm	Total sulphur (S)	<50 ppm
Lead (Pb)	<5 ppm	Zinc (Zn)	<10 ppm
Magnesium (Mg)	<10 ppm	Insoluble Acids	<0.01%
Potassium (K)	<50 ppm	Bases of indole type	Pass
Sodium (Na)	<100 ppm		

Lot analysis on the label.

4-Dimethylaminobenzaldehyde

extra pure, Ehrlich's reagent

SLR

packaging		Mfr. No
25 g	Plastic jar	D/3550/46
100 g	Plastic jar	D/3550/48
C ₉ H ₁₁ NO	EINECS: 202-819-0	
CAS: 100-10-7	GHS: Acute Tox. 4, Aquatic	
MW: 149.19	Chronic 3	
Boiling point: 176 to 177°C	H412, H302	
Melting point: 70 to 75°C	P273, P301+P312	
Flash point: 164°C		

Assay: >98%

Guaranteed analysis

Calcium (Ca)	<100 ppm	Sodium (Na)	<500 ppm
Copper (Cu)	<20 ppm	n-Butyl formate	<0.02%
Iron (Fe)	<20 ppm	Total phosphorus (P)	<100 ppm
Lead (Pb)	<20 ppm	Total silicon (Si)	<100 ppm
Magnesium (Mg)	<50 ppm	Total sulphur (S)	<200 ppm
Potassium (K)	<200 ppm	Zinc (Zn)	<50 ppm

4-Dimethylaminobenzylidenerhodanine

technical

packaging		Mfr. No
5 g	Glass jar	D/3650/44
C ₁₂ H ₁₂ N ₂ O ₅	Melting point: 285 to 288°C	
CAS: 536-17-4	EINECS: 208-625-2	
MW: 264.37		

Guaranteed analysis

Complies with standard IR spectrum: Pass



4-[4-(Dimethylamino)phenylazo]benzene-sulfonic acid, sodium salt

See Methyl orange, page 206

2-[4-(Dimethylamino)phenylazo]benzoic acid

See Methyl red, page 207

N,N-Dimethylaniline

for analysis

Certified AR

packaging		Mfr. No
100 mℓ	Coated Glass bottle	D/3802/PB05
500 mℓ	Coated Glass bottle	D/3802/PB08
C ₈ H ₁₁ N	EINECS: 204-493-5	
CAS: 121-69-7	UN: 2253	
MW: 121.18	GHS: Carc. 2, Acute Tox. 3,	
Density at 20°C: 0.950	Aquatic Chronic 2	
Boiling point: 193 to 194°C	H411, H311, H351, H331, H301	
Melting point: 1.5 to 2.5°C	P280, P312, P302+P350,	
Flash point: 63°C	P304+P340, P301+P310, P273	
ADR/RID: 6.1, II		

Assay (CG): >99%

Guaranteed analysis

Insoluble in HCl	: <0.005%	Sodium (Na)	<5 ppm
Calcium (Ca)	<2 ppm	Total phosphorus (P)	<2 ppm
Copper (Cu)	<1 ppm	Total silicon (Si)	<0.05 ppm
Iron (Fe)	<2 ppm	Total sulphur (S)	<5 ppm
Lead (Pb)	<1 ppm	Zinc (Zn)	<0.5 ppm
Magnesium (Mg)	<0.5 ppm	Aniline	<0.1%
Potassium (K)	<1 ppm		

1.4-Dimethylbenzene

See p-Xylene, page 371

1.3-Dimethylbenzene

See m-Xylene, page 371

1.2-Dimethylbenzene

See o-Xylene, page 371

5.5-Dimethyl-1.3-cyclohexanedione

See Dimedone, page 111

Dimethyldichlorosilane

pure

packaging		Mfr. No
100 mℓ	Coated Glass bottle	D/3820/PB05
250 mℓ	Coated Glass bottle	D/3820/PB07
C ₂ H ₆ Cl ₂ Si	EINECS: 200-901-0	
CAS: 75-78-5	UN: 1162	
MW: 129.06	GHS: Flam. Liq. 2, Skin Irrit.	
Density at 20°C: 1.06	2, Eye Irrit. 2, STOT SE 3	
Boiling point: 70°C	H225, H335, H315, H319	
Melting point: -76°C	P304+P340, P261,	
Flash point: -9°C	P302+P352, P280,	
ADR/RID: 3, 8, II		P305+P351+P338, P210, P240

Dimethyl digol

specified reagents for general laboratory work- extra pure

SLR

packaging		Mfr. No
500 mℓ	Coated Glass bottle	D/3830/PB08
C ₆ H ₁₄ O ₃	ADR/RID: 3, III	
CAS: 111-96-6	EINECS: 203-924-4	
MW: 134.18	UN: 3271	
Density at 20°C: 0.937	GHS: Flam. Liq. 3, Repr. 1B	
Boiling point: 162°C	H226, H360, EUH019	
Melting point: -64°C	P201, P308+P313, P233, P210	
Flash point: 51°C		

Absorbance specification

Max. Absorbance (1cm cell v HPLC water)

Refractive index 1.407 to 1.4085

Dimethylformamide

for HPLC

NEW PACK.

Certified

packaging		Mfr. No
1 ℓ	Glass bottle	D/3846/15
2.5 ℓ	Glass bottle	D/3846/17
2.5 ℓ	Coated Glass bottle	D/3846/PB17
10 ℓ	Returnable Drum	D/3846/10RSS
25 ℓ	Returnable Drum	D/3846/25RSS
30 ℓ	Returnable Drum	D/3846/30RSS

C ₃ H ₇ NO	UN: 2265	
CAS: 68-12-2	GHS: Repr. 1B, Acute Tox. 4,	
MW: 73.09	Eye Irrit. 2	
Density at 20°C: 0.945	H319, H312, H332, H360D,	
Boiling point: 153°C	H226	
Melting point: -61°C	P301+P310, P302+P352,	
Flash point: 58°C	P304+P340, P280,	
ADR/RID: 3, III	P305+P351+P338, P210	
EINECS: 200-679-5		

Guaranteed analysis

Acidity / Alkalinity	<0.001 meq/g	Residue after evaporation	<10 ppm
Water	<0.05%		

Absorbance specification

Max. Absorbance (1cm cell v HPLC water)

270 nm	<0.5 A.U.	320 nm	<0.02 A.U.
280 nm	<0.15 A.U.	340 nm	<0.015 A.U.
290 nm (sol. aqu. 1%)	: <0.05 A.U.	350 nm	<0.005 A.U.

Lot analysis on the label.

Filtered to 0.2 μm.

Dimethylformamide

for Gel Permeation Chromatography (GPC)

Certified

packaging		Mfr. No
2.5 ℓ	Coated Glass bottle	D/3847/PB17
C ₃ H ₇ NO	UN: 2265	
CAS: 68-12-2	GHS: Repr. 1B, Acute Tox. 4,	
MW: 73.09	Eye Irrit. 2	
Density at 20°C: 0.945	H319, H312, H332, H360D,	
Boiling point: 153°C	H226	
Melting point: -61°C	P301+P310, P302+P352,	
Flash point: 58°C	P304+P340, P280,	
ADR/RID: 3, III	P305+P351+P338, P210	
EINECS: 200-679-5		

Guaranteed analysis

Colour	<10 APHA	Residue after evaporation	<10 ppm
Water	<0.05%		

Fisher Chemical 5L packaging change

Our improved HDPE 5L packaging is easier to handle, more robust and preserves product quality for longer.



Dimethylformamide
for peptide synthesis

packaging		Mfr. No
2.5 ℓ	Coated Glass bottle	D/3848/PB17
C ₃ H ₇ NO	UN: 2265	
CAS: 68-12-2	GHS: Repr. 1B, Acute Tox. 4,	
MW: 73.09	Eye Irrit. 2	
Density at 20°C: 0.945	H319, H312, H332, H360D,	
Boiling point: 153°C	H226	
Melting point: -61°C	P301+P310, P302+P352,	
Flash point: 58°C	P304+P340, P280,	
ADR/RID: 3, III	P305+P351+P338, P210	
EINECS: 200-679-5		
Assay (CG): >99.8%		
Guaranteed analysis		
Acidity / Alkalinity	<0.0003 meq/g	Dimethylamine
Water	<0.03%	: <10 ppm
packed under nitrogen		

Dimethylformamide
extra dry, for synthesis

packaging		Mfr. No
1 ℓ	Glass bottle	D/3842/15
2.5 ℓ	Glass bottle	D/3842/17
2.5 ℓ	Coated Glass bottle	D/3842/PB17
C ₃ H ₇ NO	UN: 2265	
CAS: 68-12-2	GHS: Repr. 1B, Acute Tox. 4,	
MW: 73.09	Eye Irrit. 2	
Density at 20°C: 0.945	H319, H312, H332, H360D,	
Boiling point: 153°C	H226	
Melting point: -61°C	P301+P310, P302+P352,	
Flash point: 58°C	P304+P340, P280,	
ADR/RID: 3, III	P305+P351+P338, P210	
EINECS: 200-679-5		
Guaranteed analysis		
Water	<0.02%	Residue after evaporation
		<20 ppm

Dimethylformamide
for analysis

packaging		Mfr. No
500 mℓ	Glass bottle	D/3841/08
1 ℓ	Glass bottle	D/3841/15
2.5 ℓ	Glass bottle	D/3841/17
2.5 ℓ	Coated Glass bottle	D/3841/PB17
5 ℓ	HDPE Plastic bottle	D/3841/21
10 ℓ	HDPE Plastic container	D/3841/24
25 ℓ	HDPE Plastic container	D/3841/25
200 ℓ	HDPE Plastic container	D/3841/27
C ₃ H ₇ NO	UN: 2265	
CAS: 68-12-2	GHS: Repr. 1B, Acute Tox. 4,	
MW: 73.09	Eye Irrit. 2	
Density at 20°C: 0.945	H319, H312, H332, H360D,	
Boiling point: 153°C	H226	
Melting point: -61°C	P301+P310, P302+P352,	
Flash point: 58°C	P304+P340, P280,	
ADR/RID: 3, III	P305+P351+P338, P210	
EINECS: 200-679-5		
Assay (CG): >99.5%		
Guaranteed analysis		
Acidity / Alkalinity	<0.0005 meq/g	Magnesium (Mg)
Colour	<10 APHA	Potassium (K)
Water	<0.05%	Sodium (Na)
Residue after evaporation	<20 ppm	Total phosphorus (P)
Calcium (Ca)	<1 ppm	Total silicon (Si)
Copper (Cu)	<0.5 ppm	Total sulphur (S)
Iron (Fe)	<1 ppm	Zinc (Zn)
Lead (Pb)	<0.5 ppm	

Typical analysis			
Aluminium (Al)	<0.2 ppm	Manganese (Mn)	<0.05 ppm
Barium (Ba)	<0.1 ppm	Molybdenum (Mo)	<0.05 ppm
Cadmium (Cd)	<0.05 ppm	Nickel (Ni)	<0.05 ppm
Chromium (Cr)	<0.05 ppm	Strontium (Sr)	<0.05 ppm
Cobalt (Co)	<0.05 ppm	Titanium (Ti)	<0.05 ppm
Lithium (Li)	<0.05 ppm	Vanadium (V)	<0.05 ppm

Lot analysis on the label.

Dimethylformamide

specified reagents for general laboratory work- extra pure

packaging		Mfr. No
500 mℓ	Glass bottle	D/3840/08
1 ℓ	Glass bottle	D/3840/15
2.5 ℓ	Glass bottle	D/3840/17
2.5 ℓ	HDPE Plastic bottle	D/3840/PB17
5 ℓ	HDPE Plastic bottle	D/3840/21
25 ℓ	HDPE Plastic container	D/3840/25
200 ℓ	HDPE Plastic container	D/3840/27
C ₃ H ₇ NO	UN: 2265	
CAS: 68-12-2	GHS: Repr. 1B, Acute Tox. 4,	
MW: 73.09	Eye Irrit. 2	
Density at 20°C: 0.945	H319, H312, H332, H360D,	
Boiling point: 153°C	H226	
Melting point: -61°C	P301+P310, P302+P352,	
Flash point: 58°C	P304+P340, P280,	
ADR/RID: 3, III	P305+P351+P338, P210	
EINECS: 200-679-5		

Assay (CG): >99%

Guaranteed analysis			
Acidity / Alkalinity	<0.005 meq/g	Magnesium (Mg)	<2 ppm
Colour	<20 APHA	Potassium (K)	<5 ppm
Water	<0.2%	Sodium (Na)	<50 ppm
Residue after evaporation	<50 ppm	Total phosphorus (P)	<2 ppm
Calcium (Ca)	<5 ppm	Total silicon (Si)	<1 ppm
Copper (Cu)	<2 ppm	Total sulphur (S)	<5 ppm
Iron (Fe)	<5 ppm	Zinc (Zn)	<2 ppm
Lead (Pb)	<2 ppm		

Dimethylglyoxime
for analysis

Packaging			Mfr. No
100 g	Plastic jar		D/3900/48
C ₄ H ₈ N ₂ O ₂	UN: 1325	 	
CAS: 95-45-4	GHS: Flam. Sol. 2, Acute Tox. 4		
MW: 116.12	H228, H332, H302, H312		
Melting point: 239 to 241°C	P280, P302+P352,		
ADR/RID: 4.1, III	P301+P312, P304+P340,		
EINECS: 202-420-1	P233, P210		
Assay: >98%			
Guaranteed analysis			
Insoluble in alcohol	: <0.01%	Potassium (K)	<20 ppm
Calcium (Ca)	<10 ppm	Sodium (Na)	<50 ppm
Copper (Cu)	<5 ppm	Total phosphorus (P)	<20 ppm
Iron (Fe)	<5 ppm	Total silicon (Si)	<10 ppm
Lead (Pb)	<2 ppm	Total sulphur (S)	<50 ppm
Magnesium (Mg)	<5 ppm	Zinc (Zn)	<5 ppm

Lot analysis on the label.

Standard Volumetric solution

For your convenience we offer a range of NIST traceable solutions. A full list of our products can be found on page 323.

Can't see what you need?

Please email your contact details and a brief description of your requirements to begeel.sha.quotes@thermofisher.com

Dimethylglyoxime
pure

packaging		Mfr. No
25 g	Plastic jar	D/3898/46
100 g	Plastic jar	D/3898/48
C ₄ H ₈ N ₂ O ₂	UN: 1325	 
CAS: 95-45-4	GHS: Flam. Sol. 2, Acute Tox. 4	
MW: 116.12	H228, H332, H302, H312	
Melting point: 239 to 241°C	P280, P302+P352,	
ADR/RID: 4.1, III	P301+P312, P304+P340,	
EINECS: 202-420-1	P233, P210	
Guaranteed analysis		
Complies with standard IR spectrum: Pass		

Dimethylglyoxime sodium salt
pure

packaging		Mfr. No
250 g	Plastic jar	D/3950/50
C ₄ H ₈ N ₂ Na ₂ O ₂ ·8H ₂ O	GHS: Skin Irrit. 2, Eye Irrit. 2	
CAS: 75006-64-3	H315, H319	
MW: 304.22	P302+P352, P280,	
EINECS: 262-523-2	P305+P351+P338	
Guaranteed analysis		
Complies with standard IR spectrum: Pass		

Dimethyl phthalate

specified reagents for general laboratory work- extra pure

packaging		Mfr. No	
500 mℓ	Glass bottle	D/4050/08	
2.5 ℓ	Glass bottle	D/4050/17	
C ₁₀ H ₁₀ O ₄	Boiling point: 282°C		
CAS: 131-11-3	Melting point: 2°C		
MW: 194.19	Flash point: 150°C		
Density at 20°C: 1.19	EINECS: 205-011-6		
Guaranteed analysis			
Acidity	<0.001 meq/g	Residue after ignition	<0.02%
Water	<0.15%		

Dimethyl sulfoxide
for HPLC

packaging		Mfr. No
250 mℓ	Coated Glass bottle	D/4125/PB07
500 mℓ	Coated Glass bottle	D/4125/PB08
2.5 ℓ	Coated Glass bottle	D/4125/PB17
10 ℓ	Returnable Drum	D/4125/10RSS
25 ℓ	Returnable Drum	D/4125/25RSS
30 ℓ	Returnable Drum	D/4125/30RSS
(CH ₃) ₂ SO		
CAS: 67-68-5		
MW: 78.13		
Density at 20°C: 1.1		
		Boiling point: 189°C
		Melting point: 18.4°C
		Flash point: 87°C
		EINECS: 200-664-3
Assay (CG): >99%		
Guaranteed analysis		
Acidity / Alkalinity	<0.0008 meq/g	Residue after evaporation
Water	<0.05%	<10 ppm
Absorbance specification		
Max. Absorbance (1cm cell v HPLC water)		
270 nm	<0.5 A.U.	290 nm (sol. aqu. 1%)
280 nm	<0.3 A.U.	320 nm
290 nm	<0.2 A.U.	360 nm
Lot analysis on the label.		
Filtered to 0.2 μm.		

Dimethyl sulfoxide
for Gel Permeation Chromatography (GPC)

packaging		Mfr. No
2.5 ℓ	Coated Glass bottle	D/4127/PB17
(CH ₃) ₂ SO	Boiling point: 189°C	
CAS: 67-68-5	Melting point: 18.4°C	
MW: 78.13	Flash point: 87°C	
Density at 20°C: 1.1	EINECS: 200-664-3	
Guaranteed analysis		
Colour	<10 APHA	Residue after evaporation
Water	<0.05%	<10 ppm

Dimethyl sulfoxide
for amino acid analysis

packaging		Mfr. No
2.5 ℓ	Coated Glass bottle	D/4128/PB17
(CH ₃) ₂ SO	Boiling point: 189°C	
CAS: 67-68-5	Melting point: 18.4°C	
MW: 78.13	Flash point: 87°C	
Density at 20°C: 1.1	EINECS: 200-664-3	
Assay (CG): >99.5%		
Guaranteed analysis		
Colour	<5 APHA	Transmission to 550 nm : >98%
Transmission to 400 nm	: >98%	

Dimethyl sulfoxide
for analysis

packaging		Mfr. No
500 mℓ	HDPE Plastic bottle	D/4121/PB08
1 ℓ	HDPE Plastic bottle	D/4121/PB15
2.5 ℓ	HDPE Plastic bottle	D/4121/PB17
5 ℓ	HDPE Plastic bottle	D/4121/21
10 ℓ	HDPE Plastic container	D/4121/24
25 ℓ	HDPE Plastic container	D/4121/25
(CH ₃) ₂ SO		Boiling point: 189°C
CAS: 67-68-5		Melting point: 18.4°C
MW: 78.13		Flash point: 87°C
Density at 20°C: 1.1		EINECS: 200-664-3
Assay (CG): >99.9%		
Guaranteed analysis		
Acidity / Alkalinity	<0.0002 meq/g	Lead (Pb)
Colour	<10 APHA	Magnesium (Mg)
Water	<0.05%	Potassium (K)
Residue after evaporation	<50 ppm	Sodium (Na)
Calcium (Ca)	<1 ppm	Total phosphorus (P)
Copper (Cu)	<0.05 ppm	Total silicon (Si)
Iron (Fe)	<0.5 ppm	Zinc (Zn)

Typical analysis			
Aluminium (Al)	<0.1 ppm	Manganese (Mn)	<0.02 ppm
Barium (Ba)	<0.02 ppm	Molybdenum (Mo)	<0.02 ppm
Cadmium (Cd)	<0.02 ppm	Nickel (Ni)	<0.1 ppm
Chromium (Cr)	<0.05 ppm	Strontium (Sr)	<0.02 ppm
Cobalt (Co)	<0.02 ppm	Titanium (Ti)	<0.02 ppm
Lithium (Li)	<0.02 ppm	Vanadium (V)	<0.02 ppm

Lot analysis on the label.

High Purity Acids
for Elemental Analysis

Fisher Chemical offers a complete range of acids to meet the most challenging applications. Choose from our Optima®, TraceMetal® and PrimarPlus® grades to ensure that the sensitivity you require is always available.

Read more page 24

Dimethyl sulfoxide		SLR	
specified reagents for general laboratory work- extra pure			
packaging		Mfr. No	
500 mℓ	HDPE Plastic bottle	D/4120/PB08	
2.5 ℓ	HDPE Plastic bottle	D/4120/PB17	
5 ℓ	HDPE Plastic bottle	D/4120/21	
25 ℓ	HDPE Plastic container	D/4120/25	
(CH ₃) ₂ SO	Boiling point: 189°C		
CAS: 67-68-5	Melting point: 18.4°C		
MW: 78.13	Flash point: 87°C		
Density at 20°C: 1.1	EINECS: 200-664-3		
Assay (CG): >99%			
Guaranteed analysis			
Acidity / Alkalinity	<0.0005 meq/g	Lead (Pb)	<1 ppm
Colour	<20 APHA	Magnesium (Mg)	<1 ppm
Water	<0.2%	Potassium (K)	<2 ppm
Residue after evaporation	<200 ppm	Sodium (Na)	<10 ppm
Calcium (Ca)	<5 ppm	Total phosphorus (P)	<1 ppm
Copper (Cu)	<1 ppm	Total silicon (Si)	<1 ppm
Iron (Fe)	<5 ppm	Zinc (Zn)	<1 ppm

Dimidium bromide-disulfine

extra pure, indicator stock solution

packaging	Mfr. No
100 mℓ	Glass bottle
Density at 20°C: 1.0	
Boiling point: 100°C	
Melting point: 0°C	

Herring DE, Lab.Practice 1962.11.113

Diocetyl sulfosuccinate sodium salt, Aerosol® OT-100

pure

packaging	Mfr. No
500 g	Plastic jar
D/1040/53	
C ₂₀ H ₃₇ O ₇ SNa	GHS: Skin Irrit. 2, Eye Dam. 1
CAS: 577-11-7	H315, H318
MW: 444.55	P261, P301+P312,
Melting point: 153 to 157°C	P302+P352, P280,
EINECS: 209-406-4	P305+P351+P338

1.4-Dioxane stabilised with BHT

for HPLC

packaging	Mfr. No
1 ℓ	Coated Glass bottle
2.5 ℓ	Coated Glass bottle
D/4556/PB15	
D/4556/PB17	
C ₄ H ₈ O ₂	EINECS: 204-661-8
CAS: 123-91-1	UN: 1165
MW: 88.11	GHS: Flam. Liq. 2, Carc. 2,
Density at 20°C: 1.034	Eye Irrit. 2, STOT SE 3
Boiling point: 101°C	H225, H319, H351, H335,
Melting point: 12°C	EUH019, EUH066
Flash point: 12°C	P281, P261, P233, P280,
ADR/RID: 3, II	P305+P351+P338, P210

Assay (CG): >99.5%

Guaranteed analysis	
Water	<0.025%
Residue after evaporation	<2 ppm

Absorbance specification			
Max. Absorbance (1cm cell v HPLC water)			
218 nm	<1 A.U.	250 nm	<0.075 A.U.
220 nm	<0.4 A.U.	260 nm	<0.05 A.U.
230 nm	<0.22 A.U.	270 nm	<0.03 A.U.
240 nm	<0.125 A.U.	Wavelength at 1 A.U.	: <218 nM

Lot analysis on the label.
Filtered to 0.2 µm.

1.4-Dioxane stabilised with BHT

extra dry, for synthesis

packaging	Mfr. No
1 ℓ	Coated Glass bottle
2.5 ℓ	Coated Glass bottle
D/4551/PB15	
D/4551/PB17	
C ₄ H ₈ O ₂	EINECS: 204-661-8
CAS: 123-91-1	UN: 1165
MW: 88.11	GHS: Flam. Liq. 2, Carc. 2,
Density at 20°C: 1.034	Eye Irrit. 2, STOT SE 3
Boiling point: 101°C	H225, H319, H351, H335,
Melting point: 12°C	EUH019, EUH066
Flash point: 12°C	P281, P261, P233, P280,
ADR/RID: 3, II	P305+P351+P338, P210

Guaranteed analysis

Water	<0.01%	Residue after evaporation	<2 ppm
-------------	--------	---------------------------------	--------

1.4-Dioxane stabilised with BHT

for analysis

packaging	Mfr. No
500 mℓ	Glass bottle
1 ℓ	Glass bottle
1 ℓ	Coated Glass bottle
2.5 ℓ	Coated Glass bottle
25 ℓ	Metal drum
D/4550/08	
D/4550/15	
D/4550/PB15	
D/4550/PB17	
D/4550/25	
C ₄ H ₈ O ₂	EINECS: 204-661-8
CAS: 123-91-1	UN: 1165
MW: 88.11	GHS: Flam. Liq. 2, Carc. 2,
Density at 20°C: 1.034	Eye Irrit. 2, STOT SE 3
Boiling point: 101°C	H225, H319, H351, H335,
Melting point: 12°C	EUH019, EUH066
Flash point: 12°C	P281, P261, P233, P280,
ADR/RID: 3, II	P305+P351+P338, P210

Assay (CG): >99.8%

Guaranteed analysis			
Acidity / Alkalinity	<0.002 meq/g	Magnesium (Mg)	<0.05 ppm
Colour	<10 APHA	Potassium (K)	<0.1 ppm
Water	<0.05%	Sodium (Na)	<0.5 ppm
Residue after evaporation	<20 ppm	Total phosphorus (P)	<0.1 ppm
Calcium (Ca)	<0.5 ppm	Total silicon (Si)	<0.05 ppm
Copper (Cu)	<0.02 ppm	Total sulphur (S)	<0.2 ppm
Iron (Fe)	<0.2 ppm	Zinc (Zn)	<0.2 ppm
Lead (Pb)	<0.05 ppm	Acetal	<0.005%

Typical analysis			
Aluminium (Al)	<0.1 ppm	Manganese (Mn)	<0.02 ppm
Barium (Ba)	<0.02 ppm	Molybdenum (Mo)	<0.02 ppm
Cadmium (Cd)	<0.05 ppm	Nickel (Ni)	<0.02 ppm
Chromium (Cr)	<0.02 ppm	Selenium (Se)	<0.1 ppm
Cobalt (Co)	<0.02 ppm	Strontium (Sr)	<0.02 ppm
Tin (Sn)	<0.02 ppm	Titanium (Ti)	<0.02 ppm
Lithium (Li)	<0.02 ppm	Vanadium (V)	<0.02 ppm

Slow moving peroxide content depending on the length of storage
Lot analysis on the label.

1.4-Dioxane stabilised with BHT

specified reagents for general laboratory work- extra pure

packaging	Mfr. No
500 mℓ	Glass bottle
2.5 ℓ	Coated Glass bottle
5 ℓ	Fluorinated HDPE Plastic bottle
25 ℓ	Metal drum
200 ℓ	Metal drum
D/4500/08	
D/4500/PB17	
D/4500/21	
D/4500/25	
D/4500/27	
C ₄ H ₈ O ₂	EINECS: 204-661-8
CAS: 123-91-1	UN: 1165
MW: 88.11	GHS: Flam. Liq. 2, Carc. 2,
Density at 20°C: 1.034	Eye Irrit. 2, STOT SE 3
Boiling point: 101°C	H225, H319, H351, H335,
Melting point: 12°C	EUH019, EUH066
Flash point: 12°C	P281, P261, P233, P280,
ADR/RID: 3, II	P305+P351+P338, P210



Assay (CG): >99%			
Guaranteed analysis			
Acidity / Alkalinity	<0.005 meq/g	Magnesium (Mg)	<1 ppm
Colour	<20 APHA	Potassium (K)	<1 ppm
Water	<0.2%	Sodium (Na)	<5 ppm
Residue after evaporation	<200 ppm	Total phosphorus (P)	<1 ppm
Calcium (Ca)	<2 ppm	Total silicon (Si)	<1 ppm
Copper (Cu)	<1 ppm	Total sulphur (S)	<2 ppm
Iron (Fe)	<2 ppm	Zinc (Zn)	<2 ppm
Lead (Pb)	<1 ppm		

Diphenylamine

for analysis

packaging	Mfr. No
100 g	Glass jar
250 g	Glass jar
D/4700/48	
D/4700/50	
C ₁₂ H ₁₁ N	UN: 2811
CAS: 122-39-4	GHS: Acute Tox. 3, STOT RE
MW: 169.23	2, Aquatic Acute 1, Aquatic
Boiling point: 302°C	Chronic 1
Melting point: 52 to 54°C	H331, H311, H373, H410, H301
Flash point: 152°C	P280, P302+P350,
ADR/RID: 6.1, III	P304+P340, P301+P312,
EINECS: 204-539-4	P260, P273

Assay (CG): >99%

Guaranteed analysis			
Minimum sensitivity to nitrates test. : Pass		Potassium (K)	<20 ppm
Calcium (Ca)	<10 ppm	Sodium (Na)	<50 ppm
Copper (Cu)	<2 ppm	Total phosphorus (P)	<10 ppm
Iron (Fe)	<5 ppm	Total silicon (Si)	<10 ppm
Lead (Pb)	<2 ppm	Total sulphur (S)	<20 ppm
Magnesium (Mg)	<5 ppm	Zinc (Zn)	<10 ppm
Nitrate (NO ₃)	Pass	Aniline	<0.01%

Diphenylamine

pure, redox indicator

packaging	Mfr. No
100 g	Glass jar
D/4600/48	
C ₁₂ H ₁₁ N	UN: 2811
CAS: 122-39-4	GHS: Acute Tox. 3, STOT RE
MW: 169.23	2, Aquatic Acute 1, Aquatic
Boiling point: 302°C	Chronic 1
Melting point: 52 to 54°C	H331, H311, H373, H410, H301
Flash point: 152°C	P280, P302+P350,
ADR/RID: 6.1, III	P304+P340, P301+P312,
EINECS: 204-539-4	P260, P273

1.5-Diphenylcarbazide

specified reagents for general laboratory work- extra pure

packaging	Mfr. No
25 g	Plastic jar
D/4750/46	
C ₁₃ H ₁₄ N ₄ O	2, Eye Irrit. 2, STOT SE 3
CAS: 140-22-7	H335, H315, H319, H302
MW: 242.28	P261, P301+P312,
Melting point: 173 to 176°C	P302+P352, P280,
EINECS: 205-403-7	P305+P351+P338
GHS: Acute Tox. 4, Skin Irrit.	

Guaranteed analysis
Complies with standard IR spectrum: Pass

Suitable for the determination of Chromium.

1.5-Diphenylcarbazone

specified reagents for general laboratory work- extra pure

packaging	Mfr. No
10 g	Glass jar
D/4850/45	
C ₁₃ H ₁₂ N ₄ O	GHS: Flam. Sol. 2, Acute
CAS: 538-62-5	Tox. 4, Skin Irrit. 2, Eye Irrit.
MW: 240.27	2, STOT SE 3
Melting point: 153 to 157°C	H228, H319, H335, H315, H302
ADR/RID: 4.1, III	P261, P301+P312,
EINECS: 208-698-0	P302+P352, P280,
UN: 1325	P305+P351+P338, P233

Guaranteed analysis

Complies with standard IR spectrum: Pass

Conforming to the IR spectrum of reference.
Reagent for the analysis of mercury.

2.5-Diphenyloxazole

pure, Scintillation counting

packaging	Mfr. No
100 g	Plastic jar
5 kg	Plastic keg
25 kg	Plastic keg
D/4865/48	
D/4865/65	
D/4865/70	
C ₁₅ H ₁₁ NO	Boiling point: 360°C
CAS: 92-71-7	Melting point: 70 to 72°C
MW: 221.26	EINECS: 202-181-3
Density at 20°C: 1.06	P262

2.2'-Dipyridyl

for analysis

packaging	Mfr. No
5 g	Plastic jar
D/5051/44	
C ₁₀ H ₈ N ₂	UN: 2811
CAS: 366-18-7	GHS: Acute Tox. 3, Skin Irrit.
MW: 156.19	2, Eye Irrit. 2, STOT SE 3
Boiling point: 272 to 273°C	H335, H312, H319, H315,
Melting point: 69 to 72°C	H332, H301
Flash point: 121°C	P301+P310, P280,
ADR/RID: 6.1, II	P302+P352, P261,
EINECS: 206-674-4	P304+P340, P305+P351+P338

Assay: >99.5%

Lot analysis on the label.

Dithio-oxamide

pure

packaging	Mfr. No
25 g	Plastic jar
D/5150/46	
(CS.NH ₂) ₂	2, Eye Irrit. 2, STOT SE 3
CAS: 79-40-3	H302, H315, H335, H319
MW: 120.20	P261, P301+P312,
Melting point: 170°C	P302+P352, P280,
EINECS: 201-203-9	P305+P351+P338
GHS: Acute Tox. 4, Skin Irrit.	

Guaranteed analysis

Complies with standard IR spectrum: Pass

Optima LC/MS

Fisher Chemical OPTIMA® LC/MS solvents, blends, additives and modifiers have set the standard of excellence for consistent, reproducible performance in the mobile phase of LC/MS.

Find out more on page 308

Dithiothreitol
for peptide sequencing. Cleland's reagent.

packaging	Mfr. No
1 g Vial	D/P351/43
5 g Plastic jar	D/P351/44

C₄H₁₀O₃S₂
CAS: 3483-12-3
MW: 154.24
Boiling point: 125 to 130°C
Melting point: 38 to 43°C
Flash point: >110°C
EINECS: 222-468-7

GHS: Acute Tox. 4, Skin Irrit. 2, Eye Irrit. 2, STOT SE 3
H302, H335, H319, H315
P261, P301+P312, P302+P352, P280, P305+P351+P338

Absorbance specification
Max. Absorbance (1cm cell v HPLC water)
283 nm (sol. aqu. 0.02 M) ... : <0.04 A.U.

Use to disrupt disulfide bridge in order to maximise protein denaturation and facilitate deterging binding.

DMF See Dimethylformamide, page 113

DMSO See Dimethyl sulfoxide, page 116

n-Dodecane
specified reagents for general laboratory work- extra pure

SLR

packaging	Mfr. No
100 mℓ Glass bottle	D/5290/05
2.5 ℓ Glass bottle	D/5290/17

C₁₂H₂₆
CAS: 112-40-3
MW: 170.34
Density at 20°C: 0.753
Boiling point: 215 to 217°C
Melting point: -12°C
Flash point: 71°C

ADR/RID: 6.1, III
EINECS: 203-967-9
UN: 2810
GHS: Acute Tox. 4, Asp. Tox. 1
H304, H332, EUH066
P310, P304+P340, P301+P310, P331

Assay (CG): >95%

1-Dodecanesulfonic acid sodium salt
for lon pair chromatography for electrochemical detection.

Certified

packaging	Mfr. No
5 g Plastic jar	D/5297/44

C₁₂H₂₅NaO₃S
CAS: 2386-53-0
MW: 272.38

Melting point: >300°C
EINECS: 219-200-6

Assay: >98%

Guaranteed analysis
Electrochemical test Pass

Absorbance specification
Max. Absorbance (1cm cell v HPLC water)
254 nm (sol. aqu. 0.005 M) : <0.05 A.U.

Lot analysis on the label.

EC Specification:
tested by HPLC-EC to ensure low electrochemical response.
Detector - ESA Coulochem II
Mode-Oxidative
Electrode graphite

UHPLC

Filtered to 0.1 μm, Fisher Chemical UHPLC gradient grade solvents (Acetonitrile, Methanol and Water) are perfect for those working at high pressure.
See page 309 for more information

1-Dodecanesulfonic acid sodium salt
for lon pair chromatography

Certified

packaging	Mfr. No
5 g Plastic jar	D/5296/44

C₁₂H₂₅NaO₃S
CAS: 2386-53-0
MW: 272.38

Melting point: >300°C
EINECS: 219-200-6

Assay: >98%

Absorbance specification
Max. Absorbance (1cm cell v HPLC water)
254 nm (sol. aqu. 0.005 M) : <0.05 A.U.

Lot analysis on the label.

Dodecanoic acid See Lauric acid, page 177

Dodecanoyl peroxide See di-Lauroyl peroxide, page 177

Dodecyl Sodium Sulfate

See Sodium dodecyl sulfate, page 290

Dodecyl trimethyl ammonium bromide
for lon pair chromatography

Certified

packaging	Mfr. No
25 g Plastic jar	D/5311/46

C₁₅H₃₄BrN
CAS: 1119-94-4
MW: 308.35
Melting point: 246°C
EINECS: 214-290-3
GHS: Acute Tox. 4, Skin Irrit.

2, Eye Irrit. 2, STOT SE 3
H335, H315, H302, H319
P261, P301+P312, P302+P352, P280, P305+P351+P338

Assay: >98%

Absorbance specification
Max. Absorbance (1cm cell v HPLC water)
254 nm <0.05 A.U.

Lot analysis on the label.

DPX mounting medium
For microscopy

packaging	Mfr. No
100 mℓ Glass bottle	D/5319/05

Flash point: 24°C
ADR/RID: 3, III
UN: 1307
GHS: Flam. Liq. 3, Repr. 1B, Acute Tox. 4, Skin Irrit. 2,

Aquatic Acute 1
H226, H312, H332, H315, H400, H360Df
P280, P302+P350, P304+P340, P210, P201, P273

Guaranteed analysis
Identity test : Pass

Dysprosium standard metal solution 1000 ppm
for ICP solution in HNO₃ 2%, NIST certified

packaging	Mfr. No
100 mℓ Plastic bottle	SPPLDY22Y
500 mℓ Plastic bottle	SPPLDY22X

H315, H319

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurites.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

EDTA See Ethylenediaminetetraacetic acid, page 126

EDTA, disodium salt
See Ethylenediaminetetraacetic acid disodium salt solution, page 127

EDTA, Tris-Borate See Tris-Borate-EDTA buffer, page 361

EDTA, trisodium salt
See Ethylenediaminetetraacetic acid disodium salt solution, page 127

EGTA See 1, page 103

Ehrlich's Reagent
See 4-Dimethylaminobenzaldehyde, page 112

1,2-Epoxypropane
specified reagents for general laboratory work- extra pure

SLR

packaging	Mfr. No
500 mℓ Coated Glass bottle	E/0050/PB08
2.5 ℓ Coated Glass bottle	E/0050/PB17

C₃H₆O
CAS: 75-56-9
MW: 58.08
Density at 20°C: 0.830
Boiling point: 34°C
Melting point: -112°C
Flash point: -37°C
ADR/RID: 3, I
EINECS: 200-879-2
UN: 1280

GHS: Flam. Liq. 1, Carc. 1B, Muta. 1B, Acute Tox. 4, Skin Irrit. 2, Eye Irrit. 2, STOT SE 3
H224, H302, H319, H312, H332, H340, H350, H315, H335
P280, P201, P308+P313, P302+P352, P261, P301+P312, P304+P340, P305+P351+P338, P210, P240, P243, P403+P235

Assay (CG): >99%

Reagents, Solvents & Volumetric solutions comply with the European Pharmacopeia

Read more on page 234

Erbium standard metal solution 1000 ppm
for ICP solution in HNO₃ 2%, NIST certified

packaging	Mfr. No
100 mℓ Plastic bottle	SPPLER22Y
500 mℓ Plastic bottle	SPPLER22X

H315, H319

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurites.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Erbium standard metal solution 10 000 ppm
for ICP solution in HNO₃ 5%, NIST certified

packaging	Mfr. No
100 mℓ Plastic bottle	SPPLER23Y
500 mℓ Plastic bottle	SPPLER23X

H314

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurites.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Eriochrome black T C.I. 14645
pure, Metal indicator

packaging	Mfr. No
25 g Glass jar	E/0100/46
100 g Glass jar	E/0100/48

C₂₀H₁₂N₃NaO₇S
CAS: 1787-61-7
MW: 461.38
EINECS: 217-250-3

GHS: Eye Irrit. 2
H319
P280, P305+P351+P338

Guaranteed analysis
Complies with standard IR spectrum: Pass

Eriochrome blue-black R C.I. 15705
pure, Metal indicator

packaging	Mfr. No
25 g Glass jar	E/0200/46

CAS: 2538-85-4
EINECS: 219-810-2
GHS: Skin Irrit. 2, Eye Irrit. 2, STOT SE 3

H315, H335, H319
P261, P302+P352, P280, P305+P351+P338

Guaranteed analysis
Complies with standard IR spectrum: Pass

Eriochrome cyanine R C.I. 43820
pure, Metal indicator

packaging	Mfr. No
25 g Glass jar	E/0250/46

C₂₃H₁₃Na₃O₉S
CAS: 3564-18-9
MW: 536.39
EINECS: 222-641-7

Guaranteed analysis
Complies with standard IR spectrum: Pass

118

Application For analysis Specified Titration

119

Eschka’s mixture for analysis			Certified AR
packaging		Mfr. No	
500 g	Plastic jar	E/0300/53	
Guaranteed analysis			
Insoluble in HCl	: <0.005%	Potassium (K)	<50 ppm
Calcium (Ca)	<500 ppm	Total phosphorus (P)	<20 ppm
Copper (Cu)	<10 ppm	Total silicon (Si)	<50 ppm
Iron (Fe)	<100 ppm	Total sulphur (S)	<50 ppm
Lead (Pb)	<10 ppm	Zinc (Zn)	<50 ppm
Typical analysis			
Magnesium oxide	65-67%	Sodium carbonate	33-35%

Ethanal See Acetaldehyde, page 13

Ethanediol for analysis			Certified AR
packaging		Mfr. No	
500 mℓ	Glass bottle	E/0352/08	
2.5 ℓ	Glass bottle	E/0352/17	
2.5 ℓ	Coated Glass bottle	E/0352/PB17	
25 ℓ	HDPE Plastic container	E/0352/25	
200 ℓ	Metal drum	E/0352/27	
C ₂ H ₆ O ₂ CAS: 107-21-1 MW: 62.07 Density at 20°C: 1.113 Boiling point: 196 to 198°C Melting point: -13°C		Flash point: 116°C EINECS: 203-473-3 GHS: Acute Tox. 4 H302 P280, P301+P310	
Assay (CG): >99.5%			
Guaranteed analysis			
Acidity / Alkalinity	<0.0002 meq/g	Lead (Pb)	<1 ppm
Colour	<10 APHA	Magnesium (Mg)	<0.5 ppm
Water	<0.1%	Potassium (K)	<0.5 ppm
Residue after ignition	<0.003%	Sodium (Na)	<2 ppm
Calcium (Ca)	<1 ppm	Total phosphorus (P)	<1 ppm
Free Chlorine (Cl ₂)	<0.2 ppm	Total silicon (Si)	<0.2 ppm
Copper (Cu)	<1 ppm	Total sulphur (S)	<10 ppm
Iron (Fe)	<2 ppm	Zinc (Zn)	<0.5 ppm
Typical analysis			
Aluminium (Al)	<0.2 ppm	Manganese (Mn)	<0.05 ppm
Barium (Ba)	<0.1 ppm	Molybdenum (Mo)	<0.05 ppm
Cadmium (Cd)	<0.05 ppm	Nickel (Ni)	<0.05 ppm
Chromium (Cr)	<0.05 ppm	Strontium (Sr)	<0.05 ppm
Cobalt (Co)	<0.05 ppm	Titanium (Ti)	<0.05 ppm
Lithium (Li)	<0.02 ppm	Vanadium (V)	<0.05 ppm

For the glycerol determination.
Lot analysis on the label.

Ethanedionic acid See Oxalic acid dihydrated, page 222

Ethanol absolute for HPLC			Certified
packaging		Mfr. No	
1 ℓ	Glass bottle	E/0665/15	
1 ℓ	Glass bottle	E/0665DF/15	
2.5 ℓ	Glass bottle	E/0665/17	
2.5 ℓ	Glass bottle	E/0665DF/17	
10 ℓ	Returnable Drum	E/0665DF/10RS	
25 ℓ	Returnable Drum	E/0665DF/25RS	
30 ℓ	Returnable Drum	E/0665DF/30RS	
C ₂ H ₅ OH CAS: 64-17-5 MW: 46.06 Density at 20°C: 0.790 Boiling point: 78°C Melting point: -114°C Flash point: 13 to 17°C		ADR/RID: 3, II EINECS: 200-578-6 UN: 1170 GHS: Flam. Liq. 2 H225 P210, P240	

Guaranteed analysis			
Assay (GC)	: >99.8%	Water	<0.3%
Acidity / Alkalinity	<0.0002 meq/g	Residue after evaporation	<10 ppm
Absorbance specification			
Max. Absorbance (1cm cell v HPLC water)			
210 nm	<0.6 A.U.	250 nm	<0.05 A.U.
220 nm	<0.3 A.U.	270 nm	<0.025 A.U.
230 nm	<0.125 A.U.	290 nm	<0.015 A.U.
The trade of absolute ethanol is subject to local regulation and additional excise duties apply (only applicable in the UK market) to the following code: E/0665/15, E/0665/17.			
Lot analysis on the label.			
Filtered to 0.2µm.			
Trade in alcohol is regulated, for more information see Sales conditions page.			

Ethanol absolute Trace metal analysis			Primar
packaging		Mfr. No	
2.5 ℓ	Glass bottle	E/0655/17	
2.5 ℓ	Glass bottle	E/0655DF/17	
C ₂ H ₅ OH CAS: 64-17-5 MW: 46.06 Density at 20°C: 0.790 Boiling point: 78°C Melting point: -114°C Flash point: 13 to 17°C		ADR/RID: 3, II EINECS: 200-578-6 UN: 1170 GHS: Flam. Liq. 2 H225 P210, P240	
Assay (CG): >99.8%			
Guaranteed analysis			
Acidity / Alkalinity	<0.0001 meq/g	Molybdenum (Mo)	<0.01 ppm
Colour	<10 APHA	Nickel (Ni)	<0.01 ppm
Residue after evaporation	<10 ppm	Potassium (K)	<0.5 ppm
Other organics	: <0.02%	Sodium (Na)	<2 ppm
Substances		Strontium (Sr)	<0.01 ppm
reducing KMnO ₄	<0.0005%	Titanium (Ti)	<0.01 ppm
Aluminium (Al)	<0.1 ppm	Total phosphorus (P)	<0.5 ppm
Barium (Ba)	<0.01 ppm	Total silicon (Si)	<0.05 ppm
Cadmium (Cd)	<0.01 ppm	Total sulphur (S)	<1 ppm
Calcium (Ca)	<0.5 ppm	Vanadium (V)	<0.01 ppm
Chromium (Cr)	<0.01 ppm	Zinc (Zn)	<0.2 ppm
Cobalt (Co)	<0.01 ppm	Acetone	<0.001%
Copper (Cu)	<0.05 ppm	Furfuraldehyde	<0.001%
Iron (Fe)	<0.2 ppm	Substances darkened	
Lead (Pb)	<0.05 ppm	by H ₂ SO ₄	<10 APHA
Lithium (Li)	<0.01 ppm	Propan--2-ol	<0.003%
Magnesium (Mg)	<0.05 ppm	Methanol	<0.002%
Manganese (Mn)	<0.01 ppm		

The trade of absolute ethanol is subject to European excise legislation and the distributor will need a European registered excise authorisation number to purchase the following codes: E/0655DF/17.
The following codes are duty-paid items that are available in the UK market only: E/0655/17.

Lot analysis on the label.

Ethanol absolute for the determination of 17-ketosteroids		
packaging		Mfr. No
500 mℓ	Glass bottle	E/0660/08
C ₂ H ₅ OH CAS: 64-17-5 MW: 46.06 Density at 20°C: 0.790 Boiling point: 78°C Melting point: -114°C Flash point: 13 to 17°C		ADR/RID: 3, II EINECS: 200-578-6 UN: 1170 GHS: Flam. Liq. 2 H225 P210, P240

This is a duty-paid item that is available in the UK market only

Ethanol absolute for analysis			Certified AR, Eur.Ph., BP, USP
packaging		Mfr. No	
100 mℓ	Glass bottle	E/0650/05	
500 mℓ	Glass bottle	E/0650/08	
500 mℓ	Glass bottle	E/0650DF/08	
1 ℓ	Glass bottle	E/0650/15	
1 ℓ	Glass bottle	E/0650DF/15	
2.5 ℓ	Glass bottle	E/0650/17	
2.5 ℓ	HDPE Plastic bottle	E/0650/PC17	
2.5 ℓ	Glass bottle	E/0650DF/17	
2.5 ℓ	HDPE Plastic bottle	E/0650DF/C17	
2.5 ℓ	Coated Glass bottle	E/0650DF/P17	
5 ℓ	Fluorinated HDPE Plastic bottle	E/0650DF/21	
10 ℓ	HDPE Plastic container	E/0650DF/24	
25 ℓ	Metal drum	E/0650/25	
25 ℓ	Metal drum	E/0650DF/25	
C ₂ H ₅ OH CAS: 64-17-5 MW: 46.06 Density at 20°C: 0.790 Boiling point: 78°C Melting point: -114°C Flash point: 13 to 17°C		ADR/RID: 3, II EINECS: 200-578-6 UN: 1170 GHS: Flam. Liq. 2 H225 P210, P240	

Assay (GC)	: >99.8%	Iron (Fe)	<0.2 ppm
Acidity (acetic acid)	<0.002%	Lead (Pb)	<0.05 ppm
Acidity / Alkalinity	<0.0001 meq/g	Magnesium (Mg)	<0.05 ppm
Colour	<10 APHA	Potassium (K)	<0.5 ppm
Water	<0.2%	Sodium (Na)	<2 ppm
Residue after evaporation	<10 ppm	Total phosphorus (P)	<0.5 ppm
Substances darkened		Total silicon (Si)	<0.05 ppm
by (H ₂ SO ₄)	<10 APHA	Total sulphur (S)	<1 ppm
Substances		Zinc (Zn)	<0.2 ppm
reducing KMnO ₄	<0.0005%	Aldehydes	<0.002%
Acetone	<0.001%	Furfuraldehyde	<0.001%
Other organics	: <0.02%	Methanol	<0.002%
Calcium (Ca)	<0.5 ppm	Propan--2-ol	<0.003%
Copper (Cu)	<0.05 ppm		

Typical analysis			
Aluminium (Al)	<0.1 ppm	Manganese (Mn)	<0.01 ppm
Barium (Ba)	<0.01 ppm	Molybdenum (Mo)	<0.01 ppm
Cadmium (Cd)	<0.01 ppm	Nickel (Ni)	<0.01 ppm
Chromium (Cr)	<0.01 ppm	Strontium (Sr)	<0.01 ppm
Cobalt (Co)	<0.01 ppm	Titanium (Ti)	<0.01 ppm
Lithium (Li)	<0.01 ppm	Vanadium (V)	<0.01 ppm

Trade in alcohol is regulated, for more information see Sales conditions page.
The trade of absolute ethanol is subjected to the European excise legislation and the distributor will need an European registered excise authorisation number to the following codes: E/0650DF/08, E/0650DF/15, E/0650DF/17, E/0650DF/C17, E/0650DF/P17, E/0650DF/21, E/0650DF/24, E/0650DF/25.

The trade of absolute ethanol is subject to local regulation and additional excise duties apply (only applicable in the UK market) to the following codes: E/0650/05, E/0650/08, E/0650/15, E/0650/17, E/0650/PC17, E/0650/25.

Lot analysis on the label.

Ethanol absolute specified reagents for general laboratory work- extra pure			SLR
packaging		Mfr. No	
1 ℓ	Glass bottle	E/0600DF/15	
2.5 ℓ	Glass bottle	E/0600/17	
2.5 ℓ	HDPE Plastic bottle	E/0600/PC17	
2.5 ℓ	Glass bottle	E/0600DF/17	
2.5 ℓ	HDPE Plastic bottle	E/0600DF/C17	
5 ℓ	Fluorinated HDPE Plastic bottle	E/0600DF/21	
5 ℓ	Fluorinated HDPE square Drum	E/0600DF/F21	
25 ℓ	Metal drum	E/0600DF/25	
25 ℓ	HDPE Plastic container	E/0600DF/P25	
C ₂ H ₅ OH CAS: 64-17-5 MW: 46.06 Density at 20°C: 0.790 Boiling point: 78°C Melting point: -114°C Flash point: 13 to 17°C		ADR/RID: 3, II EINECS: 200-578-6 UN: 1170 GHS: Flam. Liq. 2 H225 P210, P240	

Guaranteed analysis			
Assay (GC)	: >99%	reducing KMnO ₄	<0.001%
Acidity / Alkalinity	<0.005 meq/g	Zinc (Zn)	<2 ppm
Colour	<20 APHA	Copper (Cu)	<1 ppm
Water	<0.5%	Iron (Fe)	<2 ppm
Residue after evaporation	<50 ppm	Lead (Pb)	<1 ppm
Substances		Magnesium (Mg)	<2 ppm
reducing KMnO ₄	<0.001%	Total phosphorus (P)	<2 ppm
Calcium (Ca)	<5 ppm	Total silicon (Si)	<1 ppm
Potassium (K)	<2 ppm	Total sulphur (S)	<5 ppm
Sodium (Na)	<10 ppm	Substances darkened	
Substances		by H ₂ SO ₄	<20 APHA

The trade of absolute ethanol is subjected to the European excise legislation and the distributor will need an European registered excise authorisation number to the following codes: E/0600DF/15, E/0600DF/17, E/0600DF/C17, E/0600DF/21, E/0600DF/F21, E/0600DF/P25.
The trade of absolute ethanol is subject to local regulation and additional excise duties apply (only applicable in the UK market) to the following codes: E/0600/05, E/0600/08, E/0600/17, E/0600/PC17.

Ethanol 95% v/v for analysis			Certified AR, USP
packaging		Mfr. No	
500 mℓ	Glass bottle	E/0550/08	
1 ℓ	Glass bottle	E/0550DF/15	
2.5 ℓ	Glass bottle	E/0550/17	
2.5 ℓ	HDPE Plastic bottle	E/0550/PC17	
2.5 ℓ	Glass bottle	E/0550DF/17	
5 ℓ	Fluorinated HDPE Plastic bottle	E/0550DF/21	
25 ℓ	Metal drum	E/0550DF/25	
C ₂ H ₅ OH CAS: 64-17-5 MW: 46.06 Density at 20°C: 0.790 Boiling point: 78°C Melting point: -114°C Flash point: 13 to 17°C		ADR/RID: 3, II EINECS: 200-578-6 UN: 1170 GHS: Flam. Liq. 2 H225 P210, P240	
Guaranteed analysis			
Acidity / Alkalinity	<0.0005 meq/g	Sodium (Na)	<2 ppm
Colour	<10 APHA	Substances	
Density at 15.5°C	: 0.812 to 0.816	reducing KMnO ₄	<0.0005%
Water insoluble substances test	: Pass	Total phosphorus (P)	<0.2 ppm
Residue after evaporation	<10 ppm	Total silicon (Si)	<0.5 ppm
Acetone	<0.002%	Total sulphur (S)	<2 ppm
Amyl alcohol and non volatile carbonisable substances test	: Pass	Zinc (Zn)	<0.2 ppm
Calcium (Ca)	<1 ppm	Furfuraldehyde	<0.001%
Copper (Cu)	<0.05 ppm	Methanol	<0.02%
Iron (Fe)	<0.2 ppm	Others impurities	<0.2%
Lead (Pb)	<0.05 ppm	Propan--2-ol	<0.02%
Magnesium (Mg)	<0.2 ppm	Substances darkened	
Potassium (K)	<0.1 ppm	by H ₂ SO ₄	<10 APHA

Typical analysis			
Aluminium (Al)	<0.1 ppm	Manganese (Mn)	<0.02 ppm
Barium (Ba)	<0.02 ppm	Molybdenum (Mo)	<0.02 ppm
Cadmium (Cd)	<0.05 ppm	Nickel (Ni)	<0.02 ppm
Chromium (Cr)	<0.02 ppm	Strontium (Sr)	<0.02 ppm
Cobalt (Co)	<0.02 ppm	Titanium (Ti)	<0.02 ppm
Lithium (Li)	<0.02 ppm	Vanadium (V)	<0.02 ppm

The trade of absolute ethanol is subjected to the European excise legislation and the distributor will need an European registered excise authorisation number to the following codes: E/0550DF/15, E/0550DF/17, E/0550DF/21, E/0550DF/25 .
The trade of absolute ethanol is subject to local regulation and additional excise duties apply (only applicable in the UK market) to the following codes: E/0550/08, E/0550/17, E/0550/PC17 .
Lot analysis on the label

High Volume Solvent Delivery Systems

High-volume solvent delivery systems are available in 10L to 1000L, in Europe only. Please email your contact details and a brief description of your requirements to begei.shi.quotes@thermofisher.com



Ethanol 96% v/v **NEW PACK.** Certified AR, Eur.Ph., BP, USP for analysis

packaging		Mfr. No
2.5 ℓ	Glass bottle	E/0555/17
2.5 ℓ	HDPE Plastic bottle	E/0555/PC17
2.5 ℓ	Glass bottle	E/0555DF/17
5 ℓ	Fluorinated HDPE Plastic bottle	E/0555DF/21
25 ℓ	Metal drum	E/0555/25
25 ℓ	Metal drum	E/0555DF/25
C ₂ H ₅ OH CAS: 64-17-5 MW: 46.06 Density at 20°C: 0.790 Boiling point: 78°C Melting point: -114°C Flash point: 13 to 17°C		ADR/RID: 3, II EINECS: 200-578-6 UN: 1170 GHS: Flam. Liq. 2 H225 P210, P240
Guaranteed analysis		
Acidity / Alkalinity	<0.0005 meq/g	Potassium (K) <0.1 ppm
Colour	<10 APHA	Sodium (Na) <2 ppm
Water insoluble substances test	: Pass	Substances
Residue after evaporation	<10 ppm	reducing KMnO ₄ <0.0005%
Substances darkened		Total phosphorus (P) <0.2 ppm
by (H ₂ SO ₄)	<10 APHA	Total silicon (Si) <0.05 ppm
Acetone	<0.002%	Zinc (Zn) <0.2 ppm
Amyl alcohol and non volatile carbonisable substances test	: Pass	Furfuraldehyde <0.001%
Calcium (Ca)	<1 ppm	Methanol <0.02%
Copper (Cu)	<0.05 ppm	Others impurities <0.02%
Iron (Fe)	<0.2 ppm	Propan--2-ol <0.02%
Lead (Pb)	<0.05 ppm	Substances darkened
Magnesium (Mg)	<0.2 ppm	by H ₂ SO ₄ <10 APHA
		Total sulphur (S) <2 ppm

Typical analysis	
Aluminium (Al)	<0.1 ppm
Barium (Ba)	<0.02 ppm
Cadmium (Cd)	<0.05 ppm
Chromium (Cr)	<0.02 ppm
Cobalt (Co)	<0.02 ppm
Lithium (Li)	<0.02 ppm
Manganese (Mn)	<0.02 ppm
Molybdenum (Mo)	<0.02 ppm
Nickel (Ni)	<0.02 ppm
Strontium (Sr)	<0.02 ppm
Titanium (Ti)	<0.02 ppm
Vanadium (V)	<0.02 ppm

The trade of absolute ethanol is subject to local regulation and additional excise duties apply (only applicable in the UK market) to the following codes: E/0555/17, E/0555/PC17, E/0555/25.

Trade in alcohol is regulated, for more information see Sales conditions page. The trade of absolute ethanol is subjected to the European excise legislation and the distributor will need an European registered excise authorisation number to the following codes: E/0555DF/17, E/0555DF/21, E/0555DF/25. Lot analysis on the label

Ethanol 96 % v/v **SLR** specified reagents for general laboratory work- extra pure

packaging		Mfr. No
2.5 ℓ	Glass bottle	E/0551/17
2.5 ℓ	HDPE Plastic bottle	E/0551/PC17
2.5 ℓ	Glass bottle	E/0551DF/17
2.5 ℓ	HDPE Plastic bottle	E/0551DF/C17
C ₂ H ₅ OH CAS: 64-17-5 MW: 46.06 Density at 20°C: 0.790 Boiling point: 78°C Melting point: -114°C Flash point: 13 to 17°C		ADR/RID: 3, II EINECS: 200-578-6 UN: 1170 GHS: Flam. Liq. 2 H225 P210, P240
Guaranteed analysis		
Acidity / Alkalinity	<0.005 meq/g	Lead (Pb) <1 ppm
Colour	<20 APHA	Magnesium (Mg) <2 ppm
Residue after evaporation	<100 ppm	Potassium (K) <2 ppm
Substances		Sodium (Na) <10 ppm
reducing KMnO ₄	<0.001%	Total phosphorus (P) <2 ppm
Calcium (Ca)	<5 ppm	Total silicon (Si) <1 ppm
Copper (Cu)	<1 ppm	Total sulphur (S) <5 ppm
Iron (Fe)	<2 ppm	Zinc (Zn) <2 ppm

The trade of absolute ethanol is subject to European excise legislation and the distributor will need a European registered excise authorisation number to purchase the following codes: E/0551DF/17, E/0551DF/C17 .

The following codes are duty-paid items that are available in the UK market only: E/0551/17, E/0551/PC17 .

Ethanol 95% v/v **SLR** specified reagents for general laboratory work- extra pure

packaging		Mfr. No
500 mℓ	Glass bottle	E/0500/08
500 mℓ	Glass bottle	E/0500DF/08
1 ℓ	Glass bottle	E/0500DF/15
2.5 ℓ	Glass bottle	E/0500/17
2.5 ℓ	HDPE Plastic bottle	E/0500/PC17
2.5 ℓ	Glass bottle	E/0500DF/17
2.5 ℓ	HDPE Plastic bottle	E/0500DF/C17
2.5 ℓ	Coated Glass bottle	E/0500DF/P17
5 ℓ	HDPE Plastic bottle	E/0500DF/P21
10 ℓ	HDPE Plastic container	E/0500DF/P24
25 ℓ	Plastic container	E/0500DF/P25
C ₂ H ₅ OH CAS: 64-17-5 MW: 46.06 Density at 20°C: 0.790 Boiling point: 78°C Melting point: -114°C Flash point: 13 to 17°C		ADR/RID: 3, II EINECS: 200-578-6 UN: 1170 GHS: Flam. Liq. 2 H225 P210, P240
Guaranteed analysis		
Acidity / Alkalinity	<0.005 meq/g	Zinc (Zn) <2 ppm
Colour	<20 APHA	Copper (Cu) <1 ppm
Residue after evaporation	<100 ppm	Iron (Fe) <2 ppm
Substances		Lead (Pb) <1 ppm
reducing KMnO ₄	<0.001%	Magnesium (Mg) <2 ppm
Calcium (Ca)	<5 ppm	Total phosphorus (P) <2 ppm
Potassium (K)	<2 ppm	Total silicon (Si) <1 ppm
Sodium (Na)	<10 ppm	Total sulphur (S) <5 ppm

The trade of absolute ethanol is subjected to the European excise legislation and the distributor will need an European registered excise authorisation number to the following codes:E/0500DF/08, E/0500DF/15, E/0500DF/C17, E/0500DF/P17, E/0500DF/P21, E/0500DF/P24, E/0500DF/P25 .

The trade of absolute ethanol is subject to local regulation and additional excise duties apply (only applicable in the UK market) to the following codes: E/0500/08,E/0500/17, E/0500/PC17 .

Ethanol 90 % v/v **SLR** specified reagents for general laboratory work- extra pure

packaging		Mfr. No
500 mℓ	Glass bottle	E/0400/08
2.5 ℓ	Glass bottle	E/0400/17
C ₂ H ₅ OH CAS: 64-17-5 MW: 46.06 Density at 20°C: 0.790 Boiling point: 78°C Melting point: -114.1°C Flash point: 16.66°C		ADR/RID: 3, II EINECS: 200-578-6 UN: 1170 GHS: Flam. Liq. 2 H225 P210, P240
Guaranteed analysis		
Residue after evaporation	<100 ppm	Methanol <0.1%

These are duty-paid items that are available in the UK market only

Ethanol, denatured with 2% MEK **NEW PACK.** **SLR** specified reagents for general laboratory work- extra pure

packaging		Mfr. No
5 ℓ	Fluorinated HDPE square Drum	E/0580/FP21
C ₂ H ₅ OH CAS: 64-17-5 MW: 46.06 Density at 20°C: 0.790 Boiling point: 78°C Melting point: -114°C Flash point: 13 to 17°C		ADR/RID: 3, II EINECS: 200-578-6 UN: 1170 GHS: Flam. Liq. 2 H225 P210, P240

Assay (CG): >95%

Guaranteed analysis	
Density at 20°C	: 0.788 to 0.791
Water	<0.3%

Ethanolamine **Certified AR** for analysis

packaging		Mfr. No
500 mℓ	Glass bottle	E/0701/08
2.5 ℓ	Glass bottle	E/0701/17
C ₂ H ₇ NO CAS: 141-43-5 MW: 61.08 Density at 20°C: 1.012 Boiling point: 170°C Melting point: 10.3°C Flash point: 85°C ADR/RID: 8, III		EINECS: 205-483-3 UN: 2491 GHS: Acute Tox. 4, Skin Corr. 1B, STOT SE 3 H335, H314, H302, H332, H312 P280, P305+P351+P338, P261, P308+P313
Assay (CG): >99.5%		
Guaranteed analysis		
Colour	<10 APHA	Potassium (K) <0.5 ppm
Residue after ignition	<0.005%	Sodium (Na) <5 ppm
Calcium (Ca)	<2 ppm	Total phosphorus (P) <0.5 ppm
Copper (Cu)	<1 ppm	Total silicon (Si) <0.2 ppm
Iron (Fe)	<2 ppm	Total sulphur (S) <1 ppm
Lead (Pb)	<1 ppm	Zinc (Zn) <1 ppm
Magnesium (Mg)	<0.5 ppm	

Lot analysis on the label.

Ethanolamine **SLR** specified reagents for general laboratory work- extra pure

packaging		Mfr. No
500 mℓ	Glass bottle	E/0700/08
2.5 ℓ	Glass bottle	E/0700/17
2.5 ℓ	Coated Glass bottle	E/0700/PB17
25 ℓ	Metal drum	E/0700/25
C ₂ H ₇ NO CAS: 141-43-5 MW: 61.08 Density at 20°C: 1.012 Boiling point: 170°C Melting point: 10.3°C Flash point: 85°C ADR/RID: 8, III		EINECS: 205-483-3 UN: 2491 GHS: Acute Tox. 4, Skin Corr. 1B, STOT SE 3 H335, H314, H302, H332, H312 P280, P305+P351+P338, P261, P308+P313
Assay (CG): >99%		
Guaranteed analysis		
Colour	<20 APHA	Sodium (Na) <20 ppm
Calcium (Ca)	<10 ppm	Total phosphorus (P) <2 ppm
Copper (Cu)	<5 ppm	Total silicon (Si) <1 ppm
Iron (Fe)	<10 ppm	Total sulphur (S) <5 ppm
Lead (Pb)	<5 ppm	Zinc (Zn) <2 ppm
Magnesium (Mg)	<2 ppm	Residue after ignition <0.02%
Potassium (K)	<2 ppm	

Ethenylbenzene See Styrene, page 333

Ethidium® bromide for electrophoresis

packaging		Mfr. No
5 g	Glass jar	E/P700/44
C ₂₁ H ₂₀ BrN ₃ CAS: 1239-45-8 MW: 394.30 Melting point: 260 to 262°C Flash point: >100°C ADR/RID: 6.1, I		EINECS: 214-984-6 UN: 2811 GHS: Muta. 2, Acute Tox. 2 H341, H330, H302 P310, P281, P304+P340

Assay: > 95%

Fluorescent dye for DNA/RNA in density gradients and electrophoresis gels.

Ethidium® bromide solution for electrophoresis

packaging		Mfr. No
10 mℓ	Glass bottle	E/P800/03
50 mℓ	Glass bottle	E/P800/10
C ₂₁ H ₂₀ BrN ₃ CAS: 1239-45-8 MW: 394.30 ADR/RID: 6.1, III EINECS: 214-984-6		UN: 2810 GHS: Muta. 2, Acute Tox. 3 H331, H341 P310, P281, P304+P340

2-Ethoxyethanol **Certified AR** for analysis

packaging		Mfr. No
500 mℓ	Glass bottle	E/0751/08
C ₄ H ₁₀ O ₂ CAS: 110-80-5 MW: 90.12 Density at 20°C: 0.930 Boiling point: 135°C Melting point: -100°C Flash point: 44°C ADR/RID: 3, III		EINECS: 203-804-1 UN: 1171 GHS: Flam. Liq. 3, Repr. 1B, Acute Tox. 4 H226, H360FD, H332, H312, H302 P280, P201, P308+P313, P302+P352, P304+P340, P210

Assay (CG): >99.5%

Guaranteed analysis	
Acidity / Alkalinity	<0.001 meq/g
Colour	<10 APHA
Water	<0.1%
Residue after evaporation	<50 ppm
Calcium (Ca)	<1 ppm
Copper (Cu)	<0.05 ppm
Iron (Fe)	<0.5 ppm
Lead (Pb)	<0.05 ppm
Magnesium (Mg)	<0.2 ppm
Potassium (K)	<0.5 ppm
Sodium (Na)	<2 ppm
Total phosphorus (P)	<0.5 ppm
Total silicon (Si)	<0.05 ppm
Total sulphur (S)	<2 ppm
Zinc (Zn)	<1 ppm

Typical analysis	
Aluminium (Al)	<0.1 ppm
Barium (Ba)	<0.02 ppm
Cadmium (Cd)	<0.05 ppm
Chromium (Cr)	<0.02 ppm
Cobalt (Co)	<0.02 ppm
Lithium (Li)	<0.02 ppm
Manganese (Mn)	<0.02 ppm
Molybdenum (Mo)	<0.02 ppm
Nickel (Ni)	<0.05 ppm
Strontium (Sr)	<0.02 ppm
Titanium (Ti)	<0.02 ppm
Vanadium (V)	<0.02 ppm

Lot analysis on the label

2-Ethoxyethanol **SLR** specified reagents for general laboratory work- extra pure

packaging		Mfr. No
1 ℓ	Glass bottle	E/0750/15
2.5 ℓ	Glass bottle	E/0750/17
25 ℓ	Metal drum	E/0750/25
C ₄ H ₁₀ O ₂ CAS: 110-80-5 MW: 90.12 Density at 20°C: 0.930 Boiling point: 135°C Melting point: -100°C Flash point: 44°C ADR/RID: 3, III		EINECS: 203-804-1 UN: 1171 GHS: Flam. Liq. 3, Repr. 1B, Acute Tox. 4 H226, H360FD, H332, H312, H302 P280, P201, P308+P313, P302+P352, P304+P340, P210

Assay (CG): >99%

Guaranteed analysis	
Acidity / Alkalinity	<0.005 meq/g
Colour	<20 APHA
Water	<0.2%
Residue after evaporation	<100 ppm
Calcium (Ca)	<5 ppm
Copper (Cu)	<2 ppm
Iron (Fe)	<5 ppm
Lead (Pb)	<2 ppm
Magnesium (Mg)	<2 ppm
Potassium (K)	<2 ppm
Sodium (Na)	<10 ppm
Total phosphorus (P)	<1 ppm
Total silicon (Si)	<1 ppm
Total sulphur (S)	<5 ppm
Zinc (Zn)	<2 ppm

2-Ethoxyethyl acetate
pure

packaging		Mfr. No
2.5 ℓ	Coated Glass bottle	E/0800/PB17
C ₆ H ₁₂ O ₃	UN: 1172	  
CAS: 111-15-9	GHS: Flam. Liq. 3, Repr. 1B,	
MW: 132.16	Acute Tox. 4	
Density at 20°C: 0.970	H226, H312, H360FD,	
Boiling point: 156°C	H302, H332	
Melting point: -61°C	P280, P308+P313,	
Flash point: 51°C	P304+P340, P302+P352,	
ADR/RID: 3, III	P301+P312, P210	
EINECS: 203-839-2		
Assay (CG): >98%		

Ethyl acetate **NEW PACK.**
for HPLC

Certified

packaging		Mfr. No
1 ℓ	Glass bottle	E/0906/15
2.5 ℓ	Glass bottle	E/0906/17
2.5 ℓ	Coated Glass bottle	E/0906/PB17
10 ℓ	Returnable Drum	E/0906/10RSS
25 ℓ	Returnable Drum	E/0906/25RSS
30 ℓ	Returnable Drum	E/0906/30RSS
C ₄ H ₈ O ₂	EINECS: 205-500-4	 
CAS: 141-78-6	UN: 1173	
MW: 88.11	GHS: Flam. Liq. 2, Eye Irrit.	
Density at 20°C: 0.902	2, STOT SE 3	
Boiling point: 75 to 78°C	H225, H319, H336, EUH066	
Melting point: -83.5°C	P261, P280,	
Flash point: -4°C	P305+P351+P338, P210, P240	
ADR/RID: 3, II		
Assay (CG): >99.5%		
Guaranteed analysis		
Acidity / Alkalinity	<0.0015 meq/g	Residue after evaporation
Water	<0.03%	<2 ppm
Absorbance specification		
Max. Absorbance (1cm cell v HPLC water)		
255 nm	: <1 A.U.	280 nm
260 nm	<0.12 A.U.	Wavelength at 1 A.U.
270 nm	<0.025 A.U.	: <255 nM
Lot analysis on the label.		
Filtered to 0.2 µm.		

Ethyl acetate
for HPLC, dried with water content <100pm

Certified

packaging		Mfr. No
2.5 ℓ	Glass bottle	E/0907/17
C ₄ H ₈ O ₂	EINECS: 205-500-4	 
CAS: 141-78-6	UN: 1173	
MW: 88.11	GHS: Flam. Liq. 2, Eye Irrit.	
Density at 20°C: 0.902	2, STOT SE 3	
Boiling point: 75 to 78°C	H225, H319, H336, EUH066	
Melting point: -83.5°C	P261, P280,	
Flash point: -4°C	P305+P351+P338, P210, P240	
ADR/RID: 3, II		
Assay (CG): >99.5%		
Guaranteed analysis		
Acidity / Alkalinity	<0.0015 meq/g	Residue after evaporation
Water	<0.01%	<2 ppm
Absorbance specification		
Max. Absorbance (1cm cell v HPLC water)		
255 nm	: <1 A.U.	280 nm
260 nm	<0.12 A.U.	Wavelength at 1 A.U.
270 nm	<0.025 A.U.	: <255 nM
Lot analysis on the label.		
Filtered to 0.2µm.		

Ethyl acetate
for residue and pesticide analysis

Distol

packaging		Mfr. No
1 ℓ	Glass bottle	E/0903/15
2.5 ℓ	Glass bottle	E/0903/17
C ₄ H ₈ O ₂	EINECS: 205-500-4	 
CAS: 141-78-6	UN: 1173	
MW: 88.11	GHS: Flam. Liq. 2, Eye Irrit.	
Density at 20°C: 0.902	2, STOT SE 3	
Boiling point: 75 to 78°C	H225, H319, H336, EUH066	
Melting point: -83.5°C	P261, P280,	
Flash point: -4°C	P305+P351+P338, P210, P240	
ADR/RID: 3, II		
Guaranteed analysis		
Water	<0.05%	GC (FID) Analysis
Residue after evaporation	<2 ppm	GC (NPD) analysis
GC (ECD) Analysis	Pass	Pass

GC Analysis
Analysis of pesticides: from the retention time of Lindane, after concentration of 2500 times, no impurity peak higher than that obtained with 10 pg / mL héptachlo-réoxide (ECD) and 5 pg / mL Parathion (NPD).
Analysis of hydrocarbons: no impurity chain longer than the dodecane is present at a concentration greater than 0.01 ppm (FID).
Conditions :
Column: J&W DB 5: 30m x 0.32mm x 0.1µm df.
Temperature profile: (1) 50°C, 2 min
(2) 40°C / min until 320°C
(3) 320°C for 30 min
Injection: 5µℓ cold on the column.

Ethyl acetate
for peptide sequencing

packaging		Mfr. No
2.5 ℓ	Glass bottle	E/0908/17
C ₄ H ₈ O ₂	EINECS: 205-500-4	 
CAS: 141-78-6	UN: 1173	
MW: 88.11	GHS: Flam. Liq. 2, Eye Irrit.	
Density at 20°C: 0.902	2, STOT SE 3	
Boiling point: 75 to 78°C	H225, H319, H336, EUH066	
Melting point: -83.5°C	P261, P280,	
Flash point: -4°C	P305+P351+P338, P210, P240	
ADR/RID: 3, II		
Assay (CG): >99.5%		
Guaranteed analysis		
Water	<0.03%	Acidity
Residue after evaporation	<2 ppm	<0.0015 meq/g
Absorbance specification		
Max. Absorbance (1cm cell v HPLC water)		
256 nm	: <1 A.U.	325 nm
275 nm	: <0.025 A.U.	: <0.005 A.U.

Ethyl acetate **NEW PACK.**
for analysis

Certified AR

packaging		Mfr. No
1 ℓ	Glass bottle	E/0900/15
2.5 ℓ	Glass bottle	E/0900/17
2.5 ℓ	Coated Glass bottle	E/0900/PB17
2.5 ℓ	HDPE Plastic bottle	E/0900/PC17
5 ℓ	Fluorinated HDPE Plastic bottle	E/0900/21
5 ℓ	Fluorinated HDPE square Drum	E/0900/FP21
25 ℓ	Metal drum	E/0900/25
200 ℓ	Metal drum	E/0900/27
C ₄ H ₈ O ₂	EINECS: 205-500-4	 
CAS: 141-78-6	UN: 1173	
MW: 88.11	GHS: Flam. Liq. 2, Eye Irrit.	
Density at 20°C: 0.902	2, STOT SE 3	
Boiling point: 75 to 78°C	H225, H319, H336, EUH066	
Melting point: -83.5°C	P261, P280,	
Flash point: -4°C	P305+P351+P338, P210, P240	
ADR/RID: 3, II		
Assay (CG): >99.8%		

Guaranteed analysis		
Acidity / Alkalinity	<0.0005 meq/g	Substances
Colour	<5 APHA	reducing KMnO ₄
Water	<0.03%	Total phosphorus (P)
Residue after evaporation	<10 ppm	Total silicon (Si)
iso-Propyl acetate	: <0.05%	Total sulphur (S)
n-Propyle acetate	: <0.05%	Zinc (Zn)
Copper (Cu)	<0.02 ppm	Ethanol 90 % v/v
Iron (Fe)	<0.1 ppm	Ethyl propanoate
Lead (Pb)	<0.02 ppm	Methanol
Magnesium (Mg)	<0.05 ppm	Methyl acetate
Potassium (K)	<0.5 ppm	Substances darkened
Sodium (Na)	<0.5 ppm	by H ₂ SO ₄
Typical analysis		
Aluminium (Al)	<0.1 ppm	Manganese (Mn)
Barium (Ba)	<0.01 ppm	Molybdenum (Mo)
Cadmium (Cd)	<0.01 ppm	Nickel (Ni)
Calcium (Ca)	<0.01 ppm	Strontium (Sr)
Chromium (Cr)	<0.01 ppm	Titanium (Ti)
Cobalt (Co)	<0.01 ppm	Vanadium (V)
Lithium (Li)	<0.01 ppm	

Lot analysis on the label.

Ethyl acetate **NEW PACK.**
specified reagents for general laboratory work- extra pure

SLR

packaging		Mfr. No
500 mℓ	Glass bottle	E/0850/08
1 ℓ	Glass bottle	E/0850/15
1 ℓ	Metal can	E/0850/MC15
2.5 ℓ	Glass bottle	E/0850/17
2.5 ℓ	Coated Glass bottle	E/0850/PB17
2.5 ℓ	HDPE Plastic bottle	E/0850/PC17
5 ℓ	Fluorinated HDPE Plastic bottle	E/0850/21
5 ℓ	Fluorinated HDPE square Drum	E/0850/FP21
10 ℓ	Metal can	E/0850/24
25 ℓ	Metal drum	E/0850/25
200 ℓ	Metal drum	E/0850/27
C ₄ H ₈ O ₂	EINECS: 205-500-4	 
CAS: 141-78-6	UN: 1173	
MW: 88.11	GHS: Flam. Liq. 2, Eye Irrit.	
Density at 20°C: 0.902	2, STOT SE 3	
Boiling point: 75 to 78°C	H225, H319, H336, EUH066	
Melting point: -83.5°C	P261, P280,	
Flash point: -4°C	P305+P351+P338, P210, P240	
ADR/RID: 3, II		
Assay (CG): >99%		
Guaranteed analysis		
Acidity / Alkalinity	<0.001 meq/g	Magnesium (Mg)
Colour	<10 APHA	Potassium (K)
Water	<0.05%	Sodium (Na)
Residue after evaporation	<50 ppm	Total phosphorus (P)
Substances		Total silicon (Si)
reducing KMnO ₄	<0.001%	Total sulphur (S)
Calcium (Ca)	<2 ppm	Zinc (Zn)
Copper (Cu)	<1 ppm	Substances darkened
Iron (Fe)	<2 ppm	by H ₂ SO ₄
Lead (Pb)	<1 ppm	<20 APHA

Ethyl acetoacetate
specified reagents for general laboratory work- extra pure

SLR

packaging		Mfr. No
500 mℓ	Glass bottle	E/0950/08
2.5 ℓ	Glass bottle	E/0950/17
C ₆ H ₁₀ O ₃	Boiling point: 180°C	
CAS: 141-97-9	Melting point: -45°C	
MW: 130.15	Flash point: 70°C	
Density at 20°C: 1.03	EINECS: 205-516-1	
Assay (CG): >99%		

Ethylamine
pure, solution 70% w/w in water

packaging		Mfr. No
500 mℓ	Glass bottle	E/1000/08
C ₂ H ₇ N	GHS: Flam. Liq. 2, Acute	  
CAS: 75-04-7	Tox. 3, Skin Corr. 1A, Eye	
MW: 45.90	Dam. 1, STOT SE 3	
Density at 20°C: 0.800	H225, H314, H332, H335,	
Boiling point: 39°C	H302, H311	
Melting point: -81.1°C	P301+P330+P331, P280,	
Flash point: -17°C	P305+P351+P338, P310,	
ADR/RID: 3, 8, II	P261, P301+P312,	
EINECS: 200-834-7	P304+P340, P210, P240	
UN: 2270		

Ethyl benzoate
specified reagents for general laboratory work- extra pure

SLR

Packaging		Mfr. No
500 mL	Glass bottle	E/1100/08
C ₉ H ₁₀ O ₂	Flash point: 84°C	
CAS: 93-89-0	EINECS: 202-284-3	
MW: 150.18	GHS: Skin Irrit. 2, Eye Irrit. 2	
Density at 20°C: 1.045	H315, H319	
Boiling point: 212°C	P302+P352, P280,	
Melting point: -34°C	P305+P351+P338	
Assay (CG): >99%		


Ethyl bromide See Bromoethane, page 61

Ethyl cinnamate
technical

packaging		Mfr. No
25 mℓ	Glass bottle	E/1250/04
C ₁₁ H ₁₂ O ₂	Density at 20°C: 1.045	
CAS: 103-36-6	EINECS: 203-104-6	
MW: 176.22		

Ethyl digol See Diethylene glycol monoethyl ether, page 109

Ethylencarboxamide See Acrylamide, page 26



Ethylenediaminetetraacetic acid		Certified AR
for analysis, EDTA		
packaging	Mfr. No	
100 g Plastic jar	D/0452/48	
250 g Plastic jar	D/0452/50	
C ₁₀ H ₁₆ N ₂ O ₈ CAS: 60-00-4 MW: 292.25 Density at 20°C: 0.86 Melting point: 220°C Flash point: >100°C	EINECS: 200-449-4 GHS: Acute Tox. 4, Eye Irrit. 2 H319, H332 P264, P280, P305+P351+P338, P337+P313	
Assay: >99%		

Guaranteed analysis			
Loss on drying	<0.1%	Sodium (Na)	<200 ppm
Calcium (Ca)	<50 ppm	Total chloride (Cl)	<0.02%
Copper (Cu)	<5 ppm	Total phosphorus (P)	<50 ppm
Iron (Fe)	<5 ppm	Total silicon (Si)	<50 ppm
Lead (Pb)	<5 ppm	Total sulphur (S)	<50 ppm
Magnesium (Mg)	<20 ppm	Zinc (Zn)	<20 ppm
Potassium (K)	<50 ppm		

Typical analysis			
Aluminium (Al)	<2 ppm	Lithium (Li)	<1 ppm
Silver (Ag)	<1 ppm	Manganese (Mn)	<1 ppm
Arsenic (As)	<2 ppm	Molybdenum (Mo)	<1 ppm
Barium (Ba)	<2 ppm	Nickel (Ni)	<1 ppm
Cadmium (Cd)	<1 ppm	Strontium (Sr)	<2 ppm
Chromium (Cr)	<1 ppm	Titanium (Ti)	<1 ppm
Cobalt (Co)	<1 ppm	Vanadium (V)	<1 ppm

Lot analysis on the label.

Ethylenediaminetetraacetic acid
extra pure, EDTA

SLR

packaging	Mfr. No	
250 g Plastic jar	D/0450/50	
500 g Plastic jar	D/0450/53	
10 kg Plastic keg	D/0450/68	
C ₁₀ H ₁₆ N ₂ O ₈ CAS: 60-00-4 MW: 292.25 Density at 20°C: 0.86 Melting point: 220°C Flash point: >100°C	EINECS: 200-449-4 GHS: Acute Tox. 4, Eye Irrit. 2 H319, H332 P264, P280, P305+P351+P338, P337+P313	
Assay: >99%		

Guaranteed analysis			
Calcium (Ca)	<200 ppm	Sodium (Na)	<1000 ppm
Copper (Cu)	<20 ppm	n-Butyl formate	<0.1%
Iron (Fe)	<20 ppm	Total phosphorus (P)	<200 ppm
Lead (Pb)	<20 ppm	Total silicon (Si)	<200 ppm
Magnesium (Mg)	<100 ppm	Total sulphur (S)	<200 ppm
Potassium (K)	<200 ppm	Zinc (Zn)	<100 ppm

Ethylenediaminetetraacetic acid diammonium salt
specified reagents for general laboratory work- extra pure

SLR

packaging	Mfr. No	
100 g Plastic jar	D/0500/48	
25 kg Plastic keg	D/0500/70	
(CH ₂) ₂ N(CH ₂ COOH).CH ₂ .COONH ₄) ₂ +aq. CAS: 20824-56-0	MW: 326.18 Melting point: 215 to 217°C EINECS: 244-063-4	
Guaranteed analysis		
Complies with standard IR spectrum: Pass		

Save money!
Order 4 bottles in a box



Ethylenediaminetetraacetic acid disodium salt dihydrate
for electrophoresis

packaging	Mfr. No	
100 g Plastic jar	E/P140/48	
500 g Plastic jar	E/P140/53	
5 kg Plastic keg	E/P140/65	
C ₁₀ H ₁₄ N ₂ Na ₂ O ₈ .2H ₂ O CAS: 6381-92-6 MW: 372.23 Melting point: 252°C	GHS: Acute Tox. 4 H332 P304+P312	
Assay: > 99.5%		

Ethylenediaminetetraacetic acid disodium salt dihydrate
Certified AR, Eur.Ph., BP

packaging	Mfr. No	
100 g Plastic jar	D/0700/48	
250 g Plastic jar	D/0700/50	
500 g Plastic jar	D/0700/53	
1 kg Plastic jar	D/0700/60	
5 kg Plastic jar	D/0700/65	
10 kg Plastic keg	D/0700/68	
25 kg Fibre Keg	D/0700/70	
C ₁₀ H ₁₄ N ₂ Na ₂ O ₈ .2H ₂ O CAS: 6381-92-6 MW: 372.23 Melting point: 252°C	EINECS: 205-358-3 GHS: Acute Tox. 4 H332 P304+P312	

Assay: >99%			
Guaranteed analysis			
pH (aqueous solution 5%)	4 to 5	Potassium (K)	<10 ppm
Calcium (Ca)	<20 ppm	Total chloride (Cl)	<0.005%
Copper (Cu)	<2 ppm	Total phosphorus (P)	<50 ppm
Iron (Fe)	<5 ppm	Total silicon (Si)	<50 ppm
Lead (Pb)	<2 ppm	Total sulphur (S)	<200 ppm
Magnesium (Mg)	<10 ppm	Zinc (Zn)	<20 ppm

Typical analysis			
Aluminium (Al)	<2 ppm	Lithium (Li)	<1 ppm
Silver (Ag)	<1 ppm	Manganese (Mn)	<1 ppm
Arsenic (As)	<2 ppm	Molybdenum (Mo)	<1 ppm
Barium (Ba)	<2 ppm	Nickel (Ni)	<1 ppm
Cadmium (Cd)	<1 ppm	Strontium (Sr)	<2 ppm
Chromium (Cr)	<1 ppm	Titanium (Ti)	<1 ppm
Cobalt (Co)	<1 ppm	Vanadium (V)	<1 ppm

Lot analysis on the label.

Ethylenediaminetetraacetic acid disodium salt dihydrate
specified reagents for general laboratory work- extra pure

SLR

packaging	Mfr. No	
100 g Plastic jar	D/0650/48	
250 g Plastic jar	D/0650/50	
500 g Plastic jar	D/0650/53	
1 kg Plastic jar	D/0650/60	
5 kg Plastic keg	D/0650/65	
10 kg Plastic keg	D/0650/68	
25 kg Fibre Keg	D/0650/70	
C ₁₀ H ₁₄ N ₂ Na ₂ O ₈ .2H ₂ O CAS: 6381-92-6 MW: 372.23 Melting point: 252°C	EINECS: 205-358-3 GHS: Acute Tox. 4 H332 P304+P312	

Assay: >98%			
Guaranteed analysis			
Calcium (Ca)	<100 ppm	Potassium (K)	<50 ppm
Copper (Cu)	<10 ppm	n-Butyl formate	<0.05%
Iron (Fe)	<20 ppm	Total phosphorus (P)	<200 ppm
Lead (Pb)	<10 ppm	Total silicon (Si)	<200 ppm
Magnesium (Mg)	<50 ppm	Total sulphur (S)	<1000 ppm

Ethylenediaminetetraacetic acid disodium salt solution 0.2N
HPLC for electrochemical detection

Certified

packaging	Mfr. No	
2.5 ℓ Glass bottle	J/3710/17	
C ₁₀ H ₁₄ N ₂ Na ₂ O ₈ CAS: 139-33-3 MW: 336.21	Melting point: 237 to 245°C EINECS: 205-358-3	
Guaranteed analysis		
Electrochemical test Pass		

Tested by HPLC-EC to ensure low electrochemical response.
Detector - ESA Coulochem II
Mode-oxidative
Graphite Electrode.
Lot analysis on the label.
Filtered to 0.2 μm.

Ethylenediaminetetraacetic acid disodium salt solution 0.01M (0.02N)
NIST Standard solution ready to use, For Volumetric analysis

packaging	Mfr. No	
1 ℓ HDPE Plastic bottle	J/3720/PB15	
2.5 ℓ HDPE Plastic bottle	J/3720/17	
10 ℓ Bag in box	J/3720/24	
C ₁₀ H ₁₄ N ₂ Na ₂ O ₈ CAS: 139-33-3 MW: 336.21	Melting point: 237 to 245°C EINECS: 205-358-3	

Manufactured to a tolerance of ± 0.3%.
See the full range of Standard Volumetric solution page 323.

Ethylenediaminetetraacetic acid disodium salt solution 0.01M (0.02N)
NIST Standard Concentrate, For Volumetric analysis

packaging	Mfr. No	
1 Ampoule	J/3720C/05	
6 Ampoule	J/3720C/90	
C ₁₀ H ₁₄ N ₂ Na ₂ O ₈ CAS: 139-33-3 MW: 336.21	Melting point: 237 to 245°C EINECS: 205-358-3	

Each ampoule makes 1 ℓ of solution.
See the full range of Standard Volumetric solution page 323.

Ethylenediaminetetraacetic acid disodium salt solution
0.1M (0.2N)

Eur.Ph. , BP

NIST Standard solution ready to use, For Volumetric analysis

packaging	Mfr. No	
1 ℓ HDPE Plastic bottle	J/3700/15	
2.5 ℓ HDPE Plastic bottle	J/3700/17	
10 ℓ Bag in box	J/3700/24	
C ₁₀ H ₁₄ N ₂ Na ₂ O ₈ CAS: 139-33-3 MW: 336.21	Melting point: 237 to 245°C EINECS: 205-358-3	

Manufactured to a tolerance of ± 0.1%.
See the full range of Standard Volumetric solution page 323.

Ethylenediaminetetraacetic acid trisodium salt solution 0.1M (0.2N)
NIST Standard Concentrate, For Volumetric analysis

packaging	Mfr. No	
1 Ampoule	J/3700C/05	
6 Ampoule	J/3700C/90	
C ₁₀ H ₁₄ N ₂ Na ₃ O ₈ CAS: 150-38-9 MW: 358.19	Density at 20°C: 1.0 Melting point: 237 to 245°C EINECS: 205-758-8	

Each ampoule makes 500 ml of solution.
See the full range of Standard Volumetric solution page 323.

Ethylenediaminetetraacetic acid dipotassium salt dihydrate
specified reagents for general laboratory work- extra pure

SLR

packaging	Mfr. No	
500 g Plastic jar	D/0600/53	
25 kg Fibre Keg	D/0600/70	
C ₁₀ H ₁₄ K ₂ N ₂ O ₈ .2H ₂ O CAS: 25102-12-9 MW: 404.47 Melting point: 272°C EINECS: 217-895-0	GHS: Skin Irrit. 2, Eye Irrit. 2, STOT SE 3 H319, H315, H335 P261, P302+P352, P280, P305+P351+P338	

Assay: >99%			
Guaranteed analysis			
Calcium (Ca)	<100 ppm	n-Butyl formate	<0.02%
Copper (Cu)	<10 ppm	Total phosphorus (P)	<200 ppm
Iron (Fe)	<20 ppm	Total silicon (Si)	<200 ppm
Lead (Pb)	<10 ppm	Total sulphur (S)	<1000 ppm
Magnesium (Mg)	<50 ppm	Zinc (Zn)	<100 ppm
Sodium (Na)	<2000 ppm		

Ethylenediaminetetraacetic acid tetrasodium salt
specified reagents for general laboratory work- extra pure

SLR

packaging	Mfr. No	
250 g Plastic jar	D/0750/50	
C ₁₀ H ₁₄ N ₂ O ₈ .4Na CAS: 13235-36-4 MW: 452.24 Melting point: >300°C Flash point: >100°C EINECS: 200-573-9	GHS: Acute Tox. 4, Eye Dam. 1 H302, H318, H332 P301+P330+P331, P280, P305+P351+P338, P310, P301+P312	 

Assay: >98%			
Guaranteed analysis			
pH (aqueous solution 5%)	 10 to 11.5	Lead (Pb)	<20 ppm
Iron (Fe)	<20 ppm	n-Butyl formate	<0.1%

Ethylene dibromide

See 1, page 103

Ethylene dichloride

See 1, page 105

(Ethylenedinitrilo)tetraacetic acid, disodium salt

See Ethylenediaminetetraacetic acid disodium salt solution, page 127

Ethanediol
specified reagents for general laboratory work- extra pure

SLR

packaging	Mfr. No	
500 mℓ HDPE Plastic bottle	E/0350/PB08	
1 ℓ Glass bottle	E/0350/15	
2.5 ℓ Glass bottle	E/0350/17	
2.5 ℓ HDPE Plastic bottle	E/0350/PB17	
25 ℓ HDPE Plastic container	E/0350/25	

C ₂ H ₆ O ₂ CAS: 107-21-1 MW: 62.07 Density at 20°C: 1.113 Boiling point: 196 to 198°C Melting point: -13°C	Flash point: 116°C EINECS: 203-473-3 GHS: Acute Tox. 4 H302 P280, P301+P310	
---	---	---

Assay (CG): >99%			
Guaranteed analysis			
Acidity / Alkalinity	<0.0005 meq/g	Magnesium (Mg)	<2 ppm
Colour	<20 APHA	Potassium (K)	<2 ppm
Water	<0.5%	Sodium (Na)	<10 ppm
Residue after ignition	<0.01%	Total phosphorus (P)	<5 ppm
Calcium (Ca)	<5 ppm	Total silicon (Si)	<1 ppm
Copper (Cu)	<5 ppm	Total sulphur (S)	<50 ppm
Iron (Fe)	<10 ppm	Zinc (Zn)	<2 ppm
Lead (Pb)	<5 ppm		

Ethylene glycol ethyl ether

See 2-Ethoxyethanol, page 123

Ethylene glycol monobutyl ether

See 2-n-Butoxyethanol, page 69

Ethylene glycol monoethyl ether acetate

See 2-Ethoxyethyl acetate, page 124

Ethylene glycol monomethyl ether

See 2-Methoxyethanol, page 203

Ethyl iodide

See Iodoethane, page 161

Ethyl lactate

specified reagents for general laboratory work- extra pure

packaging		Mfr. No
2.5 ℓ	Glass bottle	E/1400/17
CH ₃ CH(OH).COO.C ₂ H ₅		UN: 1192
CAS: 97-64-3		GHS: Flam. Liq. 3, Eye Dam.
MW: 118.14		1, STOT SE 3
Density at 25°C: 1.03		H226, H318, H335
Boiling point: 154°C		P261, P280,
Melting point: -23°C		P303+P361+P353,
Flash point: 46°C		P305+P351+P338,
ADR/RID: 3, III		P304+P340, P210
EINECS: 202-598-0		
Assay (CG): >98%		
Guaranteed analysis		
Residue after evaporation		<100 ppm

Typical analysis

Aluminium (Al)	<0.1 ppm	Manganese (Mn)	<0.01 ppm
Barium (Ba)	<0.01 ppm	Molybdenum (Mo)	<0.01 ppm
Cadmium (Cd)	<0.01 ppm	Nickel (Ni)	<0.01 ppm
Chromium (Cr)	<0.01 ppm	Strontium (Sr)	<0.01 ppm
Cobalt (Co)	<0.01 ppm	Titanium (Ti)	<0.01 ppm
Lithium (Li)	<0.01 ppm	Vanadium (V)	<0.01 ppm

Lot analysis on the label.

Ethyl methyl ketone

NEW PACK

SLR

specified reagents for general laboratory work- extra pure

packaging		Mfr. No
500 mℓ	Glass bottle	E/1450/08
1 ℓ	Glass bottle	E/1450/15
2.5 ℓ	Glass bottle	E/1450/17
2.5 ℓ	Coated Glass bottle	E/1450/PB17
5 ℓ	Fluorinated HDPE Plastic bottle	E/1450/21
25 ℓ	Metal drum	E/1450/25
200 ℓ	Metal drum	E/1450/27
C ₄ H ₈ O		EINECS: 201-159-0
CAS: 78-93-3		UN: 1193
MW: 72.11		GHS: Flam. Liq. 2, Eye Irrit.
Density at 20°C: 0.806		2, STOT SE 3
Boiling point: 80°C		H225, H319, H336, EUH066
Melting point: -87°C		P261, P280,
Flash point: -7°C		P305+P351+P338, P210, P240
ADR/RID: 3, II		
Assay (CG): >99%		
Guaranteed analysis		
Acidity / Alkalinity		<0.001 meq/g
Colour		<20 APHA
Water		<0.2%
Residue after evaporation		<50 ppm
Calcium (Ca)		<2 ppm
Copper (Cu)		<2 ppm
Iron (Fe)		<5 ppm
Lead (Pb)		<2 ppm
Magnesium (Mg)		<1 ppm
Potassium (K)		<1 ppm
Sodium (Na)		<2 ppm
Total phosphorus (P)		<1 ppm
Total silicon (Si)		<1 ppm
Total sulphur (S)		<2 ppm
Zinc (Zn)		<1 ppm

Ethyl methyl ketone

for analysis

Certified AR

packaging		Mfr. No
1 ℓ	Glass bottle	E/1455/15
1 ℓ	Coated Glass bottle	E/1455/PB15
2.5 ℓ	Glass bottle	E/1455/17
2.5 ℓ	Coated Glass bottle	E/1455/PB17
25 ℓ	Metal drum	E/1455/25
200 ℓ	Metal drum	E/1455/27
C ₄ H ₈ O		EINECS: 201-159-0
CAS: 78-93-3		UN: 1193
MW: 72.11		GHS: Flam. Liq. 2, Eye Irrit.
Density at 20°C: 0.806		2, STOT SE 3
Boiling point: 80°C		H225, H319, H336, EUH066
Melting point: -87°C		P261, P280,
Flash point: -7°C		P305+P351+P338, P210, P240
ADR/RID: 3, II		
Assay (CG): >99.5%		
Guaranteed analysis		
Acidity / Alkalinity		<0.0005 meq/g
Colour		<20 APHA
Water		<0.1%
Residue after evaporation		<10 ppm
Butanedione		: <1000 ppm
2-Methylpropan-2-ol		: <0.1%
Calcium (Ca)		<0.2 ppm
Copper (Cu)		<0.02 ppm
Iron (Fe)		<0.1 ppm
Lead (Pb)		<0.02 ppm
Magnesium (Mg)		<0.05 ppm
Potassium (K)		<0.5 ppm
Sodium (Na)		<0.5 ppm
Total phosphorus (P)		<0.1 ppm
Total silicon (Si)		<0.02 ppm
Total sulphur (S)		<0.2 ppm
Zinc (Zn)		<0.1 ppm
Acetone		<0.2%
Butan -2-ol		<0.02%
Butyl acetate		<0.05%
Di-n-butyl ether		<0.05%
Ethyl acetate		<0.05%
Methanol		<0.05%

Ethyl oleate

pure, for saponification

packaging		Mfr. No
500 mℓ	Glass bottle	E/1460/08
C ₂₀ H ₃₈ O ₂		Boiling point: 190 to 210°C
CAS: 111-62-6		Flash point: 160°C
MW: 310.52		EINECS: 203-889-5
Assay: 98 to 103%		
Guaranteed analysis		
Iodine value		75 to 84

Europium standard metal solution 1000 ppm

for ICP solution in HNO₃ 2%, NIST certified

packaging		Mfr. No
100 mℓ	Plastic bottle	SPPLEU22Y
500 mℓ	Plastic bottle	SPPLEU22X
H315, H319		

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Fast Green FCF C.I. 42053

pure, indicator grade

packaging		Mfr. No
25 g	Glass jar	F/P025/46
C ₃₇ H ₃₄ N ₂ Na ₂ O ₁₀ S ₃		EINECS: 219-091-5
CAS: 2353-45-9		GHS: Muta. 2
MW: 808.85		H341
Melting point: 290°C		P281
Guaranteed analysis		
Complies with standard IR spectrum: Pass		

Fehling's solution No. 1

pure, to B.P. 1988

packaging		Mfr. No
1 ℓ	HDPE Plastic bottle	J/3800/15
2.5 ℓ	HDPE Plastic bottle	J/3800/17
ADR/RID: 9, III		H411
UN: 3082		P273
GHS: Aquatic Chronic 2		
Assay: 6.85 to 7.05%		
Fehling's solution 1 and 2 are to be mixed in equal volumes before use.		
1 mℓ is equivalent to 5 mg glucose.		

Fehling's solution No. 2

pure, to B.P. 1988

packaging		Mfr. No
1 ℓ	HDPE Plastic bottle	J/3900/15
2.5 ℓ	HDPE Plastic bottle	J/3900/17
Density at 20°C: 1.24		H314
ADR/RID: 8, II		P301+P330+P331, P280,
UN: 1719		P305+P351+P338, P310
GHS: Skin Corr. 1A, Eye Dam. 1		
Fehling's solution 1 and 2 are to be mixed in equal volumes before use.		
1 mℓ is equivalent to 5 mg glucose.		

Fehling's solution No. 1

pure, to Lane and Eynon modification

packaging		Mfr. No
1 ℓ	HDPE Plastic bottle	J/3910/PB15
2.5 ℓ	Glass bottle	J/3910/17
ADR/RID: 9, III		H411
UN: 3082		P273
GHS: Aquatic Chronic 2		
Assay: 6.85 to 7.1%		
Fehling's solution 1 and 2 are to be mixed in equal volumes before use.		

Fehling's solution No. 2

pure, to Lane and Eynon modification

packaging		Mfr. No
1 ℓ	HDPE Plastic bottle	J/3920/PB15
2.5 ℓ	HDPE Plastic bottle	J/3920/17
Density at 20°C: 1.24		H314
ADR/RID: 8, II		P301+P330+P331, P280,
UN: 1719		P305+P351+P338, P310
GHS: Skin Corr. 1A, Eye Dam. 1		
Fehling's solution 1 and 2 are to be mixed in equal volumes before use.		

Ficodox® mixed COD reagent

pure

packaging		Mfr. No
2.5 ℓ	Coated Glass bottle	J/4000/PB17
Boiling point: 280°C		H314, H340, H350
ADR/RID: 8, II		P280, P301+P330+P331,
UN: 1830		P305+P351+P338,
GHS: Carc. 1B, Muta. 1B,		P302+P350, P201, P308+P313
Skin Corr. 1A		

Guaranteed analysis
Dichromate demand (mℓ FeSO₄) ... <2 mℓ

COD Test Reagents
Ficodox meets the DOE/NWC Standing Committee ot Analyst publication "Method for the Examination of Waters and Associated Materials"- Chemical Oxygen Demand (Dichromate value) of Polluted and Waste Waters- HMSO 1977.

Ficodox® Plus mixed COD reagent

pure

packaging		Mfr. No
1 ℓ	Glass bottle	J/4010/15
2.5 ℓ	Coated Glass bottle	J/4010/PB17
ADR/RID: 8, II		H314, H340, H350
UN: 1830		P280, P301+P330+P331,
GHS: Carc. 1B, Muta. 1B,		P305+P351+P338,
Skin Corr. 1A		P302+P350, P201, P308+P313

Guaranteed analysis
Dichromate demand (mℓ FeSO₄) ... <2 mℓ

COD Test reagents. Avoids the use of Mercury salt.
"Ficodox Plus "meets the DOE/NWC Standing Committee ot Analyst publication "Method for the Examination of Waters and Associated Materials"- Chemical Oxygen Demand (Dichromate value) of Polluted and Waste Waters 1986(Second Edition) available from HMSO for the estimation of COD without the use of mercury salts.

Florisil®

For chromatography, about 100-200 U.S. Mesh

packaging		Mfr. No
250 g	Plastic jar	F/1173/50
MgO ₃ Si		
CAS: 1343-88-0		
MW: 100.39		
EINECS: 215-681-1		
Guaranteed analysis		
Identity test		: Pass
Florisil® is a registered trade name of Floridin company.		

Florisil®

For chromatography, about 60-100 U.S. Mesh

packaging		Mfr. No
250 g	Plastic jar	F/1170/50
MgO ₃ Si		
CAS: 1343-88-0		
MW: 100.39		
EINECS: 215-681-1		
Guaranteed analysis		
Identity test		: Pass
Florisil® is a registered trade name of Floridin company.		

Fluorescein C.I. 45350

<i>pure</i>	
packaging	Mfr. No
25 g Glass jar	F/1200/46
100 g Glass jar	F/1200/48
C ₂₀ H ₁₂ O ₅ CAS: 2321-07-5 MW: 332.31 Melting point: 314 to 316°C	EINECS: 219-031-8 GHS: Eye Irrit. 2, Skin Sens. 1 H317, H319 P280, P305+P351+P338
Guaranteed analysis Complies with standard IR spectrum: Pass	



Fluorescein sodium

<i>pure</i>	
packaging	Mfr. No
100 g Glass jar	F/1300/48
500 g Glass jar	F/1300/53
C ₂₀ H ₁₀ Na ₂ O ₅ CAS: 518-47-8 MW: 376.28 Melting point: 320°C	EINECS: 208-253-0 GHS: Eye Irrit. 2 H319 P280, P305+P351+P338
Guaranteed analysis Complies with standard IR spectrum: Pass	



Fluorescein di-(methylimino diacetic acid) sodium salt

<i>pure, Metal indicator (Calcein)</i>	
packaging	Mfr. No
5 g Glass jar	F/1250/44
C ₃₀ H ₂₆ N ₂ O ₁₃ CAS: 1461-15-0 MW: 622.54 Melting point: 200°C EINECS: 215-957-1 GHS: Acute Tox. 4, Skin Irrit.	2, Eye Irrit. 2, STOT SE 3 H312, H332, H335, H315, H302, H319 P280, P261, P302+P352, P301+P312, P304+P340, P305+P351+P338
Guaranteed analysis Complies with standard IR spectrum: Pass	



Fluoride standard 10 ppm

<i>For Ion chromatography in H₂O, NIST certified</i>	
packaging	Mfr. No
100 mℓ Plastic bottle	SPASF91AY
500 mℓ Plastic bottle	SPASF91AX
<i>Manufactured to a tolerance of ± 0.3%. Check with your local distributor for availability of Spex product</i>	

Fluoride standard 100 ppm

<i>For Ion chromatography in H₂O, NIST certified</i>	
packaging	Mfr. No
100 mℓ Plastic bottle	SPASF91Y
500 mℓ Plastic bottle	SPASF91X
<i>Manufactured to a tolerance of ± 0.3%. Check with your local distributor for availability of Spex product</i>	

Fluoride standard solution 1000 ppm (1 mg/ml)

<i>for Ion pair chromatography, Ready-to-use</i>	
packaging	Mfr. No
100 mℓ HDPE Plastic bottle	J/4548/05
500 mℓ HDPE Plastic bottle	J/4548/08
Density at 20°C: 1.0	
<i>Manufactured to a tolerance of ± 1%.</i>	

Fluoride standard 1000 ppm

<i>For Ion chromatography in H₂O, NIST certified</i>	
packaging	Mfr. No
100 mℓ Plastic bottle	SPASF92Y
500 mℓ Plastic bottle	SPASF92X
<i>Manufactured to a tolerance of ± 0.3%. Check with your local distributor for availability of Spex product</i>	

Fluoroboric acid 50% solution in HBF₄

<i>Technical</i>	
packaging	Mfr. No
500 mℓ HDPE Plastic bottle	F/1350/08
HBF ₄ CAS: 16872-11-0 MW: 87.81 Density at 20°C: 1.41 Boiling point: 130°C Melting point: -90°C ADR/RID: 8, II	EINECS: 240-898-3 UN: 1775 GHS: Skin Corr. 1B, Eye Dam. 1 H314 P280, P301+P330+P331, P310, P305+P351+P338, P303+P361+P353



Fluorosilicic acid

<i>Technical</i>	
packaging	Mfr. No
500 mℓ HDPE Plastic bottle	F/1400/08
H ₂ F ₆ Si CAS: 16961-83-4 MW: 144.09 Density at 20°C: 1.27 ADR/RID: 8, II EINECS: 241-034-8	UN: 1778 GHS: Skin Corr. 1B, Eye Dam. 1 H314 P280, P260, P264, P301+P330+P331, P305+P351+P338, P310



Folin & Ciocalteu's phenol reagent

<i>pure</i>	
packaging	Mfr. No
500 mℓ Glass bottle	J/4100/08
ADR/RID: 8, II UN: 1760 GHS: Skin Corr. 1B, Eye Dam. 1	H314 P280, P305+P351+P338
<i>Contains 10% v / v HCl and 5% V / V H 3 PO 4.</i>	



Formaldehyde solution 37-41%

<i>specified reagents for general laboratory work- extra pure</i>		SLR
packaging	Mfr. No	
500 mℓ HDPE Plastic bottle	F/1501/PB08	
1 ℓ HDPE Plastic bottle	F/1501/PB15	
2.5 ℓ HDPE Plastic bottle	F/1501/PB17	
25 ℓ HDPE Plastic container	F/1501/25	
CH ₂ O CAS: 50-00-0 MW: 30.03 Density at 20°C: 1.083 Boiling point: 97°C Melting point: -15°C Flash point: 50°C ADR/RID: 3, 8, III EINECS: 200-001-8 UN: 1198	GHS: Flam. Liq. 3, Carc. 2, Acute Tox. 3, Skin Corr. 1B, Skin Sens. 1, STOT SE 1 H226, H314, H317, H370, H311, H331, H301, H351 P280, P301+P330+P331, P305+P351+P338, P304+P340, P302+P352, P301+P312, P210	
Assay: 37 to 41%		

Guaranteed analysis		
Acidity / Alkalinity	<0.05 meq/g	Potassium (K) <5 ppm
Colour	<20 APHA	n-Butyl formate <0.005%
Calcium (Ca)	<50 ppm	Total silicon (Si) <10 ppm
Copper (Cu)	<2 ppm	Total sulphur (S) <100 ppm
Iron (Fe)	<5 ppm	Zinc (Zn) <10 ppm
Lead (Pb)	<2 ppm	Methanol 9 to 14%
Magnesium (Mg)	<20 ppm	

Formaldehyde solution (stabilised in methanol)

<i>Technical</i>		Mfr. No
packaging		F/1451/PB17
2.5 ℓ HDPE Plastic bottle		F/1451/25
25 ℓ HDPE Plastic container		
H.CHO CAS: 50-00-0 MW: 30.03 Boiling point: 97°C ADR/RID: 3, 8, III EINECS: 200-001-8 UN: 1198 GHS: Flam. Liq. 3, Carc. 2,	Acute Tox. 3, Skin Corr. 1B, Skin Sens. 1, STOT SE 1 H226, H314, H317, H370, H311, H331, H301, H351 P280, P301+P330+P331, P305+P351+P338, P304+P340, P302+P352, P301+P312, P210	
Assay: >37 to 41 %%		
Guaranteed analysis		
Acidity / Alkalinity	<0.05 meq/g	

Formaldehyde dimethyl acetal

See Dimethoxymethane, page 112

Formaldehyde solution

<i>pure, 10% phosphate buffered</i>	
packaging	Mfr. No
5 ℓ HDPE Plastic bottle	F/1520/21
GHS: Carc. 2, Skin Sens. 1 H351, H317 P201, P261, P281, P308+P313, P302+P352	
Guaranteed analysis	
pH at 20°C	: 7.2 to 7.4

Formaldehyde solution saline

<i>pure, phosphate buffered pH 7.3</i>	
packaging	Mfr. No
25 ℓ HDPE Plastic container	F/1515/25
GHS: Carc. 2, Skin Irrit. 2, Eye Irrit. 2, Skin Sens. 1, STOT SE 3 H317, H351, H315, H319, H335 P280, P305+P351+P338, P260	



Formaldehyde solution

<i>pure, 10% w/v in 0.9% sodium chloride solution</i>	
packaging	Mfr. No
5 ℓ HDPE Plastic bottle	F/1510/21
GHS: Carc. 2, Skin Sens. 1 H317, H351 P280, P201, P302+P352, P262, P260, P308+P313	

Formamide

<i>for electrophoresis</i>	
packaging	Mfr. No
100 mℓ Coated Glass bottle	F/P565/PB05
500 mℓ Coated Glass bottle	F/P565/PB08
CH ₃ NO CAS: 75-12-7 MW: 45.04 Density at 20°C: 1.133 Boiling point: 210°C Melting point: 2 to 3°C	Flash point: 175°C EINECS: 200-842-0 GHS: Repr. 1B H360D P201, P308+P313
Assay (CG): >99.5%	
Guaranteed analysis	
Copper (Cu)	<0.1 ppm
Iron (Fe)	<1 ppm
Lead (Pb)	<0.5 ppm
Zinc (Zn)	<0.5 ppm
<i>Suitable for denaturing nucleic acids for sequencing and hybridisation.</i>	

Formamide

<i>for analysis</i>		Certified AR
packaging	Mfr. No	
500 mℓ Coated Glass bottle	F/1551/PB08	
2.5 ℓ Coated Glass bottle	F/1551/PB17	
25 ℓ HDPE Plastic container	F/1551/25	
200 ℓ HDPE Plastic container	F/1551/27	
CH ₃ NO CAS: 75-12-7 MW: 45.04 Density at 20°C: 1.133 Boiling point: 210°C Melting point: 2 to 3°C	Flash point: 175°C EINECS: 200-842-0 GHS: Repr. 1B H360D P201, P308+P313	
Assay (CG): >99.5%		
Guaranteed analysis		
Colour	<10 APHA	Magnesium (Mg) <5 ppm
Water	<0.1%	Potassium (K) <5 ppm
Ammonium formate (H.COO.NH ₄).:	<0.1 %	Sodium (Na) <2 ppm
Formic acid (H. COOH)	: <0.02%	Total phosphorus (P) <0.5 ppm
Calcium (Ca)	<10 ppm	Total silicon (Si) <0.5 ppm
Copper (Cu)	<1 ppm	Total sulphur (S) <1 ppm
Iron (Fe)	<2 ppm	Zinc (Zn) <2 ppm
Lead (Pb)	<0.5 ppm	Residue after ignition <0.5%

Typical analysis		
Formic acid (H. COOH)	<0.02 %	Manganese (Mn) <0.1 ppm
Aluminium (Al)	<0.2 ppm	Molybdenum (Mo) <0.1 ppm
Barium (Ba)	<0.2 ppm	Nickel (Ni) <0.1 ppm
Cadmium (Cd)	<0.1 ppm	Strontium (Sr) <0.1 ppm
Chromium (Cr)	<0.2 ppm	Titanium (Ti) <0.1 ppm
Cobalt (Co)	<0.1 ppm	Vanadium (V) <0.1 ppm
Lithium (Li)	<0.1 ppm	
<i>Lot analysis on the label.</i>		



Formamide

specified reagents for general laboratory work- extra pure

SLR

packaging	Mfr. No		
500 mℓ Coated Glass bottle	F/1550/PB08		
2.5 ℓ Coated Glass bottle	F/1550/PB17		
CH ₂ NO	Flash point: 175°C		
CAS: 75-12-7	EINECS: 200-842-0		
MW: 45.04	GHS: Repr. 1B		
Density at 20°C: 1.133	H360D		
Boiling point: 210°C	P201, P308+P313		
Melting point: 2 to 3°C			
Assay (CG): >99%			
Guaranteed analysis			
Colour	<20 APHA	Potassium (K)	<20 ppm
Water	<0.5%	Sodium (Na)	<50 ppm
Calcium (Ca)	<50 ppm	Total phosphorus (P)	<2 ppm
Copper (Cu)	<5 ppm	Total silicon (Si)	<2 ppm
Iron (Fe)	<10 ppm	Total sulphur (S)	<5 ppm
Lead (Pb)	<2 ppm	Zinc (Zn)	<10 ppm
Magnesium (Mg)	<20 ppm	Residue after ignition	<1000%

Formic acid

NEW

for LC/MS analysis

Optima

packaging	Mfr. No		
0.5 mℓ Ampoule	A117-05AMP		
1 mℓ Ampoule	A117-1AMP		
2 mℓ Ampoule	A117-2AMP		
10x1 mℓ Ampoule	A117-10X1AMP		
50 mℓ HDPE Plastic bottle	A117-50		
CH ₂ O ₂	UN: 1779		
CAS: 64-18-6	GHS: Flam. Liq. 3, Skin Corr.		
MW: 46.02	1A, Eye Dam. 1		
Density at 20°C: 1.22	H226, H314		
Boiling point: 101°C	P280, P260,		
Melting point: 8°C	P301+P330+P331,		
Flash point: 50°C	P305+P351+P338,		
ADR/RID: 8, 3, II	P302+P352, P210		
EINECS: 200-579-1			
Assay: <99.5%			
Guaranteed analysis			
Colour	<10 APHA	Lead (Pb)	<10 ppb
Reside after evaporation (at 0.1% aqueous)	: <1 ppm	Magnesium (Mg)	<10 ppb
Sulfite (SO ₃)	: <20 ppm	Manganese (Mn)	<10 ppb
Aluminium (Al)	<20 ppb	Nickel (Ni)	<10 ppb
Calcium (Ca)	<50 ppb	Potassium (K)	<20 ppb
Chloride (Cl)	<5 ppm	Propan--1-ol	<10 ppb
Copper (Cu)	<10 ppb	Sodium (Na)	<50 ppb
Iron (Fe)	<10 ppb	Sulfate (SO ₄)	<50 ppb
		Zinc (Zn)	<20 ppb
Absorbance specification			
Max. Absorbance (1cm cell v HPLC water)			
210 nm	<1.25 A.U.	254 nm	<0.01 A.U.
220 nm	<0.85 A.U.	LC/MS Gradient Suitability (at 0.1% aqueous)	: Pass test
230 nm	<0.55 A.U.		

This product has been tested on an LC-MS Thermo Scientific instrument to ensure suitability.

Formic acid

for analysis

Certified AR

packaging	Mfr. No
500 mℓ Coated Glass bottle	F/1900/PB08
1 ℓ Coated Glass bottle	F/1900/PB15
2.5 ℓ Coated Glass bottle	F/1900/PB17
25 ℓ HDPE Plastic container	F/1900/25
CH ₂ O ₂	UN: 1779
CAS: 64-18-6	GHS: Flam. Liq. 3, Skin Corr.
MW: 46.02	1A, Eye Dam. 1
Density at 20°C: 1.22	H226, H314
Boiling point: 101°C	P280, P260,
Melting point: 8°C	P301+P330+P331,
Flash point: 50°C	P305+P351+P338,
ADR/RID: 8, 3, II	P302+P352, P210
EINECS: 200-579-1	

Assay: >98%

Guaranteed analysis

Colour	<10 APHA	Potassium (K)	<0.5 ppm
Calcium (Ca)	<2 ppm	Sodium (Na)	<1 ppm
Chromium (Cr)	<0.05 ppm	Total chloride (Cl)	<0.0005%
Copper (Cu)	<0.05 ppm	Total nitrogen (N)	<20 ppm
Iron (Fe)	<0.5 ppm	Total phosphorus (P)	<5 ppm
Lead (Pb)	<0.05 ppm	Total silicon (Si)	<2 ppm
Magnesium (Mg)	<0.1 ppm	Total sulphur (S)	<10 ppm
Nickel (Ni)	<0.05 ppm	Zinc (Zn)	<0.1 ppm
Typical analysis			
Aluminium (Al)	<0.1 ppm	Manganese (Mn)	<0.02 ppm
Silver (Ag)	<0.02 ppm	Molybdenum (Mo)	<0.02 ppm
Barium (Ba)	<0.05 ppm	Strontium (Sr)	<0.02 ppm
Cadmium (Cd)	<0.05 ppm	Titanium (Ti)	<0.02 ppm
Cobalt (Co)	<0.02 ppm	Vanadium (V)	<0.02 ppm
Lithium (Li)	<0.02 ppm		

Lot analysis on the label.

Formic acid 90%

for analysis

Certified AR

packaging	Mfr. No		
500 mℓ HDPE Plastic bottle	F/1820/PB08		
2.5 ℓ HDPE Plastic bottle	F/1820/PB17		
CH ₂ O ₂	UN: 1779		
CAS: 64-18-6	GHS: Flam. Liq. 3, Skin Corr.		
MW: 46.03	1A, Eye Dam. 1		
Density at 20°C: 1.22	H226, H314		
Boiling point: 101°C	P280, P260,		
Melting point: 8°C	P301+P330+P331,		
Flash point: 50°C	P305+P351+P338,		
ADR/RID: 8, 3, II	P302+P352, P210		
EINECS: 200-579-1			
Assay: 89 to 91%			
Guaranteed analysis			
Colour	<10 APHA	Potassium (K)	<0.5 ppm
Calcium (Ca)	<2 ppm	Sodium (Na)	<2 ppm
Chromium (Cr)	<0.05 ppm	Total chloride (Cl)	<0.0005%
Copper (Cu)	<0.1 ppm	Total nitrogen (N)	<200 ppm
Iron (Fe)	<0.5 ppm	Total phosphorus (P)	<5 ppm
Lead (Pb)	<0.05 ppm	Total silicon (Si)	<2 ppm
Magnesium (Mg)	<0.1 ppm	Total sulphur (S)	<10 ppm
Nickel (Ni)	<0.05 ppm	Zinc (Zn)	<0.1 ppm
Typical analysis			
Aluminium (Al)	<0.1 ppm	Manganese (Mn)	<0.02 ppm
Silver (Ag)	<0.02 ppm	Molybdenum (Mo)	<0.02 ppm
Barium (Ba)	<0.05 ppm	Strontium (Sr)	<0.02 ppm
Cadmium (Cd)	<0.05 ppm	Titanium (Ti)	<0.02 ppm
Cobalt (Co)	<0.05 ppm	Vanadium (V)	<0.02 ppm
Lithium (Li)	<0.02 ppm		

Lot analysis on the label.

Formic acid

extra pure, min.98%

SLR

packaging	Mfr. No
500 mℓ Coated Glass bottle	F/1850/PB08
2.5 ℓ Coated Glass bottle	F/1850/PB17
200 ℓ HDPE Plastic container	F/1850/27
CH ₂ O ₂	UN: 1779
CAS: 64-18-6	GHS: Flam. Liq. 3, Skin Corr.
MW: 46.02	1A, Eye Dam. 1
Density at 20°C: 1.22	H226, H314
Boiling point: 101°C	P280, P260,
Melting point: 8°C	P301+P330+P331,
Flash point: 50°C	P305+P351+P338,
ADR/RID: 8, 3, II	P302+P352, P210
EINECS: 200-579-1	
Assay: >98%	

Guaranteed analysis

Colour	<20 APHA	Potassium (K)	<2 ppm
Calcium (Ca)	<10 ppm	Sodium (Na)	<10 ppm
Chromium (Cr)	<0.5 ppm	n-Butyl formate	<0.002%
Copper (Cu)	<2 ppm	Total phosphorus (P)	<20 ppm
Iron (Fe)	<5 ppm	Total silicon (Si)	<10 ppm
Lead (Pb)	<2 ppm	Total sulphur (S)	<50 ppm
Magnesium (Mg)	<1 ppm	Zinc (Zn)	<0.5 ppm
Nickel (Ni)	<0.5 ppm		

Formic acid

specified reagents for general laboratory work- extra pure

SLR

packaging	Mfr. No		
500 mℓ HDPE Plastic bottle	F/1800/PB08		
2.5 ℓ HDPE Plastic bottle	F/1800/PB17		
200 ℓ HDPE Plastic container	F/1800/27		
CH ₂ O ₂	UN: 1779		
CAS: 64-18-6	GHS: Flam. Liq. 3, Skin Corr.		
MW: 46.02	1A, Eye Dam. 1		
Density at 20°C: 1.22	H226, H314		
Boiling point: 101°C	P280, P260,		
Melting point: 8°C	P301+P330+P331,		
Flash point: 50°C	P305+P351+P338,		
ADR/RID: 8, 3, II	P302+P352, P210		
EINECS: 200-579-1			
Assay: 89 to 91%			
Guaranteed analysis			
Colour	<20 APHA	Potassium (K)	<2 ppm
Calcium (Ca)	<10 ppm	Sodium (Na)	<10 ppm
Chromium (Cr)	<0.5 ppm	n-Butyl formate	<0.002%
Copper (Cu)	<2 ppm	Total phosphorus (P)	<20 ppm
Iron (Fe)	<5 ppm	Total silicon (Si)	<10 ppm
Lead (Pb)	<2 ppm	Total sulphur (S)	<50 ppm
Magnesium (Mg)	<1 ppm	Zinc (Zn)	<0.5 ppm
Nickel (Ni)	<0.5 ppm		

Formic acid

pure, about 85%

packaging	Mfr. No
1 ℓ HDPE Plastic bottle	F/1750/PB15
2.5 ℓ HDPE Plastic bottle	F/1750/PB17
25 ℓ HDPE Plastic container	F/1750/25
CH ₂ O ₂	UN: 1779
CAS: 64-18-6	GHS: Flam. Liq. 3, Skin Corr.
MW: 46.02	1B, Eye Dam. 1
Density at 20°C: 1.22	H226, H314
Boiling point: 101°C	P280, P260,
Melting point: 8°C	P305+P351+P338,
Flash point: 50°C	P301+P330+P331,
ADR/RID: 8, 3, II	P302+P352, P210
EINECS: 200-579-1	
Assay: 84 to 86%	

Formic acid solution 0.25M

for HPLC

Certified

packaging	Mfr. No
500 mℓ Coated Glass bottle	J/4200/PB08
EINECS: 200-579-1	
Absorbance specification	
Max. Absorbance (1cm cell v HPLC water)	
Drift in gradient baseline	<40 mA.U.
Absorbance of the largest eluted peak:	<1.5 mA.U.
Molarity between 0.24 and 0.25 mol / ℓ	

Formol See Formaldehyde solution 37-41%, page 130

D(-)-Fructose

for analysis

Certified AR

packaging	Mfr. No		
100 g Plastic jar	F/1952/48		
500 g Plastic jar	F/1952/53		
1 kg Plastic jar	F/1952/60		
C ₆ H ₁₂ O ₆	Melting point: 103 to 105°C		
CAS: 57-48-7	EINECS: 200-333-3		
MW: 180.16			
Guaranteed analysis			
Loss on drying	<0.5%	Potassium (K)	<10 ppm
Specific rotation	-93.5° to -91°	Sodium (Na)	<50 ppm
Heavy metals	<5 ppm	Total chloride (Cl)	<0.001%
Limit of Hydroxymethylfurfural test	: Pass	Total phosphorus (P)	<20 ppm
Acidity	Pass	Total silicon (Si)	<5 ppm
Calcium (Ca)	<10 ppm	Total sulphur (S)	<50 ppm
Copper (Cu)	<1 ppm	Zinc (Zn)	<2 ppm
Iron (Fe)	<1 ppm	Sulfate (SO ₄)	<0.025%
Lead (Pb)	<1 ppm	Residue after ignition	<0.5%
Magnesium (Mg)	<5 ppm		
Typical analysis			
Aluminium (Al)	<2 ppm	Lithium (Li)	<1 ppm
Arsenic (As)	<1 ppm	Manganese (Mn)	<1 ppm
Barium (Ba)	<1 ppm	Nickel (Ni)	<1 ppm
Cadmium (Cd)	<1 ppm	Titanium (Ti)	<1 ppm
Chromium (Cr)	<1 ppm	Vanadium (V)	<1 ppm
Cobalt (Co)	<1 ppm		

Lot analysis on the label.

D(-)-Fructose

pure

packaging	Mfr. No
100 g Plastic jar	F/1950/48
250 g Plastic jar	F/1950/50
C ₆ H ₁₂ O ₆	Melting point: 103 to 105°C
CAS: 57-48-7	EINECS: 200-333-3
MW: 180.16	
Guaranteed analysis	
Complies with standard IR spectrum:	Pass

Fuchsin basic

technical

packaging	Mfr. No
25 g Glass jar	F/P770/46
C ₂₀ H ₁₉ N ₃ .HCl	EINECS: 211-189-6
CAS: 632-99-5	UN: 3143
MW: 337.84	GHS: Carc. 2, Acute Tox. 4
Melting point: 250°C	H351, H302
Flash point: 200°C	P281, P301+P312
ADR/RID: 6.1, III	

Pararosaniline

Fumaric acid

specified reagents for general laboratory work- extra pure

SLR

packaging	Mfr. No
500 g Plastic jar	F/2070/53
C ₄ H ₄ O ₄	Flash point: 230°C
CAS: 110-17-8	EINECS: 203-743-0
MW: 116.07	GHS: Eye Irrit. 2
Density at 20°C: 1.62	H319
Melting point: 295 to 300°C	P280, P305+P351+P338
Assay: >99%	

Furfuraldehyde
technical

packaging

2.5 l Coated Glass bottle

Mfr. No
F/2080/PB17

C₅H₄O₂
CAS: 98-01-1
MW: 96.09
Density at 20°C: 1.16
Boiling point: 159 to 161°C
Melting point: -36°C
Flash point: 60°C
ADR/RID: 6.1, 3, II
EINECS: 202-627-7

UN: 1199
GHS: Carc. 2, Acute Tox. 3,
Skin Irrit. 2, Eye Irrit. 2,
STOT SE 3
H319, H331, H312, H315,
H335, H351, H301
P260, P280, P301+P310,
P304+P340, P302+P352,
P305+P351+P338

2.5-Furandione See Maleic anhydride, page 190

Furfural See Furfuraldehyde , page 134

Fusion mixture
for analysis

Certified AR

packaging

500 g Plastic jar
1 kg Plastic jar

Mfr. No
F/2250/53
F/2250/60

K₂CO₃.Na₂CO₃
MW: 122.10
EINECS: 207-838-8
GHS: Skin Irrit. 2, Eye Irrit. 2,

STOT SE 3
H315, H319, H335
P261, P280, P301+P312,
P302+P352, P305+P351+P338

Guaranteed analysis

Loss on drying 300°C : <1%
Calcium (Ca) <100 ppm
Copper (Cu) <5 ppm
Iron (Fe) <5 ppm
Lead (Pb) <5 ppm
Magnesium (Mg) <10 ppm

Total chloride (Cl) <0.002%
Total nitrogen (N) <20 ppm
Total phosphorus (P) <10 ppm
Total silicon (Si) <20 ppm
Total sulphur (S) <10 ppm
Zinc (Zn) <5 ppm

Gadolinium standard metal solution 1000 ppm
for ICP solution in HNO₃ 2%, NIST certified

packaging

100 ml Plastic bottle
500 ml Plastic bottle
H315, H319

Mfr. No
SPPLGD22Y
SPPLGD22X

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Gadolinium standard metal solution 10 000 ppm
for ICP solution in HNO₃ 5%, NIST certified

packaging

100 ml Plastic bottle
500 ml Plastic bottle
H314

Mfr. No
SPPLGD23Y
SPPLGD23X

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

D(+)Galactose anhydrous
specified reagents for general laboratory work- extra pure

SLR

packaging

100 g Plastic jar

Mfr. No
G/0050/48

C₆H₁₂O₆
CAS: 59-23-4
MW: 180.16

Melting point: 165 to 170°C
EINECS: 200-416-4

Guaranteed analysis

Loss on drying <0.5%

Gallium standard metal solution 1000 ppm
for ICP solution in HNO₃ 2%, NIST certified

packaging

100 ml Plastic bottle
500 ml Plastic bottle
H315, H319

Mfr. No
SPPLGA22Y
SPPLGA22X

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurites.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Gelatine

SLR, Eur. Ph.

specified reagents for general laboratory work - extra pure, granular powder

packaging

250 g Plastic jar
500 g Plastic jar

Mfr. No
G/0150/50
G/0150/53

CAS: 9000-70-8
Density at 20°C: 0.68
EINECS: 232-554-6

Gentian Violet See Methyl violet, page 208

Germanium standard metal solution 1000 ppm
for ICP solution in H₂O/HF trace, NIST certified

packaging

100 ml Plastic bottle
500 ml Plastic bottle

Mfr. No
SPPLGE92Y
SPPLGE92X

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurites.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Germanium standard metal solution 10 ppm
for ICP-MS solution in H₂O/HF trace, NIST certified

packaging

100 ml Plastic bottle

Mfr. No
SPCLGE91AY

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurites.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Gibberellic acid
pure

packaging

1 g Vial

Mfr. No
G/0200/43

C₁₉H₂₂O₆
CAS: 77-06-5
MW: 346.38
Melting point: 227°C
EINECS: 201-001-0

GHS: Skin Corr. 1B, Eye Dam. 1
H314
P301+P330+P331, P280,
P305+P351+P338, P310

Guaranteed analysis

Complies with standard IR spectrum: Pass

Giemsa stain
pure

packaging

25 g Glass jar

Mfr. No
G/0275/46

CAS: 51811-82-6
Melting point: 300°C
EINECS: 257-438-2

Guaranteed analysis

Complies with standard IR spectrum: Pass

Glass balls
pure, 3.7-4.1 mm undrilled

packaging

250 g Plastic jar
500 g Plastic jar

Mfr. No
G/0300/50
G/0300/53

Guaranteed analysis

Identity test : Pass

D(+)-Glucose anhydrous
for analysis

Certified AR, Eur.Ph.,BP,USP

packaging

500 g Plastic jar
1 kg Plastic jar
2 kg Plastic jar
5 kg Plastic keg
25 kg Fibre Keg

Mfr. No
G/0500/53
G/0500/60
G/0500/61
G/0500/65
G/0500/70

C₆H₁₂O₆
CAS: 50-99-7
MW: 180.16

Density at 20°C: 15.4
Melting point: 146°C
EINECS: 200-075-1

Guaranteed analysis

pH (aqueous solution 5%) >6
Loss on drying <0.1%
Specific rotation 52.5° to 53.5°
Insoluble in alcohol : <0.005%
Calcium (Ca) <10 ppm
Copper (Cu) <1 ppm
Iron (Fe) <2 ppm
Lead (Pb) <0.5 ppm

Magnesium (Mg) <5 ppm
Potassium (K) <20 ppm
Sodium (Na) <50 ppm
Total chloride (Cl) <0.005%
Total phosphorus (P) <10 ppm
Total silicon (Si) <5 ppm
Total sulphur (S) <20 ppm
Zinc (Zn) <5 ppm

Typical analysis

Aluminium (Al) <2 ppm
Arsenic (As) <1 ppm
Barium (Ba) <1 ppm
Cadmium (Cd) <1 ppm
Chromium (Cr) <0.5 ppm
Cobalt (Co) <0.5 ppm

Lithium (Li) <0.2 ppm
Manganese (Mn) <0.2 ppm
Nickel (Ni) <0.5 ppm
Titanium (Ti) <0.5 ppm
Vanadium (V) <0.5 ppm

Lot analysis on the label.

Save money!
Order 4 bottles in a box



D(+)-Glucose anhydrous
specified reagents for general laboratory work- extra pure

SLR

packaging

500 g Plastic jar
1 kg Plastic jar
2 kg Plastic jar
5 kg Plastic keg

Mfr. No
G/0450/53
G/0450/60
G/0450/61
G/0450/65

C₆H₁₂O₆
CAS: 50-99-7
MW: 180.16

Density at 20°C: 15.4
Melting point: 146°C
EINECS: 200-075-1

Guaranteed analysis

Loss on drying <1%
Specific rotation 52° to 54°
Calcium (Ca) <50 ppm
Copper (Cu) <5 ppm
Iron (Fe) <10 ppm
Lead (Pb) <2 ppm
Magnesium (Mg) <20 ppm

Potassium (K) <100 ppm
Sodium (Na) <200 ppm
n-Butyl formate <0.02%
Total phosphorus (P) <50 ppm
Total silicon (Si) <20 ppm
Total sulphur (S) <100 ppm
Zinc (Zn) <20 ppm

D(+)-Glucose monohydrate
specified reagents for general laboratory work- extra pure

SLR, Eur.Ph.,BP

packaging

1 kg Plastic jar
5 kg Plastic keg

Mfr. No
G/0400/60
G/0400/65

C₆H₁₂O₆.H₂O
CAS: 5996-10-1
MW: 198.18

Density at 20°C: 1.54
EINECS: 200-075-1

Guaranteed analysis

Acidity / Alkalinity Pass
Loss on drying 7 to 9.5%
Specific rotation 52° to 53°
Foreign sugars/Starch/Dextrins (test) :
Pass
Sulfite (SO₃) : <15 ppm
Arsenic (As) <1 ppm

Barium (Ba) <1 ppm
Calcium (Ca) <200 ppm
Chloride (Cl) <125 ppm
Lead (Pb) <0.5 ppm
Sulfate (SO₄) <200 ppm
Residue after ignition <0.1%

Glucose-1-phosphoric acid dipotassium salt dihydrate
technical

packaging

2 g Vial

Mfr. No
G/0505/55

C₆H₁₁K₂O₈P.2H₂O
CAS: 5996-14-5
MW: 372.35
Melting point: 209°C

Guaranteed analysis

Complies with standard IR spectrum: Pass

L-Glutamic acid
specified reagents for general laboratory work- extra pure

SLR

packaging

100 g Plastic jar

Mfr. No
G/0515/48

C₅H₉NO₅
CAS: 56-86-0
MW: 147.13
Melting point: 205°C

Assay: >98.5%

Guaranteed analysis

Water <0.5%
Specific rotation 30° to 33°
Iron (Fe) <20 ppm

Lead (Pb) <10 ppm
Residue after ignition <0.5%

L(+)-Glutamic acid-5-amide
See L-Glutamine, page 136

134

Application For analysis Specified Titration

135

L-Glutamic acid monosodium salt monohydrate
pure

packaging	Mfr. No
250 g Plastic jar	G/0513/50
C ₅ H ₉ NO ₄ Na.H ₂ O	
CAS: 142-47-2	
MW: 187.14	
EINECS: 205-538-1	
Assay: >98%	

L-Glutamine
specified reagents for general laboratory work- extra pure

packaging	Mfr. No
250 g Plastic jar	G/0514/50
C ₅ H ₁₀ N ₂ O ₃	
CAS: 56-85-9	
MW: 146.15	
Assay: >99%	
Melting point: 185°C	
EINECS: 200-292-1	

Glutaraldehyde 50% solution
specified reagents for general laboratory work- extra pure

packaging	Mfr. No
500 mℓ Coated Glass bottle	G/0518/PB08
C ₅ H ₈ O ₂	
CAS: 111-30-8	
MW: 100.12	
Density at 20°C: 1.13	
Boiling point: >100°C	
Melting point: -15°C	
ADR/RID: 6.1, 8, II	
EINECS: 203-856-5	
UN: 2927	
GHS: Acute Tox. 3, Skin Corr. 1B, Eye Dam. 1, Resp. Sens. 1, Skin Sens. 1, STOT SE 3, Aquatic Acute 1	
H317, H335, H314, H334, H301, H331, H400	
P280, P305+P351+P338, P270, P260, P273, P308+P313	
Assay: 47 to 53%	

Glycerine See Glycerol, page 136

Glycerol
for electrophoresis

packaging	Mfr. No
100 mℓ Glass bottle	G/P450/05
500 mℓ Glass bottle	G/P450/08
HOCH ₂ CH(OH)CH ₂ OH	
CAS: 56-81-5	
MW: 92.09	
Density at 20°C: 1.261	
Boiling point: 290°C	
Melting point: 18°C	
Flash point: 160°C	
EINECS: 200-289-5	

For use in electrophoresis buffers and for stabilising enzyme preparations before freezing.

High Purity Acids for Elemental Analysis

Fisher Chemical offers a complete range of acids to meet the most challenging applications. Choose from our Optima®, TraceMetal® and PrimarPlus® grades to ensure that the sensitivity you require is always available.

Read more page 24

Glycerol **NEW PACK.**
for analysis Certified AR

packaging	Mfr. No
500 mℓ HDPE Plastic bottle	G/0650/08
1 ℓ HDPE Plastic bottle	G/0650/15
2.5 ℓ HDPE Plastic bottle	G/0650/17
5 ℓ Fluorinated HDPE Plastic bottle	G/0650/21
25 ℓ HDPE Plastic container	G/0650/25
200 ℓ HDPE Plastic container	G/0650/27
HOCH ₂ CH(OH)CH ₂ OH	
CAS: 56-81-5	
MW: 92.09	
Density at 20°C: 1.261	
Boiling point: 290°C	
Melting point: 18°C	
Flash point: 160°C	
EINECS: 200-289-5	
Assay: >99%	

Guaranteed analysis			
Colour	<10 APHA	Magnesium (Mg)	<2 ppm
Water	<0.5%	Potassium (K)	<5 ppm
pH (10% aq. soln.)	: 6 to 8	Sodium (Na)	<20 ppm
Residue after ignition	<0.005%	Total phosphorus (P)	<2 ppm
Identity IR test	: Pass	Total silicon (Si)	<2 ppm
Calcium (Ca)	<5 ppm	Total sulphur (S)	<5 ppm
Chloride (Cl)	<0.0005%	Zinc (Zn)	<2 ppm
Copper (Cu)	<0.5 ppm	Reducing substance	Pass
Iron (Fe)	<1 ppm	Esters of fatty acids	<0.08%
Lead (Pb)	<0.5 ppm	Fatty acids	Pass

Typical analysis			
Aluminium (Al)	<0.5 ppm	Lithium (Li)	<0.5 ppm
Arsenic (As)	<0.5 ppm	Manganese (Mn)	<0.5 ppm
Barium (Ba)	<0.5 ppm	Nickel (Ni)	<0.5 ppm
Cadmium (Cd)	<0.5 ppm	Titanium (Ti)	<0.5 ppm
Chromium (Cr)	<0.5 ppm	Vanadium (V)	<0.5 ppm
Cobalt (Co)	<0.5 ppm		

Glycerol
specified reagents for general laboratory work- extra pure

packaging	Mfr. No
500 mℓ HDPE Plastic bottle	G/0600/08
2.5 ℓ HDPE Plastic bottle	G/0600/17
25 ℓ HDPE Plastic container	G/0600/25
HOCH ₂ CH(OH)CH ₂ OH	
CAS: 56-81-5	
MW: 92.09	
Density at 20°C: 1.261	
Boiling point: 290°C	
Melting point: 18°C	
Flash point: 160°C	
EINECS: 200-289-5	
Assay (CG): >98%	

Guaranteed analysis			
pH (10% aq. soln.)	: 6 to 8	Magnesium (Mg)	<10 ppm
Residue after ignition	<0.01%	Potassium (K)	<20 ppm
Calcium (Ca)	<20 ppm	Sodium (Na)	<100 ppm
Chloride (Cl)	<0.001%	Total phosphorus (P)	<10 ppm
Copper (Cu)	<2 ppm	Total silicon (Si)	<10 ppm
Iron (Fe)	<5 ppm	Total sulphur (S)	<10 ppm
Lead (Pb)	<2 ppm	Zinc (Zn)	<10 ppm

Glycine, Buffer component
for electrophoresis and blotting

packaging	Mfr. No
500 g Plastic jar	G/P460/53
NH ₂ CH ₂ COOH	
CAS: 56-40-6	
MW: 75.07	
Density at 20°C: 1.595	
Melting point: 182 to 236°C	
EINECS: 200-272-2	
Assay: > 99.5%	

pH Buffers

NIST Certified Standard Solution & Concentrate

A full range to discover on page 63

Glycine
for analysis Certified AR

packaging	Mfr. No
100 g Plastic jar	G/0800/48
250 g Plastic jar	G/0800/50
1 kg Plastic jar	G/0800/60
NH ₂ CH ₂ COOH	
CAS: 56-40-6	
MW: 75.07	
Density at 20°C: 1.595	
Boiling point: 233°C	
Melting point: 182 to 236°C	
EINECS: 200-272-2	
Assay: >99%	

Guaranteed analysis			
pH (aqueous solution 5%)	5.5 to 6.5	Potassium (K)	<20 ppm
Ammonium (NH ₄)	<0.02%	Sodium (Na)	<50 ppm
Calcium (Ca)	<10 ppm	Total chloride (Cl)	<0.001%
Copper (Cu)	<2 ppm	Total phosphorus (P)	<10 ppm
Iron (Fe)	<5 ppm	Total silicon (Si)	<5 ppm
Lead (Pb)	<1 ppm	Total sulphur (S)	<50 ppm
Magnesium (Mg)	<5 ppm	Zinc (Zn)	<2 ppm

Typical analysis			
Aluminium (Al)	<2 ppm	Lithium (Li)	<1 ppm
Barium (Ba)	<2 ppm	Manganese (Mn)	<1 ppm
Cadmium (Cd)	<1 ppm	Nickel (Ni)	<1 ppm
Chromium (Cr)	<1 ppm	Titanium (Ti)	<1 ppm
Cobalt (Co)	<1 ppm	Vanadium (V)	<1 ppm

Lot analysis on the label.

Glycolic acid
specified reagents for general laboratory work- extra pure

packaging	Mfr. No
100 g Plastic jar	G/0803/48
C ₂ H ₃ O ₃	
CAS: 79-14-1	
MW: 76.05	
Density at 260°C: 1.36	
Boiling point: 100°C	
Melting point: 72 to 80°C	
Flash point: >300°C	
ADR/RID: 8, II	
EINECS: 201-180-5	
UN: 3261	
GHS: Acute Tox. 3, Skin Corr. 1B, Eye Dam. 1	
H302, H314, H331	
P301+P330+P331, P280, P305+P351+P338, P310, P304+P340, P301+P312	
Assay: >99%	
Guaranteed analysis	
Water <1%	

Gold standard metal solution 1000 ppm
for AAS, solution HCl 1M

packaging	Mfr. No
500 mℓ HDPE Plastic bottle	J/8028/08
GHS: Skin Irrit. 2	
H315	
P302+P352	

Manufactured to a tolerance of ± 0.5%.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.

Gold standard metal solution 1000 ppm
for ICP solution in HCL 10%, NIST certified

packaging	Mfr. No
100 mℓ Plastic bottle	SPPLAU32Y
500 mℓ Plastic bottle	SPPLAU32X
H315, H319, H335	

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Gold standard metal solution 1000 ppm
for ICP-MS solution in HCl 2%, NIST certified

packaging	Mfr. No
100 mℓ Plastic bottle	SPCLAU11Y
Manufactured to a tolerance of ± 0.5%.	
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.	
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.	
Check with your local distributor for availability of Spex product	

Gold metal
technical

packaging	Mfr. No
25 Plastic jar	G/0850/90
Au	
CAS: 7440-57-5	
MW: 196.97	
Boiling point: 2700°C	
Melting point: 1064.8°C	
EINECS: 231-165-9	

Guaranteed analysis	
Identity test	: Pass

Gram's iodine
pure

packaging	Mfr. No
250 mℓ Glass bottle	G/0910/07
CAS: 7553-56-2	
EINECS: 231-442-4	

Guaranteed analysis	
Identity test	: Pass

Graphite powder
pure

packaging	Mfr. No
1 kg Plastic jar	G/0900/60
CAS: 7782-42-5	
EINECS: 231-955-3	
GHS: Eye Irrit. 2, STOT SE 3	
H319, H335	
P261, P280, P305+P351+P338	

Guaranteed analysis	
Identity test	: Pass

Guaiacol
specified reagents for general laboratory work- extra pure

packaging	Mfr. No
100 g Glass bottle	G/0950/48
C ₇ H ₈ O ₂	
CAS: 90-05-1	
MW: 124.14	
Density at 20°C: 1.12	
Boiling point: 205°C	
Melting point: 27 to 29°C	
Flash point: 82°C	
EINECS: 201-964-7	
GHS: Acute Tox. 4, Skin Irrit. 2, Eye Irrit. 2	
H302, H315, H319	
P301+P312, P302+P352, P280, P305+P351+P338	

Guaranteed analysis	
Complies with standard IR spectrum: Pass	
Conforming to the IR spectrum of reference.	

UHPLC

Filtered to 0.1 µm, Fisher Chemical UHPLC gradient grade solvents (Acetonitrile, Methanol and Water) are perfect for those working at high pressure.

See page 309 for more information

Heptane fraction
specified reagents for general laboratory work- extra pure

SLR

packaging	Mfr. No
500 mℓ Glass bottle	H/0100/08
1 ℓ Glass bottle	H/0100/15
2.5 ℓ Glass bottle	H/0100/17
200 ℓ Metal drum	H/0100/27
C ₇ H ₁₆ CAS: 142-82-5 MW: 100.21 Density at 20°C: 0.683 Boiling point: 98°C Melting point: -91°C Flash point: -4°C ADR/RID: 3, II EINECS: 205-563-8	UN: 1206 GHS: Flam. Liq. 2, Skin Irrit. 2, STOT SE 3, Asp. Tox. 1, Aquatic Acute 1, Aquatic Chronic 1 H225, H315, H336, H304, H400, H410 P261, P301+P310, P331, P273, P302+P352, P210



1-Heptane sulfonic acid sodium salt
for lon pair chromatography for electrochemical detection

Certified

packaging	Mfr. No
25 g Plastic jar	H/0169/46
C ₇ H ₁₅ NaO ₃ S CAS: 22767-50-6 MW: 202.27 Melting point: >300°C EINECS: 245-210-5	GHS: Skin Irrit. 2, Eye Irrit. 2, STOT SE 3 H315, H335, H319 P261, P302+P352, P280, P305+P351+P338
Assay: >98%	
Guaranteed analysis	
Electrochemical test Pass	
Absorbance specification	
Max. Absorbance (1cm cell v HPLC water)	
200 nm (sol. aqu. 0.25 M) : <0.1 A.U.	240 nm (sol. aqu. 0.25 M) ... : <0.02 A.U.
210 nm (sol. aqu. 0.25 M) ... : <0.05 A.U.	250 nm (sol. aqu. 0.25 M) ... : <0.02 A.U.
220 nm (sol. aqu. 0.25 M) ... : <0.03 A.U.	260 nm (sol. aqu. 0.25 M) ... : <0.02 A.U.



Tested by HPLC-EC to show low electrochemical response.
Detector - ESA Coulochem II.
Mode-oxidative.
Electrode graphite.
Lot analysis on the label.

1-Heptane sulfonic acid sodium salt
for lon pair chromatography

Certified

packaging	Mfr. No
25 g Plastic jar	H/0168/46
100 g Plastic jar	H/0168/48
5 kg Plastic keg	H/0168/65
C ₇ H ₁₅ NaO ₃ S CAS: 22767-50-6 MW: 202.25 Melting point: >300°C EINECS: 245-210-5	GHS: Skin Irrit. 2, Eye Irrit. 2, STOT SE 3 H315, H335, H319 P261, P302+P352, P280, P305+P351+P338
Assay: >98%	
Absorbance specification	
Max. Absorbance (1cm cell v HPLC water)	
200 nm (sol. aqu. 0.25 M) : <0.1 A.U.	240 nm (sol. aqu. 0.25 M) ... : <0.02 A.U.
210 nm (sol. aqu. 0.25 M) ... : <0.05 A.U.	250 nm (sol. aqu. 0.25 M) ... : <0.02 A.U.
220 nm (sol. aqu. 0.25 M) ... : <0.03 A.U.	260 nm (sol. aqu. 0.25 M) ... : <0.02 A.U.



Lot analysis on the label.

1-Heptane sulfonic acid sodium salt solution 0.1M
for lon pair chromatography

Certified

packaging	Mfr. No
250 mℓ Glass bottle	H/0166/07
EINECS: 245-210-5	
Absorbance specification	
Max. Absorbance (1cm cell v HPLC water)	
254 nm (sol. aqu. 0.005 M) : <0.02 A.U.	
Lot analysis on the label.	

Heptan-1-ol
pure

packaging	Mfr. No
100 mℓ Glass bottle	H/0250/05
C ₇ H ₁₆ O CAS: 111-70-6 MW: 116.21 Density at 20°C: 0.822 Boiling point: 176°C Melting point: -35°C Flash point: 73°C ADR/RID: 6.1/III	EINECS: 203-897-9 UN: 2810 GHS: Acute Tox. 4, Eye Irrit. 2 H312, H319, H332 P302+P352, P301+P312, P304+P340, P280, P305+P351+P338
Assay (CG): >98%	



Hexachloroplatinic(IV) acid hydrate
specified reagents for general laboratory work- extra pure

SLR

packaging	Mfr. No
1 g Vial	H/0262/43
5 g Vial	H/0262/44
H ₂ Cl ₆ Pt.6H ₂ O CAS: 26023-84-7 MW: 517.92 Melting point: 60°C ADR/RID: 8, III EINECS: 241-010-7 UN: 2507 GHS: Acute Tox. 3, Skin	Corr. 1B, Eye Dam. 1, Resp. Sens. 1, Skin Sens. 1 H314, H301, H334, H317 P301+P310, P301+P330+P331, P305+P351+P338, P310, P261, P304+P341, P342+P311, P280



Assay: 39 to 41%

The price charged for this item is subject to alteration and will be that ruling at time of despatch.

Hexadecane
specified reagents for general laboratory work- extra pure

SLR

packaging	Mfr. No
100 mℓ Glass bottle	H/0270/05
C ₁₆ H ₃₄ CAS: 544-76-3 MW: 226.45 Density at 20°C: 0.773 Boiling point: 287°C Melting point: 18°C	Flash point: 135°C EINECS: 208-878-9 GHS: Asp. Tox. 1 H304, EUH066 P301+P310, P331, P280



Standard Volumetric solution
For your convenience we offer a range of NIST traceable solutions. A full list of our products can be found on page 323. Can't see what you need? Please email your contact details and a brief description of your requirements to begel.sha.quotes@thermofisher.com

1,1,1,3,3,3-Hexafluoropropan-2-ol
for gas chromatography

packaging	Mfr. No
25 mℓ Glass bottle	H/0285/04
C ₃ H ₂ F ₆ O CAS: 920-66-1 MW: 168.04 Density at 20°C: 1.596 Boiling point: 59°C Melting point: -4°C Flash point: >100°C ADR/RID: 8, II	EINECS: 213-059-4 UN: 3265 GHS: Acute Tox. 4, Skin Corr. 1B, Eye Dam. 1 H302, H314, H332 P280, P260, P264, P301+P330+P331, P305+P351+P338, P310
Assay (CG): >99%	



Hexafluorosilicic acid
See Fluorosilicic acid 35%, page 130

Hexamethylenediamine See 1, page 103

Hexamine
for analysis

Certified AR

packaging	Mfr. No
500 g Plastic jar	H/0320/53
1 kg Plastic jar	H/0320/60
5 kg Plastic keg	H/0320/65
C ₆ H ₁₂ N ₄ CAS: 100-97-0 MW: 140.19 Density at 20°C: 1.33 Melting point: 280°C Flash point: 250°C	ADR/RID: 4.1, III EINECS: 202-905-8 UN: 1328 GHS: Flam. Sol. 2, Skin Sens. 1 H228, H317 P280, P210, P233



Assay: >99.5%

pH (10% aq. soln.) : 8.5 to 9.5	Sodium (Na) <20 ppm
Ammonium (NH ₄) <0.005%	Total chloride (Cl) <0.002%
Calcium (Ca) <10 ppm	Total phosphorus (P) <20 ppm
Copper (Cu) <1 ppm	Total silicon (Si) <10 ppm
Iron (Fe) <2 ppm	Total sulphur (S) <50 ppm
Lead (Pb) <1 ppm	Zinc (Zn) <5 ppm
Magnesium (Mg) <5 ppm	

Typical analysis	
Aluminium (Al) <1 ppm	Lithium (Li) <1 ppm
Arsenic (As) <1 ppm	Manganese (Mn) <1 ppm
Barium (Ba) <1 ppm	Nickel (Ni) <1 ppm
Cadmium (Cd) <1 ppm	Titanium (Ti) <1 ppm
Chromium (Cr) <1 ppm	Vanadium (V) <1 ppm
Cobalt (Co) <1 ppm	

Lot analysis on the label.

Hexamine
specified reagents for general laboratory work- extra pure

SLR

packaging	Mfr. No
500 g Plastic jar	H/0300/53
1 kg Plastic jar	H/0300/60
5 kg Plastic keg	H/0300/65
25 kg Fibre Keg	H/0300/70
C ₆ H ₁₂ N ₄ CAS: 100-97-0 MW: 140.19 Density at 20°C: 1.33 Melting point: 280°C Flash point: 250°C	ADR/RID: 4.1, III EINECS: 202-905-8 UN: 1328 GHS: Flam. Sol. 2, Skin Sens. 1 H228, H317 P280, P210, P233



Assay: >99%

Guaranteed analysis		
Calcium (Ca) <50 ppm	n-Butyl formate <0.01%	
Copper (Cu) <5 ppm	Total phosphorus (P) <100 ppm	
Iron (Fe) <10 ppm	Total silicon (Si) <20 ppm	
Lead (Pb) <5 ppm	Total sulphur (S) <200 ppm	
Magnesium (Mg) <20 ppm	Zinc (Zn) <20 ppm	
Sodium (Na) <100 ppm		

n-Hexane
for analysis

Certified AR

packaging	Mfr. No
500 mℓ Coated Glass bottle	H/0421/PB08
2.5 ℓ Coated Glass bottle	H/0421/PB17
25 ℓ Metal drum	H/0421/25
200 ℓ Metal drum	H/0421/27
C ₆ H ₁₄ CAS: 110-54-3 MW: 86.18 Density at 20°C: 0.659 Boiling point: 69°C Melting point: -95°C Flash point: -22°C ADR/RID: 3, II EINECS: 203-777-6	UN: 1208 GHS: Flam. Liq. 2, Repr. 2, Skin Irrit. 2, STOT SE 3, STOT RE 2, Asp. Tox. 1, Aquatic Chronic 2 H225, H304, H336, H315, H373, H411, H361f P281, P260, P301+P310, P331, P273, P302+P352, P210




Assay (CG): >99%

Guaranteed analysis		
Acidity / Alkalinity <0.00008 meq/g	Lead (Pb) <0.02 ppm	
Colour <10 APHA	Magnesium (Mg) <0.05 ppm	
Water <0.02%	Potassium (K) <0.1 ppm	
Residue after evaporation <20 ppm	Sodium (Na) <0.1 ppm	
Substances darkened by (H ₂ SO ₄) Pass	Total phosphorus (P) <0.1 ppm	
Bromine number <0.5	Total silicon (Si) <0.05 ppm	
Calcium (Ca) <0.2 ppm	Total sulphur (S) <10 ppm	
Copper (Cu) <0.02 ppm	Zinc (Zn) <0.1 ppm	
Iron (Fe) <0.1 ppm	Aromatic Hydrocarbons <0.01%	

Typical analysis	
Aluminium (Al) <0.1 ppm	Manganese (Mn) <0.02 ppm
Barium (Ba) <0.02 ppm	Molybdenum (Mo) <0.02 ppm
Cadmium (Cd) <0.02 ppm	Nickel (Ni) <0.02 ppm
Chromium (Cr) <0.02 ppm	Strontium (Sr) <0.02 ppm
Cobalt (Co) <0.02 ppm	Titanium (Ti) <0.02 ppm
Lithium (Li) <0.02 ppm	Vanadium (V) <0.02 ppm

Lot analysis on the label.



n-Hexane **NEW PACK.** SLR
specified reagents for general laboratory work- extra pure

packaging		Mfr. No
500 mℓ	Coated Glass bottle	H/0420/PB08
2.5 ℓ	Coated Glass bottle	H/0420/PB17
5 ℓ	Fluorinated HDPE Plastic bottle	H/0420/21
25 ℓ	Metal drum	H/0420/25
C ₆ H ₁₄ CAS: 110-54-3 MW: 86.18 Density at 20°C: 0.659 Boiling point: 69°C Melting point: -95°C Flash point: -22°C ADR/RID: 3, II EINECS: 203-777-6		UN: 1208 GHS: Flam. Liq. 2, Repr. 2, Skin Irrit. 2, STOT SE 3, STOT RE 2, Asp. Tox. 1, Aquatic Chronic 2 H225, H304, H336, H315, H373, H411, H361f P281, P260, P301+P310, P331, P273, P302+P352, P210

Assay (CG): >99%			
Guaranteed analysis			
Acidity / Alkalinity	<0.0002 meq/g	Lead (Pb)	<1 ppm
Colour	<20 APHA	Magnesium (Mg)	<1 ppm
Water	<0.05%	Potassium (K)	<1 ppm
Residue after evaporation	<50 ppm	Sodium (Na)	<2 ppm
Bromine number	<0.5	Total phosphorus (P)	<1 ppm
Calcium (Ca)	<2 ppm	Total silicon (Si)	<1 ppm
Copper (Cu)	<1 ppm	Total sulphur (S)	<50 ppm
Iron (Fe)	<5 ppm	Zinc (Zn)	<2 ppm

Hexanes (95% n-Hexane approx.) **NEW PACK.** Certified
for HPLC

packaging		Mfr. No
1 ℓ	Glass bottle	H/0406/15
2.5 ℓ	Glass bottle	H/0406/17
2.5 ℓ	Coated Glass bottle	H/0406/PB17
5 ℓ	Aluminium can	H/0406/21
10 ℓ	Returnable Drum	H/0406/10RSS
25 ℓ	Returnable Drum	H/0406/25RSS
30 ℓ	Returnable Drum	H/0406/30RSS
C ₆ H ₁₄ CAS: 110-54-3 MW: 86.18 Density at 20°C: 0.659 Boiling point: 69°C Melting point: -95°C Flash point: -22°C ADR/RID: 3, II EINECS: 295-570-2 UN: 1208		GHS: Flam. Liq. 2, Repr. 2, Skin Irrit. 2, STOT SE 3, STOT RE 2, Asp. Tox. 1, Aquatic Chronic 2 H225, H304, H411, H315, H373, H336, H361f P280, P281, P302+P352, P261, P260, P301+P310, P331, P273, P210, P240

Guaranteed analysis			
Acidity / Alkalinity	<0.00015 meq/g	Residue after evaporation	<2 ppm
Water	<0.01%		

Absorbance specification			
Max. Absorbance (1cm cell v HPLC water)			
200 nm	<0.2 A.U.	230 nm	<0.01 A.U.
210 nm	<0.05 A.U.	240 nm	<0.005 A.U.
220 nm	<0.02 A.U.	250 nm	<0.005 A.U.

Lot analysis on the label.
Filtered to 0.2 µm.

UHPLC
Filtered to 0.1 µm, Fisher Chemical UHPLC
gradient grade solvents (Acetonitrile, Methanol
and Water) are perfect for those working at
high pressure.

See page 309 for more information



Hexanes (95% n-Hexane approx.) Certified
HPLC for fluorescence detection

packaging		Mfr. No
1 ℓ	Coated Glass bottle	H/0409/PB15
2.5 ℓ	Coated Glass bottle	H/0409/PB17
C ₆ H ₁₄ CAS: 92112-69-1 MW: 86.18 Density at 20°C: 0.659 Boiling point: 69°C Melting point: -95°C Flash point: -22°C ADR/RID: 3, II EINECS: 295-570-2 UN: 1208		GHS: Flam. Liq. 2, Repr. 2, Skin Irrit. 2, STOT SE 3, STOT RE 2, Asp. Tox. 1, Aquatic Chronic 2 H225, H304, H411, H315, H373, H336, H361f P280, P281, P302+P352, P261, P260, P301+P310, P331, P273, P210, P240

Assay (CG): >95%			
Guaranteed analysis			
Acidity / Alkalinity	<0.00015 meq/g	Residue after evaporation	<2 ppm
Water	<0.01%		

Absorbance specification			
Max. Absorbance (1cm cell v HPLC water)			
200 nm	<0.3 A.U.	240 nm	<0.005 A.U.
210 nm	<0.1 A.U.	250 nm	<0.005 A.U.
220 nm	<0.025 A.U.	Fluorescence	Pass
230 nm	<0.01 A.U.		

Fluorescence
Excitation in the range 250-750nm produces no fluorescence emission greater than
the solvents natural Raman signal at 350nm excitation.This is equivalent to the
emission maximum of approximately 1ppb quinine sulfate.
Lot analysis on the label (excluding fluorescence scan).
Filtered to 0.2 µm.

Hexanes (95% n-Hexane approx.) Distol
for residue and pesticide analysis

packaging		Mfr. No
1 ℓ	Glass bottle	H/0403/15
2.5 ℓ	Glass bottle	H/0403/17
C ₆ H ₁₄ CAS: 92112-69-1 MW: 86.18 Density at 20°C: 0.659 Boiling point: 69°C Melting point: -95°C Flash point: -22°C ADR/RID: 3, II EINECS: 295-570-2 UN: 1208		GHS: Flam. Liq. 2, Repr. 2, Skin Irrit. 2, STOT SE 3, STOT RE 2, Asp. Tox. 1, Aquatic Chronic 2 H225, H304, H411, H315, H373, H336, H361f P280, P281, P302+P352, P261, P260, P301+P310, P331, P273, P210, P240

Guaranteed analysis			
Water	<0.01%	GC (FID) Analysis	Pass
Residue after evaporation	<2 ppm	GC (NPD) analysis	Pass
GC (ECD) Analysis	Pass		

GC Analysis
Analysis of pesticides: from the retention time of Lindane, after concentration of
2500 times, no impurity peak higher than that obtained with 10 pg / mL heptachlo-
réoxyde (ECD) and 5 pg / mL Parathion (NPD).
Analysis of hydrocarbons: no impurity chain longer than the dodecane is present at
a concentration greater than 0.01 ppm (FID).
Conditions :
Column: J&W DB 5: 30m x 0.32mm x 0.1µm df.
Temperature profile: (1) 50°C, 2 min
(2) 40°C / min until 320°C
(3) 320°C for 30 min
Injection: 5µL cold on the column.

Hexanes (95% n-Hexane approx.)
extra dry, for synthesis

packaging		Mfr. No
1 ℓ	Coated Glass bottle	H/0402/PB15
2.5 ℓ	Coated Glass bottle	H/0402/PB17
C ₆ H ₁₄ CAS: 92112-69-1 MW: 86.18 Density at 20°C: 0.659 Boiling point: 69°C Melting point: -95°C Flash point: -22°C ADR/RID: 3, II EINECS: 295-570-2 UN: 1208		GHS: Flam. Liq. 2, Repr. 2, Skin Irrit. 2, STOT SE 3, STOT RE 2, Asp. Tox. 1, Aquatic Chronic 2 H225, H304, H411, H315, H373, H336, H361f P280, P281, P302+P352, P261, P260, P301+P310, P331, P273, P210, P240

Guaranteed analysis			
Water	<0.005%	Residue after evaporation	<2 ppm

Hexanes (95% n-Hexane approx) **NEW PACK.** Certified AR
for analysis

packaging		Mfr. No
1 ℓ	Glass bottle	H/0355/15
2.5 ℓ	Glass bottle	H/0355/17
5 ℓ	Fluorinated HDPE Plastic bottle	H/0355/21
25 ℓ	Metal drum	H/0355/25
200 ℓ	Metal drum	H/0355/27
C ₆ H ₁₄ CAS: 92112-69-1 MW: 86.18 Density at 20°C: 0.659 Boiling point: 69°C Melting point: -95°C Flash point: -22°C ADR/RID: 3, II EINECS: 295-570-2 UN: 1208		GHS: Flam. Liq. 2, Repr. 2, Skin Irrit. 2, STOT SE 3, STOT RE 2, Asp. Tox. 1, Aquatic Chronic 2 H225, H304, H411, H315, H373, H336, H361f P280, P281, P302+P352, P261, P260, P301+P310, P331, P273, P210, P240

Assay (CG): >95%			
Guaranteed analysis			
Acidity / Alkalinity	<0.00008 meq/g	Iron (Fe)	<0.1 ppm
Colour	<10 APHA	Lead (Pb)	<0.02 ppm
Water	<0.02%	Magnesium (Mg)	<0.05 ppm
Residue after evaporation	<20 ppm	Potassium (K)	<0.1 ppm
Substances darkened by (H ₂ SO ₄)	Pass	Sodium (Na)	<0.1 ppm
Bromine number	<0.5	Total phosphorus (P)	<0.1 ppm
Hexane Fractions	: >98.5%	Total silicon (Si)	<0.05 ppm
Thiophene test	: Pass	Total sulphur (S)	<10 ppm
Calcium (Ca)	<0.2 ppm	Zinc (Zn)	<0.1 ppm
Copper (Cu)	<0.02 ppm	Aromatic Hydrocarbons	<0.1%

Typical analysis			
Aluminium (Al)	<0.1 ppm	Manganese (Mn)	<0.02 ppm
Barium (Ba)	<0.02 ppm	Molybdenum (Mo)	<0.02 ppm
Cadmium (Cd)	<0.02 ppm	Nickel (Ni)	<0.02 ppm
Chromium (Cr)	<0.02 ppm	Strontium (Sr)	<0.02 ppm
Cobalt (Co)	<0.02 ppm	Titanium (Ti)	<0.02 ppm
Lithium (Li)	<0.02 ppm	Vanadium (V)	<0.02 ppm

Lot analysis on the label.

Optima LC/MS

Fisher Chemical OPTIMA® LC/MS solvents,
blends, additives and modifiers have set the
standard of excellence for consistent,
reproducible performance in the mobile
phase of LC/MS.

Find out more on page 308



Hexanes fraction **NEW PACK.** SLR
specified reagents for general laboratory work- extra pure

packaging		Mfr. No
500 mℓ	Glass bottle	H/0350/08
1 ℓ	Glass bottle	H/0350/15
2.5 ℓ	Glass bottle	H/0350/17
2.5 ℓ	Coated Glass bottle	H/0350/PB17
5 ℓ	Fluorinated HDPE Plastic bottle	H/0350/21
5 ℓ	Fluorinated HDPE square Drum	H/0350/FP21
25 ℓ	Metal drum	H/0350/25
200 ℓ	Metal drum	H/0350/27

C ₆ H ₁₄ CAS: 92112-69-1 MW: 86.18 Density at 20°C: 0.659 Boiling point: 69°C Melting point: -95°C Flash point: -22°C ADR/RID: 3, II EINECS: 295-570-2 UN: 1208		GHS: Flam. Liq. 2, Repr. 2, Skin Irrit. 2, STOT SE 3, STOT RE 2, Asp. Tox. 1, Aquatic Chronic 2 H225, H304, H411, H315, H373, H336, H361f P280, P281, P302+P352, P261, P260, P301+P310, P331, P273, P210, P240
--	--	--

Guaranteed analysis	
Residue after evaporation	<20 ppm

Hexanes fraction SLR
extra pure, low in aromatic hydrocarbons

packaging		Mfr. No
500 mℓ	Glass bottle	H/0400/08
2.5 ℓ	Glass bottle	H/0400/17

C ₆ H ₁₄ CAS: 92112-69-1 MW: 86.18 Density at 20°C: 0.659 Boiling point: 69°C Melting point: -95°C Flash point: -22°C ADR/RID: 3, II EINECS: 295-570-2 UN: 1208		GHS: Flam. Liq. 2, Repr. 2, Skin Irrit. 2, STOT SE 3, STOT RE 2, Asp. Tox. 1, Aquatic Chronic 2 H225, H304, H411, H315, H373, H336, H361f P280, P281, P302+P352, P261, P260, P301+P310, P331, P273, P210, P240
--	--	--

Guaranteed analysis			
Residue after evaporation	<20 ppm	Aromatic Hydrocarbons	<0.1%

Hexanedioic acid See Adipic acid, page 26

Label

Fisher Chemical labels show GHS pictograms clearly.
Together with actual lot analysis shown on the label,
plus a wealth of additional information, we help you
to work safely to achieve accurate and consistent
results in your analytical applications.



1-Hexane sulfonic acid sodium salt
for Ion pair chromatography for electrochemical detection

Certified

Packaging		Mfr. No
25 g	Plastic jar	H/0439/46
C ₆ H ₁₃ NaO ₃ S.xH ₂ O	GHS: Skin Irrit. 2, Eye Irrit. 2, STOT SE 3	
CAS: 2832-45-3	H335, H315, H319	
MW: 188.22	P261, P302+P352, P280, P305+P351+P338	
Melting point: >300°C		
EINECS: 220-601-3		
Assay: >98%		
Guaranteed analysis		
Electrochemical test		Pass
Absorbance specification		
Max. Absorbance (1cm cell v HPLC water)		
200 nm (sol. aqu. 0.25 M)	<0.1 A.U.	240 nm (sol. aqu. 0.25 M) ... : <0.02 A.U.
210 nm (sol. aqu. 0.25 M) ...	<0.05 A.U.	250 nm (sol. aqu. 0.25 M) ... : <0.02 A.U.
220 nm (sol. aqu. 0.25 M) ...	<0.03 A.U.	260 nm (sol. aqu. 0.25 M) ... : <0.02 A.U.

Tested by HPLC-EC to show low electrochemical response.
Detector - ESA Coulochem II
Mode-oxidative
Graphite Electrode.
Lot analysis on the label.

1-Hexane sulfonic acid sodium salt
for Ion pair chromatography

Certified

packaging		Mfr. No
25 g	Plastic jar	H/0438/46
100 g	Plastic jar	H/0438/48
C ₆ H ₁₃ NaO ₃ S.xH ₂ O	GHS: Skin Irrit. 2, Eye Irrit. 2, STOT SE 3	
CAS: 2832-45-3	H335, H315, H319	
MW: 188.22	P261, P302+P352, P280, P305+P351+P338	
Melting point: >300°C		
EINECS: 220-601-3		
Assay: >98%		
Absorbance specification		
Max. Absorbance (1cm cell v HPLC water)		
200 nm (sol. aqu. 0.25 M) : <0.1 A.U.	240 nm (sol. aqu. 0.25 M) ... : <0.02 A.U.	
210 nm (sol. aqu. 0.25 M) ... : <0.05 A.U.	250 nm (sol. aqu. 0.25 M) ... : <0.02 A.U.	
220 nm (sol. aqu. 0.25 M) ... : <0.03 A.U.	260 nm (sol. aqu. 0.25 M) ... : <0.02 A.U.	
Lot analysis on the label.		

Lot analysis on the label.

n-Hexanoic acid

SLR

packaging		Mfr. No
250 mL	Glass bottle	H/0440/07
C ₆ H ₁₂ O ₂	UN: 2829	
CAS: 142-62-1	GHS: Acute Tox. 4, Skin Corr. 1B, Eye Dam. 1	
MW: 116.16	H314, H332, H302, H312	
Density at 20°C: 0.920	P280, P301+P330+P331, P305+P351+P338, P310, P302+P352, P301+P312, P304+P340	
Boiling point: 202 to 203°C		
Melting point: -3°C		
Flash point: 104°C		
ADR/RID: 8, III		
EINECS: 205-550-7		

Guaranteed analysis
Complies with standard IR spectrum: Pass

Reagents, Solvents & Volumetric solutions comply with the European Pharmacopeia

Read more on page 234

Hexan-1-ol

SLR

packaging		Mfr. No
500 mL	Coated Glass bottle	H/0450/PB08
2.5 L	Coated Glass bottle	H/0450/PB17
25 L	Metal drum	H/0450/25
200 L	Metal drum	H/0450/27
C ₆ H ₁₄ O	ADR/RID: 3, III	
CAS: 111-27-3	EINECS: 203-852-3	
MW: 102.18	UN: 2282	
Density at 20°C: 0.814	GHS: Acute Tox. 4	
Boiling point: 156 to 157°C	H302	
Melting point: -52°C	P280, P301+P312, P210	
Flash point: 59°C		
Assay (CG): >99%		

Hex-1-ene

SLR

Packaging		Mfr. No
100 mL	Glass bottle	H/0460/05
C ₆ H ₁₂	EINECS: 209-753-1	
CAS: 592-41-6	UN: 2370	
MW: 84.16	GHS: Flam. Liq. 2, STOT SE	
Density at 20°C: 0.678	3, Asp. Tox. 1, Aquatic Chronic 2	
Boiling point: 62 to 65°C	H225, H304, H411, H336,	
Melting point: -140°C	EUH066	
Flash point: -26°C	P261, P301+P310, P331,	
ADR/RID: 3, II	P273, P210, P240	
		

L-Histidine

pure

Packaging		Mfr. No
100 g	Plastic jar	H/0467/48
3 kg	Plastic jar	H/0467/63
<chem>C6H9N3O2</chem>		Melting point: 282°C EINECS: 200-745-3
CAS: 71-00-1		
MW: 155.16		
Assay: >98%		

Holmium standard metal solution 1000 ppm

for ICP solution in HNO₃ 2%, NIST certified

packaging		Mfr. No
100 mL	Plastic bottle	SPPLHO22Y
500 mL	Plastic bottle	SPPLHO22X
H315, H319		

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Hyamine® 1622

For microscopy, solution 0.004M

packaging		Mfr. No
2.5 L	Glass bottle	H/0480/17
Guaranteed analysis		
Factor : 0.97 to 1.03		

Hyamine solution 5 % w/v

SLR

packaging		Mfr. No
1 L	Glass bottle	J/9016/15

Hyamine solution 0.05 M

SLR

packaging		Mfr. No
2.5 L	Glass bottle	J/9019/17

Hyamine solution 0.004 M

SLR

packaging		Mfr. No
2.5 L	Glass bottle	J/9018/17

Hydrazine hydrate

for analysis

Certified AR

packaging		Mfr. No
100 mL	Coated Glass bottle	H/0560/PB05
500 mL	HDPE Plastic bottle	H/0560/PB08
N ₂ H ₄ .H ₂ O	Skin Sens. 1, Aquatic Acute 1, Aquatic Chronic 1	
CAS: 7803-57-8	H311, H350, H301, H314, H331, H317, H410	
MW: 50.06	P301+P310, P280, P312, P302+P350,	
Boiling point: 118°C	P301+P330+P331,	
Melting point: -51.7°C	P305+P351+P338, P310,	
Flash point: 75°C	P201, P308+P313,	
ADR/RID: 8, 6.1, II	P304+P340, P273	
UN: 2030		
GHS: Carc. 1B, Acute Tox. 3, Skin Corr. 1B, Eye Dam. 1,		

Assay: >99%		
Guaranteed analysis		
Colour	<10 APHA	Sodium (Na) <100 ppm
Calcium (Ca)	<10 ppm	Total chloride (Cl) <0.005%
Copper (Cu)	<1 ppm	Total phosphorus (P) <20 ppm
Iron (Fe)	<2 ppm	Total silicon (Si) <50 ppm
Lead (Pb)	<1 ppm	Total sulphur (S) <50 ppm
Magnesium (Mg)	<5 ppm	Zinc (Zn) <5 ppm
Potassium (K)	<5 ppm	
Typical analysis		
Aluminium (Al)	<1 ppm	Manganese (Mn) <1 ppm
Barium (Ba)	<1 ppm	Molybdenum (Mo) <1 ppm
Cadmium (Cd)	<1 ppm	Nickel (Ni) <1 ppm
Chromium (Cr)	<1 ppm	Selenium (Se) <1 ppm
Cobalt (Co)	<1 ppm	Strontium (Sr) <1 ppm
Tin (Sn)	<1 ppm	Titanium (Ti) <1 ppm
Lithium (Li)	<1 ppm	Vanadium (V) <1 ppm

Hydrazine hydrate

SLR

packaging		Mfr. No
250 mL	Coated Glass bottle	H/0550/PB07
2.5 L	HDPE Plastic bottle	H/0550/PB17

N ₂ H ₄ .H ₂ O	Skin Sens. 1, Aquatic Acute 1, Aquatic Chronic 1	
CAS: 7803-57-8	H311, H350, H301, H314, H331, H317, H410	
MW: 50.06	P301+P310, P280, P312, P302+P350,	
Boiling point: 118°C	P301+P330+P331,	
Melting point: -51.7°C	P305+P351+P338, P310,	
Flash point: 75°C	P201, P308+P313,	
ADR/RID: 8, 6.1, II	P304+P340, P273	
UN: 2030		
GHS: Carc. 1B, Acute Tox. 3, Skin Corr. 1B, Eye Dam. 1,		

Assay: >99%		
Guaranteed analysis		
Colour	<20 APHA	Sodium (Na) <500 ppm
Calcium (Ca)	<50 ppm	n-Butyl formate <0.02%
Copper (Cu)	<5 ppm	Total phosphorus (P) <100 ppm
Iron (Fe)	<10 ppm	Total silicon (Si) <200 ppm
Lead (Pb)	<5 ppm	Total sulphur (S) <200 ppm
Magnesium (Mg)	<20 ppm	Zinc (Zn) <20 ppm
Potassium (K)	<20 ppm	

Hydrazine hydrate solution 60 %

SLR

packaging		Mfr. No
500 mL	HDPE Plastic bottle	H/0500/PB08
N ₂ H ₄ .H ₂ O	Skin Sens. 1, Aquatic Acute 1, Aquatic Chronic 1	
CAS: 7803-57-8	H311, H350, H301, H314, H331, H317, H410	
MW: 50.06	P301+P310, P280, P312, P302+P350,	
Boiling point: 118°C	P301+P330+P331,	
Melting point: -51.7°C	P305+P351+P338, P310,	
Flash point: 75°C	P201, P308+P313,	
ADR/RID: 8, 6.1, II	P304+P340, P273	
UN: 2030		
GHS: Carc. 1B, Acute Tox. 3, Skin Corr. 1B, Eye Dam. 1,		


Assay: 58 to 64%	
Guaranteed analysis	
n-Butyl formate <0.03%	

Hydrazinium chloride

SLR

packaging		Mfr. No
250 g	Plastic jar	H/0600/50
NH ₂ NH ₂ .2HCl	GHS: Carc. 1B, Acute Tox. 3,	
CAS: 5341-61-7	Skin Sens. 1, Aquatic Acute	
MW: 104.97	1, Aquatic Chronic 1	
Density at 20°C: 1.42	H311, H331, H350, H410,	
Boiling point: 200°C	H301, H317	
Melting point: 198°C	P301+P310, P280,	
ADR/RID: 6.1, III	P302+P350, P201,	
EINECS: 226-283-2	P308+P313, P304+P340, P273	
UN: 2811		

Assay: >98%



Hydrazinium sulfate		Certified AR	
for analysis			
packaging		Mfr. No	
100 g	Plastic jar	H/0700/48	
250 g	Plastic jar	H/0700/50	
H ₂ N ₂ ·H ₂ SO ₄ CAS: 10034-93-2 MW: 130.12 Melting point: 254°C ADR/RID: 6.1, III EINECS: 233-110-4 UN: 3288 GHS: Carc. 1B, Acute Tox. 3,	Skin Sens. 1, Aquatic Acute 1, Aquatic Chronic 1 H350, H331, H311, H301, H317, H410 P280, P301+P310, P302+P352, P304+P340, P308+P313, P273		
Assay: >99%			
Guaranteed analysis			
Calcium (Ca)	<10 ppm	Potassium (K)	<10 ppm
Copper (Cu)	<5 ppm	Sodium (Na)	<20 ppm
Iron (Fe)	<10 ppm	Total chloride (Cl)	<0.002%
Lead (Pb)	<5 ppm	Total silicon (Si)	<20 ppm
Magnesium (Mg)	<5 ppm	Zinc (Zn)	<5 ppm
Typical analysis			
Manganese (Mn)	<1 ppm	Cobalt (Co)	<1 ppm
Aluminium (Al)	<1 ppm	Lithium (Li)	<1 ppm
Barium (Ba)	<1 ppm	Nickel (Ni)	<1 ppm
Cadmium (Cd)	<1 ppm	Titanium (Ti)	<1 ppm
Chromium (Cr)	<1 ppm	Vanadium (V)	<1 ppm

Lot analysis on the label.

Hydrazinium sulfate		SLR	
specified reagents for general laboratory work- extra pure			
packaging		Mfr. No	
500 g	Plastic jar	H/0650/53	
25 kg	Plastic keg	H/0650/70	
H ₂ N ₂ ·H ₂ SO ₄ CAS: 10034-93-2 MW: 130.12 Melting point: 254°C ADR/RID: 6.1, III EINECS: 233-110-4 UN: 3288 GHS: Carc. 1B, Acute Tox. 3,	Skin Sens. 1, Aquatic Acute 1, Aquatic Chronic 1 H350, H331, H311, H301, H317, H410 P280, P301+P310, P302+P352, P304+P340, P308+P313, P273		
Assay: >99%			
Guaranteed analysis			
Calcium (Ca)	<50 ppm	Potassium (K)	<20 ppm
Copper (Cu)	<20 ppm	Sodium (Na)	<100 ppm
Iron (Fe)	<50 ppm	n-Butyl formate	<0.01%
Lead (Pb)	<20 ppm	Total silicon (Si)	<100 ppm
Magnesium (Mg)	<20 ppm	Zinc (Zn)	<20 ppm

Hydriodic acid 64-68%		Certified AR	
for analysis			
packaging		Mfr. No	
250 mL	Glass bottle	H/0850/07	
HI CAS: 10034-85-2 MW: 127.91 Density at 20°C: 1.7 Boiling point: 127°C Melting point: -51°C ADR/RID: 8, II EINECS: 233-109-9	UN: 1787 GHS: Skin Corr. 1A, Eye Dam. 1, STOT SE 3 H335, H314 P264, P280, P301+P330+P331, P305+P351+P338, P310, P260		
Assay: 64 to 68%			
Guaranteed analysis			
Calcium (Ca)	<10 ppm	Nickel (Ni)	<0.5 ppm
Chloride and Bromide	<0.02%	Potassium (K)	<10 ppm
Chromium (Cr)	<0.5 ppm	Sodium (Na)	<10 ppm
Copper (Cu)	<1 ppm	Total silicon (Si)	<20 ppm
Iron (Fe)	<5 ppm	Total sulphur (S)	<50 ppm
Lead (Pb)	<1 ppm	Zinc (Zn)	<2 ppm
Magnesium (Mg)	<5 ppm		

Lot analysis on the label.

Hydriodic acid 54-56%		SLR	
extra pure, 1.7 specific gravity			
packaging		Mfr. No	
100 mL	Glass bottle	H/0750/05	
HI CAS: 10034-85-2 MW: 127.91 Density at 20°C: 1.7 Boiling point: 127°C Melting point: -51°C ADR/RID: 8, II EINECS: 233-109-9	UN: 1787 GHS: Skin Corr. 1A, Eye Dam. 1, STOT SE 3 H335, H314 P264, P280, P301+P330+P331, P305+P351+P338, P310, P260		
Assay: 54 to 56%			
Guaranteed analysis			
Calcium (Ca)	<50 ppm	Nickel (Ni)	<2 ppm
Chloride and Bromide	<0.1%	Potassium (K)	<50 ppm
Chromium (Cr)	<2 ppm	Sodium (Na)	<50 ppm
Copper (Cu)	<5 ppm	Total silicon (Si)	<100 ppm
Iron (Fe)	<10 ppm	Total sulphur (S)	<200 ppm
Lead (Pb)	<5 ppm	Zinc (Zn)	<10 ppm
Magnesium (Mg)	<20 ppm		

Hydrobromic acid 44-49%		NEW		Optima	
for ultra trace elemental analysis					
packaging		Mfr. No			
500 mL	Pre-cleaned Teflon FEP bottle	A471-500			
HBr CAS: 10035-10-6 MW: 80.92 Density at 20°C: 1.48 Boiling point: 128°C Melting point: -11°C ADR/RID: 8, II EINECS: 233-113-0	UN: 1788 GHS: Skin Corr. 1B, Eye Dam. 1, STOT SE 3 H314, H335 P301+P330+P331, P280, P305+P351+P338, P261, P304+P340				
Assay: 44-49%					
Guaranteed analysis					
Zirconium (Zr)	<10 ppt	Molybdenum (Mo)	<100 ppt		
Gallium (Ga)	<100 ppt	Neodymium (Nd)	<10 ppt		
Terbium (Tb)	<10 ppt	Nickel (Ni)	<100 ppt		
Indium (In)	<10 ppt	Niobium (Nb)	<10 ppt		
Dysprosium (Dy)	<10 ppt	Palladium (Pd)	<100 ppt		
Aluminium (Al)	<100 ppt	Potassium (K)	<100 ppt		
Barium (Ba)	<100 ppt	Praseodymium (Pr)	<10 ppt		
Beryllium (Be)	<100 ppt	Rhenium (Re)	<100 ppt		
Bismuth (Bi)	<100 ppt	Rhodium (Rh)	<100 ppt		
Cadmium (Cd)	<100 ppt	Ruthenium (Ru)	<100 ppt		
Calcium (Ca)	<100 ppt	Samarium (Sm)	<10 ppt		
Cerium (Ce)	<10 ppt	Scandium (Sc)	<100 ppt		
Cesium (Cs)	<10 ppt	Propan--1-ol	<100 ppt		
Chromium (Cr)	<100 ppt	Sodium (Na)	<100 ppt		
Cobalt (Co)	<100 ppt	Strontium (Sr)	<100 ppt		
Copper (Cu)	<100 ppt	Tellurium (Te)	<10 ppt		
Erbium (Er)	<10 ppt	Thallium (Tl)	<10 ppt		
Europium (Eu)	<10 ppt	Thorium (Th)	<1 ppt		
Gadolinium (Gd)	<10 ppt	Thulium (Tm)	<10 ppt		
Gold (Au)	<100 ppt	Tin (Sn)	<100 ppt		
Hafnium (Hf)	<10 ppt	Titanium (Ti)	<100 ppt		
Holmium (Ho)	<10 ppt	Uranium (U)	<1 ppt		
Iron (Fe)	<100 ppt	Vanadium (V)	<100 ppt		
Lanthanum (La)	<10 ppt	Ytterbium (Yb)	<10 ppt		
Lead (Pb)	<1 ppt	Yttrium (Y)	<10 ppt		
Lithium (Li)	<100 ppt	Zinc (Zn)	<100 ppt		
Lutetium (Lu)	<10 ppt	Platinum (Pt)	<100 ppt		
Magnesium (Mg)	<100 ppt	Rubidium (Rb)	<100 ppt		
Manganese (Mn)	<100 ppt	Tungsten (W)	<100 ppt		

Bottled in class 10 cleanroom conditions. Individually double-bagged in class 100 cleanroom conditions. Each bottle is individually packaged. Certificate of analysis delivered with each bottle.

Hydrobromic acid 48%		Certified AR	
for analysis			
packaging		Mfr. No	
500 mL	HDPE Plastic bottle	H/0950/PB08	
2.5 L	HDPE Plastic bottle	H/0950/PB17	
HBr CAS: 10035-10-6 MW: 80.92 Density at 20°C: 1.48 Boiling point: 126 to 128°C Melting point: -11°C ADR/RID: 8, II EINECS: 233-113-0	UN: 1788 GHS: Skin Corr. 1B, Eye Dam. 1, STOT SE 3 H314, H335 P301+P330+P331, P280, P305+P351+P338, P261, P304+P340		
Assay: 47 to 49%			
Guaranteed analysis			
Calcium (Ca)	<2 ppm	Potassium (K)	<0.5 ppm
Chromium (Cr)	<0.5 ppm	Sodium (Na)	<10 ppm
Copper (Cu)	<0.1 ppm	n-Butyl formate	<0.05%
Iodide (I)	<0.0002%	Total phosphorus (P)	<2 ppm
Iron (Fe)	<0.5 ppm	Total silicon (Si)	<5 ppm
Lead (Pb)	<0.1 ppm	Total sulphur (S)	<20 ppm
Magnesium (Mg)	<0.5 ppm	Zinc (Zn)	<0.5 ppm
Nickel (Ni)	<0.5 ppm		

Hydrobromic acid 48%		SLR	
specified reagents for general laboratory work- extra pure			
packaging		Mfr. No	
500 mL	HDPE Plastic bottle	H/0900/PB08	
2.5 L	HDPE Plastic bottle	H/0900/PB17	
25 L	HDPE Plastic container	H/0950/25	
200 L	HDPE Plastic container	H/0900/27	
HBr CAS: 10034-85-2 MW: 80.92 Density at 20°C: 1.48 Boiling point: 126 to 128°C Melting point: -11°C ADR/RID: 8, II EINECS: 233-113-0	UN: 1788 GHS: Skin Corr. 1B, Eye Dam. 1, STOT SE 3 H314, H335 P301+P330+P331, P280, P305+P351+P338, P261, P304+P340		
Assay: 47 to 49%			
Guaranteed analysis			
Calcium (Ca)	<20 ppm	Potassium (K)	<10 ppm
Chromium (Cr)	<2 ppm	Sodium (Na)	<50 ppm
Copper (Cu)	<2 ppm	n-Butyl formate	<0.2%
Iodide (I)	<0.001%	Total phosphorus (P)	<20 ppm
Iron (Fe)	<5 ppm	Total silicon (Si)	<20 ppm
Lead (Pb)	<2 ppm	Total sulphur (S)	<100 ppm
Magnesium (Mg)	<10 ppm	Zinc (Zn)	<5 ppm
Nickel (Ni)	<2 ppm		

Hydrochloric acid 32-35%		NEW		Optima	
for ultra trace elemental analysis					
packaging		Mfr. No			
250 mL	Teflon PFA bottle	A466-250			
500 mL	Pre-cleaned Teflon PFA bottle	A466-500			
1 L	Pre-cleaned Teflon PFA bottle	A466-1			
HCl CAS: 7647-01-0 MW: 36.46 g/mol Density at 20°C: 1.18 Boiling point: 57°C Melting point: -35°C	ADR/RID: 8, II EINECS: 231-595-7 UN: 1789 GHS: Skin Corr. 1B, STOT SE 3 H335, H314 P305+P351+P338, P310, P280				
Assay: 32-35%					
Guaranteed analysis					
Zirconium (Zr)	<10 ppt	Arsenic (As)	<50 ppt		
Gallium (Ga)	<10 ppt	Barium (Ba)	<10 ppt		
Terbium (Tb)	<1 ppt	Beryllium (Be)	<10 ppt		
Indium (In)	<1 ppt	Bismuth (Bi)	<10 ppt		
Dysprosium (Dy)	<1 ppt	Boron (B)	<100 ppt		
Aluminium (Al)	<20 ppt	Cadmium (Cd)	<10 ppt		
Antimony (Sb)	<20 ppt				

Calcium (Ca)	<10 ppt	Niobium (Nb)	<1 ppt
Cerium (Ce)	<10 ppt	Potassium (K)	<10 ppt
Cesium (Cs)	<10 ppt	Praseodymium (Pr)	<1 ppt
Chromium (Cr)	<10 ppt	Rhenium (Re)	<10 ppt
Cobalt (Co)	<10 ppt	Rhodium (Rh)	<10 ppt
Copper (Cu)	<10 ppt	Ruthenium (Ru)	<10 ppt
Erbium (Er)	<1 ppt	Samarium (Sm)	<1 ppt
Europium (Eu)	<1 ppt	Scandium (Sc)	<10 ppt
Gadolinium (Gd)	<1 ppt	Propan--1-ol	<10 ppt
Gold (Au)	<20 ppt	Sodium (Na)	<10 ppt
Hafnium (Hf)	<10 ppt	Strontium (Sr)	<10 ppt
Holmium (Ho)	<1 ppt	Tellurium (Te)	<1 ppt
Iron (Fe)	<10 ppt	Thallium (Tl)	<10 ppt
Lanthanum (La)	<1 ppt	Thorium (Th)	<1 ppt
Lead (Pb)	<10 ppt	Thulium (Tm)	<1 ppt
Lithium (Li)	<10 ppt	Tin (Sn)	<20 ppt
Lutetium (Lu)	<10 ppt	Titanium (Ti)	<10 ppt
Magnesium (Mg)	<10 ppt	Uranium (U)	<1 ppt
Manganese (Mn)	<10 ppt	Vanadium (V)	<10 ppt
Mercury (Hg)	<50 ppt	Ytterbium (Yb)	<1 ppt
Molybdenum (Mo)	<10 ppt	Yttrium (Y)	<1 ppt
Neodymium (Nd)	<1 ppt	Zinc (Zn)	<10 ppt
Nickel (Ni)	<20 ppt	Rubidium (Rb)	<10 ppt

Bottled in class 10 cleanroom conditions. Individually double-bagged in class 100 cleanroom conditions. Each bottle is individually packaged. Certificate of analysis delivered with each bottle.

Hydrochloric acid 34-37%			NEW	Trace Metal
for trace metal analysis				
packaging			Mfr. No	
500 mL	HDPE plastic bottle		A508-P500	
1 L	HDPE plastic bottle		A508-P1	
2.5 L	HDPE plastic bottle		A508-P212	
HCl		ADR/RID: 8, II		
CAS: 7647-01-0		EINECS: 231-595-7		
MW: 36.46		UN: 1789		
Density at 20°C: 1.18		GHS: Skin Corr. 1B, STOT SE 3		
Boiling point: 57°C		H335, H314		
Melting point: -35°C		P305+P351+P338, P310, P280		
Assay: 34-37%				
Guaranteed analysis				
Colour	<10 APHA	Neodymium (Nd)	<0.1 ppb	
Aluminium (Al)	<1 ppb	Nickel (Ni)	<0.5 ppb	
Antimony (Sb)	<0.5 ppb	Niobium (Nb)	<0.1 ppb	
Arsenic (As)	<0.5 ppb	Palladium (Pd)	<0.5 ppb	
Barium (Ba)	<0.1 ppb	Platinum (Pt)	<0.5 ppb	
Beryllium (Be)	<0.1 ppb	Potassium (K)	<1 ppb	
Bismuth (Bi)	<0.1 ppb	Praseodymium (Pr)	<0.1 ppb	
Boron (B)	<1 ppb	Rhenium (Re)	<0.1 ppb	
Bromide (Br)	<10 ppm	Rhodium (Rh)	<0.1 ppb	
Cadmium (Cd)	<0.1 ppb	Rubidium (Rb)	<0.1 ppb	
Calcium (Ca)	<1 ppb	Ruthenium (Ru)	<0.1 ppb	
Cerium (Ce)	<0.1 ppb	Samarium (Sm)	<0.1 ppb	
Cesium (Cs)	<0.1 ppb	Scandium (Sc)	<0.1 ppb	
Free Chlorine (Cl ₂)	<0.5 ppm	Selenium (Se)	<1 ppb	
Chromium (Cr)	<0.5 ppb	Silver (Ag)	<1 ppb	
Cobalt (Co)	<0.1 ppb	Sodium (Na)	<1 ppb	
Copper (Cu)	<0.5 ppb	Strontium (Sr)	<0.1 ppb	
Dysprosium (Dy)	<0.1 ppb	Tantalum (Ta)	<1 ppb	
Erbium (Er)	<0.1 ppb	Tellurium (Te)	<0.1 ppb	
Europium (Eu)	<0.1 ppb	Terbium (Tb)	<0.1 ppb	
Gadolinium (Gd)	<0.1 ppb	Thallium (Tl)	<0.1 ppb	
Gallium (Ga)	<0.1 ppb	Thorium (Th)	<0.1 ppb	
Gold (Au)	<0.5 ppb	Thulium (Tm)	<0.1 ppb	
Hafnium (Hf)	<0.1 ppb	Tin (Sn)	<0.5 ppb	
Holmium (Ho)	<0.1 ppb	Titanium (Ti)	<0.5 ppb	
Indium (In)	<0.1 ppb	Total phosphorus (P)	<0.01 ppm	
Iron (Fe)	<1 ppb	Total sulphur (S)	<0.3 ppm	
Lanthanum (La)	<0.1 ppb	Tungsten (W)	<0.1 ppb	
Lead (Pb)	<0.1 ppb	Uranium (U)	<0.1 ppb	
Lithium (Li)	<0.1 ppb	Vanadium (V)	<0.5 ppb	
Lutetium (Lu)	<0.1 ppb	Ytterbium (Yb)	<0.1 ppb	
Magnesium (Mg)	<0.5 ppb	Yttrium (Y)	<0.1 ppb	
Manganese (Mn)	<0.1 ppb	Zinc (Zn)	<1 ppb	
Mercury (Hg)	<0.1 ppb	Zirconium (Zr)	<0.1 ppb	

Hydrochloric acid 37% d=1.18 **NEW PACK.** Primar Plus
for trace metal analysis

packaging		Mfr. No
250 mℓ	HDPE Plastic bottle	H/1196/PB07
500 mℓ	HDPE Plastic bottle	H/1196/PB08
1 ℓ	HDPE Plastic bottle	H/1196/PB15
2.5 ℓ	HDPE Plastic bottle	H/1196/PB17
HCl	ADR/RID: 8, II	
CAS: 7647-01-0	EINECS: 231-595-7	
MW: 36.46	UN: 1789	
Density at 20°C: 1.18	GHS: Skin Corr. 1B, STOT SE 3	
Boiling point: 57°C	H335, H314	
Melting point: -35°C	P305+P351+P338, P310, P280	


Assay: 35 to 38%			
Guaranteed analysis			
Aluminium (Al)	<0.02 ppm	Magnesium (Mg)	<0.02 ppm
Antimony (Sb)	<0.005 ppm	Manganese (Mn)	<0.002 ppm
Arsenic (As)	<0.01 ppm	Mercury (Hg)	<0.005 ppm
Barium (Ba)	<0.02 ppm	Molybdenum (Mo)	<0.002 ppm
Beryllium (Be)	<0.005 ppm	Nickel (Ni)	<0.005 ppm
Bismuth (Bi)	<0.005 ppm	Potassium (K)	<0.05 ppm
Boron (B)	<0.005 ppm	Silver (Ag)	<0.002 ppm
Cadmium (Cd)	<0.002 ppm	Sodium (Na)	<0.2 ppm
Calcium (Ca)	<0.1 ppm	Strontium (Sr)	<0.002 ppm
Free Chlorine (Cl ₂)	<1 ppm	Thallium (Tl)	<0.005 ppm
Chromium (Cr)	<0.005 ppm	Tin (Sn)	<0.005 ppm
Cobalt (Co)	<0.002 ppm	Titanium (Ti)	<0.005 ppm
Copper (Cu)	<0.005 ppm	Total nitrogen (N)	<5 ppm
Gallium (Ga)	<0.005 ppm	Total phosphorus (P)	<0.2 ppm
Germanium (Ge)	<0.005 ppm	Vanadium (V)	<0.002 ppm
Indium (In)	<0.005 ppm	Zinc (Zn)	<0.01 ppm
Iron (Fe)	<0.1 ppm	Residue after ignition	<0.0002%
Lead (Pb)	<0.002 ppm	Sulfate (SO ₄)	<1 ppm
Lithium (Li)	<0.002 ppm		

Hydrochloric acid 37% d=1.18
for AAS

packaging		Mfr. No
2.5 ℓ	Coated Glass bottle	H/1185/PB17
HCl	ADR/RID: 8, II	
CAS: 7647-01-0	EINECS: 231-595-7	
MW: 36.46	UN: 1789	
Density at 20°C: 1.18	GHS: Skin Corr. 1B, STOT SE 3	
Boiling point: 57°C	H335, H314	
Melting point: -35°C	P305+P351+P338, P310, P280	

Assay: 35 to 38%			
Guaranteed analysis			
Cadmium (Cd)	<0.02 ppm	Magnesium (Mg)	<0.2 ppm
Calcium (Ca)	<0.5 ppm	Manganese (Mn)	<0.02 ppm
Cobalt (Co)	<0.02 ppm	Nickel (Ni)	<0.02 ppm
Copper (Cu)	<0.05 ppm	Potassium (K)	<0.2 ppm
Iron (Fe)	<0.2 ppm	Sodium (Na)	<1 ppm
Lead (Pb)	<0.05 ppm	Zinc (Zn)	<0.2 ppm

Hydrochloric acid has a concentration of 11.65 M (11.65 N).



Hydrochloric acid 37% d=1.18
for electronic use (MOS)

packaging		Mfr. No
2.5 ℓ	HDPE Plastic bottle	H/1111/PB17
25 ℓ	HDPE Plastic container	H/1111/25
HCl	ADR/RID: 8, II	
CAS: 7647-01-0	EINECS: 231-595-7	
MW: 36.46	UN: 1789	
Density at 20°C: 1.18	GHS: Skin Corr. 1B, STOT SE 3	
Boiling point: 57°C	H335, H314	
Melting point: -35°C	P305+P351+P338, P310, P280	

Assay: 35.5 to 37.5%			
Guaranteed analysis			
Colour	<10 APHA	Manganese (Mn)	<0.02 ppm
Free Chlorine (Cl ₂)	<2 ppm	Molybdenum (Mo)	<0.02 ppm
Aluminium (Al)	<0.5 ppm	Nickel (Ni)	<0.01 ppm
Barium (Ba)	<0.1 ppm	Potassium (K)	<0.5 ppm
Cadmium (Cd)	<0.05 ppm	Sodium (Na)	<2 ppm
Calcium (Ca)	<1 ppm	Strontium (Sr)	<0.02 ppm
Chromium (Cr)	<0.1 ppm	Titanium (Ti)	<0.02 ppm
Cobalt (Co)	<0.02 ppm	Total phosphorus (P)	<1 ppm
Copper (Cu)	<0.1 ppm	Total silicon (Si)	<0.5 ppm
Iron (Fe)	<0.5 ppm	Total sulphur (S)	<5 ppm
Lead (Pb)	<0.1 ppm	Vanadium (V)	<0.02 ppm
Lithium (Li)	<0.02 ppm	Zinc (Zn)	<0.2 ppm
Magnesium (Mg)	<0.5 ppm		

Filtered to 0.2 μm.

Hydrochloric acid 37% **NEW SPECS** Certified, Eur.Ph.
for analysis with low Cadmium, Lead & Mercury level

packaging		Mfr. No
500 mℓ	HDPE Plastic bottle	H/1180/PB08
1 ℓ	HDPE Plastic bottle	H/1180/PB15
2.5 ℓ	HDPE Plastic bottle	H/1180/PB17
HCl	ADR/RID: 8, II	
CAS: 7647-01-0	EINECS: 231-595-7	
MW: 36.46	UN: 1789	
Density at 20°C: 1.18	GHS: Skin Corr. 1B, STOT SE 3	
Boiling point: 57°C	H335, H314	
Melting point: -35°C	P305+P351+P338, P310, P280	

Assay: 35.5 to 38%			
Guaranteed analysis			
Colour	<10 APHA	Lithium (Li)	<0.01 ppm
Free Chlorine (Cl ₂)	<1 ppm	Magnesium (Mg)	<0.1 ppm
Zirconium (Zr)	<0.02 ppm	Manganese (Mn)	<0.01 ppm
Gallium (Ga)	<0.02 ppm	Mercury (Hg)	<0.005 ppm
Aluminium (Al)	<0.05 ppm	Molybdenum (Mo)	<0.02 ppm
Arsenic (As)	<0.01 ppm	Nickel (Ni)	<0.02 ppm
Barium (Ba)	<0.02 ppm	Phosphate (PO ₄)	<0.5 ppm
Beryllium (Be)	<0.01 ppm	Potassium (K)	<0.1 ppm
Bismuth (Bi)	<0.02 ppm	Sodium (Na)	<0.5 ppm
Cadmium (Cd)	<0.01 ppm	Strontium (Sr)	<0.01 ppm
Calcium (Ca)	<0.5 ppm	Sulfate (SO ₄)	<1 ppm
Chromium (Cr)	<0.02 ppm	Thallium (Tl)	<0.02 ppm
Cobalt (Co)	<0.01 ppm	Titanium (Ti)	<0.02 ppm
Copper (Cu)	<0.02 ppm	Vanadium (V)	<0.02 ppm
Germanium (Ge)	<0.02 ppm	Zinc (Zn)	<0.05 ppm
Iron (Fe)	<0.2 ppm	Residue after ignition	<5 ppm
Lead (Pb)	<0.01 ppm		

Lot analysis on the label



Hydrochloric acid 37% d=1.18 Certified AR
for analysis

packaging		Mfr. No
500 mℓ	HDPE Plastic bottle	H/1200/PB08
1 ℓ	HDPE Plastic bottle	H/1200/PB15
1 ℓ	Coated Glass bottle	H/1200/PC15
2.5 ℓ	HDPE Plastic bottle	H/1200/PB17
2.5 ℓ	Coated Glass bottle	H/1200/PC17
25 ℓ	HDPE Plastic container	H/1200/25
HCl	ADR/RID: 8, II	
CAS: 7647-01-0	EINECS: 231-595-7	
MW: 36.4536.46 g/mol	UN: 1789	
Density at 20°C: 1.18	GHS: Skin Corr. 1B, STOT SE 3	
Boiling point: 57°C	H335, H314	
Melting point: -35°C	P305+P351+P338, P310, P280	

Assay: 35 to 38%			
Guaranteed analysis			
Colour	<10 APHA	Mercury (Hg)	<0.05 ppm
Residue after ignition	<0.0005%	Molybdenum (Mo)	<0.02 ppm
Aluminium (Al)	<0.5 ppm	Nickel (Ni)	<0.1 ppm
Arsenic (As)	<0.01 ppm	Potassium (K)	<0.5 ppm
Barium (Ba)	<0.1 ppm	Silver (Ag)	<0.02 ppm
Cadmium (Cd)	<0.02 ppm	Sodium (Na)	<2 ppm
Calcium (Ca)	<1 ppm	Strontium (Sr)	<0.02 ppm
Free Chlorine (Cl ₂)	<1 ppm	Sulfate (SO ₄)	<5 ppm
Cobalt (Co)	<0.02 ppm	Tin (Sn)	<0.1 ppm
Copper (Cu)	<0.1 ppm	Titanium (Ti)	<0.05 ppm
Iron (Fe)	<0.5 ppm	Total phosphorus (P)	<1 ppm
Lead (Pb)	<0.1 ppm	Vanadium (V)	<0.01 ppm
Lithium (Li)	<0.02 ppm	Zinc (Zn)	<0.5 ppm
Magnesium (Mg)	<0.5 ppm	Chromium (Cr)	<0.1 ppm
Manganese (Mn)	<0.02 ppm		

Lot analysis on the label.
Hydrochloric acid is equivalent to a strength of 11.65M (11.65N).

Hydrochloric acid 37% d=1.18 SLR
specified reagents for general laboratory work- extra pure

packaging		Mfr. No
1 ℓ	HDPE Plastic bottle	H/1150/PB15
2.5 ℓ	HDPE Plastic bottle	H/1150/PB17
5 ℓ	HDPE Plastic bottle	H/1150/21
25 ℓ	HDPE Plastic container	H/1150/25
45 ℓ	Plastic container	H/1150/26
HCl	ADR/RID: 8, II	
CAS: 7647-01-0	EINECS: 231-595-7	
MW: 36.46	UN: 1789	
Density at 20°C: 1.18	GHS: Skin Corr. 1B, STOT SE 3	
Boiling point: 57°C	H335, H314	
Melting point: -35°C	P305+P351+P338, P310, P280	

Assay: 35 to 38%			
Guaranteed analysis			
Colour	<20 APHA	Magnesium (Mg)	<5 ppm
Free Chlorine (Cl ₂)	<10 ppm	Potassium (K)	<5 ppm
Calcium (Ca)	<10 ppm	Sodium (Na)	<10 ppm
Copper (Cu)	<2 ppm	Total phosphorus (P)	<10 ppm
Iron (Fe)	<5 ppm	Total sulphur (S)	<20 ppm
Lead (Pb)	<2 ppm	Zinc (Zn)	<2 ppm

Hydrochloric acid is equivalent to a strength of 11.65M (11.65N).

Hydrochloric acid 32% d=1.16 Certified AR
for analysis

packaging		Mfr. No
1 ℓ	HDPE Plastic bottle	H/1100/PB15
2.5 ℓ	HDPE Plastic bottle	H/1100/PB17
45 ℓ	Plastic container	H/1100/26
HCl	ADR/RID: 8, II	
CAS: 7647-01-0	EINECS: 231-595-7	
MW: 36.46	UN: 1789	
Density at 20°C: 1.16	GHS: Skin Corr. 1B, Eye Dam. 1, STOT SE 3	
Boiling point: 57°C	H335, H314	
Melting point: -35°C	P301+P330+P331, P280, P305+P351+P338, P261, P301+P312, P304+P340	

Assay: 31.5 to 33%			
Guaranteed analysis			
Colour	<10 APHA	Nickel (Ni)	<0.05 ppm
Calcium (Ca)	<1 ppm	Potassium (K)	<0.5 ppm
Free Chlorine (Cl ₂)	<2 ppm	Sodium (Na)	<2 ppm
Chromium (Cr)	<0.05 ppm	Total phosphorus (P)	<1 ppm
Copper (Cu)	<0.05 ppm	Total silicon (Si)	<0.5 ppm
Iron (Fe)	<0.5 ppm	Total sulphur (S)	<5 ppm
Lead (Pb)	<0.05 ppm	Zinc (Zn)	<0.2 ppm
Magnesium (Mg)	<0.5 ppm		

Typical analysis			
Aluminium (Al)	<0.5 ppm	Lithium (Li)	<0.02 ppm
Silver (Ag)	<0.02 ppm	Manganese (Mn)	<0.02 ppm
Arsenic (As)	<0.05 ppm	Molybdenum (Mo)	<0.02 ppm
Barium (Ba)	<0.1 ppm	Strontium (Sr)	<0.02 ppm
Cadmium (Cd)	<0.05 ppm	Titanium (Ti)	<0.02 ppm
Cobalt (Co)	<0.02 ppm	Vanadium (V)	<0.02 ppm

Hydrochloric acid is equivalent to a strength of 10.2M (10.2N).
Lot analysis on the label.

Hydrochloric acid 32% d=1.16 SLR
specified reagents for general laboratory work- extra pure

packaging		Mfr. No
2.5 ℓ	HDPE Plastic bottle	H/1050/PB17
45 ℓ	Plastic container	H/1050/26
HCl	ADR/RID: 8, II	
CAS: 7647-01-0	EINECS: 231-595-7	
MW: 36.46	UN: 1789	
Density at 20°C: 1.16	GHS: Skin Corr. 1B, Eye Dam. 1, STOT SE 3	
Boiling point: 57°C	H335, H314	
Melting point: -35°C	P301+P330+P331, P280, P305+P351+P338, P261, P301+P312, P304+P340	

Assay: 31 to 33%			
Guaranteed analysis			
Colour	<20 APHA	Nickel (Ni)	<1 ppm
Free Chlorine (Cl ₂)	<10 ppm	Potassium (K)	<5 ppm
Calcium (Ca)	<10 ppm	Sodium (Na)	<10 ppm
Chromium (Cr)	<1 ppm	Total phosphorus (P)	<10 ppm
Copper (Cu)	<2 ppm	Total silicon (Si)	<2 ppm
Iron (Fe)	<5 ppm	Total sulphur (S)	<20 ppm
Lead (Pb)	<2 ppm	Zinc (Zn)	<2 ppm
Magnesium (Mg)	<5 ppm		

Hydrochloric acid is equivalent to a strength of 10.2M (10.2N).

Hydrochloric acid 32% d=1.16
Technical, d=1.16

packaging		Mfr. No
2.5 ℓ	HDPE Plastic bottle	H/1000/PB17
HCl	ADR/RID: 8, II	
CAS: 7647-01-0	EINECS: 231-595-7	
MW: 36.46	UN: 1789	
Density at 20°C: 1.16	GHS: Skin Corr. 1B, Eye Dam. 1, STOT SE 3	
Boiling point: 57°C	H335, H314	
Melting point: -35°C	P301+P330+P331, P280, P305+P351+P338, P261, P301+P312, P304+P340	



Hydrochloric acid 25% **NEW** Certified AR, Eur.Ph. for analysis with low Cadmium, Lead & Mercury level

packaging		Mfr. No
1 ℓ	HDPE Plastic bottle	H/1020/PB15
2.5 ℓ	HDPE Plastic bottle	H/1020/PB17
HCl		UN: 1789
CAS: 7647-01-0		GHS: Skin Corr. 1B, Eye
MW: 36.46		Dam. 1, STOT SE 3
Density at 20°C: 1.12		H335, H314
Boiling point: 107°C		P280, P301+P330+P331,
Melting point: -70°C		P310, P305+P351+P338,
ADR/RID: 8, II		P303+P361+P353, P304+P340
EINECS: 231-595-7		
Assay: >25.0%		

Guaranteed analysis			
Colour	<10 APHA	Lead (Pb)	<0.02 ppm
Free Chlorine (Cl ₂)	<1 ppm	Lithium (Li)	<0.01 ppm
Zirconium (Zr)	<0.05 ppm	Magnesium (Mg)	<0.2 ppm
Aluminium (Al)	<0.2 ppm	Manganese (Mn)	<0.01 ppm
Arsenic (As)	<0.01 ppm	Mercury (Hg)	<0.01 ppm
Barium (Ba)	<0.05 ppm	Molybdenum (Mo)	<0.02 ppm
Beryllium (Be)	<0.02 ppm	Nickel (Ni)	<0.02 ppm
Bismuth (Bi)	<0.05 ppm	Potassium (K)	<0.1 ppm
Cadmium (Cd)	<0.01 ppm	Sodium (Na)	<0.5 ppm
Calcium (Ca)	<1 ppm	Strontium (Sr)	<0.01 ppm
Chromium (Cr)	<0.02 ppm	Sulfate (SO ₄)	<1 ppm
Cobalt (Co)	<0.01 ppm	Thallium (Tl)	<0.05 ppm
Copper (Cu)	<0.02 ppm	Titanium (Ti)	<0.05 ppm
Germanium (Ge)	<0.05 ppm	Vanadium (V)	<0.01 ppm
Iron (Fe)	<0.2 ppm	Zinc (Zn)	<0.05 ppm

Lot analysis on the label

Hydrochloric acid 20% d=1.1 **NEW** Primar for trace metal analysis

packaging		Mfr. No
500 mℓ	HDPE Plastic bottle	H/1190/PB08
2.5 ℓ	HDPE Plastic bottle	H/1190/PB17
HCl		EINECS: 231-595-7
CAS: 7647-01-0		UN: 1789
MW: 36.46		GHS: Skin Irrit. 2, Eye Irrit. 2,
Density at 20°C: 1.05 to 1.12		STOT SE 3
Boiling point: 110°C		H315, H319, H335
ADR/RID: 8, II		P280, P305+P351+P338
Assay: >20%		

Guaranteed analysis			
Aluminium (Al)	<0.05 ppm	Mercury (Hg)	<0.002 ppm
Arsenic (As)	<0.02 ppm	Molybdenum (Mo)	<0.005 ppm
Barium (Ba)	<0.05 ppm	Nickel (Ni)	<0.005 ppm
Cadmium (Cd)	<0.005 ppm	Potassium (K)	<0.2 ppm
Calcium (Ca)	<0.2 ppm	Propan--1-ol	<0.005 ppm
Chromium (Cr)	<0.005 ppm	Sodium (Na)	<1 ppm
Cobalt (Co)	<0.005 ppm	Strontium (Sr)	<0.005 ppm
Copper (Cu)	<0.01 ppm	Titanium (Ti)	<0.005 ppm
Iron (Fe)	<0.1 ppm	Total phosphorus (P)	<0.2 ppm
Lead (Pb)	<0.001 ppm	Total silicon (Si)	<0.5 ppm
Lithium (Li)	<0.002 ppm	Total sulphur (S)	<0.5 ppm
Magnesium (Mg)	<0.02 ppm	Vanadium (V)	<0.005 ppm
Manganese (Mn)	<0.005 ppm	Zinc (Zn)	<0.02 ppm

Hydrochloric acid solution 2M (2N) pure, Bench Reagent

packaging		Mfr. No
2.5 ℓ	HDPE Plastic bottle	J/4250/17
Density at 20°C: 1.0 to 1.2		UN: 1789
Boiling point: 100°C		GHS: Skin Corr. 1B
Melting point: 0°C		H314
ADR/RID: 8, III		P280, P305+P351+P338,
EINECS: 231-595-7		P301+P330+P331, P310

Manufactured to a tolerance of ±5%.

Hydrochloric acid solution 5M pure, Bench Reagent

packaging		Mfr. No
2.5 ℓ	HDPE Plastic bottle	J/4270/17
25 ℓ	HDPE Plastic container	J/4270/25
CAS: 7647-01-0		UN: 1789
Density at 20°C: 1.0 to 1.2		GHS: Skin Irrit. 2, Eye Irrit. 2,
Boiling point: 81.5 to 110°C		STOT SE 3
Melting point: -74°C		H315, H319, H335
ADR/RID: 8, II		P280, P305+P351+P338
EINECS: 231-595-7		
Manufactured to a tolerance of ±2.5%.		

Hydrochloric acid solution 0.2M (0.2N) NIST Standard solution ready to use, For Volumetric analysis

packaging		Mfr. No
1 ℓ	HDPE Plastic bottle	J/4340/15
2.5 ℓ	HDPE Plastic bottle	J/4340/17
CAS: 7647-01-0		UN: 1789
Density at 20°C: 1		GHS: Skin Irrit. 2, Eye Irrit. 2,
Boiling point: 100°C		STOT SE 3
Melting point: -35°C		H315, H319, H335
ADR/RID: 8, II		P280, P305+P351+P338
EINECS: 231-595-7		

Manufactured to a tolerance of ± 0.1%. See the full range of Standard Volumetric solution page 323.

Hydrochloric acid solution 0.1M (0.1N) Eur.Ph. , BP NIST Standard solution ready to use, For Volumetric analysis

packaging		Mfr. No
1 ℓ	HDPE Plastic bottle	J/4350/15
2.5 ℓ	HDPE Plastic bottle	J/4350/17
5 ℓ	HDPE Plastic bottle	J/4350/21
10 ℓ	Bag in box	J/4350/24
25 ℓ	HDPE Plastic container	J/4350/25

CAS: 7647-01-0		UN: 1789
Density at 20°C: 1		GHS: Met. Corr. 1
Boiling point: 100°C		H290
EINECS: 231-595-7		P234, P390

Manufactured to a tolerance of ± 0.1%. See the full range of Standard Volumetric solution page 323.

Hydrochloric acid solution 0.1M (0.1N) NIST Standard Concentrate, For Volumetric analysis

packaging		Mfr. No
1	Ampoule	J/4350C/05
6	Ampoule	J/4350C/90
CAS: 7647-01-0		UN: 1789
Density at 20°C: 1.00 to 1.05		GHS: Met. Corr. 1
Boiling point: 100 to 103°C		H290
ADR/RID: 8, III		P234, P390
EINECS: 231-595-7		

Stabilised with 10 ppm of Mercury(II) chloride. Each ampoule makes 1 ℓ of solution. See the full range of Standard Volumetric solution page 323.

Hydrochloric acid solution 0.5M (0.5N) NIST Standard solution ready to use, For Volumetric analysis

packaging		Mfr. No
1 ℓ	HDPE Plastic bottle	J/4330/15
2.5 ℓ	HDPE Plastic bottle	J/4330/17
25 ℓ	HDPE Plastic container	J/4330/25
CAS: 7647-01-0		UN: 1789
Density at 20°C: 1.00 to 1.05		GHS: Met. Corr. 1
Boiling point: 100 to 103°C		H290
ADR/RID: 8, III		P234, P390
EINECS: 231-595-7		

Manufactured to a tolerance of ± 0.1%. See the full range of Standard Volumetric solution page 323.

Hydrochloric acid solution 0.02M (0.02N) NIST Standard solution ready to use, For Volumetric analysis

packaging		Mfr. No
1 ℓ	HDPE Plastic bottle	J/4360/PB15
2.5 ℓ	HDPE Plastic bottle	J/4360/17
CAS: 7647-01-0		UN: 1789
Density at 20°C: 1		GHS: Skin Irrit. 2, Eye Irrit. 2,
Boiling point: 100°C		STOT SE 3
EINECS: 231-595-7		H315, H319, H335

Manufactured to a tolerance of ± 0.3%. See the full range of Standard Volumetric solution page 323.

Hydrochloric acid solution 0.5M (0.5N) NIST Standard Concentrate, For Volumetric analysis

packaging		Mfr. No
1	Ampoule	J/4330C/05
6	Ampoule	J/4330C/90
CAS: 7647-01-0		UN: 1789
Density at 20°C: 1.18		GHS: Skin Irrit. 2, Eye Irrit. 2,
Boiling point: 57°C		STOT SE 3
Melting point: -35°C		H315, H319, H335
ADR/RID: 8, II		P280, P305+P351+P338
EINECS: 231-595-7		

Each ampoule makes 1 ℓ of solution. See the full range of Standard Volumetric solution page 323.

Hydrochloric acid solution 1M (1N) Eur.Ph. , USP, BP NIST Standard solution ready to use, For Volumetric analysis

packaging		Mfr. No
1 ℓ	HDPE Plastic bottle	J/4320/15
2.5 ℓ	HDPE Plastic bottle	J/4320/17
5 ℓ	HDPE Plastic bottle	J/4320/21
10 ℓ	Bag in box	J/4320/24

CAS: 7647-01-0		UN: 1789
Density at 20°C: 1.00 to 1.05		GHS: Met. Corr. 1
Boiling point: 100 to 103°C		H290
ADR/RID: 8, III		P234, P390
EINECS: 231-595-7		

Manufactured to a tolerance of ± 0.1%. See the full range of Standard Volumetric solution page 323.

Hydrochloric acid solution 1M (1N) NIST Standard Concentrate, For Volumetric analysis

packaging		Mfr. No
1	Ampoule	J/4320C/05
6	Ampoule	J/4320C/90
CAS: 7647-01-0		UN: 1789
Density at 20°C: 1.16		GHS: Skin Corr. 1B, Eye
Boiling point: 57°C		Dam. 1, STOT SE 3
Melting point: -35°C		H335, H314
ADR/RID: 8, II		P301+P330+P331, P280,
EINECS: 231-595-7		P305+P351+P338, P261,
UN: 1789		P301+P312, P304+P340

Each ampoule makes 1 ℓ of solution. See the full range of Standard Volumetric solution page 323.

Optima LC/MS

Fisher Chemical OPTIMA® LC/MS solvents, blends, additives and modifiers have set the standard of excellence for consistent, reproducible performance in the mobile phase of LC/MS.

Find out more on page 308

Hydrochloric acid solution 5M (5N) Standard solution ready to use for volumetric analysis

packaging		Mfr. No
1 ℓ	HDPE Plastic bottle	J/4310/15
2.5 ℓ	HDPE Plastic bottle	J/4310/17
5 ℓ	HDPE Plastic bottle	J/4310/21
CAS: 7647-01-0		UN: 1789
Density at 20°C: 1.0 to 1.2		GHS: Skin Irrit. 2, Eye Irrit. 2,
Boiling point: 81.5 to 110°C		STOT SE 3
Melting point: -74°C		H315, H319, H335
ADR/RID: 8, II		P280, P305+P351+P338
EINECS: 231-595-7		

Manufactured to a tolerance of ± 0.5%. See the full range of Standard Volumetric solution page 323.

Hydrochloric acid solution 2M (2N) Eur.Ph. , BP NIST Standard solution ready to use, For Volumetric analysis

packaging		Mfr. No
1 ℓ	HDPE Plastic bottle	J/4315/15
2.5 ℓ	HDPE Plastic bottle	J/4315/17
10 ℓ	Bag in box	J/4315/24
CAS: 7647-01-0		UN: 1789
Density at 20°C: 1.00 to 1.05		GHS: Met. Corr. 1
Boiling point: 100 to 103°C		H290
ADR/RID: 8, III		P234, P390
EINECS: 231-595-7		

Manufactured to a tolerance of ± 0.1%. See the full range of Standard Volumetric solution page 323.

Hydrofluoric acid 58-62% Technical

packaging		Mfr. No
1 ℓ	HDPE Plastic bottle	H/1300/15
2.5 ℓ	HDPE Plastic bottle	H/1300/17
HF		UN: 1790
CAS: 7664-39-3		GHS: Acute Tox. 1, Skin
MW: 20.01		Corr. 1A, Eye Dam. 1
Density at 20°C: 1.15 to 1.20		H310, H314, H300, H330
Boiling point: 105°C		P280, P303+P361+P353,
Melting point: -35°C		P301+P310,
ADR/RID: 8, 6.1, II		P301+P330+P331,
EINECS: 231-634-8		P305+P351+P338, P304+P340

Assay: 58 to 62%

Label

Fisher Chemical labels show GHS pictograms clearly. Together with actual lot analysis shown on the label, plus a wealth of additional information, we help you to work safely to achieve accurate and consistent results in your analytical applications.

Hydrofluoric acid 47-51%
for ultra trace elemental analysis

packaging		Mfr. No
250 mℓ	Pre-cleaned Teflon PFA bottle	A463-250
500 mℓ	Pre-cleaned Teflon PFA bottle	A463-500
1 ℓ	Pre-cleaned Teflon PFA bottle	A463-1
HF		UN: 1790
CAS: 7664-39-3		GHS: Acute Tox. 1, Skin
MW: 20.00		Corr. 1A, Eye Dam. 1
Density at 20°C: 1.15 to 1.20		H310, H314, H300, H330
Boiling point: 105°C		P280, P303+P361+P353,
Melting point: -35°C		P301+P310,
ADR/RID: 8, 6.1, II		P301+P330+P331,
EINECS: 231-634-8		P305+P351+P338, P304+P340

Assay: 47-51%			
Guaranteed analysis			
Zirconium (Zr)	<10 ppt	Mercury (Hg)	<50 ppt
Gallium (Ga)	<10 ppt	Molybdenum (Mo)	<10 ppt
Terbium (Tb)	<1 ppt	Neodymium (Nd)	<1 ppt
Indium (In)	<1 ppt	Nickel (Ni)	<20 ppt
Dysprosium (Dy)	<1 ppt	Niobium (Nb)	<10 ppt
Aluminium (Al)	<20 ppt	Palladium (Pd)	<20 ppt
Antimony (Sb)	<20 ppt	Potassium (K)	<10 ppt
Arsenic (As)	<50 ppt	Praseodymium (Pr)	<1 ppt
Barium (Ba)	<10 ppt	Rhenium (Re)	<10 ppt
Beryllium (Be)	<10 ppt	Rhodium (Rh)	<20 ppt
Bismuth (Bi)	<10 ppt	Ruthenium (Ru)	<20 ppt
Boron (B)	<100 ppt	Samarium (Sm)	<1 ppt
Cadmium (Cd)	<10 ppt	Scandium (Sc)	<10 ppt
Calcium (Ca)	<10 ppt	Propan--1-ol	<10 ppt
Cerium (Ce)	<10 ppt	Sodium (Na)	<10 ppt
Chromium (Cr)	<10 ppt	Strontium (Sr)	<10 ppt
Cobalt (Co)	<10 ppt	Tellurium (Te)	<1 ppt
Copper (Cu)	<10 ppt	Thallium (Tl)	<10 ppt
Erbium (Er)	<1 ppt	Thorium (Th)	<1 ppt
Europium (Eu)	<1 ppt	Thulium (Tm)	<1 ppt
Gadolinium (Gd)	<1 ppt	Tin (Sn)	<20 ppt
Germanium (Ge)	<10 ppt	Titanium (Ti)	<20 ppt
Gold (Au)	<20 ppt	Total sulphur (S)	<100 ppb
Hafnium (Hf)	<10 ppt	Uranium (U)	<1 ppt
Holmium (Ho)	<1 ppt	Vanadium (V)	<10 ppt
Iron (Fe)	<10 ppt	Ytterbium (Yb)	<1 ppt
Lanthanum (La)	<10 ppt	Yttrium (Y)	<1 ppt
Lead (Pb)	<10 ppt	Zinc (Zn)	<10 ppt
Lithium (Li)	<10 ppt	Platinum (Pt)	<20 ppt
Lutetium (Lu)	<1 ppt	Rubidium (Rb)	<20 ppt
Magnesium (Mg)	<10 ppt	Tungsten (W)	<20 ppt
Manganese (Mn)	<10 ppt		

Bottled in class 10 cleanroom conditions.Individually double-bagged in class 100 cleanroom conditions. Each bottle is individually packaged. Certificate of analysis delivered with each bottle.

Hydrofluoric acid 47-51%
for trace metal analysis

packaging		Mfr. No
500 mℓ	HDPE plastic bottle	A513-P500
HF		UN: 1790
CAS: 7664-39-3		GHS: Acute Tox. 1, Skin
MW: 20.01		Corr. 1A, Eye Dam. 1
Density at 20°C: 1.15 to 1.20		H310, H314, H300, H330
Boiling point: 105°C		P280, P303+P361+P353,
Melting point: -35°C		P301+P310,
ADR/RID: 8, 6.1, II		P301+P330+P331,
EINECS: 231-634-8		P305+P351+P338, P304+P340

Assay: 47-51%			
Guaranteed analysis			
Colour	<10 APHA	Cerium (Ce)	<0.1 ppb
Aluminium (Al)	<1ppb	Cesium (Cs)	<0.1 ppb
Antimony (Sb)	<0.2 ppb	Chromium (Cr)	<1 ppb
Arsenic (As)	<0.5 ppb	Cobalt (Co)	<0.1 ppb
Barium (Ba)	<0.1 ppb	Copper (Cu)	<0.5 ppb
Beryllium (Be)	<0.1 ppb	Dysprosium (Dy)	<0.1 ppb
Bismuth (Bi)	<0.1 ppb	Erbium (Er)	<0.1 ppb
Boron (B)	<1 ppb	Europium (Eu)	<0.1 ppb
Cadmium (Cd)	<0.1 ppb	Fluorosilicic acid (SiF ₆)	<20 ppm
Calcium (Ca)	<1ppb		

Gadolinium (Gd)	<0.1 ppb	Rubidium (Rb)	<0.1 ppb
Gallium (Ga)	<0.1 ppb	Ruthenium (Ru)	<0.1 ppb
Germanium (Ge)	<0.1 ppb	Samarium (Sm)	<0.1 ppb
Gold (Au)	<0.2 ppb	Scandium (Sc)	<0.1 ppb
Hafnium (Hf)	<0.1 ppb	Selenium (Se)	<1 ppb
Holmium (Ho)	<0.1 ppb	Silver (Ag)	<0.5 ppb
Indium (In)	<0.1 ppb	Sodium (Na)	<1 ppb
Iron (Fe)	<1 ppb	Strontium (Sr)	<0.1 ppb
Lanthanum (La)	<0.1 ppb	Tantalum (Ta)	<0.1 ppb
Lead (Pb)	<0.1 ppb	Tellurium (Te)	<0.1 ppb
Lithium (Li)	<0.1 ppb	Terbium (Tb)	<0.1 ppb
Lutetium (Lu)	<0.1 ppb	Thallium (Tl)	<0.1 ppb
Magnesium (Mg)	<1 ppb	Thorium (Th)	<0.1 ppb
Manganese (Mn)	<0.1 ppb	Thulium (Tm)	<0.1 ppb
Mercury (Hg)	<1 ppb	Tin (Sn)	<0.1 ppb
Molybdenum (Mo)	<0.1 ppb	Titanium (Ti)	<1 ppb
Neodymium (Nd)	<0.1 ppb	Total chloride (Cl)	<4 ppm
Nickel (Ni)	<0.5 ppb	Total phosphorus (P)	<0.05 ppm
Niobium (Nb)	<0.1 ppb	Total sulphur (S)	<0.5 ppm
Palladium (Pd)	<0.2 ppb	Tungsten (W)	<0.5 ppb
Platinum (Pt)	<0.2 ppb	Uranium (U)	<0.1 ppb
Potassium (K)	<1 ppb	Vanadium (V)	<0.1 ppb
Praseodymium (Pr)	<0.1 ppb	Ytterbium (Yb)	<0.1 ppb
Rhenium (Re)	<0.1 ppb	Yttrium (Y)	<0.1 ppb
Rhodium (Rh)	<0.1 ppb	Zinc (Zn)	<1 ppb

Bottled in class 10 cleanroom conditions. Lot analysis on the label.

Hydrofluoric acid 48%
for electronic use (MOS)

packaging		Mfr. No
2.5 ℓ	HDPE Plastic bottle	H/1425/17
HF		UN: 1790
CAS: 7664-39-3		GHS: Acute Tox. 1, Skin
MW: 20.01		Corr. 1A, Eye Dam. 1
Density at 20°C: 1.15 to 1.20		H310, H314, H300, H330
Boiling point: 105°C		P280, P303+P361+P353,
Melting point: -35°C		P301+P310,
ADR/RID: 8, 6.1, II		P301+P330+P331,
EINECS: 231-634-8		P305+P351+P338, P304+P340

Assay: >48%			
Guaranteed analysis			
Colour	<10 APHA	Molybdenum (Mo)	<0.05 ppm
Aluminium (Al)	<0.5 ppm	Nickel (Ni)	<0.1 ppm
Barium (Ba)	<0.02 ppm	Potassium (K)	<0.1 ppm
Cadmium (Cd)	<0.02 ppm	Sodium (Na)	<0.2 ppm
Calcium (Ca)	<0.5 ppm	Strontium (Sr)	<0.02 ppm
Chromium (Cr)	<0.05 ppm	Titanium (Ti)	<0.02 ppm
Cobalt (Co)	<0.02 ppm	n-Butyl formate	<0.0005%
Copper (Cu)	<0.1 ppm	Total phosphorus (P)	<0.5 ppm
Iron (Fe)	<0.2 ppm	Total silicon (Si)	<10 ppm
Lead (Pb)	<0.1 ppm	Total sulphur (S)	<2 ppm
Lithium (Li)	<0.02 ppm	Vanadium (V)	<0.02 ppm
Magnesium (Mg)	<0.2 ppm	Zinc (Zn)	<0.5 ppm
Manganese (Mn)	<0.05 ppm		

Hydrofluoric acid 48%
for analysis

packaging		Mfr. No
1 ℓ	HDPE Plastic bottle	H/1420/15
2.5 ℓ	HDPE Plastic bottle	H/1420/17
HF		UN: 1790
CAS: 7664-39-3		GHS: Acute Tox. 1, Skin
MW: 20.01		Corr. 1A, Eye Dam. 1
Density at 20°C: 1.15 to 1.20		H310, H314, H300, H330
Boiling point: 105°C		P280, P303+P361+P353,
Melting point: -35°C		P301+P310,
ADR/RID: 8, 6.1, II		P301+P330+P331,
EINECS: 231-634-8		P305+P351+P338, P304+P340

Assay: >48%

Guaranteed analysis			
Colour	<10 APHA	Potassium (K)	<0.1 ppm
Calcium (Ca)	<0.5 ppm	Sodium (Na)	<0.2 ppm
Chromium (Cr)	<0.05 ppm	Total chloride (Cl)	<0.0005%
Copper (Cu)	<0.1 ppm	Total phosphorus (P)	<0.5 ppm
Iron (Fe)	<0.2 ppm	Total silicon (Si)	<10 ppm
Lead (Pb)	<0.1 ppm	Total sulphur (S)	<2 ppm
Magnesium (Mg)	<0.2 ppm	Zinc (Zn)	<0.1 ppm
Manganese (Mn)	<0.05 ppm	Arsenic (As)	<0.05 ppm
Nickel (Ni)	<0.05 ppm	Residue after ignition	<0.0005%

Typical analysis			
Aluminium (Al)	<0.05 ppm	Lithium (Li)	<0.02 ppm
Silver (Ag)	<0.02 ppm	Molybdenum (Mo)	<0.05 ppm
Barium (Ba)	<0.1 ppm	Strontium (Sr)	<0.02 ppm
Cadmium (Cd)	<0.02 ppm	Titanium (Ti)	<0.02 ppm
Cobalt (Co)	<0.02 ppm	Vanadium (V)	<0.02 ppm
Tin (Sn)	<0.02 ppm		

Lot analysis on the label.

Hydrofluoric acid 40%
for analysis

packaging		Mfr. No
1 ℓ	HDPE Plastic bottle	H/1400/15
2.5 ℓ	HDPE Plastic bottle	H/1400/17
HF		UN: 1790
CAS: 7664-39-3		GHS: Acute Tox. 1, Skin
MW: 20.01		Corr. 1A, Eye Dam. 1
Density at 20°C: 1.15 to1.20		H310, H314, H300, H330
Boiling point: 105°C		P280, P303+P361+P353,
Melting point: -35°C		P301+P310,
ADR/RID: 8, 6.1, II		P301+P330+P331,
EINECS: 231-634-8		P305+P351+P338, P304+P340

Assay: >40%			
Guaranteed analysis			
Colour	<10 APHA	Potassium (K)	<0.1 ppm
Calcium (Ca)	<0.5 ppm	Sodium (Na)	<0.5 ppm
Chromium (Cr)	<0.05 ppm	Total chloride (Cl)	<0.0005%
Copper (Cu)	<0.1 ppm	Total phosphorus (P)	<0.5 ppm
Iron (Fe)	<0.2 ppm	Total silicon (Si)	<10 ppm
Lead (Pb)	<0.05 ppm	Total sulphur (S)	<2 ppm
Magnesium (Mg)	<0.2 ppm	Zinc (Zn)	<0.1 ppm
Nickel (Ni)	<0.05 ppm		
Typical analysis			
Aluminium (Al)	<0.1 ppm	Manganese (Mn)	<0.05 ppm
Silver (Ag)	<0.02 ppm	Molybdenum (Mo)	<0.05 ppm
Barium (Ba)	<0.1 ppm	Strontium (Sr)	<0.02 ppm
Cadmium (Cd)	<0.02 ppm	Titanium (Ti)	<0.02 ppm
Cobalt (Co)	<0.02 ppm	Vanadium (V)	<0.02 ppm
Lithium (Li)	<0.02 ppm		

Lot analysis on the label. The hydrofluoric acid 40% has a concentration of 22.6 M (22.6 N).

Hydrofluoric acid 40%
specified reagents for general laboratory work- extra pure

packaging		Mfr. No
1 ℓ	HDPE Plastic bottle	H/1350/15
2.5 ℓ	HDPE Plastic bottle	H/1350/17
HF		UN: 1790
CAS: 7664-39-3		GHS: Acute Tox. 1, Skin
MW: 20.01		Corr. 1A, Eye Dam. 1
Density at 20°C: 1.15 to1.20		H310, H314, H300, H330
Boiling point: 105°C		P280, P303+P361+P353,
Melting point: -35°C		P301+P310,
ADR/RID: 8, 6.1, II		P301+P330+P331,
EINECS: 231-634-8		P305+P351+P338, P304+P340

Assay: >40%			
Guaranteed analysis			
Calcium (Ca)	<10 ppm	Potassium (K)	<5 ppm
Chromium (Cr)	<1 ppm	Sodium (Na)	<10 ppm
Copper (Cu)	<5 ppm	n-Butyl formate	<0.005%
Iron (Fe)	<10 ppm	Total phosphorus (P)	<5 ppm
Lead (Pb)	<5 ppm	Total silicon (Si)	<50 ppm
Magnesium (Mg)	<5 ppm	Total sulphur (S)	<10 ppm
Nickel (Ni)	<1 ppm	Zinc (Zn)	<5 ppm

The hydrofluoric acid 40% has a concentration of 22.6 M (22.6 N).

Hydrofluoric Acid Antidote Gel
pure

packaging		Mfr. No
25 g	Pack of 6	H/1430/90
ADR/RID: /II		
Guaranteed analysis		
Identity test		: Pass

Hydrogen peroxide 60% w/v (200 volumes)
specified reagents for general laboratory work- extra pure

packaging		Mfr. No
1 ℓ	HDPE Plastic bottle	H/1862/15
H ₂ O ₂		UN: 2014
CAS: 7722-84-1		GHS: Ox. Liq. 2, Acute Tox. 4, Skin Corr. 1A, STOT SE 3
MW: 34.01		H272, H302, H332, H314, H335
Boiling point: 108°C		P280, P305+P351+P338,
Melting point: -33°C		P301+P330+P331, P310,
ADR/RID: 5.1, 8, I		P304+P340, P210
EINECS: 231-765-0		

Assay: >59%			
Guaranteed analysis			
Acidity	<0.006 meq/g	Lead (Pb)	<1 ppm
Copper (Cu)	<1 ppm	n-Butyl formate	<0.001%
Iron (Fe)	<2 ppm		

Hydrogen peroxide 60% w/v (200 volumes)
pure

packaging		Mfr. No
1 ℓ	HDPE Plastic bottle	H/1860/15
H ₂ O ₂		UN: 2014
CAS: 7722-84-1		GHS: Ox. Liq. 2, Acute Tox. 4, Skin Corr. 1A, STOT SE 3
MW: 34.01		H272, H302, H332, H314, H335
Boiling point: 108°C		P280, P305+P351+P338,
Melting point: -33°C		P301+P330+P331, P310,
ADR/RID: 5.1, 8, I		P304+P340, P210
EINECS: 231-765-0		



Hydrogen peroxyde 30-32%
for ultra trace elemental analysis

Optima

packaging		Mfr. No
500 mℓ	Pre-cleaned Teflon FEP bottle	P170-500
H ₂ O ₂	EINECS: 231-765-0	
CAS: 7722-84-1	UN: 2014	
MW: 34.01	GHS: Acute Tox. 4, Eye Dam. 1	
Density at 20°C: 1.11	H318, H302	
Boiling point: 108°C	P280, P305+P351+P338,	
Melting point: -33°C	P301+P330+P331, P301+P312	
ADR/RID: 5.1, 8, II		
Assay: 30-32%		

Guaranteed analysis			
Aluminium (Al)	<50 ppt	Molybdenum (Mo)	<10 ppt
Antimony (Sb)	<10 ppt	Neodymium (Nd)	<1 ppt
Arsenic (As)	<100 ppt	Nickel (Ni)	<20 ppt
Barium (Ba)	<10 ppt	Niobium (Nb)	<10 ppt
Beryllium (Be)	<10 ppt	Palladium (Pd)	<10 ppt
Bismuth (Bi)	<10 ppt	Potassium (K)	<20 ppt
Boron (B)	<100 ppt	Praseodymium (Pr)	<1 ppt
Cadmium (Cd)	<10 ppt	Rhenium (Re)	<10 ppt
Calcium (Ca)	<100 ppt	Rhodium (Rh)	<10 ppt
Cerium (Ce)	<1 ppt	Rubidium (Rb)	<10 ppt
Cesium (Cs)	<1 ppt	Ruthenium (Ru)	<10 ppt
Chromium (Cr)	<10 ppt	Samarium (Sm)	<1 ppt
Cobalt (Co)	<10 ppt	Scandium (Sc)	<10 ppt
Copper (Cu)	<10 ppt	Selenium (Se)	<100 ppt
Dysprosium (Dy)	<1 ppt	Silver (Ag)	<10 ppt
Erbium (Er)	<1 ppt	Sodium (Na)	<50 ppt
Europium (Eu)	<1 ppt	Strontium (Sr)	<10 ppt
Gadolinium (Gd)	<1 ppt	Tantalum (Ta)	<10 ppt
Gallium (Ga)	<10 ppt	Tellurium (Te)	<1 ppt
Germanium (Ge)	<10 ppt	Terbium (Tb)	<1 ppt
Gold (Au)	<10 ppt	Thallium (Tl)	<1 ppt
Hafnium (Hf)	<1 ppt	Thorium (Th)	<1 ppt
Holmium (Ho)	<1 ppt	Thulium (Tm)	<1 ppt
Indium (In)	<1 ppt	Tin (Sn)	<50 ppt
Iron (Fe)	<20 ppt	Titanium (Ti)	<20 ppt
Lanthanum (La)	<1 ppt	Tungsten (W)	<20 ppt
Lead (Pb)	<10 ppt	Uranium (U)	<1 ppt
Lithium (Li)	<10 ppt	Vanadium (V)	<10 ppt
Lutetium (Lu)	<1 ppt	Ytterbium (Yb)	<1 ppt
Magnesium (Mg)	<20 ppt	Yttrium (Y)	<1 ppt
Manganese (Mn)	<10 ppt	Zinc (Zn)	<50 ppt
Mercury (Hg)	<50 ppt	Zirconium (Zr)	<10 ppt

Each bottle is individually packaged
Certificate of analysis delivered with each bottle

Hydrogen peroxide 30%w/v (100 volumes)
for trace metal analysis

Primar

packaging		Mfr. No
1 ℓ	HDPE Plastic bottle	H/1820/15
H ₂ O ₂	EINECS: 231-765-0	
CAS: 7722-84-1	UN: 2014	
MW: 34.01	GHS: Acute Tox. 4, Eye Dam. 1	
Density at 20°C: 1.11	H318, H302	
Boiling point: 108°C	P280, P305+P351+P338,	
Melting point: -33°C	P301+P330+P331, P301+P312	
ADR/RID: 5.1, 8, II		
Assay: >30 %		

Guaranteed analysis			
Aluminium (Al)	<0.1 ppm	Molybdenum (Mo)	<0.02 ppm
Arsenic (As)	<0.1 ppm	Nickel (Ni)	<0.005 ppm
Barium (Ba)	<0.005 ppm	Potassium (K)	<0.2 ppm
Cadmium (Cd)	<0.005 ppm	Propan--1-ol	<0.01 ppm
Calcium (Ca)	<0.2 ppm	Sodium (Na)	<2 ppm
Chromium (Cr)	<0.05 ppm	Strontium (Sr)	<0.01 ppm
Cobalt (Co)	<0.005 ppm	Titanium (Ti)	<0.01 ppm
Copper (Cu)	<0.01 ppm	n-Butyl formate	<0.0005%
Iron (Fe)	<0.2 ppm	Total phosphorus (P)	<2 ppm
Lead (Pb)	<0.01 ppm	Total silicon (Si)	<0.1 ppm
Lithium (Li)	<0.005 ppm	Total sulphur (S)	<0.5 ppm
Magnesium (Mg)	<0.1 ppm	Vanadium (V)	<0.01 ppm
Manganese (Mn)	<0.005 ppm	Zinc (Zn)	<0.02 ppm

Hydrogen peroxyde 30%w/v (100 volumes)
for electronic use (MOS)

packaging		Mfr. No
2.5 ℓ	HDPE Plastic bottle	H/1811/17
5 ℓ	HDPE Plastic bottle	H/1811/21
H ₂ O ₂	EINECS: 231-765-0	
CAS: 7722-84-1	UN: 2014	
MW: 34.01	GHS: Acute Tox. 4, Eye Dam. 1	
Density at 20°C: 1.11	H318, H302	
Boiling point: 108°C	P280, P305+P351+P338,	
Melting point: -33°C	P301+P330+P331, P301+P312	
ADR/RID: 5.1, 8, II		
Assay: >30%		

Guaranteed analysis			
Acidity / Alkalinity	<0.005 meq/g	Molybdenum (Mo)	<0.1 ppm
Colour	<10 APHA	Nickel (Ni)	<0.1 ppm
Aluminium (Al)	<0.5 ppm	Potassium (K)	<1 ppm
Barium (Ba)	<0.1 ppm	Sodium (Na)	<5 ppm
Cadmium (Cd)	<0.1 ppm	Strontium (Sr)	<0.1 ppm
Calcium (Ca)	<5 ppm	Titanium (Ti)	<0.1 ppm
Chromium (Cr)	<0.1 ppm	n-Butyl formate	<0.0005%
Cobalt (Co)	<0.1 ppm	Total nitrogen (N)	<50 ppm
Copper (Cu)	<0.1 ppm	Total phosphorus (P)	<5 ppm
Iron (Fe)	<0.2 ppm	Total silicon (Si)	<0.5 ppm
Lead (Pb)	<0.1 ppm	Total sulphur (S)	<2 ppm
Lithium (Li)	<0.1 ppm	Vanadium (V)	<0.1 ppm
Magnesium (Mg)	<1 ppm	Zinc (Zn)	<0.5 ppm
Manganese (Mn)	<0.1 ppm		

Filtered to 0.2 µm.

Hydrogen peroxide 30%w/v (100 volumes)
for analysis

Certified AR

packaging		Mfr. No
1 ℓ	HDPE Plastic bottle	H/1800/15
2.5 ℓ	HDPE Plastic bottle	H/1800/17
25 ℓ	HDPE Plastic container	H/1800/25
H ₂ O ₂	EINECS: 231-765-0	
CAS: 7722-84-1	UN: 2014	
MW: 34.01	GHS: Acute Tox. 4, Eye Dam. 1	
Density at 20°C: 1.11	H318, H302	
Boiling point: 108°C	P280, P305+P351+P338,	
Melting point: -33°C	P301+P330+P331, P301+P312	
ADR/RID: 5.1, 8, II		
Assay: 29.5 to 31.5%		

Guaranteed analysis			
Acidity	<0.005 meq/g	Sodium (Na)	<5 ppm
Colour	<10 APHA	Total chloride (Cl)	<0.0005%
Calcium (Ca)	<5 ppm	Total nitrogen (N)	<50 ppm
Copper (Cu)	<0.1 ppm	Total phosphorus (P)	<5 ppm
Iron (Fe)	<0.2 ppm	Total silicon (Si)	<0.5 ppm
Lead (Pb)	<0.1 ppm	Total sulphur (S)	<0.5 ppm
Magnesium (Mg)	<1 ppm	Zinc (Zn)	<0.5 ppm
Potassium (K)	<1 ppm		

Typical analysis			
Aluminium (Al)	<0.5 ppm	Lithium (Li)	<0.1 ppm
Silver (Ag)	<0.1 ppm	Manganese (Mn)	<0.1 ppm
Arsenic (As)	<0.5 ppm	Molybdenum (Mo)	<0.1 ppm
Barium (Ba)	<0.1 ppm	Nickel (Ni)	<0.1 ppm
Cadmium (Cd)	<0.1 ppm	Strontium (Sr)	<0.1 ppm
Chromium (Cr)	<0.1 ppm	Titanium (Ti)	<0.1 ppm
Cobalt (Co)	<0.1 ppm	Vanadium (V)	<0.1 ppm

Lot analysis on the label.

Reagents, Solvents & Volumetric solutions comply with the European Pharmacopeia

Read more on page 234

Hydrogen peroxide 30%w/v (100 volumes)
specified reagents for general laboratory work- extra pure

SLR

packaging		Mfr. No
1 ℓ	HDPE Plastic bottle	H/1750/15
2.5 ℓ	HDPE Plastic bottle	H/1750/17
H ₂ O ₂	EINECS: 231-765-0	
CAS: 7722-84-1	UN: 2014	
MW: 34.01	GHS: Acute Tox. 4, Eye Dam. 1	
Density at 20°C: 1.11	H318, H302	
Boiling point: 108°C	P280, P305+P351+P338,	
Melting point: -33°C	P301+P330+P331, P301+P312	
ADR/RID: 5.1, 8, II		
Assay: 30%		

Guaranteed analysis			
Calcium (Ca)	<20 ppm	Sodium (Na)	<500 ppm
Copper (Cu)	<2 ppm	n-Butyl formate	<0.005%
Iron (Fe)	<5 ppm	Total phosphorus (P)	<500 ppm
Lead (Pb)	<2 ppm	Total silicon (Si)	<2 ppm
Magnesium (Mg)	<5 ppm	Total sulphur (S)	<200 ppm
Potassium (K)	<5 ppm	Zinc (Zn)	<2 ppm

Hydrogen peroxide 12%w/v (40 volumes)
pure

packaging		Mfr. No
2.5 ℓ	HDPE Plastic bottle	H/1650/17
H ₂ O ₂	UN: 2984	
CAS: 7722-84-1	GHS: Ox. Liq. 3, Acute Tox. 4, Skin Corr. 1A	
MW: 34.01	H272, H332, H314	
Density at 20°C: 1.1 to 1.2	P280, P305+P351+P338,	
Boiling point: 107.2°C	P301+P330+P331, P310,	
Melting point: -32.2°C	P304+P340, P210	
ADR/RID: 5.1, III		
EINECS: 231-765-0		

Hydrogen peroxide 6%w/v (20 volumes)
for analysis

Certified AR

packaging		Mfr. No
1 ℓ	HDPE Plastic bottle	H/1600/15
2.5 ℓ	HDPE Plastic bottle	H/1600/17
H ₂ O ₂	ADR/RID: 5.1(8)/II	
CAS: 7722-84-1	EINECS: 231-765-0	
MW: 34.01	GHS: Eye Irrit. 2	
Density at 20°C: 1.04	H319	
Boiling point: 114°C	P280, P305+P351+P338	
Melting point: -52°C		
Assay: >6 %		

Guaranteed analysis			
Acidity	<0.001 meq/g	Potassium (K)	<1 ppm
Colour	<10 APHA	Sodium (Na)	<5 ppm
Calcium (Ca)	<5 ppm	Total chloride (Cl)	<0.0001%
Copper (Cu)	<0.1 ppm	Total phosphorus (P)	<5 ppm
Iron (Fe)	<0.2 ppm	Total silicon (Si)	<0.5 ppm
Lead (Pb)	<0.1 ppm	Total sulphur (S)	<2 ppm
Magnesium (Mg)	<1 ppm	Zinc (Zn)	<0.5 ppm

Typical analysis			
Aluminium (Al)	<0.5 ppm	Lithium (Li)	<0.1 ppm
Silver (Ag)	<0.1 ppm	Manganese (Mn)	<0.1 ppm
Arsenic (As)	<0.1 ppm	Molybdenum (Mo)	<0.1 ppm
Barium (Ba)	<0.1 ppm	Nickel (Ni)	<0.1 ppm
Cadmium (Cd)	<0.1 ppm	Strontium (Sr)	<0.1 ppm
Chromium (Cr)	<0.1 ppm	Titanium (Ti)	<0.1 ppm
Cobalt (Co)	<0.1 ppm	Vanadium (V)	<0.1 ppm

Hydrogen peroxide 6%w/v (20 volumes)
specified reagents for general laboratory work- extra pure

SLR

packaging		Mfr. No
1 ℓ	HDPE Plastic bottle	H/1550/15
2.5 ℓ	HDPE Plastic bottle	H/1550/17
H ₂ O ₂	ADR/RID: 5.1(8)/II	
CAS: 7722-84-1	EINECS: 231-765-0	
MW: 34.01	GHS: Eye Irrit. 2	
Density at 20°C: 1.04	H319	
Boiling point: 114°C	P280, P305+P351+P338	
Melting point: -52°C		

Guaranteed analysis			
Calcium (Ca)	<20 ppm	Sodium (Na)	<100 ppm
Copper (Cu)	<2 ppm	n-Butyl formate	<0.001%
Iron (Fe)	<5 ppm	Total phosphorus (P)	<100 ppm
Lead (Pb)	<2 ppm	Total silicon (Si)	<2 ppm
Magnesium (Mg)	<5 ppm	Total sulphur (S)	<100 ppm
Potassium (K)	<5 ppm	Zinc (Zn)	<2 ppm

Hydroquinone See Quinol, page 268

Hydroxyacetic acid See Glycolic acid, page 137

2-Hydroxybenzaldehyde

See Salicylaldehyde, page 271

Hydroxybenzene See Phenol detached crystals, page 240

2-Hydroxybenzoic acid See Salicylic acid, page 271

3-Hydroxy-4-(1-hydroxy-2-naphthylazo)-7-nitro-1-naphthalene sulfonic acid, sodium salt

See Eriochrome black T, page 119

Label

Fisher Chemical labels show GHS pictograms clearly. Together with actual lot analysis shown on the label, plus a wealth of additional information, we help you to work safely to achieve accurate and consistent results in your analytical applications.



Hydroxylammonium chloride
for analysis

Certified AR

packaging	Mfr. No
250 g Plastic jar	H/1950/50
500 g Plastic jar	H/1950/53
1 kg Plastic jar	H/1950/60

NH₂OH.HCl

CAS: 5470-11-1

MW: 69.49

Density at 20°C: 16.7

Melting point: 155 to 158°C

ADR/RID: 8, 6.1, III

EINECS: 226-798-2

UN: 2923

GHS: Met. Corr. 1, Carc. 2,

Acute Tox. 4, Skin Irrit. 2, Eye Irrit. 2, Skin Sens. 1, STOT RE 2, Aquatic Acute 1

H319, H351, H373, H315, H317, H302, H400, H290, H312

P301+P310, P302+P352, P260, P273, P280, P305+P351+P338

!

Corrosion

!H

Health Hazard

!D

Environment

!N

Nuclear

Assay: >99%

Guaranteed analysis

pH (aqueous solution 5%)	2.5 to 3.5	Magnesium (Mg)	<5 ppm
Ammonium (NH ₄)	<0.1%	Potassium (K)	<20 ppm
Calcium (Ca)	<50 ppm	Sodium (Na)	<50 ppm
Copper (Cu)	<2 ppm	Total silicon (Si)	<5 ppm
Iron (Fe)	<5 ppm	Total sulphur (S)	<20 ppm
Lead (Pb)	<2 ppm	Zinc (Zn)	<5 ppm

Typical analysis

Aluminium (Al)	<2 ppm	Manganese (Mn)	<2 ppm
Barium (Ba)	<2 ppm	Nickel (Ni)	<2 ppm
Cadmium (Cd)	<2 ppm	Strontium (Sr)	<2 ppm
Chromium (Cr)	<2 ppm	Titanium (Ti)	<1 ppm
Cobalt (Co)	<2 ppm	Vanadium (V)	<1 ppm
Lithium (Li)	<2 ppm		

Lot analysis on the label

Hydroxylammonium chloride
specified reagents for general laboratory work- extra pure

SLR

packaging	Mfr. No
100 g Plastic jar	H/1900/48
500 g Plastic jar	H/1900/53
1 kg Plastic jar	H/1900/60

NH₂OH.HCl

CAS: 5470-11-1

MW: 69.49

Density at 20°C: 16.7

Melting point: 155 to 158°C

ADR/RID: 8, 6.1, III

EINECS: 226-798-2

UN: 2923

GHS: Met. Corr. 1, Carc. 2,

Acute Tox. 4, Skin Irrit. 2, Eye Irrit. 2, Skin Sens. 1, STOT RE 2, Aquatic Acute 1

H319, H351, H373, H315, H317, H302, H400, H290, H312

P301+P310, P302+P352, P260, P273, P280, P305+P351+P338

!

Corrosion

!H

Health Hazard

!D

Environment

!N

Nuclear

Assay: >98%

Guaranteed analysis

Calcium (Ca)	<50 ppm	Sodium (Na)	<200 ppm
Copper (Cu)	<10 ppm	Total phosphorus (P)	<100 ppm
Iron (Fe)	<20 ppm	Total silicon (Si)	<20 ppm
Lead (Pb)	<10 ppm	Total sulphur (S)	<100 ppm
Magnesium (Mg)	<20 ppm	Zinc (Zn)	<20 ppm
Potassium (K)	<100 ppm		

UHPLC

Filtered to 0.1 µm, Fisher Chemical UHPLC gradient grade solvents (Acetonitrile, Methanol and Water) are perfect for those working at high pressure.

See page 309 for more information

Hydroxylammonium sulfate
for analysis

Certified AR

packaging	Mfr. No
100 g Plastic jar	H/2001/48

(NH₂OH)₂.H₂SO₄

CAS: 10039-54-0

MW: 164.13

Melting point: 170°C

ADR/RID: 8, III

EINECS: 233-118-8

UN: 2865

GHS: Met. Corr. 1, Carc. 2,

Acute Tox. 4, Skin Irrit. 2, Eye Irrit. 2, Skin Sens. 1, STOT RE 2, Aquatic Acute 1

H351, H317, H373, H319, H315, H302, H400, H290, H312

P280, P361, P301+P312, P260, P273, P302+P352, P305+P351+P338

!

Corrosion

!H

Health Hazard

!D

Environment

!N

Nuclear

Assay: >99%

Guaranteed analysis

pH (aqueous solution 5%)	3 to 4	Potassium (K)	<10 ppm
Calcium (Ca)	<10 ppm	Sodium (Na)	<20 ppm
Copper (Cu)	<2 ppm	n-Butyl formate	<0.001%
Iron (Fe)	<5 ppm	Total phosphorus (P)	<20 ppm
Lead (Pb)	<2 ppm	Total silicon (Si)	<5 ppm
Magnesium (Mg)	<5 ppm	Zinc (Zn)	<5 ppm

N-2-Hydroxyethylpiperazine-N'-2-ethanesulfonic acid
for biochemistry, HEPES

packaging	Mfr. No
100 g Plastic jar	H/2230/48

C₈H₁₈N₂O₄S

CAS: 7365-45-9

MW: 259.18

Melting point: 234 to 238°C

EINECS: 230-907-9

Assay: >99%

Guaranteed analysis

Water	<0.5%	Sulfate (SO ₄)	<0.01%
pKa at 37°C	: 7.3 to 7.32	n-Butyl formate	<0.005%
Solution	: Pass	Total phosphorus (P)	<10 ppm
Iron (Fe)	<5 ppm	Residue after ignition	<0.1%
Lead (Pb)	<10 ppm		

Absorbance specification

Max. Absorbance (1cm cell v HPLC water)

250 nm (sol. aqu. 10%) : <0.1 A.U.

1-(1-Hydroxy-4-methyl-2-phenylazo)-2-naphthol-4-sulfonic acid
pure, metal indicator (Calmagite)

packaging	Mfr. No
5 g Glass jar	H/2350/44

C₁₇H₁₄N₂O₅S

CAS: 3147-14-6

MW: 358.37

UN: 2585

GHS: Skin Corr. 1B, Eye Dam. 1

H314

P301+P330+P331, P280, P305+P351+P338, P310

Assay: >99%

Guaranteed analysis

Calcium (Ca)	<50 ppm	Sodium (Na)	<200 ppm
Copper (Cu)	<10 ppm	Total phosphorus (P)	<100 ppm
Iron (Fe)	<20 ppm	Total silicon (Si)	<20 ppm
Lead (Pb)	<10 ppm	Total sulphur (S)	<100 ppm
Magnesium (Mg)	<20 ppm	Zinc (Zn)	<20 ppm
Potassium (K)	<100 ppm		

1-Hydroxynaphthalene

See 1-Naphthol, page 210

4-[(4-Hydroxy-1-naphthalenyl)phenylmethylene]-1(4H)-naphthalenone

See p-Naphtholbenzein, page 211

7-Hydroxy-8-(phenylazo)-1.3-naphthalenedi-sulfonic acid, disodium salt

See Orange G sodium salt, page 221

2-Hydroxy-1.2,3-propanetricarboxylic acid

See Citric acid, page 90

L-Hydroxyproline
for biochemistry

packaging	Mfr. No
5 g Plastic jar	H/2415/44

C₅H₉NO₃

CAS: 51-35-4

MW: 131.13

Melting point: 274 to 275°C

EINECS: 200-091-9

Assay: >99%

Guaranteed analysis

Loss on drying	<0.5%	Lead (Pb)	<10 ppm
Specific rotation	-77° to -75°	Residue after ignition	<0.1%
Iron (Fe)	<10 ppm		

8-Hydroxyquinoline

specified reagents for general laboratory work- extra pure

SLR

packaging	Mfr. No
100 g Glass jar	H/2500/48

C₈H₇NO

CAS: 148-24-3

MW: 145.16

GHS: Muta. 2, Acute Tox. 4, Skin Irrit. 2, Eye Irrit. 2, STOT SE 3

H319, H315, H335, H332, H312, H302, H341

P281, P260, P301+P312, P304+P340, P302+P352, P280, P305+P351+P338

Assay: >98%

Guaranteed analysis

Calcium (Ca)	<50 ppm	Sodium (Na)	<200 ppm
Copper (Cu)	<20 ppm	n-Butyl formate	<0.01%
Iron (Fe)	<50 ppm	Total phosphorus (P)	<100 ppm
Lead (Pb)	<20 ppm	Total silicon (Si)	<20 ppm
Magnesium (Mg)	<10 ppm	Total sulphur (S)	<200 ppm
Potassium (K)	<100 ppm	Zinc (Zn)	<20 ppm

DL-Hydroxysuccinic acid

See DL-Malic acid, page 190

3-Hydroxytoluene

See m-Cresol, page 98

2-Hydroxytoluene

See o-Cresol 98-100 %, page 98

Hypophosphorous acid

specified reagents for general laboratory work- extra pure

SLR

packaging	Mfr. No
500 mL HDPE Plastic bottle	H/2650/PB08

H₃PO₂

CAS: 6303-21-5

MW: 66.00

EINECS: 228-601-5

UN: 3264

GHS: Skin Corr. 1B, Eye Dam. 1

H314

P301+P330+P331, P280, P305+P351+P338, P310

Assay: >49%

Guaranteed analysis

Lead (Pb)	<20 ppm	n-Butyl formate	<0.1%
Sulfate (SO ₄)	<0.1%		

Indol-3-ylacetic acid
pure, 96+%

packaging	Mfr. No
1 g Vial	I/0370/43

C₁₀H₉NO₂

CAS: 87-51-4

MW: 175.07

STOT SE 3

H319, H315, H335

P261, P302+P352, P280, P305+P351+P338

Guaranteed analysis

Identity test : Pass

Imidazole

specified reagents for general laboratory work- extra pure

SLR

packaging	Mfr. No
100 g Plastic jar	I/0012/48
1 kg Plastic jar	I/0012/60
10 kg Plastic keg	I/0012/68

C₃H₄N₂

CAS: 288-32-4

MW: 68.08

UN: 2923

GHS: Repr. 1B, Acute Tox. 4, Skin Corr. 1B, Eye Dam. 1

H302, H314, H360D

P301+P330+P331, P280, P305+P351+P338, P310, P281, P301+P312

Assay: >99.5%

Guaranteed analysis

Copper (Cu)	<5 ppm	Lead (Pb)	<5 ppm
Iron (Fe)	<5 ppm	Mercury (Hg)	<2 ppm

Absorbance specification

Max. Absorbance (1cm cell v HPLC water)

325 nm (sol. aqu. 12%) : <0.15 A.U.

Imidazole

pure

packaging	Mfr. No
100 g Plastic jar	I/0010/48
500 g Plastic jar	I/0010/53

C₃H₄N₂

CAS: 288-32-4

MW: 68.08

UN: 2923

GHS: Repr. 1B, Acute Tox. 4, Skin Corr. 1B, Eye Dam. 1

H302, H314, H360D

P301+P330+P331, P280, P305+P351+P338, P310, P281, P301+P312

Assay: >99%

2.2'-Iminodiethanol

See Diethanolamine, page 108

High Purity Acids
for Elemental Analysis

Fisher Chemical offers a complete range of acids to meet the most challenging applications. Choose from our Optima®, TraceMetal® and PrimarPlus® grades to ensure that the sensitivity you require is always available.

Read more page 24

pH Indicator

designation	pH range	colour change	packsize		Mfr. No	page
Alizarin	11.0-13.0	pink to violet	10 g		A/1460/45	27
Alizarin red S, C.I.58005	4.0-6.0	yellow to orange-red	25 g		A/1520/46	27
Bromocresol green	3.6-5.2	yellow to blue	1 g		B/4320/43	61
	3.6-5.2	yellow to blue	5 g		B/4320/44	61
	3.6-5.2	yellow to blue	25 g		B/4320/46	61
Bromocresol purple	5.2-6.8	yellow to violet	5 g		B/4330/44	61
	5.2-6.8	yellow to violet	25 g		B/4330/46	61
Bromophenol blue	2.8-4.6	yellow to blue-violet	5 g		B/4630/44	62
	2.8-4.6	yellow to blue-violet	25 g		B/4630/46	62
Full range indicator, pH range 1-13	1.0-13.0	colour chart on bottle	100 mL		U/0025/PB05	67
	1.0-13.0	colour chart on bottle	500 mL		U/0025/PB08	67
Universal indicator, pH range 4-10	4.0-10.0	colour chart on bottle	100 mL		U/0020/PB05	67
	4.0-10.0	colour chart on bottle	500 mL		U/0020/PB08	67
	4.0-10.0	colour chart on bottle	2.5 L		U/0020/PB17	67
Congo red	3.0-5.0	violet to red-orange	25 g		C/7020/46	94
Cresol red	0.2-1.8	red to yellow	5 g		C/7460/44	98
Litmus	5.0-8.0	red to blue	25 g		L/2500/46	185
Methyl orange solid C.I. 13025	2.9-4.6	orange red to orange yellow	25 g		M/5000/46	206
	2.9-4.6	orange red to orange yellow	100 g		M/5000/48	206
Methyl red spirit soluble C.I. 13020	4.2-6.3	red to yellow	10 g		M/5150/45	207
	4.2-6.3	red to yellow	25 g		M/5150/46	207
Methyl red solid water soluble	4.2-6.3	red to yellow	5 g		M/5170/44	207
	4.2-6.3	red to yellow	100 g		M/5170/48	207
Methyl violet solid (crystal violet) C.I. 42555	0.1-3.2	yellow to violet	10 g		M/5370/45	208
	0.1-3.2	yellow to violet	25 g		M/5370/46	208
p-Naphtholbenzein	8.5-9.8	yellow to green	25 g		N/0500/46	211
Neutral red C.I 50040	6.8-8.0	red to yellow-orange	5 g		N/1270/44	212
Thymol blue	1.2-2.8	red to yellow	5 g		T/1470/44	349
Thymolphthalein solid	9.3-10.5	colourless to blue	5 g		T/1480/44	349
	9.3-10.5	colourless to blue	25 g		T/1480/46	349
Titan yellow C.I. 19540	12.0-13.0	yellow to red	25 g		T/2060/46	353

Redox Indicator

designation	transition	colour change (oxidised to reduced state)	packsize	package type		Mfr. No	page
Barium diphenylamine-4-sulfonate	+ 0.84	violet to colourless	5 g	Plastic jar		B/0620/44	54
Diphenylamine	+ 0.76	violet to colourless	100 g	Glass jar		D/4600/48	117
Methylene blue solid C.I. 52015	+ 0.53	violet to colourless	10 g	Glass jar		M/4870/45	205
	+ 0.87	violet to colourless	25 g	Glass jar		M/4870/46	205
1.10-Phenanthroline mono-hydrate	+ 1.08	faint blue to red	1 g	Vial		P/2200/43	239
	+ 1.08	faint blue to red	5 g	Glass bottle		P/2200/44	239
	+ 0.67	red to colourless	1 g	Glass jar		P/2390/43	240
	+ 0.67	red to colourless	5 g	Glass jar		P/2390/44	240
Phenosafranine	+ 0.28	blue to colourless	5 g	Glass jar		P/2450/44	241
N-Phenylanthranilic acid	+ 1.08	purple-red to colourless	5 g	Plastic jar		P/2460/44	241
Sodium diphenylamine-4-sulfonate	+ 0.84	red-violet to colourless	5 g	Glass jar		S/3780/44	289

Adsorption Indicator

designation	colour in-dex	packsize	package type		Mfr. No	page
Alizarin	C.I. 58000	10 g	Glass jar		A/1460/45	27
Bromocresol purple		5 g	Glass jar		B/4330/44	61
		25 g	Glass jar		B/4330/46	61
Bromophenol blue		5 g	Glass jar		B/4630/44	62
		25 g	Glass jar		B/4630/46	62
Congo red	C.I. 22120	25 g	Glass jar		C/7020/46	94
Dichloro (R) fluorescein		5 g	Glass jar		D/1820/44	105
1.5-Diphenylcarbazone		10 g	Glass jar		D/4850/45	117

Adsorption Indicator

designation	colour in-dex	packsize	package type		Mfr. No	page
Fluorescein sodium		100 g	Glass jar		F/1300/48	130
		500 g	Glass jar		F/1300/53	130
Mercury(II) thiocyanate		25 g	Plastic jar		M/3300/46	198
Phenosafranine	C.I. 50200	5 g	Glass jar		P/2450/44	241
Tetrabromo (R) fluorescein	C.I. 45380	10 g	Glass jar		T/0470/45	343
		50 g	Glass jar		T/0470/47	343

Indigo carmine
extra pure, C.I. 73015

SLR

packaging		Mfr. No
25 g	Glass jar	I/0300/46
C ₁₆ H ₈ N ₂ Na ₂ O ₈ S ₂	GHS: Acute Tox. 4	
CAS: 860-22-0	H302	
MW: 466.36	P301+P312	
EINECS: 212-728-8		
Guaranteed analysis		
Complies with standard IR spectrum: Pass		
Conforming to the IR spectrum of reference.		

5.5'-Indigodisulfonic acid, disodium salt

See Indigo carmine, page 160

Indium standard metal solution 1000 ppm

for ICP solution in HNO₃ 2%, NIST certified

packaging		Mfr. No
100 mL	Plastic bottle	SPPLIN22Y
500 mL	Plastic bottle	SPPLIN22X
H315, H319		

Manufactured to a tolerance of ± 0.5%.

All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.

Check with your local distributor for availability of Spex product

Indium standard metal solution 10 ppm

for ICP-MS solution in HNO₃ 2%, NIST certified

packaging		Mfr. No
100 mL	Plastic bottle	SPCLIN21AY
H315, H319		

Manufactured to a tolerance of ± 0.5%.

All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.

Check with your local distributor for availability of Spex product

myo-Inositol

pure

packaging		Mfr. No
100 g	Plastic jar	I/0385/48
C ₆ H ₁₂ O ₆	Melting point: 220 to 228°C	
CAS: 87-89-8	EINECS: 201-781-2	
MW: 180.16		

Invertase concentrate yeast

technical, Stabilized with glycerol 10%

packaging		Mfr. No
100 mL	HDPE Plastic bottle	I/0395/05
GHS: Resp. Sens. 1		
H334		
P304+P341, P260, P280, P285		
Guaranteed analysis		
Activity	11020 to 13340 U/mL	

Iodine

for analysis

Certified AR

packaging		Mfr. No
100 g	Glass bottle	I/0500/48
250 g	Glass bottle	I/0500/50
500 g	Glass bottle	I/0500/53
I ₂	UN: 3495	
CAS: 7553-56-2	GHS: Acute Tox. 4, Aquatic	
MW: 253.81	Acute 1	
Boiling point: 185°C	H332, H312, H400	
Melting point: 113°C	P280, P302+P352,	
ADR/RID: 8, 6.1, III	P273, P304+P340	
EINECS: 231-442-4		
Assay: >99.9%		

Guaranteed analysis

Insoluble matter	Pass	Magnesium (Mg)	<2 ppm
Calcium (Ca)	<5 ppm	Potassium (K)	<20 ppm
Chloride and Bromide	<0.005%	Sodium (Na)	<50 ppm
Copper (Cu)	<5 ppm	Total sulphur (S)	<50 ppm
Iron (Fe)	<10 ppm	Zinc (Zn)	<2 ppm
Lead (Pb)	<5 ppm		

Typical analysis

Aluminium (Al)	<1 ppm	Manganese (Mn)	<1 ppm
Silver (Ag)	<1 ppm	Molybdenum (Mo)	<1 ppm
Barium (Ba)	<1 ppm	Nickel (Ni)	<1 ppm
Cadmium (Cd)	<1 ppm	Strontium (Sr)	<1 ppm
Chromium (Cr)	<1 ppm	Titanium (Ti)	<1 ppm
Cobalt (Co)	<1 ppm	Vanadium (V)	<1 ppm
Lithium (Li)	<1 ppm		

Lot analysis on the label.

Iodine

specified reagents for general laboratory work- extra pure

SLR

packaging		Mfr. No
100 g	Glass bottle	I/0450/48
250 g	Glass bottle	I/0450/50
500 g	Glass bottle	I/0450/53
25 kg	Plastic keg	I/0450/70
I ₂	UN: 3495	
CAS: 7553-56-2	GHS: Acute Tox. 4, Aquatic	
MW: 253.81	Acute 1	
Boiling point: 185°C	H332, H312, H400	
Melting point: 113°C	P280, P302+P352,	
ADR/RID: 8, 6.1, III	P273, P304+P340	
EINECS: 231-442-4		
Assay: >99.5%		

Guaranteed analysis

Calcium (Ca)	<20 ppm	Magnesium (Mg)	<10 ppm
Chloride and Bromide	<0.02%	Potassium (K)	<100 ppm
Copper (Cu)	<20 ppm	Sodium (Na)	<200 ppm
Iron (Fe)	<50 ppm	Total sulphur (S)	<200 ppm
Lead (Pb)	<20 ppm	Zinc (Zn)	<10 ppm

Iodine indicator

pure, indicator grade

packaging		Mfr. No
250 g	Plastic jar	I/0505/50
1 kg	Plastic jar	I/0505/60
I ₂	ADR/RID: 8(6.1)/II	
CAS: 9005-25-8	EINECS: 232-679-6	
Density at 20°C: 1.5	UN: 2923	
Guaranteed analysis		
Identity test : Pass		
Iodine indicator is a stable white powder which is readily soluble in water. It forms a blue complex with iodine: 0.2-0.5g being the usual amount per titration.		
Sensitivity to iodine: 0.05 mL of 0.01N iodine gives a distinct blue colour in 100mL solution.		
Glassware is easily cleaned after use.		

Iodine solution - Wij's

pure, indicator grade

packaging		Mfr. No
1 L	Coated Glass bottle	J/4500/PB15
2.5 L	Coated Glass bottle	J/4500/PB17
CAS: 7782-50-5	GHS: Flam. Liq. 3, Skin Corr.	
Density at 20°C: 1.06	1A, Eye Dam. 1	
Boiling point: 117.8°C	H226, H314	
Melting point: 17.2°C	P280, P301+P330+P331,	
Flash point: 41.1°C	P305+P351+P338, P310,	
ADR/RID: 8, 3, II	P304+P340, P302+P352,	
EINECS: 200-580-7	P273, P210	
UN: 2920		

Formulation free of carbon tetrachloride.

Iodine solution 1% w/v in potassium iodide

pure

packaging		Mfr. No
1 L	Glass bottle	J/4450/15
2.5 L	Glass bottle	J/4450/17
GHS: Aquatic Chronic 3		
H412		
P273		

Iodine solution 0.05M (0.1N)

Eur.Ph. , USP, BP

NIST Standard solution ready to use, For Volumetric analysis

packaging		Mfr. No
1 L	Glass bottle	J/4410/15
2.5 L	Glass bottle	J/4410/17
CAS: 7553-56-2	EINECS: 231-442-4	
Density at 20°C: 1.02	GHS: Aquatic Chronic 3	
Boiling point: 100°C	H412	
Melting point: 0°C	P273	

Manufactured to a tolerance of ± 0.1%.

See the full range of Standard Volumetric solution page 323.

Iodine solution 0.5M (1N)

Eur.Ph. , BP

NIST Standard solution ready to use, For Volumetric analysis

packaging		Mfr. No
500 mL	Coated Glass bottle	J/4400/PB08
1 L	Coated Glass bottle	J/4400/PB15
2.5 L	Coated Glass bottle	J/4400/PB17
CAS: 7553-56-2	GHS: Aquatic Chronic 2	
Density at 20°C: 1.4	H411	
EINECS: 231-442-4	P273	

Manufactured to a tolerance of ± 0.2%.

See the full range of Standard Volumetric solution page 323.

Iodine bromide

pure

packaging		Mfr. No
100 g	Glass bottle	I/0510/48
Brl	EINECS: 232-159-9	
CAS: 7789-33-5	UN: 3260	
MW: 206.81	GHS: Skin Corr. 1B, Eye	
Density at 20°C: 4.41	Dam. 1, STOT SE 3	
Boiling point: 116°C	H335, H314	
Melting point: 42 to 50°C	P301+P330+P331, P280,	
ADR/RID: 8, II	P305+P351+P338, P310, P261	

Guaranteed analysis

Identity test : Pass

Iodine trichloride

specified reagents for general laboratory work- extra pure

SLR

packaging		Mfr. No
100 g	Coated Glass bottle	I/0650/48
ICl ₃	GHS: Skin Corr. 1B, Eye	
CAS: 865-44-1	Dam. 1, STOT SE 3	
MW: 233.26	H314, H335	
Melting point: 63°C	P305+P351+P338, P280,	
ADR/RID: 8, II	P260, P310, P264,	
EINECS: 212-739-8	P301+P330+P331	
UN: 3260		
Assay: >96%		

1-Iodobutane

specified reagents for general laboratory work- extra pure

SLR

packaging		Mfr. No
25 mL	Glass bottle	I/0697/04
C ₄ H ₉ I	Tox. 4, Skin Irrit. 2, Eye Irrit.	
CAS: 542-69-8	2, STOT SE 3	
MW: 184.02	H226, H335, H315, H319,	
Boiling point: 130 to 131°C	H302, H332	
Melting point: -103°C	P261, P301+P312,	
ADR/RID: 3, III	P302+P352, P280,	
EINECS: 208-824-4	P305+P351+P338, P210,	
UN: 1993	P304+P341	
GHS: Flam. Liq. 3, Acute		
Assay (CG): >98%		
Absorbance specification		
Max. Absorbance (1cm cell v HPLC water)		
Refractive index		1.49 to 1.51

Iodoethane

specified reagents for general laboratory work- extra pure

SLR

packaging		Mfr. No
100 mL	Coated Glass bottle	I/0700/PB05
C ₂ H ₅ I	UN: 1993	
CAS: 75-03-6	GHS: Flam. Liq. 3, Acute	
MW: 155.97	Tox. 4, Skin Irrit. 2, Eye Irrit.	
Density at 20°C: 1.95	2, STOT SE 3	
Boiling point: 69 to 73°C	H226, H319, H335, H315,	
Melting point: -108°C	H332, H302	
Flash point: 53°C	P261, P301+P312,	
ADR/RID: 3, III	P302+P352, P280,	
EINECS: 200-833-1	P305+P351+P338, P210	

Iodoform

specified reagents for general laboratory work- extra pure

SLR

packaging		Mfr. No
100 g	Plastic jar	I/0750/48
CHI ₃	GHS: Acute Tox. 4, Skin Irrit.	
CAS: 75-47-8	2, Eye Irrit. 2, STOT SE 3	
MW: 393.73	H319, H335, H332, H302,	
Melting point: 119 to 122°C	H315, H312	
Flash point: 204°C	P280, P302+P352, P261,	
ADR/RID: 6.1, III	P301+P312, P304+P340,	
EINECS: 200-874-5	P305+P351+P338	
UN: 2811		

High Volume Solvent
Delivery Systems

High-volume solvent delivery systems are available in 10L to 1000L, in Europe only. Please email your contact details and a brief description of your requirements to begeel.sha.quotes@thermofisher.com



Iodomethane		SLR
specified reagents for general laboratory work- extra pure		
packaging	Mfr. No	
100 mℓ Coated Glass bottle	I/0800/PB05	
CH ₃ I	UN: 2644	
CAS: 74-88-4	GHS: Carc. 2, Acute Tox. 3, Skin Irrit. 2, STOT SE 3	
MW: 141.94	H315, H335, H331, H301, H312, H351	
Density at 20°C: 2.28	P301+P310, P280, P304+P340, P302+P352, P261	
Boiling point: 42.5°C		
Melting point: -66°C		
ADR/RID: 6.1, I		
EINECS: 200-819-5		
Assay (CG): >98%		

Iron standard metal solution 1000 ppm		
for ICP, solution HNO ₃ 1M		
packaging	Mfr. No	
500 mℓ HDPE Plastic bottle	J/8271/08	
ADR/RID: 8, II	H314	
EINECS: 233-899-5	P280, P305+P351+P338, P301+P330+P331, P310	
UN: 2031		
GHS: Skin Corr. 1B		
Manufactured to a tolerance of ± 0.5%. All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities. See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.		

2-Iodopropane		SLR
specified reagents for general laboratory work- extra pure		
packaging	Mfr. No	
100 mℓ Coated Glass bottle	I/0831/PB05	
C ₃ H ₇ I	EINECS: 200-859-3	
CAS: 75-30-9	UN: 2392	
MW: 169.99	GHS: Flam. Liq. 3, Skin Irrit. 2, Eye Irrit. 2, STOT SE 3	
Density at 20°C: 1.7	H226, H315, H319, H335	
Boiling point: 90 to 90°C	P261, P302+P352, P280, P305+P351+P338, P210	
Melting point: -90°C		
Flash point: 42°C		
ADR/RID: 3, III		
Assay (CG): >98%		

Iron standard metal solution 10 000 ppm		
for ICP, solution HNO ₃ 1M		
packaging	Mfr. No	
100 mℓ HDPE Plastic bottle	J/8270/05	
ADR/RID: 8, II	P301+P330+P331, P280, P305+P351+P338, P310, P303+P361+P353	
UN: 2031		
GHS: Skin Corr. 1B		
H314		
Manufactured to a tolerance of ± 0.5%. All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities. See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.		

Iridium standard metal solution 1000 ppm		
for ICP solution in HCL 10%, NIST certified		
packaging	Mfr. No	
100 mℓ Plastic bottle	SPPLIR32Y	
500 mℓ Plastic bottle	SPPLIR32X	
H315, H319, H335		

Iron standard metal solution 1000 ppm		
for ICP solution in HNO ₃ 2%, NIST certified		
packaging	Mfr. No	
100 mℓ Plastic bottle	SPPLFE22Y	
250 mℓ Plastic bottle	SPPLFE22T	
500 mℓ Plastic bottle	SPPLFE22X	
H315, H319		

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Iron standard metal solution 1000 ppm		
for AAS, solution HNO ₃ 1M		
packaging	Mfr. No	
100 mℓ HDPE Plastic bottle	J/8030/05	
500 mℓ HDPE Plastic bottle	J/8030/08	
ADR/RID: 8, II	H314	
EINECS: 233-899-5	P280, P305+P351+P338, P301+P330+P331, P310	
UN: 2031		
GHS: Skin Corr. 1B		
Manufactured to a tolerance of ± 0.5%. See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.		

Iron standard metal solution 1000 ppm		
for ICP solution in HCL 2%, NIST certified		
packaging	Mfr. No	
500 mℓ Plastic bottle	SPPLFE12X	
H314		
Manufactured to a tolerance of ± 0.5%. All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities. See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314. Check with your local distributor for availability of Spex product		

Iron standard metal solution 10 000 ppm		
for ICP solution in HNO ₃ 5%, NIST certified		
packaging	Mfr. No	
100 mℓ Plastic bottle	SPPLFE23Y	
500 mℓ Plastic bottle	SPPLFE23X	
H303, H314		
Manufactured to a tolerance of ± 0.5%. All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities. See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314. Check with your local distributor for availability of Spex product		

Iron standard metal solution 10 000 ppm		
for ICP solution in HCL 5%, NIST certified		
packaging	Mfr. No	
500 mℓ Plastic bottle	SPPLFE13X	
H303		
Manufactured to a tolerance of ± 0.5%. All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities. See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314. Check with your local distributor for availability of Spex product		

Iron standard metal solution 1000 ppm		
for ICP-MS solution in HNO ₃ 2%, NIST certified		
packaging	Mfr. No	
100 mℓ Plastic bottle	SPCLFE22Y	
H315, H319		

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Iron powder		SLR
extra pure, metal reduced by hydrogen		
packaging	Mfr. No	
500 g Plastic jar	I/0950/53	
Fe	EINECS: 231-096-4	
CAS: 7439-89-6	UN: 3089	
MW: 55.85	GHS: Flam. Sol. 2	
Boiling point: 3000°C	H228	
Melting point: 1535°C	P210	
ADR/RID: 4.1, III		
Guaranteed analysis		
Identity test	: Pass	

Iron filings (coarse)		
pure		
packaging	Mfr. No	
1 kg Plastic jar	I/0855/60	
25 kg Plastic keg	I/0855/70	
Fe	Boiling point: 3000°C	
CAS: 7439-89-6	Melting point: 1535°C	
MW: 55.84	EINECS: 231-096-4	
Density at 20°C: 7.86		
Guaranteed analysis		
Identity test	: Pass	

Iron filings (fine) metal		
pure		
packaging	Mfr. No	
250 g Plastic jar	I/0850/50	
1 kg Plastic jar	I/0850/60	
3 kg Plastic jar	I/0850/63	
Fe	Boiling point: 3000°C	
CAS: 7439-89-6	Melting point: 1535°C	
MW: 55.84	EINECS: 231-096-4	
Density at 20°C: 7.86		
Guaranteed analysis		
Identity test	: Pass	

Iron metal		
technical, wire diameter ca. 0.3 mm		
packaging	Mfr. No	
100 g Plastic jar	I/1001/48	
Fe	Density at 20°C: 7.86	
CAS: 7439-89-6	Boiling point: 3000°C	
MW: 55.84	Melting point: 1535°C	
Guaranteed analysis		
Identity test	: Pass	

Iron metal strip		
technical, Strip 200 x 15 x 0.5mm		
packaging	Mfr. No	
Plastic jar	I/1021/99	
Fe	Boiling point: 3000°C	
CAS: 7439-89-6	Melting point: 1535°C	
MW: 55.84	EINECS: 231-096-4	
Density at 20°C: 7.86		
Guaranteed analysis		
Identity test	: Pass	

Iron(III) chloride hexahydrate			SLR
specified reagents for general laboratory work- extra pure			
packaging			Mfr. No
1 kg	Plastic jar		I/1020/60
5 kg	Plastic keg		I/1020/65
Cl ₃ Fe.6H ₂ O		GHS: Acute Tox. 4, Skin Irrit.	
CAS: 10025-77-1		2, Eye Dam. 1	
MW: 270.29		H318, H302, H315	
Density at 20°C: 1.82		P280, P302+P352,	
Boiling point: 280 to 285°C		P305+P351+P338, P313	
Melting point: 37°C			
Guaranteed analysis			
Calcium (Ca)	<500 ppm	Sodium (Na)	<1000 ppm
Copper (Cu)	<100 ppm	Total phosphorus (P)	<100 ppm
Iron (II)	<0.1%	Total silicon (Si)	<200 ppm
Lead (Pb)	<100 ppm	Total sulphur (S)	<500 ppm
Magnesium (Mg)	<200 ppm	Zinc (Zn)	<500 ppm
Potassium (K)	<200 ppm		

Standard Volumetric solution
For your convenience we offer a range of NIST traceable solutions. A full list of our products can be found on page 323. Can't see what you need? Please email your contact details and a brief description of your requirements to begeel.sha.quotes@thermofisher.com

Optima LC/MS

Fisher Chemical OPTIMA® LC/MS solvents, blends, additives and modifiers have set the standard of excellence for consistent, reproducible performance in the mobile phase of LC/MS.
Find out more on page 308

We Support Your Chemistry For

- Semi-bulk and bulk requirements
 - Tailored solvents and solvent blends or customized solutions
 - Testing services or customized packaging and labeling
- Please email your contact details and a brief description of your requirements to begeel.sha.quotes@thermofisher.com

Iron(III) chloride anhydrous		SLR
<i>specified reagents for general laboratory work- extra pure</i>		
packaging	Mfr. No	
250 g Plastic jar	I/1035/50	
500 g Plastic jar	I/1035/53	
3 kg Plastic jar	I/1035/63	
FeCl ₃ CAS: 7705-08-0 MW: 162.21 Density at 20°C: 2.9 Boiling point: 316°C Melting point: 300°C ADR/RID: 8, III EINECS: 231-729-4	UN: 1773 GHS: Acute Tox. 4, Skin Irrit. 2, Eye Dam. 1, Skin Sens. 1 H302, H315, H318, H317 P305+P351+P338, P261, P280, P302+P352, P301+P312, P333+P313	 
Assay: >97%		
Guaranteed analysis		
Iron (III)	<1%	

Iron(III) chloride	
<i>pure, solution about 60% w/v, FeCl₃</i>	
packaging	Mfr. No
2.5 l HDPE Plastic bottle	I/1045/PB17
FeCl ₃ CAS: 7705-08-0 MW: 162.21 ADR/RID: 8, III EINECS: 231-729-4 UN: 2582	GHS: Acute Tox. 4, Skin Irrit. 2, Eye Dam. 1, Skin Sens. 1 H302, H315, H318, H317 P305+P351+P338, P261, P280, P302+P352, P301+P312, P333+P313

Iron(III) chloride solution about 60% w/v		SLR
<i>specified reagents for general laboratory work- extra pure</i>		
packaging	Mfr. No	
500 ml Coated Glass bottle	I/1047/PB08	
FeCl ₃ CAS: 7705-08-0 MW: 162.21 ADR/RID: 8, III EINECS: 231-729-4 UN: 2582	GHS: Acute Tox. 4, Skin Irrit. 2, Eye Dam. 1, Skin Sens. 1 H302, H315, H318, H317 P305+P351+P338, P261, P280, P302+P352, P301+P312, P333+P313	 

Assay: 59 to 61%		
Guaranteed analysis		
Calcium (Ca)	<500 ppm	Total phosphorus (P)
Copper (Cu)	<100 ppm	Total silicon (Si)
Lead (Pb)	<100 ppm	Total sulphur (S)
Magnesium (Mg)	<200 ppm	Zinc (Zn)
Potassium (K)	<200 ppm	Iron (II)
Sodium (Na)	<1000 ppm	<0.1%

Iron(III) citrate trihydrate	
<i>pure</i>	
packaging	Mfr. No
500 g Plastic jar	I/1060/53
C ₆ H ₅ FeO ₇ ·3H ₂ O CAS: 17217-76-4 MW: 299.00 EINECS: 222-536-6	
Guaranteed analysis	
Identity test	



Iron(II) disulfide		SLR
<i>extra pure, powder</i>		
packaging	Mfr. No	
1 kg Glass jar	I/1165/60	
25 kg Fibre Keg	I/1165/70	
FeS ₂ CAS: 12068-85-8 MW: 119.97	EINECS: 235-106-8 EUH032 P233	
Guaranteed analysis		
Identity test		

Iron(III) nitrate nonahydrate		Certified AR
<i>for analysis</i>		
packaging	Mfr. No	
250 g Plastic jar	I/1075/50	
500 g Plastic jar	I/1075/53	
1 kg Plastic jar	I/1075/60	
Fe(NO ₃) ₃ ·9H ₂ O CAS: 7782-61-8 MW: 404.00 Density at 20°C: 1.68 Melting point: 47°C ADR/RID: 5.1, III	UN: 1466 GHS: Ox. Sol. 3, Skin Irrit. 2, Eye Irrit. 2, STOT SE 3 H272, H319, H315, H335 P261, P302+P352, P280, P305+P351+P338, P210	 

Assay: 98 to 102%		
Guaranteed analysis		
Acidity	<0.1 meq/g	Sodium (Na)
Calcium (Ca)	<50 ppm	Total chloride (Cl)
Copper (Cu)	<50 ppm	Total phosphorus (P)
Lead (Pb)	<20 ppm	Total silicon (Si)
Magnesium (Mg)	<10 ppm	Total sulphur (S)
Potassium (K)	<100 ppm	Zinc (Zn)
Typical analysis		
Silver (Ag)	<5 ppm	Lithium (Li)
Arsenic (As)	<5 ppm	Manganese (Mn)
Barium (Ba)	<5 ppm	Molybdenum (Mo)
Cadmium (Cd)	<0.005 %	Strontium (Sr)
Chromium (Cr)	<0.005 %	Titanium (Ti)
Cobalt (Co)	<0.001 %	Vanadium (V)

Lot analysis on the label.

Iron(III) nitrate nonahydrate		SLR
<i>specified reagents for general laboratory work- extra pure</i>		
packaging	Mfr. No	
250 g Plastic jar	I/1070/50	
500 g Plastic jar	I/1070/53	
1 kg Plastic jar	I/1070/60	
Fe(NO ₃) ₃ ·9H ₂ O CAS: 7782-61-8 MW: 404.00 Density at 20°C: 1.68 Melting point: 47°C ADR/RID: 5.1, III	UN: 1466 GHS: Ox. Sol. 3, Skin Irrit. 2, Eye Irrit. 2, STOT SE 3 H272, H319, H315, H335 P261, P302+P352, P280, P305+P351+P338, P210	 

Assay: >98%		
Guaranteed analysis		
Calcium (Ca)	<500 ppm	n-Butyl formate
Copper (Cu)	<100 ppm	Total phosphorus (P)
Lead (Pb)	<100 ppm	Total silicon (Si)
Magnesium (Mg)	<200 ppm	Total sulphur (S)
Potassium (K)	<200 ppm	Zinc (Zn)
Sodium (Na)	<1000 ppm	<500 ppm



Iron oxide	
<i>pure, magnetic precipitated</i>	
packaging	Mfr. No
500 g Glass jar	I/1100/53
Fe ₃ O ₄ CAS: 1317-61-9 MW: 231.50 Density at 20°C: 5.18	Melting point: 1588.9°C EINECS: 215-277-5 P262
Guaranteed analysis	
Identity test	

Iron oxide	
<i>pure, precipitated</i>	
packaging	Mfr. No
500 g Plastic jar	I/1150/53
Fe ₂ O ₃ CAS: 1309-37-1 MW: 159.69	Density at 20°C: 5.24 Melting point: 1538°C EINECS: 215-168-2
Guaranteed analysis	
Identity test	

Iron(III) oxide calcined	
<i>pure</i>	
packaging	Mfr. No
1 kg Glass jar	I/1013/60
25 kg Fibre Keg	I/1013/70
Fe ₂ O ₃ CAS: 1309-37-1 MW: 159.69	Density at 20°C: 5.24 Melting point: 1538°C EINECS: 215-168-2
Assay: >95%	

Iron(II) sulfate heptahydrate		Certified AR
<i>for analysis</i>		
packaging	Mfr. No	
250 g Plastic jar	I/1157/50	
500 g Plastic jar	I/1157/53	
1 kg Plastic jar	I/1157/60	
FeSO ₄ ·7H ₂ O CAS: 7782-63-0 MW: 278.01 Melting point: 64°C EINECS: 231-753-5	GHS: Acute Tox. 4, Skin Irrit. 2, Eye Irrit. 2 H315, H319, H302 P301+P312, P302+P352, P280, P305+P351+P338	
Assay: 98 to 102%		
Guaranteed analysis		
pH (aqueous solution 5%)	3 to 4	Potassium (K)
Identity test	: Pass	Sodium (Na)
Calcium (Ca)	<50 ppm	Total chloride (Cl)
Copper (Cu)	<10 ppm	Total nitrogen (N)
Iron (III)	<0.05%	Total phosphorus (P)
Lead (Pb)	<10 ppm	Total silicon (Si)
Magnesium (Mg)	<50 ppm	Zinc (Zn)

Typical analysis		
Silver (Ag)	<5 ppm	Lithium (Li)
Arsenic (As)	<5 ppm	Manganese (Mn)
Barium (Ba)	<5 ppm	Molybdenum (Mo)
Cadmium (Cd)	<0.01 %	Strontium (Sr)
Chromium (Cr)	<0.01 %	Titanium (Ti)
Cobalt (Co)	<0.001 %	Vanadium (V)

Lot analysis on the label.

Iron(II) sulfate heptahydrate		SLR
<i>specified reagents for general laboratory work- extra pure</i>		
packaging	Mfr. No	
1 kg Plastic jar	I/1155/60	
2.5 kg Plastic jar	I/1155/62	
5 kg Plastic jar	I/1155/65	
25 kg Plastic keg	I/1155/70	
FeSO ₄ ·7H ₂ O CAS: 7782-63-0 MW: 278.01 Melting point: 64°C EINECS: 231-753-5	GHS: Acute Tox. 4, Skin Irrit. 2, Eye Irrit. 2 H315, H319, H302 P301+P312, P302+P352, P280, P305+P351+P338	
Assay: >98%		
Guaranteed analysis		
Identity test	: Pass	Sodium (Na)
Calcium (Ca)	<200 ppm	n-Butyl formate
Copper (Cu)	<50 ppm	Total phosphorus (P)
Lead (Pb)	<50 ppm	Total silicon (Si)
Magnesium (Mg)	<200 ppm	Zinc (Zn)
Potassium (K)	<500 ppm	<500 ppm

Iron(II) sulfate heptahydrate	
<i>Technical</i>	
packaging	Mfr. No
2 kg Plastic jar	I/1154/61
FeSO ₄ ·7H ₂ O CAS: 7782-63-0 MW: 278.01 Melting point: 64°C EINECS: 231-753-5	GHS: Acute Tox. 4, Skin Irrit. 2, Eye Irrit. 2 H315, H319, H302 P301+P312, P302+P352, P280, P305+P351+P338

Iron(II) sulfate dried		SLR
<i>specified reagents for general laboratory work- extra pure</i>		
packaging	Mfr. No	
1 kg Plastic jar	I/1159/60	
FeSO ₄ CAS: 7782-63-0 MW: 278.01 Melting point: 64°C EINECS: 231-753-5	GHS: Acute Tox. 4, Skin Irrit. 2, Eye Irrit. 2 H315, H319, H302 P301+P312, P302+P352, P280, P305+P351+P338	
Assay: 77 to 90%		
Guaranteed analysis		
Copper (Cu)	<500 ppm	Lead (Pb)

Content 77 to 90 % in FeSO 4		
------------------------------	--	--

Iron(III) sulfate monohydrate	
<i>Technical</i>	
packaging	Mfr. No
500 g Plastic jar	I/1160/53
1 kg Plastic jar	I/1160/60
Fe ₂ (SO ₄) ₃ ·H ₂ O CAS: 15244-10-7 MW: 399.86 EINECS: 233-072-9	GHS: Acute Tox. 4 H302 P301+P312
Guaranteed analysis	
Identity test	

Iron(II) sulfide	
<i>Technical</i>	
packaging	Mfr. No
1 kg Plastic jar	I/1170/60
FeS CAS: 1317-37-9 MW: 87.92 Density at 20°C: 4.84	Melting point: 1194°C EINECS: 215-268-6 EUH032
Guaranteed analysis	
Identity test	

Isobutanol *See iso-Butyl alcohol, page 68*

Isohexane (contains <5% n-Hexane) **NEW PACK.**

for HPLC

Certified

packaging		Mfr. No
1 ℓ	Coated Glass bottle	H/0405/PB15
2.5 ℓ	Coated Glass bottle	H/0405/PB17
5 ℓ	Aluminium can	H/0405/21
10 ℓ	Returnable Drum	H/0405/10RSS
25 ℓ	Aluminium can	H/0405/25
25 ℓ	Returnable Drum	H/0405/25RSS
30 ℓ	Returnable Drum	H/0405/30RSS
C ₆ H ₁₄ CAS: 73513-42-5 MW: 86.18 Density at 20°C: 0.660 Boiling point: 60°C Melting point: -154°C Flash point: -23°C ADR/RID: 3, II EINECS: 203-523-4		UN: 1208 GHS: Flam. Liq. 2, Skin Irrit. 2, STOT SE 3, Asp. Tox. 1, Aquatic Chronic 2 H225, H304, H315, H336, H411 P261, P273, P301+P330+P331, P302+P352, P210
		
		
		
		
Guaranteed analysis		
Acidity / Alkalinity	<0.00015 meq/g	Residue after evaporation
Water	<0.01%	<2 ppm
Absorbance specification		
Max. Absorbance (1cm cell v HPLC water)		
200 nm	<0.3 A.U.	240 nm
210 nm	<0.15 A.U.	250 nm
220 nm	<0.07 A.U.	260 nm
230 nm	<0.03 A.U.	<0.005 A.U.

Lot analysis on the label.
Filtered to 0.2 µm.

Isohexane (contains <5% n-Hexane)

Certified

HPLC for fluorescence detection

packaging		Mfr. No
1 ℓ	Coated Glass bottle	H/0410/PB15
2.5 ℓ	Coated Glass bottle	H/0410/PB17
C ₆ H ₁₄ CAS: 73513-42-5 MW: 86.18 Density at 20°C: 0.660 Boiling point: 60°C Melting point: -154°C Flash point: -23°C ADR/RID: 3, II EINECS: 203-523-4		UN: 1208 GHS: Flam. Liq. 2, Skin Irrit. 2, STOT SE 3, Asp. Tox. 1, Aquatic Chronic 2 H225, H304, H315, H336, H411 P261, P273, P301+P330+P331, P302+P352, P210
		
		
		
		
Guaranteed analysis		
Acidity / Alkalinity	<0.00015 meq/g	Residue after evaporation
Water	<0.01%	<2 ppm
Absorbance specification		
Max. Absorbance (1cm cell v HPLC water)		
200 nm	<0.3 A.U.	240 nm
210 nm	<0.15 A.U.	250 nm
220 nm	<0.07 A.U.	260 nm
230 nm	<0.03 A.U.	Fluorescence
Excitation in the range 250-750nm produces no fluorescence emission greater than the solvents natural Raman signal at 350nm excittation.This is equivalent to the emission maximum of approximately 1ppb quinine sulfate.		
Lot analysis on the label (excluding fluorescence scan).		
Filtered to 0.2 µm.		

Isohexane (contains <5% n-Hexane)

for residue and pesticide analysis

Distol

packaging		Mfr. No
1 ℓ	Glass bottle	H/0408/15
2.5 ℓ	Glass bottle	H/0408/17
C ₆ H ₁₄ CAS: 73513-42-5 MW: 86.18 Density at 20°C: 0.660 Boiling point: 60°C Melting point: -154°C Flash point: -23°C ADR/RID: 3, II EINECS: 203-523-4		UN: 1208 GHS: Flam. Liq. 2, Skin Irrit. 2, STOT SE 3, Asp. Tox. 1, Aquatic Chronic 2 H225, H304, H315, H336, H411 P261, P273, P301+P330+P331, P302+P352, P210
		
		
		
		
Guaranteed analysis		
Water	<0.01%	GC (FID) Analysis
Residue after evaporation	<2 ppm	GC (NPD) analysis
GC (ECD) Analysis	Pass	Pass
GC Analysis		
Analysis of pesticides: from the retention time of Lindane, after concentration of 2500 times, no impurity peak higher than that obtained with 10 pg / ml héptachlo-réoxide (ECD) and 5 pg / ml Parathion (NPD).		
Analysis of hydrocarbons: no impurity chain longer than the dodecane is present at a concentration greater than 0.01 ppm (FID).		
Conditions : Column: J&W DB 5: 30m x 0.32mm x 0.1µm df. Temperature profile: (1) 50°C, 2 min (2) 40°C / min until 320°C (3) 320°C for 30 min Injection: 5µl cold on the column.		

Isohexane (contains <5% n-Hexane) **NEW PACK.**

Certified AR

for analysis

packaging		Mfr. No
1 ℓ	Coated Glass bottle	H/0411/PB15
2.5 ℓ	Coated Glass bottle	H/0411/PB17
5 ℓ	Fluorinated HDPE Plastic bottle	H/0411/21
200 ℓ	Metal drum	H/0411/27
C ₆ H ₁₄ CAS: 73513-42-5 MW: 86.18 Density at 20°C: 0.660 Boiling point: 60°C Melting point: -154°C Flash point: -23°C ADR/RID: 3, II EINECS: 203-523-4		UN: 1208 GHS: Flam. Liq. 2, Skin Irrit. 2, STOT SE 3, Asp. Tox. 1, Aquatic Chronic 2 H225, H304, H315, H336, H411 P261, P273, P301+P330+P331, P302+P352, P210
		
		
		
		
Guaranteed analysis		
Acidity / Alkalinity	<0.00008 meq/g	Lead (Pb)
Colour	<10 APHA	Magnesium (Mg)
Water	<0.02%	Potassium (K)
Residue after evaporation	<20 ppm	Sodium (Na)
Substances darkened by (H ₂ SO ₄)	Pass	Total phosphorus (P)
Bromine number	<0.5	Total silicon (Si)
Calcium (Ca)	<0.2 ppm	Total sulphur (S)
Copper (Cu)	<0.02 ppm	Zinc (Zn)
Iron (Fe)	<0.1 ppm	Aromatic Hydrocarbons
		N-Hexane
Typical analysis		
Aluminium (Al)	<0.1 ppm	Manganese (Mn)
Barium (Ba)	<0.02 ppm	Molybdenum (Mo)
Cadmium (Cd)	<0.02 ppm	Nickel (Ni)
Chromium (Cr)	<0.02 ppm	Strontium (Sr)
Cobalt (Co)	<0.02 ppm	Titanium (Ti)
Lithium (Li)	<0.02 ppm	Vanadium (V)

Lot analysis on the label.

Isohexane (contains <5% n-Hexane) **NEW PACK.**

SLR

specified reagents for general laboratory work- extra pure

packaging		Mfr. No
1 ℓ	Glass bottle	H/0401/15
2.5 ℓ	Coated Glass bottle	H/0401/PB17
5 ℓ	Fluorinated HDPE Plastic bottle	H/0401/21
10 ℓ	Metal can	H/0401/24
25 ℓ	Metal drum	H/0401/25
200 ℓ	Metal drum	H/0401/27
C ₆ H ₁₄ CAS: 73513-42-5 MW: 86.18 Density at 20°C: 0.660 Boiling point: 60°C Melting point: -154°C Flash point: -23°C ADR/RID: 3, II EINECS: 203-523-4		UN: 1208 GHS: Flam. Liq. 2, Skin Irrit. 2, STOT SE 3, Asp. Tox. 1, Aquatic Chronic 2 H225, H304, H315, H336, H411 P261, P273, P301+P330+P331, P302+P352, P210
		
		
		
		
Guaranteed analysis		
Acidity / Alkalinity	<0.0002 meq/g	Magnesium (Mg)
Colour	<20 APHA	Potassium (K)
Water	<0.05%	Sodium (Na)
Residue after evaporation	<50 ppm	Total phosphorus (P)
Bromine number	<0.5	Total silicon (Si)
Calcium (Ca)	<2 ppm	Total sulphur (S)
Copper (Cu)	<1 ppm	Zinc (Zn)
Iron (Fe)	<5 ppm	N-Hexane
Lead (Pb)	<1 ppm	<5%

Isooctane *See 2, page 360*

Isopentyl alcohol *See iso-Amyl alcohol, page 47*

Isopropanol **NEW**

Optima

for LC/MS analysis

packaging		Mfr. No
500 mℓ	Glass bottle	A461-500
1 ℓ	Glass bottle	A461-1
2.5 ℓ	Glass bottle	A461-212
4 ℓ	Glass bottle	A461-4
(CH ₃) ₂ CHOH CAS: 67-63-0 MW: 60.10 Density at 20°C: 0.785 Boiling point: 81 to 83°C Melting point: -89.5°C Flash point: 12°C ADR/RID: 3, II		EINECS: 200-661-7 UN: 1219 GHS: Flam. Liq. 2, Eye Irrit. 2, STOT SE 3 H225, H336, H319 P261, P280, P305+P351+P338, P210, P240
		
		
Assay: >99.9%		
Guaranteed analysis		
Acidity / Alkalinity	<0.0001 meq/g	Lead (Pb)
Colour	<5 APHA	Magnesium (Mg)
Water	<0.05%	Manganese (Mn)
Residue after evaporation	<1 ppm	Nickel (Ni)
Aluminium (Al)	<10 ppb	Potassium (K)
Calcium (Ca)	<10 ppb	Propan--1-ol
Copper (Cu)	<5 ppb	Sodium (Na)
Iron (Fe)	<5 ppb	Zinc (Zn)
Absorbance specification		
Max. Absorbance (1cm cell v HPLC water)		
210 nm	<0.4 A.U.	254 nm
220 nm	<0.1 A.U.	LC/MS Gradient Suitability
230 nm	<0.05 A.U.	: Pass Test
This product has been tested on an LC-MS Thermo Scientific instrument to ensure suitability. Filtered to 0.2µm.		

Isopropanol

HPLC for gradient analysis

Certified

packaging		Mfr. No
2.5 ℓ	Glass bottle	P/7508/17
(CH ₃) ₂ CHOH CAS: 67-63-0 MW: 60.10 Density at 20°C: 0.785 Boiling point: 81 to 83°C Melting point: -89.5°C Flash point: 12°C ADR/RID: 3, II		EINECS: 200-661-7 UN: 1219 GHS: Flam. Liq. 2, Eye Irrit. 2, STOT SE 3 H225, H336, H319 P261, P280, P305+P351+P338, P210, P240
		
		
Assay (CG): >99.5%		
Guaranteed analysis		
Acidity / Alkalinity	<0.0001 meq/g	Residue after evaporation
Water	<0.1%	<2 ppm
Absorbance specification		
Max. Absorbance (1cm cell v HPLC water)		
205 nm	<1 A.U.	260 nm
220 nm	<0.2 A.U.	270 nm
230 nm	<0.1 A.U.	Gradient (max eluted peak) ..
240 nm	<0.05 A.U.	Wavelength at 1 A.U.
250 nm	<0.01 A.U.	: <205 nM
HPLC analysis:		
Gradient analysis Absorbance of the largest peak <0.002A.U.		
Conditions: Column: C18 Reverse Phase Flow rate: 0.5 ml / min Wavelength 254 nm Solvent A: Water Solvent B: Propan -2-ol Gradient: Linear from 95% A/ 5% B to 100% B over 20mins. Lot analysis on the label. Filtered to 0.2 µm. Gamme d'absorbance:		

Isopropanol **NEW PACK.**

for HPLC

Certified

packaging		Mfr. No
1 ℓ	Glass bottle	P/7507/15
2.5 ℓ	Glass bottle	P/7507/17
2.5 ℓ	Coated Glass bottle	P/7507/PB17
10 ℓ	Returnable Drum	P/7507/10RSS
25 ℓ	Returnable Drum	P/7507/25RSS
30 ℓ	Returnable Drum	P/7507/30RSS
(CH ₃) ₂ CHOH CAS: 67-63-0 MW: 60.10 Density at 20°C: 0.785 Boiling point: 81 to 83°C Melting point: -89.5°C Flash point: 12°C ADR/RID: 3, II		EINECS: 200-661-7 UN: 1219 GHS: Flam. Liq. 2, Eye Irrit. 2, STOT SE 3 H225, H336, H319 P261, P280, P305+P351+P338, P210, P240
		
		
Assay (CG): >99.5%		
Guaranteed analysis		
Acidity / Alkalinity	<0.0001 meq/g	Residue after evaporation
Water	<0.1%	<2 ppm
Absorbance specification		
Max. Absorbance (1cm cell v HPLC water)		
205 nm	<1 A.U.	250 nm
220 nm	<0.2 A.U.	260 nm
230 nm	<0.1 A.U.	270 nm
240 nm	<0.05 A.U.	Wavelength at 1 A.U.
Lot analysis on the label. Filtered to 0.2 µm.		

Isopropanol
HPLC for electrochemical detection

Certified

packaging		Mfr. No
2.5 ℓ	Glass bottle	P/7509/17
(CH ₃) ₂ CHOH	EINECS: 200-661-7	
CAS: 67-63-0	UN: 1219	
MW: 60.10	GHS: Flam. Liq. 2, Eye Irrit.	
Density at 20°C: 0.785	2, STOT SE 3	
Boiling point: 81 to 83°C	H225, H336, H319	
Melting point: -89.5°C	P261, P280,	
Flash point: 12°C	P305+P351+P338, P210, P240	
ADR/RID: 3, II		
Assay (CG): >99.5%		
Guaranteed analysis		
Acidity / Alkalinity	<0.0001 meq/g	Residue after evaporation
Water	<0.05%	Electrochemical test

Absorbance specification	
Max. Absorbance (1cm cell v HPLC water)	
205 nm	<1 A.U.
220 nm	<0.2 A.U.
230 nm	<0.1 A.U.
240 nm	<0.05 A.U.
Wavelength at 1 A.U.	
: <205 nm	
Tested by HPLC-EC to ensure low electrochemical response.	
Detector - ESA Coulochem II	
Mode-oxidative	
Graphite Electrode.	
Lot analysis on the label.	
Filtered to 0.2 µm.	

Isopropanol
for residue and pesticide analysis

Distol

packaging		Mfr. No	
1 ℓ	Glass bottle	P/7501/15	
2.5 ℓ	Glass bottle	P/7501/17	
(CH ₃) ₂ CHOH	EINECS: 200-661-7		
CAS: 67-63-0	UN: 1219		
MW: 60.10	GHS: Flam. Liq. 2, Eye Irrit.		
Density at 20°C: 0.785	2, STOT SE 3		
Boiling point: 81 to 83°C	H225, H336, H319		
Melting point: -89.5°C	P261, P280,		
Flash point: 12°C	P305+P351+P338, P210, P240		
ADR/RID: 3, II			
Assay (CG): >99.8%			
Guaranteed analysis			
Water	<0.1%	GC (FID) Analysis	Pass
Residue after evaporation	<2 ppm	GC (NPD) analysis	Pass
GC (ECD) Analysis	Pass		

GC Analysis
Analysis of pesticides: from the retention time of Lindane, after concentration of 2500 times, no impurity peak higher than that obtained with 10 pg / ml heptachlor-éoxide (ECD) and 5 pg / ml Parathion (NPD).
Analysis of hydrocarbons: no impurity chain longer than the dodecane is present at a concentration greater than 0.01 ppm (FID).
Conditions :
Column: J&W DB 5; 30m x 0.32mm x 0.1µm df.
Temperature profile:
(1) 50°C for 2 min
(2) 40°C / min until 320°C
(3) 320°C for 30 min
Injection: 5µl cold on the column.

Isopropanol
extra dry, for synthesis

packaging		Mfr. No
1 ℓ	Glass bottle	P/7503/15
2.5 ℓ	Glass bottle	P/7503/17
(CH ₃) ₂ CHOH	EINECS: 200-661-7	
CAS: 67-63-0	UN: 1219	
MW: 60.10	GHS: Flam. Liq. 2, Eye Irrit.	
Density at 20°C: 0.785	2, STOT SE 3	
Boiling point: 81 to 83°C	H225, H336, H319	
Melting point: -89.5°C	P261, P280,	
Flash point: 12°C	P305+P351+P338, P210, P240	
ADR/RID: 3, II		

Assay (CG): >99.8%		
Guaranteed analysis		
Water	<0.01%	Residue after evaporation
		<2 ppm

Isopropanol
for electronic use (MOS)

packaging		Mfr. No
2.5 ℓ	Glass bottle	P/7555/17
5 ℓ	HDPE Plastic bottle	P/7555/21
25 ℓ	Metal drum	P/7555/25
(CH ₃) ₂ CHOH	EINECS: 200-661-7	
CAS: 67-63-0	UN: 1219	
MW: 60.10	GHS: Flam. Liq. 2, Eye Irrit.	
Density at 20°C: 0.785	2, STOT SE 3	
Boiling point: 81 to 83°C	H225, H336, H319	
Melting point: -89.5°C	P261, P280,	
Flash point: 12°C	P305+P351+P338, P210, P240	
ADR/RID: 3, II		

Assay (CG): >99.8%		
Guaranteed analysis		
Acidity / Alkalinity	<0.0001 meq/g	Nickel (Ni)
Colour	<5 APHA	Potassium (K)
Water	<0.1%	Sodium (Na)
Residue after evaporation	<10 ppm	Strontium (Sr)
Aluminium (Al)	<0.02 ppm	Titanium (Ti)
Barium (Ba)	<0.02 ppm	Total phosphorus (P)
Cadmium (Cd)	<0.02 ppm	Total silicon (Si)
Calcium (Ca)	<0.2 ppm	Total sulphur (S)
Chromium (Cr)	<0.02 ppm	Vanadium (V)
Cobalt (Co)	<0.02 ppm	Zinc (Zn)
Copper (Cu)	<0.02 ppm	Substances darkened
Iron (Fe)	<0.1 ppm	by H ₂ SO ₄
Lead (Pb)	<0.02 ppm	Alcohols C 4
Lithium (Li)	<0.02 ppm	Alcohols C 5
Magnesium (Mg)	<0.1 ppm	Ethanol 90 % v/v
Manganese (Mn)	<0.02 ppm	Methanol
Molybdenum (Mo)	<0.02 ppm	Propan--1-ol

Isopropanol (Propan-2-ol) NEW PACK.

Certified AR

for analysis

packaging		Mfr. No
500 mℓ	Coated Glass bottle	P/7500/PB08
1 ℓ	Glass bottle	P/7500/15
2.5 ℓ	Glass bottle	P/7500/17
2.5 ℓ	Coated Glass bottle	P/7500/PB17
2.5 ℓ	HDPE Plastic bottle	P/7500/PC17
5 ℓ	Fluorinated HDPE Plastic bottle	P/7500/21
10 ℓ	HDPE Plastic container	P/7500/24
25 ℓ	Metal drum	P/7500/25
25 ℓ	HDPE Plastic container	P/7500/PC25
200 ℓ	Metal drum	P/7500/27

(CH ₃) ₂ CHOH	EINECS: 200-661-7	
CAS: 67-63-0	UN: 1219	
MW: 60.10	GHS: Flam. Liq. 2, Eye Irrit.	
Density at 20°C: 0.785	2, STOT SE 3	
Boiling point: 81 to 83°C	H225, H336, H319	
Melting point: -89.5°C	P261, P280,	
Flash point: 12°C	P305+P351+P338, P210, P240	
ADR/RID: 3, II		

Assay (CG): >99.8%		
Guaranteed analysis		
Acidity / Alkalinity	<0.0001 meq/g	Sodium (Na)
Colour	<5 APHA	Total phosphorus (P)
Water	<0.1%	Total silicon (Si)
Residue after evaporation	<10 ppm	Total sulphur (S)
Substances		Zinc (Zn)
reducing KMnO ₄	<0.0005%	Alcohols C 4
Calcium (Ca)	<0.2 ppm	Alcohols C 5
Copper (Cu)	<0.02 ppm	Ethanol 90 % v/v
Iron (Fe)	<0.1 ppm	Methanol
Lead (Pb)	<0.02 ppm	Propan--1-ol
Magnesium (Mg)	<0.1 ppm	Substances darkened
Potassium (K)	<0.2 ppm	by H ₂ SO ₄

Typical analysis		
Aluminium (Al)	<0.02 ppm	Manganese (Mn)
Barium (Ba)	<0.02 ppm	Molybdenum (Mo)
Cadmium (Cd)	<0.02 ppm	Nickel (Ni)
Chromium (Cr)	<0.02 ppm	Strontium (Sr)
Cobalt (Co)	<0.02 ppm	Titanium (Ti)
Lithium (Li)	<0.02 ppm	Vanadium (V)

Lot analysis on the label.

Isopropanol (Propan-2-ol) NEW PACK.

specified reagents for general laboratory work- extra pure

SLR

packaging		Mfr. No
500 mℓ	Glass bottle	P/7490/08
500 mℓ	Coated Glass bottle	P/7490/PB08
1 ℓ	Glass bottle	P/7490/15
2.5 ℓ	Glass bottle	P/7490/17
2.5 ℓ	Coated Glass bottle	P/7490/PB17
5 ℓ	Fluorinated HDPE Plastic bottle	P/7490/21
5 ℓ	Fluorinated HDPE square Drum	P/7490/FP21
25 ℓ	Metal drum	P/7490/25
25 ℓ	HDPE Plastic container	P/7490/PC25
200 ℓ	Metal drum	P/7490/27

(CH ₃) ₂ CHOH	EINECS: 200-661-7	
CAS: 67-63-0	UN: 1219	
MW: 60.10	GHS: Flam. Liq. 2, Eye Irrit.	
Density at 20°C: 0.785	2, STOT SE 3	
Boiling point: 81 to 83°C	H225, H336, H319	
Melting point: -89.5°C	P261, P280,	
Flash point: 12°C	P305+P351+P338, P210, P240	
ADR/RID: 3, II		

Assay (CG): >99.5%		
Guaranteed analysis		
Acidity / Alkalinity	<0.0005 meq/g	Iron (Fe)
Colour	<10 APHA	Lead (Pb)
Water	<0.2%	Magnesium (Mg)
Residue after evaporation	<50 ppm	Potassium (K)
Substances darkened		Sodium (Na)
by (H ₂ SO ₄)	<10 APHA	Total phosphorus (P)
Substances		Total silicon (Si)
reducing KMnO ₄	<0.0005%	Total sulphur (S)
Calcium (Ca)	<5 ppm	Zinc (Zn)
Copper (Cu)	<2 ppm	

Isopropylacetone See 4-Methylpentan-2-one, page 206

Isopropyl alcohol See Isopropanol, page 168

Isopropyl iodide See 2-Iodopropane, page 162

Kaolin
technical

packaging		Mfr. No
3 kg	Plastic keg	K/0045/63
H ₂ Al ₂ Si ₂ O ₈ .H ₂ O	Density at 20°C: 1.8 to 2.	
CAS: 1332-58-7	Melting point: 1760°C	
MW: 258.20	EINECS: 310-194-1	

UHPLC
Filtered to 0.1 µm, Fisher Chemical UHPLC gradient grade solvents (Acetonitrile, Methanol and Water) are perfect for those working at high pressure.
See page 309 for more information

Kaolin light
pure, washed with acid

packaging		Mfr. No
1 kg	Plastic jar	K/0050/60
3 kg	Plastic keg	K/0050/63
25 kg	Fibre Keg	K/0050/70
H ₂ Al ₂ Si ₂ O ₈ .H ₂ O	Density at 20°C: 1.8 to 2.	
CAS: 1332-58-7	Melting point: 1750°C	
MW: 258.20	EINECS: 310-194-1	

Guaranteed analysis	
Identity test	: Pass



Karl Fischer reagents for Volumetric Analysis

For high water content.
Available as 1 or 2 component system, pyridine free.
See alphabetical entry in the main body of this catalogue for full details.

designation	packsize	package type	Mfr. No	page
Karl Fischer Aqualine® Complete 1	1 L	Glass bottle	K/1900/15	171
Karl Fischer Aqualine® Complete 2	1 L	Glass bottle	K/1950/15	171
Karl Fischer Aqualine® Complete 5	1 L	Glass bottle	K/2000/15	171
Karl Fischer Aqualine® Complete 5K	500 mℓ	Glass bottle	K/2250R/08	171
	1 L	Glass bottle	K/2250R/15	171
Karl Fischer Aqualine® Matrix K	500 mℓ	Glass bottle	K/2300R/08	171
	1 L	Glass bottle	K/2300R/15	171
Karl Fischer Aqualine® Solvent	1 L	Glass bottle	K/2100/15	171
	2.5 L	Glass bottle	K/2100/17	171
Karl Fischer Aqualine® solvent CM	1 L	Glass bottle	K/2110/15	172
	2.5 L	Glass bottle	K/2110/17	172
Karl Fischer Aqualine® Titrant 2	1 L	Glass bottle	K/2150/15	172
Karl Fischer Aqualine® Titrant 5	1 L	Glass bottle	K/2200/15	172

Karl Fischer reagents for Coulometric Analysis

For low water content at ppm level.
Anolytes & Catholyte solutions.
See alphabetical entry in the main body of this catalogue for full details.

designation	packsize	package type	Mfr. No	page
Karl Fischer Aqualine® Electrolyte A	500 mℓ	Glass bottle	K/2500/08	172
Karl Fischer Aqualine® Electrolyte AD	500 mℓ	Glass bottle	K/2510/08	172
Karl Fischer Aqualine® Electrolyte AD-G	500 mℓ	Glass bottle	K/2515/08	172
Karl Fischer Aqualine® Electrolyte AG	500 mℓ	Glass bottle	K/2520/08	173
Karl Fischer Aqualine® Electrolyte CG	25 mℓ	Vial	K/2560/04	173
Karl Fischer Aqualine® Electrolyte AG-H	500 mℓ	Glass bottle	K/2530/08	174

Standards for Volumetric & Coulometric application

designation	packsize	package type	Mfr. No	page
Karl Fischer Aqualine® sodium tartrate standard	100 g	Plastic jar	K/2770/48	173
Karl Fischer Aqualine® tri-Potassium citrate standard	10 g	Glass bottle	K/2760/45	173
Karl Fischer Aqualine® Water solution	500 mℓ	Coated Glass bottle	J/9000/PB08	173
	1 L	Glass bottle	J/9000/15	173
Karl Fischer Aqualine® Water standard 0.2 mg/g, 10 x 4mℓ ampoules	40 mℓ	Ampoule	K/2740/99	173
Karl Fischer Aqualine® Water standard 1 mg/g, 10 x 4mℓ ampoules	40 mℓ	Ampoule	K/2710/99	173
Karl Fischer Aqualine® Water standard 5 mg/mℓ	500 mℓ	Glass bottle	K/2730/08	173
Karl Fischer Aqualine® Water standard 10 mg/g, 10 x 8mℓ ampoules	80 mℓ	Ampoule	K/2720/99	174
Karl Fischer Aqualine Buffer	500 mℓ	Glass bottle	K/2350/08	174
Karl Fischer reagent solution, modified	500 mℓ	Coated Glass bottle	J/4620/PB08	174
	1 L	Coated Glass bottle	J/4620/PB15	174
	2.5 L	Coated Glass bottle	J/4620/PB17	174

Karl Fischer Aqualine® Complete 1
For Karl Fischer titration by Volumetry

packaging	Mfr. No
1 ℓ Glass bottle	K/1900/15
ADR/RID: 8, III EINECS: 203-919-7 UN: 1760 GHS: Repr. 2, Acute Tox. 2, Skin Corr. 1B, Eye Dam. 1	H314, H302, H361d, H331 P280, P301+P330+P331, P305+P351+P338, P304+P340, P302+P350, P310, P273

Guaranteed analysis
Factor : 1.3 to 1.5 mg/mℓ
Single component (composite) reagent with a water equivalence of 1mg Water/mℓ (pyridine free)

Karl Fischer Aqualine® Complete 2
For Karl Fischer titration by Volumetry

packaging	Mfr. No
1 ℓ Glass bottle	K/1950/15
2.5 ℓ Glass bottle	K/1950/17
ADR/RID: 8, III EINECS: 203-919-7 UN: 1760 GHS: Repr. 2, Acute Tox. 2, Skin Corr. 1B, Eye Dam. 1	H314, H302, H361d, H331 P280, P301+P330+P331, P305+P351+P338, P304+P340, P302+P350, P310, P273

Guaranteed analysis
Factor : 2.4 to 2.6 mg/mℓ
Single component (composite) reagent with a water equivalence of 2mg Water/mℓ (pyridine free)

Karl Fischer Aqualine® Complete 5
For Karl Fischer titration by Volumetry

packaging	Mfr. No
1 ℓ Glass bottle	K/2000/15
2.5 ℓ Glass bottle	K/2000/17
ADR/RID: 8, III EINECS: 203-919-7 UN: 1760 GHS: Repr. 2, Acute Tox. 2, Skin Corr. 1B, Eye Dam. 1	H314, H302, H361d, H331 P280, P301+P330+P331, P305+P351+P338, P304+P340, P302+P350, P310, P273

Guaranteed analysis
Factor : 5.3 to 5.6 mg/mℓ
Single component (composite) reagent with a water equivalence of 5mg Water/mℓ (pyridine free)



Karl Fischer Aqualine® Complete 5K
For Karl Fischer titration by Volumetry

packaging	Mfr. No
500 mℓ Glass bottle	K/2250R/08
1 ℓ Glass bottle	K/2250R/15
ADR/RID: 8, III EINECS: 203-919-7 UN: 1760 GHS: Repr. 2, Acute Tox. 4, Skin Corr. 1B, Aquatic Acute 1	H302, H312, H332, H314, H361d, H400 P280, P301+P330+P331, P305+P351+P338, P301+P312, P304+P340, P273

Guaranteed analysis
Factor : 5.3 to 5.5 mg/mℓ
Single component (composite) reagent for water determination in aldehydes and ketones with a water equivalence of 5mg Water/mℓ (pyridine free). Must be used in conjunction with Aqualine® Matrix K.

Karl Fischer Aqualine® Matrix K
For Karl Fischer titration by Volumetry

packaging	Mfr. No
500 mℓ Glass bottle	K/2300R/08
1 ℓ Glass bottle	K/2300R/15
Flash point: >36.1°C ADR/RID: 3, 6.1, III EINECS: 200-663-8 UN: 1992 GHS: Flam. Liq. 3, Carc. 2, Acute Tox. 4, Skin Irrit. 2,	Eye Irrit. 2, STOT RE 2 H226, H302, H312, H332, H351, H373, H319, H315 P280, P304+P340, P302+P352, P260, P305+P351+P338, P210

A mixture of solvents used as working medium for the determination of water in aldehydes and ketones. Must be used in conjunction with Aqualine®Complete 5K.

Karl Fischer Aqualine® Solvent
For Karl Fischer titration by Volumetry

packaging	Mfr. No
1 ℓ Glass bottle	K/2100/15
2.5 ℓ Glass bottle	K/2100/17
Flash point: >11°C ADR/RID: 3, 6.1, II EINECS: 200-659-6 UN: 1992 GHS: Flam. Liq. 2, Repr. 2, Acute Tox. 3, Skin Corr. 1B, Eye Dam. 1, STOT SE 1 H225, H301, H331, H370,	H311, H314, H361d P301+P310, P280, P312, P302+P350, P301+P330+P331, P305+P351+P338, P310, P281, P307+P311, P304+P340, P210, P240

Aqualine® Solvent acts as the working medium with the two-component reagent (pyridine free). Water capacity of at least 8mg water/mℓ.



Karl Fischer Aqualine® Solvent CM

For Karl Fischer titration by Volumetry

packaging		Mfr. No
1 ℓ	Glass bottle	K/2110/15
2.5 ℓ	Glass bottle	K/2110/17
Flash point: 31°C	SE 1, STOT RE 2	   
ADR/RID: 3, 6.1, II	H225, H301, H311, H331, H314, H351, H370, H373, H360D	
EINECS: 200-663-8		
UN: 1992		
GHS: Flam. Liq. 2, Carc. 2, Repr. 1B, Acute Tox. 3, Skin Corr. 1B, Eye Dam. 1, STOT	P201, P280, P302+P350, P305+P351+P338, P308+P313, P210	

Aqualine® Solvent CM acts as the working medium with the two-component reagent (pyridine free) for samples with high hydrocarbon content.

Karl Fischer Aqualine® Titrant 2

For Karl Fischer titration by Volumetry

packaging		Mfr. No
1 ℓ	Glass bottle	K/2150/15
2.5 ℓ	Glass bottle	K/2150/17
Boiling point: 64°C	GHS: Flam. Liq. 2, Acute	  
Flash point: 10°C	Tox. 3, STOT SE 1	
ADR/RID: 3, 6.1, II	H225, H331, H311, H301, H370	
EINECS: 200-659-6	P261, P280, P307+P311, P302+P352, P210, P240	
UN: 1230		

Guaranteed analysis

Factor : 2.1 to 2.3 mg/mL

Aqualine® two component reagent with a water equivalence of 2mg Water/mL (pyridine free).

Aqualine® Titrant 2 is used as the titration reagent in conjunction with Aqualine® Solvent.

Karl Fischer Aqualine® Titrant 5

For Karl Fischer titration by Volumetry

packaging		Mfr. No
1 ℓ	Glass bottle	K/2200/15
2.5 ℓ	Glass bottle	K/2200/17
Boiling point: 64°C	Tox. 3, STOT SE 1	  
Flash point: >11°C	H225, H301, H311, H331, H370	
ADR/RID: 3, 6.1, II	P301+P310, P280, P302+P350, P304+P340, P210, P240	
EINECS: 200-659-6		
UN: 1230		
GHS: Flam. Liq. 2, Acute		

Guaranteed analysis

Factor : 5.1 to 5.3 mg/mL

Aqualine® two component reagent with a water equivalence of 5mg Water/mL (pyridine free).

Aqualine® Titrant 5 is used as the titration reagent in conjunction with Aqualine® Solvent.

Karl Fischer Aqualine® Electrolyte A

For Karl Fischer titration by Coulometry

packaging		Mfr. No
500 mL	Glass bottle	K/2500/08
Flash point: 10°C	H225, H314, H330, H311, H301, H370, H373, H351	   
ADR/RID: 3, 6.1, II	P307+P311, P280, P301+P330+P331, P305+P351+P338, P304+P340, P302+P352, P260, P210	
EINECS: 200-659-6		
UN: 1992		
GHS: Flam. Liq. 2, Carc. 2, Acute Tox. 2, Skin Corr. 1B, STOT SE 1, STOT RE 2		

Anolyte solution for general use in conventional cells with a diaphragm.

Karl Fischer Aqualine® Electrolyte AD

For Karl Fischer titration by Coulometry

packaging		Mfr. No
500 mL	Glass bottle	K/2510/08
Flash point: 10°C	H225, H331, H301, H373, H311, H314, H370, H351	   
ADR/RID: 3, 6.1, II	P301+P310, P280, P312, P302+P350, P301+P330+P331, P305+P351+P338, P310, P281, P307+P311, P304+P340, P260, P210, P240	
EINECS: 200-659-6		
UN: 1992		
GHS: Flam. Liq. 2, Carc. 2, Acute Tox. 3, Skin Corr. 1B, Eye Dam. 1, STOT SE 1, STOT RE 2		

Anolyte for cells without diaphragm (fritless).

Karl Fischer Aqualine® Electrolyte AD-G

For Karl Fischer titration by Coulometry

packaging		Mfr. No
500 mL	Glass bottle	K/2515/08
Flash point: 10°C	H225, H301, H311, H314, H331, H370, H411	     
ADR/RID: 3, 6.1, II	P210, P240, P280, P301+P310, P304+P340, P302+P350, P273	
UN: 1992		
GHS: Flam. Liq. 2, Acute		
Tox. 3, Skin Corr. 1B, STOT SE 1, Aquatic Chronic 2		

Halogen-free anolyte solution for use in cells without diaphragm (fritless).

High Volume Solvent Delivery Systems

High-volume solvent delivery systems are available in 10L to 1000L, in Europe only. Please email your contact details and a brief description of your requirements to bege1.sha.quotes@thermofisher.com



Karl Fischer Aqualine® Electrolyte AG

For Karl Fischer titration by Coulometry

packaging		Mfr. No
500 mL	Glass bottle	K/2520/08
Flash point: 10°C	H225, H301, H311, H330, H314, H370, H400	    
ADR/RID: 3, 6.1, II	P280, P301+P330+P331, P305+P351+P338, P304+P340, P210, P273, P240	
UN: 1992		
GHS: Flam. Liq. 2, Acute		
Tox. 2, Skin Corr. 1B, STOT SE 1, Aquatic Acute 1		

Halogen-free anolyte solution for general use in conventional cells with diaphragm.

Karl Fischer Aqualine® Electrolyte CG

For Karl Fischer titration by Coulometry

packaging		Mfr. No
25 mL	Vial	K/2560/04
Density at 20°C: >> 1.0	Tox. 3, Skin Irrit. 2, Eye Irrit. 2, STOT SE 1, Aquatic Acute 1	   
Boiling point: 64°C	H225, H311, H331, H301, H315, H319, H370, H400	
Flash point: >11°C	P280, P304+P340, P302+P350, P302+P352, P305+P351+P338, P210, P273	
ADR/RID: 3, 6.1, II		
EINECS: 200-659-6		
UN: 1992		
GHS: Flam. Liq. 2, Acute		

Chlorine-free catholyte solution for general use in conventional cells with a diaphragm.

Karl Fischer Aqualine® Sodium tartrate-2-H₂O

For Karl Fischer titration

packaging		Mfr. No
100 g	Plastic jar	K/2770/48
(CH(OH).COONa) ₂	Density at 20°C: 1.8	
CAS: 6106-24-7	Melting point: 150°C	
MW: 230.08	EINECS: 212-773-2	
Guaranteed analysis		
Loss on drying : 15.61 to 15.71%		
Primary standard for volumetric analysis. water content 15.66±0.05%.		

Karl Fischer Aqualine® tri-Potassium citrate standard

For Karl Fischer titration

packaging		Mfr. No
10 g	Glass bottle	K/2760/45
K ₃ C ₆ H ₅ O ₇		
CAS: 6100-05-6		
EINECS: 212-755-5		
Guaranteed analysis		
Loss on drying 200°C : 5.5 to 5.61%		
Potassium citrate monohydrate for use in KF oven techniques. Actual water content by weight (as loss on drying) is specified on the label.		

Karl Fischer Aqualine® Water solution

For Karl Fischer titration

packaging		Mfr. No
500 mL	Coated Glass bottle	J/9000/PB08
1 ℓ	Glass bottle	J/9000/15
H ₂ O	UN: 1230	   
CAS: 7732-18-5	GHS: Flam. Liq. 2, Acute	
MW: 18.02	Tox. 3, STOT SE 1	
Density at 20°C: 0.79	H225, H301, H311, H331, H370	
Boiling point: 64.7°C	P301+P310, P280, P302+P350, P304+P340, P210, P240	
Flash point: 12°C		
ADR/RID: 3, 6.1, II		
EINECS: 200-659-6		

1 mL is equivalent to 5.00 + /-0.05 mg water.

Karl Fischer Aqualine® Water standard 0.2 mg

For Karl Fischer titration

packaging		Mfr. No
40 mL	Ampoule	K/2740/99
Density at 20°C: 0.865	UN: 1307	 
Boiling point: 138.5 to 142	GHS: Flam. Liq. 3, Acute	
Melting point: -47.4°C	Tox. 4, Skin Irrit. 2	
Flash point: 27°C	H226, H312, H332, H315	
ADR/RID: 3, III	P210, P261, P280, P302+P352, P304+P340	
EINECS: 215-535-7		

Guaranteed analysis

Water equivalent : 0.195 to 0.205 mg

Standard for Karl Fischer titration.

1 mL is equivalent to 0.200 + /-0.005 mg water.

Each pack contains 10 x 4mL ampoules and is supplied with a certificate of analysis.

Each ampoule has a shelf life of at least 5 years unopened.

Karl Fischer Aqualine® Water standard 1.0 mg

For Karl Fischer titration

packaging		Mfr. No
40 mL	Ampoule	K/2710/99
ADR/RID: 3, III	Dam. 1, STOT SE 3	  
EINECS: 203-273-6	H226, H332, H315, H318, H335	
UN: 2222	P304+P340, P302+P352, P210, P305+P351+P338, P261, P280	
GHS: Flam. Liq. 3, Acute		
Tox. 4, Skin Irrit. 2, Eye		

Guaranteed analysis

Water equivalent ... : 0.97 to 1.03 mg/mL

Standard for Karl Fischer titration.

1 mL is equivalent to 1.00 + /-0.03 mg water.

Each pack contains 10 x 4mL ampoules and is supplied with a certificate of analysis.

Each ampoule has a shelf life of at least 5 years unopened.

Karl Fischer Aqualine® Water standard 5 mg

For Karl Fischer titration

packaging		Mfr. No
500 mL	Glass bottle	K/2730/08
Flash point: 29°C	Tox. 4, Skin Irrit. 2, Eye	  
ADR/RID: 3, III	Dam. 1, STOT SE 3	
EINECS: 215-535-7	H226, H336, H318, H335, H315, H312, H302, H332	
UN: 1993	P210, P280, P305+P351+P338	
GHS: Flam. Liq. 3, Acute		

Guaranteed analysis

Water equivalent ... : 4.98 to 5.02 mg/mL

Standard for Karl Fischer titration.

1 mL is equivalent to 5.00 + /-0.02 mg water.

Karl Fischer Aqualine® Water standard 10 mg
For Karl Fischer titration

packaging		Mfr. No
80 mℓ	Ampoule	K/2720/99
Flash point: 32.2°C	Dam. 1, STOT SE 3	  
ADR/RID: 3, III	H226, H302, H312, H332,	
EINECS: 203-572-1	H318, H315, H335, H336	
UN: 1993	P280, P305+P351+P338,	
GHS: Flam. Liq. 3, Acute	P304+P340, P302+P352,	
Tox. 4, Skin Irrit. 2, Eye	P210, P261	

Guaranteed analysis
Water equivalent : 9.9 to 10.1 mg/mℓ

*Standard for Karl Fischer titration.
1 mℓ is equivalent to 10.00 + /-0.1 mg water.
Each pack contains 10 x 8mℓ ampoules and is supplied with a certificate of analysis.
Each ampoule has a shelf life of at least 5 years unopened.*

Karl Fischer, Aqualine buffer
For Karl Fischer titration

packaging		Mfr. No
500 mℓ	Glass bottle	K/2350/08
Flash point: 10°C	Eye Dam. 1, STOT SE 1	   
ADR/RID: 3, 8, II	H225, H370, H301, H314,	
UN: 2924	H311, H331, H361d	
GHS: Flam. Liq. 2, Repr. 2,	P210, P240,	
Acute Tox. 3, Skin Corr. 1B,	P305+P351+P338, P280	

Guaranteed analysis
Use test : Pass

Buffer for use with acidic samples

Karl Fischer reagent solution, modified
For Karl Fischer titration

packaging		Mfr. No
500 mℓ	Coated Glass bottle	J/4620/PB08
1 ℓ	Coated Glass bottle	J/4620/PB15
2.5 ℓ	Coated Glass bottle	J/4620/PB17
Density at 20°C: 0.93	GHS: Flam. Liq. 3, Repr. 1B,	  
Flash point: 38°C	Acute Tox. 3	
ADR/RID: 3, III	H226, H360, H331, H311, H301	
EINECS: 231-195-2	P301+P310, P302+P352,	
UN: 1993	P304+P312, P210, P201	

Water equivalent approximately 5 mg / mℓ when prepared.

Standard Volumetric solution

For your convenience we offer a range of NIST traceable solutions. A full list of our products can be found on page 323. Can't see what you need? Please email your contact details and a brief description of your requirements to bege.shi.quotes@thermofisher.com



Karl Fischer reagent solution
For Karl Fischer titration

packaging		Mfr. No
500 mℓ	Coated Glass bottle	J/4600/PB08
1 ℓ	Coated Glass bottle	J/4600/PB15
2.5 ℓ	Coated Glass bottle	J/4600/PB17
ADR/RID: 3, 6.1, II	H225, H301, H311, H331, H370	  
UN: 1992	P210, P240, P301+P310,	
GHS: Flam. Liq. 2, Acute	P304+P340, P302+P350,	
Tox. 3, STOT SE 1	P307+P311	

Guaranteed analysis
pH (10% aq. soln.) : 3.5 to 5.5

*Water equivalent approximately 5 mg / mℓ when prepared.
This solution deteriorates on storage even closed containers.*

Karl Fischer Aqualine® Electrolyte AG-H
For Karl Fischer titration by Coulometry

packaging		Mfr. No
500 mℓ	Glass bottle	K/2530/08
Flash point: 12°C	Tox. 3, STOT SE 1	   
ADR/RID: 3, 6.1, II	H225, H301, H311, H370, H331	
EINECS: 200-659-6	P280, P301+P310,	
UN: 1992	P302+P350, P304+P340,	
GHS: Flam. Liq. 2, Acute	P310, P210, P240	

Halogen-free anolyte solution for general use in conventional cells with a diaphragm.

Kieselguhr
pure, white

packaging		Mfr. No
1 kg	Plastic keg	K/0090/60
SiO ₂	STOT RE 2	 
CAS: 61790-53-2	H373, H319, H335	
MW: 60.08	P260, P280, P305+P351+P338	
GHS: Eye Irrit. 2, STOT SE 3,		

Guaranteed analysis
Identity test : Pass

Kieselguhr
pure, washed with acid

packaging		Mfr. No
500 g	Plastic jar	K/0100/53
1 kg	Plastic keg	K/0100/60
SiO ₂	STOT RE 2	 
CAS: 61790-53-2	H373, H319, H335	
MW: 60.08	P260, P280, P305+P351+P338	
GHS: Eye Irrit. 2, STOT SE 3,		

Guaranteed analysis
Complies with standard IR spectrum: Pass

Kjedhal Boric acid solution
pure, solution 4% alkaline with indicator

packaging		Mfr. No
1 ℓ	HDPE Plastic bottle	K/0200/15
5 ℓ	HDPE Plastic bottle	K/0200/21
Assay: 3.9 to 4.1		
Guaranteed analysis		
pH at 20°C : 3.8 to 4.2		

Kjeldahl catalyst Copper tablets
pure

Packaging		Mfr. No
500 PC	Plastic jar	K/0120/80
1000 PC	Plastic jar	K/0120/81
ADR/RID: 9, III	GHS: Aquatic Chronic 2	
EINECS: 231-820-9	H411	
UN: 3077	P273, P501	
Each tablet contains 1 g of sodium sulphate and the equivalent of 0.1 g of copper sulphate (II)		

Kjeldahl catalyst tablets Copper

packaging		Mfr. No
1000 PC	Plastic jar	K/0125/81
CAS: 7778-80-5	H412	
EINECS: 231-915-5	P273	
GHS: Aquatic Chronic 3		

Kjeldahl catalyst Mercury tablets
pure

packaging		Mfr. No
500 PC	Plastic jar	K/0130/80
ADR/RID: 6.1, III	Chronic 1	  
UN: 3288	H300, H310, H330, H373, H410	
GHS: Acute Tox. 1, STOT RE	P280, P260, P304+P340,	
2, Aquatic Acute 1, Aquatic	P301+P310, P273	

Each tablet contains 1 g of sodium sulphate and the equivalent of 0.1 g of mercury

Kjeldahl catalyst Selenium tablets
pure

Packaging		Mfr. No
500 PC	Plastic jar	K/0150/80
1 kg	Plastic jar	K/0150/81
EINECS: 231-820-9	H302, H332, H373	 
UN: 3283	P280, P260, P304+P341,	
GHS: Acute Tox. 4, STOT RE 2	P308+P313, P301+P312	

Kjeldahl catalyst Sodium sulfate tablets
pure

packaging		Mfr. No
500 PC	Plastic jar	K/0160/80
1000 PC	Plastic jar	K/0160/81
<i>Each tablet contains 1 g of sodium sulphate</i>		

Kjeldahl catalyst Titanium tablets
pure

packaging		Mfr. No
1000 PC	Plastic jar	K/0170/81
ADR/RID: 9, III	H411	
UN: 3077	P273	
GHS: Aquatic Chronic 2		
Guaranteed analysis		
Identity test : Pass		
<i>Each tablet contains 1 g of sodium sulphate and the equivalent of 0.03 g titanium dioxide and 0.03 g of copper sulphate</i>		

Kjeldahl Digestion reagent
pure

packaging		Mfr. No
2.5 ℓ	Coated Glass bottle	K/0210/PB17
ADR/RID: 8, II	P280, P301+P330+P331,	
UN: 1830	P310, P305+P351+P338,	
GHS: Skin Corr. 1A, Eye Dam. 1	P303+P361+P353	
H314		

Kjeldahl phenol reagent
pure, for determination of nitrogen

packaging		Mfr. No
1 ℓ	HDPE Plastic bottle	K/0220/15
ADR/RID: 8, 6.1, II	H341, H373	  
UN: 2922	P201, P280,	
GHS: Muta. 2, Acute Tox. 4,	P301+P330+P331,	
Skin Corr. 1A, Eye Dam. 1,	P305+P351+P338,	
STOT RE 2	P304+P340, P302+P350, P310	
H302, H312, H332, H314,		

Kjeldhal Sodium hydroxide/sodium thiosulfate solution
pure, for Kjel-Foss automatic nitrogen analyser

packaging		Mfr. No
5 ℓ	HDPE Plastic bottle	K/0310/21
ADR/RID: 8, II	P280, P301+P330+P331,	
UN: 1824	P310, P305+P351+P338,	
GHS: Skin Corr. 1A, Eye Dam. 1	P303+P361+P353	
H314		
Guaranteed analysis		
Sodium hydroxide : 39 to 41 Sodium thiosulfate : 7.7 to 8.3		
<i>Solution: Sodium Hydroxide at 400 g / ℓ and Sodium Thiosulfate 80 g / ℓ</i>		

Lacmoid indicator
pure

packaging		Mfr. No
1 g	Vial	L/0060/43
C ₁₂ H ₇ NO ₃ ·H ₃ N		
CAS: 42249-61-6		
MW: 230.22		
EINECS: 255-736-7		

pH Buffers

NIST Certified Standard Solution & Concentrate

A full range to discover on page 63



Lactic acid		Certified AR, Eur.Ph.	
for analysis			
packaging		Mfr. No	
500 mℓ	HDPE Plastic bottle	L/0150/PB08	
2.5 ℓ	HDPE Plastic bottle	L/0150/PB17	
C ₃ H ₆ O ₃	EINECS: 200-018-0		
CAS: 50-21-5	UN: 3265		
MW: 90.08	GHS: Skin Irrit. 2, Eye Dam.		
Density at 20°C: 1.2	1, STOT SE 3		
Boiling point: 122°C	H318, H335, H315		
Melting point: 18°C	P261, P280,		
Flash point: 110°C	P305+P351+P338, P302+P352		
ADR/RID: 8/III			
Assay: >88%			
Guaranteed analysis			
Colour	<10 APHA	Potassium (K)	<10 ppm
Reducing sugar	Pass	Sodium (Na)	<20 ppm
Calcium (Ca)	<10 ppm	Total chloride (Cl)	<0.001%
Copper (Cu)	<1 ppm	Total phosphorus (P)	<20 ppm
Iron (Fe)	<5 ppm	Total silicon (Si)	<20 ppm
Lead (Pb)	<1 ppm	Total sulphur (S)	<50 ppm
Magnesium (Mg)	<5 ppm	Zinc (Zn)	<10 ppm

Typical analysis			
Aluminium (Al)	<2 ppm	Manganese (Mn)	<1 ppm
Silver (Ag)	<1 ppm	Molybdenum (Mo)	<1 ppm
Barium (Ba)	<1 ppm	Nickel (Ni)	<1 ppm
Cadmium (Cd)	<1 ppm	Strontium (Sr)	<1 ppm
Chromium (Cr)	<1 ppm	Titanium (Ti)	<1 ppm
Cobalt (Co)	<1 ppm	Vanadium (V)	<1 ppm
Lithium (Li)	<1 ppm		

Lot analysis on the label.

Lactic acid		SLR	
specified reagents for general laboratory work- extra pure			
packaging		Mfr. No	
500 mℓ	HDPE Plastic bottle	L/0100/PB08	
2.5 ℓ	HDPE Plastic bottle	L/0100/PB17	
C ₃ H ₆ O ₃	EINECS: 200-018-0		
CAS: 50-21-5	UN: 3265		
MW: 90.08	GHS: Skin Irrit. 2, Eye Dam.		
Density at 20°C: 1.2	1, STOT SE 3		
Boiling point: 122°C	H318, H335, H315		
Melting point: 18°C	P261, P280,		
Flash point: 110°C	P305+P351+P338, P302+P352		
ADR/RID: 8/III			
Assay: >88%			
Guaranteed analysis			
Calcium (Ca)	<50 ppm	Potassium (K)	<50 ppm
Chromium (Cr)	<5 ppm	Sodium (Na)	<100 ppm
Copper (Cu)	<5 ppm	n-Butyl formate	<0.005%
Iron (Fe)	<10 ppm	Total phosphorus (P)	<100 ppm
Lead (Pb)	<5 ppm	Total silicon (Si)	<100 ppm
Magnesium (Mg)	<20 ppm	Total sulphur (S)	<200 ppm
Nickel (Ni)	<5 ppm	Zinc (Zn)	<50 ppm

Lactose monohydrate		Certified AR	
for analysis, fine powder			
packaging		Mfr. No	
500 g	Plastic jar	L/0250/53	
1 kg	Plastic jar	L/0250/60	
C ₁₂ H ₂₂ O ₁₁ .H ₂ O	Density at 20°C: 1.525		
CAS: 63-42-3	Melting point: 219°C		
MW: 360.32	EINECS: 200-559-2		
Guaranteed analysis			
Acidity	<0.005 meq/g	Potassium (K)	<100 ppm
Water	4.8 to 5.4%	Sodium (Na)	<100 ppm
Specific rotation	52° to 52.8°	n-Butyl formate	<0.005%
Soluble in alcohol	<0.2%	Total nitrogen (N)	<500 ppm
Calcium (Ca)	<50 ppm	Total phosphorus (P)	<200 ppm
Copper (Cu)	<1 ppm	Total silicon (Si)	<5 ppm
Iron (Fe)	<2 ppm	Total sulphur (S)	<50 ppm
Lead (Pb)	<1 ppm	Zinc (Zn)	<5 ppm
Magnesium (Mg)	<10 ppm		

Lactose monohydrate		SLR
extra pure, fine powder		
Packaging		Mfr. No
500 g	Plastic jar	L/0200/53
1 kg	Plastic jar	L/0200/60
5 kg	Plastic keg	L/0200/65
25 kg	Fibre Keg	L/0200/70
$C_{12}H_{22}O_{11} \cdot H_2O$		Density at 20°C: 1.525
CAS: 63-42-3		Melting point: 219°C
MW: 360.32		EINECS: 200-559-2
Guaranteed analysis		
Water	4.5 to 6%	Potassium (K) <500 ppm
Specific rotation	52° to 53°	Sodium (Na) <500 ppm
Calcium (Ca)	<200 ppm	n-Butyl formate <0.02%
Copper (Cu)	<5 ppm	Total phosphorus (P) <1000 ppm
Iron (Fe)	<10 ppm	Total silicon (Si) <20 ppm
Lead (Pb)	<5 ppm	Total sulphur (S) <100 ppm
Magnesium (Mg)	<50 ppm	Zinc (Zn) <20 ppm

Lanolin anhydrous		SLR
specified reagents for general laboratory work- extra pure		
Packaging		Mfr. No
500 g	Plastic jar	L/0350/53
1 kg	Plastic jar	L/0350/60
CAS: 8006-54-0		
Melting point: 38 to 40°C		
Flash point: 209°C		
EINECS: 232-348-6		
Guaranteed analysis		
Complies with standard IR spectrum: Pass		

Lanthanum chloride solution 10%w/v			
for AAS			
Packaging		Mfr. No	
50 mL	Glass bottle	L/0460/07	
1 L	Glass bottle	L/0460/15	
LaCl ₃ ·7H ₂ O		H335, H319, H315	
MW: 371.36		P305+P351+P338,	
GHS: Skin Irrit. 2, Eye Irrit. 2,		P302+P352, P261, P280	
STOT SE 3			
Assay: 9.5 to 10.5			
Guaranteed analysis			
Calcium (Ca)	<2 ppm	Lead (Pb)	<1 ppm
Copper (Cu)	<0.5 ppm	Magnesium (Mg)	<0.5 ppm
Iron (Fe)	<0.5 ppm	Zinc (Zn)	<0.5 ppm

Lanthanum standard metal solution 1000 ppm		for ICP solution in HNO ₃ 2%, NIST certified	
packaging		Mfr. No	
100 mℓ	Plastic bottle	SPPLLA22Y	
500 mℓ	Plastic bottle	SPPLLA22X	
H315, H319			

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Lanthanum standard metal solution 10 000 ppm		for ICP solution in HNO ₃ 5%, NIST certified	
packaging		Mfr. No	
100 mℓ	Plastic bottle	SPPLLA23Y	
500 mℓ	Plastic bottle	SPPLLA23X	
H314			

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Lanthanum chloride heptahydrate			
for AAS			
packaging		Mfr. No	
25 g	Plastic jar	L/0455/46	
100 g	Plastic jar	L/0455/48	
LaCl ₃ ·7H ₂ O	GHS: Skin Irrit. 2, Eye Irrit. 2, STOT SE 3		
CAS: 10025-84-0	H335, H315, H319		
MW: 371.36	P261, P302+P352, P280, P305+P351+P338		
Melting point: 91°C			
EINECS: 233-237-5			
Guaranteed analysis			
Calcium (Ca)	<10 ppm	Magnesium (Mg)	<1 ppm

Lanthanum chloride heptahydrate		Certified AR	
for analysis			
packaging		Mfr. No	
25 g	Plastic jar	L/0451/46	
100 g	Plastic jar	L/0451/48	
LaCl ₃ ·7H ₂ O	GHS: Skin Irrit. 2, Eye Irrit. 2, STOT SE 3		
CAS: 10025-84-0	H335, H315, H319		
MW: 371.36	P261, P302+P352, P280, P305+P351+P338		
Melting point: 91°C			
EINECS: 233-237-5			
Assay: >99%			
Guaranteed analysis			
Calcium (Ca)	<10 ppm	Sodium (Na)	<50 ppm
Copper (Cu)	<2 ppm	Total phosphorus (P)	<10 ppm
Iron (Fe)	<5 ppm	Total silicon (Si)	<5 ppm
Lead (Pb)	<2 ppm	Total sulphur (S)	<20 ppm
Magnesium (Mg)	<5 ppm	Zinc (Zn)	<2 ppm
Potassium (K)	<20 ppm		
Lot analysis on the label.			

Lanthanum nitrate hexahydrate		SLR
specified reagents for general laboratory work- extra pure		
packaging		Mfr. No
25 g	Plastic jar	L/0500/46
100 g	Plastic jar	L/0500/48
La(NO ₃) ₃ .6H ₂ O	EINECS: 233-238-0	
CAS: 10277-43-7	UN: 1477	
MW: 433.03	GHS: Ox. Sol. 2, Skin Irrit. 2,	
Density at 20°C: >> 1.0	Eye Irrit. 2, STOT SE 3	
Boiling point: 126°C	H272, H315, H319, H335	
Melting point: 40°C	P210, P261, P302+P352,	
ADR/RID: 5.1, II	P280, P305+P351+P338	
Assay: >99%		

Guaranteed analysis			
Identity test	: Pass	Manganese (Mn)	<=5 ppm
Cadmium (Cd)	<=5 ppm	Nickel (Ni)	<=5 ppm
Calcium (Ca)	<=20 ppm	Potassium (K)	<=10 ppm
Chromium (Cr)	<=5 ppm	Sodium (Na)	<=10 ppm
Cobalt (Co)	<=5 ppm	Total phosphorus (P)	<=50 ppm
Copper (Cu)	<=5 ppm	Total silicon (Si)	<=20 ppm
Iron (Fe)	<=10 ppm	Total sulphur (S)	<=50 ppm
Lead (Pb)	<=10 ppm	Zinc (Zn)	<=5 ppm
Magnesium (Mg)	<=5 ppm		

Lanthanum oxide		SLR
specified reagents for general laboratory work- extra pure		
packaging		Mfr. No
100 g	Plastic jar	L/0510/48
La ₂ O ₃	Melting point: 2315°C	
CAS: 1312-81-8	EINECS: 215-200-5	
MW: 325.82	GHS: Eye Irrit. 2	
Density at 20°C: 6.51	H319	
Boiling point: 4200°C	P280, P305+P351+P338	
Guaranteed analysis		
Identity test		: Pass

Lauric acid		SLR
specified reagents for general laboratory work- extra pure		
packaging		Mfr. No
500 g	Plastic jar	L/0560/53
C ₁₂ H ₂₄ O ₂	Flash point: 156°C	
CAS: 143-07-7	EINECS: 205-582-1	
MW: 200.32	GHS: Skin Irrit. 2, Eye Irrit. 2	
Density at 20°C: 0.8830	H319, H315	
Boiling point: 225°C	P280, P302+P352,	
Melting point: 44 to 46°C	P305+P351+P338	
Assay: 98 to 101%		

di-Lauroyl peroxide		
pure		
packaging		Mfr. No
100 g	Plastic jar	L/0565/48
Pack	Vial	L/0565/99
C ₂₄ H ₄₆ O ₄	EINECS: 203-326-3	
CAS: 105-74-8	UN: 3106	
MW: 398.63	GHS: Org. Perox. D	
Melting point: 53 to 57°C	H242	
ADR/RID: 5.2, II	P280, P210, P220	
Guaranteed analysis		
Complies with standard IR spectrum: Pass		

Lead standard metal solution 1000 ppm		
for AAS, solution HNO ₃ 1M		
packaging		Mfr. No
100 ml	HDPE Plastic bottle	J/8035/05
500 ml	HDPE Plastic bottle	J/8035/08
Boiling point: 100°C	H314	
ADR/RID: 8, II	P280, P305+P351+P338,	
UN: 2031	P301+P330+P331, P310	
GHS: Skin Corr. 1B		
Manufactured to a tolerance of ± 0.5%.		
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.		

Lead standard metal solution 1000 ppm
for ICP, solution HNO₃ 1M

packaging		Mfr. No
500 mℓ	HDPE Plastic bottle	J/8276/08
Boiling point: 100°C	H314	
ADR/RID: 8, II	P280, P305+P351+P338,	
UN: 2031	P301+P330+P331, P310	
GHS: Skin Corr. 1B		
<i>Manufactured to a tolerance of ± 0.5%.</i>		
<i>All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.</i>		
<i>See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.</i>		

Lead standard metal solution 10 000 ppm
for ICP, solution HNO₃ 1M

packaging		Mfr. No
100 mℓ	HDPE Plastic bottle	J/8275/05
Density at 20°C: >1.0	H360Df, H332, H302, H373,	
ADR/RID: 6.1, 8, II	H412	
UN: 3289	P280, P301+P330+P331,	
GHS: Repr. 1A, Acute Tox. 4,	P310, P305+P351+P338,	
Skin Corr. 1B, STOT RE 2	P304+P340, P308+P313	
Manufactured to a tolerance of ± 0.5%. All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities. See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.		

Lead standard metal solution 1000 ppm
for ICP solution in HNO₃ 2%, NIST certified

packaging		Mfr. No
100 mℓ	Plastic bottle	SPPLPB22Y
250 mℓ	Plastic bottle	SPPLPB22T
500 mℓ	Plastic bottle	SPPLPB22X
H315, H319		

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Lead standard metal solution 10 000 ppm
for ICP solution in HNO₃ 5%, NIST certified

packaging		Mfr. No
100 mℓ	Plastic bottle	SPPLPB23Y
500 mℓ	Plastic bottle	SPPLPB23X
H314, H360, H373, H412		

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Lead standard metal solution 1000 ppm
for ICP-MS solution in HNO₃ 2%, NIST certified

packaging		Mfr. No
100 mℓ	Plastic bottle	SPCLPB22Y
H315, H320		

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Lead metal powder
extra pure, below 200 mesh

packaging		Mfr. No
500 g	Plastic jar	L/0750/53
Pb	UN: 3077	
CAS: 7439-92-1	GHS: Repr. 1A, Acute Tox. 4,	
MW: 207.19	STOT RE 2, Aquatic Acute 1,	
Density at 20°C: 11.3	Aquatic Chronic 1	
Boiling point: 1740°C	H332, H410, H373, H302,	
Melting point: 327.4°C	H360Df	
ADR/RID: 9, III	P281, P260, P273,	
EINECS: 231-100-4	P301+P312, P304+P340	

Guaranteed analysis	
Identity test	: Pass

Lead metal foil
extra pure, about 0.13 mm

packaging		Mfr. No
500 g	Plastic jar	L/0610/53
Pb	UN: 3077	
CAS: 7439-92-1	GHS: Repr. 1A, Acute Tox. 4,	
MW: 207.19	STOT RE 2, Aquatic Acute 1,	
Density at 20°C: 11.3	Aquatic Chronic 1	
Boiling point: 1740°C	H332, H410, H373, H302,	
Melting point: 327.4°C	H360Df	
ADR/RID: 9, III	P281, P260, P273,	
EINECS: 231-100-4	P301+P312, P304+P340	

Guaranteed analysis	
Identity test	: Pass

Lead metal shot
extra pure, about 2 mm

packaging		Mfr. No
500 g	Plastic jar	L/0800/53
1 kg	Plastic jar	L/0800/60
Pb	UN: 3077	
CAS: 7439-92-1	GHS: Repr. 1A, Acute Tox. 4,	
MW: 207.19	STOT RE 2, Aquatic Acute 1,	
Density at 20°C: 11.3	Aquatic Chronic 1	
Boiling point: 1740°C	H332, H410, H373, H302,	
Melting point: 327.4°C	H360Df	
ADR/RID: 9, III	P281, P260, P273,	
EINECS: 231-100-4	P301+P312, P304+P340	

Guaranteed analysis	
Identity test	: Pass

Lead metal granulated
pure

packaging		Mfr. No
500 g	Plastic jar	L/0700/53
Pb	UN: 3077	
CAS: 7439-92-1	GHS: Repr. 1A, Acute Tox. 4,	
MW: 207.19	STOT RE 2, Aquatic Acute 1,	
Density at 20°C: 11.3	Aquatic Chronic 1	
Boiling point: 1740°C	H332, H410, H373, H302,	
Melting point: 327.4°C	H360Df	
ADR/RID: 9, III	P281, P260, P273,	
EINECS: 231-100-4	P301+P312, P304+P340	

Guaranteed analysis	
Copper (Cu)	<200 ppm
Iron (Fe)	<10 ppm

Lead metal
technical, sheet 1mm

packaging		Mfr. No
1 kg	Plastic jar	L/0605/60
Pb	UN: 3077	
CAS: 7439-92-1	GHS: Repr. 1A, Acute Tox. 4,	
MW: 207.19	STOT RE 2, Aquatic Acute 1,	
Density at 20°C: 11.3	Aquatic Chronic 1	
Boiling point: 1740°C	H332, H410, H373, H302,	
Melting point: 327.4°C	H360Df	
ADR/RID: 9, III	P281, P260, P273,	
EINECS: 231-100-4	P301+P312, P304+P340	

Guaranteed analysis	
Identity test	: Pass

Lead(II) acetate trihydrate
for analysis

packaging		Mfr. No
250 g	Plastic jar	L/1000/50
500 g	Plastic jar	L/1000/53
1 kg	Plastic jar	L/1000/60
Pb(CH ₃ CO ₂) ₂ ·3H ₂ O	GHS: Repr. 1A, STOT RE 2,	
CAS: 6080-56-4	Aquatic Acute 1, Aquatic	
MW: 379.33	Chronic 1	
Melting point: 75°C	H360Df, H373, H410	
ADR/RID: 6.1, III	P281, P260, P273,	
UN: 1616	P301+P312, P304+P340	

Assay: >99.5%	
Calcium (Ca)	<5 ppm
Copper (Cu)	<5 ppm
Iron (Fe)	<5 ppm
Magnesium (Mg)	<5 ppm
Potassium (K)	<10 ppm
Sodium (Na)	<20 ppm
Total chloride (Cl)	<0.0005%
Total nitrogen (N)	<20 ppm
Total phosphorus (P)	<20 ppm
Total silicon (Si)	<20 ppm
Total sulphur (S)	<50 ppm
Zinc (Zn)	<5 ppm

Typical analysis	
Aluminium (Al)	<1 ppm
Arsenic (As)	<2 ppm
Barium (Ba)	<1 ppm
Cadmium (Cd)	<1 ppm
Chromium (Cr)	<1 ppm
Cobalt (Co)	<1 ppm
Lithium (Li)	<1 ppm
Manganese (Mn)	<1 ppm
Nickel (Ni)	<1 ppm
Titanium (Ti)	<1 ppm
Vanadium (V)	<1 ppm

Lot analysis on the label.

Fisher Chemical 5L packaging change

Our improved HDPE 5L packaging is easier to handle, more robust and preserves product quality for longer.



Lead(II) acetate trihydrate
specified reagents for general laboratory work- extra pure

packaging		Mfr. No
250 g	Plastic jar	L/0950/50
500 g	Plastic jar	L/0950/53
1 kg	Plastic jar	L/0950/60
Pb(CH ₃ CO ₂) ₂ ·3H ₂ O	GHS: Repr. 1A, STOT RE 2,	
CAS: 6080-56-4	Aquatic Acute 1, Aquatic	
MW: 379.33	Chronic 1	
Melting point: 75°C	H360Df, H373, H410	
ADR/RID: 6.1, III	P281, P260, P273,	
UN: 1616	P301+P312, P304+P340	

Assay: >98%	
Calcium (Ca)	<100 ppm
Copper (Cu)	<50 ppm
Iron (Fe)	<50 ppm
Magnesium (Mg)	<50 ppm
Potassium (K)	<100 ppm
Sodium (Na)	<200 ppm
n-Butyl formate	<0.005%
Total phosphorus (P)	<100 ppm
Total silicon (Si)	<100 ppm
Total sulphur (S)	<200 ppm
Zinc (Zn)	<50 ppm

Lead(II) bromide
pure

packaging		Mfr. No
100 g	Plastic jar	L/1100/48
250 g	Plastic jar	L/1100/50
PbBr ₂	UN: 2291	
CAS: 10031-22-8	GHS: Repr. 1A, Acute Tox. 4,	
MW: 367.01	STOT RE 2, Aquatic Acute 1,	
Density at 20°C: 6.66	Aquatic Chronic 1	
Boiling point: 916°C	H410, H332, H373, H360Df,	
Melting point: 373°C	H302	
ADR/RID: 6.1, III	P281, P260, P273, P304+P340	
EINECS: 233-084-4		

Guaranteed analysis	
Identity test	: Pass

Lead(II) carbonate basic
extra pure, precipitated

packaging		Mfr. No
250 g	Plastic jar	L/1150/50
500 g	Plastic jar	L/1150/53
(PbCO ₃) ₂ ·Pb(OH) ₂	GHS: Repr. 1A, Acute Tox. 4,	
CAS: 1319-46-6	STOT RE 2, Aquatic Acute 1,	
MW: 775.59	Aquatic Chronic 1	
Melting point: 400°C	H332, H373, H410, H360Df,	
ADR/RID: 6.1, III	H302	
EINECS: 215-290-6	P281, P260, P273, P304+P340	
UN: 3288		

Guaranteed analysis	
Identity test	: Pass

Safe, environmentally friendly and effective packaging

Fisher Chemical outer packaging is UN compliant and 100% recyclable. It's an easy and safe way to reduce your waste management costs and contribute to a cleaner environment.



Lead(II) chloride		SLR
<i>specified reagents for general laboratory work- extra pure</i>		
packaging	Mfr. No	
250 g Plastic jar	L/1200/50	
PbCl ₂ CAS: 7758-95-4 MW: 278.09 Boiling point: 950°C Melting point: 501°C ADR/RID: 6.1, III EINECS: 231-845-5	UN: 2291 GHS: Repr. 1A, Acute Tox. 4, STOT RE 2, Aquatic Acute 1, Aquatic Chronic 1 H332, H302, H373, H410, H360Df P281, P308+P313, P260, P273	
Guaranteed analysis		
Identity test : Pass		

Lead(II) acetate basic	
<i>pure, for sugar analysis</i>	
packaging	Mfr. No
500 g Plastic jar	L/1050/53
25 kg Plastic keg	L/1050/70
(CH ₃ COO) ₂ Pb.2Pb(OH) ₂ CAS: 1335-32-6 MW: 807.69 ADR/RID: 6.1, III EINECS: 215-630-3 UN: 1616	GHS: Carc. 2, Repr. 1A, STOT RE 2, Aquatic Acute 1, Aquatic Chronic 1 H351, H373, H410, H360Df P281, P260, P273
Assay: >95%	
Guaranteed analysis	
Copper (Cu) <100 ppm Nitrate (NO ₃) <0.01% Iron (Fe) <100 ppm	

Lead(II) sulfide	
<i>pure</i>	
packaging	Mfr. No
100 g Glass jar	L/1730/48
PbS CAS: 1314-87-0 MW: 239.26 ADR/RID: 6.1, III EINECS: 215-246-6 UN: 3288	GHS: Repr. 2, Acute Tox. 4, STOT RE 2, Aquatic Acute 1, Aquatic Chronic 1 H410, H373, H360Df, H302, H332 P281, P260, P273, P301+P312
Guaranteed analysis	
Identity test : Pass	

Lead(II) iodide		SLR
<i>specified reagents for general laboratory work- extra pure</i>		
packaging	Mfr. No	
100 g Glass jar	L/1400/48	
PbI ₂ CAS: 10101-63-0 MW: 461.00 Density at 20°C: 6.16 Melting point: 402°C ADR/RID: 6.1, III EINECS: 233-256-9	UN: 2291 GHS: Repr. 1A, Acute Tox. 4, STOT RE 2, Aquatic Acute 1, Aquatic Chronic 1 H373, H410, H332, H302, H360Df P281, P260, P273, P304+P340	
Guaranteed analysis		
Identity test : Pass		

Lead(II) nitrate		Certified AR
<i>for analysis</i>		
packaging	Mfr. No	
250 g Plastic jar	L/1500/50	
500 g Plastic jar	L/1500/53	
1 kg Plastic jar	L/1500/60	
Pb(NO ₃) ₂ CAS: 10099-74-8 MW: 331.20 Density at 20°C: 4.53 Melting point: 470°C ADR/RID: 5.1, 6.1, II EINECS: 233-245-9 UN: 1469	GHS: Ox. Sol. 2, Repr. 1A, Acute Tox. 4, STOT RE 2, Aquatic Acute 1, Aquatic Chronic 1 H272, H332, H373, H410, H302, H360Df P281, P210, P301+P312, P260, P273, P304+P340	

Assay: >99.5%			
Guaranteed analysis			
Calcium (Ca)	<5 ppm	Total chloride (Cl)	<0.0005%
Copper (Cu)	<5 ppm	Total phosphorus (P)	<20 ppm
Iron (Fe)	<5 ppm	Total silicon (Si)	<20 ppm
Magnesium (Mg)	<5 ppm	Total sulphur (S)	<50 ppm
Potassium (K)	<10 ppm	Zinc (Zn)	<5 ppm
Sodium (Na)	<20 ppm		
Typical analysis			
Aluminium (Al)	<1 ppm	Lithium (Li)	<1 ppm
Arsenic (As)	<2 ppm	Manganese (Mn)	<1 ppm
Barium (Ba)	<1 ppm	Nickel (Ni)	<1 ppm
Cadmium (Cd)	<1 ppm	Titanium (Ti)	<1 ppm
Chromium (Cr)	<1 ppm	Vanadium (V)	<1 ppm
Cobalt (Co)	<1 ppm		
<i>Lot analysis on the label.</i>			

Lead(II) nitrate		SLR
<i>specified reagents for general laboratory work- extra pure</i>		
packaging	Mfr. No	
250 g Plastic jar	L/1450/50	
500 g Plastic jar	L/1450/53	
1 kg Plastic jar	L/1450/60	
Pb(NO ₃) ₂ CAS: 10099-74-8 MW: 331.20 Density at 20°C: 4.53 Melting point: 470°C ADR/RID: 5.1, 6.1, II EINECS: 233-245-9 UN: 1469	GHS: Ox. Sol. 2, Repr. 1A, Acute Tox. 4, STOT RE 2, Aquatic Acute 1, Aquatic Chronic 1 H272, H332, H373, H410, H302, H360Df P281, P210, P301+P312, P260, P273, P304+P340	

Assay: >99%			
Guaranteed analysis			
Calcium (Ca)	<100 ppm	n-Butyl formate	<0.005%
Copper (Cu)	<20 ppm	Total phosphorus (P)	<100 ppm
Iron (Fe)	<20 ppm	Total silicon (Si)	<100 ppm
Magnesium (Mg)	<20 ppm	Total sulphur (S)	<200 ppm
Potassium (K)	<100 ppm	Zinc (Zn)	<50 ppm
Sodium (Na)	<200 ppm		

UHPLC

Filtered to 0.1 µm, Fisher Chemical UHPLC gradient grade solvents (Acetonitrile, Methanol and Water) are perfect for those working at high pressure.

See page 309 for more information

Lead(II) nitrate		technical
packaging	Mfr. No	
500 g Plastic jar	L/1445/53	
Pb(NO ₃) ₃ CAS: 10099-74-8 MW: 331.20 Density at 20°C: 4.53 Melting point: 470°C ADR/RID: 5.1, 6.1, II EINECS: 233-245-9 UN: 1469	GHS: Ox. Sol. 2, Repr. 1A, Acute Tox. 4, STOT RE 2, Aquatic Acute 1, Aquatic Chronic 1 H272, H332, H373, H410, H302, H360Df P281, P210, P301+P312, P260, P273, P304+P340	
Assay: >99%		

Lead(II) oxide		SLR	
<i>specified reagents for general laboratory work- extra pure</i>			
packaging	Mfr. No		
500 g Plastic jar	L/1600/53		
PbO CAS: 1317-36-8 MW: 223.19 Boiling point: 1470°C Melting point: 886°C ADR/RID: 6.1, III EINECS: 215-267-0	UN: 2291 GHS: Repr. 1A, Acute Tox. 4, STOT RE 2, Aquatic Acute 1, Aquatic Chronic 1 H332, H373, H410, H360Df, H302 P281, P260, P273, P304+P340	 	
Assay: >98%			
Guaranteed analysis			
Calcium (Ca)	<100 ppm	n-Butyl formate	<0.02%
Copper (Cu)	<50 ppm	Total phosphorus (P)	<100 ppm
Iron (Fe)	<50 ppm	Total silicon (Si)	<100 ppm
Magnesium (Mg)	<50 ppm	Total sulphur (S)	<200 ppm
Potassium (K)	<100 ppm	Zinc (Zn)	<50 ppm
Sodium (Na)	<200 ppm		

Lead(II) oxide		Technical, canary yellow
packaging	Mfr. No	
500 g Plastic jar	L/1670/53	
PbO CAS: 1317-36-8 MW: 223.19 Boiling point: 1470°C Melting point: 886°C ADR/RID: 6.1, III EINECS: 215-267-0	UN: 2291 GHS: Repr. 1A, Acute Tox. 4, STOT RE 2, Aquatic Acute 1, Aquatic Chronic 1 H332, H373, H410, H360Df, H302 P281, P260, P273, P304+P340	
Guaranteed analysis		
Identity test : Pass		

Optima LC/MS

Fisher Chemical OPTIMA® LC/MS solvents, blends, additives and modifiers have set the standard of excellence for consistent, reproducible performance in the mobile phase of LC/MS.

Find out more on page 308

Lead (II,V)oxide		SLR
<i>extra pure, minium, red lead</i>		
packaging	Mfr. No	
500 g Plastic jar	L/1700/53	
Pb ₃ O ₄ CAS: 1314-41-6 MW: 685.57 Density at 20°C: 8.32 to 9.16 Melting point: 500 to 530°C ADR/RID: 5.1, III EINECS: 215-235-6 UN: 1479	GHS: Ox. Sol. 2, Repr. 1A, Acute Tox. 4, STOT RE 2, Aquatic Acute 1, Aquatic Chronic 1 H272, H373, H332, H302, H361, H410 P281, P308+P313, P260, P273, P210	
Guaranteed analysis		
Identity test : Pass		

Lead(IV) dioxide		SLR
<i>specified reagents for general laboratory work- extra pure</i>		
packaging	Mfr. No	
250 g Plastic jar	L/1350/50	
500 g Plastic jar	L/1350/53	
PbO ₂ CAS: 1309-60-0 MW: 239.19 Melting point: 290°C ADR/RID: 5.1, 6.1, III EINECS: 215-174-5 UN: 1872 GHS: Ox. Sol. 2, Repr. 1A,	Acute Tox. 4, STOT RE 2, Aquatic Acute 1, Aquatic Chronic 1 H271, H410, H332, H373, H360Df, H302 P281, P301+P312, P260, P273, P304+P340, P210	
Assay: >94%		
Guaranteed analysis		
Sulfate (SO ₄) <0.01%		

Lead(II) sulfate		SLR
<i>specified reagents for general laboratory work- extra pure</i>		
packaging	Mfr. No	
500 g Plastic jar	L/1720/53	
PbSO ₄ CAS: 7446-14-2 MW: 303.25 Melting point: 1170°C ADR/RID: 6.1, III EINECS: 231-198-9 UN: 3288	GHS: Repr. 1A, Acute Tox. 4, STOT RE 2, Aquatic Acute 1, Aquatic Chronic 1 H302, H410, H332, H373, H360Df P281, P260, P273, P301+P312, P304+P340	
Guaranteed analysis		
Identity test : Pass		

Lecithin	
<i>technical, egg about 90%</i>	
packaging	Mfr. No
25 g Plastic jar	L/1814/46
CAS: 8002-43-5 Density at 20°C: 1.03 Melting point: 236.1°C EINECS: 297-639-2	
Guaranteed analysis	
Complies with standard IR spectrum: Pass	

180

Application For analysis Specified Titration

181

Leishman’s stain <i>pure</i>		
packaging	Mfr. No	
10 g Glass jar	L/1815/45	
C ₆ H ₅ NO ₃ CAS: 12627-53-1 MW: 171.04 Melting point: 79 to 81°C EINECS: 235-732-1	GHS: Skin Irrit. 2, Eye Irrit. 2 H315, H319 P302+P352, P280, P305+P351+P338	
Guaranteed analysis		
Identity test : Pass		

Leishman’s stain solution <i>pure, in 100% Methanol</i>		
packaging	Mfr. No	
100 mℓ Coated Glass bottle	L/1815L/PB05	
Flash point: 10°C ADR/RID: 3, 6.1, II EINECS: 200-659-6 UN: 1992 GHS: Flam. Liq. 2, Acute	Tox. 3, STOT SE 1 H225, H301, H311, H331, H370 P210, P240, P301+P310, P304+P340, P302+P352, P308+P313	
Guaranteed analysis		
Identity test : Pass		

L-Leucine <i>specified reagents for general laboratory work- extra pure</i>		SLR
packaging	Mfr. No	
100 g Plastic jar	L/1819/48	
C ₆ H ₁₃ NO ₂ CAS: 61-90-5 MW: 131.18	Melting point: 286-288°C EINECS: 200-522-0	
Assay: >99%		

L-iso-Leucine <i>pure</i>		
packaging	Mfr. No	
25 g Plastic jar	L/1822/46	
C ₆ H ₁₃ NO ₂ CAS: 73-32-5 MW: 131.18	Melting point: 168-170°C EINECS: 200-798-2	
Assay: >99%		

Light green SF C.I 42095 <i>pure</i>		
packaging	Mfr. No	
25 g Glass jar	L/1840/46	
C ₃₇ H ₃₄ N ₂ Na ₇ O ₉ S ₃ CAS: 5141-20-8 MW: 792.84	Melting point: 288°C EINECS: 225-906-5	
Guaranteed analysis		
Complies with standard IR spectrum: Pass		

Lime water <i>pure</i>		
packaging	Mfr. No	
5 ℓ HDPE Plastic bottle	J/5000/21	
Density at 20°C: 1.0 Boiling point: 100°C Melting point: 0°C		

D(+) Limonene <i>pure</i>		
packaging	Mfr. No	
250 mℓ Glass bottle	L/1860/07	
C ₁₀ H ₁₆ CAS: 5989-27-5 MW: 136.24 Density at 20°C: 0.841 Boiling point: 175.4 to 176°C Melting point: -74°C Flash point: 48°C ADR/RID: 3, III	EINECS: 227-813-5 UN: 2052 GHS: Flam. Liq. 3, Skin Irrit. 2, Skin Sens. 1, Aquatic Acute 1, Aquatic Chronic 1 H226, H315, H317, H410 P280, P273, P302+P352, P210	

Linseed oil <i>pure, raw approx.d= 0,93</i>		
packaging	Mfr. No	
500 mℓ Glass bottle 2.5 ℓ Glass bottle	L/1950/08 L/1950/17	
CAS: 8001-26-1 Density at 20°C: 0.925 to 0.935 Boiling point: 315.6°C Melting point: -23.9°C Flash point: 206.1°C EINECS: 232-278-6	GHS: Skin Irrit. 2, Eye Irrit. 2, STOT SE 3 H319, H335, H315 P280, P261, P302+P352, P304+P340, P305+P351+P338, P312	
Guaranteed analysis		
Density at 15.5°C : 0.92 to 0.94		

Linseed oil <i>pure, boiled approx.d= 0.95</i>		
packaging	Mfr. No	
2.5 ℓ Glass bottle	L/1900/17	
CAS: 8001-26-1 EINECS: 232-278-6 GHS: Skin Irrit. 2, Eye Irrit. 2, STOT SE 3	H319, H335, H315 P280, P261, P302+P352, P304+P340, P305+P351+P338, P312	
Guaranteed analysis		
Density at 15.5°C : 0.94 to 0.96		

Lithium standard metal solution 1000 ppm <i>for AAS, solution HNO₃ 1M</i>		
packaging	Mfr. No	
100 mℓ HDPE Plastic bottle 500 mℓ HDPE Plastic bottle	J/8037/05 J/8037/08	
ADR/RID: 8, II UN: 2031 GHS: Skin Corr. 1B, Eye Dam. 1	H314 P280, P305+P351+P338, P261, P302+P352	
<i>Manufactured to a tolerance of ± 0.5%. See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.</i>		

Lithium standard metal solution 1000 ppm <i>for ICP, solution HNO₃ 1M</i>		
packaging	Mfr. No	
500 mℓ HDPE Plastic bottle	J/8281/08	
ADR/RID: 8, II UN: 2031 GHS: Skin Corr. 1B, Eye Dam. 1	H314 P280, P305+P351+P338, P261, P302+P352	
<i>Manufactured to a tolerance of ± 0.5%. All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities. See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.</i>		

Lithium standard metal solution 10 000 ppm <i>for ICP, solution HNO₃ 1M</i>		
packaging	Mfr. No	
100 mℓ HDPE Plastic bottle	J/8280/05	
ADR/RID: 8, II UN: 2031 GHS: Skin Corr. 1B, Eye Dam. 1	P280, P301+P330+P331, P310, P305+P351+P338, P303+P361+P353	
<i>Manufactured to a tolerance of ± 0.5%. All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities. See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.</i>		

Lithium standard metal solution 1000 ppm <i>for ICP solution in HNO₃ 2%, NIST certified</i>		
packaging	Mfr. No	
100 mℓ Plastic bottle 500 mℓ Plastic bottle H315, H319	SPPLLI22Y SPPLLI22X	
<i>Manufactured to a tolerance of ± 0.5%. All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities. See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314. Check with your local distributor for availability of Spex product</i>		

Lithium standard metal solution 1000 ppm <i>for ICP solution in HCL 2%, NIST certified</i>		
packaging	Mfr. No	
500 mℓ Plastic bottle H314	SPPLLI12X	
<i>Manufactured to a tolerance of ± 0.5%. All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities. See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314. Check with your local distributor for availability of Spex product</i>		

Reagents, Solvents & Volumetric solutions comply with the European Pharmacopeia

Read more on page 234

Lithium standard metal solution 10 000 ppm <i>for ICP solution in HNO₃ 5%, NIST certified</i>		
packaging	Mfr. No	
100 mℓ Plastic bottle 500 mℓ Plastic bottle H314	SPPLLI23Y SPPLLI23X	
<i>Manufactured to a tolerance of ± 0.5%. All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities. See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314. Check with your local distributor for availability of Spex product</i>		

Lithium standard metal solution 10 000 ppm <i>for ICP solution in HCL 5%, NIST certified</i>		
packaging	Mfr. No	
500 mℓ Plastic bottle H303	SPPLLI13X	
<i>Manufactured to a tolerance of ± 0.5%. All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities. See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314. Check with your local distributor for availability of Spex product</i>		

Lithium standard solution 1000 ppm (1 mg/mℓ) <i>for Ion pair chromatography, Ready-to-use</i>		
packaging	Mfr. No	
100 mℓ HDPE Plastic bottle	J/4550/05	
<i>Manufactured to a tolerance of ± 1%.</i>		

Lithium metal in liquid paraffin <i>specified reagents for general laboratory work- extra pure</i>		SLR
packaging	Mfr. No	
10 g Aluminium jar	L/2000/45	
Li CAS: 7439-93-2 MW: 6.94 Density at 20°C: 0.534 Boiling point: 1317.2°C Melting point: 180.6°C ADR/RID: 4.3, I	EINECS: 231-102-5 UN: 1415 GHS: Water-react. 1, Skin Corr. 1B, Eye Dam. 1 H260, H314, EUH014 P280, P301+P330+P331, P402+P404, P223, P310	
Guaranteed analysis		
Identity test : Pass		

Lithium acetate dihydrate <i>specified reagents for general laboratory work- extra pure</i>		SLR
packaging	Mfr. No	
250 g Plastic jar	L/2050/50	
C ₂ H ₃ LiO ₂ ·2H ₂ O CAS: 6108-17-4 MW: 102.02 Melting point: 280 to 285°C	EINECS: 208-914-3 GHS: Eye Irrit. 2 H319 P280, P305+P351+P338	
Guaranteed analysis		
Identity test : Pass		

Lithium aluminium hydride
pure

packaging		Mfr. No
10 g	Plastic jar	A/2160/45
LiAlH ₄	UN: 1410	 
CAS: 16853-85-3	GHS: Water-react. 1, Skin	
MW: 37.95	Corr. 1A, Eye Dam. 1	
Density at 20°C: 0.910	H260, H314	
Melting point: 125°C	P301+P330+P331, P280,	
ADR/RID: 4.3, I	P305+P351+P338, P310	
EINECS: 240-877-9		
Assay: >95%		

Lithium bromide anhydrous
specified reagents for general laboratory work- extra pure

packaging		Mfr. No
100 g	Plastic jar	L/2080/48
500 g	Plastic jar	L/2080/53
LiBr	Melting point: 550°C	
CAS: 7550-35-8	EINECS: 231-439-8	
MW: 86.85	GHS: Acute Tox. 4	
Density at 20°C: 3.46	H302	
Boiling point: 1265°C	P301+P312	
Assay: >98%		
Guaranteed analysis		
Iron (Fe)	<50 ppm	Sulfate (SO ₄)
Lead (Pb)	<50 ppm	<0.02%

Lithium carbonate
for analysis

Packaging		Mfr. No	
100 g	Plastic jar	L/2105/48	
Li ₂ CO ₃		 	
CAS: 554-13-2			
MW: 73.89			
Density at 20°C: 2.11			
Melting point: 720°C			
EINECS: 209-062-5			
GHS: Acute Tox. 4, Skin Irrit. 2, Eye Dam. 1, STOT SE 3 H302, H315, H318, H335 P305+P351+P338, P310, P280, P261, P302+P352			
Assay: >99%			
Guaranteed analysis			
Calcium (Ca)	<20 ppm	n-Butyl formate	<0.002%
Copper (Cu)	<5 ppm	Total nitrogen (N)	<10 ppm
Iron (Fe)	<5 ppm	Total phosphorus (P)	<20 ppm
Lead (Pb)	<5 ppm	Total silicon (Si)	<20 ppm
Magnesium (Mg)	<5 ppm	Total sulphur (S)	<50 ppm
Potassium (K)	<10 ppm	Zinc (Zn)	<10 ppm
Sodium (Na)	<50 ppm		

Lot analysis on the label.

Lithium carbonate
specified reagents for general laboratory work- extra pure

Packaging		Mfr. No
250 g	Plastic jar	L/2100/50
Li ₂ CO ₃	GHS: Acute Tox. 4, Skin Irrit.	
CAS: 554-13-2	2, Eye Dam. 1, STOT SE 3	
MW: 73.89	H302, H315, H318, H335	
Density at 20°C: 2.11	P305+P351+P338, P310,	
Melting point: 720°C	P280, P261, P302+P352	
EINECS: 209-062-5		
Assay: >99%		
Guaranteed analysis		
Calcium (Ca)	<200 ppm	Sodium (Na)
Copper (Cu)	<20 ppm	n-Butyl formate
Iron (Fe)	<20 ppm	Total phosphorus (P)
Lead (Pb)	<20 ppm	Total silicon (Si)
Magnesium (Mg)	<20 ppm	Total sulphur (S)
Potassium (K)	<200 ppm	Zinc (Zn)

Lithium chloride anhydrous
for analysis

Packaging			Mfr. No
250 g	Plastic jar		L/2201/50
5 kg	Plastic jar		L/2201/65
LiCl		EINECS: 231-212-3	
CAS: 7447-41-8		GHS: Acute Tox. 4, Skin Irrit.	
MW: 42.38		2, Eye Irrit. 2, STOT SE 3	
Density at 20°C: 2.06		H319, H315, H302, H312, H335	
Boiling point: 1382°C		P280, P302+P352,	
Melting point: 605°C		P301+P312, P305+P351+P338	
Assay: >99.5%			
Guaranteed analysis			
pH (aqueous solution 5%)	6 to 9	Potassium (K)	<50 ppm
Loss on drying	<2%	Sodium (Na)	<50 ppm
Calcium (Ca)	<20 ppm	Total nitrogen (N)	<10 ppm
Copper (Cu)	<5 ppm	Total phosphorus (P)	<20 ppm
Iron (Fe)	<5 ppm	Total silicon (Si)	<20 ppm
Lead (Pb)	<5 ppm	Total sulphur (S)	<50 ppm
Magnesium (Mg)	<5 ppm	Zinc (Zn)	<10 ppm
Typical analysis			
Aluminium (Al)	<1 ppm	Tin (Sn)	<2 ppm
Arsenic (As)	<2 ppm	Manganese (Mn)	<1 ppm
Barium (Ba)	<2 ppm	Nickel (Ni)	<2 ppm
Cadmium (Cd)	<1 ppm	Titanium (Ti)	<1 ppm
Chromium (Cr)	<1 ppm	Vanadium (V)	<1 ppm
Cobalt (Co)	<1 ppm		

Lot analysis on the label.

Lithium chloride anhydrous
specified reagents for general laboratory work- extra pure

Packaging		Mfr. No	
100 g	Plastic jar	L/2200/48	
250 g	Plastic jar	L/2200/50	
5 kg	Plastic jar	L/2200/65	
25 kg	Fibre Keg	L/2200/70	
LiCl	EINECS: 231-212-3		
CAS: 7447-41-8	GHS: Acute Tox. 4, Skin Irrit.		
MW: 42.38	2, Eye Irrit. 2, STOT SE 3		
Density at 20°C: 2.06	H319, H315, H302, H312, H335		
Boiling point: 1382°C	P280, P302+P352,		
Melting point: 605°C	P301+P312, P305+P351+P338		
Assay: >98%			
Guaranteed analysis			
Calcium (Ca)	<100 ppm	Sodium (Na)	<200 ppm
Copper (Cu)	<20 ppm	Total phosphorus (P)	<100 ppm
Iron (Fe)	<20 ppm	Total silicon (Si)	<100 ppm
Lead (Pb)	<20 ppm	Total sulphur (S)	<200 ppm
Magnesium (Mg)	<20 ppm	Zinc (Zn)	<50 ppm
Potassium (K)	<200 ppm		

Lithium hydroxide monohydrate
for analysis

packaging		Mfr. No
500 g	Plastic jar	L/2301/53
LiOH.H ₂ O	UN: 2680	 
CAS: 1310-66-3	GHS: Acute Tox. 4, Skin	
MW: 41.96	Corr. 1A, Eye Dam. 1	
Density at 20°C: 1.51	H302, H314	
Boiling point: 920°C	P301+P330+P331,	
Melting point: 462°C	P305+P351+P338, P270,	
ADR/RID: 8, II	P280, P301+P312	
EINECS: 215-183-4		
Assay: >99.5%		
Guaranteed analysis		
Calcium (Ca)	<100 ppm	Sodium (Na)
Carbonate (Na ₂ CO ₃)	<0.5%	n-Butyl formate
Copper (Cu)	<5 ppm	Total nitrogen (N)
Iron (Fe)	<5 ppm	Total phosphorus (P)
Lead (Pb)	<5 ppm	Total silicon (Si)
Magnesium (Mg)	<5 ppm	Total sulphur (S)
Potassium (K)	<10 ppm	Zinc (Zn)

Lot analysis on the label.

Lithium hydroxide monohydrate
specified reagents for general laboratory work- extra pure

packaging		Mfr. No
100 g	Plastic jar	L/2300/48
250 g	Plastic jar	L/2300/50
1 kg	Plastic jar	L/2300/60
10 kg	Plastic keg	L/2300/68
LiOH.H ₂ O	UN: 2680	 
CAS: 1310-66-3	GHS: Acute Tox. 4, Skin	
MW: 41.96	Corr. 1A, Eye Dam. 1	
Density at 20°C: 1.51	H302, H314	
Boiling point: 920°C	P301+P330+P331,	
Melting point: 462°C	P305+P351+P338, P270,	
ADR/RID: 8, II	P280, P301+P312	
EINECS: 215-183-4		
Assay: >99%		
Guaranteed analysis		
Calcium (Ca)	<200 ppm	Sodium (Na)
Carbonate (Na ₂ CO ₃)	<1%	n-Butyl formate
Copper (Cu)	<20 ppm	Total phosphorus (P)
Iron (Fe)	<20 ppm	Total silicon (Si)
Lead (Pb)	<20 ppm	Total sulphur (S)
Magnesium (Mg)	<20 ppm	Zinc (Zn)
Potassium (K)	<50 ppm	<50 ppm

Lithium iodide trihydrate
specified reagents for general laboratory work- extra pure

Packaging		Mfr. No
100 g	Glass jar	L/2320/48
LiI.3H ₂ O	GHS: Acute Tox. 4, Skin Irrit.	
CAS: 7790-22-9	2, Eye Irrit. 2, STOT SE 3	
MW: 187.89	H335, H302, H315, H319	
Boiling point: 1171°C	P261, P301+P312,	
Melting point: 446°C	P302+P352, P280,	
EINECS: 233-822-5	P305+P351+P338	
Guaranteed analysis		
Identity test		: Pass

Lithium metaborate anhydrous
specified reagents for general laboratory work- extra pure

packaging		Mfr. No
100 g	Plastic jar	L/2330/48
LiBO ₂	Melting point: 845°C	
CAS: 13453-69-5	EINECS: 236-631-5	
MW: 49.75		

Guaranteed analysis			
Calcium (Ca)	<50 ppm	Lead (Pb)	<50 ppm
Iron (Fe)	<50 ppm	Magnesium (Mg)	<20 ppm

To determine the melting point.

Lithium nitrate anhydrous
specified reagents for general laboratory work- extra pure

Packaging		Mfr. No
100 g	Plastic jar	L/2350/48
500 g	Plastic jar	L/2350/53
5 kg	Plastic jar	L/2350/65
10 kg	Plastic keg	L/2350/68
LiNO ₃	EINECS: 232-218-9	 
CAS: 7790-69-4	UN: 2722	
MW: 68.95	GHS: Ox. Sol. 3, Acute Tox.	
Density at 20°C: 2.38	4, Eye Irrit. 2	
Boiling point: 600°C	H272, H302, H319	
Melting point: 251°C	P301+P312, P210,	
ADR/RID: 5.1, III	P305+P351+P338	
Assay: >98%		

Guaranteed analysis			
Calcium (Ca)	<100 ppm	Sodium (Na)	<200 ppm
Copper (Cu)	<20 ppm	n-Butyl formate	<0.005%
Iron (Fe)	<20 ppm	Total phosphorus (P)	<100 ppm
Lead (Pb)	<20 ppm	Total silicon (Si)	<100 ppm
Magnesium (Mg)	<20 ppm	Total sulphur (S)	<200 ppm
Potassium (K)	<20 ppm	Zinc (Zn)	<20 ppm

Lithium sulfate monohydrate
for analysis

Packaging		Mfr. No	
100 g	Plastic jar	L/2450/48	
Li ₂ SO ₄ .H ₂ O	GHS: Acute Tox. 4		
CAS: 10102-25-7	H302		
MW: 127.95	P261, P301+P312,		
Melting point: 860°C	P302+P352, P280,		
EINECS: 233-820-4	P305+P351+P338		
Assay: >99%			
Guaranteed analysis			
Calcium (Ca)	<20 ppm	Sodium (Na)	<50 ppm
Copper (Cu)	<5 ppm	n-Butyl formate	<0.001%
Iron (Fe)	<5 ppm	Total nitrogen (N)	<20 ppm
Lead (Pb)	<5 ppm	Total phosphorus (P)	<20 ppm
Magnesium (Mg)	<5 ppm	Total silicon (Si)	<20 ppm
Potassium (K)	<50 ppm	Zinc (Zn)	<10 ppm

Lot analysis on the label.

Lithium sulfate monohydrate
specified reagents for general laboratory work- extra pure

Packaging		Mfr. No
500 g	Plastic jar	L/2400/53
Li ₂ SO ₄ .H ₂ O	GHS: Acute Tox. 4	
CAS: 10102-25-7	H302	
MW: 127.95	P261, P301+P312,	
Melting point: 860°C	P302+P352, P280,	
EINECS: 233-820-4	P305+P351+P338	
Assay: >98%		
Guaranteed analysis		
Calcium (Ca)	<100 ppm	Sodium (Na)
Copper (Cu)	<20 ppm	n-Butyl formate
Iron (Fe)	<20 ppm	Total phosphorus (P)
Lead (Pb)	<20 ppm	Total silicon (Si)
Magnesium (Mg)	<20 ppm	Zinc (Zn)
Potassium (K)	<200 ppm	<50 ppm

Litmus
pure, indicator grade

packaging		Mfr. No
25 g	Glass jar	L/2500/46
CAS: 1393-92-6		
EINECS: 215-739-6		

Guaranteed analysis	
pH colour change test	: Pass

Litmus solution
pure, indicator grade

packaging		Mfr. No
500 mL	Glass bottle	L/2500L/08
2.5 L	Glass bottle	L/2500L/17
Guaranteed analysis		
pH colour change test	: Pass	

Lotoxane® pure			
packaging		Mfr. No	
1 ℓ	Glass bottle	L/2550/15	
2.5 ℓ	Glass bottle	L/2550/17	
25 ℓ	Metal drum	L/2550/25	
Density at 20°C: 0.77		GHS: Asp. Tox. 1	
Boiling point: 180°C		H304, EUH066	
Flash point: 62°C		P301+P310, P331, P405, P501	

Non chlorinated aromatic solvent to replace 1.1,1-trichloroethane in cold cleaning applications.

Luminol See 3-Aminophthalhydrazide (Luminol), page 33

Lutetium standard metal solution 1000 ppm
for ICP solution in HNO₃ 2%, NIST certified

packaging		Mfr. No	
100 mℓ	Plastic bottle	SPPLLU22Y	
500 mℓ	Plastic bottle	SPPLLU22X	
H315, H319			

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

L-Lysine monohydrochloride
specified reagents for general laboratory work- extra pure

packaging		Mfr. No	
100 g	Plastic jar	L/2700/48	
1 kg	Plastic jar	L/2700/60	
C ₆ H ₁₄ N ₂ O ₂ ·HCl		Melting point: 263 to 264°C	
CAS: 657-27-2		EINECS: 211-519-9	
MW: 182.65			
Assay: >98%			
Guaranteed analysis			
Loss on drying	<0.5%	Iron (Fe)	<10 ppm
Specific rotation	19° to 22°	Lead (Pb)	<10 ppm
Residue after ignition	<0.1%		

Magnesium standard metal solution 1000 ppm
for AAS, solution HNO₃ 1M

packaging		Mfr. No	
100 mℓ	HDPE Plastic bottle	J/8040/05	
500 mℓ	HDPE Plastic bottle	J/8040/08	
ADR/RID: 8, II		H314	
UN: 2031		P280, P305+P351+P338,	
GHS: Skin Corr. 1B		P301+P330+P331, P310	

Manufactured to a tolerance of ± 0.5%.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.

UHPLC

Filtered to 0.1 µm, Fisher Chemical UHPLC gradient grade solvents (Acetonitrile, Methanol and Water) are perfect for those working at high pressure.

See page 309 for more information



Magnesium standard metal solution 1000 ppm
for ICP, solution HNO₃ 1M

packaging		Mfr. No	
500 mℓ	HDPE Plastic bottle	J/8286/08	
ADR/RID: 8, II		H314	
UN: 2031		P280, P305+P351+P338,	
GHS: Skin Corr. 1B		P301+P330+P331, P310	

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.

Magnesium standard metal solution 10 000 ppm
for ICP, solution HNO₃ 1M

packaging		Mfr. No	
100 mℓ	HDPE Plastic bottle	J/8285/05	
ADR/RID: 8, II		H314	
UN: 2031		P280, P305+P351+P338,	
GHS: Skin Corr. 1B		P301+P330+P331, P310	

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.

Magnesium standard metal solution 1000 ppm
for ICP solution in HNO₃ 2%, NIST certified

packaging		Mfr. No	
100 mℓ	Plastic bottle	SPPLMG22Y	
250 mℓ	Plastic bottle	SPPLMG22T	
H315, H320			

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Magnesium standard metal solution 1000 ppm
for ICP solution in HCL 2%, NIST certified

packaging		Mfr. No	
500 mℓ	Plastic bottle	SPPLMG12X	

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Magnesium standard metal solution 10 000 ppm
for ICP solution in HNO₃ 5%, NIST certified

packaging		Mfr. No	
100 mℓ	Plastic bottle	SPPLMG23Y	
500 mℓ	Plastic bottle	SPPLMG23X	
H314			

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Magnesium standard metal solution 10 000 ppm
for ICP solution in HCL 5%, NIST certified

packaging		Mfr. No	
500 mℓ	Plastic bottle	SPPLMG13X	
H314			

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Magnesium standard metal solution 1000 ppm
for ICP-MS solution in HNO₃ 2%, NIST certified

packaging		Mfr. No	
100 mℓ	Plastic bottle	SPCLMG22Y	
H315, H320			

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Magnesium metal powder
specified reagents for general laboratory work- extra pure

packaging		Mfr. No	
250 g	Plastic jar	M/0050/50	
500 g	Plastic jar	M/0050/53	
Mg		ADR/RID: 4.3, 4.2, II	
CAS: 7439-95-4		EINECS: 231-104-6	
MW: 24.30		UN: 1418	
Density at 20°C: 1.74		GHS: Pyr. Sol. 1, Water-react. 1	
Boiling point: 1107°C		H260, H250	
Melting point: 651°C		P223, P210, P280	
Flash point: 500°C			
Guaranteed analysis			
Identity test		: Pass	

Magnesium metal turnings
pure, for Grignard reaction

packaging		Mfr. No	
250 g	Plastic jar	M/0200/50	
500 g	Plastic jar	M/0200/53	
Mg		EINECS: 231-104-6	
CAS: 7439-95-4		UN: 1869	
MW: 24.30		GHS: Flam. Sol. 1, Self-heat.	
Boiling point: 1107°C		1, Water-react. 2	
Melting point: 651°C		H228, H261, H252	
Flash point: 500°C		P210, P223	
ADR/RID: 4.1, III			
Guaranteed analysis			
Identity test		: Pass	

High Volume Solvent
Delivery Systems

High-volume solvent delivery systems are available in 10L to 1000L, in Europe only. Please email your contact details and a brief description of your requirements to bege1.sha.quotes@thermofisher.com



Magnesium acetate tetrahydrate
for analysis

packaging		Mfr. No	
500 g	Plastic jar	M/0350/53	
C ₄ H ₆ O ₄ ·Mg·4H ₂ O		Melting point: 80°C	
CAS: 16674-78-5		EINECS: 205-554-9	
MW: 214.46			

Assay: 99.5 to 100.5%			
Guaranteed analysis			
pH (aqueous solution 5%)	7 to 9	n-Butyl formate	<0.001%
Calcium (Ca)	<20 ppm	Total nitrogen (N)	<20 ppm
Copper (Cu)	<2 ppm	Total phosphorus (P)	<20 ppm
Iron (Fe)	<5 ppm	Total silicon (Si)	<20 ppm
Lead (Pb)	<2 ppm	Total sulphur (S)	<50 ppm
Potassium (K)	<10 ppm	Zinc (Zn)	<5 ppm
Sodium (Na)	<10 ppm		

Lot analysis on the label.

Magnesium acetate tetrahydrate
specified reagents for general laboratory work- extra pure

packaging		Mfr. No	
500 g	Plastic jar	M/0300/53	
C ₄ H ₆ O ₄ ·Mg·4H ₂ O		Density at 20°C: 1.454	
CAS: 16674-78-5		Melting point: 80°C	
MW: 214.46		EINECS: 205-554-9	

Assay: >98%			
Guaranteed analysis			
Calcium (Ca)	<100 ppm	n-Butyl formate	<0.02%
Copper (Cu)	<10 ppm	Total phosphorus (P)	<100 ppm
Iron (Fe)	<20 ppm	Total silicon (Si)	<100 ppm
Lead (Pb)	<10 ppm	Total sulphur (S)	<200 ppm
Potassium (K)	<50 ppm	Zinc (Zn)	<20 ppm
Sodium (Na)	<50 ppm		

Magnesium bromide hexahydrate
specified reagents for general laboratory work- extra pure

packaging		Mfr. No	
500 g	Plastic jar	M/0400/53	
MgBr ₂ ·6H ₂ O		GHS: Skin Irrit. 2, Eye Irrit. 2,	
CAS: 13446-53-2		STOT SE 3	
MW: 292.22		H335, H319, H315	
Density at 20°C: 20		P261, P302+P352, P280,	
Melting point: 165°C		P305+P351+P338	
EINECS: 232-170-9			
Assay: >98%			

Magnesium carbonate pentahydrate
for analysis, basic light

packaging		Mfr. No	
100 g	Plastic jar	M/0515/48	
MgCO ₃ ·5H ₂ O			
CAS: 12125-28-9			
MW: 62.31			
EINECS: 235-192-7			

Assay: 99 to 103%			
Guaranteed analysis			
Calcium (Ca)	<100 ppm	n-Butyl formate	<0.005%
Copper (Cu)	<5 ppm	Total nitrogen (N)	<10 ppm
Iron (Fe)	<10 ppm	Total phosphorus (P)	<10 ppm
Lead (Pb)	<10 ppm	Total silicon (Si)	<50 ppm
Potassium (K)	<20 ppm	Total sulphur (S)	<50 ppm
Sodium (Na)	<2000 ppm	Zinc (Zn)	<10 ppm

Lot analysis on the label.

Magnesium carbonate pentahydrate		SLR, Eur.Ph.,BP	
extra pure, basic light			
packaging		Mfr. No	
500 g	Plastic jar	M/0500/53	
25 kg	Plastic keg	M/0500/70	
MgCO ₃ ·5H ₂ O			
CAS: 12125-28-9			
MW: 62.31			
EINECS: 235-192-7			
Guaranteed analysis			
Residue after ignition	42 to 45%	Sodium (Na)	<5000 ppm
Calcium (Ca)	<10000 ppm	n-Butyl formate	<0.05%
Copper (Cu)	<10 ppm	Total phosphorus (P)	<100 ppm
Iron (Fe)	<500 ppm	Total sulphur (S)	<5000 ppm
Lead (Pb)	<10 ppm	Zinc (Zn)	<100 ppm
Potassium (K)	<200 ppm		

Magnesium carbonate pentahydrate		SLR, Eur.Ph.,BP	
extra pure, basic heavy			
packaging		Mfr. No	
500 g	Plastic jar	M/0450/53	
MgCO ₃ ·5H ₂ O			
CAS: 12125-28-9			
MW: 62.31			
EINECS: 235-192-7			
Guaranteed analysis			
Residue after ignition	42 to 45%	Sodium (Na)	<5000 ppm
Calcium (Ca)	<10000 ppm	n-Butyl formate	<0.05%
Copper (Cu)	<10 ppm	Total phosphorus (P)	<100 ppm
Iron (Fe)	<500 ppm	Total silicon (Si)	<500 ppm
Lead (Pb)	<10 ppm	Total sulphur (S)	<5000 ppm
Potassium (K)	<200 ppm	Zinc (Zn)	<100 ppm

Magnesium chloride hexahydrate			Certified AR, Eur.Ph.,BP
for analysis			
packaging			Mfr. No
500 g	Plastic jar		M/0600/53
1 kg	Plastic jar		M/0600/60
2 kg	Plastic jar		M/0600/61
5 kg	Plastic keg		M/0600/65
MgCl ₂ ·6H ₂ O		EINECS: 232-094-6	
CAS: 7791-18-6		GHS: Eye Irrit. 2, STOT SE 3	
MW: 203.31		H335, H319	
Density at 20°C: 1.56		P261, P280, P305+P351+P338	
Melting point: 117°C			
Assay: 99 to 101%			
Guaranteed analysis			
Acidity / Alkalinity	Pass	Potassium (K)	<10 ppm
Water	51 to 55%	Sodium (Na)	<10 ppm
pH (aqueous solution 5%)	5 to 7	Total nitrogen (N)	<20 ppm
Bromides	: <500 ppm	Total phosphorus (P)	<5 ppm
Heavy metals	<10 ppm	Total silicon (Si)	<5 ppm
Calcium (Ca)	<20 ppm	Total sulphur (S)	<20 ppm
Copper (Cu)	<2 ppm	Zinc (Zn)	<5 ppm
Iron (Fe)	<2 ppm	Arsenic (As)	<2 ppm
Lead (Pb)	<2 ppm	Sulfate (SO ₄)	<100 ppm
Typical analysis			
Aluminium (Al)	<1 ppm	Manganese (Mn)	<5 ppm
Silver (Ag)	<1 ppm	Molybdenum (Mo)	<1 ppm
Barium (Ba)	<0.001 %	Nickel (Ni)	<1 ppm
Cadmium (Cd)	<1 ppm	Strontium (Sr)	<0.002 %
Chromium (Cr)	<1 ppm	Titanium (Ti)	<1 ppm
Cobalt (Co)	<1 ppm	Vanadium (V)	<1 ppm
Lithium (Li)	<1 ppm		
Lot analysis on the label.			

Magnesium chloride hexahydrate			SLR
specified reagents for general laboratory work- extra pure			
Packaging			Mfr. No
500 g	Plastic jar		M/0550/53
1 kg	Plastic jar		M/0550/60
5 kg	Plastic keg		M/0550/65
MgCl ₂ ·6H ₂ O		EINECS: 232-094-6	
CAS: 7791-18-6		GHS: Eye Irrit. 2, STOT SE 3	
MW: 203.31		H335, H319	
Density at 20°C: 1.56		P261, P280, P305+P351+P338	
Melting point: 117°C			
Assay: >98%			
Guaranteed analysis			
Calcium (Ca)	<100 ppm	Sodium (Na)	<50 ppm
Copper (Cu)	<10 ppm	Total phosphorus (P)	<100 ppm
Iron (Fe)	<20 ppm	Total silicon (Si)	<100 ppm
Lead (Pb)	<10 ppm	Total sulphur (S)	<200 ppm
Potassium (K)	<50 ppm	Zinc (Zn)	<20 ppm

Magnesium chloride 50% w/w solution		Certified AR	
for analysis			
Packaging		Mfr. No	
2.5 l	Glass bottle	M/0602/17	
MgCl ₂ +aq			
CAS: 7786-30-3			
MW: 203.31			
EINECS: 232-094-6			
Assay: 49 to 51%			
Guaranteed analysis			
pH (aqueous solution 5%)	5 to 5.5	Sodium (Na)	<10 ppm
Calcium (Ca)	<20 ppm	Total nitrogen (N)	<20 ppm
Copper (Cu)	<2 ppm	Total phosphorus (P)	<5 ppm
Iron (Fe)	<2 ppm	Total silicon (Si)	<5 ppm
Lead (Pb)	<2 ppm	Total sulphur (S)	<20 ppm
Potassium (K)	<10 ppm	Zinc (Zn)	<5 ppm
Lot analysis on the label.			

Magnesium hydroxide		SLR
specified reagents for general laboratory work- extra pure		
Packaging		Mfr. No
50 g	Plastic jar	M/0650/50
100 g	Plastic jar	M/0650/53
Mg(OH) ₂	GHS: Skin Irrit. 2, Eye Irrit. 2,	
CAS: 1309-42-8	STOT SE 3	
MW: 58.33	H319, H335, H315	
Density at 20°C: 2.36	P261, P302+P352, P280,	
Melting point: 350°C	P305+P351+P338	
EINECS: 215-170-3		

Magnesium nitrate hexahydrate		Certified AR	
for analysis			
Packaging		Mfr. No	
500 g	Plastic jar	M/0752/53	
1 kg	Plastic jar	M/0752/60	
MgCl ₂ ·6H ₂ O	EINECS: 233-826-7		
CAS: 13446-18-9	UN: 1474		
MW: 256.41	GHS: Skin Irrit. 2, Eye Irrit. 2,		
Boiling point: 330°C	STOT SE 3		
Melting point: 89°C	H315, H319, H335		
Flash point: 93.3°C	P261, P302+P352, P280,		
ADR/RID: 5.1/III	P305+P351+P338		
Assay: 99 to 101%			
Guaranteed analysis			
pH (aqueous solution 5%)	5 to 7	Sodium (Na)	<10 ppm
Ammonium (NH ₄)	<0.002%	Total chloride (Cl)	<0.001%
Calcium (Ca)	<20 ppm	Total phosphorus (P)	<5 ppm
Copper (Cu)	<2 ppm	Total silicon (Si)	<5 ppm
Iron (Fe)	<2 ppm	Total sulphur (S)	<20 ppm
Lead (Pb)	<2 ppm	Zinc (Zn)	<5 ppm
Potassium (K)	<10 ppm		

Typical analysis			
Aluminium (Al)	<1 ppm	Manganese (Mn)	<5 ppm
Silver (Ag)	<1 ppm	Molybdenum (Mo)	<1 ppm
Barium (Ba)	<0.001 %	Nickel (Ni)	<1 ppm
Cadmium (Cd)	<1 ppm	Strontium (Sr)	<0.002 %
Chromium (Cr)	<1 ppm	Titanium (Ti)	<1 ppm
Cobalt (Co)	<1 ppm	Vanadium (V)	<1 ppm
Lithium (Li)	<1 ppm		
Lot analysis on the label.			

Magnesium nitrate hexahydrate		SLR	
specified reagents for general laboratory work- extra pure			
packaging		Mfr. No	
500 g	Plastic jar	M/0750/53	
Mg(NO ₃) ₂ ·6H ₂ O	EINECS: 233-826-7		
CAS: 13446-18-9	UN: 1474		
MW: 256.41	GHS: Skin Irrit. 2, Eye Irrit. 2,		
Boiling point: 330°C	STOT SE 3		
Melting point: 89°C	H315, H319, H335		
Flash point: 93.3°C	P261, P302+P352, P280,		
ADR/RID: 5.1/III	P305+P351+P338		
Assay: >98%			
Guaranteed analysis			
Calcium (Ca)	<500 ppm	n-Butyl formate	<0.01%
Copper (Cu)	<10 ppm	Total phosphorus (P)	<200 ppm
Iron (Fe)	<20 ppm	Total silicon (Si)	<100 ppm
Lead (Pb)	<10 ppm	Total sulphur (S)	<500 ppm
Potassium (K)	<50 ppm	Zinc (Zn)	<20 ppm
Sodium (Na)	<50 ppm		

Magnesium oxide		SLR, Eur.Ph.,BP
extra pure, light		
packaging		Mfr. No
500 g	Plastic jar	M/0850/53
1 kg	Plastic keg	M/0850/60
MgO		Boiling point: 3600°C
CAS: 1309-48-4		Melting point: 2800°C
MW: 40.29		EINECS: 215-171-9
Density at 20°C: 3.58		
Guaranteed analysis		
Loss on ignition	<5%	Lead (Pb) <20 ppm
Arsenic (As)	<5 ppm	n-Butyl formate <0.1%
Copper (Cu)	<50 ppm	Total sulphur (S) <10000 ppm
Iron (Fe)	<500 ppm	Zinc (Zn) <50 ppm

Magnesium oxide		SLR, Eur.Ph.,USP,BP	
extra pure, heavy			
packaging		Mfr. No	
250 g	Plastic jar	M/0800/50	
500 g	Plastic jar	M/0800/53	
MgO		Boiling point: 3600°C	
CAS: 1309-48-4		Melting point: 2800°C	
MW: 40.29		EINECS: 215-171-9	
Density at 20°C: 3.58			
Guaranteed analysis			
Loss on ignition	<5%	Lead (Pb)	<20 ppm
Arsenic (As)	<5 ppm	n-Butyl formate	<0.1%
Copper (Cu)	<50 ppm	Total sulphur (S)	<10000 ppm
Iron (Fe)	<500 ppm	Zinc (Zn)	<50 ppm

UHPLC

Filtered to 0.1 µm, Fisher Chemical UHPLC gradient grade solvents (Acetonitrile, Methanol and Water) are perfect for those working at high pressure.

See page 309 for more information

Magnesium perchlorate dried		
pure, powder, Desiccating agent		
packaging		Mfr. No
250 g	Plastic jar	M/0950/50
500 g	Plastic jar	M/0950/53
25 kg	Plastic keg	M/0950/70
Mg(ClO ₄) ₂		
CAS: 10034-81-8		
MW: 223.21		
Melting point: 251°C		
ADR/RID: 5.1, II		
EINECS: 233-108-3		
UN: 1475		
GHS: Ox. Sol. 2, Skin Irrit. 2, Eye Irrit. 2, STOT SE 3		
H272, H319, H335, H315, EUH014		
P210, P261, P402+P404, P302+P352, P280, P305+P351+P338		

Identity test	: Pass
---------------	--------

Magnesium stearate		SLR, Eur.Ph.,BP,USP
specified reagents for general laboratory work- extra pure		
packaging		Mfr. No
500 g	Plastic jar	M/0955/53
C ₃₆ H ₇₀ MgO ₄		
CAS: 557-04-0		
MW: 591.27		
Melting point: 130 to 140°C		
EINECS: 209-150-3		
GHS: Skin Irrit. 2, Eye Irrit. 2, STOT SE 3 H335, H319, H315 P261, P302+P352, P280, P305+P351+P338		
Guaranteed analysis		
Complies with standard IR spectrum: Pass		

Magnesium sulfate			SLR
extra pure, dried			
packaging		Mfr. No	
500 g	Plastic jar	M/1100/53	
1 kg	Plastic jar	M/1100/60	
2 kg	Plastic jar	M/1100/61	
10 kg	Plastic keg	M/1100/68	
MgSO ₄		Density at 20°C: 2.65	
CAS: 7487-88-9		Melting point: 1124°C	
MW: 120.37		EINECS: 231-298-2	
Assay: 62 to 70%			
Guaranteed analysis			
Iron (Fe)		<10 ppm	n-Butyl formate <0.05%
Lead (Pb)		<10 ppm	

Magnesium sulfate heptahydrate		Certified AR	
for analysis			
packaging		Mfr. No	
500 g	Plastic jar	M/1050/53	
1 kg	Plastic jar	M/1050/60	
MgSO ₄ ·7H ₂ O	Density at 20°C: 1.678		
CAS: 10034-99-8	EINECS: 231-229-8		
MW: 246.48			
Assay: 99.5 to 100.5%			
Guaranteed analysis			
pH (aqueous solution 5%)	5 to 7	Sodium (Na)	<20 ppm
Calcium (Ca)	<20 ppm	Total chloride (Cl)	<0.0005%
Copper (Cu)	<2 ppm	Total nitrogen (N)	<20 ppm
Iron (Fe)	<2 ppm	Total phosphorus (P)	<5 ppm
Lead (Pb)	<2 ppm	Total silicon (Si)	<5 ppm
Potassium (K)	<20 ppm	Zinc (Zn)	<5 ppm

Magnesium sulfate heptahydrate SLR, Eur.Ph.,USP,BP
specified reagents for general laboratory work- extra pure

packaging		Mfr. No	
1 kg	Plastic jar	M/1000/60	
2.5 kg	Plastic jar	M/1000/62	
5 kg	Plastic keg	M/1000/65	
MgSO ₄ ·7H ₂ O	Density at 20°C: 1.678		
CAS: 10034-99-8	EINECS: 231-229-8		
MW: 246.48			
Assay: >99%			
Guaranteed analysis			
Calcium (Ca)	<200 ppm	Sodium (Na)	<100 ppm
Copper (Cu)	<10 ppm	n-Butyl formate	<0.05%
Iron (Fe)	<20 ppm	Total phosphorus (P)	<100 ppm
Lead (Pb)	<10 ppm	Total silicon (Si)	<500 ppm
Potassium (K)	<500 ppm	Zinc (Zn)	<20 ppm

Malachite Green (oxalate) C.I. 42000
pure

packaging		Mfr. No
25 g	Glass jar	M/1245/46
C ₂ H ₂ N ₂ ·1/2C ₂ H ₂ O ₄ ·C ₂ HO ₄ CAS: 2437-29-8 MW: 927.03 Melting point: 164°C ADR/RID: 6.1, III EINECS: 219-441-7 UN: 2811 GHS: Repr. 2, Acute Tox. 3,		Eye Dam. 1, Aquatic Acute 1, Aquatic Chronic 1 H301, H318, H410, H361d P301+P310, P301+P330+P331, P280, P305+P351+P338, P310, P281, P273

Guaranteed analysis		
pH colour change test	: Pass	Complies with standard IR spectrum: Pass

Maleic acid Certified AR
for analysis

packaging		Mfr. No
500 g	Plastic jar	M/1255/53
C ₄ H ₄ O ₄ CAS: 110-16-7 MW: 116.08 Density at 20°C: 1.59 Melting point: 134 to 138°C Flash point: 127°C EINECS: 203-742-5		GHS: Acute Tox. 4, Skin Irrit. 2, Eye Irrit. 2, Skin Sens. 1, STOT SE 3 H317, H335, H319, H315, H302 P280, P301+P330+P331, P305+P351+P338, P302+P352, P264, P261
Assay: >99.5%		
Guaranteed analysis		
Calcium (Ca)	<5 ppm	n-Butyl formate
Copper (Cu)	<5 ppm	Total nitrogen (N)
Iron (Fe)	<5 ppm	Total phosphorus (P)
Lead (Pb)	<5 ppm	Total silicon (Si)
Magnesium (Mg)	<5 ppm	Total sulphur (S)
Potassium (K)	<20 ppm	Zinc (Zn)
Sodium (Na)	<20 ppm	Fumaric acid

Lot analysis on the label.

Standard Volumetric solution

For your convenience we offer a range of NIST traceable solutions. A full list of our products can be found on page 323. Can't see what you need? Please email your contact details and a brief description of your requirements to begele.sha.quotes@thermofisher.com

Maleic acid SLR
specified reagents for general laboratory work- extra pure

Packaging		Mfr. No
500 g	Plastic jar	M/1250/53
C₄H₂O₄ CAS: 110-16-7 MW: 116.08 Density at 20°C: 1.59 Melting point: 134 to 138°C Flash point: 127°C INECS: 203-742-5		GHS: Acute Tox. 4, Skin Irrit. 2, Eye Irrit. 2, Skin Sens. 1, STOT SE 3 H317, H335, H319, H315, H302 P280, P301+P330+P331, P305+P351+P338, P302+P352, P264, P261 
Assay: >99%		

Maleic anhydride SLR
specified reagents for general laboratory work- extra pure

packaging		Mfr. No
500 g	Plastic jar	M/1300/53
10 kg	Plastic keg	M/1300/68
C ₄ H ₂ O ₃ CAS: 108-31-6 MW: 98.06 Density at 20°C: 1.48 Boiling point: 200°C Melting point: 52 to 55°C Flash point: 102°C ADR/RID: 8, III EINECS: 203-571-6		UN: 2215 GHS: Acute Tox. 4, Skin Corr. 1B, Eye Dam. 1, Resp. Sens. 1, Skin Sens. 1 H302, H314, H317, H334 P280, P301+P330+P331, P305+P351+P338, P302+P352, P310, P342+P311

Assay: >99.5%

DL-Malic acid SLR
specified reagents for general laboratory work- extra pure

Packaging		Mfr. No
500 g	Plastic jar	M/1350/53
C₄H₆O₅ CAS: 6915-15-7 MW: 134.09 Melting point: 128 to 132°C Flash point: 203°C EINECS: 230-022-8		GHS: Skin Irrit. 2, Eye Dam. 1, STOT SE 3 H318, H315, H335 P280, P305+P351+P338, P310, P261, P302+P352 
Assay: >99%		
Guaranteed analysis		
Residue after ignition	<0.03%	n-Butyl formate
Sulfate (SO ₄)	<0.02%	<0.02%

Malonic acid SLR
specified reagents for general laboratory work- extra pure

packaging		Mfr. No
100 g	Plastic jar	M/1400/48
C ₃ H ₂ O ₄ CAS: 141-82-2 MW: 104.07 Density at 20°C: 1.63 Melting point: 130 to 135°C Flash point: 172°C EINECS: 205-503-0		GHS: Acute Tox. 4, Skin Irrit. 2, Eye Dam. 1, STOT SE 3 H335, H318, H315, H302 P280, P305+P351+P338, P310, P261, P301+P312, P302+P352
Assay: >99%		

Maltose monohydrate SLR
specified reagents for general laboratory work- extra pure

packaging		Mfr. No
100 g	Plastic jar	M/1450/48
500 g	Plastic jar	M/1450/53
C ₁₂ H ₂₂ O ₁₁ ·H ₂ O CAS: 69-79-4 MW: 360.32		Melting point: 130°C EINECS: 200-716-5
Guaranteed analysis		
Complies with standard IR spectrum: Pass		

Manganese standard metal solution 1000 ppm
for AAS, solution HNO₃ 1M

packaging		Mfr. No
100 mL	HDPE Plastic bottle	J/8045/05
500 mL	HDPE Plastic bottle	J/8045/08
ADR/RID: 8, II UN: 2031 GHS: Skin Corr. 1B, Eye Dam. 1		H314 P280, P305+P351+P338, P301+P330+P331, P310
Manufactured to a tolerance of ± 0.5%. All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities. See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.		

Manganese standard metal solution 1000 ppm
for ICP, solution HNO₃ 1M

packaging		Mfr. No
500 mL	HDPE Plastic bottle	J/8291/08
ADR/RID: 8, II UN: 2031 GHS: Skin Corr. 1B, Eye Dam. 1		H314 P280, P305+P351+P338, P301+P330+P331, P310
Manufactured to a tolerance of ± 0.5%. All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities. See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.		

Manganese standard metal solution 10 000 ppm
for ICP, solution HNO₃ 1M

packaging		Mfr. No
100 mL	HDPE Plastic bottle	J/8290/05
Density at 20°C: 1.0 ADR/RID: 8, II EINECS: 231-104-6 UN: 2031		GHS: Skin Corr. 1B, Eye Dam. 1 H314 P280, P305+P351+P338, P301+P330+P331, P310
Manufactured to a tolerance of ± 0.5%. All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities. See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.		

Manganese standard metal solution 1000 ppm
for ICP solution in HNO₃ 2%, NIST certified

packaging		Mfr. No
100 mL	Plastic bottle	SPPLMN22Y
250 mL	Plastic bottle	SPPLMN22T
500 mL	Plastic bottle	SPPLMN22X
H315, H319		

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Manganese standard metal solution 10 000 ppm
for ICP solution in HNO₃ 5%, NIST certified

packaging		Mfr. No
100 mL	Plastic bottle	SPPLMN23Y
500 mL	Plastic bottle	SPPLMN23X
H314		

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Manganese standard metal solution 1000 ppm
for ICP-MS solution in HNO₃ 2%, NIST certified

packaging		Mfr. No
100 mL	Plastic bottle	SPCLMN22Y
H315, H320		

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Manganese(II) acetate tetrahydrate SLR
specified reagents for general laboratory work- extra pure

Packaging		Mfr. No
250 g	Plastic jar	M/1700/50
C ₄ H ₆ O ₄ Mn.4H ₂ O CAS: 6156-78-1 MW: 245.09 Melting point: >300°C Flash point: >130°C EINECS: 211-334-3		GHS: Repr. 2, Skin Irrit. 2, Eye Irrit. 2, STOT SE 3 H315, H319, H335, H361f P281, P261, P302+P352, P280, P305+P351+P338
Assay: >99%		 
Guaranteed analysis		
Iron (Fe)	<100 ppm	n-Butyl formate
Lead (Pb)	<100 ppm	Total sulphur (S)
		<0.005% <100 ppm

Label

Fisher Chemical labels show GHS pictograms clearly. Together with actual lot analysis shown on the label, plus a wealth of additional information, we help you to work safely to achieve accurate and consistent results in your analytical applications.

Manganese(II) chloride tetrahydrate		Certified AR
for analysis		
packaging	Mfr. No	
500 g Plastic jar	M/1850/53	
1 kg Plastic jar	M/1850/60	
25 kg Fibre Keg	M/1850/70	
MnCl ₂ ·4H ₂ O	GHS: Acute Tox. 4, STOT RE 2, Aquatic Chronic 3	
CAS: 13446-34-9	H302, H412, H373	
MW: 197.91	P280, P260, P301+P312, P273	
Melting point: 58°C		
EINECS: 231-869-6		

Assay: 98 to 102%	
Guaranteed analysis	
pH (aqueous solution 5%)	5 to 6
Substances	<20 ppm
reducing KMnO ₄	<0.0005%
Calcium (Ca)	<50 ppm
Copper (Cu)	<2 ppm
Iron (Fe)	<5 ppm
Lead (Pb)	<2 ppm
Magnesium (Mg)	<20 ppm
Potassium (K)	<20 ppm
Sodium (Na)	<50 ppm
Total nitrogen (N)	<30 ppm
Total phosphorus (P)	<50 ppm
Total silicon (Si)	<100 ppm
Total sulphur (S)	<20 ppm
Zinc (Zn)	<10 ppm
Typical analysis	
Arsenic (As)	<5 ppm
Cadmium (Cd)	<2 ppm
Chromium (Cr)	<2 ppm
Cobalt (Co)	<2 ppm
Lithium (Li)	<2 ppm
Titanium (Ti)	<2 ppm
Vanadium (V)	<2 ppm

Lot analysis on the label.

Manganese(II) chloride tetrahydrate		SLR
specified reagents for general laboratory work- extra pure		
packaging	Mfr. No	
250 g Plastic jar	M/1800/50	
500 g Plastic jar	M/1800/53	
MnCl ₂ ·4H ₂ O	GHS: Acute Tox. 4, STOT RE 2, Aquatic Chronic 3	
CAS: 13446-34-9	H302, H412, H373	
MW: 197.91	P280, P260, P301+P312, P273	
Melting point: 58°C		
EINECS: 231-869-6		

Assay: >98%	
Guaranteed analysis	
Calcium (Ca)	<200 ppm
Copper (Cu)	<20 ppm
Iron (Fe)	<20 ppm
Lead (Pb)	<20 ppm
Magnesium (Mg)	<50 ppm
Potassium (K)	<50 ppm
Sodium (Na)	<200 ppm
Total phosphorus (P)	<200 ppm
Total silicon (Si)	<200 ppm
Total sulphur (S)	<500 ppm
Zinc (Zn)	<200 ppm

Manganese (II) nitrate tetrahydrate

NEW

Certified AR

for analysis

packaging	Mfr. No	
500 g Plastic jar	M/2052/53	
Mn(NO ₃) ₂ ·4H ₂ O	GHS: Ox. Sol. 2, Skin Irrit. 2, Eye Irrit. 2, STOT SE 3	
CAS: 10377-66-9	H315, H319, H335, H272	
ADR/RID: 5.1, III	P210, P261, P280, P305+P351+P338	
EINECS: 231-104-6		
UN: 2724		

Assay: >98%	
Guaranteed analysis	
Calcium (Ca)	<50 ppm
Copper (Cu)	<10 ppm
Iron (Fe)	<10 ppm
Lead (Pb)	<10 ppm
Magnesium (Mg)	<20 ppm
Potassium (K)	<20 ppm
Sodium (Na)	<500 ppm
n-Butyl formate	<0.002%
Total phosphorus (P)	<100 ppm
Total sulphur (S)	<50 ppm
Zinc (Zn)	<50 ppm

Typical analysis	
Cadmium (Cd)	<1 ppm
Chromium (Cr)	<1 ppm
Cobalt (Co)	<1 ppm
Lithium (Li)	<1 ppm
Titanium (Ti)	<1 ppm
Vanadium (V)	<1 ppm

Lot analysis on the label

Manganese(IV) oxide

extra pure, precipitated

SLR

packaging	Mfr. No	
250 g Glass jar	M/1665/50	
500 g Glass jar	M/1665/53	
MnO ₂	EINECS: 215-202-6	
CAS: 1313-13-9	UN: 1479	
MW: 86.94	GHS: Acute Tox. 4	
Boiling point: DC°C	H332, H302	
Melting point: 535°C	P210, P270, P264, P301+P312, P304+P340	
ADR/RID: 5.1, II		

Guaranteed analysis	
Identity test	: Pass

Manganese(IV) oxide

Technical, powder

packaging	Mfr. No	
500 g Glass jar	M/1660/53	
MnO ₂	EINECS: 215-202-6	
CAS: 1313-13-9	UN: 1479	
MW: 86.94	GHS: Acute Tox. 4	
Boiling point: DC°C	H332, H302	
Melting point: 535°C	P210, P270, P264, P301+P312, P304+P340	
ADR/RID: 5.1, II		

Guaranteed analysis	
Identity test	: Pass

Manganese(II) sulfate tetrahydrate

for analysis

Certified AR

packaging	Mfr. No	
250 g Plastic jar	M/2300/50	
500 g Plastic jar	M/2300/53	
MnSO ₄ ·4H ₂ O	GHS: STOT RE 2, Aquatic Chronic 2	
CAS: 10101-68-5	H373, H411	
MW: 223.05	P260, P273	
ADR/RID: 9, III		
UN: 3077		

Assay: >99%	
Guaranteed analysis	
pH (aqueous solution 5%)	5 to 6
Substances	<20 ppm
reducing KMnO ₄	<0.0005%
Calcium (Ca)	<50 ppm
Copper (Cu)	<2 ppm
Iron (Fe)	<5 ppm
Lead (Pb)	<2 ppm
Magnesium (Mg)	<20 ppm
Potassium (K)	<20 ppm
Sodium (Na)	<50 ppm
n-Butyl formate	<0.001%
Total nitrogen (N)	<10 ppm
Total phosphorus (P)	<50 ppm
Total silicon (Si)	<100 ppm
Zinc (Zn)	<10 ppm

Typical analysis	
Arsenic (As)	<5 ppm
Cadmium (Cd)	<2 ppm
Chromium (Cr)	<2 ppm
Cobalt (Co)	<2 ppm
Lithium (Li)	<2 ppm
Titanium (Ti)	<2 ppm
Vanadium (V)	<2 ppm

Lot analysis on the label

pH Buffers

NIST Certified Standard Solution & Concentrate

A full range to discover on page 63



Manganese(II) sulfate tetrahydrate

specified reagents for general laboratory work- extra pure

SLR

packaging	Mfr. No	
500 g Plastic jar	M/2250/53	
1 kg Plastic jar	M/2250/60	
MnSO ₄ ·4H ₂ O	EINECS: 232-089-9	
CAS: 10101-68-5	UN: 3077	
MW: 223.05	GHS: STOT RE 2, Aquatic Chronic 2	
Boiling point: 850°C	H373, H411	
Melting point: -700°C	P260, P273	
ADR/RID: 9, III		

Assay: >98%	
Guaranteed analysis	
Calcium (Ca)	<200 ppm
Copper (Cu)	<20 ppm
Iron (Fe)	<20 ppm
Lead (Pb)	<20 ppm
Magnesium (Mg)	<50 ppm
Potassium (K)	<50 ppm
Sodium (Na)	<200 ppm
n-Butyl formate	<0.05%
Total phosphorus (P)	<200 ppm
Total silicon (Si)	<200 ppm
Zinc (Zn)	<100 ppm

Manganese(II) sulfate monohydrate

pure, dried powder

packaging	Mfr. No	
1 kg Plastic jar	M/2350/60	
MnSO ₄ ·H ₂ O	EINECS: 232-089-9	
CAS: 10034-96-5	UN: 3077	
MW: 169.02	GHS: STOT RE 2, Aquatic Chronic 2	
Boiling point: 850°C	H373, H411	
Melting point: 700°C	P260, P273	
ADR/RID: 9, III		

Guaranteed analysis	
Identity test	: Pass

Mannitol

for analysis

Certified AR, Eur.Ph., BP

packaging	Mfr. No	
500 g Plastic jar	M/2405/53	
C ₆ H ₁₄ O ₆	Melting point: 165 to 170°C	
CAS: 69-65-8	EINECS: 200-711-8	
MW: 182.17		

Assay (CG): 99 to 102%	
Guaranteed analysis	
Water	<0.5%
Specific rotation	23° to 24°
Reducing sugars	<0.05%
Conductivity	<20 µS
Isomalt	<2%
D-Maltitol	<2%
D-Sorbitol	<2%
Acidity	<0.0008 meq/g
Calcium (Ca)	<10 ppm
Copper (Cu)	<5 ppm
Iron (Fe)	<5 ppm
Lead (Pb)	<0.5 ppm
Magnesium (Mg)	<5 ppm
Potassium (K)	<10 ppm
Sodium (Na)	<50 ppm
Total chloride (Cl)	<0.001%
Total phosphorus (P)	<20 ppm
Total silicon (Si)	<20 ppm
Total sulphur (S)	<50 ppm
Zinc (Zn)	<10 ppm

Typical analysis	
Aluminium (Al)	<1 ppm
Arsenic (As)	<5 ppm
Barium (Ba)	<1 ppm
Cadmium (Cd)	<1 ppm
Chromium (Cr)	<1 ppm
Cobalt (Co)	<1 ppm
Tin (Sn)	<2 ppm
Lithium (Li)	<1 ppm
Manganese (Mn)	<1 ppm
Nickel (Ni)	<1 ppm
Titanium (Ti)	<1 ppm
Vanadium (V)	<1 ppm

Lot analysis on the label.

Save money!

Order 4 bottles in a box



Mannitol

specified reagents for general laboratory work- extra pure

SLR

packaging	Mfr. No	
500 g Plastic jar	M/2400/53	
C ₆ H ₁₄ O ₆	Melting point: 165 to 170°C	
CAS: 69-65-8	EINECS: 200-711-8	
MW: 182.17		

Guaranteed analysis	
Water	<0.5%
Reducing sugars	: <0.2%
Calcium (Ca)	<100 ppm
Copper (Cu)	<50 ppm
Iron (Fe)	<50 ppm
Lead (Pb)	<50 ppm
Magnesium (Mg)	<50 ppm
Potassium (K)	<50 ppm
Sodium (Na)	<200 ppm
n-Butyl formate	<0.005%
Total phosphorus (P)	<100 ppm
Total silicon (Si)	<100 ppm
Total sulphur (S)	<200 ppm
Zinc (Zn)	<10 ppm

D(+)-Mannose

specified reagents for general laboratory work- extra pure

SLR

packaging	Mfr. No	
10 g Plastic jar	M/2420/45	
C ₆ H ₁₂ O ₆	Melting point: 133°C	
CAS: 3458-28-4	EINECS: 222-392-4	
MW: 180.16		

Guaranteed analysis	
Loss on drying	<0.5%
Specific rotation	13.3° to 14.3°
Residue after ignition	<0.1%
Lead (Pb)	<10 ppm

Marble chips

pure

packaging	Mfr. No	
1 kg Plastic jar	M/2450/60	
3 kg Plastic jar	M/2450/63	
CCaO ₃	GHS: Skin Irrit. 2, Eye Irrit. 2	
CAS: 471-34-1	H315, H319	
MW: 100.09	P302+P352, P280, P305+P351+P338	
Melting point: 825°C		
EINECS: 207-439-9		

Guaranteed analysis	
Identity test	: Pass

Marble chips

technical, granule size 2-3 mm

packaging	Mfr. No	
1 kg Plastic jar	M/2451/60	
3 kg Plastic jar	M/2451/63	
CCaO ₃	EINECS: 207-439-9	
CAS: 471-34-1	GHS: Skin Irrit. 2, Eye Irrit. 2	
MW: 100.09	H319, H315	
Density at 20°C: 2.7 to 2.9	P261, P280, P302+P352, P305+P351+P338	
Melting point: 825°C		

Guaranteed analysis	
Identity test	: Pass

MEK See Ethyl methyl ketone, page 128

D(+)-Melibiose monohydrate

pure

packaging	Mfr. No	
10 g Plastic jar	M/2510/45	
C ₁₂ H ₂₂ O ₁₁ ·H ₂ O	Melting point: 182°C	
CAS: 66009-10-7	EINECS: 209-568-6	
MW: 360.32		

Guaranteed analysis	
Loss on drying	<6%
Specific rotation	133° to 137°
Residue after ignition	<0.2%
Lead (Pb)	<10 ppm

Menthol		SLR
extra pure, natural		
packaging	Mfr. No	
100 g Glass jar	M/2550/48	
C ₁₀ H ₂₀ O CAS: 2216-51-5 MW: 156.28 Density at 20°C: 0.890 Boiling point: 212°C Melting point: 41 to 43°C Flash point: 93°C	EINECS: 218-690-9 GHS: Skin Irrit. 2, Eye Dam. 1, STOT SE 3 H315, H318, H335 P280, P305+P351+P338, P310, P261, P302+P352	
Guaranteed analysis		
Specific rotation -50° to -49°		

Mercaptoethanol		
technical		
packaging	Mfr. No	
100 mℓ Glass bottle	M/P200/05	
C ₂ H ₆ OS CAS: 60-24-2 MW: 78.13 Density at 20°C: 1.11 Boiling point: 157 to 158°C Melting point: -100°C Flash point: 73°C ADR/RID: 6.1, II EINECS: 200-464-6	UN: 2966 GHS: Acute Tox. 3, Skin Irrit. 2, Eye Irrit. 2, STOT SE 3, Aquatic Chronic 2 H302, H411, H315, H335, H311, H319 P280, P312, P302+P350, P261, P273, P301+P312, P305+P351+P338	 
Assay (CG): >98%		
Guaranteed analysis		
Colour <30 APHA		

Mercury standard metal solution 1000 ppm		
for AAS, solution HNO ₃ 1M		
packaging	Mfr. No	
100 mℓ HDPE Plastic bottle	J/8047/05	
500 mℓ HDPE Plastic bottle	J/8047/08	
Density at 20°C: 1.0 ADR/RID: 6.1, 8, II UN: 3289 GHS: Skin Corr. 1B, Eye	Dam. 1, STOT RE 2 H373, H314 P280, P305+P351+P338, P301+P330+P331, P310, P260	 

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.

Mercury standard metal solution 1000 ppm		
for ICP, solution HNO ₃ 1M		
packaging	Mfr. No	
500 mℓ HDPE Plastic bottle	J/8293/08	
Density at 20°C: 1.0 ADR/RID: 6.1, 8, II UN: 3289 GHS: Skin Corr. 1B, Eye	Dam. 1, STOT RE 2 H373, H314 P280, P305+P351+P338, P301+P330+P331, P310, P260	 

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.

Mercury standard metal solution 10 000 ppm		
for ICP, solution HNO ₃ 1M		
packaging	Mfr. No	
100 mℓ HDPE Plastic bottle	J/8292/05	
ADR/RID: 6.1, 8, II UN: 3289 GHS: Repr. 1B, Acute Tox. 3, Skin Corr. 1B, Eye Dam. 1, STOT RE 2 H314, H331, H372, H360D	P202, P280, P260, P304+P340, P301+P330+P331, P305+P351+P338, P302+P352, P310	  

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.

Mercury standard metal solution 1000 ppm		
for ICP solution in HNO ₃ 10%, NIST certified		
packaging	Mfr. No	
100 mℓ Plastic bottle	SPPLHG42Y	
250 mℓ Plastic bottle	SPPLHG42T	
500 mℓ Plastic bottle	SPPLHG42X	
H314		

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Mercury standard metal solution 10 ppm		
for ICP solution in HNO ₃ 5%, NIST certified		
packaging	Mfr. No	
100 mℓ Plastic bottle	SPPLHG21AY	
500 mℓ Plastic bottle	SPPLHG21AX	
500 mℓ Plastic bottle	SPPLHG21X	
H314		

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314
Check with your local distributor for availability of Spex product

Mercury standard metal solution 10 000 ppm		
for ICP solution in HNO ₃ 10%, NIST certified		
packaging	Mfr. No	
100 mℓ Plastic bottle	SPPLHG43Y	
500 mℓ Plastic bottle	SPPLHG43X	
H314, H412		

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Mercury standard metal solution 10 ppm		
for ICP-MS solution in HNO ₃ 5%, NIST certified		
packaging	Mfr. No	
100 mℓ Plastic bottle	SPCLHG21AY	
H314		

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Mercury (metal)		Certified AR
for analysis,		
packaging	Mfr. No	
100 g Plastic jar	M/3750/48	
250 g Plastic jar	M/3750/50	
500 g Plastic jar	M/3750/53	
1 kg HDPE Plastic bottle	M/3750/60	
Hg CAS: 7439-97-6 MW: 200.59 Density at 20°C: 13.54 Boiling point: 356.5°C Melting point: -38.9°C ADR/RID: 8, III EINECS: 231-106-7	UN: 2809 GHS: Repr. 1B, Acute Tox. 2, STOT RE 1, Aquatic Acute 1, Aquatic Chronic 1 H410, H330, H360D, H372 P280, P308+P313, P304+P340, P260, P273	  
Assay: >99.99%		
Guaranteed analysis		
Cadmium (Cd) <0.2 ppm Magnesium (Mg) <0.5 ppm		
Calcium (Ca) <1 ppm Nickel (Ni) <0.2 ppm		
Copper (Cu) <0.5 ppm Potassium (K) <1 ppm		
Iron (Fe) <0.5 ppm Sodium (Na) <2 ppm		
Lead (Pb) <0.5 ppm Total sulphur (S) <1 ppm		
Lithium (Li) <0.2 ppm Zinc (Zn) <0.5 ppm		

The prices applied to mercury and derivatives vary with the price of the raw material.
Lot analysis on the label.

Mercury (metal)		SLR
specified reagents for general laboratory work- extra pure		
packaging	Mfr. No	
250 g HDPE Plastic bottle	M/3650/50	
1 kg Plastic jar	M/3650/60	
Hg CAS: 7439-97-6 MW: 200.59 Density at 20°C: 13.54 Boiling point: 356.5°C Melting point: -38.9°C ADR/RID: 8, III EINECS: 231-106-7	UN: 2809 GHS: Repr. 1B, Acute Tox. 2, STOT RE 1, Aquatic Acute 1, Aquatic Chronic 1 H410, H330, H360D, H372 P280, P308+P313, P304+P340, P260, P273	  

Assay: > 99%
Guaranteed analysis
Cadmium (Cd) <0.2 ppm Lead (Pb) <0.5 ppm
The price charged for this item is subject to alteration and will be that ruling at time of despatch..

Mercury (metal)		
pure, for polarography		
packaging	Mfr. No	
100 g Plastic jar	M/3751/48	
Hg CAS: 7439-97-6 MW: 200.59 Density at 20°C: 13.54 Boiling point: 356.5°C Melting point: -38.9°C ADR/RID: 8, III EINECS: 231-106-7	UN: 2809 GHS: Repr. 1B, Acute Tox. 2, STOT RE 1, Aquatic Acute 1, Aquatic Chronic 1 H410, H330, H360D, H372 P280, P308+P313, P304+P340, P260, P273	  

Guaranteed analysis
Identity test : Pass
The price charged for this item is subject to alteration and will be that ruling at time of despatch.

Mercury clean up kit		
pure		
packaging	Mfr. No	
Pack Box	M/3754/90	
ADR/RID: 8, III UN: 3260 GHS: Acute Tox. 4, Skin Corr. 1B, Eye Dam. 1, Aquatic Acute 1, Aquatic Chronic 1	H314, H302, H410 P301+P330+P331, P280, P305+P351+P338, P310, P273, P301+P312	  

The Mercury Clean-up kit is a compact safety pack used for cleaning up split mercury. It contains the special mercury absorber, safety gloves, brush, and scoop, and a spare container for holding the collected mercury prior to disposal. Full instructions are only included with each kit.

Mercury absorber spillage treatment		
technical		
packaging	Mfr. No	
100 g Plastic jar	M/3755/48	
ADR/RID: 8, III UN: 3260 GHS: Acute Tox. 4, Skin Corr. 1B, Eye Dam. 1, Aquatic Acute 1, Aquatic Chronic 1	H314, H302, H410 P301+P330+P331, P280, P305+P351+P338, P310, P273, P301+P312	  



Mercury(II) acetate
for analysis

Certified AR

packaging		Mfr. No
50 g	Plastic jar	M/2652/47
250 g	Plastic jar	M/2652/50
C ₄ H ₆ HgO ₄ CAS: 1600-27-7 MW: 318.68 Density at 20°C: 32.5 Melting point: 179 to 182°C ADR/RID: 6.1, II EINECS: 216-491-1 UN: 1629		GHS: Acute Tox. 1, STOT RE 2, Aquatic Acute 1, Aquatic Chronic 1 H330, H300, H310, H373, H410 P301+P310, P280, P361, P302+P350, P310, P304+P340, P260, P273



Assay: >99%			
Guaranteed analysis			
Calcium (Ca)	<20 ppm	Total chloride (Cl)	<0.005%
Copper (Cu)	<5 ppm	Total nitrogen (N)	<50 ppm
Iron (Fe)	<10 ppm	Total phosphorus (P)	<20 ppm
Lead (Pb)	<10 ppm	Total silicon (Si)	<20 ppm
Magnesium (Mg)	<10 ppm	Total sulphur (S)	<50 ppm
Potassium (K)	<10 ppm	Zinc (Zn)	<5 ppm
Sodium (Na)	<20 ppm		
Typical analysis			
Aluminium (Al)	<2 ppm	Lithium (Li)	<1 ppm
Barium (Ba)	<1 ppm	Manganese (Mn)	<1 ppm
Cadmium (Cd)	<1 ppm	Nickel (Ni)	<2 ppm
Chromium (Cr)	<1 ppm	Titanium (Ti)	<1 ppm
Cobalt (Co)	<1 ppm	Vanadium (V)	<1 ppm

Lot analysis on the label.

Mercury(II) acetate
specified reagents for general laboratory work- extra pure

SLR

packaging		Mfr. No
250 g	Plastic jar	M/2650/50
C ₄ H ₆ HgO ₄ CAS: 1600-27-7 MW: 318.68 Density at 20°C: 32.5 Melting point: 179 to 182°C ADR/RID: 6.1, II EINECS: 216-491-1 UN: 1629		GHS: Acute Tox. 1, STOT RE 2, Aquatic Acute 1, Aquatic Chronic 1 H330, H300, H310, H373, H410 P301+P310, P280, P361, P302+P350, P310, P304+P340, P260, P273



Assay: >98%			
Guaranteed analysis			
Calcium (Ca)	<100 ppm	Sodium (Na)	<100 ppm
Copper (Cu)	<50 ppm	n-Butyl formate	<0.005%
Iron (Fe)	<50 ppm	Total phosphorus (P)	<50 ppm
Lead (Pb)	<50 ppm	Total silicon (Si)	<50 ppm
Magnesium (Mg)	<50 ppm	Total sulphur (S)	<100 ppm
Potassium (K)	<50 ppm	Zinc (Zn)	<50 ppm

Mercury(I) chloride
for analysis

Certified AR

packaging		Mfr. No
100 g	Plastic jar	M/3450/48
Hg ₂ Cl ₂ CAS: 10112-91-1 MW: 472.09 Density at 20°C: 7.15 Melting point: 400°C ADR/RID: 6.1, III EINECS: 233-307-5 UN: 2025		GHS: Acute Tox. 4, Skin Irrit. 2, Eye Irrit. 2, STOT SE 3, Aquatic Acute 1, Aquatic Chronic 1 H302, H319, H335, H315, H410 P280, P304+P340, P260, P273, P305+P351+P338, P302+P352



Assay: >99.5%

Guaranteed analysis			
Calcium (Ca)	<20 ppm	Sodium (Na)	<20 ppm
Copper (Cu)	<5 ppm	Total nitrogen (N)	<20 ppm
Iron (Fe)	<5 ppm	Total phosphorus (P)	<20 ppm
Lead (Pb)	<5 ppm	Total silicon (Si)	<20 ppm
Magnesium (Mg)	<5 ppm	Total sulphur (S)	<100 ppm
Potassium (K)	<10 ppm	Zinc (Zn)	<10 ppm

Lot analysis on the label.

Mercury (II) chloride
for analysis

Certified AR, BP & USP

packaging		Mfr. No
100 g	Plastic jar	M/2850/48
250 g	Plastic jar	M/2850/50
500 g	Plastic jar	M/2850/53
Cl ₂ Hg CAS: 7487-94-7 MW: 271.50 Density at 25°C: 5.44 Boiling point: 302°C Melting point: 277°C ADR/RID: 6.1, II EINECS: 231-299-8 UN: 1624 GHS: Muta. 2, Repr. 2,		Acute Tox. 2, Skin Corr. 1B, STOT RE 1, Aquatic Acute 1, Aquatic Chronic 1 H300, H314, H410, H372, H341, H361f P301+P310, P280, P361, P302+P350, P301+P330+P331, P305+P351+P338, P310, P273



Assay: >99.5%			
Guaranteed analysis			
Acidity / Alkalinity	Pass	Potassium (K)	<10 ppm
Loss on drying	<1%	Sodium (Na)	<20 ppm
Mercurous chloride test	: Pass	Total nitrogen (N)	<10 ppm
Calcium (Ca)	<20 ppm	Total phosphorus (P)	<20 ppm
Copper (Cu)	<5 ppm	Total silicon (Si)	<20 ppm
Iron (Fe)	<10 ppm	Total sulphur (S)	<50 ppm
Lead (Pb)	<5 ppm	Zinc (Zn)	<10 ppm
Magnesium (Mg)	<5 ppm		

Typical analysis			
Aluminium (Al)	<2 ppm	Lithium (Li)	<1 ppm
Barium (Ba)	<1 ppm	Manganese (Mn)	<1 ppm
Cadmium (Cd)	<1 ppm	Nickel (Ni)	<1 ppm
Chromium (Cr)	<1 ppm	Titanium (Ti)	<1 ppm
Cobalt (Co)	<1 ppm	Vanadium (V)	<1 ppm

Lot analysis on the label

Mercury(II) chloride
specified reagents for general laboratory work- extra pure

SLR

packaging		Mfr. No
100 g	Plastic jar	M/2800/48
250 g	Plastic jar	M/2800/50
500 g	Plastic jar	M/2800/53
HgCl ₂ CAS: 7487-94-7 MW: 271.50 Density at 25°C: 5.44 Boiling point: 302°C Melting point: 277°C ADR/RID: 6.1, II EINECS: 231-299-8 UN: 1624 GHS: Muta. 2, Repr. 2, Acute Tox. 2, Skin Corr. 1B,		Eye Dam. 1, STOT RE 1, Aquatic Acute 1, Aquatic Chronic 1 H300, H314, H341, H410, H361f, H372 P301+P310, P280, P361, P302+P350, P301+P330+P331, P305+P351+P338, P310, P281, P304+P340, P260, P273



Assay: >98%			
Guaranteed analysis			
Calcium (Ca)	<100 ppm	Sodium (Na)	<100 ppm
Copper (Cu)	<50 ppm	Total phosphorus (P)	<50 ppm
Iron (Fe)	<50 ppm	Total silicon (Si)	<50 ppm
Lead (Pb)	<50 ppm	Total sulphur (S)	<100 ppm
Magnesium (Mg)	<50 ppm	Zinc (Zn)	<50 ppm
Potassium (K)	<50 ppm		

Mercury(II) iodide
specified reagents for general laboratory work- extra pure

SLR

packaging		Mfr. No
25 g	Plastic jar	M/2950/46
100 g	Plastic jar	M/2950/48
HgI ₂ CAS: 7774-29-0 MW: 454.39 Boiling point: 354°C Melting point: 259°C ADR/RID: 6.1, II EINECS: 231-873-8 UN: 1638		GHS: Acute Tox. 1, STOT RE 2, Aquatic Acute 1, Aquatic Chronic 1 H410, H373, H300, H310, H330 P301+P310, P280, P302+P350, P310, P304+P340, P273



Assay: >99%

Mercury(II) nitrate monohydrate
for analysis

Certified AR

packaging		Mfr. No
100 g	Plastic jar	M/3012/48
HgN ₂ O ₆ .H ₂ O CAS: 7783-34-8 MW: 342.60 ADR/RID: 6.1, II UN: 1625 GHS: Acute Tox. 1, STOT RE		2, Aquatic Acute 1, Aquatic Chronic 1 H373, H300, H410, H310, H330 P280, P302+P350, P310, P304+P340, P309+P311, P260, P273



Assay: >98%			
Guaranteed analysis			
Iron (Fe)	<10 ppm	n-Butyl formate	<0.002%
Sulfate (SO ₄)	<0.002%	Residue after ignition	<0.01%

Lot analysis on the label

Mercury(I) nitrate dihydrate
for analysis

Certified AR

packaging		Mfr. No
100 g	Plastic jar	M/3550/48
Hg ₂ (NO ₃) ₂ .2H ₂ O CAS: 7782-86-7 MW: 561.22 Density at 20°C: 4.78 Melting point: 70°C ADR/RID: 6.1, II UN: 1627		GHS: Acute Tox. 1, STOT RE 2, Aquatic Acute 1, Aquatic Chronic 1 H300, H310, H330, H410, H373 P260, P280, P301+P310, P302+P350, P304+P340, P273



Assay: >98%			
Guaranteed analysis			
Calcium (Ca)	<20 ppm	Sodium (Na)	<20 ppm
Copper (Cu)	<5 ppm	n-Butyl formate	<0.005%
Iron (Fe)	<5 ppm	Total phosphorus (P)	<20 ppm
Lead (Pb)	<5 ppm	Total silicon (Si)	<20 ppm
Magnesium (Mg)	<5 ppm	Total sulphur (S)	<50 ppm
Potassium (K)	<10 ppm	Zinc (Zn)	<10 ppm

Lot analysis on the label.

Optima LC/MS

Fisher Chemical OPTIMA® LC/MS solvents, blends, additives and modifiers have set the standard of excellence for consistent, reproducible performance in the mobile phase of LC/MS.

Find out more on page 308

Mercury(I) nitrate dihydrate
specified reagents for general laboratory work- extra pure

SLR

packaging		Mfr. No
500 g	Plastic jar	M/3500/53
Hg ₂ (NO ₃) ₂ .2H ₂ O CAS: 7782-86-7 MW: 561.22 Density at 20°C: 4.78 Melting point: 70°C ADR/RID: 6.1, II UN: 1627		GHS: Acute Tox. 1, STOT RE 2, Aquatic Acute 1, Aquatic Chronic 1 H300, H310, H330, H410, H373 P260, P280, P301+P310, P302+P350, P304+P340, P273



Assay: >95%			
Guaranteed analysis			
Calcium (Ca)	<100 ppm	Sodium (Na)	<100 ppm
Copper (Cu)	<50 ppm	n-Butyl formate	<0.01%
Iron (Fe)	<50 ppm	Total phosphorus (P)	<50 ppm
Lead (Pb)	<50 ppm	Total silicon (Si)	<50 ppm
Magnesium (Mg)	<50 ppm	Total sulphur (S)	<100 ppm
Potassium (K)	<50 ppm	Zinc (Zn)	<50 ppm

Mercury(II) oxide yellow
pure

packaging		Mfr. No
100 g	Plastic jar	M/3100/48
HgO CAS: 21908-53-2 MW: 216.59 Melting point: 500°C ADR/RID: 6.1, II EINECS: 244-654-7 UN: 1641		GHS: Acute Tox. 1, STOT RE 2, Aquatic Acute 1, Aquatic Chronic 1 H310, H410, H373, H330, H300 P301+P310, P280, P361, P302+P350, P310, P304+P340, P260, P273



Assay: >99%	
Guaranteed analysis	
n-Butyl formate	<0.2%

Mercury(II) oxide red
specified reagents for general laboratory work- extra pure

SLR

packaging		Mfr. No
100 g	Plastic jar	M/3050/48
HgO CAS: 21908-53-2 MW: 216.59 Melting point: 500°C ADR/RID: 6.1, II EINECS: 244-654-7 UN: 1641		GHS: Acute Tox. 1, STOT RE 2, Aquatic Acute 1, Aquatic Chronic 1 H310, H410, H373, H330, H300 P301+P310, P280, P361, P302+P350, P310, P304+P340, P260, P273



Assay: >99%

High Purity Acids
for Elemental Analysis

Fisher Chemical offers a complete range of acids to meet the most challenging applications. Choose from our Optima®, TraceMetal® and PrimarPlus® grades to ensure that the sensitivity you require is always available.

Read more page 24

Mercury(II) sulfate
for analysis

Certified AR

packaging	Mfr. No
250 g Plastic jar	M/3220/50
HgSO ₄ CAS: 7783-35-9 MW: 296.64 Density at 20°C: 6.47 ADR/RID: 6.1, II EINECS: 231-992-5 UN: 1645	GHS: Acute Tox. 1, STOT RE 2, Aquatic Acute 1, Aquatic Chronic 1 H410, H330, H300, H373, H310 P260, P280, P301+P310, P302+P350, P304+P340, P273
Assay: >99%	
Guaranteed analysis	
Calcium (Ca)<50 ppm Copper (Cu)<5 ppm Iron (Fe)<10 ppm Lead (Pb)<5 ppm Magnesium (Mg)<5 ppm Potassium (K)<10 ppm	Sodium (Na)<20 ppm n-Butyl formate<0.002% Total nitrogen (N)<50 ppm Total phosphorus (P)<20 ppm Total silicon (Si)<20 ppm Zinc (Zn)<10 ppm
Lot analysis on the label	

Mercury(II) sulfate

specified reagents for general laboratory work- extra pure

SLR

packaging	Mfr. No
100 g Plastic jar 250 g Plastic jar	M/3200/48 M/3200/50
HgSO ₄ CAS: 7783-35-9 MW: 296.64 Density at 20°C: 6.47 ADR/RID: 6.1, II EINECS: 231-992-5 UN: 1645	GHS: Acute Tox. 1, STOT RE 2, Aquatic Acute 1, Aquatic Chronic 1 H410, H330, H300, H373, H310 P260, P280, P301+P310, P302+P350, P304+P340, P273
Assay: >98%	
Guaranteed analysis	
Calcium (Ca)<100 ppm Copper (Cu)<50 ppm Iron (Fe)<50 ppm Lead (Pb)<50 ppm Magnesium (Mg)<50 ppm Potassium (K)<50 ppm	Sodium (Na)<100 ppm n-Butyl formate<0.01% Total phosphorus (P)<50 ppm Total silicon (Si)<50 ppm Zinc (Zn)<50 ppm

Mercury(II) sulfate solution 20% w/v

pure

packaging	Mfr. No
500 mL Glass bottle	J/5100/08
Density at 20°C: 1.3 ADR/RID: 8, 6.1, II UN: 2922 GHS: Acute Tox. 2, Skin Irrit. 2, Eye Irrit. 2, STOT RE 2, Aquatic Acute 1, Aquatic Chronic 1	H301, H330, H373, H410, H319, H315 P280, P260, P304+P340, P301+P310, P305+P351+P338, P302+P352, P310, P273
Guaranteed analysis	
Dichromate demand (mL FeSO ₄) ... <2 mL	

Mercury(II) thiocyanate

specified reagents for general laboratory work- extra pure

SLR

packaging	Mfr. No
25 g Plastic jar	M/3300/46
C ₂ HgN ₂ S ₂ CAS: 592-85-8 MW: 316.75 Melting point: 165°C ADR/RID: 6.1, II EINECS: 209-773-0 UN: 1646	GHS: Acute Tox. 1, STOT RE 2, Aquatic Acute 1, Aquatic Chronic 1 H330, H300, H310, H373, H410 P301+P310, P280, P302+P350, P310, P304+P340, P260, P273
Guaranteed analysis	
Identity test : Pass	

Metaphosphoric acid

specified reagents for general laboratory work- extra pure

SLR

packaging	Mfr. No
500 g Plastic jar	M/3851/53
HPO ₃ CAS: 37267-86-0 ADR/RID: 8, III EINECS: 253-433-4 UN: 3260 GHS: Skin Corr. 1B, Eye	Dam. 1, STOT SE 3 H335, H314 P280, P301+P330+P331, P310, P305+P351+P338, P303+P361+P353, P261
Guaranteed analysis	
Arsenic (As)<5 ppm Iron (Fe)<100 ppm	Lead (Pb)<10 ppm n-Butyl formate<0.002%

Methanesulfonic acid

For protein hydrolysis with reduced TRP losses

packaging	Mfr. No
100 g Coated Glass bottle	M/3874/48
CH ₄ O ₃ S CAS: 75-75-2 MW: 96.10 Density at 20°C: 1.481 Boiling point: 167°C Melting point: 19°C Flash point: 189°C ADR/RID: 8, 6.1, II	EINECS: 200-898-6 UN: 2922 GHS: Skin Corr. 1B H314 P280, P361, P302+P350, P301+P330+P331, P305+P351+P338
Assay: >98%	

Methanesulfonic acid, sodium salt

Ion pair chromatography

Certified

packaging	Mfr. No
25 g Plastic jar	M/3876/46
CH ₃ SO ₃ Na CAS: 2386-57-4 MW: 118.09 EINECS: 219-203-2 GHS: Skin Irrit. 2, Eye Irrit. 2,	STOT SE 3 H315, H319, H335 P302+P352, P280, P305+P351+P338
Assay: >95%	
Lot analysis on the label.	

Reagents, Solvents & Volumetric solutions comply with the European Pharmacopeia

Read more on page 234

Methanol

NEW

for LC/MS analysis

Optima

packaging	Mfr. No
500 mL Glass bottle 1 L Glass bottle 2.5 L Glass bottle 4 L Glass bottle	A456-500 A456-1 A456-212 A456-4
CH ₃ OH CAS: 67-56-1 MW: 32.04 Density at 20°C: 0.791 Boiling point: 64.7°C Melting point: -98°C Flash point: 12°C ADR/RID: 3, 6.1, II	EINECS: 200-659-6 UN: 1230 GHS: Flam. Liq. 2, Acute Tox. 3, STOT SE 1 H225, H301, H311, H331, H370 P301+P310, P280, P302+P350, P304+P340, P210, P240

Assay (CG): <99.9%

Guaranteed analysis

Acidity 0.0003 meq/g Alkalinity<0.0002 meq/g Water<0.02% Residue after evaporation<1 ppm Aluminium (Al)<10 ppb Barium (Ba)<10 ppb Cadmium (Cd)<10 ppb Calcium (Ca)<20 ppb Chromium (Cr)<10 ppb Cobalt (Co)<10 ppb Copper (Cu)<10 ppb	Iron (Fe)<10 ppb Lead (Pb)<10 ppb Magnesium (Mg)<10 ppb Manganese (Mn)<10 ppb Nickel (Ni)<10 ppb Potassium (K)<10 ppb Propan--1-ol<10 ppb Sodium (Na)<50 ppb Tin (Sn)<10 ppb Zinc (Zn)<10 ppb
---	--

Absorbance specification
Max. Absorbance (1cm cell v HPLC water)
210 nm<0.5 A.U.
214 nm: <0.4 A.U.
220 nm<0.2 A.U.
230 nm<0.1 A.U.
254 nm<0.01 A.U.
260 nm<0.005 A.U.
280 nm<0.005 A.U.
LC/MS at Negative Mode as Aldicarb<

50 ppb
LC/MS at Positive Mode as Reserpine<
50 ppb
LC/UV Gradient (single max peak) at 220nm<0.005 A.U.
LC/UV Gradient (single max peak) at 254nm<0.001 A.U.

This product has been tested on an LC-MS Thermo Scientific instrument to ensure suitability.
Filtered to 0.2µm.

Methanol

for HPLC-MS

Certified

packaging	Mfr. No
1 L Glass bottle 2.5 L Glass bottle	M/4062/15 M/4062/17
CH ₃ OH CAS: 67-56-1 MW: 32.04 Density at 20°C: 0.791 Boiling point: 64.7°C Melting point: -98°C Flash point: 12°C ADR/RID: 3, 6.1, II	EINECS: 200-659-6 UN: 1230 GHS: Flam. Liq. 2, Acute Tox. 3, STOT SE 1 H225, H301, H311, H331, H370 P301+P310, P280, P302+P350, P304+P340, P210, P240

Assay (CG): >99.9%

Guaranteed analysis

Acidity<0.0002 meq/g Alkalinity<0.0002 meq/g Water<0.02% Residue after evaporation<5 ppm Aluminium (Al)<0.2 ppm Barium (Ba)<0.01 ppm Cadmium (Cd)<0.01 ppm Calcium (Ca)<0.05 ppm Chromium (Cr)<0.01 ppm Cobalt (Co)<0.01 ppm Copper (Cu)<0.01 ppm	Iron (Fe)<0.01 ppm Lead (Pb)<0.01 ppm Magnesium (Mg)<0.01 ppm Manganese (Mn)<0.01 ppm Nickel (Ni)<0.01 ppm Potassium (K)<0.01 ppm Propan--1-ol<0.01 ppm Sodium (Na)<0.5 ppm Tin (Sn)<0.01 ppm Zinc (Zn)<0.01 ppm
---	---

Absorbance specification

Max. Absorbance (1cm cell v HPLC water)
230 nm<0.05 A.U.
260 nm<0.005 A.U.

Gradient (max eluted peak).<0.5 mA.U.
LC/MS Gradient Suitability: Pass test

Gradient analysis:
Absorbance of the largest eluted peak: <0.5 mA.U.
Conditions :
Column: C18 Reverse Phase.
Flow rate: 2 mL / min
Wavelength 254 nm
Absorbance range: 0-0.08 AU.
Solvent A: Water
Solvent B: Methanol
Gradient: linear 95% /5% B 100% B during 20 min.
Lot analysis on the label.
This product has been tested on an LC-MS instrument to ensure suitability.
Filtered to 0.2 µm.

Methanol

NEW

for UHPLC gradient grade analysis

Certified UHPLC

packaging	Mfr. No
1 L Coated Glass bottle 2.5 L Coated Glass bottle	M/4070/PB15 M/4070/PB17
CH ₃ OH CAS: 67-56-1 MW: 32.04 Density at 20°C: 0.791 Boiling point: 64.7°C Melting point: -98°C ADR/RID: 3, 6.1, II EINECS: 200-659-6	UN: 1230 GHS: Flam. Liq. 2, Acute Tox. 3, STOT SE 1 H225, H301, H311, H331, H370 P301+P310, P280, P302+P350, P304+P340, P210, P240

Assay (CG): >99.9%

Guaranteed analysis

Acidity 0.0002 meq/g Alkalinity<0.0001 meq/g	Water<0.02% Residue after evaporation<1 ppm
--	--

Absorbance specification
Max. Absorbance (1cm cell v HPLC water)
210 nm<0.2 A.U.
220 nm<0.1 A.U.
230 nm<0.05 A.U.

240 nm<0.01 A.U.
250 nm<0.005 A.U.
260 nm<0.005 A.U.

Gradient analysis

Absorbance of the largest eluted peak at 220nm: <0.1mA.U.
Absorbance of the largest eluted peak at 254nm: <0.1mA.U.
Baseline Drift at 220nm<150mA.U.
Conditions:
Column:C18 reverse phase
Wavelength: 220nm
Wavelength: 254nm
Solvent A: Water
Solvent B: Methanol
Gradient: Linear gradient from 95%A/5%B to 100%B
Suitable for UHPLC instrument.
Lot analysis on the label.
Filtered to 0.1µm

Label

Fisher Chemical labels show GHS pictograms clearly. Together with actual lot analysis shown on the label, plus a wealth of additional information, we help you to work safely to achieve accurate and consistent results in your analytical applications.

Methanol

HPLC for gradient analysis

Certified Advanced HPLC

packaging		Mfr. No	
1 ℓ	Coated Glass bottle	M/4060/PB15	
2.5 ℓ	Coated Glass bottle	M/4060/PB17	
CH ₃ OH	EINECS: 200-659-6	  	
CAS: 67-56-1	UN: 1230		
MW: 32.04	GHS: Flam. Liq. 2, Acute		
Density at 20°C: 0.791	Tox. 3, STOT SE 1		
Boiling point: 64.7°C	H225, H301, H311, H331, H370		
Melting point: -98°C	P301+P310, P280,	  	
Flash point: 12°C	P302+P350, P304+P340,		
ADR/RID: 3, 6.1, II	P210, P240		
Assay (CG): >99.9%			
Guaranteed analysis			
Acidity / Alkalinity	0.0002 meq/g	Residue after evaporation	
Water	<0.05%	<5 ppm	
Absorbance specification			
Max. Absorbance (1cm cell v HPLC water)			
210 nm	<0.2 A.U.	250 nm	
220 nm	<0.1 A.U.	260 nm	
230 nm	<0.05 A.U.	Drift in gradient baseline	
240 nm	<0.01 A.U.	<150 mA.U.	
Gradient analysis(30% to 95% of the solvent)			
Absorbance of the largest eluted peak: <0.1 mA.U.			
Conditions :			
Column: C18 Reverse Phase.			
Flow rate: 3 mL / min			
Wavelength 220 nm			
Solvent A: Water			
Solvent B: Methanol			
Gradient: 10 min at 30% solvent B followed by a linear gradient to 95% solvent B over 15 min.			
Lot analysis on the label.			
Filtered to 0.2 µm.			

Methanol HPLC gradient grade

NEW PACK.

Certified

packaging		Mfr. No	
1 ℓ	Glass bottle	M/4058/15	
2.5 ℓ	Glass bottle	M/4058/17	
2.5 ℓ	Coated Glass bottle	M/4058/PB17	
5 ℓ	Aluminium can	M/4058/21	
10 ℓ	Returnable Drum	M/4058/10RSS	
25 ℓ	Returnable Drum	M/4058/25RSS	
30 ℓ	Returnable Drum	M/4058/30RSS	
CH ₃ OH	EINECS: 200-659-6	  	
CAS: 67-56-1	UN: 1230		
MW: 32.04	GHS: Flam. Liq. 2, Acute		
Density at 20°C: 0.791	Tox. 3, STOT SE 1		
Boiling point: 64.7°C	H225, H301, H311, H331, H370		
Melting point: -98°C	P301+P310, P280,	  	
Flash point: 12°C	P302+P350, P304+P340,		
ADR/RID: 3, 6.1, II	P210, P240		
Assay (CG): >99.9%			
Guaranteed analysis			
Acidity / Alkalinity	0.0002 meq/g	Residue after evaporation	
Water	<0.05<0.03%	<5<5 ppm	
Absorbance specification			
Max. Absorbance (1cm cell v HPLC water)			
210 nm	<0.2 A.U.	250 nm	
220 nm	<0.1 A.U.	260 nm	
230 nm	<0.05 A.U.	Gradient (max eluted peak). ..	
240 nm	<0.01 A.U.	<0.001 A.U.	
Gradient analysis			
Absorbance of the largest eluted peak: <0.001A.U.			
Conditions :			
Column: C18 Reverse Phase.			
Flow rate: 2 mL / min			
Wavelength 254 nm			
Absorbance range: 0-0.08 AU.			
Solvent A: Water			

Methanol

HPLC for fluorescence detection

Certified

packaging		Mfr. No	
2.5 ℓ	Coated Glass bottle	M/4057/PB17	
CH ₃ OH	EINECS: 200-659-6	  	
CAS: 67-56-1	UN: 1230		
MW: 32.04	GHS: Flam. Liq. 2, Acute		
Density at 20°C: 0.791	Tox. 3, STOT SE 1		
Boiling point: 64.7°C	H225, H301, H311, H331, H370		
Melting point: -98°C	P301+P310, P280,	  	
Flash point: 12°C	P302+P350, P304+P340,		
ADR/RID: 3, 6.1, II	P210, P240		
Assay (CG): >99.8%			
Guaranteed analysis			
Acidity / Alkalinity	0.0002 meq/g	Residue after evaporation	
Water	<0.02%	<5 ppm	
Absorbance specification			
Max. Absorbance (1cm cell v HPLC water)			
210 nm	<0.2 A.U.	250 nm	
220 nm	<0.1 A.U.	260 nm	
230 nm	<0.05 A.U.	Fluorescence	
240 nm	<0.01 A.U.	Pass	
Excitation in the range 250-750nm produces no fluorescence emission greater than the solvents natural Raman signal at 350nm excitation.This is equivalent to the emission maximum of approximately 1ppb quinine sulfate.			
Lot analysis on the label (excluding fluorescence scan).			
Filtered to 0.2 µm.			

Solvent B: Methanol
Gradient: linear 95% /5% B 100% B during 20 min.
Lot analysis on the label (except gradient analysis).
Filtered to 0.2 µm.

Methanol

NEW PACK.

Certified

packaging		Mfr. No	
1 ℓ	Glass bottle	M/4056/15	
2.5 ℓ	Glass bottle	M/4056/17	
2.5 ℓ	Coated Glass bottle	M/4056/PB17	
5 ℓ	Aluminium can	M/4056/21	
10 ℓ	Returnable Drum	M/4056/10RSS	
25 ℓ	Returnable Drum	M/4056/25RSS	
30 ℓ	Returnable Drum	M/4056/30RSS	
CH ₃ OH	EINECS: 200-659-6	  	
CAS: 67-56-1	UN: 1230		
MW: 32.04	GHS: Flam. Liq. 2, Acute		
Density at 20°C: 0.791	Tox. 3, STOT SE 1		
Boiling point: 64.7°C	H225, H301, H311, H331, H370		
Melting point: -98°C	P301+P310, P280,	  	
Flash point: 12°C	P302+P350, P304+P340,		
ADR/RID: 3, 6.1, II	P210, P240		
Assay (CG): >99.8%			
Guaranteed analysis			
Acidity / Alkalinity	<0.0002 meq/g	Residue after evaporation	
Water	<0.02%	<5 ppm	
Absorbance specification			
Max. Absorbance (1cm cell v HPLC water)			
205 nm	<1 A.U.	240 nm	
210 nm	<0.2 A.U.	250 nm	
220 nm	<0.1 A.U.	260 nm	
230 nm	<0.05 A.U.	Wavelength at 1 A.U.	
Lot analysis on the label.		<205 nM	
Filtered to 0.2 µm.			

Methanol

HPLC for fluorescence detection

Certified

packaging		Mfr. No	
2.5 ℓ	Coated Glass bottle	M/4057/PB17	
CH ₃ OH	EINECS: 200-659-6	  	
CAS: 67-56-1	UN: 1230		
MW: 32.04	GHS: Flam. Liq. 2, Acute		
Density at 20°C: 0.791	Tox. 3, STOT SE 1		
Boiling point: 64.7°C	H225, H301, H311, H331, H370		
Melting point: -98°C	P301+P310, P280,	  	
Flash point: 12°C	P302+P350, P304+P340,		
ADR/RID: 3, 6.1, II	P210, P240		
Assay (CG): >99.8%			
Guaranteed analysis			
Acidity / Alkalinity	0.0002 meq/g	Residue after evaporation	
Water	<0.02%	<5 ppm	
Absorbance specification			
Max. Absorbance (1cm cell v HPLC water)			
210 nm	<0.2 A.U.	250 nm	
220 nm	<0.1 A.U.	260 nm	
230 nm	<0.05 A.U.	Fluorescence	
240 nm	<0.01 A.U.	Pass	
Excitation in the range 250-750nm produces no fluorescence emission greater than the solvents natural Raman signal at 350nm excitation.This is equivalent to the emission maximum of approximately 1ppb quinine sulfate.			
Lot analysis on the label (excluding fluorescence scan).			
Filtered to 0.2 µm.			

Excitation in the range 250-750nm produces no fluorescence emission greater than the solvents natural Raman signal at 350nm excitation.This is equivalent to the emission maximum of approximately 1ppb quinine sulfate.
Lot analysis on the label (excluding fluorescence scan).
Filtered to 0.2 µm.

Methanol

HPLC for electrochemical detection

Certified

packaging		Mfr. No	
2.5 ℓ	Coated Glass bottle	M/4059/PB17	
CH ₃ OH	EINECS: 200-659-6	  	
CAS: 67-56-1	UN: 1230		
MW: 32.04	GHS: Flam. Liq. 2, Acute		
Density at 20°C: 0.791	Tox. 3, STOT SE 1		
Boiling point: 64.7°C	H225, H301, H311, H331, H370		
Melting point: -98°C	P301+P310, P280,	  	
Flash point: 12°C	P302+P350, P304+P340,		
ADR/RID: 3, 6.1, II	P210, P240		
Assay (CG): >99.8%			
Guaranteed analysis			
Acidity / Alkalinity	0.0002 meq/g	Residue after evaporation	
Water	<0.02%	Electrochemical test	
Tested by HPLC-EC to ensure low electrochemical response.		Pass	
Detector - ESA Coulochem II			
Mode-oxidative			
Graphite Electrode.			
Lot analysis on the label.			
Filtered to 0.2 µm.			

Methanol

for trace metal analysis

Primar

packaging		Mfr. No	
2.5 ℓ	Glass bottle	M/4020/17	
CH ₃ OH	EINECS: 200-659-6	  	
CAS: 67-56-1	UN: 1230		
MW: 32.04	GHS: Flam. Liq. 2, Acute		
Density at 20°C: 0.791	Tox. 3, STOT SE 1		
Boiling point: 64.7°C	H225, H301, H311, H331, H370		
Melting point: -98°C	P301+P310, P280,	  	
Flash point: 12°C	P302+P350, P304+P340,		
ADR/RID: 3, 6.1, II	P210, P240		
Assay (CG): >99.9%			
Guaranteed analysis			
Acidity / Alkalinity	<0.0002 meq/g	Manganese (Mn)	
Colour	<5 APHA	Molybdenum (Mo)	
Water	<0.05%	Nickel (Ni)	
Residue after evaporation	<2 ppm	Potassium (K)	
Aluminium (Al)	<0.02 ppm	Propan--1-ol	
Barium (Ba)	<0.005 ppm	Sodium (Na)	
Cadmium (Cd)	<0.005 ppm	Strontium (Sr)	
Calcium (Ca)	<0.05 ppm	Titanium (Ti)	
Chromium (Cr)	<0.005 ppm	Total phosphorus (P)	
Cobalt (Co)	<0.005 ppm	Total silicon (Si)	
Copper (Cu)	<0.01 ppm	Total sulphur (S)	
Iron (Fe)	<0.05 ppm	Vanadium (V)	
Lead (Pb)	<0.005 ppm	Zinc (Zn)	
Lithium (Li)	<0.005 ppm	Propan--2-ol	
Magnesium (Mg)	<0.02 ppm	Ethanol 90 % v/v	

Methanol

for residue and pesticide analysis

Distol

packaging		Mfr. No	
1 ℓ	Glass bottle	M/4053/15	
2.5 ℓ	Glass bottle	M/4053/17	
CH ₃ OH	EINECS: 200-659-6	  	
CAS: 67-56-1	UN: 1230		
MW: 32.04	GHS: Flam. Liq. 2, Acute		
Density at 20°C: 0.791	Tox. 3, STOT SE 1		
Boiling point: 64.7°C	H225, H301, H311, H331, H370		
Melting point: -98°C	P301+P310, P280,	 	
Flash point: 12°C	P302+P350, P304+P340,		
ADR/RID: 3, 6.1, II	P210, P240		
Assay (CG): >99.9%			
Guaranteed analysis			
Acidity / Alkalinity	<0.0002 meq/g	Manganese (Mn)	
Colour	<5 APHA	Molybdenum (Mo)	
Water	<0.05%	Nickel (Ni)	
Residue after evaporation	<2 ppm	Potassium (K)	
Aluminium (Al)	<0.02 ppm	Propan--1-ol	
Barium (Ba)	<0.005 ppm	Sodium (Na)	
Cadmium (Cd)	<0.005 ppm	Strontium (Sr)	
Calcium (Ca)	<0.05 ppm	Titanium (Ti)	
Chromium (Cr)	<0.005 ppm	Total phosphorus (P)	
Cobalt (Co)	<0.005 ppm	Total silicon (Si)	
Copper (Cu)	<0.01 ppm	Total sulphur (S)	
Iron (Fe)	<0.05 ppm	Vanadium (V)	
Lead (Pb)	<0.005 ppm	Zinc (Zn)	
Lithium (Li)	<0.005 ppm	Propan--2-ol	
Magnesium (Mg)	<0.02 ppm	Ethanol 90 % v/v	

Guaranteed analysis

Water	<0.1%	GC (FID) Analysis	Pass
Residue after evaporation	<2 ppm	GC (NPD) analysis	Pass
GC (ECD) Analysis	Pass		

GC Analysis

Analysis of pesticides: from the retention time of Lindane, after concentration of 2500 times, no impurity peak higher than that obtained with 10 pg / mL heptachlo-reoxide (ECD) and 5 pg / mL Parathion (NPD).

Analysis of hydrocarbons: no impurity chain longer than the dodecane is present at a concentration greater than 0.01 ppm (FID).

Conditions :

Column: J&W DB 5: 30m x 0.32mm x 0.1µm df.

Temperature profile: (1) 50°C, 2 min

(2) 40°C / min until 320°C

(3) 320°C for 30 min

Injection: 5µL cold on the column.

Methanol

extra dry, for synthesis

packaging		Mfr. No
500 mL	Glass bottle	M/4050/08
1 ℓ	Glass bottle	M/4050/15
2.5 ℓ	Glass bottle	M/4050/17
2.5 ℓ	Coated Glass bottle	M/4050/PB17
25 ℓ	Metal drum	M/4050/25
CH ₃ OH	EINECS: 200-659-6	  
CAS: 67-56-1	UN: 1230	
MW: 32.04	GHS: Flam. Liq. 2, Acute	
Density at 20°C: 0.791	Tox. 3, STOT SE 1	
Boiling point: 64.7°C	H225, H301, H311, H331, H370	
Melting point: -98°C	P301+P310, P280,	  
Flash point: 12°C	P302+P350, P304+P340,	

Methanol **NEW PACK.** Certified AR
for analysis

packaging		Mfr. No
500 mℓ	Coated Glass bottle	M/4000/PB08
1 ℓ	Glass bottle	M/4000/15
1 ℓ	Coated Glass bottle	M/4000/PB15
2.5 ℓ	Glass bottle	M/4000/17
2.5 ℓ	Coated Glass bottle	M/4000/PB17
2.5 ℓ	HDPE Plastic bottle	M/4000/PC17
5 ℓ	Fluorinated HDPE Plastic bottle	M/4000/21
5 ℓ	Fluorinated HDPE square Drum	M/4000/FP21
25 ℓ	Metal drum	M/4000/25
25 ℓ	HDPE Plastic container	M/4000/PC25
200 ℓ	Metal drum	M/4000/27

CH ₃ OH	EINECS: 200-659-6	
CAS: 67-56-1	UN: 1230	
MW: 32.04	GHS: Flam. Liq. 2, Acute	
Density at 20°C: 0.791	Tox. 3, STOT SE 1	
Boiling point: 64.7°C	H225, H301, H311, H331, H370	
Melting point: -98°C	P301+P310, P280,	
Flash point: 12°C	P302+P350, P304+P340,	
ADR/RID: 3, 6.1, II	P210, P240	

Assay (CG): >99.9%

Guaranteed analysis			
Acidity / Alkalinity	<0.0005 meq/g	Iron (Fe)	<0.1 ppm
Colour	<5 APHA	Lead (Pb)	<0.02 ppm
Water	<0.05%	Magnesium (Mg)	<0.05 ppm
Residue after evaporation	<10 ppm	Potassium (K)	<0.5 ppm
Substances darkened		Sodium (Na)	<0.5 ppm
by (H ₂ SO ₄)	<5 APHA	Total phosphorus (P)	<0.05 ppm
Substances		Total silicon (Si)	<0.05 ppm
reducing KMnO ₄	<0.00025%	Total sulphur (S)	<0.1 ppm
Carbonyl compounds	: <0.001%	Zinc (Zn)	<0.05 ppm
Solution	: Pass	Ethanol 90 % v/v	<0.05%
Calcium (Ca)	<0.2 ppm	Propan--2-ol	<0.005%
Copper (Cu)	<0.02 ppm		

Typical analysis			
Aluminium (Al)	<0.05 ppm	Molybdenum (Mo)	<0.01 ppm
Barium (Ba)	<0.01 ppm	Nickel (Ni)	<0.01 ppm
Cadmium (Cd)	<0.02 ppm	Selenium (Se)	<0.01 ppm
Chromium (Cr)	<0.01 ppm	Strontium (Sr)	<0.01 ppm
Cobalt (Co)	<0.01 ppm	Titanium (Ti)	<0.01 ppm
Lithium (Li)	<0.01 ppm	Vanadium (V)	<0.01 ppm
Manganese (Mn)	<0.01 ppm		

Lot analysis on the label.

Methanol **NEW PACK.** SLR
specified reagents for general laboratory work- extra pure

packaging		Mfr. No
1 ℓ	Glass bottle	M/3950/15
2.5 ℓ	Glass bottle	M/3950/17
2.5 ℓ	Coated Glass bottle	M/3950/PB17
2.5 ℓ	HDPE Plastic bottle	M/3950/PC17
5 ℓ	Fluorinated HDPE Plastic bottle	M/3950/21
5 ℓ	Fluorinated HDPE square Drum	M/3950/FP21
25 ℓ	Metal drum	M/3950/25
25 ℓ	HDPE Plastic container	M/3950/PC25
200 ℓ	Metal drum	M/3950/27

CH ₃ OH	EINECS: 200-659-6	
CAS: 67-56-1	UN: 1230	
MW: 32.04	GHS: Flam. Liq. 2, Acute	
Density at 20°C: 0.791	Tox. 3, STOT SE 1	
Boiling point: 64.7°C	H225, H301, H311, H331, H370	
Melting point: -98°C	P301+P310, P280,	
Flash point: 12°C	P302+P350, P304+P340,	
ADR/RID: 3, 6.1, II	P210, P240	

Assay (CG): >99.5%

Guaranteed analysis			
Acidity / Alkalinity	<0.0005 meq/g	Magnesium (Mg)	<1 ppm
Colour	<10 APHA	Potassium (K)	<2 ppm
Water	<0.2%	Sodium (Na)	<2 ppm
Residue after evaporation	<50 ppm	Total phosphorus (P)	<1 ppm
Substances		Total silicon (Si)	<1 ppm
reducing KMnO ₄	<0.001%	Total sulphur (S)	<1 ppm
Calcium (Ca)	<2 ppm	Zinc (Zn)	<1 ppm
Copper (Cu)	<2 ppm	Substances darkened	
Iron (Fe)	<5 ppm	by H ₂ SO ₄	<20 APHA
Lead (Pb)	<2 ppm		

Methanol
Technical

packaging		Mfr. No
2.5 ℓ	Glass bottle	M/3900/17
2.5 ℓ	Coated Glass bottle	M/3900/PB17
25 ℓ	Metal drum	M/3900/25
200 ℓ	Metal drum	M/3900/27

CH ₃ OH	EINECS: 200-659-6	
CAS: 67-56-1	UN: 1230	
MW: 32.04	GHS: Flam. Liq. 2, Acute	
Density at 20°C: 0.791	Tox. 3, STOT SE 1	
Boiling point: 64.7°C	H225, H301, H311, H331, H370	
Melting point: -98°C	P301+P310, P280,	
Flash point: 12°C	P302+P350, P304+P340,	
ADR/RID: 3, 6.1, II	P210, P240	

Assay: > 98%

Methanol with 0.1% formic acid (v/v) **NEW** Certified
Methanol with 0.1% formic and HPLC gradient grade

packaging		Mfr. No
2.5 ℓ	Glass bottle	FE/1300/17
Density at 20°C: 0.781	GHS: Flam. Liq. 2, Acute	
Boiling point: 64.7°C	Tox. 3, STOT SE 1	
Melting point: -98°C	H225, H301, H311, H331, H370	
Flash point: 12°C	P301+P310, P280,	
ADR/RID: 3, 6.1, II	P302+P350, P304+P340,	
UN: 1230	P210, P240	

Guaranteed analysis		
Residue after evaporation	<=1 ppm	Gradient analysis at 254 nm. : <0.01 A.U.
Formic acid (% v/v)	: 0.093 to 0.107%	

Absorbance specification			
Max. Absorbance (1cm cell v HPLC water)			
230 nm	<1 A.U.	254 nm	<0.02 A.U.
240 nm	<0.27 A.U.	260 nm	<0.01 A.U.
250 nm	<0.05 A.U.		

Lot analysis on the label.
Filtered to 0.2µm.

Methenamine *See Hexamine, page 141*

L-Methionine
pure

packaging		Mfr. No
100 g	Plastic jar	M/4065/48

CH ₃ S(CH ₂) ₂ CHNH ₂ COOH	
CAS: 63-68-3	
MW: 149.21	
EINECS: 200-562-9	

Assay: >95%

Methoxybenzene *See Anisole, page 49*

2-Methoxyethanol Certified AR
for analysis

packaging		Mfr. No
500 mℓ	Glass bottle	M/4101/08
2.5 ℓ	Coated Glass bottle	M/4101/PB17
25 ℓ	Metal drum	M/4101/25

C ₃ H ₈ O ₂	EINECS: 203-713-7	
CAS: 109-86-4	UN: 1188	
MW: 76.09	GHS: Flam. Liq. 3, Repr. 1B,	
Density at 20°C: 0.960	Acute Tox. 4	
Boiling point: 124°C	H226, H312, H332,	
Melting point: -85°C	H360FD, H302	
Flash point: 38°C	P280, P201, P308+P313,	
ADR/RID: 3, III	P302+P352, P304+P340, P210	

Assay (CG): >99.5%

Guaranteed analysis			
Acidity / Alkalinity	<0.0008 meq/g	Magnesium (Mg)	<0.1 ppm
Colour	<10 APHA	Potassium (K)	<1 ppm
Water	<0.1%	Sodium (Na)	<1 ppm
Residue after evaporation	<20 ppm	Total phosphorus (P)	<0.1 ppm
Iodine consuming substances. :	<0.005%	Total silicon (Si)	<0.05 ppm
Calcium (Ca)	<0.5 ppm	Total sulphur (S)	<0.2 ppm
Copper (Cu)	<0.5 ppm	Zinc (Zn)	<1 ppm
Iron (Fe)	<1 ppm	Peroxide	<0.001%
Lead (Pb)	<0.5 ppm		

Typical analysis			
Aluminium (Al)	<0.2 ppm	Manganese (Mn)	<0.02 ppm
Barium (Ba)	<0.05 ppm	Molybdenum (Mo)	<0.02 ppm
Cadmium (Cd)	<0.05 ppm	Nickel (Ni)	<0.05 ppm
Chromium (Cr)	<0.05 ppm	Strontium (Sr)	<0.02 ppm
Cobalt (Co)	<0.05 ppm	Titanium (Ti)	<0.02 ppm
Lithium (Li)	<0.02 ppm	Vanadium (V)	<0.02 ppm

Lot analysis on the label.

2-Methoxyethanol SLR
specified reagents for general laboratory work- extra pure

packaging		Mfr. No
1 ℓ	Glass bottle	M/4100/15
2.5 ℓ	Coated Glass bottle	M/4100/PB17

C ₃ H ₈ O ₂	EINECS: 203-713-7	
CAS: 109-86-4	UN: 1188	
MW: 76.09	GHS: Flam. Liq. 3, Repr. 1B,	
Density at 20°C: 0.960	Acute Tox. 4	
Boiling point: 124°C	H226, H312, H332,	
Melting point: -85°C	H360FD, H302	
Flash point: 38°C	P280, P201, P308+P313,	
ADR/RID: 3, III	P302+P352, P304+P340, P210	

Assay (CG): >99%

Guaranteed analysis			
Acidity / Alkalinity	<0.005 meq/g	Magnesium (Mg)	<2 ppm
Colour	<20 APHA	Potassium (K)	<2 ppm
Water	<0.3%	Sodium (Na)	<5 ppm
Residue after evaporation	<50 ppm	Total phosphorus (P)	<1 ppm
Calcium (Ca)	<5 ppm	Total silicon (Si)	<1 ppm
Copper (Cu)	<5 ppm	Total sulphur (S)	<1 ppm
Iron (Fe)	<10 ppm	Zinc (Zn)	<2 ppm
Lead (Pb)	<5 ppm		

2-Methoxyphenol *See Guaiacol, page 137*

Methyl acetate SLR
specified reagents for general laboratory work- extra pure

packaging		Mfr. No
500 mℓ	Glass bottle	M/4150/08
2.5 ℓ	Glass bottle	M/4150/17

C ₃ H ₆ O ₂	EINECS: 201-185-2	
CAS: 79-20-9	UN: 1231	
MW: 74.08	GHS: Flam. Liq. 2, Eye Irrit.	
Density at 20°C: 0.930	2, STOT SE 3	
Boiling point: 57.4°C	H225, H336, H319, EUH066	
Melting point: -98°C	P261, P280,	
Flash point: -10°C	P305+P351+P338, P210,	
ADR/RID: 3, II	P240, P304+P340	

Assay: >98%

Guaranteed analysis		
Acidity	<0.0003 meq/g	Residue after evaporation
		<100 ppm

Methylamine solution 25/30 % w/v in water
technical

packaging		Mfr. No
500 mℓ	Glass bottle	M/4200/08

Density at 20°C: 0.900	Tox. 3, Skin Corr. 1B, Eye	
Boiling point: 48°C	Dam. 1, STOT SE 3	
Melting point: -38°C	H225, H314, H335, H301, H331	
Flash point: -18°C	P301+P310,	
ADR/RID: 3, 8, II	P301+P330+P331, P280,	
UN: 1235	P305+P351+P338, P310,	
GHS: Flam. Liq. 2, Acute	P304+P340, P261, P210, P240	

4-Methylamino phenol sulphate SLR
specified reagents for general laboratory work- extra pure

packaging		Mfr. No
250 g	Plastic jar	M/5400/50

C ₁₄ H ₁₈ N ₂ O ₂ ·H ₂ SO ₄	UN: 2811	
CAS: 55-55-0	GHS: Acute Tox. 4, Skin	
MW: 344.38	Sens. 1, STOT RE 2, Aquatic	
Density at 20°C: 1.577	Acute 1, Aquatic Chronic 1	
Melting point: 260°C	H410, H317, H373, H302	
Flash point: 260°C	P280, P260, P273, P301+P312	
ADR/RID: 6.1, III		
EINECS: 200-237-1		

Assay: >99%

Guaranteed analysis			
Calcium (Ca)	<200 ppm	Potassium (K)	<100 ppm
Copper (Cu)	<20 ppm	Sodium (Na)	<200 ppm
Iron (Fe)	<50 ppm	Total phosphorus (P)	<100 ppm
Lead (Pb)	<20 ppm	Total silicon (Si)	<50 ppm
Magnesium (Mg)	<100 ppm	Zinc (Zn)	<50 ppm

4-Methylaniline *See p-Toluidine, page 355*

We Support Your Chemistry For

- Semi-bulk and bulk requirements
- Tailored solvents and solvent blends or customized solutions
- Testing services or customized packaging and labeling

Please email your contact details and a brief description of your requirements to begeel.sha.quotes@thermofisher.com



Methyl orange solid C.I. 13025
pure, indicator grade

packaging		Mfr. No
25 g	Glass jar	M/5000/46
100 g	Glass jar	M/5000/48
C ₁₄ H ₁₄ N ₃ NaO ₃ S	EINECS: 208-925-3	
CAS: 547-58-0	UN: 3143	
MW: 327.34	GHS: Acute Tox. 3	
Melting point: >300°C	H301	
ADR/RID: 6.1, III	P301+P310, P270, P264	
Guaranteed analysis		
pH colour change test : Pass		Complies with standard IR spectrum: Pass
pH 2.9-4.6 red-orange to yellow-orange.		

Methyl orange solution 0.04%
pure, indicator grade

packaging			Mfr. No	
500 mL Glass bottle			M/5000L/08	
Flash point: 37°C		GHS: Flam. Liq. 3		
ADR/RID: 3, III		H226		
UN: 1170		P210		
Guaranteed analysis				
pH at 20°C		: 6.5 to 8	Sensitivity test (ACS test)	: Pass
pH colour change test		: Pass		
pH 2.9-4.6 red-orange to yellow-orange.				

Methyl orange-xylene cyanol FF
pure, indicator grade

packaging		Mfr. No
500 mL	Glass bottle	M/5020L/08
Flash point: 37°C	GHS: Flam. Liq. 2	
ADR/RID: 3, III	H225	
UN: 1170	P210, P240	
Guaranteed analysis		
pH colour change test		: Pass
pH 2.9-4.6 red to green.		
Solution methyl orange 0.2% / xylene cyanol 0.28%		

2-Methylpentan-2.4-diol
specified reagents for general laboratory work- extra pure

packaging		Mfr. No
500 mℓ	Glass bottle	M/5050/08
C ₆ H ₁₄ O ₂	Flash point: 93°C	
CAS: 107-41-5	EINECS: 203-489-0	
MW: 118.18	GHS: Skin Irrit. 2, Eye Irrit. 2	
Density at 20°C: 0.922	H315, H319	
Boiling point: 197°C	P302+P352, P280,	
Melting point: -40°C	P305+P351+P338	
Assay (CG): >98%		
Absorbance specification		
Max. Absorbance (1cm cell v HPLC water)		
Refractive index 1.427 to 1.428		

4-Methylpentan-2-one
for AAS

packaging		Mfr. No
500 mℓ	Glass bottle	M/5080/08
2.5 ℓ	Glass bottle	M/5080/17
C ₆ H ₁₂ O	EINECS: 203-550-1	
CAS: 108-10-1	UN: 1245	
MW: 100.16	GHS: Flam. Liq. 2, Acute	
Density at 20°C: 0.800	Tox. 4, Eye Irrit. 2, STOT SE 3	
Boiling point: 117.4°C	H225, H319, H335, H332,	
Melting point: -84°C	EUH066	
Flash point: 14°C	P304+P340, P261, P280,	
ADR/RID: 3, II	P305+P351+P338, P210, P240	

Guaranteed analysis

Calcium (Ca)	<0.05 ppm	Magnesium (Mg)	<0.01 ppm
Cobalt (Co)	<0.02 ppm	Nickel (Ni)	<0.02 ppm
Copper (Cu)	<0.02 ppm	Potassium (K)	<0.05 ppm
Iron (Fe)	<0.05 ppm	Sodium (Na)	<0.1 ppm
Lead (Pb)	<0.02 ppm	Zinc (Zn)	<0.1 ppm

4-Methylpentan-2-one
for analysis

packaging		Mfr. No
1 ℓ	Glass bottle	M/5085/15
1 ℓ	Coated Glass bottle	M/5085/PB15
2.5 ℓ	Glass bottle	M/5085/17
2.5 ℓ	Coated Glass bottle	M/5085/PB17
200 ℓ	Metal drum	M/5085/27
C ₆ H ₁₂ O		EINECS: 203-550-1
CAS: 108-10-1		UN: 1245
MW: 100.16		GHS: Flam. Liq. 2, Acute
Density at 20°C: 0.800		Tox. 4, Eye Irrit. 2, STOT SE 3
Boiling point: 117.4°C		H225, H319, H335, H332,
Melting point: -84°C		EUH066
Flash point: 14°C		P304+P340, P261, P280,
ADR/RID: 3, II		P305+P351+P338, P210, P240

Assay (CG): >99%

Guaranteed analysis

Acidity / Alkalinity	<0.0002 meq/g	Potassium (K)	<0.1 ppm
Colour	<10 APHA	Sodium (Na)	<0.5 ppm
Water	<0.1%	Total phosphorus (P)	<0.05 ppm
Residue after evaporation	<10 ppm	Total silicon (Si)	<0.02 ppm
Butan -one	<0.2%	Total sulphur (S)	<0.1 ppm
Calcium (Ca)	<0.5 ppm	Zinc (Zn)	<0.05 ppm
Copper (Cu)	<0.02 ppm	Acetone	<0.1%
Iron (Fe)	<0.1 ppm	2-Hexanol	<0.2%
Lead (Pb)	<0.02 ppm	Pentan--2-one	<0.2<0.05%
Magnesium (Mg)	<0.05 ppm		

Typical analysis

Aluminium (Al)	<0.1 ppm	Manganese (Mn)	<0.01 ppm
Barium (Ba)	<0.01 ppm	Molybdenum (Mo)	<0.01 ppm
Cadmium (Cd)	<0.01 ppm	Nickel (Ni)	<0.01 ppm
Chromium (Cr)	<0.01 ppm	Strontium (Sr)	<0.01 ppm
Cobalt (Co)	<0.01 ppm	Titanium (Ti)	<0.01 ppm
Lithium (Li)	<0.01 ppm	Vanadium (V)	<0.01 ppm

Lot analysis on the label.

4-Methylpentan-2-one
specified reagents for general laboratory work- extra pure

packaging		Mfr. No
500 mℓ	Glass bottle	M/5078/08
500 mℓ	Coated Glass bottle	M/5078/PB08
1 ℓ	Glass bottle	M/5078/15
2.5 ℓ	Glass bottle	M/5078/17
2.5 ℓ	Coated Glass bottle	M/5078/PB17
5 ℓ	Fluorinated HDPE Plastic bottle	M/5078/21
25 ℓ	Metal drum	M/5078/25
C ₆ H ₁₂ O	EINECS: 203-550-1	
CAS: 108-10-1	UN: 1245	
MW: 100.16	GHS: Flam. Liq. 2, Acute	
Density at 20°C: 0.800	Tox. 4, Eye Irrit. 2, STOT SE 3	
Boiling point: 117.4°C	H225, H319, H335, H332,	
Melting point: -84°C	EUH066	
Flash point: 14°C	P304+P340, P261, P280,	
ADR/RID: 3, II	P305+P351+P338, P210, P240	

Assay (CG): >99%

Guaranteed analysis

Acidity / Alkalinity	<0.0005 meq/g	Magnesium (Mg)	<2 ppm
Colour	<10 APHA	Potassium (K)	<2 ppm
Water	<0.2%	Sodium (Na)	<2 ppm
Residue after evaporation	<50 ppm	Total phosphorus (P)	<1 ppm
Calcium (Ca)	<5 ppm	Total silicon (Si)	<1 ppm
Copper (Cu)	<2 ppm	Total sulphur (S)	<2 ppm
Iron (Fe)	<5 ppm	Zinc (Zn)	<1 ppm
Lead (Pb)	<2 ppm		

2-Methylphenol
See o-Cresol 98-100 %, page 98

3-Methylphenol
See m-Cresol, page 98

2-Methyl-2-propanol
See tert-Butyl alcohol, page 69

2-Methylpropyl acetate
See iso-Butyl acetate, page 69

Methyl propyl carbinol
See Pentan-2-ol, page 227

N-Methylpyrrolidone
See 1-Methyl-2-pyrrolidinone, page 208

Methyl red spirit soluble C.I. 13020
pure, indicator grade

Packaging		Mfr. No
10 g	Glass jar	M/5150/45
25 g	Glass jar	M/5150/46
C ₁₅ H ₁₅ N ₃ O ₂	EINECS: 207-776-1	
CAS: 493-52-7	GHS: Acute Tox. 4	
MW: 269.30	H302	
Melting point: 178 to 182°C	P301+P312	
Guaranteed analysis		
pH colour change test		: Pass Complies with standard IR spectrum: Pass
pH 4.2-6.3 red to yellow.		

Methyl red solid water soluble
pure, indicator grade

Packaging			Mfr. No
5 g	Glass jar		M/5170/44
100 g	Glass jar		M/5170/48
$C_{15}H_{15}N_3O_2$		EINECS: 207-776-1	
CAS: 493-52-7		GHS: Acute Tox. 4	
MW: 269.30		H302	
Melting point: 178 to 182°C		P301+P312	
Guaranteed analysis			
pH colour change test		: Pass	Complies with standard IR spectrum; Pass
pH 4.2-6.3 red to yellow.			

Methyl red solution 0.01% water soluble
pure, indicator grade

Packaging		Mfr. No
500 mℓ	Glass bottle	M/5171L/08
Flash point: 37°C	GHS: Flam. Liq. 3, STOT SE 2	
ADR/RID: 3, III	H226, H371	
UN: 1170	P260, P210	
Guaranteed analysis		
pH at 20°C		: 6 to 8
pH 4.2-6.3 red to yellow.		
In solution in water 81%, ethanol 15% and methanol 4%.		

Methyl red solution 0.01% spirit soluble
pure, indicator grade

packaging		Mfr. No
500 mℓ	Glass bottle	M/5150L/08
Flash point: 16.66°C	GHS: Flam. Liq. 2	
ADR/RID: 3, II	H225	
UN: 1170	P210, P240	
Guaranteed analysis		
pH colour change test : Pass		
pH 4.2-6.3 red to yellow.		
In solution 50% ethanol 45% water and 5% methanol.		

Methyl red solution 0.025% water soluble
pure, indicator grade

packaging		Mfr. No
500 mℓ	Glass bottle	M/5170L/08
Flash point: 37°C	GHS: Flam. Liq. 3, STOT SE 2	
ADR/RID: 3, III	H226, H371	
UN: 1170	P260, P210	
Guaranteed analysis		
pH colour change test : Pass		
pH 4.2-6.3 red to yellow.		
In solution 81% water, 15% ethanol and 4% methanol.		

Methyl violet 2B (basic violet) C.I. 42535
pure, indicator grade

packaging		Mfr. No
10 g	Glass jar	M/5345/45
C ₂₄ H ₂₈ ClN ₃	UN: 2811	
CAS: 8004-87-3	GHS: Acute Tox. 4, Eye Irrit. 2	
MW: 393.95	H302, H319	
Melting point: 137°C	P301+P312, P280,	
ADR/RID: 6.1, III	P305+P351+P338	
Guaranteed analysis		
Complies with standard IR spectrum: Pass		

Methyl violet (crystal violet) C.I. 42555
pure, indicator grade

packaging		Mfr. No
100 g	Glass jar	M/5350/48
C ₂₅ H ₃₀ ClN ₄	GHS: Carc. 2, Acute Tox. 4,	
CAS: 548-62-9	Eye Dam. 1, Aquatic Acute	
MW: 407.97	1, Aquatic Chronic 1	
Melting point: 215°C	H410, H302, H318, H351	
ADR/RID: 9, III	P280, P305+P351+P338,	
EINECS: 208-953-6	P310, P201, P308+P313,	
UN: 3077	P273, P301+P312	

Reagents, Solvents & Volumetric solutions comply with the European Pharmacopeia

Read more on page 234



Methyl violet solid (crystal violet) C.I. 42555
pure, indicator grade

packaging		Mfr. No
10 g	Glass jar	M/5370/45
25 g	Glass jar	M/5370/46
C ₂₅ H ₃₀ ClN ₃	GHS: Carc. 1B, Acute Tox. 4, Eye Dam. 1, Aquatic Acute 1, Aquatic Chronic 1	
CAS: 90-94-8	H410, H302, H318, H350	
MW: 407.97	P280, P270,	
Melting point: 215°C	P305+P351+P338, P201,	
ADR/RID: 9, III	P308+P313, P273	
EINECS: 208-953-6		
UN: 3077		
Guaranteed analysis		
pH colour change test		: Pass
pH 0.1-3.2 yellow to purple.		

Methyl violet 1% aqueous solution (crystal violet) C.I 42555
pure, indicator grade

packaging		Mfr. No
500 mℓ	Glass bottle	M/5370L/08
GHS: Carc. 2, Eye Irrit. 2, Aquatic Chronic 3		
H351, H412, H319		
P281, P308+P313, P305+P351+P338, P273		
Guaranteed analysis		
pH at 20°C		: 3 to 5
pH colour change test		: Pass
<i>pH 0.1-3.2 yellow to purple.</i>		

1-Methyl-2-pyrrolidinone Certified AR
for analysis

packaging		Mfr. No	
500 mℓ	Glass bottle	M/5125/08	
25 ℓ	Metal drum	M/5125/25	
C ₅ H ₉ NO	EINECS: 212-828-1		
CAS: 872-50-4	GHS: Repr. 1B, Skin Irrit. 2, Eye Irrit. 2, STOT SE 3		
MW: 99.13	H319, H315, H335, H360D		
Density at 20°C: 1.03	P201, P308+P313, P261,		
Boiling point: 202°C	P302+P352, P280, P305+P351+P338		
Melting point: -24°C			
Flash point: 91°C			
Assay (CG): >99.7%			
Guaranteed analysis			
Water	<0.05%	Sodium (Na)	<5 ppm
Calcium (Ca)	<1 ppm	Total phosphorus (P)	<0.5 ppm
Copper (Cu)	<1 ppm	Total silicon (Si)	<0.1 ppm
Iron (Fe)	<1 ppm	Total sulphur (S)	<1 ppm
Lead (Pb)	<1 ppm	Zinc (Zn)	<0.2 ppm
Magnesium (Mg)	<0.5 ppm	Residue after ignition	<0.01%
Potassium (K)	<1 ppm		
Typical analysis			
Aluminium (Al)	<0.5 ppm	Manganese (Mn)	<0.1 ppm
Barium (Ba)	<0.1 ppm	Molybdenum (Mo)	<0.1 ppm
Cadmium (Cd)	<0.1 ppm	Nickel (Ni)	<0.1 ppm
Chromium (Cr)	<0.1 ppm	Strontium (Sr)	<0.1 ppm
Cobalt (Co)	<0.1 ppm	Titanium (Ti)	<0.1 ppm
Lithium (Li)	<0.1 ppm	Vanadium (V)	<0.1 ppm
Lot analysis on the label.			

1-Methyl-2-pyrrolidinone SLR
specified reagents for general laboratory work- extra pure

Packaging		Mfr. No	
1000 ml	Glass bottle	M/5120/05	
5000 ml	Glass bottle	M/5120/08	
25 l	Metal drum	M/5120/25	
C₅H₉NO CAS: 872-50-4 MW: 99.13 Density at 20°C: 1.03 Boiling point: 202°C Melting point: -24°C Flash point: 91°C		EINECS: 212-828-1 GHS: Repr. 1B, Skin Irrit. 2, Eye Irrit. 2, STOT SE 3 H319, H315, H335, H360D P201, P308+P313, P261, P302+P352, P280, P305+P351+P338 	
Assay (CG): >99%			
Guaranteed analysis			
Water	<0.2%	Sodium (Na)	<20 ppm
Calcium (Ca)	<5 ppm	Total phosphorus (P)	<2 ppm
Copper (Cu)	<5 ppm	Total silicon (Si)	<2 ppm
Iron (Fe)	<5 ppm	Total sulphur (S)	<5 ppm
Lead (Pb)	<5 ppm	Zinc (Zn)	<2 ppm
Magnesium (Mg)	<2 ppm	Residue after ignition	<0.1%
Potassium (K)	<5 ppm		

Methyl salicylate SLR
specified reagents for general laboratory work- extra pure

packaging		Mfr. No
500 mℓ	Coated Glass bottle	M/5200/PB08
2.5 ℓ	Coated Glass bottle	M/5200/PB17
C ₈ H ₈ O ₃	EINECS: 204-317-7	
CAS: 119-36-8	GHS: Skin Irrit. 2, Eye Irrit. 2, STOT SE 3	
MW: 152.15	H319, H315, H335	
Density at 20°C: 1.183	P261, P302+P352, P280,	
Boiling point: 222°C	P305+P351+P338	
Melting point: -8--7°C		
Flash point: 110°C		
Assay (CG): >99%		
Guaranteed analysis		
Acidity		<0.01 meq/g

MIBK See 4-Methylpentan-2-one, page 206

Millon's reagent
pure

Packaging		Mfr. No
1000 mL	Coated Glass bottle	J/5300/PB05
ADR/RID: 6.1, II		
UN: 2024		
GHS: Acute Tox. 1, Skin		
Corr. 1B, Eye Dam. 1, Resp.		
Sens. 1, STOT RE 2, Aquatic		
Acute 1, Aquatic Chronic 1		
H300, H310, H330, H410,		
H373, H314		
P260, P280, P301+P310,		
P302+P350, P304+P340, P273		

Mohr's salt

See Ammonium Iron (II) Sulfate hexahydrated, page 42

pH Buffers

NIST Certified Standard
Solution & Concentrate

A full range to discover on page 63

Molecular sieve type 13X (10A)
pure, particle size 1.6-2.5mm

packaging		Mfr. No
500 g	Plastic jar	MS/1070/53
CAS: 63231-69-6		
Guaranteed analysis		
Identity test		: Pass
Porosity 10 Angstroms.		
Density 0.71 g / mℓ.		

Molecular sieve type 4A
pure, particle size 1.6-2.5mm

packaging		Mfr. No
500 g	Plastic jar	MS/1030/53
2.5 kg	Plastic jar	MS/1030/62
CAS: 70955-01-0		
Guaranteed analysis		
Identity test		: Pass
Porosity 4 Angstroms.		
Density 0.74 g / mℓ.		

Molecular sieve type 4A
pure, particle size 2.5-5.0mm

packaging		Mfr. No
500 g	Plastic jar	MS/1020/53
2.5 kg	Plastic jar	MS/1020/62
CAS: 70955-01-0		
Guaranteed analysis		
Identity test		: Pass
Porosity 4 Angstroms.		
Density 0.72 g / mℓ.		

Molecular sieve type 3A
pure, particle size 1.6-2.5mm

packaging		Mfr. No
250 g	Plastic jar	MS/1010/50
500 g	Plastic jar	MS/1010/53
CAS: 308080-99-1		
Guaranteed analysis		
Identity test		: Pass
Porosity 3 Angstroms.		
Density 0.74 g / mℓ.		

Molybdenum standard metal solution 1000 ppm
for AAS, solution HNO₃ 0.1M

packaging		Mfr. No
100 mℓ	HDPE Plastic bottle	J/8050/05
500 mℓ	HDPE Plastic bottle	J/8050/08
Manufactured to a tolerance of ± 0.5%.		
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.		

Molybdenum standard metal solution 1000 ppm
for ICP, solution NH3 0.1M

packaging		Mfr. No
500 mℓ	HDPE Plastic bottle	J/8296/08
Manufactured to a tolerance of ± 0.5%.		
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.		
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.		

Molybdenum standard metal solution 10 000 ppm
for ICP, solution NH3 1%

Packaging		Mfr. No
100 ml	HDPE Plastic bottle	J/8294/05
GHS: Carc. 2 H351		
P281, P201, P202, P308+P313		
Manufactured to a tolerance of $\pm 0.5\%$.		
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.		
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.		

Molybdenum standard metal solution 1000 ppm
for ICP solution in H₂O, NIST certified

packaging		Mfr. No
100 mℓ	Plastic bottle	SPPLMO92Y
250 mℓ	Plastic bottle	SPPLMO92T
500 mℓ	Plastic bottle	SPPLMO92X
Manufactured to a tolerance of ± 0.5%.		
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.		
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.		
Check with your local distributor for availability of Spex product		

Molybdenum standard metal solution 10 000 ppm
for ICP solution in H₂O, NIST certified

packaging		Mfr. No
100 mℓ	Plastic bottle	SPPLMO93Y
500 mℓ	Plastic bottle	SPPLMO93X
Manufactured to a tolerance of ± 0.5%.		
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.		
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.		
Check with your local distributor for availability of Spex product		

Molybdenum standard metal solution 1000 ppm
for ICP-MS solution in H₂O, NIST certified

packaging		Mfr. No
100 mℓ	Plastic bottle	SPCLMO92Y
Manufactured to a tolerance of ± 0.5%.		
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.		
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.		
Check with your local distributor for availability of Spex product		

dodeca-Molybdophosphoric acid
pure

Packaging		Mfr. No
100 g	Plastic jar	M/5650/48
H ₃ Mo ₁₂ O ₄₀ P.xH ₂ O	GHS: Ox. Sol. 2, Skin Corr.	
CAS: 51429-74-4	1B, Eye Dam. 1	
MW: 1825.25	H272, H314	
ADR/RID: 5.1, 8, III	P301+P330+P331, P280,	
UN: 3085	P305+P351+P338, P310, P210	
Guaranteed analysis		
Iron (Fe).....	<200 ppm	Total sulphur (S) <1000 ppm
Lead (Pb).....	<100 ppm	

MonoFax 057 surfactant
pure, (mono and di-phosphate esters)

packaging	Mfr. No
500 mℓ Coated Glass bottle	M/5715/PB08
ADR/RID: 8/III UN: 1760	
Guaranteed analysis	
Neutralisation value	: 5 to 6

Mordant Black 11 See Eriochrome black T, page 119

Mordant Blue 3 See Eriochrome cyanine R, page 119

Mordant Red 11 See Alizarin, page 27

Morpholine
pure

packaging	Mfr. No
500 mℓ Coated Glass bottle	M/5750/PB08
C ₄ H ₉ NO CAS: 110-91-8 MW: 87.12 Density at 20°C: 0.990 Boiling point: 126 to 130°C Melting point: -5°C Flash point: 32°C ADR/RID: 8, 3, I EINECS: 203-815-1	UN: 2054 GHS: Flam. Liq. 3, Acute Tox. 4, Skin Corr. 1B, Eye Dam. 1 H226, H332, H302, H312, H314 P280, P312, P302+P350, P301+P330+P331, P305+P351+P338, P310, P304+P340, P210
Assay (CG): >99%	

2-(N-Morpholino)ethane-sulfonic acid
specified reagents for general laboratory work- extra pure

packaging	Mfr. No
100 g Plastic jar	M/5760/48
C ₆ H ₁₃ NO ₄ S.H ₂ O CAS: 145224-94-8 MW: 213.25 Melting point: 300°C GHS: Skin Irrit. 2, Eye Irrit. 2,	STOT SE 3 H315, H335, H319 P280, P305+P351+P338, P261, P302+P352, P264, P271
Assay: >99%	
Guaranteed analysis	
Residue after ignition	<0.2% Lead (Pb)
Iron (Fe)	<5 ppm Sulfate (SO ₄)
Absorbance specification	
Max. Absorbance (1cm cell v HPLC water)	
250 nm (sol. aqu. 10%)	

Naphtha solvent
Technical

packaging	Mfr. No
2.5 ℓ Coated Glass bottle 25 ℓ Metal drum	N/0050/PB17 N/0050/25
CAS: 8030-30-6 Density at 150°C: 0.85 Boiling point: 158 to 190°C Flash point: 42°C ADR/RID: 3, III EINECS: 232-443-2 UN: 1268	GHS: Flam. Liq. 3, STOT SE 3, Asp. Tox. 1, Aquatic Chronic 2 H226, H304, H335, H336, H411, EUH066 P210, P261, P273, P301+P310, P331, P305+P351+P338

Naphthalene
specified reagents for general laboratory work- extra pure

packaging	Mfr. No
500 g Glass jar	N/0150/53
C ₁₀ H ₈ CAS: 91-20-3 MW: 128.17 Density at 20°C: 0.990 Boiling point: 218°C Melting point: 79 to 82°C Flash point: 78°C ADR/RID: 4.1, III	EINECS: 202-049-5 UN: 1334 GHS: Carc. 2, Acute Tox. 4, Aquatic Acute 1, Aquatic Chronic 1 H351, H410, H302 P281, P273, P301+P312, P233
Assay (CG): >99%	
Guaranteed analysis	
Residue after evaporation	

Naphthalene black 10B
Technical

packaging	Mfr. No
10 g Glass jar	N/0211/45
C ₂₂ H ₁₄ N ₆ Na ₂ O ₉ S ₂ CAS: 1064-48-8 MW: 616.48 EINECS: 213-903-1 GHS: Skin Irrit. 2, Eye Irrit. 2,	STOT SE 3 H315, H319, H335 P261, P302+P352, P280, P305+P351+P338
Guaranteed analysis	
Identity test	

Naphthalene-2-sulfonic acid sodium salt
for lon pair chromatography

packaging	Mfr. No
25 g Plastic jar	N/0216/46
C ₁₀ H ₇ NaO ₃ S CAS: 532-02-5 MW: 230.22 EINECS: 208-523-8 GHS: Skin Irrit. 2, Eye Irrit. 2,	STOT SE 3 H315, H319, H335 P261, P302+P352, P280, P305+P351+P338
Guaranteed analysis	
Complies with standard IR spectrum: Pass	

1-Naphthol
for analysis

packaging	Mfr. No
100 g Glass jar	N/0350/48
C ₁₀ H ₈ O CAS: 90-15-3 MW: 144.17 Boiling point: 278 to 280°C Melting point: 95 to 97°C Flash point: 125°C ADR/RID: 6.1, III EINECS: 201-969-4	UN: 2811 GHS: Acute Tox. 4, Skin Irrit. 2, Eye Dam. 1, STOT SE 3 H315, H335, H302, H318, H312 P280, P302+P350, P305+P351+P338, P310, P261, P301+P312
Assay (CG): >99%	
Guaranteed analysis	
Insoluble in alcohol	: <0.01%
Insoluble matter in alcohol	<3 4
Aluminium (Al)	<2 ppm
Barium (Ba)	<2 ppm
Cadmium (Cd)	<2 ppm
Calcium (Ca)	<20 ppm
Chromium (Cr)	<2 ppm
Cobalt (Co)	<2 ppm
Copper (Cu)	<5 ppm
Iron (Fe)	<5 ppm
Lead (Pb)	<5 ppm
Lithium (Li)	<2 ppm
Magnesium (Mg)	<10 ppm
Manganese (Mn)	<2 ppm
Nickel (Ni)	<2 ppm
Potassium (K)	<10 ppm
Sodium (Na)	<20 ppm
Titanium (Ti)	<2 ppm
n-Butyl formate	<0.005%
Total phosphorus (P)	<5 ppm
Total silicon (Si)	<2 ppm
Total sulphur (S)	<10 ppm
Vanadium (V)	<5 ppm
Zinc (Zn)	<5 ppm
Naphthalene	<0.1%
2-Naphthol	<0.5%

1-Naphthol
specified reagents for general laboratory work- extra pure

packaging	Mfr. No
250 g Glass jar	N/0300/50
C ₁₀ H ₈ O CAS: 90-15-3 MW: 144.17 Boiling point: 278 to 280°C Melting point: 95 to 97°C Flash point: 125°C ADR/RID: 6.1, III EINECS: 201-969-4	UN: 2811 GHS: Acute Tox. 4, Skin Irrit. 2, Eye Dam. 1, STOT SE 3 H315, H335, H302, H318, H312 P280, P302+P350, P305+P351+P338, P310, P261, P301+P312
Assay (CG): >97%	
Guaranteed analysis	
Calcium (Ca)	<100 ppm
Copper (Cu)	<50 ppm
Iron (Fe)	<50 ppm
Lead (Pb)	<50 ppm
Magnesium (Mg)	<50 ppm
Potassium (K)	<50 ppm
Sodium (Na)	<100 ppm
n-Butyl formate	<0.02%
Total phosphorus (P)	<20 ppm
Total silicon (Si)	<10 ppm
Total sulphur (S)	<50 ppm
Zinc (Zn)	<20 ppm

2-Naphthol
specified reagents for general laboratory work- extra pure

packaging	Mfr. No
500 g Plastic jar	N/0400/53
C ₁₀ H ₈ O CAS: 135-19-3 MW: 144.17 Boiling point: 285 to 286°C Melting point: 120 to 124°C Flash point: 153°C ADR/RID: 9, III	EINECS: 205-182-7 UN: 3077 GHS: Acute Tox. 4, Aquatic Acute 1 H332, H302, H400 P304+P340, P301+P312, P273
Assay (CG): >97%	
Guaranteed analysis	
Calcium (Ca)	<100 ppm
Copper (Cu)	<50 ppm
Iron (Fe)	<50 ppm
Lead (Pb)	<50 ppm
Magnesium (Mg)	<50 ppm
Potassium (K)	<50 ppm
Sodium (Na)	<100 ppm
n-Butyl formate	<0.02%
Total phosphorus (P)	<20 ppm
Total silicon (Si)	<10 ppm
Total sulphur (S)	<50 ppm
Zinc (Zn)	<20 ppm

p-Naphtholbenzein
pure, indicator grade

packaging	Mfr. No
25 g Glass jar	N/0500/46
C ₂₇ H ₁₈ O ₂ CAS: 145-50-6 MW: 374.44 Melting point: 230 to 235°C EINECS: 205-656-3	GHS: Skin Irrit. 2, Eye Irrit. 2, STOT SE 3 H335, H315, H319 P261, P302+P352, P280, P305+P351+P338
Guaranteed analysis	
pH colour change test	
pH 8.5-9.8 yellow to green.	

1-Naphthylamine-7-sulfonic acid
specified reagents for general laboratory work- extra pure

packaging	Mfr. No
100 g Plastic jar	N/0910/48
C ₁₁ H ₈ (NH ₂).SO ₃ H CAS: 119-28-8 MW: 223.25 ADR/RID: 8, III EINECS: 204-311-4	UN: 2585 GHS: Skin Irrit. 2, Eye Irrit. 2 H319, H315 P280, P302+P352, P305+P351+P338
Guaranteed analysis	
Complies with standard IR spectrum: Pass	

N-1-Naphthylethylene diamine dihydrochloride
for analysis

packaging	Mfr. No
5 g Glass jar 10 g Glass jar	N/1252/44 N/1252/45
C ₁₂ H ₁₄ N ₂ .2HCl CAS: 1465-25-4 MW: 259.18 Melting point: 200°C EINECS: 215-981-2 GHS: Acute Tox. 4, Skin Irrit.	2, Eye Irrit. 2, STOT SE 3 H315, H302, H319, H335 P261, P301+P312, P302+P352, P280, P305+P351+P338
Assay: >99%	
Guaranteed analysis	
Sensitivity test (sulfanilamide)	
Residue after ignition	
Lot analysis on the label.	

N-1-Naphthylethylene diamine dihydrochloride
specified reagents for general laboratory work- extra pure

packaging	Mfr. No
10 g Glass jar 100 g Glass jar	N/1250/45 N/1250/48
C ₁₂ H ₁₄ N ₂ .2HCl CAS: 1465-25-4 MW: 259.18 Melting point: 200°C EINECS: 215-981-2 GHS: Acute Tox. 4, Skin Irrit.	2, Eye Irrit. 2, STOT SE 3 H315, H302, H319, H335 P261, P301+P312, P302+P352, P280, P305+P351+P338
Assay: >97%	
Guaranteed analysis	
Sensitivity test (sulfanilamide)	

Neodymium standard metal solution 1000 ppm
for ICP solution in HNO₃ 2%, NIST certified

packaging	Mfr. No
100 mℓ Plastic bottle 500 mℓ Plastic bottle	SPPLND22Y SPPLND22X
H315, H320	

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Nessler’s solution
pure, for detection and determination of ammonia and ammonium salts

packaging	Mfr. No
500 mℓ Coated Glass bottle 1 ℓ Coated Glass bottle	J/5500/PB08 J/5500/PB15
Density at 20°C: 1.1 to 1.3 ADR/RID: 8, 6.1, II UN: 2922 GHS: Acute Tox. 4, Skin Corr. 1A, Eye Dam. 1, STOT	RE 2, Aquatic Chronic 3 H412, H302, H373, H314 P301+P330+P331, P280, P305+P351+P338, P310, P260, P273, P301+P312

Neutral red C.I 50040
pure, indicator grade

packaging		Mfr. No
5 g	Glass jar	N/1270/44
C ₁₅ H ₁₇ ClN ₄	EINECS: 209-035-8	
CAS: 553-24-2	GHS: Muta. 2	
MW: 288.77	H341	
Melting point: 290°C	P281	
Guaranteed analysis		
pH colour change test : Pass Complies with standard IR spectrum: Pass		

Nickel standard metal solution 1000 ppm
for AAS, solution HNO₃ 1M

packaging		Mfr. No
100 mℓ	HDPE Plastic bottle	J/8055/05
500 mℓ	HDPE Plastic bottle	J/8055/08
Density at 20°C: 1.0	H314	
ADR/RID: 8, II	P280, P305+P351+P338,	
UN: 2031	P301+P330+P331, P310	
GHS: Skin Corr. 1B, Eye Dam. 1		
Manufactured to a tolerance of ± 0.5%.		
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.		

Nickel standard metal solution 1000 ppm
for ICP, solution HNO₃ 1M

packaging		Mfr. No
500 mℓ	HDPE Plastic bottle	J/8297/08
Density at 20°C: 1.0	H314	
ADR/RID: 8, II	P280, P305+P351+P338,	
UN: 2031	P301+P330+P331, P310	
GHS: Skin Corr. 1B, Eye Dam. 1		
Manufactured to a tolerance of ± 0.5%.		
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.		
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.		

Nickel standard metal solution 10 000 ppm
for ICP, solution HNO₃ 1M

packaging		Mfr. No
100 mℓ	HDPE Plastic bottle	J/8295/05
Density at 20°C: 1.0	H314, H372, H317	
ADR/RID: 8, II	P305+P351+P338,	
UN: 2031	P301+P330+P331, P310,	
GHS: Skin Corr. 1B, Skin	P280, P260	
Sens. 1, STOT RE 1		

Manufactured to a tolerance of ± 0.5%.

All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.

See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.

High Volume Solvent Delivery Systems

High-volume solvent delivery systems are available in 10L to 1000L, in Europe only. Please email your contact details and a brief description of your requirements to begal.sha.quotes@thermofisher.com



Nickel standard metal solution 1000 ppm
for ICP solution in HNO₃ 2%, NIST certified

Packaging		Mfr. No
100 mL	Plastic bottle	SPPLNI22Y
250 mL	Plastic bottle	SPPLNI22T
500 mL	Plastic bottle	SPPLNI22X
H315, H319		

Nickel standard metal solution 10 000 ppm
for ICP solution in HNO₃ 5%, NIST certified

packaging		Mfr. No
100 mL	Plastic bottle	SPPLNI23Y
500 mL	Plastic bottle	SPPLNI23X
H314, H317, H351		  

Manufactured to a tolerance of ± 0.5%.

All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.

See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.

Check with your local distributor for availability of Spex product

Nickel standard metal solution 1000 ppm
for ICP-MS solution in HNO₃ 2%, NIST certified

Packaging		Mfr. No
100 mL	Plastic bottle	SPCLNI22Y
H315, H319		
Manufactured to a tolerance of ± 0.5%.		
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.		
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.		
Check with your local distributor for availability of Spex product		

Nickel metal powder
specified reagents for general laboratory work- extra pure

Packaging		Mfr. No
250 g	Plastic jar	N/1300/50
Ni	EINECS: 231-111-4	
CAS: 7440-02-0	UN: 3089	
MW: 58.70	GHS: Carc. 2, Skin Sens. 1, STOT RE 1, Aquatic Chronic 3	
Boiling point: 2730°C	H351, H372, H412, H317	
Melting point: 1455°C	P280, P260, P273	
ADR/RID: 4.1/II		
Guaranteed analysis		
Identity test : Pass		

Nickel metal sheet
pure, about 0.5 mm

Packaging		Mfr. No
100 g	Plastic jar	N/1370/48
250 g	Plastic jar	N/1370/50
Ni	EINECS: 231-111-4	
CAS: 7440-02-0	GHS: Carc. 2, Skin Sens. 1, STOT RE 1	
MW: 58.70	H351, H372, H317	
Boiling point: 2732 to 2837°C	P280, P281, P260, P273	
Melting point: 1455°C		
Guaranteed analysis		
Identity test : Pass		

Nickel-aluminium alloy
pure, powder for preparing catalytic nickel

Packaging		Mfr. No
250 g	Plastic jar	N/1460/50
AlNi	UN: 3208	
CAS: 12635-29-9	GHS: Water-react. 3, Carc.	
MW: 85.68	2, Resp. Sens. 1, Skin Sens. 1	
Density at 20°C: 3.46	H261, H351, H317, H334	
Melting point: 1350°C	P261, P304+P341,	
ADR/RID: 4.3, III	P342+P311, P280, P281	
Guaranteed analysis		
Identity test : Pass		

Nickel(II) carbonate tetrahydrate
specified reagents for general laboratory work- extra pure

Packaging		Mfr. No
250 g	Plastic jar	N/1550/50
500 g	Plastic jar	N/1550/53
approx. $\text{CH}_4\text{Ni}_3\text{O}_7 \cdot 4\text{H}_2\text{O}$		
CAS: 12607-70-4		
ADR/RID: 9, III		
EINECS: 235-715-9		
UN: 3077		
GHS: Muta. 2, Repr. 1B,		
Acute Tox. 4, Skin Irrit. 2,		
Resp. Sens. 1, Skin Sens. 1,		
Guaranteed analysis		
Identity test : Pass		

Nickel(II) chloride hexahydrate
for analysis

Packaging		Mfr. No
250 g	Plastic jar	N/1650/50
1 kg	Plastic jar	N/1650/60
NiCl ₂ ·6H ₂ O	Resp. Sens. 1, Skin Sens. 1,	
CAS: 7791-20-0	STOT RE 1, Aquatic Acute 1,	
MW: 237.71	Aquatic Chronic 1	
Density at 20°C: 3.55	H301, H331, H334, H341,	
ADR/RID: 6.1, III	H410, H317, H315, H372,	
EINECS: 231-743-0	H350i, H360D	
UN: 3288	P301+P310, P308+P313,	
GHS: Muta. 2, Repr. 1B,	P304+P340, P280, P260,	
Acute Tox. 3, Skin Irrit. 2,	P273, P302+P352	
Assay: 99 to 101%		

Guaranteed analysis			
pH (aqueous solution 5%)	3 to 7	Magnesium (Mg)	<20 ppm
Aluminium (Al)	<0.001%	Potassium (K)	<50 ppm
Barium (Ba)	<0.001%	Propan--1-ol	<0.001%
Cadmium (Cd)	<0.001%	Sodium (Na)	<100 ppm
Calcium (Ca)	<50 ppm	Strontium (Sr)	<0.001%
Chromium (Cr)	<0.001%	Titanium (Ti)	<0.001%
Cobalt (Co)	<50 ppm	Total phosphorus (P)	<20 ppm
Copper (Cu)	<10 ppm	Total silicon (Si)	<10 ppm
Iron (Fe)	<20 ppm	Total sulphur (S)	<50 ppm
Lead (Pb)	<20 ppm	Vanadium (V)	<0.001%
Lithium (Li)	<0.001%	Zinc (Zn)	<10 ppm

Lot analysis on the label.

Nickel(II) chloride hexahydrate
specified reagents for general laboratory work- extra pure

packaging		Mfr. No
500 g	Plastic jar	N/1600/53
NiCl ₂ ·6H ₂ O	Resp. Sens. 1, Skin Sens. 1,	
CAS: 7791-20-0	STOT RE 1, Aquatic Acute 1,	
MW: 237.71	Aquatic Chronic 1	
Density at 20°C: 3.55	H301, H331, H334, H341,	
ADR/RID: 6.1, III	H410, H317, H315, H372,	
EINECS: 231-743-0	H350i, H360D	
UN: 3288	P301+P310, P308+P313,	
GHS: Muta. 2, Repr. 1B,	P304+P340, P280, P260,	
Acute Tox. 3, Skin Irrit. 2,	P273, P302+P352	

Guaranteed analysis			
Calcium (Ca)	<500 ppm	Potassium (K)	<500 ppm
Cobalt (Co)	<2000 ppm	Sodium (Na)	<2000 ppm
Copper (Cu)	<100 ppm	Total phosphorus (P)	<200 ppm
Iron (Fe)	<200 ppm	Total silicon (Si)	<200 ppm
Lead (Pb)	<200 ppm	Total sulphur (S)	<500 ppm
Magnesium (Mg)	<200 ppm	Zinc (Zn)	<100 ppm

Nickel(II) nitrate hexahydrate
specified reagents for general laboratory work- extra pure

packaging		Mfr. No
250 g	Plastic jar	N/1750/50
500 g	Plastic jar	N/1750/53
Ni(NO ₃) ₂ ·6H ₂ O	Repr. 1B, Acute Tox. 4, Skin	
CAS: 13478-00-7	Irrit. 2, Eye Dam. 1, Resp.	
MW: 290.79	Sens. 1, Skin Sens. 1, STOT	
Boiling point: 137°C	RE 1, Aquatic Acute 1, Aquatic Chronic 1	
Melting point: 56.7°C		
ADR/RID: 5.1, III	H272, H350i, H341, H360D,	
EINECS: 236-068-5	H372, H332, H302, H318,	
UN: 2725	H315, H317, H410	
GHS: Ox. Sol. 2, Muta. 2,	P281, P210	

Guaranteed analysis			
Calcium (Ca)	<500 ppm	Sodium (Na)	<1000 ppm
Cobalt (Co)	<2000 ppm	n-Butyl formate	<0.05%
Copper (Cu)	<100 ppm	Total phosphorus (P)	<200 ppm
Iron (Fe)	<200 ppm	Total silicon (Si)	<200 ppm
Lead (Pb)	<200 ppm	Total sulphur (S)	<500 ppm
Magnesium (Mg)	<200 ppm	Zinc (Zn)	<100 ppm
Potassium (K)	<500 ppm		

Nickel(II) oxide (black) <i>pure</i>		
packaging	Mfr. No	
500 g Plastic jar	N/1850/53	
NiO CAS: 1313-99-1 MW: 74.71 Density at 20°C: 6.67 Melting point: 1960°C ADR/RID: 6.1/III	EINECS: 215-215-7 UN: 3288 GHS: Skin Sens. 1, STOT RE 1, Aquatic Chronic 4 H350i, H372, H317, H413 P280, P201, P308+P313, P260	
Guaranteed analysis		
Identity test : Pass		

Nicotinic acid <i>specified reagents for general laboratory work- extra pure</i>			SLR
packaging	Mfr. No		
250 g Plastic jar	N/2110/50		
C ₆ H ₅ NO ₂ CAS: 59-67-6 MW: 123.11 Melting point: 236 to 239°C	EINECS: 200-441-0 GHS: Eye Irrit. 2 H319 P280, P305+P351+P338		
Assay: >99.5%			
Guaranteed analysis			
Water	<1%	Melting Point tested	234 to 237°C
Residue after ignition	<0.2%		

Nickel(II) sulfate hexahydrate <i>for analysis</i>		Certified AR
packaging	Mfr. No	
250 g Plastic jar	N/2051/50	
NiSO ₄ ·6H ₂ O CAS: 10101-97-0 MW: 262.85 Melting point: 53°C ADR/RID: 6.1, III UN: 3288 GHS: Muta. 2, Repr. 1B, Acute Tox. 4, Skin Irrit. 2,	Resp. Sens. 1, Skin Sens. 1, STOT RE 1, Aquatic Acute 1, Aquatic Chronic 1 H350i, H341, H360D, H315, H332, H302, H317, H334, H372, H410 P308+P313, P280, P281, P201, P260, P273	
Assay: 99 to 101%		
Guaranteed analysis		
pH (aqueous solution 5%)	4 to 6	Potassium (K)<50 ppm
Calcium (Ca)	<50 ppm	Sodium (Na)<100 ppm
Cobalt (Co)	<50 ppm	Total chloride (Cl)<0.001%
Copper (Cu)	<10 ppm	Total nitrogen (N)<50 ppm
Iron (Fe)	<20 ppm	Total phosphorus (P)<20 ppm
Lead (Pb)	<20 ppm	Total silicon (Si)<10 ppm
Magnesium (Mg)	<100 ppm	Zinc (Zn)<10 ppm
Typical analysis		
Aluminium (Al)	<0.001 %	Silver (Ag)<0.001 %
Barium (Ba)	<0.001 %	Strontium (Sr)<0.001 %
Cadmium (Cd)	<0.001 %	Titanium (Ti)<0.001 %
Chromium (Cr)	<0.001 %	Vanadium (V)<0.001 %
Lithium (Li)	<0.001 %	

Lot analysis on the label.

Nickel(II) sulfate hexahydrate <i>specified reagents for general laboratory work- extra pure</i>			SLR
packaging	Mfr. No		
250 g Plastic jar 500 g Plastic jar 1 kg Plastic jar	N/2000/50 N/2000/53 N/2000/60		
NiSO ₄ ·6H ₂ O CAS: 10101-97-0 MW: 262.85 Melting point: 53°C ADR/RID: 6.1, III UN: 3288 GHS: Muta. 2, Repr. 1B, Acute Tox. 4, Skin Irrit. 2,	Resp. Sens. 1, Skin Sens. 1, STOT RE 1, Aquatic Acute 1, Aquatic Chronic 1 H350i, H341, H360D, H315, H332, H302, H317, H334, H372, H410 P308+P313, P280, P281, P201, P260, P273	 	
Assay: >98%			
Guaranteed analysis			
Calcium (Ca)	<500 ppm	Potassium (K)<500 ppm	
Cobalt (Co)	<2000 ppm	Sodium (Na)<1000 ppm	
Copper (Cu)	<100 ppm	n-Butyl formate<0.01%	
Iron (Fe)	<200 ppm	Total phosphorus (P)<200 ppm	
Lead (Pb)	<200 ppm	Total silicon (Si)<200 ppm	
Magnesium (Mg)	<200 ppm	Zinc (Zn)<100 ppm	

Ninhydrin monohydrate <i>specified reagents for general laboratory work- extra pure</i>			SLR
packaging	Mfr. No		
100 g Glass jar 500 g Glass jar	N/2158/48 N/2158/53		
C ₉ H ₆ O ₄ CAS: 485-47-2 MW: 178.15 Melting point: 250°C EINECS: 207-618-1 GHS: Acute Tox. 4, Skin Irrit.	2, Eye Irrit. 2, STOT SE 3 H319, H315, H302, H335 P261, P301+P312, P302+P352, P280, P305+P351+P338		
Assay: >98%			
Guaranteed analysis			
Residue after ignition	<0.05%	Pass	
Suitability for detection of amino acid ... :			

Ninhydrin spray <i>pure</i>		
packaging	Mfr. No	
120 ml Aerosol	N/2160/90	
CAS: 71-36-3 Density at 20°C: 0.81 Boiling point: 117.4°C Melting point: -90°C Flash point: 37°C ADR/RID: 3, III EINECS: 207-618-1 UN: 1120	GHS: Flam. Liq. 3, Acute Tox. 4, Skin Irrit. 2, Eye Dam. 1, STOT SE 3 H226, H315, H302, H336, H318, H335 P280, P305+P351+P338, P310, P261, P301+P312, P304+P340, P302+P352, P210	
Guaranteed analysis		
Residue after evaporation	0.4 to 0.6%	

Pressurise container carefullyby pumping the removable piston several times.

Ninhydrin spray refill <i>pure, 0.5% ninhydrin in butan-1-ol</i>		
packaging	Mfr. No	
500 ml Glass bottle	N/2161/90	
CAS: 71-36-3 Density at 20°C: 0.81 Boiling point: 117.4°C Melting point: -90°C Flash point: 37°C ADR/RID: 3, III EINECS: 207-618-1 UN: 1120	GHS: Flam. Liq. 3, Acute Tox. 4, Skin Irrit. 2, Eye Dam. 1, STOT SE 3 H226, H315, H302, H336, H318, H335 P280, P305+P351+P338, P310, P261, P301+P312, P304+P340, P302+P352, P210	
Guaranteed analysis		
Residue after evaporation	0.4 to 0.6%	

Niobium standard metal solution 1000 ppm <i>for ICP solution in H₂O/HF trace, NIST certified</i>		
packaging	Mfr. No	
100 ml Plastic bottle 500 ml Plastic bottle H319	SPPLNB92Y SPPLNB92X	

*Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurites.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.*

Niobium standard metal solution 10 000 ppm <i>for ICP solution in H₂O/HF trace, NIST certified</i>		
packaging	Mfr. No	
100 ml Plastic bottle 500 ml Plastic bottle H319	SPPLNB93Y SPPLNB93X	

*Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurites.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product*

Nitrate standard solution 1000 ppm (1 mg/ml) <i>for Ion pair chromatography, Ready-to-use</i>		
packaging	Mfr. No	
100 ml HDPE Plastic bottle	J/4556/05	
CAS: 7757-79-1 EINECS: 231-818-8		

Manufactured to a tolerance of ± 1%.

Nitrate standard 1000 ppm <i>For Ion chromatography in H₂O, NIST certified</i>		
packaging	Mfr. No	
100 ml Plastic bottle 500 ml Plastic bottle CAS: 7757-79-1	SPASNO392Y SPASNO392X	

*Manufactured to a tolerance of ± 0.3%.
Check with your local distributor for availability of Spex product*

Nitrate-Nitrogen standard 1000 ppm <i>For Ion chromatography in H₂O, NIST certified</i>		
packaging	Mfr. No	
100 ml Plastic bottle CAS: 7757-79-1	SPASNO3N92Y	

*Manufactured to a tolerance of ± 0.3%.
Check with your local distributor for availability of Spex product*

Nitrite standard 1000 ppm <i>For Ion chromatography in H₂O, NIST certified</i>		
packaging	Mfr. No	
100 ml Plastic bottle 500 ml Plastic bottle CAS: 7632-00-0	SPASNO292Y SPASNO292X	

*Manufactured to a tolerance of ± 0.3%.
Check with your local distributor for availability of Spex product*

Nitrite-Nitrogen standard 1000 ppm <i>For Ion chromatography in H₂O, NIST certified</i>		
packaging	Mfr. No	
100 ml Plastic bottle 500 ml Plastic bottle CAS: 7632-00-0	SPASNO2N92Y SPASNO2N92X	

*Manufactured to a tolerance of ± 0.3%.
Check with your local distributor for availability of Spex product*

Nitric acid 67-69% <i>for ultra trace elemental analysis</i>			Optima
packaging	Mfr. No		
250 ml Pre-cleaned Teflon FEP bottle 500 ml Pre-cleaned Teflon FEP bottle 1 l Pre-cleaned Teflon FEP bottle	A467-250 A467-500 A467-1		
HNO ₃ CAS: 7697-37-2 MW: 63.01 Density at 20°C: 1.42 Boiling point: 122°C Melting point: -42°C ADR/RID: 8, 5.1, II	EINECS: 231-714-2 UN: 2031 GHS: Ox. Liq. 3, Skin Corr. 1A H314, H272 P301+P330+P331, P305+P351+P338, P310, P304+P340, P280, P210	 	
Assay: 67-69%			

Guaranteed analysis			
Zirconium (Zr)	<10 ppt	Manganese (Mn)	<10 ppt
Gallium (Ga)	<10 ppt	Mercury (Hg)	<50 ppt
Terbium (Tb)	<1 ppt	Molybdenum (Mo)	<10 ppt
Indium (In)	<1 ppt	Neodymium (Nd)	<1 ppt
Dysprosium (Dy)	<1 ppt	Nickel (Ni)	<20 ppt
Aluminium (Al)	<20 ppt	Niobium (Nb)	<1 ppt
Antimony (Sb)	<10 ppt	Palladium (Pd)	<20 ppt
Arsenic (As)	<20 ppt	Potassium (K)	<10 ppt
Barium (Ba)	<10 ppt	Praseodymium (Pr)	<1 ppt
Beryllium (Be)	<10 ppt	Rhenium (Re)	<10 ppt
Bismuth (Bi)	<10 ppt	Rhodium (Rh)	<10 ppt
Boron (B)	<10 ppt	Ruthenium (Ru)	<20 ppt
Cadmium (Cd)	<10 ppt	Samarium (Sm)	<1 ppt
Calcium (Ca)	<10 ppt	Scandium (Sc)	<10 ppt
Cerium (Ce)	<10 ppt	Propan--1-ol	<10 ppt
Cesium (Cs)	<10 ppt	Sodium (Na)	<10 ppt
Chromium (Cr)	<10 ppt	Strontium (Sr)	<10 ppt
Cobalt (Co)	<10 ppt	Tantalum (Ta)	<10 ppt
Copper (Cu)	<10 ppt	Tellurium (Te)	<1 ppt
Erbium (Er)	<1 ppt	Thallium (Tl)	<10 ppt
Europium (Eu)	<1 ppt	Thorium (Th)	<1 ppt
Gadolinium (Gd)	<1 ppt	Thulium (Tm)	<1 ppt
Germanium (Ge)	<10 ppt	Tin (Sn)	<20 ppt
Gold (Au)	<20 ppt	Titanium (Ti)	<10 ppt
Hafnium (Hf)	<10 ppt	Uranium (U)	<1 ppt
Holmium (Ho)	<1 ppt	Vanadium (V)	<10 ppt
Iron (Fe)	<10 ppt	Ytterbium (Yb)	<1 ppt
Lanthanum (La)	<1 ppt	Yttrium (Y)	<1 ppt
Lead (Pb)	<10 ppt	Zinc (Zn)	<10 ppt
Lithium (Li)	<10 ppt	Platinum (Pt)	<20 ppt
Lutetium (Lu)	<1 ppt	Rubidium (Rb)	<10 ppt
Magnesium (Mg)	<10 ppt	Tungsten (W)	<10 ppt

*Bottled in class 10 cleanroom conditions.Individually double-bagged in class 100 cleanroom conditions. Each bottle is individually packaged.
Certificate of analysis delivered with each bottle.*

Safe, environmentally friendly and effective packaging

Fisher Chemical outer packaging is UN compliant and 100% recyclable. It's an easy and safe way to reduce your waste management costs and contribute to a cleaner environment.

Nitric acid 67-69% for trace metal analysis		Trace Metal
packaging		Mfr. No
500 mℓ	HDPE plastic bottle	A509-P500
1 ℓ	HDPE plastic bottle	A509-P1
2.5 ℓ	HDPE plastic bottle	A509-P212
HNO ₃	EINECS: 231-714-2	
CAS: 7697-37-2	UN: 2031	
MW: 63.01	GHS: Ox. Liq. 3, Skin Corr. 1A	
Density at 20°C: 1.42	H314, H272	
Boiling point: 122°C	P301+P330+P331,	
Melting point: -42°C	P305+P351+P338, P310,	
ADR/RID: 8, 5.1, II	P304+P340, P280, P210	
Assay: 67-69%		

Guaranteed analysis			
Colour	<10 APHA	Nickel (Ni)	<0.5 ppb
Aluminium (Al)	<1 ppb	Niobium (Nb)	<0.1 ppb
Antimony (Sb)	<0.5 ppb	Palladium (Pd)	<0.5 ppb
Arsenic (As)	<0.5 ppb	Platinum (Pt)	<0.5 ppb
Barium (Ba)	<0.1 ppb	Potassium (K)	<1 ppb
Beryllium (Be)	<0.1 ppb	Praseodymium (Pr)	<0.1 ppb
Bismuth (Bi)	<0.1 ppb	Rhenium (Re)	<0.1 ppb
Boron (B)	<1 ppb	Rhodium (Rh)	<0.5 ppb
Cadmium (Cd)	<0.5 ppb	Rubidium (Rb)	<0.1 ppb
Calcium (Ca)	<1 ppb	Ruthenium (Ru)	<0.5 ppb
Cerium (Ce)	<0.1 ppb	Samarium (Sm)	<0.1 ppb
Cesium (Cs)	<0.1 ppb	Scandium (Sc)	<0.1 ppb
Chromium (Cr)	<1 ppb	Selenium (Se)	<1 ppb
Cobalt (Co)	<0.5 ppb	Silver (Ag)	<0.1 ppb
Copper (Cu)	<0.5 ppb	Sodium (Na)	<1 ppb
Dysprosium (Dy)	<0.1 ppb	Strontium (Sr)	<0.1 ppb
Erbium (Er)	<0.1 ppb	Tantalum (Ta)	<0.1 ppb
Europium (Eu)	<0.1 ppb	Tellurium (Te)	<0.1 ppb
Gadolinium (Gd)	<0.1 ppb	Terbium (Tb)	<0.1 ppb
Gallium (Ga)	<0.1 ppb	Thallium (Tl)	<0.1 ppb
Germanium (Ge)	<0.1 ppb	Thorium (Th)	<0.1 ppb
Gold (Au)	<0.1 ppb	Thulium (Tm)	<0.1 ppb
Hafnium (Hf)	<0.1 ppb	Tin (Sn)	<0.5 ppb
Holmium (Ho)	<0.1 ppb	Titanium (Ti)	<0.5 ppb
Indium (In)	<0.1 ppb	Total chloride (Cl)	<0.2 ppm
Iron (Fe)	<1 ppb	Total phosphorus (P)	<0.01 ppm
Lanthanum (La)	<0.1 ppb	Total sulphur (S)	<0.3 ppm
Lead (Pb)	<0.1 ppb	Tungsten (W)	<0.1 ppb
Lithium (Li)	<0.1 ppb	Uranium (U)	<0.1 ppb
Lutetium (Lu)	<0.1 ppb	Vanadium (V)	<0.5 ppb
Magnesium (Mg)	<1 ppb	Ytterbium (Yb)	<0.1 ppb
Manganese (Mn)	<0.1 ppb	Yttrium (Y)	<0.1 ppb
Mercury (Hg)	<0.1 ppb	Zinc (Zn)	<0.5 ppb
Molybdenum (Mo)	<0.1 ppb	Zirconium (Zr)	<0.1 ppb
Neodymium (Nd)	<0.1 ppb		

Bottled in class 10 cleanroom conditions.
Lot analysis on the label.

Nitric acid 68% d =1.42		NEW PACK.	Primar Plus
for trace metal analysis			
packaging			Mfr. No
500 mℓ	HDPE Plastic bottle		N/2272/PB08
1 ℓ	HDPE Plastic bottle		N/2272/PB15
2.5 ℓ	HDPE Plastic bottle		N/2272/PB17
HNO ₃	EINECS: 231-714-2		
CAS: 7697-37-2	UN: 2031		
MW: 63.01	GHS: Ox. Liq. 3, Skin Corr. 1A		
Density at 20°C: 1.42	H314, H272		
Boiling point: 122°C	P301+P330+P331,		
Melting point: -42°C	P305+P351+P338, P310,		
ADR/RID: 8, 5.1, II	P304+P340, P280, P210		
Assay: >68%			
Guaranteed analysis			
Aluminium (Al)	<0.1 ppm	Calcium (Ca)	<0.05 ppm
Antimony (Sb)	<0.005 ppm	Chromium (Cr)	<0.005 ppm
Arsenic (As)	<0.005 ppm	Copper (Cu)	<0.005 ppm
Barium (Ba)	<0.05 ppm	Gallium (Ga)	<0.005 ppm
Beryllium (Be)	<0.005 ppm	Germanium (Ge)	<0.005 ppm
Bismuth (Bi)	<0.005 ppm	Indium (In)	<0.005 ppm
Boron (B)	<0.005 ppm	Iron (Fe)	<0.05 ppm
Cadmium (Cd)	<0.002 ppm		

Lead (Pb)	<0.002 ppm	Strontium (Sr)	<0.002 ppm
Lithium (Li)	<0.002 ppm	Thallium (Tl)	<0.005 ppm
Magnesium (Mg)	<0.02 ppm	Tin (Sn)	<0.005 ppm
Manganese (Mn)	<0.002 ppm	Titanium (Ti)	<0.002 ppm
Mercury (Hg)	<0.005 ppm	Total chloride (Cl)	<0.00002%
Molybdenum (Mo)	<0.002 ppm	Total phosphorus (P)	<0.05 ppm
Nickel (Ni)	<0.002 ppm	Total sulphur (S)	<0.2 ppm
Potassium (K)	<0.02 ppm	Vanadium (V)	<0.002 ppm
Selenium (Se)	<0.002 ppm	Zinc (Zn)	<0.01 ppm
Silver (Ag)	<0.005 ppm	Zirconium (Zr)	<0.005 ppm
Sodium (Na)	<0.1 ppm		

Nitric acid 70% d=1.42
for electronic use (MOS)

packaging		Mfr. No
2.5 ℓ	HDPE Plastic bottle	N/2222/PB17
HNO ₃	EINECS: 231-714-2	
CAS: 7697-37-2	UN: 2031	
MW: 63.01	GHS: Ox. Liq. 3, Skin Corr. 1A	
Density at 20°C: 1.42	H272, H314	
Boiling point: 122°C	P301+P330+P331, P280,	
Melting point: -42°C	P305+P351+P338, P310,	
ADR/RID: 8, 5.1, II	P221, P210	
Assay: 69 to 71%		
Guaranteed analysis		
Colour	<10 APHA	Molybdenum (Mo) <0.02 ppm
Aluminium (Al)	<0.1 ppm	Nickel (Ni) <0.05 ppm
Barium (Ba)	<0.02 ppm	Potassium (K) <0.2 ppm
Cadmium (Cd)	<0.02 ppm	Sodium (Na) <2 ppm
Calcium (Ca)	<1 ppm	Strontium (Sr) <0.02 ppm
Chromium (Cr)	<0.05 ppm	Titanium (Ti) <0.02 ppm
Cobalt (Co)	<0.02 ppm	n-Butyl formate <0.00005%
Copper (Cu)	<0.1 ppm	Total phosphorus (P) <0.5 ppm
Iron (Fe)	<0.2 ppm	Total silicon (Si) <0.5 ppm
Lead (Pb)	<0.1 ppm	Total sulphur (S) <1 ppm
Lithium (Li)	<0.02 ppm	Vanadium (V) <0.02 ppm
Magnesium (Mg)	<0.2 ppm	Zinc (Zn) <0.1 ppm
Manganese (Mn)	<0.02 ppm	

The concentration of nitric acid solution is 15.8 M (15.8 N).
Filtered to 10 µm.

Nitric acid 69%

for analysis with low Cadmium, Lead & Mercury level

Packaging		Mfr. No	
1 ℓ	HDPE Plastic bottle	N/2320/PB15	
2.5 ℓ	HDPE Plastic bottle	N/2320/PB17	
HNO ₃			
CAS: 7697-37-2			
MW: 63.01			
Density at 20°C: 1.42			
Boiling point: 122°C			
Melting point: -42°C			
ADR/RID: 8, 5.1, II			
Assay: 68-71%			
Guaranteed analysis			
Colour	<10 APHA	Lithium (Li)	<0.01 ppm
Zirconium (Zr)	<0.02 ppm	Magnesium (Mg)	<0.05 ppm
Aluminium (Al)	<0.05 ppm	Manganese (Mn)	<0.01 ppm
Arsenic (As)	<0.01 ppm	Mercury (Hg)	<0.005 ppm
Barium (Ba)	<0.01 ppm	Molybdenum (Mo)	<0.01 ppm
Beryllium (Be)	<0.01 ppm	Nickel (Ni)	<0.02 ppm
Bismuth (Bi)	<0.05 ppm	Phosphate (PO ₄)	<0.5 ppm
Cadmium (Cd)	<0.005 ppm	Potassium (K)	<0.1 ppm
Calcium (Ca)	<0.1 ppm	Strontium (Sr)	<0.01 ppm
Chloride (Cl)	<0.2 ppm	Sulfate (SO ₄)	<1 ppm
Chromium (Cr)	<0.02 ppm	Thallium (Tl)	<0.02 ppm
Cobalt (Co)	<0.01 ppm	Titanium (Ti)	<0.02 ppm
Copper (Cu)	<0.02 ppm	Vanadium (V)	<0.01 ppm
Germanium (Ge)	<0.02 ppm	Zinc (Zn)	<0.02 ppm
Iron (Fe)	<0.2 ppm	Residue after ignition	<5 ppm
Lead (Pb)	<0.01 ppm		

Lot analysis on the label

Nitric acid 68 % d= 1.42 for analysis		Certified AR
packaging		Mfr. No
500 mℓ	Coated Glass bottle	N/2300/PB08
1 ℓ	HDPE Plastic bottle	N/2300/PB15
2.5 ℓ	HDPE Plastic bottle	N/2300/PB17
2.5 ℓ	Coated Glass bottle	N/2300/PC17
HNO ₃	EINECS: 231-714-2	
CAS: 7697-37-2	UN: 2031	
MW: 63.01	GHS: Ox. Liq. 3, Skin Corr. 1A	
Density at 20°C: 1.42	H272, H314	
Boiling point: 122°C	P301+P330+P331, P280,	
Melting point: -42°C	P305+P351+P338, P310,	
ADR/RID: 8, 5.1, II	P221, P210	
Assay: 68 to 70%		

Guaranteed analysis			
Colour	<10 APHA	Nickel (Ni)	<0.02 ppm
Residue after ignition	<0.0005%	Potassium (K)	<0.1 ppm
Calcium (Ca)	<0.5 ppm	Residue after ignition	<0.0005%
Chromium (Cr)	<0.05 ppm	Total chloride (Cl)	<0.00005%
Copper (Cu)	<0.05 ppm	Total phosphorus (P)	<0.2 ppm
Iron (Fe)	<0.2 ppm	Total sulphur (S)	<1 ppm
Lead (Pb)	<0.05 ppm	Zinc (Zn)	<0.05 ppm
Magnesium (Mg)	<0.2 ppm	Sulfate (SO ₄)	<1 ppm
Mercury (Hg)	<0.05 ppm		

Typical analysis			
Aluminium (Al)	<0.05 ppm	Silver (Ag)	<0.01 ppm
Arsenic (As)	<0.01 ppm	Sodium (Na)	<0.5 ppm
Barium (Ba)	<0.02 ppm	Strontium (Sr)	<0.01 ppm
Cadmium (Cd)	<0.01 ppm	Tin (Sn)	<0.01 ppm
Cobalt (Co)	<0.01 ppm	Titanium (Ti)	<0.01 ppm
Lithium (Li)	<0.01 ppm	Total silicon (Si)	<0.5 ppm
Manganese (Mn)	<0.02 ppm	Vanadium (V)	<0.01 ppm
Molybdenum (Mo)	<0.01 ppm		

Lot analysis on the label.
The concentration of nitric acid solution is 15.8 M (15.8 N).

Nitric acid 70% d=1.42
specified reagents for general laboratory work- extra pure

packaging		Mfr. No	
1 ℓ	HDPE Plastic bottle	N/2250/PB15	
2.5 ℓ	HDPE Plastic bottle	N/2250/PB17	
25 ℓ	Plastic container	N/2250/25	
HNO ₃	EINECS: 231-714-2		
CAS: 7697-37-2	UN: 2031		
MW: 63.01	GHS: Ox. Liq. 3, Skin Corr. 1A		
Density at 20°C: 1.42	H272, H314		
Boiling point: 122°C	P301+P330+P331, P280,		
Melting point: -42°C	P305+P351+P338, P310,		
ADR/RID: 8, 5.1, II	P221, P210		
Assay: 68 to 72%			
Guaranteed analysis			
Calcium (Ca)	<20 ppm	Potassium (K)	<10 ppm
Chromium (Cr)	<2 ppm	Sodium (Na)	<50 ppm
Copper (Cu)	<2 ppm	n-Butyl formate	<0.005%
Iron (Fe)	<10 ppm	Total phosphorus (P)	<10 ppm
Lead (Pb)	<2 ppm	Total silicon (Si)	<10 ppm
Magnesium (Mg)	<5 ppm	Total sulphur (S)	<20 ppm
Nickel (Ni)	<2 ppm	Zinc (Zn)	<2 ppm

The concentration of nitric acid solution is 15.8 M (15.8 N).

Nitric acid 70% d=1.42
Technical

packaging		Mfr. No
2.5 ℓ	HDPE Plastic bottle	N/2200/PB17
HNO ₃	EINECS: 231-714-2	
CAS: 7697-37-2	UN: 2031	
MW: 63.01	GHS: Ox. Liq. 3, Skin Corr. 1A	
Density at 20°C: 1.42	H272, H314	
Boiling point: 122°C	P301+P330+P331, P280,	
Melting point: -42°C	P305+P351+P338, P310,	
ADR/RID: 8, 5.1, II	P221, P210	
Assay: >70%		

Nitric Acid 65% d = 1.40 for analysis		Certified AR
packaging		Mfr. No
1 ℓ	HDPE Plastic bottle	N/2185/PB15
2.5 ℓ	HDPE Plastic bottle	N/2185/PB17
HNO ₃	EINECS: 231-714-2	
CAS: 7697-37-2	UN: 2031	
MW: 63.01	GHS: Ox. Liq. 3, Skin Corr. 1A	
Density at 20°C: 1.40	H314, H272	
Boiling point: 122°C	P301+P330+P331,	
Melting point: -41°C	P305+P351+P338, P310,	
ADR/RID: 8, 5.1, II	P304+P340, P280, P210	
Assay: <66%		

Guaranteed analysis			
Colour	<10 APHA	Nickel (Ni)	<0.05 ppm
Calcium (Ca)	<1 ppm	Potassium (K)	<0.2 ppm
Chromium (Cr)	<0.05 ppm	Sodium (Na)	<2 ppm
Copper (Cu)	<0.1 ppm	Total chloride (Cl)	<0.00005%
Iron (Fe)	<0.2 ppm	Total phosphorus (P)	<0.5 ppm
Lead (Pb)	<0.1 ppm	Total sulphur (S)	<1 ppm
Magnesium (Mg)	<0.2 ppm	Zinc (Zn)	<0.1 ppm
Mercury (Hg)	<0.05 ppm		

Typical analysis			
Arsenic (As)	<0.02 ppm	Silver (Ag)	<0.02 ppm
Barium (Ba)	<0.02 ppm	Strontium (Sr)	<0.02 ppm
Cobalt (Co)	<0.02 ppm	Tin (Sn)	<0.02 ppm
Copper (Li)	<0.02 ppm	Titanium (Ti)	<0.02 ppm
Manganese (Mn)	<0.02 ppm	Total silicon (Si)	<0.5 ppm
Molybdenum (Mo)	<0.02 ppm	Vanadium (V)	<0.02 ppm

Lot analysis on the label.

Nitric Acid 65% d = 1.40
specified reagents for general laboratory work- extra pure

Packaging		Mfr. No	
1 ℓ	HDPE Plastic bottle	N/2180/PB15	
2.5 ℓ	HDPE Plastic bottle	N/2180/PB17	
HNO ₃	EINECS: 231-714-2		
CAS: 7697-37-2	UN: 2031		
MW: 63.01	GHS: Ox. Liq. 3, Skin Corr. 1A		
Density at 20°C: 1.40	H314, H272		
Boiling point: 122°C	P301+P330+P331,		
Melting point: -41°C	P305+P351+P338, P310,		
ADR/RID: 8, 5.1, II	P304+P340, P280, P210		
Assay: <66%			
Guaranteed analysis			
Colour	<10 APHA	Potassium (K)	<0.2 ppm
Calcium (Ca)	<1 ppm	Sodium (Na)	<2 ppm
Chromium (Cr)	<0.05 ppm	n-Butyl formate	<0.00005%
Copper (Cu)	<0.1 ppm	Total phosphorus (P)	<0.5 ppm
Iron (Fe)	<0.2 ppm	Total silicon (Si)	<0.5 ppm
Lead (Pb)	<0.1 ppm	Total sulphur (S)	<1 ppm
Magnesium (Mg)	<0.2 ppm	Zinc (Zn)	<0.1 ppm
Nickel (Ni)	<0.05 ppm		

Nitric Acid 53% d = 1.33
for analysis

Certified AR

packaging		Mfr. No
1 ℓ	HDPE Plastic bottle	N/2175/PB15
HNO ₃	UN: 2031	
CAS: 7697-37-2	GHS: Skin Corr. 1A, Eye Dam. 1	
MW: 63.01	H314	
Density at 20°C: 1.33	P280, P305+P351+P338,	
Boiling point: 122°C	P303+P361+P353,	
Melting point: -42°C	P301+P330+P331, P310,	
ADR/RID: 8, II	P304+P340	
EINECS: 231-714-2		
Assay: 52.5 to 55%		

Guaranteed analysis		
Colour	<10 APHA	Mercury (Hg) <0.05 ppm
Density at 20°C	: 1.325 to 1.335	Nickel (Ni) <0.05 ppm
Arsenic (As)	<0.05 ppm	Sodium (Na) <2 ppm
Calcium (Ca)	<1 ppm	Total chloride (Cl) <0.00005%
Chromium (Cr)	<0.05 ppm	Total phosphorus (P) <0.5 ppm
Copper (Cu)	<0.1 ppm	Total silicon (Si) <0.5 ppm
Lead (Pb)	<0.1 ppm	Total sulphur (S) <1 ppm
Magnesium (Mg)	<0.2 ppm	Zinc (Zn) <0.1 ppm

Typical analysis		
Aluminium (Al)	<0.1 ppm	Molybdenum (Mo) <0.02 ppm
Barium (Ba)	<0.02 ppm	Potassium (K) <0.2 ppm
Cadmium (Cd)	<0.02 ppm	Silver (Ag) <0.02 ppm
Cobalt (Co)	<0.02 ppm	Strontium (Sr) <0.02 ppm
Iron (Fe)	<0.2 ppm	Tin (Sn) <0.02 ppm
Lithium (Li)	<0.02 ppm	Titanium (Ti) <0.02 ppm
Manganese (Mn)	<0.02 ppm	Vanadium (V) <0.02 ppm

Lot analysis on the label.

Nitric Acid

packaging		Mfr. No
5 ℓ	HDPE Plastic bottle	N/2170/PB21
HNO ₃	UN: 2031	
CAS: 7697-37-2	GHS: Skin Corr. 1A, Eye Dam. 1	
MW: 63.01	H314	
Density at 20°C: 1.33	P280, P305+P351+P338,	
Boiling point: 122°C	P303+P361+P353,	
Melting point: -42°C	P301+P330+P331, P310,	
ADR/RID: 8, II	P304+P340	
EINECS: 231-714-2		
Assay: 52.5 to 55%		

Nitric acid solution 2M (2N)
NIST Standard solution ready to use, For Volumetric analysis

packaging		Mfr. No
2.5 ℓ	HDPE Plastic bottle	J/5545/PB17
CAS: 7697-37-2	UN: 2031	
Density at 20°C: 1.03 to 1.12	GHS: Skin Corr. 1B	
Boiling point: 100°C	H314	
ADR/RID: 8, II	P280, P305+P351+P338,	
EINECS: 231-714-2	P301+P330+P331, P310	

Manufactured to a tolerance of ± 0.5%.

Nitric acid solution 1M (1N)
NIST Standard solution ready to use, For Volumetric analysis

Eur.Ph. , USP, BP

packaging		Mfr. No
1 ℓ	HDPE Plastic bottle	J/5550/PB15
2.5 ℓ	HDPE Plastic bottle	J/5550/PB17
10 ℓ	Bag in box	J/5550/24
CAS: 7697-37-2	UN: 2031	
Density at 20°C: 1.03 to 1.12	GHS: Skin Corr. 1B	
Boiling point: 100°C	H314	
ADR/RID: 8, II	P280, P305+P351+P338,	
EINECS: 231-714-2	P301+P330+P331, P310	

Manufactured to a tolerance of ± 0.1%.
See the full range of Standard Volumetric solution page 323.

Nitric acid solution 1M (1N)
NIST Standard Concentrate, For Volumetric analysis

packaging		Mfr. No
1	Ampoule	J/5550C/05
6	Ampoule	J/5550C/90
HNO ₃	UN: 2031	
CAS: 7697-37-2	GHS: Skin Corr. 1A, Eye Dam. 1	
MW: 63.01	H314	
Density at 20°C: 1.33	P280, P305+P351+P338,	
Boiling point: 122°C	P303+P361+P353,	
Melting point: -42°C	P301+P330+P331, P310,	
ADR/RID: 8, II	P304+P340	
EINECS: 231-714-2		

Each ampoule makes 1 ℓ of solution.
See the full range of Standard Volumetric solution page 323.

Nitric acid solution 0.1M (0.1N)
NIST Standard solution ready to use, For Volumetric analysis

packaging		Mfr. No
1 ℓ	HDPE Plastic bottle	J/5600/PB15
2.5 ℓ	HDPE Plastic bottle	J/5600/17
10 ℓ	Bag in box	J/5600/24
CAS: 7697-37-2		
Density at 20°C: 1		
EINECS: 231-714-2		

Manufactured to a tolerance of ± 0.1%.
See the full range of Standard Volumetric solution page 323.

2.2',2''-Nitrilotriethanol See Triethanolamine, page 357

p-Nitroaniline See 4-Nitroaniline, page 218

4-Nitroaniline
pure

packaging		Mfr. No
100 g	Glass jar	N/2600/48
C ₆ H ₅ N ₂ O ₂	UN: 1661	
CAS: 100-01-6	GHS: Acute Tox. 3, STOT RE	
MW: 138.13	2, Aquatic Chronic 3	
Boiling point: 332°C	H412, H311, H373, H331, H301	
Melting point: 146 to 151°C	P301+P310, P310,	
Flash point: 199°C	P304+P340, P280,	
ADR/RID: 6.1, II	P302+P350, P260, P273	
EINECS: 202-810-1		

Nitromethane
specified reagents for general laboratory work- extra pure

SLR

packaging		Mfr. No
100 mℓ	Coated Glass bottle	N/3000/PB05
500 mℓ	Coated Glass bottle	N/3000/PB08
200 ℓ	Metal drum	N/3000/27
CH ₃ NO ₂	ADR/RID: 3, II	
CAS: 75-52-5	EINECS: 200-876-6	
MW: 61.04	UN: 1261	
Density at 20°C: 1.12	GHS: Flam. Liq. 3, Acute Tox. 4	
Boiling point: 100 to 102°C	H226, H302	
Melting point: -29°C	P301+P312, P304+P340, P210	
Flash point: 35°C		

3-Nitrophenol
specified reagents for general laboratory work- extra pure

SLR

packaging		Mfr. No
25 g	Plastic jar	N/3150/46
C ₆ H ₅ NO ₃	GHS: Acute Tox. 4, Skin Irrit.	
CAS: 554-84-7	2, Eye Irrit. 2, STOT SE 3,	
MW: 139.11	STOT RE 2	
Boiling point: 194°C	H373, H319, H315, H335,	
Melting point: 94 to 98°C	H302, H312, H332	
ADR/RID: 6.1, III	P280, P302+P352, P261,	
EINECS: 209-073-5	P260, P301+P312,	
UN: 1663	P304+P340, P305+P351+P338	
Guaranteed analysis		
Melting Point tested 96 to 98°C		

4-Nitrophenol
extra pure, pH indicator

SLR

packaging		Mfr. No
100 g	Glass jar	N/3200/48
C ₆ H ₅ NO ₃	EINECS: 202-811-7	
CAS: 100-02-7	UN: 1663	
MW: 139.11	GHS: Acute Tox. 4, STOT RE 2	
Boiling point: 279°C	H373, H332, H312, H302	
Melting point: >112°C	P280, P304+P340,	
Flash point: 169°C	P302+P352, P260	
ADR/RID: 6.1, III		
Guaranteed analysis		
Melting Point tested 111 to 114°C		
For the determination of pH. pH 5.6-7.6 colorless to yellow.		

p-Nitrophenylamine See 4-Nitroaniline, page 218

4-Nitrophenyl disodium orthophosphate
pure

Packaging		Mfr. No
5 g	Plastic jar	N/3390/44
25 g	Plastic jar	N/3390/46
C ₆ H ₄ NNa ₂ PO ₆ ·6H ₂ O	Melting point: >300°C	
CAS: 4264-83-9	EINECS: 224-246-5	
MW: 371.15		
Guaranteed analysis		
Complies with standard IR spectrum: Pass		

Nitroso-R-salt
extra pure, determination of cobalt

SLR

Packaging		Mfr. No
25 g	Plastic jar	N/3650/46
$C_{10}H_5NNa_2O_8S_2 \cdot xH_2O$	2, Eye Irrit. 2, STOT SE 3	
CAS: 525-05-3	H315, H302, H332, H319, H312, H335	
MW: 377.26	P280, P302+P352, P261, P301+P312, P304+P340, P305+P351+P338	
Melting point: 300°C		
EINECS: 208-369-1		
GHS: Acute Tox. 4, Skin Irrit.		
Assay: >90%		

UHPLC

Filtered to 0.1 μm, Fisher Chemical UHPLC gradient grade solvents (Acetonitrile, Methanol and Water) are perfect for those working at high pressure.

See page 309 for more information



n-Nonane
pure

packaging		Mfr. No
100 mℓ	Glass bottle	N/3900/05
C ₉ H ₂₀	UN: 1920	
CAS: 111-84-2	GHS: Flam. Liq. 3, Acute	
MW: 128.26	Tox. 4, Skin Irrit. 2, STOT SE	
Density at 20°C: 0.718	3, Asp. Tox. 1, Aquatic Chronic 4	
Boiling point: 151°C	H226, H332, H336, H304,	
Melting point: -53°C	H315, H413	
Flash point: 31°C	P261, P301+P310, P331,	
ADR/RID: 3, III	P302+P352, P210, P304+P340	
EINECS: 203-913-4		
Assay (CG): >99%		

1-Nonanol
specified reagents for general laboratory work- extra pure

SLR

Packaging		Mfr. No
500 ml	Glass bottle	N/4020/08
C ₉ H ₂₀ O	EINECS: 205-583-7	
CAS: 143-08-8	UN: 3082	
MW: 144.26	GHS: Acute Tox. 4, Skin Irrit.	
Density at 20°C: 0.827	2, Eye Irrit. 2, Aquatic Chronic 2	
Boiling point: 213 to 215°C	H319, H411, H315, H332	
Melting point: -5--8°C	P273, P304+P340,	
Flash point: 74°C	P302+P352, P280,	
ADR/RID: 9, III	P305+P351+P338	
Assay (CG): >98%		

Nonyl alcohol See 1-Nonanol, page 219

n-Octane
specified reagents for general laboratory work- extra pure

SLR

packaging		Mfr. No
100 mℓ	Glass bottle	O/0020/05
500 mℓ	Glass bottle	O/0020/08
C ₈ H ₁₈	UN: 1262	
CAS: 111-65-9	GHS: Flam. Liq. 2, Skin Irrit.	
MW: 114.23	2, STOT SE 3, Asp. Tox. 1,	
Density at 20°C: 0.708	Aquatic Acute 1, Aquatic	
Boiling point: 125 to 127°C	Chronic 1	
Melting point: -57°C	H225, H315, H410, H336, H304	
Flash point: 13°C	P261, P301+P310, P331,	
ADR/RID: 3, II	P273, P302+P352, P210, P240	
EINECS: 203-892-1		
Assay (CG): >99%		

High Purity Acids for Elemental Analysis

Fisher Chemical offers a complete range of acids to meet the most challenging applications. Choose from our Optima®, TraceMetal® and PrimarPlus® grades to ensure that the sensitivity you require is always available.

Read more page 24



1-Octane sulfonic acid sodium salt
for Ion pair chromatography

packaging		Mfr. No
5 g	Glass jar	O/0028/44
25 g	Plastic jar	O/0028/46
100 g	Plastic jar	O/0028/48
C ₈ H ₁₇ NaSO ₃ CAS: 5324-84-5 MW: 216.28 Melting point: >300°C EINECS: 226-195-4		GHS: Skin Irrit. 2, Eye Irrit. 2, STOT SE 3 H335, H315, H319 P261, P302+P352, P280, P305+P351+P338
Assay: >98%		
Guaranteed analysis		
Solution (5% in methanol)		Pass
Absorbance specification		
Max. Absorbance (1cm cell v HPLC water)		
200 nm (sol. aqu. 0.25 M)		<0.1 A.U.
210 nm (sol. aqu. 0.25 M)		<0.05 A.U.
220 nm (sol. aqu. 0.25 M)		<0.03 A.U.
240 nm (sol. aqu. 0.25 M) ...		<0.02 A.U.
250 nm (sol. aqu. 0.25 M) ...		<0.02 A.U.
260 nm (sol. aqu. 0.25 M) ...		<0.02 A.U.
Lot analysis on the label.		

1-Octane sulfonic acid sodium salt
for Ion pair chromatography in electrochemical detection

packaging		Mfr. No
25 g	Plastic jar	O/0029/46
C ₈ H ₁₇ NaSO ₃ CAS: 5324-84-5 MW: 216.28 Melting point: >300°C EINECS: 226-195-4		GHS: Skin Irrit. 2, Eye Irrit. 2, STOT SE 3 H335, H315, H319 P261, P302+P352, P280, P305+P351+P338
Assay: >98%		
Guaranteed analysis		
Electrochemical test		Pass
Absorbance specification		
Max. Absorbance (1cm cell v HPLC water)		
200 nm (sol. aqu. 0.25 M)		<0.1 A.U.
210 nm (sol. aqu. 0.25 M) ...		<0.05 A.U.
220 nm (sol. aqu. 0.25 M) ...		<0.02 A.U.
240 nm (sol. aqu. 0.25 M) ...		<0.02 A.U.
250 nm (sol. aqu. 0.25 M) ...		<0.02 A.U.
260 nm (sol. aqu. 0.25 M) ...		<0.02 A.U.
Tested by HPLC-EC to ensure low electrochemical response.		
Detector - ESA Coulochem II		
Mode-oxidative		
Graphite Electrode.		
Lot analysis on the label.		

1-Octanesulfonic acid sodium salt solution 0.1M
for Ion pair chromatography

packaging		Mfr. No
250 mL	Glass bottle	O/0026/07
Absorbance specification		
Max. Absorbance (1cm cell v HPLC water)		
200 nm (sol. aqu. 0.005 M) :		<0.04 A.U.
210 nm (sol. aqu. 0.005 M) :		<0.02 A.U.
220 nm (sol. aqu. 0.005 M):		<0.016 A.U.
230 nm (sol. aqu. 0.005 M):		<0.012 A.U.
240 nm (sol. aqu. 0.005 M):		<0.008 A.U.
250 nm (sol. aqu. 0.005 M):		<0.008 A.U.
260 nm (sol. aqu. 0.005 M):		<0.008 A.U.
Lot analysis on the label.		

We Support Your Chemistry For

• Semi-bulk and bulk requirements

• Tailored solvents and solvent blends or customized solutions

• Testing services or customized packaging and labeling

Please email your contact details and a brief description of your requirements to begei.sha.quotes@thermofisher.com

n-Octanoic acid
specified reagents for general laboratory work- extra pure

packaging		Mfr. No
500 mL	Glass bottle	O/0030/08
2.5 L	HDPE Plastic bottle	O/0030/PB17
25 L	HDPE Plastic container	O/0030/25
C ₈ H ₁₆ O ₂ CAS: 124-07-2 MW: 144.21 Density at 20°C: 0.910 Boiling point: 237°C Melting point: 16 to 16.5°C Flash point: 130°C		ADR/RID: 8, III EINECS: 204-677-5 UN: 3265 GHS: Skin Corr. 1B H314 P301+P330+P331, P280, P305+P351+P338, P310
Assay: >98%		

Octan-1-ol
pure

packaging		Mfr. No
500 mL	Glass bottle	O/0050/08
2.5 L	Glass bottle	O/0050/17
2.5 L	Coated Glass bottle	O/0050/PB17
C ₈ H ₁₈ O CAS: 111-87-5 MW: 130.23 Density at 20°C: 0.824 Boiling point: 195°C Melting point: -16°C		Flash point: 81°C EINECS: 203-917-6 GHS: Eye Irrit. 2 H319 P280, P305+P351+P338
Assay (CG): >98%		

Octan-2-ol
pure

packaging		Mfr. No
500 mL	Glass bottle	O/0100/08
2.5 L	Glass bottle	O/0100/17
C ₈ H ₁₈ O CAS: 123-96-6 MW: 130.23 Density at 20°C: 0.810 Boiling point: 174 to 181°C Melting point: -38°C Flash point: 72°C ADR/RID: 6.1, III		EINECS: 204-667-0 UN: 2810 GHS: Acute Tox. 3, Skin Irrit. 2, Eye Irrit. 2, STOT SE 3 H315, H319, H335, H301 P261, P302+P352, P280, P305+P351+P338
Assay (CG): >98%		

Oct-1-ene
pure

packaging		Mfr. No
100 mL	Coated Glass bottle	O/0105/PB05
C ₈ H ₁₆ CAS: 111-66-0 MW: 112.21 Density at 20°C: 0.715 Boiling point: 121 to 123°C Melting point: -102°C Flash point: 13°C ADR/RID: 3, II EINECS: 203-893-7 UN: 1993		GHS: Flam. Liq. 2, Skin Irrit. 2, Eye Irrit. 2, STOT SE 3, Asp. Tox. 1, Aquatic Chronic 2 H225, H335, H319, H315, H304, H411 P261, P301+P310, P331, P273, P302+P352, P280, P305+P351+P338, P210, P240, P304+P340

Assay (CG): >95%

Oleic acid
specified reagents for general laboratory work- extra pure

packaging		Mfr. No
1 L	Glass bottle	O/0200/15
2.5 L	Glass bottle	O/0200/17
C ₁₈ H ₃₄ O ₂ CAS: 112-80-1 MW: 282.47 Density at 20°C: 0.890 Boiling point: 360°C Melting point: 13°C Flash point: 189°C		EINECS: 204-007-1 GHS: Skin Irrit. 2, Eye Irrit. 2, STOT SE 3 H335, H319, H315 P261, P280, P302+P352, P305+P351+P338
Guaranteed analysis		
Residue after ignition		<0.1% Iodine value
Acid value (Mg KOH/g) . .		197 to 202 mg 85 to 92

Oleic acid
pure

packaging		Mfr. No
2.5 L	Glass bottle	O/0150/17
25 L	HDPE Plastic container	O/0150/25
200 L	Metal drum	O/0150/27
C ₁₈ H ₃₄ O ₂ CAS: 112-80-1 MW: 282.46 Density at 20°C: 0.890 Boiling point: 360°C Melting point: 13°C Flash point: 189°C		EINECS: 204-007-1 GHS: Skin Irrit. 2, Eye Irrit. 2, STOT SE 3 H335, H319, H315 P261, P280, P302+P352, P305+P351+P338
Guaranteed analysis		
Complies with standard IR spectrum:		Pass

Oleum See Sulfuric acid, page 339

Olive oil
pure

packaging		Mfr. No
500 mL	Glass bottle	O/0250/08
2.5 L	Glass bottle	O/0250/17
CAS: 8001-25-0 Density at 20°C: 0.914 Melting point: 0°C		Flash point: 225°C EINECS: 232-277-0
Guaranteed analysis		
Identity test		: Pass

Orange G sodium salt
pure

packaging		Mfr. No
25 g	Glass jar	O/0275/46
C ₁₆ H ₁₀ N ₂ Na ₂ O ₇ S ₂ CAS: 1936-15-8 MW: 452.37 Melting point: 141°C EINECS: 217-705-6		GHS: Skin Irrit. 2, Eye Irrit. 2, STOT SE 3 H335, H315, H319 P261, P302+P352, P280, P305+P351+P338
Guaranteed analysis		
Complies with standard IR spectrum:		Pass

Fisher Chemical 5L packaging change

Our improved HDPE 5L packaging is easier to handle, more robust and preserves product quality for longer.

Orcein staining
pure

packaging		Mfr. No
5 g	Glass jar	O/0300/44
C ₂₈ H ₂₄ N ₂ O ₇ CAS: 1400-62-0 MW: 500.51 ADR/RID: 6.1, III EINECS: 215-750-6 UN: 3143		GHS: Acute Tox. 4, Skin Irrit. 2, Eye Irrit. 2 H302, H319, H315 P301+P312, P302+P352, P280, P305+P351+P338
Guaranteed analysis		
Complies with standard IR spectrum:		Pass

Orcinol monohydrate
pure

packaging		Mfr. No
25 g	Plastic jar	O/0351/46
C ₇ H ₈ O ₂ ·H ₂ O CAS: 6153-39-5 MW: 142.15 Boiling point: 290°C Melting point: 56 to 61°C EINECS: 207-984-2		GHS: Acute Tox. 4, Skin Irrit. 2, Eye Irrit. 2, STOT SE 3 H319, H302, H315, H335 P261, P301+P312, P302+P352, P280, P305+P351+P338
Guaranteed analysis		
Complies with standard IR spectrum:		Pass

Orthophosphoric acid 85%
HPLC for electrochemical detection

packaging		Mfr. No
500 mL	Coated Glass bottle	O/0515/PB08
H ₃ PO ₄ CAS: 7664-38-2 MW: 98.00 Density at 20°C: 1.68 Boiling point: 158°C Melting point: 21°C ADR/RID: 8, III		EINECS: 231-633-2 UN: 1805 GHS: Skin Corr. 1B, Eye Dam. 1 H314 P280, P301+P330+P331, P310, P305+P351+P338, P303+P361+P353
Assay: >85%		
Guaranteed analysis		
Electrochemical test		Pass
Iron (Fe)		<10 ppm
Lead (Pb)		<5 ppm
n-Butyl formate		<0.001%
Absorbance specification		
Max. Absorbance (1cm cell v HPLC water)		
254 nm (sol. aqu. 1 M)		<0.02 A.U.
Tested by HPLC-EC to ensure low electrochemical response.		
Detector - ESA Coulochem II		
Mode-oxidative		
Graphite Electrode.		
Orthophosphoric acid 85% is equivalent to a strength of 15.2 M (45.6 N).		
Lot analysis on the label.		
Filtered to 0.2 µm.		

Orthophosphoric acid solution 1M
for HPLC

packaging		Mfr. No
500 mL	Coated Glass bottle	J/5655/PB08
H ₃ PO ₄ CAS: 7664-38-2 MW: 98.00 ADR/RID: 8, III EINECS: 231-633-2		UN: 1805 GHS: Skin Irrit. 2, Eye Irrit. 2 H315, H319 P280, P302+P350, P305+P351+P338
Absorbance specification		
Max. Absorbance (1cm cell v HPLC water)		
Drift in gradient baseline		<40 mA.U.
Absorbance of the largest eluted peak:		<1.5 mA.U.

Orthophosphoric acid 85%
for electronic use (MOS)

packaging		Mfr. No
2.5 ℓ	HDPE Plastic bottle	O/0505/PB17
200 ℓ	Metal drum	O/0505/27
H ₃ PO ₄	EINECS: 231-633-2	
CAS: 7664-38-2	UN: 1805	
MW: 98.00	GHS: Skin Corr. 1B, Eye Dam. 1	
Density at 20°C: 1.68	H314	
Boiling point: 158°C	P280, P301+P330+P331,	
Melting point: 21°C	P310, P305+P351+P338,	
ADR/RID: 8, III	P303+P361+P353	
Assay: >85%		
Guaranteed analysis		
Colour	<10 APHA	Molybdenum (Mo)
Aluminium (Al)	<2 ppm	Nickel (Ni)
Barium (Ba)	<1 ppm	Potassium (K)
Cadmium (Cd)	<1 ppm	Sodium (Na)
Calcium (Ca)	<50 ppm	Strontium (Sr)
Cobalt (Co)	<1 ppm	Titanium (Ti)
Copper (Cu)	<5 ppm	n-Butyl formate
Iron (Fe)	<10 ppm	Total nitrogen (N)
Lead (Pb)	<5 ppm	Total silicon (Si)
Lithium (Li)	<1 ppm	Total sulphur (S)
Magnesium (Mg)	<10 ppm	Vanadium (V)
Manganese (Mn)	<1 ppm	Zinc (Zn)

Orthophosphoric acid 85% is equivalent to a strength of 15.2 M (45.6 N).
Filtered to 10 µm.

Orthophosphoric acid 85+%
for analysis

packaging		Mfr. No
500 mℓ	HDPE Plastic bottle	O/0500/PB08
1 ℓ	HDPE Plastic bottle	O/0500/PB15
2.5 ℓ	HDPE Plastic bottle	O/0500/PB17
2.5 ℓ	Coated Glass bottle	O/0500/PC17
25 ℓ	HDPE Plastic container	O/0500/25
200 ℓ	HDPE Plastic container	O/0500/27
H ₃ PO ₄	EINECS: 231-633-2	
CAS: 7664-38-2	UN: 1805	
MW: 98.00	GHS: Skin Corr. 1B, Eye Dam. 1	
Density at 20°C: 1.68	H314	
Boiling point: 158°C	P280, P301+P330+P331,	
Melting point: 21°C	P310, P305+P351+P338,	
ADR/RID: 8, III	P303+P361+P353	
Assay: >85%		
Guaranteed analysis		
Colour	<10 APHA	Potassium (K)
Calcium (Ca)	<50 ppm	Sodium (Na)
Copper (Cu)	<5 ppm	Total chloride (Cl)
Iron (Fe)	<10 ppm	Total nitrogen (N)
Lead (Pb)	<5 ppm	Total silicon (Si)
Magnesium (Mg)	<10 ppm	Total sulphur (S)
Nickel (Ni)	<5 ppm	Zinc (Zn)
Typical analysis		
Aluminium (Al)	<2 ppm	Manganese (Mn)
Arsenic (As)	<2 ppm	Molybdenum (Mo)
Barium (Ba)	<1 ppm	Silver (Ag)
Cadmium (Cd)	<1 ppm	Strontium (Sr)
Cobalt (Co)	<1 ppm	Titanium (Ti)
Lithium (Li)	<1 ppm	Vanadium (V)

Orthophosphoric acid is equivalent to a strength of 15.2 M or 45.6 N.
Lot analysis on the label.

UHPLC

Filtered to 0.1 µm, Fisher Chemical UHPLC gradient grade solvents (Acetonitrile, Methanol and Water) are perfect for those working at high pressure.

See page 309 for more information



Orthophosphoric acid 85+%
specified reagents for general laboratory work- extra pure

packaging		Mfr. No
500 mℓ	HDPE Plastic bottle	O/0450/PB08
1 ℓ	HDPE Plastic bottle	O/0450/PB15
2.5 ℓ	HDPE Plastic bottle	O/0450/PB17
5 ℓ	HDPE Plastic bottle	O/0450/PB21
25 ℓ	HDPE Plastic container	O/0450/25
H ₃ PO ₄	EINECS: 231-633-2	
CAS: 7664-38-2	UN: 1805	
MW: 98.00	GHS: Skin Corr. 1B, Eye Dam. 1	
Density at 20°C: 1.68	H314	
Boiling point: 158°C	P280, P301+P330+P331,	
Melting point: 21°C	P310, P305+P351+P338,	
ADR/RID: 8, III	P303+P361+P353	
Assay: >85%		
Guaranteed analysis		
Calcium (Ca)	<200 ppm	Potassium (K)
Copper (Cu)	<20 ppm	Sodium (Na)
Iron (Fe)	<50 ppm	n-Butyl formate
Lead (Pb)	<20 ppm	Total silicon (Si)
Magnesium (Mg)	<50 ppm	Total sulphur (S)
Nickel (Ni)	<20 ppm	Zinc (Zn)
Orthophosphoric acid is equivalent to a strength of 15.2 M or 45.6 N.		

Orthophosphoric acid
pure, about 100% H3PO4 (contains some H4P2O7)

packaging		Mfr. No
500 g	Plastic jar	O/0510/53
H ₃ PO ₄	EINECS: 231-633-2	
CAS: 7664-38-2	UN: 1805	
MW: 98.00	GHS: Skin Corr. 1B, Eye Dam. 1	
Density at 20°C: 1.68	H314	
Boiling point: 158°C	P280, P301+P330+P331,	
Melting point: 21°C	P310, P305+P351+P338,	
ADR/RID: 8, III	P303+P361+P353	
Guaranteed analysis		
Identity test	: Pass	

Oxalic acid anhydrous
specified reagents for general laboratory work- extra pure

packaging		Mfr. No
500 g	Plastic jar	O/0580/53
C ₂ H ₂ O ₄	ADR/RID: 8, 6.1, II	
CAS: 144-62-7	EINECS: 205-634-3	
MW: 90.04	UN: 2923	
Density at 20°C: 1.9	GHS: Acute Tox. 4	
Boiling point: DC°C	H312, H302	
Melting point: 189 to 191°C	P280, P301+P330+P331,	
Flash point: >93.4°C	P302+P352, P301+P312	
Guaranteed analysis		
Complies with standard IR spectrum: Pass		

Oxalic acid dihydrate
for analysis

packaging		Mfr. No
500 g	Plastic jar	O/0650/53
1 kg	Plastic jar	O/0650/60
5 kg	Plastic keg	O/0650/65
25 kg	Plastic keg	O/0650/70
C ₂ H ₂ O ₄ ·2H ₂ O	UN: 2923	
CAS: 6153-56-6	GHS: Acute Tox. 4	
MW: 126.07	H302, H312	
Melting point: 98 to 102°C	P280, P305+P351+P338,	
ADR/RID: 8, 6.1, II	P302+P352, P301+P312	
EINECS: 205-634-3		
Assay: >99.5%		

Guaranteed analysis		
Calcium (Ca)	<10 ppm	Total chloride (Cl)
Copper (Cu)	<2 ppm	Total nitrogen (N)
Iron (Fe)	<2 ppm	Total phosphorus (P)
Lead (Pb)	<2 ppm	Total silicon (Si)
Magnesium (Mg)	<10 ppm	Total sulphur (S)
Potassium (K)	<20 ppm	Zinc (Zn)
Sodium (Na)	<50 ppm	

Typical analysis		
Aluminium (Al)	<1 ppm	Molybdenum (Mo)
Barium (Ba)	<1 ppm	Nickel (Ni)
Cadmium (Cd)	<1 ppm	Silver (Ag)
Chromium (Cr)	<1 ppm	Strontium (Sr)
Cobalt (Co)	<1 ppm	Titanium (Ti)
Lithium (Li)	<1 ppm	Vanadium (V)
Manganese (Mn)	<1 ppm	

Lot analysis on the label.

Oxalic acid dihydrate
specified reagents for general laboratory work- extra pure

packaging		Mfr. No
500 g	Plastic jar	O/0600/53
1 kg	Plastic jar	O/0600/60
2.5 kg	Plastic jar	O/0600/62
5 kg	Plastic keg	O/0600/65
C ₂ H ₂ O ₄ ·2H ₂ O	UN: 2923	
CAS: 6153-56-6	GHS: Acute Tox. 4	
MW: 126.07	H302, H312	
Melting point: 98 to 102°C	P280, P305+P351+P338,	
ADR/RID: 8, 6.1, II	P302+P352, P301+P312	
EINECS: 205-634-3		
Assay: >99%		
Guaranteed analysis		
Calcium (Ca)	<50 ppm	Sodium (Na)
Copper (Cu)	<10 ppm	n-Butyl formate
Iron (Fe)	<10 ppm	Total phosphorus (P)
Lead (Pb)	<10 ppm	Total silicon (Si)
Magnesium (Mg)	<50 ppm	Total sulphur (S)
Potassium (K)	<100 ppm	Zinc (Zn)

Oxalic acid solution 0.05M (0.1N)
NIST Standard Concentrate, For Volumetric analysis

packaging		Mfr. No
1	Ampoule	J/5650C/05
6	Ampoule	J/5650C/90

Each ampoule makes 1 ℓ of solution.
See the full range of Standard Volumetric solution page 323.

Oxalic acid, sodium salt See Sodium oxalate, page 301

Oxygen mixture
pure, (Potassium(V)chlorate & Manganèse(IV)oxide)

packaging		Mfr. No
500 g	Glass jar	O/0850/53
ADR/RID: 5.1, II	H271, H332, H302, H411	  
UN: 1479	P261, P280, P301+P312,	
GHS: Ox. Sol. 1, Acute Tox.	P306+P360,	
4, Aquatic Chronic 2	P371+P380+P375, P273	
Guaranteed analysis		
Identity test : Pass		

Palladium standard metal solution 1000 ppm
for ICP solution in HCL 10%, NIST certified

packaging		Mfr. No
100 mℓ	Plastic bottle	SPPLPD32Y
500 mℓ	Plastic bottle	SPPLPD32X
H315, H319, H335		

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Palladium chloride(II)
specified reagents for general laboratory work- extra pure

packaging		Mfr. No
1 g	Vial	P/0120/43
PdCl ₂	GHS: Carc. 2, Acute Tox. 1, Skin Corr. 1B, Eye Dam. 1, Skin Sens. 1	  
CAS: 7647-10-1		
MW: 177.31		
Density at 180°C: 4.0	H300, H314, H317, H351	
Melting point: 678 to 680°C	P301+P310,	
ADR/RID: 6.1, 8, II	P301+P330+P331,	
EINECS: 231-596-2	P305+P351+P338, P310,	
UN: 3290	P280, P281	
Guaranteed analysis		
Identity test	: Pass	

The palladium prices vary depending on the course of the raw material.

Palmitic acid
pure, 98+%

packaging		Mfr. No
500 g	Plastic jar	P/0280/53
C ₁₆ H ₃₂ O ₂	Flash point: 206°C	
CAS: 57-10-3	EINECS: 200-312-9	
MW: 256.42	GHS: Skin Irrit. 2, Eye Irrit. 2	
Density at 20°C: 0.850	H315, H319	
Boiling point: 351.5°C	P280, P302+P352,	
Melting point: 59 to 63°C	P305+P351+P338	
Assay (CG): >98%		

(+)-Pantothenic acid calcium salt
specified reagents for general laboratory work- extra pure

Packaging			Mfr. No
100 g	Plastic jar		P/0295/48
C ₈ H ₃₂ Ca(NO ₃) ₂	Melting point: 190°C		
CAS: 137-08-6	EINECS: 205-278-9		
MW: 476.54			
Guaranteed analysis			
Loss on drying	<5%	Calcium (Ca)	82000 to 86000 ppm
Specific rotation	25° to 27°	Lead (Pb)	<30 ppm

Paraffin liquid d=0.86
For IR Spectroscopy

Packaging		Mfr. No
500 ml	Glass bottle	P/0405/08
CAS: 8012-95-1	EINECS: 232-384-2	
Density at 20°C: 0.880	GHS: Asp. Tox. 1	
Boiling point: 300°C	H304	
Melting point: -24°C	P301+P310, P331	
Flash point: 218°C		
Guaranteed analysis		
Density at 15.5°C	: 0.83 to 0.86	IR Analysis Conform

Paraffin soft white <i>pure</i>	
packaging	Mfr. No
500 g Plastic jar	P/0480/53
CAS: 64771-71-7 EINECS: 265-232-9	
Guaranteed analysis	
Complies with standard IR spectrum; Pass	

Paraffin soft yellow <i>pure</i>	
packaging	Mfr. No
500 g Plastic jar	P/0440/53
CAS: 64771-71-7 EINECS: 265-232-9	
Guaranteed analysis	
Complies with standard IR spectrum; Pass	

Paraffin liquid water white <i>pure, d=0.86-0.89</i>	
packaging	Mfr. No
2.5 ℓ Glass bottle	P/0400/17
25 ℓ HDPE Plastic container	P/0400/25
CAS: 8012-95-1 Density at 20°C: 0.880 Boiling point: 300°C	Melting point: -24°C Flash point: 218°C EINECS: 232-384-2
Guaranteed analysis	
Density 15.5°C 0.86 to 0.89	

Paraffin liquid water white <i>pure, d=0.83-0.86</i>	
packaging	Mfr. No
500 mℓ Glass bottle	P/0360/08
2.5 ℓ Glass bottle	P/0360/17
25 ℓ HDPE Plastic container	P/0360/25
200 ℓ HDPE Plastic container	P/0360/27
CAS: 8012-95-1 Density at 20°C: 0.880 Boiling point: 300°C Melting point: -24°C Flash point: 218°C	EINECS: 232-384-2 GHS: Asp. Tox. 1 H304 P301+P310, P331
Guaranteed analysis	
Density at 15.5°C : 0.83 to 0.86	

Paraffin liquid <i>Technical, d=0.88 for oil baths</i>	
packaging	Mfr. No
2.5 ℓ Glass bottle	P/0320/17
25 ℓ HDPE Plastic container	P/0320/25
CAS: 8012-95-1 Density at 20°C: 0.880 Boiling point: 300°C	Melting point: -24°C Flash point: 218°C EINECS: 232-384-2
Guaranteed analysis	
Density at 15.5°C : 0.86 to 0.89	

Paraffin wax pastillated <i>pure, melting point about 60°C</i>	
packaging	Mfr. No
450 g Plastic jar	P/0680/90
CAS: 8002-74-2 Boiling point: 322°C Melting point: 58 to 62°C	Flash point: 198°C EINECS: 232-315-6

Paraffin wax pastillated <i>pure, melting point about 57°C</i>	
packaging	Mfr. No
450 g Plastic jar	P/0650/90
CAS: 8002-74-2 Boiling point: 322°C Melting point: 56 to 58°C	Flash point: 198°C EINECS: 232-315-6

Paraffin wax pastillated <i>pure, melting point about 52°C</i>	
packaging	Mfr. No
450 g Plastic jar	P/0600/90
CAS: 8002-74-2 Boiling point: 322°C Melting point: 51 to 53°C	Flash point: 198°C EINECS: 232-315-6

Paraffin wax <i>pure, clearing point about 65-71°C</i>	
packaging	Mfr. No
450 g Plastic jar	P/0720/90
CAS: 8002-74-2 Boiling point: 322°C Melting point: 58 to 62°C	Flash point: 198°C EINECS: 232-315-6
Guaranteed analysis	
Clearing point : 65 to 71°C	

Paraformaldehyde <i>specified reagents for general laboratory work- extra pure</i>		SLR
packaging	Mfr. No	
500 g Plastic jar	P/0840/53	
(CH2 O) _n CAS: 30525-89-4 Density at 20°C: 1.46 Melting point: 120 to 170°C Flash point: 71°C ADR/RID: 4.1, III EINECS: 200-001-8 UN: 2213 GHS: Flam. Sol. 2, Carc. 2,	Acute Tox. 4, Skin Irrit. 2, Eye Irrit. 2, Skin Sens. 1, STOT SE 3 H228, H302, H315, H335, H319, H351, H317, H332 P280, P281, P261, P301+P312, P302+P352, P305+P351+P338, P233, P210	
Guaranteed analysis		
Complies with standard IR spectrum: Pass		

Pectin <i>pure, citrus - 121 grade slow setting 121</i>	
packaging	Mfr. No
100 g Plastic jar	P/0920/48
CAS: 9000-69-5 EINECS: 232-553-0	
Guaranteed analysis	
Identity test : Pass	

Pelargonic alcohol *See 1-Nonanol, page 219*

Save money!

Order 4 bottles in a box



iso-Pentane NEW PACK.		Certified AR
<i>for analysis</i>		
packaging	Mfr. No	
500 mℓ Glass bottle	P/1031/08	
2.5 ℓ Glass bottle	P/1031/17	
5 ℓ Mild steel drum	P/1031/MS21	
C ₅ H ₁₂ CAS: 78-78-4 MW: 72.15 Density at 20°C: 0.620 Boiling point: 30°C Melting point: -160°C Flash point: -51°C ADR/RID: 3, I	EINECS: 201-142-8 UN: 1265 GHS: Flam. Liq. 1, STOT SE 3, Asp. Tox. 1, Aquatic Chronic 2 H224, H336, H411, H304, EUH066 P261, P301+P310, P331, P273, P210, P243	

Assay (CG): >99%		
Guaranteed analysis		
Acidity / Alkalinity	<0.0002 meq/g	Magnesium (Mg) <0.1 ppm
Colour	<10 APHA	Manganese (Mn) <0.02 ppm
Water	<0.01%	Molybdenum (Mo) <0.02 ppm
Residue after evaporation	<50 ppm	Nickel (Ni) <0.02 ppm
Aluminium (Al)	<0.1 ppm	Potassium (K) <0.2 ppm
Barium (Ba)	<0.02 ppm	Sodium (Na) <1 ppm
Cadmium (Cd)	<0.02 ppm	Strontium (Sr) <0.02 ppm
Calcium (Ca)	<0.5 ppm	Titanium (Ti) <0.02 ppm
Chromium (Cr)	<0.02 ppm	Total phosphorus (P) <0.5 ppm
Cobalt (Co)	<0.02 ppm	Total silicon (Si) <0.05 ppm
Copper (Cu)	<0.05 ppm	Total sulphur (S) <5 ppm
Iron (Fe)	<0.2 ppm	Vanadium (V) <0.02 ppm
Lead (Pb)	<0.05 ppm	Zinc (Zn) <0.1 ppm
Lithium (Li)	<0.02 ppm	Pentane (mixed isomers) <1%
<i>Lot analysis on the label.</i>		

iso-Pentane NEW PACK.		SLR
<i>specified reagents for general laboratory work- extra pure</i>		
packaging	Mfr. No	
500 mℓ Glass bottle	P/1030/08	
2.5 ℓ Glass bottle	P/1030/17	
5 ℓ Mild steel drum	P/1030/MS21	
C ₅ H ₁₂ CAS: 78-78-4 MW: 72.15 Density at 20°C: 0.620 Boiling point: 30°C Melting point: -160°C Flash point: -51°C ADR/RID: 3, I	EINECS: 201-142-8 UN: 1265 GHS: Flam. Liq. 1, STOT SE 3, Asp. Tox. 1, Aquatic Chronic 2 H224, H336, H411, H304, EUH066 P261, P301+P310, P331, P273, P210, P243	

Assay (CG): >99%		
Guaranteed analysis		
Acidity / Alkalinity	<0.002 meq/g	Magnesium (Mg) <2 ppm
Colour	<10 APHA	Potassium (K) <2 ppm
Water	<0.05%	Sodium (Na) <10 ppm
Residue after evaporation	<50 ppm	Total phosphorus (P) <2 ppm
Calcium (Ca)	<5 ppm	Total silicon (Si) <2 ppm
Copper (Cu)	<2 ppm	Total sulphur (S) <20 ppm
Iron (Fe)	<5 ppm	Zinc (Zn) <2 ppm
Lead (Pb)	<2 ppm	

pH Buffers

NIST Certified Standard Solution & Concentrate

A full range to discover on page 63



Pentane (mixed isomers) NEW PACK.	
for HPLC	
packaging	Mfr. No
1 ℓ Glass bottle	P/1006/15
2.5 ℓ Glass bottle	P/1006/17
10 ℓ Returnable Drum	P/1006/10RSS
25 ℓ Returnable Drum	P/1006/25RSS
30 ℓ Returnable Drum	P/1006/30RSS
C ₅ H ₁₂ CAS: 109-66-0 MW: 72.15 Density at 20°C: 0.626 Boiling point: 36°C Melting point: -130°C Flash point: -49°C ADR/RID: 3, II	EINECS: 203-692-4 UN: 1265 GHS: Flam. Liq. 2, STOT SE 3, Asp. Tox. 1, Aquatic Chronic 2 H225, H304, H336, H411, EUH066 P261, P301+P310, P331, P273, P210, P240



Assay (CG): >95%		
Guaranteed analysis		
Acidity / Alkalinity	<0.0002 meq/g	Residue after evaporation <2 ppm
Water	<0.01%	
Absorbance specification		
Max. Absorbance (1cm cell v HPLC water)		
200 nm	<0.25 A.U.	230 nm <0.01 A.U.
210 nm	<0.2 A.U.	240 nm <0.005 A.U.
220 nm	<0.05 A.U.	250 nm <0.005 A.U.
<i>Lot analysis on the label.</i>		
<i>Filtered to 0.2 µm.</i>		

Pentane, (95% n-Pentane approx.) <i>for residue and pesticide analysis</i>		Distol
packaging	Mfr. No	
1 ℓ Glass bottle	P/1024/15	
2.5 ℓ Glass bottle	P/1024/17	
C ₅ H ₁₂ CAS: 109-66-0 MW: 72.15 Density at 20°C: 0.626 Boiling point: 36°C Melting point: -130°C Flash point: -49°C ADR/RID: 3, II	EINECS: 203-692-4 UN: 1265 GHS: Flam. Liq. 2, STOT SE 3, Asp. Tox. 1, Aquatic Chronic 2 H225, H304, H336, H411, EUH066 P261, P301+P310, P331, P273, P210, P240	

Assay (CG): >95%		
Guaranteed analysis		
Water	<0.02%	GC (ECD) Analysis Pass
Residue after evaporation	<2 ppm	GC (FID) Analysis Pass
GC (NPD) analysis	Pass	

GC Analysis
Analysis of pesticides: *from the retention time of Lindane, after concentration of 2500 times, no impurity peak higher than that obtained with 10 pg / mL héptachlo-réoxide (ECD) and 5 pg / mL Parathion (NPD).*
Analysis of hydrocarbons: *no impurity chain longer than the dodecane is present at a concentration greater than 0.01 ppm (FID).*
Conditions :
Column: J&W DB 5: 30m x 0.32mm x 0.1µm df.
Temperature profile: (1) 50°C, 2 min
(2) 40°C / min until 320°C
(3) 320°C for 30 min
Injection: 5µℓ cold on the column.

n-Pentane
for analysis

Certified AR

packaging		Mfr. No
1 ℓ	Glass bottle	P/1021/15
2.5 ℓ	Glass bottle	P/1021/17
2.5 ℓ	Coated Glass bottle	P/1021/PB17
C ₅ H ₁₂ CAS: 109-66-0 MW: 72.15 Density at 20°C: 0.626 Boiling point: 36°C Melting point: -130°C Flash point: -49°C ADR/RID: 3, II	EINECS: 203-692-4 UN: 1265 GHS: Flam. Liq. 2, STOT SE 3, Asp. Tox. 1, Aquatic Chronic 2 H225, H304, H336, H411, EUH066 P261, P301+P310, P331, P273, P210, P240	   

Assay (CG): >99%		
Guaranteed analysis		
Acidity / Alkalinity	<0.0001 meq/g	Magnesium (Mg) <0.1 ppm
Colour	<10 APHA	Potassium (K) <0.2 ppm
Water	<0.01%	Sodium (Na) <1 ppm
Residue after evaporation	<20 ppm	Total phosphorus (P) <0.5 ppm
Calcium (Ca)	<0.5 ppm	Total silicon (Si) <0.02 ppm
Copper (Cu)	<0.05 ppm	Total sulphur (S) <5 ppm
Iron (Fe)	<0.2 ppm	Zinc (Zn) <0.1 ppm
Lead (Pb)	<0.05 ppm	

Typical analysis			
Aluminium (Al)	<0.1 ppm	Manganese (Mn)	<0.02 ppm
Barium (Ba)	<0.02 ppm	Molybdenum (Mo)	<0.02 ppm
Cadmium (Cd)	<0.02 ppm	Nickel (Ni)	<0.02 ppm
Chromium (Cr)	<0.02 ppm	Strontium (Sr)	<0.02 ppm
Cobalt (Co)	<0.02 ppm	Titanium (Ti)	<0.02 ppm
Lithium (Li)	<0.02 ppm	Vanadium (V)	<0.02 ppm

Lot analysis on the label.

n-Pentane **NEW PACK.**
specified reagents for general laboratory work- extra pure

SLR

packaging		Mfr. No
500 mℓ	Glass bottle	P/1020/08
2.5 ℓ	Glass bottle	P/1020/17
2.5 ℓ	Coated Glass bottle	P/1020/PB17
5 ℓ	Fluorinated HDPE Plastic bottle	P/1020/21
25 ℓ	Metal drum	P/1020/25
C ₅ H ₁₂ CAS: 109-66-0 MW: 72.15 Density at 20°C: 0.626 Boiling point: 36°C Melting point: -130°C Flash point: -49°C ADR/RID: 3, II	EINECS: 203-692-4 UN: 1265 GHS: Flam. Liq. 2, STOT SE 3, Asp. Tox. 1, Aquatic Chronic 2 H225, H304, H336, H411, EUH066 P261, P301+P310, P331, P273, P210, P240	   

Assay (CG): >99%		
Guaranteed analysis		
Acidity / Alkalinity	<0.001 meq/g	Magnesium (Mg) <2 ppm
Colour	<10 APHA	Potassium (K) <2 ppm
Water	<0.05%	Sodium (Na) <10 ppm
Residue after evaporation	<50 ppm	Total phosphorus (P) <2 ppm
Calcium (Ca)	<5 ppm	Total silicon (Si) <2 ppm
Copper (Cu)	<2 ppm	Total sulphur (S) <20 ppm
Iron (Fe)	<5 ppm	Zinc (Zn) <2 ppm
Lead (Pb)	<2 ppm	

Pentanedial See Glutaraldehyde 50% solution, page 136

1-Pentanesulfonic acid sodium salt
for lon pair chromatography in electrochemical detection

Certified

packaging		Mfr. No
25 g	Plastic jar	P/1039/46
C ₅ H ₁₁ NaO ₃ S CAS: 22767-49-3 MW: 174.21 Melting point: 300°C EINECS: 245-208-4	GHS: Skin Irrit. 2, Eye Irrit. 2, STOT SE 3 H319, H335, H315 P261, P302+P352, P280, P305+P351+P338	
Assay: >98%		
Guaranteed analysis		
Electrochemical test Pass		
Absorbance specification		
Max. Absorbance (1cm cell v HPLC water)		
200 nm (sol. aqu. 0.25 M)	<0.1 A.U.	240 nm (sol. aqu. 0.25 M) ... <0.02 A.U.
210 nm (sol. aqu. 0.25 M) ...	<0.05 A.U.	250 nm (sol. aqu. 0.25 M) ... <0.02 A.U.
220 nm (sol. aqu. 0.25 M) ...	<0.03 A.U.	260 nm (sol. aqu. 0.25 M) ... <0.02 A.U.

tested by HPLC-EC to ensure low electrochemical response.
Detector - ESA Coulochem II
Mode-oxidative
Graphite Electrode.
Lot analysis on the label.

1-Pentanesulfonic acid sodium salt
for lon pair chromatography

packaging		Mfr. No
25 g	Plastic jar	P/1038/46
100 g	Plastic jar	P/1038/48
C ₅ H ₁₁ NaO ₃ S CAS: 22767-49-3 MW: 174.21 Melting point: 300°C EINECS: 245-208-4	GHS: Skin Irrit. 2, Eye Irrit. 2, STOT SE 3 H319, H335, H315 P261, P302+P352, P280, P305+P351+P338	
Assay: >98%		
Absorbance specification		
Max. Absorbance (1cm cell v HPLC water)		
200 nm (sol. aqu. 0.25 M)	<0.1 A.U.	240 nm (sol. aqu. 0.25 M) ... <0.02 A.U.
210 nm (sol. aqu. 0.25 M) ...	<0.05 A.U.	250 nm (sol. aqu. 0.25 M) ... <0.02 A.U.
220 nm (sol. aqu. 0.25 M) ...	<0.03 A.U.	260 nm (sol. aqu. 0.25 M) ... <0.02 A.U.

Lot analysis on the label.

1-Pentanesulfonic acid sodium salt solution 0.1M
for lon pair chromatography

packaging		Mfr. No
250 mℓ	Glass bottle	P/1036/07
Absorbance specification		
Max. Absorbance (1cm cell v HPLC water)		
200 nm	<0.04 A.U.	240 nm <0.008 A.U.
210 nm	<0.02 A.U.	250 nm <0.008 A.U.
220 nm	<0.016 A.U.	260 nm <0.008 A.U.
230 nm	<0.012 A.U.	

Lot analysis on the label.

Pentan-1-ol
for analysis

Certified AR

packaging		Mfr. No
500 mℓ	Coated Glass bottle	P/1041/PB08
2.5 ℓ	Coated Glass bottle	P/1041/PB17
C ₅ H ₁₂ O CAS: 71-41-0 MW: 88.15 Density at 20°C: 0.811 Boiling point: 137 to 139°C Melting point: -78°C Flash point: 49°C ADR/RID: 3, III	EINECS: 200-752-1 UN: 1105 GHS: Flam. Liq. 3, Acute Tox. 4, Skin Irrit. 2, STOT SE 3 H226, H315, H335, H332 P261, P304+P340, P302+P352, P210	 
Assay (CG): >98.5%		

Guaranteed analysis		
Acidity / Alkalinity	<0.002 meq/g	Manganese (Mn) <0.02 ppm
Colour	<10 APHA	Molybdenum (Mo) <0.02 ppm
Water	<0.2%	Nickel (Ni) <0.05 ppm
Residue after evaporation	<10 ppm	Potassium (K) <0.1 ppm
Aluminium (Al)	<0.05 ppm	Sodium (Na) <1 ppm
Barium (Ba)	<0.02 ppm	Strontium (Sr) <0.02 ppm
Cadmium (Cd)	<0.02 ppm	Titanium (Ti) <0.02 ppm
Calcium (Ca)	<0.5 ppm	Total phosphorus (P) <0.2 ppm
Chromium (Cr)	<0.05 ppm	Total silicon (Si) <0.05 ppm
Cobalt (Co)	<0.02 ppm	Total sulphur (S) <0.5 ppm
Copper (Cu)	<0.05 ppm	Vanadium (V) <0.02 ppm
Iron (Fe)	<0.2 ppm	Zinc (Zn) <0.5 ppm
Lead (Pb)	<0.05 ppm	iso-Amyl acetate <0.1%
Lithium (Li)	<0.02 ppm	iso-Amyl alcohol <0.5%
Magnesium (Mg)	<0.1 ppm	

Lot analysis on the label.

Pentan-1-ol
specified reagents for general laboratory work- extra pure

SLR

packaging		Mfr. No
500 mℓ	Coated Glass bottle	P/1040/PB08
2.5 ℓ	Coated Glass bottle	P/1040/PB17
C ₅ H ₁₂ O CAS: 71-41-0 MW: 88.15 Density at 20°C: 0.811 Boiling point: 137 to 139°C Melting point: -78°C Flash point: 49°C ADR/RID: 3, III	EINECS: 200-752-1 UN: 1105 GHS: Flam. Liq. 3, Acute Tox. 4, Skin Irrit. 2, STOT SE 3 H226, H315, H335, H332 P261, P304+P340, P302+P352, P210	 
Assay (CG): >98%		
Guaranteed analysis		
Acidity / Alkalinity	<0.005 meq/g	Magnesium (Mg) <2 ppm
Colour	<10 APHA	Potassium (K) <2 ppm
Water	<0.3%	Sodium (Na) <10 ppm
Residue after evaporation	<100 ppm	Total phosphorus (P) <2 ppm
Calcium (Ca)	<5 ppm	Total silicon (Si) <2 ppm
Copper (Cu)	<2 ppm	Total sulphur (S) <5 ppm
Iron (Fe)	<5 ppm	Zinc (Zn) <2 ppm
Lead (Pb)	<2 ppm	

Pentan-2-ol (sec-Amyl alcohol)
specified reagents for general laboratory work- extra pure

SLR

packaging		Mfr. No
500 mℓ	Coated Glass bottle	P/1080/PB08
C ₅ H ₁₂ O CAS: 6032-29-7 MW: 88.15 Density at 20°C: 0.812 Boiling point: 118 to 119°C Melting point: -50°C Flash point: 34°C ADR/RID: 3, III	EINECS: 227-907-6 UN: 1105 GHS: Flam. Liq. 3, Acute Tox. 4, STOT SE 3 H226, H332, H335, H302, EUH066 P261, P301+P312, P304+P340, P210	 
Assay (CG): >98%		

Pentan-3-one
technical

packaging		Mfr. No
250 mℓ	Glass bottle	P/1100/07
C ₅ H ₁₀ O CAS: 96-22-0 MW: 86.13 Density at 20°C: 0.815 Boiling point: 100 to 102°C Melting point: -40°C Flash point: 7°C	ADR/RID: 3, II EINECS: 202-490-3 UN: 1156 GHS: Flam. Liq. 2, STOT SE 3 H225, H336, H335, EUH066 P261, P210, P240	 
Guaranteed analysis		
Complies with standard IR spectrum: Pass		
Diethyl ketone		

1-Pentyl acetate See n-Amyl acetate, page 47

Pentyl alcohol See Pentan-1-ol, page 226

Pepsin powder
specified reagents for general laboratory work- extra pure


SLR

packaging		Mfr. No
25 g	Plastic jar	P/1120/46
100 g	Plastic jar	P/1120/48
500 g	Plastic jar	P/1120/53
CAS: 9001-75-6 EINECS: 232-629-3 GHS: Skin Irrit. 2, Eye Irrit. 2, Resp. Sens. 1, STOT SE 3	H315, H319, H335, H334 P304+P341, P342+P311, P261, P302+P352, P280, P305+P351+P338	 

Guaranteed analysis	
Complies with standard IR spectrum: Pass	

Peptone
pure, For bacteriological work

packaging		Mfr. No
100 g	Plastic jar	P/1160/48
250 g	Plastic jar	P/1160/50
CAS: 73049-73-7		
Guaranteed analysis		
Complies with standard IR spectrum: Pass		



Perchloric acid 65-71%
for ultra trace elemental analysis

packaging		Mfr. No
250 mℓ	Pre-cleaned Teflon FEP bottle	A469-250
500 mℓ	Pre-cleaned Teflon FEP bottle	A469-500
1 ℓ	Pre-cleaned Teflon FEP bottle	A469-1
HCℓO ₄	EINECS: 231-512-4	
CAS: 7601-90-3	UN: 1873	
MW: 100.46	GHS: Ox. Liq. 1, Skin Corr.	
Density at 20°C: 1.66	1A, Eye Dam. 1	
Boiling point: 203°C	H271, H314	
Melting point: -18°C	P301+P330+P331, P280,	
Flash point: 113°C	P305+P351+P338, P310, P210	
ADR/RID: 5.1, 8, I		

Assay: 65-71%			
Guaranteed analysis			
Aluminium (Al)	<100 ppt	Nickel (Ni)	<100 ppt
Antimony (Sb)	<100 ppt	Niobium (Nb)	<10 ppt
Arsenic (As)	<100 ppt	Palladium (Pd)	<10 ppt
Barium (Ba)	<100 ppt	Platinum (Pt)	<100 ppt
Beryllium (Be)	<100 ppt	Potassium (K)	<100 ppt
Bismuth (Bi)	<10 ppt	Praseodymium (Pr)	<10 ppt
Cadmium (Cd)	<10 ppt	Rhodium (Rh)	<10 ppt
Calcium (Ca)	<100 ppt	Rubidium (Rb)	<10 ppt
Cerium (Ce)	<10 ppt	Samarium (Sm)	<10 ppt
Cesium (Cs)	<10 ppt	Scandium (Sc)	<100 ppt
Cobalt (Co)	<100 ppt	Silver (Ag)	<100 ppt
Copper (Cu)	<100 ppt	Sodium (Na)	<100 ppt
Dysprosium (Dy)	<10 ppt	Strontium (Sr)	<100 ppt
Erbium (Er)	<10 ppt	Tantalum (Ta)	<30 ppt
Europium (Eu)	<10 ppt	Tellurium (Te)	<10 ppt
Gadolinium (Gd)	<10 ppt	Terbium (Tb)	<10 ppt
Gallium (Ga)	<10 ppt	Thallium (Tl)	<100 ppt
Hafnium (Hf)	<10 ppt	Thorium (Th)	<10 ppt
Holmium (Ho)	<10 ppt	Thulium (Tm)	<10 ppt
Indium (In)	<10 ppt	Tin (Sn)	<100 ppt
Iron (Fe)	<100 ppt	Titanium (Ti)	<100 ppt
Lanthanum (La)	<10 ppt	Tungsten (W)	<20 ppt
Lead (Pb)	<10 ppt	Uranium (U)	<10 ppt
Lithium (Li)	<100 ppt	Vanadium (V)	<100 ppt
Lutetium (Lu)	<10 ppt	Ytterbium (Yb)	<10 ppt
Magnesium (Mg)	<100 ppt	Yttrium (Y)	<10 ppt
Manganese (Mn)	<100 ppt	Zinc (Zn)	<100 ppt
Molybdenum (Mo)	<100 ppt	Zirconium (Zr)	<100 ppt
Neodymium (Nd)	<10 ppt		

Bottled in class 10 cleanroom conditions. Individually double-bagged in class 100 cleanroom conditions. Each bottle is individually packaged. Certificate of analysis delivered with each bottle.

Perchloric acid 65-71%
for trace metal analysis

packaging		Mfr. No
500 mℓ	HDPE plastic bottle	A511-P500
500 mℓ	Coated Glass bottle	P/1292/PB08
1 ℓ	HDPE plastic bottle	A511-P1
1 ℓ	Coated Glass bottle	P/1292/PB15
2.5 ℓ	HDPE plastic bottle	A511-P212
HCℓO ₄	EINECS: 231-512-4	
CAS: 7601-90-3	UN: 1873	
MW: 100.46 g/mol	GHS: Ox. Liq. 1, Skin Corr.	
Density at 20°C: 1.66°C	1A, Eye Dam. 1	
Boiling point: 203°C	H271, H314	
Melting point: -18°C	P301+P330+P331, P280,	
Flash point: 113°C	P305+P351+P338, P310, P210	
ADR/RID: 5.1, 8, I		

Assay: 65-71%			
Guaranteed analysis			
Colour	<10 APHA	Cerium (Ce)	<0.5 ppb
Aluminium (Al)	<1 ppb	Cesium (Cs)	<0.5 ppb
Antimony (Sb)	<0.5 ppb	Cobalt (Co)	<0.5 ppb
Arsenic (As)	<0.5 ppb	Copper (Cu)	<0.5 ppb
Barium (Ba)	<1 ppb	Dysprosium (Dy)	<0.5 ppb
Beryllium (Be)	<0.5 ppb	Erbium (Er)	<0.5 ppb
Bismuth (Bi)	<0.5 ppb	Europium (Eu)	<0.5 ppb
Cadmium (Cd)	<1 ppb	Gadolinium (Gd)	<0.5 ppb
Calcium (Ca)	<1 ppb		

Gallium (Ga)	<0.5 ppb	Rubidium (Rb)	<0.5 ppb
Hafnium (Hf)	<0.5 ppb	Samarium (Sm)	<0.5 ppb
Holmium (Ho)	<0.5 ppb	Scandium (Sc)	<0.5 ppb
Indium (In)	<0.5 ppb	Silver (Ag)	<1 ppb
Iron (Fe)	<1 ppb	Sodium (Na)	<1 ppb
Lanthanum (La)	<0.5 ppb	Strontium (Sr)	<0.5 ppb
Lead (Pb)	<1 ppb	Tantalum (Ta)	<1 ppb
Lithium (Li)	<0.5 ppb	Tellurium (Te)	<0.5 ppb
Lutetium (Lu)	<0.5 ppb	Terbium (Tb)	<0.5 ppb
Magnesium (Mg)	<1 ppb	Thallium (Tl)	<0.5 ppb
Manganese (Mn)	<1 ppb	Thorium (Th)	<1 ppb
Molybdenum (Mo)	<0.5 ppb	Thulium (Tm)	<0.5 ppb
Neodymium (Nd)	<0.5 ppb	Tin (Sn)	<1 ppb
Nickel (Ni)	<1 ppb	Titanium (Ti)	<1 ppb
Niobium (Nb)	<10 ppb	Tungsten (W)	<0.5 ppb
Palladium (Pd)	<0.5 ppb	Uranium (U)	<0.5 ppb
Platinum (Pt)	<0.5 ppb	Vanadium (V)	<0.5 ppb
Potassium (K)	<1 ppb	Ytterbium (Yb)	<0.5 ppb
Praseodymium (Pr)	<0.5 ppb	Yttrium (Y)	<0.5 ppb
Rhodium (Rh)	<0.5 ppb	Zinc (Zn)	<0.5 ppb

Bottled in class 10 cleanroom conditions. Lot analysis on the label.

Perchloric acid 70% d=1.67
for analysis

packaging		Mfr. No
500 mℓ	Coated Glass bottle	P/1280/PB08
1 ℓ	Coated Glass bottle	P/1280/PB15
2.5 ℓ	Coated Glass bottle	P/1280/PB17

HCℓO ₄	EINECS: 231-512-4	
CAS: 7601-90-3	UN: 1873	
MW: 100.46	GHS: Ox. Liq. 1, Skin Corr.	
Density at 20°C: 1.66	1A, Eye Dam. 1	
Boiling point: 203°C	H271, H314	
Melting point: -18°C	P301+P330+P331, P280,	
Flash point: 113°C	P305+P351+P338, P310, P210	
ADR/RID: 5.1, 8, I		

Assay: 68 to 72%			
Guaranteed analysis			
Colour	<10 APHA	Potassium (K)	<0.2 ppm
Calcium (Ca)	<1 ppm	Sodium (Na)	<2 ppm
Chlorate	<0.001%	Total chloride (Cl)	<0.0001%
Chromium (Cr)	<0.5 ppm	Total nitrogen (N)	<20 ppm
Copper (Cu)	<0.2 ppm	Total phosphorus (P)	<2 ppm
Iron (Fe)	<1 ppm	Total silicon (Si)	<1 ppm
Lead (Pb)	<0.2 ppm	Total sulphur (S)	<5 ppm
Magnesium (Mg)	<0.5 ppm	Zinc (Zn)	<0.5 ppm
Nickel (Ni)	<0.1 ppm		

Typical analysis			
Aluminium (Al)	<0.5 ppm	Manganese (Mn)	<0.02 ppm
Arsenic (As)	<0.2 ppm	Molybdenum (Mo)	<0.02 ppm
Barium (Ba)	<0.05 ppm	Silver (Ag)	<0.02 ppm
Cadmium (Cd)	<0.02 ppm	Strontium (Sr)	<0.02 ppm
Cobalt (Co)	<0.05 ppm	Titanium (Ti)	<0.02 ppm
Lithium (Li)	<0.02 ppm	Vanadium (V)	<0.02 ppm

Perchloric acid 70% is equivalent to a strength of 11.6 M (11.6 N). Lot analysis on the label.

Perchloric acid 60%
for HPLC

packaging		Mfr. No
100 mℓ	Coated Glass bottle	P/1290/PB05

HCℓO ₄	EINECS: 231-512-4	
CAS: 7601-90-3	UN: 1873	
MW: 100.46	GHS: Ox. Liq. 1, Skin Corr.	
Density at 20°C: 1.66	1A, Eye Dam. 1	
Boiling point: 203°C	H271, H314	
Melting point: -18°C	P301+P330+P331, P280,	
Flash point: 113°C	P305+P351+P338, P310, P210	
ADR/RID: 5.1, 8, I		

Assay: 60 to 62%			
Guaranteed analysis			
Residue after evaporation	<30 ppm	Lead (Pb)	<0.5 ppm
Iron (Fe)	<1 ppm		

Absorbance specification	
Max. Absorbance (1 cm cell v HPLC water)	
254 nm (sol. aqu. 1 M)	< 0.01 A.U.

Perchloric acid 60% is equivalent to a strength of 9.2 M (9.2 N). Lot analysis on the label. Filtered to 0.2 μm.

Perchloric acid 60% d=1.53
for analysis

packaging		Mfr. No
500 mℓ	Coated Glass bottle	P/1240/PB08
2.5 ℓ	Coated Glass bottle	P/1240/PB17
25 ℓ	Glass bottle	P/1240/25

HCℓO ₄	EINECS: 231-512-4	
CAS: 7601-90-3	UN: 1873	
MW: 100.46	GHS: Ox. Liq. 1, Skin Corr.	
Density at 20°C: 1.66	1A, Eye Dam. 1	
Boiling point: 203°C	H271, H314	
Melting point: -18°C	P301+P330+P331, P280,	
Flash point: 113°C	P305+P351+P338, P310, P210	
ADR/RID: 5.1, 8, I		

Assay: 60 to 62%			
Guaranteed analysis			
Colour	<10 APHA	Potassium (K)	<0.5 ppm
Calcium (Ca)	<1 ppm	Sodium (Na)	<5 ppm
Chlorate	<0.002%	Total chloride (Cl)	<0.0005%
Chromium (Cr)	<0.2 ppm	Total nitrogen (N)	<20 ppm
Copper (Cu)	<0.2 ppm	Total phosphorus (P)	<2 ppm
Iron (Fe)	<1 ppm	Total silicon (Si)	<1 ppm
Lead (Pb)	<0.2 ppm	Total sulphur (S)	<5 ppm
Magnesium (Mg)	<0.5 ppm	Zinc (Zn)	<0.1 ppm
Nickel (Ni)	<0.1 ppm		

Typical analysis			
Aluminium (Al)	<0.5 ppm	Manganese (Mn)	<0.02 ppm
Arsenic (As)	<0.1 ppm	Molybdenum (Mo)	<0.05 ppm
Barium (Ba)	<0.05 ppm	Silver (Ag)	<0.02 ppm
Cadmium (Cd)	<0.05 ppm	Strontium (Sr)	<0.02 ppm
Cobalt (Co)	<0.05 ppm	Titanium (Ti)	<0.02 ppm
Lithium (Li)	<0.02 ppm	Vanadium (V)	<0.02 ppm

Perchloric acid 60% is equivalent to a strength of 9.2 M (9.2 N). Lot analysis on the label.

Perchloric acid 60%
specified reagents for general laboratory work- extra pure

packaging		Mfr. No
500 mℓ	Coated Glass bottle	P/1200/PB08
2.5 ℓ	Coated Glass bottle	P/1200/PB17

HCℓO ₄	EINECS: 231-512-4	
CAS: 7601-90-3	UN: 1873	
MW: 100.46	GHS: Ox. Liq. 1, Skin Corr.	
Density at 20°C: 1.66	1A, Eye Dam. 1	
Boiling point: 203°C	H271, H314	
Melting point: -18°C	P301+P330+P331, P280,	
Flash point: 113°C	P305+P351+P338, P310, P210	
ADR/RID: 5.1, 8, I		

Assay: >60%			
Guaranteed analysis			
Calcium (Ca)	<10 ppm	Potassium (K)	<5 ppm
Chromium (Cr)	<2 ppm	Sodium (Na)	<50 ppm
Copper (Cu)	<5 ppm	n-Butyl formate	<0.001%
Iron (Fe)	<10 ppm	Total phosphorus (P)	<20 ppm
Lead (Pb)	<5 ppm	Total silicon (Si)	<5 ppm
Magnesium (Mg)	<5 ppm	Total sulphur (S)	<50 ppm
Nickel (Ni)	<2 ppm	Zinc (Zn)	<5 ppm

Perchloric acid 60% is equivalent to a strength of 9.2 M (9.2 N).

Perchloric acid solution 0.1M (0.1N)
NIST Standard solution ready to use, For Volumetric analysis

packaging		Mfr. No
1 ℓ	HDPE Plastic bottle	J/5820/PB15
2.5 ℓ	HDPE Plastic bottle	J/5820/PB17

CAS: 7601-90-3	1A, Eye Dam. 1	
Flash point: -40°C	H226, H314	
ADR/RID: 8, 3, II	P280, P301+P330+P331,	
EINECS: 231-512-4	P305+P351+P338, P310,	
UN: 2920	P302+P352, P304+P340, P210	
GHS: Flam. Liq. 3, Skin Corr.		

Manufactured to a tolerance of ± 0.2%. See the full range of Standard Volumetric solution page 323.

Perchloroethylene See Tetrachloroethylene, page 344


Periodic acid
for analysis

packaging		Mfr. No
100 g	Plastic jar	P/1322/48

H ₅ IO ₆	EINECS: 233-937-0	
CAS: 10450-60-9	UN: 3085	
MW: 227.94	GHS: Ox. Sol. 2, Skin Corr. 1B	
Density at 20°C: 1.4	H272, H314	
Melting point: 122°C	P301+P330+P331, P280,	
ADR/RID: 5.1, 8, II	P305+P351+P338, P310, P210	

Assay: >99.5%			
Guaranteed analysis			
Aluminium (Al)	<2 ppm	Molybdenum (Mo)	<2 ppm
Barium (Ba)	<2 ppm	Nickel (Ni)	<5 ppm
Cadmium (Cd)	<2 ppm	Potassium (K)	<20 ppm
Calcium (Ca)	<10 ppm	Propan--1-ol	<1 ppm
Chromium (Cr)	<5 ppm	Sodium (Na)	<50 ppm
Cobalt (Co)	<2 ppm	Strontium (Sr)	<2 ppm
Copper (Cu)	<5 ppm	Titanium (Ti)	<1 ppm
Iodide (I)	<0.001%	Total nitrogen (N)	<50 ppm
Iron (Fe)	<5 ppm	Total silicon (Si)	<20 ppm
Lead (Pb)	<5 ppm	Total sulphur (S)	<50 ppm
Lithium (Li)	<1 ppm	Vanadium (V)	<1 ppm
Magnesium (Mg)	<5 ppm	Zinc (Zn)	<10 ppm
Manganese (Mn)	<2 ppm		

Lot analysis on the label.



Petroleum ether 30-40°C

NEW PACK.

Certified AR

for analysis, n-hexane <0.1%

packaging	Mfr. No
1 ℓ Glass bottle	P/1720/15
2.5 ℓ Glass bottle	P/1720/17
5 ℓ Mild steel drum	P/1720/MS21

CAS: 64742-49-0
Density at 20°C: 0.62
Boiling point: 30 to 40°C
Melting point: -130°C
Flash point: -49°C
ADR/RID: 3, I
EINECS: 265-151-9

UN: 1268
GHS: Flam. Liq. 1, STOT SE 3, Asp. Tox. 1, Aquatic Chronic 2
H224, H304, H336, H411, EUH066
P210, P273, P301+P310, P331, P261, P271



Guaranteed analysis			
Acidity / Alkalinity	<0.0001 meq/g	Lead (Pb)	<0.05 ppm
Colour	<10 APHA	Magnesium (Mg)	<0.05 ppm
Water	<0.01%	Potassium (K)	<0.2 ppm
Residue after evaporation	<10 ppm	Sodium (Na)	<0.5 ppm
Bromine number	<2	Total phosphorus (P)	<0.5 ppm
Calcium (Ca)	<0.5 ppm	Total silicon (Si)	<0.05 ppm
Copper (Cu)	<0.05 ppm	Total sulphur (S)	<10 ppm
Iron (Fe)	<0.2 ppm	Zinc (Zn)	<0.1 ppm

Typical analysis			
Aluminium (Al)	<0.1 ppm	Molybdenum (Mo)	<0.05 ppm
Barium (Ba)	<0.02 ppm	Nickel (Ni)	<0.02 ppm
Cadmium (Cd)	<0.02 ppm	Selenium (Se)	<0.05 ppm
Chromium (Cr)	<0.02 ppm	Strontium (Sr)	<0.01 ppm
Cobalt (Co)	<0.02 ppm	Tin (Sn)	<0.02 ppm
Lithium (Li)	<0.02 ppm	Titanium (Ti)	<0.01 ppm
Manganese (Mn)	<0.02 ppm	Vanadium (V)	<0.01 ppm

Petroleum ether 30-40°C

NEW PACK.

SLR

specified reagents for general laboratory work- extra pure

packaging

Mfr. No

500 mℓ Glass bottle	P/1400/08
2.5 ℓ Glass bottle	P/1400/17
5 ℓ Mild steel drum	P/1400/MS21
25 ℓ Metal drum	P/1400/25

CAS: 64742-49-0
Density at 20°C: 0.62
Boiling point: 30 to 40°C
Melting point: -130°C
Flash point: -49°C
ADR/RID: 3, I
EINECS: 265-151-9

UN: 1268
GHS: Flam. Liq. 1, STOT SE 3, Asp. Tox. 1, Aquatic Chronic 2
H224, H304, H336, H411, EUH066
P210, P273, P301+P310, P331, P261, P271



Guaranteed analysis			
Acidity / Alkalinity	<0.0005 meq/g	Lead (Pb)	<2 ppm
Colour	<10 APHA	Magnesium (Mg)	<2 ppm
Water	<0.03%	Potassium (K)	<2 ppm
Residue after evaporation	<20 ppm	Sodium (Na)	<5 ppm
Bromine number	<5	Total phosphorus (P)	<1 ppm
Calcium (Ca)	<5 ppm	Total silicon (Si)	<1 ppm
Copper (Cu)	<2 ppm	Total sulphur (S)	<50 ppm
Iron (Fe)	<5 ppm	Zinc (Zn)	<2 ppm

Fisher Chemical 5L packaging change

Our improved HDPE 5L packaging is easier to handle, more robust and preserves product quality for longer.



Petroleum ether 40-60°C

Distol

for residue and pesticide analysis

packaging	Mfr. No
1 ℓ Glass bottle	P/2092/15
2.5 ℓ Glass bottle	P/2092/17

CAS: 64742-49-0
Density at 20°C: 0.648
Boiling point: 36 to 70°C
Melting point: -40°C
Flash point: -45°C
ADR/RID: 3, II
EINECS: 265-151-9

UN: 1268
GHS: Flam. Liq. 2, Skin Irrit. 2, STOT SE 3, Asp. Tox. 1, Aquatic Chronic 2
H225, H315, H411, H304, H336
P210, P240, P260, P262, P273, P301+P330+P331



Guaranteed analysis			
Water	<0.02%	GC (FID) Analysis	Pass
Residue after evaporation	<2 ppm	GC (NPD) analysis	Pass
GC (ECD) Analysis	Pass		

Petroleum ether 40-60°C

extra dry, for synthesis

packaging

Mfr. No

2.5 ℓ Glass bottle	P/2082/17
--------------------	-----------

CAS: 64742-49-0
Density at 20°C: 0.648
Boiling point: 36 to 70°C
Melting point: -40°C
Flash point: -45°C
ADR/RID: 3, II
EINECS: 265-151-9

UN: 1268
GHS: Flam. Liq. 2, Skin Irrit. 2, STOT SE 3, Asp. Tox. 1, Aquatic Chronic 2
H225, H315, H411, H304, H336
P210, P240, P260, P262, P273, P301+P330+P331



Guaranteed analysis			
Water	<0.01%	Residue after evaporation	<2 ppm

Petroleum ether 40-60°C

NEW PACK.

Certified AR

for analysis, n-hexane <2%

packaging

Mfr. No

1 ℓ Glass bottle	P/1760/15
2.5 ℓ Glass bottle	P/1760/17
2.5 ℓ Coated Glass bottle	P/1760/PB17
5 ℓ Fluorinated HDPE Plastic bottle	P/1760/21
25 ℓ Metal drum	P/1760/25
200 ℓ Metal drum	P/1760/27

CAS: 64742-49-0
Density at 20°C: 0.648
Boiling point: 36 to 36 to 70°C
Melting point: -40°C
Flash point: -45°C
ADR/RID: 3, II
EINECS: 265-151-9

UN: 1268
GHS: Flam. Liq. 2, Skin Irrit. 2, STOT SE 3, Asp. Tox. 1, Aquatic Chronic 2
H225, H315, H411, H304, H336
P210, P240, P260, P262, P273, P301+P330+P331



Guaranteed analysis			
Acidity / Alkalinity	<0.0001 meq/g	Magnesium (Mg)	<0.05 ppm
Colour	<10 APHA	Potassium (K)	<0.2 ppm
Water	<0.01%	Sodium (Na)	<0.5 ppm
Residue after evaporation	<10 ppm	Tin (Sn)	<0.02 ppm
Bromine number	<2	Total phosphorus (P)	<0.5 ppm
Calcium (Ca)	<0.5 ppm	Total silicon (Si)	<0.05 ppm
Copper (Cu)	<0.05 ppm	Total sulphur (S)	<10 ppm
Iron (Fe)	<0.2 ppm	Zinc (Zn)	<0.1 ppm
Lead (Pb)	<0.05 ppm		

Typical analysis			
Aluminium (Al)	<0.1 ppm	Molybdenum (Mo)	<0.02 ppm
Barium (Ba)	<0.02 ppm	Nickel (Ni)	<0.02 ppm
Cadmium (Cd)	<0.02 ppm	Selenium (Se)	<0.02 ppm
Chromium (Cr)	<0.02 ppm	Strontium (Sr)	<0.01 ppm
Cobalt (Co)	<0.02 ppm	Titanium (Ti)	<0.01 ppm
Lithium (Li)	<0.02 ppm	Vanadium (V)	<0.01 ppm
Manganese (Mn)	<0.05 ppm		

Petroleum ether 40-60°C

Certified AR

for analysis, low in aromatic

packaging

Mfr. No

2.5 ℓ Glass bottle	P/2000/17
--------------------	-----------

CAS: 64742-49-0
Density at 20°C: 0.648
Boiling point: 36 to 70°C
Melting point: -40°C
Flash point: -45°C
ADR/RID: 3, II
EINECS: 265-151-9

UN: 1268
GHS: Flam. Liq. 2, Skin Irrit. 2, STOT SE 3, Asp. Tox. 1, Aquatic Chronic 2
H225, H315, H411, H304, H336
P210, P240, P260, P262, P273, P301+P330+P331



Guaranteed analysis			
Acidity / Alkalinity	<0.00003 meq/g	Residue after evaporation	<10 ppm
Water	<0.02%	Total sulphur (S)	<5 ppm

Petroleum ether 40-60°C

NEW PACK.

SLR

specified reagents for general laboratory work- extra pure

packaging

Mfr. No

500 mℓ Glass bottle	P/1440/08
1 ℓ Glass bottle	P/1440/15
2.5 ℓ Glass bottle	P/1440/17
2.5 ℓ Coated Glass bottle	P/1440/PB17
5 ℓ Fluorinated HDPE Plastic bottle	P/1440/21
25 ℓ Metal drum	P/1440/25
200 ℓ Metal drum	P/1440/27

CAS: 64742-49-0
Density at 20°C: 0.648
Boiling point: 36 to 70°C
Melting point: -40°C
Flash point: -45°C
ADR/RID: 3, II
EINECS: 265-151-9

UN: 1268
GHS: Flam. Liq. 2, Skin Irrit. 2, STOT SE 3, Asp. Tox. 1, Aquatic Chronic 2
H225, H315, H411, H304, H336
P210, P240, P260, P262, P273, P301+P330+P331



Guaranteed analysis			
Acidity / Alkalinity	<0.0005 meq/g	Lead (Pb)	<2 ppm
Colour	<10 APHA	Magnesium (Mg)	<2 ppm
Water	<0.03%	Potassium (K)	<2 ppm
Residue after evaporation	<20 ppm	Sodium (Na)	<5 ppm
Bromine number	<5	Total phosphorus (P)	<1 ppm
Calcium (Ca)	<5 ppm	Total silicon (Si)	<1 ppm
Copper (Cu)	<2 ppm	Total sulphur (S)	<50 ppm
Iron (Fe)	<5 ppm	Zinc (Zn)	<2 ppm

Optima LC/MS

Fisher Chemical OPTIMA® LC/MS solvents, blends, additives and modifiers have set the standard of excellence for consistent, reproducible performance in the mobile phase of LC/MS.

Find out more on page 308



Petroleum ether 60-80°C

Distol

for residue and pesticide analysis

packaging	Mfr. No
2.5 ℓ Glass bottle	P/2094/17

CAS: 64742-49-0
Density at 20°C: 0.670
Boiling point: 60 to 80°C
Melting point: -40°C
Flash point: -30°C
ADR/RID: 3, II
EINECS: 265-151-9
UN: 1268

GHS: Flam. Liq. 2, Repr. 2, Skin Irrit. 2, STOT SE 3, STOT RE 2, Asp. Tox. 1, Aquatic Chronic 2
H225, H336, H315, H411, H304, H373, H361f
P201, P308+P313, P281, P260, P273, P210



Guaranteed analysis			
Water	<0.02%	GC (FID) Analysis	Conform
Residue after evaporation	<2 ppm	GC (NPD) analysis	Conform
GC (ECD) Analysis	Conform		

Petroleum ether 60-80°C

Certified AR

for analysis, n-hexane <50%

packaging

Mfr. No

1 ℓ Glass bottle	P/1800/15
2.5 ℓ Glass bottle	P/1800/17
2.5 ℓ Coated Glass bottle	P/1800/PB17
25 ℓ Metal drum	P/1800/25
200 ℓ Metal drum	P/1800/27

CAS: 64742-49-0
Density at 20°C: 0.670
Boiling point: 60 to 80°C
Melting point: -40°C
Flash point: -30°C
ADR/RID: 3, II
EINECS: 265-151-9
UN: 1268

GHS: Flam. Liq. 2, Repr. 2, Skin Irrit. 2, STOT SE 3, STOT RE 2, Asp. Tox. 1, Aquatic Chronic 2
H225, H336, H315, H411, H304, H373, H361f
P201, P308+P313, P281, P260, P273, P210



Guaranteed analysis			
Acidity / Alkalinity	<0.0001 meq/g	Magnesium (Mg)	<0.05 ppm
Colour	<10 APHA	Potassium (K)	<0.2 ppm
Water	<0.01%	Sodium (Na)	<0.5 ppm
Residue after evaporation	<10 ppm	Tin (Sn)	<0.02 ppm
Bromine number	<2	Total phosphorus (P)	<0.5 ppm
Calcium (Ca)	<0.5 ppm	Total silicon (Si)	<0.05 ppm
Copper (Cu)	<0.05 ppm	Total sulphur (S)	<10 ppm
Iron (Fe)	<0.2 ppm	Zinc (Zn)	<0.1 ppm
Lead (Pb)	<0.05 ppm		

Typical analysis			
Aluminium (Al)	<0.1 ppm	Molybdenum (Mo)	<0.02 ppm
Barium (Ba)	<0.02 ppm	Nickel (Ni)	<0.02 ppm
Cadmium (Cd)	<0.02 ppm	Selenium (Se)	<0.02 ppm
Chromium (Cr)	<0.02 ppm	Strontium (Sr)	<0.01 ppm
Cobalt (Co)	<0.02 ppm	Titanium (Ti)	<0.01 ppm
Lithium (Li)	<0.02 ppm	Vanadium (V)	<0.01 ppm
Manganese (Mn)	<0.05 ppm		

Standard Volumetric solution

For your convenience we offer a range of NIST traceable solutions. A full list of our products can be found on page 323. Can't see what you need? Please email your contact details and a brief description of your requirements to begei.sha.quotes@thermofisher.com



Petroleum ether 60-80°C		Certified AR
for analysis, low in aromatic		
packaging	Mfr. No	
2.5 ℓ Glass bottle	P/2040/17	
CAS: 64742-49-0	GHS: Flam. Liq. 2, Repr. 2, Skin Irrit. 2, STOT SE 3,	
Density at 20°C: 0.670	STOT RE 2, Asp. Tox. 1,	
Boiling point: 60 to 80°C	Aquatic Chronic 2	
Melting point: -40°C	H225, H336, H315, H411,	
Flash point: -30°C	H304, H373, H361f	
ADR/RID: 3, II	P201, P308+P313, P281,	
EINECS: 265-151-9	P260, P273, P210	
UN: 1268		
Guaranteed analysis		
Acidity / Alkalinity	<0.00003 meq/g	Total sulphur (S)
Water	<0.02%	Aromatic Hydrocarbons
Residue after evaporation	<10 ppm	<5 ppm
		<0.05%

Petroleum ether 60-80°C		SLR
Extra Dry, for synthesis		
packaging	Mfr. No	
500 mℓ Coated Glass bottle	P/2084/PB08	
2.5 ℓ Glass bottle	P/2084/17	
CAS: 64742-49-0	GHS: Flam. Liq. 2, Repr. 2, Skin Irrit. 2, STOT SE 3,	
Density at 20°C: 0.670	STOT RE 2, Asp. Tox. 1,	
Boiling point: 60 to 80°C	Aquatic Chronic 2	
Melting point: -40°C	H225, H336, H315, H411,	
Flash point: -30°C	H304, H373, H361f	
ADR/RID: 3, II	P201, P308+P313, P281,	
EINECS: 265-151-9	P260, P273, P210	
UN: 1268		
Guaranteed analysis		
Water	<0.01%	Residue after evaporation
		<2 ppm

Petroleum ether 60-80°C		NEW PACK.	SLR
specified reagents for general laboratory work- extra pure			
packaging	Mfr. No		
500 mℓ Glass bottle	P/1480/08		
1 ℓ Glass bottle	P/1480/15		
2.5 ℓ Glass bottle	P/1480/17		
5 ℓ Fluorinated HDPE Plastic bottle	P/1480/21		
25 ℓ Metal drum	P/1480/25		
200 ℓ Metal drum	P/1480/27		
CAS: 64742-49-0	GHS: Flam. Liq. 2, Repr. 2, Skin Irrit. 2, STOT SE 3,		
Density at 20°C: 0.670	STOT RE 2, Asp. Tox. 1,		
Boiling point: 60 to 80°C	Aquatic Chronic 2		
Melting point: -40°C	H225, H336, H315, H411,		
Flash point: -30°C	H304, H373, H361f		
ADR/RID: 3, II	P201, P308+P313, P281,		
EINECS: 265-151-9	P260, P273, P210		
UN: 1268			

Guaranteed analysis			
Acidity / Alkalinity	<0.0005 meq/g	Lead (Pb)	<2 ppm
Colour	<10 APHA	Magnesium (Mg)	<2 ppm
Water	<0.03%	Potassium (K)	<2 ppm
Residue after evaporation	<20 ppm	Sodium (Na)	<5 ppm
Bromine number	<5	Total phosphorus (P)	<1 ppm
Calcium (Ca)	<5 ppm	Total silicon (Si)	<1 ppm
Copper (Cu)	<2 ppm	Total sulphur (S)	<50 ppm
Iron (Fe)	<5 ppm	Zinc (Zn)	<2 ppm

Petroleum ether 80-100°C		Certified AR
for analysis, n-hexane <20%		
packaging	Mfr. No	
2.5 ℓ Glass bottle	P/1840/17	
25 ℓ Metal drum	P/1840/25	
200 ℓ Metal drum	P/1840/27	
CAS: 64742-49-0	GHS: Flam. Liq. 2, Skin Irrit. 2, STOT SE 3, Asp. Tox. 1,	
Density at 20°C: 0.700	Aquatic Acute 1, Aquatic Chronic 1	
Boiling point: 80 to 100°C	H225, H304, H336, H315, H410	
Melting point: -40°C	P201, P308+P313,	
Flash point: -20°C	P301+P310, P331, P210, P240	
ADR/RID: 3, II		
EINECS: 265-151-9		
UN: 1268		
Guaranteed analysis		
Acidity / Alkalinity	<0.0001 meq/g	Magnesium (Mg)
Colour	<10 APHA	Potassium (K)
Water	<0.01%	Sodium (Na)
Residue after evaporation	<10 ppm	Total phosphorus (P)
Bromine number	<2	Total silicon (Si)
Calcium (Ca)	<0.5 ppm	Total sulphur (S)
Copper (Cu)	<0.05 ppm	Vanadium (V)
Iron (Fe)	<0.2 ppm	Zinc (Zn)
Lead (Pb)	<0.05 ppm	

Typical analysis			
Aluminium (Al)	<0.1 ppm	Molybdenum (Mo)	<0.02 ppm
Barium (Ba)	<0.02 ppm	Nickel (Ni)	<0.02 ppm
Cadmium (Cd)	<0.02 ppm	Selenium (Se)	<0.05 ppm
Chromium (Cr)	<0.02 ppm	Strontium (Sr)	<0.01 ppm
Cobalt (Co)	<0.02 ppm	Tin (Sn)	<0.02 ppm
Lithium (Li)	<0.02 ppm	Titanium (Ti)	<0.01 ppm
Manganese (Mn)	<0.05 ppm		

Lot analysis on the label.

Petroleum ether 80-100°C		NEW PACK.	SLR
specified reagents for general laboratory work- extra pure			
packaging	Mfr. No		
500 mℓ Glass bottle	P/1520/08		
2.5 ℓ Glass bottle	P/1520/17		
5 ℓ Fluorinated HDPE Plastic bottle	P/1520/21		
25 ℓ Metal drum	P/1520/25		
CAS: 64742-49-0	GHS: Flam. Liq. 2, Skin Irrit. 2, STOT SE 3, Asp. Tox. 1,		
Density at 20°C: 0.700	Aquatic Acute 1, Aquatic Chronic 1		
Boiling point: 80 to 100°C	H225, H304, H336, H315, H410		
Melting point: -40°C	P201, P308+P313,		
Flash point: -20°C	P301+P310, P331, P210, P240		
ADR/RID: 3, II			
EINECS: 265-151-9			
UN: 1268			

Guaranteed analysis			
Acidity / Alkalinity	<0.0005 meq/g	Lead (Pb)	<2 ppm
Colour	<10 APHA	Magnesium (Mg)	<2 ppm
Water	<0.03%	Potassium (K)	<2 ppm
Residue after evaporation	<20 ppm	Sodium (Na)	<5 ppm
Bromine number	<5	Total phosphorus (P)	<1 ppm
Calcium (Ca)	<5 ppm	Total silicon (Si)	<1 ppm
Copper (Cu)	<2 ppm	Total sulphur (S)	<50 ppm
Iron (Fe)	<5 ppm	Zinc (Zn)	<2 ppm

Petroleum ether 100-120°C		Certified AR
for analysis, n-hexane <2%		
packaging	Mfr. No	
2.5 ℓ Glass bottle	P/1880/17	
25 ℓ Metal drum	P/1880/25	
200 ℓ Metal drum	P/1880/27	
CAS: 64742-49-0	GHS: Flam. Liq. 2, Skin Irrit. 2, STOT SE 3, Asp. Tox. 1,	
Density at 20°C: 0.735	Aquatic Acute 1, Aquatic Chronic 1	
Boiling point: 100 to 120°C	H225, H304, H336, H315, H410	
Melting point: -30°C	P201, P308+P313,	
Flash point: 2°C	P301+P310, P331, P210, P240	
ADR/RID: 3, II		
EINECS: 265-151-9		
UN: 1268		
Guaranteed analysis		
Acidity / Alkalinity	<0.0001 meq/g	Lead (Pb)
Colour	<10 APHA	Magnesium (Mg)
Water	<0.01%	Potassium (K)
Residue after evaporation	<10 ppm	Sodium (Na)
Bromine number	<2	Total phosphorus (P)
Calcium (Ca)	<0.5 ppm	Total silicon (Si)
Copper (Cu)	<0.05 ppm	Total sulphur (S)
Iron (Fe)	<0.2 ppm	Zinc (Zn)

Typical analysis			
Aluminium (Al)	<0.1 ppm	Molybdenum (Mo)	<0.02 ppm
Barium (Ba)	<0.02 ppm	Nickel (Ni)	<0.02 ppm
Cadmium (Cd)	<0.02 ppm	Selenium (Se)	<0.05 ppm
Chromium (Cr)	<0.02 ppm	Strontium (Sr)	<0.01 ppm
Cobalt (Co)	<0.02 ppm	Tin (Sn)	<0.02 ppm
Lithium (Li)	<0.02 ppm	Titanium (Ti)	<0.01 ppm
Manganese (Mn)	<0.05 ppm	Vanadium (V)	<0.01 ppm

Lot analysis on the label.

Petroleum ether 100-120°C		SLR
specified reagents for general laboratory work- extra pure		
packaging	Mfr. No	
500 mℓ Glass bottle	P/1560/08	
1 ℓ Glass bottle	P/1560/15	
2.5 ℓ Glass bottle	P/1560/17	
25 ℓ Metal drum	P/1560/25	
CAS: 64742-49-0	GHS: Flam. Liq. 2, Skin Irrit. 2, STOT SE 3, Asp. Tox. 1,	
Density at 20°C: 0.735	Aquatic Acute 1, Aquatic Chronic 1	
Boiling point: 100 to 120°C	H225, H304, H336, H315, H410	
Melting point: -30°C	P201, P308+P313,	
Flash point: 2°C	P301+P310, P331, P210, P240	
ADR/RID: 3, II		
EINECS: 265-151-9		
UN: 1268		

Guaranteed analysis			
Acidity / Alkalinity	<0.0005 meq/g	Lead (Pb)	<2 ppm
Colour	<10 APHA	Magnesium (Mg)	<2 ppm
Water	<0.03%	Potassium (K)	<2 ppm
Residue after evaporation	<20 ppm	Sodium (Na)	<5 ppm
Bromine number	<5	Total phosphorus (P)	<1 ppm
Calcium (Ca)	<5 ppm	Total silicon (Si)	<1 ppm
Copper (Cu)	<2 ppm	Total sulphur (S)	<50 ppm
Iron (Fe)	<5 ppm	Zinc (Zn)	<2 ppm

Petroleum ether 100-120°C		
technical, water determination in petroleum products		
packaging	Mfr. No	
2.5 ℓ Coated Glass bottle	P/2100/PB17	
CAS: 64742-49-0	GHS: Flam. Liq. 2, Skin Irrit. 2, STOT SE 3, Asp. Tox. 1,	
Density at 20°C: 0.735	Aquatic Acute 1, Aquatic Chronic 1	
Boiling point: 100 to 120°C	H225, H304, H336, H315, H410	
Flash point: 2°C	P201, P308+P313,	
ADR/RID: 3, II	P301+P310, P331, P210, P240	
EINECS: 265-151-9		
UN: 1268		
Guaranteed analysis		
Residue after evaporation	<20 ppm	

Petroleum ether > 120°C		SLR
specified reagents for general laboratory work- extra pure		
packaging	Mfr. No	
2.5 ℓ Glass bottle	P/1600/17	
CAS: 64742-49-0	GHS: Flam. Liq. 2, Skin Irrit. 2, STOT SE 3, Asp. Tox. 1,	
Density at 20°C: 0.755	Aquatic Acute 1, Aquatic Chronic 1	
Boiling point: 120°C	H225, H304, H336, H315, H410	
Melting point: -40°C	P210, P240, P260, P262,	
Flash point: 8°C	P273, P301+P330+P331	
ADR/RID: 3, II		
EINECS: 265-151-9		
UN: 1268		

Guaranteed analysis			
Acidity / Alkalinity	<0.0005 meq/g	Lead (Pb)	<2 ppm
Colour	<10 APHA	Magnesium (Mg)	<2 ppm
Water	<0.03%	Potassium (K)	<2 ppm
Residue after evaporation	<20 ppm	Sodium (Na)	<5 ppm
Bromine number	<5	Total phosphorus (P)	<1 ppm
Calcium (Ca)	<5 ppm	Total silicon (Si)	<1 ppm
Copper (Cu)	<2 ppm	Total sulphur (S)	<50 ppm
Iron (Fe)	<5 ppm	Zinc (Zn)	<2 ppm

Petroleum ether (petroleum spirit with <15% n-hexane)		
pure		
packaging	Mfr. No	
2.5 ℓ Glass bottle	P/1360/17	
25 ℓ Metal drum	P/1360/25	
CAS: 8032-32-4	GHS: Flam. Liq. 2, Repr. 2, Skin Irrit. 2, STOT SE 3,	
Density at 20°C: 0.702	STOT RE 2, Asp. Tox. 1,	
Boiling point: 80 to 110°C	Aquatic Chronic 2	
Melting point: -40°C	H225, H315, H411, H304,	
Flash point: -20°C	H336, H373, H361f	
ADR/RID: 3, II	P201, P308+P313, P281,	
EINECS: 232-453-7	P210, P273, P260	
UN: 1268		

Reagents, Solvents & Volumetric solutions conform to the European Pharmacopeia (Eur.Ph), British (BP) & US (USP) **NEW PACK.** **NEW**

designation	EP	BP	USP	packsize	package type	Mfr. No	page
Acetic add glacial				500 mℓ	HDPE Plastic bottle	A/0400/PB08	15
				1 L	Coated Glass bottle	A/0400/PC15	15
				2.5 L	HDPE Plastic bottle	A/0400/PB17	15
					Coated Glass bottle	A/0400/PC17	15
				25 L	HDPE Plastic container	A/0400/25	15
				45 L	Plastic container	A/0400/26	15
Acetic acid glacial 99-100%				1 L	HDPE Plastic bottle	A/0360/PB15	15
				2.5 L	HDPE Plastic bottle	A/0360/PB17	15
				25 L	HDPE Plastic container	A/0360/25	15
Acetone				1 L	Glass bottle	A/0600/15	17
				2.5 L	Glass bottle	A/0600/17	17
					HDPE Plastic bottle	A/0600/PC17	17
				5 L	Fluorinated HDPE Plastic bottle	A/0600/21	17
					HDPE Plastic bottle	A/0600/PC21	17
				25 L	Metal drum	A/0600/25	17
				200 L	Metal drum	A/0600/27	17
Acetonitrile HPLC gradient grade				1 L	Glass bottle	A/0627/15	19
				2.5 L	Glass bottle	A/0627/17	19
					Coated Glass bottle	A/0627/PB17	19
				5 L	Aluminium can	A/0627/21	19
				10 L	Returnable Drum	A/0627/10RSS	19
				25 L	Returnable Drum	A/0627/25RSS	19
				30 L	Returnable Drum	A/0627/30RSS	19
Aluminium chloride hexahydrate				500 g	Plastic jar	A/1960/53	29
Ammonium chloride				500 g	Plastic jar	A/3920/53	39
				1 kg	Plastic jar	A/3920/60	39
				2 kg	Plastic jar	A/3920/61	39
				5 kg	Plastic keg	A/3920/65	39
				25 kg	Fibre Keg	A/3920/70	39
Ammonium thiocyanate solution 0.1M (0.1N)				1 L	HDPE Plastic bottle	J/1500/PB15	46
				2.5 L	HDPE Plastic bottle	J/1500/PB17	46
Benzyl alcohol				500 mℓ	Glass bottle	B/2550/08	57
				2.5 L	Coated Glass bottle	B/2550/PB17	57
Boric acid				500 g	Plastic jar	B/3750/53	59
				1 kg	Plastic jar	B/3750/60	59
				5 kg	Plastic keg	B/3750/65	59
Boric acid				500 g	Plastic jar	B/3760/53	59
				1 kg	Plastic jar	B/3760/60	59
				5 kg	Plastic keg	B/3760/65	59
				25 kg	Fibre Keg	B/3760/70	59
Calcium carbonate				1 kg	Plastic jar	C/1040/60	75
				3 kg	Plastic jar	C/1040/63	75
				5 kg	Plastic keg	C/1040/65	75
Calcium gluconate				500 g	Plastic jar	C/1610/53	75
Calcium lactate pentahydrate				500 g	Plastic jar	C/1820/53	76
				5 kg	Plastic keg	C/1820/65	76
Cerium(IV) sulfate solution 0.1M (0.1N)				1 L	Glass bottle	J/3000/15	81
				2.5 L	Glass bottle	J/3000/17	81
Chloral hydrate				500 g	Plastic jar	C/4280/53	83
				10 kg	Plastic keg	C/4280/68	83
Chloramine-T Sodium salt				250 g	Plastic jar	C/4320/50	83
				500 g	Plastic jar	C/4320/53	83
Citric acid, anhydrous				500 g	Plastic jar	C/6240/53	90
Citric acid monohydrate				500 g	Plastic jar	C/6200/53	90
					Plastic jar	C/6200/53	90
					Plastic jar	C/6200/60	90
				1 kg	Plastic jar	C/6200/60	90
					Plastic jar	C/6200/62	90
				2.5 kg	Plastic jar	C/6200/62	90
					Plastic jar	C/6200/62	90
				5 kg	Plastic keg	C/6200/65	90
					Plastic keg	C/6200/65	90
				25 kg	Fibre Keg	C/6200/70	90
					Fibre Keg	C/6200/70	90
					Fibre Keg	C/6200/70	90

Reagents, Solvents & Volumetric solutions conform to the European Pharmacopeia (Eur.Ph), British (BP) & US (USP) **NEW PACK.** **NEW**

designation	EP	BP	USP	packsize	package type	Mfr. No	page
Diethyl ether stabilized with BHT				500 mℓ	Coated Glass bottle	D/2450/PB08	110
					Coated Glass bottle	D/2450/PB08	110
				1 L	Glass bottle	D/2450/15	110
				2.5 L	Glass bottle	D/2450/17	110
				5 L	Mild steel drum	D/2450/MS21	110
				25 L	Metal drum	D/2450/25	110
				200 L	Metal drum	D/2450/27	110
Ethanol absolute				500 mℓ	Glass bottle	E/0650DF/08	121
				1 L	Glass bottle	E/0650DF/15	121
				205 L	Metal drum	E/0650DF/28	121
				2.5 L	Glass bottle	E/0650DF/17	121
					HDPE Plastic bottle	E/0650DF/C17	121
				5 L	Fluorinated HDPE Plastic bottle	E/0650DF/21	121
				10 L	HDPE Plastic container	E/0650DF/24	121
Ethanol				25 L	Metal drum	E/0650DF/25	121
				500 mℓ	Glass bottle	E/0550/08	121
				1 L	Glass bottle	E/0550DF/15	121
				2.5 L	Glass bottle	E/0550DF/17	121
					Glass bottle	E/0550DF/17	121
				5 L	Fluorinated HDPE Plastic bottle	E/0550DF/21	121
				25 L	Metal drum	E/0550DF/25	121
Ethanol 96% v/v				2.5 L	HDPE Plastic bottle	E/0555/PC17	122
					Glass bottle	E/0555DF/17	122
				5 L	Fluorinated HDPE Plastic bottle	E/0555DF/21	122
				25 L	Metal drum	E/0555DF/25	122
Ethylenediaminetetraacetic acid				100 g	Plastic jar	D/0700/48	126
				250 g	Plastic jar	D/0700/50	126
				500 g	Plastic jar	D/0700/53	126
				1 kg	Plastic jar	D/0700/60	126
				10 kg	Plastic keg	D/0700/68	126
				25 kg	Fibre Keg	D/0700/70	126
Ethylenediaminetetraacetic acid disodium salt solution 0.2N				1 L	HDPE Plastic bottle	J/3700/15	127
				2.5 L	HDPE Plastic bottle	J/3700/17	127
				10 L	Bag in box	J/3700/24	127
Gelatine				500 g	Plastic jar	G/0150/53	134
D(+)-Glucose anhydrous				500 g	Plastic jar	G/0500/53	135
				1 kg	Plastic jar	G/0500/60	135
				2 kg	Plastic jar	G/0500/61	135
				5 kg	Plastic keg	G/0500/65	135
				25 kg	Fibre Keg	G/0500/70	135
D(+)-Glucose monohydrate				1 kg	Plastic jar	G/0400/60	135
				5 kg	Plastic keg	G/0400/65	135
Hydrochloric acid solution 0.1M (0.1N)				1 L	HDPE Plastic bottle	J/4350/15	150
				2.5 L	HDPE Plastic bottle	J/4350/17	150
				5 L	HDPE Plastic bottle	J/4350/21	150
				10 L	Bag in box	J/4350/24	150
Hydrochloric acid solution 1M (1N)				1 L	HDPE Plastic bottle	J/4320/15	151
				2.5 L	HDPE Plastic bottle	J/4320/17	151
				5 L	HDPE Plastic bottle	J/4320/21	151
				10 L	Bag in box	J/4320/24	151
Hydrochloric acid solution 2M (2N)				1 L	HDPE Plastic bottle	J/4315/15	151
				2.5 L	HDPE Plastic bottle	J/4315/17	151
				10 L	Bag in box	J/4315/24	151
Iodine solution 0.05M (0.1N)				1 L	Glass bottle	J/4410/15	161
				2.5 L	Glass bottle	J/4410/17	161
Iodine solution 0.5M (1N)				500 mℓ	Coated Glass bottle	J/4400/PB08	161
				1 L	Coated Glass bottle	J/4400/PB15	161
				2.5 L	Coated Glass bottle	J/4400/PB17	161
Lactic acid				500 mℓ	HDPE Plastic bottle	L/0150/PB08	176
				2.5 L	HDPE Plastic bottle	L/0150/PB17	176
Magnesium carbonate pentahydrate				500 g	Plastic jar	M/0500/53	188
Magnesium carbonate pentahydrate				25 kg	Plastic keg	M/0500/70	188
				500 g	Plastic jar	M/0450/53	188

Reagents, Solvents & Volumetric solutions conform to the European Pharmacopeia (Eur.Ph), British (BP) & US (USP) NEW PACK. NEW

designation	EP	BP	USP	packsize	package type	Mfr. No	page
Magnesium chloride hexahydrate				500 g	Plastic jar	M/0600/53	188
				1 kg	Plastic jar	M/0600/60	188
				2 kg	Plastic jar	M/0600/61	188
				5 kg	Plastic keg	M/0600/65	188
Magnesium oxide				500 g	Plastic jar	M/0850/53	189
				1 kg	Plastic keg	M/0850/60	189
Magnesium oxide				500 g	Plastic jar	M/0800/53	189
Magnesium stearate				500 g	Plastic jar	M/0955/53	189
Magnesium sulfate heptahydrate				1 kg	Plastic jar	M/1000/60	190
				2.5 kg	Plastic jar	M/1000/62	190
				5 kg	Plastic keg	M/1000/65	190
				500 g	Plastic jar	M/2405/53	193
Nitric acid solution 1M (1N)				1 L	HDPE Plastic bottle	J/5550/PB15	218
				2.5 L	HDPE Plastic bottle	J/5550/PB17	218
				10 L	Bag in box	J/5550/24	218
Phenol detached crystals				500 g	Plastic jar	P/2360/53	240
				1 kg	Plastic jar	P/2360/60	240
				2 kg	Plastic jar	P/2360/61	240
				25 kg	Fibre Keg	P/2360/70	240
Potassium carbonate anhydrous				250 g	Plastic jar	P/4120/50	248
				500 g	Plastic jar	P/4120/53	248
				1 kg	Plastic jar	P/4120/60	248
				2.5 kg	Plastic jar	P/4120/62	248
Potassium chloride				500 g	Plastic jar	P/4240/53	249
				1 kg	Plastic jar	P/4240/60	249
				2.5 kg	Plastic jar	P/4240/62	249
				5 kg	Plastic jar	P/4240/65	249
tri-Potassium citrate				500 g	Plastic jar	P/4442/53	250
				3 kg	Plastic jar	P/4442/63	250
Potassium dichromate solution 0.0167M (0.1N)				1 L	Glass bottle	J/6212/15	251
				2.5 L	Glass bottle	J/6212/17	251
Potassium hydrogen carbonate				500 g	Plastic jar	P/5080/53	254
				2.5 kg	Plastic jar	P/5080/62	254
Potassium hydroxide				500 g	Plastic jar	P/5640/53	256
				1 kg	Plastic jar	P/5640/60	256
				5 kg	Plastic jar	P/5640/65	256
				10 kg	Plastic keg	P/5640/68	256
				25 kg	Plastic keg	P/5640/70	256
				1 L	HDPE Plastic bottle	J/6625/15	256
Potassium hydroxide solution 0.1M (0.1N), (in methanol)				1 L	HDPE Plastic bottle	J/6630/PB15	257
				2.5 L	HDPE Plastic bottle	J/6630/17	257
Potassium iodate 0.05M solution stabilised with 10 ppm mercury(II) chloride				1 L	Glass bottle	J/6690/15	258
				2.5 L	Glass bottle	J/6690/17	258
Potassium nitrate				250 g	Plastic jar	P/6120/50	259
				500 g	Plastic jar	P/6120/53	259
				1 kg	Plastic jar	P/6120/60	259
				5 kg	Plastic jar	P/6120/65	259
Potassium permanganate				250 g	Glass jar	P/6560/50	261
				500 g	Glass jar	P/6560/53	261
				5 kg	Plastic keg	P/6560/65	261
				25 kg	Plastic keg	P/6560/70	261
Propane-1,2-diol				500 mL	Glass bottle	P/7440/08	264
				2.5 L	Glass bottle	P/7440/17	264
				25 L	Metal drum	P/7440/25	264
				200 L	Metal drum	P/7440/27	264
Silver nitrate solution 0.1M (0.1N)				1 L	HDPE Plastic bottle	J/7330/15	277
				2.5 L	HDPE Plastic bottle	J/7330/17	277
				10 L	Bag in box	J/7330/24	277
Sodium benzoate				500 g	Plastic jar	S/2440/53	283
				1 kg	Plastic jar	S/2440/60	283
				25 kg	Plastic keg	S/2440/70	283

Reagents, Solvents & Volumetric solutions conform to the European Pharmacopeia (Eur.Ph), British (BP) & US (USP) NEW PACK. NEW

designation	EP	BP	USP	packsize	package type	Mfr. No	page
Sodium chloride 99.5%				500 g	Plastic jar	S/3161/53	285
				1 kg	Plastic jar	S/3161/60	285
				3 kg	Plastic jar	S/3161/63	285
				5 kg	Plastic jar	S/3161/65	285
				10 kg	Plastic keg	S/3161/68	285
				25 kg	Plastic keg	S/3161/70	285
Sodium chloride				500 g	Plastic jar	S/3160/53	286
				1 kg	Plastic jar	S/3160/60	286
				3 kg	Plastic jar	S/3160/63	286
				5 kg	Plastic jar	S/3160/65	286
				10 kg	Plastic keg	S/3160/68	286
				25 kg	Plastic keg	S/3160/70	286
tri-Sodium citrate dihydrate				500 g	Plastic jar	S/3320/53	287
				1 kg	Plastic jar	S/3320/60	287
				5 kg	Plastic keg	S/3320/65	287
				25 kg	Fibre Keg	S/3320/70	287
Sodium dihydrogen orthophosphate dihydrate				500 g	Plastic jar	S/3760/53	289
				1 kg	Plastic jar	S/3760/60	289
				5 kg	Plastic keg	S/3760/65	289
				25 kg	Plastic keg	S/3760/70	289
Sodium hydrogen carbonate				500 g	Plastic jar	S/4240/53	291
				1 kg	Plastic jar	S/4240/60	291
				3 kg	Plastic jar	S/4240/63	291
				5 kg	Plastic keg	S/4240/65	291
di-Sodium hydrogen orthophosphate anhydrous				500 g	Plastic jar	S/4520/53	292
				1 kg	Plastic jar	S/4520/60	292
				3 kg	Plastic jar	S/4520/63	292
				5 kg	Plastic keg	S/4520/65	292
				25 kg	Fibre Keg	S/4520/70	292
				500 g	Plastic jar	S/4400/53	293
di-Sodium hydrogen orthophosphate dodecahydrate				1 kg	Plastic jar	S/4400/60	293
				10 kg	Plastic keg	S/4400/68	293
				500 g	Plastic jar	S/4920/53	294
Sodium hydroxide pellets				1 kg	Plastic jar	S/4920/60	294
				3 kg	Plastic jar	S/4920/63	294
				5 kg	Plastic jar	S/4920/65	294
				10 kg	Plastic keg	S/4920/68	294
				25 kg	Plastic keg	S/4920/70	294
				1 L	HDPE Plastic bottle	J/7610/PB15	296
Sodium hydroxide solution 2M (2N)				2.5 L	HDPE Plastic bottle	J/7610/17	296
				5 L	HDPE Plastic bottle	J/7610/21	296
				1 L	HDPE Plastic bottle	J/7620/15	296
Sodium hydroxide solution 1M (1N)				2.5 L	HDPE Plastic bottle	J/7620/17	296
				5 L	HDPE Plastic bottle	J/7620/21	296
				10 L	Bag in box	J/7620/24	296
				1 L	HDPE Plastic bottle	J/7660/15	297
Sodium hydroxide solution 0.1M (0.1N)				2.5 L	HDPE Plastic bottle	J/7660/17	297
				5 L	HDPE Plastic bottle	J/7660/21	297
				10 L	Bag in box	J/7660/24	297
Sodium metabisulfite				500 g	Plastic jar	S/5240/53	298
				1 kg	Plastic jar	S/5240/60	298
				2.5 kg	Plastic jar	S/5240/62	298
di-Sodium tetraborate decahydrate				1 kg	Plastic jar	S/7010/60	305
				3 kg	Plastic jar	S/7010/63	305
				25 kg	Plastic keg	S/7010/70	305
				1 L	HDPE Plastic bottle	J/7950/15	306
Sodium thiosulfate solution 0.1M (0.1N)				2.5 L	HDPE Plastic bottle	J/7950/17	306
				5 L	HDPE Plastic bottle	J/7950/21	306
				10 L	Bag in box	J/7950/24	306
				500 g	Plastic jar	S/7400/53	313
Sorbitol				3 kg	Plastic jar	S/7400/63	313
				25 kg	Plastic keg	S/7400/70	313
				500 g	Plastic jar	S/8600/53	334
Sucrose				1 kg	Plastic jar	S/8600/60	334
				3 kg	Plastic jar	S/8600/63	334

Reagents, Solvents & Volumetric solutions conform to the European Pharmacopeia (Eur.Ph), British (BP) & US (USP) NEW PACK. NEW

designation	EP	BP	USP	packsize	package type	Mfr. No	page
Sucrose				5 kg	Plastic keg	S/8600/65	334
				10 kg	Plastic keg	S/8600/68	334
				25 kg	Fibre Keg	S/8600/70	334
Sulfuric acid solution 0.5M (1N)				1 L	HDPE Plastic bottle	J/8430/15	339
				2.5 L	HDPE Plastic bottle	J/8430/17	339
				10 L	Bag in box	J/8430/24	339
Sulfuric acid solution 0.05M (0.1N)				1 L	HDPE Plastic bottle	J/8460/15	340
					HDPE Plastic bottle	J/8460/15	340
				2.5 L	HDPE Plastic bottle	J/8460/17	340
					HDPE Plastic bottle	J/8460/17	340
				10 L	Bag in box	J/8460/24	340
					Bag in box	J/8460/24	340
				25 L	HDPE Drum	J/8460/25	340
(+) - Tartaric acid				500 g	Plastic jar	T/0250/53	340
				1 kg	Plastic jar	T/0250/60	340
Zinc oxide				500 g	Plastic jar	Z/1300/53	376
				1 kg	Plastic jar	Z/1300/60	376

1.10-Phenanthroline monohydrate

Certified AR

for analysis

packaging	Mfr. No
1 g Vial	P/2200/43
5 g Glass bottle	P/2200/44
100 g Plastic jar	P/2200/48
C ₁₂ H ₈ N ₂ ·H ₂ O	UN: 2811
CAS: 5144-89-8	GHS: Acute Tox. 3, Aquatic
MW: 198.23	Acute 1, Aquatic Chronic 1
Melting point: 97 to 101°C	H301, H410
ADR/RID: 6.1, III	P301+P310, P273

Assay: >99.5%

Guaranteed analysis

Residue after ignition <0.05% Sensitivity test (iron) : Pass

Indicator Redox

Lot analysis on the label.

1.10-Phenanthroline monohydrate

pure

packaging	Mfr. No
5 g Plastic jar	P/2190/44
C ₁₂ H ₈ N ₂ ·H ₂ O	UN: 2811
CAS: 5144-89-8	GHS: Acute Tox. 3, Aquatic
MW: 198.23	Acute 1, Aquatic Chronic 1
Melting point: 97 to 101°C	H301, H410
ADR/RID: 6.1, III	P301+P310, P273

Assay: >98%

1.10-Phenanthroline ferrous solution 0.025M

pure

packaging	Mfr. No
100 mL Coated Glass bottle	P/2220/PB05
C ₁₂ H ₁₄ FeNSO ₄	Density at 20°C: 0.9999
CAS: 14634-91-4	Melting point: 0°C
MW: 692.53	EINECS: 238-676-6

Guaranteed analysis

Identity test : Pass

Suitable for use in the Chemical Oxygen Demand test (DCO).

2-Phenylethanol

SLR

specified reagents for general laboratory work- extra pure

packaging	Mfr. No
500 mL Glass bottle	P/2526/08
C ₈ H ₁₀ O	EINECS: 200-456-2
CAS: 60-12-8	UN: 2810
MW: 122.17	GHS: Acute Tox. 3, Skin Irrit.
Density at 20°C: 1.02	2, Eye Irrit. 2
Boiling point: 219°C	H302, H311, H319, H315
Melting point: -27°C	P280, P302+P352,
Flash point: 102°C	P301+P312, P305+P351+P338
ADR/RID: 6.1, III	

Assay (CG): >99%

2-Phenylethylamine

SLR

specified reagents for general laboratory work- extra pure

packaging	Mfr. No
100 mL Coated Glass bottle	P/2530/PB05
C ₈ H ₁₁ N	EINECS: 200-574-4
CAS: 64-04-0	UN: 2735
MW: 121.18	GHS: Acute Tox. 4, Skin
Density at 20°C: 0.962	Corr. 1B, Eye Dam. 1
Boiling point: 194 to 195°C	H302, H314
Melting point: -60°C	P301+P330+P331, P280,
Flash point: 90°C	P305+P351+P338, P310,
ADR/RID: 8, I	P301+P312

Assay: >98.5%

Phenol liquified 80% w/w in water

Certified AR

for analysis

packaging	Mfr. No
2.5 L Coated Glass bottle	P/2316/PB17
C ₆ H ₆ O	GHS: Muta. 2, Acute Tox. 3,
CAS: 108-95-2	Skin Corr. 1B, Eye Dam. 1,
MW: 94.11	STOT RE 2
Density at 20°C: 1.06	H341, H311, H314, H373,
Boiling point: 182°C	H331, H301
Melting point: 40.5°C	P280, P312, P302+P350,
Flash point: 79°C	P301+P330+P331,
ADR/RID: 6.1, II	P305+P351+P338, P310,
EINECS: 203-632-7	P281, P260, P304+P340
UN: 2821	

Assay: 77 to 81.5%

Guaranteed analysis

Acidity / Alkalinity <0.0005 meq/g Sodium (Na) <10 ppm

Other organic impurities test : Pass Total chloride (Cl) <0.0005%

Calcium (Ca) <5 ppm Total phosphorus (P) <20 ppm

Copper (Cu) <0.5 ppm Total silicon (Si) <5 ppm

Iron (Fe) <1 ppm Total sulphur (S) <50 ppm

Lead (Pb) <0.5 ppm Zinc (Zn) <0.5 ppm

Magnesium (Mg) <0.5 ppm Others organic impurities No reaction

Potassium (K) <5 ppm

Lot analysis on the label.

Phenol liquified 80% w/w in water

pure

packaging	Mfr. No
2.5 L Coated Glass bottle	P/2315/PB17
C ₆ H ₆ O	GHS: Muta. 2, Acute Tox. 3,
CAS: 108-95-2	Skin Corr. 1B, Eye Dam. 1,
MW: 94.11	STOT RE 2
Density at 20°C: 1.06	H341, H311, H314, H373,
Boiling point: 182°C	H331, H301
Melting point: 40.5°C	P280, P312, P302+P350,
Flash point: 79°C	P301+P330+P331,
ADR/RID: 6.1, II	P305+P351+P338, P310,
EINECS: 203-632-7	P281, P260, P304+P340
UN: 2821	



Phenol detached crystals
for analysis

Certified AR, Eur.Ph.,BP,USP

packaging	Mfr. No
500 g Plastic jar	P/2360/53
1 kg Plastic jar	P/2360/60
2 kg Plastic jar	P/2360/61
25 kg Fibre Keg	P/2360/70
C ₆ H ₆ O	UN: 1671
CAS: 108-95-2	GHS: Muta. 2, Acute Tox. 3,
MW: 94.11	Skin Corr. 1B, STOT RE 2
Density at 20°C: 1.07	H311, H331, H341, H314,
Boiling point: 182°C	H373, H301
Melting point: 39 to 42°C	P280, P302+P350,
Flash point: 79°C	P301+P330+P331,
ADR/RID: 6.1, II	P305+P351+P338,
EINECS: 203-632-7	P304+P340, P301+P310, P260



Assay (CG): >99%

Guaranteed analysis			
Acidity / Alkalinity	<0.0005 meq/g	Sodium (Na)	<10 ppm
Water insoluble matter	<0.005%	Total chloride (Cl)	<0.0005%
Calcium (Ca)	<5 ppm	Total phosphorus (P)	<20 ppm
Copper (Cu)	<0.5 ppm	Total silicon (Si)	<5 ppm
Iron (Fe)	<1 ppm	Total sulphur (S)	<50 ppm
Lead (Pb)	<0.5 ppm	Zinc (Zn)	<0.5 ppm
Magnesium (Mg)	<0.5 ppm	Others organic impurities	No reaction
Potassium (K)	<5 ppm		

Lot analysis on the label.

Phenol detached crystals

SLR

specified reagents for general laboratory work- extra pure

packaging	Mfr. No
500 g Plastic jar	P/2320/53
1 kg Plastic jar	P/2320/60
C ₆ H ₆ O	UN: 1671
CAS: 108-95-2	GHS: Muta. 2, Acute Tox. 3,
MW: 94.11	Skin Corr. 1B, STOT RE 2
Density at 20°C: 1.07	H311, H331, H341, H314,
Boiling point: 182°C	H373, H301
Melting point: 39 to 42°C	P280, P302+P350,
Flash point: 79°C	P301+P330+P331,
ADR/RID: 6.1, II	P305+P351+P338,
EINECS: 203-632-7	P304+P340, P301+P310, P260



Assay (CG): >98%

Guaranteed analysis			
Calcium (Ca)	<20 ppm	Sodium (Na)	<50 ppm
Copper (Cu)	<2 ppm	n-Butyl formate	<0.005%
Iron (Fe)	<5 ppm	Total phosphorus (P)	<100 ppm
Lead (Pb)	<2 ppm	Total silicon (Si)	<20 ppm
Magnesium (Mg)	<2 ppm	Total sulphur (S)	<200 ppm
Potassium (K)	<20 ppm	Zinc (Zn)	<2 ppm

Phenol red

pure, indicator grade

packaging	Mfr. No
5 g Glass jar	P/2420/44
25 g Glass jar	P/2420/46
500 g Metal can	P/2420/53
C ₁₉ H ₁₄ O ₅ S	GHS: Skin Irrit. 2, Eye Irrit. 2,
CAS: 143-74-8	STOT SE 3
MW: 354.38	H319, H315, H335
Melting point: >300°C	P261, P302+P352, P280,
EINECS: 205-609-7	P305+P351+P338



Guaranteed analysis			
pH colour change test	: Pass	Complies with standard IR spectrum:	Pass
pH 6.8-8.4 yellow to red.			

Phenol red solution 0.02% in water

pure, indicator grade

packaging	Mfr. No
500 mℓ Glass bottle	P/2420L/08
CAS: 34487-61-1	UN: 1170
Flash point: 32°C	GHS: Flam. Liq. 3
ADR/RID: 3, III	H226
EINECS: 252-057-8	P280, P210, P240



Guaranteed analysis

pH colour change test

: Pass

pH 6.8-8.4 yellow to red.

In solution 75% water and 25% ethanol.

Phenol disulfonic acid solution 25%w/v in sulfuric acid

pure

packaging	Mfr. No
500 mℓ Coated Glass bottle	J/5900/PB08
Density at 20°C: 1.82	H314
ADR/RID: 8, II	P280, P301+P330+P331,
EINECS: 231-639-5	P310, P305+P351+P338,
UN: 1803	P303+P361+P353
GHS: Skin Corr. 1A	



Phenol-indo-2.6-dichlorophenol sodium salt

pure

packaging	Mfr. No
1 g Glass jar	P/2390/43
5 g Glass jar	P/2390/44
C ₁₂ H ₆ Cl ₂ NNaO ₂	2, Eye Irrit. 2, STOT SE 3
CAS: 620-45-1	H302, H315, H319, H335
MW: 290.08	P280, P305+P351+P338,
EINECS: 210-640-4	P301+P312, P302+P350,
GHS: Acute Tox. 4, Skin Irrit.	P304+P341



Guaranteed analysis

Complies with standard IR spectrum: Pass

Phenol-indo-2.6-dichlorophenol sodium salt

pure, tablets

packaging	Mfr. No
20 PC Glass bottle	P/2392/82
C ₁₂ H ₆ Cl ₂ NNaO ₂	2, Eye Irrit. 2, STOT SE 3
CAS: 620-45-1	H302, H315, H319, H335
MW: 290.07	P280, P305+P351+P338,
EINECS: 210-640-4	P301+P312, P302+P350,
GHS: Acute Tox. 4, Skin Irrit.	P304+P341



Guaranteed analysis

Identity test

: Pass

Phenolphthalein solid

pure, indicator grade

packaging	Mfr. No
100 g Plastic jar	P/2400/48
500 g Plastic jar	P/2400/53
10 kg Plastic keg	P/2400/68
C ₂₀ H ₁₄ O ₄	EINECS: 201-004-7
GHS: 77-09-8	GHS: Carc. 1B, Muta. 2, Repr. 2
MW: 318.33	H350, H341, H361f
Melting point: 258 to 263°C	P201, P308+P313, P281



Guaranteed analysis

pH colour change test

: Pass

pH 8.3-10.0 colourless to red.

Complies with standard IR spectrum: Pass

Phenolphthalein solution 0.2% in water

pure, indicator grade

packaging	Mfr. No
500 mℓ Glass bottle	P/2400L/08
2.5 ℓ Glass bottle	P/2400L/17
CAS: 77-09-8	Tox. 4, STOT SE 2
Flash point: 13°C	H225, H302, H312, H332, H371
ADR/RID: 3, 6.1, II	P210, P240, P280,
EINECS: 201-004-7	P301+P312, P304+P340,
UN: 1992	P303+P361+P353, P308+P313
GHS: Flam. Liq. 2, Acute	



Guaranteed analysis

pH colour change test

: Pass

For non-aqueous titration.

Phenolphthalein solution 1% in water

pure, indicator grade

packaging	Mfr. No
500 mℓ Glass bottle	P/2401L/08
CAS: 77-09-8	GHS: Flam. Liq. 2, Muta. 2,
Density at 20°C: 0.79	Acute Tox. 4, STOT SE 2
Flash point: 13°C	H225, H371, H350, H341,
ADR/RID: 3, 6.1, II	H332, H312, H302
EINECS: 201-004-7	P280, P312, P302+P350,
UN: 1992	P308+P313, P260, P210



Guaranteed analysis

pH colour change test

: Pass

For non-aqueous titration.

Phenosafarine C.I. 50200

pure

packaging	Mfr. No
5 g Glass jar	P/2450/44
C ₁₈ H ₁₅ ClN ₄	GHS: Skin Irrit. 2, Eye Irrit. 2,
CAS: 81-93-6	STOT SE 3
MW: 322.80	H319, H335, H315
Melting point: 300°C	P261, P302+P352, P280,
EINECS: 201-387-0	P305+P351+P338



Guaranteed analysis

Complies with standard IR spectrum: Pass

DL-beta-Phenylalanine

specified reagents for general laboratory work- extra pure

SLR

packaging	Mfr. No
100 g Plastic jar	P/2461/48
C ₉ H ₁₁ NO ₂	Melting point: 266°C
CAS: 150-30-1	EINECS: 205-756-7
MW: 165.19	

Assay: >98.5%

Guaranteed analysis

Loss on drying	<0.5%	Iron (Fe)	<10 ppm
Residue after ignition	<0.1%	Lead (Pb)	<10 ppm
Melting Point tested	265 to 275°C		

N-Phenylanthranilic acid

pure

packaging	Mfr. No
5 g Plastic jar	P/2460/44
C ₁₃ H ₁₁ NO ₂	STOT SE 3
CAS: 91-40-7	H319, H335, H315
EINECS: 202-066-8	P261, P302+P352, P280,
GHS: Skin Irrit. 2, Eye Irrit. 2,	P305+P351+P338



Guaranteed analysis

Complies with standard IR spectrum: Pass

Phenylbromide See Bromobenzene 99+ %, page 60

p-Phenylenediamine

specified reagents for general laboratory work- extra pure

SLR

packaging	Mfr. No
100 g Glass jar	P/2520/48
C ₆ H ₈ N ₂	GHS: Acute Tox. 3, Eye Irrit.
CAS: 106-50-3	2, Skin Sens. 1, Aquatic
MW: 108.14	Acute 1, Aquatic Chronic 1
Boiling point: 267°C	H410, H311, H317, H331,
Melting point: 139 to 141°C	H319, H301
Flash point: 156°C	P280, P312, P302+P350,
ADR/RID: 6.1, III	P304+P340, P273,
EINECS: 203-404-7	P305+P351+P338
UN: 1673	



Assay: >97%

Phenylhydrazine

specified reagents for general laboratory work- extra pure

SLR

packaging	Mfr. No
250 mℓ Coated Glass bottle	P/2560/PB07
C ₆ H ₈ N ₂	GHS: Carc. 1B, Muta. 2,
CAS: 100-63-0	Acute Tox. 3, Skin Irrit. 2,
MW: 108.15	Eye Irrit. 2, Skin Sens. 1,
Density at 20°C: 1.09	STOT RE 1, Aquatic Acute 1
Boiling point: 238 to 241°C	H372, H319, H341, H350,
Freezing point: >18°C	H331, H311, H315, H317,
Melting point: 19°C	H400, H301
Flash point: 89°C	P280, P302+P350,
ADR/RID: 6.1, II	P304+P340, P301+P312,
EINECS: 202-873-5	P260, P273, P305+P351+P338
UN: 2572	



Phenylhydrazinium chloride

specified reagents for general laboratory work- extra pure

SLR

packaging	Mfr. No
250 g Glass jar	P/2600/50
C ₆ H ₈ N ₂ .HCl	Eye Irrit. 2, Skin Sens. 1,
CAS: 59-88-1	STOT RE 1, Aquatic Acute 1
MW: 144.61	H331, H301, H311, H341,
Melting point: 250 to 254°C	H350, H317, H315, H319,
ADR/RID: 6.1, II	H372, H400
EINECS: 200-444-7	P280, P301+P310,
UN: 2811	P302+P350, P304+P340,
GHS: Carc. 1B, Muta. 2,	P305+P351+P338,
Acute Tox. 3, Skin Irrit. 2,	P308+P313, P273



Assay: >97%

Guaranteed analysis

Calcium (Ca)	<200 ppm	Sodium (Na)	<1000 ppm
Copper (Cu)	<50 ppm	Total phosphorus (P)	<500 ppm
Iron (Fe)	<50 ppm	Total silicon (Si)	<200 ppm
Lead (Pb)	<50 ppm	Total sulphur (S)	<500 ppm
Magnesium (Mg)	<200 ppm	Zinc (Zn)	<200 ppm
Potassium (K)	<500 ppm		

Phenylamine See Aniline, page 48

Phenylmethysulphonyl fluoride		SLR
<i>specified reagents for general laboratory work- extra pure</i>		
packaging	Mfr. No	
5 g Glass jar	P/2805/44	
C ₇ H ₇ FO ₂ S CAS: 329-98-6 MW: 174.19 Melting point: 90 to 94°C ADR/RID: 8, 6.1, II EINECS: 206-350-2 UN: 2923	GHS: Acute Tox. 3, Skin Corr. 1B, Eye Dam. 1 H301, H314 P301+P310, P301+P330+P331, P280, P305+P351+P338, P310	 
Assay: >99%		
Guaranteed analysis		
IR Analysis Conform		

3-Phenyl-2-propenal *See Cinnamaldehyde, page 90*

Phenyl salicylate		SLR
<i>specified reagents for general laboratory work- extra pure</i>		
packaging	Mfr. No	
250 g Plastic jar	P/2840/50	
C ₁₃ H ₁₀ O ₃ CAS: 118-55-8 MW: 214.22 Density at 20°C: 1.25	Boiling point: 172 to 173°C Melting point: 40 to 44°C Flash point: >110°C EINECS: 204-259-2	
Guaranteed analysis		
Complies with standard IR spectrum: Pass		

Phloroglucinol dihydrate		Certified AR
<i>for analysis</i>		
packaging	Mfr. No	
25 g Plastic jar	P/2922/46	
C ₆ H ₆ O ₃ ·2H ₂ O CAS: 6099-90-7 MW: 162.14 Melting point: 218 to 221°C EINECS: 203-611-2	GHS: Skin Irrit. 2, Eye Irrit. 2, STOT SE 3 H335, H315, H319 P261, P302+P352, P280, P305+P351+P338	
Assay: >98%		
Guaranteed analysis		
Acidity / Alkalinity	<0.01 meq/g	Potassium (K) <10 ppm
Calcium (Ca)	<5 ppm	Sodium (Na) <50 ppm
Copper (Cu)	<10 ppm	Total phosphorus (P) <20 ppm
Iron (Fe)	<5 ppm	Total silicon (Si) <10 ppm
Lead (Pb)	<10 ppm	Total sulphur (S) <50 ppm
Magnesium (Mg)	<10 ppm	Zinc (Zn) <10 ppm
<i>Lot analysis on the label.</i>		

Phloroglucinol dihydrate		SLR
<i>specified reagents for general laboratory work- extra pure</i>		
packaging	Mfr. No	
10 g Plastic jar 25 g Plastic jar	P/2920/45 P/2920/46	
C ₆ H ₆ O ₃ ·2H ₂ O CAS: 6099-90-7 MW: 162.14 Melting point: 218 to 221°C EINECS: 203-611-2	GHS: Skin Irrit. 2, Eye Irrit. 2, STOT SE 3 H335, H315, H319 P261, P302+P352, P280, P305+P351+P338	
Assay: >97%		
Guaranteed analysis		
Calcium (Ca)	<100 ppm	Sodium (Na) <200 ppm
Copper (Cu)	<20 ppm	Total phosphorus (P) <100 ppm
Iron (Fe)	<50 ppm	Total silicon (Si) <50 ppm
Lead (Pb)	<20 ppm	Total sulphur (S) <200 ppm
Magnesium (Mg)	<50 ppm	Zinc (Zn) <50 ppm
Potassium (K)	<50 ppm	

Phosphate standard solution 1000 ppm (1 mg/ml)	
<i>for Ion pair chromatography, Ready-to-use</i>	
packaging	Mfr. No
100 ml HDPE Plastic bottle	J/4558/05
<i>Manufactured to a tolerance of ± 1%.</i>	
Phosphate standard 1000 ppm	
<i>For Ion chromatography in H₂O, NIST certified</i>	
packaging	Mfr. No
100 ml Plastic bottle	SPASPO492Y
500 ml Plastic bottle	SPASPO492X
<i>Manufactured to a tolerance of ± 0.3%.</i>	
<i>Check with your local distributor for availability of Spex product</i>	

Phosphate-phosphorus standard 1000 ppm	
<i>For Ion chromatography in H₂O, NIST certified</i>	
packaging	Mfr. No
100 ml Plastic bottle	SPASPO4P92Y
500 ml Plastic bottle	SPASPO4P92X
<i>Manufactured to a tolerance of ± 0.3%.</i>	
<i>Check with your local distributor for availability of Spex product</i>	

Phosphorous acid		SLR
<i>extra pure, crystalline</i>		
packaging	Mfr. No	
500 g Plastic jar	P/3040/53	
H ₃ PO ₃ CAS: 13598-36-2 MW: 82.00 Density at 20°C: 1.65 Boiling point: 200°C Melting point: 73°C Flash point: >100°C	ADR/RID: 8, III EINECS: 237-066-7 UN: 2834 GHS: Skin Corr. 1A, Eye Dam. 1 H314 P301+P330+P331, P280, P305+P351+P338, P310	
Assay: >98%		
Guaranteed analysis		
Iron (Fe)	<50 ppm	n-Butyl formate <0.005%
Sulfate (SO ₄)	<0.03%	

Phosphorus standard metal solution 1000 ppm	
<i>for ICP, solution aqueous</i>	
packaging	Mfr. No
500 ml HDPE Plastic bottle	J/8299/08
<i>Manufactured to a tolerance of ± 0.5%.</i>	
<i>All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.</i>	
<i>See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.</i>	

Phosphorus standard metal solution 10 000 ppm	
<i>for ICP, aqueous solution</i>	
packaging	Mfr. No
100 ml HDPE Plastic bottle	J/8298/05
<i>Manufactured to a tolerance of ± 0.5%.</i>	
<i>All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.</i>	
<i>See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.</i>	

Phosphorus standard metal solution 1000 ppm	
<i>for ICP solution in H₂O, NIST certified</i>	
packaging	Mfr. No
100 ml Plastic bottle	SPPLP92Y
250 ml Plastic bottle	SPPLP92T
<i>Manufactured to a tolerance of ± 0.5%.</i>	
<i>All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurites.</i>	
<i>See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.</i>	
<i>Check with your local distributor for availability of Spex product</i>	

Phosphorus standard metal solution 10 000 ppm	
<i>for ICP solution in H₂O, NIST certified</i>	
packaging	Mfr. No
100 ml Plastic bottle	SPPLP93Y
500 ml Plastic bottle	SPPLP93X
<i>Manufactured to a tolerance of ± 0.5%.</i>	
<i>All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurites.</i>	
<i>See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.</i>	
<i>Check with your local distributor for availability of Spex product</i>	

Phosphorus pentachloride		SLR
<i>specified reagents for general laboratory work- extra pure</i>		
packaging	Mfr. No	
250 g Coated Glass bottle	P/3160/50	
PCl ₅ CAS: 10026-13-8 MW: 208.24 Density at 20°C: 1.6 Boiling point: 160°C Melting point: 167°C ADR/RID: 8, II EINECS: 233-060-3	UN: 1806 GHS: Acute Tox. 2, Skin Corr. 1B, STOT RE 2 H330, H302, H373, H314, EUH029, EUH014 P280, P301+P330+P331, P305+P351+P338, P304+P340, P310, P402+P404	  
Assay: >98%		
Guaranteed analysis		
Residue after ignition <0.1%		

Phosphorus pentoxide		Certified AR
<i>for analysis</i>		
packaging	Mfr. No	
250 g Plastic jar 500 g Plastic jar	P/3000/50 P/3000/53	
P ₂ O ₅ CAS: 1314-56-3 MW: 141.94 Melting point: 340 to 360°C ADR/RID: 8, II EINECS: 215-236-1	UN: 1807 GHS: Skin Corr. 1A, Eye Dam. 1 H314 P301+P330+P331, P280, P305+P351+P338, P310	
Assay: >99%		
Guaranteed analysis		
Lower oxides	<0.01%	Potassium (K) <20 ppm
Calcium (Ca)	<20 ppm	Sodium (Na) <50 ppm
Copper (Cu)	<5 ppm	Total silicon (Si) <20 ppm
Iron (Fe)	<10 ppm	Total sulphur (S) <10 ppm
Lead (Pb)	<5 ppm	Zinc (Zn) <5 ppm
Magnesium (Mg)	<10 ppm	
Typical analysis		
Aluminium (Al)	<2 ppm	Molybdenum (Mo) <2 ppm
Barium (Ba)	<2 ppm	Nickel (Ni) <2 ppm
Cadmium (Cd)	<2 ppm	Silver (Ag) <2 ppm
Cobalt (Co)	<2 ppm	Strontium (Sr) <2 ppm
Lithium (Li)	<2 ppm	Titanium (Ti) <2 ppm
Manganese (Mn)	<2 ppm	Vanadium (V) <2 ppm
<i>Lot analysis on the label.</i>		

Phosphorus pentoxide		SLR
<i>specified reagents for general laboratory work- extra pure</i>		
packaging	Mfr. No	
100 g Plastic jar 250 g Plastic jar 500 g Plastic jar	P/2960/48 P/2960/50 P/2960/53	
P ₂ O ₅ CAS: 1314-56-3 MW: 141.94 Melting point: 340 to 360°C ADR/RID: 8, II EINECS: 215-236-1	UN: 1807 GHS: Skin Corr. 1A, Eye Dam. 1 H314 P301+P330+P331, P280, P305+P351+P338, P310	
Assay: >98%		
Guaranteed analysis		
Lower oxides	<0.05%	Potassium (K) <100 ppm
Calcium (Ca)	<100 ppm	Sodium (Na) <200 ppm
Copper (Cu)	<20 ppm	Total silicon (Si) <100 ppm
Iron (Fe)	<50 ppm	Total sulphur (S) <50 ppm
Lead (Pb)	<20 ppm	Zinc (Zn) <20 ppm
Magnesium (Mg)	<50 ppm	

Phosphorus tribromide		SLR
<i>specified reagents for general laboratory work- extra pure</i>		
packaging	Mfr. No	
100 ml Coated Glass bottle	P/3200/PB05	
PBr ₃ CAS: 7789-60-8 MW: 270.69 Density at 20°C: 2.88 Boiling point: 175°C Melting point: -40°C ADR/RID: 8, II EINECS: 232-178-2	UN: 1808 GHS: Skin Corr. 1B, Eye Dam. 1, STOT SE 3 H314, H335, EUH014 P301+P330+P331, P280, P305+P351+P338, P310, P261, P402+P404	 
Assay: >99%		

o-Phthalaldehyde		SLR
<i>specified reagents for general laboratory work- extra pure</i>		
packaging	Mfr. No	
5 g Glass jar	P/3298/44	
C ₈ H ₆ O ₂ CAS: 643-79-8 MW: 134.14 Density at 20°C: 1.13 Boiling point: 83 to 84°C Melting point: 54 to 57°C Flash point: 132°C ADR/RID: 6.1, 8, II EINECS: 211-402-2	UN: 2928 GHS: Acute Tox. 3, Skin Corr. 1B, Eye Dam. 1, Skin Sens. 1 H314, H317, H301 P301+P310, P301+P330+P331, P305+P351+P338, P310, P280	 
Guaranteed analysis		
Complies with standard IR spectrum: Pass		

For the automated fluorescence analysis of amino acids and peptides.

Phthalein Complexon
See o-Cresolphthaleindi-(methyliminodi-acetic acid) metal indicator, page 99

High Volume Solvent Delivery Systems

High-volume solvent delivery systems are available in 10L to 1000L, in Europe only. Please email your contact details and a brief description of your requirements to beigel.sha.quotes@thermofisher.com



Phthalic anhydride		SLR
<i>specified reagents for general laboratory work- extra pure</i>		
packaging	Mfr. No	
500 g Plastic jar	P/3360/53	
1 kg Plastic jar	P/3360/60	
C ₈ H ₄ O ₃ CAS: 85-44-9 MW: 148.12 Density at 20°C: 15.3 Boiling point: 284°C Melting point: 129 to 132°C Flash point: 152°C ADR/RID: 8, III EINECS: 201-607-5	UN: 2214 GHS: Acute Tox. 4, Skin Irrit. 2, Eye Dam. 1, Resp. Sens. 1, Skin Sens. 1, STOT SE 3 H317, H302, H334, H335, H318, H315 P280, P305+P351+P338, P260, P309+P311	  
Assay: >99%		
Guaranteed analysis		
Residue after ignition <0.005%		

Piperazine hexahydrate		SLR
<i>specified reagents for general laboratory work- extra pure</i>		
packaging	Mfr. No	
250 g Plastic jar	P/3500/50	
C ₄ H ₁₀ N ₂ ·6H ₂ O CAS: 142-63-2 MW: 194.24 Boiling point: 145 to 156°C Melting point: 41 to 45°C Flash point: 87°C ADR/RID: 8, III EINECS: 203-808-3 UN: 2579	GHS: Skin Corr. 1B, Eye Dam. 1, Resp. Sens. 1, Skin Sens. 1, Aquatic Chronic 3 H412, H314, H317, H334 P301+P330+P331, P305+P351+P338, P310, P261, P304+P341, P342+P311, P280, P273	 
Assay: >98%		

Piperidine		
<i>pure, 99+%</i>		
packaging	Mfr. No	
100 mL Coated Glass bottle	P/3520/PB05	
250 mL Coated Glass bottle	P/3520/PB07	
C ₅ H ₁₁ N CAS: 110-89-4 MW: 85.15 Density at 20°C: 0.862 Boiling point: 106°C Melting point: -11°C Flash point: 16°C ADR/RID: 8, 3, I EINECS: 203-813-0	UN: 2401 GHS: Flam. Liq. 2, Acute Tox. 3, Skin Corr. 1B, Eye Dam. 1 H225, H311, H331, H314 P280, P302+P350, P301+P330+P331, P305+P351+P338, P310, P304+P340, P210, P240	  
Assay (CG): >99%		

Platinum standard metal solution 1000 ppm		
<i>for ICP solution in HCL 10%, NIST certified</i>		
packaging	Mfr. No	
100 mL Plastic bottle	SPPLPT32Y	
500 mL Plastic bottle	SPPLPT32X	
H315, H319, H335		

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Platinum		SLR
<i>extra pure, metal wire 0.127mm</i>		
packaging	Mfr. No	
50 cm Plastic jar	P/3610/88	
Pt CAS: 7440-06-4 MW: 195.09	Boiling point: 3827°C Melting point: 1772°C EINECS: 231-116-1	
Assay: 99.9%		
Guaranteed analysis		
Identity test : Pass		
<i>The prices applied to platinum and derivatives vary with the price of the raw material.</i>		

Platinum		
<i>pure, metal wire 0.3 mm</i>		
packaging	Mfr. No	
50 cm Plastic jar	P/3600/88	
100 cm Plastic jar	P/3600/89	
1 M Plastic jar	P/3600/87	
Pt CAS: 7440-06-4 MW: 195.09	Boiling point: 3827°C Melting point: 1772°C EINECS: 231-116-1	
Assay: 99.9%		
Guaranteed analysis		
Identity test : Pass		

Platinum		SLR
<i>extra pure, metal wire 0.5mm</i>		
packaging	Mfr. No	
50 cm Plastic jar	P/3640/88	
100 cm Plastic jar	P/3640/89	
1 M Plastic jar	P/3640/87	
Pt CAS: 7440-06-4 MW: 195.09	Boiling point: 3827°C Melting point: 1772°C EINECS: 231-116-1	
Assay: 99.9%		
Guaranteed analysis		
Identity test : Pass		
<i>The prices applied to platinum and derivatives vary with the price of the raw material.</i>		

Platinum(IV) dioxide		
<i>pure</i>		
packaging	Mfr. No	
1 g Vial	P/3652/43	
PtO ₂ CAS: 1314-15-4 MW: 227.09	UN: 1479 GHS: Ox. Sol. 2, Eye Irrit. 2 H272, H319 P280, P305+P351+P338, P310, P210	 
Melting point: 450°C ADR/RID: 5-1, II EINECS: 215-223-0		
Guaranteed analysis		
Identity test : Pass		

Polyethylene glycol 6000 grade		SLR
<i>specified reagents for general laboratory work- extra pure</i>		
packaging	Mfr. No	
500 g Plastic jar	P/3681/53	
25 kg Fibre Keg	P/3681/70	
HO(C ₂ H ₄ O) _n H CAS: 25322-68-3 MW: 6000.00	Melting point: 55 to 60°C Flash point: 270°C EINECS: 500-038-2	
Guaranteed analysis		
Identity test : Pass		
<i>Grade number indicates approximate molecular weight.</i>		

Polyethylene glycol 4000 grade		SLR
<i>specified reagents for general laboratory work- extra pure</i>		
packaging	Mfr. No	
500 g Plastic jar	P/3680/53	
HO(C ₂ H ₄ O) _n H CAS: 25322-68-3 MW: 4000.00 Density at 20°C: 1.15	Melting point: 54 to 58°C Flash point: 270°C EINECS: 500-038-2	
Guaranteed analysis		
Acidity <0.008 meq/g		
<i>Grade number indicates approximate molecular weight</i>		

Polypropylene glycol 2000 grade		SLR
<i>specified reagents for general laboratory work- extra pure</i>		
packaging	Mfr. No	
1 L Glass bottle	P/3694/15	
CAS: 25322-69-4 Melting point: -40°C Flash point: 230°C EINECS: 500-039-8		
Absorbance specification		
Max. Absorbance (1cm cell v HPLC water)		
Refractive index 1.450 to 1.452		
<i>Grade number indicates approximate weight.</i>		

Polyethylene glycol 1500 grade		SLR
<i>specified reagents for general laboratory work- extra pure</i>		
packaging	Mfr. No	
500 g Plastic jar	P/3679/53	
25 kg Fibre Keg	P/3679/70	
HO(C ₂ H ₄ O) _n H CAS: 25322-68-3 MW: 1500.00	Melting point: -65°C Flash point: >250°C EINECS: 500-038-2	
Guaranteed analysis		
Acidity <0.008 meq/g		
<i>Grade number indicates approximate molecular weight</i>		

Polyethylene glycol 1000 grade		SLR
<i>specified reagents for general laboratory work- extra pure</i>		
packaging	Mfr. No	
500 g Plastic jar	P/3678/53	
25 kg Metal drum	P/3678/70	
HO(C ₂ H ₄ O) _n H CAS: 25322-68-3 MW: 1000.00	Melting point: 35 to 40°C Flash point: >260°C EINECS: 500-038-2	
Guaranteed analysis		
Acidity <0.008 meq/g		
<i>Grade number indicates approximate molecular weight</i>		

Polyethylene glycol 600 grade		SLR
<i>specified reagents for general laboratory work- extra pure</i>		
packaging	Mfr. No	
500 mL Glass bottle	P/3677/08	
HO(C ₂ H ₄ O) _n H CAS: 25322-68-3 MW: 600.00 Density at 20°C: 1.127	Melting point: 17 to 22°C Flash point: 252°C EINECS: 500-038-2	
Guaranteed analysis		
Acidity <0.008 meq/g		
<i>Grade number indicates approximate molecular weight.</i>		

Polyethylene glycol 400 grade		SLR
<i>specified reagents for general laboratory work- extra pure</i>		
packaging	Mfr. No	
500 mL Glass bottle	P/3676/08	
25 L HDPE Plastic container	P/3676/25	
HO(C ₂ H ₄ O) _n H CAS: 25322-68-3 MW: 400.00 Density at 20°C: 1.128	Boiling point: >200°C Melting point: 4 to 8°C Flash point: 235°C EINECS: 500-038-2	
Guaranteed analysis		
Acidity <0.008 meq/g		
<i>Grade number indicates approximate molecular weight</i>		

Polyethylene glycol 200 grade		SLR
<i>specified reagents for general laboratory work- extra pure</i>		
packaging	Mfr. No	
500 mL Glass bottle	P/3675/08	
25 L HDPE Plastic container	P/3675/25	
HO(C ₂ H ₄ O) _n H CAS: 25322-68-3 MW: 200.00 Density at 20°C: 1.125	Melting point: -65°C Flash point: 171°C EINECS: 500-038-2	
Guaranteed analysis		
Acidity <0.008 meq/g		
<i>Grade number indicates approximate molecular weight</i>		

Polyvinyl pyrrolidone		
<i>pure</i>		
packaging	Mfr. No	
100 g Plastic jar	P/3700/48	
(C ₆ H ₉ NO) _n CAS: 9003-39-8 MW: approx. 44.000	Density at 20°C: 1.1 Melting point: 100°C Flash point: >215°C	
Guaranteed analysis		
Complies with standard IR spectrum: Pass		

Potassium standard solution 1000 ppm (1 mg/mL)		
<i>for Ion pair chromatography, Ready-to-use</i>		
packaging	Mfr. No	
100 mL HDPE Plastic bottle	J/4560/05	
Density at 20°C: 1.0 Boiling point: 100°C Melting point: 0°C		
<i>Manufactured to a tolerance of ± 1%.</i>		

Potassium standard metal solution 1000 ppm		
<i>for AAS, solution HNO₃ 1M</i>		
packaging	Mfr. No	
100 mL HDPE Plastic bottle	J/8060/05	
500 mL HDPE Plastic bottle	J/8060/08	
ADR/RID: 8, II UN: 2031 GHS: Skin Corr. 1B, Eye Dam. 1	H314 P280, P305+P351+P338, P301+P330+P331, P310	
<i>Manufactured to a tolerance of ± 0.5%. See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.</i>		

Potassium standard metal solution 1000 ppm

for ICP, solution HNO₃ 1M

packaging	Mfr. No
500 mℓ HDPE Plastic bottle	J/8301/08
ADR/RID: 8, II	H314
UN: 2031	P280, P305+P351+P338,
GHS: Skin Corr. 1B, Eye Dam. 1	P301+P330+P331, P310

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.

Potassium standard metal solution 10 000 ppm

for ICP, solution HNO₃ 1M

packaging	Mfr. No
100 mℓ HDPE Plastic bottle	J/8300/05
ADR/RID: 8, II	H314
UN: 2031	P280, P305+P351+P338,
GHS: Skin Corr. 1B	P301+P330+P331, P310

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.

Potassium standard metal solution 1000 ppm

for ICP solution in HNO₃ 2%, NIST certified

packaging	Mfr. No
100 mℓ Plastic bottle	SPPLK22Y
250 mℓ Plastic bottle	SPPLK22T
500 mℓ Plastic bottle	SPPLK22X
H315, H319	

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Potassium standard metal solution 1000 ppm

for ICP solution in HCL 2%, NIST certified

packaging	Mfr. No
500 mℓ Plastic bottle	SPPLK12X

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Potassium standard metal solution 10 000 ppm

for ICP solution in HNO₃ 5%, NIST certified

packaging	Mfr. No
100 mℓ Plastic bottle	SPPLK23Y
500 mℓ Plastic bottle	SPPLK23X
H314	

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Potassium standard metal solution 10 000 ppm

for ICP solution in HCL 2%, NIST certified

packaging	Mfr. No
500 mℓ Plastic bottle	SPPLK13X
H314	

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Potassium standard metal solution 1000 ppm

for ICP-MS solution in HNO₃ 2%, NIST certified

packaging	Mfr. No
100 mℓ Plastic bottle	SPCLK22Y
H315, H319	

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Potassium metal in liquid paraffin

pure

packaging	Mfr. No
10 g Aluminium jar	P/3720/45
100 g Aluminium jar	P/3720/48
K	UN: 2257
CAS: 7440-09-7	GHS: Water-react. 1, Skin
MW: 39.10	Corr. 1B, Eye Dam. 1
Boiling point: 770°C	H260, H314, EUH014
Melting point: 63°C	P301+P330+P331, P280,
ADR/RID: 4.3, I	P305+P351+P338, P310,
EINECS: 231-119-8	P402+P404

Guaranteed analysis

Identity test : Pass

Potassium acetate

for analysis

Certified AR

packaging	Mfr. No
500 g Plastic jar	P/3761/53
1 kg Plastic jar	P/3761/60
2.5 kg Plastic keg	P/3761/62

KC ₂ H ₃ O ₂	Melting point: 292°C
CAS: 127-08-2	Flash point: >250°C
MW: 98.14	EINECS: 204-822-2

Assay: >99%

Guaranteed analysis

pH (aqueous solution 5%)	7 to 9	Total chloride (Cl)	<0.001%
Calcium (Ca)	<20 ppm	Total nitrogen (N)	<5 ppm
Copper (Cu)	<2 ppm	Total phosphorus (P)	<5 ppm
Iron (Fe)	<5 ppm	Total silicon (Si)	<10 ppm
Lead (Pb)	<2 ppm	Total sulphur (S)	<20 ppm
Magnesium (Mg)	<10 ppm	Zinc (Zn)	<2 ppm
Sodium (Na)	<500 ppm		

Typical analysis

Aluminium (Al)	<2 ppm	Manganese (Mn)	<1 ppm
Barium (Ba)	<1 ppm	Molybdenum (Mo)	<1 ppm
Cadmium (Cd)	<1 ppm	Nickel (Ni)	<1 ppm
Chromium (Cr)	<1 ppm	Silver (Ag)	<1 ppm
Cobalt (Co)	<1 ppm	Titanium (Ti)	<1 ppm
Lithium (Li)	<1 ppm	Vanadium (V)	<1 ppm

Lot analysis on the label.

Potassium acetate

specified reagents for general laboratory work- extra pure

Potassium bromate

specified reagents for general laboratory work- extra pure

packaging	Mfr. No
500 g Plastic jar	P/3760/53
1 kg Plastic jar	P/3760/60
2.5 kg Plastic keg	P/3760/62
5 kg Plastic keg	P/3760/65

KC ₂ H ₃ O ₂	Melting point: 292°C
CAS: 127-08-2	Flash point: >250°C
MW: 98.14	EINECS: 204-822-2

Assay: >99%

Guaranteed analysis

pH (aqueous solution 5%)	7 to 9	Sodium (Na)	<2000 ppm
Calcium (Ca)	<100 ppm	n-Butyl formate	<0.005%
Copper (Cu)	<10 ppm	Total phosphorus (P)	<20 ppm
Iron (Fe)	<20 ppm	Total silicon (Si)	<50 ppm
Lead (Pb)	<10 ppm	Total sulphur (S)	<100 ppm
Magnesium (Mg)	<50 ppm	Zinc (Zn)	<10 ppm

Potassium bicarbonate

See Potassium hydrogen carbonate, page 254

Potassium biphtalate

See Potassium hydrogen phthalate, page 255

Potassium bitartrate

See Potassium hydrogen tartrate, page 255

Potassium bromate

for analysis

Certified AR

packaging	Mfr. No
250 g Plastic jar	P/3920/50
500 g Plastic jar	P/3920/53
1 kg Plastic jar	P/3920/60

KBrO ₃	UN: 1484
CAS: 7758-01-2	GHS: Ox. Sol. 1, Carc. 1B,
MW: 167.01	Acute Tox. 3
Density at 20°C: 3.27	H271, H301, H350
Melting point: 350°C	P301+P310, P201,
ADR/RID: 5.1, II	P308+P313, P210
EINECS: 231-829-8	

Assay: >99.9%

Guaranteed analysis

pH (aqueous solution 5%)	5 to 9	Lead (Pb)	<5 ppm
Loss on drying	<0.05%	Magnesium (Mg)	<5 ppm
Bromide	<0.01%	Sodium (Na)	<100 ppm
Bromide (Br)	<0.01 %	Total phosphorus (P)	<20 ppm
Calcium (Ca)	<20 ppm	Total silicon (Si)	<10 ppm
Copper (Cu)	<5 ppm	Total sulphur (S)	<20 ppm
Iron (Fe)	<5 ppm	Zinc (Zn)	<5 ppm

Typical analysis

Aluminium (Al)	<2 ppm	Molybdenum (Mo)	<2 ppm
Barium (Ba)	<2 ppm	Nickel (Ni)	<2 ppm
Cadmium (Cd)	<2 ppm	Silver (Ag)	<2 ppm
Chromium (Cr)	<2 ppm	Strontium (Sr)	<2 ppm
Cobalt (Co)	<2 ppm	Titanium (Ti)	<2 ppm
Lithium (Li)	<2 ppm	Vanadium (V)	<2 ppm
Manganese (Mn)	<2 ppm		

Lot analysis on the label.

packaging	Mfr. No
250 g Plastic jar	P/3880/50
500 g Plastic jar	P/3880/53

KBrO ₃	UN: 1484
CAS: 7758-01-2	GHS: Ox. Sol. 1, Carc. 1B,
MW: 167.01	Acute Tox. 3
Density at 20°C: 3.27	H271, H301, H350
Melting point: 350°C	P301+P310, P201,
ADR/RID: 5.1, II	P308+P313, P210
EINECS: 231-829-8	

Assay: >99.5%

Guaranteed analysis

pH (aqueous solution 5%)	5 to 9	Magnesium (Mg)	<20 ppm
Bromide	: <0.05%	Sodium (Na)	<500 ppm
Bromide (Br)	<0.05%	Total phosphorus (P)	<100 ppm
Calcium (Ca)	<100 ppm	Total silicon (Si)	<50 ppm
Copper (Cu)	<20 ppm	Total sulphur (S)	<100 ppm
Iron (Fe)	<20 ppm	Zinc (Zn)	<20 ppm
Lead (Pb)	<20 ppm		

Potassium bromate-bromide solution 0.0167 M (0.1N)

NIST Standard solution ready to use, For Volumetric analysis

packaging	Mfr. No
1 ℓ Glass bottle	J/6100/15
2.5 ℓ Glass bottle	J/6100/17

GHS: Carc. 1B, Skin Irrit. 2, Eye Irrit. 2, STOT SE 3
H350, H315, H319, H335
P261, P280, P308+P313, P302+P352, P305+P351+P338

Manufactured to a tolerance of ± 0.1%.
See the full range of Standard Volumetric solution page 323.

Potassium bromate-bromide solution 0.0167 M (0.1N)

NIST Standard Concentrate, For Volumetric analysis

packaging	Mfr. No
1 Ampoule	J/6100C/05
6 Ampoule	J/6100C/90

GHS: Carc. 1B, Skin Irrit. 2, Eye Irrit. 2, STOT SE 3
H350, H315, H319, H335
P261, P280, P308+P313, P302+P352, P305+P351+P338

Each ampoule makes 1 ℓ of solution.
See the full range of Standard Volumetric solution page 323.

Potassium bromide

For IR Spectroscopy

packaging	Mfr. No
100 g Plastic jar	P/4020/48

KBr	EINECS: 231-830-3
CAS: 7758-02-3	GHS: Eye Irrit. 2
MW: 119.01	H319
Boiling point: 1435°C	P280, P305+P351+P338
Melting point: 734°C	

Guaranteed analysis

IR Analysis : Pass

Lot analysis on the label.

Potassium bromide		Certified AR
for analysis		
packaging		Mfr. No
250 g	Plastic jar	P/4000/50
500 g	Plastic jar	P/4000/53
1 kg	Plastic jar	P/4000/60
KBr	EINECS: 231-830-3	
CAS: 7758-02-3	GHS: Eye Irrit. 2	
MW: 119.01	H319	
Boiling point: 1435°C	P280, P305+P351+P338	
Melting point: 734°C		
Assay: >99.5%		
Guaranteed analysis		
pH (aqueous solution 5%)	5 to 8	Sodium (Na) <1000 ppm
Loss on drying	<0.5%	Total chloride (Cl) <0.2%
Bromate (BrO ₃)	<0.0005%	Total nitrogen (N) <10 ppm
Calcium (Ca)	<10 ppm	Total phosphorus (P) <10 ppm
Copper (Cu)	<2 ppm	Total silicon (Si) <20 ppm
Iron (Fe)	<5 ppm	Total sulphur (S) <20 ppm
Lead (Pb)	<2 ppm	Zinc (Zn) <2 ppm
Magnesium (Mg)	<5 ppm	
Typical analysis		
Aluminium (Al)	<2 ppm	Molybdenum (Mo) <1 ppm
Barium (Ba)	<2 ppm	Nickel (Ni) <1 ppm
Cadmium (Cd)	<1 ppm	Silver (Ag) <2 ppm
Chromium (Cr)	<1 ppm	Strontium (Sr) <2 ppm
Cobalt (Co)	<1 ppm	Titanium (Ti) <1 ppm
Lithium (Li)	<1 ppm	Vanadium (V) <1 ppm
Manganese (Mn)	<1 ppm	
Lot analysis on the label.		

Potassium bromide		SLR
specified reagents for general laboratory work- extra pure		
packaging		Mfr. No
250 g	Plastic jar	P/3960/50
500 g	Plastic jar	P/3960/53
1 kg	Plastic jar	P/3960/60
25 kg	Plastic keg	P/3960/70
KBr	EINECS: 231-830-3	
CAS: 7758-02-3	GHS: Eye Irrit. 2	
MW: 119.01	H319	
Boiling point: 1435°C	P280, P305+P351+P338	
Melting point: 734°C		
Assay: >99%		
Guaranteed analysis		
pH (aqueous solution 5%)	5 to 8	Sodium (Na) <5000 ppm
Calcium (Ca)	<50 ppm	n-Butyl formate <0.5%
Copper (Cu)	<10 ppm	Total phosphorus (P) <50 ppm
Iron (Fe)	<20 ppm	Total silicon (Si) <100 ppm
Lead (Pb)	<10 ppm	Total sulphur (S) <100 ppm
Magnesium (Mg)	<20 ppm	Zinc (Zn) <10 ppm

Potassium carbonate anhydrous		Certified AR, Eur.Ph.
for analysis		
packaging		Mfr. No
250 g	Plastic jar	P/4120/50
500 g	Plastic jar	P/4120/53
1 kg	Plastic jar	P/4120/60
2.5 kg	Plastic jar	P/4120/62
K ₂ CO ₃	2, Eye Irrit. 2, STOT SE 3	
CAS: 584-08-7	H302, H335, H315, H319	
MW: 138.21	P280, P261,	
Density at 190°C: 2.43	P305+P351+P338,	
Melting point: 891°C	P301+P312, P302+P352,	
EINECS: 209-529-3	P304+P341	
GHS: Acute Tox. 4, Skin Irrit.		
Assay: 99.5 to 101%		

Guaranteed analysis			
Loss at 300°C	<1%	Sodium (Na)	<1000 ppm
Loss on drying	<5%	Total chloride (Cl)	<0.002%
Cadmium (Cd)	<1 ppm	Total nitrogen (N)	<20 ppm
Calcium (Ca)	<10 ppm	Total phosphorus (P)	<5 ppm
Copper (Cu)	<5 ppm	Total silicon (Si)	<20 ppm
Iron (Fe)	<5 ppm	Total sulphur (S)	<10 ppm
Lead (Pb)	<5 ppm	Zinc (Zn)	<2 ppm
Magnesium (Mg)	<5 ppm		
Typical analysis			
Aluminium (Al)	<2 ppm	Molybdenum (Mo)	<1 ppm
Barium (Ba)	<1 ppm	Nickel (Ni)	<1 ppm
Chromium (Cr)	<1 ppm	Silver (Ag)	<1 ppm
Cobalt (Co)	<1 ppm	Strontium (Sr)	<5 ppm
Lithium (Li)	<1 ppm	Titanium (Ti)	<1 ppm
Manganese (Mn)	<1 ppm	Vanadium (V)	<1 ppm

Lot analysis on the label.

Potassium carbonate anhydrous		SLR
specified reagents for general laboratory work- extra pure		
packaging		Mfr. No
500 g	Plastic jar	P/4080/53
1 kg	Plastic jar	P/4080/60
2.5 kg	Plastic jar	P/4080/62
5 kg	Plastic keg	P/4080/65
25 kg	Plastic keg	P/4080/70
K ₂ CO ₃		
CAS: 584-08-7		
MW: 138.21		
Density at 190°C: 2.43		
Melting point: 891°C		
EINECS: 209-529-3		
GHS: Acute Tox. 4, Skin Irrit.		
Assay: >99%		

Potassium chlorate		SLR
specified reagents for general laboratory work- extra pure		
packaging	Mfr. No	
500 g Plastic jar	P/4160/53	
KClO ₃	UN: 1485	
CAS: 3811-04-9	GHS: Ox. Sol. 1, Acute Tox.	
MW: 122.55	4, Aquatic Chronic 2	
Density at 20°C: 2.52	H271, H411, H302, H332	
Melting point: 356 to 368°C	P210, P273,	
ADR/RID: 5.1, II	P301+P312, P304+P340	
EINECS: 223-289-7		

Assay: >99%			
Guaranteed analysis			
Calcium (Ca)	<100 ppm	n-Butyl formate	<0.05%
Copper (Cu)	<10 ppm	Total phosphorus (P)	<200 ppm
Iron (Fe)	<20 ppm	Total silicon (Si)	<100 ppm
Lead (Pb)	<10 ppm	Total sulphur (S)	<500 ppm
Magnesium (Mg)	<50 ppm	Zinc (Zn)	<50 ppm
Sodium (Na)	<500 ppm		

Save money!

Order 4 bottles in a box



Potassium chloride		Certified AR	
for analysis			
packaging		Mfr. No	
500 g	Plastic jar	P/4280/53	
1 kg	Plastic jar	P/4280/60	
5 kg	Plastic keg	P/4280/65	
10 kg	Plastic keg	P/4280/68	
25 kg	Plastic keg	P/4280/70	
KCl	Boiling point: 1420°C		
CAS: 7447-40-7	Melting point: 770°C		
MW: 74.54	EINECS: 231-211-8		
Density at 20°C: 1.987			
Assay: >99.5%			
Guaranteed analysis			
pH (aqueous solution 5%)	5 to 8	Manganese (Mn)	<1 ppm
Loss on drying	<0.5%	Sodium (Na)	<500 ppm
Calcium (Ca)	<10 ppm	Total nitrogen (N)	<10 ppm
Copper (Cu)	<2 ppm	Total phosphorus (P)	<5 ppm
Iron (Fe)	<5 ppm	Total silicon (Si)	<10 ppm
Lead (Pb)	<2 ppm	Total sulphur (S)	<10 ppm
Magnesium (Mg)	<5 ppm	Zinc (Zn)	<2 ppm
Typical analysis			
Aluminium (Al)	<5 ppm	Molybdenum (Mo)	<1 ppm
Barium (Ba)	<5 ppm	Nickel (Ni)	<1 ppm
Cadmium (Cd)	<1 ppm	Silver (Ag)	<1 ppm
Chromium (Cr)	<1 ppm	Strontium (Sr)	<5 ppm
Cobalt (Co)	<1 ppm	Titanium (Ti)	<1 ppm
Lithium (Li)	<1 ppm	Vanadium (V)	<1 ppm
Lot analysis on the label.			

Potassium chloride		SLR, Eur. Ph.	
specified reagents for general laboratory work- extra pure			
packaging		Mfr. No	
500 g	Plastic jar	P/4240/53	
1 kg	Plastic jar	P/4240/60	
2.5 kg	Plastic jar	P/4240/62	
5 kg	Plastic jar	P/4240/65	
25 kg	Plastic keg	P/4240/70	
KCl	Boiling point: 1420°C		
CAS: 7447-40-7	Melting point: 770°C		
MW: 74.54	EINECS: 231-211-8		
Density at 20°C: 1.987			
Assay: 99 to 100.5%			
Guaranteed analysis			
Acidity / Alkalinity	Pass	Iron (Fe)	<20 ppm
pH (aqueous solution 5%)	5 to 8	Lead (Pb)	<10 ppm
Loss on drying	<1%	Magnesium (Mg)	<20 ppm
Bromide	: <0.1%	Sodium (Na)	<1000 ppm
Heavy metals	<10 ppm	Sulfate (SO ₄)	<300 ppm
Mg and alkaline earth metals.:	<200 ppm	Total phosphorus (P)	<20 ppm
Aluminium (Al)	<1 ppm	Total silicon (Si)	<50 ppm
Calcium (Ca)	<50 ppm	Total sulphur (S)	<50 ppm
Copper (Cu)	<10 ppm	Zinc (Zn)	<10 ppm
Iodide (I)	Pass		

Potassium chloride solution 4M	
pure, Electrode filling solution	
packaging	Mfr. No
100 mℓ Glass bottle	J/2986/05
KCl	
CAS: 7447-40-7	
MW: 74.54	
EINECS: 231-211-8	

Potassium chloride solution 4M	
pure, Electrode filling solution with silver chloride	
packaging	Mfr. No
100 mℓ Glass bottle	J/2984/05
KCl	
CAS: 7447-40-7	
MW: 74.54	
EINECS: 231-211-8	

Potassium chloride satured solution	
pure, Electrode filling solution	
packaging	Mfr. No
100 mℓ Glass bottle	J/2983/05
CAS: 7447-40-7	
EINECS: 231-211-8	

Potassium Chloride solution 3M	
pure	
packaging	Mfr. No
1 ℓ HDPE Plastic bottle	J/2892/15
CAS: 7447-40-7	
EINECS: 231-211-8	
Guaranteed analysis	
Value between 2.99 to 3.01	

Potassium chromate		Certified AR	
for analysis			
packaging	Mfr. No		
250 g Plastic jar	P/4400/50		
500 g Plastic jar	P/4400/53		
1 kg Plastic jar	P/4400/60		
25 kg Fibre Keg	P/4400/70		
K ₂ CrO ₄	Skin Irrit. 2, Eye Irrit. 2, Skin		
CAS: 7789-00-6	Sens. 1, STOT SE 3, Aquatic		
MW: 194.20	Acute 1, Aquatic Chronic 1		
Melting point: 975°C	H317, H315, H335, H340,		
ADR/RID: 6.1, 5.1, II	H319, H350i, H400, H410		
EINECS: 232-140-5	P308+P313, P280, P261,		
UN: 3086	P273, P302+P352,		
GHS: Carc. 1B, Muta. 1B,	P305+P351+P338		
Assay: >99.5%			
Guaranteed analysis			
pH (aqueous solution 5%)	9 to 10	Sodium (Na)	<2000 ppm
Calcium (Ca)	<50 ppm	Total chloride (Cl)	<0.001%
Copper (Cu)	<5 ppm	Total phosphorus (P)	<100 ppm
Iron (Fe)	<5 ppm	Total silicon (Si)	<50 ppm
Lead (Pb)	<20 ppm	Total sulphur (S)	<100 ppm
Magnesium (Mg)	<5 ppm	Zinc (Zn)	<5 ppm
Typical analysis			
Aluminium (Al)	<0.002 %	Manganese (Mn)	<2 ppm
Barium (Ba)	<2 ppm	Nickel (Ni)	<2 ppm
Cadmium (Cd)	<2 ppm	Titanium (Ti)	<2 ppm
Cobalt (Co)	<2 ppm	Vanadium (V)	<2 ppm
Lithium (Li)	<2 ppm		
Lot analysis on the label.			

UHPLC

Filtered to 0.1 μm, Fisher Chemical UHPLC gradient grade solvents (Acetonitrile, Methanol and Water) are perfect for those working at high pressure.

See page 309 for more information



Potassium chromate		SLR
specified reagents for general laboratory work- extra pure		
packaging	Mfr. No	
250 g Plastic jar	P/4360/50	
500 g Plastic jar	P/4360/53	
1 kg Plastic jar	P/4360/60	
K ₂ CrO ₄ CAS: 7789-00-6 MW: 194.20 Melting point: 975°C ADR/RID: 6.1, 5.1, II EINECS: 232-140-5 UN: 3086 GHS: Carc. 1B, Muta. 1B,	Skin Irrit. 2, Eye Irrit. 2, Skin Sens. 1, STOT SE 3, Aquatic Acute 1, Aquatic Chronic 1 H317, H315, H335, H340, H319, H350i, H400, H410 P308+P313, P280, P261, P273, P302+P352, P305+P351+P338	

Assay: >99%		
Guaranteed analysis		
pH (aqueous solution 5%)	9 to 10	Sodium (Na) <5000 ppm
Calcium (Ca)	<200 ppm	n-Butyl formate <0.005%
Copper (Cu)	<20 ppm	Total phosphorus (P) <500 ppm
Iron (Fe)	<20 ppm	Total silicon (Si) <200 ppm
Lead (Pb)	<100 ppm	Total sulphur (S) <500 ppm
Magnesium (Mg)	<20 ppm	Zinc (Zn) <20 ppm

tri-Potassium citrate Certified AR, Eur.Ph.,BP,USP for analysis

Packaging		Mfr. No	
500 g	Plastic jar	P/4442/53	
3 kg	Plastic jar	P/4442/63	
C ₆ H ₅ K ₃ O ₇ ·H ₂ O		Density at 20°C: 1.98 Melting point: 230°C EINECS: 212-755-5	
CAS: 6100-05-6			
MW: 324.42			
Assay: 99 to 101%			
Guaranteed analysis			
pH (aqueous solution 5%)	8 to 9	Sodium (Na)	<1500 ppm
Calcium (Ca)	<20 ppm	Total chloride (Cl)	<0.001%
Copper (Cu)	<2 ppm	Total nitrogen (N)	<10 ppm
Iron (Fe)	<5 ppm	Total phosphorus (P)	<10 ppm
Lead (Pb)	<2 ppm	Total silicon (Si)	<20 ppm
Magnesium (Mg)	<5 ppm	Total sulphur (S)	<20 ppm
Oxalate (C ₂ O ₄)	<0.01%	Zinc (Zn)	<5 ppm
Typical analysis			
Aluminium (Al)	<2 ppm	Manganese (Mn)	<1 ppm
Arsenic (As)	<2 ppm	Molybdenum (Mo)	<1 ppm
Barium (Ba)	<1 ppm	Nickel (Ni)	<1 ppm
Cadmium (Cd)	<1 ppm	Silver (Ag)	<1 ppm
Chromium (Cr)	<1 ppm	Strontium (Sr)	<5 ppm
Cobalt (Co)	<1 ppm	Titanium (Ti)	<1 ppm
Lithium (Li)	<1 ppm	Vanadium (V)	<1 ppm

Lot analysis on the label.

tri-Potassium citrate SLR specified reagents for general laboratory work- extra pure

Packaging		Mfr. No
500 g	Plastic jar	P/4440/53
C ₆ H ₅ K ₃ O ₇ ·H ₂ O	Density at 20°C: 1.98	
CAS: 6100-05-6	Melting point: 230°C	
MW: 324.42	EINECS: 212-755-5	
Assay: >99%		
Guaranteed analysis		
pH (aqueous solution 5%)	>8	Sodium (Na) <2000 ppm
Calcium (Ca)	<100 ppm	n-Butyl formate <0.005%
Copper (Cu)	<10 ppm	Total phosphorus (P) <50 ppm
Iron (Fe)	<20 ppm	Total silicon (Si) <100 ppm
Lead (Pb)	<10 ppm	Total sulphur (S) <100 ppm
Magnesium (Mg)	<20 ppm	Zinc (Zn) <20 ppm

Potassium cyanate		SLR
specified reagents for general laboratory work- extra pure		
packaging		Mfr. No
600 g	Plastic jar	P/4520/53
KOCN	EINECS: 209-676-3	
CAS: 590-28-3	GHS: Acute Tox. 4	
MW: 81.12	H302	
Melting point: 315°C	P301+P312, P280	
Assay: >98%		
Guaranteed analysis		
Loss on drying	<2%	Potassium cyanide (KCN) <0.5%

Potassium cyanide Certified AR for analysis

packaging	Mfr. No	
250 g Plastic jar	P/4600/50	
25 kg Metal drum	P/4600/70	
KCN CAS: 151-50-8 MW: 65.12 Density at 160°C: 1.52 Boiling point: 1625°C Melting point: 634°C ADR/RID: 6.1, I EINECS: 205-792-3 UN: 1680	GHS: Met. Corr. 1, Acute Tox. 1, Aquatic Acute 1, Aquatic Chronic 1 H310, H410, H330, H300, H290, EUH032 P301+P310, P280, P302+P350, P305+P351+P338, P304+P340, P273	

Assay: >99%			
Guaranteed analysis			
Aluminium (Al)	<2 ppm	Molybdenum (Mo)	<1 ppm
Barium (Ba)	<1 ppm	Nickel (Ni)	<1 ppm
Cadmium (Cd)	<1 ppm	Propan--1-ol	<2 ppm
Calcium (Ca)	<20 ppm	Sodium (Na)	<5000 ppm
Chromium (Cr)	<1 ppm	Strontium (Sr)	<5 ppm
Cobalt (Co)	<1 ppm	Titanium (Ti)	<1 ppm
Copper (Cu)	<5 ppm	n-Butyl formate	<0.02%
Iron (Fe)	<50 ppm	Total phosphorus (P)	<200 ppm
Lead (Pb)	<5 ppm	Total silicon (Si)	<50 ppm
Lithium (Li)	<1 ppm	Total sulphur (S)	<100 ppm
Magnesium (Mg)	<10 ppm	Vanadium (V)	<1 ppm
Manganese (Mn)	<1 ppm	Zinc (Zn)	<20 ppm

Lot analysis on the label.

Potassium cyanide SLR specified reagents for general laboratory work- extra pure

packaging	Mfr. No	
250 g Plastic jar	P/4560/50	
KCN CAS: 151-50-8 MW: 65.12 Density at 160°C: 1.52 Boiling point: 1625°C Melting point: 634°C ADR/RID: 6.1, I EINECS: 205-792-3 UN: 1680	GHS: Met. Corr. 1, Acute Tox. 1, Aquatic Acute 1, Aquatic Chronic 1 H310, H410, H330, H300, H290, EUH032 P301+P310, P280, P302+P350, P305+P351+P338, P304+P340, P273	

Assay: >97%			
Guaranteed analysis			
Calcium (Ca)	<100 ppm	Sodium (Na)	<10000 ppm
Copper (Cu)	<20 ppm	Total phosphorus (P)	<500 ppm
Iron (Fe)	<200 ppm	Total silicon (Si)	<200 ppm
Lead (Pb)	<20 ppm	Total sulphur (S)	<500 ppm
Magnesium (Mg)	<50 ppm	Zinc (Zn)	<50 ppm

Potassium dichromate		Certified AR
for analysis		
packaging	Mfr. No	
250 g Plastic jar	P/4720/50	
500 g Plastic jar	P/4720/53	
1 kg Plastic jar	P/4720/60	
25 kg Plastic keg	P/4720/70	
K ₂ Cr ₂ O ₇ CAS: 7778-50-9 MW: 294.19 Density at 20°C: 2.676 Boiling point: 500°C Melting point: 398°C ADR/RID: 5.1, 6.1, II EINECS: 231-906-6 UN: 3086 GHS: Ox. Sol. 2, Carc. 1B, Muta. 1B, Acute Tox. 2, Skin	Corr. 1B, Eye Dam. 1, Resp. Sens. 1, Skin Sens. 1, Aquatic Acute 1, Aquatic Chronic 1 H272, H301, H312, H314, H317, H330, H334, H335, H340, H350, H360FD, H372, H410 P221, P280, P301+P310, P302+P350, P304+P341, P305+P351+P338, P273	

Assay: >99.9%			
Guaranteed analysis			
Loss on drying	<0.05%	Sodium (Na)	<100 ppm
Calcium (Ca)	<20 ppm	Total chloride (Cl)	<0.001%
Copper (Cu)	<10 ppm	Total phosphorus (P)	<50 ppm
Iron (Fe)	<20 ppm	Total silicon (Si)	<20 ppm
Lead (Pb)	<20 ppm	Total sulphur (S)	<50 ppm
Magnesium (Mg)	<10 ppm	Zinc (Zn)	<10 ppm
Typical analysis			
Aluminium (Al)	<5 ppm	Manganese (Mn)	<1 ppm
Barium (Ba)	<1 ppm	Nickel (Ni)	<1 ppm
Cadmium (Cd)	<1 ppm	Titanium (Ti)	<1 ppm
Cobalt (Co)	<1 ppm	Vanadium (V)	<1 ppm
Lithium (Li)	<1 ppm		

Lot analysis on the label.

Potassium dichromate SLR extra pure, powder

packaging	Mfr. No	
500 g Plastic jar	P/4690/53	
K ₂ Cr ₂ O ₇ CAS: 7778-50-9 MW: 294.19 Density at 20°C: 2.676 Boiling point: 500°C Melting point: 398°C ADR/RID: 5.1, 6.1, II EINECS: 231-906-6 UN: 3086 GHS: Ox. Sol. 2, Carc. 1B, Muta. 1B, Acute Tox. 2, Skin	Corr. 1B, Eye Dam. 1, Resp. Sens. 1, Skin Sens. 1, Aquatic Acute 1, Aquatic Chronic 1 H272, H301, H312, H314, H317, H330, H334, H335, H340, H350, H360FD, H372, H410 P221, P280, P301+P310, P302+P350, P304+P341, P305+P351+P338, P273	

Assay: >99.5%			
Guaranteed analysis			
Loss on drying	<0.1%	Sodium (Na)	<500 ppm
Calcium (Ca)	<100 ppm	n-Butyl formate	<0.02%
Copper (Cu)	<50 ppm	Total phosphorus (P)	<500 ppm
Iron (Fe)	<100 ppm	Total silicon (Si)	<100 ppm
Lead (Pb)	<100 ppm	Total sulphur (S)	<500 ppm
Magnesium (Mg)	<50 ppm	Zinc (Zn)	<50 ppm

Potassium dichromate		SLR
extra pure, crystals		
packaging	Mfr. No	
250 g Plastic jar	P/4680/50	
500 g Plastic jar	P/4680/53	
1 kg Plastic jar	P/4680/60	
5 kg Plastic jar	P/4680/65	
K ₂ Cr ₂ O ₇ CAS: 7778-50-9 MW: 294.19 Density at 20°C: 2.676 Boiling point: 500°C Melting point: 398°C ADR/RID: 5.1, 6.1, II EINECS: 231-906-6 UN: 3086 GHS: Ox. Sol. 2, Carc. 1B, Muta. 1B, Acute Tox. 2, Skin	Corr. 1B, Eye Dam. 1, Resp. Sens. 1, Skin Sens. 1, Aquatic Acute 1, Aquatic Chronic 1 H272, H301, H312, H314, H317, H330, H334, H335, H340, H350, H360FD, H372, H410 P221, P280, P301+P310, P302+P350, P304+P341, P305+P351+P338, P273	

Assay: >99.5%			
Guaranteed analysis			
Loss on drying	<0.1%	Sodium (Na)	<500 ppm
Calcium (Ca)	<100 ppm	n-Butyl formate	<0.02%
Copper (Cu)	<50 ppm	Total phosphorus (P)	<500 ppm
Iron (Fe)	<100 ppm	Total silicon (Si)	<100 ppm
Lead (Pb)	<100 ppm	Total sulphur (S)	<500 ppm
Magnesium (Mg)	<50 ppm	Zinc (Zn)	<50 ppm

Potassium dichromate pure

packaging	Mfr. No	
1 kg Plastic jar	P/4640/60	
3 kg Plastic jar	P/4640/63	
K ₂ Cr ₂ O ₇ CAS: 7778-50-9 MW: 294.19 Density at 20°C: 2.676 Boiling point: 500°C Melting point: 398°C ADR/RID: 5.1, 6.1, II EINECS: 231-906-6 UN: 3086 GHS: Ox. Sol. 2, Carc. 1B, Muta. 1B, Acute Tox. 2, Skin	Corr. 1B, Eye Dam. 1, Resp. Sens. 1, Skin Sens. 1, Aquatic Acute 1, Aquatic Chronic 1 H272, H301, H312, H314, H317, H330, H334, H335, H340, H350, H360FD, H372, H410 P221, P280, P301+P310, P302+P350, P304+P341, P305+P351+P338, P273	

Guaranteed analysis	
Identity test	: Pass

Potassium dichromate solution 0.0167M (0.1N) Eur.Ph. , USP, BP NIST Standard solution ready to use, For Volumetric analysis

packaging	Mfr. No	
1 ℓ Glass bottle	J/6212/15	
2.5 ℓ Glass bottle	J/6212/17	
CAS: 7778-50-9 Density at 20°C: 1 ADR/RID: /II EINECS: 231-906-6	GHS: Carc. 1B, Muta. 1B, Repr. 1B, Aquatic Chronic 3 H340, H350, H412, H360FD P201, P308+P313, P273	

Manufactured to a tolerance of ± 0.1%.
See the full range of Standard Volumetric solution page 323.

Potassium dichromate solution 0.0167M (0.1N)

NIST Standard Concentrate, For Volumetric analysis

packaging		Mfr. No
1	Ampoule	J/6212C/05
6	Ampoule	J/6212C/90
CAS: 7778-50-9		Sens. 1, Skin Sens. 1, STOT
ADR/RID: 6.1, III		RE 2, Aquatic Chronic 2
EINECS: 231-906-6		H350, H340, H360FD,
UN: 3287		H331, H373, H315, H319,
GHS: Carc. 1B, Muta. 1B,		H334, H317, H411
Repr. 1B, Acute Tox. 3, Skin		P280, P305+P351+P338,
Irrit. 2, Eye Irrit. 2, Resp.		P201, P308+P313, P260, P273



Each ampoule makes 1 l of solution.

See the full range of Standard Volumetric solution page 323.

Potassium dichromate solution 0.021 M (0.125N)

NIST For COD and volumetric analysis

packaging		Mfr. No
1 l	Glass bottle	J/6220/15
2.5 l	Glass bottle	J/6220/17
CAS: 7778-50-9		Repr. 1B
EINECS: 231-906-6		H350, H340, H360FD, EUH203
GHS: Carc. 1B, Muta. 1B,		P201, P281, P308+P313



Manufactured to a tolerance of ± 0.2%.

See the full range of Standard Volumetric solution page 323.

Potassium dihydrogen orthophosphate

for HPLC

Certified

packaging		Mfr. No
250 g	Plastic jar	P/4806/50
1 kg	Plastic jar	P/4806/60
25 kg	Plastic keg	P/4806/70
KH ₂ PO ₄		Melting point: 253°C
CAS: 7778-77-0		EINECS: 231-913-4
MW: 136.09		
Assay: >99%		

Guaranteed analysis		
Reaction (5% solution) at 25°C ... : 4.1 to 4.5	Lead (Pb)	<5 ppm
	n-Butyl formate	<0.0005%
Iron (Fe)	<5 ppm	

Absorbance specification		
Max. Absorbance (1cm cell v HPLC water)		
254 nm (sol. aqu. 1 M)		
: <0.05 A.U.		

Lot analysis on the label.

Filtered to 0.2 µm.



Potassium dihydrogen orthophosphate

HPLC for electrochemical detection

Certified

packaging		Mfr. No
250 g	Plastic jar	P/4810/50
KH ₂ PO ₄		Melting point: 253°C
CAS: 7778-77-0		EINECS: 231-913-4
MW: 136.09		
Assay: >99%		
Guaranteed analysis		
Electrochemical test	Pass	Iron (Fe)
Reaction (5% solution) at 25°C ... : 4.1 to 4.5	Lead (Pb)	<5 ppm
	n-Butyl formate	<0.0005%

Absorbance specification		
Max. Absorbance (1cm cell v HPLC water)		
254 nm (sol. aqu. 1 M)		
: <0.05 A.U.		

Tested by HPLC-EC to ensure low electrochemical response.

Detector - ESA Coulochem II

Mode-oxidative

Graphite Electrode.

Lot analysis on the label.

Filtered to 0.2 µm.

Potassium dihydrogen orthophosphate

for analysis

Certified AR

packaging		Mfr. No
250 g	Plastic jar	P/4800/50
500 g	Plastic jar	P/4800/53
1 kg	Plastic jar	P/4800/60
3 kg	Plastic jar	P/4800/63
25 kg	Plastic keg	P/4800/70

KH ₂ PO ₄		Melting point: 253°C
CAS: 7778-77-0		EINECS: 231-913-4
MW: 136.09		

Assay: >99.5%		
Guaranteed analysis		
pH (aqueous solution 5%)	4.3 to 4.4	Magnesium (Mg)
Loss on drying	<0.05%	<20 ppm
Heavy metals	<10 ppm	Sodium (Na)
Insoluble matter	<0.01%	<50 ppm
Calcium (Ca)	<50 ppm	Total chloride (Cl)
Copper (Cu)	<5 ppm	<0.0005%
Iron (Fe)	<10 ppm	Total nitrogen (N)
Lead (Pb)	<5 ppm	<10 ppm
		Total silicon (Si)
		<50 ppm
		Total sulphur (S)
		<100 ppm
		Zinc (Zn)
		<10 ppm
		Sulfate (SO ₄)
		<0.003%

Typical analysis		
Aluminium (Al)	<2 ppm	Molybdenum (Mo)
Barium (Ba)	<2 ppm	<2 ppm
Cadmium (Cd)	<2 ppm	Nickel (Ni)
Chromium (Cr)	<2 ppm	<2 ppm
Cobalt (Co)	<2 ppm	Silver (Ag)
Lithium (Li)	<2 ppm	<2 ppm
Manganese (Mn)	<2 ppm	Strontium (Sr)
		<5 ppm
		Titanium (Ti)
		<2 ppm
		Vanadium (V)
		<2 ppm

Lot analysis on the label.

Potassium dihydrogen orthophosphate

specified reagents for general laboratory work- extra pure

SLR

packaging		Mfr. No
500 g	Plastic jar	P/4760/53
1 kg	Plastic jar	P/4760/60
3 kg	Plastic jar	P/4760/63
5 kg	Plastic jar	P/4760/65
25 kg	Plastic keg	P/4760/70

KH ₂ PO ₄		Melting point: 253°C
CAS: 7778-77-0		EINECS: 231-913-4
MW: 136.09		

Assay: >99%		
Guaranteed analysis		
pH (aqueous solution 5%)	4.2 to 4.5	Magnesium (Mg)
Loss on drying	<0.1%	<100 ppm
Calcium (Ca)	<200 ppm	Sodium (Na)
Copper (Cu)	<20 ppm	<1000 ppm
Iron (Fe)	<50 ppm	n-Butyl formate
Lead (Pb)	<20 ppm	<0.005%
		Total silicon (Si)
		<200 ppm
		Total sulphur (S)
		<500 ppm
		Zinc (Zn)
		<50 ppm

Potassium dihydrogen phosphate

See Potassium dihydrogen orthophosphate, page 252

Potassium fluoride anhydrous

for analysis

Certified AR

packaging		Mfr. No
500 g	Plastic jar	P/5010/53
KF		EINECS: 232-151-5
CAS: 7789-23-3		UN: 1812
MW: 58.09		GHS: Acute Tox. 3
Boiling point: 1505°C		H331, H311, H301
Melting point: 860°C		P301+P310, P280, P312,
ADR/RID: 6.1, III		P304+P340, P302+P352, P261



Assay: >99%		
Guaranteed analysis		
Calcium (Ca)	<10 ppm	n-Butyl formate
Copper (Cu)	<10 ppm	<0.005%
Iron (Fe)	<10 ppm	Total phosphorus (P)
Lead (Pb)	<10 ppm	<300 ppm
Magnesium (Mg)	<10 ppm	Total silicon (Si)
Sodium (Na)	<2000 ppm	<50 ppm
		Total sulphur (S)
		<50 ppm
		Zinc (Zn)
		<20 ppm

Typical analysis		
Aluminium (Al)	<5 ppm	Molybdenum (Mo)
Barium (Ba)	<5 ppm	<5 ppm
Cadmium (Cd)	<5 ppm	Nickel (Ni)
Chromium (Cr)	<5 ppm	<5 ppm
Cobalt (Co)	<5 ppm	Silver (Ag)
Lithium (Li)	<5 ppm	<5 ppm
Manganese (Mn)	<5 ppm	Strontium (Sr)
		<5 ppm
		Titanium (Ti)
		<5 ppm
		Vanadium (V)
		<5 ppm

Lot analysis on the label

Potassium fluoride anhydrous

specified reagents for general laboratory work- extra pure

SLR

packaging		Mfr. No
500 g	Plastic jar	P/5000/53
KF		EINECS: 232-151-5
CAS: 7789-23-3		UN: 1812
MW: 58.09		GHS: Acute Tox. 3
Boiling point: 1505°C		H331, H311, H301
Melting point: 860°C		P301+P310, P280, P312,
ADR/RID: 6.1, III		P304+P340, P302+P352, P261



Assay: >97%		
Guaranteed analysis		
Calcium (Ca)	<200 ppm	n-Butyl formate
Copper (Cu)	<200 ppm	<0.05%
Iron (Fe)	<100 ppm	Total phosphorus (P)
Lead (Pb)	<100 ppm	<500 ppm
Magnesium (Mg)	<100 ppm	Total silicon (Si)
Sodium (Na)	<5000 ppm	<5000 ppm
		Total sulphur (S)
		<500 ppm
		Zinc (Zn)
		<200 ppm

Potassium fluoride dihydrate

specified reagents for general laboratory work- extra pure

SLR

packaging		Mfr. No
500 g	Plastic jar	P/4980/53
KF.2H ₂ O		UN: 1812
CAS: 13455-21-5		GHS: Acute Tox. 3
MW: 94.13		H301, H331, H311
Boiling point: 156°C		P304+P340, P301+P310,
Melting point: 41°C		P280, P312, P302+P350
ADR/RID: 6.1, III		



Assay: 99 to 102%		
Guaranteed analysis		
Iron (Fe)	<200 ppm	Sulfate (SO ₄)
Lead (Pb)	<50 ppm	<0.2%
Sodium (Na)	<5000 ppm	n-Butyl formate
		<0.02%

Potassium hexachloroplatinate(IV)

specified reagents for general laboratory work- extra pure

SLR

packaging		Mfr. No
1 g	Vial	P/4320/43
K ₂ PtCl ₆		GHS: Acute Tox. 3, Eye
CAS: 16921-30-5		Dam. 1, Resp. Sens. 1, Skin
MW: 486.01		Sens. 1
Melting point: 250°C		H334, H301, H318, H317
ADR/RID: 8, 6.1, III		P301+P310,
EINECS: 240-979-3		P305+P351+P338, P261,
UN: 2923		P304+P341, P342+P311, P280



Guaranteed analysis	
Identity test	: Pass

Potassium hexacyanoferrate(II)

for analysis

Certified AR

packaging		Mfr. No
500 g	Plastic jar	P/4960/53
1 kg	Plastic jar	P/4960/60
25 kg	Plastic keg	P/4960/70
K ₄ Fe(CN) ₆ .3H ₂ O		EINECS: 237-722-2
CAS: 14459-95-1		GHS: Aquatic Chronic 3
MW: 422.41		H412
Density at 20°C: 1.85		P273
Melting point: 70°C		

Assay: 99 to 101%		
Guaranteed analysis		
Aluminium (Al)	<5 ppm	Nickel (Ni)
Barium (Ba)	<5 ppm	<5 ppm
Cadmium (Cd)	<5 ppm	Propan--1-ol
Calcium (Ca)	<50 ppm	<5 ppm
Chromium (Cr)	<5 ppm	Sodium (Na)
Cobalt (Co)	<5 ppm	<200 ppm
Copper (Cu)	<20 ppm	Strontium (Sr)
Lead (Pb)	<20 ppm	<5 ppm
Lithium (Li)	<5 ppm	Titanium (Ti)
Magnesium (Mg)	<20 ppm	<5 ppm
Manganese (Mn)	<5 ppm	n-Butyl formate
Molybdenum (Mo)	<5 ppm	<0.02%
		Total phosphorus (P)
		<50 ppm
		Total silicon (Si)
		<50 ppm
		Total sulphur (S)
		<50 ppm
		Vanadium (V)
		<5 ppm
		Zinc (Zn)
		<20 ppm

Potassium hexacyanoferrate(II)

specified reagents for general laboratory work- extra pure

SLR

packaging		Mfr. No
500 g	Plastic jar	P/4920/53
25 kg	Plastic keg	P/4920/70
C ₆ FeK ₄ N ₆ .3H ₂ O		EINECS: 237-722-2
CAS: 14459-95-1		GHS: Aquatic Chronic 3
MW: 422.41		H412
Density at 20°C: 1.85		P273

Potassium hexacyanoferrate(III) for analysis		Certified AR	
packaging		Mfr. No	
250 g	Plastic jar	P/4880/50	
500 g	Plastic jar	P/4880/53	
1 kg	Plastic jar	P/4880/60	
25 kg	Plastic keg	P/4880/70	
K ₃ Fe(CN) ₆		EINECS: 237-323-3	
CAS: 13746-66-2		EUH032	
MW: 329.26		P233	
Assay: >99%			
Guaranteed analysis			
Insoluble matter	Pass	Nickel (Ni)	<5 ppm
Aluminium (Al)	<5 ppm	Propan--1-ol	<5 ppm
Barium (Ba)	<5 ppm	Sodium (Na)	<200 ppm
Cadmium (Cd)	<5 ppm	Strontium (Sr)	<5 ppm
Calcium (Ca)	<50 ppm	Sulfate (SO ₄)	<0.01%
Chromium (Cr)	<5 ppm	Titanium (Ti)	<5 ppm
Cobalt (Co)	<5 ppm	n-Butyl formate	<0.01%
Copper (Cu)	<20 ppm	Total phosphorus (P)	<50 ppm
Iron cyanide	<0.02%	Total silicon (Si)	<50 ppm
Lead (Pb)	<20 ppm	Total sulphur (S)	<50 ppm
Lithium (Li)	<5 ppm	Vanadium (V)	<5 ppm
Magnesium (Mg)	<20 ppm	Zinc (Zn)	<50 ppm
Manganese (Mn)	<5 ppm	Ferrocyanide	<0.02%
Molybdenum (Mo)	<5 ppm		

Lot analysis on the label.

Potassium hexacyanoferrate(III) <i>specified reagents for general laboratory work- extra pure</i>			SLR
packaging			Mfr. No
250 g	Plastic jar		P/4840/50
500 g	Plastic jar		P/4840/53
C ₆ FeK ₃ N ₆		EINECS: 237-323-3	
CAS: 13746-66-2		EUH032	
MW: 329.26		P233	
Assay: >99%			
Guaranteed analysis			
Calcium (Ca)	<500 ppm	n-Butyl formate	<0.2%
Copper (Cu)	<200 ppm	Total phosphorus (P)	<500 ppm
Lead (Pb)	<200 ppm	Total silicon (Si)	<500 ppm
Magnesium (Mg)	<200 ppm	Total sulphur (S)	<500 ppm
Sodium (Na)	<2000 ppm	Zinc (Zn)	<200 ppm

Potassium hydrogen carbonate		Certified AR	
for analysis			
packaging		Mfr. No	
500 g	Plastic jar	P/5120/53	
1 kg	Plastic jar	P/5120/60	
KHCO ₃		Melting point: 292°C	
CAS: 298-14-6		EINECS: 206-059-0	
MW: 100.12			
Assay: >99.8%			
Guaranteed analysis			
Calcium (Ca)	<20 ppm	Total chloride (Cl)	<0.001%
Copper (Cu)	<2 ppm	Total nitrogen (N)	<10 ppm
Iron (Fe)	<5 ppm	Total phosphorus (P)	<10 ppm
Lead (Pb)	<2 ppm	Total silicon (Si)	<10 ppm
Magnesium (Mg)	<10 ppm	Total sulphur (S)	<50 ppm
Sodium (Na)	<200 ppm	Zinc (Zn)	<20 ppm
Typical analysis			
Aluminium (Al)	<5 ppm	Molybdenum (Mo)	<5 ppm
Barium (Ba)	<5 ppm	Nickel (Ni)	<5 ppm
Cadmium (Cd)	<5 ppm	Silver (Ag)	<5 ppm
Chromium (Cr)	<5 ppm	Strontium (Sr)	<5 ppm
Cobalt (Co)	<5 ppm	Titanium (Ti)	<5 ppm
Lithium (Li)	<5 ppm	Vanadium (V)	<5 ppm
Manganese (Mn)	<5 ppm		

Lot analysis on the label.

Potassium hydrogen carbonate		SLR, Eur.Ph., BP, USP	
specified reagents for general laboratory work- extra pure			
Packaging		Mfr. No.	
500 g	Plastic jar	P/5080/53	
2.5 kg	Plastic jar	P/5080/62	
KHCO ₃		Melting point: 292°C	
CAS: 298-14-6		EINECS: 206-059-0	
MW: 100.12			
Assay: >99%			
Guaranteed analysis			
Calcium (Ca)	<100 ppm	n-Butyl formate	<0.01%
Copper (Cu)	<10 ppm	Total phosphorus (P)	<50 ppm
Iron (Fe)	<20 ppm	Total silicon (Si)	<50 ppm
Lead (Pb)	<20 ppm	Total sulphur (S)	<200 ppm
Magnesium (Mg)	<50 ppm	Zinc (Zn)	<100 ppm
Sodium (Na)	<2000 ppm		

di-Potassium hydrogen orthophosphate anhydrous for analysis		Certified AR
Packaging		Mfr. No
500 g	Plastic jar	P/5245/53
1 kg	Plastic jar	P/5245/60
25 kg	Fibre Keg	P/5245/70

Lot analysis on the label.

di-Potassium hydrogen orthophosphate anhydrous specified reagents for general laboratory work- extra pure		SLR	
--	--	-----	--

Packaging		Mfr. No	
500 g	Plastic jar	P/5240/53	
1 kg	Plastic jar	P/5240/60	
5 kg	Plastic keg	P/5240/65	
10 kg	Plastic keg	P/5240/68	
K ₂ HPO ₄		Melting point: 340°C	
CAS: 7758-11-4		EINECS: 231-834-5	
MW: 174.18			
Assay: >99%			
Guaranteed analysis			
Loss on drying	<0.5%	Magnesium (Mg)	<200 ppm
Calcium (Ca)	<200 ppm	n-Butyl formate	<0.02%
Copper (Cu)	<20 ppm	Total silicon (Si)	<200 ppm
Iron (Fe)	<20 ppm	Total sulphur (S)	<200 ppm
Lead (Pb)	<20 ppm	Zinc (Zn)	<100 ppm

pH Buffers

NIST Certified Standard
Solution & Concentrate

A full range to discover on page 63



di-Potassium hydrogen orthophosphate trihydrate for analysis		Certified AR
packaging		Mfr. No
500 g	Plastic jar	P/5250/53
1 kg	Plastic jar	P/5250/60
25 kg	Plastic keg	P/5250/70
K ₂ HPO ₄ ·3H ₂ O		
CAS: 16788-57-1		
MW: 228.22		
EINECS: 231-834-5		
Assay: 99 to 101%		
Guaranteed analysis		
pH (aqueous solution 5%)	9 to 9.2	Sodium (Na) <500 ppm
Calcium (Ca)	<50 ppm	Total chloride (Cl) <0.001%
Copper (Cu)	<5 ppm	Total nitrogen (N) <10 ppm
Iron (Fe)	<5 ppm	Total silicon (Si) <50 ppm
Lead (Pb)	<5 ppm	Total sulphur (S) <50 ppm
Magnesium (Mg)	<20 ppm	Zinc (Zn) <20 ppm
Typical analysis		
Aluminium (Al)	<5 ppm	Manganese (Mn) <2 ppm
Arsenic (As)	<2 ppm	Molybdenum (Mo) <2 ppm
Barium (Ba)	<2 ppm	Nickel (Ni) <2 ppm
Cadmium (Cd)	<2 ppm	Silver (Ag) <2 ppm
Chromium (Cr)	<2 ppm	Strontium (Sr) <5 ppm
Cobalt (Co)	<2 ppm	Titanium (Ti) <2 ppm
Lithium (Li)	<2 ppm	Vanadium (V) <2 ppm

Lot analysis on the label.

di-Potassium hydrogen orthophosphate trihydrate specified reagents for general laboratory work- extra pure		SLR	
---	--	-----	--

packaging		Mfr. No
500 g	Plastic jar	P/5248/53
1 kg	Plastic jar	P/5248/60
K ₂ HPO ₄ ·3H ₂ O		
CAS: 16788-57-1		
MW: 228.22		
EINECS: 231-834-5		
Assay: >99%		
Guaranteed analysis		
Calcium (Ca)	<200 ppm	n-Butyl formate <0.02%
Copper (Cu)	<20 ppm	Total silicon (Si) <200 ppm
Iron (Fe)	<20 ppm	Total sulphur (S) <200 ppm
Lead (Pb)	<20 ppm	Zinc (Zn) <100 ppm
Magnesium (Mg)	<200 ppm	

Potassium hydrogen phthalate for analysis		Certified AR	
packaging		Mfr. No	
100 g	Plastic jar	P/5320/48	
250 g	Plastic jar	P/5320/50	
500 g	Plastic jar	P/5320/53	
1 kg	Plastic jar	P/5320/60	
25 kg	Fibre Keg	P/5320/70	
C ₈ H ₅ KO ₄		Melting point: 295 to 300°C	
CAS: 877-24-7		EINECS: 212-889-4	
MW: 204.22			
Assay: 99.9 to 100.1%			
Guaranteed analysis			
Loss on drying	<0.1%	Sodium (Na)	<200 ppm
Calcium (Ca)	<20 ppm	Total chloride (Cl)	<0.002%
Copper (Cu)	<5 ppm	Total phosphorus (P)	<50 ppm
Iron (Fe)	<5 ppm	Total silicon (Si)	<20 ppm
Lead (Pb)	<5 ppm	Total sulphur (S)	<50 ppm
Magnesium (Mg)	<10 ppm	Zinc (Zn)	<10 ppm
Typical analysis			
Aluminium (Al)	<2 ppm	Molybdenum (Mo)	<2 ppm
Barium (Ba)	<2 ppm	Nickel (Ni)	<2 ppm
Cadmium (Cd)	<2 ppm	Silver (Ag)	<2 ppm
Chromium (Cr)	<2 ppm	Strontium (Sr)	<2 ppm
Cobalt (Co)	<2 ppm	Titanium (Ti)	<2 ppm
Lithium (Li)	<2 ppm	Vanadium (V)	<2 ppm
Manganese (Mn)	<2 ppm		

Lot analysis on the label.

Potassium hydrogen sulfate for analysis			Certified AR
Packaging		Mfr. No	
250 g	Plastic jar	P/5440/50	
500 g	Plastic jar	P/5440/53	
1 kg	Plastic jar	P/5440/60	
KHSO ₄		EINECS: 231-594-1	
CAS: 7646-93-7		UN: 2509	
MW: 136.17		GHS: Skin Corr. 1B, STOT SE 3	
Density at 20°C: 2.24		H314, H335	
Melting point: 214°C		P301+P330+P331, P280,	
ADR/RID: 8, II		P305+P351+P338, P310, P261	
			
Assay: >99%			
Guaranteed analysis			
Calcium (Ca)	<50 ppm	Total chloride (Cl)	<0.0005%
Copper (Cu)	<5 ppm	Total nitrogen (N)	<10 ppm
Iron (Fe)	<5 ppm	Total phosphorus (P)	<10 ppm
Lead (Pb)	<5 ppm	Total silicon (Si)	<20 ppm
Magnesium (Mg)	<20 ppm	Zinc (Zn)	<20 ppm
Sodium (Na)	<100 ppm		
Typical analysis			
Aluminium (Al)	<5 ppm	Molybdenum (Mo)	<5 ppm
Barium (Ba)	<5 ppm	Nickel (Ni)	<5 ppm
Cadmium (Cd)	<5 ppm	Silver (Ag)	<5 ppm
Chromium (Cr)	<5 ppm	Strontium (Sr)	<5 ppm
Cobalt (Co)	<5 ppm	Titanium (Ti)	<5 ppm
Lithium (Li)	<5 ppm	Vanadium (V)	<5 ppm
Manganese (Mn)	<5 ppm		

Lot analysis on the label.

Potassium hydrogen sulfate		SLR	
specified reagents for general laboratory work- extra pure			
Packaging		Mfr. No	
500 g	Plastic jar	P/5420/53	
KHSO ₄	EINECS: 231-594-1		
CAS: 7646-93-7	UN: 2509		
MW: 136.17	GHS: Skin Corr. 1B, STOT SE 3		
Density at 20°C: 2.24	H314, H335		
Melting point: 214°C	P301+P330+P331, P280,		
ADR/RID: 8, II	P305+P351+P338, P310, P261		
Assay: >98%			
Guaranteed analysis			
Calcium (Ca)	<200 ppm	Sodium (Na)	<5000 ppm
Copper (Cu)	<20 ppm	n-Butyl formate	<0.002%
Iron (Fe)	<20 ppm	Total phosphorus (P)	<50 ppm
Lead (Pb)	<20 ppm	Total silicon (Si)	<100 ppm
Magnesium (Mg)	<100 ppm	Zinc (Zn)	<100 ppm

Potassium hydrogen tartrate		SLR
specified reagents for general laboratory work- extra pure		
Packaging		Mfr. No
1 kg	Plastic jar	P/5480/60
C ₄ H ₅ KO ₆	Density at 20°C: 1.95	
CAS: 868-14-4	EINECS: 212-769-1	
MW: 188.18		
Assay: > 97%		
Guaranteed analysis		
Complies with standard IR spectrum: Pass		

Potassium hydroxide pellets
for analysis

Certified AR, Eur.Ph., BP

packaging		Mfr. No
500 g	Plastic jar	P/5640/53
1 kg	Plastic jar	P/5640/60
5 kg	Plastic jar	P/5640/65
10 kg	Plastic keg	P/5640/68
25 kg	Plastic keg	P/5640/70
KOH		UN: 1813
CAS: 1310-58-3		GHS: Acute Tox. 4, Skin
MW: 56.11		Corr. 1A, Eye Dam. 1
Boiling point: 1320°C		H314, H302
Melting point: 360°C		P301+P330+P331, P280,
ADR/RID: 8, II		P305+P351+P338, P301+P312
EINECS: 215-181-3		

Assay: >86%

Guaranteed analysis			
Calcium (Ca)	<10 ppm	Total chloride (Cl)	<0.002%
Carbonate (Na ₂ CO ₃)	<1%	Total nitrogen (N)	<10 ppm
Copper (Cu)	<5 ppm	Total phosphorus (P)	<5 ppm
Iron (Fe)	<5 ppm	Total silicon (Si)	<50 ppm
Lead (Pb)	<5 ppm	Total sulphur (S)	<20 ppm
Magnesium (Mg)	<5 ppm	Zinc (Zn)	<10 ppm
Sodium (Na)	<5000 ppm		

Typical analysis			
Aluminium (Al)	<5 ppm	Manganese (Mn)	<5 ppm
Arsenic (As)	<0.5 ppm	Molybdenum (Mo)	<5 ppm
Barium (Ba)	<5 ppm	Nickel (Ni)	<5 ppm
Cadmium (Cd)	<5 ppm	Silver (Ag)	<5 ppm
Chromium (Cr)	<5 ppm	Strontium (Sr)	<5 ppm
Cobalt (Co)	<5 ppm	Titanium (Ti)	<5 ppm
Lithium (Li)	<5 ppm	Vanadium (V)	<5 ppm

Lot analysis on the label.

Potassium hydroxide pellets
specified reagents for general laboratory work- extra pure

SLR

packaging		Mfr. No
500 g	Plastic jar	P/5600/53
1 kg	Plastic jar	P/5600/60
2 kg	Plastic jar	P/5600/61
5 kg	Plastic jar	P/5600/65
10 kg	Plastic keg	P/5600/68
KOH		UN: 1813
CAS: 1310-58-3		GHS: Acute Tox. 4, Skin
MW: 56.11		Corr. 1A, Eye Dam. 1
Boiling point: 1320°C		H314, H302
Melting point: 360°C		P301+P330+P331, P280,
ADR/RID: 8, II		P305+P351+P338, P301+P312
EINECS: 215-181-3		

Assay: >85%

Guaranteed analysis			
Calcium (Ca)	<50 ppm	Sodium (Na)	<10000 ppm
Carbonate (Na ₂ CO ₃)	<2%	n-Butyl formate	<0.01%
Copper (Cu)	<20 ppm	Total phosphorus (P)	<200 ppm
Iron (Fe)	<20 ppm	Total silicon (Si)	<200 ppm
Lead (Pb)	<20 ppm	Total sulphur (S)	<100 ppm
Magnesium (Mg)	<20 ppm	Zinc (Zn)	<50 ppm

High Purity Acids
for Elemental Analysis

Fisher Chemical offers a complete range of acids to meet the most challenging applications. Choose from our Optima®, TraceMetal® and PrimarPlus® grades to ensure that the sensitivity you require is always available.

Read more page 24

Potassium hydroxide flake
specified reagents for general laboratory work- extra pure

SLR

packaging		Mfr. No
1 kg	Plastic jar	P/5560/60
KOH		UN: 1813
CAS: 1310-58-3		GHS: Acute Tox. 4, Skin
MW: 56.11		Corr. 1A, Eye Dam. 1
Boiling point: 1320°C		H314, H302
Melting point: 360°C		P301+P330+P331, P280,
ADR/RID: 8, II		P305+P351+P338, P301+P312
EINECS: 215-181-3		

Assay: >85%

Guaranteed analysis			
Calcium (Ca)	<50 ppm	Sodium (Na)	<10000 ppm
Carbonate (Na ₂ CO ₃)	<2%	n-Butyl formate	<0.01%
Copper (Cu)	<20 ppm	Total phosphorus (P)	<200 ppm
Iron (Fe)	<20 ppm	Total silicon (Si)	<200 ppm
Lead (Pb)	<20 ppm	Total sulphur (S)	<100 ppm
Magnesium (Mg)	<20 ppm	Zinc (Zn)	<50 ppm

Potassium hydroxide solution 48-50 % liquor
pure

packaging		Mfr. No
2.5 ℓ	HDPE Plastic bottle	P/5645/17
KOH		UN: 1814
CAS: 1310-58-3		GHS: Acute Tox. 4, Skin Corr. 1A
MW: 56.11		H318, H302
Density at 20°C: 1.51		P280, P301+P330+P331,
Melting point: -28.8°C		P310, P305+P351+P338,
ADR/RID: 8, II		P303+P361+P353
EINECS: 215-181-3		

Assay: 48 to 50%

Potassium hydroxide solution 0.1M (0.1N), (in methanol)

Eur.Ph. , BP

NIST Standard solution ready to use, For Volumetric analysis

packaging		Mfr. No
1 ℓ	HDPE Plastic bottle	J/6625/15
CAS: 1310-58-3		Tox. 3, STOT SE 1
Density at 20°C: .7910		H225, H311, H370, H301, H331
Flash point: 12°C		P301+P310, P280, P312,
ADR/RID: 3, 6.1, II		P302+P350, P307+P311,
EINECS: 215-181-3		P304+P340,
UN: 1992		P305+P351+P338, P210, P240
GHS: Flam. Liq. 2, Acute		

Manufactured to a tolerance of ± 0.5%.
See the full range of Standard Volumetric solution page 323.

Potassium hydroxide solution 0.1M (0.1N), (in ethanol)

NIST Standard solution ready to use, For Volumetric analysis

packaging		Mfr. No
1 ℓ	HDPE Plastic bottle	J/6620/15
2.5 ℓ	HDPE Plastic bottle	J/6620/17
Density at 20°C: 0.8		2, STOT SE 2
Flash point: 16.66°C		H225, H302, H312, H315,
ADR/RID: 3, 8, II		H319, H332, H371
EINECS: 200-578-6		P280, P301+P312,
UN: 2924		P305+P351+P338,
GHS: Flam. Liq. 2, Acute		P304+P341, P261, P210
Tox. 4, Skin Irrit. 2, Eye Irrit.		

Manufactured to a tolerance of ± 0.2%.
See the full range of Standard Volumetric solution page 323.

Potassium hydroxide solution 0.1M (0.1N)
NIST Standard Concentrate, For Volumetric analysis

packaging		Mfr. No
1	Ampoule	J/6640C/05
6	Ampoule	J/6640C/90
CAS: 1310-58-3		GHS: Skin Corr. 1A, Eye Dam. 1
ADR/RID: 8, II		H314
EINECS: 215-181-3		P280, P301+P330+P331,
UN: 1814		P305+P351+P338, P310

Each ampoule makes 1 ℓ of solution.
See the full range of Standard Volumetric solution page 323.

Potassium hydroxide solution 0.5M (0.5N), (in ethanol)

NIST Standard solution ready to use, For Volumetric analysis

packaging		Mfr. No
1 ℓ	HDPE Plastic bottle	J/6612/15
2.5 ℓ	HDPE Plastic bottle	J/6612/17
Flash point: 16.66°C		H225, H314, H371, H332,
ADR/RID: 3, 8, II		H312, H302
EINECS: 200-578-6		P301+P330+P331, P280,
UN: 2924		P305+P351+P338, P310,
GHS: Flam. Liq. 2, Acute		P260, P304+P340, P210, P240
Tox. 4, Skin Corr. 1B, STOT SE 2		

Manufactured to a tolerance of ± 0.2%.
see also the full range of "Standard solution" in the Standard solution section.
See the full range of Standard Volumetric solution page 323.

Potassium hydroxide solution 1M (1N)

Eur.Ph. , BP

NIST Standard solution ready to use, For Volumetric analysis

packaging		Mfr. No
1 ℓ	HDPE Plastic bottle	J/6630/PB15
2.5 ℓ	HDPE Plastic bottle	J/6630/17
CAS: 1310-58-3		GHS: Skin Corr. 1A, Eye Dam. 1
Density at 20°C: 1.05		H314
ADR/RID: 8, II		P280, P301+P330+P331,
EINECS: 215-181-3		P305+P351+P338, P310
UN: 1814		

Manufactured to a tolerance of ± 0.1%.
See the full range of Standard Volumetric solution page 323.

Potassium hydroxide solution 1M (1N)
NIST Standard Concentrate, For Volumetric analysis

packaging		Mfr. No
1	Ampoule	J/6630C/05
6	Ampoule	J/6630C/90
ADR/RID: 8, II		H314
EINECS: 215-181-3		P301+P330+P331, P280,
UN: 1814		P305+P351+P338, P310
GHS: Skin Corr. 1A, Eye Dam. 1		

Each ampoule makes 1 ℓ of solution.
See the full range of Standard Volumetric solution page 323.

UHPLC

Filtered to 0.1 µm, Fisher Chemical UHPLC gradient grade solvents (Acetonitrile, Methanol and Water) are perfect for those working at high pressure.

See page 309 for more information

Potassium hydroxide solution 1M (1N), (alcoholic)
NIST Standard solution ready to use, For Volumetric analysis

packaging		Mfr. No
1 ℓ	HDPE Plastic bottle	J/6610/15
2.5 ℓ	HDPE Plastic bottle	J/6610/17
Density at 20°C: 0.8		Dam. 1, STOT SE 2
Flash point: 19°C		H225, H314, H371, H332,
ADR/RID: 3, 8, II		H312, H302
EINECS: 200-578-6		P301+P330+P331, P280,
UN: 2924		P305+P351+P338, P310,
GHS: Flam. Liq. 2, Acute		P260, P304+P340, P210, P240
Tox. 4, Skin Corr. 1A, Eye		

Manufactured to a tolerance of ± 0.1%.
See the full range of Standard Volumetric solution page 323.

Potassium iodate
for analysis

Certified AR

packaging		Mfr. No
100 g	Plastic jar	P/5800/48
250 g	Plastic jar	P/5800/50
500 g	Plastic jar	P/5800/53
1 kg	Plastic jar	P/5800/60

KIO ₃	GHS: Ox. Sol. 2, Acute Tox.
CAS: 7758-05-6	4, Skin Irrit. 2, Eye Irrit. 2,
MW: 214.00	STOT SE 3
Density at 20°C: 3.93	H272, H315, H302, H335, H319
Melting point: 560°C	P210, P261, P301+P312,
ADR/RID: 5.1, II	P302+P352, P280,
EINECS: 231-831-9	P305+P351+P338
UN: 1479	

Assay: >99.9%

Guaranteed analysis			
pH (aqueous solution 5%)	5 to 7	Lead (Pb)	<10 ppm
Loss on drying	<0.05%	Magnesium (Mg)	<5 ppm
Bromate, bromide, chlorate & chloride: < 0.02%		Sodium (Na)	<100 ppm
Calcium (Ca)	<200 ppm	Total nitrogen (N)	<20 ppm
Copper (Cu)	<10 ppm	Total silicon (Si)	<10 ppm
Iodide (I)	<0.001%	Total sulphur (S)	<100 ppm
Iron (Fe)	<10 ppm	Zinc (Zn)	<10 ppm

Typical analysis			
Aluminium (Al)	<5 ppm	Molybdenum (Mo)	<5 ppm
Barium (Ba)	<5 ppm	Nickel (Ni)	<5 ppm
Cadmium (Cd)	<5 ppm	Silver (Ag)	<5 ppm
Chromium (Cr)	<5 ppm	Strontium (Sr)	<5 ppm
Cobalt (Co)	<5 ppm	Titanium (Ti)	<2 ppm
Lithium (Li)	<5 ppm	Vanadium (V)	<2 ppm
Manganese (Mn)	<5 ppm		

Lot analysis on the label.

Label

Fisher Chemical labels show GHS pictograms clearly. Together with actual lot analysis shown on the label, plus a wealth of additional information, we help you to work safely to achieve accurate and consistent results in your analytical applications.



Potassium iodate
specified reagents for general laboratory work- extra pure

SLR

packaging		Mfr. No
100 g	Plastic jar	P/5760/48
500 g	Plastic jar	P/5760/53
25 kg	Plastic keg	P/5760/70
KIO ₃ CAS: 7758-05-6 MW: 214.00 Density at 20°C: 3.93 Melting point: 560°C ADR/RID: 5.1, II EINECS: 231-831-9 UN: 1479		GHS: Ox. Sol. 2, Acute Tox. 4, Skin Irrit. 2, Eye Irrit. 2, STOT SE 3 H272, H315, H302, H335, H319 P210, P261, P301+P312, P302+P352, P280, P305+P351+P338
Assay: >99.5%		
Guaranteed analysis		
Loss on drying	<0.2%	Lead (Pb) <50 ppm
Bromate, bromide, chlorate and chloride	<0.1%	Magnesium (Mg) <20 ppm
Calcium (Ca)	<500 ppm	Sodium (Na) <500 ppm
Copper (Cu)	<50 ppm	Total silicon (Si) <50 ppm
Iodide (I)	<0.005%	Total sulphur (S) <500 ppm
Iron (Fe)	<50 ppm	Zinc (Zn) <50 ppm



Potassium iodate 0.025M solution stabilised with 10 ppm mercury(II) chloride
NIST Standard solution ready to use, For Volumetric analysis

packaging		Mfr. No
2.5 ℓ	Glass bottle	J/6700/17
Guaranteed analysis		
Value between	0.024	

Manufactured to a tolerance of ± 0.1%.
See the full range of Standard Volumetric solution page 323.

Potassium iodate 0.05M solution stabilised with 10 ppm mercury(II) chloride
Eur.Ph. , USP, BP
NIST Standard solution ready to use, For Volumetric analysis

packaging		Mfr. No
1 ℓ	Glass bottle	J/6690/15
2.5 ℓ	Glass bottle	J/6690/17
P262		

Manufactured to a tolerance of ± 0.1%.
See the full range of Standard Volumetric solution page 323.

Potassium iodide
for analysis

Certified AR

packaging		Mfr. No
100 g	Plastic jar	P/5880/48
250 g	Plastic jar	P/5880/50
500 g	Plastic jar	P/5880/53
1 kg	Plastic jar	P/5880/60
3 kg	Plastic jar	P/5880/63
25 kg	Plastic keg	P/5880/70
KI CAS: 7681-11-0 MW: 166.01		Boiling point: 1330°C Melting point: 680°C EINECS: 231-659-4
Assay: >99.9%		
Guaranteed analysis		
Water insoluble matter	<0.005%	Lead (Pb) <2 ppm
pH (aqueous solution 5%)	6 to 8	Magnesium (Mg) <5 ppm
Loss on drying	<0.2%	Sodium (Na) <500 ppm
Heavy metals	<10 ppm	Total chloride (Cl) <0.1%
Iodine reducing substances	<0.001%	Total nitrogen (N) <10 ppm
Calcium (Ca)	<5 ppm	Total phosphorus (P) <10 ppm
Chloride and Bromide	<0.01%	Total silicon (Si) <10 ppm
Copper (Cu)	<2 ppm	Total sulphur (S) <20 ppm
Iodate (IO ₃)	<2 ppm	Zinc (Zn) <10 ppm
Iron (Fe)	<2 ppm	

Typical analysis		
Aluminium (Al)	<2 ppm	Manganese (Mn) <2 ppm
Arsenic (As)	<0.1 ppm	Molybdenum (Mo) <2 ppm
Barium (Ba)	<2 ppm	Nickel (Ni) <2 ppm
Cadmium (Cd)	<2 ppm	Silver (Ag) <2 ppm
Chromium (Cr)	<2 ppm	Strontium (Sr) <5 ppm
Cobalt (Co)	<2 ppm	Titanium (Ti) <2 ppm
Lithium (Li)	<2 ppm	Vanadium (V) <2 ppm

Lot analysis on the label.

Potassium iodide
specified reagents for general laboratory work- extra pure

SLR

packaging		Mfr. No
100 g	Plastic jar	P/5840/48
250 g	Plastic jar	P/5840/50
500 g	Plastic jar	P/5840/53
1 kg	Plastic jar	P/5840/60
5 kg	Plastic jar	P/5840/65
10 kg	Plastic keg	P/5840/68
KI CAS: 7681-11-0 MW: 166.01		Boiling point: 1330°C Melting point: 680°C EINECS: 231-659-4
Assay: >99.5%		
Guaranteed analysis		
Loss on drying	<0.5%	Magnesium (Mg) <20 ppm
Calcium (Ca)	<20 ppm	Sodium (Na) <2000 ppm
Copper (Cu)	<10 ppm	Total phosphorus (P) <50 ppm
Iodate (IO ₃)	<10 ppm	Total silicon (Si) <50 ppm
Iron (Fe)	<10 ppm	Total sulphur (S) <100 ppm
Lead (Pb)	<10 ppm	Zinc (Zn) <50 ppm

Potassium iron(III)cyanide

See Potassium hexacyanoferrate (III), page 254

Potassium metabisulfite
for analysis

Certified AR

packaging		Mfr. No
500 g	Plastic jar	P/5960/53
1 kg	Plastic jar	P/5960/60
K ₂ O ₅ S ₂ CAS: 16731-55-8 MW: 222.33 Melting point: 190°C EINECS: 240-795-3		GHS: Eye Dam. 1, STOT SE 3 H318, H335, EUH031 P280, P305+P351+P338, P310, P261, P233
Assay: >96%		
Guaranteed analysis		
Calcium (Ca)	<50 ppm	Sodium (Na) <1000 ppm
Copper (Cu)	<10 ppm	Total chloride (Cl) <0.02%
Iron (Fe)	<10 ppm	Total phosphorus (P) <20 ppm
Lead (Pb)	<10 ppm	Total silicon (Si) <10 ppm
Magnesium (Mg)	<20 ppm	Zinc (Zn) <10 ppm
Typical analysis		
Aluminium (Al)	<5 ppm	Manganese (Mn) <5 ppm
Arsenic (As)	<0.2 ppm	Molybdenum (Mo) <5 ppm
Barium (Ba)	<5 ppm	Nickel (Ni) <5 ppm
Cadmium (Cd)	<5 ppm	Silver (Ag) <2 ppm
Chromium (Cr)	<5 ppm	Strontium (Sr) <5 ppm
Cobalt (Co)	<5 ppm	Titanium (Ti) <5 ppm
Lithium (Li)	<5 ppm	Vanadium (V) <5 ppm

Lot analysis on the label.

Potassium metabisulfite
specified reagents for general laboratory work- extra pure

SLR

packaging		Mfr. No
500 g	Plastic jar	P/5920/53
1 kg	Plastic jar	P/5920/60
K ₂ O ₅ S ₂ CAS: 16731-55-8 MW: 222.33 Melting point: 190°C EINECS: 240-795-3		GHS: Eye Dam. 1, STOT SE 3 H318, H335, EUH031 P280, P305+P351+P338, P310, P261, P233
Assay: >95%		
Guaranteed analysis		
Calcium (Ca)	<200 ppm	Sodium (Na) <5000 ppm
Copper (Cu)	<50 ppm	n-Butyl formate <0.1%
Iron (Fe)	<50 ppm	Total phosphorus (P) <50 ppm
Lead (Pb)	<50 ppm	Total silicon (Si) <50 ppm
Magnesium (Mg)	<100 ppm	Zinc (Zn) <50 ppm



Potassium nitrate
for analysis

Certified AR, Eur.Ph., BP

packaging		Mfr. No
250 g	Plastic jar	P/6120/50
500 g	Plastic jar	P/6120/53
1 kg	Plastic jar	P/6120/60
5 kg	Plastic jar	P/6120/65
KNO ₃ CAS: 7757-79-1 MW: 101.11 Boiling point: 400°C Melting point: 334°C ADR/RID: 5.1, III		EINECS: 231-818-8 UN: 1486 GHS: Ox. Sol. 3 H272 P210
Assay: 99.5 to 101%		
Guaranteed analysis		
Acidity / Alkalinity	Pass	Nitrite (NO ₂) <0.0001%
pH (aqueous solution 5%)	5 to 8	Sodium (Na) <500 ppm
Heavy metals	<10 ppm	Total chloride (Cl) <0.0005%
Ammonium (NH ₄)	<0.001%	Total phosphorus (P) <5 ppm
Calcium (Ca)	<20 ppm	Total silicon (Si) <20 ppm
Copper (Cu)	<1 ppm	Total sulphur (S) <20 ppm
Iodate (IO ₃)	<0.5 ppm	Zinc (Zn) <10 ppm
Iron (Fe)	<1 ppm	Sulfate (SO ₄) <150 ppm
Lead (Pb)	<1 ppm	Reducing substance Pass
Magnesium (Mg)	<10 ppm	



Typical analysis		
Aluminium (Al)	<5 ppm	Manganese (Mn) <2 ppm
Arsenic (As)	<2 ppm	Molybdenum (Mo) <2 ppm
Barium (Ba)	<5 ppm	Nickel (Ni) <2 ppm
Cadmium (Cd)	<2 ppm	Silver (Ag) <2 ppm
Chromium (Cr)	<2 ppm	Strontium (Sr) <2 ppm
Cobalt (Co)	<2 ppm	Titanium (Ti) <2 ppm
Lithium (Li)	<2 ppm	Vanadium (V) <2 ppm

Potassium nitrate
specified reagents for general laboratory work- extra pure

SLR

packaging		Mfr. No
500 g	Plastic jar	P/6080/53
1 kg	Plastic jar	P/6080/60
2 kg	Plastic jar	P/6080/61
5 kg	Plastic jar	P/6080/65
KNO ₃ CAS: 7757-79-1 MW: 101.11 Boiling point: 400°C Melting point: 334°C ADR/RID: 5.1, III		EINECS: 231-818-8 UN: 1486 GHS: Ox. Sol. 3 H272 P210
Assay: >99%		



Guaranteed analysis		
Calcium (Ca)	<100 ppm	n-Butyl formate <0.01%
Copper (Cu)	<20 ppm	Total phosphorus (P) <20 ppm
Iron (Fe)	<20 ppm	Total silicon (Si) <100 ppm
Lead (Pb)	<20 ppm	Total sulphur (S) <100 ppm
Magnesium (Mg)	<50 ppm	Zinc (Zn) <50 ppm
Sodium (Na)	<2000 ppm	

Potassium nitrate

Technical

packaging		Mfr. No
1 kg	Plastic jar	P/6040/60
3 kg	Plastic jar	P/6040/63
KNO ₃ CAS: 7757-79-1 MW: 101.11 Boiling point: 400°C Melting point: 334°C ADR/RID: 5.1, III		EINECS: 231-818-8 UN: 1486 GHS: Ox. Sol. 3 H272 P210
Guaranteed analysis		
Identity test		: Pass



Potassium nitrite
for analysis

Certified AR

packaging		Mfr. No
250 g	Plastic jar	P/6162/50
KNO ₂ CAS: 7758-09-0 MW: 85.11 Density at 20°C: 19.15 Melting point: 387°C ADR/RID: 5.1, II		EINECS: 231-832-4 UN: 1488 GHS: Ox. Sol. 2, Acute Tox. 3, Aquatic Acute 1 H272, H301, H400 P301+P310, P210, P273
Assay: >99%		
Guaranteed analysis		
Aluminium (Al)	<5 ppm	Molybdenum (Mo) <5 ppm
Barium (Ba)	<5 ppm	Nickel (Ni) <5 ppm
Cadmium (Cd)	<5 ppm	Propan--1-ol <5 ppm
Calcium (Ca)	<50 ppm	Sodium (Na) <5000 ppm
Chromium (Cr)	<5 ppm	Strontium (Sr) <5 ppm
Cobalt (Co)	<5 ppm	Titanium (Ti) <5 ppm
Copper (Cu)	<5 ppm	n-Butyl formate <0.005%
Iron (Fe)	<5 ppm	Total phosphorus (P) <50 ppm
Lead (Pb)	<5 ppm	Total silicon (Si) <50 ppm
Lithium (Li)	<5 ppm	Total sulphur (S) <50 ppm
Magnesium (Mg)	<20 ppm	Vanadium (V) <5 ppm
Manganese (Mn)	<5 ppm	Zinc (Zn) <10 ppm

Lot analysis on the label.



Potassium nitrite		SLR
<i>specified reagents for general laboratory work- extra pure</i>		
packaging	Mfr. No	
500 g Plastic jar	P/6160/53	
KNO ₂ CAS: 7758-09-0 MW: 85.11 Density at 20°C: 19.15 Melting point: 387°C ADR/RID: 5.1, II	EINECS: 231-832-4 UN: 1488 GHS: Ox. Sol. 2, Acute Tox. 3, Aquatic Acute 1 H272, H301, H400 P301+P310, P210, P273	  
Assay: >97%		
Guaranteed analysis		
Calcium (Ca)	<500 ppm	n-Butyl formate
Copper (Cu)	<50 ppm	Total phosphorus (P)
Iron (Fe)	<50 ppm	Total silicon (Si)
Lead (Pb)	<50 ppm	Total sulphur (S)
Magnesium (Mg)	<200 ppm	Zinc (Zn)
Sodium (Na)	<10000 ppm	

tri-Potassium orthophosphate monohydrate		SLR
<i>specified reagents for general laboratory work- extra pure</i>		
packaging	Mfr. No	
500 g Plastic jar	P/6280/53	
K ₃ PO ₄ .H ₂ O CAS: 27176-10-9 MW: 230.29 Density at 20°C: 2.564 Melting point: 1340°C ADR/RID: 8, II	EINECS: 231-907-1 UN: 3262 GHS: Skin Corr. 1A, Eye Dam. 1 H314 P301+P330+P331, P280, P305+P351+P338, P310	
Assay: >95%		
Guaranteed analysis		
Lead (Pb)	<50 ppm	n-Butyl formate
Sulfate (SO ₄)	<0.005%	

Potassium oxalate monohydrate		Certified AR
<i>for analysis</i>		
packaging	Mfr. No	
250 g Plastic jar	P/6360/50	
500 g Plastic jar	P/6360/53	
25 kg Fibre Keg	P/6360/70	
C ₂ K ₂ O ₄ .H ₂ O CAS: 6487-48-5 MW: 184.23 Density at 20°C: 2.13 Melting point: 356°C ADR/RID: 6.1, III	EINECS: 209-506-8 UN: 2811 GHS: Acute Tox. 4 H312, H302 P280, P302+P352, P301+P312	
Assay: 99.5 to 100.5%		
Guaranteed analysis		
pH (aqueous solution 5%)	7 to 8.5	Molybdenum (Mo)
Aluminium (Al)	<5 ppm	Nickel (Ni)
Barium (Ba)	<5 ppm	Propan--1-ol
Cadmium (Cd)	<5 ppm	Sodium (Na)
Calcium (Ca)	<50 ppm	Strontium (Sr)
Chromium (Cr)	<5 ppm	Titanium (Ti)
Cobalt (Co)	<5 ppm	n-Butyl formate
Copper (Cu)	<5 ppm	Total nitrogen (N)
Iron (Fe)	<5 ppm	Total phosphorus (P)
Lead (Pb)	<5 ppm	Total silicon (Si)
Lithium (Li)	<5 ppm	Total sulphur (S)
Magnesium (Mg)	<20 ppm	Vanadium (V)
Manganese (Mn)	<5 ppm	Zinc (Zn)
<i>Lot analysis on the label.</i>		

Potassium oxalate monohydrate		SLR
<i>specified reagents for general laboratory work- extra pure</i>		
packaging	Mfr. No	
500 g Plastic jar	P/6320/53	
1 kg Plastic jar	P/6320/60	
C ₂ K ₂ O ₄ .H ₂ O CAS: 6487-48-5 MW: 184.23 Density at 20°C: 2.13 Melting point: 356°C ADR/RID: 6.1, III	EINECS: 209-506-8 UN: 2811 GHS: Acute Tox. 4 H312, H302 P280, P302+P352, P301+P312	
Assay: >99%		
Guaranteed analysis		
Calcium (Ca)	<200 ppm	n-Butyl formate
Copper (Cu)	<20 ppm	Total phosphorus (P)
Iron (Fe)	<20 ppm	Total silicon (Si)
Lead (Pb)	<20 ppm	Total sulphur (S)
Magnesium (Mg)	<100 ppm	Zinc (Zn)
Sodium (Na)	<1000 ppm	

Potassium perchlorate		SLR
<i>specified reagents for general laboratory work- extra pure</i>		
packaging	Mfr. No	
500 g Plastic jar	P/6400/53	
KClO ₄ CAS: 7778-74-7 MW: 138.55 Density at 20°C: 2.52 Melting point: 400°C ADR/RID: 5.1, II	EINECS: 231-912-9 UN: 1489 GHS: Ox. Sol. 1, Acute Tox. 4 H271, H302 P210, P220, P301+P312	 
Assay: >99%		
Guaranteed analysis		
Sodium (Na)	<100 ppm	n-Butyl formate
Sulfate (SO ₄)	<0.005%	

Potassium periodate		Certified AR
<i>for analysis</i>		
packaging	Mfr. No	
100 g Plastic jar	P/6480/48	
KIO ₄ CAS: 7790-21-8 MW: 230.00 Density at 20°C: 3.62 Melting point: 582°C ADR/RID: 5.1, II	GHS: Ox. Sol. 2, Acute Tox. 4, Skin Irrit. 2, Eye Irrit. 2, STOT SE 3 H272, H319, H302, H315, H312, H332, H335 P280, P210, P302+P352, P261, P301+P312, P304+P340, P305+P351+P338	 
Assay: >99.8%		
Guaranteed analysis		
Loss on drying	<0.1%	Lead (Pb)
Bromate, bromide, chlorate & chloride: < 0.01%		Magnesium (Mg)
Calcium (Ca)	<30 ppm	Sodium (Na)
Copper (Cu)	<5 ppm	Total silicon (Si)
Iodide (I)	<0.001%	Total sulphur (S)
Iron (Fe)	<20 ppm	Zinc (Zn)
<i>Lot analysis on the label</i>		

Potassium periodate		SLR
<i>specified reagents for general laboratory work- extra pure</i>		
packaging	Mfr. No	
500 g Plastic jar	P/6440/53	
KIO ₄ CAS: 7790-21-8 MW: 230.00 Density at 20°C: 3.62 Melting point: 582°C ADR/RID: 5.1, II	EINECS: 232-196-0 UN: 1479	 



Assay: >99%		
Guaranteed analysis		
Loss on drying	<0.5%	Magnesium (Mg)
Calcium (Ca)	<100 ppm	Sodium (Na)
Copper (Cu)	<20 ppm	Total silicon (Si)
Iodide (I)	<0.005%	Total sulphur (S)
Iron (Fe)	<54 ppm	Zinc (Zn)
Lead (Pb)	<20 ppm	

Potassium permanganate		Certified AR, Eur.Ph., USP, BP
<i>for analysis</i>		
packaging	Mfr. No	
250 g Glass jar	P/6560/50	
500 g Glass jar	P/6560/53	
5 kg Plastic keg	P/6560/65	
25 kg Plastic keg	P/6560/70	
KMnO ₄ CAS: 7722-64-7 MW: 158.04 Melting point: 240°C ADR/RID: 5.1, II	UN: 1490 GHS: Ox. Sol. 2, Acute Tox. 4, Aquatic Acute 1, Aquatic Chronic 1 H272, H410, H302 P210, P273, P301+P312	  
EINECS: 231-760-3		

Assay: >99%		
Guaranteed analysis		
Calcium (Ca)	<100 ppm	Sodium (Na)
Copper (Cu)	<20 ppm	Total chloride (Cl)
Iron (Fe)	<20 ppm	Total nitrogen (N)
Lead (Pb)	<20 ppm	Total sulphur (S)
Magnesium (Mg)	<20 ppm	Zinc (Zn)
Typical analysis		
Arsenic (As)	<2 ppm	Molybdenum (Mo)
Barium (Ba)	<0.001 %	Strontium (Sr)
Cadmium (Cd)	<5 ppm	Titanium (Ti)
Cobalt (Co)	<5 ppm	Vanadium (V)
Lithium (Li)	<5 ppm	

Lot analysis on the label.

Potassium permanganate		SLR
<i>specified reagents for general laboratory work- extra pure</i>		
packaging	Mfr. No	
500 g Glass jar	P/6520/53	
5 kg Glass jar	P/6520/65	
KMnO ₄ CAS: 7722-64-7 MW: 158.04 Melting point: 240°C ADR/RID: 5.1, II	UN: 1490 GHS: Ox. Sol. 2, Acute Tox. 4, Aquatic Acute 1, Aquatic Chronic 1 H272, H410, H302 P210, P273, P301+P312	  
EINECS: 231-760-3		

Assay: >99%		
Guaranteed analysis		
Calcium (Ca)	<500 ppm	Sodium (Na)
Copper (Cu)	<100 ppm	n-Butyl formate
Iron (Fe)	<100 ppm	Total sulphur (S)
Lead (Pb)	<100 ppm	Zinc (Zn)
Magnesium (Mg)	<100 ppm	



Potassium permanganate solution 0.02M (0.1N)	
<i>NIST Standard solution ready to use, For Volumetric analysis</i>	
packaging	Mfr. No
1 ℓ Glass bottle	J/6900/15
2.5 ℓ Coated Glass bottle	J/6900/PB17
KMnO ₄ CAS: 7722-64-7 MW: 158.04 Density at 20°C: 1	EINECS: 231-760-3 GHS: Aquatic Chronic 3 H412 P273

Manufactured to a tolerance of ± 0.1%.
See the full range of Standard Volumetric solution page 323.

Potassium persulfate		Certified AR
<i>for analysis</i>		
packaging	Mfr. No	
250 g Plastic jar	P/6640/50	
500 g Plastic jar	P/6640/53	
K ₂ S ₂ O ₈ CAS: 7727-21-1 MW: 270.31 Melting point: 100°C ADR/RID: 5.1, III	Resp. Sens. 1, Skin Sens. 1, STOT SE 3 H272, H334, H302, H319, H335, H315, H317 P304+P341, P342+P311, P261, P301+P312, P302+P352, P280, P305+P351+P338, P210	   
EINECS: 231-781-8 UN: 1492 GHS: Ox. Sol. 3, Acute Tox. 4, Skin Irrit. 2, Eye Irrit. 2,		
Assay: >99%		
Guaranteed analysis		
Acidity	<0.01 meq/g	Sodium (Na)
Calcium (Ca)	<50 ppm	Total chloride (Cl)
Copper (Cu)	<5 ppm	Total nitrogen (N)
Iron (Fe)	<5 ppm	Total silicon (Si)
Lead (Pb)	<5 ppm	Zinc (Zn)
Magnesium (Mg)	<20 ppm	
Typical analysis		
Aluminium (Al)	<5 ppm	Molybdenum (Mo)
Barium (Ba)	<5 ppm	Nickel (Ni)
Cadmium (Cd)	<5 ppm	Silver (Ag)
Chromium (Cr)	<5 ppm	Strontium (Sr)
Cobalt (Co)	<5 ppm	Titanium (Ti)
Lithium (Li)	<5 ppm	Vanadium (V)
Manganese (Mn)	<1 ppm	

Assay: >99%		
Guaranteed analysis		
Iron (Fe)	<20 ppm	n-Butyl formate
<0.0005%		

Potassium persulfate	
<i>for electrophoresis</i>	
packaging	Mfr. No
10 g Glass jar	P/P560/45
K ₂ S ₂ O ₈ CAS: 7727-21-1 MW: 270.31 Melting point: 100°C ADR/RID: 5.1, III	Resp. Sens. 1, Skin Sens. 1, STOT SE 3 H272, H334, H302, H319, H335, H315, H317 P304+P341, P342+P311, P261, P301+P312, P302+P352, P280, P305+P351+P338, P210
EINECS: 231-781-8 UN: 1492 GHS: Ox. Sol. 3, Acute Tox. 4, Skin Irrit. 2, Eye Irrit. 2,	
Assay: >99%	
Guaranteed analysis	
Iron (Fe)	

Assay: >99%	
Guaranteed analysis	
Iron (Fe)	



Potassium persulfate			SLR
specified reagents for general laboratory work- extra pure			
packaging		Mfr. No	
250 g	Plastic jar	P/6600/50	
500 g	Plastic jar	P/6600/53	
K ₂ S ₂ O ₈ CAS: 7727-21-1 MW: 270.31 Melting point: 100°C ADR/RID: 5.1, III EINECS: 231-781-8 UN: 1492 GHS: Ox. Sol. 3, Acute Tox. 4, Skin Irrit. 2, Eye Irrit. 2,		Resp. Sens. 1, Skin Sens. 1, STOT SE 3 H272, H334, H302, H319, H335, H315, H317 P304+P341, P342+P311, P261, P301+P312, P302+P352, P280, P305+P351+P338, P210	  
Assay: >98%			
Guaranteed analysis			
Calcium (Ca)	<200 ppm	Sodium (Na)	<500 ppm
Copper (Cu)	<20 ppm	n-Butyl formate	<0.05%
Iron (Fe)	<20 ppm	Total silicon (Si)	<100 ppm
Lead (Pb)	<20 ppm	Zinc (Zn)	<50 ppm
Magnesium (Mg)	<100 ppm		

Potassium pyrogallate solution			
pure			
packaging		Mfr. No	
2.5 ℓ	HDPE Plastic bottle	J/7000/PB17	
ADR/RID: 8, II UN: 1760 GHS: Skin Corr. 1A, Eye Dam. 1 H314		P280, P301+P330+P331, P310, P305+P351+P338, P303+P361+P353	

Potassium pyrosulfate			SLR
specified reagents for general laboratory work- extra pure			
packaging		Mfr. No	
500 g	Plastic jar	P/6720/53	
K ₂ S ₂ O ₇ CAS: 7790-62-7 MW: 254.31 Density at 20°C: 2.51 Melting point: 325°C ADR/RID: 8, III EINECS: 232-216-8		UN: 3260 GHS: Skin Corr. 1C, Eye Dam. 1 H314 P305+P351+P338, P310, P280, P260, P264, P301+P330+P331	
Assay: >97%			
Guaranteed analysis			
Iron (Fe)	<20 ppm	n-Butyl formate	<0.003%
Lead (Pb)	<10 ppm		

Potassium sodium tartrate tetrahydrate		Certified AR	
for analysis, crystals			
packaging		Mfr. No	
500 g	Plastic jar	P/6880/53	
1 kg	Plastic jar	P/6880/60	
2.5 kg	Plastic jar	P/6880/62	
C ₄ H ₄ KNaO ₆ ·4H ₂ O		Melting point: 70 to 80°C EINECS: 206-156-8	
CAS: 6381-59-5			
MW: 282.23			
Assay: 99 to 101%			
Guaranteed analysis			
pH (aqueous solution 5%)	7 to 8.5	Total nitrogen (N)	<50 ppm
Calcium (Ca)	<20 ppm	Total phosphorus (P)	<10 ppm
Copper (Cu)	<2 ppm	Total silicon (Si)	<20 ppm
Iron (Fe)	<5 ppm	Total sulphur (S)	<50 ppm
Lead (Pb)	<2 ppm	Zinc (Zn)	<2 ppm
Magnesium (Mg)	<10 ppm	Reducing substance	Pass
Total chloride (Cl)	<0.005%		

Typical analysis			
Aluminium (Al)	<2 ppm	Molybdenum (Mo)	<1 ppm
Barium (Ba)	<1 ppm	Nickel (Ni)	<1 ppm
Cadmium (Cd)	<1 ppm	Silver (Ag)	<1 ppm
Chromium (Cr)	<1 ppm	Strontium (Sr)	<5 ppm
Cobalt (Co)	<1 ppm	Titanium (Ti)	<1 ppm
Lithium (Li)	<1 ppm	Vanadium (V)	<1 ppm
Manganese (Mn)	<1 ppm		

Lot analysis on the label.

Potassium sodium tartrate tetrahydrate		SLR
extra pure, crystals		
packaging		Mfr. No
500 g	Plastic jar	P/6800/53
1 kg	Plastic jar	P/6800/60
5 kg	Plastic keg	P/6800/65
25 kg	Fibre Keg	P/6800/70
C ₄ H ₄ KNaO ₆ ·4H ₂ O CAS: 6381-59-5 MW: 282.23		Melting point: 70 to 80°C EINECS: 206-156-8

Assay: >99%			
Guaranteed analysis			
pH (aqueous solution 5%)	7 to 9	n-Butyl formate	<0.02%
Calcium (Ca)	<100 ppm	Total phosphorus (P)	<500 ppm
Copper (Cu)	<10 ppm	Total silicon (Si)	<100 ppm
Iron (Fe)	<20 ppm	Total sulphur (S)	<200 ppm
Lead (Pb)	<10 ppm	Zinc (Zn)	<10 ppm
Magnesium (Mg)	<50 ppm		

Potassium sodium tartrate tetrahydrate		SLR
extra pure, powder		
packaging		Mfr. No
500 g	Plastic jar	P/6840/53
C ₄ H ₄ KNaO ₆ ·4H ₂ O CAS: 6381-59-5 MW: 282.23		Melting point: 70 to 80°C EINECS: 206-156-8

Assay: >99%			
Guaranteed analysis			
Water insoluble matter	<0.05%	Magnesium (Mg)	<50 ppm
pH (aqueous solution 5%)	7 to 9	n-Butyl formate	<0.02%
Heavy metals	<5 ppm	Total phosphorus (P)	<50 ppm
Calcium (Ca)	<100 ppm	Total silicon (Si)	<100 ppm
Copper (Cu)	<10 ppm	Total sulphur (S)	<200 ppm
Iron (Fe)	<20 ppm	Zinc (Zn)	<10 ppm
Lead (Pb)	<10 ppm		

Potassium sulfate		Certified AR	
for analysis			
packaging		Mfr. No	
500 g	Plastic jar	P/7000/53	
1 kg	Plastic jar	P/7000/60	
5 kg	Plastic jar	P/7000/65	
K ₂ SO ₄		Boiling point: 1689°C	
CAS: 7778-80-5		Melting point: 1067°C	
MW: 174.27		EINECS: 231-915-5	
Assay: >99.5%			
Guaranteed analysis			
pH (aqueous solution 5%)	5.5 to 7.5	Sodium (Na)	<500 ppm
Calcium (Ca)	<20 ppm	Total chloride (Cl)	<0.001%
Copper (Cu)	<2 ppm	Total nitrogen (N)	<10 ppm
Iron (Fe)	<5 ppm	Total silicon (Si)	<50 ppm
Lead (Pb)	<2 ppm	Zinc (Zn)	<2 ppm
Magnesium (Mg)	<10 ppm		
Typical analysis			
Aluminium (Al)	<2 ppm	Manganese (Mn)	<1 ppm
Barium (Ba)	<1 ppm	Molybdenum (Mo)	<1 ppm
Cadmium (Cd)	<1 ppm	Nickel (Ni)	<1 ppm
Chromium (Cr)	<1 ppm	Silver (Ag)	<1 ppm
Cobalt (Co)	<1 ppm	Titanium (Ti)	<1 ppm
Lithium (Li)	<1 ppm	Vanadium (V)	<1 ppm
Lot analysis on the label.			

Lot analysis on the label.

Potassium sulfate		SLR
specified reagents for general laboratory work- extra pure		
packaging		Mfr. No
500 g	Plastic jar	P/6960/53
1 kg	Plastic jar	P/6960/60
5 kg	Plastic jar	P/6960/65
K ₂ SO ₄ CAS: 7778-80-5 MW: 174.27		Boiling point: 1689°C Melting point: 1067°C EINECS: 231-915-5
Assay: >99%		
Guaranteed analysis		
Calcium (Ca)	<200 ppm	Sodium (Na) <2000 ppm
Copper (Cu)	<20 ppm	n-Butyl formate <0.02%
Iron (Fe)	<20 ppm	Total silicon (Si) <200 ppm
Lead (Pb)	<20 ppm	Zinc (Zn) <20 ppm
Magnesium (Mg)	<100 ppm	

Potassium (+)-tartrate		SLR	
specified reagents for general laboratory work- extra pure			
packaging		Mfr. No	
500 g	Plastic jar	P/7080/53	
C ₄ H ₄ K ₂ O ₆ CAS: 921-53-9 MW: 235.28 EINECS: 213-067-8			
Assay: 99 to 102%			
Guaranteed analysis			
Calcium (Ca)	<100 ppm	Sodium (Na)	<5000 ppm
Copper (Cu)	<50 ppm	Sulfate (SO ₄)	<1000 ppm
Iron (Fe)	<50 ppm	n-Butyl formate	<0.05%
Lead (Pb)	<50 ppm	Total sulphur (S)	<200 ppm
Magnesium (Mg)	<50 ppm	Zinc (Zn)	<50 ppm

Potassium tetraborate tetrahydrate			SLR	
specified reagents for general laboratory work- extra pure				
packaging			Mfr. No	
500 g	Plastic jar		P/7140/53	
K ₂ B ₄ O ₇ ·4H ₂ O		EINECS: 215-575-5		
CAS: 12045-78-2		GHS: Repr. 2		
MW: 305.51		H361f		
Melting point: 850°C		P281		
Assay: >98.5%				
Guaranteed analysis				
Sulfate (SO ₄)		<0.02%	n-Butyl formate	<0.005%

Potassium tetrachloroplatinate(II)		SLR
specified reagents for general laboratory work- extra pure		
packaging		Mfr. No
1 g	Plastic jar	P/7170/43
K ₂ PtCl ₄ CAS: 10025-99-7 MW: 415.11 Density at 20°C: 3.499 Melting point: 250°C ADR/RID: 6.1, III EINECS: 233-050-9 UN: 3288		GHS: Acute Tox. 3, Skin Irrit. 2, Eye Dam. 1, Resp. Sens. 1, Skin Sens. 1 H301, H315, H318, H334, H317 P280, P305+P351+P338, P304+P341, P342+P311, P310, P261
Guaranteed analysis		
Identity test		: Pass


Potassium tetraoxalate dihydrate		Certified AR	
for analysis			
packaging		Mfr. No	
250 g	Plastic jar	P/7200/50	
C ₄ H ₂ KO ₈ ·2H ₂ O	UN: 2811		
CAS: 6100-20-5	GHS: Acute Tox. 4		
MW: 254.20	H312, H302		
ADR/RID: 6.1, III	P280, P302+P352, P301+P312		
EINECS: 204-874-6			
Assay: 99.5 to 100.2%			
Guaranteed analysis			
Aluminium (Al)	<5 ppm	Molybdenum (Mo)	<5 ppm
Barium (Ba)	<5 ppm	Nickel (Ni)	<5 ppm
Cadmium (Cd)	<5 ppm	Propan--1-ol	<5 ppm
Calcium (Ca)	<20 ppm	Sodium (Na)	<100 ppm
Chromium (Cr)	<5 ppm	Strontium (Sr)	<5 ppm
Cobalt (Co)	<5 ppm	Titanium (Ti)	<5 ppm
Copper (Cu)	<10 ppm	n-Butyl formate	<0.0005%
Iron (Fe)	<10 ppm	Total nitrogen (N)	<100 ppm
Lead (Pb)	<10 ppm	Total silicon (Si)	<20 ppm
Lithium (Li)	<5 ppm	Total sulphur (S)	<50 ppm
Magnesium (Mg)	<10 ppm	Vanadium (V)	<5 ppm
Manganese (Mn)	<5 ppm	Zinc (Zn)	<5 ppm
Lot analysis on the label.			

Lot analysis on the label.

Potassium thiocyanate		Certified AR
for analysis		
packaging		Mfr. No
250 g	Plastic jar	P/7280/50
500 g	Plastic jar	P/7280/53
KSCN		
CAS: 333-20-0		
MW: 97.18		
Density at 20°C: 1.886		
Boiling point: D500°C		
Melting point: 170 to 179°C		
EINECS: 206-370-1		GHS: Acute Tox. 4, Aquatic Chronic 3 H412, H302, H332, H312, EUH032 P280, P302+P352, P273, P233, P301+P312, P304+P340
Assay: >99.5%		
Guaranteed analysis		
Iodine consuming substances ...	: <0.01%	Magnesium (Mg) <10 ppm
Ammonium (NH ₄)	<0.002%	Sodium (Na) <100 ppm
Cadmium (Cd)	<5 ppm	Sulfate (SO ₄) <0.002%
Calcium (Ca)	<20 ppm	Sulfide (S) <0.001%
Copper (Cu)	<2 ppm	Total chloride (Cl) <0.005%
Iron (Fe)	<1 ppm	Total silicon (Si) <20 ppm
Lead (Pb)	<2 ppm	Zinc (Zn) <10 ppm

Typical analysis			
Aluminium (Al)	<5 ppm	Molybdenum (Mo)	<5 ppm
Barium (Ba)	<5 ppm	Nickel (Ni)	<5 ppm
Chromium (Cr)	<5 ppm	Silver (Ag)	<5 ppm
Cobalt (Co)	<5 ppm	Strontium (Sr)	<5 ppm
Lithium (Li)	<5 ppm	Titanium (Ti)	<5 ppm
Manganese (Mn)	<5 ppm	Vanadium (V)	<5 ppm

Lot analysis on the label.



Potassium thiocyanate
specified reagents for general laboratory work- extra pure

SLR

packaging	Mfr. No
250 g Plastic jar	P/7240/50
500 g Plastic jar	P/7240/53
KSCN	GHS: Acute Tox. 4, Aquatic Chronic 3
CAS: 333-20-0	H412, H302, H332, H312, EUH032
MW: 97.18	P280, P302+P352, P273, P233, P301+P312, P304+P340
Density at 20°C: 1.886	
Boiling point: D500°C	
Melting point: 170 to 179°C	
EINECS: 206-370-1	
Assay: >98%	
Guaranteed analysis	
Calcium (Ca)	<100 ppm
Copper (Cu)	<20 ppm
Iron (Fe)	<10 ppm
Lead (Pb)	<20 ppm
Magnesium (Mg)	<50 ppm
Sodium (Na)	<2000 ppm
Sulfate (SO ₄)	<0.05%
n-Butyl formate	<0.01%
Total silicon (Si)	<100 ppm
Zinc (Zn)	<50 ppm

Potassium thiocyanate solution 0.1M (0.1N)
NIST Standard solution ready to use, For Volumetric analysis

packaging	Mfr. No
1 ℓ Glass bottle	J/7050/15
2.5 ℓ HDPE Plastic bottle	J/7050/PB17
GHS: Acute Tox. 4	
H302, H312, EUH032	
P280, P301+P312, P302+P352	

Manufactured to a tolerance of ± 0.2%.
See the full range of Standard Volumetric solution page 323.

Praseodymium standard metal solution 1000 ppm
for ICP solution in HNO₃ 2%, NIST certified

packaging	Mfr. No
100 mℓ Plastic bottle	SPPLPR22Y
500 mℓ Plastic bottle	SPPLPR22X
H315, H319	

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Preservative
For microscopy, odourless

packaging	Mfr. No
500 mℓ Coated Glass bottle	P/7300/PB08
GHS: Acute Tox. 4	
H302	
P264, P270, P280, P301+P312, P330	

L-Proline
pure

packaging	Mfr. No
100 g Glass jar	P/7425/48
C ₃ H ₅ NO ₂	Density at 20°C: 1.35
CAS: 147-85-3	Melting point: 228 to 233°C
MW: 115.13	EINECS: 205-702-2
Assay: >99%	
Guaranteed analysis	
Water	<0.5%
Specific rotation	-85.5° to -84°
Residue after ignition	<0.1%
Iron (Fe)	<10 ppm
Lead (Pb)	<10 ppm

Propane-1.2-diol
specified reagents for general laboratory work- extra pure

SLR, Eur.Ph., BP, USP

packaging	Mfr. No
500 mℓ Glass bottle	P/7440/08
2.5 ℓ Glass bottle	P/7440/17
25 ℓ Metal drum	P/7440/25
200 ℓ Metal drum	P/7440/27
C ₃ H ₈ O ₂	Boiling point: 187°C
CAS: 57-55-6	Melting point: -60°C
MW: 76.09	Flash point: 99°C
Density at 20°C: 1.03	EINECS: 200-338-0
Assay (CG): >99%	
Guaranteed analysis	
Acidity	<0.0005 meq/g

Propane-1.3-diol
pure

packaging	Mfr. No
500 mℓ Glass bottle	P/7450/08
C ₃ H ₈ O ₂	Melting point: -32°C
CAS: 504-63-2	Flash point: 140°C
MW: 76.09	EINECS: 207-997-3
Boiling point: 214°C	

2-Propanol See Isopropanol, page 168

Propan-1-ol
for HPLC

Certified

packaging	Mfr. No
1 ℓ Glass bottle	P/7486/15
2.5 ℓ Glass bottle	P/7486/17
C ₃ H ₇ O	EINECS: 200-746-9
CAS: 71-23-8	UN: 1274
MW: 60.10	GHS: Flam. Liq. 2, Eye Dam. 1, STOT SE 3
Density at 20°C: 0.800	H225, H336, H318
Boiling point: 97°C	P280, P305+P351+P338, P310, P261, P301+P312, P210, P240, P304+P340
Melting point: -127°C	
Flash point: 15°C	
ADR/RID: 3, II	
Assay (CG): >99.5%	
Guaranteed analysis	
Acidity / Alkalinity	<0.0001 meq/g
Residue after evaporation	<10 ppm
Water	<0.05%
Absorbance specification	
Max. Absorbance (1cm cell v HPLC water)	
220 nm	<0.7 A.U.
230 nm	<0.25 A.U.
240 nm	<0.1 A.U.
250 nm	<0.05 A.U.
260 nm	<0.05 A.U.
270 nm	<0.01 A.U.
Lot analysis on the label.	
Filtered to 0.2 μm.	

Assay (CG): >99.5%

Guaranteed analysis

Acidity / Alkalinity

Residue after evaporation

Water

Absorbance specification

Max. Absorbance (1cm cell v HPLC water)

220 nm

230 nm

240 nm

250 nm

260 nm

270 nm

Lot analysis on the label.

Filtered to 0.2 μm.

We Support Your Chemistry For

Semi-bulk and bulk requirements

Tailored solvents and solvent blends or customized solutions

Testing services or customized packaging and labeling

Please email your contact details and a brief description of your requirements to begel.sha.quotes@thermofisher.com

Propan-1-ol
for analysis

NEW PACK.

Certified AR

packaging	Mfr. No
1 ℓ Glass bottle	P/7481/15
1 ℓ Coated Glass bottle	P/7481/PB15
2.5 ℓ Glass bottle	P/7481/17
5 ℓ Fluorinated HDPE Plastic bottle	P/7481/21
25 ℓ Metal drum	P/7481/25
200 ℓ Metal drum	P/7481/27
C ₃ H ₇ O	EINECS: 200-746-9
CAS: 71-23-8	UN: 1274
MW: 60.10	GHS: Flam. Liq. 2, Eye Dam. 1, STOT SE 3
Density at 20°C: 0.800	H225, H336, H318
Boiling point: 97°C	P280, P305+P351+P338, P310, P261, P301+P312, P210, P240, P304+P340
Melting point: -127°C	
Flash point: 15°C	
ADR/RID: 3, II	
Assay (CG): >99.5%	
Guaranteed analysis	
Acidity / Alkalinity	<0.0002 meq/g
Colour	<10 APHA
Water	<0.1%
Residue after evaporation	<10 ppm
Calcium (Ca)	<0.2 ppm
Copper (Cu)	<0.02 ppm
Iron (Fe)	<0.1 ppm
Lead (Pb)	<0.02 ppm
Magnesium (Mg)	<0.1 ppm
Potassium (K)	<0.2 ppm
Sodium (Na)	<1 ppm
Substances reducing KMnO ₄	<0.0005%
Total phosphorus (P)	<0.1 ppm
Total silicon (Si)	<0.05 ppm
Total sulphur (S)	<0.5 ppm
Zinc (Zn)	<0.1 ppm
Ethanol 90 % v/v	<0.01%
Methanol	<0.005%
Propan--2-ol	<0.005%
Substances darkened by H ₂ SO ₄	<10 APHA
Typical analysis	
Aluminium (Al)	<0.02 ppm
Barium (Ba)	<0.02 ppm
Cadmium (Cd)	<0.02 ppm
Chromium (Cr)	<0.02 ppm
Cobalt (Co)	<0.02 ppm
Lithium (Li)	<0.02 ppm
Manganese (Mn)	<0.02 ppm
Molybdenum (Mo)	<0.02 ppm
Nickel (Ni)	<0.02 ppm
Strontium (Sr)	<0.02 ppm
Titanium (Ti)	<0.02 ppm
Vanadium (V)	<0.02 ppm
Lot analysis on the label.	

Propan-1-ol
specified reagents for general laboratory work- extra pure

NEW PACK.

SLR

packaging	Mfr. No
500 mℓ Glass bottle	P/7480/08
2.5 ℓ Glass bottle	P/7480/17
5 ℓ Fluorinated HDPE Plastic bottle	P/7480/21
25 ℓ Metal drum	P/7480/25
200 ℓ Metal drum	P/7480/27
C ₃ H ₇ O	EINECS: 200-746-9
CAS: 71-23-8	UN: 1274
MW: 60.10	GHS: Flam. Liq. 2, Eye Dam. 1, STOT SE 3
Density at 20°C: 0.800	H225, H336, H318
Boiling point: 97°C	P280, P305+P351+P338, P310, P261, P301+P312, P210, P240, P304+P340
Melting point: -127°C	
Flash point: 15°C	
ADR/RID: 3, II	
Assay (CG): >99%	
Guaranteed analysis	
Acidity / Alkalinity	<0.0005 meq/g
Colour	<10 APHA
Water	<0.5%
Residue after evaporation	<50 ppm
Substances reducing KMnO ₄	<0.01%
Calcium (Ca)	<5 ppm
Copper (Cu)	<2 ppm
Iron (Fe)	<5 ppm
Lead (Pb)	<2 ppm
Magnesium (Mg)	<2 ppm
Potassium (K)	<2 ppm
Sodium (Na)	<5 ppm
Total phosphorus (P)	<2 ppm
Total silicon (Si)	<2 ppm
Total sulphur (S)	<2 ppm
Zinc (Zn)	<2 ppm
Substances darkened by H ₂ SO ₄	<20 APHA
Lot analysis on the label.	

Assay (CG): >99%

Guaranteed analysis

Acidity / Alkalinity

Colour

Water

Residue after evaporation

Substances reducing KMnO₄

Calcium (Ca)

Copper (Cu)

Iron (Fe)

Lead (Pb)

Magnesium (Mg)

Potassium (K)

Sodium (Na)

Total phosphorus (P)

Total silicon (Si)

Total sulphur (S)

Zinc (Zn)

Substances darkened by H₂SO₄

Propionaldehyde
specified reagents for general laboratory work- extra pure

SLR

packaging	Mfr. No
250 mℓ Coated Glass bottle	P/7520/PB07
C ₃ H ₆ O	EINECS: 204-623-0
CAS: 123-38-6	UN: 1275
MW: 58.08	GHS: Flam. Liq. 2, Skin Irrit. 2, Eye Irrit. 2, STOT SE 3
Density at 20°C: 0.798	H225, H315, H319, H335
Boiling point: 48°C	P261, P302+P352, P280, P305+P351+P338, P210, P240
Melting point: -81°C	
Flash point: -40°C	
ADR/RID: 3, II	
Assay (CG): >95%	
Guaranteed analysis	
Acidity	<0.1 meq/g

Propionic acid
for analysis

Certified AR

packaging	Mfr. No
500 mℓ HDPE Plastic bottle	P/7564/PB08
C ₃ H ₆ O ₂	EINECS: 201-176-3
CAS: 79-09-4	UN: 3463
MW: 74.08	GHS: Flam. Liq. 3, Skin Corr. 1B, Eye Dam. 1
Density at 20°C: 0.990	H314, H226
Boiling point: 141°C	P301+P330+P331, P280, P305+P351+P338, P310
Melting point: -22°C	
Flash point: 51°C	
ADR/RID: 8, 3, II	
Assay: >99%	
Guaranteed analysis	
Aluminium (Al)	<0.5 ppm
Barium (Ba)	<0.5 ppm
Cadmium (Cd)	<0.2 ppm
Calcium (Ca)	<5 ppm
Chromium (Cr)	<0.5 ppm
Cobalt (Co)	<0.2 ppm
Copper (Cu)	<2 ppm
Iron (Fe)	<5 ppm
Lead (Pb)	<2 ppm
Lithium (Li)	<0.2 ppm
Magnesium (Mg)	<2 ppm
Manganese (Mn)	<0.2 ppm
Molybdenum (Mo)	<0.2 ppm
Nickel (Ni)	<0.2 ppm
Potassium (K)	<1 ppm
Sodium (Na)	<5 ppm
Strontium (Sr)	<0.2 ppm
Titanium (Ti)	<0.2 ppm
Total phosphorus (P)	<2 ppm
Total silicon (Si)	<1 ppm
Total sulphur (S)	<5 ppm
Vanadium (V)	<0.2 ppm
Zinc (Zn)	<0.2 ppm
Absorbance specification	
Max. Absorbance (1cm cell v HPLC water)	
280 nm	<0.1 A.U.
Lot analysis on the label.	

Propionic acid
specified reagents for general laboratory work- extra pure

SLR

packaging	Mfr. No
500 mℓ HDPE Plastic bottle	P/7560/PB08
2.5 ℓ HDPE Plastic bottle	P/7560/PB17
25 ℓ HDPE Plastic container	P/7560/25
C ₃ H ₆ O ₂	EINECS: 201-176-3
CAS: 79-09-4	UN: 3463
MW: 74.08	GHS: Flam. Liq. 3, Skin Corr. 1B, Eye Dam. 1
Density at 20°C: 0.990	H314, H226
Boiling point: 141°C	P301+P330+P331, P280, P305+P351+P338, P310
Melting point: -22°C	
Flash point: 51°C	
ADR/RID: 8, 3, II	
Assay: >98%	
Guaranteed analysis	
Calcium (Ca)	<20 ppm
Copper (Cu)	<10 ppm
Iron (Fe)	<20 ppm
Lead (Pb)	<10 ppm
Magnesium (Mg)	<20 ppm
Potassium (K)	<10 ppm
Sodium (Na)	<20 ppm
Total phosphorus (P)	<10 ppm
Total silicon (Si)	<10 ppm
Total sulphur (S)	<20 ppm
Zinc (Zn)	<10 ppm

iso-Propyl acetate			
pure			
packaging		Mfr. No	
500 mℓ	Glass bottle	P/7640/08	
2.5 ℓ	Glass bottle	P/7640/17	
2.5 ℓ	Coated Glass bottle	P/7640/PB17	
25 ℓ	Metal drum	P/7640/25	
C ₅ H ₁₀ O ₂	EINECS: 203-561-1	 	
CAS: 108-21-4	UN: 1220		
MW: 102.13	GHS: Flam. Liq. 2, Eye Irrit. 2, STOT SE 3		
Density at 20°C: 0.872	H225, H319, H336, EUH066		
Boiling point: 88.8°C	P261, P280,		
Melting point: -73°C	P305+P351+P338, P210, P240		
Flash point: 4°C			
ADR/RID: 3, II			
Assay (CG): >98%			
Guaranteed analysis			
Residue after evaporation	<100 ppm	Acidity	<0.002 meq/g

Pyridine			
For IR Spectroscopy			
packaging		Mfr. No	
100 mℓ	Coated Glass bottle	P/7965/PB05	
500 mℓ	Coated Glass bottle	P/7965/PB08	
C ₅ H ₅ N	ADR/RID: 3, II	 	
CAS: 110-86-1	EINECS: 203-809-9		
MW: 79.10	UN: 1282		
Density at 20°C: 0.978	GHS: Flam. Liq. 2, Acute Tox. 4		
Boiling point: 115 to 116°C	H225, H312, H332, H302		
Melting point: -42°C	P280, P302+P352,		
Flash point: 17°C	P304+P340, P210, P240		
Guaranteed analysis			
Complies with standard IR spectrum: Pass			
Lot analysis on the label.			

n-Propyl acetate		SLR
specified reagents for general laboratory work- extra pure		
packaging		Mfr. No
500 mℓ	Glass bottle	P/7600/08
2.5 ℓ	Glass bottle	P/7600/17
C ₅ H ₁₀ O ₂	EINECS: 203-686-1	 
CAS: 109-60-4	UN: 1276	
MW: 102.13	GHS: Flam. Liq. 2, Eye Irrit. 2, STOT SE 3	
Density at 20°C: 0.880	H225, H336, H319, EUH066	
Boiling point: 102°C	P261, P280,	
Melting point: -95°C	P305+P351+P338, P210, P240	
Flash point: 10°C		
ADR/RID: 3, II		
Assay (CG): >98%		
Guaranteed analysis		
Residue after evaporation	<100 ppm	

Propyl alcohol See Propan-1-ol, page 264

Propylene glycol See Propane-1, page 264

Pumice stone	
Technical, granular about 7-9mm	
packaging	
500 g	Plastic jar
CAS: 1332-09-8	
Guaranteed analysis	
Identity test	: Pass

Pumice stone	
Technical, powder 240 mesh	
packaging	
3 kg	Plastic jar
CAS: 1332-09-8	
Guaranteed analysis	
Identity test	: Pass

Pyridine		Certified AR
for analysis		
packaging		Mfr. No
500 mℓ	Glass bottle	P/7960/08
500 mℓ	Coated Glass bottle	P/7960/PB08
1 ℓ	Glass bottle	P/7960/15
1 ℓ	Coated Glass bottle	P/7960/PB15
2.5 ℓ	Coated Glass bottle	P/7960/PB17
25 ℓ	Metal drum	P/7960/25
C ₅ H ₅ N	ADR/RID: 3, II	 
CAS: 110-86-1	EINECS: 203-809-9	
MW: 79.10	UN: 1282	
Density at 20°C: 0.978	GHS: Flam. Liq. 2, Acute Tox. 4	
Boiling point: 115 to 116°C	H225, H312, H332, H302	
Melting point: -42°C	P280, P302+P352,	
Flash point: 17°C	P304+P340, P210, P240	
Assay: >99.5%		
Assay (CG): >99.5%		
Guaranteed analysis		
Water	<0.05%	Potassium (K)
Residue after evaporation	<20 ppm	Sodium (Na)
Ammonium (NH ₄)	<0.002 %	Substances
Calcium (Ca)	<0.5 ppm	reducing KMnO ₄
Chlorine	<5 ppm	Total phosphorus (P)
Free Chlorine (Cl ₂)	<5 ppm	Total silicon (Si)
Copper (Cu)	<0.05 ppm	Total sulphur (S)
Iron (Fe)	<0.2 ppm	Zinc (Zn)
Lead (Pb)	<0.05 ppm	Ammonia (NH ₃)
Magnesium (Mg)	<0.1 ppm	
Typical analysis		
Aluminium (Al)	0.02 ppm	Manganese (Mn)
Barium (Ba)	0.02 ppm	Molybdenum (Mo)
Cadmium (Cd)	0.02 ppm	Nickel (Ni)
Chromium (Cr)	0.02 ppm	Strontium (Sr)
Cobalt (Co)	0.02 ppm	Titanium (Ti)
Lithium (Li)	0.02 ppm	Vanadium (V)

Pyridine NEW PACK.		SLR
specified reagents for general laboratory work- extra pure		
packaging		Mfr. No
500 mℓ	Glass bottle	P/7920/08
2.5 ℓ	Coated Glass bottle	P/7920/PB17
5 ℓ	Fluorinated HDPE Plastic bottle	P/7920/21
25 ℓ	Metal drum	P/7920/25
C ₅ H ₅ N	ADR/RID: 3, II	 
CAS: 110-86-1	EINECS: 203-809-9	
MW: 79.10	UN: 1282	
Density at 20°C: 0.978	GHS: Flam. Liq. 2, Acute Tox. 4	
Boiling point: 115 to 116°C	H225, H312, H332, H302	
Melting point: -42°C	P280, P302+P352,	
Flash point: 17°C	P304+P340, P210, P240	
Assay (CG): >99%		

Guaranteed analysis			
Colour	>0 APHA	Potassium (K)	<2 ppm
Water	<0.1%	Sodium (Na)	<5 ppm
Residue after evaporation	<200 ppm	Substances	
Calcium (Ca)	<5 ppm	reducing KMnO ₄	Pass
Copper (Cu)	<2 ppm	Total phosphorus (P)	<2 ppm
Iron (Fe)	<5 ppm	Total silicon (Si)	<2 ppm
Lead (Pb)	<2 ppm	Total sulphur (S)	<2 ppm
Magnesium (Mg)	<2 ppm	Zinc (Zn)	<2 ppm

Pyridine anhydrous	
for DNA synthesis	
packaging	
1 ℓ	Coated Glass bottle
C ₅ H ₅ N	ADR/RID: 3, II
CAS: 110-86-1	EINECS: 203-809-9
MW: 79.10	UN: 1282
Density at 20°C: 0.978	GHS: Flam. Liq. 2, Acute Tox. 4
Boiling point: 115 to 116°C	H225, H312, H332, H302
Melting point: -42°C	P280, P302+P352,
Flash point: 17°C	P304+P340, P210, P240
Guaranteed analysis	
Water	<0.01%
Containing activated 3A molecular sieve.	

Pyridine, anhydrous	
specified reagents for general laboratory work- extra pure	

packaging		Mfr. No
500 mℓ	Coated Glass bottle	P/7940/PB08
2.5 ℓ	Coated Glass bottle	P/7940/PB17
C ₅ H ₅ N	ADR/RID: 3, II	 
CAS: 110-86-1	EINECS: 203-809-9	
MW: 79.10	UN: 1282	
Density at 20°C: 0.978	GHS: Flam. Liq. 2, Acute Tox. 4	
Boiling point: 115 to 116°C	H225, H312, H332, H302	
Melting point: -42°C	P280, P302+P352,	
Flash point: 17°C	P304+P340, P210, P240	
Guaranteed analysis		
Water	<0.02%	
<i>For the analysis by Karl Fischer titration in non-aqueous.</i>		

Pyridine solution 65 % v/v in water	
pure	
packaging	
500 mℓ	Glass bottle
ADR/RID: 3, III	P210, P280, P301+P312,
UN: 1993	P304+P340,
GHS: Flam. Liq. 2, Acute Tox. 4	P305+P351+P338,
H226, H302, H312, H332	P303+P361+P353

3-Pyridinecarboxylic acid See Nicotinic acid, page 214

Pyridinium chloride	
specified reagents for general laboratory work- extra pure	
packaging	
100 g	Glass jar
C ₅ H ₅ N.HCl	GHS: Acute Tox. 4, Skin Irrit. 2, Eye Irrit. 2, STOT SE 3
CAS: 628-13-7	H319, H302, H312, H332, H315, H335
MW: 115.56	P280, P302+P352, P261,
Boiling point: 222 to 224°C	P301+P312, P304+P340,
Melting point: 140 to 146°C	P305+P351+P338
ADR/RID: 6.1, III	
EINECS: 211-027-4	
UN: 2811	
Assay: >97%	

1-(2-Pyridylazo)-2-naphthol	
pure	
packaging	
1 g	Vial
5 g	Glass jar
C ₁₅ H ₁₁ N ₃ O	GHS: Skin Irrit. 2, Eye Irrit. 2, STOT SE 3
CAS: 85-85-8	H315, H335, H319
MW: 249.27	P261, P302+P352, P280,
Melting point: 138 to 141°C	P305+P351+P338
EINECS: 201-637-9	
Guaranteed analysis	
Complies with standard IR spectrum: Pass	

4-(2-Pyridylazo)-resorcinol mono sodium salt	
pure	
packaging	
1 g	Vial
C ₁₁ H ₈ N ₂ NaO ₂ ·H ₂ O	GHS: Skin Irrit. 2, Eye Irrit. 2, STOT SE 3
CAS: 1141-59-9	H335, H319, H315
MW: 215.21	P261, P302+P352, P280,
Melting point: 188°C	P305+P351+P338
EINECS: 214-528-6	
Guaranteed analysis	
Complies with standard IR spectrum: Pass	

Pyrocatechol See Catechol, page 80

Pyrogallol		Certified AR
for analysis		
packaging		Mfr. No
100 g	Glass jar	P/8122/48
C ₆ H ₆ O ₃	EINECS: 201-762-9	 
CAS: 87-66-1	UN: 2811	
MW: 126.11	GHS: Muta. 2, Acute Tox. 4, Aquatic Chronic 3	
Density at 20°C: 1.453	H412, H302, H341, H332, H312	
Boiling point: 309°C	P280, P302+P352, P273,	
Melting point: 131 to 135°C	P301+P312, P304+P340	
ADR/RID: 6.1, II		
Assay (CG): >99%		
Guaranteed analysis		
pH (aqueous solution 5%)	4 to 5	Sodium (Na)
Calcium (Ca)	<50 ppm	n-Butyl formate
Copper (Cu)	<10 ppm	Total phosphorus (P)
Iron (Fe)	<10 ppm	Total silicon (Si)
Lead (Pb)	<10 ppm	Total sulphur (S)
Magnesium (Mg)	<20 ppm	Zinc (Zn)
Potassium (K)	<20 ppm	
Lot analysis on the label.		

Pyrogallol		SLR
specified reagents for general laboratory work- extra pure		
packaging		Mfr. No
100 g	Glass jar	P/8120/48
C ₆ H ₆ O ₃	EINECS: 201-762-9	 
CAS: 87-66-1	UN: 2811	
MW: 126.11	GHS: Muta. 2, Acute Tox. 4, Aquatic Chronic 3	
Density at 20°C: 1.453	H412, H302, H341, H332, H312	
Boiling point: 309°C	P280, P302+P352, P273,	
Melting point: 131 to 135°C	P301+P312, P304+P340	
ADR/RID: 6.1, II		
Assay (CG): >98%		
Guaranteed analysis		
Calcium (Ca)	<200 ppm	Sodium (Na)
Copper (Cu)	<50 ppm	n-Butyl formate
Iron (Fe)	<50 ppm	Total phosphorus (P)
Lead (Pb)	<50 ppm	Total silicon (Si)
Magnesium (Mg)	<100 ppm	Total sulphur (S)
Potassium (K)	<100 ppm	Zinc (Zn)

Pyrogallol <i>pure</i>		
packaging	Mfr. No	
100 g Glass jar	P/8080/48	
C ₆ H ₆ O ₃ CAS: 87-66-1 MW: 126.11 Density at 20°C: 1.453 Boiling point: 309°C Melting point: 131 to 135°C ADR/RID: 6.1, II	EINECS: 201-762-9 UN: 2811 GHS: Muta. 2, Acute Tox. 4, Aquatic Chronic 3 H412, H302, H341, H332, H312 P280, P302+P352, P273, P301+P312, P304+P340	 
Assay (CG): >98%		

Pyrrolidine <i>specified reagents for general laboratory work- extra pure</i>		SLR
packaging	Mfr. No	
100 mL Coated Glass bottle	P/8140/PB05	
C ₄ H ₉ N CAS: 123-75-1 MW: 71.12 Density at 20°C: 0.866 Boiling point: 86 to 88°C Melting point: -63°C Flash point: 3°C ADR/RID: 3, 8, II	EINECS: 204-648-7 UN: 1922 GHS: Flam. Liq. 2, Acute Tox. 4, Skin Corr. 1A, Eye Dam. 1 H225, H314, H332, H302 P301+P330+P331, P280, P305+P351+P338, P310, P304+P340, P210	  
Assay (CG): >97%		
Guaranteed analysis		
Water	<3%	

Quinhydrone <i>specified reagents for general laboratory work- extra pure</i>		SLR
packaging	Mfr. No	
100 g Plastic jar	Q/0120/48	
C ₁₂ H ₁₀ O ₄ CAS: 106-34-3 MW: 218.21 Density at 20°C: 14 Melting point: 168 to 172°C ADR/RID: 6.1, III EINECS: 203-387-6	UN: 2811 GHS: Acute Tox. 3, Skin Irrit. 2, Eye Irrit. 2, STOT SE 3, Aquatic Acute 1 H301, H319, H335, H315, H400 P261, P273, P302+P352, P280, P305+P351+P338	 
Guaranteed analysis		
Complies with standard IR spectrum: Pass	Potassium (K)	<50 ppm
Calcium (Ca)	Sodium (Na)	<500 ppm
Copper (Cu)	Total phosphorus (P)	<50 ppm
Iron (Fe)	Total silicon (Si)	<20 ppm
Lead (Pb)	Total sulphur (S)	<100 ppm
Magnesium (Mg)	Zinc (Zn)	<20 ppm

Quinine sulfate dihydrate <i>pure</i>		
packaging	Mfr. No	
25 g Glass jar	Q/0136/46	
C ₂₀ H ₄₈ N ₄ O ₄ ·H ₂ SO ₄ ·2H ₂ O CAS: 6119-70-6 MW: 782.96 Melting point: 233 to 235°C ADR/RID: 6.1, III EINECS: 212-359-2	UN: 1544 GHS: Acute Tox. 4, Resp. Sens. 1, Skin Sens. 1 H334, H302, H317 P261, P304+P341, P342+P311, P280, P301+P312	
Guaranteed analysis		
Complies with standard IR spectrum: Pass		

Quinol <i>for analysis</i>		Certified AR
packaging	Mfr. No	
500 g Plastic jar	Q/0152/53	
C ₆ H ₆ O ₂ CAS: 123-31-9 MW: 110.11 Density at 20°C: 1.32 Boiling point: 285 to 286°C Melting point: 170 to 174°C Flash point: 165°C ADR/RID: 9, III EINECS: 204-617-8	UN: 3077 GHS: Carc. 2, Muta. 2, Acute Tox. 4, Eye Dam. 1, Skin Sens. 1, Aquatic Acute 1 H318, H341, H317, H351, H302, H400 P280, P312, P302+P350, P305+P351+P338, P301+P312, P273	   

Assay: >99%			
Guaranteed analysis			
Calcium (Ca)	<50 ppm	Sodium (Na)	<100 ppm
Copper (Cu)	<10 ppm	Total phosphorus (P)	<100 ppm
Iron (Fe)	<10 ppm	Total silicon (Si)	<20 ppm
Lead (Pb)	<10 ppm	Total sulphur (S)	<200 ppm
Magnesium (Mg)	<20 ppm	Zinc (Zn)	<10 ppm
Potassium (K)	<50 ppm	Pyrocatechol	<0.02%

Quinol <i>specified reagents for general laboratory work- extra pure</i>		SLR
packaging	Mfr. No	
250 g Plastic jar 500 g Plastic jar	Q/0150/50 Q/0150/53	
C ₆ H ₆ O ₂ CAS: 123-31-9 MW: 110.11 Density at 20°C: 1.32 Boiling point: 285 to 286°C Melting point: 170 to 174°C Flash point: 165°C ADR/RID: 9, III EINECS: 204-617-8	UN: 3077 GHS: Carc. 2, Muta. 2, Acute Tox. 4, Eye Dam. 1, Skin Sens. 1, Aquatic Acute 1 H318, H341, H317, H351, H302, H400 P280, P312, P302+P350, P305+P351+P338, P301+P312, P273	   

Assay: >98%			
Guaranteed analysis			
Calcium (Ca)	<20 ppm	Sodium (Na)	<500 ppm
Copper (Cu)	<50 ppm	Total phosphorus (P)	<200 ppm
Iron (Fe)	<50 ppm	Total silicon (Si)	<100 ppm
Lead (Pb)	<50 ppm	Total sulphur (S)	<500 ppm
Magnesium (Mg)	<100 ppm	Zinc (Zn)	<50 ppm
Potassium (K)	<200 ppm		

Quinone *See p-Benzquinone (dried) , page 56*

Raffinose hydrate		SLR
specified reagents for general laboratory work- extra pure		
packaging	Mfr. No	
25 g Plastic jar	R/0020/46	
C ₁₈ H ₃₂ O ₁₆ ·5H ₂ O	Melting point: 80°C	
CAS: 512-69-6	EINECS: 208-146-9	
MW: 594.52		
Guaranteed analysis		
Water	15 to 15.5%	Lead (Pb) <10 ppm
Specific rotation	103° to 106°	Residue after ignition <0.1%
Reducing sugars	: <0.3%	

Rennet <i>pure</i>		
packaging	Mfr. No	
25 g Glass jar	R/0030/46	
CAS: 9042-08-4 EINECS: 232-926-8		
Guaranteed analysis		
Identity test	: Pass	

Resazurin tablets <i>pure</i>		
packaging	Mfr. No	
10 PC Glass jar 1000 PC Glass jar	R/0040/99 R/0040/79	
CAS: 550-82-3 EINECS: 208-987-1		
Guaranteed analysis		
Complies with standard IR spectrum: Pass		

Resorcinol <i>for analysis</i>		Certified AR	
packaging		Mfr. No	
100 g	Glass jar	R/0200/48	
250 g	Glass jar	R/0200/50	
C ₆ H ₆ O ₂	EINECS: 203-585-2	 	
CAS: 108-46-3	UN: 2876		
MW: 110.11	GHS: Acute Tox. 4, Skin Irrit.		
Density at 20°C: 1.272	2, Eye Irrit. 2, Aquatic Acute 1		
Boiling point: 281°C	H319, H315, H410, H302		
Melting point: 109 to 111°C	P273, P302+P352, P280,		
Flash point: 127°C	P305+P351+P338		
ADR/RID: 6.1, III			
Assay (CG): >99%			
Guaranteed analysis			
Acidity	<0.0005 meq/g	Potassium (K)	<20 ppm
Calcium (Ca)	<50 ppm	Sodium (Na)	<100 ppm
Copper (Cu)	<10 ppm	Total phosphorus (P)	<100 ppm
Iron (Fe)	<10 ppm	Total silicon (Si)	<20 ppm
Lead (Pb)	<10 ppm	Total sulphur (S)	<200 ppm
Magnesium (Mg)	<20 ppm	Zinc (Zn)	<10 ppm
Lot analysis on the label.			

Resorcinol specified reagents for general laboratory work- extra pure		SLR	
packaging		Mfr. No	
100 g	Glass jar	R/0150/48	
250 g	Glass jar	R/0150/50	
C₆H₆O₂ CAS: 108-46-3 MW: 110.11 Density at 20°C: 1.272 Boiling point: 281°C Melting point: 109 to 111°C Flash point: 127°C ADR/RID: 6.1, III		EINECS: 203-585-2 UN: 2876 GHS: Acute Tox. 4, Skin Irrit. 2, Eye Irrit. 2, Aquatic Acute 1 H319, H315, H410, H302 P273, P302+P352, P280, P305+P351+P338  	
Assay (CG): >98%			
Guaranteed analysis			
Calcium (Ca)	<500 ppm	Sodium (Na)	<500 ppm
Copper (Cu)	<50 ppm	Total phosphorus (P)	<200 ppm
Iron (Fe)	<50 ppm	Total silicon (Si)	<100 ppm
Lead (Pb)	<50 ppm	Total sulphur (S)	<500 ppm
Magnesium (Mg)	<100 ppm	Zinc (Zn)	<50 ppm
Potassium (K)	<100 ppm		

Rhenium standard metal solution 1000 ppm <i>for ICP solution in H₂O, NIST certified</i>		
packaging	Mfr. No	
100 mL Plastic bottle 500 mL Plastic bottle	SPPLRE92Y SPPLRE92X	

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Rhodamine B C.I. 45170 <i>pure</i>		
packaging	Mfr. No	
100 g Glass jar	R/0250/48	
C ₂₈ H ₃₁ ClN ₂ O ₃ CAS: 81-88-9 MW: 479.02 Density at 20°C: >1.0 Melting point: 210 to 211°C ADR/RID: 6.1, III EINECS: 201-383-9	UN: 2811 GHS: Muta. 2, Acute Tox. 4, Eye Dam. 1 H318, H341, H302 P280, P305+P351+P338, P301+P312, P260	  

Guaranteed analysis	
Complies with standard IR spectrum: Pass	
Conforming to the IR spectrum of reference.	

Rhodium standard metal solution 1000 ppm <i>for ICP solution in HCL 10%, NIST certified</i>		
packaging	Mfr. No	
100 mL Plastic bottle 500 mL Plastic bottle	SPPLRH32Y SPPLRH32X	
H315, H319, H335		

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Rhodium standard metal solution 10 ppm <i>for ICP-MS solution in HCl 2%, NIST certified</i>		
packaging	Mfr. No	
100 mL Plastic bottle	SPCLRH11AY	

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Rhodizonic acid sodium salt <i>pure</i>		
packaging	Mfr. No	
5 g Glass jar	R/0300/44	
C ₆ Na ₂ O ₆ CAS: 523-21-7 MW: 214.04	Melting point: >300°C EINECS: 208-340-3	
Guaranteed analysis		
Identity test	: Pass	

Riboflavin
for electrophoresis

packaging	Mfr. No
25 g Plastic jar	R/P300/46
C ₁₇ H ₂₀ N ₄ O ₆ CAS: 83-88-5 MW: 376.36	Melting point: 290°C EINECS: 201-507-1
Assay: 98 to 102%	

Ricinus oil

See Castor oil, page 80

Ringer's solution
pure, tablets

packaging	Mfr. No
1000 PC Plastic jar	R/0350/79
EINECS: 233-140-8	
Guaranteed analysis	
Complies with standard IR spectrum: Pass	
Quarter strength for milk testing.One tablet dissolved in 500ml of glass distilled water yields quarter strength solution.	

Rose Bengal C.I. 45440
pure

packaging	Mfr. No
25 g Glass jar	R/0400/46
C ₂₀ H ₂ Cl ₄ I ₄ Na ₂ O ₃ CAS: 632-69-9 MW: 101.65 Flash point: >92°C	EINECS: 211-183-3 GHS: Eye Irrit. 2 H319 P280, P305+P351+P338
Guaranteed analysis	
Complies with standard IR spectrum: Pass	

Rose water
pure

packaging	Mfr. No
100 ml Glass bottle	R/0370/05
CAS: 7732-18-5 Density at 20°C: 1 Boiling point: 100°C	Melting point: 0°C EINECS: 231-791-2

Rubidium standard metal solution 1000 ppm
for ICP solution in HNO₃ 2%, NIST certified

packaging	Mfr. No
100 ml Plastic bottle	SPPLRB22Y
500 ml Plastic bottle	SPPLRB22X
H315, H320	
Manufactured to a tolerance of ± 0.5%. All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities. See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314. Check with your local distributor for availability of Spex product	

Rubidium standard metal solution 10 000 ppm
for ICP solution in HNO₃ 5%, NIST certified

packaging	Mfr. No
100 ml Plastic bottle	PLRB23Y
500 ml Plastic bottle	PLRB23X
H314	
Manufactured to a tolerance of ± 0.5%. All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities. See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314. Check with your local distributor for availability of Spex product	

Rubidium chloride
specified reagents for general laboratory work- extra pure

packaging	Mfr. No
5 g Plastic jar	R/0500/44
RbCl CAS: 7791-11-9 MW: 120.92	Boiling point: 1390°C Melting point: 718°C EINECS: 232-240-9
Guaranteed analysis	
Identity test : Pass	

Ruthenium standard metal solution 1000 ppm
for ICP solution in HCL 10%, NIST certified

packaging	Mfr. No
100 ml Plastic bottle	SPPLRU32Y
500 ml Plastic bottle	SPPLRU32X
H315, H319, H335	
Manufactured to a tolerance of ± 0.5%. All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities. See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314. Check with your local distributor for availability of Spex product	

Ruthenium trichloride
pure, powder

packaging	Mfr. No
1 g Vial	R/0600/43
RuCl ₃ CAS: 10049-08-8 MW: 207.43 Density at 20°C: 3.11 Melting point: 500°C ADR/RID: 6.1, III EINECS: 233-167-5 UN: 3260	GHS: Acute Tox. 4, Skin Corr. 1B, Eye Dam. 1, Aqua- tic Chronic 3 H412, H314, H302 P301+P330+P331, P280, P305+P351+P338, P310, P273, P301+P312
Guaranteed analysis	
Identity test : Pass	

Saccharin
specified reagents for general laboratory work- extra pure

packaging	Mfr. No		
100 g Plastic jar	S/0040/48		
C ₇ H ₅ NO ₂ S CAS: 81-07-2 MW: 183.19 Density at 20°C: 0.828 Melting point: 226 to 230°C	EINECS: 201-321-0 GHS: Skin Corr. 1B, Eye Dam. 1 H314 P301+P330+P331, P280, P305+P351+P338, P310		
Assay: >98%			
Guaranteed analysis			
Water	<1%	Arsenic (As)	<2 ppm
Residue after ignition	<0.2%	Lead (Pb)	<10 ppm

Saccharin sodium dihydrate
specified reagents for general laboratory work- extra pure

packaging	Mfr. No		
100 g Plastic jar	S/0080/48		
C ₇ H ₄ NNaO ₃ S.2H ₂ O CAS: 6155-57-3 MW: 241.21 EINECS: 204-886-1	H314 P301+P330+P331, P280, P305+P351+P338, P310		
Assay: >98%			
Guaranteed analysis			
Loss on drying	12 to 16%	Lead (Pb)	<10 ppm
Arsenic (As)	<2 ppm		

Safranine C.I 50240
pure

packaging	Mfr. No
25 g Glass jar	S/0120/46
C ₂₀ H ₁₉ CIN ₄ CAS: 477-73-6 MW: 350.84 EINECS: 207-518-8	GHS: Skin Irrit. 2, Eye Irrit. 2 H319, H315 P302+P352, P280, P305+P351+P338
Guaranteed analysis	
Complies with standard IR spectrum: Pass	

Salicylaldehyde
specified reagents for general laboratory work- extra pure

packaging	Mfr. No
100 ml Coated Glass bottle	S/0200/PB05
C ₇ H ₆ O ₂ CAS: 90-02-8 MW: 122.12 Density at 20°C: 1.16 Boiling point: 197°C Melting point: -7°C Flash point: 76°C ADR/RID: 6.1, III	EINECS: 201-961-0 UN: 2810 GHS: Acute Tox. 3, Skin Irrit. 2, Eye Irrit. 2 H315, H319, H311, H302 P280, P302+P352, P264, P270, P305+P351+P338, P301+P312

Salicylic acid
for analysis

packaging	Mfr. No		
500 g Plastic jar	S/0282/53		
C ₇ H ₆ O ₃ CAS: 69-72-7 MW: 138.12 Boiling point: 211°C Melting point: 158 to 161°C Flash point: 157°C EINECS: 200-712-3	GHS: Acute Tox. 4, Skin Irrit. 2, Eye Dam. 1, STOT SE 3 H302, H318, H315, H335 P280, P305+P351+P338, P304+P340, P261, P301+P312, P302+P352		
Assay: >99.8%			
Guaranteed analysis			
Calcium (Ca)	<20 ppm	Sodium (Na)	<50 ppm
Copper (Cu)	<1 ppm	n-Butyl formate	<0.001%
Iron (Fe)	<1 ppm	Total phosphorus (P)	<20 ppm
Lead (Pb)	<1 ppm	Total silicon (Si)	<10 ppm
Magnesium (Mg)	<10 ppm	Total sulphur (S)	<100 ppm
Potassium (K)	<20 ppm	Zinc (Zn)	<5 ppm

Lot analysis on the label

UHPLC

Filtered to 0.1 µm, Fisher Chemical UHPLC gradient grade solvents (Acetonitrile, Methanol and Water) are perfect for those working at high pressure.

See page 309 for more information

Salicylic acid
specified reagents for general laboratory work- extra pure

packaging	Mfr. No		
250 g Plastic jar	S/0280/50		
500 g Plastic jar	S/0280/53		
1 kg Plastic jar	S/0280/60		
C ₇ H ₆ O ₃ CAS: 69-72-7 MW: 138.12 Boiling point: 211°C Melting point: 158 to 161°C Flash point: 157°C EINECS: 200-712-3	GHS: Acute Tox. 4, Skin Irrit. 2, Eye Dam. 1, STOT SE 3 H302, H318, H315, H335 P280, P305+P351+P338, P304+P340, P261, P301+P312, P302+P352		
Assay: >99.5%			
Guaranteed analysis			
Calcium (Ca)	<100 ppm	Sodium (Na)	<200 ppm
Copper (Cu)	<5 ppm	n-Butyl formate	<0.005%
Iron (Fe)	<5 ppm	Total phosphorus (P)	<100 ppm
Lead (Pb)	<5 ppm	Total silicon (Si)	<50 ppm
Magnesium (Mg)	<50 ppm	Total sulphur (S)	<500 ppm
Potassium (K)	<100 ppm	Zinc (Zn)	<20 ppm

Samarium standard metal solution 1000 ppm
for ICP solution in HNO₃ 2%, NIST certified

packaging	Mfr. No
100 ml Plastic bottle	SPPLSM22Y
500 ml Plastic bottle	SPPLSM22X
H315, H319	

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Sand Ottawa
extra pure, 20-30 mesh

packaging	Mfr. No
1 kg Plastic jar	S/0365/60
3 kg Plastic jar	S/0365/63
SiO ₂ CAS: 14808-60-7 MW: 60.08 Density at 20°C: 2.65	Boiling point: 2230°C Melting point: 1610°C EINECS: 238-878-4
Guaranteed analysis	
Identity test : Pass	

Sand low iron
extra pure, 40-100 mesh

packaging	Mfr. No
1 kg Plastic jar	S/0360/60
3 kg Plastic jar	S/0360/63
SiO ₂ CAS: 14808-60-7 MW: 60.08 Density at 20°C: 2.65	Boiling point: 2230°C Melting point: 1610°C EINECS: 238-878-4
Guaranteed analysis	
Identity test : Pass	

Semicarbazide hydrochloride
specified reagents for general laboratory work- extra pure

SLR

packaging		Mfr. No
100 g	Plastic jar	S/0560/48
25 kg	Plastic keg	S/0560/70
CH ₃ N ₃ O.HCl	GHS: Carc. 2, Acute Tox. 3, Skin Corr. 1B, Eye Dam. 1	
CAS: 563-41-7	H351, H314, H301	
MW: 111.53	P301+P310,	
Melting point: 170 to 178°C	P301+P330+P331, P280,	
ADR/RID: 6.1, III	P305+P351+P338, P310, P281	
EINECS: 209-247-0		
UN: 2811		

Assay: >98%

Guaranteed analysis			
Calcium (Ca)	<200 ppm	Sodium (Na)	<200 ppm
Copper (Cu)	<50 ppm	Total phosphorus (P)	<100 ppm
Iron (Fe)	<50 ppm	Total silicon (Si)	<50 ppm
Lead (Pb)	<50 ppm	Total sulphur (S)	<100 ppm
Magnesium (Mg)	<100 ppm	Zinc (Zn)	<50 ppm
Potassium (K)	<100 ppm	Hydrazine	<0.5%

Shellac orange
pure, flake

packaging		Mfr. No
250 g	Plastic jar	S/0600/50
500 g	Plastic jar	S/0600/53
CAS: 9000-59-3		
EINECS: 232-549-9		
Guaranteed analysis		
Complies with standard IR spectrum: Pass		

Silica gel incorporating 12-13% calcium sulfate as a binder
For Thin layer chromatography

packaging		Mfr. No
500 g	Plastic jar	S/0790/53
SiO ₂		
MW: 60.08		
Guaranteed analysis		
Identity test : Pass		

Silica gel 40 Angstrom pore size, 0.2-0.5mm
For chromatography column

packaging		Mfr. No
500 g	Plastic jar	S/0780/53
1 kg	Plastic jar	S/0780/60
SiO ₂		
CAS: 112926-00-8		
MW: 60.08		

Silica gel 60 Angstrom pore size, 6-35 micron
For chromatography column, Davisil® Silicagel

packaging		Mfr. No
500 g	Plastic jar	S/0691/53
1 kg	Plastic jar	S/0691/60
25 kg	Fibre Keg	S/0691/70
SiO ₂	Boiling point: 2230°C	
CAS: 7631-86-9	Melting point: 1710°C	
MW: 60.08	EINECS: 231-545-4	

Silica gel 60 Angstrom pore size, 20-45 micron
For chromatography column, Davisil® Silicagel

packaging		Mfr. No
500 g	Plastic jar	S/0692/53
1 kg	Plastic jar	S/0692/60
25 kg	Fibre Keg	S/0692/70
SiO ₂	Boiling point: 2230°C	
CAS: 7631-86-9	Melting point: 1710°C	
MW: 60.08	EINECS: 231-545-4	

Silica gel 60 Angstrom pore size, 40-63 micron
For chromatography column, Davisil® Silicagel

packaging		Mfr. No
500 g	Plastic jar	S/0694/53
1 kg	Plastic jar	S/0694/60
25 kg	Fibre Keg	S/0694/70
SiO ₂	Boiling point: 2230°C	
CAS: 7631-86-9	Melting point: 1710°C	
MW: 60.08	EINECS: 231-545-4	

Silica gel 60 Angstrom pore size, 60-200 micron
For chromatography column, Davisil® Silicagel

packaging		Mfr. No
500 g	Plastic jar	S/0694/53
1 kg	Plastic jar	S/0694/60
25 kg	Fibre Keg	S/0694/70
SiO ₂	Boiling point: 2230°C	
CAS: 7631-86-9	Melting point: 1710°C	
MW: 60.08	EINECS: 231-545-4	

Silica gel 60 Angstrom pore size, 35-70 micron
For chromatography column, Davisil® Silicagel

packaging		Mfr. No
500 g	Plastic jar	S/0693/53
1 kg	Plastic jar	S/0693/60
5 kg	Plastic keg	S/0693/65
25 kg	Fibre Keg	S/0693/70
SiO ₂	Boiling point: 2230°C	
CAS: 7631-86-9	Melting point: 1710°C	
MW: 60.08	EINECS: 231-545-4	

Silica gel 60 Angstrom pore size, 70-200 micron
For chromatography column, Davisil® Silicagel

packaging		Mfr. No
500 g	Plastic jar	S/0695/53
1 kg	Plastic jar	S/0695/60
25 kg	Fibre Keg	S/0695/70
SiO ₂	Boiling point: 2230°C	
CAS: 7631-86-9	Melting point: 1710°C	
MW: 60.08	EINECS: 231-545-4	

Silica gel 100 Angstrom pore size, 35-70 micron
For chromatography column, Davisil® Silicagel

packaging		Mfr. No
25 kg	Fibre Keg	S/0696/70
SiO ₂	Melting point: 1710°C	
CAS: 7631-86-9	EINECS: 231-545-4	
Boiling point: 2230°C		

Silica gel 150 Angstrom pore size, 35-70 micron
for chromatography column, Davisil Silicagel

packaging		Mfr. No
500 g	Plastic jar	S/0699/53
1 kg	Plastic jar	S/0699/60
SiO ₂		
CAS: 7631-86-9		
MW: 60.08		
EINECS: 231-545-4		

Silica gel grade 923 Davisil®
For chromatography column

packaging		Mfr. No
100 g	Plastic jar	S/0697/48
250 g	Plastic jar	S/0697/50
500 g	Plastic jar	S/0697/53
SiO ₂	Boiling point: 2230°C	
CAS: 7631-86-9	Melting point: 1710°C	
MW: 60.08	EINECS: 231-545-4	

Silica gel (2.5-6.0 mm) self indicating
dessicant agent

packaging		Mfr. No
500 g	Plastic jar	S/0762/53
1 kg	Plastic jar	S/0762/60
5 kg	Plastic keg	S/0762/65
SiO ₂	Boiling point: 2230°C	
CAS: 7631-86-9	Melting point: 1710°C	
MW: 60.08	EINECS: 231-545-4	
Density at 20°C: 2.20		
Colour change orange to green as water absorbed		

Silica gel (2.5-6.0 mm) self indicating (iron salt)
pure, Dessicant agent

packaging		Mfr. No
500 g	Plastic jar	S/0761/53
1 kg	Plastic jar	S/0761/60
5 kg	Plastic keg	S/0761/65
SiO ₂	Boiling point: 2230°C	
CAS: 7631-86-9	Melting point: 1710°C	
MW: 60.08	EINECS: 231-545-4	
Density at 20°C: 2.20		
Contains iron salts, Colour change from orange to colorless as water is absorbed.		

Silica gel (2.5-6.0 mm)
pure, Dessicant agent

packaging		Mfr. No
500 g	Plastic jar	S/0720/53
1 kg	Plastic jar	S/0720/60
SiO ₂		
CAS: 7631-86-9		
MW: 60.08		
EINECS: 231-545-4		

Silica gel (1.0-3.0 mm)
pure, Dessicant agent

packaging		Mfr. No
500 g	Plastic jar	S/0730/53
SiO ₂		
CAS: 7631-86-9		
MW: 60.08		
EINECS: 231-545-4		

Silica precipitated
pure

packaging		Mfr. No
500 g	Plastic jar	S/0680/53
SiO ₂	Boiling point: 2230°C	
CAS: 7631-86-9	Melting point: 1710°C	
MW: 60.09	EINECS: 231-545-4	

Silicon powder
technical

packaging		Mfr. No
500 g	Plastic jar	S/0850/53
CAS: 7440-21-3	UN: 1346	
Density at 20°C: 0.6	GHS: Flam. Sol. 2	
Melting point: 1410°C	H228	
ADR/RID: 4.1, III	P280, P210	
EINECS: 231-130-8		
Guaranteed analysis		
Identity test : Pass		

Silicon standard metal solution 1000 ppm
for AAS, solution aqueous

packaging		Mfr. No
100 mℓ	HDPE Plastic bottle	J/8062/05
500 mℓ	HDPE Plastic bottle	J/8062/08
Density at 20°C: 1		
EINECS: 240-968-3		
Manufactured to a tolerance of ± 0.5%. See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.		

Silicon standard metal solution 1000 ppm
for ICP solution in H₂O/HF trace, NIST certified

packaging		Mfr. No
100 mℓ	Plastic bottle	SPPLSI92Y
250 mℓ	Plastic bottle	SPPLSI92T
500 mℓ	Plastic bottle	SPPLSI92X
Manufactured to a tolerance of ± 0.5%. All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurites. See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314. Check with your local distributor for availability of Spex product		

Silicon standard metal solution 10 000 ppm
for ICP solution in H₂O/HF trace, NIST certified

packaging		Mfr. No
100 mℓ	Plastic bottle	SPPLSI93Y
500 mℓ	Plastic bottle	SPPLSI93X
Manufactured to a tolerance of ± 0.5%. All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurites. See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314. Check with your local distributor for availability of Spex product		

UHPLC

Filtered to 0.1 µm, Fisher Chemical UHPLC gradient grade solvents (Acetonitrile, Methanol and Water) are perfect for those working at high pressure.
See page 309 for more information



Silicon standard metal solution 10 000 ppm
for ICP solution in H₂O, NIST certified

packaging	Mfr. No
500 mℓ Plastic bottle	SPPLS19A3X
H314	



Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Silicon tetrachloride
specified reagents for general laboratory work- extra pure

SLR

packaging	Mfr. No
500 mℓ Coated Glass bottle	S/0844/PB08
SiCl ₄ CAS: 10026-04-7 MW: 169.90 Density at 20°C: 1.48 Boiling point: 58°C Melting point: -70°C ADR/RID: 8, II EINECS: 233-054-0	UN: 1818 GHS: Acute Tox. 3, Skin Corr. 1A, Eye Dam. 1, STOT SE 3 H335, H301, H331, H314, EUH014 P280, P305+P351+P338, P261, P402+P404



Silver standard metal solution 1000 ppm
for AAS, solution HNO₃ 1M

packaging	Mfr. No
100 mℓ HDPE Plastic bottle	J/8064/05
500 mℓ HDPE Plastic bottle	J/8064/08
ADR/RID: 8, III UN: 3264 GHS: Skin Corr. 1B H314	P301+P330+P331, P280, P305+P351+P338, P310, P303+P361+P353



Manufactured to a tolerance of ± 0.5%.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.

Silver standard metal solution 1000 ppm
for ICP, solution HNO₃ 1M

packaging	Mfr. No
500 mℓ HDPE Plastic bottle	J/8310/08
ADR/RID: 8, III UN: 3264 GHS: Skin Corr. 1B H314	P301+P330+P331, P280, P305+P351+P338, P310, P303+P361+P353



Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.

Silver standard metal solution 10 000 ppm
for ICP, solution HNO₃ 1M

packaging	Mfr. No
100 mℓ HDPE Plastic bottle	J/8309/05
ADR/RID: 8, III UN: 3264	

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.

Silver standard metal solution 1000 ppm
for ICP solution in HNO₃ 2%, NIST certified

packaging	Mfr. No
100 mℓ Plastic bottle	SPPLAG22Y
250 mℓ Plastic bottle	SPPLAG22T
500 mℓ Plastic bottle	SPPLAG22X
H315, H319	



Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Silver standard metal solution 10 000 ppm
for ICP solution in HNO₃ 5%, NIST certified

packaging	Mfr. No
100 mℓ Plastic bottle	SPPLAG23Y
500 mℓ Plastic bottle	SPPLAG23X
H314	



Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Silver standard metal solution 1000 ppm
for ICP-MS solution in HNO₃ 2%, NIST certified

packaging	Mfr. No
100 mℓ Plastic bottle	SPCLAG22Y
H315, H319	



Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Silver metal powder

specified reagents for general laboratory work- extra pure

SLR

packaging	Mfr. No
10 g Vial	S/0960/45
Ag CAS: 7440-22-4 MW: 107.87	Boiling point: 2212°C Melting point: 961°C EINECS: 231-131-3

Guaranteed analysis
Identity test : Pass

The prices apply to the silver and derivatives vary with the price of the raw material.

Silver metal foil

extra pure, 0.1 mm thick, approx 10 g

SLR

packaging	Mfr. No
10 g Plastic jar	S/0880/90
Ag CAS: 7440-22-4 MW: 107.87 Density at 20°C: 10.5	Boiling point: 2212°C Melting point: 961°C EINECS: 231-131-3

Guaranteed analysis
Identity test : Pass

The prices apply to the silver and derivatives vary with the price of the raw material.

Silver metal wire 1.2 mm

specified reagents for general laboratory work- extra pure

SLR

packaging	Mfr. No
10 g Glass jar	S/1010/45
Ag CAS: 7440-22-4 MW: 107.87	Boiling point: 2212°C Melting point: 961°C EINECS: 231-131-3

Guaranteed analysis
Identity test : Pass

The prices apply to the silver and derivatives vary with the price of the raw material.

Silver chloride

specified reagents for general laboratory work- extra pure

SLR

packaging	Mfr. No
25 g Glass jar	S/1120/46
AgCl CAS: 7783-90-6 MW: 143.32	Boiling point: 1150°C Melting point: 455°C EINECS: 232-033-3

Assay: >99%

Silver diethyldithiocarbamate

specified reagents for general laboratory work- extra pure

SLR

packaging	Mfr. No
5 g Glass jar	S/1200/44
C ₅ H ₁₀ AgNS ₂ CAS: 1470-61-7 MW: 256.14 Melting point: 172 to 175°C EINECS: 216-003-7 GHS: Acute Tox. 4, Skin Irrit.	2, Eye Irrit. 2, STOT SE 3 H335, H302, H315, H319 P261, P301+P312, P302+P352, P280, P305+P351+P338



Guaranteed analysis
Identity test : Pass

Silver nitrate
for analysis

Certified AR

packaging	Mfr. No
25 g Glass jar	S/1280/46
100 g Glass jar	S/1280/48
250 g Glass jar	S/1280/50
500 g Glass jar	S/1280/53
10 kg Plastic keg	S/1280/68

AgNO ₃ CAS: 7761-88-8 MW: 169.87 Boiling point: 444°C Melting point: 212°C ADR/RID: 5.1, II EINECS: 231-853-9 UN: 1493	GHS: Ox. Sol. 2, Skin Corr. 1B, Eye Dam. 1, Aquatic Acute 1, Aquatic Chronic 1 H272, H410, H314 P301+P330+P331, P280, P305+P351+P338, P310, P210, P273
--	--



Assay: >99.9%

Guaranteed analysis		
Water insoluble matter	<0.005%	Potassium (K) <20 ppm
Not precipitated by HCl	: <0.01%	Sodium (Na) <20 ppm
Calcium (Ca)	<10 ppm	Total chloride (Cl) <0.0002%
Copper (Cu)	<2 ppm	Total phosphorus (P) <10 ppm
Iron (Fe)	<2 ppm	Total silicon (Si) <5 ppm
Lead (Pb)	<2 ppm	Total sulphur (S) <20 ppm
Magnesium (Mg)	<10 ppm	Zinc (Zn) <5 ppm

Typical analysis		
Aluminium (Al)	<2 ppm	Manganese (Mn) <1 ppm
Barium (Ba)	<1 ppm	Molybdenum (Mo) <1 ppm
Cadmium (Cd)	<1 ppm	Nickel (Ni) <1 ppm
Chromium (Cr)	<1 ppm	Strontium (Sr) <1 ppm
Cobalt (Co)	<1 ppm	Titanium (Ti) <1 ppm
Lithium (Li)	<1 ppm	Vanadium (V) <1 ppm

Lot analysis on the label.

Silver nitrate

specified reagents for general laboratory work- extra pure

SLR

packaging	Mfr. No
25 g Glass jar	S/1240/46
100 g Glass jar	S/1240/48
250 g Glass jar	S/1240/50
500 g Glass jar	S/1240/53

AgNO ₃ CAS: 7761-88-8 MW: 169.87 Boiling point: 444°C Melting point: 212°C ADR/RID: 5.1, II EINECS: 231-853-9 UN: 1493	GHS: Ox. Sol. 2, Skin Corr. 1B, Eye Dam. 1, Aquatic Acute 1, Aquatic Chronic 1 H272, H410, H314 P301+P330+P331, P280, P305+P351+P338, P310, P210, P273
--	--



Assay: >99.8%

Guaranteed analysis

Calcium (Ca)	<100 ppm	Sodium (Na)	<200 ppm
Copper (Cu)	<20 ppm	n-Butyl formate	<0.002%
Iron (Fe)	<20 ppm	Total phosphorus (P)	<100 ppm
Lead (Pb)	<20 ppm	Total silicon (Si)	<50 ppm
Magnesium (Mg)	<100 ppm	Total sulphur (S)	<200 ppm
Potassium (K)	<200 ppm	Zinc (Zn)	<50 ppm

Silver nitrate solution 1000 g/ℓ

For COD

packaging	Mfr. No
500 mℓ Coated Glass bottle	J/7380/PB08
ADR/RID: 8, II UN: 1760 GHS: Skin Corr. 1B, Aquatic Chronic 2	H314, H411 P301+P330+P331, P280, P305+P351+P338, P310, P273



Silver nitrate solution 1M (1N)

Standard solution ready to use for volumetric analysis

packaging	Mfr. No
1 ℓ HDPE Plastic bottle	J/7320/PB15
2.5 ℓ HDPE Plastic bottle	J/7320/PB17

CAS: 7761-88-8 Density at 20°C: 1.13 ADR/RID: 8, II EINECS: 231-853-9 UN: 1760	GHS: Skin Corr. 1B, Aquatic Chronic 2 H314, H411 P301+P330+P331, P280, P305+P351+P338, P310, P273
--	---



Manufactured to a tolerance of ± 0.1%.
See the full range of Standard Volumetric solution page 323.

Silver nitrate solution 0.1M (0.1N)

NIST Standard solution ready to use, For Volumetric analysis

Eur.Ph. , USP, BP

packaging	Mfr. No
1 ℓ HDPE Plastic bottle	J/7330/15
2.5 ℓ HDPE Plastic bottle	J/7330/17
10 ℓ Bag in box	J/7330/24
25 ℓ HDPE Drum	J/7330/25

CAS: 7761-88-8 Density at 20°C: 1.015 EINECS: 231-853-9 GHS: Skin Irrit. 2, Eye Irrit. 2,	Aquatic Chronic 3 H319, H412, H315 P273, P302+P352, P280, P305+P351+P338
--	--



Manufactured to a tolerance of ± 0.1%.
See the full range of Standard Volumetric solution page 323.

Silver nitrate solution 0.1M (0.1N)			
NIST Standard Concentrate, For Volumetric analysis			
packaging		Mfr. No	
1	Ampoule	J/7330C/05	
6	Ampoule	J/7330C/90	
CAS: 7761-88-8		GHS: Skin Corr. 1B, Aquatic Chronic 2	
Density at 20°C: 1.13		H314, H411	
ADR/RID: 8, II		P301+P330+P331, P280,	
EINECS: 231-853-9		P305+P351+P338, P310, P273	
UN: 1760			

Each ampoule makes 1 ℓ of solution.
See the full range of Standard Volumetric solution page 323.

Silver nitrate solution 0.05M (0.05N)			
NIST Standard solution ready to use, For Volumetric analysis			
packaging		Mfr. No	
1 ℓ	HDPE Plastic bottle	J/7340/15	
2.5 ℓ	HDPE Plastic bottle	J/7340/17	
CAS: 7761-88-8		Aquatic Chronic 3	
Density at 20°C: 1.015		H319, H412, H315	
EINECS: 231-853-9		P273, P302+P352, P280,	
GHS: Skin Irrit. 2, Eye Irrit. 2,		P305+P351+P338	

Manufactured to a tolerance of ± 0.2%.
See the full range of Standard Volumetric solution page 323.

Silver nitrate solution 0.0282M (0.0282N)			
NIST Standard solution ready to use, For Volumetric analysis			
packaging		Mfr. No	
250 mℓ	Glass bottle	J/7370/07	
1 ℓ	HDPE Plastic bottle	J/7370/15	
2.5 ℓ	HDPE Plastic bottle	J/7370/17	
CAS: 7761-88-8		Aquatic Chronic 3	
Density at 20°C: 1.015		H319, H412, H315	
EINECS: 231-853-9		P273, P302+P352, P280,	
GHS: Skin Irrit. 2, Eye Irrit. 2,		P305+P351+P338	

Manufactured to a tolerance of ± 0.3%.
See the full range of Standard Volumetric solution page 323.

Silver nitrate solution 0.025M (0.025N)			
Standard solution ready to use, For volumetric analysis			
packaging		Mfr. No	
1 ℓ	Glass bottle	J/7310/15	
2.5 ℓ	HDPE Plastic bottle	J/7310/17	
CAS: 7761-88-8		Aquatic Chronic 3	
Density at 20°C: 1.015		H319, H412, H315	
ADR/RID: /II		P273, P302+P352, P280,	
EINECS: 231-853-9		P305+P351+P338	

See the full range of Standard Volumetric solution page 323.
Manufactured to a tolerance of ± 0.1%

Silver nitrate solution 0.02M (0.02N)			
NIST Standard solution ready to use, For Volumetric analysis			
packaging		Mfr. No	
1 ℓ	HDPE Plastic bottle	J/7350/15	
2.5 ℓ	HDPE Plastic bottle	J/7350/17	
CAS: 7761-88-8		Aquatic Chronic 3	
Density at 20°C: 1.015		H319, H412, H315	
EINECS: 231-853-9		P273, P302+P352, P280,	
GHS: Skin Irrit. 2, Eye Irrit. 2,		P305+P351+P338	

Manufactured to a tolerance of ± 0.3%.
See the full range of Standard Volumetric solution page 323.

Silver oxide			
specified reagents for general laboratory work- extra pure			
packaging		Mfr. No	
10 g	Glass jar	S/1320/45	
25 g	Glass jar	S/1320/46	
Ag ₂ O		UN: 3085	
CAS: 20667-12-3		GHS: Ox. Sol. 2, Skin Corr. 1B, Eye Dam. 1	
MW: 231.74		H272, H314	
Density at 20°C: 7.14		P301+P330+P331, P280,	
ADR/RID: 5.1, 8, II		P305+P351+P338, P310, P210	
EINECS: 243-957-1			

Guaranteed analysis
Identity test : Pass

Silver sulfate			
for analysis			
packaging		Mfr. No	
25 g	Glass jar	S/1400/46	
100 g	Glass jar	S/1400/48	
500 g	Glass jar	S/1400/53	
Ag ₂ SO ₄		EINECS: 233-653-7	
CAS: 10294-26-5		GHS: Eye Dam. 1	
MW: 311.79		H318	
Boiling point: 1085°C		P280, P305+P351+P338, P310	
Melting point: 652°C			
Assay: >99%			
Guaranteed analysis			
Calcium (Ca)	<20 ppm	Total nitrogen (N)	<20 ppm
Copper (Cu)	<5 ppm	Zinc (Zn)	<10 ppm
Iron (Fe)	<5 ppm	Copper (Cu)	<5 ppm
Lead (Pb)	<5 ppm	Iron (Fe)	<5 ppm
Magnesium (Mg)	<20 ppm	Lead (Pb)	<5 ppm
Potassium (K)	<50 ppm	n-Butyl formate	<0.001%
Sodium (Na)	<50 ppm	Total nitrogen (N)	<20 ppm
Total chloride (Cl)	<0.001%		
Typical analysis			
Aluminium (Al)	<5 ppm	Manganese (Mn)	<5 ppm
Barium (Ba)	<5 ppm	Molybdenum (Mo)	<5 ppm
Cadmium (Cd)	<5 ppm	Nickel (Ni)	<5 ppm
Chromium (Cr)	<5 ppm	Strontium (Sr)	<5 ppm
CO	<5 ppm	Titanium (Ti)	<5 ppm
Lithium (Li)	<5 ppm	Vanadium (V)	<5 ppm
Lot analysis on the label.			

Silver sulfate			
specified reagents for general laboratory work- extra pure			
packaging		Mfr. No	
25 g	Glass jar	S/1360/46	
100 g	Glass jar	S/1360/48	
Ag ₂ SO ₄		EINECS: 233-653-7	
CAS: 10294-26-5		GHS: Eye Dam. 1	
MW: 311.79		H318	
Boiling point: 1085°C		P280, P305+P351+P338, P310	
Melting point: 652°C			
Assay: >98.5%			

Guaranteed analysis			
Calcium (Ca)	<100 ppm	Potassium (K)	<200 ppm
Copper (Cu)	<20 ppm	Sodium (Na)	<200 ppm
Iron (Fe)	<20 ppm	n-Butyl formate	<0.002%
Lead (Pb)	<20 ppm	Total silicon (Si)	<50 ppm
Magnesium (Mg)	<100 ppm	Zinc (Zn)	<50 ppm

pH Buffers

NIST Certified Standard Solution & Concentrate

A full range to discover on page 63



Silver sulfate solution 5%w/v (d = 1.84) in sulfuric acid		
For COD		
Packaging		Mfr. No
2.5 ℓ	Coated Glass bottle	J/7410/PB17
Ag ₂ SO ₄		
MW: 311.79		
ADR/RID: 8,		
I UN: 1830		
Guaranteed analysis		
Assay (w/w)	: 4.99 to 5.01 %	Dichromate demand (mℓ FeSO ₄) ... <2 mℓ

Silver sulfate solution 1%w/v (d = 1.84) in sulfuric acid			
For COD			
packaging		Mfr. No	
1 ℓ	Coated Glass bottle	J/7400/PC15	
2.5 ℓ	Coated Glass bottle	J/7400/PB17	
Ag ₂ SO ₄		UN: 1830	
CAS: 10294-26-5		GHS: Skin Corr. 1A, Eye Dam. 1	
MW: 311.79		H314	
Density at 20°C: 1.84		P280, P301+P330+P331, P310, P305+P351+P338, P303+P361+P353	
ADR/RID: 8, II			
EINECS: 233-653-7			
Guaranteed analysis			
Assay (silver)		: 0.9 to 1.1%	Dichromate demand (mℓ FeSO ₄) ... <2 mℓ

Soap powder			
pure, hard			
packaging		Mfr. No	
500 g	Plastic jar	S/1450/53	
Guaranteed analysis			
Identity test : Pass			

Soap solution Clarks'			
pure			
packaging		Mfr. No	
2.5 ℓ	Glass bottle	J/7490/17	
Flash point: 25°C		GHS: Flam. Liq. 3	
ADR/RID: 3, III		H226	
UN: 1170		P210	

Soap solution Wanklyn			
pure			
packaging		Mfr. No	
2.5 ℓ	Glass bottle	J/7500/17	
Boiling point: 78°C		GHS: Flam. Liq. 3	
Flash point: 27°C		H226	
ADR/RID: 3, III		P210	
UN: 1170			

Soda lime, granules			
extra pure, 1.0-2.5mm not hygroscopic self indicating			
packaging		Mfr. No	
500 g	Plastic jar	S/1700/53	
CAS: 8006-28-8		H314	
ADR/RID: 8, III		P301+P330+P331, P280,	
UN: 1907		P305+P351+P338, P310	
GHS: Skin Corr. 1B, Eye Dam. 1			
Guaranteed analysis			
Identity test : Pass			
Virage from green to brown.			

Soda lime			
extra pure, 1.0-2.5mm not hygroscopic granules			
packaging		Mfr. No	
500 g	Plastic jar	S/1600/53	
CAS: 8006-28-8		H314	
ADR/RID: 8, III		P301+P330+P331, P280,	
UN: 1907		P305+P351+P338, P310	
GHS: Skin Corr. 1B, Eye Dam. 1			
Guaranteed analysis			
Identity test : Pass			

Soda lime			
extra pure, 2.5-5mm not hygroscopic granules			
packaging		Mfr. No	
250 g	Plastic jar	S/1560/50	
500 g	Plastic jar	S/1560/53	
1 kg	Plastic jar	S/1560/60	
5 kg	Plastic keg	S/1560/65	
CAS: 8006-28-8		H314	
ADR/RID: 8, III		P301+P330+P331, P280,	
UN: 1907		P305+P351+P338, P310	
GHS: Skin Corr. 1B, Eye Dam. 1			
Guaranteed analysis			
Identity test : Pass			

Sodium standard solution 1000 ppm (1 mg/mℓ)			
for Ion pair chromatography, Ready-to-use			
packaging		Mfr. No	
100 mℓ	HDPE Plastic bottle	J/4562/05	
Manufactured to a tolerance of ± 1%.			

Sodium standard metal solution 1000 ppm			
for AAS, solution HNO ₃ 1M			
packaging			Mfr. No
100 mℓ	HDPE Plastic bottle		J/8065/05
500 mℓ	HDPE Plastic bottle		J/8065/08
ADR/RID: 8, II		P280, P301+P330+P331,	
UN: 2031		P310, P303+P361+P353,	
GHS: Skin Corr. 1B, Eye Dam. 1		P305+P351+P338	
H314			
Manufactured to a tolerance of ± 0.5%.			
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.			

Sodium standard metal solution 10 000 ppm

for ICP, solution HNO₃ 1M

packaging		Mfr. No
100 mℓ	HDPE Plastic bottle	J/8315/05
ADR/RID: 8, II	H314	
UN: 2031	P280, P305+P351+P338,	
GHS: Skin Corr. 1B	P301+P330+P331, P310	

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.

Sodium standard metal solution 1000 ppm

for ICP solution in HNO₃ 2%, NIST certified

packaging		Mfr. No
100 mL	Plastic bottle	SPPLNA22Y
250 mL	Plastic bottle	SPPLNA22T
500 mL	Plastic bottle	SPPLNA22X
H315, H319		

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Sodium standard metal solution 1000 ppm

for ICP solution in HCL 2%, NIST certified

packaging		Mfr. No
500 mℓ	Plastic bottle	SPPLNA12X

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Sodium standard metal solution 10 000 ppm

for ICP solution in HNO₃ 5%, NIST certified

packaging		Mfr. No
100 mL	Plastic bottle	SPPLNA23Y
500 mL	Plastic bottle	SPPLNA23X
H314		

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Sodium standard metal solution 10 000 ppm

for ICP solution in HCL 5%, NIST certified

packaging		Mfr. No
500 mℓ	Plastic bottle	SPPLNA13X

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Sodium standard metal solution 1000 ppm

for ICP-MS solution in HNO₃ 2%, NIST certified

Packaging		Mfr. No
100 mL	Plastic bottle	SPCLNA22Y
H315, H319		

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Sodium

extra pure, metal sticks in liquid paraffin

packaging		Mfr. No
100 g	Aluminium jar	S/1920/48
Na	GHS: Water-react. 1, Skin Corr. 1B, Eye Dam. 1	
CAS: 7440-23-5	H260, H314, EUH014	
MW: 23.99	P301+P310,	
Boiling point: 883°C	P301+P330+P331, P280,	
Melting point: 98°C	P305+P351+P338, P310,	
ADR/RID: 4.3, I	P402+P404	
EINECS: 231-132-9		
UN: 1428		

Guaranteed analysis
Calcium (Ca) <1000 ppm n-Butyl formate <0.01%
Iron (Fe) <10 ppm

Sodium acetate anhydrous

for analysis

packaging		Mfr. No
250 g	Plastic jar	S/2120/50
500 g	Plastic jar	S/2120/53
1 kg	Plastic jar	S/2120/60
25 kg	Plastic keg	S/2120/70

C₂H₃NaO₂
CAS: 127-09-3
MW: 82.04
Melting point: 324°C
Flash point: >250°C
EINECS: 204-823-8

Assay: >99%
Guaranteed analysis
Acidity <0.0005 meq/g Magnesium (Mg) <5 ppm
Alkalinity <0.001 meq/g Potassium (K) <50 ppm
pH (aqueous solution 5%) 7 to 9.2 Total chloride (Cl) <0.0005%
Loss on drying <1% Total nitrogen (N) <20 ppm
Insoluble matter <0.01% Total phosphorus (P) <5 ppm
Calcium (Ca) <50 ppm Total silicon (Si) <20 ppm
Carbonate (Na₂CO₃) <0.02% Total sulphur (S) <20 ppm
Copper (Cu) <5 ppm Zinc (Zn) <10 ppm
Iron (Fe) <5 ppm Phosphate (PO₄) <0.001%
Lead (Pb) <5 ppm Sulfate (SO₄) <0.003%

Typical analysis
Aluminium (Al) <2 ppm Manganese (Mn) <2 ppm
Arsenic (As) <2 ppm Molybdenum (Mo) <2 ppm
Barium (Ba) <2 ppm Nickel (Ni) <2 ppm
Cadmium (Cd) <2 ppm Silver (Ag) <2 ppm
Chromium (Cr) <2 ppm Strontium (Sr) <2 ppm
Cobalt (Co) <2 ppm Titanium (Ti) <2 ppm
Lithium (Li) <2 ppm Vanadium (V) <2 ppm

Lot analysis on the label.

UHPLC

Filtered to 0.1 µm, Fisher Chemical UHPLC gradient grade solvents (Acetonitrile, Methanol and Water) are perfect for those working at high pressure.

See page 309 for more information



Sodium acetate anhydrous

specified reagents for general laboratory work- extra pure

packaging		Mfr. No
500 g	Plastic jar	S/2080/53
1 kg	Plastic jar	S/2080/60
5 kg	Plastic keg	S/2080/65

C₂H₃NaO₂
CAS: 127-09-3
MW: 82.04
Melting point: 324°C
Flash point: >250°C
EINECS: 204-823-8

Assay: >98%
Guaranteed analysis
Calcium (Ca) <200 ppm n-Butyl formate <0.01%
Copper (Cu) <20 ppm Total phosphorus (P) <10 ppm
Iron (Fe) <20 ppm Total silicon (Si) <100 ppm
Lead (Pb) <20 ppm Total sulphur (S) <100 ppm
Magnesium (Mg) <20 ppm Zinc (Zn) <50 ppm
Potassium (K) <1000 ppm

Sodium acetate, concentrated solution 0.1M, (0.5N)

Standard concentrate, For volumetric analysis

packaging		Mfr. No
1	Ampoule	J/7440C/05

Each ampoule makes up 500mℓ of solution

Sodium acetate trihydrate

for HPLC

packaging		Mfr. No
250 g	Plastic jar	S/2046/50
1 kg	Plastic jar	S/2046/60

C₂H₃NaO₂·3H₂O
CAS: 6131-90-4
MW: 136.08
Density at 20°C: 1.45
Melting point: 58°C
EINECS: 204-814-9

Assay: >99%
Guaranteed analysis
Reaction (5% solution) at 25°C ... : 7.5 to 9.2 Lead (Pb) <5 ppm
Iron (Fe) <3 ppm n-Butyl formate <0.0005%

Absorbance specification
Max. Absorbance (1cm cell v HPLC water)
254 nm (sol. aqu. 1 M) : <0.02 A.U.

Lot analysis on the label.

Filtered to 0.2 µm.

Sodium acetate trihydrate

HPLC for electrochemical detection

packaging		Mfr. No
250 g	Plastic jar	S/2052/50
C ₂ H ₃ NaO ₂ ·3H ₂ O	Density at 20°C: 1.45	
CAS: 6131-90-4	Melting point: 58°C	
MW: 136.08	EINECS: 204-814-9	

Assay: >99%
Guaranteed analysis
Electrochemical test Pass Iron (Fe) <3 ppm
Reaction (5% solution) at 25°C ... : 7.5 to 9.2 Lead (Pb) <5 ppm
n-Butyl formate <0.0005%

Absorbance specification
Max. Absorbance (1cm cell v HPLC water)
254 nm (sol. aqu. 1 M) : <0.02 A.U.

Tested by HPLC-EC to ensure low electrochemical response.

Detector - ESA Coulochem II

Mode-oxidative

Electrode- Graphite.

Lot analysis on the label.

Filtered to 0.2 µm.

Sodium acetate trihydrate

for analysis, crystals

packaging		Mfr. No
500 g	Plastic jar	S/2040/53
1 kg	Plastic jar	S/2040/60
2.5 kg	Plastic jar	S/2040/62
5 kg	Plastic keg	S/2040/65

C₂H₃NaO₂·3H₂O
CAS: 6131-90-4
MW: 136.08
Density at 20°C: 1.45
Melting point: 58°C
EINECS: 204-814-9

Assay: 99 to 101%
Guaranteed analysis
pH (aqueous solution 5%) 8 to 9 Potassium (K) <50 ppm
Substances <0.0005% Total chloride (Cl) <0.0005%
reducing KMnO₄ <0.005% Total nitrogen (N) <10 ppm
Calcium (Ca) <10 ppm Total phosphorus (P) <5 ppm
Copper (Cu) <5 ppm Total silicon (Si) <20 ppm
Iron (Fe) <5 ppm Total sulphur (S) <10 ppm
Lead (Pb) <5 ppm Zinc (Zn) <5 ppm
Magnesium (Mg) <5 ppm

Typical analysis
Aluminium (Al) <2 ppm Manganese (Mn) <1 ppm
Arsenic (As) <2 ppm Molybdenum (Mo) <2 ppm
Barium (Ba) <1 ppm Nickel (Ni) <1 ppm
Cadmium (Cd) <1 ppm Silver (Ag) <2 ppm
Chromium (Cr) <1 ppm Strontium (Sr) <5 ppm
Cobalt (Co) <1 ppm Titanium (Ti) <1 ppm
Lithium (Li) <1 ppm Vanadium (V) <1 ppm

Lot analysis on the label.

Sodium acetate trihydrate

extra pure, crystals

packaging		Mfr. No
500 g	Plastic jar	S/2000/53
1 kg	Plastic jar	S/2000/60
2.5 kg	Plastic jar	S/2000/62
5 kg	Plastic keg	S/2000/65

C₂H₃NaO₂·3H₂O
CAS: 6131-90-4
MW: 136.08
Density at 20°C: 1.45
Melting point: 58°C
EINECS: 204-814-9

Assay: 99 to 101%
Guaranteed analysis
pH (aqueous solution 5%) 7.5 to 9 Magnesium (Mg) <20 ppm
Loss on drying 39 to 40.5% Potassium (K) <200 ppm
Heavy metals <10 ppm Sulfate (SO₄) <200 ppm
Arsenic (As) <2 ppm n-Butyl formate <0.01%
Calcium (Ca) <50 ppm Total phosphorus (P) <20 ppm
Copper (Cu) <20 ppm Total silicon (Si) <100 ppm
Iron (Fe) <10 ppm Total sulphur (S) <50 ppm
Lead (Pb) <20 ppm Zinc (Zn) <20 ppm

Sodium alginate

specified reagents for general laboratory work- extra pure

Packaging		Mfr. No.
500 g	Plastic jar	S/2160/53
CAS: 9005-38-3		
Melting point: >300°C		

Guaranteed analysis
Complies with standard IR spectrum: Pass

UHPLC

Filtered to 0.1 µm, Fisher Chemical UHPLC gradient grade solvents (Acetonitrile, Methanol and Water) are perfect for those working at high pressure.

See page 309 for more information



Sodium aluminate		
Technical		
packaging		Mfr. No
1 kg	Plastic jar	S/2200/60
AlO ₂ .Na	UN: 2812	
CAS: 1302-42-7	GHS: Skin Corr. 1B, Eye Dam. 1	
MW: 81.97	H314	
Melting point: 1650°C	P301+P330+P331, P280,	
ADR/RID: 8, III	P305+P351+P338, P310	
EINECS: 215-100-1		
Guaranteed analysis		
Identity test		: Pass

Sodium arsenate hydrate		Certified AR
for analysis		
packaging		Mfr. No
250 g	Plastic jar	S/2272/50
Na ₂ HAsO ₄ ·7H ₂ O	Aquatic Acute 1, Aquatic Chronic 1	
CAS: 10048-95-0	H301, H331, H350, H400	
MW: 312.01	P201, P280, P301+P310, P304+P340, P310, P273	
ADR/RID: 6.1, II		
UN: 1685		
GHS: Carc. 1B, Acute Tox. 3,		
Assay: 99 to 102%		

Sodium arsenite		SLR
specified reagents for general laboratory work- extra pure		
packaging		Mfr. No
100 g	Plastic jar	S/2330/48
NaAsO ₂	EINECS: 232-070-5	
CAS: 7784-46-5	UN: 2027	
MW: 129.91	GHS: Acute Tox. 2, Aquatic	
Density at 20°C: 1.87	Acute 1, Aquatic Chronic 1	
Melting point: 615°C	H350, H331, H300, H310, H410	
ADR/RID: 6.1, II	P201, P313, P273, P280	
		

Sodium arsenite solution 0.05M (0.1N)			Mfr. No
NIST Standard solution ready to use, For Volumetric analysis			
packaging			
1 ℓ	HDPE Plastic bottle		J/7510/PB15
2.5 ℓ	HDPE Plastic bottle		J/7510/PB17
ADR/RID: 6.1, III	H350, H412		
UN: 1686	P281, P308+P313, P273		
GHS: Aquatic Chronic 3			
Manufactured to a tolerance of ± 0.2%.			
See the full range of Standard Volumetric solution page 323.			

Sodium azide		Certified AR	
for analysis			
packaging		Mfr. No	
100 g	Plastic jar	S/2380/48	
500 g	Plastic jar	S/2380/53	
NaN ₃		EINECS: 247-852-1	
CAS: 26628-22-8		UN: 1687	
MW: 65.01		GHS: Acute Tox. 2, Aquatic	
Density at 20°C: 1.85		Acute 1, Aquatic Chronic 1	
Boiling point: 300°C		H300, H410, EUH032	
Melting point: 275°C		P301+P310, P280,	
ADR/RID: 6.1, II		P302+P350, P273	
Assay: >99.5%			
Guaranteed analysis			
Calcium (Ca)	<20 ppm	Potassium (K)	<100 ppm
Copper (Cu)	<5 ppm	Total chloride (Cl)	<0.001%
Iron (Fe)	<5 ppm	Total silicon (Si)	<20 ppm
Lead (Pb)	<5 ppm	Total sulphur (S)	<50 ppm
Magnesium (Mg)	<10 ppm	Zinc (Zn)	<10 ppm
Typical analysis			
Aluminium (Al)	<5 ppm	Molybdenum (Mo)	<5 ppm
Barium (Ba)	<5 ppm	Nickel (Ni)	<5 ppm
Cadmium (Cd)	<5 ppm	Silver (Ag)	<5 ppm
Chromium (Cr)	<5 ppm	Strontium (Sr)	<5 ppm
Cobalt (Co)	<5 ppm	Titanium (Ti)	<5 ppm
Lithium (Li)	<5 ppm	Vanadium (V)	<5 ppm
Manganese (Mn)	<5 ppm		
Lot analysis on the label.			

Sodium azide			SLR
specified reagents for general laboratory work- extra pure			
packaging			Mfr. No
100 g	Plastic jar		S/2360/48
500 g	Plastic jar		S/2360/53
10 kg	Plastic keg		S/2360/68
NaN ₃		EINECS: 247-852-1	
CAS: 26628-22-8		UN: 1687	
MW: 65.01		GHS: Acute Tox. 2, Aquatic	
Density at 20°C: 1.85		Acute 1, Aquatic Chronic 1	
Boiling point: 300°C		H300, H410, EUH032	
Melting point: 275°C		P301+P310, P280,	
ADR/RID: 6.1, II		P302+P350, P273	
Assay: >99%			
Guaranteed analysis			
Calcium (Ca)	<100 ppm	Potassium (K)	<500 ppm
Copper (Cu)	<20 ppm	n-Butyl formate	<0.02%
Iron (Fe)	<20 ppm	Total silicon (Si)	<100 ppm
Lead (Pb)	<20 ppm	Total sulphur (S)	<200 ppm
Magnesium (Mg)	<50 ppm	Zinc (Zn)	<50 ppm

Sodium benzoate		Certified AR	
for analysis			
packaging		Mfr. No	
500 g	Plastic jar	S/2441/53	
C ₇ H ₅ NaO ₂	EINECS: 208-534-8		
CAS: 532-32-1	GHS: Eye Irrit. 2		
MW: 144.11	H319		
Melting point: 300°C	P280, P305+P351+P338		
Assay: >99.5%			
Guaranteed analysis			
pH (aqueous solution 5%)	7 to 8	Iron (Fe)	<1 ppm
Loss on drying	<1%	Lead (Pb)	<1 ppm
Aluminium (Al)	<2 ppm	Lithium (Li)	<2 ppm
Barium (Ba)	<2 ppm	Magnesium (Mg)	<50 ppm
Cadmium (Cd)	<2 ppm	Manganese (Mn)	<2 ppm
Calcium (Ca)	<100 ppm	Molybdenum (Mo)	<2 ppm
Chlorine (Cl) compounds	<0.01%	Nickel (Ni)	<2 ppm
Chromium (Cr)	<2 ppm	Potassium (K)	<500 ppm
Cobalt (Co)	<2 ppm	Propan--1-ol	<2 ppm
Copper (Cu)	<1 ppm		

Strontium (Sr)	<2 ppm	Total silicon (Si)	<50 ppm
Titanium (Ti)	<2 ppm	Total sulphur (S)	<100 ppm
n-Butyl formate	<0.005%	Vanadium (V)	<2 ppm
Total phosphorus (P)			
<50 ppm			
Lot analysis on the label.			

Sodium benzoate		SLR, Eur.Ph., BP	
specified reagents for general laboratory work- extra pure			
packaging		Mfr. No	
500 g	Plastic jar	S/2440/53	
1 kg	Plastic jar	S/2440/60	
25 kg	Plastic keg	S/2440/70	
C ₇ H ₅ NaO ₂	EINECS: 208-534-8		
CAS: 532-32-1	GHS: Eye Irrit. 2		
MW: 144.11	H319		
Melting point: 300°C	P280, P305+P351+P338		
Assay: >99%			
Guaranteed analysis			
Loss on drying	<2%	Potassium (K)	<1000 ppm
Calcium (Ca)	<500 ppm	n-Butyl formate	<0.05%
Copper (Cu)	<10 ppm	Total phosphorus (P)	<200 ppm
Iron (Fe)	<10 ppm	Total silicon (Si)	<200 ppm
Lead (Pb)	<10 ppm	Total sulphur (S)	<500 ppm
Magnesium (Mg)	<200 ppm	Zinc (Zn)	<100 ppm

Sodium bichromate	
See Sodium dichromate dihydrated, page 288	

Sodium bismuthate		SLR
specified reagents for general laboratory work- extra pure		
packaging		Mfr. No
100 g	Plastic jar	S/2480/48
NaBiO ₃	GHS: Acute Tox. 4, Skin Irrit.	
CAS: 12232-99-4	2, Eye Irrit. 2, STOT SE 3	
MW: 279.96	H302, H315, H319, H335	
EINECS: 235-455-6	P264, P270, P301+P330+P331	
Guaranteed analysis		
Identity test		: Pass

Sodium borate decahydrate	
See di-Sodium tetraborate decahydrated, page 305	

Sodium borohydride		SLR
specified reagents for general laboratory work- extra pure		
packaging		Mfr. No
10 g	Plastic jar	S/2560/45
25 g	Plastic jar	S/2560/46
100 g	Plastic jar	S/2560/48
NaBH ₄	UN: 1426	  
CAS: 16940-66-2	GHS: Water-react. 1, Acute	
MW: 37.83	Tox. 3, Skin Corr. 1C, Eye Dam. 1	
Melting point: 400°C	H260, H314, H301, H332	
Flash point: 69°C	P280, P305+P351+P338,	
ADR/RID: 4.3, I	P310, P223, P260, P264	
EINECS: 241-004-4		
Assay: >98%		
Guaranteed analysis		
Identity test : Pass		

Sodium bromate			SLR
specified reagents for general laboratory work- extra pure			
packaging			Mfr. No
500 g	Plastic jar		S/2600/53
NaBrO ₃	UN: 1494		
CAS: 7789-38-0	GHS: Ox. Sol. 1, Eye Irrit. 2,		
MW: 150.90	STOT SE 3		
Density at 20°C: 3.33	H271, H335, H319		
Melting point: 381°C	P210, P261, P280,		
ADR/RID: 5.1, II			
EINECS: 232-160-4			
Assay: >95%			

Sodium bromide			SLR
specified reagents for general laboratory work- extra pure			
packaging			Mfr. No
500 g	Plastic jar		S/2640/53
3 kg	Plastic jar		S/2640/63
NaBr		Boiling point: 1390°C	
CAS: 7647-15-6		Melting point: 755°C	
MW: 102.89		EINECS: 231-599-9	
Density at 20°C: 3.208			
Assay: >99%			
Guaranteed analysis			
Alkalinity	<0.002 meq/g	Lead (Pb)	<5 ppm
Loss on drying	<3%	Sulfate (SO ₄)	<0.02%
Iron (Fe)	<10 ppm	n-Butyl formate	<0.4%

Sodium cacodylate trihydrate			SLR
specified reagents for general laboratory work- extra pure			
packaging			Mfr. No
100 g	Plastic jar		S/2660/48
Na(CH ₃) ₂ AsO ₂ ·3H ₂ O	UN: 1688		
CAS: 124-65-2	GHS: Acute Tox. 3, Aquatic		
MW: 214.03	Acute 1, Aquatic Chronic 1		
Melting point: -61°C	H331, H410, H301, H350		
ADR/RID: 6.1, II	P301+P310, P201,		
EINECS: 204-708-2	P308+P313, P304+P340, P273		
Assay: >97%			
Guaranteed analysis			
Acidity	<0.07 meq/g	Sulfate (SO ₄)	<0.3%
Alkalinity	<0.025 meq/g	n-Butyl formate	<0.1%

Sodium carbonate anhydrous		Certified AR
for analysis		
packaging	Mfr. No	
500 g Plastic jar	S/2920/53	
1 kg Plastic jar	S/2920/60	
2.5 kg Plastic jar	S/2920/62	
5 kg Plastic keg	S/2920/65	
CNa ₂ O ₃ CAS: 497-19-8 MW: 105.99 Density at 20°C: 2.53 Boiling point: 1600°C	Melting point: 851°C EINECS: 207-838-8 GHS: Eye Irrit. 2 H319 P280, P305+P351+P338	
Assay: >99.9%		

Guaranteed analysis			
Water insoluble matter	<0.005%	Total chloride (Cl)	<0.001%
Loss on drying 300°C	: <1%	Total nitrogen (N)	<10 ppm
Calcium (Ca)	<100 ppm	Total phosphorus (P)	<5 ppm
Copper (Cu)	<5 ppm	Total silicon (Si)	<20 ppm
Iron (Fe)	<5 ppm	Total sulphur (S)	<10 ppm
Lead (Pb)	<5 ppm	Zinc (Zn)	<5 ppm
Magnesium (Mg)	<20 ppm	Phosphate (PO ₄)	<0.001%
Potassium (K)	<50 ppm		
Typical analysis			
Aluminium (Al)	<2 ppm	Manganese (Mn)	<1 ppm
Arsenic (As)	<2 ppm	Molybdenum (Mo)	<1 ppm
Barium (Ba)	<1 ppm	Nickel (Ni)	<1 ppm
Cadmium (Cd)	<1 ppm	Silver (Ag)	<1 ppm
Chromium (Cr)	<1 ppm	Strontium (Sr)	<2 ppm
Cobalt (Co)	<1 ppm	Titanium (Ti)	<1 ppm
Lithium (Li)	<1 ppm	Vanadium (V)	<1 ppm

Lot analysis on the label.

Sodium carbonate anhydrous
specified reagents for general laboratory work- extra pure

packaging		Mfr. No
500 g	Plastic jar	S/2880/53
1 kg	Plastic jar	S/2880/60
2.5 kg	Plastic jar	S/2880/62
5 kg	Plastic keg	S/2880/65
25 kg	Fibre Keg	S/2880/70
Na ₂ CO ₃ CAS: 497-19-8 MW: 105.99 Density at 20°C: 2.53 Boiling point: 1600°C		Melting point: 851°C EINECS: 207-838-8 GHS: Eye Irrit. 2 H319 P280, P305+P351+P338
		
Assay: >99.5%		

Guaranteed analysis			
Loss on drying	<2%	Potassium (K)	<200 ppm
Calcium (Ca)	<500 ppm	n-Butyl formate	<0.02%
Copper (Cu)	<20 ppm	Total phosphorus (P)	<50 ppm
Iron (Fe)	<20 ppm	Total silicon (Si)	<200 ppm
Lead (Pb)	<20 ppm	Total sulphur (S)	<200 ppm
Magnesium (Mg)	<100 ppm	Zinc (Zn)	<100 ppm

Sodium carbonate anhydrous
Technical

packaging		Mfr. No
2.5 kg	Plastic jar	S/2840/62
25 kg	Fibre Keg	S/2840/70
Na ₂ CO ₃ CAS: 497-19-8 MW: 105.99 Density at 20°C: 2.53 Boiling point: 1600°C		Melting point: 851°C EINECS: 207-838-8 GHS: Eye Irrit. 2 H319 P280, P305+P351+P338
		
Guaranteed analysis		
Identity test : Pass		

Sodium carbonate decahydrate		Certified AR	
for analysis			
Packaging		Mfr. No	
500 g	Plastic jar	S/2760/53	
1 kg	Plastic jar	S/2760/60	
5 kg	Plastic keg	S/2760/65	
Na ₂ CO ₃ ·10H ₂ O			
CAS: 6132-02-1			
GHS: Eye Irrit. 2			
H319			
P280, P305+P351+P338			
MW: 286.14			
Density at 20°C: 1.46			
Melting point: 34°C			
Assay: >99%			
Guaranteed analysis			
Aluminium (Al)	<5 ppm	Total chloride (Cl)	<0.0015%
Calcium (Ca)	<20 ppm	Total nitrogen (N)	<10 ppm
Iron (Fe)	<2 ppm	Total phosphorus (P)	<5 ppm
Lead (Pb)	<5 ppm	Total silicon (Si)	<50 ppm
Magnesium (Mg)	<2 ppm	Total sulphur (S)	<20 ppm
Potassium (K)	<50 ppm		

Typical analysis			
Arsenic (As)	<2 ppm	Molybdenum (Mo)	<2 ppm
Barium (Ba)	<2 ppm	Nickel (Ni)	<2 ppm
Cadmium (Cd)	<2 ppm	Silver (Ag)	<2 ppm
Chromium (Cr)	<2 ppm	Strontium (Sr)	<2 ppm
Cobalt (Co)	<2 ppm	Titanium (Ti)	<2 ppm
Lithium (Li)	<2 ppm	Vanadium (V)	<2 ppm
Manganese (Mn)	<2 ppm		

Lot analysis on the label.

Sodium carbonate decahydrate
specified reagents for general laboratory work- extra pure

Packaging		Mfr. No	
1 kg	Plastic jar	S/2720/60	
2.5 kg	Plastic jar	S/2720/62	
5 kg	Plastic keg	S/2720/65	
Na ₂ CO ₃ ·10H ₂ O CAS: 6132-02-1 MW: 286.14 Density at 20°C: 1.46 Melting point: 34°C		EINECS: 207-838-8 GHS: Eye Irrit. 2 H319 P280, P305+P351+P338	
			
Assay: >98%			
Guaranteed analysis			
Calcium (Ca)	<200 ppm	n-Butyl formate	<0.02%
Iron (Fe)	<20 ppm	Total phosphorus (P)	<50 ppm
Lead (Pb)	<20 ppm	Total silicon (Si)	<200 ppm
Magnesium (Mg)	<100 ppm	Total sulphur (S)	<200 ppm
Potassium (K)	<200 ppm		

Sodium carbonate decahydrate
technical

Packaging		Mfr. No
3 kg	Plastic jar	S/2700/63
Na ₂ CO ₃ ·10H ₂ O CAS: 6132-02-1 MW: 286.14 Density at 20°C: 1.46	Melting point: 34°C GHS: Eye Irrit. 2 H319 P280, P305+P351+P338	
Guaranteed analysis		
Identity test : Pass		

Sodium carbonate solution 0.05M (0.1N)
NIST Standard solution ready to use, For Volumetric analysis

Packaging		Mfr. No
1 ℓ	HDPE Plastic bottle	J/7580/PB15
2.5 ℓ	HDPE Plastic bottle	J/7580/17
262		
Manufactured to a tolerance of ± 0.1%.		
See the full range of Standard Volumetric solution page 323.		

Sodium carbonate solution 0.05M (0.1N)
NIST Standard Concentrate, For Volumetric analysis

packaging		Mfr. No
1	Ampoule	J/7580C/05
6	Ampoule	J/7580C/90
P262		
<i>Each ampoule makes 1 l of solution.</i>		
<i>See the full range of Standard Volumetric solution page 323.</i>		

Sodium carbonate solution 0.5M (1N)
NIST Standard solution ready to use, For Volumetric analysis

packaging		Mfr. No
1 ℓ	HDPE Plastic bottle	J/7520/PB15
2.5 ℓ	HDPE Plastic bottle	J/7520/17
P262		
Manufactured to a tolerance of ± 0.1%.		
See the full range of Standard Volumetric solution page 323.		

Sodium chlorate
for analysis

packaging		Mfr. No
250 g	Plastic jar	S/3045/50
NaClO ₃	EINECS: 231-887-4	
CAS: 7775-09-9	UN: 1495	
MW: 106.44	GHS: Ox. Sol. 1, Acute Tox. 4, Aquatic Chronic 2	
Density at 20°C: 2.49	H271, H302, H411	
Melting point: 248 to 261°C	P210, P273, P301+P312	
ADR/RID: 5.1, II		

Guaranteed analysis			
Aluminium (Al)	<5 ppm	Nickel (Ni)	<5 ppm
Barium (Ba)	<5 ppm	Potassium (K)	<100 ppm
Cadmium (Cd)	<5 ppm	Propan--1-ol	<5 ppm
Calcium (Ca)	<20 ppm	Strontium (Sr)	<5 ppm
Chromium (Cr)	<5 ppm	Titanium (Ti)	<5 ppm
Cobalt (Co)	<5 ppm	n-Butyl formate	<0.01%
Copper (Cu)	<5 ppm	Total nitrogen (N)	<100 ppm
Iron (Fe)	<5 ppm	Total phosphorus (P)	<50 ppm
Lead (Pb)	<5 ppm	Total silicon (Si)	<50 ppm
Lithium (Li)	<5 ppm	Total sulphur (S)	<50 ppm
Magnesium (Mg)	<10 ppm	Vanadium (V)	<5 ppm
Manganese (Mn)	<5 ppm	Zinc (Zn)	<20 ppm
Molybdenum (Mo)	<5 ppm		

Sodium chlorate
specified reagents for general laboratory work- extra pure

packaging		Mfr. No
500 g	Plastic jar	S/3040/53
1 kg	Plastic jar	S/3040/60
NaClO ₃		EINECS: 231-887-4
CAS: 7775-09-9		UN: 1495
MW: 106.44		GHS: Ox. Sol. 1, Acute Tox.
Density at 20°C: 2.49		4, Aquatic Chronic 2
Melting point: 248 to 261°C		H271, H302, H411
ADR/RID: 5.1, II		P210, P273, P301+P312
		
		
		

Guaranteed analysis			
Calcium (Ca)	<200 ppm	n-Butyl formate	<0.05%
Copper (Cu)	<50 ppm	Total phosphorus (P)	<200 ppm
Iron (Fe)	<50 ppm	Total silicon (Si)	<200 ppm
Lead (Pb)	<50 ppm	Total sulphur (S)	<500 ppm
Magnesium (Mg)	<100 ppm	Zinc (Zn)	<200 ppm

Sodium chlorate
Technical

packaging	Mfr. No	
1 kg Plastic jar	S/3000/60	
3 kg Plastic jar	S/3000/63	
NaClO ₃ CAS: 7775-09-9 MW: 106.44 Density at 20°C: 2.49 Melting point: 248 to 261°C ADR/RID: 5.1, II	EINECS: 231-887-4 UN: 1495 GHS: Ox. Sol. 1, Acute Tox. 4, Aquatic Chronic 2 H271, H302, H411 P210, P273, P301+P312	

Sodium chloride
for HPLC

Packaging		Mfr. No	
500 g	Plastic jar	S/3165/53	
NaCl		Boiling point: 1461°C	
CAS: 7647-14-5		Melting point: 801°C	
MW: 58.44		EINECS: 231-598-3	
Density at 20°C: 2.165			
Assay: 0 to 150%			
Guaranteed analysis			
Bromide and iodide	<0.005%	Potassium (K)	<50 ppm
Iron (Fe)	<1 ppm	Total phosphorus (P)	<2 ppm
Lead (Pb)	<2 ppm	Total sulphur (S)	<10 ppm
Absorbance specification			
Max. Absorbance (1cm cell v HPLC water)			
206 nm	: <1 A.U.	254 nm	<0.02 A.U.
Lot analysis on the label.			
Filtered to 0.2 µm.			

Sodium chloride
for analysis

Packaging		Mfr. No	
500 g	Plastic jar	S/3161/53	
1 kg	Plastic jar	S/3161/60	
3 kg	Plastic jar	S/3161/63	
5 kg	Plastic jar	S/3161/65	
10 kg	Plastic keg	S/3161/68	
25 kg	Plastic keg	S/3161/70	
NaCl		Boiling point: 1461°C Melting point: 801°C EINECS: 231-598-3	
CAS: 7647-14-5			
MW: 58.44			
Assay: 99.5 to 100.5%			
Guaranteed analysis			
Alkalinity	Pass	Arsenic (As)	<1 ppm
Loss on drying	<0.5%	Barium (Ba)	Pass
Bacterial endotoxins	: <5	Bromide (Br)	<100 ppm
Ferrocyanide (Pass/Fail)	: Pass	Iodide (I)	Pass
Heavy metals	<5 ppm	Iron (Fe)	<2 ppm
Mg and alkaline earth metals	: <100 ppm	Phosphate (PO ₄)	<25 ppm
Nitrite at 354 nm	: <0.01 A.U.	Potassium (K)	<500 ppm
Acidity	Pass	Sulfate (SO ₄)	<200 ppm
Aluminium (Al)	<0.2 ppm		

Lot analysis on the label

High Volume Solvent Delivery Systems

High-volume solvent delivery systems are available in 10L to 1000L, in Europe only. Please email your contact details and a brief description of your requirements to bege1.sha.quotes@thermofisher.com

Sodium chloride
for analysis

Certified AR, Eur.Ph.

packaging		Mfr. No
500 g	Plastic jar	S/3160/53
1 kg	Plastic jar	S/3160/60
3 kg	Plastic jar	S/3160/63
5 kg	Plastic jar	S/3160/65
10 kg	Plastic keg	S/3160/68
25 kg	Plastic keg	S/3160/70

NaCl	Boiling point: 1461°C
CAS: 7647-14-5	Melting point: 801°C
MW: 58.44	EINECS: 231-598-3
Density at 20°C: 2.165	

Assay: >99.9%

Guaranteed analysis			
Water insoluble matter	<0.005%	Potassium (K)	<50 ppm
pH (aqueous solution 5%)	5 to 9	Total nitrogen (N)	<5 ppm
Insoluble matter	<0.005%	Total phosphorus (P)	<2 ppm
Bromide and iodide	<0.002%	Total silicon (Si)	<20 ppm
Calcium (Ca)	<20 ppm	Total sulphur (S)	<20 ppm
Copper (Cu)	<2 ppm	Zinc (Zn)	<5 ppm
Iron (Fe)	<2 ppm	Phosphate (PO ₄)	<5 ppm
Iron cyanide	<0.0001 %	Sulfate (SO ₄)	<0.004%
Lead (Pb)	<2 ppm	Ferrocyanide	<0.0001%
Magnesium (Mg)	<5 ppm		

Typical analysis			
Aluminium (Al)	<2 ppm	Manganese (Mn)	<1 ppm
Arsenic (As)	<2 ppm	Molybdenum (Mo)	<1 ppm
Barium (Ba)	<2 ppm	Nickel (Ni)	<1 ppm
Cadmium (Cd)	<1 ppm	Silver (Ag)	<1 ppm
Chromium (Cr)	<1 ppm	Strontium (Sr)	<1 ppm
Cobalt (Co)	<1 ppm	Titanium (Ti)	<1 ppm
Lithium (Li)	<1 ppm	Vanadium (V)	<1 ppm

Lot analysis on the label.

Sodium chloride
specified reagents for general laboratory work- extra pure

SLR

packaging		Mfr. No
500 g	Plastic jar	S/3120/53
1 kg	Plastic jar	S/3120/60
3 kg	Plastic jar	S/3120/63
5 kg	Plastic jar	S/3120/65
25 kg	Plastic keg	S/3120/70

NaCl	Boiling point: 1461°C
CAS: 7647-14-5	Melting point: 801°C
MW: 58.44	EINECS: 231-598-3
Density at 20°C: 2.165	

Assay: >99.5%

Guaranteed analysis			
Heavy metals	<10 ppm	Magnesium (Mg)	<20 ppm
Bromide (Br)	<0.005%	Potassium (K)	<200 ppm
Calcium (Ca)	<100 ppm	Total phosphorus (P)	<20 ppm
Copper (Cu)	<10 ppm	Total silicon (Si)	<100 ppm
Iron (Fe)	<10 ppm	Total sulphur (S)	<100 ppm
Lead (Pb)	<10 ppm	Zinc (Zn)	<20 ppm

Sodium chloride (rock salt)
Technical

packaging		Mfr. No
3 kg	Plastic jar	S/3100/63

NaCl	Boiling point: 1461°C
CAS: 7647-14-5	Melting point: 801°C
MW: 58.44	EINECS: 231-598-3
Density at 20°C: 2.165	

Guaranteed analysis
Identity test : Pass

Sodium chloride
Technical

packaging		Mfr. No
3 kg	Plastic jar	S/3105/63
25 kg	Plastic keg	S/3105/70

NaCl	Boiling point: 1461°C
CAS: 7647-14-5	Melting point: 801°C
MW: 58.44	EINECS: 231-598-3
Density at 20°C: 2.165	

Sodium chloride solution 0.1M (0.1N)
NIST Standard Concentrate, For Volumetric analysis

packaging		Mfr. No
1	Ampoule	J/7590C/05
6	Ampoule	J/7590C/90

NaCl	
CAS: 7647-14-5	
MW: 58.44	
EINECS: 231-598-3	

Each ampoule makes 1 ℓ of solution.
See the full range of Standard Volumetric solution page 323.

Sodium chloroacetate
specified reagents for general laboratory work- extra pure

SLR

packaging		Mfr. No
250 g	Plastic jar	S/3200/50

C ₂ H ₃ ClNaO ₂	UN: 2659	
CAS: 3926-62-3	GHS: Acute Tox. 3, Skin Irrit.	
MW: 116.48	2, Aquatic Acute 1	
Melting point: 199°C	H301, H315, H400	
ADR/RID: 6.1, III	P301+P310, P273, P302+P352	
EINECS: 223-498-3		

Guaranteed analysis
Identity test : Pass

Sodium chromate anhydrous
specified reagents for general laboratory work- extra pure

SLR

packaging		Mfr. No
500 g	Plastic jar	S/3240/53

Na ₂ CrO ₄	RE 1, Aquatic Acute 1, Aquatic Chronic 1	
CAS: 7775-11-3	H301, H312, H372, H314, H330, H334, H340, H410, H317, H350, H360FD	
MW: 161.97	P280, P301+P330+P331, P305+P351+P338, P304+P341, P201, P308+P313, P304+P340, P302+P352, P260, P273	
Melting point: 792°C		
ADR/RID: 6.1, II		
EINECS: 231-889-5		
UN: 3288		
GHS: Carc. 1B, Muta. 1B, Repr. 1B, Acute Tox. 2, Skin Corr. 1B, Eye Dam. 1, Resp. Sens. 1, Skin Sens. 1, STOT		

Guaranteed analysis
Identity test : Pass

Standard Volumetric solution

For your convenience we offer a range of NIST traceable solutions. A full list of our products can be found on page 323. Can't see what you need? Please email your contact details and a brief description of your requirements to beigel.sha.quotes@thermofisher.com

tri-Sodium citrate dihydrate
HPLC for electrochemical detection

Certified

packaging		Mfr. No
500 g	Plastic jar	S/3380/53

C ₆ H ₅ Na ₃ O ₇ ·2H ₂ O	Melting point: 300°C
CAS: 6132-04-3	EINECS: 200-675-3
MW: 294.09	

Assay: 99.5 to 100.5%

Guaranteed analysis			
Electrochemical test	Pass	Lead (Pb)	<5 ppm
Reaction (5% solution) at 25°C	: 8 to 9	n-Butyl formate	<0.005%
Iron (Fe)	<5 ppm		

Absorbance specification
Max. Absorbance (1 cm cell v HPLC water)
254 nm (sol. aqu. 1 M) : <0.1 A.U. 270 nm (sol. aqu. 1 M) : <0.05 A.U.

Purified then tested by HPLC-EC to show low electrochemical response.
Detector - ESA Coulochem II
Mode-oxydative
Electrode- Graphite.
Lot analysis on the label.
Filtered to 0.2 µm.

tri-Sodium citrate dihydrate
for electrophoresis

packaging		Mfr. No
500 g	Plastic jar	S/P500/53
5 kg	Plastic keg	S/P500/65

C ₆ H ₅ Na ₃ O ₇ ·2H ₂ O	Melting point: 300°C
CAS: 6132-04-3	EINECS: 200-675-3
MW: 294.09	
Melting point: 300°C	

Assay: 99.5 to 100.5%

Suitable as a buffer component in electrophoresis.

tri-Sodium citrate dihydrate
for analysis

Certified AR, Eur.Ph., BP, USP

packaging		Mfr. No
500 g	Plastic jar	S/3320/53
1 kg	Plastic jar	S/3320/60
5 kg	Plastic keg	S/3320/65
25 kg	Fibre Keg	S/3320/70

C ₆ H ₅ Na ₃ O ₇ ·2H ₂ O	Melting point: 300°C
CAS: 6132-04-3	EINECS: 200-675-3
MW: 294.09	

Assay: 99.5 to 100.5%

Guaranteed analysis			
pH (aqueous solution 5%)	8 to 9	Total chloride (Cl)	<0.005%
Calcium (Ca)	<20 ppm	Total nitrogen (N)	<20 ppm
Copper (Cu)	<5 ppm	Total phosphorus (P)	<10 ppm
Iron (Fe)	<5 ppm	Total silicon (Si)	<20 ppm
Lead (Pb)	<5 ppm	Total sulphur (S)	<20 ppm
Magnesium (Mg)	<5 ppm	Zinc (Zn)	<2 ppm
Potassium (K)	<100 ppm	Readily carbonisable substances test:	Pass

Typical analysis			
Aluminium (Al)	<2 ppm	Manganese (Mn)	<1 ppm
Arsenic (As)	<2 ppm	Molybdenum (Mo)	<1 ppm
Barium (Ba)	<1 ppm	Nickel (Ni)	<1 ppm
Cadmium (Cd)	<1 ppm	Silver (Ag)	<1 ppm
Chromium (Cr)	<1 ppm	Strontium (Sr)	<1 ppm
Cobalt (Co)	<1 ppm	Titanium (Ti)	<1 ppm
Lithium (Li)	<1 ppm	Vanadium (V)	<1 ppm

Lot analysis on the label.

tri-Sodium citrate dihydrate
specified reagents for general laboratory work- extra pure

SLR

packaging		Mfr. No
500 g	Plastic jar	S/3280/53
1 kg	Plastic jar	S/3280/60
5 kg	Plastic keg	S/3280/65

C ₆ H ₅ Na ₃ O ₇ ·2H ₂ O	Melting point: 300°C
CAS: 6132-04-3	EINECS: 200-675-3
MW: 294.09	

Assay: >99%

Guaranteed analysis			
pH (aqueous solution 5%)	8 to 9	Potassium (K)	<200 ppm
Calcium (Ca)	<100 ppm	n-Butyl formate	<0.01%
Copper (Cu)	<20 ppm	Total phosphorus (P)	<50 ppm
Iron (Fe)	<20 ppm	Total silicon (Si)	<100 ppm
Lead (Pb)	<20 ppm	Total sulphur (S)	<100 ppm
Magnesium (Mg)	<20 ppm	Zinc (Zn)	<10 ppm

Sodium cyanide
for analysis

Certified AR

packaging		Mfr. No
250 g	Plastic jar	S/3481/50

NaCN	GHS: Acute Tox. 1, Aquatic Acute 1, Aquatic Chronic 1	
CAS: 143-33-9	H330, H310, H410, H300, EUH032	
MW: 49.01	P301+P310, P280, P361, P302+P350,	
Boiling point: 1497°C	P305+P351+P338, P310, P304+P340, P273	
Melting point: 562°C		
ADR/RID: 6.1, I		
EINECS: 205-599-4		
UN: 1689		

Assay: >99%

Guaranteed analysis			
Calcium (Ca)	<50 ppm	n-Butyl formate	<0.01%
Copper (Cu)	<5 ppm	Total phosphorus (P)	<50 ppm
Iron (Fe)	<50 ppm	Total silicon (Si)	<50 ppm
Lead (Pb)	<5 ppm	Total sulphur (S)	<100 ppm
Magnesium (Mg)	<10 ppm	Zinc (Zn)	<5 ppm
Potassium (K)	>0 ppm		

Typical analysis			
Aluminium (Al)	<5 ppm	Molybdenum (Mo)	<2 ppm
Barium (Ba)	<2 ppm	Nickel (Ni)	<2 ppm
Cadmium (Cd)	<2 ppm	Silver (Ag)	<2 ppm
Chromium (Cr)	<2 ppm	Strontium (Sr)	<2 ppm
Cobalt (Co)	<2 ppm	Titanium (Ti)	<2 ppm
Lithium (Li)	<2 ppm	Vanadium (V)	<2 ppm
Manganese (Mn)	<2 ppm		

Lot analysis on the label

Sodium cyanide
specified reagents for general laboratory work- extra pure

SLR

packaging		Mfr. No
500 g	Plastic jar	S/3480/53

NaCN	GHS: Acute Tox. 1, Aquatic Acute 1, Aquatic Chronic 1	
CAS: 143-33-9	H330, H310, H410, H300, EUH032	
MW: 49.01	P301+P310, P280, P361, P302+P350,	
Density at 20°C: 16	P305+P351+P338, P310, P304+P340, P273	
Boiling point: 1497°C		
Melting point: 562°C		
ADR/RID: 6.1, I		
EINECS: 205-599-4		
UN: 1689		

Assay: >98%

Guaranteed analysis			
Calcium (Ca)	<200 ppm	n-Butyl formate	<0.05%
Copper (Cu)	<20 ppm	Total phosphorus (P)	<200 ppm
Iron (Fe)	<200 ppm	Total silicon (Si)	<200 ppm
Lead (Pb)	<20 ppm	Total sulphur (S)	<500 ppm
Magnesium (Mg)	<50 ppm	Zinc (Zn)	<20 ppm
Potassium (K)	<1000 ppm		

Sodium dichromate dihydrate Certified AR for analysis

packaging		Mfr. No
500 g	Plastic jar	S/3561/53
Na ₂ Cr ₂ O ₇ ·2H ₂ O	Dam. 1, Resp. Sens. 1, Skin Sens. 1, STOT SE 3, STOT RE 1, Aquatic Acute 1, Aquatic Chronic 1	
CAS: 7789-12-0	H272, H340, H317, H330, H334, H314, H350, H360FD, H372, H301, H312, H335, H410	
MW: 298.00	P280, P305+P351+P338, P310, P201, P210, P273	
Boiling point: 400°C		
Melting point: 357°C		
ADR/RID: 5.1, 6.1, II		
EINECS: 234-190-3		
UN: 3087		
GHS: Ox. Sol. 2, Carc. 1B, Muta. 1B, Repr. 1B, Acute Tox. 2, Skin Corr. 1B, Eye		

Assay: 99.5 to 101%			
Guaranteed analysis			
Aluminium (Al)	<5 ppm	Magnesium (Mg)	<20 ppm
Barium (Ba)	<1 ppm	Nickel (Ni)	<1 ppm
Cadmium (Cd)	<1 ppm	Potassium (K)	<200 ppm
Calcium (Ca)	<20 ppm	Titanium (Ti)	<1 ppm
Cobalt (Co)	<1 ppm	n-Butyl formate	<0.002%
Copper (Cu)	<10 ppm	Total silicon (Si)	<50 ppm
Iron (Fe)	<20 ppm	Total sulphur (S)	<100 ppm
Lead (Pb)	<20 ppm	Vanadium (V)	<1 ppm
Lithium (Li)	<1 ppm	Zinc (Zn)	<20 ppm

Lot analysis on the label.

Sodium dichromate dihydrate SLR specified reagents for general laboratory work- extra pure

packaging		Mfr. No
500 g	Plastic jar	S/3560/53
1 kg	Plastic jar	S/3560/60
3 kg	Plastic jar	S/3560/63
Na ₂ Cr ₂ O ₇ ·2H ₂ O	Dam. 1, Resp. Sens. 1, Skin Sens. 1, STOT SE 3, STOT RE 1, Aquatic Acute 1, Aquatic Chronic 1	
CAS: 7789-12-0	H272, H340, H317, H330, H334, H314, H350, H360FD, H372, H301, H312, H335, H410	
MW: 298.00	P280, P305+P351+P338, P310, P201, P210, P273	
Boiling point: 400°C		
Melting point: 357°C		
ADR/RID: 5.1, 6.1, II		
EINECS: 234-190-3		
UN: 3087		
GHS: Ox. Sol. 2, Carc. 1B, Muta. 1B, Repr. 1B, Acute Tox. 2, Skin Corr. 1B, Eye		

Assay: >99%			
Guaranteed analysis			
Calcium (Ca)	<100 ppm	Potassium (K)	<1000 ppm
Copper (Cu)	<50 ppm	n-Butyl formate	<0.02%
Iron (Fe)	<100 ppm	Total silicon (Si)	<200 ppm
Lead (Pb)	<100 ppm	Total sulphur (S)	<500 ppm
Magnesium (Mg)	<200 ppm	Zinc (Zn)	<100 ppm

Reagents, Solvents & Volumetric solutions comply with the European Pharmacopeia

Read more on page 234

Sodium dichromate dihydrate Technical

packaging		Mfr. No
1 kg	Plastic jar	S/3520/60
3 kg	Plastic jar	S/3520/63
Na ₂ Cr ₂ O ₇ ·2H ₂ O	Dam. 1, Resp. Sens. 1, Skin Sens. 1, STOT SE 3, STOT RE 1, Aquatic Acute 1, Aquatic Chronic 1	
CAS: 7789-12-0	H272, H340, H317, H330, H334, H314, H350, H360FD, H372, H301, H312, H335, H410	
MW: 298.00	P280, P305+P351+P338, P310, P201, P210, P273	
Boiling point: 400°C		
Melting point: 357°C		
ADR/RID: 5.1, 6.1, II		
EINECS: 234-190-3		
UN: 3087		
GHS: Ox. Sol. 2, Carc. 1B, Muta. 1B, Repr. 1B, Acute Tox. 2, Skin Corr. 1B, Eye		

Assay: >98%

Sodium diethyldithiocarbamate trihydrate SLR specified reagents for general laboratory work- extra pure

packaging		Mfr. No
25 g	Plastic jar	S/3640/46
100 g	Plastic jar	S/3640/48
C ₃ H ₁₀ NNaS ₂ ·3H ₂ O	GHS: Acute Tox. 4, Skin Irrit. 2, Eye Irrit. 2, STOT SE 3, Aquatic Acute 1	
CAS: 20624-25-3	H302, H335, H315, H319, H400	
MW: 225.31	P280, P270, P305+P351+P338, P310, P273, P261	
Melting point: 95 to 98.5°C		
ADR/RID: 9, III		
EINECS: 205-710-6		
UN: 3077		

Assay: >98.5%
Guaranteed analysis
Sensitivity to Cu molar absorptivity : > 12000

Sodium dihydrogen orthophosphate dihydrate for electrophoresis, monobasic

packaging		Mfr. No
500 g	Plastic jar	S/P520/53
3 kg	Plastic jar	S/P520/63
NaH ₂ PO ₄ ·2H ₂ O	Density at 20°C: ~ 2.0	
CAS: 13472-35-0	Melting point: 60°C	
MW: 156.01		

Assay: 99.5 to 100.5%			
Guaranteed analysis			
pH (aqueous solution 5%)	4.3 to 4.5	Lead (Pb)	<1 ppm
Copper (Cu)	<1 ppm	n-Butyl formate	<0.001%
Iron (Fe)	<5 ppm	Total sulphur (S)	<10 ppm

Suitable for use in in electrophoresis.

Safe, environmentally friendly and effective packaging

Fisher Chemical outer packaging is UN compliant and 100% recyclable. It's an easy and safe way to reduce your waste management costs and contribute to a cleaner environment.



Sodium dihydrogen orthophosphate dihydrate Certified HPLC for electrochemical detection

packaging		Mfr. No
500 g	Plastic jar	S/3770/53
NaH ₂ PO ₄ ·2H ₂ O	Density at 20°C: ~ 2.0	
CAS: 13472-35-0	Melting point: 60°C	
MW: 156.01	EINECS: 231-449-2	
Assay: 99 to 101%		
Guaranteed analysis		
Electrochemical test	Pass	Iron (Fe) <10 ppm
Reaction (5% solution) at 25°C ...	: 4.3 to 4.5	Lead (Pb) <2 ppm
		n-Butyl formate <0.001%

Absorbance specification	
Max. Absorbance (1 cm cell v HPLC water)	
254 nm (sol. aqu. 1 M) : <0.05 A.U.	
Tested by HPLC-EC to ensure low electrochemical response.	
Detector - ESA Coulochem II	
Mode-oxidative	
Graphite Electrode.	
Lot analysis on the label.	
Filtered to 0.2 µm.	

Sodium dihydrogen orthophosphate dihydrate Certified AR, Eur.Ph., BP for analysis

packaging		Mfr. No
500 g	Plastic jar	S/3760/53
1 kg	Plastic jar	S/3760/60
5 kg	Plastic keg	S/3760/65
25 kg	Plastic keg	S/3760/70
NaH ₂ PO ₄ ·2H ₂ O	Density at 20°C: ~ 2.0	
CAS: 13472-35-0	Melting point: 60°C	
MW: 156.01	EINECS: 231-449-2	

Assay: 99 to 101%			
Guaranteed analysis			
pH (aqueous solution 5%)	4.3 to 4.5	Potassium (K)	<100 ppm
Heavy metals	<10 ppm	Total chloride (Cl)	<0.001%
Calcium (Ca)	<20 ppm	Total silicon (Si)	<50 ppm
Copper (Cu)	<2 ppm	Total sulphur (S)	<20 ppm
Iron (Fe)	<10 ppm	Zinc (Zn)	<2 ppm
Lead (Pb)	<2 ppm	Reducing substance	Pass
Magnesium (Mg)	<10 ppm		

Typical analysis			
Aluminium (Al)	<2 ppm	Molybdenum (Mo)	<1 ppm
Arsenic (As)	<2 ppm	Nickel (Ni)	<1 ppm
Barium (Ba)	<1 ppm	Silver (Ag)	<1 ppm
Cadmium (Cd)	<1 ppm	Strontium (Sr)	<0.001 %
Cobalt (Co)	<1 ppm	Titanium (Ti)	<1 ppm
Lithium (Li)	<1 ppm	Vanadium (V)	<1 ppm
Manganese (Mn)	<1 ppm		

Lot analysis on the label.

Sodium dihydrogen orthophosphate dihydrate SLR extra pure, crystals

packaging		Mfr. No
500 g	Plastic jar	S/3720/53
1 kg	Plastic jar	S/3720/60
5 kg	Plastic keg	S/3720/65
25 kg	Plastic keg	S/3720/70
NaH ₂ PO ₄ ·2H ₂ O	Density at 20°C: ~ 2.0	
CAS: 13472-35-0	Melting point: 60°C	
MW: 156.01	EINECS: 231-449-2	

Assay: >98%			
Guaranteed analysis			
Calcium (Ca)	<100 ppm	Potassium (K)	<500 ppm
Copper (Cu)	<10 ppm	n-Butyl formate	<0.005%
Iron (Fe)	<50 ppm	Total silicon (Si)	<100 ppm
Lead (Pb)	<10 ppm	Total sulphur (S)	<100 ppm
Magnesium (Mg)	<50 ppm	Zinc (Zn)	<10 ppm

Sodium dihydrogen orthophosphate 0.25M solution Certified HPLC for gradient analysis

packaging		Mfr. No
500 mL	Coated Glass bottle	J/7600/PB08
EINECS: 231-449-2		
Absorbance specification		
Max. Absorbance (1cm cell v HPLC water)		
Drift in gradient baseline <40 mA.U.		

Absorbance of the largest eluted peak: <1.5 mA.U.	
Conditions:	
Column: C18 Reverse Phase	
Flow Rate: 3ml/min	
Wavelength: 254nm	
Solvent A: Water + 10%Buffer	
Solvent B: Acetonitrile + 10% Buffer	
Gradient: 10 min. At 5% solvent B followed by a linear gradient to 95% solvent B over 15min.	
Sodium phosphate buffers should be adjusted to the appropriate pH with 1M Orthophosphoric acid	

Sodium diphenylamine-4-sulfonate pure

packaging		Mfr. No
5 g	Glass jar	S/3780/44
C ₁₂ H ₁₀ NNaO ₃ S	GHS: Skin Irrit. 2, Eye Irrit. 2, STOT SE 3	
CAS: 6152-67-6	H319, H315, H335	
MW: 271.26	P261, P302+P352, P280, P305+P351+P338	
Melting point: 200°C		
EINECS: 228-165-6		

Guaranteed analysis
Complies with standard IR spectrum: Pass

Sodium dithionite SLR specified reagents for general laboratory work- extra pure

packaging		Mfr. No
500 g	Plastic jar	S/3800/53
3 kg	Plastic jar	S/3800/63
Na ₂ S ₂ O ₄	EINECS: 231-890-0	
CAS: 7775-14-6	UN: 1384	
MW: 174.10	GHS: Self-heat. 1, Acute Tox. 4	
Density at 20°C: 1.4	H251, H302, EUH031	
Melting point: 300°C	P235+P410, P233, P301+P312	
ADR/RID: 4.2, II		

Assay: >85%

Sodium dodecyl sulfate for electrophoresis

packaging		Mfr. No
100 g	Plastic jar	S/P530/48
500 g	Plastic jar	S/P530/53
25 kg	Fibre Keg	S/P530/70

C ₁₂ H ₂₅ NaO ₄ S	GHS: Flam. Sol. 2, Acute Tox. 3, Skin Irrit. 2, Eye Irrit. 2, STOT SE 3	
CAS: 151-21-3	H228, H311, H302, H335, H315, H319	
Melting point: 206°C	P280, P302+P352, P210, P270, P301+P312, P305+P351+P338	
Flash point: >100°C		
ADR/RID: 4.1, 6.1, III		
EINECS: 205-788-1		
UN: 2926		

This grade of SDS has been purified for use as a denaturant for proteins and also acts as an excellent solubiliser.

Sodium dodecyl sulfate		
Technical		
packaging	Mfr. No	
500 g Plastic jar	S/5200/53	
C ₁₂ H ₂₅ NaO ₄ S	GHS: Flam. Sol. 2, Acute	
CAS: 151-21-3	Tox. 3, Skin Irrit. 2, Eye Irrit.	
MW: 288.38	2, STOT SE 3	
Melting point: 206°C	H228, H311, H302, H335,	
Flash point: >100°C	H315, H319	
ADR/RID: 4.1, 6.1, III	P280, P302+P352, P210,	
EINECS: 205-788-1	P270, P301+P312,	
UN: 2926	P305+P351+P338	
Assay (CG): >97.5%		

Sodium sulfate alkyl contains mainly sodium lauryl sulfate.

Sodium fluoride		Certified AR
for analysis		
packaging	Mfr. No	
250 g Plastic jar	S/3920/50	
500 g Plastic jar	S/3920/53	
1 kg Plastic jar	S/3920/60	
3 kg Plastic jar	S/3920/63	
25 kg Plastic keg	S/3920/70	
NaF	EINECS: 231-667-8	
CAS: 7681-49-4	UN: 1690	
MW: 41.99	GHS: Acute Tox. 3, Skin Irrit.	
Density at 20°C: 2.78	2, Eye Irrit. 2	
Boiling point: 1700°C	H315, H301, H319, EUH032	
Melting point: 993°C	P301+P310, P280,	
ADR/RID: 6.1, III	P302+P350, P305+P351+P338	
Assay: >99%		

Lot analysis on the label.

Sodium fluoride		SLR
specified reagents for general laboratory work- extra pure		
packaging		Mfr. No
250 g	Plastic jar	S/3880/50
500 g	Plastic jar	S/3880/53
NaF		EINECS: 231-667-8
CAS: 7681-49-4		UN: 1690
MW: 41.99		GHS: Acute Tox. 3, Skin Irrit.
Density at 20°C: 2.78		2, Eye Irrit. 2
Boiling point: 1700°C		H315, H301, H319, EUH032
Melting point: 993°C		P301+P310, P280,
ADR/RID: 6.1, III		P302+P350, P305+P351+P338
Assay: >98%		
Guaranteed analysis		
Calcium (Ca)	<200 ppm	n-Butyl formate <0.02%
Copper (Cu)	<20 ppm	Total phosphorus (P) <200 ppm
Iron (Fe)	<50 ppm	Total silicon (Si) <1000 ppm
Lead (Pb)	<50 ppm	Total sulphur (S) <200 ppm
Magnesium (Mg)	<100 ppm	Zinc (Zn) <20 ppm
Potassium (K)	<500 ppm	

Sodium formaldehydesulfoxylate dihydrate		
pure		
packaging	Mfr. No	
500 g Plastic jar	S/4040/53	
25 kg Fibre Keg	S/4040/70	
CH ₃ NaO ₂ S.2H ₂ O	UN: 3288	
CAS: 6035-47-8	GHS: Muta. 2, Acute Tox. 4	
MW: 154.12	H341, EUH031	
Melting point: 63 to 64°C	P281, P308+P313	
ADR/RID: 6.1, III		
Guaranteed analysis		
Complies with standard IR spectrum: Pass		

Sodium formate		Certified
for HPLC		
packaging		Mfr. No
250 g	Plastic jar	S/4086/50
CHNaO ₂	Boiling point: 360°C	
CAS: 141-53-7	Melting point: 258 to 261°C	
MW: 68.01	EINECS: 205-488-0	
Density at 20°C: 1.92		
Assay: >98%		
Absorbance specification		
Max. Absorbance (1cm cell v HPLC water)		
254 nm (sol. aqu. 1 M) : <0.02 A.U.		

Sodium formate		Certified AR	
for analysis			
packaging	Mfr. No		
500 g Plastic jar	S/4082/53		
CHNaO ₂	Boiling point: 360°C		
CAS: 141-53-7	Melting point: 258 to 261°C		
MW: 68.01	EINECS: 205-488-0		
Density at 20°C: 1.92			
Assay: >99%			
Guaranteed analysis			
Aluminium (Al)	<5 ppm	Manganese (Mn)	<1 ppm
Barium (Ba)	<2 ppm	Nickel (Ni)	<1 ppm
Cadmium (Cd)	<1 ppm	Potassium (K)	<100 ppm
Calcium (Ca)	<50 ppm	Propan--1-ol	<1 ppm
Chromium (Cr)	<1 ppm	Titanium (Ti)	<1 ppm
Cobalt (Co)	<1 ppm	n-Butyl formate	<0.005%
Copper (Cu)	<10 ppm	Total silicon (Si)	<50 ppm
Iron (Fe)	<10 ppm	Total sulphur (S)	<50 ppm
Lead (Pb)	<10 ppm	Vanadium (V)	<1 ppm
Lithium (Li)	<1 ppm	Zinc (Zn)	<10 ppm
Magnesium (Mg)	<5 ppm		
Lot analysis on the label.			

Sodium formate		SLR	
specified reagents for general laboratory work- extra pure			
packaging		Mfr. No	
500 g	Plastic jar	S/4080/53	
CHNaO ₂	Boiling point: 360°C		
CAS: 141-53-7	Melting point: 258 to 261°C		
MW: 68.01	EINECS: 205-488-0		
Density at 20°C: 1.92			
Assay: >98%			
Guaranteed analysis			
Calcium (Ca)	<200 ppm	Potassium (K)	<500 ppm
Copper (Cu)	<50 ppm	n-Butyl formate	<0.02%
Iron (Fe)	<50 ppm	Total silicon (Si)	<200 ppm
Lead (Pb)	<50 ppm	Total sulphur (S)	<200 ppm
Magnesium (Mg)	<20 ppm	Zinc (Zn)	<50 ppm

Sodium gluconate			SLR
specified reagents for general laboratory work- extra pure			
packaging			Mfr. No
500 g Plastic jar			S/4100/53
C ₆ H ₁₁ NaO ₇		Melting point: 206 to 209°C	
CAS: 527-07-1		EINECS: 208-407-7	
MW: 218.14			
Assay: >99%			
Guaranteed analysis			
pH (aqueous solution 5%)	5.5 to 7.5	Sulfate (SO ₄)	<0.02%
Iron (Fe)	<10 ppm	n-Butyl formate	<0.005%
Lead (Pb)	<20 ppm		

Sodium β-glycerophosphate pentahydrate		SLR
specified reagents for general laboratory work- extra pure		
packaging		Mfr. No
100 g	Plastic jar	S/4110/48
C ₃ H ₇ Na ₃ O ₈ P ₅ ·5H ₂ O	Melting point: 102 to 104°C	
CAS: 819-83-0	EINECS: 212-464-3	
MW: 315.11		
Assay: >97%		
Guaranteed analysis		
pH (aqueous solution 2%)	8 to 10.5	Lead (Pb) <20 ppm
Inorganic Phosphate	<0.2%	

Sodium hexametaphosphate		SLR
extra pure, powder mainly (NaPO ₃) ₆		
packaging		Mfr. No
500 g	Plastic jar	S/4160/53
1 kg	Plastic jar	S/4160/60
Na ₆ O ₁₈ P ₆	Density at 20°C: 2.181	
CAS: 10124-56-8	Melting point: 640°C	
MW: 611.52	EINECS: 233-243-1	
Guaranteed analysis		
Identity test : Pass		

Sodium hexametaphosphate		
pure, Flake, mainly (NaPO ₃) ₆		
packaging		Mfr. No
500 g	Plastic jar	S/4120/53
1 kg	Plastic jar	S/4120/60
Na ₆ O ₁₈ P ₆	Density at 20°C: 2.181	
CAS: 10124-56-8	Melting point: 640°C	
MW: 611.52	EINECS: 233-243-1	
Guaranteed analysis		
Identity test : Pass		

Sodium hydrogen carbonate		Certified
for HPLC		
packaging	Mfr. No	
250 g Plastic jar	S/4246/50	
NaHCO ₃	Melting point: 270°C	
CAS: 144-55-8	EINECS: 205-633-8	
MW: 84.01		
Assay: >99.5%		
Guaranteed analysis		
Iron (Fe)	<5 ppm	Ammonium (NH ₄) <0.0005%
Lead (Pb)	<2 ppm	pH (aqueous solution 1%) <8.4
Absorbance specification		
Max. Absorbance (1 cm cell v HPLC water)		
254 nm (sol. aqu. 1 M) : <0.03 A.U.		

Sodium hydrogen carbonate		Certified AR, Eur.Ph.	
for analysis			
packaging	Mfr. No		
500 g Plastic jar	S/4240/53		
1 kg Plastic jar	S/4240/60		
3 kg Plastic jar	S/4240/63		
5 kg Plastic keg	S/4240/65		
25 kg Plastic keg	S/4240/70		
NaHCO ₃	Melting point: 270°C		
CAS: 144-55-8	EINECS: 205-633-8		
MW: 84.01			
Assay: 99.7 to 100.3%			
Guaranteed analysis			
pH (aqueous solution 1%)	8 to 8.5	Total chloride (Cl)	<0.01%
Insoluble matter	<0.015%	Total nitrogen (N)	<20 ppm
Calcium (Ca)	<50 ppm	Total phosphorus (P)	<2 ppm
Copper (Cu)	<2 ppm	Total silicon (Si)	<20 ppm
Iron (Fe)	<5 ppm	Total sulphur (S)	<10 ppm
Lead (Pb)	<2 ppm	Zinc (Zn)	<5 ppm
Magnesium (Mg)	<50 ppm	Ammonium (NH ₄)	<0.0005%
Potassium (K)	<50 ppm	Phosphate (PO ₄)	<0.001%
Typical analysis			
Aluminium (Al)	<2 ppm	Manganese (Mn)	<1 ppm
Arsenic (As)	<2 ppm	Molybdenum (Mo)	<1 ppm
Barium (Ba)	<1 ppm	Nickel (Ni)	<1 ppm
Cadmium (Cd)	<1 ppm	Silver (Ag)	<1 ppm
Chromium (Cr)	<1 ppm	Strontium (Sr)	<5 ppm
Cobalt (Co)	<1 ppm	Titanium (Ti)	<1 ppm
Lithium (Li)	<1 ppm	Vanadium (V)	<1 ppm

Lot analysis on the label.

Sodium hydrogen carbonate		SLR
specified reagents for general laboratory work- extra pure		
packaging		Mfr. No
1 kg	Plastic jar	S/4200/60
3 kg	Plastic jar	S/4200/63
5 kg	Plastic keg	S/4200/65
25 kg	Plastic keg	S/4200/70
NaHCO ₃		Melting point: 270°C
CAS: 144-55-8		EINECS: 205-633-8
MW: 84.01		
Assay: >99%		
Guaranteed analysis		
Calcium (Ca)	<200 ppm	n-Butyl formate <0.02%
Copper (Cu)	<10 ppm	Total phosphorus (P) <10 ppm
Iron (Fe)	<20 ppm	Total silicon (Si) <100 ppm
Lead (Pb)	<10 ppm	Total sulphur (S) <50 ppm
Magnesium (Mg)	<200 ppm	Zinc (Zn) <20 ppm
Potassium (K)	<500 ppm	

di-Sodium hydrogen citrate sesquihydrate			SLR	
specified reagents for general laboratory work- extra pure				
packaging			Mfr. No	
500 g Plastic jar			S/4280/53	
C ₆ H ₆ Na ₂ O ₇ ·1.5H ₂ O		Melting point: 149°C		
CAS: 6132-05-4		EINECS: 205-623-3		
MW: 263.11				
Assay: 99 to 101%				
Guaranteed analysis				
pH (solution aq. 3%)		: 4.9 to 5.2	Sulfate (SO ₄)	<0.02%
Lead (Pb)		<10 ppm	n-Butyl formate	<0.01%
Potassium (K)		<200 ppm		

di-Sodium hydrogen orthophosphate anhydrous
for electrophoresis, dibasic

packaging	Mfr. No
500 g Plastic jar	S/P525/53
3 kg Plastic keg	S/P525/63
Na ₂ HPO ₄ CAS: 7558-79-4 MW: 141.96 EINECS: 231-448-7	
Assay: >99%	
Guaranteed analysis	
n-Butyl formate	<0.01%

di-Sodium hydrogen orthophosphate anhydrous
Certified AR, Eur.Ph., BP - for analysis

packaging	Mfr. No
500 g Plastic jar	S/4520/53
1 kg Plastic jar	S/4520/60
3 kg Plastic jar	S/4520/63
5 kg Plastic keg	S/4520/65
25 kg Fibre Keg	S/4520/70
Na ₂ HPO ₄ CAS: 7558-79-4 MW: 141.96 EINECS: 231-448-7	
Assay: 99.5 to 101%	
Guaranteed analysis	
pH (aqueous solution 2%)	9 to 9.3
pH(aq. Solution 5%) at 25°C.	8.7 to 9.3
Loss on drying	<0.2%
Heavy metals	<10 ppm
Insoluble matter	<0.01%
Monosodium phosphate	: <2.5%
Calcium (Ca)	<20 ppm
Copper (Cu)	<5 ppm
Iron (Fe)	<10 ppm
Lead (Pb)	<5 ppm
Typical analysis	
Aluminium (Al)	<2 ppm
Arsenic (As)	<2 ppm
Barium (Ba)	<1 ppm
Cadmium (Cd)	<1 ppm
Cobalt (Co)	<1 ppm
Lithium (Li)	<1 ppm
Manganese (Mn)	<1 ppm
Lot analysis on the label.	

di-Sodium hydrogen orthophosphate anhydrous SLR
specified reagents for general laboratory work- extra pure

packaging	Mfr. No
500 g Plastic jar	S/4480/53
1 kg Plastic jar	S/4480/60
5 kg Plastic keg	S/4480/65
Na ₂ HPO ₄ CAS: 7558-79-4 MW: 141.96 EINECS: 231-448-7	
Assay: >99%	
Guaranteed analysis	
Loss on drying	<0.5%
Calcium (Ca)	<100 ppm
Copper (Cu)	<50 ppm
Iron (Fe)	<50 ppm
Lead (Pb)	<20 ppm
Magnesium (Mg)	<50 ppm
Potassium (K)	<500 ppm
n-Butyl formate	<0.01%
Total silicon (Si)	<200 ppm
Total sulphur (S)	<100 ppm
Zinc (Zn)	<200 ppm

di-Sodium hydrogen orthophosphate dihydrate Certified
HPLC for electrochemical detection

packaging	Mfr. No
500 g Plastic jar	S/4460/53
Na ₂ HPO ₄ ·2H ₂ O CAS: 10028-24-7 MW: 177.99	
Melting point: 243 to 245°C EINECS: 231-448-7	
Assay: >99%	
Guaranteed analysis	
Electrochemical test	Pass
Reaction (5% solution) at 25°C	: 7.5 to 9.2
Iron (Fe)	<10 ppm
Lead (Pb)	<5 ppm
n-Butyl formate	<0.0005%
Absorbance specification	
Max. Absorbance (1cm cell v HPLC water)	
254 nm (sol. aqu. 0.5 M)	
: <0.05 A.U.	
Tested by HPLC-EC to ensure low electrochemical response.	
Detector - ESA Coulochem II	
Mode-oxidative	
Graphite Electrode.	
Lot analysis on the label.	
Filtered to 0.2 µm.	

di-Sodium hctrogen orthophosphate dihydrate Certified AR
for analysis

packaging	Mfr. No
500 g Plastic jar	S/4450/53
1 kg Plastic jar	S/4450/60
5 kg Plastic keg	S/4450/65
10 kg Plastic keg	S/4450/68
25 kg Plastic keg	S/4450/70
Na ₂ HPO ₄ ·2H ₂ O CAS: 10028-24-7 MW: 177.99	
Melting point: 243 to 245°C EINECS: 231-448-7	
Assay: 99 to 101%	
Guaranteed analysis	
pH (aqueous solution 2%)	9 to 9.3
Calcium (Ca)	<20 ppm
Copper (Cu)	<5 ppm
Iron (Fe)	<10 ppm
Lead (Pb)	<5 ppm
Magnesium (Mg)	<10 ppm
Potassium (K)	<100 ppm
Total chloride (Cl)	<0.002%
Total nitrogen (N)	<20 ppm
Total silicon (Si)	<50 ppm
Total sulphur (S)	<20 ppm
Zinc (Zn)	<5 ppm
Typical analysis	
Aluminium (Al)	<2 ppm
Arsenic (As)	<2 ppm
Barium (Ba)	<1 ppm
Cadmium (Cd)	<1 ppm
Cobalt (Co)	<1 ppm
Lithium (Li)	<1 ppm
Manganese (Mn)	<1 ppm
Molybdenum (Mo)	<1 ppm
Nickel (Ni)	<1 ppm
Silver (Ag)	<1 ppm
Strontium (Sr)	<2 ppm
Titanium (Ti)	<1 ppm
Vanadium (V)	<1 ppm
Lot analysis on the label.	

di-Sodium hydrogen orthophosphate dihydrate SLR
specified reagents for general laboratory work- extra pure

packaging	Mfr. No
500 g Plastic jar	S/4440/53
1 kg Plastic jar	S/4440/60
5 kg Plastic keg	S/4440/65
Na ₂ HPO ₄ ·2H ₂ O CAS: 10028-24-7 MW: 177.99	
Melting point: 243 to 245°C EINECS: 231-448-7	
Assay: >98%	
Guaranteed analysis	
Calcium (Ca)	<100 ppm
Copper (Cu)	<20 ppm
Iron (Fe)	<50 ppm
Lead (Pb)	<20 ppm
Magnesium (Mg)	<50 ppm
Potassium (K)	<500 ppm
n-Butyl formate	<0.01%
Total nitrogen (N)	<20 ppm
Total silicon (Si)	<200 ppm
Total sulphur (S)	<100 ppm
Zinc (Zn)	<20 ppm

di-Sodium hydrogen orthophosphate dodecahydrate
Certified AR, Eur.Ph., BP
for analysis

packaging	Mfr. No
500 g Plastic jar	S/4400/53
1 kg Plastic jar	S/4400/60
10 kg Plastic keg	S/4400/68
Na ₂ HPO ₄ ·12H ₂ O CAS: 10039-32-4 MW: 358.14	
Melting point: 35°C EINECS: 231-448-7	
Assay: 99 to 102%	
Guaranteed analysis	
pH (aqueous solution 5%)	9.1 to 9.4
Calcium (Ca)	<20 ppm
Copper (Cu)	<2 ppm
Iron (Fe)	<5 ppm
Lead (Pb)	<2 ppm
Magnesium (Mg)	<5 ppm
Potassium (K)	<50 ppm
Total chloride (Cl)	<0.002%
Total silicon (Si)	<50 ppm
Total sulphur (S)	<20 ppm
Zinc (Zn)	<5 ppm
Typical analysis	
Aluminium (Al)	<2 ppm
Arsenic (As)	<2 ppm
Barium (Ba)	<1 ppm
Cadmium (Cd)	<1 ppm
Cobalt (Co)	<1 ppm
Lithium (Li)	<1 ppm
Manganese (Mn)	<1 ppm
Molybdenum (Mo)	<1 ppm
Nickel (Ni)	<1 ppm
Silver (Ag)	<1 ppm
Strontium (Sr)	<2 ppm
Titanium (Ti)	<1 ppm
Vanadium (V)	<1 ppm
Lot analysis on the label.	

di-Sodium hydrogen orthophosphate dodecahydrate SLR
specified reagents for general laboratory work- extra pure

packaging	Mfr. No
500 g Plastic jar	S/4360/53
1 kg Plastic jar	S/4360/60
2.5 kg Plastic jar	S/4360/62
Na ₂ HPO ₄ ·12H ₂ O CAS: 10039-32-4 MW: 358.14	
Melting point: 35°C EINECS: 231-448-7	
Assay: >99%	
Guaranteed analysis	
Calcium (Ca)	<100 ppm
Copper (Cu)	<20 ppm
Iron (Fe)	<50 ppm
Lead (Pb)	<20 ppm
Magnesium (Mg)	<50 ppm
Potassium (K)	<500 ppm
n-Butyl formate	<0.01%
Total silicon (Si)	<200 ppm
Total sulphur (S)	<100 ppm
Zinc (Zn)	<20 ppm

di-Sodium hydrogen orthophosphate 0.25M solution Certified
for HPLC

packaging	Mfr. No
500 mL Coated Glass bottle	J/7650/PB08
Na ₂ HPO ₆ MW: 141.98 EINECS: 231-448-7	
Absorbance specification	
Max. Absorbance (1cm cell v HPLC water)	
Drift in gradient baseline	
<40 mA.U.	
Absorbance of the largest eluted peak: <1.5 mA.U.	

Optima LC/MS

Fisher Chemical OPTIMA® LC/MS solvents, blends, additives and modifiers have set the standard of excellence for consistent, reproducible performance in the mobile phase of LC/MS.

Find out more on page 308

Sodium hydrogen sulfate monohydrate Certified AR
for analysis

packaging	Mfr. No
500 g Plastic jar	S/4690/53
NaHSO ₄ ·H ₂ O CAS: 10034-88-5 MW: 138.08 Density at 20°C: 2.12 Melting point: 58°C ADR/RID: 8, III	
EINECS: 231-665-7 UN: 3260 GHS: Eye Dam. 1 H318 P305+P351+P338, P310, P280	
Assay: >99%	
Guaranteed analysis	
Aluminium (Al)	<5 ppm
Barium (Ba)	<5 ppm
Cadmium (Cd)	<5 ppm
Calcium (Ca)	<20 ppm
Chromium (Cr)	<5 ppm
Cobalt (Co)	<5 ppm
Copper (Cu)	<5 ppm
Iron (Fe)	<5 ppm
Lead (Pb)	<5 ppm
Lithium (Li)	<5 ppm
Magnesium (Mg)	<10 ppm
Manganese (Mn)	<5 ppm
Molybdenum (Mo)	<5 ppm
Nickel (Ni)	<5 ppm
Potassium (K)	<20 ppm
Propan--1-ol	<5 ppm
Strontium (Sr)	<5 ppm
n-Butyl formate	<0.001%
Total nitrogen (N)	<10 ppm
Total phosphorus (P)	<5 ppm
Total silicon (Si)	<50 ppm
Zinc (Zn)	<10 ppm
Lot analysis on the label.	

Sodium hydrogen tartrate monohydrate SLR
specified reagents for general laboratory work- extra pure

packaging	Mfr. No
500 g Plastic jar	S/4720/53
NaHCH ₄ O ₆ ·H ₂ O CAS: 6131-98-2 MW: 190.08 EINECS: 208-400-9	
Assay: > 98%	
Guaranteed analysis	
pH (aqueous solution 5%)	
3.3 - 3.6	

Sodium hydrogen tartrate decahydrate Certified AR
for analysis

packaging	Mfr. No
500 g Plastic jar	S/4760/53
NaHCH ₄ O ₆ ·H ₂ O CAS: 6131-98-2 MW: 190.08	
Assay: >99%	
Guaranteed analysis	
pH (aqueous solution 5%)	3.3 to 3.6
Loss on drying	<10%
Calcium (Ca)	<50 ppm
Copper (Cu)	<5 ppm
Iron (Fe)	<10 ppm
Lead (Pb)	<5 ppm
Magnesium (Mg)	<20 ppm
Potassium (K)	<100 ppm
n-Butyl formate	<0.001%
Total phosphorus (P)	<50 ppm
Total silicon (Si)	<20 ppm
Total sulphur (S)	<50 ppm
Zinc (Zn)	<10 ppm
Typical analysis	
Aluminium (Al)	<5 ppm
Barium (Ba)	<0.5 ppm
Cadmium (Cd)	<2 ppm
Chromium (Cr)	<2 ppm
Cobalt (Co)	<2 ppm
Lithium (Li)	<2 ppm
Manganese (Mn)	<2 ppm
Molybdenum (Mo)	<2 ppm
Nickel (Ni)	<2 ppm
Silver (Ag)	<2 ppm
Strontium (Sr)	<2 ppm
Titanium (Ti)	<2 ppm
Vanadium (V)	<2 ppm
Lot analysis on the label	

Sodium hydroxide pellets
for analysis

Certified AR, Eur.Ph., BP

packaging		Mfr. No	
500 g	Plastic jar	S/4920/53	
1 kg	Plastic jar	S/4920/60	
3 kg	Plastic jar	S/4920/63	
5 kg	Plastic jar	S/4920/65	
10 kg	Plastic keg	S/4920/68	
25 kg	Plastic keg	S/4920/70	
NaOH		EINECS: 215-185-5	
CAS: 1310-73-2		UN: 1823	
MW: 40.00		GHS: Skin Corr. 1A, Eye Dam. 1	
Boiling point: 1390°C		H314	
Melting point: 318°C		P280, P301+P330+P331,	
ADR/RID: 8, II		P305+P351+P338, P302+P352	
Assay: >98%			
Guaranteed analysis			
Heavy metals	<20 ppm	Total nitrogen (N)	<10 ppm
Calcium (Ca)	<20 ppm	Total phosphorus (P)	<5 ppm
Carbonate (Na ₂ CO ₃)	<1%	Total silicon (Si)	<50 ppm
Copper (Cu)	<5 ppm	Total sulphur (S)	<20 ppm
Iron (Fe)	<5 ppm	Zinc (Zn)	<5 ppm
Lead (Pb)	<5 ppm	Mercury (Hg)	<0.1 ppm
Magnesium (Mg)	<5 ppm	Phosphate (PO ₄)	<0.001%
Potassium (K)	<200 ppm	Sulfate (SO ₄)	<0.003%
Total chloride (Cl)	<0.005%		
Typical analysis			
Aluminium (Al)	<2 ppm	Manganese (Mn)	<1 ppm
Arsenic (As)	<2 ppm	Molybdenum (Mo)	<1 ppm
Barium (Ba)	<1 ppm	Nickel (Ni)	<10 ppm
Cadmium (Cd)	<1 ppm	Silver (Ag)	<1 ppm
Chromium (Cr)	<1 ppm	Strontium (Sr)	<10 ppm
Cobalt (Co)	<1 ppm	Titanium (Ti)	<1 ppm
Lithium (Li)	<1 ppm	Vanadium (V)	<1 ppm

Lot analysis on the label.

Sodium hydroxide pellets
specified reagents for general laboratory work- extra pure

SLR

Packaging		Mfr. No	
250 g	Plastic jar	S/4880/50	
500 g	Plastic jar	S/4880/53	
1 kg	Plastic jar	S/4880/60	
3 kg	Plastic jar	S/4880/63	
5 kg	Plastic jar	S/4880/65	
10 kg	Plastic keg	S/4880/68	
NaOH			
CAS: 1310-73-2			
UN: 1823			
GHS: Skin Corr. 1A, Eye Dam. 1			
H314			
Boiling point: 1390°C		P280, P301+P330+P331,	
Melting point: 318°C		P305+P351+P338, P302+P352	
ADR/RID: 8, II			
Assay: >97%			
Guaranteed analysis			
Calcium (Ca)	<100 ppm	Potassium (K)	<2000 ppm
Carbonate (Na ₂ CO ₃)	<2%	n-Butyl formate	<0.02%
Copper (Cu)	<20 ppm	Total phosphorus (P)	<50 ppm
Iron (Fe)	<20 ppm	Total silicon (Si)	<500 ppm
Lead (Pb)	<20 ppm	Total sulphur (S)	<200 ppm
Magnesium (Mg)	<100 ppm	Zinc (Zn)	<500 ppm

High Purity Acids for Elemental Analysis

Fisher Chemical offers a complete range of acids to meet the most challenging applications. Choose from our Optima®, TraceMetal® and PrimarPlus® grades to ensure that the sensitivity you require is always available.

Read more page 24



Sodium hydroxide pearl
for analysis

Certified AR

Packaging			Mfr. No
1 kg	Plastic jar		S/4845/60
25 kg	Plastic keg		S/4845/70
NaOH		EINECS: 215-185-5	
CAS: 1310-73-2		UN: 1823	
MW: 40.00		GHS: Skin Corr. 1A, Eye Dam. 1	
Boiling point: 1390°C		H314	
Melting point: 318°C		P280, P301+P330+P331,	
ADR/RID: 8, II		P305+P351+P338, P302+P352	
Assay: >98%			
Guaranteed analysis			
Calcium (Ca)	<20 ppm	Total chloride (Cl)	<0.005%
Carbonate (Na ₂ CO ₃)	<1%	Total nitrogen (N)	<5 ppm
Copper (Cu)	<5 ppm	Total phosphorus (P)	<10 ppm
Iron (Fe)	<5 ppm	Total silicon (Si)	<100 ppm
Lead (Pb)	<5 ppm	Total sulphur (S)	<50 ppm
Magnesium (Mg)	<20 ppm	Zinc (Zn)	<10 ppm
Potassium (K)	<500 ppm		
Typical analysis			
Aluminium (Al)	<2 ppm	Manganese (Mn)	<2 ppm
Arsenic (As)	<2 ppm	Molybdenum (Mo)	<2 ppm
Barium (Ba)	<2 ppm	Nickel (Ni)	<2 ppm
Cadmium (Cd)	<2 ppm	Silver (Ag)	<2 ppm
Chromium (Cr)	<2 ppm	Strontium (Sr)	<10 ppm
Cobalt (Co)	<2 ppm	Titanium (Ti)	<2 ppm
Lithium (Li)	<2 ppm	Vanadium (V)	<2 ppm

Sodium hydroxide pearl
specified reagents for general laboratory work- extra pure

SLR

Packaging			Mfr. No
1 kg	Plastic jar		S/4840/60
3 kg	Plastic jar		S/4840/63
NaOH		EINECS: 215-185-5	
CAS: 1310-73-2		UN: 1823	
MW: 40.00		GHS: Skin Corr. 1A, Eye Dam. 1	
Boiling point: 1390°C		H314	
Melting point: 318°C		P280, P301+P330+P331,	
ADR/RID: 8, II		P305+P351+P338, P302+P352	
Assay: >97%			
Guaranteed analysis			
Calcium (Ca)	<100 ppm	Potassium (K)	<2000 ppm
Carbonate (Na ₂ CO ₃)	<2%	n-Butyl formate	<0.02%
Copper (Cu)	<20 ppm	Total phosphorus (P)	<50 ppm
Iron (Fe)	<20 ppm	Total silicon (Si)	<500 ppm
Lead (Pb)	<20 ppm	Total sulphur (S)	<200 ppm
Magnesium (Mg)	<100 ppm	Zinc (Zn)	<50 ppm

Sodium hydroxide pearl
Technical

packaging		Mfr. No	
1 kg	Plastic jar	S/4800/60	
3 kg	Plastic jar	S/4800/63	
5 kg	Plastic jar	S/4800/65	
10 kg	Plastic keg	S/4800/68	
25 kg	Plastic keg	S/4800/70	
NaOH	EINECS: 215-185-5		
CAS: 1310-73-2	UN: 1823		
MW: 40.00	GHS: Skin Corr. 1A, Eye Dam. 1		
Boiling point: 1390°C	H314		
Melting point: 318°C	P280, P301+P330+P331,		
ADR/RID: 8, II	P305+P351+P338, P302+P352		

Sodium hydroxide 46-48% solution
HPLC for electrochemical detection

Certified

packaging		Mfr. No	
2.5 ℓ	HDPE Plastic bottle	S/4940/17	
NaOH	EINECS: 215-185-5		
CAS: 1310-73-2	UN: 1824		
MW: 40.00	GHS: Skin Corr. 1A, Eye Dam. 1		
Density at 20°C: 1.5	H314		
Boiling point: 145°C	P280, P301+P330+P331,		
Melting point: >12°C	P310, P305+P351+P338,		
ADR/RID: 8, II	P303+P361+P353		
Assay: 46 to 51%			
Guaranteed analysis			
Electrochemical test	Pass	Lead (Pb)	<5 ppm
Iron (Fe)	<5 ppm	n-Butyl formate	<0.01%
Absorbance specification			
Max. Absorbance (1cm cell v HPLC water)			
254 nm (sol. aqu. 1 M) : <0.05 A.U.			
Purified then tested by HPLC-EC to show low electrochemical response.			
Detector - ESA Coulochem II			
Mode-oxydative			
Graphite Electrode.			
The concentration of a solution of sodium hydroxide 46/48% is 17.6 M (17.6 N).			
Lot analysis on the label.			
Filtered to 0.2 µm.			

Sodium hydroxide 46-48% solution d=1.50
for analysis

Certified AR

packaging		Mfr. No	
100 mℓ	HDPE Plastic bottle	S/4930/05	
1 ℓ	HDPE Plastic bottle	S/4930/15	
2.5 ℓ	HDPE Plastic bottle	S/4930/17	
5 ℓ	HDPE Plastic bottle	S/4930/21	
25 ℓ	HDPE Plastic container	S/4930/25	
200 ℓ	HDPE Plastic container	S/4930/27	
NaOH	EINECS: 215-185-5		
CAS: 1310-73-2	UN: 1824		
MW: 40.00	GHS: Skin Corr. 1A, Eye Dam. 1		
Density at 20°C: 1.5	H314		
Boiling point: 145°C	P280, P301+P330+P331,		
Melting point: >12°C	P310, P305+P351+P338,		
ADR/RID: 8, II	P303+P361+P353		
Assay: 46 to 51%			
Guaranteed analysis			
Calcium (Ca)	<20 ppm	Total chloride (Cl)	<0.01%
Carbonate (Na ₂ CO ₃)	<0.1%	Total nitrogen (N)	<10 ppm
Copper (Cu)	<2 ppm	Total phosphorus (P)	<2 ppm
Iron (Fe)	<5 ppm	Total silicon (Si)	<20 ppm
Lead (Pb)	<5 ppm	Total sulphur (S)	<20 ppm
Magnesium (Mg)	<5 ppm	Zinc (Zn)	<5 ppm
Potassium (K)	<200 ppm		
Typical analysis			
Aluminium (Al)	<2 ppm	Manganese (Mn)	<1 ppm
Arsenic (As)	<2 ppm	Molybdenum (Mo)	<1 ppm
Barium (Ba)	<1 ppm	Nickel (Ni)	<1 ppm
Cadmium (Cd)	<1 ppm	Silver (Ag)	<1 ppm
Chromium (Cr)	<1 ppm	Strontium (Sr)	<10 ppm
Cobalt (Co)	<1 ppm	Titanium (Ti)	<1 ppm
Lithium (Li)	<1 ppm	Vanadium (V)	<1 ppm

Advantages:

(I) The use of a solution ready to use saving the time required for preparation and to limit the risks arising from manipulation of solid sodium hydroxide.

(II) An extremely low carbonate content for the preparation of carbonate-free sodium hydroxide solutions.

(III) Low nitrogen content of the solution permitting use in the determination of nitrogen by Kjeldahl method.

The concentration of a solution of sodium hydroxide 46/48% is 17.6 M (17.6 N).

Lot analysis on the label.

Sodium hydroxyde solution 32%
for analysis, for Nitrogen determination

NEW

Certified AR

Packaging		Mfr. No	
1 ℓ	HDPE Plastic bottle	S/4955/PB15	
2.5 ℓ	HDPE Plastic bottle	S/4955/PB17	
NaOH	ADR/RID: 8, II		
CAS: 1310-73-2	EINECS: 215-185-5		
MW: 40.00	UN: 1824		
Density at 20°C: 1.33 to1.36	H314		
Boiling point: 120°C	P301+P330+P331, P280,		
Melting point: 8°C	P305+P351+P338, P310		
Assay: 32 to 33%			
Guaranteed analysis			
Heavy metals	5 ppm	Iron (Fe)	5 ppm
Silicate	: 10 ppm	Phosphate (PO ₄)	5 ppm
Aluminium (Al)	5 ppm	Sulfate (SO ₄)	5 ppm
Carbonate (Na ₂ CO ₃)	<=0.2 %	Total nitrogen (N)	1 ppm
Chloride (Cl)	5 ppm		
Not analysis on the label			

Lot analysis on the label

Sodium hydroxide 30% solution
for analysis

Certified AR

Packaging		Mfr. No	
1 ℓ	HDPE Plastic bottle	S/4950/PB15	
2.5 ℓ	HDPE Plastic bottle	S/4950/PB17	
NaOH		EINECS: 215-185-5	
CAS: 1310-73-2		UN: 1824	
MW: 40.00		GHS: Skin Corr. 1A, Eye Dam. 1	
Density at 20°C: 1.33 to1.36		H314	
Boiling point: 120°C		P301+P330+P331, P280,	
Melting point: 8°C		P305+P351+P338, P310	
ADR/RID: 8, II			
Assay: >30%			
Guaranteed analysis			
Calcium (Ca)	<10 ppm	Potassium (K)	<200 ppm
Carbonate (Na ₂ CO ₃)	<0.1%	Total chloride (Cl)	<0.1%
Copper (Cu)	<2 ppm	Total nitrogen (N)	<10 ppm
Iron (Fe)	<5 ppm	Total phosphorus (P)	<2 ppm
Lead (Pb)	<1 ppm	Total silicon (Si)	<20 ppm
Magnesium (Mg)	<5 ppm	Total sulphur (S)	<10 ppm
Manganese (Mn)	<1 ppm	Zinc (Zn)	<5 ppm
Nickel (Ni)	<1 ppm		

Sodium hydroxide 30% solution
specified reagents for general laboratory work- extra pure

SLR

packaging		Mfr. No	
1 ℓ	HDPE Plastic bottle	S/4945/PB15	
5 ℓ	HDPE Plastic bottle	S/4945/PB21	
20 ℓ	HDPE Plastic container	S/4945/31	
NaOH	EINECS: 215-185-5		
CAS: 1310-73-2	UN: 1824		
MW: 40.00	GHS: Skin Corr. 1A, Eye Dam. 1		
Density at 20°C: 1.33 to1.36	H314		
Boiling point: 120°C	P301+P330+P331, P280,		
Melting point: 8°C	P305+P351+P338, P310		
ADR/RID: 8, II			
Assay: >30%			
Guaranteed analysis			
Calcium (Ca)	<10 ppm	Potassium (K)	<200 ppm
Carbonate (Na ₂ CO ₃)	<0.1%	n-Butyl formate	<0.1%
Copper (Cu)	<2 ppm	Total nitrogen (N)	<10 ppm
Iron (Fe)	<5 ppm	Total phosphorus (P)	<2 ppm
Lead (Pb)	<1 ppm	Total silicon (Si)	<20 ppm
Magnesium (Mg)	<5 ppm	Total sulphur (S)	<10 ppm
Manganese (Mn)	<1 ppm	Zinc (Zn)	<5 ppm
Nickel (Ni)	<1 ppm		

Sodium hydroxide solution 2 % v/v
pure

packaging	Mfr. No
500 mℓ HDPE Plastic bottle	J/9013/08
Assay: 1.9 to 2.1%	

Sodium hydroxide solution 2M
pure, Bench Reagent

packaging	Mfr. No
2.5 ℓ HDPE Plastic bottle	J/7810/17
25 ℓ HDPE Plastic container	J/7810/25

ADR/RID: 8, II	H314
EINECS: 215-185-5	P301+P330+P331, P280,
UN: 1824	P305+P351+P338, P310
GHS: Skin Corr. 1A, Eye Dam. 1	

Sodium hydroxide solution 10M (10N) ca. 40w/v
Standard solution ready to use, For volumetric analysis

packaging	Mfr. No
1 ℓ HDPE Plastic bottle	J/7800/15
2.5 ℓ HDPE Plastic bottle	J/7800/17
5 ℓ HDPE Plastic bottle	J/7800/21
10 ℓ Bag in box	J/7800/24
25 ℓ HDPE Plastic container	J/7800/25

Density at 20°C: 1.33 to 1.36	UN: 1824
Boiling point: 120°C	GHS: Skin Corr. 1A, Eye Dam. 1
Melting point: 8°C	H314
ADR/RID: 8, II	P301+P330+P331, P280,
EINECS: 215-185-5	P305+P351+P338, P310

Assay: 39.8 to 40.2%

Manufactured to a tolerance of ± 0.5%.
See the full range of Standard Volumetric solution page 323.

Sodium hydroxide solution 2M (2N) BP

packaging	Mfr. No
1 ℓ HDPE Plastic bottle	J/7610/PB15
2.5 ℓ HDPE Plastic bottle	J/7610/17
5 ℓ HDPE Plastic bottle	J/7610/21

Density at 20°C: 1.36	UN: 1824
Boiling point: 120°C	GHS: Skin Corr. 1A, Eye Dam. 1
Melting point: 8°C	H314
ADR/RID: 8, II	P301+P330+P331, P280,
EINECS: 215-185-5	P305+P351+P338, P310

Manufactured to a tolerance of ± 0.1%.
See the full range of Standard Volumetric solution page 323.

Sodium hydroxide solution 1M (1N) Eur.Ph. , USP, BP

packaging	Mfr. No
1 ℓ HDPE Plastic bottle	J/7620/15
2.5 ℓ HDPE Plastic bottle	J/7620/17
5 ℓ HDPE Plastic bottle	J/7620/21
10 ℓ Bag in box	J/7620/24

Density at 20°C: 1.040	UN: 1824
Boiling point: 120°C	GHS: Skin Corr. 1B, Eye Dam. 1
Melting point: 8°C	H314
ADR/RID: 8, II	P280, P305+P351+P338,
EINECS: 215-185-5	P301+P330+P331, P310

Manufactured to a tolerance of ± 0.1%.
See the full range of Standard Volumetric solution page 323.

Sodium hydroxide solution 1M (1N)
NIST Standard Concentrate, For Volumetric analysis

packaging	Mfr. No
1 Ampoule	J/7620C/05
6 Ampoule	J/7620C/90

Boiling point: 120°C	GHS: Skin Corr. 1A, Eye Dam. 1
ADR/RID: 8, II	H314
EINECS: 215-185-5	P301+P330+P331, P280,
UN: 1824	P305+P351+P338, P310

Each ampoule makes 1 ℓ of solution.
See the full range of Standard Volumetric solution page 323.

Sodium hydroxide solution methanolic 1M (1N)
NIST Standard solution ready to use, For Volumetric analysis

packaging	Mfr. No
1 ℓ HDPE Plastic bottle	J/7625/PB15
2.5 ℓ HDPE Plastic bottle	J/7625/17

Density at 20°C: 1.040	Tox. 3, Skin Corr. 1B, Eye
Boiling point: 120°C	Dam. 1, STOT SE 1
Melting point: 8°C	H225, H331, H311, H301,
ADR/RID: 3, 6.1. 8, II	H370, H314
EINECS: 215-185-5	P280, P260, P264, P270,
UN: 3286	P210, P307+P311
GHS: Flam. Liq. 2, Acute	

Manufactured to a tolerance of ± 0.2%.
See the full range of Standard Volumetric solution page 323.

Sodium hydroxide solution 0.5M (0.5N)
NIST Standard solution ready to use, For Volumetric analysis

packaging	Mfr. No
1 ℓ HDPE Plastic bottle	J/7630/15
2.5 ℓ HDPE Plastic bottle	J/7630/17
5 ℓ HDPE Plastic bottle	J/7630/21
10 ℓ Bag in box	J/7630/24
200 ℓ HDPE Plastic container	J/7630/27

Density at 20°C: 1.020	GHS: Skin Corr. 1A, Eye Dam. 1
Boiling point: 100°C	H314
ADR/RID: 8, II	P280, P305+P351+P338,
EINECS: 215-185-5	P301+P330+P331, P310
UN: 1824	

Manufactured to a tolerance of ± 0.1%.
See the full range of Standard Volumetric solution page 323.

Sodium hydroxide solution 0.5M (0.5N)
NIST Standard Concentrate, For Volumetric analysis

packaging	Mfr. No
1 Ampoule	J/7630C/05
6 Ampoule	J/7630C/90

Boiling point: 100°C	H314
ADR/RID: 8, II	P280, P302+P350,
EINECS: 215-185-5	P305+P351+P338,
UN: 1824	P301+P330+P331, P315
GHS: Skin Corr. 1A, Eye Dam. 1	

Each ampoule makes 1 ℓ of solution.
See the full range of Standard Volumetric solution page 323.

Sodium hydroxide solution 0.2M (0.2N)
NIST Standard solution ready to use, For Volumetric analysis

packaging	Mfr. No
1 ℓ HDPE Plastic bottle	J/7640/15
2.5 ℓ HDPE Plastic bottle	J/7640/17
5 ℓ HDPE Plastic bottle	J/7640/21
10 ℓ Bag in box	J/7640/24

Density at 20°C: 1.010	GHS: Skin Irrit. 2, Eye Irrit. 2
ADR/RID: 8, III	H315, H319
EINECS: 215-185-5	P280, P305+P351+P338,
UN: 1824	P302+P350

Manufactured to a tolerance of ± 0.1%.
See the full range of Standard Volumetric solution page 323.

Sodium hydroxide solution 0.1M (0.1N) Eur.Ph. , BP

packaging	Mfr. No
1 ℓ HDPE Plastic bottle	J/7660/15
2.5 ℓ HDPE Plastic bottle	J/7660/17
5 ℓ HDPE Plastic bottle	J/7660/21
10 ℓ Bag in box	J/7660/24

Density at 20°C: 1.000	EINECS: 215-185-5
Boiling point: 100°C	UN: 1824
ADR/RID: 8/III	

Manufactured to a tolerance of ± 0.1%.
See the full range of Standard Volumetric solution page 323.

Sodium hydroxide solution 0.1M (0.1N)
NIST Standard Concentrate, For Volumetric analysis

packaging	Mfr. No
1 Ampoule	J/7660C/05
6 Ampoule	J/7660C/90

Density at 20°C: 1.040	GHS: Skin Corr. 1B, Eye Dam. 1
ADR/RID: 8, II	H314
EINECS: 215-185-5	P280, P305+P351+P338,
UN: 1824	P301+P330+P331, P310

Each ampoule makes 1 ℓ of solution.
See the full range of Standard Volumetric solution page 323.

Sodium hydroxide solution 0.01M (0.01N)
NIST Standard solution ready to use, For volumetric analysis

packaging	Mfr. No
2.5 ℓ HDPE Plastic bottle	J/9010/17

Manufactured to a tolerance of ± 0.1%.
See the full range of Standard Volumetric solution page...

Sodium hypochlorite solution
Technical, min. 14% free chlorine

packaging	Mfr. No
1 ℓ HDPE Plastic bottle	S/5042/15

NaOCl	UN: 1791
CAS: 7681-52-9	GHS: Skin Corr. 1B, Aquatic
MW: 74.44	Acute 1
Density at 20°C: 1.200	H314, H400, EUH031
Melting point: -16°C	P301+P330+P331, P280,
ADR/RID: 8, III	P305+P351+P338, P310,
EINECS: 231-668-3	P273, P233

Guaranteed analysis
Chlorine >14%

Sodium hypochlorite solution self indicating
Technical

packaging	Mfr. No
5 ℓ HDPE Plastic bottle	S/5041/21

NaClO	EINECS: 231-668-3
CAS: 7681-52-9	UN: 1791
MW: 74.44	GHS: Skin Corr. 1B, Eye Dam. 1
Density at 20°C: 1.1	H314, EUH031
Melting point: -16°C	P280, P305+P351+P338,
ADR/RID: 8, III	P310, P273, P233, P302+P352

Guaranteed analysis
Chlorine >8%

Contains about 20 ppm Potassium permanganate.

Sodium hypochlorite solution
Technical

packaging	Mfr. No
1 ℓ HDPE Plastic bottle	S/5040/PB15
2.5 ℓ HDPE Plastic bottle	S/5040/PB17
5 ℓ HDPE Plastic bottle	S/5040/21

NaClO	UN: 1791
CAS: 7681-52-9	GHS: Skin Corr. 1B, Aquatic
MW: 74.44	Acute 1
Density at 20°C: 1.200	H314, H400, EUH031
Melting point: -16°C	P301+P330+P331, P280,
ADR/RID: 8, III	P305+P351+P338, P310,
EINECS: 231-668-3	P273, P233

Guaranteed analysis
Alkalinity <2% Chlorine >14%

Sodium hypophosphite monohydrate Certified AR
for analysis

packaging	Mfr. No
500 g Plastic jar	S/5081/53

NaH ₂ PO ₂ ·H ₂ O	Melting point: 90°C
CAS: 10039-56-2	EINECS: 216-948-5
MW: 105.99	

Assay: >99%

Guaranteed analysis	
pH (aqueous solution 5%)	5 to 7
Aluminium (Al)	<2 ppm
Ammonia (NH ₃)	<0.0005%
Arsenic (As)	<2 ppm
Barium (Ba)	<2 ppm
Cadmium (Cd)	<2 ppm
Calcium (Ca)	<100 ppm
Chromium (Cr)	<2 ppm
Cobalt (Co)	<2 ppm
Copper (Cu)	<5 ppm
Iron (Fe)	<5 ppm
Lead (Pb)	<5 ppm
Lithium (Li)	<2 ppm
Magnesium (Mg)	<5 ppm
Manganese (Mn)	<2 ppm
Molybdenum (Mo)	<2 ppm
Nickel (Ni)	<2 ppm
Potassium (K)	<100 ppm
Propan--1-ol	<2 ppm
Strontium (Sr)	<2 ppm
Titanium (Ti)	<2 ppm
n-Butyl formate	<0.05%
Total silicon (Si)	<100 ppm
Total sulphur (S)	<200 ppm
Vanadium (V)	<2 ppm
Zinc (Zn)	<10 ppm

Lot analysis on the label.

Sodium iodate SLR
specified reagents for general laboratory work- extra pure

packaging	Mfr. No
100 g Plastic jar	S/5100/48

NaIO ₃	GHS: Ox. Sol. 2, Acute Tox.
CAS: 7681-55-2	4, Skin Irrit. 2, Eye Irrit. 2,
MW: 197.89	STOT SE 3
Density at 20°C: 4.28	H272, H335, H302, H315, H319
ADR/RID: 5.1, II	P210, P261, P301+P312,
EINECS: 231-672-5	P302+P352, P280,
UN: 1479	P305+P351+P338

Assay: >99.5%

Guaranteed analysis
Chloride and chlorate <0.04% Sulfate (SO₄) <0.03%
Iodide (I) <0.002%

Sodium iodide		Certified AR
for analysis		
packaging		Mfr. No
100 g	Plastic jar	S/5130/48
500 g	Plastic jar	S/5130/53
10 kg	Plastic keg	S/5130/68
NaI		Boiling point: 1300°C
CAS: 7681-82-5		Melting point: 661°C
MW: 149.89		EINECS: 231-679-3
Density at 20°C: 3.66		
Assay: >99.5%		

Guaranteed analysis			
Acidity / Alkalinity	<0.002 meq/g	Lead (Pb)	<2 ppm
Loss on drying	<0.5%	Magnesium (Mg)	<5 ppm
Iodine reducing substances	: <0.005%	Potassium (K)	<100 ppm
Calcium (Ca)	<10 ppm	Total nitrogen (N)	<10 ppm
Chloride and Bromide	<0.01%	Total phosphorus (P)	<10 ppm
Copper (Cu)	<2 ppm	Total silicon (Si)	<20 ppm
Iodate (IO ₃)	<2 ppm	Total sulphur (S)	<100 ppm
Iron (Fe)	<2 ppm	Zinc (Zn)	<10 ppm

Typical analysis			
Aluminium (Al)	<2 ppm	Manganese (Mn)	<2 ppm
Arsenic (As)	<0.1 ppm	Molybdenum (Mo)	<2 ppm
Barium (Ba)	<2 ppm	Nickel (Ni)	<2 ppm
Cadmium (Cd)	<2 ppm	Silver (Ag)	<2 ppm
Chromium (Cr)	<2 ppm	Strontium (Sr)	<2 ppm
Cobalt (Co)	<2 ppm	Titanium (Ti)	<2 ppm
Lithium (Li)	<2 ppm	Vanadium (V)	<2 ppm

Lot analysis on the label.

Sodium iodide
specified reagents for general laboratory work- extra pure

packaging		Mfr. No
100 g	Plastic jar	S/5120/48
250 g	Plastic jar	S/5120/50
500 g	Plastic jar	S/5120/53
NaI		Boiling point: 1300°C
CAS: 7681-82-5		Melting point: 661°C
MW: 149.89		EINECS: 231-679-3
Density at 20°C: 3.66		
Assay: >99%		

Guaranteed analysis			
Loss on drying	<1%	Magnesium (Mg)	<20 ppm
Calcium (Ca)	<50 ppm	Potassium (K)	<500 ppm
Chloride and Bromide	<0.05%	Total phosphorus (P)	<50 ppm
Copper (Cu)	<20 ppm	Total silicon (Si)	<100 ppm
Iodate (IO ₃)	<10 ppm	Total sulphur (S)	<100 ppm
Iron (Fe)	<20 ppm	Zinc (Zn)	<50 ppm
Lead (Pb)	<20 ppm		

Sodium lactate solution about 60% w/w
specified reagents for general laboratory work- extra pure

packaging	Mfr. No	
500 mL Glass bottle	S/5161/08	
C ₃ H ₅ NaO ₃	Boiling point: 140°C	
CAS: 72-17-3	Melting point: 17.2°C	
MW: 112.06	EINECS: 200-772-0	
Density at 20°C: 1.3	P262	

Sodium lauryl sulfate

See Sodium dodecyl sulfate, page 290

Sodium lauryl sulfate solution 0.05 M
Standard solution ready to use, For volumetric analysis

Packaging		Mfr. No
2.5 l	Glass bottle	J/9021/17
CAS: 151-21-3		
Density at 20°C: >> 1.0		
EINECS: 205-788-1		

Manufactured to a tolerance of ± 0.1%.
See the full range of Standard Volumetric solution page...

Sodium metabisulfite
for analysis

packaging	Mfr. No	
500 g Plastic jar	S/5250/53	
1 kg Plastic jar	S/5250/60	
Na ₂ S ₂ O ₃	EINECS: 231-673-0	
CAS: 7681-57-4	GHS: Acute Tox. 4, Eye Dam. 1	
MW: 190.10	H302, H318, EUH031	
Density at 20°C: 1.4	P280, P305+P351+P338,	
Melting point: 150°C	P310, P233, P301+P312	

Assay: 99 to 100.5%		
---------------------	--	--

Guaranteed analysis			
Heavy metals	<20 ppm	Magnesium (Mg)	<20 ppm
Thiosulfates test	: Pass	Potassium (K)	<500 ppm
Calcium (Ca)	<50 ppm	Total chloride (Cl)	<0.005%
Copper (Cu)	<5 ppm	Total silicon (Si)	<20 ppm
Iron (Fe)	<5 ppm	Zinc (Zn)	<5 ppm
Lead (Pb)	<5 ppm		

Typical analysis			
Aluminium (Al)	<2 ppm	Molybdenum (Mo)	<2 ppm
Arsenic (As)	<0.5 ppm	Nickel (Ni)	<2 ppm
Barium (Ba)	<2 ppm	Silver (Ag)	<2 ppm
Cadmium (Cd)	<2 ppm	Strontium (Sr)	<2 ppm
Chromium (Cr)	<2 ppm	Titanium (Ti)	<2 ppm
Cobalt (Co)	<2 ppm	Vanadium (V)	<2 ppm
Lithium (Li)	<2 ppm		

Lot analysis on the label.

Sodium metabisulfite
specified reagents for general laboratory work- extra pure

packaging	Mfr. No	
500 g Plastic jar	S/5240/53	
1 kg Plastic jar	S/5240/60	
2.5 kg Plastic jar	S/5240/62	
Na ₂ S ₂ O ₃	EINECS: 231-673-0	
CAS: 7681-57-4	GHS: Acute Tox. 4, Eye Dam. 1	
MW: 190.10	H302, H318, EUH031	
Density at 20°C: 1.4	P280, P305+P351+P338,	
Melting point: 150°C	P310, P233, P301+P312	

Assay: >93%		
-------------	--	--

Guaranteed analysis			
Calcium (Ca)	<200 ppm	Potassium (K)	<500 ppm
Copper (Cu)	<20 ppm	n-Butyl formate	<0.5%
Iron (Fe)	<20 ppm	Total silicon (Si)	<100 ppm
Lead (Pb)	<20 ppm	Zinc (Zn)	<20 ppm
Magnesium (Mg)	<100 ppm		

High Volume Solvent Delivery Systems

High-volume solvent delivery systems are available in 10L to 1000L, in Europe only. Please email your contact details and a brief description of your requirements to bege1.sha.quotes@thermofisher.com

Sodium metaborate
specified reagents for general laboratory work- extra pure

packaging	Mfr. No	
500 g Plastic jar	S/5280/53	
NaBO ₂ ·4H ₂ O	GHS: Repr. 2, Skin Irrit. 2,	
CAS: 10555-76-7	Eye Irrit. 2, STOT SE 3	
MW: 137.87	H315, H319, H335, H361fd	
Melting point: 57°C	P280, P305+P351+P338,	
EINECS: 231-891-6	P201, P308+P313	

Assay: >98%			
Guaranteed analysis			
Sulfate (SO ₄)	<0.03%	n-Butyl formate	<0.01%

Sodium metasilicate pentahydrate
Technical

packaging	Mfr. No	
1 kg Plastic jar	S/5360/60	
3 kg Plastic jar	S/5360/63	
Na ₂ O ₃ ·Si ₂ ·5H ₂ O	UN: 3253	
CAS: 10213-79-3	GHS: Skin Corr. 1B, STOT SE 3	
MW: 212.14	H314, H335	
Melting point: 72°C	P280, P305+P351+P338,	
ADR/RID: 8, III	P310, P301+P330+P331, P261	
EINECS: 229-912-9		

Guaranteed analysis	
Identity test	: Pass

Sodium molybdate dihydrate
for analysis

packaging	Mfr. No	
100 g Plastic jar	S/5480/48	
250 g Plastic jar	S/5480/50	
Na ₂ MoO ₄ ·2H ₂ O	Density at 20°C: 3.28	
CAS: 10102-40-6	EINECS: 231-551-7	
MW: 241.95	P262	

Assay: >99%			
Guaranteed analysis			
pH (aqueous solution 5%)	7 to 10	Magnesium (Mg)	<20 ppm
Solution	: Pass	Potassium (K)	<500 ppm
Calcium (Ca)	<50 ppm	Total nitrogen (N)	<50 ppm
Total chloride (Cl)	<0.001 %	Total sulphur (S)	<200 ppm
Copper (Cu)	<20 ppm	Zinc (Zn)	<10 ppm
Iron (Fe)	<20 ppm	n-Butyl formate	<0.001%
Lead (Pb)	<20 ppm		

Typical analysis			
Barium (Ba)	<2 ppm	Nickel (Ni)	<2 ppm
Chromium (Cr)	<2 ppm	Silver (Ag)	<2 ppm
Cobalt (Co)	<2 ppm	Strontium (Sr)	<2 ppm
Lithium (Li)	<2 ppm	Titanium (Ti)	<2 ppm
Manganese (Mn)	<2 ppm	Vanadium (V)	<2 ppm

Lot analysis on the label.

Sodium molybdate dihydrate
specified reagents for general laboratory work- extra pure

packaging		Mfr. No
250 g	Plastic jar	S/5440/50
500 g	Plastic jar	S/5440/53
Na ₂ MoO ₄ ·2H ₂ O		Density at 20°C: 3.28
CAS: 10102-40-6		EINECS: 231-551-7
MW: 241.95		P262
Assay: >95%		

Sodium nitrate
for analysis

packaging	Mfr. No	
500 g Plastic jar	S/5560/53	
3 kg Plastic jar	S/5560/63	
NaNO ₃	GHS: Ox. Sol. 2, Acute Tox.	
CAS: 7631-99-4	4, Skin Irrit. 2, Eye Irrit. 2,	
MW: 84.99	STOT SE 3	
Boiling point: 380°C	H272, H315, H335, H302, H319	
Melting point: 306°C	P261, P301+P312,	
ADR/RID: 5.1, III	P302+P352, P280,	
EINECS: 231-554-3	P305+P351+P338, P210	
UN: 1498		

Assay: >99.5%			
Guaranteed analysis			
pH (aqueous solution 5%)	5 to 7	Nitrite (NO ₂)	<0.001%
Calcium (Ca)	<20 ppm	Potassium (K)	<100 ppm
Copper (Cu)	<2 ppm	Total chloride (Cl)	<0.0005%
Iodate (IO ₃)	<5 ppm	Total phosphorus (P)	<2 ppm
Iron (Fe)	<2 ppm	Total silicon (Si)	<20 ppm
Lead (Pb)	<2 ppm	Total sulphur (S)	<20 ppm
Magnesium (Mg)	<10 ppm	Zinc (Zn)	<5 ppm

Typical analysis			
Aluminium (Al)	<2 ppm	Manganese (Mn)	<1 ppm
Arsenic (As)	<2 ppm	Molybdenum (Mo)	<1 ppm
Barium (Ba)	<1 ppm	Nickel (Ni)	<1 ppm
Cadmium (Cd)	<1 ppm	Silver (Ag)	<1 ppm
Chromium (Cr)	<1 ppm	Strontium (Sr)	<1 ppm
Cobalt (Co)	<1 ppm	Titanium (Ti)	<1 ppm
Lithium (Li)	<1 ppm	Vanadium (V)	<1 ppm

Lot analysis on the label.

Sodium nitrate
specified reagents for general laboratory work- extra pure

packaging	Mfr. No	
500 g Plastic jar	S/5520/53	
1 kg Plastic jar	S/5520/60	
3 kg Plastic jar	S/5520/63	
NaNO ₃	GHS: Ox. Sol. 2, Acute Tox.	
CAS: 7631-99-4	4, Skin Irrit. 2, Eye Irrit. 2,	
MW: 84.99	STOT SE 3	
Boiling point: 380°C	H272, H315, H335, H302, H319	
Melting point: 306°C	P261, P301+P312,	
ADR/RID: 5.1, III	P302+P352, P280,	
EINECS: 231-554-3	P305+P351+P338, P210	
UN: 1498		

Assay: >98%			
Guaranteed analysis			
Calcium (Ca)	<100 ppm	n-Butyl formate	<0.005%
Copper (Cu)	<10 ppm	Total phosphorus (P)	<10 ppm
Iron (Fe)	<10 ppm	Total silicon (Si)	<100 ppm
Lead (Pb)	<10 ppm	Total sulphur (S)	<100 ppm
Magnesium (Mg)	<50 ppm	Zinc (Zn)	<20 ppm
Potassium (K)	<500 ppm		

Label

Fisher Chemical labels show GHS pictograms clearly. Together with actual lot analysis shown on the label, plus a wealth of additional information, we help you to work safely to achieve accurate and consistent results in your analytical applications.

Sodium nitrite		Certified AR	
for analysis			
packaging		Mfr. No	
500 g	Plastic jar	S/5640/53	
2.5 kg	Plastic jar	S/5640/62	
NaNO ₂		UN: 1500	
CAS: 7632-00-0		GHS: Ox. Sol. 3, Acute Tox. 3, Aquatic Acute 1	
MW: 69.00		H272, H301, H400	
Boiling point: 320°C		P301+P310, P304+P340, P273, P210	
Melting point: 271°C			
ADR/RID: 5.1, 6.1, III			
EINECS: 231-555-9			
Assay: >99%			
Guaranteed analysis			
Calcium (Ca)	<10 ppm	Total chloride (Cl)	<0.002%
Copper (Cu)	<5 ppm	Total phosphorus (P)	<50 ppm
Iron (Fe)	<10 ppm	Total silicon (Si)	<50 ppm
Lead (Pb)	<5 ppm	Total sulphur (S)	<50 ppm
Magnesium (Mg)	<10 ppm	Zinc (Zn)	<10 ppm
Potassium (K)	<100 ppm		
Typical analysis			
Aluminium (Al)	<2 ppm	Manganese (Mn)	<2 ppm
Arsenic (As)	<0.5 ppm	Molybdenum (Mo)	<2 ppm
Barium (Ba)	<2 ppm	Nickel (Ni)	<2 ppm
Cadmium (Cd)	<2 ppm	Silver (Ag)	<2 ppm
Chromium (Cr)	<2 ppm	Strontium (Sr)	<2 ppm
Cobalt (Co)	<2 ppm	Titanium (Ti)	<2 ppm
Lithium (Li)	<2 ppm	Vanadium (V)	<2 ppm
Lot analysis on the label.			

Sodium nitrite		SLR
specified reagents for general laboratory work- extra pure		
packaging		Mfr. No
500 g	Plastic jar	S/5600/53
1 kg	Plastic jar	S/5600/60
2.5 kg	Plastic jar	S/5600/62
NaNO ₂		UN: 1500
CAS: 7632-00-0		GHS: Ox. Sol. 3, Acute Tox.
MW: 69.00		3, Aquatic Acute 1
Boiling point: 320°C		H272, H301, H400
Melting point: 271°C		P301+P310, P304+P340,
ADR/RID: 5.1, 6.1, III		P273, P210
EINECS: 231-555-9		
Assay: >97%		
Guaranteed analysis		
Calcium (Ca)	<50 ppm	n-Butyl formate
Copper (Cu)	<20 ppm	Total phosphorus (P)
Iron (Fe)	<50 ppm	Total silicon (Si)
Lead (Pb)	<20 ppm	Total sulphur (S)
Magnesium (Mg)	<50 ppm	Zinc (Zn)
Potassium (K)	<500 ppm	

Sodium nitrite solution 0.5M (0.5N)	
NIST Standard solution ready to use, For Volumetric analysis	
packaging	Mfr. No
1 ℓ Glass bottle	J/7900/15
GHS: Acute Tox. 4	
H302	
P261, P280, P301+P312	
Manufactured to a tolerance of ± 0.3%.	
See the full range of Standard Volumetric solution page 323.	

Sodium nitroprusside		Certified AR	
for analysis			
packaging		Mfr. No	
100 g	Plastic jar	S/5720/48	
C ₅ FeN ₆ Na ₂ O.2H ₂ O	GHS: Acute Tox. 3		
CAS: 13755-38-9	H301, H311, H331, EUH032		
MW: 297.94	P301+P310, P280, P312,		
Density at 20°C: 1.72	P302+P350, P310,		
ADR/RID: 6.1, II	P304+P340, P233		
UN: 1588			
Assay: >99%			
Guaranteed analysis			
Sulfate (SO ₄)	<0.01%	Ferricyanide	<0.01%
n-Butyl formate	<0.01%	Ferrocyanide	<0.02%
Lot analysis on the label			

Sodium nitroprusside		SLR
specified reagents for general laboratory work- extra pure		
packaging		Mfr. No
100 g	Plastic jar	S/5680/48
C ₅ FeN ₆ Na ₂ O.2H ₂ O	UN: 1588	
CAS: 13755-38-9	GHS: Acute Tox. 3	
MW: 297.94	H301, H311, H331, EUH032	
Density at 20°C: 1.72	P301+P310, P280, P312,	
ADR/RID: 6.1, II	P302+P350, P310,	
EINECS: 238-373-9	P304+P340, P233	
Assay: >98%		
Guaranteed analysis		
Ferricyanide	<0.02%	Sulfate (SO ₄)
Ferrocyanide	<0.1%	<0.05%

tri-Sodium orthophosphate dodecahydrate		Certified AR	
for analysis			
packaging	Mfr. No		
500 g Plastic jar	S/5761/53		
Na ₃ PO ₄ ·12H ₂ O	UN: 3262		
CAS: 10101-89-0	GHS: Skin Corr. 1A, Eye Dam. 1		
MW: 380.13	H314		
Density at 20°C: 1.62	P301+P330+P331, P305+P351+P338, P310, P260, P280, P264		
Melting point: 75°C			
ADR/RID: 8, II			
EINECS: 231-509-8			
Assay: 99 to 101%			
Guaranteed analysis			
Alkalinity	<0.7 meq/g	Manganese (Mn)	<2 ppm
Aluminium (Al)	<2 ppm	Molybdenum (Mo)	<2 ppm
Arsenic (As)	<2 ppm	Nickel (Ni)	<2 ppm
Barium (Ba)	<2 ppm	Potassium (K)	<50 ppm
Cadmium (Cd)	<2 ppm	Propan--1-ol	<2 ppm
Calcium (Ca)	<50 ppm	Strontium (Sr)	<2 ppm
Chromium (Cr)	<2 ppm	Titanium (Ti)	<2 ppm
Cobalt (Co)	<2 ppm	n-Butyl formate	<0.001%
Copper (Cu)	<5 ppm	Total nitrogen (N)	<10 ppm
Iron (Fe)	<10 ppm	Total silicon (Si)	<50 ppm
Lead (Pb)	<5 ppm	Total sulphur (S)	<50 ppm
Lithium (Li)	<2 ppm	Vanadium (V)	<2 ppm
Magnesium (Mg)	<20 ppm	Zinc (Zn)	<10 ppm
Lot analysis on the label.			

Standard Volumetric solution

For your convenience we offer a range of NIST traceable solutions. A full list of our products can be found on page 323.

Can't see what you need?

Please email your contact details and a brief description of your requirements to begeel.sha.quotes@thermofisher.com



tri-Sodium orthophosphate dodecahydrate		SLR
specified reagents for general laboratory work- extra pure		
packaging		Mfr. No
250 g	Plastic jar	S/5760/50
500 g	Plastic jar	S/5760/53
1 kg	Plastic jar	S/5760/60
5 kg	Plastic keg	S/5760/65
25 kg	Plastic keg	S/5760/70
Na ₃ PO ₄ .12H ₂ O		
CAS: 10101-89-0		
MW: 380.13		
Density at 20°C: 1.62		
Melting point: 75°C		
ADR/RID: 8, II		
EINECS: 231-509-8		
UN: 3262		
GHS: Skin Corr. 1A, Eye Dam. 1		
H314		
P301+P330+P331,		
P305+P351+P338, P310,		
P260, P280, P264		

Sodium oxalate		Certified AR	
for analysis			
packaging	Mfr. No		
250 g Plastic jar	S/5840/50		
500 g Plastic jar	S/5840/53		
C ₂ Na ₂ O ₄	EINECS: 200-550-3		
CAS: 62-76-0	GHS: Acute Tox. 4		
MW: 134.00	H312, H302		
Density at 20°C: 2.34	P280, P302+P352, P301+P312		
Melting point: 250 to 270°C			
Assay: >99.9%			
Guaranteed analysis			
pH (solution aq. 3%)	: 7 to 8.5	Potassium (K)	<50 ppm
Loss on drying	<0.05%	n-Butyl formate	<0.001%
Calcium (Ca)	<50 ppm	Total nitrogen (N)	<50 ppm
Copper (Cu)	<5 ppm	Total phosphorus (P)	<10 ppm
Iron (Fe)	<10 ppm	Total silicon (Si)	<10 ppm
Lead (Pb)	<5 ppm	Total sulphur (S)	<20 ppm
Magnesium (Mg)	<10 ppm	Zinc (Zn)	<10 ppm
Lot analysis on the label			

Sodium oxalate		pure
packaging	Mfr. No	
100 g Plastic jar	S/5800/48	
500 g Plastic jar	S/5800/53	
C ₂ Na ₂ O ₄	EINECS: 200-550-3	
CAS: 62-76-0	UN: 2811	
MW: 134.00	GHS: Acute Tox. 4	
Density at 20°C: 2.3421	H312, H302	
Melting point: 250 to 270°C	P280, P302+P352, P301+P312	
ADR/RID: 6.1/III		

Sodium perborate trihydrate		SLR
specified reagents for general laboratory work- extra pure		
packaging	Mfr. No	
500 g Plastic jar	S/5880/53	
NaBO ₃ .4H ₂ O	UN: 1479	
CAS: 10486-00-7	GHS: Repr. 1B, Eye Dam. 1, STOT SE 3	
MW: 153.87	H360Df, H319, H335	
Melting point: 60°C	P280, P305+P351+P338, P201, P308+P313, P261	
ADR/RID: 5.1/III		
EINECS: 231-556-4		
Guaranteed analysis		
Identity test	: Pass	

Sodium perchlorate		Certified
for HPLC		
packaging	Mfr. No	
250 g Plastic jar	S/5966/50	
NaClO ₄ ·H ₂ O	EINECS: 231-511-9	
CAS: 7791-07-3	UN: 1502	
MW: 140.46	GHS: Ox. Sol. 1, Acute Tox. 4	
Density at 20°C: 2.02	H271, H302	
Melting point: 130°C	P210, P301+P312	
ADR/RID: 5.1, II		
Assay: >99%		
Guaranteed analysis		
Iron (Fe)	<2 ppm	
Lead (Pb)	<2 ppm	
Absorbance specification		
Max. Absorbance (1cm cell v HPLC water)		
254 nm (sol. aqu. 1 M)		
Lot analysis on the label.		
Filtered to 0.2 µm.		

Sodium perchlorate		Certified AR	
for analysis			
packaging	Mfr. No		
500 g Plastic jar	S/5965/53		
NaClO ₄ ·H ₂ O	EINECS: 231-511-9		
CAS: 7791-07-3	UN: 1502		
MW: 140.46	GHS: Ox. Sol. 1, Acute Tox. 4		
Density at 20°C: 2.02	H271, H302		
Melting point: 130°C	P210, P301+P312		
ADR/RID: 5.1, II			
Assay: >99%			
Guaranteed analysis			
pH (aqueous solution 5%)	4.5 to 7	Magnesium (Mg)	<10 ppm
Calcium (Ca)	<20 ppm	Potassium (K)	<50 ppm
Chloride and chlorate	<0.002%	Total phosphorus (P)	<10 ppm
Copper (Cu)	<5 ppm	Total silicon (Si)	<10 ppm
Iron (Fe)	<5 ppm	Total sulphur (S)	<20 ppm
Lead (Pb)	<5 ppm	Zinc (Zn)	<5 ppm
Typical analysis			
Aluminium (Al)	<2 ppm	Molybdenum (Mo)	<2 ppm
Cadmium (Cd)	<2 ppm	Nickel (Ni)	<2 ppm
Chromium (Cr)	<2 ppm	Silver (Ag)	<2 ppm
Cobalt (Co)	<2 ppm	Strontium (Sr)	<2 ppm
Lithium (Li)	<2 ppm	Titanium (Ti)	<2 ppm
Manganese (Mn)	<2 ppm	Vanadium (V)	<2 ppm
Lot analysis on the label.			

Sodium perchlorate		SLR	
specified reagents for general laboratory work- extra pure			
packaging	Mfr. No		
250 g Plastic jar	S/5960/50		
500 g Plastic jar	S/5960/53		
NaClO ₄ ·H ₂ O	EINECS: 231-511-9		
CAS: 7791-07-3	UN: 1502		
MW: 140.46	GHS: Ox. Sol. 1, Acute Tox. 4		
Density at 20°C: 2.02	H271, H302		
Melting point: 130°C	P210, P301+P312		
ADR/RID: 5.1, II			
Assay: >98%			
Guaranteed analysis			
Calcium (Ca)	<200 ppm	Potassium (K)	<500 ppm
Chloride and chlorate	<0.02%	Total phosphorus (P)	<100 ppm
Copper (Cu)	<50 ppm	Total silicon (Si)	<100 ppm
Iron (Fe)	<50 ppm	Total sulphur (S)	<200 ppm
Lead (Pb)	<50 ppm	Zinc (Zn)	<50 ppm
Magnesium (Mg)	<100 ppm		

Sodium periodate
for analysis

Certified AR

packaging		Mfr. No
100 g	Plastic jar	S/5980/48
250 g	Plastic jar	S/5980/50
5 kg	Plastic jar	S/5980/65
NaO ₄ CAS: 7790-28-5 MW: 213.89 Melting point: 300°C ADR/RID: 5.1, 6.1, II EINECS: 232-197-6 UN: 3087		GHS: Ox. Sol. 2, Acute Tox. 3, Skin Irrit. 2, Eye Irrit. 2, STOT SE 3 H272, H315, H319, H335, H301 P210, P302+P352, P280, P261, P301+P312, P305+P351+P338



Assay: >99%		
Guaranteed analysis		
pH (aqueous solution 5%)	3.5 to 4.5	Iron (Fe) <5 ppm
Loss on drying	<0.05%	Lead (Pb) <5 ppm
Bromate, bromide, chlorate & chloride: < 0.01%		Magnesium (Mg) <10 ppm
Calcium (Ca)	<20 ppm	Potassium (K) <200 ppm
Copper (Cu)	<5 ppm	Total silicon (Si) <20 ppm
Iodide (I)	<0.001%	Total sulphur (S) <100 ppm
		Zinc (Zn) <10 ppm

Typical analysis		
Aluminium (Al)	<5 ppm	Molybdenum (Mo) <5 ppm
Cadmium (Cd)	<5 ppm	Nickel (Ni) <5 ppm
Chromium (Cr)	<5 ppm	Silver (Ag) <5 ppm
Cobalt (Co)	<5 ppm	Titanium (Ti) <5 ppm
Lithium (Li)	<5 ppm	Vanadium (V) <5 ppm
Manganese (Mn)	<5 ppm	

Lot analysis on the label

Sodium peroxide
for analysis

Certified AR

packaging		Mfr. No
500 g	Metal can	S/6041/53
Na ₂ O ₂ CAS: 1313-60-6 MW: 77.98 Density at 20°C: 2.8 Melting point: 460°C ADR/RID: 5.1, I		EINECS: 215-209-4 UN: 1504 GHS: Ox. Sol. 1, Skin Corr. 1A H271, H314 P301+P330+P331, P280, P305+P351+P338, P310, P210



Assay: >95%		
Guaranteed analysis		
Calcium (Ca)	<500 ppm	Total chloride (Cl) <0.005%
Copper (Cu)	<20 ppm	Total phosphorus (P) <20 ppm
Iron (Fe)	<20 ppm	Total silicon (Si) <100 ppm
Lead (Pb)	<20 ppm	Total sulphur (S) <20 ppm
Magnesium (Mg)	<500 ppm	Zinc (Zn) <50 ppm

Typical analysis		
Aluminium (Al)	<5 ppm	Molybdenum (Mo) <5 ppm
Cadmium (Cd)	<5 ppm	Nickel (Ni) <5 ppm
Chromium (Cr)	<5 ppm	Silver (Ag) <5 ppm
Cobalt (Co)	<5 ppm	Titanium (Ti) <5 ppm
Lithium (Li)	<5 ppm	Vanadium (V) <5 ppm
Manganese (Mn)	<5 ppm	

Sodium peroxide
extra pure, granulated

SLR

packaging		Mfr. No
500 g	Metal can	S/6040/53
Na ₂ O ₂ CAS: 1313-60-6 MW: 77.98 Density at 20°C: 2.8 Melting point: 460°C ADR/RID: 5.1, I		EINECS: 215-209-4 UN: 1504 GHS: Ox. Sol. 1, Skin Corr. 1A H271, H314 P301+P330+P331, P280, P305+P351+P338, P310, P210



Assay: > 94%		
Guaranteed analysis		
Total chloride (Cl)	<0.005%	Magnesium (Mg) <500 ppm
Lead (Pb)	<20 ppm	

Sodium peroxydisulfate
See Sodium persulfate, page 302

Sodium persulfate
specified reagents for general laboratory work- extra pure

SLR

packaging		Mfr. No
500 g	Plastic jar	S/6080/53
Na ₂ S ₂ O ₈ CAS: 7775-27-1 MW: 238.09 Density at 20°C: 2.6 Melting point: 100°C ADR/RID: 5.1, III EINECS: 231-892-1 UN: 1505 GHS: Ox. Sol. 3, Acute Tox.		4, Skin Irrit. 2, Eye Irrit. 2, Resp. Sens. 1, Skin Sens. 1, STOT SE 3 H272, H315, H302, H317, H334, H319, H335 P261, P280, P308+P313, P302+P352, P210, P305+P351+P338



Assay: >96%		
Guaranteed analysis		
Ammonia (NH ₃)	<0.2%	n-Butyl formate <0.005%

Sodium propionate
specified reagents for general laboratory work- extra pure

SLR

packaging		Mfr. No
500 g	Plastic jar	S/6100/53
C ₃ H ₇ NaO ₂ CAS: 137-40-6 MW: 96.06 Melting point: 285 to 286°C Flash point: >250°C		EINECS: 205-290-4 GHS: Acute Tox. 4, Skin Irrit. 2 H315, H312 P280, P302+P352



tetra-Sodium pyrophosphate decahydrate
for analysis

Certified AR

packaging		Mfr. No
500 g	Plastic jar	S/6140/53
Na ₄ P ₂ O ₇ ·12H ₂ O CAS: 13472-36-1 MW: 446.05 Boiling point: 93.8°C Melting point: 80°C EINECS: 205-290-4		GHS: Skin Irrit. 2, Eye Irrit. 2, STOT SE 3 H315, H319, H335 P261, P302+P352, P280, P305+P351+P338



Assay: 99 to 101%		
Guaranteed analysis		
pH (aqueous solution 5%)	10.2 to 10.7	Magnesium (Mg) <10 ppm
Cadmium (Cd)	<5 ppm	Potassium (K) <50 ppm
Calcium (Ca)	<20 ppm	Total chloride (Cl) <0.002%
Copper (Cu)	<2 ppm	Total silicon (Si) <50 ppm
Iron (Fe)	<2 ppm	Total sulphur (S) <50 ppm
Lead (Pb)	<2 ppm	Zinc (Zn) <10 ppm

Typical analysis		
Aluminium (Al)	<5 ppm	Molybdenum (Mo) <5 ppm
Barium (Ba)	<5 ppm	Nickel (Ni) <5 ppm
Chromium (Cr)	<5 ppm	Silver (Ag) <5 ppm
Cobalt (Co)	<5 ppm	Strontium (Sr) <5 ppm
Lithium (Li)	<5 ppm	Titanium (Ti) <5 ppm
Manganese (Mn)	<5 ppm	Vanadium (V) <5 ppm

Reagents, Solvents & Volumetric solutions comply with the European Pharmacopeia

Read more on page 234

tetra-Sodium pyrophosphate decahydrate
specified reagents for general laboratory work- extra pure

SLR

packaging		Mfr. No
1 kg	Plastic jar	S/6120/60
Na ₄ P ₂ O ₇ ·12H ₂ O CAS: 13472-36-1 MW: 446.05 Density at 20°C: 1.82 Boiling point: 93.8°C Melting point: 80°C		EINECS: 205-290-4 GHS: Skin Irrit. 2, Eye Irrit. 2, STOT SE 3 H315, H319, H335 P261, P302+P352, P280, P305+P351+P338
Assay: >98%		
Guaranteed analysis		
Calcium (Ca)	<100 ppm	Potassium (K) <200 ppm
Copper (Cu)	<10 ppm	n-Butyl formate <0.01%
Iron (Fe)	<10 ppm	Total silicon (Si) <200 ppm
Lead (Pb)	<10 ppm	Total sulphur (S) <200 ppm
Magnesium (Mg)	<50 ppm	Zinc (Zn) <50 ppm



Sodium pyrosulfite
See Sodium metabisulfite, page 298

Sodium salicylate
specified reagents for general laboratory work- extra pure

SLR

packaging		Mfr. No
500 g	Plastic jar	S/6200/53
25 kg	Fibre Keg	S/6200/70
C ₇ H ₅ NaO ₃ CAS: 54-21-7 MW: 160.11 Melting point: 200°C EINECS: 200-198-0 GHS: Acute Tox. 4, Skin Irrit.		2, Eye Irrit. 2, STOT SE 3 H315, H335, H302, H319 P261, P301+P312, P302+P352, P280, P305+P351+P338



Assay: >99.5%		
Guaranteed analysis		
Water	<1%	Sulfate (SO ₄) <0.1%
Arsenic (As)	<2 ppm	n-Butyl formate <0.04%
Lead (Pb)	<10 ppm	

Sodium selenite
specified reagents for general laboratory work- extra pure

SLR

packaging		Mfr. No
25 g	Plastic jar	S/6280/46
Na ₂ SeO ₃ CAS: 10102-18-8 MW: 172.94 Density at 20°C: 3.121 Melting point: 350°C ADR/RID: 6.1, I EINECS: 233-267-9		UN: 2630 GHS: Acute Tox. 2, Skin Sens. 1, Aquatic Chronic 2 H317, H331, H411, H300, EUH031 P301+P310, P304+P340, P280, P260, P273



Assay: >99%		
Guaranteed analysis		
pH (aqueous solution 2%)	9 to 11	

Sodium sesquicarbonate dihydrate
pure

packaging		Mfr. No
3 kg	Plastic jar	S/6310/63
Na ₂ CO ₃ ·NaHCO ₃ ·2H ₂ O CAS: 6106-20-3 MW: 226.03 EINECS: 208-580-9		GHS: Eye Irrit. 2 H319 P280, P305+P351+P338



Guaranteed analysis		
Identity test	: Pass	

Sodium silicate powder
Technical

packaging		Mfr. No
1 kg	Plastic jar	S/6400/60
25 kg	Fibre Keg	S/6400/70
SiO ₂ ·Na ₂ O CAS: 6834-92-0 MW: 122.06 Melting point: 72°C ADR/RID: 8, III EINECS: 229-912-9		UN: 3253 GHS: Skin Corr. 1B, STOT SE 3 H314, H335 P280, P301+P330+P331, P310, P305+P351+P338, P304+P340, P302+P352



Guaranteed analysis		
Identity test	: Pass	

Wt ratio SiO₂: Na₂O approx.2:1

Sodium silicate solution d=1.5
Technical

packaging		Mfr. No
2.5 ℓ	HDPE Plastic bottle	S/6340/17
25 ℓ	HDPE Plastic container	S/6340/25
SiO ₂ ·Na ₂ O CAS: 1344-09-8 Density at 20°C: 1.4 Boiling point: 100.6°C Melting point: -0.6°C ADR/RID: 8/II		EINECS: 215-687-4 UN: 1760 GHS: Skin Irrit. 2, Eye Dam. 1, STOT SE 3 H318, H315, H335 P280, P305+P351+P338



Guaranteed analysis		
Identity test	: Pass	

Wt ratio SiO₂: Na₂ O approx. 2:1

Sodium succinate disodium salt hexahydrate
for analysis

Certified AR

packaging		Mfr. No
100 g	Plastic jar	S/6490/48
10 kg	Plastic keg	S/6490/68
C ₄ H ₄ Na ₂ O ₄ ·6H ₂ O CAS: 6106-21-4 MW: 270.15 Melting point: >300°C		

Assay: 99.5 to 100.5%		
Guaranteed analysis		
Calcium (Ca)	<20 ppm	n-Butyl formate <0.001%
Copper (Cu)	<5 ppm	Total nitrogen (N) <20 ppm
Iron (Fe)	<5 ppm	Total silicon (Si) <20 ppm
Lead (Pb)	<5 ppm	Total sulphur (S) <20 ppm
Magnesium (Mg)	<10 ppm	Zinc (Zn) <10 ppm
Potassium (K)	<50 ppm	

Lot analysis on the label

Sodium succinate disodium salt hexahydrate
specified reagents for general laboratory work- extra pure

SLR

packaging		Mfr. No
500 g	Plastic jar	S/6480/53
C ₄ H ₄ Na ₂ O ₄ ·6H ₂ O CAS: 6106-21-4 MW: 270.15		Melting point: >300°C EINECS: 205-778-7

Assay: >99%		
Guaranteed analysis		
Calcium (Ca)	<100 ppm	Potassium (K) <200 ppm
Copper (Cu)	<20 ppm	n-Butyl formate <0.005%
Iron (Fe)	<20 ppm	Total silicon (Si) <100 ppm
Lead (Pb)	<20 ppm	Total sulphur (S) <100 ppm
Magnesium (Mg)	<50 ppm	Zinc (Zn) <50 ppm

Sodium sulfate anhydrous <i>for residue and pesticide analysis</i>		Distol
packaging		Mfr. No
1 kg	Plastic jar	S/6660/60
Na ₂ SO ₄	Boiling point: 1700°C	
CAS: 7757-82-6	Melting point: 884 to 888°C	
MW: 142.04	EINECS: 231-820-9	
Assay: > 99.5%		
Guaranteed analysis		
Loss on drying	<0.2%	Potassium (K) <100 ppm

Sodium sulfate anhydrous for analysis, granular		Certified AR
packaging		Mfr. No
500 g	Plastic jar	S/6640/53
1 kg	Plastic jar	S/6640/60
3 kg	Plastic jar	S/6640/63
5 kg	Plastic jar	S/6640/65
10 kg	Plastic keg	S/6640/68
Na ₂ SO ₄		Boiling point: 1700°C
CAS: 7757-82-6		Melting point: 884 to 888°C
MW: 142.04		EINECS: 231-820-9

Lot analysis on the label.

Sodium sulfate anhydrous for analysis, fine powder		Certified AR
packaging		Mfr. No
500 g	Plastic jar	S/6650/53
1 kg	Plastic jar	S/6650/60
5 kg	Plastic jar	S/6650/65

Lot analysis on the label.

Sodium sulfate anhydrous		SLR
specified reagents for general laboratory work- extra pure		
packaging		Mfr. No
500 g	Plastic jar	S/6600/53
1 kg	Plastic jar	S/6600/60
3 kg	Plastic jar	S/6600/63
5 kg	Plastic jar	S/6600/65
25 kg	Plastic keg	S/6600/70
Na ₂ SO ₄	Boiling point: 1700°C	
CAS: 7757-82-6	Melting point: 884 to 888°C	
MW: 142.04	EINECS: 231-820-9	
Density at 20°C: 2.68		

Assay: >99%

Guaranteed analysis			
Loss on drying	<0.5%	Potassium (K)	<500 ppm
Calcium (Ca)	<200 ppm	n-Butyl formate	<0.01%
Copper (Cu)	<10 ppm	Total phosphorus (P)	<100 ppm
Iron (Fe)	<20 ppm	Total silicon (Si)	<100 ppm
Lead (Pb)	<10 ppm	Zinc (Zn)	<50 ppm
Magnesium (Mg)	<100 ppm		

Sodium sulfate decahydrate		Certified AR
for analysis		
packaging		Mfr. No
500 g	Plastic jar	S/6560/53
1 kg	Plastic jar	S/6560/60
Na ₂ SO ₄ ·10H ₂ O	Density at 20°C: 1.46	
CAS: 7727-73-3	Melting point: 884°C	
MW: 322.19	EINECS: 231-820-9	

Assay: 99 to 101%

Guaranteed analysis			
pH (aqueous solution 5%)	6 to 7	Magnesium (Mg)	<5 ppm
Reducing substances	<0.001%	Potassium (K)	<50 ppm
Oxidizing substance	Pass	Total chloride (Cl)	<0.0005%
Calcium (Ca)	<50 ppm	Total nitrogen (N)	<5 ppm
Copper (Cu)	<1 ppm	Total phosphorus (P)	<20 ppm
Iron (Fe)	<2 ppm	Total silicon (Si)	<20 ppm
Lead (Pb)	<1 ppm	Zinc (Zn)	<5 ppm
Typical analysis			
Aluminium (Al)	<2 ppm	Manganese (Mn)	<1 ppm
Arsenic (As)	<2 ppm	Molybdenum (Mo)	<1 ppm
Barium (Ba)	<1 ppm	Nickel (Ni)	<1 ppm
Cadmium (Cd)	<1 ppm	Silver (Ag)	<1 ppm
Chromium (Cr)	<1 ppm	Strontium (Sr)	<2 ppm
Cobalt (Co)	<1 ppm	Titanium (Ti)	<1 ppm
Lithium (Li)	<1 ppm	Vanadium (V)	<1 ppm

Lot analysis on the label.

Sodium sulfate decahydrate specified reagents for general laboratory work- extra pure		SLR
--	--	-----

packaging		Mfr. No	
1 kg	Plastic jar	S/6520/60	
Na ₂ SO ₄ ·10H ₂ O	Density at 20°C: 1.46		
CAS: 7727-73-3	Melting point: 884°C		
MW: 322.19	EINECS: 231-820-9		
Assay: >99%			
Guaranteed analysis			
Calcium (Ca)	<200 ppm	Potassium (K)	<200 ppm
Copper (Cu)	<5 ppm	n-Butyl formate	<0.02%
Iron (Fe)	<10 ppm	Total phosphorus (P)	<50 ppm
Lead (Pb)	<5 ppm	Total silicon (Si)	<100 ppm
Magnesium (Mg)	<20 ppm	Zinc (Zn)	<20 ppm

Sodium sulfite anhydrous for analysis		Certified AR
packaging		Mfr. No
500 g	Plastic jar	S/6850/53
1 kg	Plastic jar	S/6850/60
Na ₂ SO ₃	Density at 20°C: 2.63	
CAS: 7757-83-7	Melting point: >500°C	
MW: 126.05	EINECS: 231-821-4	
Assay: >98%		
Guaranteed analysis		
Calcium (Ca)	<50 ppm	Potassium (K) <500 ppm
Copper (Cu)	<10 ppm	Total chloride (Cl) <0.05%
Iron (Fe)	<10 ppm	Total silicon (Si) <50 ppm
Lead (Pb)	<10 ppm	Zinc (Zn) <10 ppm
Magnesium (Mg)	<20 ppm	
Typical analysis		
Aluminium (Al)	<5 ppm	Manganese (Mn) <5 ppm
Arsenic (As)	<0.5 ppm	Molybdenum (Mo) <5 ppm
Barium (Ba)	<5 ppm	Nickel (Ni) <5 ppm
Cadmium (Cd)	<5 ppm	Silver (Ag) <5 ppm
Chromium (Cr)	<5 ppm	Strontium (Sr) <5 ppm
Cobalt (Co)	<5 ppm	Titanium (Ti) <5 ppm
Lithium (Li)	<5 ppm	Vanadium (V) <5 ppm

Lot analysis on the label.

Sodium sulfite anhydrous specified reagents for general laboratory work- extra pure		SLR
--	--	-----

packaging		Mfr. No	
500 g	Plastic jar	S/6840/53	
1 kg	Plastic jar	S/6840/60	
3 kg	Plastic jar	S/6840/63	
Na ₂ SO ₃		Density at 20°C: 2.63	
CAS: 7757-83-7		Melting point: >500°C	
MW: 126.05		EINECS: 231-821-4	
Assay: >97%			
Guaranteed analysis			
Calcium (Ca)	<200 ppm	Potassium (K)	<1000 ppm
Copper (Cu)	<50 ppm	n-Butyl formate	<0.1%
Iron (Fe)	<50 ppm	Total silicon (Si)	<200 ppm
Lead (Pb)	<50 ppm	Zinc (Zn)	<50 ppm
Magnesium (Mg)	<100 ppm		

Sodium tauroglycocholate technical	
---------------------------------------	--

packaging	Mfr. No
100 g Plastic jar	S/6960/48
C ₂₈ H ₄₇ N ₂ NaO ₈ S	
CAS: 41945-48-6	
MW: 230.08	
EINECS: 255-596-7	
Guaranteed analysis	
Identity test : Pass	

di-Sodium tetraborate decahydrate for analysis		Certified AR
packaging		Mfr. No
250 g	Plastic jar	S/7040/50
500 g	Plastic jar	S/7040/53
1 kg	Plastic jar	S/7040/60
5 kg	Plastic keg	S/7040/65
Na ₂ B ₄ O ₇ ·10H ₂ O		EINECS: 215-540-4
CAS: 1303-96-4		GHS: Repr. 1B
MW: 381.36		H360FD
Density at 20°C: 17.3		P201, P308+P313
Melting point: 75°C		
Assay: 99.5 to 100.5%		

Guaranteed analysis			
pH (aqueous solution 5%)	9.2 to 9.4	Total chloride (Cl)	<0.0005%
Insoluble matter	<0.005%	Total phosphorus (P)	<10 ppm
Calcium (Ca)	<20 ppm	Total silicon (Si)	<50 ppm
Copper (Cu)	<5 ppm	Total sulphur (S)	<50 ppm
Iron (Fe)	<5 ppm	Zinc (Zn)	<10 ppm
Lead (Pb)	<5 ppm	Phosphate (PO ₄)	<10 ppm
Magnesium (Mg)	<10 ppm	Sulfate (SO ₄)	<50 ppm
Potassium (K)	<100 ppm		

Typical analysis			
Aluminium (Al)	<2 ppm	Manganese (Mn)	<2 ppm
Arsenic (As)	<2 ppm	Molybdenum (Mo)	<2 ppm
Barium (Ba)	<2 ppm	Nickel (Ni)	<2 ppm
Cadmium (Cd)	<2 ppm	Silver (Ag)	<2 ppm
Chromium (Cr)	<2 ppm	Strontium (Sr)	<2 ppm
Cobalt (Co)	<2 ppm	Titanium (Ti)	<2 ppm
Lithium (Li)	<2 ppm	Vanadium (V)	<2 ppm

Lot analysis on the label.

di-Sodium tetraborate decahydrate extra pure, crystals		SLR
---	--	-----

packaging		Mfr. No
1 kg	Plastic jar	S/7000/60
3 kg	Plastic jar	S/7000/63
Na ₂ B ₄ O ₇ ·10H ₂ O		EINECS: 215-540-4
CAS: 1303-96-4		GHS: Repr. 1B
MW: 381.36		H360FD
Density at 20°C: 17.3		P201, P308+P313
Melting point: 75°C		

Assay: >99%

Calcium (Ca)	<200 ppm	Potassium (K)	<500 ppm
Copper (Cu)	<50 ppm	n-Butyl formate	<0.02%
Iron (Fe)	<50 ppm	Total phosphorus (P)	<100 ppm
Lead (Pb)	<50 ppm	Total sulphur (S)	<500 ppm
Magnesium (Mg)	<100 ppm	Zinc (Zn)	<100 ppm

di-Sodium tetraborate decahydrate specified reagents for general laboratory work- extra pure		SLR, Eur. Ph.
---	--	---------------

packaging		Mfr. No
1 kg	Plastic jar	S/7010/60
3 kg	Plastic jar	S/7010/63
25 kg	Plastic keg	S/7010/70
Na ₂ B ₄ O ₇ ·10H ₂ O		EINECS: 215-540-4
CAS: 1303-96-4		GHS: Repr. 1B
MW: 381.36		H360FD
Density at 20°C: 17.3		P201, P308+P313
Melting point: 75°C		

Assay: >99%

Guaranteed analysis			
Calcium (Ca)	<200 ppm	Potassium (K)	<500 ppm
Copper (Cu)	<50 ppm	n-Butyl formate	<0.02%
Iron (Fe)	<50 ppm	Total phosphorus (P)	<100 ppm
Lead (Pb)	<50 ppm	Total sulphur (S)	<500 ppm
Magnesium (Mg)	<100 ppm	Zinc (Zn)	<100 ppm

Sodium tetrahydroborate

See Sodium borohydride, page 283

Sodium thiocyanate, anhydrous

specified reagents for general laboratory work- extra pure

SLR

packaging		Mfr. No
500 g	Plastic jar	S/7120/53
NaSCN	Chronic 3	
CAS: 540-72-7	H312, H412, H302, H332, EUH032	
MW: 81.07	P280, P302+P352, P273, P233, P301+P312, P304+P340	
Melting point: 287°C		
EINECS: 208-754-4		
GHS: Acute Tox. 4, Aquatic		
Assay: >98%		
Guaranteed analysis		
Calcium (Ca)	<200 ppm	Magnesium (Mg) <50 ppm
Copper (Cu)	<10 ppm	Potassium (K) <500 ppm
Iron (Fe)	<20 ppm	n-Butyl formate <0.05%
Lead (Pb)	<10 ppm	Zinc (Zn) <100 ppm

Sodium thiosulfate anhydrous

specified reagents for general laboratory work- extra pure

SLR

packaging		Mfr. No
500 g	Plastic jar	S/7250/53
1 kg	Plastic jar	S/7250/60
3 kg	Plastic jar	S/7250/63
Na ₂ S ₂ O ₃	Boiling point: 100°C	
CAS: 7772-98-7	Melting point: 48°C	
MW: 158.11	EINECS: 231-867-5	
Density at 20°C: 1.66		
Guaranteed analysis		
Identity test : Pass		

Sodium thiosulfate pentahydrate

for analysis

Certified AR

packaging		Mfr. No
500 g	Plastic jar	S/7240/53
1 kg	Plastic jar	S/7240/60
3 kg	Plastic jar	S/7240/63
5 kg	Plastic keg	S/7240/65
25 kg	Plastic keg	S/7240/70
Na ₂ S ₂ O ₃ ·5H ₂ O	Melting point: 48.5°C	
CAS: 10102-17-7	EINECS: 231-867-5	
MW: 248.18		
Assay: 99.5 to 100.5%		
Guaranteed analysis		
pH (aqueous solution 5%)	5.5 to 7.5	Potassium (K) <50 ppm
Sulfate & sulfite	<0.02%	Sulfide (S) <0.0002%
Calcium (Ca)	<20 ppm	Sulphate and sulphite >0.02 %
Copper (Cu)	<5 ppm	Total chloride (Cl) <0.02%
Iron (Fe)	<5 ppm	Total nitrogen (N) <50 ppm
Lead (Pb)	<5 ppm	Total silicon (Si) <20 ppm
Magnesium (Mg)	<10 ppm	Zinc (Zn) <10 ppm
Typical analysis		
Aluminium (Al)	<2 ppm	Molybdenum (Mo) <2 ppm
Barium (Ba)	<2 ppm	Nickel (Ni) <2 ppm
Cadmium (Cd)	<2 ppm	Silver (Ag) <2 ppm
Chromium (Cr)	<2 ppm	Strontium (Sr) <2 ppm
Cobalt (Co)	<2 ppm	Titanium (Ti) <2 ppm
Lithium (Li)	<2 ppm	Vanadium (V) <2 ppm
Manganese (Mn)	<2 ppm	

Lot analysis on the label.

Sodium thiosulfate pentahydrate

specified reagents for general laboratory work- extra pure

SLR

packaging		Mfr. No
500 g	Plastic jar	S/7200/53
1 kg	Plastic jar	S/7200/60
2 kg	Plastic jar	S/7200/61
3 kg	Plastic jar	S/7200/63
Na ₂ S ₂ O ₃ ·5H ₂ O	Density at 20°C: 1.72	
CAS: 10102-17-7	Melting point: 48.5°C	
MW: 248.18		EINECS: 231-867-5
Assay: >99%		
Guaranteed analysis		
Calcium (Ca)	<100 ppm	Potassium (K) <500 ppm
Copper (Cu)	<20 ppm	n-Butyl formate <0.1%
Iron (Fe)	<20 ppm	Total silicon (Si) <100 ppm
Lead (Pb)	<20 ppm	Zinc (Zn) <50 ppm
Magnesium (Mg)	<50 ppm	

Sodium thiosulfate pentahydrate

Technical

packaging		Mfr. No
3 kg	Plastic jar	S/7160/63
25 kg	Plastic keg	S/7160/70
Na ₂ S ₂ O ₃ ·5H ₂ O	Density at 20°C: 1.72	
CAS: 10102-17-7	Melting point: 48.5°C	
MW: 248.18		EINECS: 231-867-5

Sodium thiosulfate solution 0.1M (0.1N)

Eur.Ph. , USP, BP
NIST Standard solution ready to use, For Volumetric analysis

packaging		Mfr. No
1 ℓ	HDPE Plastic bottle	J/7950/15
2.5 ℓ	HDPE Plastic bottle	J/7950/17
5 ℓ	HDPE Plastic bottle	J/7950/21
10 ℓ	Bag in box	J/7950/24
Density at 20°C: 1.12		

Manufactured to a tolerance of ± 0.1%.

See the full range of Standard Volumetric solution page 323.

Sodium thiosulfate solution 0.1M (0.1N)

NIST Standard Concentrate, For Volumetric analysis

packaging		Mfr. No
1	Ampoule	J/7950C/05
6	Ampoule	J/7950C/90

Each ampoule makes 1 ℓ of solution.

See the full range of Standard Volumetric solution page 323.

Sodium thiosulfate solution 0.5M (0.5N)

NIST Standard Concentrate, For Volumetric analysis

packaging		Mfr. No
1	Ampoule	J/7951C/05
6	Ampoule	J/7951C/90

Density at 20°C: 1.0 to 1.1

Each ampoule makes 500 mL of solution.

See the full range of Standard Volumetric solution page 323.

Sodium thiosulfate solution 0.01N

Nist Standard solution ready to use, For volumetric analysis

Eur.Ph. BP

packaging		Mfr. No
2.5 ℓ	HDPE Plastic bottle	J/9011/17

Manufactured to a tolerance of ± 0.1%.

See the full range of Standard Volumetric solution page...

Sodium tungstate dihydrate

specified reagents for general laboratory work- extra pure

SLR

packaging		Mfr. No	
100 g	Plastic jar	S/7320/48	
Na ₂ WO ₄ ·2H ₂ O	GHS: Acute Tox. 4, Eye Irrit. 2		
CAS: 10213-10-2	H319, H302		
MW: 329.86	P301+P312, P280,		
Density at 20°C: 4.18	P305+P351+P338		
EINECS: 236-743-4			
Assay: 98 to 102%			
Guaranteed analysis			
pH (aqueous solution 2%)	7 to 10	Potassium (K)	<500 ppm
Calcium (Ca)	<50 ppm	n-Butyl formate	<0.01%
Copper (Cu)	<20 ppm	Total silicon (Si)	<500 ppm
Iron (Fe)	<20 ppm	Total sulphur (S)	<500 ppm
Lead (Pb)	<20 ppm	Zinc (Zn)	<20 ppm
Magnesium (Mg)	<20 ppm		



LC/MS application –including Optima LC/MS Solvents, Blends and Reagents

NEW

Solvents tested on LC-MS Thermo Scientific instrument to ensure suitability. Filtered to 0.2µm.

designation	packsize	package type	Mfr. No	page
Acetic acid	1 mL ampoule	Ampoule	A113-1AMP	13
	10x1 mL ampoule	Ampoule	A113-10X1AMP	13
	50 mL ampoule	HDPE Plastic bottle	A113-50	13
Acetonitrile	500 mL	Glass bottle	A955-500	18
	1 L	Glass bottle	A955-1	18
	2.5 L	Glass bottle	A955-212	18
	4 L	Glass bottle	A955-4	18
LCMSKIT	Box	Glass bottle	LCMSKIT	18
Acetonitrile	1 L	Glass bottle	A/0638/15	18
	2.5 L	Glass bottle	A/0638/17	18
Acetonitrile with 0.1% formic acid (v/v)	500 mL	Glass bottle	LS120-500	21
	1 L	Glass bottle	LS120-1	21
	2.5 L	Glass bottle	LS120-212	21
	4 L	Glass bottle	LS120-4	21
Acetonitrile with 0.1% trifluoroacetic acid (v/v)	500 mL	Glass bottle	LS121-500	22
	1 L	Glass bottle	LS121-1	22
	2.5 L	Glass bottle	LS121-212	22
	4 L	Glass bottle	LS121-4	22
Ammonium acetate	50 g	Glass bottle	A114-50	36
Ammonium formate	50 g	Glass bottle	A115-50	40
Formic acid	0.5 mL ampoule	Ampoule	A117-05AMP	132
	1 mL ampoule	Ampoule	A117-1AMP	132
	2 mL ampoule	Ampoule	A117-2AMP	132
	10x1 mL ampoule	Ampoule	A117-10X1AMP	132
	50 mL	HDPE Plastic bottle	A117-50	132
Isopropanol	500 mL	Glass bottle	A461-500	167
	1 L	Glass bottle	A461-1	167
	2.5 L	Glass bottle	A461-212	167
	4 L	Glass bottle	A461-4	167
Methanol	500 mL	Glass bottle	A456-500	199
	1 L	Glass bottle	A456-1	199
	2.5 L	Glass bottle	A456-212	199
	4 L	Glass bottle	A456-4	199
Methanol	1 L	Glass bottle	M/4062/15	199
	2.5 L	Glass bottle	M/4062/17	199
Trifluoroacetic acid	0.5 mL ampoule	Ampoule	A116-05AMP	358
	1 mL ampoule	Ampoule	A116-1AMP	358
	2 mL ampoule	Ampoule	A116-2AMP	358
	10x1 mL ampoule	Ampoule	A116-10X1AMP	358
	50 mL	HDPE bottle	A116-50	358
Water	500 mL	Glass bottle	W6-500	367
	1 L	Glass bottle	W6-1	367
	2.5 L	Glass bottle	W6-212	367
	4 L	Glass bottle	W6-4	367
Water	1 L	Glass bottle	W/0112/15	367
	2.5 L	Glass bottle	W/0112/17	367
Water with 0.1% formic acid (v/v)	500 mL	Glass bottle	LS118-500	369
	1 L	Glass bottle	LS118-1	369
	2.5 L	Glass bottle	LS118-212	369
	4 L	Glass bottle	LS118-4	369
Water with 0.1% trifluoroacetic acid (v/v)	500 mL	Glass bottle	LS119-500	370
	1 L	Glass bottle	LS119-1	370
	2.5 L	Glass bottle	LS119-212	370
	4 L	Glass bottle	LS119-4	370

UHPLC application - UHPLC gradient grade Solvents

NEW

Solvents tested on UHPLC Thermo Scientific instrument to ensure suitability.

Filtered to 0.1µm. Lot analysis on the label.

packsize	package type	Mfr. No	page
1 L	Coated Glass bottle	A/0650/PB15	18
2.5 L	Coated Glass bottle	A/0650/PB17	18
Box	Coated Glass bottle	UHPLCKIT	19
1 L	Coated Glass bottle	M/4070/PB15	199
2.5 L	Coated Glass bottle	M/4070/PB17	199
1 L	Coated Glass bottle	W/0120/PB15	367
2.5 L	Coated Glass bottle	W/0120/PB17	367

HPLC Gradient grade- Solvents, Blends and Reagents

NEW PACK.

NEW

Solvents tested on gradient analysis conditions. Filtered to 0.2µm. Lot analysis on the label.

designation	packsize	package type	Mfr. No	page
Acetonitrile Advanced gradient grade	1 L	Coated Glass bottle	A/0632/PB15	19
	2.5 L	Coated Glass bottle	A/0632/PB17	19
Acetonitrile HPLC gradient grade	1 L	Glass bottle	A/0627/15	19
	2.5 L	Glass bottle	A/0627/17	19
		Coated Glass bottle	A/0627/PB17	19
	5 L	Aluminium can	A/0627/21	19
	10 L	Returnable Drum	A/0627/10RSS	19
	25 L	Returnable Drum	A/0627/25RSS	19
Acetonitrile HPLC grade	30 L	Returnable Drum	A/0627/30RSS	19
	1 L	Glass bottle	A/0636/15	19
Acetonitrile with 0.1% formic acid (v/v)	2.5 L	Glass bottle	A/0636/17	19
	2.5 L	Glass bottle	FE/1100/17	22
Acetonitrile with 0.1% trifluoroacetic acid (v/v)	2.5 L	Glass bottle	FE/0200/17	22
Isopropanol HPLC gradient grade	2.5 L	Glass bottle	P/7508/17	167
Methanol Advanced gradient grade	1 L	Coated Glass bottle	M/4060/PB15	200
	2.5 L	Coated Glass bottle	M/4060/PB17	200
Methanol HPLC gradient grade	1 L	Glass bottle	M/4058/15	200
	2.5 L	Glass bottle	M/4058/17	200
		Coated Glass bottle	M/4058/PB17	200
	5 L	Aluminium can	M/4058/21	200
	10 L	Returnable Drum	M/4058/10RSS	200
	25 L	Returnable Drum	M/4058/25RSS	200
Methanol with 0.1% formic acid (v/v)	30 L	Returnable Drum	M/4058/30RSS	200
	2.5 L	Glass bottle	FE/1300/17	202
Water Advanced gradient grade	1 L	Coated Glass bottle	W/0110/PB15	368
	2.5 L	Coated Glass bottle	W/0110/PB17	368
Water HPLC gradient grade	1 L	Glass bottle	W/0106/15	368
	2.5 L	Glass bottle	W/0106/17	368
		Coated Glass bottle	W/0106/PB17	368
	10 L	Returnable Drum	W/0106/10RSS	368
	25 L	Returnable Drum	W/0106/25RSS	368
	30 L	Returnable Drum	W/0106/30RSS	368
Water with 0.1% Formic Acid HPLC gradient grade	2.5 L	Glass bottle	FE/1200/17	370
Acetonitrile with 0.1% Formic acid HPLC gradient grade	2.5 L	Glass bottle	FE/0100/17	370

HPLC grade- Solvents and Reagents

NEW PACK.

Solvents filtered to 0.2µm. Lot analysis on the label.

designation	packsize	package type	Mfr. No	page
Acetic acid glacial	500 mL	Coated Glass bottle	A/0406/PB08	13
	1 L	Coated Glass bottle	A/0406/PB15	13
Acetone	1 L	Glass bottle	A/0606/15	16
	2.5 L	Glass bottle	A/0606/17	16
	10 L	Returnable Drum	A/0606/10RSS	16
	25 L	Returnable Drum	A/0606/25RSS	16
	30 L	Returnable Drum	A/0606/30RSS	16

HPLC grade- Solvents and Reagents **NEW PACK.**

designation	packsize	package type	Mfr. No	page
Acetonitrile	1 L	Glass bottle	A/0626/15	20
	2.5 L	Glass bottle	A/0626/17	20
		Coated Glass bottle	A/0626/PB17	20
	5 L	Aluminium can	A/0626/21	20
	10 L	Returnable Drum	A/0626/10RSS	20
	25 L	Returnable Drum	A/0626/25RSS	20
	30 L	Returnable Drum	A/0626/30RSS	20
Ammonia solution 35% d=0.88	100 mℓ	Coated Glass bottle	A/3295/PB05	34
Ammonia solution 0.25M	500 mℓ	Coated Glass bottle	J/0100/PB08	35
Ammonium acetate	250 g	Plastic jar	A/3446/50	36
Ammonium acetate solution 0.25M	500 mℓ	HDPE plastic bottle	J/0200/PB08	37
Ammonium carbonate	250 g	Plastic jar	A/3686/50	38
1-Butane sulfonic acid sodium salt	25 g	Plastic jar	B/4791/46	67
tert-Butyl methyl ether	2.5 L	Glass bottle	M/4496/17	71
	10 L	Returnable Drum	M/4496/10RSS	71
	25 L	Returnable Drum	M/4496/25RSS	71
	30 L	Returnable Drum	M/4496/30RSS	71
1-Chlorobutane	1 L	Glass bottle	C/4756/15	85
	2.5 L	Glass bottle	C/4756/17	85
Chloroform stabilized with amylene	1 L	Glass bottle	C/4966/15	85
	2.5 L	Glass bottle	C/4966/17	85
		Coated Glass bottle	C/4966/PB17	85
	10 L	Returnable Drum	C/4966/10RSS	85
	25 L	Returnable Drum	C/4966/25RSS	85
	30 L	Returnable Drum	C/4966/30RSS	85
Cyclohexane	1 L	Glass bottle	C/8936/15	99
	2.5 L	Glass bottle	C/8936/17	99
		Coated Glass bottle	C/8936/PB17	99
	10 L	Returnable Drum	C/8936/10RSS	99
	25 L	Returnable Drum	C/8936/25RSS	99
	30 L	Returnable Drum	C/8936/30RSS	99
Di-n-butylamine	250 mℓ	Coated Glass bottle	D/1338/PB07	103
Di-n-butylamine phosphate	100 mℓ	Glass bottle	D/1341/05	104
1,2-Dichloroethane	1 L	Glass bottle	D/1756/15	105
	2.5 L	Glass bottle	D/1756/17	105
Dichloromethane unstabilized	1 L	Glass bottle	D/1857/15	106
	2.5 L	Glass bottle	D/1857/17	106
		Coated Glass bottle	D/1857/PB17	106
	10 L	Returnable Drum	D/1857/10RSS	106
	25 L	Returnable Drum	D/1857/25RSS	106
	30 L	Returnable Drum	D/1857/30RSS	106
Dichloromethane stabilized with amylene	1 L	Glass bottle	D/1856/15	106
	2.5 L	Glass bottle	D/1856/17	106
		Coated Glass bottle	D/1856/PB17	106
	5 L	Aluminium can	D/1856/21	106
	10 L	Returnable Drum	D/1856/10RSS	106
	25 L	Returnable Drum	D/1856/25RSS	106
	30 L	Returnable Drum	D/1856/30RSS	106
Dichloromethane stabilised with methanol	1 L	Glass bottle	D/1859/15	106
	2.5 L	Glass bottle	D/1859/17	106
Diethyl ether stabilized with ethanol	1 L	Glass bottle	D/2506/15	109
	2.5 L	Glass bottle	D/2506/17	109
		Coated Glass bottle	D/2506/PB17	109
Dimethylformamide	1 L	Glass bottle	D/3846/15	113
	2.5 L	Glass bottle	D/3846/17	113
		Coated Glass bottle	D/3846/PB17	113
	10 L	Returnable Drum	D/3846/10RSS	113
	25 L	Returnable Drum	D/3846/25RSS	113
	30 L	Returnable Drum	D/3846/30RSS	113
Dimethyl sulfoxide	250 mℓ	Coated Glass bottle	D/4125/PB07	115
	500 mℓ	Coated Glass bottle	D/4125/PB08	115
	2.5 L	Coated Glass bottle	D/4125/PB17	115
	10 L	Returnable Drum	D/4125/10RSS	115
	25 L	Returnable Drum	D/4125/25RSS	115
	30 L	Returnable Drum	D/4125/30RSS	115

HPLC grade- Solvents and Reagents **NEW PACK.**

designation	packsize	package type	Mfr. No	page
1,4-Dioxane stabilised with BHT	1 L	Coated Glass bottle	D/4556/PB15	116
	2.5 L	Coated Glass bottle	D/4556/PB17	116
Ethanol absolute	1 L	Glass bottle	E/0665/15	120
		Glass bottle	E/0665DF/15	120
	2.5 L	Glass bottle	E/0665/17	120
		Glass bottle	E/0665DF/17	120
	10 L	Returnable Drum	E/0665DF/10RS	120
	25 L	Returnable Drum	E/0665DF/25RS	120
	30 L	Returnable Drum	E/0665DF/30RS	120
Ethyl acetate	1 L	Glass bottle	E/0906/15	124
	2.5 L	Glass bottle	E/0906/17	124
		Coated Glass bottle	E/0906/PB17	124
	10 L	Returnable Drum	E/0906/10RSS	124
	25 L	Returnable Drum	E/0906/25RSS	124
	30 L	Returnable Drum	E/0906/30RSS	124
Ethyl acetate	2.5 L	Glass bottle	E/0907/17	124
Formic acid solution 0.25M	500 mℓ	Coated Glass bottle	J/4200/PB08	133
Heptane (approx. 99% n-Heptane)	1 L	Glass bottle	H/0106/15	138
	2.5 L	Glass bottle	H/0106/17	138
		Coated Glass bottle	H/0106/PB17	138
	10 L	Returnable Drum	H/0106/10RSS	138
	25 L	Returnable Drum	H/0106/25RSS	138
	30 L	Returnable Drum	H/0106/30RSS	138
Hexanes (95% n-Hexane approx.)	1 L	Glass bottle	H/0406/15	142
	2.5 L	Glass bottle	H/0406/17	142
		Coated Glass bottle	H/0406/PB17	142
	5 L	Aluminium can	H/0406/21	142
	10 L	Returnable Drum	H/0406/10RSS	142
	25 L	Returnable Drum	H/0406/25RSS	142
	30 L	Returnable Drum	H/0406/30RSS	142
Isohexane (contains <5% n-Hexane)	1 L	Coated Glass bottle	H/0405/PB15	166
	2.5 L	Coated Glass bottle	H/0405/PB17	166
	5 L	Aluminium can	H/0405/21	166
	10 L	Returnable Drum	H/0405/10RSS	166
	25 L	Aluminium can	H/0405/25	166
		Returnable Drum	H/0405/25RSS	166
	30 L	Returnable Drum	H/0405/30RSS	166
Isopropanol	1 L	Glass bottle	P/7507/15	167
	2.5 L	Glass bottle	P/7507/17	167
		Coated Glass bottle	P/7507/PB17	167
	10 L	Returnable Drum	P/7507/10RSS	167
	25 L	Returnable Drum	P/7507/25RSS	167
	30 L	Returnable Drum	P/7507/30RSS	167
Methanol	1 L	Glass bottle	M/4056/15	200
	2.5 L	Glass bottle	M/4056/17	200
		Coated Glass bottle	M/4056/PB17	200
	5 L	Aluminium can	M/4056/21	200
	10 L	Returnable Drum	M/4056/10RSS	200
	25 L	Returnable Drum	M/4056/25RSS	200
	30 L	Returnable Drum	M/4056/30RSS	200
Orthophosphoric acid solution 1M	500 mℓ	Coated Glass bottle	J/5655/PB08	221
Pentane (mixed isomers)	1 L	Glass bottle	P/1006/15	225
	2.5 L	Glass bottle	P/1006/17	225
	10 L	Returnable Drum	P/1006/10RSS	225
	25 L	Returnable Drum	P/1006/25RSS	225
	30 L	Returnable Drum	P/1006/30RSS	225
Perchloric acid 60%	100 mℓ	Coated Glass bottle	P/1290/PB05	228
Potassium dihydrogen orthophosphate	250 g	Plastic jar	P/4806/50	252
	1 kg	Plastic jar	P/4806/60	252
	25 kg	Plastic keg	P/4806/70	252
Propan-1-ol	1 L	Glass bottle	P/7486/15	264
	2.5 L	Glass bottle	P/7486/17	264
Sodium acetate trihydrate	250 g	Plastic jar	S/2046/50	281
	1 kg	Plastic jar	S/2046/60	281
Sodium chloride	500 g	Plastic jar	S/3165/53	285
Sodium dihydrogen orthophosphate 0.25M solution	500 mℓ	Coated Glass bottle	J/7600/PB08	289

HPLC grade- Solvents and Reagents **NEW PACK.**

designation	packsize	package type	Mfr. No	page
Sodium formate	250 g	Plastic jar	S/4086/50	290
Sodium hydrogen carbonate	250 g	Plastic jar	S/4246/50	291
di-Sodium hydrogen orthophosphate 0.25M solution	500 mL	Coated Glass bottle	J/7650/PB08	293
Sodium perchlorate	250 g	Plastic jar	S/5966/50	301
Tetrabutyl ammonium bromide	25 g	Plastic jar	T/0476/46	343
Tetra-n-butyl ammonium hydroxide	250 mL	HDPE Plastic bottle	T/0486/07	343
Tetra-n-butyl ammonium hydroxide	250 mL	HDPE Plastic bottle	T/0481/07	343
Tetradecyltrimethylammonium bromide	25 g	Plastic jar	T/0611/46	345
Tetraethylammonium bromide	25 g	Plastic jar	T/0621/46	345
Tetrahydrofuran unstabilized	1 L	Coated Glass bottle	T/0706/PB15	345
	2.5 L	Coated Glass bottle	T/0706/PB17	345
	10 L	Returnable Drum	T/0706/10RSS	345
	25 L	Returnable Drum	T/0706/25RSS	345
	30 L	Returnable Drum	T/0706/30RSS	345
Tetramethyl ammonium chloride	25 g	Plastic jar	T/0931/46	346
Tetramethyl ammonium hydrogen sulfate	5 g	Glass jar	T/0935/44	347
Toluene	1 L	Glass bottle	T/2306/15	353
	2.5 L	Glass bottle	T/2306/17	353
	10 L	Returnable Drum	T/2306/10RSS	353
	25 L	Returnable Drum	T/2306/25RSS	353
	30 L	Returnable Drum	T/2306/30RSS	353
Triethylamine	250 mL	Glass bottle	T/3203/07	357
Trifluoroacetic acid	25 mL	Coated Glass bottle	T/3258/04	358
	100 mL	Coated Glass bottle	T/3258/PB05	358
Trifluoroacetic acid solution 1%	500 mL	Coated Glass bottle	J/8800/PB08	359
Tri-n-hexylamine	100 mL	Coated Glass bottle	T/3296/PB05	359
2,2,4-Trimethylpentane	1 L	Glass bottle	T/3606/15	360
	2.5 L	Glass bottle	T/3606/17	360
	5 L	Aluminium can	T/3606/21	360
	10 L	Returnable Drum	T/3606/10RSS	360
	25 L	Returnable Drum	T/3606/25RSS	360
	30 L	Returnable Drum	T/3606/30RSS	360

Other Applications & Grades- Solvents and Reagents

designation	packsize	package type	Mfr. No	page
Acetic acid glacial	500 mL	Coated Glass bottle	A/0407/PB08	13
Acetonitrile HPLC for fluorescence detection	1 L	Coated Glass bottle	A/0630/PB15	20
	2.5 L	Coated Glass bottle	A/0630/PB17	20
Acetonitrile HPLC for electrochemical detection	1 L	Coated Glass bottle	A/0635/PB15	20
	2.5 L	Coated Glass bottle	A/0635/PB17	20
Ammonium acetate	250 g	Plastic jar	A/3450/50	36
Citric acid monohydrate	500 g	Plastic jar	C/6230/53	90
Ethylenediaminetetraacetic acid disodium salt solution 0.2N	2.5 L	Glass bottle	J/3710/17	127
Heptane (approx. 99% n-Heptane)	1 L	Glass bottle	H/0107/15	139
	2.5 L	Glass bottle	H/0107/17	139
Hexanes (95% n-Hexane approx.)	1 L	Coated Glass bottle	H/0409/PB15	142
	2.5 L	Coated Glass bottle	H/0409/PB17	142
Isohexane (contains <5% n-Hexane)	1 L	Coated Glass bottle	H/0410/PB15	166
	2.5 L	Coated Glass bottle	H/0410/PB17	166
Isopropanol	2.5 L	Glass bottle	P/7509/17	168
Methanol HPLC for fluorescence detection	2.5 L	Coated Glass bottle	M/4057/PB17	200
Methanol HPLC for electrochemical detection	2.5 L	Coated Glass bottle	M/4059/PB17	201
Orthophosphoric acid 85%	500 mL	Coated Glass bottle	O/0515/PB08	221
Potassium dihydrogen orthophosphate	250 g	Plastic jar	P/4810/50	252
Sodium acetate trihydrate	250 g	Plastic jar	S/2052/50	281
tri-Sodium citrate dihydrate	500 g	Plastic jar	S/3380/53	287
Sodium dihydrogen orthophosphate dihydrate	500 g	Plastic jar	S/3770/53	289
di-Sodium hydrogen orthophosphate dihydrate	500 g	Plastic jar	S/4460/53	292
Sodium hydroxide 46-48% solution	2.5 L	HDPE Plastic bottle	S/4940/17	295
Water HPLC for fluorescence detection	2.5 L	Glass bottle	W/0107/17	368
Water HPLC for electrochemical detection	2.5 L	Glass bottle	W/0108/17	368


Sorbic acid SLR
specified reagents for general laboratory work- extra pure

packaging	Mfr. No
100 g Plastic jar	S/7390/48
C ₆ H ₈ O ₂ CAS: 110-44-1 MW: 112.13 Density at 20°C: 1.2 Melting point: 132 to 135°C Flash point: 127°C	EINECS: 203-768-7 GHS: Skin Irrit. 2, Eye Irrit. 2, STOT SE 3 H315, H335, H319 P261, P302+P352, P280, P305+P351+P338
Assay: >98.5%	



Sorbitol SLR, Eur.Ph., USP
specified reagents for general laboratory work- extra pure

packaging	Mfr. No
500 g Plastic jar	S/7400/53
3 kg Plastic jar	S/7400/63
25 kg Plastic keg	S/7400/70
C ₆ H ₁₄ O ₆ CAS: 50-70-4 MW: 182.18	Density at 20°C: 1.489 Melting point: 95 to 99°C EINECS: 200-061-5
Guaranteed analysis	
Water	<2% Melting Point tested 95 to 100°C
Residue after ignition	<0.1%



Fisher Chemical Standard mono-element solution for AAS

Manufactured to a tolerance of ± 0.5%
See alphabetical entry in the main body of this catalogue for full details

designation	matrix	packsize	package type	Mfr. No	page
Aluminium standard metal solution 1000 ppm	HNO ₃ 1 mol/l	100 mℓ	HDPE Plastic bottle	J/8000/05	27
	HNO ₃ 1 mol/l	500 mℓ	HDPE Plastic bottle	J/8000/08	27
Antimony standard metal solution 1000 ppm	HCl 1 mol/l	100 mℓ	HDPE Plastic bottle	J/8002/05	49
	HCl 1 mol/l	500 mℓ	HDPE Plastic bottle	J/8002/08	49
Arsenic standard metal solution 1000 ppm	HCl 1 mol/l	100 mℓ	HDPE Plastic bottle	J/8003/05	50
	HCl 1 mol/l	500 mℓ	HDPE Plastic bottle	J/8003/08	50
Barium standard metal solution 1000 ppm	HNO ₃ 1 mol/l	100 mℓ	HDPE Plastic bottle	J/8004/05	52
	HNO ₃ 1 mol/l	500 mℓ	HDPE Plastic bottle	J/8004/08	52
Boron standard metal solution 1000 ppm	H ₂ O	100 mℓ	HDPE Plastic bottle	J/8006/05	59
	H ₂ O	500 mℓ	HDPE Plastic bottle	J/8006/08	59
Cadmium standard metal solution 1000 ppm	HNO ₃ 1 mol/l	100 mℓ	HDPE Plastic bottle	J/8005/05	71
	HNO ₃ 1 mol/l	500 mℓ	HDPE Plastic bottle	J/8005/08	71
Calcium standard metal solution 1000 ppm	HNO ₃ 1 mol/l	100 mℓ	HDPE Plastic bottle	J/8010/05	73
	HNO ₃ 1 mol/l	500 mℓ	HDPE Plastic bottle	J/8010/08	73
Chromium standard metal solution 1000 ppm	HNO ₃ 1 mol/l	100 mℓ	HDPE Plastic bottle	J/8015/05	87
	HNO ₃ 1 mol/l	500 mℓ	HDPE Plastic bottle	J/8015/08	87
Cobalt standard metal solution 1000 ppm	HNO ₃ 1 mol/l	100 mℓ	HDPE Plastic bottle	J/8020/05	91
	HNO ₃ 1 mol/l	500 mℓ	HDPE Plastic bottle	J/8020/08	91
Copper standard metal solution 1000 ppm	HNO ₃ 1 mol/l	100 mℓ	HDPE Plastic bottle	J/8025/05	94
	HNO ₃ 1 mol/l	500 mℓ	HDPE Plastic bottle	J/8025/08	94
Gold standard metal solution 1000 ppm	HCl 1 mol/l	500 mℓ	HDPE Plastic bottle	J/8028/08	137
Iron standard metal solution 1000 ppm	HNO ₃ 1 mol/l	100 mℓ	HDPE Plastic bottle	J/8030/05	162
	HNO ₃ 1 mol/l	500 mℓ	HDPE Plastic bottle	J/8030/08	162
Lead standard metal solution 1000 ppm	HNO ₃ 1 mol/l	100 mℓ	HDPE Plastic bottle	J/8035/05	177
	HNO ₃ 1 mol/l	500 mℓ	HDPE Plastic bottle	J/8035/08	177
Lithium standard metal solution 1000 ppm	HNO ₃ 1 mol/l	100 mℓ	HDPE Plastic bottle	J/8037/05	182
	HNO ₃ 1 mol/l	500 mℓ	HDPE Plastic bottle	J/8037/08	182
Magnesium standard metal solution 1000 ppm	HNO ₃ 1 mol/l	100 mℓ	HDPE Plastic bottle	J/8040/05	186
	HNO ₃ 1 mol/l	500 mℓ	HDPE Plastic bottle	J/8040/08	186
Manganese standard metal solution 1000 ppm	HNO ₃ 1 mol/l	100 mℓ	HDPE Plastic bottle	J/8045/05	191
	HNO ₃ 1 mol/l	500 mℓ	HDPE Plastic bottle	J/8045/08	191
Mercury standard metal solution 1000 ppm	HNO ₃ 1 mol/l	100 mℓ	HDPE Plastic bottle	J/8047/05	194
	HNO ₃ 1 mol/l	500 mℓ	HDPE Plastic bottle	J/8047/08	194
Molybdenum standard metal solution 1000 ppm	HNO ₃ 0.1 mol/l	100 mℓ	HDPE Plastic bottle	J/8050/05	209
	HNO ₃ 0.1 mol/l	500 mℓ	HDPE Plastic bottle	J/8050/08	209
Nickel standard metal solution 1000 ppm	HNO ₃ 1 mol/l	100 mℓ	HDPE Plastic bottle	J/8055/05	212
	HNO ₃ 1 mol/l	500 mℓ	HDPE Plastic bottle	J/8055/08	212
Potassium standard metal solution 1000 ppm	HNO ₃ 1 mol/l	100 mℓ	HDPE Plastic bottle	J/8060/05	245
	HNO ₃ 1 mol/l	500 mℓ	HDPE Plastic bottle	J/8060/08	245
Selenium standard metal solution 1000 ppm	HNO ₃ 1 mol/l	100 mℓ	HDPE Plastic bottle	J/8061/05	273
	HNO ₃ 1 mol/l	500 mℓ	HDPE Plastic bottle	J/8061/08	273
Silicon standard metal solution 1000 ppm	H ₂ O	100 mℓ	HDPE Plastic bottle	J/8062/05	275
	H ₂ O	500 mℓ	HDPE Plastic bottle	J/8062/08	275
Silver standard metal solution 1000 ppm	HNO ₃ 1 mol/l	100 mℓ	HDPE Plastic bottle	J/8064/05	276
	HNO ₃ 1 mol/l	500 mℓ	HDPE Plastic bottle	J/8064/08	276
Sodium standard metal solution 1000 ppm	HNO ₃ 1 mol/l	100 mℓ	HDPE Plastic bottle	J/8065/05	279
	HNO ₃ 1 mol/l	500 mℓ	HDPE Plastic bottle	J/8065/08	279
Strontium standard metal solution 1000 ppm	HNO ₃ 1 mol/l	500 mℓ	HDPE Plastic bottle	J/8166/08	331
Tellurium standard metal solution 1000 ppm	HCl 20%	100 mℓ	HDPE Plastic bottle	J/8100/05	342
	HCl 20%	500 mℓ	HDPE Plastic bottle	J/8100/08	342
Tin standard metal solution 1000 ppm	HCl 1 mol/l	100 mℓ	HDPE Plastic bottle	J/8067/05	349
	HCl 1 mol/l	500 mℓ	HDPE Plastic bottle	J/8067/08	349
Titanium standard metal solution 1000 ppm	HCl 2 mol/l	100 mℓ	HDPE Plastic bottle	J/8110/05	351
	HCl 2 mol/l	500 mℓ	HDPE Plastic bottle	J/8110/08	351
Tungsten standard metal solution 1000 ppm	NH3 0.1 mol/l	100 mℓ	HDPE Plastic bottle	J/8120/05	363
	NH3 0.1 mol/l	500 mℓ	HDPE Plastic bottle	J/8120/08	363
Vanadium standard metal solution 1000 ppm	HNO ₃ 1 mol/l	100 mℓ	HDPE Plastic bottle	J/8150/05	366
	HNO ₃ 1 mol/l	500 mℓ	HDPE Plastic bottle	J/8150/08	366
Zinc standard metal solution 1000 ppm	HNO ₃ 1 mol/l	100 mℓ	HDPE Plastic bottle	J/8070/05	373
	HNO ₃ 1 mol/l	500 mℓ	HDPE Plastic bottle	J/8070/08	373
Zirconium standard metal solution 1000 ppm	HCl 1 mol/l	500 mℓ	HDPE Plastic bottle	J/8200/08	377

Fisher Chemical Standard mono-element solution for ICP

Manufactured to a tolerance of ± 0.5%
All standards are manufactured from Ultra-pure metals/acids
All standards are guaranteed to contain <10 ppm total metal impurities
See alphabetical entry in the main body of this catalogue for full details

designation	packsize	package type	Mfr. No	page
Aluminium standard metal solution 1000 ppm	500 mℓ	HDPE Plastic bottle	J/8206/08	27
Aluminium standard metal solution 10 000 ppm	100 mℓ	HDPE Plastic bottle	J/8205/05	27
Antimony standard metal solution 1000 ppm	500 mℓ	HDPE Plastic bottle	J/8215/08	49
Antimony standard metal solution 10 000 ppm	100 mℓ	HDPE Plastic bottle	J/8210/05	49
Barium standard metal solution 1000 ppm	500 mℓ	HDPE Plastic bottle	J/8221/08	52
Barium standard metal solution 10 000 ppm	100 mℓ	HDPE Plastic bottle	J/8220/05	52
Boron standard metal solution 1000 ppm	500 mℓ	HDPE Plastic bottle	J/8231/08	59
Boron standard metal solution 10 000 ppm	100 mℓ	HDPE Plastic bottle	J/8230/05	59
Cadmium standard metal solution 1000 ppm	500 mℓ	HDPE Plastic bottle	J/8236/08	71
Cadmium standard metal solution 10 000 ppm	100 mℓ	HDPE Plastic bottle	J/8235/05	72
Calcium standard metal solution 1000 ppm	500 mℓ	HDPE Plastic bottle	J/8241/08	73
Calcium standard metal solution 10 000 ppm	100 mℓ	HDPE Plastic bottle	J/8240/05	73
Chromium standard metal solution 1000 ppm	500 mℓ	HDPE Plastic bottle	J/8246/08	87
Chromium standard metal solution 10 000 ppm	100 mℓ	HDPE Plastic bottle	J/8245/05	88
Cobalt standard metal solution 1000 ppm	500 mℓ	HDPE Plastic bottle	J/8251/08	91
Cobalt standard metal solution 10 000 ppm	100 mℓ	HDPE Plastic bottle	J/8250/05	91
Copper standard metal solution 1000 ppm	500 mℓ	HDPE Plastic bottle	J/8256/08	94
Copper standard metal solution 10 000 ppm	100 mℓ	HDPE Plastic bottle	J/8255/05	94
Iron standard metal solution 1000 ppm	500 mℓ	HDPE Plastic bottle	J/8271/08	162
Iron standard metal solution 10000 ppm	100 mℓ	HDPE Plastic bottle	J/8270/05	162
Lead standard metal solution 1000 ppm	500 mℓ	HDPE Plastic bottle	J/8276/08	178
Lead standard metal solution 10 000 ppm	100 mℓ	HDPE Plastic bottle	J/8275/05	178
Lithium standard metal solution 1000 ppm	500 mℓ	HDPE Plastic bottle	J/8281/08	183
Lithium standard metal solution 10 000 ppm	100 mℓ	HDPE Plastic bottle	J/8280/05	183
Magnesium standard metal solution 1000 ppm	500 mℓ	HDPE Plastic bottle	J/8286/08	186
Magnesium standard metal solution 10 000 ppm	100 mℓ	HDPE Plastic bottle	J/8285/05	186
Manganese standard metal solution 1000 ppm	500 mℓ	HDPE Plastic bottle	J/8291/08	191
Manganese standard metal solution 10 000 ppm	100 mℓ	HDPE Plastic bottle	J/8290/05	191
Mercury standard metal solution 1000 ppm	500 mℓ	HDPE Plastic bottle	J/8293/08	194
Mercury standard metal solution 10 000 ppm	100 mℓ	HDPE Plastic bottle	J/8292/05	194
Molybdenum standard metal solution 1000 ppm	500 mℓ	HDPE Plastic bottle	J/8296/08	209
Molybdenum standard metal solution 10000 ppm	100 mℓ	HDPE Plastic bottle	J/8294/05	209
Nickel standard metal solution 1000 ppm	500 mℓ	HDPE Plastic bottle	J/8297/08	212
Nickel standard metal solution 10 000 ppm	100 mℓ	HDPE Plastic bottle	J/8295/05	212
Phosphorus standard metal solution 1000 ppm	500 mℓ	HDPE Plastic bottle	J/8299/08	242
Phosphorus standard metal solution 10 000 ppm	100 mℓ	HDPE Plastic bottle	J/8298/05	242
Potassium standard metal solution 1000 ppm	500 mℓ	HDPE Plastic bottle	J/8301/08	246
Potassium standard metal solution 10 000 ppm	100 mℓ	HDPE Plastic bottle	J/8300/05	246
Selenium standard metal solution 1000 ppm	500 mℓ	HDPE Plastic bottle	J/8306/08	273
Selenium standard metal solution 10 000 ppm	100 mℓ	HDPE Plastic bottle	J/8305/05	273
Silver standard metal solution 1000 ppm	500 mℓ	HDPE Plastic bottle	J/8310/08	276
Silver standard metal solution 10 000 ppm	100 mℓ	HDPE Plastic bottle	J/8309/05	276
Sodium standard metal solution 1000 ppm	500 mℓ	HDPE Plastic bottle	J/8316/08	279
Sodium standard metal solution 10 000 ppm	100 mℓ	HDPE Plastic bottle	J/8315/05	280
Strontium standard metal solution 1000 ppm	500 mℓ	HDPE Plastic bottle	J/8321/08	331
Strontium standard metal solution 10 000 ppm	100 mℓ	HDPE Plastic bottle	J/8320/05	332
Sulfur standard standard metal solution 1000 ppm	500 mℓ	HDPE Plastic bottle	J/8324/08	335
Sulfur standard standard metal solution 10 000 ppm	100 mℓ	HDPE Plastic bottle	J/8323/05	335
Tellurium standard metal solution 1000 ppm	500 mℓ	HDPE Plastic bottle	J/8327/08	342
Tin standard metal solution 1000 ppm	500 mℓ	HDPE Plastic bottle	J/8326/08	349
Tin standard metal solution 10 000 ppm	100 mℓ	HDPE Plastic bottle	J/8325/05	349
Titanium standard metal solution 1000 ppm	500 mℓ	HDPE Plastic bottle	J/8331/08	352
Titanium standard metal solution 10 000 ppm	100 mℓ	HDPE Plastic bottle	J/8330/05	352
Tungsten standard metal solution 1000 ppm	500 mℓ	HDPE Plastic bottle	J/8336/08	363
Zinc standard metal solution 1000 ppm	500 mℓ	HDPE Plastic bottle	J/8341/08	373
Zinc standard metal solution 10 000 ppm	100 mℓ	HDPE Plastic bottle	J/8340/05	373

SPEX Standard mono-element solution for ICP & ICP-MS

NIST Certified, manufactured to a tolerance of ± 0.5%
All standards are manufactured from Ultra-pure metals/salts
Actual measured impurities of up to 30 elements
See alphabetical entry in the main body of this catalogue for full details
High purity acids for elemental analysis

designation	matrix	packsize	package type	Mfr. No	page
Aluminium standard metal solution 1000 ppm	for ICP solution in HNO ₃ 2%	100 mL	Plastic bottle	SPPLAL22Y	28
	for ICP solution in HNO ₃ 2%	250 mL	Plastic bottle	SPPLAL22T	28
	for ICP solution in HNO ₃ 2%	500 mL	Plastic bottle	SPPLAL22X	28
Aluminium standard metal solution 1000 ppm	for ICP solution in HCL 2%	500 mL	Plastic bottle	SPPLAL12X	28
Aluminium standard metal solution 10 000 ppm	for ICP solution in HNO ₃ 5%	100 mL	Plastic bottle	SPPLAL23Y	28
	for ICP solution in HNO ₃ 5%	500 mL	Plastic bottle	SPPLAL23X	28
Aluminium standard metal solution 10 000 ppm	for ICP solution in HCL 5%	500 mL	Plastic bottle	SPPLAL13X	28
Aluminium standard metal solution 1000 ppm	for ICP-MS solution in HNO ₃ 2%	100 mL	Plastic bottle	SPCLAL22Y	28
Antimony standard metal solution 1000 ppm	for ICP-MS solution in H ₂ O/0.6% Tartaric Acid/HNO ₃ trace	100 mL	Plastic bottle	SPCLSB72Y	49
Antimony standard metal solution 1000 ppm	for ICP solution in H ₂ O/0.6% Tartaric Acid/HNO ₃ trace	100 mL	Plastic bottle	SPPLSB72Y	49
	for ICP solution in H ₂ O/0.6% Tartaric Acid/HNO ₃ trace	250 mL	Plastic bottle	SPPLSB72T	49
	for ICP solution in H ₂ O/0.6% Tartaric Acid/HNO ₃ trace	500 mL	Plastic bottle	SPPLSB72X	49
Antimony standard metal solution 10 000 ppm	for ICP solution in H ₂ O/0.6% Tartaric Acid/HNO ₃ trace	100 mL	Plastic bottle	SPPLSB73Y	49
	for ICP solution in H ₂ O/0.6% Tartaric Acid/HNO ₃ trace	500 mL	Plastic bottle	SPPLSB73X	49
Antimony standard metal solution 10 000 ppm	for ICP solution in HCL 20%	500 mL	Plastic bottle	PLSB53X	50
Arsenic standard metal solution 1000 ppm	for ICP solution in HNO ₃ 2%	100 mL	Plastic bottle	SPPLAS22Y	50
	for ICP solution in HNO ₃ 2%	250 mL	Plastic bottle	SPPLAS22T	50
	for ICP solution in HNO ₃ 2%	500 mL	Plastic bottle	SPPLAS22X	50
Arsenic standard metal solution 1000 ppm	for ICP solution in HCL 2%	500 mL	Plastic bottle	SPPLAS12X	51
Arsenic standard metal solution 10 000 ppm	for ICP solution in HNO ₃ 5%	100 mL	Plastic bottle	SPPLAS23Y	51
	for ICP solution in HNO ₃ 5%	500 mL	Plastic bottle	SPPLAS23X	51
Arsenic standard metal solution 1000 ppm	for ICP-MS solution in HNO ₃ 2%	100 mL	Plastic bottle	SPCLAS22Y	51
Barium standard metal solution 1000 ppm	for ICP solution in HNO ₃ 2%	100 mL	Plastic bottle	SPPLBA22Y	52
	for ICP solution in HNO ₃ 2%	250 mL	Plastic bottle	SPPLBA22T	52
	for ICP solution in HNO ₃ 2%	500 mL	Plastic bottle	SPPLBA22X	52
Barium standard metal solution 1000 ppm	for ICP solution in HCL 2%	500 mL	Plastic bottle	PLBA12X	53
Barium standard metal solution 10 000 ppm	for ICP solution in HNO ₃ 5%	100 mL	Plastic bottle	SPPLBA23Y	53
	for ICP solution in HNO ₃ 5%	500 mL	Plastic bottle	SPPLBA23X	53
Barium standard metal solution 1000 ppm	for ICP-MS solution in HNO ₃ 2%	100 mL	Plastic bottle	SPCLBA22Y	53
Beryllium standard metal solution 1000 ppm	for ICP solution in HNO ₃ 2%	100 mL	Plastic bottle	SPPLBE22Y	57
	for ICP solution in HNO ₃ 2%	250 mL	Plastic bottle	SPPLBE22T	57
	for ICP solution in HNO ₃ 2%	500 mL	Plastic bottle	SPPLBE22X	57

SPEX Standard mono-element solution for ICP & ICP-MS

designation	matrix	packsize	package type	Mfr. No	page
Beryllium standard metal solution 10 000 ppm	for ICP solution in HNO ₃ 5%	100 mL	Plastic bottle	SPPLBE23Y	57
	for ICP solution in HNO ₃ 5%	500 mL	Plastic bottle	SPPLBE23X	57
Beryllium standard metal solution 1000 ppm	for ICP-MS solution in HNO ₃ 2%	100 mL	Plastic bottle	SPCLBE22Y	57
Bismuth standard metal solution 1000 ppm	for ICP solution in HNO ₃ 10%	100 mL	Plastic bottle	SPPLBI42Y	58
	for ICP solution in HNO ₃ 10%	500 mL	Plastic bottle	SPPLBI42X	58
Bismuth standard metal solution 10 ppm	for ICP-MS solution in HNO ₃ 2%	100 mL	Plastic bottle	SPCLBI21AY	58
Boron standard metal solution 1000 ppm	for ICP solution in H ₂ O	100 mL	Plastic bottle	SPPLB92Y	59
	for ICP solution in H ₂ O	250 mL	Plastic bottle	SPPLB92T	59
	for ICP solution in H ₂ O	500 mL	Plastic bottle	SPPLB92X	59
Boron standard metal solution 10 000 ppm	for ICP solution in H ₂ O	100 mL	Plastic bottle	SPPLB93Y	60
	for ICP solution in H ₂ O	500 mL	Plastic bottle	SPPLB93X	60
Cadmium standard metal solution 1000 ppm	for ICP solution in HNO ₃ 2%	100 mL	Plastic bottle	SPPLCD22Y	72
	for ICP solution in HNO ₃ 2%	250 mL	Plastic bottle	SPPLCD22T	72
	for ICP solution in HNO ₃ 2%	500 mL	Plastic bottle	SPPLCD22X	72
Cadmium standard metal solution 1000 ppm	for ICP solution in HCL 2%	500 mL	Plastic bottle	PLCD12X	72
Cadmium standard metal solution 10 000 ppm	for ICP solution in HNO ₃ 5%	100 mL	Plastic bottle	SPPLCD23Y	72
	for ICP solution in HNO ₃ 5%	500 mL	Plastic bottle	SPPLCD23X	72
Cadmium standard metal solution 1000 ppm	for ICP-MS solution in HNO ₃ 2%	100 mL	Plastic bottle	SPCLCD22Y	72
Calcium standard metal solution 1000 ppm	for ICP solution in HNO ₃ 2%	100 mL	Plastic bottle	SPPLCA22Y	73
	for ICP solution in HNO ₃ 2%	250 mL	Plastic bottle	SPPLCA22T	73
Calcium standard metal solution 1000 ppm	for ICP solution in HCL 2%	500 mL	Plastic bottle	SPPLCA12X	74
Calcium standard metal solution 10 000 ppm	for ICP solution in HNO ₃ 5%	100 mL	Plastic bottle	SPPLCA23Y	74
	for ICP solution in HNO ₃ 5%	250 mL	Plastic bottle	SPPLCA23T	74
	for ICP solution in HNO ₃ 5%	500 mL	Plastic bottle	SPPLCA23X	74
Calcium standard metal solution 10 000 ppm	for ICP solution in HCL 5%	500 mL	Plastic bottle	SPPLCA13X	74
Calcium standard metal solution 1000 ppm	for ICP-MS solution in HNO ₃ 2%	100 mL	Plastic bottle	SPCLCA22Y	74
Carbon standard solution 1000 ppm	for ICP solution in H ₂ O	100 mL	Plastic bottle	SPPLC92Y	78
	for ICP solution in H ₂ O	500 mL	Plastic bottle	SPPLC92X	78
Cerium standard metal solution 1000 ppm	for ICP solution in HNO ₃ 2%	100 mL	Plastic bottle	SPPLCE22Y	80
	for ICP solution in HNO ₃ 2%	500 mL	Plastic bottle	SPPLCE22X	80
Cerium standard metal solution 10 000 ppm	for ICP solution in HNO ₃ 5%	100 mL	Plastic bottle	SPPLCE23Y	80
	for ICP solution in HNO ₃ 5%	500 mL	Plastic bottle	SPPLCE23X	80
Cesium standard metal solution 1000 ppm	for ICP solution in HNO ₃ 2%	100 mL	Plastic bottle	SPPLCS22Y	81
	for ICP solution in HNO ₃ 2%	500 mL	Plastic bottle	SPPLCS22X	81
Cesium standard metal solution 10 000 ppm	for ICP solution in HNO ₃ 5%	100 mL	Plastic bottle	SPPLCS23Y	81
	for ICP solution in HNO ₃ 5%	500 mL	Plastic bottle	SPPLCS23X	81
Chromium standard metal solution 1000 ppm	for ICP solution in HNO ₃ 2%	100 mL	Plastic bottle	SPPLCR22Y	88
	for ICP solution in HNO ₃ 2%	250 mL	Plastic bottle	SPPLCR22T	88
	for ICP solution in HNO ₃ 2%	500 mL	Plastic bottle	SPPLCR22X	88
Chromium standard metal solution 1000 ppm	for ICP solution in HCL 2%	500 mL	Plastic bottle	SPPLCR12X	88
Chromium standard metal solution 1000 ppm	for ICP solution in H ₂ O	500 mL	Plastic bottle	SPPLCR92X	88
Chromium standard metal solution 10 000 ppm	for ICP solution in HNO ₃ 2%	100 mL	Plastic bottle	SPPLCR23Y	88
	for ICP solution in HNO ₃ 5%	500 mL	Plastic bottle	SPPLCR23X	88
Chromium standard metal solution 10 000 ppm	for ICP solution in HCL 5%	500 mL	Plastic bottle	PLCR13X	88
Chromium standard metal solution 10 000 ppm	for ICP solution in H ₂ O	500 mL	Plastic bottle	SPPLCR93X	88
Chromium standard metal solution 1000 ppm	for ICP-MS solution in HNO ₃ 2%	100 mL	Plastic bottle	SPCLCR22Y	88
Cobalt standard metal solution 1000 ppm	for ICP solution in HNO ₃ 2%	100 mL	Plastic bottle	SPPLCO22Y	91
	for ICP solution in HNO ₃ 2%	250 mL	Plastic bottle	SPPLCO22T	91
	for ICP solution in HNO ₃ 2%	500 mL	Plastic bottle	SPPLCO22X	91
Cobalt standard metal solution 1000 ppm	for ICP solution in HCL 2%	500 mL	Plastic bottle	SPPLCO12X	91
Cobalt standard metal solution 10 000 ppm	for ICP solution in HNO ₃ 5%	100 mL	Plastic bottle	SPPLCO23Y	92
	for ICP solution in HNO ₃ 5%	500 mL	Plastic bottle	SPPLCO23X	92
Cobalt standard metal solution 1000 ppm	for ICP-MS solution in HNO ₃ 2%	100 mL	Plastic bottle	SPCLCO22Y	92
Copper standard metal solution 1000 ppm	for ICP solution in HNO ₃ 2%	100 mL	Plastic bottle	SPPLCU22Y	95
	for ICP solution in HNO ₃ 2%	250 mL	Plastic bottle	SPPLCU22T	95
	for ICP solution in HNO ₃ 2%	500 mL	Plastic bottle	SPPLCU22X	95
Copper standard metal solution 1000 ppm	for ICP solution in HCL 2%	500 mL	Plastic bottle	SPPLCU12X	95
Copper standard metal solution 10 000 ppm	for ICP solution in HNO ₃ 5%	100 mL	Plastic bottle	SPPLCU23Y	95
	for ICP solution in HNO ₃ 5%	500 mL	Plastic bottle	SPPLCU23X	95
Copper standard metal solution 10 000 ppm	for ICP solution in HCL 5%	500 mL	Plastic bottle	SPPLCU13X	95
Copper standard metal solution 1000 ppm	for ICP-MS solution in HNO ₃ 2%	100 mL	Plastic bottle	SPCLCU22Y	95
Dysprosium standard metal solution 1000 ppm	for ICP solution in HNO ₃ 2%	100 mL	Plastic bottle	SPPLDY22Y	119
	for ICP solution in HNO ₃ 2%	500 mL	Plastic bottle	SPPLDY22X	119
Erbium standard metal solution 1000 ppm	for ICP solution in HNO ₃ 2%	100 mL	Plastic bottle	SPPLER22Y	119
	for ICP solution in HNO ₃ 2%	500 mL	Plastic bottle	SPPLER22X	119

SPEX Standard mono-element solution for ICP & ICP-MS

designation	matrix	packsize	package type	Mfr. No	page
Erbium standard metal solution 10 000 ppm	for ICP solution in HNO ₃ 5%	100 mL	Plastic bottle	SPPLER23Y	119
	for ICP solution in HNO ₃ 5%	500 mL	Plastic bottle	SPPLER23X	119
Europium standard metal solution 1000 ppm	for ICP solution in HNO ₃ 2%	100 mL	Plastic bottle	SPPLEU22Y	128
	for ICP solution in HNO ₃ 2%	500 mL	Plastic bottle	SPPLEU22X	128
Gadolinium standard metal solution 1000 ppm	for ICP solution in HNO ₃ 2%	100 mL	Plastic bottle	SPPLGD22Y	134
	for ICP solution in HNO ₃ 2%	500 mL	Plastic bottle	SPPLGD22X	134
Gadolinium standard metal solution 10 000 ppm	for ICP solution in HNO ₃ 5%	100 mL	Plastic bottle	SPPLGD23Y	134
	for ICP solution in HNO ₃ 5%	500 mL	Plastic bottle	SPPLGD23X	134
Gallium standard metal solution 1000 ppm	for ICP solution in HNO ₃ 2%	100 mL	Plastic bottle	SPPLGA22Y	134
	for ICP solution in HNO ₃ 2%	500 mL	Plastic bottle	SPPLGA22X	134
Germanium standard metal solution 1000 ppm	for ICP solution in H ₂ O/HF trace	100 mL	Plastic bottle	SPPLGE92Y	134
	for ICP solution in H ₂ O/HF trace	500 mL	Plastic bottle	SPPLGE92X	134
Germanium standard metal solution 10 ppm	for ICP-MS solution in H ₂ O/HF trace	100 mL	Plastic bottle	SPCLGE91AY	134
Gold standard metal solution 1000 ppm	for ICP solution in HCL 10%	100 mL	Plastic bottle	SPPLAU32Y	137
	for ICP solution in HCL 10%	500 mL	Plastic bottle	SPPLAU32X	137
Gold standard metal solution 1000 ppm	for ICP-MS solution in HCL 2%	100 mL	Plastic bottle	SPCLAU11Y	137
	for ICP solution in HCL 2%	100 mL	Plastic bottle	SPPLHF12Y	138
Hafnium standard metal solution 1000 ppm	for ICP solution in HCL 2%	500 mL	Plastic bottle	SPPLHF12X	138
Hafnium standard metal solution 10 000 ppm	for ICP solution in HCL 5%	500 mL	Plastic bottle	SPPLHF13X	
Holmium standard metal solution 1000 ppm	for ICP solution in HNO ₃ 2%	100 mL	Plastic bottle	SPPLHO22Y	144
	for ICP solution in HNO ₃ 2%	500 mL	Plastic bottle	SPPLHO22X	144
Indium standard metal solution 1000 ppm	for ICP solution in HNO ₃ 2%	100 mL	Plastic bottle	SPPLIN22Y	160
	for ICP solution in HNO ₃ 2%	500 mL	Plastic bottle	SPPLIN22X	160
Indium standard metal solution 10 ppm	for ICP-MS solution in HNO ₃ 2%	100 mL	Plastic bottle	SPCLIN21AY	160
Iridium standard metal solution 1000 ppm	for ICP solution in HCL 10%	100 mL	Plastic bottle	SPPLIR32Y	162
	for ICP solution in HCL 10%	500 mL	Plastic bottle	SPPLIR32X	162
Iron standard metal solution 1000 ppm	for ICP solution in HNO ₃ 2%	100 mL	Plastic bottle	SPPLFE22Y	162
	for ICP solution in HNO ₃ 2%	250 mL	Plastic bottle	SPPLFE22T	162
	for ICP solution in HNO ₃ 2%	500 mL	Plastic bottle	SPPLFE22X	162
Iron standard metal solution 1000 ppm	for ICP solution in HCL 2%	500 mL	Plastic bottle	SPPLFE12X	162
Iron standard metal solution 10 000 ppm	for ICP solution in HNO ₃ 5%	100 mL	Plastic bottle	SPPLFE23Y	163
	for ICP solution in HNO ₃ 5%	500 mL	Plastic bottle	SPPLFE23X	163
Iron standard metal solution 10 000 ppm	for ICP solution in HCL 5%	500 mL	Plastic bottle	SPPLFE13X	163
Iron standard metal solution 1000 ppm	for ICP-MS solution in HNO ₃ 2%	100 mL	Plastic bottle	SPCLFE22Y	163
Lanthanum standard metal solution 1000 ppm	for ICP solution in HNO ₃ 2%	100 mL	Plastic bottle	SPPLLA22Y	176
	for ICP solution in HNO ₃ 2%	500 mL	Plastic bottle	SPPLLA22X	176
Lanthanum standard metal solution 10 000 ppm	for ICP solution in HNO ₃ 5%	100 mL	Plastic bottle	SPPLLA23Y	177
	for ICP solution in HNO ₃ 5%	500 mL	Plastic bottle	SPPLLA23X	177
Lead standard metal solution 1000 ppm	for ICP solution in HNO ₃ 2%	100 mL	Plastic bottle	SPPLPB22Y	178
	for ICP solution in HNO ₃ 2%	250 mL	Plastic bottle	SPPLPB22T	178
	for ICP solution in HNO ₃ 2%	500 mL	Plastic bottle	SPPLPB22X	178
Lead standard metal solution 10 000 ppm	for ICP solution in HNO ₃ 5%	100 mL	Plastic bottle	SPPLPB23Y	178
	for ICP solution in HNO ₃ 5%	500 mL	Plastic bottle	SPPLPB23X	178
Lead standard metal solution 1000 ppm	for ICP-MS solution in HNO ₃ 2%	100 mL	Plastic bottle	SPCLPB22Y	178
Lithium standard metal solution 1000 ppm	for ICP solution in HNO ₃ 2%	100 mL	Plastic bottle	SPPLLI22Y	183
	for ICP solution in HNO ₃ 2%	500 mL	Plastic bottle	SPPLLI22X	183
Lithium standard metal solution 1000 ppm	for ICP solution in HCL 2%	500 mL	Plastic bottle	SPPLLI12X	183
Lithium standard metal solution 10 000 ppm	for ICP solution in HNO ₃ 5%	100 mL	Plastic bottle	SPPLLI23Y	183
	for ICP solution in HNO ₃ 5%	500 mL	Plastic bottle	SPPLLI23X	183
Lithium standard metal solution 10 000 ppm	for ICP solution in HCL 5%	500 mL	Plastic bottle	SPPLLI13X	183
Lutetium standard metal solution 1000 ppm	for ICP solution in HNO ₃ 2%	100 mL	Plastic bottle	SPPLLU22Y	186
	for ICP solution in HNO ₃ 2%	500 mL	Plastic bottle	SPPLLU22X	186
Magnesium standard metal solution 1000 ppm	for ICP solution in HNO ₃ 2%	100 mL	Plastic bottle	SPPLMG22Y	186
	for ICP solution in HNO ₃ 2%	250 mL	Plastic bottle	SPPLMG22T	186
Magnesium standard metal solution 1000 ppm	for ICP solution in HCL 2%	500 mL	Plastic bottle	SPPLMG12X	186
Magnesium standard metal solution 10 000 ppm	for ICP solution in HNO ₃ 5%	100 mL	Plastic bottle	SPPLMG23Y	186
	for ICP solution in HNO ₃ 5%	500 mL	Plastic bottle	SPPLMG23X	186
Magnesium standard metal solution 10 000 ppm	for ICP solution in HCL 5%	500 mL	Plastic bottle	SPPLMG13X	187
Magnesium standard metal solution 1000 ppm	for ICP-MS solution in HNO ₃ 2%	100 mL	Plastic bottle	SPCLMG22Y	187

SPEX Standard mono-element solution for ICP & ICP-MS

designation	matrix	packsize	package type	Mfr. No	page
Manganese standard metal solution 1000 ppm	for ICP solution in HNO ₃ 2%	100 mL	Plastic bottle	SPPLMN22Y	191
	for ICP solution in HNO ₃ 2%	250 mL	Plastic bottle	SPPLMN22T	191
	for ICP solution in HNO ₃ 2%	500 mL	Plastic bottle	SPPLMN22X	191
Manganese standard metal solution 10 000 ppm	for ICP solution in HNO ₃ 5%	100 mL	Plastic bottle	SPPLMN23Y	191
	for ICP solution in HNO ₃ 5%	500 mL	Plastic bottle	SPPLMN23X	191
Manganese standard metal solution 1000 ppm	for ICP-MS solution in HNO ₃ 2%	100 mL	Plastic bottle	SPCLMN22Y	191
Mercury standard metal solution 1000 ppm	for ICP solution in HNO ₃ 10%	100 mL	Plastic bottle	SPPLHG42Y	194
	for ICP solution in HNO ₃ 10%	250 mL	Plastic bottle	SPPLHG42T	194
	for ICP solution in HNO ₃ 10%	500 mL	Plastic bottle	SPPLHG42X	194
Mercury standard metal solution 10 ppm	for ICP solution in HNO ₃ 5%	100 mL	Plastic bottle	SPPLHG21AY	194
	for ICP solution in HNO ₃ 5%	500 mL	Plastic bottle	SPPLHG21AX	194
Mercury standard metal solution 100 ppm	for ICP solution in HNO ₃ 2%		Plastic bottle	SPPLHG21X	194
Mercury standard metal solution 10 000 ppm	for ICP solution in HNO ₃ 10%	100 mL	Plastic bottle	SPPLHG43Y	194
	for ICP solution in HNO ₃ 10%	500 mL	Plastic bottle	SPPLHG43X	194
Mercury standard metal solution 10 ppm	for ICP-MS solution in HNO ₃ 5%	100 mL	Plastic bottle	SPCLHG21AY	195
Molybdenum standard metal solution 1000 ppm	for ICP solution in H ₂ O	100 mL	Plastic bottle	SPPLMO92Y	209
	for ICP solution in H ₂ O	250 mL	Plastic bottle	SPPLMO92T	209
	for ICP solution in H ₂ O	500 mL	Plastic bottle	SPPLMO92X	209
Molybdenum standard metal solution 10 000 ppm	for ICP solution in H ₂ O	100 mL	Plastic bottle	SPPLMO93Y	209
	for ICP solution in H ₂ O	500 mL	Plastic bottle	SPPLMO93X	209
Molybdenum standard metal solution 1000 ppm	for ICP-MS solution in H ₂ O	100 mL	Plastic bottle	SPCLMO92Y	209
Neodymium standard metal solution 1000 ppm	for ICP solution in HNO ₃ 2%	100 mL	Plastic bottle	SPPLND22Y	211
	for ICP solution in HNO ₃ 2%	500 mL	Plastic bottle	SPPLND22X	211
Nickel standard metal solution 1000 ppm	for ICP solution in HNO ₃ 2%	100 mL	Plastic bottle	SPPLNI22Y	212
	for ICP solution in HNO ₃ 2%	250 mL	Plastic bottle	SPPLNI22T	212
	for ICP solution in HNO ₃ 2%	500 mL	Plastic bottle	SPPLNI22X	212
Nickel standard metal solution 10 000 ppm	for ICP solution in HNO ₃ 5%	100 mL	Plastic bottle	SPPLNI23Y	212
	for ICP solution in HNO ₃ 5%	500 mL	Plastic bottle	SPPLNI23X	212
Nickel standard metal solution 1000 ppm	for ICP-MS solution in HNO ₃ 2%	100 mL	Plastic bottle	SPCLNI22Y	212
Niobium standard metal solution 1000 ppm	for ICP solution in H ₂ O/HF trace	100 mL	Plastic bottle	SPPLNB92Y	215
	for ICP solution in H ₂ O/HF trace	500 mL	Plastic bottle	SPPLNB92X	215
Niobium standard metal solution 10 000 ppm	for ICP solution in H ₂ O/HF trace	100 mL	Plastic bottle	SPPLNB93Y	215
	for ICP solution in H ₂ O/HF trace	500 mL	Plastic bottle	SPPLNB93X	215
Palladium standard metal solution 1000 ppm	for ICP solution in HCL 10%	100 mL	Plastic bottle	SPPLPD32Y	223
	for ICP solution in HCL 10%	500 mL	Plastic bottle	SPPLPD32X	223
Phosphorus standard metal solution 1000 ppm	for ICP solution in H ₂ O	100 mL	Plastic bottle	SPPLP92Y	243
	for ICP solution in H ₂ O	250 mL	Plastic bottle	SPPLP92T	243
Phosphorus standard metal solution 10 000 ppm	for ICP solution in H ₂ O	100 mL	Plastic bottle	SPPLP93Y	243
	for ICP solution in H ₂ O	500 mL	Plastic bottle	SPPLP93X	243
Platinum standard metal solution 1000 ppm	for ICP solution in HCL 10%	100 mL	Plastic bottle	SPPLPT32Y	244
	for ICP solution in HCL 10%	500 mL	Plastic bottle	SPPLPT32X	244
Potassium standard metal solution 1000 ppm	for ICP solution in HNO ₃ 2%	100 mL	Plastic bottle	SPPLK22Y	246
	for ICP solution in HNO ₃ 2%	250 mL	Plastic bottle	SPPLK22T	246
	for ICP solution in HNO ₃ 2%	500 mL	Plastic bottle	SPPLK22X	246
Potassium standard metal solution 1000 ppm	for ICP solution in HCL 2%	500 mL	Plastic bottle	SPPLK12X	246
Potassium standard metal solution 10 000 ppm	for ICP solution in HNO ₃ 5%	100 mL	Plastic bottle	SPPLK23Y	246
	for ICP solution in HNO ₃ 5%	500 mL	Plastic bottle	SPPLK23X	246
Potassium standard metal solution 10 000 ppm	for ICP solution in HCL 2%	500 mL	Plastic bottle	SPPLK13X	246
Potassium standard metal solution 1000 ppm	for ICP-MS solution in HNO ₃ 2%	100 mL	Plastic bottle	SPCLK22Y	246
Praseodymium standard metal solution 1000 ppm	for ICP solution in HNO ₃ 2%	100 mL	Plastic bottle	SPPLPR22Y	264
	for ICP solution in HNO ₃ 2%	500 mL	Plastic bottle	SPPLPR22X	264
Rhenium standard metal solution 1000 ppm	for ICP solution in H ₂ O	100 mL	Plastic bottle	SPPLRE92Y	269
	for ICP solution in H ₂ O	500 mL	Plastic bottle	SPPLRE92X	269
Rhodium standard metal solution 1000 ppm	for ICP solution in HCL 10%	100 mL	Plastic bottle	SPPLRH32Y	269
	for ICP solution in HCL 10%	500 mL	Plastic bottle	SPPLRH32X	269
Rhodium standard metal solution 10 ppm	for ICP-MS solution in HCL 2%	100 mL	Plastic bottle	SPCLRH11AY	269
Rubidium standard metal solution 1000 ppm	for ICP solution in HNO ₃ 2%	100 mL	Plastic bottle	SPPLRB22Y	270
	for ICP solution in HNO ₃ 2%	500 mL	Plastic bottle	SPPLRB22X	270

SPEX Standard mono-element solution for ICP & ICP-MS

designation	matrix	packsize	package type	Mfr. No	page
Rubidium standard metal solution 10 000 ppm	for ICP solution in HNO ₃ 5%	100 mL	Plastic bottle	PLRB23Y	270
	for ICP solution in HNO ₃ 5%	500 mL	Plastic bottle	PLRB23X	270
Ruthenium standard metal solution 1000 ppm	for ICP solution in HCL 10%	100 mL	Plastic bottle	SPPLRU32Y	270
	for ICP solution in HCL 10%	500 mL	Plastic bottle	SPPLRU32X	270
Samarium standard metal solution 1000 ppm	for ICP solution in HNO ₃ 2%	100 mL	Plastic bottle	SPPLSM22Y	271
	for ICP solution in HNO ₃ 2%	500 mL	Plastic bottle	SPPLSM22X	271
Scandium standard metal solution 1000 ppm	for ICP solution in HNO ₃ 2%	100 mL	Plastic bottle	SPPLSC22Y	272
	for ICP solution in HNO ₃ 2%	250 mL	Plastic bottle	SPPLSC22T	272
	for ICP solution in HNO ₃ 2%	500 mL	Plastic bottle	SPPLSC22X	272
Scandium standard metal solution 10 000 ppm	for ICP solution in HNO ₃ 5%	100 mL	Plastic bottle	SPPLSC23Y	272
	for ICP solution in HNO ₃ 5%	500 mL	Plastic bottle	SPPLSC23X	272
Scandium standard metal solution 10 ppm	for ICP-MS solution in HNO ₃ 2%	100 mL	Plastic bottle	SPCLSC21AY	272
	for ICP solution in HNO ₃ 2%	100 mL	Plastic bottle	SPPLSE22Y	273
Selenium standard metal solution 1000 ppm	for ICP solution in HNO ₃ 2%	250 mL	Plastic bottle	SPPLSE22T	273
	for ICP solution in HNO ₃ 2%	500 mL	Plastic bottle	SPPLSE22X	273
	for ICP solution in HNO ₃ 5%	100 mL	Plastic bottle	SPPLSE23Y	273
Selenium standard metal solution 10 000 ppm	for ICP solution in HNO ₃ 5%	500 mL	Plastic bottle	SPPLSE23X	273
Silicon standard metal solution 1000 ppm	for ICP-MS solution in HNO ₃ 2%	100 mL	Plastic bottle	SPCLSE22Y	273
	for ICP solution in H ₂ O/HF trace	100 mL	Plastic bottle	SPPLSI92Y	275
	for ICP solution in H ₂ O/HF trace	250 mL	Plastic bottle	SPPLSI92T	275
Silicon standard metal solution 10 000 ppm	for ICP solution in H ₂ O/HF trace	500 mL	Plastic bottle	SPPLSI92X	275
	for ICP solution in H ₂ O/HF trace	100 mL	Plastic bottle	SPPLSI93Y	275
	for ICP solution in H ₂ O/HF trace	500 mL	Plastic bottle	SPPLSI93X	275
Silicon standard metal solution 10 000 ppm	for ICP solution in H ₂ O	500 mL	Plastic bottle	SPPLSI9A3X	276
	for ICP solution in HNO ₃ 2%	100 mL	Plastic bottle	SPPLAG22Y	276
	for ICP solution in HNO ₃ 2%	250 mL	Plastic bottle	SPPLAG22T	276
Silver standard metal solution 1000 ppm	for ICP solution in HNO ₃ 2%	500 mL	Plastic bottle	SPPLAG22X	276
	for ICP solution in HNO ₃ 5%	100 mL	Plastic bottle	SPPLAG23Y	276
	for ICP solution in HNO ₃ 5%	500 mL	Plastic bottle	SPPLAG23X	276
Silver standard metal solution 1000 ppm	for ICP-MS solution in HNO ₃ 2%	100 mL	Plastic bottle	SPCLAG22Y	276
	for ICP solution in HNO ₃ 2%	100 mL	Plastic bottle	SPPLNA22Y	280
Sodium standard metal solution 1000 ppm	for ICP solution in HNO ₃ 2%	250 mL	Plastic bottle	SPPLNA22T	280
	for ICP solution in HCL 10%	500 mL	Plastic bottle	SPPLNA22X	280
Sodium standard metal solution 1000 ppm	for ICP solution in HCL 2%	500 mL	Plastic bottle	SPPLNA12X	280
Sodium standard metal solution 10 000 ppm	for ICP solution in HNO ₃ 5%	100 mL	Plastic bottle	SPPLNA23Y	280
	for ICP solution in HNO ₃ 5%	500 mL	Plastic bottle	SPPLNA23X	280
Sodium standard metal solution 10 000 ppm	for ICP solution in HCL 5%	500 mL	Plastic bottle	SPPLNA13X	280
Sodium standard metal solution 1000 ppm	for ICP-MS solution in HNO ₃ 2%	100 mL	Plastic bottle	SPCLNA22Y	280
	for ICP solution in HNO ₃ 2%	100 mL	Plastic bottle	SPPLSR22Y	332
Strontium standard metal solution 1000 ppm	for ICP solution in HNO ₃ 2%	250 mL	Plastic bottle	SPPLSR22T	332
	for ICP solution in HNO ₃ 2%	500 mL	Plastic bottle	SPPLSR22X	332
Strontium standard metal solution 1000 ppm	for ICP solution in HCL 2%	500 mL	Plastic bottle	SPPLSR12X	332
Strontium standard metal solution 10 000 ppm	for ICP solution in HNO ₃ 5%	100 mL	Plastic bottle	SPPLSR23Y	332
	for ICP solution in HNO ₃ 5%	500 mL	Plastic bottle	SPPLSR23X	332
Strontium standard metal solution 10 000 ppm	for ICP solution in HCL 5%	500 mL	Plastic bottle	PLSR13X	332
	for ICP solution in H ₂ O	100 mL	Plastic bottle	SPPLS92Y	336
Sulfur standard metal solution 1000 ppm	for ICP solution in H ₂ O	250 mL	Plastic bottle	SPPLS92T	336
	for ICP solution in H ₂ O	500 mL	Plastic bottle	SPPLS92X	336
Sulfur standard metal solution 10 000 ppm	for ICP solution in H ₂ O	100 mL	Plastic bottle	SPPLS93Y	336
	for ICP solution in H ₂ O	500 mL	Plastic bottle	SPPLS93X	336
Tantalum standard metal solution 1000 ppm	for ICP solution in H ₂ O/HF trace	100 mL	Plastic bottle	SPPLTA92Y	341
	for ICP solution in H ₂ O/HF trace	500 mL	Plastic bottle	SPPLTA92X	341
Tantalum standard metal solution 10 000 ppm	for ICP solution in H ₂ O/HF trace	100 mL	Plastic bottle	SPPLTA93Y	341
	for ICP solution in H ₂ O/HF trace	500 mL	Plastic bottle	SPPLTA93X	341
Tellurium standard metal solution 1000 ppm	for ICP solution in HNO ₃ 2%	100 mL	Plastic bottle	SPPLTE22Y	342
	for ICP solution in HNO ₃ 5%	500 mL	Plastic bottle	SPPLTE22X	342
Tellurium standard metal solution 1000 ppm	for ICP solution in HCL 10%	500 mL	Plastic bottle	PLTE32X	342
Tellurium standard metal solution 10 000 ppm	for ICP solution in HCL 40%	500 mL	Plastic bottle	PLTE53X	342
Terbium standard metal solution 1000 ppm	for ICP solution in HNO ₃ 2%	100 mL	Plastic bottle	SPPLTB22Y	343
	for ICP solution in HNO ₃ 2%	500 mL	Plastic bottle	SPPLTB22X	343
Terbium standard metal solution 10 000 ppm	for ICP solution in HNO ₃ 5%	100 mL	Plastic bottle	PLTB23Y	343
	for ICP solution in HNO ₃ 5%	500 mL	Plastic bottle	PLTB23X	343
Terbium standard metal solution 10 ppm	for ICP-MS solution in HNO ₃ 2%	100 mL	Plastic bottle	SPCLTB21AY	343

SPEX Standard mono-element solution for ICP & ICP-MS

designation	matrix	packsize	package type	Mfr. No	page
Thallium standard metal solution 1000 ppm	for ICP solution in HNO ₃ 2%	100 mL	Plastic bottle	SPPLTL22Y	347
	for ICP solution in HNO ₃ 2%	250 mL	Plastic bottle	SPPLTL22T	347
	for ICP solution in HNO ₃ 2%	500 mL	Plastic bottle	SPPLTL22X	347
Thallium standard metal solution 10 000 ppm	for ICP solution in HNO ₃ 5%	100 mL	Plastic bottle	PLTL23Y	347
	for ICP solution in HNO ₃ 5%	500 mL	Plastic bottle	PLTL23X	347
Thallium standard metal solution 1000 ppm	for ICP-MS solution in HNO ₃ 2%	100 mL	Plastic bottle	SPCLTL22Y	347
Thorium standard metal solution 1000 ppm	for ICP solution in HNO ₃ 2%	100 mL	Plastic bottle	SPPLTH22Y	348
	for ICP solution in HNO ₃ 2%	500 mL	Plastic bottle	SPPLTH22X	348
Thorium standard metal solution 1000 ppm	for ICP-MS solution in HNO ₃ 2%	100 mL	Plastic bottle	SPCLTH22Y	348
Thullium standard metal solution 1000 ppm	for ICP solution in HNO ₃ 2%	100 mL	Plastic bottle	SPPLTM22Y	348
	for ICP solution in HNO ₃ 2%	500 mL	Plastic bottle	SPPLTM22X	348
Tin standard metal solution 1000 ppm	for ICP solution in HNO ₃ 1%/HF 1%	500 mL	Plastic bottle	SPPLSN22X	350
	for ICP solution in HCL 20%	100 mL	Plastic bottle	SPPLSN52Y	350
	for ICP solution in HCL 20%	250 mL	Plastic bottle	SPPLSN52T	350
Tin standard metal solution 1000 ppm	for ICP solution in HCL 20%	500 mL	Plastic bottle	SPPLSN52X	350
	for ICP solution in HNO ₃ 2%/HF 2%	500 mL	Plastic bottle	SPPLSN23X	350
Tin standard metal solution 10 000 ppm	for ICP solution in HCL 20%	100 mL	Plastic bottle	SPPLSN53Y	350
	for ICP solution in HCL 20%	500 mL	Plastic bottle	SPPLSN53X	350
Tin standard metal solution 1000 ppm	for ICP-MS solution in HNO ₃ 1%/HF 1%	100 mL	Plastic bottle	SPCLSN22Y	350
Titanium standard metal solution 1000 ppm	for ICP solution in HCL 20%	500 mL	Plastic bottle	SPPLTI52X	352
	for ICP solution in H ₂ O/HF trace	100 mL	Plastic bottle	SPPLTI92Y	352
	for ICP solution in H ₂ O/HF trace	250 mL	Plastic bottle	SPPLTI92T	352
Titanium standard metal solution 1000 ppm	for ICP solution in H ₂ O/HF trace	500 mL	Plastic bottle	SPPLTI92X	352
	for ICP solution in HCL 40%	500 mL	Plastic bottle	SPPLTI53X	352
Titanium standard metal solution 10 000 ppm	for ICP solution in H ₂ O/HF trace	100 mL	Plastic bottle	SPPLTI93Y	352
	for ICP solution in H ₂ O/HF trace	500 mL	Plastic bottle	SPPLTI93X	352
Titanium standard metal solution 1000 ppm	for ICP-MS solution in H ₂ O/HF trace	100 mL	Plastic bottle	SPCLTI92Y	352
Tungsten standard metal solution 1000 ppm	for ICP solution in HNO ₃ 1%/HF 2%	500 mL	Plastic bottle	SPPLW22X	363
Tungsten standard metal solution 1000 ppm	for ICP solution in H ₂ O/HF trace	100 mL	Plastic bottle	SPPLW92Y	363
Tungsten standard metal solution 10 000 ppm	for ICP solution in HNO ₃ 2%/HF 5%	500 mL	Plastic bottle	SPPLW23X	363
Tungsten standard metal solution 10 000 ppm	for ICP solution in H ₂ O/HF trace	100 mL	Plastic bottle	SPPLW93Y	363
	for ICP solution in H ₂ O/HF trace	500 mL	Plastic bottle	SPPLW93X	363
	for ICP solution in HNO ₃ 2%	100 mL	Plastic bottle	SPPLU22Y	365
Uranium standard metal solution 1000 ppm	for ICP solution in HNO ₃ 2%	500 mL	Plastic bottle	SPPLU22X	365
	for ICP solution in HNO ₃ 5%	100 mL	Plastic bottle	SPPLU23Y	365
Uranium standard metal solution 10 000 ppm	for ICP solution in HNO ₃ 5%	500 mL	Plastic bottle	SPPLU23X	365
Uranium standard metal solution 1000 ppm	for ICP-MS solution in HNO ₃ 2%	100 mL	Plastic bottle	SPCLU22Y	365
Vanadium standard metal solution 1000 ppm	for ICP solution in HNO ₃ 2%	100 mL	Plastic bottle	SPPLV22Y	366
	for ICP solution in HNO ₃ 2%	500 mL	Plastic bottle	SPPLV22X	366
Vanadium standard metal solution 1000 ppm	for ICP solution in HCL 2%	500 mL	Plastic bottle	SPPLV12X	366
Vanadium standard metal solution 10 000 ppm	for ICP solution in HNO ₃ 15%	100 mL	Plastic bottle	SPPLV43Y	366
	for ICP solution in HNO ₃ 15%	500 mL	Plastic bottle	SPPLV43X	366
Vanadium standard metal solution 10 000 ppm	for ICP solution in HCL 15%	500 mL	Plastic bottle	SPPLV33X	366
Vanadium standard metal solution 1000 ppm	for ICP-MS solution in HNO ₃ 2%	100 mL	Plastic bottle	SPCLV22Y	366
Water blank in HCL 5%	for ICP	500 mL	Plastic bottle	SPPLBLKHCL	369
Water blank in HNO ₃ 5%	for ICP	500 mL	Plastic bottle	SPPLBLKHNO3	369
Water blank in HCL 2%	for ICP-MS	100 mL	Plastic bottle	SPCLBLKHCL	369
Water blank in HNO ₃ 2%	for ICP-MS	100 mL	Plastic bottle	SPCLBLKHNO3	369
Water Blank	for ICP-MS	100 mL	Plastic bottle	SPCLBLKH2O	369
Ytterbium standard metal solution 1000 ppm	for ICP solution in HNO ₃ 2%	100 mL	Plastic bottle	SPPLY22Y	372
	for ICP solution in HNO ₃ 2%	250 mL	Plastic bottle	SPPLY22T	372
	for ICP solution in HNO ₃ 2%	500 mL	Plastic bottle	SPPLY22X	372
Ytterbium standard metal solution 1000 ppm	for ICP solution in HNO ₃ 2%	100 mL	Plastic bottle	SPPLYB22Y	372
	for ICP solution in HNO ₃ 2%	500 mL	Plastic bottle	SPPLYB22X	372
Ytterbium standard metal solution 10 ppm	for ICP-MS solution in HNO ₃ 2%	100 mL	Plastic bottle	SPCLY21AY	372
Yttrium standard metal solution 10 000 ppm	for ICP solution in HNO ₃ 5%	100 mL	Plastic bottle	SPPLY23Y	372
	for ICP solution in HNO ₃ 5%	500 mL	Plastic bottle	SPPLY23X	372
Zinc standard metal solution 1000 ppm	for ICP solution in HNO ₃ 2%	100 mL	Plastic bottle	SPPLZN22Y	373
	for ICP solution in HNO ₃ 2%	250 mL	Plastic bottle	SPPLZN22T	373
Zinc standard metal solution 1000 ppm	for ICP solution in HCL 2%	500 mL	Plastic bottle	SPPLZN12X	373
	for ICP solution in HNO ₃ 5%	100 mL	Plastic bottle	SPPLZN23Y	373
Zinc standard metal solution 10 000 ppm	for ICP solution in HNO ₃ 5%	500 mL	Plastic bottle	SPPLZN23X	373
	for ICP solution in HCL 5%	500 mL	Plastic bottle	SPPLZN13X	373
Zinc standard metal solution 1000 ppm	for ICP-MS solution in HNO ₃ 2%	100 mL	Plastic bottle	SPCLZN22Y	373

SPEX Standard mono-element solution for ICP & ICP-MS

designation	matrix	packsize	package type	Mfr. No	page
Zirconium standard metal solution 1000 ppm	for ICP solution in HNO ₃ 2%	100 mL	Plastic bottle	SPPLZR22Y	377
	for ICP solution in HNO ₃ 2%	250 mL	Plastic bottle	SPPLZR22T	377
	for ICP solution in HNO ₃ 2%	500 mL	Plastic bottle	SPPLZR22X	377
Zirconium standard metal solution 1000 ppm	for ICP solution in HCL 10%	500 mL	Plastic bottle	SPPLZR32X	377
Zirconium standard metal solution 10 000 ppm	for ICP solution in HNO ₃ 5%	100 mL	Plastic bottle	SPPLZR23Y	377
	for ICP solution in HNO ₃ 5%	500 mL	Plastic bottle	SPPLZR23X	377
Zirconium standard metal solution 10 000 ppm	for ICP solution in HCL 10%	500 mL	Plastic bottle	SPPLZR33X	377

SPEX Standard Instrument check, Calibration & Multi-element solution for ICP & ICP-MS

NIST Certified

designation	packsize	package type	Mfr. No	page
Standard Instrument Calibration 1	100 mL	Plastic bottle	SPCLCAL1	327
Standard Instrument Calibration 2	100 mL	Plastic bottle	SPCLCAL2	327
Standard Instrument Calibration 3	100 mL	Plastic bottle	SPCLCAL3	327
Standard Instrument check 3	100 mL	Plastic bottle	SPCALMIX3100	327
Standard Instrument check 10	500 mL	Plastic bottle	SPCALMIX10500	327
Standard Instrument check 3	500 mL	Plastic bottle	SPCALMIX3500	327
Standard Instrument check 4	100 mL	Plastic bottle	SPCALMIX4100	327
Standard Instrument check 4	500 mL	Plastic bottle	SPCALMIX4500	327
Standard Instrument check 7	100 mL	Plastic bottle	SPCALMIX7100	327
Standard Instrument check 7	500 mL	Plastic bottle	SPCALMIX7500	327
Standard Instrument check 8	100 mL	Plastic bottle	SPCALMIX8100	327
Standard Instrument check 8	500 mL	Plastic bottle	SPCALMIX8500	327
Standard Instrument check 10	100 mL	Plastic bottle	SPCALMIX10100	328
Standard 1 solution Multi-element	100 mL	Plastic bottle	SPCLMS1	328
Standard Set of 7 Multi-element solution (no Hg)	100 mL	Plastic bottle	SPCLMSSETN	328
Standard 2 solution Multi-element	100 mL	Plastic bottle	SPCLMS2	328
Standard 2A solution Multi-element	100 mL	Plastic bottle	SPCLMS2A	328
Standard 2A solution Multi-element (no Hg)	100 mL	Plastic bottle	SPCLMS2AN	328
Standard 2 solution Multi-element (no Hg)	100 mL	Plastic bottle	SPCLMS2N	328
Standard 3 solution Multi-element	100 mL	Plastic bottle	SPCLMS3	328
Standard 4 solution Multi-element	100 mL	Plastic bottle	SPCLMS4	328
Standard 5 solution Multi-element	100 mL	Plastic bottle	SPCLMS5	328
Standard Set of 7 Multi-element solution	100 mL	Plastic bottle	SPCLMSSET	328
Standard 1 Mixed Multi-element calibration	100 mL	Plastic bottle	SPMIXSTD1100	328
Standard 3 Mixed Multi-element calibration	500 mL	Plastic bottle	SPMIXSTD3500	329
Standard 3A Mixed Multi-element calibration	100 mL	Plastic bottle	SPMIXSTD3A100	329
Standard 3A Mixed Multi-element calibration	500 mL	Plastic bottle	SPMIXSTD3A500	329
Standard 4 Mixed Multi-element calibration	100 mL	Plastic bottle	SPMIXSTD4100	329
Standard 4 Mixed Multi-element calibration	500 mL	Plastic bottle	SPMIXSTD4500	329
Standard 4A Mixed Multi-element calibration	100 mL	Plastic bottle	SPMIXSTD4A100	329
Standard 4A Mixed Multi-element calibration	500 mL	Plastic bottle	SPMIXSTD4A500	329
Standard 5 Mixed Multi-element calibration	100 mL	Plastic bottle	SPMIXSTD5100	329
Standard 5 Mixed Multi-element calibration	500 mL	Plastic bottle	SPMIXSTD5500	329
Standard 5A Mixed Multi-element calibration	100 mL	Plastic bottle	SPMIXSTD5A100	329
Standard 1 Mixed Multi-element calibration	500 mL	Plastic bottle	SPMIXSTD1500	329
Standard 5A Mixed Multi-element calibration	500 mL	Plastic bottle	SPMIXSTD5A500	329
Standard 1A Mixed Multi-element calibration	100 mL	Plastic bottle	SPMIXSTD1A100	330
Standard 1A Mixed Multi-element calibration	500 mL	Plastic bottle	SPMIXSTD1A500	330
Standard 2 Mixed Multi-element calibration	100 mL	Plastic bottle	SPMIXSTD2100	330
Standard 2 Mixed Multi-element calibration	500 mL	Plastic bottle	SPMIXSTD2500	330
Standard 2A Mixed Multi-element calibration	100 mL	Plastic bottle	SPMIXSTD2A100	330
Standard 2A Mixed Multi-element calibration	500 mL	Plastic bottle	SPMIXSTD2A500	330
Standard 3 Mixed Multi-element calibration	100 mL	Plastic bottle	SPMIXSTD3100	330
Standard Plasma emission Kit 1 (38 elements)	100 mL	Plastic bottle	SPICPKIT1	330
Standard Plasma emission Kit 2 (13 elements)	100 mL	Plastic bottle	SPICPKIT2	330
Standard Plasma emission Kit 3 (10 elements)	100 mL	Plastic bottle	SPICPKIT3	330
Standard Plasma emission Kit 3 (15 elements)	100 mL	Plastic bottle	SPICPKIT4	330
Standard Tuning solution 1	100 mL	Plastic bottle	SPCLTUNE1	330

Standard solutions Ready to use for volumetric analysis

designation	comment	NIST	tolerance	packsize	package type	Mfr. No	page
Ammonium iron (II) sulfate solution 0.1M (0.1N)		●	± 0.005	1 L	Glass bottle	J/1101/15	43
		●	± 0.005	2.5 L	HDPE Plastic bottle	J/1101/PB17	43
Ammonium iron(II) sulfate solution 0.025M			± 0.005	2.5 L	Glass bottle	J/1100/17	43
Ammonium thiocyanate solution 0.1M (0.1N)	Eur.Ph. , USP, BP	●	± 0.002	1 L	HDPE Plastic bottle	J/1500/PB15	46
	Eur.Ph. , USP, BP	●	± 0.002	2.5 L	HDPE Plastic bottle	J/1500/PB17	46
Cerium(IV) sulfate solution 0.1M (0.1N)	Eur.Ph. , BP	●	± 0.002	1 L	Glass bottle	J/3000/15	81
	Eur.Ph. , BP	●	± 0.002	2.5 L	Glass bottle	J/3000/17	81
Ethylenediaminetetraacetic acid disodium salt solution 0.2N		●	± 0.003	1 L	HDPE Plastic bottle	J/3720/PB15	127
		●	± 0.003	2.5 L	HDPE Plastic bottle	J/3720/17	127
Ethylenediaminetetraacetic acid disodium salt solution 0.2N	Eur.Ph. , BP	●	± 0.001	1 L	HDPE Plastic bottle	J/3700/15	127
	Eur.Ph. , BP	●	± 0.001	2.5 L	HDPE Plastic bottle	J/3700/17	127
	Eur.Ph. , BP	●	± 0.001	10 L	Bag in box	J/3700/24	127
Hydrochloric acid solution 0.2M (0.2N)		●	± 0.001	1 L	HDPE Plastic bottle	J/4340/15	150
		●	± 0.001	2.5 L	HDPE Plastic bottle	J/4340/17	150
Hydrochloric acid solution 0.1M (0.1N)	Eur.Ph. , BP	●	± 0.001	1 L	HDPE Plastic bottle	J/4350/15	150
	Eur.Ph. , BP	●	± 0.001	2.5 L	HDPE Plastic bottle	J/4350/17	150
	Eur.Ph. , BP	●	± 0.001	5 L	HDPE Plastic bottle	J/4350/21	150
	Eur.Ph. , BP	●	± 0.001	10 L	Bag in box	J/4350/24	150
Hydrochloric acid solution 0.5M (0.5N)		●	± 0.001	1 L	HDPE Plastic bottle	J/4330/15	150
		●	± 0.001	2.5 L	HDPE Plastic bottle	J/4330/17	150
		●	± 0.001	25 L	HDPE Plastic container	J/4330/25	150
Hydrochloric acid solution 0.02M (0.02N)		●	± 0.003	1 L	HDPE Plastic bottle	J/4360/PB15	151
		●	± 0.003	2.5 L	HDPE Plastic bottle	J/4360/17	151
Hydrochloric acid solution 1M (1N)	Eur.Ph. , USP, BP	●	± 0.001	1 L	HDPE Plastic bottle	J/4320/15	151
	Eur.Ph. , USP, BP	●	± 0.001	2.5 L	HDPE Plastic bottle	J/4320/17	151
	Eur.Ph. , USP, BP	●	± 0.001	5 L	HDPE Plastic bottle	J/4320/21	151
	Eur.Ph. , USP, BP	●	± 0.001	10 L	Bag in box	J/4320/24	151
Hydrochloric acid solution 5M (5N)			± 0.005	1 L	HDPE Plastic bottle	J/4310/15	151
			± 0.005	2.5 L	HDPE Plastic bottle	J/4310/17	151
Hydrochloric acid solution 2M (2N)	Eur.Ph. , BP	●	± 0.001	1 L	HDPE Plastic bottle	J/4315/15	151
	Eur.Ph. , BP	●	± 0.001	2.5 L	HDPE Plastic bottle	J/4315/17	151
Iodine solution 0.05M (0.1N)	Eur.Ph. , USP, BP	●	± 0.001	1 L	Glass bottle	J/4410/15	161
	Eur.Ph. , USP, BP	●	± 0.001	2.5 L	Glass bottle	J/4410/17	161
Iodine solution 0.5M (1N)	Eur.Ph. , BP	●	± 0.002	500 mL	Coated Glass bottle	J/4400/PB08	161
	Eur.Ph. , BP	●	± 0.002	1 L	Coated Glass bottle	J/4400/PB15	161
	Eur.Ph. , BP	●	± 0.002	2.5 L	Coated Glass bottle	J/4400/PB17	161
Nitric acid solution 2M (2N)		●	± 0.005	2.5 L	HDPE Plastic bottle	J/5545/PB17	218
Nitric acid solution 1M (1N)	Eur.Ph. , USP, BP	●	± 0.001	1 L	HDPE Plastic bottle	J/5550/PB15	218
	Eur.Ph. , USP, BP	●	± 0.001	2.5 L	HDPE Plastic bottle	J/5550/PB17	218
	Eur.Ph. , USP, BP	●	± 0.001	10 L	Bag in box	J/5550/24	218
Nitric acid solution 0.1M (0.1N)		●	± 0.001	1 L	HDPE Plastic bottle	J/5600/PB15	218
		●	± 0.001	2.5 L	HDPE Plastic bottle	J/5600/17	218
		●	± 0.001	10 L	Bag in box	J/5600/24	218
Perchloric acid solution 0.1M (0.1N)		●	± 0.002	1 L	HDPE Plastic bottle	J/5820/PB15	229
		●	± 0.002	2.5 L	HDPE Plastic bottle	J/5820/PB17	229
Potassium bromate-bromide solution 0.0167 M (0.1N)		●	± 0.001	1 L	Glass bottle	J/6100/15	247
		●	± 0.001	2.5 L	Glass bottle	J/6100/17	247
Potassium dichromate solution 0.0167M (0.1N)	Eur.Ph. , USP, BP	●	± 0.001	1 L	Glass bottle	J/6212/15	251
	Eur.Ph. , USP, BP	●	± 0.001	2.5 L	Glass bottle	J/6212/17	251
Potassium hydroxide solution 0.1M (0.1N), (in methanol)	Eur.Ph. , BP	●	± 0.005	1 L	HDPE Plastic bottle	J/6625/15	256

Standard solutions Ready to use for volumetric analysis

designation	comment	NIST	tolerance	packsize	package type	Mfr. No	page
Potassium hydroxide solution 0.1M (0.1N), (in ethanol)		●	± 0.002	1 L	HDPE Plastic bottle	J/6620/15	256
		●	± 0.002	2.5 L	HDPE Plastic bottle	J/6620/17	256
Potassium hydroxide solution 0.5M (0.5N), (in ethanol)		●	± 0.002	1 L	HDPE Plastic bottle	J/6612/15	257
		●	± 0.002	2.5 L	HDPE Plastic bottle	J/6612/17	257
Potassium hydroxide solution 1M (1N)	Eur.Ph. , BP	●	± 0.001	1 L	HDPE Plastic bottle	J/6630/PB15	257
	Eur.Ph. , BP	●	± 0.001	2.5 L	HDPE Plastic bottle	J/6630/17	257
Potassium hydroxide solution 1M (1N), (alcoholic)		●	± 0.001	1 L	HDPE Plastic bottle	J/6610/15	257
		●	± 0.001	2.5 L	HDPE Plastic bottle	J/6610/17	257
Potassium iodate 0.025M solution stabilised with 10 ppm mercury(II) chloride		●	± 0.001	2.5 L	Glass bottle	J/6700/17	258
Potassium iodate 0.05M solution stabilised with 10 ppm mercury(II) chloride	Eur.Ph. , USP, BP	●	± 0.001	1 L	Glass bottle	J/6690/15	258
	Eur.Ph. , USP, BP	●	± 0.001	2.5 L	Glass bottle	J/6690/17	258
Potassium permanganate solution 0.02M (0.1N)		●	± 0.001	1 L	Glass bottle	J/6900/15	261
		●	± 0.001	2.5 L	Coated Glass bottle	J/6900/PB17	261
Potassium thiocyanate solution 0.1M (0.1N)		●	± 0.002	1 L	Glass bottle	J/7050/15	264
		●	± 0.002	2.5 L	HDPE Plastic bottle	J/7050/PB17	264
Silver nitrate solution 1M (1N)			± 0.001	1 L	HDPE Plastic bottle	J/7320/PB15	277
			± 0.001	2.5 L	HDPE Plastic bottle	J/7320/PB17	277
Silver nitrate solution 0.1M (0.1N)	Eur.Ph. , USP, BP	●	± 0.001	1 L	HDPE Plastic bottle	J/7330/15	277
	Eur.Ph. , USP, BP	●	± 0.001	2.5 L	HDPE Plastic bottle	J/7330/17	277
	Eur.Ph. , USP, BP	●	± 0.001	10 L	Bag in box	J/7330/24	277
	Eur.Ph. , USP, BP	●	± 0.001	25 L	HDPE Drum	J/7330/25	277
Silver nitrate solution 0.05M (0.05N)		●	± 0.002	1 L	HDPE Plastic bottle	J/7340/15	278
		●	± 0.002	2.5 L	HDPE Plastic bottle	J/7340/17	278
Silver nitrate solution 0.0282M (0.0282N)		●	± 0.003	250 mL	Glass bottle	J/7370/07	278
		●	± 0.003	1 L	HDPE Plastic bottle	J/7370/15	278
		●	± 0.003	2.5 L	HDPE Plastic bottle	J/7370/17	278
Silver nitrate solution 0.02M (0.02N)		●	± 0.003	1 L	HDPE Plastic bottle	J/7350/15	278
		●	± 0.003	2.5 L	HDPE Plastic bottle	J/7350/17	278
Sodium arsenite solution 0.05M (0.1N)		●	± 0.002	1 L	HDPE Plastic bottle	J/7510/PB15	282
		●	± 0.002	2.5 L	HDPE Plastic bottle	J/7510/PB17	282
Sodium carbonate solution 0.05M (0.1N)		●	± 0.001	1 L	HDPE Plastic bottle	J/7580/PB15	284
		●	± 0.001	2.5 L	HDPE Plastic bottle	J/7580/17	284
Sodium carbonate solution 0.5M (1N)		●	± 0.001	1 L	HDPE Plastic bottle	J/7520/PB15	285
		●	± 0.001	2.5 L	HDPE Plastic bottle	J/7520/17	285
Sodium hydroxide solution 10M (10N) ca. 40w/v			± 0.005	1 L	HDPE Plastic bottle	J/7800/15	296
			± 0.005	2.5 L	HDPE Plastic bottle	J/7800/17	296
			± 0.005	5 L	HDPE Plastic bottle	J/7800/21	296
			± 0.005	10 L	Bag in box	J/7800/24	296
			± 0.005	25 L	HDPE Plastic container	J/7800/25	296
Sodium hydroxide solution 2M (2N)	BP	●	± 0.001	1 L	HDPE Plastic bottle	J/7610/PB15	296
	BP	●	± 0.001	2.5 L	HDPE Plastic bottle	J/7610/17	296
	BP	●	± 0.001	5 L	HDPE Plastic bottle	J/7610/21	296
Sodium hydroxide solution 1M (1N)	Eur.Ph. , USP, BP	●	± 0.001	1 L	HDPE Plastic bottle	J/7620/15	296
	Eur.Ph. , USP, BP	●	± 0.001	2.5 L	HDPE Plastic bottle	J/7620/17	296
	Eur.Ph. , USP, BP	●	± 0.001	5 L	HDPE Plastic bottle	J/7620/21	296
	Eur.Ph. , USP, BP	●	± 0.001	10 L	Bag in box	J/7620/24	296

Standard solutions Ready to use for volumetric analysis

designation	comment	NIST	tolerance	packsize	package type	Mfr. No	page
Sodium hydroxide solution methanolic 1M (1N)		●	± 0.002	1 L	HDPE Plastic bottle	J/7625/PB15	296
		●	± 0.002	2.5 L	HDPE Plastic bottle	J/7625/17	296
Sodium hydroxide solution 0.5M (0.5N)		●	± 0.001	1 L	HDPE Plastic bottle	J/7630/15	296
		●	± 0.001	2.5 L	HDPE Plastic bottle	J/7630/17	296
		●	± 0.001	5 L	HDPE Plastic bottle	J/7630/21	296
		●	± 0.001	10 L	Bag in box	J/7630/24	296
		●	± 0.001	200 L	HDPE Plastic container	J/7630/27	296
Sodium hydroxide solution 0.2M (0.2N)		●	± 0.001	1 L	HDPE Plastic bottle	J/7640/15	297
		●	± 0.001	2.5 L	HDPE Plastic bottle	J/7640/17	297
		●	± 0.001	5 L	HDPE Plastic bottle	J/7640/21	297
		●	± 0.001	10 L	Bag in box	J/7640/24	297
Sodium hydroxide solution 0.1M (0.1N)	Eur.Ph. , BP	●	± 0.001	1 L	HDPE Plastic bottle	J/7660/15	297
	Eur.Ph. , BP	●	± 0.001	2.5 L	HDPE Plastic bottle	J/7660/17	297
	Eur.Ph. , BP	●	± 0.001	5 L	HDPE Plastic bottle	J/7660/21	297
	Eur.Ph. , BP	●	± 0.001	10 L	Bag in box	J/7660/24	297
Sodium hydroxide solution (0.01M) 0.01N		●	± 0.001	2.5 L	HDPE Plastic bottle	J/9010/17	297
Sodium nitrite solution 0.5M (0.5N)		●	± 0.003	1 L	Glass bottle	J/7900/15	300
Sodium thiosulfate solution 0.1M (0.1N)	Eur.Ph. , USP, BP	●	± 0.001	1 L	HDPE Plastic bottle	J/7950/15	306
	Eur.Ph. , USP, BP	●	± 0.001	2.5 L	HDPE Plastic bottle	J/7950/17	306
	Eur.Ph. , USP, BP	●	± 0.001	5 L	HDPE Plastic bottle	J/7950/21	306
	Eur.Ph. , USP, BP	●	± 0.001	10 L	Bag in box	J/7950/24	306
Sodium thiosulfate solution 0.01M (0.01N)	Eur.Ph, BP	●	± 0.001	2.5 L	HDPE Plastic bottle	J/9011/17	306
Sulfuric acid solution 1M (2N)		●	± 0.001	1 L	HDPE Plastic bottle	J/8420/PB15	339
		●	± 0.001	2.5 L	HDPE Plastic bottle	J/8420/17	339
		●	± 0.001	10 L	Bag in box	J/8420/24	339
		●	± 0.001	25 L	HDPE Drum	J/8420/25	339
Sulfuric acid solution 0.5M (1N)	Eur.Ph. , USP, BP	●	± 0.001	1 L	HDPE Plastic bottle	J/8430/15	339
	Eur.Ph. , USP, BP	●	± 0.001	2.5 L	HDPE Plastic bottle	J/8430/17	339
	Eur.Ph. , USP, BP	●	± 0.001	10 L	Bag in box	J/8430/24	339
	Eur.Ph. , USP, BP	●	± 0.001	25 L	HDPE Drum	J/8430/25	339
Sulfuric acid solution 0.25M (0.5N)		●	± 0.001	1 L	HDPE Plastic bottle	J/8440/PB15	339
		●	± 0.001	2.5 L	HDPE Plastic bottle	J/8440/17	339
Sulfuric acid solution 0.1M (0.2N)		●	± 0.001	1 L	HDPE Plastic bottle	J/8450/15	339
		●	± 0.001	2.5 L	HDPE Plastic bottle	J/8450/17	339
		●	± 0.001	10 L	Bag in box	J/8450/24	339
		●	± 0.001	25 L	HDPE Drum	J/8450/25	339
Sulfuric acid solution 0.05M (0.1N)	Eur.Ph. , BP	●	± 0.001	1 L	HDPE Plastic bottle	J/8460/15	340
	Eur.Ph. , BP	●	± 0.001	2.5 L	HDPE Plastic bottle	J/8460/17	340
	Eur.Ph. , BP	●	± 0.001	10 L	Bag in box	J/8460/24	340
	Eur.Ph. , BP	●	± 0.001	25 L	HDPE Drum	J/8460/25	340
Sulfuric acid solution 0.01M (0.2N)		●	± 0.003	1 L	HDPE Plastic bottle	J/8470/PB15	340
		●	± 0.003	2.5 L	HDPE Plastic bottle	J/8470/17	340
		●	± 0.003	10 L	Bag in box	J/8470/24	340
Tetra-n-butyl ammonium hydroxide solution 0.1 M (0.1N) in toluene/methanol			± 0.05	500 mL	Glass bottle	J/8750/08	344
			± 0.05	1 L	Glass bottle	J/8750/15	344
			± 0.05	2.5 L	Glass bottle	J/8750/17	344
			± 0.05	10 L	Bag in box	J/8750/24	344
Titanium(III) chloride solution 0.2M (0.2N)			± 0.1	2.5 L	Glass bottle	J/8780/17	353

Standard solutions concentrated for volumetric analysis

designation	each ampoule ma- kes	packsize		Mfr. No	page
Acetic acid solution 0.1M (0.1N)	1ℓ	1 Ampoule		J/0551C/05	15
		6 Ampoules		J/0551C/90	15
Ammonium thiocyanate solution 0.1M (0.1N)	1ℓ	1 Ampoule		J/1500C/05	46
		6 Ampoules		J/1500C/90	46
Cerium(IV) sulfate solution 0.05M (0.05N)		6 Ampoules		J/3010C/90	81
Ethylenediaminetetraacetic acid disodium salt solution 0.01M (0.02N)	1ℓ	1 Ampoule		J/3720C/05	127
		6 Ampoule		J/3720C/90	127
Ethylenediaminetetraacetic acid trisodium salt solution 0.1M (0.2N)	500 mℓ	1 Ampoule		J/3700C/05	127
		6 Ampoules		J/3700C/90	127
Hydrochloric acid solution 0.1M (0.1N)	1ℓ	1 Ampoule		J/4350C/05	150
		6 Ampoules		J/4350C/90	150
Hydrochloric acid solution 0.5M (0.5N)	1ℓ	1 Ampoule		J/4330C/05	151
		6 Ampoules		J/4330C/90	151
Hydrochloric acid solution 1M (1N)	1ℓ	1 Ampoule		J/4320C/05	151
		6 Ampoules		J/4320C/90	151
Nitric acid solution 1M (1N)	1ℓ	1 Ampoule		J/5550C/05	218
		6 Ampoules		J/5550C/90	218
Oxalic acid solution 0.05M (0.1N)	1ℓ	1 Ampoule		J/5650C/05	223
		6 Ampoules		J/5650C/90	223
Potassium bromate-bromide solution 0.0167 M (0.1N)	1ℓ	1 Ampoule		J/6100C/05	247
		6 Ampoules		J/6100C/90	247
Potassium dichromate solution 0.0167M (0.1N)	1ℓ	1 Ampoule		J/6212C/05	252
		6 Ampoules		J/6212C/90	252
Potassium hydroxide solution 0.1M (0.1N)	1ℓ	1 Ampoule		J/6640C/05	257
		6 Ampoules		J/6640C/90	257
Potassium hydroxide solution 1M (1N)	1ℓ	1 Ampoule		J/6630C/05	257
		6 Ampoules		J/6630C/90	257
Silver nitrate solution 0.1M (0.1N)	1ℓ	1 Ampoule		J/7330C/05	278
		6 Ampoules		J/7330C/90	278
Sodium acetate solution 0.1M	1 ℓ	1 Ampoule		J/7440C/05	281
Sodium carbonate solution 0.05M (0.1N)	1ℓ	1 Ampoule		J/7580C/05	285
		6 Ampoules		J/7580C/90	285
Sodium chloride solution 0.1M (0.1N)	1ℓ	1 Ampoule		J/7590C/05	286
		6 Ampoules		J/7590C/90	286
Sodium hydroxide solution 1M (1N)	1ℓ	1 Ampoule		J/7620C/05	296
		6 Ampoules		J/7620C/90	296
Sodium hydroxide solution 0.5M (0.5N)	1ℓ	1 Ampoule		J/7630C/05	296
		6 Ampoules		J/7630C/90	296
Sodium hydroxide solution 0.1M (0.1N)	1ℓ	1 Ampoule		J/7660C/05	297
		6 Ampoules		J/7660C/90	297
Sodium thiosulfate solution 0.1M (0.1N)	1ℓ	1 Ampoule		J/7950C/05	306
		6 Ampoules		J/7950C/90	306
Sodium thiosulfate solution 0.5M (0.5N)	500 mℓ	1 Ampoule		J/7951C/05	306
		6 Ampoules		J/7951C/90	306
Sulfuric acid solution 0.05M (0.1N)	1ℓ	1 Ampoule		J/8460C/05	340
		6 Ampoules		J/8460C/90	340
Sulfuric acid solution 0.5M (1N)	1ℓ	1 Ampoule		J/8430C/05	340
		6 Ampoules		J/8430C/90	340
Zinc sulfate solution 0.05M (0.05N)	1ℓ	1 Ampoule		J/8490C/05	377
		6 Ampoules		J/8490C/90	377

Standard Instrument Calibration 1
for ICP-MS, NIST certified

packaging	Mfr. No
100 mℓ Plastic bottle	SPCLCAL1
H314	

see the complete information (Application, Content & Matrix) in the Standard solution section page 322
Check with your local distributor for availability of Spex product

Standard Instrument Calibration 2
for ICP-MS, NIST certified

packaging	Mfr. No
100 mℓ Plastic bottle	SPCLCAL2
H314	

see the complete information (Application, Content & Matrix) in the Standard solution section page 322
Check with your local distributor for availability of Spex product

Standard Instrument Calibration 3
for ICP-MS, NIST certified

packaging	Mfr. No
100 mℓ Plastic bottle	SPCLCAL3
H314	

see the complete information (Application, Content & Matrix) in the Standard solution section page 322
Check with your local distributor for availability of Spex product

Standard Instrument check 3
for ICP, NIST certified

packaging	Mfr. No
100 mℓ Plastic bottle	SPCALMIX3100
H314, H318	

see the complete information (Application, Content & Matrix) in the Standard solution section page 322
Check with your local distributor for availability of Spex product

Standard Instrument check 10
for ICP, NIST certified

packaging	Mfr. No
500 mℓ Plastic bottle	SPCALMIX10500
H315, H319	

see the complete information (Application, Content & Matrix) in the Standard solution section page 322
Check with your local distributor for availability of Spex product

Standard Instrument check 3
for ICP, NIST certified

packaging	Mfr. No
500 mℓ Plastic bottle	SPCALMIX3500
H314, H318	

see the complete information (Application, Content & Matrix) in the Standard solution section page 322
Check with your local distributor for availability of Spex product

Standard Instrument check 4
for ICP, NIST certified

packaging	Mfr. No
100 mℓ Plastic bottle	SPCALMIX4100
H315, H320	

see the complete information (Application, Content & Matrix) in the Standard solution section page 322
Check with your local distributor for availability of Spex product

Standard Instrument check 4
for ICP, NIST certified

packaging	Mfr. No
500 mℓ Plastic bottle	SPCALMIX4500
H315, H320	

see the complete information (Application, Content & Matrix) in the Standard solution section page 322
Check with your local distributor for availability of Spex product

Standard Instrument check 7
for ICP, NIST certified

packaging	Mfr. No
100 mℓ Plastic bottle	SPCALMIX7100
H315, H319	

see the complete information (Application, Content & Matrix) in the Standard solution section page 322
Check with your local distributor for availability of Spex product

Standard Instrument check 7
for ICP, NIST certified

packaging	Mfr. No
500 mℓ Plastic bottle	SPCALMIX7500
H315, H319	

see the complete information (Application, Content & Matrix) in the Standard solution section page 322
Check with your local distributor for availability of Spex product

Standard Instrument check 8
for ICP, NIST certified

packaging	Mfr. No
100 mℓ Plastic bottle	SPCALMIX8100
H315, H320	

see the complete information (Application, Content & Matrix) in the Standard solution section page 322
Check with your local distributor for availability of Spex product

Standard Instrument check 8
for ICP, NIST certified

packaging	Mfr. No
500 mℓ Plastic bottle	SPCALMIX8500
H315, H320	

see the complete information (Application, Content & Matrix) in the Standard solution section page 322
Check with your local distributor for availability of Spex product

Standard 4A Mixed Multi-element calibration

for ICP, NIST certified

packaging	Mfr. No.
500 mL Plastic bottle	SPMIXSTD4A500
H314	

see the complete information (Application, Content & Matrix) in the Standard solution section page 322
Check with your local distributor for availability of Spex product

Standard 5 Mixed Multi-element calibration

for ICP, NIST certified

packaging	Mfr. No.
100 mL Plastic bottle	SPMIXSTD5100
H314	

see the complete information (Application, Content & Matrix) in the Standard solution section page 322
Check with your local distributor for availability of Spex product

Standard 5 Mixed Multi-element calibration

for ICP, NIST certified

500 ml	Plastic bottle	SPMIXSTD5500
H314		

see the complete information (Application, Content & Matrix) in the Standard solution section page 322
Check with your local distributor for availability of Spex product

Standard 5A Mixed Multi-element calibration

for ICP, NIST certified

for ICP, NIST certified	
packaging	Mfr. No.
100 mL Plastic bottle	SPMIXSTD5A100

see the complete information (Application, Content & Matrix) in the Standard solution section page 322
Check with your local distributor for availability of Spex product

Standard 1 Mixed Multi-element calibration

for ICP, NIST certified

packaging	Mfr. No
500 mL Plastic bottle	SPMIXST1500

see the complete information (Application, Content & Matrix) in the Standard solution section page 322
Check with your local distributor for availability of Spex product

Standard 5A Mixed Multi-element calibration

for ICP, NIST certified

packaging	Mfr. No
500 mL Plastic bottle	SPMIXSTD5A500
H314	

see the complete information (Application, Content & Matrix) in the Standard solution section page 322

Check with your local distributor for availability of Spex product

Standard 1A Mixed Multi-element calibration
for ICP, NIST certified

packaging	Mfr. No
100 mℓ Plastic bottle	SPMIXSTD1A100
H314	



see the complete information (Application, Content & Matrix) in the Standard solution section page 322

Check with your local distributor for availability of Spex product

Standard 1A Mixed Multi-element calibration
for ICP, NIST certified

packaging	Mfr. No
500 mℓ Plastic bottle	SPMIXSTD1A500
H314	



see the complete information (Application, Content & Matrix) in the Standard solution section page 322

Check with your local distributor for availability of Spex product

Standard 2 Mixed Multi-element calibration
for ICP, NIST certified

packaging	Mfr. No
100 mℓ Plastic bottle	SPMIXSTD2100
H314	



see the complete information (Application, Content & Matrix) in the Standard solution section page 322

Check with your local distributor for availability of Spex product

Standard 2 Mixed Multi-element calibration
for ICP, NIST certified

packaging	Mfr. No
500 mℓ Plastic bottle	SPMIXSTD2500
H314	



see the complete information (Application, Content & Matrix) in the Standard solution section page 322

Check with your local distributor for availability of Spex product

Standard 2A Mixed Multi-element calibration
for ICP, NIST certified

packaging	Mfr. No
100 mℓ Plastic bottle	SPMIXSTD2A100
H314	



see the complete information (Application, Content & Matrix) in the Standard solution section page 322

Check with your local distributor for availability of Spex product

Standard 2A Mixed Multi-element calibration
for ICP, NIST certified

packaging	Mfr. No
500 mℓ Plastic bottle	SPMIXSTD2A500
H314	



see the complete information (Application, Content & Matrix) in the Standard solution section page 322

Check with your local distributor for availability of Spex product

Standard 3 Mixed Multi-element calibration
for ICP, NIST certified

packaging	Mfr. No
100 mℓ Plastic bottle	SPMIXSTD3100
H315, H319	



see the complete information (Application, Content & Matrix) in the Standard solution section page 322

Check with your local distributor for availability of Spex product

Standard Plasma emission Kit 1 (38 elements)
for ICP, NIST certified

packaging	Mfr. No
100 mℓ Plastic bottle	SPICPKIT1

see the complete information (Application, Content & Matrix) in the Standard solution section page 322

Check with your local distributor for availability of Spex product

Standard Plasma emission Kit 2 (13 elements)
for ICP, NIST certified

packaging	Mfr. No
100 mℓ Plastic bottle	SPICPKIT2

see the complete information (Application, Content & Matrix) in the Standard solution section page 322

Check with your local distributor for availability of Spex product

Standard Plasma emission Kit 3 (10 elements)
for ICP, NIST certified

packaging	Mfr. No
100 mℓ Plastic bottle	SPICPKIT3

see the complete information (Application, Content & Matrix) in the Standard solution section page 322

Check with your local distributor for availability of Spex product

Standard Plasma emission Kit 3 (15 elements)
for ICP, NIST certified

packaging	Mfr. No
100 mℓ Plastic bottle	SPICPKIT4

see the complete information (Application, Content & Matrix) in the Standard solution section page 322

Check with your local distributor for availability of Spex product

Standard Tuning solution 1
for ICP-MS, NIST certified

packaging	Mfr. No
100 mℓ Plastic bottle	SPCLTUNE1
H315, H319	



see the complete information (Application, Content & Matrix) in the Standard solution section page 322

Check with your local distributor for availability of Spex product

Starch
for analysis, soluble

Certified AR

packaging	Mfr. No
250 g Plastic jar	S/8000/50
500 g Plastic jar	S/8000/53
1 kg Plastic jar	S/8000/60

(C₆H₁₀O₅)_n
CAS: 9005-25-8
Density at 20°C: 1.5
EINECS: 232-679-6

Guaranteed analysis

Loss on drying	<15%	Sensitivity to iodine	: Pass
Reducing substances	: <0.2%	Acidity	<0.01 meq/g
Sulfated ash	<0.5%	Total chloride (Cl)	<0.005%

Lot analysis on the label.

Starch, soluble

SLR

specified reagents for general laboratory work- extra pure

packaging	Mfr. No
250 g Plastic jar	S/7960/50
500 g Plastic jar	S/7960/53
1 kg Plastic jar	S/7960/60

(C₆H₁₀O₅)_n
CAS: 9005-25-8
Density at 20°C: 1.5
EINECS: 232-679-6

Guaranteed analysis

Identity test	: Pass
---------------	--------

Starch, maize

SLR

specified reagents for general laboratory work- extra pure

packaging	Mfr. No
1 kg Plastic jar	S/7880/60

(C₆H₁₀O₅)_n
CAS: 9005-25-8
Density at 20°C: 1.5
EINECS: 232-679-6

Guaranteed analysis

Identity test	: Pass
---------------	--------

Starch, potato

SLR

specified reagents for general laboratory work- extra pure

packaging	Mfr. No
1 kg Plastic jar	S/7920/60

(C₆H₁₀O₅)_n
CAS: 9005-25-8
Density at 20°C: 1.5
EINECS: 232-679-6

Guaranteed analysis

Loss on drying	<20%	Residue after ignition	<0.5%
----------------	------	------------------------	-------

Starch, wheat

SLR

specified reagents for general laboratory work- extra pure

packaging	Mfr. No
1 kg Plastic jar	S/8040/60

(C₆H₁₀O₅)_n
CAS: 9005-25-8
Density at 20°C: 1.5
EINECS: 232-679-6

Guaranteed analysis

Residue after ignition	<0.3%	Lead (Pb)	<10 ppm
------------------------	-------	-----------	---------

Stearic acid

SLR

specified reagents for general laboratory work- extra pure

packaging	Mfr. No
100 g Plastic jar	S/8130/48
500 g Plastic jar	S/8130/53

C₁₈H₃₆O₂
CAS: 57-11-4
MW: 284.47
Density at 20°C: 0.840
Boiling point: 361°C
Melting point: 67 to 69°C
Flash point: 196°C
ADR/RID: 4.1, III

EINECS: 200-313-4
UN: 1325
GHS: Flam. Sol. 2, Skin Irrit. 2, Eye Irrit. 2, STOT SE 3
H228, H319, H335, H315
P261, P302+P352, P280, P305+P351+P338, P233, P210



Assay: >97%

The product may contain other chain length fatty acids.

Stearic anhydride

SLR

specified reagents for general laboratory work- extra pure

packaging	Mfr. No
100 g Plastic jar	S/8144/48

C₃₆H₇₀O₃
CAS: 638-08-4
MW: 550.95
Density at 20°C: 0.840
Melting point: 69 to 71°C
ADR/RID: 8, III

EINECS: 211-318-6
UN: 3261
GHS: Skin Corr. 1B, Eye Dam. 1
H314
P301+P330+P331, P280, P305+P351+P338, P310



Assay: >97%

In determining the hydroxyl groups in the alcohol and phenols.
Sully, Analyst, 1962.87.

Strontium standard metal solution 1000 ppm
for AAS, solution HNO₃ 1M

packaging	Mfr. No
100 mℓ HDPE Plastic bottle	J/8166/05
500 mℓ HDPE Plastic bottle	J/8166/08

ADR/RID: 8, II
UN: 2031
GHS: Skin Corr. 1B, Eye Dam. 1

H314
P280, P305+P351+P338, P301+P330+P331, P310



Manufactured to a tolerance of ± 0.5%.

See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.

Strontium standard metal solution 1000 ppm
for ICP, solution HNO₃ 1M

packaging	Mfr. No
500 mℓ HDPE Plastic bottle	J/8321/08

ADR/RID: 8, II
UN: 2031
GHS: Skin Corr. 1B, Eye Dam. 1

H314
P280, P305+P351+P338, P301+P330+P331, P310



Manufactured to a tolerance of ± 0.5%.

All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.

See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.

All ICP standards are manufactured from acids, metal & metal salts Ultra pure quality and contain less than <10 ppm total metal impurities.
see also the full range of "Standard solution" in the Standard solution section.

UHPLC

Filtered to 0.1 µm, Fisher Chemical UHPLC gradient grade solvents (Acetonitrile, Methanol and Water) are perfect for those working at high pressure.

See page 309 for more information



Strontium standard metal solution 10 000 ppm

for ICP, solution HNO₃ 1M

packaging	Mfr. No
100 mℓ HDPE Plastic bottle	J/8320/05
ADR/RID: 8, II	H314
UN: 2031	P280, P305+P351+P338,
GHS: Skin Corr. 1B, Eye Dam. 1	P301+P330+P331, P310

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.

Strontium standard metal solution 1000 ppm

for ICP solution in HNO₃ 2%, NIST certified

packaging	Mfr. No
100 mℓ Plastic bottle	SPPLSR22Y
250 mℓ Plastic bottle	SPPLSR22T
500 mℓ Plastic bottle	SPPLSR22X
H315, H320	

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Strontium standard metal solution 1000 ppm

for ICP solution in HCL 2%, NIST certified

packaging	Mfr. No
500 mℓ Plastic bottle	SPPLSR12X

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Strontium standard metal solution 10 000 ppm

for ICP solution in HNO₃ 5%, NIST certified

packaging	Mfr. No
100 mℓ Plastic bottle	SPPLSR23Y
500 mℓ Plastic bottle	SPPLSR23X
H314	

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Strontium standard metal solution 10 000 ppm

for ICP solution in HCL 5%, NIST certified

packaging	Mfr. No
500 mℓ Plastic bottle	PLSR13X

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Strontium carbonate

specified reagents for general laboratory work- extra pure

SLR

packaging	Mfr. No
100 g Plastic jar	S/8160/48
500 g Plastic jar	S/8160/53
SrCO ₃	Density at 20°C: 3.7
CAS: 1633-05-2	EINECS: 216-643-7
MW: 147.63	

Assay: >98%

Guaranteed analysis

Iron (Fe)	<20 ppm	n-Butyl formate	<0.02%
-----------	---------	-----------------	--------

Strontium chloride

for analysis

Certified AR

packaging	Mfr. No
250 g Plastic jar	S/8320/50
500 g Plastic jar	S/8320/53
Cl ₂ Sr.6H ₂ O	GHS: Skin Irrit. 2, Eye Dam.
CAS: 10025-70-4	1, STOT SE 3
MW: 266.62	H315, H335, H318
Density at 20°C: 1.9	P280, P305+P351+P338
Melting point: 115°C	

Assay: 99 to 102%

Guaranteed analysis

Acidity	<0.0005 meq/g	Potassium (K)	<100 ppm
Calcium (Ca)	<500 ppm	Sodium (Na)	<200 ppm
Copper (Cu)	<2 ppm	Total nitrogen (N)	<50 ppm
Iron (Fe)	<5 ppm	Total silicon (Si)	<50 ppm
Lead (Pb)	<2 ppm	Total sulphur (S)	<20 ppm
Magnesium (Mg)	<50 ppm	Zinc (Zn)	<20 ppm

Typical analysis

Aluminium (Al)	<5 ppm	Manganese (Mn)	<5 ppm
Barium (Ba)	<0.002%	Molybdenum (Mo)	<5 ppm
Cadmium (Cd)	<5 ppm	Nickel (Ni)	<5 ppm
Chromium (Cr)	<5 ppm	Silver (Ag)	<5 ppm
Cobalt (Co)	<5 ppm	Titanium (Ti)	<5 ppm
Lithium (Li)	<5 ppm	Vanadium (V)	<5 ppm

Lot analysis on the label

Strontium chloride hexahydrate

specified reagents for general laboratory work- extra pure

SLR

packaging	Mfr. No
250 g Plastic jar	S/8280/50
500 g Plastic jar	S/8280/53
Cl ₂ Sr.6H ₂ O	EINECS: 233-971-6
CAS: 10025-70-4	GHS: Skin Irrit. 2, Eye Dam.
MW: 266.62	1, STOT SE 3
Density at 20°C: 1.92	H315, H335, H318
Melting point: 115°C	P280, P305+P351+P338

Assay: >98%

Guaranteed analysis

Calcium (Ca)	<1000 ppm	Potassium (K)	<200 ppm
Copper (Cu)	<10 ppm	Sodium (Na)	<500 ppm
Iron (Fe)	<20 ppm	Total silicon (Si)	<200 ppm
Lead (Pb)	<10 ppm	Total sulphur (S)	<100 ppm
Magnesium (Mg)	<200 ppm	Zinc (Zn)	<100 ppm

UHPLC
Filtered to 0.1 µm, Fisher Chemical UHPLC gradient grade solvents (Acetonitrile, Methanol and Water) are perfect for those working at high pressure.
See page 309 for more information

Strontium nitrate anhydrous

for analysis

Certified AR

packaging	Mfr. No
500 g Plastic jar	S/8402/53
N ₂ O ₆ Sr	UN: 1507
CAS: 10042-76-9	GHS: Ox. Sol. 3, Acute Tox.
MW: 211.63	4, Skin Irrit. 2, Eye Irrit. 2,
Boiling point: 645°C	STOT SE 3
Melting point: 570°C	H272, H315, H335, H319, H302
ADR/RID: 5.1, III	P220, P261, P280,
EINECS: 233-131-9	P305+P351+P338

Assay: >99.5%

Guaranteed analysis

pH (aqueous solution 5%)	5 to 7	Manganese (Mn)	<5 ppm
Aluminium (Al)	<5 ppm	Molybdenum (Mo)	<5 ppm
Barium (Ba)	<0.2%	Nickel (Ni)	<5 ppm
Cadmium (Cd)	<5 ppm	Potassium (K)	<54 ppm
Calcium (Ca)	<500 ppm	Propan--1-ol	<5 ppm
Chromium (Cr)	<5 ppm	Sodium (Na)	<50 ppm
Cobalt (Co)	<5 ppm	Titanium (Ti)	<5 ppm
Copper (Cu)	<2 ppm	n-Butyl formate	<0.001%
Iron (Fe)	<5 ppm	Total silicon (Si)	<50 ppm
Lead (Pb)	<2.4 ppm	Total sulphur (S)	<50 ppm
Lithium (Li)	<5 ppm	Vanadium (V)	<5 ppm
Magnesium (Mg)	<10 ppm	Zinc (Zn)	<20 ppm

Strontium nitrate anhydrous

specified reagents for general laboratory work- extra pure

SLR

packaging	Mfr. No
500 g Plastic jar	S/8400/53
N ₂ O ₆ Sr	UN: 1507
CAS: 10042-76-9	GHS: Ox. Sol. 3, Acute Tox.
MW: 211.63	4, Skin Irrit. 2, Eye Irrit. 2,
Density at 20°C: 2.99	STOT SE 3
Boiling point: 645°C	H272, H315, H335, H319, H302
Melting point: 570°C	P220, P261, P280,
ADR/RID: 5.1, III	P305+P351+P338
EINECS: 233-131-9	

Assay: >99%

Guaranteed analysis

Calcium (Ca)	<1000 ppm	Sodium (Na)	<200 ppm
Copper (Cu)	<10 ppm	n-Butyl formate	<0.01%
Iron (Fe)	<20 ppm	Total silicon (Si)	<200 ppm
Lead (Pb)	<10 ppm	Total sulphur (S)	<200 ppm
Magnesium (Mg)	<50 ppm	Zinc (Zn)	<100 ppm
Potassium (K)	<200 ppm		

Styrene stabilized with 10-15 ppm tert-butylcatechol pure

packaging	Mfr. No
250 mℓ Glass bottle	S/8480/07
1 ℓ Glass bottle	S/8480/15
2.5 ℓ Glass bottle	S/8480/17
25 ℓ Metal drum	S/8480/25
C ₈ H ₈	EINECS: 202-851-5
CAS: 100-42-5	UN: 2055
MW: 104.15	GHS: Flam. Liq. 3, Acute
Density at 20°C: 0.906	Tox. 4, Skin Irrit. 2, Eye Irrit. 2
Boiling point: 145 to 146°C	H226, H332, H319, H315
Melting point: -31°C	P304+P340, P302+P352,
Flash point: 31°C	P280, P305+P351+P338, P210
ADR/RID: 3, III	

Fisher Chemical 5L packaging change

Our improved HDPE 5L packaging is easier to handle, more robust and preserves product quality for longer.

Succinic acid

for analysis

Certified AR

packaging	Mfr. No
250 g Plastic jar	S/8522/50
500 g Plastic jar	S/8522/53
10 kg Plastic keg	S/8522/68
C ₄ H ₆ O ₄	EINECS: 203-740-4
CAS: 110-15-6	GHS: Skin Irrit. 2, Eye Irrit. 2,
MW: 118.09	STOT SE 3
Density at 20°C: 1.56	H315, H319, H335
Boiling point: 235°C	P261, P302+P352, P280,
Melting point: 185 to 190°C	P305+P351+P338
Flash point: 206°C	

Assay: >99.5%

Guaranteed analysis

Aluminium (Al)	<2 ppm	Manganese (Mn)	<2 ppm
Barium (Ba)	<2 ppm	Nickel (Ni)	<2 ppm
Cadmium (Cd)	<2 ppm	Potassium (K)	<50 ppm
Calcium (Ca)	<50 ppm	Sodium (Na)	<50 ppm
Chromium (Cr)	<2 ppm	Titanium (Ti)	<2 ppm
Cobalt (Co)	<2 ppm	n-Butyl formate	<0.001%
Copper (Cu)	<5 ppm	Total nitrogen (N)	<20 ppm
Iron (Fe)	<5 ppm	Total silicon (Si)	<20 ppm
Lead (Pb)	<5 ppm	Total sulphur (S)	<20 ppm
Lithium (Li)	<2 ppm	Vanadium (V)	<2 ppm
Magnesium (Mg)	<20 ppm	Zinc (Zn)	<10 ppm

Lot analysis on the label.

Succinic acid

specified reagents for general laboratory work- extra pure

SLR

packaging	Mfr. No
500 g Plastic jar	S/8520/53
C ₄ H ₆ O ₄	EINECS: 203-740-4
CAS: 110-15-6	GHS: Skin Irrit. 2, Eye Irrit. 2,
MW: 118.09	STOT SE 3
Density at 20°C: 1.56	H315, H319, H335
Boiling point: 235°C	P261, P302+P352, P280,
Melting point: 185 to 190°C	P305+P351+P338
Flash point: 206°C	

Assay: >99%

Guaranteed analysis

Calcium (Ca)	<200 ppm	Sodium (Na)	<200 ppm
Copper (Cu)	<20 ppm	n-Butyl formate	<0.005%
Iron (Fe)	<20 ppm	Total silicon (Si)	<100 ppm
Lead (Pb)	<20 ppm	Total sulphur (S)	<100 ppm
Magnesium (Mg)	<100 ppm	Zinc (Zn)	<50 ppm
Potassium (K)	<200 ppm		

Sucrose

for electrophoresis

packaging	Mfr. No
500 g Plastic jar	S/P770/53
3 kg Plastic jar	S/P770/63

C ₁₂ H ₂₂ O ₁₁	Melting point: 190 to 192°C
CAS: 57-50-1	EINECS: 200-334-9
MW: 342.30	

Guaranteed analysis

Reducing sugars	: <0.005%	Lead (Pb)	<0.5 ppm
Copper (Cu)	<0.5 ppm	Total chloride (Cl)	<5 ppm
Iron (Fe)	<0.5 ppm		

UHPLC

Filtered to 0.1 µm, Fisher Chemical UHPLC gradient grade solvents (Acetonitrile, Methanol and Water) are perfect for those working at high pressure.
See page 309 for more information

Sucrose for analysis		Certified AR, Eur.Ph., BP	
packaging		Mfr. No	
500 g	Plastic jar	S/8600/53	
1 kg	Plastic jar	S/8600/60	
3 kg	Plastic jar	S/8600/63	
5 kg	Plastic keg	S/8600/65	
10 kg	Plastic keg	S/8600/68	
25 kg	Fibre Keg	S/8600/70	
C ₁₂ H ₂₂ O ₁₁ CAS: 57-50-1 MW: 342.30		Melting point: 190 to 192°C EINECS: 200-334-9	
Guaranteed analysis			
Acidity / Alkalinity	<0.0008 meq/g	Magnesium (Mg)	<5 ppm
Loss on drying	<0.03%	Potassium (K)	<20 ppm
Specific rotation	66.4° to 66.6°	Sodium (Na)	<50 ppm
Reducing sugars	: <0.005%	Total chloride (Cl)	<0.001%
Insoluble matter	<0.005%	Total nitrogen (N)	<20 ppm
Invert sugar	: <0.05%	Total phosphorus (P)	<20 ppm
Calcium (Ca)	<10 ppm	Total silicon (Si)	<5 ppm
Copper (Cu)	<1 ppm	Total sulphur (S)	<20 ppm
Iron (Fe)	<1 ppm	Zinc (Zn)	<10 ppm
Lead (Pb)	<1 ppm	Residue after ignition	<0.01%
Typical analysis			
Aluminium (Al)	<2 ppm	Lithium (Li)	<2 ppm
Arsenic (As)	<2 ppm	Manganese (Mn)	<2 ppm
Barium (Ba)	<2 ppm	Nickel (Ni)	<2 ppm
Cadmium (Cd)	<2 ppm	Titanium (Ti)	<2 ppm
Chromium (Cr)	<2 ppm	Vanadium (V)	<2 ppm
Cobalt (Co)	<2 ppm		

Lot analysis on the label.

Sucrose		SLR
specified reagents for general laboratory work- extra pure		
packaging		Mfr. No
500 g	Plastic jar	S/8560/53
1 kg	Plastic jar	S/8560/60
3 kg	Plastic jar	S/8560/63
5 kg	Plastic keg	S/8560/65
C ₁₂ H ₂₂ O ₁₁ CAS: 57-50-1 MW: 342.30		Melting point: 190 to 192°C EINECS: 200-334-9
Guaranteed analysis		
Calcium (Ca)	<50 ppm	Sodium (Na) <200 ppm
Copper (Cu)	<5 ppm	n-Butyl formate <0.002%
Iron (Fe)	<5 ppm	Total phosphorus (P) <100 ppm
Lead (Pb)	<5 ppm	Total silicon (Si) <20 ppm
Magnesium (Mg)	<20 ppm	Total sulphur (S) <100 ppm
Potassium (K)	<100 ppm	Zinc (Zn) <500 ppm

Sudan black B C.I. 26150	
pure	
packaging	Mfr. No
25 g Glass jar	S/8630/46
C ₂₆ H ₂₄ N ₆ CAS: 4197-25-5 MW: 456.54	
Melting point: 120 to 124°C EINECS: 224-087-1	
Guaranteed analysis	
Complies with standard IR spectrum: Pass	

Sudan III, C.I. 26100	
technical	
packaging	Mfr. No
25 g Glass jar	S/8620/46
C ₂₂ H ₁₆ N ₄ O CAS: 85-86-9 MW: 352.38 EINECS: 201-638-4	
Guaranteed analysis	
Complies with standard IR spectrum: Pass	

Sudan IV, C.I. 26105	
technical	
packaging	Mfr. No
25 g Glass jar	S/8625/46
C ₂₄ H ₂₀ N ₄ O CAS: 85-83-6 MW: 380.45 Melting point: 181.1°C EINECS: 201-635-8	
GHS: Skin Irrit. 2, Eye Irrit. 2, STOT SE 3 H319, H315, H335 P261, P280, P302+P352, P305+P351+P338	
Guaranteed analysis	
Complies with standard IR spectrum: Pass	

Sulfate standard solution 1000 ppm (1 mg/ml)	
for Ion pair chromatography, Ready-to-use	
packaging	Mfr. No
100 ml HDPE Plastic bottle	J/4564/05
Manufactured to a tolerance of ± 1%.	

Sulfate standard 1000 ppm	
For Ion chromatography in H ₂ O, NIST certified	
packaging	Mfr. No
100 ml Plastic bottle	SPASSO492Y
500 ml Plastic bottle	SPASSO492X
Manufactured to a tolerance of ± 0.3%.	
Check with your local distributor for availability of Spex product	

Sulfate-sulfur standard 1000 ppm	
For Ion chromatography in H ₂ O, NIST certified	
packaging	Mfr. No
100 ml Plastic bottle	SPASSO4592Y
500 ml Plastic bottle	SPASSO4592X
Manufactured to a tolerance of ± 0.3%.	
Check with your local distributor for availability of Spex product	

Sulfamic acid		Certified AR
for analysis		
Packaging		Mfr. No
500 g	Plastic jar	S/8641/53
10 kg	Plastic keg	S/8641/68
$\text{H}_2\text{SO}_3\text{N}$ CAS: 5329-14-6 MW: 97.09 Density at 20°C: 2.15 Melting point: 200 to 208°C ADR/RID: 8, III EINECS: 226-218-8		UN: 2967 GHS: Skin Irrit. 2, Eye Irrit. 2, Aquatic Chronic 3 H412, H315, H319 P280, P305+P351+P338, P273, P302+P350

Cadmium (Cd)	<5 ppm	Potassium (K)	<50 ppm
Calcium (Ca)	<20 ppm	Sodium (Na)	<50 ppm
Copper (Cu)	<5 ppm	Sulfate (SO ₄)	<0.02%
Iron (Fe)	<5 ppm	Total chloride (Cl)	<0.001%
Lead (Pb)	<5 ppm	Total silicon (Si)	<20 ppm
Magnesium (Mg)	<10 ppm	Zinc (Zn)	<10 ppm

Typical analysis			
Aluminium (Al)	<5 ppm	Molybdenum (Mo)	<5 ppm
Barium (Ba)	<5 ppm	Nickel (Ni)	<5 ppm
Chromium (Cr)	<5 ppm	Silver (Ag)	<5 ppm
Cobalt (Co)	<5 ppm	Strontium (Sr)	<5 ppm
Lithium (Li)	<5 ppm	Titanium (Ti)	<5 ppm
Manganese (Mn)	<5 ppm	Vanadium (V)	<5 ppm

Lot analysis on the label.

Sulfamic acid		SLR	
specified reagents for general laboratory work- extra pure			
packaging		Mfr. No	
500 g	Plastic jar	S/8640/53	
25 kg	Fibre Keg	S/8640/70	
H ₂ SO ₃ N			
CAS: 5329-14-6			
UN: 2967			
GHS: Skin Irrit. 2, Eye Irrit. 2, Aquatic Chronic 3			
H412, H315, H319			
P280, P305+P351+P338, P273, P302+P350			
MW: 97.09			
Density at 20°C: 2.15			
Melting point: 200 to 208°C			
ADR/RID: 8, III			
EINECS: 226-218-8			
Assay: >99%			
Guaranteed analysis			
Calcium (Ca)	<100 ppm	Potassium (K)	<200 ppm
Copper (Cu)	<20 ppm	Sodium (Na)	<200 ppm
Iron (Fe)	<20 ppm	Sulfate (SO ₄)	<0.05%
Lead (Pb)	<20 ppm	Total silicon (Si)	<100 ppm
Magnesium (Mg)	<50 ppm	Zinc (Zn)	<50 ppm

Sulfanilamide		SLR
specified reagents for general laboratory work- extra pure		
packaging	Mfr. No	
250 g Plastic jar	S/8660/50	
C ₆ H ₄ N ₂ O ₂ S	GHS: Skin Irrit. 2, Eye Irrit. 2,	
CAS: 63-74-1	STOT SE 3	
MW: 172.21	H315, H335, H319	
Melting point: 163 to 167°C	P261, P302+P352, P280,	
EINECS: 200-563-4	P305+P351+P338	
Assay: >99%		
Guaranteed analysis		
Loss on drying	<0.5%	Iron (Fe) <40 ppm
Residue after ignition	<0.1%	Lead (Pb) <10 ppm
Melting Point tested	164.5 to 166.5°C	n-Butyl formate <0.04%

Sulfanilic acid		Certified AR	
for analysis			
packaging		Mfr. No	
100 g	Plastic jar	S/8720/48	
C ₆ H ₇ NO ₃ S		GHS: Skin Irrit. 2, Eye Irrit. 2,	
CAS: 121-57-3		Skin Sens. 1	
MW: 173.19		H319, H317, H315	
Density at 250°C: 1.485		P280, P302+P352,	
Melting point: 365°C		P305+P351+P338	
EINECS: 204-482-5			
Assay: >99%			
Guaranteed analysis			
Calcium (Ca)	<20 ppm	Potassium (K)	<50 ppm
Copper (Cu)	<5 ppm	Sodium (Na)	<50 ppm
Iron (Fe)	<5 ppm	Sulfate (SO ₄)	<0.005%
Lead (Pb)	<5 ppm	n-Butyl formate	<0.002%
Magnesium (Mg)	<10 ppm	Total silicon (Si)	<20 ppm
Nitrite (NO ₂)	<0.00005%	Zinc (Zn)	<10 ppm

Lot analysis on the label.

We Support Your Chemistry For

- Semi-bulk and bulk requirements
- Tailored solvents and solvent blends or customized solutions
- Testing services or customized packaging and labeling

Please email your contact details and a brief description of your requirements to begei.sha.quotes@thermofisher.com

Sulfanilic acid		SLR	
specified reagents for general laboratory work- extra pure			
Packaging		Mfr. No	
500 g	Plastic jar	S/8680/53	
C ₆ H ₇ NO ₃ S		GHS: Skin Irrit. 2, Eye Irrit. 2,	
CAS: 121-57-3		Skin Sens. 1	
MW: 173.19		H319, H317, H315	
Density at 20°C: 1.485 to 250 C		P280, P302+P352,	
Melting point: 365°C		P305+P351+P338	
EINECS: 204-482-5			
Assay: >99%			
Guaranteed analysis			
Calcium (Ca)	<100 ppm	Sodium (Na)	<200 ppm
Copper (Cu)	<20 ppm	Sulfate (SO ₄)	<0.05%
Iron (Fe)	<20 ppm	n-Butyl formate	<0.02%
Lead (Pb)	<20 ppm	Total silicon (Si)	<100 ppm
Magnesium (Mg)	<50 ppm	Zinc (Zn)	<50 ppm
Potassium (K)	<200 ppm		

5-Sulfosalicylic acid dihydrate		Certified AR	
for analysis			
packaging		Mfr. No	
100 g	Plastic jar	S/8850/48	
500 g	Plastic jar	S/8850/53	
<chem>C7H6O6S.2H2O</chem> CAS: 5965-83-3 MW: 254.21 Melting point: 106 to 110°C ADR/RID: 8, III EINECS: 202-555-6 UN: 2585		GHS: Acute Tox. 4, Skin Corr. 1B, Eye Dam. 1 H314, H302 P301+P330+P331, P280, P305+P351+P338, P310, P303+P361+P353	
Assay: 99.5 to 101.5%		 	
Guaranteed analysis			
Calcium (Ca)	<50 ppm	Sodium (Na)	<50 ppm
Copper (Cu)	<10 ppm	n-Butyl formate	<0.001%
Iron (Fe)	<10 ppm	Total silicon (Si)	<20 ppm
Lead (Pb)	<10 ppm	Zinc (Zn)	<10 ppm
Magnesium (Mg)	<20 ppm	Salicylic acid	<0.01%
Potassium (K)	<50 ppm		

Sulfur standard standard metal solution 1000 ppm	
for ICP, solution aqueous	
packaging	Mfr. No
500 ml HDPE Plastic bottle	J/8324/08
Manufactured to a tolerance of ± 0.5%.	
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.	
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.	

Sulfur standard standard metal solution 10 000 ppm	
for ICP, solution aqueous	
packaging	Mfr. No
100 ml HDPE Plastic bottle	J/8323/05
Manufactured to a tolerance of ± 0.5%.	
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.	
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.	

Save money!

Order 4 bottles in a box

Sulfur standard metal solution 1000 ppm
for ICP solution in H₂O, NIST certified

packaging	Mfr. No
100 mℓ Plastic bottle	SPPLS92Y
250 mℓ Plastic bottle	SPPLS92T
500 mℓ Plastic bottle	SPPLS92X

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Sulfur standard metal solution 10 000 ppm
for ICP solution in H₂O, NIST certified

packaging	Mfr. No
100 mℓ Plastic bottle	SPPLS93Y
500 mℓ Plastic bottle	SPPLS93X

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Sulfur precipitated
specified reagents for general laboratory work- extra pure

packaging	Mfr. No
500 g Plastic jar	S/8990/53
S	ADR/RID: 4.1, III
CAS: 7704-34-9	EINECS: 231-722-6
MW: 32.05	UN: 1350
Density at 20°C: 2.07	GHS: Skin Irrit. 2
Boiling point: 445°C	H315
Melting point: 113°C	P280, P302+P352,
Flash point: 168°C	P321, P332+P313
Guaranteed analysis	
Acidity	<0.01 meq/g Sulfated ash <0.025%
Water	<0.2% Arsenic (As) <2 ppm

Sulfur flowers
specified reagents for general laboratory work- extra pure

packaging	Mfr. No
1 kg Plastic jar	S/8920/60
3 kg Plastic jar	S/8920/63
S	ADR/RID: 4.1, III
CAS: 7704-34-9	EINECS: 231-722-6
MW: 32.05	UN: 1350
Density at 20°C: 2.07	GHS: Skin Irrit. 2
Boiling point: 445°C	H315
Melting point: 113°C	P264, P280, P302+P352,
Flash point: 207.2°C	P321, P332+P313, P362
Guaranteed analysis	
Identity test	: Pass

Sulfur roll
pure

packaging	Mfr. No
1 kg Plastic jar	S/8970/60
S	ADR/RID: 4.1/III
CAS: 7704-34-9	EINECS: 231-722-6
MW: 32.05	UN: 1350
Density at 20°C: 2.07	GHS: Skin Irrit. 2
Boiling point: 445°C	H315
Melting point: 113°C	P264, P280, P302+P352,
Flash point: 207.2°C	P321, P332+P313, P362
Guaranteed analysis	
Identity test	: Pass

Sulfuric acid 93-98%
for ultra trace elemental analysis

packaging	Mfr. No
250 mℓ Pre-cleaned Teflon FEP bottle	A468-250
500 mℓ Pre-cleaned Teflon FEP bottle	A468-500
1 ℓ Pre-cleaned Teflon FEP bottle	A468-1
H ₂ SO ₄	EINECS: 231-639-5
CAS: 7664-93-9	UN: 1830
MW: 98.07	GHS: Skin Corr. 1A, Eye Dam. 1
Density at 20°C: 1.84	H314
Boiling point: 290°C	P301+P330+P331,
Melting point: 10°C	P305+P351+P338, P280
ADR/RID: 8, II	
Assay: 93-98%	

Guaranteed analysis		
Zirconium (Zr)	<10 ppt	Magnesium (Mg) <50 ppt
Gallium (Ga)	<10 ppt	Manganese (Mn) <10 ppt
Terbium (Tb)	<10 ppt	Mercury (Hg) <100 ppt
Indium (In)	<10 ppt	Molybdenum (Mo) <10 ppt
Dysprosium (Dy)	<10 ppt	Neodymium (Nd) <50 ppt
Aluminium (Al)	<50 ppt	Nickel (Ni) <50 ppt
Antimony (Sb)	<50 ppt	Potassium (K) <50 ppt
Arsenic (As)	<500 ppt	Praseodymium (Pr) <10 ppt
Barium (Ba)	<10 ppt	Rhodium (Rh) <50 ppt
Beryllium (Be)	<10 ppt	Samarium (Sm) <10 ppt
Bismuth (Bi)	<10 ppt	Scandium (Sc) <10 ppt
Cadmium (Cd)	<10 ppt	Selenium (Se) <500 ppt
Calcium (Ca)	<50 ppt	Propan--1-ol <50 ppt
Cerium (Ce)	<10 ppt	Sodium (Na) <50 ppt
Cesium (Cs)	<10 ppt	Tantalum (Ta) <100 ppt
Chromium (Cr)	<10 ppt	Tellurium (Te) <50 ppt
Cobalt (Co)	<10 ppt	Thallium (Tl) <10 ppt
Copper (Cu)	<10 ppt	Thorium (Th) <10 ppt
Erbium (Er)	<10 ppt	Thulium (Tm) <10 ppt
Europium (Eu)	<10 ppt	Tin (Sn) <50 ppt
Gadolinium (Gd)	<10 ppt	Titanium (Ti) <50 ppt
Germanium (Ge)	<100 ppt	Uranium (U) <10 ppt
Hafnium (Hf)	<10 ppt	Vanadium (V) <10 ppt
Holmium (Ho)	<10 ppt	Ytterbium (Yb) <10 ppt
Iron (Fe)	<50 ppt	Yttrium (Y) <10 ppt
Lanthanum (La)	<10 ppt	Zinc (Zn) <50 ppt
Lead (Pb)	<10 ppt	Platinum (Pt) <50 ppt
Lithium (Li)	<10 ppt	Rubidium (Rb) <10 ppt
Lutetium (Lu)	<10 ppt	Tungsten (W) <10 ppt

Bottled in class 10 cleanroom conditions. Individually double-bagged in class 100 cleanroom conditions. Each bottle is individually packaged.
Certificate of analysis delivered with each bottle.

Sulfuric acid 93-98%
for trace metal analysis

packaging	Mfr. No
500 mℓ HDPE plastic bottle	A510-P500
1 ℓ HDPE plastic bottle	A510-P1
2.5 ℓ HDPE plastic bottle	A510-P212
H ₂ SO ₄	EINECS: 231-639-5
CAS: 7664-93-9	UN: 1830
MW: 98.07	GHS: Skin Corr. 1A, Eye Dam. 1
Density at 20°C: 1.84	H314
Boiling point: 290°C	P301+P330+P331,
Melting point: 10°C	P305+P351+P338, P280
ADR/RID: 8, II	
Assay: 93-98%	

Guaranteed analysis		
Colour	<10 APHA	Cobalt (Co) <0.5 ppm
Aluminium (Al)	<1 ppm	Copper (Cu) <0.5 ppm
Antimony (Sb)	<1 ppm	Dysprosium (Dy) <0.1 ppm
Arsenic (As)	<0.5 ppm	Erbium (Er) <0.1 ppm
Barium (Ba)	<0.1 ppm	Europium (Eu) <0.1 ppm
Beryllium (Be)	<0.1 ppm	Gadolinium (Gd) <0.1 ppm
Bismuth (Bi)	<0.1 ppm	Gallium (Ga) <0.1 ppm
Cadmium (Cd)	<0.5 ppm	Germanium (Ge) <1 ppm
Calcium (Ca)	<1 ppm	Gold (Au) <0.5 ppb
Cerium (Ce)	<0.1 ppm	Hafnium (Hf) <0.1 ppm
Cesium (Cs)	<0.1 ppm	Holmium (Ho) <0.1 ppm
Chromium (Cr)	<0.5 ppm	



Indium (In)	<0.1 ppm	Silver (Ag)	<1 ppm
Iron (Fe)	<1 ppm	Sodium (Na)	<1 ppm
Lanthanum (La)	<0.1 ppm	Strontium (Sr)	<0.5 ppm
Lead (Pb)	<0.1 ppm	Substances	
Lithium (Li)	<0.5 ppm	reducing KMnO ₄	<20 ppm
Lutetium (Lu)	<0.1 ppm	Tellurium (Te)	<0.1 ppm
Magnesium (Mg)	<1 ppm	Terbium (Tb)	<0.1 ppm
Manganese (Mn)	<0.5 ppm	Thallium (Tl)	<0.1 ppm
Mercury (Hg)	<0.1 ppm	Thorium (Th)	<0.1 ppm
Molybdenum (Mo)	<0.5 ppm	Thulium (Tm)	<0.1 ppm
Neodymium (Nd)	<0.1 ppm	Tin (Sn)	<1 ppm
Nickel (Ni)	0.5	Titanium (Ti)	<1 ppm
Niobium (Nb)	<0.1 ppm	Total chloride (Cl)	<0.7 ppm
Nitrate (NO ₃)	<0.2 ppm	Total phosphorus (P)	<0.05 ppm
Potassium (K)	<1 ppm	Tungsten (W)	<0.5 ppm
Praseodymium (Pr)	<0.1 ppm	Uranium (U)	<0.1 ppm
Rhodium (Rh)	<0.5 ppm	Vanadium (V)	<0.5 ppm
Rubidium (Rb)	<0.5 ppm	Ytterbium (Yb)	<0.1 ppm
Samarium (Sm)	<0.1 ppm	Yttrium (Y)	<0.1 ppm
Scandium (Sc)	<0.1 ppm	Zinc (Zn)	<1 ppm
Selenium (Se)	<10 ppm		

Bottled in class 10 cleanroom conditions.
Lot analysis on the label.

Sulfuric acid 95-98% d=1.83
for trace metal analysis

packaging	Mfr. No
250 mℓ HDPE Plastic bottle	S/9231/PB07
500 mℓ HDPE Plastic bottle	S/9231/PB08
1 ℓ HDPE Plastic bottle	S/9231/PB15
2.5 ℓ HDPE Plastic bottle	S/9231/PB17
H ₂ SO ₄	EINECS: 231-639-5
CAS: 7664-93-9	UN: 1830
MW: 98.07	GHS: Skin Corr. 1A, Eye Dam. 1
Density at 20°C: 1.84	H314
Boiling point: 290°C	P301+P330+P331,
Melting point: 10°C	P305+P351+P338, P280
ADR/RID: 8, II	
Assay: >95%	

Guaranteed analysis			
Aluminium (Al)	<0.02 ppm	Mercury (Hg)	<0.005 ppm
Antimony (Sb)	<0.005 ppm	Molybdenum (Mo)	<0.002 ppm
Arsenic (As)	<0.005 ppm	Nickel (Ni)	<0.002 ppm
Barium (Ba)	<0.02 ppm	Potassium (K)	<0.02 ppm
Beryllium (Be)	<0.005 ppm	Selenium (Se)	<0.005 ppm
Bismuth (Bi)	<0.005 ppm	Silver (Ag)	<0.002 ppm
Boron (B)	<0.005 ppm	Sodium (Na)	<0.2 ppm
Cadmium (Cd)	<0.002 ppm	Strontium (Sr)	<0.002 ppm
Calcium (Ca)	<0.05 ppm	Thallium (Tl)	<0.005 ppm
Chromium (Cr)	<0.005 ppm	Tin (Sn)	<0.005 ppm
Cobalt (Co)	<0.002 ppm	Titanium (Ti)	<0.002 ppm
Copper (Cu)	<0.005 ppm	Total chloride (Cl)	<0.00001%
Gallium (Ga)	<0.005 ppm	Total nitrogen (N)	<5 ppm
Germanium (Ge)	<0.005 ppm	Total phosphorus (P)	<0.05 ppm
Indium (In)	<0.005 ppm	Vanadium (V)	<0.002 ppm
Iron (Fe)	<0.05 ppm	Zinc (Zn)	<0.01 ppm
Lead (Pb)	<0.002 ppm	Zirconium (Zr)	<0.005 ppm
Lithium (Li)	<0.002 ppm	Residue after ignition	<0.0002%
Magnesium (Mg)	<0.02 ppm	Substances	
Manganese (Mn)	<0.002 ppm	reducing KMnO ₄	<0.0001%

Sulfuric acid min 95%
for electronic use (MOS)

packaging	Mfr. No
2.5 ℓ HDPE Plastic bottle	S/9222/PB17
2.5 ℓ Coated Glass bottle	S/9222/PC17
H ₂ SO ₄	EINECS: 231-639-5
CAS: 7664-93-9	UN: 1830
MW: 98.07	GHS: Skin Corr. 1A, Eye Dam. 1
Density at 20°C: 1.84	H314
Boiling point: 290°C	P301+P330+P331,
Melting point: 10°C	P305+P351+P338, P280
ADR/RID: 8, II	
Assay: >95%	

Guaranteed analysis		
Colour	<10 APHA	Manganese (Mn) <0.02 ppm
Substances		Molybdenum (Mo) <0.02 ppm
reducing KMnO ₄	<0.0001%	Nickel (Ni) <0.1 ppm
Aluminium (Al)	<0.1 ppm	Nitrate (NO ₃) <0.00002%
Barium (Ba)	<0.05 ppm	Potassium (K) <1 ppm
Cadmium (Cd)	<0.02 ppm	Sodium (Na) <5 ppm
Calcium (Ca)	<1 ppm	Strontium (Sr) <0.02 ppm
Chromium (Cr)	<0.1 ppm	Titanium (Ti) <0.02 ppm
Cobalt (Co)	<0.02 ppm	n-Butyl formate <0.00002%
Copper (Cu)	<0.1 ppm	Total nitrogen (N) <5 ppm
Iron (Fe)	<0.5 ppm	Total silicon (Si) <0.5 ppm
Lead (Pb)	<0.1 ppm	Vanadium (V) <0.02 ppm
Lithium (Li)	<0.02 ppm	Zinc (Zn) <0.5 ppm
Magnesium (Mg)	<0.5 ppm	

Filtered to 10 µm.

Sulfuric acid 95%
for analysis with low Cadmium, Lead & Mercury level

packaging	Mfr. No
500 mℓ HDPE Plastic bottle	S/9220/PB08
1 ℓ HDPE Plastic bottle	S/9220/PB15
2.5 ℓ HDPE Plastic bottle	S/9220/PB17
H ₂ SO ₄	EINECS: 231-639-5
CAS: 7664-93-9	UN: 1830
MW: 98.07	GHS: Skin Corr. 1A, Eye Dam. 1
Density at 20°C: 1.84	H314
Boiling point: 290°C	P301+P330+P331,
Melting point: 10°C	P305+P351+P338, P280
ADR/RID: 8, II	
Assay: >95.0%	

Guaranteed analysis		
Colour	<10 APHA	Lead (Pb) <0.005 ppm
Substances		Lithium (Li) <0.01 ppm
reducing KMnO ₄	<2 ppm	Magnesium (Mg) <0.05 ppm
Zirconium (Zr)	<0.05 ppm	Manganese (Mn) <0.01 ppm
Aluminium (Al)	<0.05 ppm	Mercury (Hg) <0.005 ppm
Arsenic (As)	<0.01 ppm	Molybdenum (Mo) <0.005 ppm
Barium (Ba)	<0.05 ppm	Nickel (Ni) <0.02 ppm
Beryllium (Be)	<0.01 ppm	Nitrate (NO ₃) <0.2 ppm
Bismuth (Bi)	<0.1 ppm	Phosphate (PO ₄) <1 ppm
Cadmium (Cd)	<0.002 ppm	Potassium (K) <0.1 ppm
Calcium (Ca)	<0.2 ppm	Strontium (Sr) <0.02 ppm
Chloride (Cl)	<0.1 ppm	Thallium (Tl) <0.05 ppm
Chromium (Cr)	<0.05 ppm	Titanium (Ti) <0.05 ppm
Cobalt (Co)	<0.01 ppm	Vanadium (V) <0.01 ppm
Copper (Cu)	<0.02 ppm	Zinc (Zn) <0.05 ppm
Germanium (Ge)	<0.05 ppm	Residue after ignition <5 ppm
Iron (Fe)	<0.2 ppm	

Lot analysis on the label



Sulfuric acid min 95% d=1.83
for analysis

Certified AR

packaging		Mfr. No
1 ℓ	HDPE Plastic bottle	S/9240/PB15
1 ℓ	Coated Glass bottle	S/9240/PC15
2.5 ℓ	Glass bottle	S/9240/17
2.5 ℓ	HDPE Plastic bottle	S/9240/PB17
2.5 ℓ	Coated Glass bottle	S/9240/PC17
25 ℓ	HDPE Plastic container	S/9240/25
45 ℓ	Plastic container	S/9240/26
H ₂ SO ₄ CAS: 7664-93-9 MW: 98.07 Density at 20°C: 1.84 Boiling point: 290°C Melting point: 10°C ADR/RID: 8, II		EINECS: 231-639-5 UN: 1830 GHS: Skin Corr. 1A, Eye Dam. 1 H314 P301+P330+P331, P305+P351+P338, P280
Assay: >95%		

Guaranteed analysis			
Colour	<10 APHA	Manganese (Mn)	<0.01 ppm
Residue after ignition	<0.0005%	Molybdenum (Mo)	<0.02 ppm
Substances		Nickel (Ni)	<0.05 ppm
reducing KMnO ₄	<0.0001%	Nitrate (NO ₃)	<0.00002%
Aluminium (Al)	<0.05 ppm	Potassium (K)	<0.5 ppm
Arsenic (As)	<0.01 ppm	Selenium (Se)	<0.5 ppm
Barium (Ba)	<0.02 ppm	Sodium (Na)	<1 ppm
Cadmium (Cd)	<0.01 ppm	Strontium (Sr)	<0.02 ppm
Calcium (Ca)	<0.5 ppm	Titanium (Ti)	<0.02 ppm
Chromium (Cr)	<0.1 ppm	Total chloride (Cl)	<0.00002%
Cobalt (Co)	<0.01 ppm	Total nitrogen (N)	<5 ppm
Copper (Cu)	<0.05 ppm	Total silicon (Si)	<0.5 ppm
Iron (Fe)	<0.5 ppm	Vanadium (V)	<0.01 ppm
Lead (Pb)	<0.02 ppm	Zinc (Zn)	<0.2 ppm
Lithium (Li)	<0.01 ppm	Propan--1-ol	<0.05 ppm
Magnesium (Mg)	<0.2 ppm		

Sulfuric acid 95% is equivalent to a strength of 18.0 M (36.80N).
Lot analysis on the label.

Sulfuric acid min 95% d=1.83
extra pure , low in nitrogen

SLR

packaging		Mfr. No
500 mℓ	HDPE Plastic bottle	S/9200/PB08
2.5 ℓ	HDPE Plastic bottle	S/9200/PB17
H ₂ SO ₄ CAS: 7664-93-9 MW: 98.07 Density at 20°C: 1.840 Boiling point: 290°C Melting point: 10°C ADR/RID: 8, II		EINECS: 231-639-5 UN: 1830 GHS: Skin Corr. 1A, Eye Dam. 1 H314 P301+P330+P331, P305+P351+P338, P280
		
Assay: >95%		
Guaranteed analysis		
Colour	<20 APHA	Nitrate (NO ₃) <0.00002%
Calcium (Ca)	<10 ppm	Potassium (K) <10 ppm
Copper (Cu)	<5 ppm	Sodium (Na) <50 ppm
Iron (Fe)	<10 ppm	n-Butyl formate <0.005%
Lead (Pb)	<5 ppm	Total nitrogen (N) <2 ppm
Magnesium (Mg)	<5 ppm	Zinc (Zn) <2 ppm
Sulfuric acid 95% is equivalent to a strength of 18.0 M (36.80N).		

Sulfuric acid 95% is equivalent to a strength of 18.0 M (36.80N).

Safe, environmentally friendly and effective packaging

Fisher Chemical outer packaging is UN compliant and 100% recyclable. It's an easy and safe way to reduce your waste management costs and contribute to a cleaner environment.

Sulfuric acid min 95% d=1.83
specified reagents for general laboratory work- extra pure

SLR

packaging		Mfr. No
1 ℓ	HDPE Plastic bottle	S/9160/PB15
2.5 ℓ	HDPE Plastic bottle	S/9160/PB17
5 ℓ	HDPE Plastic bottle	S/9160/PB21
25 ℓ	HDPE Plastic container	S/9160/25
45 ℓ	Plastic container	S/9160/26
H ₂ SO ₄ CAS: 7664-93-9 MW: 98.07 Density at 20°C: 1.840 Boiling point: 290°C Melting point: 10°C ADR/RID: 8, II		EINECS: 231-639-5 UN: 1830 GHS: Skin Corr. 1A, Eye Dam. 1 H314 P301+P330+P331, P305+P351+P338, P280
Assay: >95%		

Guaranteed analysis			
Colour	<20 APHA	Potassium (K)	<10 ppm
Calcium (Ca)	<10 ppm	Sodium (Na)	<50 ppm
Copper (Cu)	<5 ppm	n-Butyl formate	<0.005%
Iron (Fe)	<10 ppm	Total nitrogen (N)	<20 ppm
Lead (Pb)	<5 ppm	Zinc (Zn)	<2 ppm
Magnesium (Mg)	<5 ppm		

Sulfuric acid 95% is equivalent to a strength of 18.0 M (36.80N).

Sulfuric acid min 95% d=1.83
extra pure, for milk testing, (Gerber test)

SLR

Packaging		Mfr. No
2.5 ℓ	HDPE Plastic bottle	S/9360/PB17
H₂SO₄ CAS: 7664-93-9 MW: 98.07 Density at 20°C: 1.840 Boiling point: 290°C Melting point: 10°C ADR/RID: 8, II		EINECS: 231-639-5 UN: 1830 GHS: Skin Corr. 1A, Eye Dam. 1 H314 P301+P330+P331, P305+P351+P338, P280 
Assay: 89.5 to 91%		
Guaranteed analysis		
Colour	<20 APHA	

Sulfuric acid min 95% d=1.83
Technical

packaging		Mfr. No
2.5 ℓ	HDPE Plastic bottle	S/9120/PB17
H ₂ SO ₄ CAS: 7664-93-9 MW: 98.07 Density at 20°C: 1.840 Boiling point: 290°C Melting point: 10°C ADR/RID: 8, II		EINECS: 231-639-5 UN: 1830 GHS: Skin Corr. 1A, Eye Dam. 1 H314 P301+P330+P331, P305+P351+P338, P280
Assay: >95%		

Sulfuric acid 95% is equivalent to a strength of 18.0 M (36.80N).

Sulfuric acid 50% w/v
pure, for formaldehyde determination

packaging		Mfr. No
2.5 ℓ	HDPE Plastic bottle	S/8350/PB17
H ₂ SO ₄ CAS: 7664-93-9 MW: 98.07 Density at 20°C: 1.10 to 1.30 ADR/RID: 8, II EINECS: 231-639-5		UN: 2796 GHS: Skin Corr. 1A H314 P301+P330+P331, P280, P305+P351+P338, P310

Sulfuric acid 25% w/v
pure

packaging		Mfr. No
5 ℓ	HDPE Plastic bottle	H/0564/21
H ₂ SO ₄ CAS: 7664-93-9 MW: 98.07 Density at 20°C: 1.10 to 1.30 ADR/RID: 8, II EINECS: 231-639-5		UN: 2796 GHS: Skin Corr. 1A H314 P301+P330+P331, P280, P305+P351+P338, P310

Sulfuric acid

packaging		Mfr. No
2.5 ℓ	HDPE Plastic bottle	S/9500/PB17
H ₂ SO ₄ CAS: 7664-93-9 MW: 98.07 Density at 20°C: 1.10 to 1.30 ADR/RID: 8, II EINECS: 231-639-5		UN: 2796 GHS: Skin Corr. 1A H314 P301+P330+P331, P280, P305+P351+P338, P310

Sulfuric acid fuming
pure

packaging		Mfr. No
500 mℓ	Coated Glass bottle	S/9440/PB08
H ₂ SO ₄ .SO ₃ CAS: 8014-95-7 MW: 178.14 Density at 20°C: 1.92 Boiling point: 138°C Melting point: 2°C ADR/RID: 8, 6.1, I		UN: 1831 GHS: Skin Corr. 1A, Eye Dam. 1 H314, EUH014 P301+P330+P331, P305+P351+P338, P310, P304+P340, P280, P402+P404
Assay: >104.5%		

Sulfuric acid solution 5M approx. (10N)
pure

packaging		Mfr. No
1 ℓ	HDPE Plastic bottle	J/8400/PB15
2.5 ℓ	HDPE Plastic bottle	J/8400/PB17
10 ℓ	HDPE Plastic container	J/8400/24
H ₂ SO ₄ CAS: 7664-93-9 MW: 98.07 Density at 20°C: 1.10 to 1.30 ADR/RID: 8, II EINECS: 231-639-5		UN: 2796 GHS: Skin Corr. 1A H314 P301+P330+P331, P280, P305+P351+P338, P310
Guaranteed analysis		
Molarity	5.3 M	

Sulfuric acid solution 2M (4N)
pure, Bench Reagent

packaging		Mfr. No
1 ℓ	HDPE Plastic bottle	J/8410/PB15
2.5 ℓ	HDPE Plastic bottle	J/8410/17
10 ℓ	Bag in box	J/8410/24
H ₂ SO ₄ CAS: 7664-93-9 MW: 98.07 Density at 20°C: 1.10 to 1.30 ADR/RID: 8, II EINECS: 231-639-5		UN: 2796 GHS: Skin Corr. 1A H314 P301+P330+P331, P280, P305+P351+P338, P310

Sulfuric acid solution 1M (2N)
NIST Standard solution ready to use, For Volumetric analysis

packaging		Mfr. No
1 ℓ	HDPE Plastic bottle	J/8420/PB15
2.5 ℓ	HDPE Plastic bottle	J/8420/17
10 ℓ	Bag in box	J/8420/24
25 ℓ	HDPE Drum	J/8420/25
CAS: 7664-93-9 Density at 20°C: 1.03 to 1.07 Boiling point: 100°C ADR/RID: 8, II EINECS: 231-639-5		UN: 2796 GHS: Skin Irrit. 2, Eye Irrit. 2 H319, H315 P305+P351+P338, P280

Manufactured to a tolerance of ± 0.1%.
See the full range of Standard Volumetric solution page 323.

Sulfuric acid solution 0.5M (1N)
NIST Standard solution ready to use, For Volumetric analysis

Eur.Ph. , USP, BP

packaging		Mfr. No
1 ℓ	HDPE Plastic bottle	J/8430/15
2.5 ℓ	HDPE Plastic bottle	J/8430/17
10 ℓ	Bag in box	J/8430/24
25 ℓ	HDPE Drum	J/8430/25
CAS: 7664-93-9 Density at 20°C: 1.03 to 1.07 Boiling point: 100°C ADR/RID: 8, II EINECS: 231-639-5		UN: 2796 GHS: Skin Irrit. 2, Eye Irrit. 2 H319, H315 P305+P351+P338, P280

Manufactured to a tolerance of ± 0.1%.
See the full range of Standard Volumetric solution page 323.

Sulfuric acid solution 0.25M (0.5N)
NIST Standard solution ready to use, For Volumetric analysis

packaging		Mfr. No
1 ℓ	HDPE Plastic bottle	J/8440/PB15
2.5 ℓ	HDPE Plastic bottle	J/8440/17
CAS: 7664-93-9 Density at 20°C: 1.01 to 1.03 Boiling point: 100°C ADR/RID: 8, III EINECS: 231-639-5		UN: 3264 GHS: Met. Corr. 1 H290 P234, P390

Manufactured to a tolerance of ± 0.1%.
See the full range of Standard Volumetric solution page 323.

Sulfuric acid solution 0.1M (0.2N)
NIST Standard solution ready to use, For Volumetric analysis

packaging		Mfr. No
1 ℓ	HDPE Plastic bottle	J/8450/15
2.5 ℓ	HDPE Plastic bottle	J/8450/17
10 ℓ	Bag in box	J/8450/24
25 ℓ	HDPE Drum	J/8450/25
CAS: 7664-93-9 Density at 20°C: 1.01 to 1.03 Boiling point: 100°C ADR/RID: 8, III EINECS: 231-639-5		UN: 3264 GHS: Met. Corr. 1 H290 P234, P390

Manufactured to a tolerance of ± 0.1%.
See the full range of Standard Volumetric solution page 323.

UHPLC

Filtered to 0.1 µm, Fisher Chemical UHPLC gradient grade solvents (Acetonitrile, Methanol and Water) are perfect for those working at high pressure.

See page 309 for more information

Sulfuric acid solution 0.05M (0.1N)
NIST Standard solution ready to use, For Volumetric analysis

packaging

Mfr. No

1 ℓ HDPE Plastic bottle
2.5 ℓ HDPE Plastic bottle
10 ℓ Bag in box
25 ℓ HDPE Drum

J/8460/15
J/8460/17
J/8460/24
J/8460/25

CAS: 7664-93-9
Density at 20°C: 1.01 to 1.03
Boiling point: 100°C
ADR/RID: 8, III
EINECS: 231-639-5

UN: 3264
GHS: Met. Corr. 1
H290
P234, P390



Manufactured to a tolerance of ± 0.1%.
See the full range of Standard Volumetric solution page 323.

Sulfuric acid solution 0.05M (0.1N)
NIST Standard Concentrate, For Volumetric analysis

packaging

Mfr. No

1 Ampoule
6 Ampoule

J/8460C/05
J/8460C/90

CAS: 7664-93-9
Density at 20°C: 1.03 to 1.07
ADR/RID: 8, II
EINECS: 231-639-5

UN: 2796
GHS: Skin Irrit. 2, Eye Irrit. 2
H319, H315
P305+P351+P338, P280



Each ampoule makes 1 ℓ of solution.
See the full range of Standard Volumetric solution page 323.

Sulfuric acid solution 0.01M (0.2N)
NIST Standard solution ready to use, For Volumetric analysis

packaging

Mfr. No

1 ℓ HDPE Plastic bottle
2.5 ℓ HDPE Plastic bottle
10 ℓ Bag in box

J/8470/PB15
J/8470/17
J/8470/24

CAS: 7664-93-9
Density at 20°C: 1.00
Boiling point: 100°C
EINECS: 231-639-5



Manufactured to a tolerance of ± 0.3%.
See the full range of Standard Volumetric solution page 323.

Sulfuric acid solution 0.5M (1N)
NIST Standard Concentrate, For Volumetric analysis

packaging

Mfr. No

1 Ampoule
6 Ampoule

J/8430C/05
J/8430C/90

CAS: 7664-93-9
Density at 20°C: 1.10 to 1.30
ADR/RID: 8, II
EINECS: 231-639-5
UN: 2796

GHS: Skin Corr. 1A
H314
P301+P330+P331, P280,
P305+P351+P338, P310



Each ampoule makes 1 ℓ of solution.
See the full range of Standard Volumetric solution page 323.

Sulfurous Acid
pure, approx.6% w/v SO2

packaging

Mfr. No

2.5 ℓ Glass bottle

J/8500/17

H2SO3
CAS: 7782-99-2
MW: 82.07
ADR/RID: 8, II
EINECS: 231-973-1
UN: 1833

GHS: Acute Tox. 4, Skin Corr. 1B
H332, H314
P280, P301+P330+P331,
P305+P351+P338, P310,
P302+P352, P304+P340




Sulfuryl chloride
Technical

packaging

Mfr. No

500 mℓ Coated Glass bottle

S/9560/PB08

SO2Cl2
CAS: 7791-25-5
MW: 134.97
Density at 20°C: 1.667
Boiling point: 69.1°C
Melting point: -54.1°C
ADR/RID: 6.1, 8, I

EINECS: 232-245-6
UN: 1834
GHS: Skin Corr. 1B, STOT SE 3
H335, H314, EUH014
P280, P305+P351+P338,
P304+P312, P260



Talc fine powder
pure

packaging

Mfr. No

1 kg Plastic jar
2.5 kg Plastic jar
10 kg Plastic keg

T/0050/60
T/0050/62
T/0050/68

4MgO4SiO2·H2O
CAS: 14807-96-6
MW: 379.28

Density at 20°C: 2.5 to 2.8
Melting point: 800°C
EINECS: 238-877-9

Guaranteed analysis
Identity test : Pass

Talc fine powder purified by acids
pure

packaging

Mfr. No

500 g Plastic jar

T/0100/53

6MgO4SiO2·H2O
CAS: 14807-96-6
MW: 379.28

Density at 20°C: 2.5 to 2.8
Melting point: 800°C
EINECS: 238-877-9

Guaranteed analysis
Identity test : Pass

Tannic acid
specified reagents for general laboratory work- extra pure

packaging

Mfr. No

500 g Plastic jar

T/0150/53

C76H52O46
CAS: 1401-55-4
MW: 170.22

Melting point: 218°C
Flash point: 198°C
EINECS: 215-753-2

Guaranteed analysis
Loss on drying <9% Residue after ignition <0.2%

(+)-Tartaric acid
for analysis

packaging

Mfr. No

500 g Plastic jar
1 kg Plastic jar

T/0250/53
T/0250/60

C4H6O6
CAS: 87-69-4
MW: 150.09
Density at 20°C: 1.76
Melting point: 168 to 172°C
Flash point: 210°C

EINECS: 201-766-0
GHS: Skin Irrit. 2, Eye Irrit. 2,
STOT SE 3
H335, H315, H319
P280, P305+P351+P338,
P310, P302+P352, P261



Assay: >99.5%
Guaranteed analysis
Calcium (Ca) <20 ppm Sodium (Na) <50 ppm
Copper (Cu) <2 ppm Total chloride (Cl) <0.0005%
Iron (Fe) <2 ppm Total phosphorus (P) <20 ppm
Lead (Pb) <2 ppm Total silicon (Si) <20 ppm
Magnesium (Mg) <10 ppm Total sulphur (S) <100 ppm
Oxalate (C2O4) <0.01% Zinc (Zn) <10 ppm
Potassium (K) <20 ppm



Typical analysis			
Arsenic (As)	<2 ppm	Molybdenum (Mo)	<2 ppm
Barium (Ba)	<2 ppm	Nickel (Ni)	<2 ppm
Cadmium (Cd)	<2 ppm	Silver (Ag)	<2 ppm
Chromium (Cr)	<2 ppm	Strontium (Sr)	<2 ppm
Cobalt (Co)	<2 ppm	Titanium (Ti)	<2 ppm
Lithium (Li)	<2 ppm	Vanadium (V)	<2 ppm
Manganese (Mn)	<2 ppm		

Lot analysis on the label.

(+)-Tartaric acid
extra pure, crystals

packaging

Mfr. No

250 g Plastic jar
500 g Plastic jar
1 kg Plastic jar

T/0200/50
T/0200/53
T/0200/60

C4H6O6
CAS: 87-69-4
MW: 150.09
Density at 20°C: 1.76 C
Melting point: 168 to 172°C
Flash point: 210°C

EINECS: 201-766-0
GHS: Skin Irrit. 2, Eye Irrit. 2,
STOT SE 3
H335, H315, H319
P280, P305+P351+P338,
P310, P302+P352, P261



Assay: >99%
Guaranteed analysis
Calcium (Ca) <100 ppm Sodium (Na) <200 ppm
Copper (Cu) <10 ppm n-Butyl formate <0.002%
Iron (Fe) <10 ppm Total phosphorus (P) <50 ppm
Lead (Pb) <10 ppm Total silicon (Si) <100 ppm
Magnesium (Mg) <10 ppm Total sulphur (S) <500 ppm
Potassium (K) <100 ppm Zinc (Zn) <50 ppm

Tantalum standard metal solution 1000 ppm
for ICP solution in H2O/HF trace, NIST certified

packaging

Mfr. No

100 mℓ Plastic bottle
500 mℓ Plastic bottle

SPPLTA92Y
SPPLTA92X

H319



Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Tantalum standard metal solution 10 000 ppm
for ICP solution in H2O/HF trace, NIST certified

packaging

Mfr. No

100 mℓ Plastic bottle
500 mℓ Plastic bottle

SPPLTA93Y
SPPLTA93X

H319



Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Tap, nickel plated brass for 25ℓ metal drum
3/4" BSP, nickel plated brasss body, chrome ball valve PTFE seat

packaging

Mfr. No

Pack Box

TAP-105-K

(+)-Tartaric acid
extra pure, powder

packaging

Mfr. No

500 g Plastic jar
1 kg Plastic jar
5 kg Plastic keg

T/0210/53
T/0210/60
T/0210/65

C4H6O6
CAS: 87-69-4
MW: 150.09
Density at 20°C: 1.76 C
Melting point: 168 to 172°C
Flash point: 210°C

EINECS: 201-766-0
GHS: Skin Irrit. 2, Eye Irrit. 2,
STOT SE 3
H335, H315, H319
P280, P305+P351+P338,
P310, P302+P352, P261



(+)-Tartaric acid diammonium salt
specified reagents for general laboratory work- extra pure

packaging

Mfr. No

500 g Plastic jar

A/6590/53

(NH4)2C4H4O6
CAS: 3164-29-2
MW: 184.15

Density at 20°C: 1.6
EINECS: 221-618-9

Assay: >98%
Guaranteed analysis
Calcium (Ca) <50 ppm Potassium (K) <100 ppm
Copper (Cu) <50 ppm Sodium (Na) <100 ppm
Iron (Fe) <50 ppm n-Butyl formate <0.005%
Lead (Pb) <50 ppm Zinc (Zn) <50 ppm
Magnesium (Mg) <100 ppm

Sodium tartrate dihydrate
for analysis

packaging

Mfr. No

500 g Plastic jar

S/6922/53

C4H4Na2O6·2H2O
CAS: 6106-24-7
MW: 230.08
EINECS: 212-773-2

Assay: 99 to 101%
Guaranteed analysis
pH (aqueous solution 5%) 7.5 to 8.5 Molybdenum (Mo) <5 ppm
Aluminium (Al) <5 ppm Nickel (Ni) <5 ppm
Arsenic (As) <5 ppm Potassium (K) <20 ppm
Barium (Ba) <5 ppm Propan--1-ol <5 ppm
Cadmium (Cd) <5 ppm Strontium (Sr) <5 ppm
Calcium (Ca) <50 ppm Titanium (Ti) <5 ppm
Chromium (Cr) <5 ppm n-Butyl formate <0.0005%
Cobalt (Co) <5 ppm Total nitrogen (N) <20 ppm
Copper (Cu) <5 ppm Total phosphorus (P) <50 ppm
Iron (Fe) <5 ppm Total silicon (Si) <20 ppm
Lead (Pb) <5 ppm Total sulphur (S) <50 ppm
Lithium (Li) <5 ppm Vanadium (V) <5 ppm
Magnesium (Mg) <20 ppm Zinc (Zn) <10 ppm
Manganese (Mn) <5 ppm

Lot analysis on the label.

pH Buffers
NIST Certified Standard
Solution & Concentrate
A full range to discover on page 63



Tetra-n-butyl ammonium hydroxide solution		SLR
<i>specified reagents for general laboratory work- extra pure</i>		
packaging	Mfr. No	
100 mℓ HDPE Plastic bottle	T/0480/05	
250 mℓ HDPE Plastic bottle	T/0480/07	
C ₄ H ₃₃ NO CAS: 2052-49-5 MW: 259.46 Density at 20°C: 0.995 Boiling point: >100°C ADR/RID: 8, II EINECS: 218-147-6	UN: 3267 GHS: Skin Corr. 1A, Eye Dam. 1 H314 P280, P301+P330+P331, P310, P305+P351+P338, P303+P361+P353	
Assay: 39.5 to 41%		
Guaranteed analysis		
Residue after ignition	<0.05%	Bromide (Br) <0.1%
Bromide	: <0.1%	Sulfate (SO ₄) <0.02%

Tetra-n-butyl ammonium hydroxide solution 0.1 M (0.1N) in toluene/methanol

Standard solution ready to use, For volumetric analysis

packaging	Mfr. No	
500 mℓ Glass bottle	J/8750/08	
1 ℓ Glass bottle	J/8750/15	
2.5 ℓ Glass bottle	J/8750/17	
10 ℓ Bag in box	J/8750/24	
CAS: 2052-49-5 Density at 20°C: 0.850 Flash point: 4°C ADR/RID: 3, 6.1, II EINECS: 218-147-6 UN: 1992 GHS: Flam. Liq. 2, Repr. 2, Acute Tox. 3, Skin Irrit. 2,	Eye Irrit. 2, STOT SE 1, STOT RE 2, Asp. Tox. 1 H225, H331, H311, H301, H370, H361d, H304, H373, H315, H319, H336 P301+P330+P331, P280, P305+P351+P338, P310, P210, P240	  

Manufactured to a tolerance of 5%.
See the full range of Standard Volumetric solution page 323.

Tetra-n-butylammonium iodide		SLR
<i>specified reagents for general laboratory work- extra pure</i>		
packaging	Mfr. No	
100 g Plastic jar	T/0490/48	
C ₁₆ H ₃₆ IN CAS: 311-28-4 MW: 369.38 Melting point: 141 to 150°C EINECS: 206-220-5	GHS: Acute Tox. 4, Skin Irrit. 2, Eye Irrit. 2, STOT SE 3 H302, H335, H315, H319 P261, P302+P352, P280, P305+P351+P338	
Assay: >98%		
Guaranteed analysis		
Residue after ignition	<0.05%	

Tetrabutyl ammonium phosphate		SLR
<i>specified reagents for general laboratory work- extra pure</i>		
packaging	Mfr. No	
5 g Glass jar	T/0505/44	
C ₁₆ H ₃₈ NPO ₄ CAS: 5574-97-0 MW: 339.45 Melting point: 148 to 151°C EINECS: 226-947-1	GHS: Skin Irrit. 2, Eye Irrit. 2 H319, H315 P302+P352, P280, P305+P351+P338	
Assay: >97%		

Tetrachloroethylene		Certified
<i>for UV spectroscopy</i>		
packaging	Mfr. No	
2.5 ℓ Coated Glass bottle	T/0605/PB17	
C ₂ Cl ₄ CAS: 127-18-4 MW: 165.83 Density at 20°C: 1.625 Boiling point: 120 to 122°C Melting point: -22°C	ADR/RID: 6.1, III EINECS: 204-825-9 UN: 1897 GHS: Carc. 2, Aquatic Chronic 2 H351, H411 P281, P273, P260	 
Assay (CG): >99.9%		
Guaranteed analysis		
Alkalinity	<0.00006 meq/g	Residue after evaporation <5 ppm
Water	<0.01%	Acidity <0.0003 meq/g
Absorbance specification		
Max. Absorbance (1cm cell v HPLC water)		
290 nm (sol. aqu. 1%)	: <0.125 A.U.	340 nm <0.05 A.U.
310 nm	<0.08 A.U.	360 nm <0.05 A.U.
320 nm	<0.05 A.U.	400 nm <0.01 A.U.

Tetrachloroethylene		Certified AR
<i>for analysis</i>		
packaging	Mfr. No	
2.5 ℓ Glass bottle	T/0601/17	
2.5 ℓ Coated Glass bottle	T/0601/PB17	
25 ℓ Metal drum	T/0601/25	
C ₂ Cl ₄ CAS: 127-18-4 MW: 165.83 Density at 20°C: 1.625 Boiling point: 120 to 122°C Melting point: -22°C	ADR/RID: 6.1, III EINECS: 204-825-9 UN: 1897 GHS: Carc. 2, Aquatic Chronic 2 H351, H411 P281, P273, P260	 
Assay (CG): >99.5%		
Guaranteed analysis		
Acidity / Alkalinity	<0.0003 meq/g	Iron (Fe) <0.1 ppm
Colour	<10 APHA	Lead (Pb) <0.05 ppm
Water	<0.01%	Magnesium (Mg) <0.2 ppm
Residue after evaporation	<20 ppm	Potassium (K) <0.5 ppm
Calcium (Ca)	<0.5 ppm	Sodium (Na) <0.5 ppm
Chloride (Cl)	<0.0001%	Total phosphorus (P) <0.1 ppm
Chlorine	<0.3 ppm	Total sulphur (S) <0.5 ppm
Free Chlorine (Cl ₂)	<0.3 ppm	Zinc (Zn) <0.1 ppm
Copper (Cu)	<0.05 ppm	Acetone <0.02%
Typical analysis		
Aluminium (Al)	<0.02 ppm	Molybdenum (Mo) <0.02 ppm
Barium (Ba)	<0.02 ppm	Nickel (Ni) <0.02 ppm
Cadmium (Cd)	<0.02 ppm	Strontium (Sr) <0.02 ppm
Chromium (Cr)	<0.02 ppm	Titanium (Ti) <0.02 ppm
Cobalt (Co)	<0.02 ppm	Total silicon (Si) <0.02 ppm
Lithium (Li)	<0.02 ppm	Vanadium (V) <0.02 ppm
Manganese (Mn)	<0.02 ppm	

Lot analysis on the label.

Tetrachloroethylene		NEW PACK.	SLR
<i>specified reagents for general laboratory work- extra pure</i>			
packaging	Mfr. No		
2.5 ℓ Glass bottle	T/0600/17		
2.5 ℓ Coated Glass bottle	T/0600/PB17		
5 ℓ Fluorinated HDPE Plastic bottle	T/0600/21		
25 ℓ Metal drum	T/0600/25		
C ₂ Cl ₄ CAS: 127-18-4 MW: 165.83 Density at 20°C: 1.625 Boiling point: 120 to 122°C Melting point: -22°C	ADR/RID: 6.1, III EINECS: 204-825-9 UN: 1897 GHS: Carc. 2, Aquatic Chronic 2 H351, H411 P281, P273, P260	 	
Assay (CG): >99%			

Guaranteed analysis			
Acidity / Alkalinity	<0.003 meq/g	Magnesium (Mg)	<2 ppm
Colour	<10 APHA	Potassium (K)	<2 ppm
Water	<0.1%	Sodium (Na)	<5 ppm
Residue after evaporation	<50 ppm	Total phosphorus (P)	<2 ppm
Calcium (Ca)	<5 ppm	Total silicon (Si)	<2 ppm
Copper (Cu)	<2 ppm	Total sulphur (S)	<2 ppm
Iron (Fe)	<5 ppm	Zinc (Zn)	<2 ppm
Lead (Pb)	<2 ppm		

Tetrachloroethylene

Technical

packaging	Mfr. No	
2.5 ℓ Glass bottle	T/0550/17	
2.5 ℓ Coated Glass bottle	T/0550/PB17	
25 ℓ Metal drum	T/0550/25	
C ₂ Cl ₄ CAS: 127-18-4 MW: 165.83 Density at 20°C: 1.625 Boiling point: 120 to 122°C Melting point: -22°C	ADR/RID: 6.1, III EINECS: 204-825-9 UN: 1897 GHS: Carc. 2, Aquatic Chronic 2 H351, H411 P281, P273, P260	 

Tetrachlorosilane

See Silicon tetrachloride, page 276

Tetradecyltrimethylammonium bromide		Certified
<i>for HPLC</i>		
packaging	Mfr. No	
25 g Plastic jar	T/0611/46	
C ₁₇ H ₃₈ NBr CAS: 1119-97-7 MW: 336.39 Melting point: 245 to 250°C ADR/RID: 6.1, II EINECS: 214-291-9	UN: 1759 GHS: Skin Corr. 1B, Eye Dam. 1 H314 P301+P330+P331, P280, P305+P351+P338, P310	
Assay: >98%		
Absorbance specification		
Max. Absorbance (1cm cell v HPLC water)		
254 nm	<0.05 A.U.	

Lot analysis on the label.
Filtered to 0.2 µm.

Tetraethylammonium bromide		Certified
<i>for HPLC</i>		
packaging	Mfr. No	
25 g Plastic jar	T/0621/46	
C ₈ H ₂₀ NBr CAS: 71-91-0 MW: 210.16 Melting point: 285°C EINECS: 200-769-4	GHS: Skin Irrit. 2, Eye Irrit. 2, STOT SE 3 H315, H335, H319 P261, P302+P352, P280, P305+P351+P338	
Assay: >98%		
Absorbance specification		
Max. Absorbance (1cm cell v HPLC water)		
254 nm (sol. aqu. 0.1 M)	: <0.05 A.U.	

Lot analysis on the label.
Filtered to 0.2 µm.




Tetraethylammonium bromide		SLR
<i>specified reagents for general laboratory work- extra pure</i>		
packaging	Mfr. No	
100 g Plastic jar	T/0620/48	
C ₈ H ₂₀ NBr CAS: 71-91-0 MW: 210.16 Melting point: 285°C EINECS: 200-769-4	GHS: Skin Irrit. 2, Eye Irrit. 2, STOT SE 3 H315, H335, H319 P261, P302+P352, P280, P305+P351+P338	
Assay: >98%		

Tetrahydrofuran unstabilised		NEW PACK.	Certified
For HPLC			
Packaging		Mfr. No	
1 ℓ	Coated Glass bottle	T/0706/PB15	
2.5 ℓ	Coated Glass bottle	T/0706/PB17	
10 ℓ	Returnable Drum	T/0706/10RSS	
25 ℓ	Returnable Drum	T/0706/25RSS	
30 ℓ	Returnable Drum	T/0706/30RSS	
C ₄ H ₈ O		EINECS: 203-726-8	
CAS: 109-99-9		UN: 2056	
MW: 72.11		GHS: Flam. Liq. 2, Eye Irrit.	
Density at 20°C: 0.880		2, STOT SE 3	
Boiling point: 66°C		H225, H319, H335, EUH019	
Melting point: -108.4°C		P261, P233, P280,	
Flash point: -21°C		P305+P351+P338, P210, P240	
ADR/RID: 3, II			
Assay (CG): >99.8%			

Tetrahydrofuran stabilised with 0.025% BHT

for Gel Permeation Chromatography (GPC)

packaging	Mfr. No	
2.5 ℓ Coated Glass bottle	T/0709/PB17	
C ₄ H ₈ O CAS: 109-99-9 MW: 72.11 Density at 20°C: 0.880 Boiling point: 66°C Melting point: -108.4°C Flash point: -21°C ADR/RID: 3, II EINECS: 203-726-8	UN: 2056 GHS: Flam. Liq. 2, Eye Irrit. 2, STOT SE 3 H225, H335, H319, EUH019 P310, P304+P340, P261, P273, P233, P302+P352, P280, P305+P351+P338, P210, P240	 
Guaranteed analysis		
Colour	<10 APHA	Residue after evaporation <20 ppm
Water	<0.03%	Residue after ignition <20%
Stabilised with 0.025% BHT.		



Tetrahydrofuran unstabilised **NEW PACK.**
extra dry, for synthesis

packaging		Mfr. No
1 ℓ	Coated Glass bottle	T/0702/PB15
2.5 ℓ	Coated Glass bottle	T/0702/PB17
C ₄ H ₈ O	EINECS: 203-726-8	
CAS: 109-99-9	UN: 2056	
MW: 72.11	GHS: Flam. Liq. 2, Eye Irrit.	
Density at 20°C: 0.880	2, STOT SE 3	
Boiling point: 66°C	H225, H319, H335, EUH019	
Melting point: -108.4°C	P261, P233, P280,	
Flash point: -21°C	P305+P351+P338, P210, P240	
ADR/RID: 3, II		
Assay (CG): >99.8		
Guaranteed analysis		
Water	<0.01%	Residue after evaporation <2 ppm
Recommended for the synthesis of DNA and peptide sequencing.		

Tetrahydrofuran stabilised with 0.025% BHT **NEW PACK.** Certified AR
for analysis

packaging		Mfr. No
1 ℓ	Glass bottle	T/0701/15
1 ℓ	Coated Glass bottle	T/0701/PB15
2.5 ℓ	Glass bottle	T/0701/17
2.5 ℓ	Coated Glass bottle	T/0701/PB17
5 ℓ	Fluorinated HDPE Plastic bottle	T/0701/21
25 ℓ	Metal drum	T/0701/25
C ₄ H ₈ O	UN: 2056	
CAS: 109-99-9	GHS: Flam. Liq. 2, Eye Irrit.	
MW: 72.11	2, STOT SE 3	
Density at 20°C: 0.880	H225, H335, H319, EUH019	
Boiling point: 66°C	P310, P304+P340, P261,	
Melting point: -108.4°C	P273, P233, P302+P352,	
Flash point: -21°C	P280, P305+P351+P338,	
ADR/RID: 3, II	P210, P240	
EINECS: 203-726-8		
Assay (CG): >99.8%		
Guaranteed analysis		
Acidity / Alkalinity	<0.0002 meq/g	Magnesium (Mg) <0.1 ppm
Colour	<10 APHA	Potassium (K) <0.5 ppm
Water	<0.03%	Sodium (Na) <1 ppm
Residue after ignition	<0.002%	Total phosphorus (P) <0.1 ppm
Calcium (Ca)	<0.5 ppm	Total sulphur (S) <0.5 ppm
Copper (Cu)	<0.05 ppm	Zinc (Zn) <0.2 ppm
Iron (Fe)	<0.1 ppm	Methyl tetrahydrofuran <0.1%
Lead (Pb)	<0.05 ppm	Peroxide <0.005%
Typical analysis		
Aluminium (Al)	<0.02 ppm	Molybdenum (Mo) <0.02 ppm
Barium (Ba)	<0.02 ppm	Nickel (Ni) <0.02 ppm
Cadmium (Cd)	<0.02 ppm	Strontium (Sr) <0.02 ppm
Chromium (Cr)	<0.02 ppm	Titanium (Ti) <0.02 ppm
Cobalt (Co)	<0.02 ppm	Total silicon (Si) <0.02 ppm
Lithium (Li)	<0.02 ppm	Vanadium (V) <0.02 ppm
Manganese (Mn)	<0.02 ppm	
Lot analysis on the label.		
Stabilised with 0.025% BHT.		

We Support Your Chemistry For

- Semi-bulk and bulk requirements
- Tailored solvents and solvent blends or customized solutions
- Testing services or customized packaging and labeling

Please email your contact details and a brief description of your requirements to begei.sha.quotes@thermofisher.com



Tetrahydrofuran stabilised with 0.025% BHT **NEW PACK.** SLR
specified reagents for general laboratory work- extra pure

packaging		Mfr. No
500 mℓ	Glass bottle	T/0700/08
2.5 ℓ	Glass bottle	T/0700/17
2.5 ℓ	Coated Glass bottle	T/0700/PB17
5 ℓ	HDPE Plastic bottle	T/0700/21
5 ℓ	HDPE Plastic bottle	T/0700/PC21
25 ℓ	Metal drum	T/0700/25
200 ℓ	Metal drum	T/0700/27
C ₄ H ₈ O	UN: 2056	
CAS: 109-99-9	GHS: Flam. Liq. 2, Eye Irrit.	
MW: 72.11	2, STOT SE 3	
Density at 20°C: 0.880	H225, H335, H319, EUH019	
Boiling point: 66°C	P310, P304+P340, P261,	
Melting point: -108.4°C	P273, P233, P302+P352,	
Flash point: -21°C	P280, P305+P351+P338,	
ADR/RID: 3, II	P210, P240	
EINECS: 203-726-8		
Assay (CG): >99.5%		
Guaranteed analysis		
Acidity / Alkalinity	<0.0005 meq/g	Potassium (K) <2 ppm
Colour	<10 APHA	Sodium (Na) <5 ppm
Water	<0.05%	Total phosphorus (P) <2 ppm
Calcium (Ca)	<5 ppm	Total silicon (Si) <2 ppm
Copper (Cu)	<2 ppm	Total sulphur (S) <2 ppm
Iron (Fe)	<5 ppm	Zinc (Zn) <2 ppm
Lead (Pb)	<2 ppm	Residue after ignition <0.005%
Magnesium (Mg)	<2 ppm	Peroxide <0.01%
Stabilised with 0.025% BHT.		

Tetrahydronaphthalene SLR
specified reagents for general laboratory work- extra pure

packaging		Mfr. No
500 mℓ	Glass bottle	T/0850/08
C ₁₀ H ₁₂	EINECS: 204-340-2	
CAS: 119-64-2	UN: 3082	
MW: 132.21	GHS: Skin Irrit. 2, Eye Irrit. 2,	
Density at 20°C: 0.973	Aquatic Chronic 2	
Boiling point: 207°C	H319, H411, H315, EUH019	
Melting point: -35°C	P273, P233, P302+P352,	
Flash point: 77°C	P280, P305+P351+P338	
ADR/RID: 9, III		
Assay (CG): >97%		

Tetramethyl ammonium chloride Certified
for HPLC

Packaging		Mfr. No
25 g	Plastic jar	T/0931/46
C ₄ H ₁₂ NCl	GHS: Acute Tox. 3, Skin Irrit.	
CAS: 75-57-0	2, Eye Irrit. 2, STOT SE 3	
MW: 109.60	H335, H312, H315, H301, H319	
Melting point: 420°C	P301+P310, P280,	
ADR/RID: 6.1, II	P302+P352, P261,	
EINECS: 200-880-8	P305+P351+P338	
UN: 2811		
Assay: >98%		
Absorbance specification		
Max. Absorbance (1cm cell v HPLC water)		
254 nm (sol. aqu. 0.1 M) : <0.05 A.U.		
Lot analysis on the label.		
Filtered to 0.2 µm.		

pH Buffers

NIST Certified Standard Solution & Concentrate

A full range to discover on page 63



Tetramethyl ammonium hydrogen sulfate Certified
for HPLC

packaging		Mfr. No
5 g	Glass jar	T/0935/44
C ₄ H ₁₃ NSO ₄	GHS: Skin Irrit. 2, Eye Irrit. 2, STOT SE 3	
CAS: 80526-82-5	H319, H315, H335	
MW: 171.21	P261, P302+P352, P280, P305+P351+P338	
Melting point: 295 to 297°C		
EINECS: 279-490-5		
Assay: >99%		
Guaranteed analysis		
Loss on drying 100°C		<2%
Absorbance specification		
Max. Absorbance (1cm cell v HPLC water)		
200 nm	<0.3 A.U.	220 nm <0.05 A.U.
Lot analysis on the label.		
Filtered to 0.2 µm.		

Tetramethyl ammonium hydroxide solution 25% in water SLR
specified reagents for general laboratory work- extra pure

packaging		Mfr. No
25 mℓ	HDPE Plastic bottle	T/0940/04
100 mℓ	HDPE Plastic bottle	T/0940/05
C ₄ H ₁₃ NO	UN: 1835	
CAS: 75-59-2	GHS: Acute Tox. 3, Skin Corr. 1B	
MW: 91.15	H314, H301, H311	
Density at 20°C: 1.02	P280, P301+P330+P331,	
Boiling point: 110°C	P305+P351+P338, P310,	
ADR/RID: 8, II	P302+P352	
EINECS: 200-882-9		
Assay: >24.5%		
Guaranteed analysis		
Colour	<20 APHA	Halogenated compounds <0.2%
Halides (Cl)	<0.2%	Residue after ignition <0.3%

Polyethylene glycol mono [4-(1.1,3.3-tetra-methylbutyl)phenyl] ether

See Triton X-100®, page 363



N,N,N',N'-Tetramethyl-p-phenylene-diamine dihydrochloride SLR
specified reagents for general laboratory work- extra pure

packaging		Mfr. No
5 g	Plastic jar	T/0960/44
10 g	Plastic jar	T/0960/45
C ₁₀ H ₁₆ N ₂ ·2HCl	GHS: Acute Tox. 4, Skin Irrit.	
CAS: 637-01-4	2, Eye Irrit. 2, STOT SE 3	
MW: 237.17	H312, H315, H302, H335,	
Melting point: 220 to 224°C	H319, H332	
ADR/RID: 6.1, III	P280, P302+P352, P261,	
EINECS: 211-274-8	P301+P312, P304+P340,	
UN: 2811	P305+P351+P338	
Assay: >97%		

TFA

See Trifluoroacetic acid, page 358



UHPLC

Filtered to 0.1 μm, Fisher Chemical UHPLC gradient grade solvents (Acetonitrile, Methanol and Water) are perfect for those working at high pressure.

See page 309 for more information



Thallium standard metal solution 1000 ppm
for ICP solution in HNO₃ 2%, NIST certified

packaging		Mfr. No
100 mℓ	Plastic bottle	SPPLTL22Y
250 mℓ	Plastic bottle	SPPLTL22T
500 mℓ	Plastic bottle	SPPLTL22X
H315, H319		
Manufactured to a tolerance of ± 0.5%.		
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.		
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.		
Check with your local distributor for availability of Spex product		

Thallium standard metal solution 10 000 ppm
for ICP solution in HNO₃ 5%, NIST certified

packaging		Mfr. No
100 mℓ	Plastic bottle	PLTL23Y
500 mℓ	Plastic bottle	PLTL23X
H302, H314		
Manufactured to a tolerance of ± 0.5%.		
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.		
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.		
Check with your local distributor for availability of Spex product		

Thallium standard metal solution 1000 ppm
for ICP-MS solution in HNO₃ 2%, NIST certified

packaging		Mfr. No
100 mℓ	Plastic bottle	SPCLTL22Y
H315, H320		
Manufactured to a tolerance of ± 0.5%.		
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.		
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.		
Check with your local distributor for availability of Spex product		

Thermit ignition powder
pure

Packaging		Mfr. No
250 g	Glass jar	T/0978/50
ADR/RID: 5.1, II	react. 2, Ox. Sol. 2, Acute Tox. 4	
EINECS: 215-128-4	H228, H261, H272, H302, H332	
UN: 1479	P210, P220, P221,	
GHS: Flam. Sol. 1, Water-	P301+P312, P304+P340, P240	
Guaranteed analysis		
Identity test : Pass		

Thermit powder		
pure		
packaging	Mfr. No	
500 g Glass jar	T/0977/53	
ADR/RID: 4.1, III	H228, H261	
EINECS: 215-168-2	P233, P210, P240, P280,	
UN: 3178	P223, P231+P232	
GHS: Flam. Sol. 1, Water-react. 2		
Guaranteed analysis		
Identity test : Pass		

Thioacetamide		SLR
specified reagents for general laboratory work- extra pure		
packaging	Mfr. No	
25 g Plastic jar	T/0980/46	
C ₂ H ₅ NS	Skin Irrit. 2, Eye Irrit. 2,	
CAS: 62-55-5	Aquatic Chronic 3	
MW: 75.13	H315, H350, H412, H319, H302	
Density at 20°C: 1.37	P201, P308+P313, P273,	
Melting point: 112 to 114°C	P301+P312, P302+P352,	
EINECS: 200-541-4	P280, P305+P351+P338	
GHS: Carc. 1B, Acute Tox. 4,		
Assay: >98%		
Guaranteed analysis		
Residue after ignition <0.05%		

2,2'-Thiodiethanol			
for amino acid analysis, (thioglycolic acid)			
packaging		Mfr. No	
100 g	Coated Glass bottle	T/0998/48	
C ₄ H ₁₀ O ₂ S	EINECS: 203-874-3		
CAS: 111-48-8	GHS: Eye Irrit. 2		
MW: 122.19	H319		
Density at 20°C: 1.18	P280, P305+P351+P338		
Boiling point: 165°C			
Guaranteed analysis			
Colour	<5 APHA Ninhydrin reaction test		: Pass
packed under nitrogen			

Thioglycolic acid		SLR
specified reagents for general laboratory work- extra pure		
packaging		Mfr. No
100 mℓ	Coated Glass bottle	T/1000/PB05
250 mℓ	Coated Glass bottle	T/1000/PB07
2.5 ℓ	Coated Glass bottle	T/1000/PB17
C ₂ H ₄ O ₂ S	UN: 1940	
CAS: 68-11-1	GHS: Acute Tox. 3, Skin	
MW: 92.12	Corr. 1B, Eye Dam. 1	
Density at 20°C: 1.325	H314, H301, H311, H331	
Boiling point: 220°C	P301+P310, P280,	
Melting point: -16°C	P302+P350,	
Flash point: 126°C	P301+P330+P331,	
ADR/RID: 8, II	P305+P351+P338, P304+P340	
EINECS: 200-677-4		
Assay: >97%		
Guaranteed analysis		
Sensitivity test (iron) : Pass		

Save money!

Order 4 bottles in a box



Thiourea			
<i>specified reagents for general laboratory work- extra pure</i>			SLR
packaging	Mfr. No		
250 g Plastic jar	T/1150/50		
500 g Plastic jar	T/1150/53		
10 kg Plastic keg	T/1150/68		
CH ₄ N ₂ S	UN: 2811		
CAS: 62-56-6	GHS: Carc. 2, Repr. 2, Acute		
MW: 76.12	Tox. 4, Aquatic Chronic 2		
Density at 20°C: 1.405	H351, H361d, H302, H411		
Melting point: 171 to 175°C	P301+P310, P310, P281,		
ADR/RID: 6.1, III	P304+P340, P273		
EINECS: 200-543-5			

Assay: >99%		
Guaranteed analysis		
Calcium (Ca)	<200 ppm	Sodium (Na)
Copper (Cu)	<20 ppm	n-Butyl formate
Iron (Fe)	<20 ppm	Total phosphorus (P)
Lead (Pb)	<20 ppm	Total silicon (Si)
Magnesium (Mg)	<100 ppm	Zinc (Zn)
Potassium (K)	<100 ppm	

Thorium standard metal solution 1000 ppm			
<i>for ICP solution in HNO₃ 2%, NIST certified</i>			
packaging	Mfr. No		
100 mL Plastic bottle	SPPLTH22Y		
500 mL Plastic bottle	SPPLTH22X		
H315, H319			

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Thorium standard metal solution 1000 ppm			
<i>for ICP-MS solution in HNO₃ 2%, NIST certified</i>			
packaging	Mfr. No		
100 mL Plastic bottle	SPCLTH22Y		
H315, H319			

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Thullium standard metal solution 1000 ppm			
<i>for ICP solution in HNO₃ 2%, NIST certified</i>			
packaging	Mfr. No		
100 mL Plastic bottle	SPPLTM22Y		
500 mL Plastic bottle	SPPLTM22X		
H315, H319			

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Thymol		SLR
specified reagents for general laboratory work- extra pure		
packaging		Mfr. No
100 g	Plastic jar	T/1450/48
C ₁₀ H ₁₄ O	EINECS: 201-944-8	
CAS: 89-83-8	UN: 2430	
MW: 150.22	GHS: Acute Tox. 4, Skin	
Density at 20°C: 0.960	Corr. 1B, Aquatic Chronic 2	
Boiling point: 232°C	H411, H302, H314	
Freezing point: <49°C	P301+P330+P331, P280,	
Melting point: 48 to 52°C	P305+P351+P338, P310,	
Flash point: 102°C	P273, P301+P312	
ADR/RID: 8, III		
Guaranteed analysis		
Residue after evaporation <0.05%		

Thymol blue solution 0.04%		
pure, indicator grade		
packaging		Mfr. No
500 mL	Glass bottle	T/1470L/08
Flash point: 32°C	GHS: Flam. Liq. 3	
ADR/RID: 3, III	H226	
UN: 1170	P210	
Guaranteed analysis		
pH colour change test : Pass		

Thymol blue	
pure, indicator grade	
packaging	Mfr. No
5 g Glass jar	T/1470/44
C ₂₇ H ₃₀ O ₅ S	Melting point: 221°C
CAS: 76-61-9	EINECS: 200-973-3
MW: 466.58	
Guaranteed analysis	
pH colour change test : Pass Complies with standard IR spectrum: Pass	
pH 1.2-2.8 red to yellow.	
pH 8.0-9.6 yellow to blue.	

Thymolphthalein solid		
pure, indicator grade		
packaging		Mfr. No
5 g	Glass jar	T/1480/44
25 g	Glass jar	T/1480/46
C ₂₈ H ₃₀ O ₄	Melting point: 251 to 253°C	
CAS: 125-20-2	EINECS: 204-729-7	
MW: 430.53		
Guaranteed analysis		
pH colour change test : Pass Complies with standard IR spectrum: Pass		
pH 9.3-10.5 colorless to blue.		
Conforming to the IR spectrum of reference.		

Thymolphthalein solution 0.2%	
pure, indicator grade	
packaging	Mfr. No
500 ml Glass bottle	T/1480L/08
Flash point: 22°C	GHS: Flam. Liq. 2
ADR/RID: 3, II	H225
EINECS: 200-578-6	P210, P240, P241, P280,
UN: 1170	P303+P361+P353
	
Guaranteed analysis	
pH colour change test : Pass	

Thymolsulfonephthaleindi-(methyl-iminodiacetic acid) sodium salt		
pure		
packaging	Mfr. No	
1 g Glass jar	T/1500/43	
C₃₇H₄₁N₂Na₃O₁₃S CAS: 1945-77-3 MW: 822.78 EINECS: 217-743-3		
Guaranteed analysis		
Complies with standard IR spectrum: Pass		

Tillman’s reagent *See Phenol-indo-2, page 240*

Tin standard metal solution 1000 ppm			
<i>for AAS, solution HCl 1M</i>			
packaging	Mfr. No		
100 mL HDPE Plastic bottle	J/8067/05		
500 mL HDPE Plastic bottle	J/8067/08		
GHS: Skin Irrit. 2			
H315			
P302+P352			

Manufactured to a tolerance of ± 0.5%.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.

Tin standard metal solution 1000 ppm			
<i>for ICP, solution HCl 1M</i>			
packaging	Mfr. No		
500 mL HDPE Plastic bottle	J/8326/08		
GHS: Skin Irrit. 2			
H315			
P302+P352			

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.

Tin standard metal solution 10 000 ppm			
<i>for ICP, solution HCl 1M</i>			
packaging	Mfr. No		
100 mL HDPE Plastic bottle	J/8325/05		
GHS: Skin Irrit. 2			
H315			
P302+P352			

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.

Safe, environmentally friendly and effective packaging

Fisher Chemical outer packaging is UN compliant and 100% recyclable. It's an easy and safe way to reduce your waste management costs and contribute to a cleaner environment.



Tin standard metal solution 1000 ppm
for ICP solution in HNO₃ 1%/HF 1%, NIST certified

packaging	Mfr. No
500 mℓ Plastic bottle	SPPLSN22X
H311, H319	

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Tin standard metal solution 1000 ppm
for ICP solution in HCL 20%, NIST certified

packaging	Mfr. No
100 mℓ Plastic bottle	SPPLSN52Y
250 mℓ Plastic bottle	SPPLSN52T
500 mℓ Plastic bottle	SPPLSN52X
H315, H319, H335	

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Tin standard metal solution 10 000 ppm
for ICP solution in HNO₃ 2%/HF 2%, NIST certified

packaging	Mfr. No
500 mℓ Plastic bottle	SPPLSN23X
H311, H314, H333	

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Tin standard metal solution 10 000 ppm
for ICP solution in HCL 20%, NIST certified

packaging	Mfr. No
100 mℓ Plastic bottle	SPPLSN53Y
500 mℓ Plastic bottle	SPPLSN53X
H315, H319, H335	

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Tin standard metal solution 1000 ppm
for ICP-MS solution in HNO₃ 1%/HF 1%, NIST certified

packaging	Mfr. No
100 mℓ Plastic bottle	SPCLSN22Y
H311, H319	

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurites.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Tin (metal)
for analysis

packaging	Mfr. No
250 g Plastic jar	T/1640/50
Sn	Boiling point: 2270°C
CAS: 7440-31-5	Melting point: 231.9°C
MW: 118.69	EINECS: 231-141-8
Density at 20°C: 7.31	
Assay: >99.5%	
Guaranteed analysis	
Antimony (Sb) <200 ppm	Iron (Fe) <100 ppm
Bismuth (Bi) <40 ppm	Lead (Pb) <50 ppm
Copper (Cu) <10 ppm	Zinc (Zn) <20 ppm

The prices applied to the tin and derivatives vary with the price of the raw material.

Tin metal powder fine
pure

packaging	Mfr. No
100 g Plastic jar	T/1650/48
Sn	Melting point: 231.9°C
CAS: 7440-31-5	EINECS: 231-141-8
MW: 118.69	GHS: Eye Irrit. 2, STOT SE 3
Density at 20°C: 7.31	H319, H335
Boiling point: 2270°C	P280, P305+P351+P338
Guaranteed analysis	
Identity test	: Pass

The prices applied to the tin and derivatives vary with the price of the raw material.

Tin metal foil (0.10 mm)
pure

packaging	Mfr. No
100 g Plastic jar	T/1560/48
500 g Plastic jar	T/1560/53
Sn	Boiling point: 2270°C
CAS: 7440-31-5	Melting point: 231.9°C
MW: 118.69	EINECS: 231-141-8
Density at 20°C: 7.31	
Guaranteed analysis	
Identity test	: Pass

The prices applied to the tin and derivatives vary with the price of the raw material.

Tin metal
pure

packaging	Mfr. No
100 g Plastic jar	T/1600/48
500 g Plastic jar	T/1600/53
Sn	Boiling point: 2270°C
CAS: 7440-31-5	Melting point: 231.9°C
MW: 118.69	EINECS: 231-141-8
Density at 20°C: 7.31	
Guaranteed analysis	
Identity test	: Pass

The prices applied to the tin and derivatives vary with the price of the raw material.

Tin(II) chloride anhydrous
specified reagents for general laboratory work- extra pure

packaging	Mfr. No
100 g Plastic jar	T/1645/48
SnCl ₂	UN: 3260
CAS: 7772-99-8	GHS: Acute Tox. 4, Skin
MW: 189.60	Corr. 1B, Eye Dam. 1, Skin
Density at 20°C: 3.95	Sens. 1
Boiling point: 652°C	H302, H314, H317
Melting point: 246°C	P301+P330+P331,
ADR/RID: 8, III	P305+P351+P338, P310,
EINECS: 231-868-0	P280, P301+P312
Guaranteed analysis	
Identity test	: Pass

Tin(II) chloride dihydrate
for analysis

packaging	Mfr. No
100 g Plastic jar	T/1654/48
250 g Plastic jar	T/1654/50
500 g Plastic jar	T/1654/53
1 kg Plastic jar	T/1654/60
SnCl ₂ ·2H ₂ O	GHS: Acute Tox. 4, Skin
CAS: 10025-69-1	Corr. 1B, Eye Dam. 1, Skin
MW: 225.63	Sens. 1
Boiling point: 652°C	H302, H314, H317
Melting point: 37 to 38°C	P280, P301+P330+P331,
ADR/RID: 8, III	P305+P351+P338,
EINECS: 231-868-0	P302+P350, P310
UN: 3260	
Assay: >97%	
Guaranteed analysis	
Calcium (Ca) <50 ppm	Potassium (K) <20 ppm
Copper (Cu) <50 ppm	Sodium (Na) <50 ppm
Iron (Fe) <50 ppm	Total silicon (Si) <50 ppm
Lead (Pb) <50 ppm	Total sulphur (S) <50 ppm
Magnesium (Mg) <50 ppm	Zinc (Zn) <50 ppm
Typical analysis	
Aluminium (Al) <1 ppm	Manganese (Mn) <1 ppm
Arsenic (As) <2 ppm	Molybdenum (Mo) <1 ppm
Barium (Ba) <1 ppm	Nickel (Ni) <1 ppm
Cadmium (Cd) <1 ppm	Silver (Ag) <1 ppm
Chromium (Cr) <1 ppm	Strontium (Sr) <1 ppm
Cobalt (Co) <1 ppm	Titanium (Ti) <1 ppm
Lithium (Li) <1 ppm	Vanadium (V) <1 ppm
Corrosive when wet.	
Lot analysis on the label.	

Tin(II) chloride dihydrate
specified reagents for general laboratory work- extra pure

packaging	Mfr. No
100 g Plastic jar	T/1652/48
250 g Plastic jar	T/1652/50
500 g Plastic jar	T/1652/53
1 kg Plastic jar	T/1652/60
SnCl ₂ ·2H ₂ O	UN: 3260
CAS: 10025-69-1	GHS: Acute Tox. 4, Skin
MW: 225.63	Corr. 1B, Eye Dam. 1, Skin
Density at 20°C: 2.71	Sens. 1
Boiling point: 652°C	H302, H314, H317
Melting point: 37 to 38°C	P280, P301+P330+P331,
ADR/RID: 8, III	P305+P351+P338,
EINECS: 231-868-0	P302+P350, P310
Assay: >96%	
Guaranteed analysis	
Calcium (Ca) <200 ppm	Potassium (K) <100 ppm
Copper (Cu) <200 ppm	Sodium (Na) <200 ppm
Iron (Fe) <200 ppm	Total silicon (Si) <200 ppm
Lead (Pb) <200 ppm	Total sulphur (S) <200 ppm
Magnesium (Mg) <200 ppm	Zinc (Zn) <200 ppm
Corrosive when wet.	

Tin(IV) chloride anhydrous
pure

packaging	Mfr. No
100 mℓ Coated Glass bottle	T/1665/PB05
SnCl ₄	UN: 1827
CAS: 7646-78-8	GHS: Skin Corr. 1B, Eye
MW: 260.52	Dam. 1, STOT SE 3, Aquatic
Density at 20°C: 2.22	Chronic 3
Boiling point: 114°C	H412, H314
Melting point: -33°C	P301+P330+P331, P280,
ADR/RID: 8, II	P305+P351+P338,
EINECS: 231-588-9	P342+P311, P261, P273

Tin(IV) chloride pentahydrate
pure

packaging	Mfr. No
100 g Plastic jar	T/1661/48
500 g Plastic jar	T/1661/53
SnCl ₄ ·5H ₂ O	UN: 2440
CAS: 10026-06-9	GHS: Skin Corr. 1B, Eye
MW: 350.58	Dam. 1, Aquatic Chronic 3
Density at 20°C: 2.2(H ₂ O=1)	H314, H412
Melting point: 56°C	P301+P330+P331, P280,
ADR/RID: 8, III	P305+P351+P338, P310, P273
EINECS: 231-588-9	
Guaranteed analysis	
Identity test	: Pass

Tin(IV) oxide
pure

packaging	Mfr. No
100 g Plastic jar	T/1670/48
SnO ₂	Density at 20°C: 6.95
CAS: 18282-10-5	Melting point: 1127°C
MW: 150.69	EINECS: 242-159-0
Guaranteed analysis	
Identity test	: Pass

TISAB solution
pure

packaging	Mfr. No
5 ℓ HDPE Plastic bottle	F/1325/21
Density at 20°C: 1.04	
Guaranteed analysis	
pH at 20°C	: 5 to 7
Contents per Litre:	
Acetic acid glacial: 12.6 mℓ	
trans-1.2-Diaminocyclohexane-N,N,N',N'-tetra acetic acid: 0.1 g	
Sodium chloride: 58 g	
Sodium hydroxyde AR: 8 g	

Titanium standard metal solution 1000 ppm
for AAS, solution HCl 2M

packaging	Mfr. No
100 mℓ HDPE Plastic bottle	J/8110/05
500 mℓ HDPE Plastic bottle	J/8110/08
Manufactured to a tolerance of ± 0.5%.	
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.	

Titanium standard metal solution 1000 ppm
for ICP, solution HCl 2M

packaging	Mfr. No
500 mℓ HDPE Plastic bottle	J/8331/08
Manufactured to a tolerance of ± 0.5%. All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities. See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.	

Titanium standard metal solution 10 000 ppm
for ICP, solution HCl 2M

packaging	Mfr. No
100 mℓ HDPE Plastic bottle	J/8330/05
Manufactured to a tolerance of ± 0.5%. All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities. See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.	

Titanium standard metal solution 1000 ppm
for ICP solution in HCL 20%, NIST certified

packaging	Mfr. No
500 mℓ Plastic bottle	SPPLTI52X
H315, H319, H335	

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Titanium standard metal solution 1000 ppm
for ICP solution in H₂O/HF trace, NIST certified

packaging	Mfr. No
100 mℓ Plastic bottle	SPPLTI92Y
250 mℓ Plastic bottle	SPPLTI92T
500 mℓ Plastic bottle	SPPLTI92X

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Titanium standard metal solution 10 000 ppm
for ICP solution in HCL 40%, NIST certified

packaging	Mfr. No
500 mℓ Plastic bottle	SPPLTI53X
H314, H335	

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Titanium standard metal solution 10 000 ppm
for ICP solution in H₂O/HF trace, NIST certified

packaging	Mfr. No
100 mℓ Plastic bottle	SPPLTI93Y
500 mℓ Plastic bottle	SPPLTI93X

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Titanium standard metal solution 1000 ppm
for ICP-MS solution in H₂O/HF trace, NIST certified

packaging	Mfr. No
100 mℓ Plastic bottle	SPCLTI92Y

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Titanium
pure, metal, crushed sponge

packaging	Mfr. No
100 g Plastic jar	T/1840/48
Ti	
CAS: 7440-32-6	EINECS: 231-142-3
MW: 47.88	UN: 2546
Density at 20°C: 4.51	GHS: Flam. Sol. 2, Pyr. Sol. 1
Boiling point: 3287°C	H228, H250
Melting point: 1660°C	P210, P222, P280,
ADR/RID: 4.2, II	P335+P334, P370+P378
Guaranteed analysis	
Identity test : Pass	

Titanium(III) chloride
for amino acid analysis, 15% solution w/v

packaging	Mfr. No
500 mℓ Coated Glass bottle	T/1958/PB08
TiCl ₃	
CAS: 7705-07-9	GHS: Skin Corr. 1B, Eye
MW: 154.26	Dam. 1, STOT SE 3
Density at 20°C: 1.21	H314, H335
ADR/RID: 8, II	P301+P330+P331,
EINECS: 231-728-9	P305+P351+P338, P310,
UN: 2869	P280, P402+P404, P302+P350

Titanium(III) chloride solution 12.5 % w/w TiCl₃
Technical

packaging	Mfr. No
500 mℓ Coated Glass bottle	T/1950/PB08
TiCl ₃	
CAS: 7705-07-9	UN: 2869
MW: 154.26	GHS: Skin Corr. 1B, Eye
Density at 20°C: 1.14 to 1.16	Dam. 1, STOT SE 3
Melting point: 440°C	H314, H335
ADR/RID: 8, II	P301+P330+P331,
EINECS: 231-728-9	P305+P351+P338, P310,
	P280, P402+P404, P302+P350

Titanium(III) chloride solution 0.2M (0.2N)
Standard solution ready to use, For volumetric analysis

packaging	Mfr. No
2.5 ℓ Glass bottle	J/8780/17
ADR/RID: 8, III	
UN: 3264	H319, H335, H315
GHS: Skin Irrit. 2, Eye Irrit. 2,	P261, P280, P302+P352,
STOT SE 3	P304+P340, P305+P351+P338
Manufactured to a tolerance of ± 0.1%. See the full range of Standard Volumetric solution page 323.	

Titanium(IV) dioxide
specified reagents for general laboratory work- extra pure

packaging	Mfr. No
500 g Plastic jar	T/1900/53
TiO ₂	
CAS: 13463-67-7	Boiling point: 2900°C
MW: 79.88	Melting point: 1855°C
Density at 20°C: 4.23	EINECS: 236-675-5
Assay: >98%	
Guaranteed analysis	
Loss on ignition	<0.5% Iron (Fe) <500 ppm

Titanium(IV) sulfate solution 15 % w/v approx
Technical

packaging	Mfr. No
1 ℓ Coated Glass bottle	T/1800/PB15
KOH	
CAS: 13693-11-3	UN: 2796
Density at 20°C: 1.12	GHS: Skin Corr. 1A
ADR/RID: 8, II	H314
EINECS: 237-215-6	P301+P330+P331, P280,
	P305+P351+P338, P310
Assay: 14 to 16%	

Titanium(III) sulfate solution approx.15 % w/v Ti₂ (SO₄)₃
Technical

packaging	Mfr. No
500 mℓ Coated Glass bottle	T/2000/PB08
CAS: 10343-61-0	
Density at 20°C: 1.1 to 1.3	GHS: Met. Corr. 1, Skin Irrit.
ADR/RID: 8, III	2, Eye Irrit. 2
EINECS: 233-749-9	H290, H315, H319
UN: 1760	P261, P280, P302+P352,
	P305+P351+P338, P390

Titan yellow C.I. 19540
pure, indicator grade

packaging	Mfr. No
25 g Glass jar	T/2060/46
C ₂₈ H ₁₉ N ₅ Na ₂ O ₆ S ₄	
CAS: 1829-00-1	GHS: Acute Tox. 4
MW: 695.71	H312, H332, H302
EINECS: 217-377-4	P280, P302+P352,
	P301+P312, P304+P340
Guaranteed analysis	
Complies with standard IR spectrum: Pass	
pH 12.0-13.0 yellow to red. Suitable for magnesium determination.	

Toluene
for HPLC

packaging	Mfr. No
1 ℓ Glass bottle	T/2306/15
2.5 ℓ Glass bottle	T/2306/17
10 ℓ Returnable Drum	T/2306/10RSS
25 ℓ Returnable Drum	T/2306/25RSS
30 ℓ Returnable Drum	T/2306/30RSS
C ₇ H ₈	
CAS: 108-88-3	UN: 1294
MW: 92.14	GHS: Flam. Liq. 2, Repr. 2,
Density at 20°C: 0.866	Skin Irrit. 2, STOT SE 3,
Boiling point: 111°C	STOT RE 2, Asp. Tox. 1
Melting point: -95°C	H225, H336, H373, H304,
Flash point: 4°C	H315, H361d
ADR/RID: 3, II	P281, P260, P331,
EINECS: 203-625-9	P301+P310, P304+P340,
	P302+P352, P210

Assay (CG): >99.8%

Guaranteed analysis

Acidity / Alkalinity	<0.0015 meq/g	Residue after evaporation	<2 ppm
Water	<0.03%		

Absorbance specification

Max. Absorbance (1cm cell v HPLC water)

290 nm	<0.2 A.U.	330 nm	<0.005 A.U.
290 nm (sol. aqu. 1%)	: <0.05 A.U.	340 nm	<0.005 A.U.
310 nm	<0.025 A.U.	350 nm	<0.005 A.U.
320 nm	<0.01 A.U.		


Lot analysis on the label.
Filtered to 0.2 µm.

Toluene
for Gel Permeation Chromatography (GPC)

packaging	Mfr. No
2.5 ℓ Coated Glass bottle	T/2310/PB17
C ₇ H ₈	
CAS: 108-88-3	UN: 1294
MW: 92.14	GHS: Flam. Liq. 2, Repr. 2,
Density at 20°C: 0.866	Skin Irrit. 2, STOT SE 3,
Boiling point: 111°C	STOT RE 2, Asp. Tox. 1
Melting point: -95°C	H225, H336, H373, H304,
Flash point: 4°C	H315, H361d
ADR/RID: 3, II	P281, P260, P331,
EINECS: 203-625-9	P301+P310, P304+P340,
	P302+P352, P210

Guaranteed analysis

Colour	<10 APHA	Residue after evaporation	<2 ppm
Water	<0.03%		



Toluene

for trace metal analysis

Primar

packaging

Mfr. No

2.5 ℓ Glass bottle

T/2301/17

C₇H₈

CAS: 108-88-3

MW: 92.14

Density at 20°C: 0.866

Boiling point: 111°C

Melting point: -95°C

Flash point: 4°C

ADR/RID: 3, II

EINECS: 203-625-9

UN: 1294

GHS: Flam. Liq. 2, Repr. 2,
Skin Irrit. 2, STOT SE 3,
STOT RE 2, Asp. Tox. 1

H225, H336, H373, H304,
H315, H361d

P281, P260, P331,
P301+P310, P304+P340,
P302+P352, P210







Assay (CG): >99.9%

Guaranteed analysis

Acidity / Alkalinity	<0.0002 meq/g	Manganese (Mn)	<0.005 ppm
Colour	<5 APHA	Molybdenum (Mo)	<0.005 ppm
Water	<0.05%	Nickel (Ni)	<0.005 ppm
Residue after evaporation	<2 ppm	Potassium (K)	<0.05 ppm
Aluminium (Al)	<0.02 ppm	Propan--1-ol	<0.005 ppm
Barium (Ba)	<0.005 ppm	Sodium (Na)	<0.2 ppm
Cadmium (Cd)	<0.005 ppm	Strontium (Sr)	<0.005 ppm
Calcium (Ca)	<0.05 ppm	Titanium (Ti)	<0.005 ppm
Chromium (Cr)	<0.005 ppm	Total phosphorus (P)	<0.02 ppm
Cobalt (Co)	<0.005 ppm	Total silicon (Si)	<0.02 ppm
Copper (Cu)	<0.01 ppm	Total sulphur (S)	<0.1 ppm
Iron (Fe)	<0.05 ppm	Vanadium (V)	<0.005 ppm
Lead (Pb)	<0.005 ppm	Zinc (Zn)	<0.01 ppm
Lithium (Li)	<0.005 ppm	Benzene	<0.02%
Magnesium (Mg)	<0.02 ppm	Xylenes	<0.02%

Toluene		Distol
for residue and pesticide analysis		
packaging	Mfr. No	
1 ℓ Glass bottle	T/2323/15	
2.5 ℓ Glass bottle	T/2323/17	
C ₇ H ₈	UN: 1294	
CAS: 108-88-3	GHS: Flam. Liq. 2, Repr. 2,	
MW: 92.14	Skin Irrit. 2, STOT SE 3,	
Density at 20°C: 0.866	STOT RE 2, Asp. Tox. 1	
Boiling point: 111°C	H225, H336, H373, H304,	
Melting point: -95°C	H315, H361d	
Flash point: 4°C	P281, P260, P331,	
ADR/RID: 3, II	P301+P310, P304+P340,	
EINECS: 203-625-9	P302+P352, P210	
Assay (CG): >99.8%		
Guaranteed analysis		
Water	<0.05%	GC (FID) Analysis Pass
Residue after evaporation	<2 ppm	GC (NPD) analysis Pass
GC (ECD) Analysis	Pass	
GC Analysis		
Analysis of pesticides: from the retention time of Lindane, after concentration of		

Optima LC/MS

Fisher Chemical OPTIMA® LC/MS solvents, blends, additives and modifiers have set the standard of excellence for consistent, reproducible performance in the mobile phase of LC/MS.

Find out more on page 308



Toluene		
extra dry, for synthesis		
packaging		Mfr. No
1 ℓ	Coated Glass bottle	T/2302/PB15
2.5 ℓ	Coated Glass bottle	T/2302/PB17
C ₇ H ₈	UN: 1294	
CAS: 108-88-3	GHS: Flam. Liq. 2, Repr. 2,	
MW: 92.14	Skin Irrit. 2, STOT SE 3,	
Density at 20°C: 0.866	STOT RE 2, Asp. Tox. 1	
Boiling point: 111°C	H225, H336, H373, H304,	
Melting point: -95°C	H315, H361d	
Flash point: 4°C	P281, P260, P331,	
ADR/RID: 3, II	P301+P310, P304+P340,	
EINECS: 203-625-9	P302+P352, P210	
Assay (CG): >99.8%		
Guaranteed analysis		
Water	<0.01%	Residue after evaporation <2 ppm

Toluene		NEW PACK.	Certified AR
for analysis			
packaging		Mfr. No	
1 ℓ	Glass bottle	T/2300/15	
1 ℓ	Coated Glass bottle	T/2300/PB15	
2.5 ℓ	Glass bottle	T/2300/17	
2.5 ℓ	Coated Glass bottle	T/2300/PB17	
5 ℓ	Fluorinated HDPE Plastic bottle	T/2300/21	
25 ℓ	Metal drum	T/2300/25	
200 ℓ	Metal drum	T/2300/27	
C ₇ H ₈		UN: 1294	  
CAS: 108-88-3		GHS: Flam. Liq. 2, Repr. 2,	
MW: 92.14		Skin Irrit. 2, STOT SE 3,	
Density at 20°C: 0.866		STOT RE 2, Asp. Tox. 1	
Boiling point: 111°C		H225, H336, H373, H304,	
Melting point: -95°C		H315, H361d	
Flash point: 4°C		P281, P260, P331,	
ADR/RID: 3, II		P301+P310, P304+P340,	
EINECS: 203-625-9		P302+P352, P210	
Assay (CG): >99.8%			
Guaranteed analysis			
Acidity / Alkalinity	<0.0001 meq/g	Lead (Pb)	<0.02 ppm
Colour	<10 APHA	Magnesium (Mg)	<0.05 ppm
Water	<0.03%	Potassium (K)	<0.5 ppm
Residue after evaporation	<10 ppm	Sodium (Na)	<1 ppm
Substances darkened		Total phosphorus (P)	<0.1 ppm
by (H ₂ SO ₄)	<50 APHA	Total silicon (Si)	<0.05 ppm
Substances		Total sulphur (S)	<10 ppm
reducing KMnO ₄	<0.0005%	Zinc (Zn)	<0.05 ppm
Calcium (Ca)	<0.1 ppm	Benzene	<0.05%
Copper (Cu)	<0.02 ppm	Xylenes	<0.05%
Iron (Fe)	<0.1 ppm		
Typical analysis			
Aluminium (Al)	<0.02 ppm	Manganese (Mn)	<0.02 ppm
Barium (Ba)	<0.02 ppm	Molybdenum (Mo)	<0.02 ppm
Cadmium (Cd)	<0.02 ppm	Nickel (Ni)	<0.02 ppm
Chromium (Cr)	<0.02 ppm	Strontium (Sr)	<0.02 ppm
Cobalt (Co)	<0.02 ppm	Titanium (Ti)	<0.02 ppm
Lithium (Li)	<0.02 ppm	Vanadium (V)	<0.02 ppm

Lot analysis on the label.

High Volume Solvent Delivery Systems

High-volume solvent delivery systems are available in 10L to 1000L, in Europe only. Please email your contact details and a brief description of your requirements to bege1.sha.quotes@thermofisher.com



Toluene		SLR
specified reagents for general laboratory work- extra pure, low in sulfur		
packaging		Mfr. No
2.5 ℓ	Glass bottle	T/2250/17
2.5 ℓ	Coated Glass bottle	T/2250/PB17
25 ℓ	Metal drum	T/2250/25
C ₇ H ₈	UN: 1294	
CAS: 108-88-3	GHS: Flam. Liq. 2, Repr. 2,	
MW: 92.14	Skin Irrit. 2, STOT SE 3,	
Density at 20°C: 0.866	STOT RE 2, Asp. Tox. 1	
Boiling point: 111°C	H225, H336, H373, H304,	
Melting point: -95°C	H315, H361d	
Flash point: 4°C	P281, P260, P331,	
ADR/RID: 3, II	P301+P310, P304+P340,	
EINECS: 203-625-9	P302+P352, P210	
Guaranteed analysis		
Residue after evaporation	<50 ppm	Sulfur compounds (SC ₂)
		<0.0005%

Toluene NEW PACK.		SLR
specified reagents for general laboratory work- extra pure		
packaging		Mfr. No
500 ml	Glass bottle	T/2200/08
2.5 l	Glass bottle	T/2200/17
2.5 l	Coated Glass bottle	T/2200/PB17
5 l	Fluorinated HDPE Plastic bottle	T/2200/21
5 l	Fluorinated HDPE square Drum	T/2200/FP21
25 l	Metal drum	T/2200/25
200 l	Metal drum	T/2200/27
C ₇ H ₈	UN: 1294	
CAS: 108-88-3	GHS: Flam. Liq. 2, Repr. 2,	
MW: 92.14	Skin Irrit. 2, STOT SE 3,	
Density at 20°C: 0.866	STOT RE 2, Asp. Tox. 1	
Boiling point: 111°C	H225, H336, H373, H304,	
Melting point: -95°C	H315, H361d	
Flash point: 4°C	P281, P260, P331,	
ADR/RID: 3, II	P301+P310, P304+P340,	
EINECS: 203-625-9	P302+P352, P210	
Assay (CG): >99%		
Guaranteed analysis		
Acidity / Alkalinity	<0.0005 meq/g	Magnesium (Mg)
Colour	<10 APHA	Potassium (K)
Water	<0.05%	Sodium (Na)
Residue after evaporation	<30 ppm	Total phosphorus (P)
Substances		Total silicon (Si)
reducing KMnO ₄	<0.001%	Total sulphur (S)
Calcium (Ca)	<5 ppm	Zinc (Zn)
Copper (Cu)	<2 ppm	Substances darkened
Iron (Fe)	<5 ppm	by H ₂ SO ₄
Lead (Pb)	<2 ppm	

Toluene		
Technical		
packaging	Mfr. No	
1 ℓ Coated Glass bottle	T/2150/PB15	
2.5 ℓ Glass bottle	T/2150/17	
2.5 ℓ Coated Glass bottle	T/2150/PB17	
25 ℓ Metal drum	T/2150/25	
200 ℓ Metal drum	T/2150/27	
C ₇ H ₈	UN: 1294	
CAS: 108-88-3	GHS: Flam. Liq. 2, Repr. 2,	
MW: 92.14	Skin Irrit. 2, STOT SE 3,	
Density at 20°C: 0.866	STOT RE 2, Asp. Tox. 1	
Boiling point: 111°C	H225, H336, H373, H304,	
Melting point: -95°C	H315, H361d	
Flash point: 4°C	P281, P260, P331,	
ADR/RID: 3, II	P301+P310, P304+P340,	
EINECS: 203-625-9	P302+P352, P210	

Toluene-4-sulfonic acid		SLR	
specified reagents for general laboratory work- extra pure			
packaging		Mfr. No	
100 g	Plastic jar	T/2500/48	
500 g	Plastic jar	T/2500/53	
C ₇ H ₆ SO ₃ ·H ₂ O	UN: 2585		
CAS: 6192-52-5	GHS: Skin Irrit. 2, Eye Irrit. 2,		
MW: 190.22	STOT SE 3		
Melting point: 100 to 108°C	H315, H335, H319		
Flash point: 180°C	P261, P302+P352, P280,		
ADR/RID: 8, III	P305+P351+P338		
EINECS: 203-180-0			
Assay: >98%			
Guaranteed analysis			
Calcium (Ca)	<500 ppm	Potassium (K)	<2000 ppm
Copper (Cu)	<50 ppm	Sodium (Na)	<500 ppm
Iron (Fe)	<100 ppm	n-Butyl formate	<0.02%
Lead (Pb)	<50 ppm	Total silicon (Si)	<200 ppm
Magnesium (Mg)	<200 ppm	Zinc (Zn)	<100 ppm

p-Toluidine		SLR
specified reagents for general laboratory work- extra pure		
packaging	Mfr. No	
500 g Glass jar	T/2650/53	
C ₇ H ₉ N	UN: 3451	
CAS: 106-49-0	GHS: Carc. 2, Acute Tox. 3,	
MW: 107.16	Eye Irrit. 2, Skin Sens. 1,	
Density at 20°C: 1.05	Aquatic Acute 1	
Boiling point: 200°C	H311, H351, H319, H317,	
Melting point: 41 to 46°C	H331, H400, H301	
Flash point: 87°C	P280, P312, P302+P350,	
ADR/RID: 6.1, II	P304+P340, P273,	
EINECS: 203-403-1	P305+P351+P338	
Assay (CG): >98%		

Toluidine blue C.I. 52040		
pure		
packaging	Mfr. No	
25 g Glass jar	T/2660/46	
C ₁₅ H ₁₆ N ₃ SCI	GHS: Acute Tox. 4	
CAS: 92-31-9	H302	
MW: 305.83	P301+P312	
EINECS: 202-146-2		
Guaranteed analysis		
Complies with standard IR spectrum: Pass		

Total hardness indicator tablets

pure, indicator grade

packaging	Mfr. No
1000 PC Glass jar	T/2750/79
Guaranteed analysis	
Identity test	: Pass
For use as an indicator in conjunction with EDTA solution.	

Standard Volumetric solution

For your convenience we offer a range of NIST traceable solutions. A full list of our products can be found on page 323. Can't see what you need? Please email your contact details and a brief description of your requirements to bege1.sha.quotes@thermofisher.com



Trichloroacetic acid
for electrophoresis

packaging	Mfr. No
250 g Plastic jar	T/P500/50
1 kg Plastic jar	T/P500/60

C₂HCl₃O₂
CAS: 76-03-9
MW: 163.39
Density at 20°C: 1.62
Boiling point: 196°C
Melting point: 52 to 58°C
ADR/RID: 8, II
EINECS: 200-927-2

UN: 1839
GHS: Skin Corr. 1A, STOT SE 3, Aquatic Acute 1, Aquatic Chronic 1
H314, H335, H410
P301+P330+P331, P280, P305+P351+P338, P310, P261, P273

Assay: >99.8%			
Guaranteed analysis			
Iron (Fe)	<5 ppm	Total phosphorus (P)	<10 ppm
Lead (Pb)	<2 ppm	Total sulphur (S)	<200 ppm
n-Butyl formate	<0.001%		
<i>Suitable for use in fixing solutions and as a protein precipitation reagent.</i>			

Suitable for use in fixing solutions and as a protein precipitation reagent.

Trichloroacetic acid
for analysis

Certified AR

packaging	Mfr. No
250 g Plastic jar	T/3000/50
500 g Plastic jar	T/3000/53
1 kg Plastic jar	T/3000/60

C₂HCl₃O₂
CAS: 76-03-9
MW: 163.39
Density at 20°C: 1.62
Boiling point: 196°C
Melting point: 52 to 58°C
ADR/RID: 8, II
EINECS: 200-927-2

UN: 1839
GHS: Skin Corr. 1A, STOT SE 3, Aquatic Acute 1, Aquatic Chronic 1
H314, H335, H410
P301+P330+P331, P280, P305+P351+P338, P310, P261, P273

Assay: >99.5%			
Guaranteed analysis			
Substances reducing KMnO ₄	<0.0002%	Sodium (Na)	<20 ppm
Calcium (Ca)	<10 ppm	Total chloride (Cl)	<0.001%
Copper (Cu)	<2 ppm	Total nitrogen (N)	<50 ppm
Iron (Fe)	<5 ppm	Total phosphorus (P)	<10 ppm
Lead (Pb)	<2 ppm	Total silicon (Si)	<20 ppm
Magnesium (Mg)	<10 ppm	Total sulphur (S)	<200 ppm
Potassium (K)	<10 ppm	Zinc (Zn)	<5 ppm
Typical analysis			
Aluminium (Al)	<2 ppm	Molybdenum (Mo)	<2 ppm
Barium (Ba)	<2 ppm	Nickel (Ni)	<2 ppm
Cadmium (Cd)	<2 ppm	Silver (Ag)	<2 ppm
Chromium (Cr)	<2 ppm	Strontium (Sr)	<2 ppm
Cobalt (Co)	<2 ppm	Titanium (Ti)	<2 ppm
Lithium (Li)	<2 ppm	Vanadium (V)	<2 ppm
Manganese (Mn)	<2 ppm		

Lot analysis on the label.

Trichloroacetic acid
specified reagents for general laboratory work- extra pure

SLR

packaging	Mfr. No
250 g Plastic jar	T/2950/50
500 g Plastic jar	T/2950/53
1 kg Plastic jar	T/2950/60

C₂HCl₃O₂
CAS: 76-03-9
MW: 163.39
Density at 20°C: 1.62
Boiling point: 196°C
Melting point: 52 to 58°C
ADR/RID: 8, II
EINECS: 200-927-2

UN: 1839
GHS: Skin Corr. 1A, STOT SE 3, Aquatic Acute 1, Aquatic Chronic 1
H314, H335, H410
P301+P330+P331, P280, P305+P351+P338, P310, P261, P273

Assay: >99%			
Guaranteed analysis			
Calcium (Ca)	<50 ppm	Sodium (Na)	<100 ppm
Copper (Cu)	<10 ppm	n-Butyl formate	<0.005%
Iron (Fe)	<20 ppm	Total phosphorus (P)	<50 ppm
Lead (Pb)	<10 ppm	Total silicon (Si)	<100 ppm
Magnesium (Mg)	<50 ppm	Total sulphur (S)	<500 ppm
Potassium (K)	<50 ppm	Zinc (Zn)	<20 ppm

Trichloroethylene

for analysis, stabilized with triethylamine

packaging	Mfr. No
500 mℓ Glass bottle	T/3101/08
2.5 ℓ Coated Glass bottle	T/3101/PB17

C₂HCl₃
CAS: 79-01-6
MW: 131.39
Density at 20°C: 1.46
Boiling point: 87°C
Melting point: -86°C
ADR/RID: 6.1, III
EINECS: 201-167-4
UN: 1710

GHS: Carc. 1B, Muta. 2, Skin Irrit. 2, Eye Irrit. 2, STOT SE 3, Aquatic Chronic 3
H336, H412, H350, H341, H319, H315
P201, P261, P273, P302+P352, P280, P305+P351+P338

Assay (CG): >99.8%			
Guaranteed analysis			
Acidity / Alkalinity	<0.0003 meq/g	Lithium (Li)	<0.02 ppm
Colour	<10 APHA	Magnesium (Mg)	<0.2 ppm
Water	<0.02%	Manganese (Mn)	<0.02 ppm
Residue after evaporation	<10 ppm	Molybdenum (Mo)	<0.02 ppm
Free Chlorine (Cl ₂)	<5 ppm	Nickel (Ni)	<0.02 ppm
Aluminium (Al)	<0.02 ppm	Potassium (K)	<0.5 ppm
Barium (Ba)	<0.02 ppm	Sodium (Na)	<0.5 ppm
Cadmium (Cd)	<0.02 ppm	Strontium (Sr)	<0.02 ppm
Calcium (Ca)	<0.5 ppm	Titanium (Ti)	<0.02 ppm
Chlorine	<5 ppm	Total phosphorus (P)	<0.1 ppm
Chromium (Cr)	<0.02 ppm	Total silicon (Si)	<0.05 ppm
Cobalt (Co)	<0.02 ppm	Total sulphur (S)	<0.5 ppm
Copper (Cu)	<0.05 ppm	Vanadium (V)	<0.02 ppm
Iron (Fe)	<0.05 ppm	Zinc (Zn)	<0.1 ppm
Lead (Pb)	<0.05 ppm	Triethylamine	0.02 to 0.04%

Lot analysis on the label.
Stabilised with 0.04% Triethylamine

Trichloroethylene

extra pure, stabilized with triethylamine

NEW PACK.

SLR

packaging	Mfr. No
500 mℓ Glass bottle	T/3100/08
2.5 ℓ Glass bottle	T/3100/17
2.5 ℓ Coated Glass bottle	T/3100/PB17
5 ℓ Fluorinated HDPE Plastic bottle	T/3100/21
25 ℓ Metal drum	T/3100/25
200 ℓ Metal drum	T/3100/27

C ₂ HCl ₃ CAS: 79-01-6 MW: 131.39 Density at 20°C: 1.46 Boiling point: 87°C Melting point: -86°C ADR/RID: 6.1, III EINECS: 201-167-4 UN: 1710	GHS: Carc. 1B, Muta. 2, Skin Irrit. 2, Eye Irrit. 2, STOT SE 3, Aquatic Chronic 3 H336, H412, H350, H341, H319, H315 P201, P261, P273, P302+P352, P280, P305+P351+P338	
---	--	--

Assay (CG): >99%			
Guaranteed analysis			
Acidity / Alkalinity	<0.002 meq/g	Lead (Pb)	<2 ppm
Colour	<10 APHA	Magnesium (Mg)	<2 ppm
Water	<0.05%	Potassium (K)	<2 ppm
Residue after evaporation	<50 ppm	Sodium (Na)	<5 ppm
Free Chlorine (Cl ₂)	<10 ppm	Total phosphorus (P)	<2 ppm
Calcium (Ca)	<5 ppm	Total silicon (Si)	<2 ppm
Chlorine	<10 ppm	Total sulphur (S)	<2 ppm
Copper (Cu)	<2 ppm	Zinc (Zn)	<2 ppm
Iron (Fe)	<2 ppm	Triethylamine	0.02 to 0.03%

Stabilised with 0.04% Triethylamine.

Trichloroethylene

Technical, stabilized

packaging	Mfr. No
500 mℓ Glass bottle	T/3050/08
2.5 ℓ Coated Glass bottle	T/3050/PB17
25 ℓ Metal drum	T/3050/25

C ₂ HCl ₃ CAS: 79-01-6 MW: 131.39 Density at 20°C: 1.46 Boiling point: 87°C Melting point: -86°C ADR/RID: 6.1, III EINECS: 201-167-4 UN: 1710	GHS: Carc. 1B, Muta. 2, Skin Irrit. 2, Eye Irrit. 2, STOT SE 3, Aquatic Chronic 3 H336, H412, H350, H341, H319, H315 P201, P261, P273, P302+P352, P280, P305+P351+P338	
---	--	--

Tricine

for biology ,N-Tris(hydroxymethyl) methyl glycine , 99+%

packaging	Mfr. No
100 g Plastic jar	T/P505/48

C₆H₁₃NO₃
CAS: 5704-04-1
MW: 179.17
Melting point: 182 to 184°C
EINECS: 227-193-6

GHS: Skin Irrit. 2, Eye Irrit. 2, STOT SE 3
H315, H319, H335
P261, P302+P352, P280, P305+P351+P338

Assay: >99%			
Guaranteed analysis			
Loss on drying	<0.2%	Lead (Pb)	<10 ppm
Iron (Fe)	<5 ppm	Sulfate (SO ₄)	<0.01%

Absorbance specification	
Max. Absorbance (1cm cell v HPLC water)	
250 nm (sol. aqu. 10%)	: <0.1 A.U.

Triethanolamine

for analysis

Certified AR

packaging	Mfr. No
500 mℓ Glass bottle	T/3160/08
2.5 ℓ Glass bottle	T/3160/17

C₆H₁₅NO₃
CAS: 102-71-6
MW: 149.19
Density at 20°C: 1.125
Boiling point: 360°C
Melting point: 21°C
Flash point: 190°C

ADR/RID: 8/III
EINECS: 203-049-8
UN: 3259
GHS: Eye Irrit. 2
H319
P280, P305+P351+P338

Assay (CG): >99%			
Guaranteed analysis			
Water	<0.2%	Sodium (Na)	<10 ppm
Calcium (Ca)	<2 ppm	Total phosphorus (P)	<0.5 ppm
Chloride (Cl)	<0.0001%	Total silicon (Si)	<5 ppm
Copper (Cu)	<0.5 ppm	Total sulphur (S)	<1 ppm
Iron (Fe)	<1 ppm	Zinc (Zn)	<0.5 ppm
Lead (Pb)	<0.5 ppm	Diethanolamine	<0.5%
Magnesium (Mg)	<1 ppm	Ethanolamine	<0.1%
Potassium (K)	<5 ppm	Residue after ignition	<0.005%

Typical analysis			
Aluminium (Al)	<0.5 ppm	Manganese (Mn)	<0.5 ppm
Barium (Ba)	<0.5 ppm	Molybdenum (Mo)	<0.5 ppm
Cadmium (Cd)	<0.5 ppm	Nickel (Ni)	<0.5 ppm
Chromium (Cr)	<0.5 ppm	Strontium (Sr)	<0.5 ppm
Cobalt (Co)	<0.5 ppm	Titanium (Ti)	<0.5 ppm
Lithium (Li)	<0.5 ppm	Vanadium (V)	<0.5 ppm

Lot analysis on the label.

Triethanolamine

specified reagents for general laboratory work- extra pure

SLR

packaging	Mfr. No
500 mℓ Glass bottle	T/3150/08
2.5 ℓ Glass bottle	T/3150/17
25 ℓ Metal drum	T/3150/25

C ₆ H ₁₅ NO ₃ CAS: 102-71-6 MW: 149.19 Density at 20°C: 1.125 Boiling point: 360°C Melting point: 21°C Flash point: 190°C	ADR/RID: 8/III EINECS: 203-049-8 UN: 3259 GHS: Eye Irrit. 2 H319 P280, P305+P351+P338	
--	--	--

Assay (CG): >97%			
Guaranteed analysis			
Water	<1%	Potassium (K)	<20 ppm
Calcium (Ca)	<10 ppm	Sodium (Na)	<50 ppm
Chloride (Cl)	<0.0005%	Total phosphorus (P)	<2 ppm
Copper (Cu)	<2 ppm	Total silicon (Si)	<2 ppm
Iron (Fe)	<5 ppm	Total sulphur (S)	<5 ppm
Lead (Pb)	<2 ppm	Zinc (Zn)	<2 ppm
Magnesium (Mg)	<5 ppm	Residue after ignition	<0.05%

Triethylamine

for HPLC

Certified

packaging	Mfr. No
250 mℓ Glass bottle	T/3203/07


C₆H₁₅N
CAS: 121-44-8
MW: 101.19
Density at 20°C: 0.728
Boiling point: 90°C
Melting point: -115°C
Flash point: -11°C
ADR/RID: 3, 8, II
EINECS: 204-469-4
UN: 1296

GHS: Flam. Liq. 2, Acute Tox. 3, Skin Corr. 1A, Eye Dam. 1, STOT SE 3
H225, H311, H302, H332, H314, H335
P280, P302+P350, P305+P351+P338, P304+P340, P308+P313, P261, P210

Assay (CG): >99.5%			
Guaranteed analysis			
Water	<0.05%	Iron (Fe)	<0.2 ppm
Residue after evaporation	<100 ppm	Lead (Pb)	<0.02 ppm

Absorbance specification	
Max. Absorbance (1cm cell v HPLC water)	
254 nm (sol. aqu. 0.1 M) ...	: <0.125 A.U.

Lot analysis on the label.
Filtered to 0.2 μm.



Triethylamine		Certified AR
for analysis		
packaging	Mfr. No	
1 ℓ Glass bottle	T/3201/15	
2.5 ℓ Metal drum	T/3201/25	
2.5 ℓ Coated Glass bottle	T/3201/PB17	
C ₆ H ₁₅ N CAS: 121-44-8 MW: 101.19 Density at 20°C: 0.728 Boiling point: 90°C Melting point: -115°C Flash point: -11°C ADR/RID: 3, 8, II EINECS: 204-469-4 UN: 1296	GHS: Flam. Liq. 2, Acute Tox. 3, Skin Corr. 1A, Eye Dam. 1, STOT SE 3 H225, H311, H302, H332, H314, H335 P280, P302+P350, P305+P351+P338, P304+P340, P308+P313, P261, P210	

Assay (CG): >99.5%			
Guaranteed analysis			
Colour	<10 APHA	Potassium (K)	<0.5 ppm
Water	<0.05%	Sodium (Na)	<1 ppm
Residue after evaporation	<50 ppm	Total phosphorus (P)	<0.2 ppm
Calcium (Ca)	<0.5 ppm	Total silicon (Si)	<0.1 ppm
Copper (Cu)	<0.05 ppm	Total sulphur (S)	<0.5 ppm
Iron (Fe)	<0.2 ppm	Zinc (Zn)	<0.2 ppm
Lead (Pb)	<0.05 ppm	Diethylamine	<0.1%
Magnesium (Mg)	<0.2 ppm		
Typical analysis			
Aluminium (Al)	<0.05 ppm	Manganese (Mn)	<0.05 ppm
Barium (Ba)	<0.05 ppm	Molybdenum (Mo)	<0.05 ppm
Cadmium (Cd)	<0.05 ppm	Nickel (Ni)	<0.05 ppm
Chromium (Cr)	<0.05 ppm	Strontium (Sr)	<0.05 ppm
Cobalt (Co)	<0.05 ppm	Titanium (Ti)	<0.05 ppm
Lithium (Li)	<0.05 ppm	Vanadium (V)	<0.05 ppm

Lot analysis on the label.

Triethylamine		SLR
specified reagents for general laboratory work- extra pure		
packaging	Mfr. No	
500 mℓ Glass bottle	T/3200/08	
500 mℓ Coated Glass bottle	T/3200/PB08	
2.5 ℓ Glass bottle	T/3200/17	
2.5 ℓ Coated Glass bottle	T/3200/PB17	
25 ℓ Metal drum	T/3200/25	
C ₆ H ₁₅ N CAS: 121-44-8 MW: 101.19 Density at 20°C: 0.728 Boiling point: 90°C Melting point: -115°C Flash point: -11°C ADR/RID: 3, 8, II EINECS: 204-469-4 UN: 1296	GHS: Flam. Liq. 2, Acute Tox. 3, Skin Corr. 1A, Eye Dam. 1, STOT SE 3 H225, H311, H302, H332, H314, H335 P280, P302+P350, P305+P351+P338, P304+P340, P308+P313, P261, P210	

Assay (CG): >99%			
Guaranteed analysis			
Colour	<20 APHA	Potassium (K)	<2 ppm
Water	<0.2%	Sodium (Na)	<5 ppm
Residue after evaporation	<100 ppm	Total phosphorus (P)	<2 ppm
Calcium (Ca)	<5 ppm	Total silicon (Si)	<2 ppm
Copper (Cu)	<2 ppm	Total sulphur (S)	<5 ppm
Iron (Fe)	<5 ppm	Zinc (Zn)	<2 ppm
Lead (Pb)	<2 ppm	Diethylamine	<0.5%
Magnesium (Mg)	<2 ppm		

UHPLC

Filtered to 0.1 μm, Fisher Chemical UHPLC gradient grade solvents (Acetonitrile, Methanol and Water) are perfect for those working at high pressure.

See page 309 for more information

Triethylene glycol		SLR
specified reagents for general laboratory work- extra pure		
packaging	Mfr. No	
500 mℓ Glass bottle	T/3270/08	
2.5 ℓ Glass bottle	T/3270/17	
C ₆ H ₁₄ O ₄ CAS: 112-27-6 MW: 150.18 Density at 20°C: 1.12	Boiling point: 285°C Melting point: -7°C Flash point: 165°C EINECS: 203-953-2	
Assay (CG): >95%		
Guaranteed analysis		
Acidity	<0.005 meq/g	

Trifluoroacetic acid		NEW	Optima
for LC/MS analysis			
packaging	Mfr. No		
0.5 mℓ Ampoule	A116-05AMP		
1 mℓ Ampoule	A116-1AMP		
2 mℓ Ampoule	A116-2AMP		
10x1 mℓ Ampoule	A116-10X1AMP		
50 mℓ HDPE bottle	A116-50		
C ₂ HF ₃ O ₂ CAS: 76-05-1 MW: 114.02 Density at 20°C: 1.535 Boiling point: 72°C Melting point: -15°C ADR/RID: 8, I EINECS: 200-929-3	UN: 2699 GHS: Acute Tox. 4, Skin Corr. 1A, Eye Dam. 1, Aquatic Chronic 3 H314, H332, H412 P301+P330+P331, P280, P305+P351+P338, P310, P273, P304+P340	 	
Assay: <99.5%			
Guaranteed analysis			
Colour	<10 APHA	Lead (Pb)	<10 ppb
Water	<5%	Magnesium (Mg)	<10 ppb
Reside after evaporation (at 0.1% aqueous)	: <1 ppm	Manganese (Mn)	<10 ppb
Fluoride (F)	<10 ppm	Nickel (Ni)	<10 ppb
Aluminium (Al)	<20 ppb	Potassium (K)	<20 ppb
Calcium (Ca)	<50 ppb	Propan--1-ol	<10 ppb
Chloride (Cl)	<5 ppm	Sodium (Na)	<50 ppb
Copper (Cu)	<10 ppb	Sulfate (SO ₄)	<2 ppm
Iron (Fe)	<10 ppb	Zinc (Zn)	<20 ppb

Absorbance specification			
Max. Absorbance (1cm cell v HPLC water)			
210 nm	<0.55 A.U.	275 nm	: <0.06 A.U.
220 nm	<0.18 A.U.	300 nm	<0.03 A.U.
230 nm	<0.06 A.U.	LC/MS Gradient Suitability (at 0.1% aqueous)	
254 nm	<0.003 A.U.		

This product has been tested on an LC-MS Thermo Scientific instrument to ensure suitability.

Trifluoroacetic acid		Certified
for HPLC		
packaging	Mfr. No	
25 mℓ Coated Glass bottle	T/3258/04	
100 mℓ Coated Glass bottle	T/3258/PB05	
C ₂ HF ₃ O ₂ CAS: 76-05-1 MW: 114.02 Density at 20°C: 1.535 Boiling point: 72°C Melting point: -15°C ADR/RID: 8, I EINECS: 200-929-3	UN: 2699 GHS: Acute Tox. 4, Skin Corr. 1A, Eye Dam. 1, Aqua- tic Chronic 3 H314, H332, H412 P301+P330+P331, P280, P305+P351+P338, P310, P273, P304+P340	
Assay: >99%		

Lot analysis on the label.
Filtered to 0.2 μm.

Trifluoroacetic acid solution 1%		Certified
HPLC for gradient analysis		
packaging	Mfr. No	
500 mℓ Coated Glass bottle	J/8800/PB08	
C ₂ HF ₃ O ₂ CAS: 76-05-1 MW: 114.02 EINECS: 200-929-3	GHS: Skin Irrit. 2, Eye Irrit. 2 H319, H315 P280, P261, P302+P352, P305+P351+P338	
Assay: 0.98 to 1.02%		
Conditions:		
Column: C18 Reverse Phase		
Flow Rate: 3ml/min		
Wavelength: 254nm		
Solvent A: Water + 10%Buffer		
Solvent B: Acetonitrile + 10% Buffer		
Gradient: 10 min. At 5% solvent B followed by a linear gradient to 95% solvent B over 15min.		

Trifluoroacetic acid			
for peptide sequencing			
packaging		Mfr. No	
100 ml	Coated Glass bottle	T/3255/PB05	
500 ml	Coated Glass bottle	T/3255/PB08	
C ₂ HF ₃ O ₂	UN: 2699		
CAS: 76-05-1	GHS: Acute Tox. 4, Skin		
MW: 114.02	Corr. 1A, Eye Dam. 1, Aqua- tic Chronic 3		
Density at 20°C: 1.535	H314, H332, H412		
Boiling point: 72°C	P301+P330+P331, P280,		
Melting point: -15°C	P305+P351+P338, P310,		
ADR/RID: 8, I	P273, P304+P340		
EINECS: 200-929-3			
Assay: >99.9%			
Guaranteed analysis			
Water	<0.05%	Chloride (Cl)	<0.001%
Residue after evaporation	<0.0002%	Iron (Fe)	<1 ppm
Trifluoroacetic anhydride test	: Pass	n-Butyl formate	<0.001%
Fluoride (F)	<10 ppm		
Absorbance specification			
Max. Absorbance (1cm cell v HPLC water)			
230 nm (sol. aqu. 0.1 M)	: <0.15 A.U.	254 nm (sol. aqu. 0.1%)	: <0.01 A.U.

Trifluoroacetic acid		SLR
specified reagents for general laboratory work- extra pure		
packaging	Mfr. No	
100 mℓ Coated Glass bottle	T/3256/PB05	
C ₂ HF ₃ O ₂ CAS: 76-05-1 MW: 114.02 Density at 20°C: 1.535 Boiling point: 72°C Melting point: -15°C ADR/RID: 8, I EINECS: 200-929-3	UN: 2699 GHS: Acute Tox. 4, Skin Corr. 1A, Eye Dam. 1, Aquatic Chronic 3 H314, H332, H412 P301+P330+P331, P280, P305+P351+P338, P310, P273, P304+P340	
Assay: >99.5%		

Trifluoroacetic anhydride		SLR
specified reagents for general laboratory work- extra pure		
packaging	Mfr. No	
100 mℓ Coated Glass bottle	T/3257/PB05	
C ₄ F ₆ O ₃ CAS: 407-25-0 MW: 210.03 Density at 20°C: 1.51 Boiling point: 39.1°C Melting point: -63.5°C ADR/RID: 8, I EINECS: 206-982-9	UN: 3265 GHS: Acute Tox. 4, Skin Corr. 1A, Aquatic Chronic 3 H332, H412, H314, EUH014 P301+P330+P331, P305+P351+P338, P280, P260, P273, P304+P340	
Assay: >99%		

Tri-n-hexylamine		Certified
for HPLC		
packaging	Mfr. No	
100 mℓ Coated Glass bottle	T/3296/PB05	
(CH ₃ (CH ₂) ₅) ₃ N CAS: 102-86-3 MW: 269.52 ADR/RID: 9, III EINECS: 203-062-9 UN: 3082	GHS: Acute Tox. 4, Skin Irrit. 2, Eye Irrit. 2, Aquatic Chronic 2 H319, H315, H302, H411 P280, P301+P312, P302+P352, P305+P351+P338, P273	

Assay (CG): > 97%
HPLC reagent suitable for HPLC analysis.

1,2,3-Trihydroxybenzene See Pyrogallol, page 267

1,3,5-Trihydroxybenzene dihydrate

See Phloroglucinol dihydrated, page 242

3,4,5-Trihydroxybenzoic acid
pure

packaging	Mfr. No	
100 g Plastic jar	T/3300/48	
C ₇ H ₆ O ₅ ·H ₂ O CAS: 149-91-7 MW: 188.14 Density at 20°C: 1.694 Melting point: 253°C EINECS: 205-749-9	GHS: Skin Irrit. 2, Eye Irrit. 2, STOT SE 3 H315, H335, H319 P261, P302+P352, P280, P305+P351+P338	

Guaranteed analysis
Complies with standard IR spectrum: Pass

Triiodomethane See Iodoform, page 161

3-Trimethoxysilyl-propyl-methacrylate

pure

packaging	Mfr. No	
100 g Plastic jar	T/3350/48	
C ₁₆ H ₂₆ O ₅ Si CAS: 2530-85-0 MW: 248.35 Density at 20°C: 1.04 Boiling point: 190°C Flash point: 92°C ADR/RID: 3/III	EINECS: 219-785-8 GHS: Skin Irrit. 2, Eye Irrit. 2, STOT SE 3 H319, H335, H315 P280, P261, P302+P352, P304+P340, P305+P351+P338, P312	

Guaranteed analysis
Complies with standard IR spectrum: Pass

Trimethylanilinium hydroxide solution 0.1 M
for gas chromatography

packaging	Mfr. No	
50 mℓ Glass bottle	T/3510/10	
C ₆ H ₁₃ NO CAS: 1899-02-1 MW: 153.23 Density at 20°C: 0.8 Flash point: 10°C ADR/RID: 3, 6.1, II EINECS: 217-592-3 UN: 1230	GHS: Flam. Liq. 2, Acute Tox. 3, Skin Irrit. 2, Eye Irrit. 2, STOT SE 1 H225, H315, H319, H370, H301, H311, H331 P307+P311, P304+P340, P302+P352, P280, P305+P351+P338, P210	

Guaranteed analysis
Complies with standard IR spectrum: Pass

Trimethylchlorosilane

See Chlorotrimethylsilane, page 87

Trimethylene glycol See Propane-1, page 264

2,2,4-Trimethylpentane

NEW PACK.

Certified

for HPLC

packaging		Mfr. No
1 ℓ	Glass bottle	T/3606/15
2.5 ℓ	Glass bottle	T/3606/17
5 ℓ	Aluminium can	T/3606/21
10 ℓ	Returnable Drum	T/3606/10RSS
25 ℓ	Returnable Drum	T/3606/25RSS
30 ℓ	Returnable Drum	T/3606/30RSS

C₈H₁₈
CAS: 540-84-1
MW: 114.23
Density at 20°C: 0.690
Boiling point: 98 to 99°C
Melting point: -107°C
Flash point: -12°C
ADR/RID: 3, II
EINECS: 208-759-1

UN: 1262
GHS: Flam. Liq. 2, Skin Irrit. 2, STOT SE 3, Asp. Tox. 1, Aquatic Acute 1, Aquatic Chronic 1
H225, H304, H410, H336, H315
P261, P301+P310, P331, P273, P302+P352, P210, P240

Assay (CG): >99.5%

Guaranteed analysis

Acidity / Alkalinity <0.0002 meq/g
Residue after evaporation <2 ppm
Water <0.01%

Absorbance specification

Max. Absorbance (1cm cell v HPLC water)

200 nm <0.4 A.U. 230 nm <0.02 A.U.
210 nm <0.2 A.U. 240 nm <0.005 A.U.
220 nm <0.1 A.U. 250 nm <0.005 A.U.

Lot analysis on the label.
Filtered to 0.2 µm.

(2) 40°C / min until 320°C
(3) 320°C for 30 min
Injection: 5µℓ cold on the column.

2,2,4-Trimethylpentane

Certified AR

for analysis

packaging		Mfr. No
1 ℓ	Glass bottle	T/3601/15
2.5 ℓ	Glass bottle	T/3601/17
25 ℓ	Metal drum	T/3601/25

C₈H₁₈
CAS: 540-84-1
MW: 114.23
Density at 20°C: 0.690
Boiling point: 98 to 99 to 99°C
Melting point: -107°C
Flash point: -12°C
ADR/RID: 3, II
EINECS: 208-759-1

UN: 1262
GHS: Flam. Liq. 2, Skin Irrit. 2, STOT SE 3, Asp. Tox. 1, Aquatic Acute 1, Aquatic Chronic 1
H225, H304, H410, H336, H315
P261, P301+P310, P331, P273, P302+P352, P210, P240

Assay (CG): >99.5%

Guaranteed analysis

Acidity / Alkalinity <0.0001 meq/g
Colour <10 APHA
Water <0.02%
Residue after evaporation ... <20<20 ppm
Calcium (Ca) <0.5 ppm
Copper (Cu) <0.1 ppm
Iron (Fe) <0.2 ppm
Lead (Pb) <0.1 ppm
Magnesium (Mg) <0.2 ppm
Potassium (K) <0.2 ppm
Sodium (Na) <0.5 ppm
Total phosphorus (P) <0.1 ppm
Total silicon (Si) <0.05 ppm
Total sulphur (S) <1 ppm
Zinc (Zn) <0.2 ppm

Typical analysis

Aluminium (Al) <0.05 ppm
Barium (Ba) <0.05 ppm
Cadmium (Cd) <0.05 ppm
Chromium (Cr) <0.05 ppm
Cobalt (Co) <0.05 ppm
Lithium (Li) <0.05 ppm
Manganese (Mn) <0.05 ppm
Molybdenum (Mo) <0.05 ppm
Nickel (Ni) <0.05 ppm
Strontium (Sr) <0.05 ppm
Titanium (Ti) <0.05 ppm
Vanadium (V) <0.05 ppm

Lot analysis on the label.

2,2,4-Trimethylpentane

for residue and pesticide analysis

Distol

packaging		Mfr. No
1 ℓ	Glass bottle	T/3608/15
2.5 ℓ	Glass bottle	T/3608/17

C₈H₁₈
CAS: 540-84-1
MW: 114.23
Density at 20°C: 0.690
Boiling point: 98 to 99°C
Melting point: -107°C
Flash point: -12°C
ADR/RID: 3, II
EINECS: 208-759-1

UN: 1262
GHS: Flam. Liq. 2, Skin Irrit. 2, STOT SE 3, Asp. Tox. 1, Aquatic Acute 1, Aquatic Chronic 1
H225, H304, H410, H336, H315
P261, P301+P310, P331, P273, P302+P352, P210, P240

Assay (CG): >95%

Guaranteed analysis

Water <0.02% GC (ECD) Analysis Pass
Residue after evaporation <2 ppm GC (FID) Analysis Pass
GC (NPD) analysis Pass

Analysis of pesticides: from the retention time of Lindane, after concentration of 2500 times, no impurity peak higher than that obtained with 10 pg / mL heptachloro-éoxide (ECD) and 5 pg / mL Parathion (NPD).

Analysis of hydrocarbons: no impurity chain longer than the dodecane is present at a concentration greater than 0.01 ppm (FID).

Conditions :
Column: J&W DB 5: 30m x 0.32mm x 0.1µm df.
Temperature profile:
(1) 50°C for 2 min

2,2,4-Trimethylpentane

NEW PACK.

SLR

specified reagents for general laboratory work- extra pure

packaging		Mfr. No
500 mℓ	Glass bottle	T/3600/08
2.5 ℓ	Glass bottle	T/3600/17
2.5 ℓ	Coated Glass bottle	T/3600/PB17
5 ℓ	Fluorinated HDPE Plastic bottle	T/3600/21
25 ℓ	Metal drum	T/3600/25
200 ℓ	Metal drum	T/3600/27

C₈H₁₈
CAS: 540-84-1
MW: 114.23
Density at 20°C: 0.690
Boiling point: 98 to 99°C
Melting point: -107°C
Flash point: -12°C
ADR/RID: 3, II
EINECS: 208-759-1

UN: 1262
GHS: Flam. Liq. 2, Skin Irrit. 2, STOT SE 3, Asp. Tox. 1, Aquatic Acute 1, Aquatic Chronic 1
H225, H304, H410, H336, H315
P261, P301+P310, P331, P273, P302+P352, P210, P240

Assay (CG): >99%

Guaranteed analysis

Acidity / Alkalinity <0.001 meq/g
Colour <10 APHA
Water <0.1%
Residue after evaporation <100 ppm
Calcium (Ca) <5 ppm
Copper (Cu) <2 ppm
Iron (Fe) <5 ppm
Lead (Pb) <2 ppm
Magnesium (Mg) <2 ppm
Potassium (K) <2 ppm
Sodium (Na) <5 ppm
Total silicon (Si) <2 ppm
Total sulphur (S) <10 ppm
Zinc (Zn) <2 ppm

NO-Bis (Trimethylsilyl) acetamide (BSA)

for gas chromatography

packaging		Mfr. No
10 mℓ	Glass bottle	B/3640/03

C₈H₂₁NOSi₂
CAS: 10416-59-8
MW: 203.43
Density at 20°C: 0.820
Boiling point: 71 to 73°C
Melting point: 24°C
Flash point: 42°C
ADR/RID: 4.1, 8, II
EINECS: 233-892-7

UN: 2925
GHS: Flam. Liq. 3, Acute Tox. 4, Skin Corr. 1B, Eye Dam. 1
H226, H314, H302, EUH014
P301+P330+P331, P280, P305+P351+P338, P310, P402+P404, P301+P312, P210, P240

Assay (CG): >90%

Trimethylsilyl chloride

See Chlorotrimethylsilane, page 87

N-Trimethylsilylimidazole

for gas chromatography

packaging		Mfr. No
25 g	Glass bottle	T/3622/46

C₆H₁₂N₂Si
CAS: 18156-74-6
MW: 140.26
Density at 20°C: 0.950
Boiling point: 93 to 94°C
Melting point: -42°C
Flash point: 5°C
ADR/RID: 3, II

EINECS: 242-040-3
UN: 1993
GHS: Flam. Liq. 2, Skin Irrit. 2, Eye Irrit. 2, STOT SE 3
H225, H319, H335, H315, EUH014
P210, P261, P264, P280, P305+P351+P338

Assay (CG): >98%

N,N-Bis(Trimethylsilyl) trifluoroacetamide

for gas chromatography

packaging		Mfr. No
10 g	Glass bottle	B/3647/45

C₈H₁₈F₃NOSi₂
CAS: 21149-38-2
MW: 257.40
ADR/RID: 3, III
EINECS: 244-242-7

UN: 1993
GHS: Flam. Liq. 3
H226
P210

Assay (CG): >98%

Triphenyl phosphate

pure

packaging		Mfr. No
500 g	Plastic jar	T/3650/53

C₁₈H₁₅PO₄
CAS: 115-86-6
MW: 326.29
Density at 20°C: 1.3
Boiling point: 370°C
Melting point: 47 to 53°C
Flash point: 220°C

ADR/RID: 6.1, III
EINECS: 204-112-2
UN: 2811
GHS: Aquatic Acute 1, Aquatic Chronic 1
H410
P273

Guaranteed analysis
Complies with standard IR spectrum: Pass

Triphenyl phosphine

pure

packaging		Mfr. No
100 g	Plastic jar	T/3652/48

C₁₈H₁₅P
CAS: 603-35-0
MW: 262.29
Boiling point: 377°C
Melting point: 78.5 to 81.5°C
Flash point: 181°C

ADR/RID: 6.1/II
EINECS: 210-036-0
GHS: Acute Tox. 4, Skin Sens. 1, STOT RE 2
H317, H302, H373
P280, P301+P312

Guaranteed analysis
Melting Point tested 78 to 82°C

2,3,5-Triphenyltetrazolium chloride

technical, (tetrazolium salt)

packaging		Mfr. No
10 g	Plastic jar	T/3700/45
25 g	Plastic jar	T/3700/46
100 g	Plastic jar	T/3700/48

C₁₉H₁₅CIN₄
CAS: 298-96-4
MW: 334.81
Melting point: >235°C
EINECS: 206-071-6

GHS: Skin Irrit. 2, Eye Irrit. 2, STOT SE 3
H319, H315, H335
P261, P280, P302+P352, P305+P351+P338

Guaranteed analysis
Complies with standard IR spectrum: Pass

Tris-Borate-EDTA buffer

for electrophoresis,(10x concentration)

packaging		Mfr. No
1 ℓ	Glass bottle	T/P050/15
10 ℓ	Bag in box	T/P050/24

CAS: 610769-35-2
Boiling point: 219 to 220°C
Melting point: 170°C
GHS: Repr. 1B, Skin Irrit. 2,

Eye Irrit. 2
H360FD, H315, H319
P280, P305+P351+P338, P308+P313, P201

Guaranteed analysis
Identity test : Pass

A pre-mixed concentrate prepared from high purity reagents.Each 1ℓ contains 121.1gr Tris buffer, 61.8g boric acid and 7.4g EDTA in filtered deionised water.

Tris-Glycine-SDS buffer

for electrophoresis,(10x concentration)

packaging		Mfr. No
1 ℓ	HDPE Plastic bottle	T/P010/15

Guaranteed analysis

Identity test : Pass

To use - dilute 1:1

High Purity Acids for Elemental Analysis

Fisher Chemical offers a complete range of acids to meet the most challenging applications. Choose from our Optima®, TraceMetal® and PrimarPlus® grades to ensure that the sensitivity you require is always available.

Read more page 24

360

Application For analysis Specified Titration

361

T | Tris(hydroxymethyl)methylamine

Tris(hydroxymethyl)methylamine
Tris Buffer, biological buffer

packaging	Mfr. No
100 g Plastic jar	T/P630/48
500 g Plastic jar	T/P630/53
1 kg Plastic jar	T/P630/60

C₄H₁₁NO₃
CAS: 77-86-1
MW: 121.14
Boiling point: 219 to 220°C
Melting point: 168.5°C
EINECS: 201-064-4

GHS: Skin Irrit. 2, Eye Irrit. 2, STOT SE 3
H335, H315, H319
P280, P305+P351+P338, P261, P302+P352

Assay: >99.8%

Guaranteed analysis

Loss on drying	<0.1%	Sodium (Na)	<50 ppm
Calcium (Ca)	<20 ppm	n-Butyl formate	<0.001%
Copper (Cu)	<1 ppm	Total phosphorus (P)	<20 ppm
Iron (Fe)	<1 ppm	Total silicon (Si)	<20 ppm
Lead (Pb)	<1 ppm	Total sulphur (S)	<50 ppm
Magnesium (Mg)	<10 ppm	Zinc (Zn)	<10 ppm
Potassium (K)	<50 ppm		

Absorbance specification

Max. Absorbance (1 cm cell v HPLC water)

280 nm (sol. aqu. 1 M) : <0.05 A.U.

Tris(hydroxymethyl)methylamine hydrochloride
Tris-HCl, for biochemistry

packaging	Mfr. No
100 g Plastic jar	T/P631/48
500 g Plastic jar	T/P631/53
1 kg Plastic jar	T/P631/60

C₄H₁₁NO₃.HCl
CAS: 1185-53-1
MW: 157.60
Melting point: 150 to 153°C
EINECS: 214-684-5

GHS: Skin Irrit. 2, Eye Irrit. 2, STOT SE 3
H335, H319, H315
P261, P302+P352, P280, P305+P351+P338

Assay: >99%

Guaranteed analysis

Arsenic (As)	<0.1 ppm	Lead (Pb)	<5 ppm
Copper (Cu)	<5 ppm	Zinc (Zn)	<5 ppm
Iron (Fe)	<5 ppm		

Absorbance specification

Max. Absorbance (1 cm cell v HPLC water)

260 nm (sol. aqu. 0.5 M) : <0.015 A.U. 280 nm (sol. aqu. 0.5 M) : <0.01 A.U.

N-Tris(hydroxymethyl)methyl-2-aminoethanesulfonic acid
for HPLC, biological buffer

packaging	Mfr. No
100 g Plastic jar	T/3720/48

C₆H₁₅NO₆S
CAS: 7365-44-8
MW: 229.25
Melting point: 223 to 225°C
EINECS: 230-906-3

GHS: Skin Irrit. 2, Eye Irrit. 2, STOT SE 3
H335, H315, H319
P261, P302+P352, P280, P305+P351+P338

Assay: >99%

Guaranteed analysis

Water	<0.2%	Sulfate (SO ₄)	<0.01%
pKa at 37°C	: 7.15 to 7.17	n-Butyl formate	<0.02%
Iron (Fe)	<5 ppm	Residue after ignition	<0.2%
Lead (Pb)	<10 ppm		

Absorbance specification

Max. Absorbance (1 cm cell v HPLC water)

250 nm (sol. aqu. 10%) : <0.1 A.U.

iso-Valeric acid
specified reagents for general laboratory work- extra pure

packaging	Mfr. No
100 mℓ Coated Glass bottle	V/0011/PB05

C₅H₁₀O₂
CAS: 503-74-2
MW: 102.13
Density at 20°C: 0.926
Boiling point: 176°C
Melting point: -35°C
Flash point: 70°C
ADR/RID: 8, 6.1, II
EINECS: 207-975-3

GHS: Acute Tox. 3, Skin Corr. 1B, Eye Dam. 1
H302, H311, H314
P280, P312, P302+P350, P301+P330+P331, P305+P351+P338, P310, P301+P312

Assay: >98%

Tris(hydroxymethyl) methylamine, 'Tris buffer'
specified reagents for general laboratory work- extra pure

packaging	Mfr. No
100 g Plastic jar	T/3710/48
250 g Plastic jar	T/3710/50
1 kg Plastic jar	T/3710/60

C₄H₁₁NO₃
CAS: 77-86-1
MW: 121.14
Boiling point: 219 to 220°C
Melting point: 168.5°C
EINECS: 201-064-4

GHS: Skin Irrit. 2, Eye Irrit. 2, STOT SE 3
H335, H315, H319
P280, P305+P351+P338, P261, P302+P352

Assay: >99%

Guaranteed analysis

Calcium (Ca)	<100 ppm	Sodium (Na)	<200 ppm
Copper (Cu)	<5 ppm	n-Butyl formate	<0.005%
Iron (Fe)	<5 ppm	Total phosphorus (P)	<100 ppm
Lead (Pb)	<5 ppm	Total silicon (Si)	<100 ppm
Magnesium (Mg)	<50 ppm	Total sulphur (S)	<200 ppm
Potassium (K)	<200 ppm	Zinc (Zn)	<50 ppm

Triton Kaowool® epsilon ceramic fibre
pure

packaging	Mfr. No
100 g Plastic jar	T/3740/48
500 g Plastic jar	T/3740/53

Melting point: .1650°C
GHS: Carc. 2, Acute Tox. 4, Skin Irrit. 2, Eye Irrit. 2, STOT SE 3
H351, H332, H315, H319, H335
P261, P280, P308+P313, P302+P352, P305+P351+P338, P304+P340

Guaranteed analysis

Identity test : Pass

Triton X-100® reduced
for biochemistry

packaging	Mfr. No
5 g Glass bottle	T/3753/44

C₃₄H₆₂O₁₁
CAS: 9002-93-1
MW: 646.85
Density at 20°C: 1.067
Boiling point: 270°C
Melting point: 6°C
Flash point: 274°C

ADR/RID: 9, III
UN: 3082
GHS: Acute Tox. 4, Eye Dam. 1, Aquatic Chronic 3
H412, H318, H302
P280, P305+P351+P338, P310, P273, P301+P312

Guaranteed analysis

Refractive index 1.465 to 1.478

Low fluorescence and UV absorption.



Triton X-100®
specified reagents for general laboratory work- extra pure

packaging	Mfr. No
500 mℓ HDPE Plastic bottle	T/3751/08
2.5 ℓ HDPE Plastic bottle	T/3751/17

C₃₄H₆₂O₁₁
CAS: 9002-93-1
MW: 646.85
Density at 20°C: 1.067
Boiling point: 270°C
Melting point: 6°C
Flash point: 274°C

ADR/RID: 9, III
UN: 3082
GHS: Acute Tox. 4, Eye Dam. 1, Aquatic Chronic 3
H412, H318, H302
P280, P305+P351+P338, P310, P273, P301+P312

Triton X-114®
pure

packaging	Mfr. No
500 mℓ HDPE Plastic bottle	T/3757/08

C₃₀H₅₄O₉
CAS: 9036-19-5
MW: 558.75

Density at 20°C: <1.05 UN: 3082

Trypsin
pure, from beef pancreas

packaging	Mfr. No
100 g Plastic jar	T/3760/48

CAS: 9002-07-7
EINECS: 232-650-8
GHS: Skin Irrit. 2, Eye Irrit. 2, Resp. Sens. 1, STOT SE 3

H319, H334, H315, H335
P304+P341, P342+P311, P261, P302+P352, P280, P305+P351+P338

Guaranteed analysis

Activity 0.4 to 0.6 U Anson

L-Tryptophan
specified reagents for general laboratory work- extra pure

packaging	Mfr. No
25 g Plastic jar	T/3770/46

C₁₁H₁₂N₂O₂
CAS: 73-22-3
MW: 204.23

Melting point: 280 to 285°C
EINECS: 200-795-6

Assay: >99%

Guaranteed analysis

Loss on drying	<0.5%	Iron (Fe)	<10 ppm
Specific rotation	-33° to -30°	Lead (Pb)	<10 ppm
Residue after ignition	<0.1%		

Tungsten standard metal solution 1000 ppm
for AAS, solution NH3 0.1M

packaging	Mfr. No
100 mℓ HDPE Plastic bottle	J/8120/05
500 mℓ HDPE Plastic bottle	J/8120/08

Manufactured to a tolerance of ± 0.5%.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.

pH Buffers
NIST Certified Standard Solution & Concentrate
A full range to discover on page 63

Tungsten standard metal solution 10 000 ppm | T

Tungsten standard metal solution 1000 ppm
for ICP, solution NH3 0.1M

packaging	Mfr. No
500 mℓ HDPE Plastic bottle	J/8336/08

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
All ICP standards are manufactured from acids, metal & metal salts Ultra pure quality and contain less than <10 ppm total metal impurities
.

Tungsten standard metal solution 1000 ppm
for ICP solution in HNO₃ 1%/HF 2%, NIST certified

packaging	Mfr. No
500 mℓ Plastic bottle	SPPLW22X

H311, H314

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Tungsten standard metal solution 1000 ppm
for ICP solution in H₂O/HF trace, NIST certified

packaging	Mfr. No
100 mℓ Plastic bottle	SPPLW92Y

Density at 20°C: 1

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.

Tungsten standard metal solution 10 000 ppm
for ICP solution in HNO₃ 2%/HF 5%, NIST certified

packaging	Mfr. No
500 mℓ Plastic bottle	SPPLW23X

H314

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Tungsten standard metal solution 10 000 ppm
for ICP solution in H₂O/HF trace, NIST certified

packaging	Mfr. No
100 mℓ Plastic bottle	SPPLW93Y
500 mℓ Plastic bottle	SPPLW93X

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Tungstophosphoric acid hydrate
for analysis

Certified AR

packaging	Mfr. No
100 g Plastic jar	T/3902/48
H ₃ PO ₄ ·12WO ₃ +aq CAS: 12501-23-4 MW: 2880.17 ADR/RID: 8, II EINECS: 235-087-6	UN: 3260 GHS: Skin Corr. 1B H314 P301+P330+P331, P305+P351+P338, P310, P280

Guaranteed analysis		
Aluminium (Al)	<5 ppm	Potassium (K) <500 ppm
Barium (Ba)	<5 ppm	Propan--1-ol <5 ppm
Cadmium (Cd)	<5 ppm	Sodium (Na) <500 ppm
Calcium (Ca)	<50 ppm	Strontium (Sr) <5 ppm
Copper (Cu)	<10 ppm	Titanium (Ti) <5 ppm
Iron (Fe)	<20 ppm	n-Butyl formate <0.001%
Lead (Pb)	<10 ppm	Total nitrogen (N) <20 ppm
Lithium (Li)	<5 ppm	Total sulphur (S) <100 ppm
Magnesium (Mg)	<50 ppm	Zinc (Zn) <50 ppm
Nickel (Ni)	<5 ppm	

Lot analysis on the label.

Tungstosilicic acid hydrate
specified reagents for general laboratory work- extra pure

SLR

packaging	Mfr. No
100 g Plastic jar	T/3940/48
SiO ₂ ·12WO ₃ ·XH ₂ O CAS: 12027-43-9 EINECS: 234-719-8	

Guaranteed analysis
Identity test : Pass

Turpentine oil
extra pure, rectified

SLR

packaging	Mfr. No
1 ℓ Glass bottle	T/4050/15
C ₁₀ H ₁₆ CAS: 8006-64-2 MW: 136112.00 Density at 20°C: 0.86 Boiling point: 154 to 170°C Melting point: -50--60°C Flash point: 35°C ADR/RID: 3, III EINECS: 232-350-7 UN: 1299	GHS: Flam. Liq. 3, Acute Tox. 4, Skin Irrit. 2, Eye Irrit. 2, Skin Sens. 1, Asp. Tox. 1, Aquatic Chronic 2 H226, H315, H411, H304, H302, H317, H319, H312, H332 P280, P302+P352, P301+P310, P331, P273, P301+P312, P304+P340, P305+P351+P338, P210

Turpentine substitute
specified reagents for general laboratory work- extra pure

SLR

packaging	Mfr. No
1 ℓ Coated Glass bottle	T/4200/PB15
2.5 ℓ Glass bottle	T/4200/17
2.5 ℓ HDPE Plastic bottle	T/4200/PB17
25 ℓ Metal drum	T/4200/25
CAS: 64742-88-7 Density at 20°C: 0.770 to 0.806 Boiling point: 149 to 204°C Flash point: 35°C ADR/RID: 3, III EINECS: 265-191-7 UN: 1300	GHS: Flam. Liq. 3, Acute Tox. 4, Asp. Tox. 1, Aquatic Chronic 2 H226, H304, H302, H411 P304+P340, P301+P310, P331, P210

Tween® 20
pure

packaging	Mfr. No
1 kg Glass bottle	T/4206/60
C ₅₈ H ₁₁₄ O ₂₆ CAS: 9005-64-5 Density at 20°C: 1.1 Boiling point: 100°C	Melting point: 0.1°C Flash point: >150°C EINECS: 500-018-3

Guaranteed analysis
Complies with standard IR spectrum: Pass

Tween® 40
pure

packaging	Mfr. No
1 kg Glass bottle	T/4205/60
CAS: 9005-66-7 Density at 20°C: 1.05 Melting point: 0.1°C Flash point: >200°C	

Guaranteed analysis
Complies with standard IR spectrum: Pass

Tween® 80
pure

packaging	Mfr. No
1 kg Glass bottle	T/4203/60
CAS: 9005-65-6 MW: approx. 1300 Density at 20°C: 1.073 to 1080	Flash point: >110°C EINECS: 500-019-9

L-Tyrosine
specified reagents for general laboratory work- extra pure

SLR

packaging	Mfr. No
100 g Plastic jar	T/4210/48
C ₉ H ₁₁ NO ₃ CAS: 60-18-4 MW: 181.19 Melting point: 290°C Flash point: 176°C EINECS: 200-460-4	GHS: Skin Irrit. 2, Eye Irrit. 2, STOT SE 3 H335, H319, H315 P261, P302+P352, P280, P305+P351+P338

Assay: >98.5%			
Guaranteed analysis			
Water	<0.5%	Residue after ignition	<0.1%
Loss on drying	<0.5%	Iron (Fe)	<50 ppm
Specific rotation	-11.2° to -9.8°	Lead (Pb)	<10 ppm

n-Undecane
pure

packaging	Mfr. No
100 mℓ Glass bottle	U/0005/05
C ₁₁ H ₂₄ CAS: 1120-21-4 MW: 156.31 Density at 20°C: 0.740 Boiling point: 196°C Melting point: -26°C Flash point: 60°C	ADR/RID: 3, III EINECS: 214-300-6 UN: 2330 GHS: Flam. Liq. 3, Asp. Tox. 1 H226, H304, EUH066 P310, P304+P340, P301+P310, P331, P210

Assay (CG): >99%

Uracil
specified reagents for general laboratory work- extra pure

SLR

packaging	Mfr. No
25 g Plastic jar	U/0030/46
C ₄ H ₄ N ₂ O ₂ CAS: 66-22-8 MW: 112.09	Melting point: 335°C EINECS: 200-621-9

Guaranteed analysis
Complies with standard IR spectrum: Pass

Uranium standard metal solution 1000 ppm
for ICP solution in HNO₃ 2%, NIST certified

packaging	Mfr. No
100 mℓ Plastic bottle	SPPLU22Y
500 mℓ Plastic bottle	SPPLU22X
H315, H319	

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurites.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Uranium standard metal solution 10 000 ppm
for ICP solution in HNO₃ 5%, NIST certified

packaging	Mfr. No
100 mℓ Plastic bottle	SPPLU23Y
500 mℓ Plastic bottle	SPPLU23X
H302, H314	

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurites.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Uranium standard metal solution 1000 ppm
for ICP-MS solution in HNO₃ 2%, NIST certified

packaging	Mfr. No
100 mℓ Plastic bottle	SPCLU22Y
H315, H319	

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurites.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

We Support Your Chemistry For

- Semi-bulk and bulk requirements
- Tailored solvents and solvent blends or customized solutions
- Testing services or customized packaging and labeling

Please email your contact details and a brief description of your requirements to beigel.sha.quotes@thermofisher.com

Urea
for biochemistry

packaging	Mfr. No
250 g Plastic jar	U/P610/50
1 kg Plastic jar	U/P610/60
CH ₄ N ₂ O CAS: 57-13-6 MW: 60.06	Melting point: 131 to 135°C EINECS: 200-315-5

Assay: >99.5%			
Guaranteed analysis			
Insolubles (5M aq. Soln)	: <0.005%	Lead (Pb)	<0.5 ppm
Cadmium (Cd)	<0.5 ppm	Magnesium (Mg)	<0.2 ppm
Cobalt (Co)	<0.5 ppm	Potassium (K)	<5 ppm
Copper (Cu)	<0.2 ppm	n-Butyl formate	<0.0005%
Iron (Fe)	<0.2 ppm		

Absorbance specification
Max. Absorbance (1cm cell v HPLC water) 260 nm (sol. aqu. 5 M) <0.06 A.U. 280 nm (sol. aqu. 5 M) <0.04 A.U.
For denaturation and solubilisation of proteins or nucleic acids.

Urea
for analysis

Certified AR

packaging	Mfr. No
500 g Plastic jar	U/0500/53
3 kg Plastic jar	U/0500/63
5 kg Plastic keg	U/0500/65
25 kg Fibre Keg	U/0500/70

CH ₄ N ₂ O CAS: 57-13-6 MW: 60.06 Melting point: 131 to 135°C	EINECS: 200-315-5 GHS: Skin Irrit. 2 H315 P302+P352
--	--

Assay: >99.5%			
Guaranteed analysis			
Acidity / Alkalinity	<0.0005 meq/g	Sodium (Na)	<20 ppm
Calcium (Ca)	<20 ppm	Total chloride (Cl)	<0.0005%
Copper (Cu)	<1 ppm	Total phosphorus (P)	<10 ppm
Iron (Fe)	<2 ppm	Total silicon (Si)	<5 ppm
Lead (Pb)	<1 ppm	Total sulphur (S)	<10 ppm
Magnesium (Mg)	<10 ppm	Zinc (Zn)	<5 ppm
Potassium (K)	<10 ppm	Biuret	<0.1%

Typical analysis			
Aluminium (Al)	<2 ppm	Molybdenum (Mo)	<2 ppm
Barium (Ba)	<2 ppm	Nickel (Ni)	<2 ppm
Cadmium (Cd)	<2 ppm	Silver (Ag)	<2 ppm
Chromium (Cr)	<2 ppm	Strontium (Sr)	<2 ppm
Cobalt (Co)	<2 ppm	Titanium (Ti)	<2 ppm
Lithium (Li)	<2 ppm	Vanadium (V)	<2 ppm
Manganese (Mn)	<2 ppm		

Lot analysis on the label.

Urea
specified reagents for general laboratory work- extra pure

SLR

packaging	Mfr. No
500 g Plastic jar	U/0450/53
1 kg Plastic jar	U/0450/60
3 kg Plastic jar	U/0450/63

CH ₄ N ₂ O CAS: 57-13-6 MW: 60.06 Melting point: 131 to 135°C	EINECS: 200-315-5 GHS: Skin Irrit. 2 H315 P302+P352
--	--

Assay: >99%			
Guaranteed analysis			
Calcium (Ca)	<100 ppm	Sodium (Na)	<100 ppm
Copper (Cu)	<5 ppm	n-Butyl formate	<0.002%
Iron (Fe)	<10 ppm	Total phosphorus (P)	<50 ppm
Lead (Pb)	<5 ppm	Total silicon (Si)	<20 ppm
Magnesium (Mg)	<50 ppm	Total sulphur (S)	<50 ppm
Potassium (K)	<50 ppm	Zinc (Zn)	<20 ppm

Urease-active meal
pure, from Jack beans

packaging	Mfr. No
25 g Plastic jar	U/0655/46
Guaranteed analysis	
Identity test	: Pass

Urease tablets

packaging	Mfr. No
60 PC Glass jar	U/0650/78
Guaranteed analysis	
Identity test	: Pass

Each tablet hydrolyses approx. 400mg in 3 hours at 37°C.

Vanadium standard metal solution 1000 ppm
for AAS, solution HNO₃ 1M

packaging	Mfr. No
100 mℓ HDPE Plastic bottle	J/8150/05
500 mℓ HDPE Plastic bottle	J/8150/08
Density at 20°C: 1.0	H314
ADR/RID: 8, II	P280, P301+P330+P331,
UN: 2031	P310, P305+P351+P338,
GHS: Skin Corr. 1B	P303+P361+P353

Manufactured to a tolerance of ± 0.5%.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.

Vanadium standard metal solution 1000 ppm
for ICP solution in HNO₃ 2%, NIST certified

packaging	Mfr. No
100 mℓ Plastic bottle	SPPLV22Y
500 mℓ Plastic bottle	SPPLV22X
H315, H319	

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Vanadium standard metal solution 1000 ppm
for ICP solution in HCL 2%, NIST certified

packaging	Mfr. No
500 mℓ Plastic bottle	SPPLV12X

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Reagents, Solvents & Volumetric solutions comply with the European Pharmacopeia

Read more on page 234

Vanadium standard metal solution 10 000 ppm
for ICP solution in HNO₃ 15%, NIST certified

packaging	Mfr. No
100 mℓ Plastic bottle	SPPLV43Y
500 mℓ Plastic bottle	SPPLV43X
H314	

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Vanadium standard metal solution 10 000 ppm
for ICP solution in HCL 15%, NIST certified

packaging	Mfr. No
500 mℓ Plastic bottle	SPPLV33X
H315, H319, H335	

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Vanadium standard metal solution 1000 ppm
for ICP-MS solution in HNO₃ 2%, NIST certified

packaging	Mfr. No
100 mℓ Plastic bottle	SPCLV22Y
H315, H319	

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Vanadium pentoxide
pure

packaging	Mfr. No
100 g Plastic jar	V/0100/48
500 g Plastic jar	V/0100/53
25 kg Fibre Keg	V/0100/70

V ₂ O ₅	UN: 2862
CAS: 1314-62-1	GHS: Muta. 2, Repr. 2,
MW: 181.88	Acute Tox. 4, STOT SE 3,
Density at 20°C: 3.35	STOT RE 1, Aquatic Chronic 2
Boiling point: 1750°C	H361d, H341, H302, H332,
Melting point: 690°C	H335, H372, H411
ADR/RID: 6.1, III	P281, P304+P340, P260, P273
EINECS: 215-239-8	

Guaranteed analysis

Identity test

Vanadyl sulfate hydrate
specified reagents for general laboratory work- extra pure

packaging	Mfr. No
100 g Plastic jar	V/0105/48
VOSO ₄ ·H ₂ O	GHS: Skin Irrit. 2, Eye Irrit. 2,
CAS: 12439-96-2	STOT SE 3
MW: 163.00	H335, H319, H315
ADR/RID: 6.1, II	P261, P302+P352, P280,
UN: 2931	P305+P351+P338

Guaranteed analysis

Identity test

Vinyl alcohol
pure, polymerized

packaging	Mfr. No
500 g Plastic jar	V/0200/53
CAS: 9002-89-5	Melting point: 200°C
MW: approx. 25 000	Flash point: >113°C
Density at 20°C: 1.26	

Guaranteed analysis

Degree of hydrolysis>84% Viscosity 4.5 to 6.5 cP

Water
for LC/MS analysis

packaging	Mfr. No
500 mℓ Glass bottle	W6-500
1 ℓ Glass bottle	W6-1
2.5 ℓ Glass bottle	W6-212
4 ℓ Glass bottle	W6-4

H ₂ O	Density at 20°C: 1
CAS: 7732-18-5	Boiling point: 100°C
MW: 18.02	EINECS: 231-791-2

Guaranteed analysis

Total Halogens (as chloride) : Not detected	Lead (Pb)<10 ppb
Aluminium (Al)<10 ppb	Magnesium (Mg)<10 ppb
Barium (Ba)<10 ppb	Manganese (Mn)<10 ppb
Cadmium (Cd)<10 ppb	Nickel (Ni)<10 ppb
Calcium (Ca)<20 ppb	Potassium (K)<10 ppb
Chromium (Cr)<10 ppb	Propan--1-ol<10 ppb
Cobalt (Co)<10 ppb	Sodium (Na)<20
Copper (Cu)<10 ppb	Tin (Sn)<10 ppm
Iron (Fe)<10 ppb	Zinc (Zn)<10 ppm

Absorbance specification

Max. Absorbance (1cm cell v HPLC water)	
210 nm<0.01 A.U.	50 ppb
220 nm<0.01 A.U.	LC/MS at Positive Mode as Reserpine : <
230 nm<0.01 A.U.	50 ppb
240 nm<0.01 A.U.	LC/UV Gradient (single max peak) at
254 nm<0.005 A.U.	210nm<0.005 A.U.
260 nm<0.005 A.U.	LC/UV Gradient (single max peak) at
280 nm<0.005 A.U.	254nm<0.0005 A.U.
LC/MS at Negative Mode as Aldicarb : <	

This product has been tested on an LC-MS Thermo Scientific instrument to ensure suitability.
Filtered to 0.2µm.

Safe, environmentally friendly and effective packaging

Fisher Chemical outer packaging is UN compliant and 100% recyclable. It's an easy and safe way to reduce your waste management costs and contribute to a cleaner environment.



Water
for HPLC-MS

packaging	Mfr. No
1 ℓ Glass bottle	W/0112/15
2.5 ℓ Glass bottle	W/0112/17

H ₂ O	Density at 20°C: 1
CAS: 7732-18-5	Boiling point: 100°C
MW: 18.02	EINECS: 231-791-2

Guaranteed analysis

Residue after evaporation<1 ppm	Lead (Pb)<0.01 ppm
LC-MS analysis Conform	Magnesium (Mg)<0.01 ppm
Aluminium (Al)<0.05 ppm	Manganese (Mn)<0.01 ppm
Barium (Ba)<0.01 ppm	Nickel (Ni)<0.01 ppm
Cadmium (Cd)<0.01 ppm	Potassium (K)<0.1 ppm
Calcium (Ca)<0.05 ppm	Propan--1-ol<0.01 ppm
Chromium (Cr)<0.01 ppm	Sodium (Na)<0.01 ppm
Cobalt (Co)<0.01 ppm	Tin (Sn)<0.01 ppm
Copper (Cu)<0.01 ppm	Zinc (Zn)<0.01 ppm
Iron (Fe)<0.01 ppm	

Absorbance specification

Max. Absorbance (1cm cell v HPLC water)	
200 nm<0.005 A.U.	Gradient (max eluted peak). <0.2 mA.U.
230 nm<0.005 A.U.	LC/MS Gradient Suitability: Pass test

Gradient analysis:

Absorbance of the largest eluted peak: <0.2 mA.U.
Wavelength 254 nm
This product has been tested on an LC-MS instrument to ensure suitability.
Lot analysis on the label.
Filtered to 0.2 µm.
All water by Fisher Chemical is deionised not distilled.

Water
for UHPLC gradient grade analysis

packaging	Mfr. No
1 ℓ Coated Glass bottle	W/0120/PB15
2.5 ℓ Coated Glass bottle	W/0120/PB17

H ₂ O	Density at 20°C: 1
CAS: 7732-18-5	Boiling point: 100°C
MW: 18.02	EINECS: 231-791-2

Guaranteed analysis

Residue after evaporation<1 ppm	Residue after ignition<0.0001%
Conductivity<1 µS	

Gradient analysis

Absorbance of the largest eluted peak at 220nm: <0.1mA.U.
Absorbance of the largest eluted peak at 254nm: <0.1mA.U.
Conditions
Column: C18 reverse phase
Wavelength: 220nm
Wavelength: 254nm
Solvent A: Water
Solvent B: Acetonitrile
Gradient: 10 minutes at 100% solvent A followed by a linear gradient 95% solvent B
Suitable for UHPLC instrument.
Lot analysis on the label.
Filtered to 0.1µm.

High Purity Acids
for Elemental Analysis

Fisher Chemical offers a complete range of acids to meet the most challenging applications. Choose from our Optima®, TraceMetal® and PrimarPlus® grades to ensure that the sensitivity you require is always available.
Read more page 24

Water		Certified Advanced HPLC
HPLC for gradient analysis		
packaging		Mfr. No
1 ℓ	Coated Glass bottle	W/0110/PB15
2.5 ℓ	Coated Glass bottle	W/0110/PB17
H ₂ O	Density at 20°C: 1	
CAS: 7732-18-5	Boiling point: 100°C	
MW: 18.02	EINECS: 231-791-2	
Guaranteed analysis		
Residue after evaporation	<5 ppm	Residue after ignition <0.002%
Absorbance specification		
Max. Absorbance (1cm cell v HPLC water)		
Gradient (max eluted peak). <0.1 mA.U.		

Gradient analysis (5% to 95% of the solvent)
Absorbance of the largest eluted peak: <0.1 mA.U.
Absorbance of the largest eluted peak <0.1 mA.U.
Conditions :
Column: C18 Reverse Phase.
Flow rate: 3 mL / min
Wavelength 220 nm
Solvent A: Water
Solvent B: Acetonitrile
10 min solvent A to 100% linear gradient until then 95% solvent for B 15 min.
All water by Fisher Chemical is deionised not distilled.
Lot analysis on the label.
Filtered to 0.2µm

Water NEW PACK.		Certified
HPLC for gradient analysis		
packaging		Mfr. No
1 ℓ	Glass bottle	W/0106/15
2.5 ℓ	Glass bottle	W/0106/17
2.5 ℓ	Coated Glass bottle	W/0106/PB17
10 ℓ	Returnable Drum	W/0106/10RSS
25 ℓ	Returnable Drum	W/0106/25RSS
30 ℓ	Returnable Drum	W/0106/30RSS
H ₂ O		Density at 20°C: 1
CAS: 7732-18-5		Boiling point: 100°C
MW: 18.02		EINECS: 231-791-2
Guaranteed analysis		
Residue after evaporation		<5 ppm Conductivity
Residue after ignition		<5 µS
		<0.0003%

Gradient analysis
Absorbance of the largest eluted peak:<0.002A.U.
Conditions :
Column: C18 Reverse Phase.
Flow rate: 2 mL / min
Wavelength 254 nm
Absorbance range: 0-0.04 AU.
Solvent A: Water
Solvent B: Acetonitrile
20 min solvent was followed by a linear gradient 100% A solvent 100%
All water by Fisher Chemical is deionised not distilled.
Lot analysis on the label.
Filtered to 0.2 µm.

We Support Your Chemistry For

- Semi-bulk and bulk requirements
- Tailored solvents and solvent blends or customized solutions
- Testing services or customized packaging and labeling

Please email your contact details and a brief description of your requirements to begel.sha.quotes@thermofisher.com

Water		Certified
HPLC for fluorescence detection		
packaging		Mfr. No
2.5 l	Glass bottle	W/0107/17
H ₂ O	Density at 20°C: 1	
CAS: 7732-18-5	Boiling point: 100°C	
MW: 18.02	EINECS: 231-791-2	
Guaranteed analysis		
Residue after evaporation		<5 ppm Residue after ignition <0.0003%
Absorbance specification		
Max. Absorbance (1 cm cell v HPLC water)		
Fluorescence		Pass

Excitation in the range 250-750nm produces no fluorescence emission greater than the solvents natural Raman signal at 350nm excittation.This is equivalent to the emission maximum of approximately 1ppb quinine sulfate.
All water by Fisher Chemical is deionised not distilled.
Lot analysis on the label (excluding fluorescence scan).
Filtered to 0.2 µm.

Water		Certified
HPLC for electrochemical detection		
packaging		Mfr. No
2.5 ℓ	Glass bottle	W/0108/17
H ₂ O	Density at 20°C: 1	
CAS: 7732-18-5	Boiling point: 100°C	
MW: 18.02	EINECS: 231-791-2	
Guaranteed analysis		
Residue after evaporation	<5 ppm	Electrochemical test
Residue after ignition	<0.0003%	Pass

Tested by HPLC-EC to ensure low electrochemical response.
Detector - ESA Coulochem II
Mode-oxidative
Graphite Electrode.
All water by Fisher Chemical is deionised not distilled.
Lot analysis on the label.
Filtered to 0.2 µm.

Water		Optima
for ultra trace elemental analysis		
packaging	Mfr. No	
500 mℓ LDPE bottle	W9-500	
1 ℓ LDPE bottle	W9-1	
2 ℓ LDPE bottle	W9-2	
H ₂ O	Density at 20°C: 1	
CAS: 7732-18-5	Boiling point: 100°C	
MW: 18.02	EINECS: 231-791-2	

Guaranteed analysis		
Colour	<10 APHA	Gadolinium (Gd) <1 ppt
Sulfate (SO ₄)	<1 ppb	Germanium (Ge) <10 ppt
Total chloride (Cl)	<1 ppb	Gold (Au) <10 ppt
Total phosphorus (P)	<1 ppb	Hafnium (Hf) <1 ppt
Zirconium (Zr)	<10 ppt	Holmium (Ho) <1 ppt
Gallium (Ga)	<10 ppt	Iron (Fe) <20 ppt
Terbium (Tb)	<10 ppt	Lanthanum (La) <1 ppt
Indium (In)	<1 ppt	Lead (Pb) <10 ppt
Dysprosium (Dy)	<1 ppt	Lithium (Li) <10 ppt
Aluminium (Al)	<20 ppt	Lutetium (Lu) <1 ppt
Antimony (Sb)	<10 ppt	Magnesium (Mg) <10 ppt
Arsenic (As)	<10 ppt	Manganese (Mn) <10 ppt
Barium (Ba)	<10 ppt	Mercury (Hg) <20 ppt
Beryllium (Be)	<10 ppt	Molybdenum (Mo) <10 ppt
Bismuth (Bi)	<10 ppt	Neodymium (Nd) <1 ppt
Boron (B)	<50 ppt	Nickel (Ni) <10 ppt
Cadmium (Cd)	<10 ppt	Niobium (Nb) <10 ppt
Calcium (Ca)	<20 ppt	Palladium (Pd) <10 ppt
Cerium (Ce)	<10 ppt	Potassium (K) <10 ppt
Cesium (Cs)	<10 ppt	Praseodymium (Pr) <10 ppt
Chromium (Cr)	<10 ppt	Rhenium (Re) <10 ppt
Cobalt (Co)	<10 ppt	Rhodium (Rh) <10 ppt
Copper (Cu)	<10 ppt	Ruthenium (Ru) <10 ppt
Erbium (Er)	<1 ppt	Samarium (Sm) <10 ppt
Europium (Eu)	<1 ppt	

Scandium (Sc)	<10 ppt	Tin (Sn)	<10 ppt
Selenium (Se)	<50 ppt	Titanium (Ti)	<10 ppt
Propan--1-ol	<10 ppt	Uranium (U)	<1 ppt
Sodium (Na)	<10 ppt	Vanadium (V)	<10 ppt
Strontium (Sr)	<10 ppt	Ytterbium (Yb)	<10 ppt
Tantalum (Ta)	<10 ppt	Yttrium (Y)	<1 ppt
Tellurium (Te)	<1 ppt	Zinc (Zn)	<10 ppt
Thallium (Tl)	<10 ppt	Platinum (Pt)	<10 ppt
Thorium (Th)	<1 ppt	Rubidium (Rb)	<10 ppt
Thulium (Tm)	<10 ppt		

Bottled in class 10 cleanroom conditions.Individually double-bagged in class 100 cleanroom conditions. Each bottle is individually packaged.
Certificate of analysis delivered with each bottle.

Water blank in HCl 5%	
for ICP, NIST certified	
packaging	Mfr. No
500 ml Plastic bottle	SPPLBLKHCL
Manufactured to a tolerance of ± 0.5%.	
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.	
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.	
Check with your local distributor for availability of Spex product	

Water blank in HNO ₃ 5%		Certified
for ICP, NIST certified		
packaging	Mfr. No	
500 mℓ Plastic bottle	SPPLBLKHNO3	
H314		

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Water blank in HCl 2%		Certified
for ICP-MS, NIST certified		
packaging	Mfr. No	
100 mℓ Plastic bottle	SPCLBLKHCL	
H314, H318		

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Water blank in HNO ₃ 2%		Certified
for ICP-MS, NIST certified		
packaging	Mfr. No	
100 mℓ Plastic bottle	SPCLBLKHNO3	
H315, H319		

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Water Blank	
for ICP-MS, NIST certified	
packaging	Mfr. No
100 mL Plastic bottle	SPCLBLKH20
Manufactured to a tolerance of ± 0.5%.	
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.	
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.	
Check with your local distributor for availability of Spex product	

Water		Certified AR
for analysis		
packaging	Mfr. No	
5 ℓ HDPE Plastic bottle	W/0100/21	
5 ℓ Fluorinated HDPE square Drum	W/0100/FP21	
10 ℓ Bag in box	W/0100/24	
25 ℓ HDPE Plastic container	W/0100/25	
200 ℓ HDPE Plastic container	W/0100/27	
H ₂ O	Density at 20°C: 1	
CAS: 7732-18-5	Boiling point: 100°C	
MW: 18.02	EINECS: 231-791-2	

Guaranteed analysis		
pH at 20°C	: 4.5 to 7.5	Substances
Conductivity	><5 µS	reducing KMnO ₄ <0.00001%
Calcium (Ca)	<0.05 ppm	Total chloride (Cl) <0.00001%
Copper (Cu)	<0.01 ppm	Total nitrogen (N) <0.1 ppm
Iron (Fe)	<0.01 ppm	Total phosphorus (P) <0.01 ppm
Lead (Pb)	<0.005 ppm	Total silicon (Si) <0.2 ppm
Magnesium (Mg)	<0.05 ppm	Total sulphur (S) <0.5 ppm
Potassium (K)	<0.05 ppm	Zinc (Zn) <0.01 ppm
Sodium (Na)	<0.2 ppm	

Typical analysis		
Aluminium (Al)	<0.01 ppm	Molybdenum (Mo) <0.01 ppm
Barium (Ba)	<0.01 ppm	Nickel (Ni) <0.01 ppm
Cadmium (Cd)	<0.01 ppm	Silver (Ag) <0.01 ppm
Chromium (Cr)	<0.01 ppm	Strontium (Sr) <0.01 ppm
Cobalt (Co)	<0.01 ppm	Tin (Sn) <0.01 ppm
Lithium (Li)	<0.01 ppm	Titanium (Ti) <0.01 ppm
Manganese (Mn)	<0.01 ppm	Vanadium (V) <0.01 ppm

All water by Fisher Chemical is deionised not distilled.
Lot analysis on the label.

Water with 0.1% formic acid (v/v)		Optima
for LC/MS analysis		
packaging	Mfr. No	
500 mℓ Glass bottle	LS118-500	
1 ℓ Glass bottle	LS118-1	
2.5 ℓ Glass bottle	LS118-212	
4 ℓ Glass bottle	LS118-4	

Density at 20°C: 1		
Boiling point: 100°C		
Guaranteed analysis		
Assay (v/v)	: 0.095 to0.105%	Magnesium (Mg) <10 ppb
Colour	<10 APHA	Manganese (Mn) <10 ppb
Residue after evaporation	<1 ppm	Nickel (Ni) <10 ppb
Aluminium (Al)	<20 ppb	Potassium (K) <20 ppb
Calcium (Ca)	<50 ppb	Propan--1-ol <10 ppb
Copper (Cu)	<10 ppb	Sodium (Na) <50 ppb
Iron (Fe)	<10 ppb	Zinc (Zn) <20 ppb
Lead (Pb)	<10 ppb	

Absorbance specification		
Max. Absorbance (1cm cell v HPLC water)		
210 nm	<1.25 A.U.	254 nm <0.01 A.U.
220 nm	<0.85 A.U.	LC/MS Gradient Suitability : Pass test
230 nm	<0.55 A.U.	

This product has been tested on an LC-MS Thermo Scientific instrument to ensure suitability.
Filtered to 0.2µm.

Water with 0.1% formic acid (v/v) NEW Certified
HPLC for gradient analysis

packaging		Mfr. No	
2.5 ℓ	Glass bottle	FE/1200/17	
Density at 20°C: 1.00			
Boiling point: 100°C			
Guaranteed analysis			
Residue after evaporation	<=1 ppm	Gradient analysis at 254 nm . : <0.01 A.U.	
Formic acid (% v/v)	: 0.095 to 0.105%		
Absorbance specification			
Max. Absorbance (1cm cell v HPLC water)			
210 nm	<1.25 A.U.	250 nm	<0.015 A.U.
220 nm	<0.85 A.U.	254 nm	<0.005 A.U.
230 nm	<0.35 A.U.	260 nm	<0.005 A.U.
240 nm	<0.085 A.U.		
Lot analysis on the label.			
Filtered to 0.2µm.			

Water with 0.1% trifluoroacetic acid (v/v) NEW Optima
for LC/MS analysis

packaging		Mfr. No	
500 mℓ	Glass bottle	LS119-500	
1 ℓ	Glass bottle	LS119-1	
2.5 ℓ	Glass bottle	LS119-212	
4 ℓ	Glass bottle	LS119-4	
CAS: 76-05-1			
Density at 20°C: 1			
EINECS: 200-929-3			
Guaranteed analysis			
Assay (v/v)	: 0.095 to0.105%	Magnesium (Mg)	<10 ppb
Colour	<10 APHA	Manganese (Mn)	<10 ppb
Residue after evaporation ...	<1 ppm ppm	Nickel (Ni)	<10 ppb
Aluminium (Al)	<20 ppb	Potassium (K)	<20 ppb
Calcium (Ca)	<50 ppb	Propan--1-ol	<10 ppb
Copper (Cu)	<10 ppb	Sodium (Na)	<50 ppb
Iron (Fe)	<10 ppb	Zinc (Zn)	<20 ppb
Lead (Pb)	<10 ppb		
Absorbance specification			
Max. Absorbance (1cm cell v HPLC water)			
210 nm	<0.55 A.U.	254 nm	<0.003 A.U.
220 nm	<0.18 A.U.	LC/MS Gradient Suitability	: Pass test
230 nm	<0.06 A.U.		

Water with 0.1% trifluoroacetic acid (v/v) NEW Certified
HPLC for gradient analysis

packaging		Mfr. No
2.5 ℓ	Glass bottle	FE/0100/17
CAS: 76-05-1 Density at 20°C: 1 EINECS: 200-929-3		
Guaranteed analysis		
Residue after evaporation	<=1 ppm	Trifluoroacetic acid (% v/v) : 0.095 to 0.105%
Gradient analysis at 254 nm. :	<0.01 A.U.	
Absorbance specification		
Max. Absorbance (1cm cell v HPLC water)		
205 nm	<1 A.U.	240 nm <0.03 A.U.
210 nm	<0.7 A.U.	250 nm <0.005 A.U.
220 nm	<0.2 A.U.	254 nm <0.005 A.U.
230 nm	<0.07 A.U.	260 nm <0.005 A.U.
Lot analysis on the label. Filtered to 0.2µm.		

Wax Beeswax White
pure

Packaging		Mfr. No
500 g	Plastic jar	W/0150/53
CAS: 8012-89-3		
Melting point: 61 to 65°C		
EINECS: 232-383-7		
Guaranteed analysis		
Identity test		: Pass

Wax Beeswax
pure

Packaging		Mfr. No
50 g	Plastic jar	W/0200/50
CAS: 8012-89-3		
Melting point: 61 to 65°C		
INECS: 232-383-7		
Guaranteed analysis		
Complies with standard IR spectrum: Pass		

Wax Carnauba yellow
pure

Packaging		Mfr. No
100 g	Plastic jar	W/0255/48
CAS: 8015-86-9		
Melting point: 81 to 86°C		
EINECS: 232-399-4		
Guaranteed analysis		
Identity test		: Pass

Wax Ceresin white
technical

Packaging		Mfr. No
50 g	Plastic jar	W/0300/50
CAS: 8001-75-0		
Melting point: 61 to 78°C		
INECS: 232-290-1		
Guaranteed analysis		
Complies with standard IR spectrum: Pass		

Wax lanette emulsifying
pure

Packaging		Mfr. No
100 g	Plastic jar	W/0320/53
Melting point: 52°C		
Guaranteed analysis		
Identity test : Pass		

Wood's metal
pure, m.p. ca. 70°C

Packaging		Mfr. No
100 g	Plastic container	W/0350/48
500 g	Plastic jar	W/0350/53
CAS: 76093-98-6		Aquatic Chronic 1
Melting point: 70°C		H410, H302, H360Df, H350,
ADR/RID: 6.1, II		H330, H372, H341
UN: 2570		P280, P260, P273,
GHS: Carc. 1B, Muta. 2,		P301+P312, P304+P340,
Repr. 1A, Acute Tox. 2,		P201, P307+P311
STOT RE 1, Aquatic Acute 1,		
		

Guaranteed analysis

Identity test : Pass

o-Xylene
pure

packaging		Mfr. No	
500 mℓ	Glass bottle	X/0300/08	
2.5 ℓ	Glass bottle	X/0300/17	
200 ℓ	Metal drum	X/0300/27	
C ₈ H ₁₀ CAS: 95-47-6 MW: 106.17 Density at 20°C: 0.878 Boiling point: 143 to 145°C Melting point: -25°C Flash point: 31°C ADR/RID: 3, III		EINECS: 202-422-2 UN: 1307 GHS: Flam. Liq. 3, Acute Tox. 4, Skin Irrit. 2, Asp. Tox. 1 H226, H304, H312, H315, H332 P280, P301+P310, P331, P305+P351+P338, P210	
			

Assay (CG): >98%

m-Xylene
pure

packaging		Mfr. No	
500 mℓ	Glass bottle	X/0350/08	
C ₈ H ₁₀ CAS: 108-38-3 MW: 106.17 Density at 20°C: 0.864 Boiling point: 139 to 139°C Melting point: -48°C Flash point: 25°C		ADR/RID: 3, III EINECS: 203-576-3 UN: 1307 GHS: Flam. Liq. 3, Acute Tox. 4, Skin Irrit. 2 H226, H315, H312, H332 P280, P305+P351+P338, P210	
			

Assay (CG): >98%

p-Xylene
pure

packaging		Mfr. No	
500 mℓ	Glass bottle	X/0400/08	
C ₈ H ₁₀ CAS: 106-42-3 MW: 106.17 Density at 20°C: 0.866 Boiling point: 138°C Melting point: 13°C Flash point: 25°C		ADR/RID: 3, III EINECS: 203-396-5 UN: 1307 GHS: Flam. Liq. 3, Acute Tox. 4, Skin Irrit. 2 H226, H315, H312, H332 P280, P210, P305+P351+P338	
			

Assay (CG): >98%

Xylene
for electronic use (MOS)

packaging		Mfr. No
2.5 ℓ	Glass bottle	X/0255/17
5 ℓ	Aluminium can	X/0255/21
200 ℓ	Metal drum	X/0255/27
C ₈ H ₁₀ CAS: 1330-20-7 MW: 106.17 Density at 20°C: 0.865 Boiling point: 136 to 140°C Melting point: -34°C Flash point: 21°C ADR/RID: 3, II		EINECS: 215-535-7 UN: 1307 GHS: Flam. Liq. 2, Acute Tox. 4, Skin Irrit. 2 H225, H332, H315, H312 P280, P302+P352, P304+P340, P210, P240
		
Assay (CG): >98%		

Guaranteed analysis			
Acidity / Alkalinity		<0.0001 meq/g	Molybdenum (Mo)
Colour		<10 APHA	Nickel (Ni)
Water		<0.05%	Potassium (K)
Aluminium (Al)		<0.02 ppm	Sodium (Na)
Barium (Ba)		<0.02 ppm	Strontium (Sr)
Cadmium (Cd)		<0.02 ppm	Tin (Sn)
Calcium (Ca)		<0.5 ppm	Titanium (Ti)
Chromium (Cr)		<0.02 ppm	Total phosphorus (P)
Cobalt (Co)		<0.02 ppm	Total silicon (Si)
Copper (Cu)		<0.05 ppm	Total sulphur (S)
Iron (Fe)		<0.1 ppm	Vanadium (V)
Lead (Pb)		<0.05 ppm	Zinc (Zn)
Lithium (Li)		<0.02 ppm	Benzene
Magnesium (Mg)		<0.2 ppm	n-Propyl benzene
Manganese (Mn)		<0.02 ppm	Toluene
			<1%
			<1%

Filtered to 0.2 µm.

Xylene NEW PACK. Certified AR
for analysis

packaging		Mfr. No	
1 ℓ	Glass bottle	X/0250/15	
2.5 ℓ	Glass bottle	X/0250/17	
2.5 ℓ	Coated Glass bottle	X/0250/PB17	
5 ℓ	Fluorinated HDPE Plastic bottle	X/0250/21	
25 ℓ	Metal drum	X/0250/25	
200 ℓ	Metal drum	X/0250/27	
C ₈ H ₁₀ CAS: 1330-20-7 MW: 106.17 Density at 20°C: 0.865 Boiling point: 136 to 136 to 140 °C Melting point: -34°C Flash point: 21°C		ADR/RID: 3, II EINECS: 215-535-7 UN: 1307 GHS: Flam. Liq. 2, Acute Tox. 4, Skin Irrit. 2 H225, H332, H315, H312 P280, P302+P352, P304+P340, P210, P240	
			

Assay (CG): >98%

Guaranteed analysis			
Acidity / Alkalinity		<0.0001 meq/g	Nickel (Ni)
Colour		<10 APHA	Potassium (K)
Water		<0.05%	Sodium (Na)
Residue after evaporation		<20 ppm	Total phosphorus (P)
Substances			Total silicon (Si)
reducing KMnO ₄		<0.0005%	Total sulphur (S)
Calcium (Ca)		<0.5 ppm	Zinc (Zn)
Copper (Cu)		<0.05 ppm	Benzene
Iron (Fe)		<0.1 ppm	n-Propyl benzene
Lead (Pb)		<0.05 ppm	Toluene
Magnesium (Mg)		<0.2 ppm	<1%
			<1%

Typical analysis			
Aluminium (Al)		<0.02 ppm	Molybdenum (Mo)
Barium (Ba)		<0.02 ppm	Selenium (Se)
Cadmium (Cd)		<0.02 ppm	Strontium (Sr)
Chromium (Cr)		<0.02 ppm	Tin (Sn)
Cobalt (Co)		<0.02 ppm	Titanium (Ti)
Lithium (Li)		<0.02 ppm	Vanadium (V)
Manganese (Mn)		<0.02 ppm	<0.02 ppm

Lot analysis on the label.



Xylene **NEW PACK** SLR
specified reagents for general laboratory work- extra pure

packaging		Mfr. No
500 mℓ	Glass bottle	X/0200/08
2.5 ℓ	Glass bottle	X/0200/17
2.5 ℓ	Coated Glass bottle	X/0200/PB17
5 ℓ	Fluorinated HDPE Plastic bottle	X/0200/21
25 ℓ	Metal drum	X/0200/25
200 ℓ	Metal drum	X/0200/27
C ₈ H ₁₀ CAS: 1330-20-7 MW: 106.17 Density at 20°C: 0.865 Boiling point: 136 to 140°C Melting point: -34°C Flash point: 21°C ADR/RID: 3, II		EINECS: 215-535-7 UN: 1307 GHS: Flam. Liq. 2, Acute Tox. 4, Skin Irrit. 2 H225, H332, H315, H312 P280, P302+P352, P304+P340, P210, P240
Assay (CG): >97%		
Guaranteed analysis		
Acidity / Alkalinity	<0.003 meq/g	Lead (Pb) <2 ppm
Colour	<10 APHA	Magnesium (Mg) <2 ppm
Water	<0.05%	Potassium (K) <2 ppm
Residue after evaporation	<100 ppm	Sodium (Na) <5 ppm
Substances		Total phosphorus (P) <2 ppm
reducing KMnO ₄	<0.001%	Total silicon (Si) <2 ppm
Calcium (Ca)	<5 ppm	Total sulphur (S) <10 ppm
Copper (Cu)	<2 ppm	Zinc (Zn) <2 ppm
Iron (Fe)	<5 ppm	

Xylene
Technical

packaging		Mfr. No
1 ℓ	Metal can	X/0100/MC15
2.5 ℓ	Glass bottle	X/0100/17
2.5 ℓ	Coated Glass bottle	X/0100/PB17
25 ℓ	Metal drum	X/0100/25
200 ℓ	Metal drum	X/0100/27
C ₈ H ₁₀ CAS: 1330-20-7 MW: 106.17 Density at 20°C: 0.865 Boiling point: 136 to 140°C Melting point: -34°C Flash point: 21°C ADR/RID: 3, II		EINECS: 215-535-7 UN: 1307 GHS: Flam. Liq. 2, Acute Tox. 4, Skin Irrit. 2 H225, H332, H315, H312 P280, P302+P352, P304+P340, P210, P240
Assay: > 97%		
Assay (CG): >97%		

Xylenol Orange
See o-Cresolsulfonphthaleindi-(methyl-iminodiacetic acid) sodium salt, page 99

Yeast extract powder
pure

packaging		Mfr. No
100 g	Plastic jar	Y/0200/48
Guaranteed analysis		
Identity test	: Pass	

Optima LC/MS

Fisher Chemical OPTIMA® LC/MS solvents, blends, additives and modifiers have set the standard of excellence for consistent, reproducible performance in the mobile phase of LC/MS.

Find out more on page 308

Ytterbium standard metal solution 1000 ppm
for ICP solution in HNO₃ 2%, NIST certified

packaging		Mfr. No
100 mℓ	Plastic bottle	SPPLY22Y
250 mℓ	Plastic bottle	SPPLY22T
500 mℓ	Plastic bottle	SPPLY22X
H315, H319		

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Ytterbium standard metal solution 1000 ppm
for ICP solution in HNO₃ 2%, NIST certified

packaging		Mfr. No
100 mℓ	Plastic bottle	SPPLYB22Y
500 mℓ	Plastic bottle	SPPLYB22X
H315, H319		

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Ytterbium standard metal solution 10 ppm
for ICP-MS solution in HNO₃ 2%, NIST certified

packaging		Mfr. No
100 mℓ	Plastic bottle	SPCLY21AY
H315, H319		

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Yttrium standard metal solution 10 000 ppm
for ICP solution in HNO₃ 5%, NIST certified

packaging		Mfr. No
100 mℓ	Plastic bottle	SPPLY23Y
500 mℓ	Plastic bottle	SPPLY23X
H314		

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Save money!

Order 4 bottles in a box

Zinc standard metal solution 1000 ppm
for AAS, solution HNO₃ 1M

packaging		Mfr. No
100 mℓ	HDPE Plastic bottle	J/8070/05
500 mℓ	HDPE Plastic bottle	J/8070/08
ZnO MW: 81.38 ADR/RID: 8, II EINECS: 231-943-8 UN: 2031		GHS: Skin Corr. 1B, Eye Dam. 1 H314 P280, P305+P351+P338, P301+P330+P331, P310

Manufactured to a tolerance of ± 0.5%.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.

Zinc standard metal solution 1000 ppm
for ICP, solution HNO₃ 1M

packaging		Mfr. No
500 mℓ	HDPE Plastic bottle	J/8341/08
ZnO CAS: 7779-88-6 MW: 81.38 ADR/RID: 8, II EINECS: 231-943-8		UN: 2031 GHS: Skin Corr. 1B, Eye Dam. 1 H314 P280, P305+P351+P338, P301+P330+P331, P310

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.

Zinc standard metal solution 10 000 ppm
for ICP, solution HNO₃ 1M

packaging		Mfr. No
100 mℓ	HDPE Plastic bottle	J/8340/05
ZnO MW: 81.38 ADR/RID: 8, III UN: 3264 GHS: Skin Corr. 1B, Eye Dam. 1		H314 P280, P301+P330+P331, P310, P305+P351+P338, P303+P361+P353

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.

Zinc standard metal solution 1000 ppm
for ICP solution in HNO₃ 2%, NIST certified

packaging		Mfr. No
100 mℓ	Plastic bottle	SPPLZN22Y
250 mℓ	Plastic bottle	SPPLZN22T
H315, H319		

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Zinc standard metal solution 1000 ppm
for ICP solution in HCL 2%, NIST certified

packaging		Mfr. No
500 mℓ	Plastic bottle	SPPLZN12X

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Zinc standard metal solution 10 000 ppm
for ICP solution in HNO₃ 5%, NIST certified

packaging		Mfr. No
100 mℓ	Plastic bottle	SPPLZN23Y
500 mℓ	Plastic bottle	SPPLZN23X
H314		

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Zinc standard metal solution 10 000 ppm
for ICP solution in HCL 5%, NIST certified

packaging		Mfr. No
500 mℓ	Plastic bottle	SPPLZN13X

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Zinc standard metal solution 1000 ppm
for ICP-MS solution in HNO₃ 2%, NIST certified

packaging		Mfr. No
100 mℓ	Plastic bottle	SPCLZN22Y
H315, H320		

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Zinc granulated (metal)
for analysis

packaging		Mfr. No
500 g	Plastic jar	Z/0300/53
1 kg	Plastic jar	Z/0300/60
25 kg	Plastic keg	Z/0300/70
Zn CAS: 7440-66-6 MW: 65.36 Density at 20°C: 7.14 Boiling point: 907°C Melting point: 419°C ADR/RID: 9, III		EINECS: 231-175-3 UN: 3077 GHS: Aquatic Acute 1, Aquatic Chronic 1 H410 P273

Assay: >99.9%

Guaranteed analysis

Arsenic (As)	<0.05 ppm	Lead (Pb)	<40 ppm
Cadmium (Cd)	<20 ppm	Sulfide (S)	<0.00001%
Copper (Cu)	<30 ppm	Tin (Sn)	<10 ppm
Iron (Fe)	<20 ppm		

Lot analysis on the label.

Zinc powder (metal) <i>for analysis</i>		Certified AR
packaging	Mfr. No	
500 g Plastic jar	Z/0500/53	
1 kg Plastic jar	Z/0500/60	
Zn	EINECS: 231-175-3	
CAS: 7440-66-6	UN: 1436	
MW: 65.36	GHS: Pyr. Sol. 1, Water-react. 1, Aquatic Acute 1, Aquatic Chronic 1	
Density at 20°C: 7.14	H260, H250, H410 P273	
Boiling point: 907°C		
Melting point: 420°C		
ADR/RID: 4.3, 4.2, II		
Assay: >90%		
Guaranteed analysis		
Substances	Iron (Fe)	<50 ppm
reducing KMnO ₄	Total nitrogen (N)	<100 ppm
<i>Lot analysis on the label.</i>		

Zinc metal granulated <i>extra pure, activated with copper</i>		SLR
packaging	Mfr. No	
500 g Plastic jar	Z/0400/53	
Zn	EINECS: 231-175-3	
CAS: 7440-66-6	UN: 3077	
MW: 65.36	GHS: Aquatic Acute 1, Aquatic Chronic 1	
Density at 20°C: 7.14	H410	
Boiling point: 907°C	P273	
Melting point: 419°C		
ADR/RID: 9, III		
Guaranteed analysis		
Arsenic (As)	<1 ppm	Insoluble Acids
		<0.05%
<i>For the determination of Arsenic.</i>		

Zinc metal granulated <i>specified reagents for general laboratory work- extra pure</i>		SLR
packaging	Mfr. No	
500 g Plastic jar	Z/0250/53	
2 kg Plastic jar	Z/0250/61	
Zn	EINECS: 231-175-3	
CAS: 7440-66-6	UN: 3077	
MW: 65.36	GHS: Aquatic Acute 1, Aquatic Chronic 1	
Density at 20°C: 7.14	H410	
Boiling point: 907°C	P273	
Melting point: 419°C		
ADR/RID: 9, III		
Guaranteed analysis		
Identity test	: Pass	

Zinc metal powder <i>pure</i>		
packaging	Mfr. No	
500 g Plastic jar	Z/0450/53	
1 kg Plastic jar	Z/0450/60	
2 kg Plastic jar	Z/0450/61	
Zn	EINECS: 231-175-3	
CAS: 7440-66-6	UN: 1436	
MW: 65.36	GHS: Pyr. Sol. 1, Water-react. 1, Aquatic Acute 1, Aquatic Chronic 1	
Density at 20°C: 7.14	H260, H250, H410 P273	
Boiling point: 907°C		
Melting point: 420°C		
ADR/RID: 4.3, 4.2, II		
Assay: >90%		

Zinc metal stick <i>pure, low in arsenic</i>		
packaging	Mfr. No	
500 g Plastic jar	Z/0550/53	
Zn	EINECS: 231-175-3	
CAS: 7440-66-6	UN: 3077	
MW: 65.36	GHS: Aquatic Acute 1, Aquatic Chronic 1	
Density at 20°C: 7.14	H410	
Boiling point: 907°C	P273	
Melting point: 419°C		
ADR/RID: 9, III		
Guaranteed analysis		
Identity test	: Pass	

Zinc metal foil, about 0.3 mm <i>pure</i>		
packaging	Mfr. No	
500 g Plastic jar	Z/0160/53	
Zn	EINECS: 231-175-3	
CAS: 7440-66-6	UN: 3077	
MW: 65.36	GHS: Aquatic Acute 1, Aquatic Chronic 1	
Density at 20°C: 7.14	H410	
Boiling point: 907°C	P273	
Melting point: 419°C		
ADR/RID: 9, III		
Guaranteed analysis		
Identity test	: Pass	

Zinc acetate dihydrate <i>for analysis</i>		Certified AR
packaging	Mfr. No	
250 g Plastic jar	Z/0700/50	
500 g Plastic jar	Z/0700/53	
1 kg Plastic jar	Z/0700/60	
C ₄ H ₆ O ₄ Zn.2H ₂ O	UN: 3077	
CAS: 5970-45-6	GHS: Acute Tox. 4, Eye Irrit. 2, Aquatic Acute 1	
MW: 219.49	H319, H302, H400	
Density at 20°C: 1.84	P273, P301+P312, P280, P305+P351+P338	
Boiling point: 100°C		
Melting point: 100°C		
ADR/RID: 9, III		
EINECS: 209-170-2		
Assay: >99.5%		
Guaranteed analysis		
Calcium (Ca)	<10 ppm	Sodium (Na)
Copper (Cu)	<5 ppm	Total chloride (Cl)
Iron (Fe)	<5 ppm	Total nitrogen (N)
Lead (Pb)	<5 ppm	Total phosphorus (P)
Magnesium (Mg)	<10 ppm	Total silicon (Si)
Potassium (K)	<10 ppm	Total sulphur (S)
Typical analysis		
Aluminium (Al)	<5 ppm	Molybdenum (Mo)
Barium (Ba)	<5 ppm	Nickel (Ni)
Cadmium (Cd)	<5 ppm	Silver (Ag)
Chromium (Cr)	<5 ppm	Strontium (Sr)
Cobalt (Co)	<5 ppm	Titanium (Ti)
Lithium (Li)	<5 ppm	Vanadium (V)
Manganese (Mn)	<5 ppm	

Lot analysis on the label.

High Volume Solvent Delivery Systems

High-volume solvent delivery systems are available in 10L to 1000L, in Europe only. Please email your contact details and a brief description of your requirements to bege1.sha.quotes@thermofisher.com





Zinc acetate dihydrate <i>specified reagents for general laboratory work- extra pure</i>		SLR
packaging	Mfr. No	
500 g Plastic jar	Z/0650/53	
C ₄ H ₆ O ₄ Zn.2H ₂ O	UN: 3077	
CAS: 5970-45-6	GHS: Acute Tox. 4, Eye Irrit. 2, Aquatic Acute 1	
MW: 219.49	H319, H302, H400	
Density at 20°C: 1.84	P273, P301+P312, P280, P305+P351+P338	
Melting point: 100°C		
ADR/RID: 9, III		
EINECS: 209-170-2		
Assay: >99%		
Guaranteed analysis		
Calcium (Ca)	<100 ppm	Sodium (Na)
Copper (Cu)	<20 ppm	n-Butyl formate
Iron (Fe)	<20 ppm	Total phosphorus (P)
Lead (Pb)	<20 ppm	Total silicon (Si)
Magnesium (Mg)	<100 ppm	Total sulphur (S)
Potassium (K)	<100 ppm	

Zinc bromide <i>specified reagents for general laboratory work- extra pure</i>		SLR
packaging	Mfr. No	
100 g Plastic jar	Z/0750/48	
ZnBr ₂	EINECS: 231-718-4	
CAS: 7699-45-8	UN: 3260	
MW: 225.19	GHS: Skin Corr. 1B, Aquatic Acute 1	
Density at 20°C: 4.22	H314, H400	
Boiling point: 697°C	P301+P330+P331, P280, P305+P351+P338, P310, P273	
Melting point: 394°C		
ADR/RID: 8, II		
Assay: >98%		
Guaranteed analysis		
Iron (Fe)	<10 ppm	Sulfate (SO ₄)
		<0.02%

Zinc carbonate <i>pure, basic</i>		
packaging	Mfr. No	
500 g Plastic jar	Z/0800/53	
Zn(CO ₃)	Density at 20°C: 4.39	
CAS: 5263-02-5	EINECS: 226-076-7	
MW: 549.01		
Guaranteed analysis		
Identity test	: Pass	

Zinc chloride <i>for analysis</i>		Certified AR
packaging	Mfr. No	
500 g Plastic jar	Z/0870/53	
1 kg Plastic jar	Z/0870/60	
ZnCl ₂	GHS: Acute Tox. 4, Skin Corr. 1B, STOT SE 3, Aquatic Acute 1, Aquatic Chronic 1	
CAS: 7646-85-7	H314, H335, H302, H410	
MW: 136.29	P301+P330+P331, P280, P305+P351+P338, P261, P273, P301+P312	
Boiling point: 732°C		
Melting point: 293°C		
ADR/RID: 8, III		
EINECS: 231-592-0		
UN: 2331		
Assay: >98%		
Guaranteed analysis		
Basic salt	: <1%	Potassium (K)
Calcium (Ca)	<10 ppm	Sodium (Na)
Copper (Cu)	<5 ppm	Total phosphorus (P)
Iron (Fe)	<10 ppm	Total silicon (Si)
Lead (Pb)	<5 ppm	Total sulphur (S)
Magnesium (Mg)	<10 ppm	

Typical analysis			
Aluminium (Al)	<5 ppm	Molybdenum (Mo)	<5 ppm
Barium (Ba)	<5 ppm	Nickel (Ni)	<5 ppm
Cadmium (Cd)	<5 ppm	Silver (Ag)	<5 ppm
Chromium (Cr)	<5 ppm	Strontium (Sr)	<5 ppm
Cobalt (Co)	<5 ppm	Titanium (Ti)	<5 ppm
Lithium (Li)	<5 ppm	Vanadium (V)	<5 ppm
Manganese (Mn)	<5 ppm		

Lot analysis on the label.

Zinc chloride <i>specified reagents for general laboratory work- extra pure</i>		SLR
packaging	Mfr. No	
500 g Plastic jar	Z/0850/53	
1 kg Plastic jar	Z/0850/60	
10 kg Plastic keg	Z/0850/68	
ZnCl ₂	GHS: Acute Tox. 4, Skin Corr. 1B, STOT SE 3, Aquatic Acute 1, Aquatic Chronic 1	
CAS: 7646-85-7	H314, H335, H302, H410	
MW: 136.29	P301+P330+P331, P280, P305+P351+P338, P261, P273, P301+P312	
Boiling point: 732°C		
Melting point: 293°C		
ADR/RID: 8, III		
EINECS: 231-592-0		
UN: 2331		
Assay: >97%		
Guaranteed analysis		
Basic salt	: <5%	Potassium (K)
Calcium (Ca)	<2000 ppm	Sodium (Na)
Copper (Cu)	<20 ppm	Total phosphorus (P)
Iron (Fe)	<50 ppm	Total silicon (Si)
Lead (Pb)	<20 ppm	Total sulphur (S)
Magnesium (Mg)	<1000 ppm	

Zinc iodide <i>specified reagents for general laboratory work- extra pure</i>		SLR
packaging	Mfr. No	
100 g Plastic jar	Z/1100/48	
ZnI ₂	EINECS: 233-396-0	
CAS: 10139-47-6	UN: 3260	
MW: 319.18	GHS: Skin Corr. 1B, Eye Dam. 1	
Density at 20°C: 4.74	H314	
Melting point: 446°C	P301+P330+P331, P280, P305+P351+P338, P310	
ADR/RID: 8, III		
Assay: >98%		
Guaranteed analysis		
Iron (Fe)	<100 ppm	Sulfate (SO ₄)
		<0.05%

Zinc nitrate hexahydrate <i>specified reagents for general laboratory work- extra pure</i>		SLR
packaging	Mfr. No	
250 g Plastic jar	Z/1150/50	
500 g Plastic jar	Z/1150/53	
Zn(NO ₃) ₂ .6H ₂ O	UN: 1514	
CAS: 10196-18-6	GHS: Ox. Sol. 2, Acute Tox. 4, Skin Irrit. 2, Eye Irrit. 2, STOT SE 3	
MW: 297.47	H272, H335, H315, H319, H302	
Density at 20°C: 2.06	P210, P261, P302+P352, P280, P305+P351+P338	
Melting point: 36°C		
ADR/RID: 5.1, II		
EINECS: 231-943-8		
Assay: >97%		
Guaranteed analysis		
Calcium (Ca)	<50 ppm	Sodium (Na)
Copper (Cu)	<20 ppm	n-Butyl formate
Iron (Fe)	<20 ppm	Total phosphorus (P)
Lead (Pb)	<50 ppm	Total silicon (Si)
Magnesium (Mg)	<100 ppm	Total sulphur (S)
Potassium (K)	<50 ppm	

Zinc nitrate hexahydrate		
Technical		
packaging	Mfr. No	
250 g Plastic jar	Z/1145/50	
500 g Plastic jar	Z/1145/53	
N ₂ O ₆ Zn.6H ₂ O	UN: 1514	
CAS: 10196-18-6	GHS: Ox. Sol. 2, Acute Tox. 4, Skin Irrit. 2, STOT SE 3	
MW: 297.47	H272, H335, H315, H319, H302	
Density at 20°C: 2.06	P210, P261, P302+P352, P280, P305+P351+P338	
Melting point: 36°C		
ADR/RID: 5.1, II		
EINECS: 231-943-8		
Guaranteed analysis		
Identity test	: Pass	

Zinc oxide		
for analysis		
packaging	Mfr. No	Certified AR
500 g Plastic jar	Z/1350/53	
1 kg Plastic jar	Z/1350/60	
25 kg Fibre Keg	Z/1350/70	
ZnO	EINECS: 215-222-5	
CAS: 1314-13-2	UN: 3077	
MW: 81.38	GHS: Aquatic Acute 1, Aquatic Chronic 1	
Density at 20°C: 5.6	H410	
Melting point: 1975°C	P273	
ADR/RID: 9, III		
Assay: >99.5%		
Guaranteed analysis		
Loss on drying 600°C	: <0.5%	Manganese (Mn) <5 ppm
Reducing substances	: <0.002%	Molybdenum (Mo) <5 ppm
Aluminium (Al)	<5 ppm	Nickel (Ni) <5 ppm
Barium (Ba)	<5 ppm	Potassium (K) <50 ppm
Cadmium (Cd)	<5 ppm	Propan--1-ol <5 ppm
Calcium (Ca)	<50 ppm	Sodium (Na) <50 ppm
Carbonate (Na ₂ CO ₃)	<0.2%	Strontium (Sr) <5 ppm
Chromium (Cr)	<5 ppm	Titanium (Ti) <5 ppm
Cobalt (Co)	<5 ppm	n-Butyl formate <0.002%
Copper (Cu)	<5 ppm	Total phosphorus (P) <50 ppm
Iron (Fe)	<5 ppm	Total silicon (Si) <200 ppm
Lead (Pb)	<50 ppm	Total sulphur (S) <100 ppm
Lithium (Li)	<5 ppm	Vanadium (V) <5 ppm
Magnesium (Mg)	<50 ppm	

Lot analysis on the label.

Zinc oxide		
specified reagents for general laboratory work- extra pure		
packaging	Mfr. No	SLR, Eur.Ph.,BP, USP
500 g Plastic jar	Z/1300/53	
1 kg Plastic jar	Z/1300/60	
ZnO	EINECS: 215-222-5	
CAS: 1314-13-2	UN: 3077	
MW: 81.38	GHS: Aquatic Acute 1, Aquatic Chronic 1	
Density at 20°C: 5.6	H410	
Melting point: 1975°C	P273	
ADR/RID: 9, III		
Assay: >99%		
Guaranteed analysis		
Loss on ignition	<1%	Potassium (K) <200 ppm
Calcium (Ca)	<200 ppm	Sodium (Na) <200 ppm
Copper (Cu)	<20 ppm	n-Butyl formate <0.02%
Iron (Fe)	<20 ppm	Total phosphorus (P) <200 ppm
Lead (Pb)	<200 ppm	Total silicon (Si) <500 ppm
Magnesium (Mg)	<200 ppm	Total sulphur (S) <500 ppm

Zinc stearate		
Technical		
packaging	Mfr. No	
500 g Plastic jar	Z/1400/53	
C ₃₆ H ₇₀ O ₄ Zn	Melting point: 118 to 128°C	
CAS: 557-05-1	Flash point: 100°C	
MW: 632.34	EINECS: 209-151-9	
Guaranteed analysis		
Identity test	: Pass	

Zinc sulfate heptahydrate		
for analysis		
packaging	Mfr. No	Certified AR
500 g Plastic jar	Z/1600/53	
5 kg Plastic keg	Z/1600/65	
25 kg Plastic keg	Z/1600/70	
ZnSO ₄ .7H ₂ O	UN: 3077	
CAS: 7446-20-0	GHS: Acute Tox. 4, Eye Dam. 1, Aquatic Acute 1, Aquatic Chronic 1	
MW: 287.53	H318, H410, H302	
Density at 25°C: 3.54	P280, P305+P351+P338, P310, P273, P301+P312	
Melting point: 100°C		
ADR/RID: 9, III		
EINECS: 231-793-3		
Assay: 99.5 to 104%		
Guaranteed analysis		
pH (aqueous solution 5%)	4.5 to 6	Potassium (K) <50 ppm
Calcium (Ca)	<50 ppm	Sodium (Na) <50 ppm
Copper (Cu)	<5 ppm	Total chloride (Cl) <0.001%
Iron (Fe)	<5 ppm	Total nitrogen (N) <20 ppm
Lead (Pb)	<5 ppm	Total silicon (Si) <50 ppm
Magnesium (Mg)	<50 ppm	
Typical analysis		
Aluminium (Al)	<5 ppm	Molybdenum (Mo) <5 ppm
Barium (Ba)	<5 ppm	Nickel (Ni) <5 ppm
Cadmium (Cd)	<5 ppm	Silver (Ag) <5 ppm
Chromium (Cr)	<5 ppm	Strontium (Sr) <5 ppm
Cobalt (Co)	<5 ppm	Titanium (Ti) <5 ppm
Total phosphorus (P)	<50 ppm	Vanadium (V) <5 ppm
Total silicon (Si)	<200 ppm	
Total sulphur (S)	<100 ppm	
Vanadium (V)	<5 ppm	

Zinc sulfate heptahydrate		
specified reagents for general laboratory work- extra pure		
packaging	Mfr. No	SLR
500 g Plastic jar	Z/1550/53	
1 kg Plastic jar	Z/1550/60	
2 kg Plastic jar	Z/1550/61	
5 kg Plastic keg	Z/1550/65	
ZnSO ₄ .7H ₂ O	UN: 3077	
CAS: 7446-20-0	GHS: Acute Tox. 4, Eye Dam. 1, Aquatic Acute 1, Aquatic Chronic 1	
MW: 287.53	H318, H410, H302	
Density at 20°C: 3.54 to 25°C	P280, P305+P351+P338, P310, P273, P301+P312	
Melting point: 100°C		
ADR/RID: 9, III		
EINECS: 231-793-3		
Assay: >99%		
Guaranteed analysis		
Calcium (Ca)	<200 ppm	Potassium (K) <200 ppm
Copper (Cu)	<20 ppm	Sodium (Na) <200 ppm
Iron (Fe)	<20 ppm	n-Butyl formate <0.005%
Lead (Pb)	<20 ppm	Total silicon (Si) <200 ppm
Magnesium (Mg)	<200 ppm	

Zinc sulfate solution 0.05M (0.05N)		
NIST Standard Concentrate, For Volumetric analysis		
packaging	Mfr. No	
1 Ampoule	J/8490C/05	
6 Ampoule	J/8490C/90	
Boiling point: 100°C	H411, H318	
Melting point: 0	P280, P305+P351+P338, P310, P273	
GHS: Eye Dam. 1, Aquatic Chronic 2		

Each ampoule makes 1 ℓ of solution.
See the full range of Standard Volumetric solution page 323.

Zinc sulfide		
pure		
packaging	Mfr. No	
500 g Plastic jar	Z/1700/53	
ZnS	EINECS: 215-251-3	
CAS: 1314-98-3	EUH031	
MW: 97.44	P233	
Density at 20°C: 4.1		
Guaranteed analysis		
Identity test	: Pass	

Zirconium standard metal solution 1000 ppm		
for AAS, solution HCl 1M		
packaging	Mfr. No	
500 mℓ HDPE Plastic bottle	J/8200/08	
ADR/RID: 8,		
III UN: 1789		
Manufactured to a tolerance of ± 0.5%. See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.		

Zirconium standard metal solution 1000 ppm		
for ICP solution in HNO ₃ 2%, NIST certified		
packaging	Mfr. No	
100 mℓ Plastic bottle	SPPLZR22Y	
250 mℓ Plastic bottle	SPPLZR22T	
500 mℓ Plastic bottle	SPPLZR22X	
H315, H319		
Manufactured to a tolerance of ± 0.5%. All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities. See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314. Check with your local distributor for availability of Spex product		

Zirconium standard metal solution 1000 ppm		
for ICP solution in HCL 10%, NIST certified		
packaging	Mfr. No	
500 mℓ Plastic bottle	SPPLZR32X	
H315, H319, H335		
Manufactured to a tolerance of ± 0.5%. All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities. See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314. Check with your local distributor for availability of Spex product		

Zirconium standard metal solution 10 000 ppm		
for ICP solution in HNO ₃ 5%, NIST certified		
packaging	Mfr. No	
100 mℓ Plastic bottle	SPPLZR23Y	
500 mℓ Plastic bottle	SPPLZR23X	
H314		

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Zirconium standard metal solution 10 000 ppm		
for ICP solution in HCL 10%, NIST certified		
packaging	Mfr. No	
500 mℓ Plastic bottle	SPPLZR33X	
H315, H319, H335		

Manufactured to a tolerance of ± 0.5%.
All standards are manufactured from acids, metals & metal salts and contain less than <10 ppm total metal impurities.
See also the full range in the Standard Elemental solution for Trace metal analysis section see page 314.
Check with your local distributor for availability of Spex product

Zirconium dioxide		
specified reagents for general laboratory work- extra pure		
packaging	Mfr. No	SLR
250 g Plastic jar	Z/1850/50	
ZrO ₂	Melting point: >2700°C	
CAS: 1314-23-4	EINECS: 215-227-2	
MW: 123.22	GHS: Eye Irrit. 2, STOT SE 3	
Density at 20°C: 5.6	H335, H319	
Boiling point: 4300°C	P261, P280, P305+P351+P338	
Guaranteed analysis		
Identity test	: Pass	

Zirconyl chloride octahydrate		
specified reagents for general laboratory work- extra pure		
packaging	Mfr. No	SLR
500 g Plastic jar	Z/1900/53	
ZrCl ₂ O.8H ₂ O	GHS: Acute Tox. 4, Skin Corr. 1B, Eye Dam. 1	
CAS: 13520-92-8	H314, H302	
MW: 322.25	P301+P330+P331, P280, P305+P351+P338, P310, P301+P312	
Density at 20°C: 1.91		
ADR/RID: 8, II		
EINECS: 231-717-9		
UN: 3260		
Guaranteed analysis		
Identity test	: Pass	

Standard Volumetric solution

For your convenience we offer a range of NIST traceable solutions. A full list of our products can be found on page 323. Can't see what you need? Please email your contact details and a brief description of your requirements to begeel.sha.quotes@thermofisher.com





- **Molecular Biology** p. 385
- **Protein Chemistry** p. 438
- **Cell Biology** p. 452
- **Core Bioreagents & complementary products** p. 479
- **Index by product name** p. 499
- **Index by manufacturer number** p. 501

Fisher BioReagents®

End-to-end solutions for molecular biology, protein chemistry and cell biology

Agar
Agarose
Amino Acids
Antibiotics: Bacterial Selection
Antibiotics: Cell Culture
Biological Buffers
Cell Regulation Compounds
DNA Electrophoresis Buffers
DNA Electrophoresis Products
DNA Ladders
DNA Modifying Enzymes
DNA Precipitation Products
dNTPs
Microbial Media
Northern Blot Products
Nucleic Acid Purification
Nucleic Acid Sequencing
Peptide Sequencing Products
Peptide Synthesis Products
PCR Enzymes
Protein Electrophoresis Buffers
Protein Purification
Protein Stains
Restriction Enzymes
Reverse Transcriptases
RNA Purification Kits
RNA Purification
Southern Blot
Sugar Additives
Tissue Stains
Vitamins
Water
Western Blot

Vital reagents for life science research

- Ultra-high purity
- Prequalified for the life science application
- Fisher BioReagents are among the finest in the industry

Over 1,000 products

Our extensive line includes products for:

- Nucleic Acid Electrophoresis, Purification, Hybridisation and Sequencing
- Polymerase Chain Reaction (PCR)
- Protein Purification and Electrophoresis
- Buffers and Detergents
- Cell Biology/Cell Culture
- Microbiology/Immunology

Purity grades to meet your specific needs

Reagents are prequalified and guaranteed to be suitable for the designated technique.

Stringent specifications

Address critical factors such as purity, water content, levels of contaminants and absence of DNase, RNase or protease activity.

Innovative packaging design

Packaging designed for safety, convenient handling and storage, and preservation of product integrity.

Manufacturing quality

Rigorous QA standards, supplier certification process, customer site audit, cGMP-certified manufacturing facilities, ISO 9001:2008 Quality Management System.

Bench to batch

Fisher BioReagents are packaged in sizes to meet your needs.



Fisher BioReagents®: Purity Grades for Every Application

Material Grade	Definition
Certified	Reagent chemicals for which the purity standard is established by Fisher Chemical. Purity is guaranteed to meet published maximum limits of impurities.
DNA Grade	Designates reagents suitable for use in molecular biology applications involving the manipulation of DNA. Tested for specific contaminants, such as DNase and protease.
DNA Synthesis	Designates reagents suitable for use with automated DNA synthesis instrumentation.
Electrophoresis	Material used specifically for electrophoresis applications.
Genetic Analysis Grade	Material that is specially prepared for various molecular cloning applications. Tested for specific contaminants, such as DNase and RNase.
IEF Grade	Material suitable for use with isoelectric focusing of proteins.
Islet Isolation Grade	Material suitable for isolation of pancreatic islets.
Molecular Biology Grade	Designates reagents suitable for use in molecular biology applications. Tested for specific contaminants, such as nucleases and bacteria, where appropriate.
Molecular Genetics	Reagent chemicals that have been specifically purified and assayed for molecular genetics applications.
PCR Grade	Material suitable for use in Polymerase Chain Reaction (PCR).
Peptide Synthesis	Designates reagents suitable for use with protein synthesis instrumentation.
Protein Electrophoresis Grade	Material used specifically for protein electrophoresis applications.
Sequencing	Material designed for use with automated DNA or protein sequencing equipment.
Super Pure	Material with a purity level exceeding the various monograph grades.
Tissue Culture Grade	Materials of superior quality where there are no published standards and that are suitable for use in tissue culture applications.

Packaging for safety, convenience and product quality

Fisher BioReagents come in a wide variety of innovative packaging designed for safety, environmental protection, convenient handling and storage, and preservation of product integrity. The primary container is included in the product description of most chemicals in this catalogue.

Primary containers include:

- Plastic and glass bottles, jars
- Specialised acid containers
- Square poly bottles
- Sterility proof sachets
- Poly pails
- Polypac™ containers
- Compact, laminated boxes



How to use this catalogue

Refer to the index page: “BioReagents by Product Name” to look for the page number of the product you’re interested in. Products are also Categorised by application: Molecular Biology, Protein Chemistry, Cell Biology and Core BioReagents (used for multiple applications). Refer to the appropriate chapter to look for the product you’re interested in.

- 1

Product name
- 2

Quality grade
- 3

CAS number
- 4

Packaging and ordering information
- 5

Product specifications
- 6

GHS code
- 7

Application information
- 8

Recommended storage conditions

1

2

3

4

5

6

Agarose

Low-Melting, <1kb DNA/RNA

Genetic Analysis Grade

packaging

100g Polybottle

Cat. No.

1058-3355

Mfr. No.

BP1360-100

CAS

9012-36-6

MW

334.12

DNase

Not detected

PRO (A-M)

<0.15

Cell Strength

>500g/cm²

Color at Temperature

Colorless

Melting Temperature

135°C

Melting Concentration

<0.1%

Residue

Not detected

Sulfate Content

<0.1%

Applications: Low melting temperature agarose. Certified recovery of small nucleic acid fragments. Outstanding resolution of PCR and RT-PCR fragments from 50 to 1000 bp. In-gel PCR, In-gel gel electrophoresis, transformations. Recommended Storage: R

Fisher BioReagents®

Your source for high purity products for nucleic acid electrophoresis

All Fisher BioReagents agaroses are DNase- and RNase-free to ensure optimal results for your nucleic acid application.

Fisher BioReagents offers 3 different grades of agarose that are functionally tested and pre-qualified for specific applications:

- **Genetic Analysis Grade** – Agarose that yields biologically active DNA or RNA. Testing includes enzymatic performance measurements.
- **Molecular Biology Grade** – Agarose that is suitable for analytical separation of DNA or RNA.
- **PCR Grade** – Agarose that is suitable for the analytical separation of PCR amplicons (<1kb).

Agarose Selection Guide					
Type of Agarose	Low EEO	Low Melting (>200bp)	Low Melting (<1kb)	Wide Separation Range	PCR Grade
Cat. No	BP160	BP165	BP1360	BP1356	BP2410
Recovery of DNA or RNA	•	•	•	•	•
Southern and Northern Blots	•				
DNA/RNA separation 50bp to 1kb			•		•
DNA/RNA separation >1kb	•	•		•	
PCR fragment analysis	•	•		•	•
In-gel reactions (ligation, transformations, PCR)			•		
Colony lifts	•				
Available pack sizes	100g and 500g	25g	100g	100g and 500g	100g
Agarose grade	Molecular Biology	Molecular Biology	Genetic Analysis	Genetic Analysis	PCR Grade

Buffers for DNA Electrophoresis Applications

Two buffers commonly used for DNA electrophoresis are Tris-acetate with EDTA and Tris-borate with EDTA. Because the pH of these buffers is neutral, the phosphate backbone of DNA has a net negative charge and migrates to the anode. TAE and TBE have different properties which makes one more suitable than the other for a specific purpose.

TAE: DNase-, RNase- and Protease-free			TBE: DNase- and RNase-free		
Cat. No.	Concentration	Size	Cat. No.	Concentration	Size
BP2434-4	1X	4L	BP2430-1	1X	1L
BP2434-20	1X	20L	BP2430-4	1X	4L
BP1335-500	10X	500mL	BP2430-20	1X	20L
BP1335-1	10X	1L	BP1396-86	5X	1L*
BP1335-4	10X	4L	BP1333-1	10X	1L
BP1335-20	10X	20L	BP1333-4	10X	4L
BP1330-1	25X	1L	BP1333-20	10X	20L
BP1332-500	50X	500mL	BP1334-1	10X	1L**
BP1332-1	50X	1L			
BP1332-4	50X	4L			
BP1332-20	50X	20L			
BP1331-1	25X	1L**			

*Pre-weighed powder in poly bottle. Dissolve in water.
**Pre-weighed powder in foil pack. Dissolve in water.

Buffer Components for DNA Electrophoresis

Buffer Component		
	Size	Cat. No.
Tris Base	500g	BP152-500
	1kg	BP152-1
	5kg	BP152-5
	10kg	BP152-10
	25kg	BP152-25
Boric Acid	500g	BP168-500
	1kg	BP168-1
EDTA Disodium Salt	500g	BP120-500
	1kg	BP120-1

Buffers for RNA Electrophoresis Applications


MOPS is a commonly used buffer system for RNA electrophoresis using formaldehyde or formamide-denatured RNA. It is important to use RNase-free chemicals, water and containers when preparing the buffer solution. The typical formulation of 10X MOPS running buffer is 0.4M MOPS (pH 7.0), 0.1M sodium acetate and 0.01M EDTA.

MOPS: DNase-, RNase- and Protease-free		
Cat. No.	Description	Size
BP308-100	Powder	100g
BP308-500	Powder	500g
BP2900-500	10X Buffer Solution	500mL
BP2900-1	10X Buffer Solution	1L
Water		
BP2484-50	Nuclease-Free	50mL
BP2484-100	Nuclease-Free	100mL
BP2470-1	DNA Grade	1L
BP561-1	RNA Grade	1L
Formaldehyde		
BP531-25	37% by weight	25mL
BP531-500	37% by weight	500mL



NEW!
Ethanol, Molecular Biology Grade, is an ultrapure molecular biology grade ethanol used for the purification and precipitation of biomolecules such as nucleic acids and proteins.

Cat. No.	Size
BP2818-100	100mL
BP2818-500	500mL
BP2818-4	4 L



exACTGene® and routine DNA ladders

Ready-to-use (pre-mixed with loading dye), room temperature stable DNA ladders are available for all common electrophoresis applications.

exACTGene DNA ladders are ideal for qualitative analysis, quantitative estimation and size assessment				
Cat. No.	Application	Size Range	Number of Bands	Number of Loadings
BP2570-100	PCR fragment analysis	25 to 650bp	14	100/10µL
BP2571-100	PCR fragment analysis, small DNA digests	25 to 1,000bp	12	100/10µL
BP2572-100	Quick check of PCR or enzyme digestion results	50 to 2,000bp	8	100/10µL
BP2573-100	General purpose, small DNA fragments	100 to 1,000bp	10	100/10µL
BP2574-100	Fast run times, small DNA fragments	100 to 2,000bp	11	100/10µL
BP2575-100	Clone identification	100 to 2,686bp	14	100/10µL
BP2576-100	Large size PCR or cloning	300 to 5,000bp	10	100/10µL
BP2577-100	Small and large cloning application	100 to 5,000bp	16	100/10µL
BP2578-100	General purpose, large digested DNA	300 to 10,000bp	13	100/10µL
BP2579-100	General purpose, wide separation range	100 to 10,000bp	19	100/10µL
BP2580-100	General purpose, extra large DNA fragments	300 to 24,000bp	15	100/10µL



Fisher BioReagents®
Proteomics Products
Convenience, Quality and Consistency

EZ-Run™ Protein Gel Solution

- Ready to use
- Superior resolution
- Wide separation range on same mini-gel
- No stacking gel required
- Proprietary gel chemistry
- Stable for two years at room temperature
- Compatible with all conventional staining methods
- Suitable for post-electrophoresis applications such as Western Blot transfer and MALDI analysis

EZ-Run™ Protein Gel Solutions		
Description	Quantity	Cat. No.
10% EZ-Run Protein Gel Solution	100mL	BP7710-100
	500mL	BP7710-500
12.5% EZ-Run Protein Gel Solution	100mL	BP7712-100
	500mL	BP7712-500
15% EZ-Run Protein Gel Solution	100mL	BP7715-100
	500mL	BP7715-500
20X Running Buffer for EZ-Run Protein Gel Solution	500mL	BP7700-500

Acrylamide, Bis-Acrylamide and Catalysts		
Description	Quantity	Cat. No.
Acrylamide	100g	BP170-100
	500g	BP170-500
	5kg	BP170-500
Acrylamide Solution, 40%	1L	BP1402-1
Bis-Acrylamide	25g	BP171-25
	100g	BP171-100
Bis-Acrylamide Solution, 2%	250mL	BP1404-250
Acrylamide:Bis-Acrylamide, 19:1, Powder	100g	BP1364-100
Acrylamide:Bis-Acrylamide, 29:1, Powder	100g	BP1366-100
Acrylamide:Bis-Acrylamide, 37.5:1, Powder	100g	BP1368-100
Acrylamide:Bis-Acrylamide, 19:1, 40% Solution	1L	BP1406-1
Acrylamide:Bis-Acrylamide, 29:1, 40% Solution	1L	BP1408-1
Acrylamide:Bis-Acrylamide, 37.5:1, 40% Solution	1L	BP1410-1
Ammonium Persulfate	25g	BP179-25
	100g	BP179-100
Sodium Persulfate	1kg	BP2637-1
TEMED	20g	BP150-20
	100g	BP150-100

Buffers for Protein Electrophoresis		
Description	Quantity	Cat. No.
Tris-Glycine, 10X Solution	1L	BP1306-1
	4L	BP1306-4
	1L*	BP1307-1
Tris-Glycine-SDS, 10X Solution	1L	BP1341-1
	4L	BP1341-4
Tris-Glycine-SDS, 5X Solution, Pre-Weighed Powder*	1L	BP1398-92
Tris-Glycine-SDS, 10X Solution, Pre-Weighed Powder*	1L	BP1342-1
PBS, 10X Solution	500mL	BP399-500
PBS, 10X Solution	1L	BP399-1
PBS, 10X Solution	4L	BP399-4
PBS, 10X Solution	20L	BP399-20
PBS Tablets, 1X Solution	100 tabs**	BP2944-100
PBS with Tween 20, 1X Solution, Pre-Weighed Powder*	10 foil pouches	BP2938-10
Tris-Buffered Saline, 10X (pH 7.4)	100mL	BP2471-100
	500mL	BP2471-500
	1L	BP2471-1
Tris Base	500g	BP152-500
	1kg	BP152-1
	5kg	BP152-5
	10kg	BP152-10
Glycine	25kg	BP152-25
	500g	BP381-500
	1kg	BP381-1
	5kg	BP381-5

*Pre-weighed powder to make 1L. Dissolve in water.
**One tablet dissolved in 200mL water yields 0.01M phosphate buffer, 0.0027M KCl, and 0.137M NaCl, pH 7.4 at 25°C.






Detergents and Denaturing Agents		
Description	Quantity	Cat. No.
Brij 35	500g	BP345-500
CHAPS	1g	BP571-1
	5g	BP571-5
SDS	100g	BP166-100
	500g	BP166-500
	5kg	BP166-5
SDS, 10% Solution	200mL	BP2436-200
	1L	BP2436-1
SDS, 20% Solution	200mL	BP1311-200
	1L	BP1311-1
Triton X-100	100mL	BP151-100
	500mL	BP151-500
Tween 20	100mL	BP337-100
	500mL	BP337-500
Tween 80	500mL	BP338-500

EZ-Run™ Protein Standards

- Highly purified markers and ladders provide compact and clear bands
- Prestained standards are indispensable in monitoring protein separation and transfer efficiency
- Reference bands allow quick gel progress assessment
- Unstained standards are most suitable for precise sizing of proteins
- All standards are supplied in loading buffer and are ready to use



Description	MW Range	No. of Bands	Reference Band	Source	Quantity	Cat. No.
Unstained Protein Standards	14.4 to 118.0 kDa	7		Native Proteins	500uL	BP3600-500
					2 x 500uL	BP3600-1
	10.0 to 200.0 kDa	14	50 kDa	Recombinant Proteins	500uL	BP3602-500
					2 x 500uL	BP3602-1
Prestained Protein Standards	20.0 to 118.0 kDa	6		Native Proteins	500uL	BP3601-500
					2 x 500uL	BP3601-1
	11.0 to 170.0 kDa	10	72 kDa	Recombinant Proteins	500uL	BP3603-500
					2 x 500uL	BP3603-1

Agarose Broad Separation Range for DNA/RNA		Genetic Analysis Grade
packaging		Mfr. No
100 g	PolyBottle	BP1356-100
500 g	PolyBottle	BP1356-500
C ₁₂ H ₁₈ O ₉ CAS: 9012-36-6 MW: 306.12		
DNase	Not detected	
EEO (-Mr)	<0.15	
Gel Strength	>1200g/cm2	
Gelation Temperature	34.5°-37.5°C	
Moisture Content	<10%	
RNase	Not detected	
Sulfate Content	<0.15%	
Applications: Routine nucleic acid electrophoresis; broad resolution range 500 bp to 25kb. Certified for recovery of DNA and RNA. Reliable digestion and ligation of recovered DNA or RNA fragments. Cloning.		
Recommended Storage: RT		

Agarose Intermediate Melting		PCR Grade
packaging		Mfr. No
100 g	AmberGlass	BP2410-100
C ₁₂ H ₁₈ O ₉ CAS: 9012-36-6 MW: 306.12		
DNase	Not detected	
EEO (-Mr)	<=0.12	
Gel Strength	>=1200g/cm2	
Gelation Temperature	<35.5°C	
Melting Temperature	80°C	
RNase	Not detected	
Sulfate Content	<=0.11%	
Applications: Separation of low MW DNA <100 bp. PCR analysis applications.		
Recommended Storage: RT		

Agarose High-EEO		Protein Electrophoresis Grade
packaging		Mfr. No
100 g	PolyBottle	BP162-100
C ₁₂ H ₁₈ O ₉ CAS: 9012-36-6 MW: 306.12		
DNase	Not detected	
EEO (-Mr)	0.23-0.26	
Gel Strength	>650g/cm2	
Gelation Temperature	34.5°-37.5°C	
Moisture Content	<10%	
RNase	Not detected	
Sulfate Content	<0.15%	
Applications: Immunoelectrophoresis (IEP). Crossed-IEP. Counter-immunoelectrophoresis (CIEP). Serum protein electrophoresis.		
Recommended Storage: RT		

Agarose Low-Melting, Nucleic Acid Recovery		Molecular Biology Grade
packaging		Mfr. No
25 g	PolyBottle	BP165-25
C ₁₂ H ₁₈ O ₉ CAS: 9012-36-6 MW: 306.12		
DNase	Not detected	
EEO (-Mr)	<0.10	
Gel Strength	>200g/cm2	
Gelation Temperature	26°-30°C	
Melting Temperature	65°C	
Moisture Content	<10%	
RNase	Not detected	
Sulfate Content	<0.10%	
Applications: Original low melting temperature agarose. Nucleic acid separation from 200 bp to 25kb. Preparative DNA and RNA electrophoresis. Ideally suited for DNA and RNA recovery. Cloning of tissue culture cells and viral plaque assays.		
Recommended Storage: RT		

Agarose High Gelling Temperature		DNA Grade
packaging		Mfr. No
25 g	PolyBottle	BP164-25
100 g	PolyBottle	BP164-100
500 g	PolyBottle	BP164-500
C ₁₂ H ₁₈ O ₉ CAS: 9012-36-6 MW: 306.12		
DNase	Not detected	
EEO (-Mr)	<0.1	
Gel Strength	>800g/cm2	
Gelation Temperature	<40.5°-43.5°C	
Moisture Content	<10%	
RNase	Not detected	
Sulfate Content	<0.15%	
Recommended Storage: RT		

Agarose Low-EEO/Multi-Purpose		Molecular Biology Grade
packaging		Mfr. No
100 g	PolyBottle	BP160-100
500 g	PolyBottle	BP160-500
C ₁₂ H ₁₈ O ₉ CAS: 9012-36-6 MW: 306.12		
DNase	Not detected	
EEO (-Mr)	0.09-0.13	
Gel Strength	>1200g/cm2	
Gelation Temperature	34°-45°C	
Moisture Content	<10%	
RNase	Not detected	
Sulfate Content	<0.15%	
Applications: Routine electrophoresis of DNA and RNA; wide resolution range 500 bp to 23kb. High gel strength ideal for Southern and Northern blotting. PCR >1kb. Immunoprecipitation techniques. Baculovirus screening. Colony lifts.		
Recommended Storage: RT		

Agarose Low-Melting, <1kb DNA/RNA		Genetic Analysis Grade
packaging	Mfr. No	
100 g PolyBottle	BP1360-100	
C ₁₂ H ₁₈ O ₉ CAS: 9012-36-6 MW: 306.12		
DNase Not detected EEO (-Mr) <0.15 Gel Strength >500g/cm2 Gelation Temperature <35°C Melting Temperature 65°C Moisture Content <10% RNase Not detected Sulfate Content <0.15%		

Applications: Low melting temperature agarose. Certified recovery of small nucleic acid fragments. Outstanding resolution of PCR and RT-PCR fragments from 50 to 1000 bp. In-gel PCR. In-gel ligations/transformations.
Recommended Storage: RT

Agarose Medium-EEO		Protein Electrophoresis Grade
packaging	Mfr. No	
100 g PolyBottle	BP161-100	
C ₁₂ H ₁₈ O ₉ CAS: 9012-36-6 MW: 306.12		
DNase Not detected EEO (-Mr) 0.16-0.19 Gel Strength >1000g/cm2 Gelation Temperature 34.5°-37.5°C Moisture Content <10% RNase Not detected Sulfate Content <0.15%		

Applications: Ideal alternative to polyacrylamide for serum protein electrophoresis. Immunoelectrophoresis. Electrophoresis of nucleic acids.
Recommended Storage: RT

Agarose Agarose Gel-Loading Dye, 6X		Molecular Biology Grade
packaging	Mfr. No	
5 mL PolyBottle	BP633-5	
description Dark purple slightly viscous liquid DNase Not detected Optical Absorbance (dilution 1:500 with deionized water) at 525nm 0.3-0.6 Optical Absorbance (dilution 1:500 with deionized water) at 588-594nm 0.5-0.8 Optical Absorbance (dilution 1:500 with deionized water) at 635-641nm 0.5-0.7 RNase Not detected		

Applications: This tracking dye is added to DNA and RNA samples prior to electrophoresis on agarose gels.
Components: Amaranth (<1.0%), Ficoll* (15%), Proprietary Component (<1.0%), and Water.
[915-67-3 (Amaranth)] ; [26873-85-8 (Ficoll*)]
Recommended Storage: RT

Bromophenol Blue Free Acid		
packaging	Mfr. No	
25 g AmberGlass	BP115-25	
C ₁₉ H ₁₀ Br ₄ O ₅ S CAS: 115-39-9 MW: 669.96 EINECS: 204-086-2		
E ₁ % _{1cm} >=1000g-1cm-1 Lambda Max. 591nm ±2nm Loss on Drying (at 105°C) <=5% Solubility To pass test		

Applications: Used as a tracking dye in electrophoresis.
Recommended storage: RT

Bromophenol Blue Sodium Salt		
packaging	Mfr. No	
25 g AmberGlass	BP114-25	
C ₁₉ H ₉ Br ₄ NaO ₅ S CAS: 62625-28-9 MW: 691.94 H335, H319, H315		
P261, P302+P352, P280, P305+P351+P338		
Clarity of Solution To pass test E ₁ % _{1cm} >=800g-1cm-1 Lambda Max. in 0.05M Sodium Acetate 585-595nm Visual Transition Interval To pass test		

Applications: Commonly used as a tracking dye in electrophoresis.
Recommended storage: RT

Boric Acid Crystalline		Electrophoresis
packaging	Mfr. No	
500 g PolyBottle 1 kg PolyBottle	BP168-500 BP168-1	
H ₃ BO ₃ CAS: 10043-35-3 MW: 61.83 EINECS: 233-139-2		
H360FD P310, P201, P308+P313, P304+P340		
Assay >=99.0% Calcium (Ca) <=0.005% Chloride (Cl) <=0.001% DNase Not detected Electrophoresis To pass test Heavy Metals (Pb) <=0.001% Insoluble in Methanol <=0.005% Iron <=0.001% Nonvolatile with Methanol <=0.05% Phosphate (PO ₄) <=0.001% Sulfate (SO ₄) <=0.010%		

Applications: This high-purity electrophoresis-grade reagent is used in the preparation of Tris/Borate/EDTA electrophoresis buffers.
Recommended Storage: RT

Ethidium Bromide		
packaging	Mfr. No	
1 g AmberGlass 5 g AmberGlass	BP102-1 BP102-5	
C ₂₁ H ₂₀ BrN ₃ CAS: 1239-45-8 MW: 394.3 EINECS: 214-984-6		H341, H330, H302 P310, P281, P304+P340
E ₁ % _{1cm} >=100ℓ g-1cm-1 Lambda Max. 480nm ±2nm Loss on Drying (at 105°C) <=5% Solubility To pass test		

Applications: Fluorometric detection of double stranded nucleic acids. Also acts as an RNA polymerase inhibitor, and in separation of high molecular weight DNAs.
Recommended storage: RT

Ethidium Bromide 1% Solution		Molecular Biology
packaging	Mfr. No	
10 mL AmberPolyBottle	BP1302-10	
C ₂₁ H ₂₀ BrN ₃ CAS: 1239-45-8 MW: 394.3 EINECS: 214-984-6		H331, H341 P310, P281, P304+P340
Appearance Clear, reddish-brown liquid DNase Not detected Electrophoresis To pass test E ₁ % _{1cm} >100 absorbance units Lambda Max. 480nm ±2nm Protease Not detected RNase Not detected		

Filtered through a 0.2-micron filter.
Recommended storage: RT

Formaldehyde 37% by Weight		Molecular Biology
packaging	Mfr. No	
25 mL AmberGlass 500 mL AmberGlass,EcoSafPak*	BP531-25 BP531-500	
CH ₂ O CAS: 50-00-0 MW: 30.02 EINECS: 200-001-8 H226, H314, H317, H370, H311, H331, H301, H351		P280, P301+P330+P331, P305+P351+P338, P304+P340, P302+P352, P301+P312, P210

Assay (HCHO)	>=36.5-38%
Chloride (Cl)	<=5ppm
Color (APHA)	<=10
Heavy Metals (Pb)	<=5ppm
Iron	<=2ppm
Methanol	10-15%
pH (at 25°C)	3 to 4
Residue after ignition	<=0.005%
Sulfate (SO ₄)	<=0.002%
Titrateable Acid	<=0.006mEq/g

Applications: Molecular Biology-grade Formaldehyde is used for denaturing RNA.
Components: 37% Formaldehyde, 15% Methyl Alcohol and 48% Water.
[50-00-0 (Formaldehyde)] ; [67-56-1 (Methyl Alcohol)]
Recommended Storage: RT
UN 1198; DOT Class 3:Flammable Liquid
EcoSafPak* is an environmentally friendly packaging system made of 100% recyclable material by an SFI certified supplier.

Glycerol Gel-Loading Dye, 5X Contains 30% Glycerol		Molecular Biology
packaging	Mfr. No	
1 mL PolyTube 5 mL PolyBottle	BP645-1 BP645-5	
DNase Not detected Optical Absorbance (dilution 1:500 with deionized water) at 525nm 0.25-0.45 Optical Absorbance (dilution 1:500 with deionized water) at 588-594nm 0.30-0.45 Optical Absorbance (dilution 1:500 with deionized water) at 635-641nm 0.25-0.40 RNase Not detected		

Applications: This dye is added to DNA and RNA samples prior to electrophoresis on agarose gels.
Components: Glycerol (30%), Proprietary Components (<6%), and Water [56-81-5 (Glycerol)]
Recommended Storage: RT

Glyoxal 40% w/v		Electrophoresis
packaging	Mfr. No	
500 mL AmberGlass	BP1370-500	
C ₂ H ₂ O ₂ CAS: 107-22-2 MW: 58.04 EINECS: 203-474-9 H341, H332, H319, H315, H317		P280, P261, P308+P313, P305+P351+P338, P302+P352, P333+P313
Density (at 20°C) 1.26-1.29g/mL Electrophoresis of RNA To pass test pH (at 25°C) 2.0 to 3.5		

Applications: Glyoxal is used to denature RNA prior to agarose gel electrophoresis. Since Glyoxal is readily oxidized by air, it should be deionized with a mixed-bed resin, dispensed into aliquots, and then stored at -20°C.
Recommended Storage: 4°C
UN 1760; DOT Class 8:Corrosive



Methanol Peroxide-free		Sequencing
packaging		Mfr. No
1 ℓ	AmberGlass, EcoSafPak*	BP1105-1
4 ℓ	AmberGlass	BP1105-4
CH ₃ O	P301+P310, P280,	
CAS: 67-56-1	P302+P350, P304+P340,	
MW: 32.04	P210, P240	
EINECS: 200-659-6		
H225, H301, H311, H331, H370		
Acetone		To pass test
Assay		>=99.9%
Color (APHA)		<=5
Fluorescence Background (as Quinine Sulfate)		To pass test
Fluorescent Derivatizable Amines		To pass test
IR		Conforms to standard
LC Gradient Suitability		To pass test
Optical Absorbance at 205nm		<=1.00
Optical Absorbance at 220nm		<=0.30
Optical Absorbance at 230nm		<=0.15
Optical Absorbance at 254nm		<=0.025
Peroxides		<=0.001%
Reactive Impurities with H ₂ SO ₄		Not detected
Refractive Index (at 25°C)		1.3260-1.3300
Residue after evaporation		<=3 ppm
Solubility in H ₂ O		To pass test
Substances Reducing Permanganate		To pass test
Titratable Acid		<=0.0003mEq/g
Titratable Base		<=0.0002mEq/g
Water		<=0.1%

Applications: Used in staining and destaining protein electrophoresis gels, in HPLC, and in other biological applications.
Recommended storage: RT
Also available in recyclable FisherPak* and NOWPak* containers. EcoSafPak* is an environmentally friendly packaging system made of 100% recyclable material by an SFI certified supplier.

Sequencing Gel-Loading Dye, 3X
Contains 98% Formamide

packaging		Mfr. No
1 mL	PolyTube	BP639-1
DNase		Not detected
Optical Absorbance (dilution 1250 with deionized water) at 525nm		0.40-0.50
Optical Absorbance (dilution 1250 with deionized water) at 588-594nm		0.60-0.75
Optical Absorbance (dilution 1250 with deionized water) at 635-641nm		0.50-0.67
RNase		Not detected

Applications: This loading dye also denatures nucleic acid samples prior to electrophoresis on sequencing gels.
Components: 98% Formamide (<90%), Proprietary Component (<3%), and Water.
[75-12-7 (98% Formamide)]

TBE Buffer Mix		Electrophoresis
Dry Powder Mix of Tris/Boric Acid/EDTA		
Packaging		Mfr. No
86.1 g	PolyBottle	BP1396-86
CAS: 610769-35-2		
H360FD, H315, H319, H335		
P280, P305+P351+P338, P201, P308+P313		
DNase	Not detected	
Electrophoresis of DNA with 5X diluted buffer	Pass test	
pH (at 25°C)	8.3±0.2	
RNase	Not detected	
Solubility (in 1ℓ final volume)	Clear, colorless with no suspending impurities	

Recommended Storage: RT
Applications: Electrophoresis buffer.
Makes 1 liter of 5X buffer. Final concentration of components in 1X buffer: 0.089M Tris, 0.089M Boric Acid, 0.002M EDTA.
[77-86-1 (Tris)] ; [6381-92-6 (EDTA)] ; [10043-35-3 (Boric Acid)]

Tris-Acetate-EDTA 25X Powder		Electrophoresis
packaging		Mfr. No
1	FoilPack	BP1331-1
CAS: 135852-26-5		

Conductivity of a 1X Solution (at 25°C)	1300±200µmhos/cm
DNase	Not detected
Electrophoresis	To pass test
pH (1X solution) (at 25°C)	8.3±0.1
Protease	Not detected
RNase	Not detected
Solubility	To pass test

Applications: Tris-Acetate-EDTA (TAE) is commonly used as a buffer for nucleic acid electrophoresis.
Each pack contains preweighed powder to make 1ℓ of a 25X solution (1M Tris-Acetate and 0.025M EDTA).
Components: Tris (64.36%), Acetic Acid, Sodium Salt (31.91%) and Ethylenediamine Tetraacetic Acid (3.72%).
[77-86-1 (Tris)] ; [127-09-3 (Acetic Acid, Sodium Salt)] ; [60-00-4 (EDTA)]
Recommended Storage: RT

Tris-Borate-EDTA 10X Powder		Electrophoresis
Packaging		Mfr. No
1		BP1334-1
CAS: 610769-35-2		
H360FD, H315, H319, H335		
P280, P305+P351+P338, P201, P308+P313		
DNase	Not detected	
Electrophoresis	To pass test	
pH (1X solution) (at 25°C)	8.3±0.1	
Protease	Not detected	
RNase	Not detected	
Solubility (1X concentration 1g/60mℓ H ₂ O)	Clear, colorless with no suspending impurities	

Applications: Tris-Borate-EDTA (TBE) is commonly used as a buffer for nucleic acid electrophoresis.
Each pack contains preweighed powder to make 1ℓ of a 10X solution (0.89M Tris Base, 0.89M Boric Acid, and 0.02M EDTA).
Components: Tris (63.91%), Boric Acid (32.54%), and EDTA (3.55%).
[77-86-1 (Tris)] ; [10043-35-3 (Boric Acid)] ; [60-00-4 (EDTA)]
Recommended Storage: RT

Tris Buffer 2M Solution		Molecular Biology
packaging		Mfr. No
100 mL	PolyBottle	BP1759-100
500 mL	PolyBottle	BP1759-500
CAS: 77-86-1	P280, P305+P351+P338,	
EINECS: 201-064-4	P261, P302+P352, P264, P271	
H315, H319, H335		
Arsenic	<0.0005%	
Copper	<0.0001%	
DNase	Not detected	
Iron	<0.0001%	
Lead	<0.0001%	
Magnesium (Mg)	<0.0005%	
Optical Absorbance of a 1M Solution at 280nm	<0.05	
pH (at 25°C)	10.1-11.1	
Protease	Not detected	
RNase	Not detected	
Tris Concentration	2±0.05M	

Applications: Tris is a buffer component in molecular biology, tissue culture, and electrophoresis procedures.
Tris Solution, 2M, provides a convenient stock solution for preparing Tris buffers.
Recommended Storage: RT

Tris Buffer 0.3M Solution		Molecular Biology
packaging		Mfr. No
500 mL	PolyBottle	BP1761-500
1 ℓ	PolyBottle	BP1761-1
CAS: 77-86-1		
EINECS: 201-064-4		
DNase		Not detected
Molarity		0.3±0.010
pH (at 25°C)		10.1 to 11.1
Protease		Not detected
RNase		Not detected

Applications: Tris is a buffer component in molecular biology, tissue culture, and electrophoresis procedures.
Tris Solution, 0.3M, provides a convenient stock solution for these applications.
Recommended Storage: RT

Tris-EDTA 100X Powder		Molecular Biology	
packaging		Mfr. No	
1	FoilPack	BP1339-1	
CAS: 38641-82-6			
H315, H319, H335			
P280, P261, P305+P351+P338, P302+P352, P264, P271			
Conductivity of a 1X Solution (at 25°C)			560 to 740µmhos/cm
DNase			Not detected
pH (1X solution) (at 23°C)		8.0±0.1	
Protease		Not detected	
RNase		Not detected	
Solubility		To pass test	

Applications: Tris-EDTA (TE) is routinely used for suspending nucleic acid samples.
Each pack contains preweighed powder to make 1ℓ of a 100X solution (1.0M Tris Base and 0.1M EDTA).
[77-86-1 (Tris)] ; [60-00-4 (EDTA)]
Recommended Storage: RT

Tris Hydrochloride Solution pH 7.0 1M solution, High Purity, Low Metal		Molecular Biology
packaging		Mfr. No
100 mL	PolyBottle	BP1756-100
500 mL	PolyBottle	BP1756-500
CAS: 1185-53-1	P280, P305+P351+P338,	
EINECS: 214-684-5	P302+P352	
H315, H319		
Arsenic	<0.0005%	
Calcium (Ca)	<0.0002%	
Copper	<0.0001%	
DNase	Not detected	
Iron	<0.0001%	
Lead	<0.0001%	
Magnesium (Mg)	<0.0005%	
Molar Concentration	1±0.04M	
Protease	Not detected	
RNase	Not detected	
Zinc	<0.0001%	

Applications: Tris is a buffer component in molecular biology, tissue culture, and electrophoresis procedures.
[77-86-1 (Tris)] ; [1185-53-1 (Tris HCl)]
Recommended Storage: RT
Filtered through a 0.2-micron filter.

Tris Hydrochloride Solution pH 7.5 1M solution, High Purity, Low Metal		Molecular Biology
packaging		Mfr. No
100 mL	PolyBottle	BP1757-100
500 mL	PolyBottle	BP1757-500
CAS: 1185-53-1	P280, P305+P351+P338,	
EINECS: 214-684-5	P302+P352	
H315, H319		
Arsenic	<0.0005%	
Calcium (Ca)	<0.0002%	
Copper	<0.0001%	
DNase	Not detected	
Iron	<0.0001%	
Lead	<0.0001%	
Magnesium (Mg)	<0.0005%	
Molar Concentration	1±0.04M	
Protease	Not detected	
RNase	Not detected	
Zinc	<0.0001%	

Applications: Tris is a buffer component in molecular biology, tissue culture, and electrophoresis procedures.
[77-86-1 (Tris)] ; [1185-53-1 (Tris HCl)]
Recommended Storage: RT
Filtered through a 0.2-micron filter.

Tris Hydrochloride Solution pH 8.0 1M solution, High Purity, Low Metal		Molecular Biology
packaging	Mfr. No	
100 mℓ PolyBottle	BP1758-100	
500 mℓ PolyBottle	BP1758-500	
CAS: 1185-53-1	P280, P305+P351+P338,	
EINECS: 214-684-5	P302+P352	
H315, H319		
Arsenic	<0.0005%	
Calcium (Ca)	<0.0002%	
Copper	<0.0001%	
DNase	Not detected	
Iron	<0.0001%	
Lead	<0.0001%	
Magnesium (Mg)	<0.0005%	
Molar Concentration	1±0.04M	
Protease	Not detected	
RNase	Not detected	
Zinc	<0.0001%	

Applications: Tris is a buffer component in molecular biology, tissue culture, and electrophoresis procedures.
[77-86-1 (Tris)] ; [1185-53-1 (Tris HCl)]
Recommended Storage: RT
Filtered through a 0.2-micron filter.

Tris-Acetate-EDTA 10X Solution		Electrophoresis
packaging	Mfr. No	
500 mℓ PolyBottle	BP1335-500	
1 ℓ PolyBottle	BP1335-1	
4 ℓ PolyPac*	BP1335-4	
20 ℓ PolyPac*	BP1335-20	
CAS: 135852-26-5		
H315, H319		
P264, P280, P305+P351+P338,		
P302+P352		
DNase	Not detected	
Electrophoresis	To pass test	
pH (1X solution) (at 25°C)	8.2-8.4	
Protease	Not detected	
RNase	Not detected	

Applications: Tris-Acetate-EDTA (TAE) is commonly used as a buffer for nucleic acid electrophoresis.
0.4M Tris-Acetate and 0.01M EDTA.
Components: Tris (<5.0%), Acetic Acid (1.0%), and EDTA (<1%).
[77-86-1 (Tris)] ; [64-19-7 (Acetic Acid)] ; [60-00-4 (EDTA)]
Recommended Storage: RT
Filtered through a 0.2-micron filter.

Tris-Acetate-EDTA 25X Solution		Electrophoresis
packaging	Mfr. No	
1 ℓ PolyBottle	BP1330-1	
CAS: 135852-26-5		
H315, H319		
P264, P280, P305+P351+P338,		
P302+P352		
Conductivity of a 1X Solution (at 25°C)	Report	
DNase	Not detected	
Electrophoresis	To pass test	
pH (1X solution) (at 25°C)	8.3±0.1	
Protease	Not detected	
RNase	Not detected	

Applications: Tris-Acetate-EDTA (TAE) is commonly used as a buffer for nucleic acid electrophoresis.
1M Tris-Acetate and 0.025M EDTA.
Components: Tris (11.31%), EDTA Disodium Salt (0.87%), Acetic Acid, Sodium Salt (7.66%), and Hydrogen Chloride (3.74%).
[77-86-1 (Tris)] ; [6381-92-6 (EDTA Disodium Salt)] ; [127-09-3 (Acetic Acid, Sodium Salt)] ; [7647-01-0 (Hydrogen Chloride)]
Recommended Storage: RT
Filtered through a 0.2-micron filter.

Tris-Acetate-EDTA 50X Solution		Electrophoresis
packaging	Mfr. No	
500 mℓ PolyBottle	BP1332-500	
1 ℓ PolyBottle	BP1332-1	
4 ℓ PolyPac*	BP1332-4	
20 ℓ PolyPac*	BP1332-20	
CAS: 135852-26-5		
H315, H319		
P264, P280, P305+P351+P338,		
P302+P352		
Conductivity of a 1X Solution (at 25°C)	Report	
DNase	Not detected	
Electrophoresis	To pass test	
pH (1X solution) (at 25°C)	8.2-8.4	
Protease	Not detected	
RNase	Not detected	

Applications: Tris-Acetate-EDTA (TAE) is commonly used as a buffer for nucleic acid electrophoresis.
2M Tris-Acetate and 0.050M EDTA.
Components: Tris (24%), Acetic Acid (5.0%), and EDTA (<2%).
[77-86-1 (Tris)] ; [64-19-7 (Acetic Acid)] ; [60-00-4 (EDTA)]
Recommended Storage: RT
Filtered through a 0.2-micron filter.

Tris-Acetate-EDTA 1X Solution		Electrophoresis
packaging	Mfr. No	
4 ℓ PolyPac*	BP2434-4	
20 ℓ PolyPac*	BP2434-20	
CAS: 135852-26-5		
Conductivity of a 1X Solution (at 25°C)	Report	
DNase	Not detected	
Electrophoresis	To pass test	
pH (1X solution) (at 25°C)	8.3±0.1	
Protease	Not detected	
RNase	Not detected	

Applications: Tris-Acetate-EDTA (TAE) is commonly used as a buffer for nucleic acid electrophoresis.
0.04M Tris-Acetate and 0.001M EDTA.
[77-86-1 (Tris)] ; [64-19-7 (Acetic Acid)] ; [60-00-4 (EDTA)]
Recommended Storage: RT
Filtered through a 0.2-micron filter.

Tris-Borate-EDTA 10X Solution		Electrophoresis
packaging	Mfr. No	
1 ℓ PolyBottle	BP1333-1	
4 ℓ PolyPac*	BP1333-4	
20 ℓ PolyPac*	BP1333-20	
CAS: 610769-35-2		
H360FD, H315, H319		
P280, P305+P351+P338, P308+P313,		
P201		
Conductivity of 1X solution (at 25°C)	750-1250μmhos/cm	
DNase	Not detected	
Electrophoresis	To pass test	
pH (1X solution) (at 25°C)	8.3±0.1	
Protease	Not detected	
RNase	Not detected	

Applications: Tris-Borate-EDTA (TBE) is commonly used as a buffer for nucleic acid electrophoresis.
0.89M Tris Base, 0.89M Boric Acid, and 0.02M EDTA.
Components: Tris (10.8%), Boric Acid (5.5%), and EDTA (0.6%) [77-86-1 (Tris)] ; [10043-35-3 (Boric Acid)] ; [60-00-4 (EDTA)]
Recommended Storage: RT
Filtered through a 0.2-micron filter.

Tris-Borate-EDTA 1X Solution		Electrophoresis
packaging	Mfr. No	
1 ℓ PolyBottle	BP2430-1	
4 ℓ PolyPac*	BP2430-4	
20 ℓ PolyPac*	BP2430-20	
CAS: 610769-35-2		
DNase	Not detected	
Electrophoresis	To pass test	
pH of 1X solution (at 25°C)	8.3±0.1	
Protease	Not detected	
RNase	Not detected	

Applications: Tris-Borate-EDTA (TBE) is commonly used as a buffer for nucleic acid electrophoresis.
0.089M Tris Base, 0.089M Boric Acid, and 0.002M EDTA.
[77-86-1 (Tris)] ; [6381-92-6 (EDTA)] ; [10043-35-3 (Boric Acid)]
Recommended Storage: RT
Filtered through a 0.2-micron filter.

Tris-EDTA 100X Solution		Molecular Biology
packaging	Mfr. No	
1 ℓ PolyBottle	BP1338-1	
4 ℓ PolyPac*	BP1338-4	
CAS: 38641-82-6		
H315, H319		
P280, P305+P351+P338,		
P302+P352		
DNase	Not detected	
pH (1X solution) (at 23°C)	8.0±0.1	
Protease	Not detected	
RNase	Not detected	

Applications: Tris-EDTA (TE) is routinely used for suspending nucleic acid samples.
1.0M Tris Base and 0.1M EDTA.
[77-86-1 (Tris)] ; [60-00-4 (EDTA)]
Recommended Storage: RT
Filtered through a 0.2-micron filter.

Xylene Cyanole FF Bluish-green Powder		Electrophoresis
packaging	Mfr. No	
10 g AmberGlass	BP565-10	
C ₂₅ H ₂₇ N ₂ NaO ₆ S ₂	H315, H319, H335	
CAS: 2650-17-1	P280, P305+P351+P338	
MW: 538.61		
EINECS: 220-167-5		
Absorption Maxima	614±2nm	
DNase	Not detected	
Loss on Drying (at 105°C)	<=10.0%	
Molar Absorbance (Em, H2O)	>30.000	
RNase	Not detected	

Applications: This dye is used in electrophoretic applications to mark the migration position of nucleic acid fragments.
Recommended Storage: RT
C.I. 42135

100 Base-Pair Ladder	
packaging	Mfr. No
50 μg PolyMicroTube	BP2551-50

Applications: For precise sizing of double-stranded DNA fragments from 100-2000 bp on agarose gels.
Supplied at 1μg/μℓ in TE buffer.
Consists of multiple repeats of a 100 bp fragment.
Not designed for quantifying DNA content in a sample.
Can be visualized by ethidium bromide staining or by autoradiography (after radiolabeling).
Recommended Gel: 1.5% agarose with loading amount of 2.0μg/lane.
Do not heat before loading.
Recommended Storage: -20°C
Not on TSCA inventory: for R and D use only; not for manufacturing or commercial purposes.

Fisher BioReagents* *exACTGene**
50 bp Mini DNA Ladder


- Ready-to-use DNA size standards are stable at room temperature—no need to thaw and refreeze
- ▶ 25-650 bp range
 - ▶ **Number of Bands:** 14
 - ▶ Higher intensity reference band at 200 and 400 bp
 - ▶ Ideal for size determination and quantitative determination of DNA mass
 - ▶ Contain loading dye to reduce pipetting steps and save time
 - ▶ Provided in sufficient quantity to load 100 lanes

APPLICATIONS: PCR Size Confirmation

INCLUDES: 1mL of premixed ladder (0.5µg/10µL) in loading dye (10mM EDTA, 10% glycerol, 0.015% bromophenol blue, 0.17% SDS).

STORAGE CONDITIONS: Stable at room temperature for two years.

Quantity	Packaging	MFR No.
1000µL	Poly Micro Tube	BP2570-100



BP2570-100

Fisher BioReagents* *exACTGene**
Low Range DNA Ladder


- Ready-to-use DNA size standards are stable at room temperature—no need to thaw and refreeze
- ▶ 50-2000 bp range
 - ▶ **Number of Bands:** 8
 - ▶ Higher intensity reference band at 500 bp
 - ▶ Ideal for size assessment and quantitative determination of DNA mass
 - ▶ Contain loading dye to reduce pipetting steps and save time
 - ▶ Provided in sufficient quantity to load 100 lanes

APPLICATIONS: Quick Check of PCR Amplicons or Enzyme Digestion results

INCLUDES: 1mL of premixed ladder (0.5µg/10µL) in loading dye (10mM EDTA, 10% glycerol, 0.015% bromophenol blue, 0.17% SDS).

STORAGE CONDITIONS: Stable at room temperature for two years.

Quantity	Packaging	MFR No.
1000µL	Poly Micro Tube	BP2572-100



BP2572-100

Fisher BioReagents* *exACTGene**
100 bp DNA Ladder


- Ready-to-use DNA size standards are stable at room temperature—no need to thaw and refreeze
- ▶ 100-1000 bp range
 - ▶ **Number of Bands:** 10
 - ▶ Higher intensity reference band at 500 bp
 - ▶ Ideal for size assessment and quantitative determination of DNA mass
 - ▶ Contain loading dye to reduce pipetting steps and save time
 - ▶ Provided in sufficient quantity to load 100 lanes

APPLICATIONS: General Purpose and Small Size DNA

INCLUDES: 1mL of premixed ladder (0.5µg/10µL) in loading dye (10mM EDTA, 10% glycerol, 0.015% bromophenol blue, 0.17% SDS).

STORAGE CONDITIONS: Stable at room temperature for two years.

Quantity	Packaging	MFR No.
1000µL	Poly Micro Tube	BP2573-100



BP2573-100

Fisher BioReagents* *exACTGene**
Mid Range DNA Ladder


- Ready-to-use DNA size standards are stable at room temperature—no need to thaw and refreeze
- ▶ 300-5000 bp range
 - ▶ **Number of Bands:** 10
 - ▶ Ideal for size assessment and quantitative determination of DNA mass
 - ▶ Contain loading dye to reduce pipetting steps and save time
 - ▶ Provided in sufficient quantity to load 100 lanes

APPLICATIONS: Large Size PCR confirmation or cloning application

INCLUDES: 1mL of premixed ladder (0.5µg/10µL) in loading dye (10mM EDTA, 10% glycerol, 0.015% bromophenol blue, 0.17% SDS).

STORAGE CONDITIONS: Stable at room temperature for two years.

Quantity	Packaging	MFR No.
1000µL	Poly Micro Tube	BP2576-100



BP2576-100

Fisher BioReagents* *exACTGene**
100 bp PCR DNA Ladder


- Ready-to-use DNA size standards are stable at room temperature—no need to thaw and refreeze
- ▶ 25-1000 bp range
 - ▶ **Number of Bands:** 12
 - ▶ Higher intensity reference band at 500 bp
 - ▶ Ideal for size assessment and quantitative determination of DNA mass
 - ▶ Contain loading dye to reduce pipetting steps and save time
 - ▶ Provided in sufficient quantity to load 100 lanes

APPLICATIONS: PCR Size Confirmation and Small DNA Digests

INCLUDES: 1mL of premixed ladder (0.5µg/10µL) in loading dye (10mM EDTA, 10% glycerol, 0.015% bromophenol blue, 0.17% SDS).

STORAGE CONDITIONS: Stable at room temperature for two years.

Quantity	Packaging	MFR No.
1000µL	Poly Micro Tube	BP2571-100



BP2571-100

Fisher BioReagents* *exACTGene**
Low Range Plus DNA Ladder


- Ready-to-use DNA size standards are stable at room temperature—no need to thaw and refreeze
- ▶ 100-2000 bp range
 - ▶ **Number of Bands:** 11
 - ▶ Higher intensity reference band at 500 bp
 - ▶ Ideal for size assessment and quantitative determination of DNA mass
 - ▶ Contain loading dye to reduce pipetting steps and save time
 - ▶ Provided in sufficient quantity to load 100 lanes

APPLICATIONS: Fast Run Times, Small Size DNA

INCLUDES: 1mL of premixed ladder (0.5µg/10µL) in loading dye (10mM EDTA, 10% glycerol, 0.015% bromophenol blue, 0.17% SDS).

STORAGE CONDITIONS: Stable at room temperature for two years.

Quantity	Packaging	MFR No.
1000µL	Poly Micro Tube	BP2574-100



BP2574-100

Fisher BioReagents* *exACTGene**
Cloning DNA Ladder


- Ready-to-use DNA size standards are stable at room temperature—no need to thaw and refreeze
- ▶ 100-2686 bp range
 - ▶ **Number of Bands:** 14
 - ▶ Higher intensity reference band at 500 and 1000 bp
 - ▶ Ideal for size assessment and quantitative determination of DNA mass
 - ▶ A 2686 bp (pUC 19) reference band for clone identification
 - ▶ Contain loading dye to reduce pipetting steps and save time
 - ▶ Provided in sufficient quantity to load 100 lanes

APPLICATIONS: Clone Identification

INCLUDES: 1mL of premixed ladder (0.5µg/10µL) in loading dye (10mM EDTA, 10% glycerol, 0.015% bromophenol blue, 0.17% SDS).

STORAGE CONDITIONS: Stable at room temperature for two years.

Quantity	Packaging	MFR No.
1000µL	Poly Micro Tube	BP2575-100



BP2575-100

Fisher BioReagents* *exACTGene**
Mid Range Plus DNA Ladder


- Ready-to-use DNA size standards are stable at room temperature—no need to thaw and refreeze
- ▶ 100-5000 bp range
 - ▶ **Number of Bands:** 16
 - ▶ Higher intensity reference band at 500 and 1000 bp
 - ▶ Ideal for size assessment and quantitative determination of DNA mass
 - ▶ Contain loading dye to reduce pipetting steps and save time
 - ▶ Provided in sufficient quantity to load 100 lanes

APPLICATIONS: Small and Large Cloning Applications

INCLUDES: 1mL of premixed ladder (0.5µg/10µL) in loading dye (10mM EDTA, 10% glycerol, 0.015% bromophenol blue, 0.17% SDS).

STORAGE CONDITIONS: Stable at room temperature for two years.

Quantity	Packaging	MFR No.
1000µL	Poly Micro Tube	BP2577-100



BP2577-100

Fisher BioReagents* exACTGene* 1kb DNA Ladder

Ready-to-use DNA size standards are stable at room temperature—no need to thaw and refreeze


- ▶ 300-10,000 bp range
- ▶ **Number of Bands:** 13
- ▶ Higher intensity reference band at 5000 bp
- ▶ Ideal for size assessment and quantitative determination of DNA mass
- ▶ Contain loading dye to reduce pipetting steps and save time
- ▶ Provided in sufficient quantity to load 100 lanes

APPLICATIONS: General Purpose, Large Digested DNA

INCLUDES: 1mL of premixed ladder (0.5µg/10µL) in loading dye (10mM EDTA, 10% glycerol, 0.015% bromophenol blue, 0.17% SDS).

STORAGE CONDITIONS: Stable at room temperature for two years.

Quantity	Packaging	MFR No.
1000µL	Poly Micro Tube	BP2578-100



BP2578-100

Fisher BioReagents* exACTGene* 24kb Max DNA Ladder

Ready-to-use DNA size standards are stable at room temperature—no need to thaw and refreeze


- ▶ 300-24,000 bp range
- ▶ **Number of Bands:** 15
- ▶ Higher intensity reference band at 10,000 bp
- ▶ Ideal for wide range size assessment and quantitative determination of DNA mass
- ▶ Contain loading dye to reduce pipetting steps and save time
- ▶ Provided in sufficient quantity to load 100 lanes

APPLICATIONS: General Purpose, Extra Large Size DNA

INCLUDES: 1mL of premixed ladder (0.5µg/10µL) in loading dye (10mM EDTA, 10% glycerol, 0.015% bromophenol blue, 0.17% SDS).

STORAGE CONDITIONS: Stable at room temperature for two years.

Quantity	Packaging	MFR No.
1000µL	Poly Micro Tube	BP2580-100



BP2580-100

Fisher BioReagents* exACTGene* 1kb Plus DNA Ladder

Ready-to-use DNA size standards are stable at room temperature—no need to thaw and refreeze


- ▶ 100-10,000 bp range
- ▶ **Number of Bands:** 19
- ▶ Higher intensity reference band at 500, 1000 and 5000 bp
- ▶ Ideal for wide range size assessment and quantitative determination of DNA mass
- ▶ Contain loading dye to reduce pipetting steps and save time
- ▶ Provided in sufficient quantity to load 100 lanes

APPLICATIONS: General Purpose, Wide DNA Size Range

INCLUDES: 1mL of premixed ladder (0.5µg/10µL) in loading dye (10mM EDTA, 10% glycerol, 0.015% bromophenol blue, 0.17% SDS).

STORAGE CONDITIONS: Stable at room temperature for two years.

Quantity	Packaging	MFR No.
1000µL	Poly Micro Tube	BP2579-100




BP2579-100

Fisher BioReagents* Routine DNA Ladders


Ready-to-use DNA size standards are room-temperature stable for up to two years

- ▶ Both ladders include a landmark doublet at 1400 and 1550 bp

Marker (Base Pair Range)	No. of Bands	Application	Size	MFR No.
100 bp Low Scale DNA Ladder (50 to 2000 bp)	11	PCR Size Confirmation	1000µL (200 loadings)	BP2581-200
1kb Full Scale DNA Ladder (50 to 10,000 bp)	16	General Purpose, Wide DNA Size Range	1000µL (200 loadings)	BP2582-200



BP2582-200



BP2581-200

250 Base-Pair Ladder

packaging	Mfr. No
200 µg PolyMicroTube	BP2552-100

Applications: For precise sizing of double-stranded DNA fragments from 250-2000 bp on agarose gels. Supplied at 1µg/µL in TE buffer (pH 7.5). Consists of multiple repeats of a 250 bp fragment. Not designed for quantifying DNA content in a sample. **Recommended Gel:** 1% agarose with loading amount of 2.0µg/lane. Do not heat before loading. **Recommended Storage:** -20°C **Not on TSCA inventory:** for R and D use only; not for manufacturing or commercial purposes.

50 Base-Pair Ladder

packaging	Mfr. No
200 µg PolyMicroTube	BP2550-100

Applications: For precise sizing of average-sized PCR products and other DNAs from 50-500 bp on agarose gels. Supplied at 1µg/µL in TE buffer (pH 7.5). Consists of multiple repeats of a 50 bp fragment. Can be visualized by ethidium bromide staining or by autoradiography (after radiolabeling). **Recommended Gel:** 2% agarose/ 6% polyacrylamide with loading amount of 2.0µg/lane. Do not heat before loading. **Recommended Storage:** -20°C **Not on TSCA inventory:** for R and D use only; not for manufacturing or commercial purposes.

Kilobase DNA Marker

packaging	Mfr. No
50 µg PolyMicroTube	BP2553-100

Applications: For quick and precise sizing of double-stranded DNA fragments from 500 bp to 10.0kb on agarose gels. Supplied at 0.5µg/µL in TE buffer (pH 8.0). Consists of 11 easily identifiable bands - 0.5, 1.0, 1.5, 2.0, 2.5, 3.0, 4.0, 5.0, 6.0, 8.0, and 10.0kb - when run on a 0.8% agarose gel. Not designed for quantifying DNA content in a sample. **Recommended Gel:** 0.8% agarose with loading amount of 0.5µg/lane. Do not heat before loading. **Recommended Storage:** -20°C **Not on TSCA inventory:** for R and D use only; not for manufacturing or commercial purposes.

DNA-Hind III Digest

packaging	Mfr. No
100 µg PolyMicroTube	BP2556-200

Applications: For sizing DNA fragments from 125 to 23.130 bp on agarose gels. **Recommended Gel:** 1% agarose with loading amount of 0.5µg/lane. Contains lambda DNA that has been digested with restriction enzyme. **Recommended Storage:** -20°C **Not on TSCA inventory:** for R and D use only; not for manufacturing or commercial purposes.

DNA-Hind III/phi X-174 RF DNA-Hae III Digest

packaging	Mfr. No
25 µg PolyMicroTube	BP2555-25

Applications: For sizing DNA fragments from 72 to 23.130 bp on agarose gels. Contains a combination of lambda DNA and phi X-174 RF DNA that has been digested with restriction enzyme. **Recommended Gel:** 1% agarose with loading amount of 0.5-1.0µg/lane. Heat before loading at 60°-65°C for 2 minutes. **Recommended Storage:** -20°C **Not on TSCA inventory:** for R and D use only; not for manufacturing or commercial purposes.

DNA-Hind III/phi X-174 RF DNA-Hae III Digest Lyophilized Powder

packaging	Mfr. No
25 µg PolyMicroTube	BP2554-25

Applications: For sizing DNA fragments from 72 to 23.130 bp on agarose gels. Contains a combination of lambda DNA and phi X-174 RF DNA that has been digested with restriction enzyme. **Recommended Gel:** 1% agarose with loading amount of 0.5-1.0µg/lane. Heat before loading at 60°-65°C for 2 minutes. **Recommended Storage:** -20°C **Not on TSCA inventory:** for R and D use only; not for manufacturing or commercial purposes.

Ø X-174 RF DNA-Hae III Digest

packaging	Mfr. No
10 µg PolyMicroTube	BP2558-10
100 µg PolyMicroTube	BP2558-100

Applications: For sizing DNA fragments from 72 to 1.353 bp on agarose gels. **Recommended Gel:** 1% agarose with loading amount of 0.5µg/lane. Do not heat before loading. Contains phi X-174 RF DNA that has been digested with restriction enzyme. **Recommended Storage:** -20°C **Not on TSCA inventory:** for R and D use only; not for manufacturing or commercial purposes.

Collagenase From Clostridium histolyticum Islet Isolation Grade Yellowish Brown Lyophilized Powder

packaging	Mfr. No
1 g AmberGlass	BP2649-1

CAS: 9001-12-1

Activity >=3.500units/g

Applications: Collagenase is most suitable for pancreatic islet isolation. **Unit Definition:** One unit will liberate peptides from collagen equivalent in ninhydrin color to 3.0µmole of L-leucine in 18 hr. at pH 7.4 and 37°C. Inhibitors: Hg²⁺, Pb²⁺, Cd²⁺, Cu²⁺, Zn²⁺, EDTA, o-phenanthroline. **Recommended Storage:** 0° to 5°C, desiccate



Cytosine-β-D-Arabinofuranoside Hydrochloride
White Powder

packaging	Mfr. No
100 mg AmberGlass	BP2512-100
500 mg AmberGlass	BP2512-500
5 g AmberGlass	BP2512-5
10 g AmberGlass	BP2512-10

C₉H₁₃N₅O₅·HCl

CAS: 69-74-9

MW: 279.67

H317, H319, H361d

P280, P281,

P305+P351+P338

FTIR

Conforms to standard

Melting Point

185°-195°C

Optical Rotation α_D²⁵

+130° ± 5° (c=1, water)

Tested for

RNase activity, protease activity, and specific performance tests.

UV/VIS: lambda max (H₂O)

280 ± 6nm

Applications: Cytosine-β-D-Arabinofuranoside Hydrochloride inhibits incorporation of labeled thymidine into DNA.
Recommended Storage: 4°C
Not on TSCA inventory: for R and D use only; not for manufacturing or commercial purposes.

Lysozyme, Egg White
White Crystalline Powder

packaging	Mfr. No
1 g AmberGlass	BP535-1
5 g AmberGlass	BP535-5
10 g AmberGlass	BP535-10
25 g AmberGlass	BP535-25
100 g AmberGlass	BP535-100

CAS: 12650-88-3

EINECS: 235-747-3

Activity	Approx. 20.000 units/mg protein
Salmonella	None detected
Solubility	Pass test
Tested for	RNase activity, protease activity, and specific performance tests.

Applications: This enzyme is used to lyse E. coli for preparation of plasmid DNA from transformed cell cultures.
Recommended Storage: <0° C
Not on TSCA inventory: for R and D use only; not for manufacturing or commercial purposes.
UN 1198

OPTIZYME™ Alkaline Phosphatase
Source: Recombinant E. Coli Strain

packaging	Mfr. No
1000 units Poly Tube	BP8097-1
3000 units Poly Tube	BP8097-5

Tested for

RNase activity, protease activity, and specific performance tests.

Applications: Dephosphorylation of cloning vector DNA to prevent recircularization during ligation; PCR product clean-up: nucleotide degradation prior to sequencing of PCR product; Dephosphorylation of nucleic acid 5'-termini prior to labeling with T4 Polynucleotide Kinase; Other applications where dephosphorylation of DNA and RNA substrates is necessary; Protein dephosphorylation.
Description: OPTIZYME™ Alkaline Phosphatase catalyzes the release of 5'- and 3'- phosphate groups from DNA, RNA and both ribo- and deoxyribonucleoside triphosphates. OPTIZYME Alkaline Phosphatase also removes phosphate groups from proteins.
Supplied With: 10X OPTIZYME™ AP Buffer (100mM Tris-HCl (pH 8.0 at 37°C), 50mM MgCl2, 1M KCl, 0.2% (v/v) Triton X-100, 1 mg/ml BSA)
Concentration: 1 u/μl
Storage Buffer Components: 20 mM HEPES-NaOH (pH 7.4), 1 mM MgCl2, 0.1 mM ZnCl2, 0.1% (v/v) Triton X-100 and 50% (v/v) glycerol.
Unit Definition: One unit is the amount of the enzyme required to dephosphorylate 5'-termini of 1 μg of linearized pUC57 DNA in 10 min at 37°C in OPTIZYME AP Buffer.
Tested for: Endo-, exodeoxyribonuclease and ribonuclease-free; Functionally tested for dephosphorylation of 5'-termini of overhanging, recessed and blunt DNA.
Recommended storage: -20°C

OPTIZYME DNA Polymerase I
Source: Recombinant E. Coli Strain

packaging	Mfr. No
500 units Poly Tube	BP8109-1
2500 units Poly Tube	BP8109-5

Tested for

RNase activity, protease activity, and specific performance tests.

Applications: DNA labeling by nick-translation in conjunction with DNase I; Second-strand cDNA synthesis in conjunction with RNase H; Filling-in of 5'-overhangs
Description: OPTIZYME™ DNA Polymerase I, a template-dependent DNA polymerase, catalyzes 5'>3' synthesis of DNA. The enzyme also exhibits 3'>5' exonuclease (proofreading) activity, 5'>3' exonuclease activity and ribonuclease H activity.
Supplied With: 10X OPTIZYME™ DNA Pol I Buffer (500mM Tris-HCl (pH 7.5 at 25°C), 100mM MgCl2, 10mM DTT).
Concentration: 10u/μl
Storage Buffer Components: 25 mM Tris-HCl (pH 7.5), 0.1 mM EDTA, 1 mM DTT and 50% (v/v) glycerol.
Unit Definition: One unit of the enzyme catalyzes the incorporation of 10 nmol of deoxyribonucleotides into a polynucleotide fraction (adsorbed on DE-81) in 30 min at 37°C, using poly(dA-dT)-poly(dA-dT) as a template-primer. Enzyme activity is assayed in the following mixture: 67 mM potassium phosphate (pH 7.4), 6.7 mM MgCl2, 1 mM 2-mercaptoethanol, 0.033 mM dATP, 0.033 mM dTTP, 0.4 MBq/ml -dTTP and 62.5 μg/ml poly(dA-dT) poly(dA-dT).
Tested for: Endodeoxyribonuclease-free
Recommended storage: -20°C

OPTIZYME™ DNase I, RNase-Free
Source: Recombinant E. Coli Strain

packaging	Mfr. No
1000 units Poly Tube	BP8107-1

Tested for

RNase activity, protease activity, and specific performance tests.

Applications: Preparation of DNA-free RNA; Removal of template DNA following in vitro transcription; Preparation of DNA-free RNA prior to RT-PCR and RT-qPCR; DNA labeling by nick-translation in conjunction with DNA Polymerase I; Studies of DNA-protein interactions by DNase footprinting.
Description: OPTIZYME™ DNase I, RNase-Free is an endonuclease that digests single- and double-stranded DNA. It hydrolyzes phosphodiester bonds producing mono- and oligodeoxyribonucleotides with phosphate and OH groups. The enzyme activity is strictly dependent on Ca2+ and is activated by Mg2+ or Mn2+ ions: (1) In the presence of Mg2+, DNase I cleaves each strand of dsDNA independently, in a statistically random fashion; (2) In the presence of Mn2+, the enzyme cleaves both DNA strands at approximately the same site, producing DNA fragments with blunt ends or with one or two nucleotide overhangs.
Supplied With: 10X OPTIZYME™ DNase I Buffer (100mM Tris-HCl (pH 7.5 at 25°C), 25mM MgCl2, 1mM CaCl2); 50mM OPTIZYME™ EDTA
Concentration: 1u/μl
Storage Buffer Components: 50 mM Tris-HCl (pH 7.5), 10 mM CaCl2 and 50% (v/v) glycerol.
Unit Definition: One unit of the enzyme completely degrades 1 μg of plasmid DNA in 10 min at 37°C. Enzyme activity is assayed in the following mixture: 10 mM Tris-HCl (pH 7.5 at 25°C), 2.5 mM MgCl2, 0.1 mM CaCl2, 1 μg of pUC19 DNA. One DNase I, RNase-Free unit is equivalent to 0.3 Kunitz unit.
Tested for: Ribonuclease-free; Functionally tested for the digestion of template DNA after in vitro transcription.
Recommended storage: -20°C

OPTIZYME Exonuclease III
Source: Recombinant E. Coli Strain

packaging	Mfr. No
4000 units Poly Tube	BP8108-1

Tested for

RNase activity, protease activity, and specific performance tests.

Applications: Creation of unidirectional deletions in DNA fragments in conjunction with S1 Nuclease; Generation of a single-stranded template for dideoxy sequencing of DNA; Site-directed mutagenesis; Cloning of PCR products; Preparation of strand-specific probes.
Description: OPTIZYME™ Exonuclease III exhibits four different catalytic activities: (1) 3'>5' exodeoxyribonuclease activity specific for double-stranded DNA. Exo III degrades dsDNA from blunt ends, 5'-overhangs or nicks, releasing 5'-mononucleotides from the 3'-ends of DNA strands and producing stretches of single stranded DNA. It is not active on 3'-overhang ends of DNA that are at least four bases long and do not carry a 3'-terminal C-residue (on single-stranded DNA, or on phosphorothioate-linked nucleotides); (2) 3'-phosphatase activity: Exo III removes the 3'-terminal phosphate and generates a 3'-OH group; (3) RNase H activity: Exo III exonucleolytically degrades the RNA strand in DNA-RNA hybrids; (4) Apurinic/apyrimidinic-endonuclease activity: Exo III cleaves phosphodiester bonds at apurinic or apyrimidinic sites to produce 5'-termini that are base free deoxyribose 5'-phosphate residues.
Supplied With: 10X OPTIZYME™ ExoIII Buffer (660mM Tris-HCl (pH 8.0 at 30°C), 6.6 mM MgCl2).
Concentration: 200u/μl
Storage Buffer Components: 50 mM Tris-HCl (pH 8.0), 50 mM KCl, 1 mM DTT and 50% (v/v) glycerol.
Unit Definition: One unit of the enzyme catalyzes the release of 1 nmol of acid soluble reaction products from E.coli -DNA in 30 min at 37°C. Enzyme activity is assayed in the following mixture: 50 mM Tris-HCl (pH 8.0), 5 mM MgCl2, 1 mM DTT and 0.05 mM sonicated E.coli -DNA.
Tested for: Endodeoxyribonuclease-free; Functionally tested for the creation of unidirectional deletions in DNA fragments.
Recommended storage: -20°C

OPTIZYME™ Klenow Fragment
Source: Recombinant E. Coli Strain

packaging	Mfr. No
300 units Poly Tube	BP8106-1
1500 units Poly Tube	BP8106-5

Tested for

RNase activity, protease activity, and specific performance tests.

Applications: DNA blunting by fill-in of 5'-overhangs; Random-primed DNA labeling; Labeling by fill-in 5'-overhangs of dsDNA; DNA sequencing by the Sanger method; Site-specific mutagenesis of DNA with synthetic oligonucleotides; Second strand synthesis of cDNA.
Description: OPTIZYME™ Klenow Fragment is the large fragment of DNA Polymerase I (E.coli). It exhibits 5'>3' polymerase activity and 3'> 5' exonuclease (proofreading) activity, but lacks the 5'> 3' exonuclease activity of DNA Polymerase I.
Supplied With: 10X OPTIZYME™ Klenow Fragment Buffer (500mM Tris-HCl (pH 8.0 at 25°C), 50mM MgCl2, 10mM DTT).
Concentration: 10u/μl
Storage Buffer Components: 25 mM Tris-HCl (pH 7.5), 0.1 mM EDTA, 1 mM DTT and 50% (v/v) glycerol.
Unit Definition: One unit of the enzyme catalyzes the incorporation of 10 nmol of deoxyribonucleotides into a polynucleotide fraction (adsorbed on DE-81) in 30 min at 37°C. Enzyme activity is assayed in the following mixture: 50 mM Tris-HCl (pH 8.0 at 25°C), 5 mM MgCl2, 1 mM DTT, 0.033 mM dNTP, 0.4 M Bq/ml -dTTP and 62.5 μg/ml activated salmon milt DNA.
Tested for: Endodeoxyribonuclease-free; Functionally tested for fill in of 5'-overhanging DNA termini and for random primed DNA labeling
Recommended storage: -20°C

OPTIZYME™ M-MLV Reverse Transcriptase
Source: Recombinant E. Coli Strain

packaging	Mfr. No
5000 units Poly Tube	BP8104-1
25000 units Poly Tube	BP8104-5

Tested for

RNase activity, protease activity, and specific performance tests.

Applications: First strand cDNA synthesis for RT-PCR and real-time RT-PCR; Synthesis of cDNA for cloning and expression; Generation of labeled cDNA probes for microarrays; DNA labeling; Analysis of RNA by primer extension.
Description: OPTIZYME™ M-MLV Reverse Transcriptase possesses RNA-dependent and DNA-dependent polymerase activity and RNase H activity specific to RNA in RNA-DNA hybrids, which is significantly lower than that of Avian Myeloblastis Virus (AMV) reverse transcriptase. M-MLV RT activity is optimal at 42°C (active up to 50°C). The enzyme is capable of both first strand cDNA synthesis (<=13 kb) and incorporation of modified nucleotides.
Supplied With: 5X OPTIZYME™ M-MLV RT Buffer (250mM Tris-HCl (pH 8.3 at 25°C), 250mM KCl, 20mM MgCl2, 50mM DTT).
Concentration: 200 u/μl
Storage Buffer Components: 50 mM Tris-HCl (pH 7.5), 0.1 M NaCl, 1 mM EDTA, 5 mM DTT, 0.1% (v/v) Triton X-100 and 50% (v/v) glycerol.
Unit Definition: One unit of the enzyme incorporates 1 nmol of dTMP into a polynucleotide fraction (adsorbed on DE-81) in 10 min at 37°C. Enzyme activity is assayed in the following mixture: 50 mM Tris-HCl (pH 8.3), 4 mM MgCl2, 10 mM DTT, 50 mM KCl, 0.5 mM dTTP, 0.4 MBq/ml -dTTP, 0.4 mM polyA-oligo (dT)12-18.
Tested for: Endo-, exodeoxyribonuclease, phosphatase and ribonuclease-free; Functionally tested in first strand cDNA and RT-PCR.
Recommended storage: -20°C



rDNase I (RNase-free) Source: Recombinant E. Coli Stain		OPTIZYME™
packaging		Mfr. No
1000 units PolyTube		BP3226-1
2000 units PolyTube		BP3226-2

Tested for RNase activity, protease activity, and specific performance tests.

Applications: Robust ability to remove genomic and template DNA from RNA preparations; meets or exceeds the performance of native bovine DNase I. Developed by recombinant methods; non-animal source (yeast) eliminates potential RNase contamination-safest choice to address environmental concerns.
Certified RNase free. Highly purified for excellent lot-to-lot consistency and optimal performance. Can even be used in the characterization of DNA:Protein interactions.
Ideal for use in RT-PCR, cDNA synthesis, and in vitro RNA transcription.
Description: Recombinant DNase I (rDNase I) efficiently hydrolyzes phosphodiester linkages of double- or single-stranded DNA. Requires magnesium and calcium for maximum activity.
Concentration: 2 units/ μ l
Storage Buffer Components: 20mM HEPES (pH 7.5), 10mM CaCl₂, 10mM MgCl₂, 1mM DTT, 50% (v/v) Glycerol.
Provided with 10X Reaction Buffer: 100mM Tris-HCl (pH 7.5), 25mM MgCl₂, 5mM CaCl₂.
Unit Definition: One unit is defined as the amount of enzyme required to completely degrade 1 μ g DNA in 10 minutes at 37°C.
Recommended Storage: -20°C

OPTIZYME™ recombinant DNase I		10X Reaction Buffer
packaging		Mfr. No
1 mL PolyTube		BP3227-1

Applications: Provides optimal pH and ionic conditions for use with Optizyme rDNase I.
10X Reaction Buffer Components:100mM Tris-HCl (pH 7.5), 25mM MgCl₂, 5mM CaCl₂.
Recommended Storage: -20°C

Ribonuclease Inhibitor Source: Recombinant		OPTIZYME™
packaging		Mfr. No
2500 units PolyTube		BP3222-1
10000 units PolyTube		BP3222-5

Tested for DNase and Nickase contamination, absence of RNase activity, and specific performance tests

Applications: Recombinant Ribonuclease Inhibitor displays a broad spectrum of inhibitory activity against RNases, and does not have activity against other polymerases and reverse transcriptases. Suitable for use in common molecular biology applications such as isolation and purification of RNA, cDNA synthesis, RT-PCR, in vitro RNA transcription/translation, ribonuclease protection assay, and preparation of RNase-free antibodies.
Description: Supplied as a recombinant product from E. coli (originally isolated from rat lung) and provides superior protection of RNA from degradation by RNases. This RNase inhibitor is commonly added to all solutions used during the isolation of RNA, and will not interfere with the performance of most enzymes in downstream applications. Active over a broad temperature range, and even provides some RNase inhibition at 60°C which is useful when performing reverse transcription reactions at elevated temperatures to overcome secondary structure in RNA.
Storage Buffer Components: 20mM HEPES-KOH (pH 7.6 at 4°C), 50 mM KCl, 8 mM DTT, and 50% glycerol (v/v).
Note: Inactivates RNase A, RNase B, and RNase T2. No activity against SP6, T7, or T3 RNA polymerases, AMV or MMLV reverse transcriptases, and Taq DNA polymerase.
Unit Activity: Defined as the amount of inhibitor required to inhibit 50% of the activity of 5ng of RNase A.
Recommended Storage: -20°C

OPTIZYME™ SP6 RNA Polymerase Source: Recombinant E. Coli Strain	
packaging	Mfr. No
2000 units Poly Tube	BP8100-1
10000 units Poly Tube	BP8100-5

Tested for RNase activity, protease activity, and specific performance tests.

Applications: Synthesis of unlabeled & labeled RNA that can be used: (1) for hybridization and for in vitro RNA translation, (2) as RNA or siRNA, (3) as a substrate in RNase protection assays or as template for genomic DNA sequencing and (4) in studies of RNA secondary structure and RNA-protein interactions.
Description: OPTIZYME™ SP6 RNA Polymerase is a DNA-dependent RNA polymerase with strict specificity for its respective double-stranded promoter. The enzyme catalyzes the 5'->3' synthesis of RNA on either single-stranded DNA or double-stranded DNA downstream from the SP6 promoter. SP6 RNA Polymerase accepts modified nucleotides (e.g., biotin-, digoxigenin-, fluorescein-labeled nucleotides) as substrates for RNA synthesis.
Supplied With: 5X OPTIZYME™ Transcription Buffer (200mM Tris-HCl (pH 7.9 at 25°C), 30mM MgCl₂, 50mM DTT, 50mM NaCl, 10mM spermidine).
Concentration: 20u/ μ l
Storage Buffer Components: 50 mM Tris-HCl (pH 8.0), 150 mM NaCl, 5 mM DTT, 0.1 mg/mL BSA, 0.03% (v/v) ELUGENT Detergent and 50% (v/v) glycerol.
Unit Definition: One unit of the enzyme incorporates 1 nmol of AMP into a polynucleotide fraction (adsorbed on DE-81) in 60 minutes at 37°C. Enzyme activity is assayed in the following mixture: 40 mM Tris-HCl (pH 8.0), 6 mM MgCl₂, 10 mM DTT, 2 mM spermidine, 0.5 mM of each NTP, 0.6 MBq/mL -ATP, 20 μ g/mL plasmid DNA containing the SP6 promoter sequences.
Tested for: Endo-, exodeoxyribonuclease and ribonuclease-free; Functionally tested in in vitro transcription reaction.
Recommended storage: -20°C

OPTIZYME™ T3 RNA Polymerase Source: Recombinant E. Coli Strain	
packaging	Mfr. No
1000 units Poly Tube	BP8101-1
5000 units Poly Tube	BP8101-5

Tested for RNase activity, protease activity, and specific performance tests.

Applications: Synthesis of unlabeled & labeled RNA that can be used: (1) for hybridization and for in vitro RNA translation, (2) as RNA or siRNA, (3) as a substrate in RNase protection assays or as template for genomic DNA sequencing and (4) in studies of RNA secondary structure and RNA-protein interactions.
Description: OPTIZYME™ T3 RNA Polymerase is a DNA-dependent RNA polymerase with strict specificity for its respective double-stranded promoter. The enzyme catalyzes the 5'->3' synthesis of RNA on either single-stranded DNA or double-stranded DNA downstream from the T3 promoter. Bacteriophage T3 RNA Polymerase accepts modified nucleotides (e.g. biotin-, digoxigenin-, fluorescein-labeled nucleotides) as substrates for RNA synthesis.
Supplied With: 5X OPTIZYME™ Transcription Buffer (200mM Tris-HCl (pH 7.9 at 25°C), 30mM MgCl₂, 50mM DTT, 50mM NaCl, 10mM spermidine).
Concentration: 20u/ μ l
Storage Buffer Components: 50 mM Tris-HCl (pH 8.0), 150 mM NaCl, 5 mM DTT, 0.1 mg/mL BSA, 0.03% (v/v) ELUGENT Detergent and 50% (v/v) glycerol.
Unit Definition: One unit of the enzyme incorporates 1 nmol of AMP into a polynucleotide fraction (adsorbed on DE-81) in 60 minutes at 37°C. Enzyme activity is assayed in the following mixture: 40 mM Tris-HCl (pH 8.0), 6 mM MgCl₂, 10 mM DTT, 2 mM spermidine, 0.5 mM of each NTP, 0.6 MBq/mL -ATP, 20 μ g/mL plasmid DNA containing the T3 promoter sequence.
Tested for: Endo-, exodeoxyribonuclease and ribonuclease-free; Functionally tested in in vitro transcription reaction.
Recommended storage: -20°C

OPTIZYME™ T4 DNA Ligase Source: Recombinant E. Coli Strain	
packaging	Mfr. No
100 units Poly Tube	BP8099-1
1500 units Poly Tube	BP8099-15
500 units Poly Tube	BP8099-5

Tested for RNase activity, protease activity, and specific performance tests.

Applications: Molecular cloning; Joining of double-stranded oligonucleotide linkers or adaptors to DNA; Site-directed mutagenesis; Amplified fragment length polymorphism (AFLP); Ligase-mediated RNA detection; Nick repair in duplex DNA, RNA or DNA/RNA hybrids.
Description: OPTIZYME™ T4 DNA Ligase catalyzes the formation of a phosphodiester bond between juxtaposed 5'-phosphate and 3'-hydroxyl termini in duplex DNA or RNA. The enzyme repairs single-strand nicks in duplex DNA, RNA or DNA/RNA hybrids and joins DNA fragments with either cohesive or blunt termini. The T4 DNA Ligase requires ATP as a cofactor.
Supplied With: 10X OPTIZYME™ T4 DNA Ligase Buffer (400mM Tris-HCl, 100mM MgCl₂, 100mM DTT, 5mM ATP (pH 7.8 at 25°C); 50% OPTIZYME™ PEG 4000 Solution (50% (w/v) polyethylene glycol).
Concentration: 5u/ μ l
Storage Buffer Components: 20 mM Tris-HCl (pH 7.5), 50 mM KCl, 1 mM DTT, 0.1 mM EDTA and 50% (v/v) glycerol.
Unit Definition: One Weiss unit of the enzyme catalyzes the conversion of 1 nmol of into Norit-adsorbable form in 20 min at 37°C. One Weiss unit is equivalent to approximately 200 cohesive end ligation units (CEU)*. Enzyme activity is assayed in the following mixture: 66 mM Tris-HCl (pH 7.6), 6.6 mM MgCl₂, 0.066 mM ATP, 10 mM DTT, 3.3 μ M . (*One CEU is defined as the amount of enzyme required to give 50% ligation of HindIII fragments of lambda DNA in 30 min at 16°C).
Tested for: Endo-, exodeoxyribonuclease, phosphatase and ribonuclease-free; Functionally tested for ligation of cohesive and blunt end DNA fragments.
Recommended storage: -20°C

OPTIZYME™ T4 DNA Polymerase Source: Recombinant E. Coli Strain	
packaging	Mfr. No
100 units Poly Tube	BP8105-1
500 units Poly Tube	BP8105-5

Tested for RNase activity, protease activity, and specific performance tests.

Applications: Blunting of DNA ends: fill-in 5'-overhangs and/or removal of 3'-overhangs; Blunting of PCR products with 3'-dA overhangs; Synthesis of labeled DNA probes by the replacement reaction; Oligonucleotide-directed site-specific mutagenesis; Ligation-independent cloning of PCR products.
Description: OPTIZYME™ T4 DNA Polymerase is a template -dependent DNA polymerase which catalyzes 5'->3' synthesis from primed single-stranded DNA. The enzyme possesses 3'->5' exonuclease activity, but lacks 5'->3' exonuclease activity.
Supplied With: 5X OPTIZYME™ T4 DNA Pol Buffer (335mM Tris-HCl (pH 8.8 at 25°C), 33mM MgCl₂, 5mM DTT, 84mM (NH₄)₂SO₄)
Concentration: 5 u/ μ l
Storage Buffer Components: 20 mM potassium phosphate (pH 7.5), 200 mM KCl, 2 mM DTT and 50% (v/v) glycerol.
Unit Definition: One unit of the enzyme catalyzes the incorporation of 10 nmol of deoxyribonucleotides into a polynucleotide fraction (adsorbed on DE-81) in 30 min at 37°C. Enzyme activity is assayed in the following mixture: 67 mM Tris-HCl (pH 8.8), 6.7 mM MgCl₂, 1 mM DTT, 16.7 mM (NH₄)₂SO₄, 0.2 mg/mL BSA, 0.033 mM of each dNTP, 0.4 MBq/mL -dTTP and 0.2 mM heat-denatured and nuclease-digested calf thymus DNA.
Tested for: Endodeoxyribonuclease-free
Recommended storage: -20°C

OPTIZYME™ T4 Polynucleotide Kinase Source: Recombinant E. Coli Strain	
packaging	Mfr. No
500 units Poly Tube	BP8098-1
2500 units Poly Tube	BP8098-5

Tested for RNase activity, protease activity, and specific performance tests.

Applications: Labeling 5'-termini of nucleic acids to be used as: probes for hybridization, probes for transcript mapping markers for gel-electrophoresis primers for DNA sequencing, and primers for PCR; 5'-phosphorylation of oligonucleotides, PCR products, and other DNA or RNA prior to ligation; Phosphorylation of PCR primers; Detection of DNA modification by the post-labeling assay; Removal of 3'-phosphate groups.
Description: OPTIZYME™ T4 PNK catalyzes the transfer of the Upsilon gamma-phosphate from ATP to the 5'-OH group of single- and double-stranded DNAs and RNAs, oligonucleotides or nucleoside 3'-monophosphates (forward reaction). The reaction is reversible. In the presence of ADP, T4 PNK exhibits 5'-phosphatase activity and catalyzes the exchange of phosphate groups between 5'-P-oligo-/polynucleotides and ATP (exchange reaction). The enzyme is also a 3'-phosphatase.
Supplied With: 10X OPTIZYME T4 PNK Buffer 1 (500 mM Tris-HCl (pH 7.6 at 25°C), 100 mM MgCl₂, 50 mM DTT, 1 mM spermidine); 10X OPTIZYME T4 PNK Buffer 2 (2500 mM imidazole-HCl (pH 6.4 at 25°C), 180 mM MgCl₂, 50 mM DTT, 1 mM spermidine and 1 mM ADP); 24% PEG 6000 Sol for T4 PNK (24% (w/v) polyethylene glycol 6000)
Concentration: 10 u/ μ l
Storage Buffer Components: 20 mM Tris-HCl (pH 7.5), 25 mM KCl, 0.1 mM EDTA, 2 mM DTT and 50% (v/v) glycerol.
Unit Definition: One unit of the enzyme transfers 1 nmol of g-phosphate from ATP to 5'-OH DNA in 30 min at 37°C. Enzyme activity is assayed in the following mixture: 100 mM Tris-HCl (pH 8.0), 10 mM MgCl₂, 5 mM DTT, 0.5 mM 5'-OH DNA, 0.05 mM ATP and 0.1 MBq/mL [g -33P]-ATP.
Tested for: Endo-, exodeoxyribonuclease and ribonuclease-free; Functionally tested for labeling of 5'-termini of DNA.
Recommended storage: -20°C



OPTIZYME™ T7 RNA Polymerase

Source: Recombinant E. Coli Strain

packaging	Mfr. No
5000 units Poly Tube	BP8102-1
25000 units Poly Tube	BP8102-5

Tested for RNase activity, protease activity, and specific performance tests.

Applications: Synthesis of unlabeled & labeled RNA that can be used: (1) for hybridization and for in vitro RNA translation, (2) as RNA or siRNA, (3) as a substrate in RNase protection assays or as template for genomic DNA sequencing and (4) in studies of RNA secondary structure and RNA-protein interactions.

Description: OPTIZYME™ T7 RNA Polymerase is a DNA dependent RNA polymerase with strict specificity for its respective double-stranded promoter. The enzyme catalyzes the 5 -> 3’ synthesis of RNA on either single-stranded DNA or double-stranded DNA downstream from the T7 promoter. Bacteriophage T7 RNA Polymerase accepts modified nucleotides (e.g., biotin-, digoxigenin-, fluorescein-labeled nucleotides) as substrates for RNA synthesis.

Supplied With: 5X OPTIZYME™ Transcription Buffer (200mM Tris-HCl (pH 7.9 at 25°C), 30mM MgCl2, 50mM DTT, 50mM NaCl, 10mM spermidine).

Concentration: 20 u/μl

Storage Buffer Components: 50 mM Tris-HCl (pH 8.0), 150 mM NaCl, 5 mM DTT, 0.1 mg/ml BSA, 0.03% (v/v) ELUGENT Detergent and 50% (v/v) glycerol.

Unit Definition: One unit of the enzyme incorporates 1 nmol of AMP into a polynucleotide fraction (adsorbed on DE-81) in 60 minutes at 37°C. Enzyme activity is assayed in the following mixture: 40 mM Tris-HCl (pH 8.0), 6 mM MgCl2, 10 mM DTT, 2 mM spermidine, 0.5 mM of each NTP, 0.6 MBq/ml -ATP, 20 μg/ml plasmid DNA containing the T7 promoter sequence.

Recommended storage: -20°C

OPTIZYME™ Terminal Deoxynucleotidyl Transferase

Source: Recombinant E. Coli Strain

packaging	Mfr. No
300 units Poly Tube	BP8103-1
1500 units Poly Tube	BP8103-5

Tested for RNase activity, protease activity, and specific performance tests.

Applications: Production of synthetic homo- and heteropolymers; Homopolymeric tailing of linear duplex DNA with any type of 3’-OH terminus; Oligodeoxyribonucleotide and DNA labeling; 5’-RACE (Rapid Amplification of cDNA Ends); In situ localization of apoptosis.

Description: OPTIZYME™ TdT, a template-independent DNA polymerase, catalyzes the repetitive addition of deoxyribonucleotides to the 3’-OH of oligodeoxyribonucleotides and single-stranded, or double-stranded DNA. TdT requires an oligonucleotide of at least three nucleotides to serve as a primer. With RNA as template, TdT shows variable performance which strongly depends upon the tertiary structure of acceptor RNA 3’-end and the nature of nucleotide. Generally, it is lower than using DNA as a template.

Supplied With: 5X OPTIZYME™ TdT Buffer (1M potassium cacodylate, 125mM Tris, 0.05% (v/v) Triton X-100, 5mM CoCl2 (pH 7.2 at 25°C))

Concentration: 20 u/μl

Storage Buffer Components: 100 mM potassium acetate (pH 6.8), 2 mM 2 -mercaptoethanol, 0.01% (v/v) Triton X-100 and 50% (v/v) glycerol.

Unit Definition: One unit of the enzyme catalyzes the incorporation of 1 nmol of deoxythymidylate into a polynucleotide fraction (adsorbed on DE-81) in 60 min at 37°C. Enzyme activity is assayed in the following mixture: 200 mM potassium cacodylate (pH 7.2), 1 mM CoCl2, 0.01% (v/v) Triton X 100, 10 μM oligo(dT)10, 1 mM dTTP and 0.4 MBq/ml -dTTP.

Tested for: Endo-, exodeoxyribonuclease, phosphatase and ribonuclease-free

Recommended storage: -20°C

Ribonuclease H

packaging	Mfr. No
50 units PolyTube	BP3215-1

CAS: 9050-76-4

H315, H319, H335

P280, P305+P351+P338

Tested for RNase activity, protease activity, and specific performance tests.

Applications: Ribonuclease H is used in removal of the RNA strand prior to second strand cDNA synthesis and analysis of in vitro polyadenylation reaction products.

Description: Ribonuclease H (RNase H) from E. coli is an endoribonuclease that specifically hydrolyzes the phosphodiester bonds of RNA hybridized to DNA to produce 3’-OH and 5’-P terminated products. It will not degrade single-stranded nucleic acids, duplex DNA or double-stranded RNA.

Concentration: 0.5-2u/μl

Storage Buffer Components: 20mm HEPES-KOH (pH 7.8), 1mM DTT, 50mM KCl and 50% (v/v) glycerol, and 0.2mg/ml BSA.

Unit Definition: One unit is defined as the amount of enzyme required to produce 1nmol acid-soluble ribonucleotides from radiolabeled poly(rA):poly(dT) in 20 minutes at 37°C in 20mM HEPES-KOH (pH 7.8), 50mM KCl, 10mM MgCl2, 1mM DTT, 20μM radiolabeled poly(rA):poly(dT).

Recommended Storage: -20°C

OPTIZYME™ Aarl

Source: E. coli

packaging	Mfr. No
125 units Poly Tube	BP8069-1

Buffer 1	No Reaction%
Buffer 2	No Reaction%
Buffer 3	0-20%
Buffer 4	No Reaction%
Buffer 5	0-20%
Buffer Aarl, Scal, Pasl	100%
Tested for	RNase activity, protease activity, and specific performance tests.

Applications: OPTIZYME™ Aarl digests dsDNA with the recognition sequence indicated below.

Conditions for 100% Activity: 1X OPTIZYME Buffer Aarl, Scal, Pasl: 10mM Bis-Tris Propane-HCl (pH 6.5 at 37°C), 10mM MgCl2, 100mM KCl, 0.1mg/ml BSA; 0.5μM of oligonucleotide; Incubate at 37°C.

Recognition Sequence: 5’...C A C C T G C (N)4 ^...3’

Supplied With: 10X OPTIZYME Buffer Aarl, Scal, Pasl and 50X Oligonucleotide

Concentration: 2u/μl

Storage Buffer Components: 10mM potassium phosphate (pH 7.4 at 25°C), 100mM KCl, 1mM EDTA, 1mM DTT, 0.2mg/ml BSA and 50% (v/v) glycerol

Unit Definition: One unit is defined as the amount of Aarl at which no change in the fragmentation pattern is observed with further increase of enzyme. 1 μg of lambda DNA is incubated with the enzyme for 1 hour at 37°C in 50 μl of recommended reaction buffer. The cleavage of DNA by Aarl is never complete.

Tested for: Overdigestion Assay, Ligation/Re-Cutting Assay, Blue/White Colony Certified, Genome Qualified

Recommended storage: -20°C

OPTIZYME™ AatII

Source: E. coli

packaging	Mfr. No
300 units Poly Tube	BP8041-1

Buffer 1	50-100%
Buffer 2	20-50%
Buffer 3	0-20%
Buffer 4	100%
Buffer 5	0-20%

Applications: OPTIZYME™ AatII digests dsDNA with the recognition sequence indicated below.

Recognition Sequence: 5’...G A C G T ^C...3’

Supplied With: 10X OPTIZYME Buffer 4

Conditions for 100% Activity: 1X OPTIZYME Buffer 4: 33mM Tris-acetate (pH 7.9 at 37°C), 10mM Mg-acetate, 66mM K-acetate and 0.1mg/ ml BSA; Incubate at 37°C

Concentration: 10u/μl

Storage Buffer Components: 10mM Tris-HCl (pH 7.4 at 25°C), 100mM KCl, 1mM DTT, 1mM EDTA, 0.2mg/ml BSA and 50% (v/v) glycerol

Unit Definition: One unit is defined as the amount of AatII required to digest 1 μg of lambda DNA in 1 hour at 37°C in 50 μl of recommended reaction buffer.

Tested for: Overdigestion Assay, Ligation/Re-Cutting Assay

Recommended storage: -20°C

OPTIZYME™ Alol

Source: E. coli

packaging	Mfr. No
100 units Poly Tube	BP8075-1

Buffer 1	0-20%
Buffer 2	0-20%
Buffer 3	0-20%
Buffer 4	20-50%
Buffer 5	100%
Tested for	RNase activity, protease activity, and specific performance tests.

Applications: OPTIZYME™ Alol digests dsDNA with the recognition sequence indicated below.

Recognition Sequence: 5’...^ 7(N) G A A C (N)6 T C C (N)12-13^...3’

Supplied With: 10X OPTIZYME Buffer 5

Conditions for 100% Activity: 1X OPTIZYME Buffer 5: 10mM Tris-HCl (pH 8.5 at 37°C), 10mM MgCl2, 100mM KCl and 0.1mg/ml BSA; Incubate at 30°C.

Concentration: 1-3u/μl

Storage Buffer Components: 10mM Tris-HCl (pH 7.4 at 25°C), 100mM KCl, 1mM DTT, 1mM EDTA, 0.2mg/ml BSA and 50% (v/v) glycerol

Unit Definition: One unit is defined as the amount of Alol required to digest 1 μg of lambda DNA in 1 hour at 30°C in 50 μl of recommended reaction buffer.

Tested for: Overdigestion Assay, Ligation/Re-Cutting Assay, Blue/White Colony Certified

Recommended storage: -20°C

OPTIZYME™ Alul

Source: Arthrobacter luteus

packaging	Mfr. No
600 units Poly Tube	BP8015-1

Buffer 1	50-100%
Buffer 2	0-20%
Buffer 3	0-20%
Buffer 4	100%
Buffer 5	0-20%
Tested for	RNase activity, protease activity, and specific performance tests.

Applications: OPTIZYME™ Alul digests dsDNA with the recognition sequence indicated below.

Recognition Sequence: 5’...A G ^C T...3’

Supplied With: 10X OPTIZYME Buffer 4

Conditions for 100% Activity: 1X OPTIZYME Buffer 4: 33mM Tris-acetate (pH 7.9 at 37°C), 10mM Mg-acetate, 66mM K-acetate and 0.1mg/ml BSA; Incubate at 37°C

Concentration: 10u/μl

Storage Buffer Components: 10mM Tris-HCl (pH 7.5 at 25°C), 50mM KCl, 1 mM DTT, 0.1mM EDTA, 0.15% Triton X-100, 0.2mg/ml BSA and 50% (v/v) glycerol

Unit Definition: One unit is defined as the amount of Alul required to digest 1 μg of lambda DNA in 1 hour at 37°C in 50 μl of recommended reaction buffer.

Tested for: Overdigestion Assay, Ligation/Re-Cutting Assay

Recommended storage: -20°C

OPTIZYME™ Alw44I

Source: Acinetobacter lwoffii RFL44

packaging	Mfr. No
1000 units Poly Tube	BP8059-1

Buffer 1	50-100%
Buffer 2	100%
Buffer 3	0-20%
Buffer 4	100%
Buffer 5	50-100%

Applications: OPTIZYME™ Alw44I digests dsDNA with the recognition sequence indicated below.

Recognition Sequence: 5’...G ^T G C A C...3’

Supplied With: 10X OPTIZYME Buffer 4

Conditions for 100% Activity: 1X OPTIZYME Buffer 4: 33mM Tris-acetate (pH 7.9 at 37°C), 10mM Mg-acetate, 66mM K-acetate and 0.1mg/ ml BSA; Incubate at 37°C

Concentration: 10u/μl

Storage Buffer Components: 10mM Tris-HCl (pH 7.5 at 25°C), 50mM KCl, 1mM DTT, 0.1mM EDTA, 0.2mg/ml BSA and 50% (v/v) glycerol

Unit Definition: One unit is defined as the amount of Alw44I required to digest 1 μg of lambda DNA-Smal fragments in 1 hour at 37°C in 50 μl of recommended reaction buffer.

Tested for: Overdigestion Assay, Ligation/Re-Cutting Assay, Genome Qualified

Recommended storage: -20°C

OPTIZYME™ Apal Source: E. coli	
packaging	Mfr. No
5000 units Poly Tube	BP8025-1
Buffer 1	100%
Buffer 2	20-50%
Buffer 3	0-20%
Buffer 4	20-50%
Buffer 5	0-20%

Applications: OPTIZYME™ Apal digests dsDNA with the recognition sequence indicated below.
Recognition Sequence: 5’...G G G C C^A C...3’
Supplied With: 10X OPTIZYME Buffer 1
Conditions for 100% Activity: 1X OPTIZYME Buffer 1: 10mM Tris-HCl (pH 7.5 at 37°C), 10mM MgCl2 and 0.1mg/ml BSA; Incubate at 37°C.
Concentration: 10u/μl
Storage Buffer Components: 10mM Tris-HCl (pH 7.4 at 25°C), 50mM NaCl, 1 mM DTT, 1mM EDTA, 0.2mg/ml BSA and 50% (v/v) glycerol
Unit Definition: One unit is defined as the amount of Apal required to digest 1 μg of lambda DNA-Cpol fragments in 1 hour at 37°C in 50 μl of recommended reaction buffer.
Tested for: Overdigestion Assay, Ligation/Re-Cutting Assay, Blue/White Colony Certified, Genome Qualified
Recommended storage: -20°C

OPTIZYME™ Aval Source: E. coli	
packaging	Mfr. No
1000 units Poly Tube	BP8035-1
Buffer 1	100%
Buffer 2	50-100%
Buffer 3	0-20%
Buffer 4	100%
Buffer 5	0-20%

Applications: OPTIZYME™ Aval digests dsDNA with the recognition sequence indicated below.
Recognition Sequence: 5’...C^Y C G R G...3’
Supplied With: 10X OPTIZYME Buffer 4
Conditions for 100% Activity: 1X OPTIZYME Buffer 4: 33mM Tris-acetate (pH 7.9 at 37°C), 10mM Mg-acetate, 66mM K-acetate and 0.1mg/ml BSA; Incubate at 37°C
Concentration: 10u/μl
Storage Buffer Components: 10mM Tris-HCl (pH 7.4 at 25°C), 100mM KCl, 1mM DTT, 1mM EDTA, 0.2mg/ml BSA and 50% (v/v) glycerol
Unit Definition: One unit is defined as the amount of Aval required to digest 1 μg of lambda DNA in 1 hour at 37°C in 50 μl of recommended reaction buffer.
Tested for: Overdigestion Assay, Ligation/Re-Cutting Assay, Blue/White Colony Certified, Genome Qualified
Recommended storage: -20°C

OPTIZYME™ Avall Source: E. coli	
packaging	Mfr. No
800 units Poly Tube	BP8043-1
Buffer 1	0-20%
Buffer 2	50-100%
Buffer 3	50-100%
Buffer 4	50-100%
Buffer 5	100%

Applications: OPTIZYME™ Avall digests dsDNA with the recognition sequence indicated below.
Recognition Sequence: 5’...G^G W C C...3’
Supplied With: 10X OPTIZYME Buffer 5
Conditions for 100% Activity: 1X OPTIZYME Buffer 5: 10mM Tris-HCl (pH 8.5 at 37°C), 10mM MgCl2, 100mM KCl and 0.1mg/ml BSA; Incubate at 37°C
Concentration: 10u/μl
Storage Buffer Components: 10mM Tris-HCl (pH 7.4 at 25°C), 100mM KCl, 1mM DTT, 1mM EDTA, 0.2mg/ml BSA and 50% (v/v) glycerol
Unit Definition: One unit is defined as the amount of Avall required to digest 1 μg of lambda DNA in 1 hour at 37°C in 50 μl of recommended reaction buffer.
Tested for: Overdigestion Assay, Ligation/Re-Cutting Assay
Recommended storage: -20°C

OPTIZYME™ Ball Source: Micrococcus luteus Ng 16-122	
packaging	Mfr. No
200 units Poly Tube	BP8039-1
Buffer 1	0-20%
Buffer 2	20-50%
Buffer 3	0-20%
Buffer 4	20-50%
Buffer 5	100%

Applications: OPTIZYME™ Ball digests dsDNA with the recognition sequence indicated below.
Recognition Sequence: 5’...T G G^C C A...3’
Supplied With: 10X OPTIZYME Buffer 5
Conditions for 100% Activity: 1X OPTIZYME Buffer 5: 10mM Tris-HCl (pH 8.5 at 37°C), 10mM MgCl2, 100mM KCl and 0.1mg/ml BSA; Incubate at 37°C
Concentration: 5u/μl
Storage Buffer Components: 10mM Tris-HCl (pH 7.4 at 25°C), 100mM KCl, 1mM DTT, 1mM EDTA, 0.2mg/ml BSA and 50% (v/v) glycerol
Unit Definition: One unit is defined as the amount of Ball required to digest 1 μg of lambda DNA dcm- in 1 hour at 37°C in 50 μl of recommended reaction buffer.
Tested for: Overdigestion Assay, Ligation/Re-Cutting Assay, Blue/White Colony Certified, Genome Qualified
Recommended storage: -20°C

OPTIZYME™ BamHI Source: E. coli	
packaging	Mfr. No
12500 units Poly Tube	BP8005-5
Buffer 1	20-50%
Buffer 2	100%
Buffer 3	20-50%
Buffer 4	100%
Buffer 5	50-100%
Buffer BamHI	100%

Applications: OPTIZYME™ BamHI digests dsDNA with the recognition sequence indicated below.
Recognition Sequence: 5’...G^G A T C C...3’
Supplied With: 10X OPTIZYME BamHI Buffer
Conditions for 100% Activity: 1X OPTIZYME Buffer BamHI: 10mM Tris-HCl (pH 8.0 at 37°C), 5mM MgCl2, 100mM KCl, 0.02% Triton X-100 and 0.1mg/ml BSA; Incubate at 37°C.
Concentration: 10u/μl
Storage Buffer Components: 10mM Tris-HCl (pH 7.4 at 25°C), 200mM NaCl, 1mM DTT, 0.1mM EDTA, 0.15% Triton X-100, 0.2mg/ml BSA and 50% (v/v) glycerol.
Unit Definition: One unit is defined as the amount of BamHI required to digest 1 μg of lambda DNA-Bsp120I fragments in 1 hour at 37°C in 50 μl of recommended reaction buffer.
Tested for: Overdigestion Assay, Ligation/Re-Cutting Assay, Blue/White Colony Certified, Genome Qualified
Recommended storage: -20°C

OPTIZYME™ BclI Source: Bacillus caldolyticus	
packaging	Mfr. No
3000 units Poly Tube	BP8053-1
Buffer 1	20-50%
Buffer 2	100%
Buffer 3	20-50%
Buffer 4	100%
Buffer 5	20-50%

Applications: OPTIZYME™ BclI digests dsDNA with the recognition sequence indicated below.
Recognition Sequence: 5’...T^A G A T C A...3’
Supplied With: 10X OPTIZYME Buffer 2
Conditions for 100% Activity: 1X OPTIZYME Buffer 2: 10mM Tris-HCl (pH 7.5 at 37°C), 10mM MgCl2, 50mM NaCl and 0.1mg/ml BSA; Incubate at 55°C.
Concentration: 10u/μl
Storage Buffer Components: 10mM Tris-HCl (pH 7.4 at 25°C), 100mM KCl, 1mM DTT, 1mM EDTA, 0.2mg/ml BSA and 50% (v/v) glycerol
Unit Definition: One unit is defined as the amount of BclI required to digest 1 μg of lambda DNA dam- in 1 hour at 55°C in 50 μl of recommended reaction buffer.
Tested for: Overdigestion Assay, Ligation/Re-Cutting Assay, Blue/White Colony Certified
Recommended storage: -20°C

OPTIZYME™ BglI Bacillus globigii	
packaging	Mfr. No
2000 units Poly Tube	BP8046-1
Buffer 1	0-20%
Buffer 2	50-100%
Buffer 3	100%
Buffer 4	0-20%
Buffer 5	100%

Applications: OPTIZYME™ BglI digests dsDNA with the recognition sequence indicated below.
Recognition Sequence: 5’...G C C N N N N^N G G C...3’
Supplied With: 10X OPTIZYME Buffer 3
Conditions for 100% Activity: 1X OPTIZYME Buffer 3: 50mM Tris-HCl (pH 7.5 at 37°C), 10mM MgCl2, 100mM NaCl and 0.1mg/ml BSA; Incubate at 37°C
Concentration: 10u/μl
Storage Buffer Components: 10mM Tris-HCl (pH 7.5 at 25°C), 300mM NaCl, 1mM DTT, 0.1mM EDTA, 0.2mg/ml BSA and 50% (v/v) glycerol
Unit Definition: One unit is defined as the amount of BglI required to digest 1 μg of lambda DNA in 1 hour at 37°C in 50 μl of recommended reaction buffer.
Tested for: Overdigestion Assay, Ligation/Re-Cutting Assay, Genome Qualified
Recommended storage: -20°C

OPTIZYME™ BglII Source: E. coli	
packaging	Mfr. No
1000 units Poly Tube	BP8014-1
2500 units Poly Tube	BP8014-5
Buffer 1	0-20%
Buffer 2	20-50%
Buffer 3	100%
Buffer 4	0-20%
Buffer 5	50-100%

Applications: OPTIZYME™ BglII digests dsDNA with the recognition sequence indicated below.
Recognition Sequence: 5’...A^A G A T C T...3’
Supplied With: 10X OPTIZYME Buffer 3
Conditions for 100% Activity: 1X OPTIZYME Buffer 3: 50mM Tris-HCl (pH 7.5 at 37°C), 10mM MgCl2, 100mM NaCl and 0.1mg/ml BSA; Incubate at 37°C
Concentration: 10u/μl
Storage Buffer Components: 10mM Tris-HCl (pH 8.2 at 25°C), 200 mM KCl, 1mM DTT, 0.1mM EDTA, 0.5mg/ml BSA and 50% (v/v) glycerol
Unit Definition: One unit is defined as the amount of BglII required to digest 1 μg of lambda DNA in 1 hour at 37°C in 50 μl of recommended reaction buffer.
Tested for: Overdigestion Assay, Ligation/Re-Cutting Assay, Blue/White Colony Certified, Genome Qualified
Recommended storage: -20°C



OPTIZYME™ Bpil
Source: Bacillus pumillus Sw 4-3

packaging	Mfr. No
1000 units Poly Tube	BP8072-1

Buffer 1	20-50%
Buffer 2	100%
Buffer 3	50-100%
Buffer 4	50-100%
Buffer 5	50-100%

Applications: OPTIZYME™ Bpil digests dsDNA with the recognition sequence indicated below.
Recognition Sequence: 5’...G A A G A C (N)2^...3’
Supplied With: 10X OPTIZYME Buffer 2
Conditions for 100% Activity: 1X OPTIZYME Buffer 2: 10mM Tris-HCl (pH 7.5 at 37°C), 10mM MgCl2, 50mM NaCl and 0.1mg/ml BSA; Incubate at 37°C.
Concentration: 10u/μl
Storage Buffer Components: 10mM Tris-HCl (pH 7.4 at 25°C), 100mM KCl, 1mM DTT, 1mM EDTA, 0.2mg/ml BSA and 50% (v/v) glycerol
Unit Definition: One unit is defined as the amount of Bpil required to digest 1 μg of lambda DNA in 1 hour at 37°C in 50 μl of recommended reaction buffer.
Tested for: Overdigestion Assay, Ligation/Re-Cutting Assay, Blue/White Colony Certified, Genome Qualified
Recommended storage: -20°C

OPTIZYME™ Bsh1236I
Source: Bacillus sphaericus RFL1236

packaging	Mfr. No
2500 units Poly Tube	BP8071-1

Buffer 1	0-20%
Buffer 2	0-20%
Buffer 3	50-100%
Buffer 4	20-50%
Buffer 5	100%

Applications: OPTIZYME™ Bsh1236I digests dsDNA with the recognition sequence indicated below.
Recognition Sequence: 5’...C G^C G...3’
Supplied With: 10X OPTIZYME Buffer 5
Conditions for 100% Activity: 1X OPTIZYME Buffer 5: 10mM Tris-HCl (pH 8.5 at 37°C), 10mM MgCl2, 100mM KCl and 0.1mg/ml BSA; Incubate at 37°C.
Concentration: 10u/μl
Storage Buffer Components: 10mM Tris-HCl (pH 7.5 at 25°C), 50mM KCl, 1mM DTT, 0.1mM EDTA, 0.2mg/ml BSA and 50% (v/v) glycerol
Unit Definition: One unit is defined as the amount of Bsh1236I required to digest 1 μg lambda DNA in 1 hour at 37°C in 50 μl of reaction buffer.
Tested for: Overdigestion Assay, Ligation/Re-Cutting Assay
Recommended storage: -20°C

OPTIZYME™ BshTI
Source: Bacillus sphaericus Jo 22-024

packaging	Mfr. No
1000 units Poly Tube	BP8078-1

Buffer 1	0-20%
Buffer 2	20-50%
Buffer 3	100%
Buffer 4	20-50%
Buffer 5	50-100%

Applications: OPTIZYME™ BshTI digests dsDNA with the recognition sequence indicated below.
Recognition Sequence: 5’...A^C C G G T...3’
Supplied With: 10X OPTIZYME Buffer 3
Conditions for 100% Activity: 1X OPTIZYME Buffer 3: 50mM Tris-HCl (pH 7.5 at 37°C), 10mM MgCl2, 100mM NaCl and 0.1mg/ml BSA; Incubate at 37°C
Concentration: 10u/μl
Storage Buffer Components: 10mM Tris-HCl (pH 7.4 at 25°C), 100mM KCl, 1mM DTT, 1mM EDTA, 0.2mg/ml BSA and 50% (v/v) glycerol
Unit Definition: One unit is defined as the amount of BshTI required to digest 1 μg of lambda DNA in 1 hour at 37°C in 50 μl of recommended reaction buffer.
Tested for: Overdigestion Assay, Ligation/Re-Cutting Assay, Blue/White Colony Certified, Genome Qualified
Recommended storage: -20°C

OPTIZYME™ BssHII
Source: Paracoccus alcaliphilus ZVK3-3

packaging	Mfr. No
200 units Poly Tube	BP8036-1

Buffer 1	0-20%
Buffer 2	0-20%
Buffer 3	100%
Buffer 4	0-20%
Buffer 5	100%

Applications: OPTIZYME™ BssHII digests dsDNA with the recognition sequence indicated below.
Recognition Sequence: 5’...G^C G C G C...3’
Supplied With: 10X OPTIZYME Buffer 5
Conditions for 100% Activity: 1X OPTIZYME Buffer 5: 10mM Tris-HCl (pH 8.5 at 37°C), 10mM MgCl2, 100mM KCl and 0.1mg/ml BSA; Incubate at 37°C
Concentration: 10u/μl
Storage Buffer Components: 10mM Tris-HCl (pH 7.5 at 25°C), 50mM KCl, 1mM DTT, 0.1mM EDTA, 0.2mg/ml BSA and 50% (v/v) glycerol
Unit Definition: One unit is defined as the amount of BssHII required to digest 1 μg lambda DNA in 1 hour at 37°C in 50 μl of recommended reaction buffer.
Tested for: Overdigestion Assay, Ligation/Re-Cutting Assay, Blue/White Colony Certified, Genome Qualified
Recommended storage: -20°C

OPTIZYME™ BstEII
Source: Escherichia coli RFL91

packaging	Mfr. No
2000 units Poly Tube	BP8038-1

Buffer 1	20-50%
Buffer 2	20-50%
Buffer 3	100%
Buffer 4	No Reaction%
Buffer 5	50-100%

Applications: OPTIZYME™ BstEII digests dsDNA with the recognition sequence indicated below.
Recognition Sequence: 5’...G^G T N A C C...3’
Supplied With: 10X OPTIZYME Buffer 3
Conditions for 100% Activity: 1X OPTIZYME Buffer 3: 50mM Tris-HCl (pH 7.5 at 37°C), 10mM MgCl)2 , 100mM NaCl and 0.1mg/ml BSA; Incubate at 37°C
Concentration: 10u/μl
Storage Buffer Components: 10mM Tris-HCl (pH 7.4 at 25°C), 100mM KCl, 1mM DTT, 1mM EDTA, 0.2mg/ml BSA and 50% (v/v) glycerol
Unit Definition: One unit is defined as the amount of BstEII required to digest 1 μg of lambda DNA in 1 hour at 37°C in 50 μl of recommended reaction buffer.
Tested for: Overdigestion Assay, Ligation/Re-Cutting Assay, Blue/White Colony Certified, Genome Qualified
Recommended storage: -20°C

OPTIZYME™ Cfr9I
Source: E. coli

packaging	Mfr. No
1500 units Poly Tube	BP8081-1

Buffer 1	0-20%
Buffer 2	0-20%
Buffer 3	0-20%
Buffer 4	20-50%
Buffer 5	0-20%
Buffer Cfr9I	100%

Applications: OPTIZYME™ Cfr9I digests dsDNA with the recognition sequence indicated below.
Recognition Sequence: 5’...C^C C G G G...3’
Supplied With: 10X OPTIZYME Cfr9I
Conditions for 100% Activity: 1X OPTIZYME Buffer Cfr9I: 10mM Tris-HCl (pH 7.2 at 37°C), 5mM MgCl2, 200mM sodium glutamate and 0.1mg/ml BSA; Incubate at 37°C.
Concentration: 10u/μl
Storage Buffer Components: 10mM Tris-HCl (pH 7.5 at 25°C), 250mM KCl, 1mM DTT, 0.1mM EDTA, 0.2mg/ml BSA and 50% (v/v) glycerol
Unit Definition: One unit is defined as the amount of Cfr9I required to digest 1 μg of lambda DNA-HindIII fragments in 1 hour at 37°C in 50 μl of recommended reaction buffer (containing 2 μg DNA fragments).
Tested for: Overdigestion Assay, Ligation/Re-Cutting Assay, Genome Qualified
Recommended storage: -20°C

OPTIZYME™ Clal
Source: Bacillus subtilis 15

packaging	Mfr. No
600 units Poly Tube	BP8024-1
2500 units Poly Tube	BP8024-5

Buffer 1	20-50%
Buffer 2	20-50%
Buffer 3	20-50%
Buffer 4	100%
Buffer 5	20-50%

Applications: OPTIZYME™ Clal digests dsDNA with the recognition sequence indicated below.
Recognition Sequence: 5’...A T^C G A T...3’
Supplied With: 10X OPTIZYME Buffer 4
Conditions for 100% Activity: 1X OPTIZYME Buffer 4: 33mM Tris-acetate (pH 7.9 at 37°C), 10mM Mg-acetate, 66mM K-acetate and 0.1mg/ml BSA; Incubate at 37°C
Concentration: 10u/μl
Storage Buffer Components: 10mM Tris-HCl (pH 8.0 at 25°C), 50mM KCl, 1mM DTT, 0.1mM EDTA, 0.5mg/ml BSA and 50% (v/v) glycerol
Unit Definition: One unit is defined as the amount of Clal required to digest 1 μg of lambda DNA in 1 hour at 37°C in 50 μl of recommended reaction buffer.
Tested for: Overdigestion Assay, Ligation/Re-Cutting Assay, Blue/White Colony Certified, Genome Qualified
Recommended storage: -20°C

OPTIZYME™ Csp6I
Source: Corynebacterium species RFL6

packaging	Mfr. No
1500 units Poly Tube	BP8068-1

Buffer 1	100%
Buffer 2	50-100%
Buffer 3	0-20%
Buffer 4	50-100%
Buffer 5	0-20%

Applications: OPTIZYME™ Csp6I digests dsDNA with the recognition sequence indicated below.
Recognition Sequence: 5’...G^T A C...3’
Supplied With: 10X OPTIZYME Buffer 1
Conditions for 100% Activity: 1X OPTIZYME Buffer 1: 10mM Tris-HCl (pH 7.5 at 37°C), 10mM MgCl)2 and 0.1mg/ml BSA; Incubate at 37°C.
Concentration: 10u/μl
Storage Buffer Components: 10mM Tris-HCl (pH 7.4 at 25°C), 100mM KCl, 1mM DTT, 1mM EDTA, 0.2mg/ml BSA and 50% (v/v) glycerol
Unit Definition: One unit is defined as the amount of Csp6I required to digest 1 μg of lambda DNA in 1 hour at 37°C in 50 μl of recommended reaction buffer.
Tested for: Overdigestion Assay, Ligation/Re-Cutting Assay
Recommended storage: -20°C



OPTIZYME™ Ddel Source: E. coli	
packaging	Mfr. No
500 units Poly Tube	BP8060-1
1000 units Poly Tube	BP8060-5

Buffer 1	20-50%
Buffer 2	20-50%
Buffer 3	20-50%
Buffer 4	100%
Buffer 5	20-50%

Applications: OPTIZYME™ Ddel digests dsDNA with the recognition sequence indicated below.
Recognition Sequence: 5'...C^T N A G...3'
Supplied With: 10X OPTIZYME Buffer 4
Conditions for 100% Activity: 1X OPTIZYME Buffer 4: 33mM Tris-acetate (pH 7.9 at 37°C), 10mM Mg-acetate, 66mM K-acetate and 0.1mg/ml BSA; Incubate at 37°C
Concentration: 10u/μl
Storage Buffer Components: 10mM Tris-HCl (pH 7.4 at 25°C), 100mM KCl, 1mM DTT, 1mM EDTA, 0.2mg/ml BSA and 50% (v/v) glycerol
Unit Definition: One unit is defined as the amount of Ddel required to digest 1 μg of lambda DNA in 1 hour at 37°C in 50 μl of recommended reaction buffer.
Tested for: Overdigestion Assay, Ligation/Re-Cutting Assay, Genome Qualified
Recommended storage: -20°C

OPTIZYME™ DpnI Source: E. coli	
packaging	Mfr. No
500 units Poly Tube	BP8009-1

Buffer 1	100%
Buffer 2	100%
Buffer 3	50-100%
Buffer 4	100%
Buffer 5	50-100%

Applications: OPTIZYME™ DpnI digests dsDNA with the recognition sequence indicated below.
Recognition Sequence: 5'...G m6A^T C...3'
Supplied With: 10X OPTIZYME Buffer 4
Conditions for 100% Activity: 1X OPTIZYME Buffer 4: 33mM Tris-acetate (pH 7.9 at 37°C), 10mM Mg-acetate, 66mM K-acetate and 0.1mg/ml BSA; Incubate at 37°C
Concentration: 10u/μl
Storage Buffer Components: 10mM Tris-HCl (pH 7.4 at 25°C), 400mM KCl, 1mM DTT, 0.1mM EDTA, 0.2mg/ml BSA and 50% (v/v) glycerol
Unit Definition: One unit is defined as the amount of DpnI required to digest 1 μg of pBR322 DNA (dam methylated) in 1 hour at 37°C in 50 μl of recommended reaction buffer.
Tested for: Overdigestion Assay, Ligation/Re-Cutting Assay
Recommended storage: -20°C

OPTIZYME™ DraI Source: Deinococcus radiophilus	
packaging	Mfr. No
2000 units Poly Tube	BP8026-1

Buffer 1	50-100%
Buffer 2	50-100%
Buffer 3	20-50%
Buffer 4	100%
Buffer 5	20-50%

Applications: OPTIZYME™ DraI digests dsDNA with the recognition sequence indicated below.
Recognition Sequence: 5'...T T T^A A A...3'
Supplied With: 10X OPTIZYME Buffer 4
Conditions for 100% Activity: 1X OPTIZYME Buffer 4: 33mM Tris-acetate (pH 7.9 at 37°C), 10mM Mg-acetate, 66mM K-acetate and 0.1mg/ml BSA; Incubate at 37°C
Concentration: 10u/μl
Storage Buffer Components: 10mM Tris-HCl (pH 7.5 at 25°C), 50mM KCl, 1mM DTT, 0.1mM EDTA, 0.15% Triton X-100, 0.2mg/ml BSA and 50% (v/v) glycerol
Unit Definition: One unit is defined as the amount of DraI required to digest 1 μg lambda DNA in 1 hour at 37°C in 50 μl of recommended reaction buffer.
Tested for: Overdigestion Assay, Ligation/Re-Cutting Assay
Recommended storage: -20°C

OPTIZYME™ Ecl136II Source: Enterobacter cloacae RFL136	
packaging	Mfr. No
1500 units Poly Tube	BP8080-1

Buffer 1	50-100%
Buffer 2	20-50%
Buffer 3	0-20%
Buffer 4	50-100%
Buffer 5	0-20%
Buffer Ecl136II, SacI	100%

Applications: OPTIZYME™ Ecl136II digests dsDNA with the recognition sequence indicated below.
Recognition Sequence: 5'...G A G^C T C...3'
Supplied With: 10X OPTIZYME Ecl136II, SacI
Conditions for 100% Activity: 1X Buffer Ecl136II, SacI: 10mM Bis-Tris Propane-HCl (pH 6.5 at 37°C), 10mM MgCl2 and 0.1mg/ml BSA; Incubate at 37°C.
Concentration: 10u/μl
Storage Buffer Components: 10mM Tris-HCl (pH 7.4 at 25°C), 100mM KCl, 1mM DTT, 1mM EDTA, 0.2mg/ml BSA and 50% (v/v) glycerol
Unit Definition: One unit is defined as the amount of Ecl136II required to digest 1 μg of lambda DNA-HindIII in 1 hour at 37°C in 50 μl of recommended reaction buffer.
Tested for: Overdigestion Assay, Ligation/Re-Cutting Assay, Blue/White Colony Certified, Genome Qualified
Recommended storage: -20°C

OPTIZYME™ Eco57I Source: E. coli	
packaging	Mfr. No
1000 units Poly Tube	BP8066-1

Buffer 1	100%
Buffer 2	100%
Buffer 3	20-50%
Buffer 4	50-100%
Buffer 5	20-50%

Applications: OPTIZYME™ Eco57I digests dsDNA with the recognition sequence indicated below.
Recognition Sequence: 5'...C T G A A G (N)16^...3'
Supplied With: 10X OPTIZYME Buffer 2 and 50X S-adenosylmethionine.
Conditions for 100% Activity: 1X OPTIZYME Buffer 2: 10mM Tris-HCl (pH 7.5 at 37°C), 10mM MgCl2 , 50mM NaCl and 0.1mg/ml BSA; 0.01mM S-adenosylmethionine; Incubate at 37°C
Concentration: 5u/μl
Storage Buffer Components: 10mM potassium phosphate (pH 7.4 at 25°C), 100mM NaCl, 1mM EDTA, 1mM DTT and 50% (v/v) glycerol
Unit Definition: One unit is defined as the amount of Eco57I at which no change in the fragmentation pattern is observed with further increase of enzyme. 1 μg of lambda DNA is incubated with Eco57I for 1 hour at 37°C in 50 μl of recommended reaction buffer.
Tested for: Overdigestion Assay, Ligation/Re-Cutting Assay
Recommended storage: -20°C

OPTIZYME™ EcoRI Source: E. coli	
packaging	Mfr. No
15000 units Poly Tube	BP8003-5

Buffer 1	0-20%
Buffer 2	No Reaction%
Buffer 3	100%
Buffer 4	No Reaction%
Buffer 5	100%
Buffer EcoRI	100%

Applications: OPTIZYME™ EcoRI digests dsDNA with the recognition sequence indicated below.
Recognition Sequence: 5'...G^A A T T C...3'
Supplied With: 10X OPTIZYME EcoRI Buffer
Conditions for 100% Activity: 1X OPTIZYME Buffer EcoRI: 50mM Tris-HCl (pH 7.5 at 37°C), 10mM MgCl2, 100mM NaCl, 0.02% Triton X-100 and 0.1mg/ml BSA; Incubate at 37°C.
Concentration: 10u/μl
Storage Buffer Components: 10mM potassium phosphate (pH 7.4 at 25°C), 300mM NaCl, 1mM EDTA, 1mM DTT, 0.2mg/ml BSA, 0.15% Triton X-100 and 50% (v/v) glycerol.
Unit Definition: One unit is defined as the amount of EcoRI required to digest 1 μg of lambda DNA in 1 hour at 37°C in 50 μl of recommended reaction buffer.
Tested for: Overdigestion Assay, Ligation/Re-Cutting Assay, Blue/White Colony Certified, Genome Qualified
Recommended storage: -20°C

OPTIZYME™ EcoRI Source: E. coli	
packaging	Mfr. No
25000 units Poly Tube	BP8054-1

Buffer 1	0-20%
Buffer 2	No Reaction%
Buffer 3	100%
Buffer 4	No Reaction%
Buffer 5	100%
Buffer EcoRI	100%

Applications: OPTIZYME™ EcoRI digests dsDNA with the recognition sequence indicated below.
Recognition Sequence: 5'...G^A A T T C...3'
Supplied With: 10X OPTIZYME EcoRI Buffer
Conditions for 100% Activity: 1X OPTIZYME Buffer EcoRI: 50mM Tris-HCl (pH 7.5 at 37°C), 10mM MgCl2, 100mM NaCl, 0.02% Triton X-100 and 0.1mg/ml BSA; Incubate at 37°C.
Concentration: 50u/μl
Storage Buffer Components: 10mM potassium phosphate (pH 7.4 at 25°C), 300mM NaCl, 1mM EDTA, 1mM DTT, 0.2mg/ml BSA, 0.15% Triton X-100 and 50% (v/v) glycerol.
Unit Definition: One unit is defined as the amount of EcoRI required to digest 1 μg of lambda DNA in 1 hour at 37°C in 50 μl of recommended reaction buffer.
Tested for: Overdigestion Assay, Ligation/Re-Cutting Assay, Blue/White Colony Certified, Genome Qualified
Recommended storage: -20°C

OPTIZYME™ EcoRV Source: Escherichia coli RFL32	
packaging	Mfr. No
3000 units Poly Tube	BP8012-1

Buffer 1	0-20%
Buffer 2	50-100%
Buffer 3	50-100%
Buffer 4	20-50%
Buffer 5	100%

Applications: OPTIZYME™ EcoRV digests dsDNA with the recognition sequence indicated below.
Recognition Sequence: 5'...G A T^A T C...3'
Supplied With: 10X OPTIZYME Buffer 5
Conditions for 100% Activity: 1X OPTIZYME Buffer 5: 10mM Tris-HCl (pH 8.5 at 37°C), 10mM MgCl2, 100mM KCl and 0.1mg/ml BSA; Incubate at 37°C
Concentration: 10u/μl
Storage Buffer Components: 25mM Tris-HCl (pH 7.5 at 25°C), 200mM NaCl, 1mM DTT, 1mM EDTA, 0.2mg/ml BSA and 50% (v/v) glycerol
Unit Definition: One unit is defined as the amount of EcoRV required to digest 1 μg lambda DNA in 1 hour at 37°C in 50 μl of recommended reaction buffer.
Tested for: Overdigestion Assay, Ligation/Re-Cutting Assay, Blue/White Colony Certified, Genome Qualified
Recommended storage: -20°C



OPTIZYME™ Esp3I Source: E. coli	
packaging	Mfr. No
1000 units Poly Tube	BP8070-1
Buffer 1	100%
Buffer 2	20-50%
Buffer 3	0-20%
Buffer 4	100%
Buffer 5	0-20%

Applications: OPTIZYME™ Esp3I digests dsDNA with the recognition sequence indicated below.
Recognition Sequence: 5’...C G T C T C (N)1^...3’
Supplied With: 10X OPTIZYME Buffer 4
Conditions for 100% Activity: 1X OPTIZYME Buffer 4 + DTT: 33mM Tris-acetate (pH 7.9 at 37°C), 10mM Mg-acetate, 66mM K-acetate and 0.1mg/ml BSA; 1.0mM DTT; Incubate at 37°C.
Concentration: 10u/μl
Storage Buffer Components: 10mM potassium phosphate (pH 7.4 at 25°C), 100mM KCl, 1mM EDTA, 7mM 2-mercaptoethanol, 0.5mg/ml BSA and 50% (v/v) glycerol
Unit Definition: One unit is defined as the amount of Esp3I required to digest 1 μg of lambda DNA in 1 hour at 37°C in 50 μl of recommended reaction buffer.
Tested for: Overdigestion Assay, Ligation/Re-Cutting Assay, Genome Qualified
Recommended storage: -20°C

OPTIZYME™ HaeIII Source: Bacillus subtilis R	
packaging	Mfr. No
2500 units Poly Tube	BP8002-1
10000 units Poly Tube	BP8002-5
Buffer 1	0-20%
Buffer 2	0-20%
Buffer 3	0-20%
Buffer 4	50-100%
Buffer 5	100%

Applications: OPTIZYME™ HaeIII digests dsDNA with the recognition sequence indicated below.
Recognition Sequence: 5’...G G^C C...3’
Supplied With: 10X OPTIZYME Buffer 5
Conditions for 100% Activity: 1X OPTIZYME Buffer 5: 10mM Tris-HCl (pH 8.5 at 37°C), 10mM MgCl2 , 100mM KCl and 0.1mg/ml BSA; Incubate at 37°C.
Concentration: 10u/μl
Storage Buffer Components: 10mM Tris-HCl (pH 7.5 at 25°C), 50mM KCl, 1mM DTT, 0.1mM EDTA, 0.2mg/ml BSA and 50% (v/v) glycerol.
Unit Definition: One unit is defined as the amount of HaeIII required to digest 1 μg lambda DNA in 1 hour at 37°C in 50 μl of recommended reaction buffer.
Tested for: Overdigestion Assay, Ligation/Re-Cutting Assay
Recommended storage: -20°C

OPTIZYME™ HincII Source: E. coli	
packaging	Mfr. No
500 units Poly Tube	BP8034-1
Buffer 1	50-100%
Buffer 2	50-100%
Buffer 3	20-50%
Buffer 4	100%
Buffer 5	50-100%

Applications: OPTIZYME™ HincII digests dsDNA with the recognition sequence indicated below.
Recognition Sequence: 5’...G T Y^R A C...3’
Supplied With: 10X OPTIZYME Buffer 4
Conditions for 100% Activity: 1X OPTIZYME Buffer 4: 33mM Tris-acetate (pH 7.9 at 37°C), 10mM Mg-acetate, 66mM K-acetate and 0.1mg/ml BSA; Incubate at 37°C
Concentration: 10u/μl
Storage Buffer Components: 10mM Tris-HCl (pH 7.5 at 25°C), 200mM KCl, 1mM DTT, 0.1mM EDTA, 0.5mg/ml BSA and 50% (v/v) glycerol
Unit Definition: One unit is defined as the amount of HincII required to digest 1 μg of lambda DNA in 1 hour at 37°C in 50 μl of recommended reaction buffer.
Tested for: Overdigestion Assay, Ligation/Re-Cutting Assay
Recommended storage: -20°C

OPTIZYME™ HindIII Source: Haemophilus influenzae Rd	
packaging	Mfr. No
15000 units Poly Tube	BP8006-5
Buffer 1	0-20%
Buffer 2	20-50%
Buffer 3	0-20%
Buffer 4	50-100%
Buffer 5	100%

Applications: OPTIZYME™ HindIII digests dsDNA with the recognition sequence indicated below.
Recognition Sequence: 5’...A^A G C T T...3’
Supplied With: 10X OPTIZYME Buffer 5
Conditions for 100% Activity: 1X OPTIZYME Buffer 5: 10mM Tris-HCl (pH 8.5 at 37°C), 10mM MgCl2, 100mM KCl and 0.1mg/ml BSA; Incubate at 37°C
Concentration: 10u/μl
Storage Buffer Components: 10mM Tris-HCl (pH 7.5 at 25°C), 250mM KCl, 1mM DTT, 0.1mM EDTA, 0.2mg/ml BSA and 50% (v/v) glycerol
Unit Definition: One unit is defined as the amount of HindIII required to digest 1 μg lambda DNA in 1 hour at 37°C in 50 μl of recommended reaction buffer.
Tested for: Overdigestion Assay, Ligation/Re-Cutting Assay, Blue/White Colony Certified, Genome Qualified
Recommended storage: -20°C

OPTIZYME™ HinfI Source: Haemophilus influenzae Rf	
packaging	Mfr. No
2000 units Poly Tube	BP8051-1
Buffer 1	0-20%
Buffer 2	20-50%
Buffer 3	50-100%
Buffer 4	50-100%
Buffer 5	100%

Applications: OPTIZYME™ HinfI digests dsDNA with the recognition sequence indicated below.
Recognition Sequence: 5’...G^A N T C...3’
Supplied With: 10X OPTIZYME Buffer 5
Conditions for 100% Activity: 1X OPTIZYME Buffer 5: 10mM Tris-HCl (pH 8.5 at 37°C), 10mM MgCl2, 100mM KCl and 0.1mg/ml BSA; Incubate at 37°C
Concentration: 10u/μl
Storage Buffer Components: 10mM Tris-HCl (pH 7.4 at 25°C), 100mM KCl, 1mM DTT, 1mM EDTA, 0.2mg/ml BSA and 50% (v/v) glycerol
Unit Definition: One unit is defined as the amount of HinfI required to digest 1 μg lambda DNA in 1 hour at 37°C in 50 μl of recommended reaction buffer.
Tested for: Overdigestion Assay, Ligation/Re-Cutting Assay
Recommended storage: -20°C

OPTIZYME™ HpaI Source: Kurthia species N88	
packaging	Mfr. No
500 units Poly Tube	BP8049-1
Buffer 1	100%
Buffer 2	50-100%
Buffer 3	20-50%
Buffer 4	100%
Buffer 5	20-50%

Applications: OPTIZYME™ HpaI digests dsDNA with the recognition sequence indicated below.
Recognition Sequence: 5’...G T T^A A C...3’
Supplied With: 10X OPTIZYME Buffer 1
Conditions for 100% Activity: 1X OPTIZYME Buffer 1: 10mM Tris-HCl (pH 7.5 at 37°C), 10mM MgCl2 and 0.1mg/ml BSA; Incubate at 37°C.
Concentration: 10u/μl
Storage Buffer Components: 10mM Tris-HCl (pH 7.5 at 25°C), 50mM KCl, 1mM DTT, 0.1mM EDTA, 0.2mg/ml BSA and 50% (v/v) glycerol
Unit Definition: One unit is defined as the amount of HpaI required to digest 1 μg of lambda DNA in 1 hour at 37°C in 50 μl of recommended reaction buffer.
Tested for: Overdigestion Assay, Ligation/Re-Cutting Assay, Genome Qualified
Recommended storage: -20°C

OPTIZYME™ HpaII Source: Haemophilus parainfluenzae	
packaging	Mfr. No
1000 units Poly Tube	BP8032-1
Buffer 1	50-100%
Buffer 2	50-100%
Buffer 3	0-20%
Buffer 4	100%
Buffer 5	20-50%

Applications: OPTIZYME™ HpaII digests dsDNA with the recognition sequence indicated below.
Recognition Sequence: 5’...C^C G G...3’
Supplied With: 10X OPTIZYME Buffer 4
Conditions for 100% Activity: 1X OPTIZYME Buffer 4: 33mM Tris-acetate (pH 7.9 at 37°C), 10mM Mg-acetate, 66mM K-acetate and 0.1mg/ml BSA; Incubate at 37°C
Concentration: 10u/μl
Storage Buffer Components: 10mM Tris-HCl (pH 7.4 at 25°C), 100mM NaCl, 1mM DTT, 1mM EDTA, 0.2mg/ml BSA and 50% (v/v) glycerol
Unit Definition: One unit is defined as the amount of HpaII required to digest 1 μg lambda DNA in 1 hour at 37°C in 50 μl of recommended reaction buffer.
Tested for: Overdigestion Assay, Ligation/Re-Cutting Assay
Recommended storage: -20°C

OPTIZYME™ Hpy8I Source: E. coli	
packaging	Mfr. No
1000 units Poly Tube	BP8079-1
Buffer 1	50-100%
Buffer 2	50-100%
Buffer 3	0-20%
Buffer 4	100%
Buffer 5	20-50%

Applications: OPTIZYME™ Hpy8I digests dsDNA with the recognition sequence indicated below.
Recognition Sequence: 5’...G T N^N A C...3’
Supplied With: 10X OPTIZYME Buffer 4
Conditions for 100% Activity: 1X OPTIZYME Buffer 4: 33mM Tris-acetate (pH 7.9 at 37°C), 10mM Mg-acetate, 66mM K-acetate and 0.1mg/ml BSA; Incubate at 37°C
Concentration: 10u/μl
Storage Buffer Components: 10mM Tris-HCl (pH 7.4 at 25°C), 100mM KCl, 1mM DTT, 1mM EDTA, 0.2mg/ml BSA and 50% (v/v) glycerol
Unit Definition: One unit is defined as the amount of Hpy8I required to digest 1 μg of lambda DNA in 1 hour at 37°C in 50 μl of recommended reaction buffer.
Tested for: Overdigestion Assay, Ligation/Re-Cutting Assay
Recommended storage: -20°C



OPTIZYME™ KpnI
Source: Klebsiella pneumoniae OK8

packaging	Mfr. No
4000 units Poly Tube	BP8083-1

Buffer 1	20-50%
Buffer 2	0-20%
Buffer 3	0-20%
Buffer 4	20-50%
Buffer 5	0-20%
Buffer KpnI	100%

Applications: OPTIZYME™ KpnI digests dsDNA with the recognition sequence indicated below.
Recognition Sequence: 5’...G G T A C^C...3’
Supplied With: 10X OPTIZYME KpnI
Conditions for 100% Activity: 1X OPTIZYME Buffer KpnI: 10mM Tris-HCl (pH 7.5 at 37°C), 10mM MgCl2, 0.02% Triton X-100 and 0.1mg/ml BSA; Incubate at 37°C.
Concentration: 10u/μl
Storage Buffer Components: 10mM Tris-HCl (pH 7.5 at 25°C), 50mM KCl, 1mM DTT, 0.1mM EDTA, 0.2mg/ml BSA and 50% (v/v) glycerol
Unit Definition: One unit is defined as the amount of KpnI required to digest 1 μg of lambda DNA-BamHI fragments in 1 hour at 37°C in 50 μl of recommended reaction buffer.
Tested for: Overdigestion Assay, Ligation/Re-Cutting Assay, Blue/White Colony Certified, Genome Qualified
Recommended storage: -20°C

OPTIZYME™ LglI
Source: Lysobacter gummosus RFL1

packaging	Mfr. No
500 units Poly Tube	BP8067-1

Buffer 1	20-50%
Buffer 2	50-100%
Buffer 3	20-50%
Buffer 4	100%
Buffer 5	20-50%

Applications: OPTIZYME™ LglI digests dsDNA with the recognition sequence indicated below.
Recognition Sequence: 5’...G C T C T T C (N)1^...3’
Supplied With: 10X OPTIZYME Buffer 4
Conditions for 100% Activity: 1X OPTIZYME Buffer 4: 33mM Tris-acetate (pH 7.9 at 37°C), 10mM Mg-acetate, 66mM K-acetate and 0.1mg/ml BSA; Incubate at 37°C
Concentration: 5U/μl
Storage Buffer Components: 10mM Tris-HCl (pH 7.4 at 25°C), 100mM KCl, 1mM DTT, 1mM EDTA, 0.2mg/ml BSA and 50% (v/v) glycerol
Unit Definition: One unit is defined as the amount of LglI required to digest 1 μg of lambda DNA in 1 hour at 37°C in 50 μl of recommended reaction buffer.
Tested for: Overdigestion Assay, Ligation/Re-Cutting Assay, Genome Qualified
Recommended storage: -20°C

OPTIZYME™ MluI
Source: Micrococcus luteus

packaging	Mfr. No
1000 units Poly Tube	BP8021-1

Buffer 1	0-20%
Buffer 2	20-50%
Buffer 3	50-100%
Buffer 4	20-50%
Buffer 5	100%

Applications: OPTIZYME™ MluI digests dsDNA with the recognition sequence indicated below.
Recognition Sequence: 5’...A^C G C G T...3’
Supplied With: 10X OPTIZYME Buffer 5
Conditions for 100% Activity: 1X OPTIZYME Buffer 5: 10mM Tris-HCl (pH 8.5 at 37°C), 10mM MgCl2, 100mM KCl and 0.1mg/ml BSA; Incubate at 37°C
Concentration: 10u/μl
Storage Buffer Components: 10mM Tris-HCl (pH 7.5 at 25°C), 50mM KCl, 1mM DTT, 0.1mM EDTA, 0.2mg/ml BSA and 50% (v/v) glycerol
Unit Definition: One unit is defined as the amount of MluI required to digest 1 μg of lambda DNA in 1 hour at 37°C in 50 μl of recommended reaction buffer.
Tested for: Overdigestion Assay, Ligation/Re-Cutting Assay, Blue/White Colony Certified, Genome Qualified
Recommended storage: -20°C

OPTIZYME™ MspI
Source: Moraxella species

packaging	Mfr. No
3000 units Poly Tube	BP8048-1

Buffer 1	50-100%
Buffer 2	50-100%
Buffer 3	0-20%
Buffer 4	100%
Buffer 5	0-20%

Applications: OPTIZYME™ MspI digests dsDNA with the recognition sequence indicated below.
Recognition Sequence: 5’...C^C G G...3’
Supplied With: 10X OPTIZYME Buffer 4
Conditions for 100% Activity: 1X OPTIZYME Buffer 4: 33mM Tris-acetate (pH 7.9 at 37°C), 10mM Mg-acetate, 66mM K-acetate and 0.1mg/ml BSA; Incubate at 37°C
Concentration: 10u/μl
Storage Buffer Components: 10mM potassium phosphate (pH 7.5 at 25°C), 200mM NaCl, 1mM DTT, 1mM EDTA, 0.2mg/ml BSA and 50% (v/v) glycerol
Unit Definition: One unit is defined as the amount of MspI required to digest 1 μg lambda DNA in 1 hour at 37°C in 50 μl of recommended reaction buffer.
Tested for: Overdigestion Assay, Ligation/Re-Cutting Assay
Recommended storage: -20°C

OPTIZYME™ NaeI
Source: E. coli

packaging	Mfr. No
250 units Poly Tube	BP8057-1

Buffer 1	50-100%
Buffer 2	20-50%
Buffer 3	0-20%
Buffer 4	100%
Buffer 5	0-20%

Applications: OPTIZYME™ NaeI digests dsDNA with the recognition sequence indicated below.
Recognition Sequence: 5’...G C C^G G C...3’
Supplied With: 10X OPTIZYME Buffer 4
Conditions for 100% Activity: 1X OPTIZYME Buffer 4: 33mM Tris-acetate (pH 7.9 at 37°C), 10mM Mg-acetate, 66mM K-acetate and 0.1mg/ml BSA; Incubate at 37°C
Concentration: 10u/μl
Storage Buffer Components: 10mM Tris-HCl (pH 7.5 at 25°C), 500mM KCl, 1mM DTT, 0.1mM EDTA, 0.2mg/ml BSA, 0.15% Triton X-100 and 50% (v/v) glycerol
Unit Definition: One unit is defined as the amount of NaeI required to digest 1 μg of pBR322 DNA-NdeI fragments in 1 hour at 37°C in 50 μl of recommended reaction buffer.
Tested for: Overdigestion Assay, Ligation/Re-Cutting Assay, Blue/White Colony Certified, Genome Qualified
Recommended storage: -20°C

OPTIZYME™ NdeI
Source: E. coli

packaging	Mfr. No
1500 units Poly Tube	BP8020-1

Buffer 1	0-20%
Buffer 2	0-20%
Buffer 3	100%
Buffer 4	0-20%
Buffer 5	50-100%

Applications: OPTIZYME™ NdeI digests dsDNA with the recognition sequence indicated below.
Recognition Sequence: 5’...C A^T A T G...3’
Supplied With: 10X OPTIZYME Buffer 3
Conditions for 100% Activity: 1X OPTIZYME Buffer 3: 50mM Tris-HCl (pH 7.5 at 37°C), 10mM MgCl2, 100mM NaCl and 0.1mg/ml BSA; Incubate at 37°C
Concentration: 10u/μl
Storage Buffer Components: 10mM Tris-HCl (pH 7.4 at 25°C), 100mM KCl, 1mM DTT, 1mM EDTA, 0.2mg/ml BSA and 50% (v/v) glycerol
Unit Definition: One unit is defined as the amount of NdeI required to digest 1 μg of lambda DNA in 1 hour at 37°C in 50 μl of recommended reaction buffer.
Tested for: Overdigestion Assay, Ligation/Re-Cutting Assay
Recommended storage: -20°C

OPTIZYME™ NheI
Source: Neisseria mucosa heidelbergensis

packaging	Mfr. No
500 units Poly Tube	BP8019-1

Buffer 1	100%
Buffer 2	20-50%
Buffer 3	0-20%
Buffer 4	100%
Buffer 5	0-20%

Applications: OPTIZYME™ NheI digests dsDNA with the recognition sequence indicated below.
Recognition Sequence: 5’...G^C T A G C...3’
Supplied With: 10X OPTIZYME Buffer 4
Conditions for 100% Activity: 1X OPTIZYME Buffer 4: 33mM Tris-acetate (pH 7.9 at 37°C), 10mM Mg-acetate, 66mM K-acetate and 0.1mg/ml BSA; Incubate at 37°C
Concentration: 10u/μl
Storage Buffer Components: 10mM Tris-HCl (pH 8.0 at 25°C), 50mM KCl, 1mM DTT, 0.1mM EDTA, 0.2mg/ml BSA and 50% (v/v) glycerol
Unit Definition: One unit is defined as the amount of NheI required to digest 1 μg of lambda DNA-HindIII fragments in 1 hour at 37°C in 50 μl of recommended reaction buffer.
Tested for: Overdigestion Assay, Ligation/Re-Cutting Assay, Blue/White Colony Certified, Genome Qualified
Recommended storage: -20°C

OPTIZYME™ NcoI
Source: E. coli

packaging	Mfr. No
500 units Poly Tube	BP8017-1
2500 units Poly Tube	BP8017-5

Buffer 1	20-50%
Buffer 2	20-50%
Buffer 3	20-50%
Buffer 4	100%
Buffer 5	50-100%

Applications: OPTIZYME™ NcoI digests dsDNA with the recognition sequence indicated below.
Recognition Sequence: 5’...C^C A T G...3’
Supplied With: 10X OPTIZYME Buffer 4
Conditions for 100% Activity: 1X OPTIZYME Buffer 4: 33mM Tris-acetate (pH 7.9 at 37°C), 10mM Mg-acetate, 66mM K-acetate and 0.1mg/ml BSA; Incubate at 37°C
Concentration: 10u/μl
Storage Buffer Components: 10mM Tris-HCl (pH 7.5 at 25°C), 50mM KCl, 1mM DTT, 0.1mM EDTA, 0.2mg/ml BSA and 50% (v/v) glycerol
Unit Definition: One unit is defined as the amount of NcoI required to digest 1 μg of lambda DNA in 1 hour at 37°C in 50 μl of recommended reaction buffer.
Tested for: Overdigestion Assay, Ligation/Re-Cutting Assay, Blue/White Colony Certified, Genome Qualified
Recommended storage: -20°C



OPTIZYME™ NotI Source: Nocardia otitidis-caviarum	
packaging	Mfr. No
300 units Poly Tube	BP8004-1
1000 units Poly Tube	BP8004-5

Buffer 1	0-20%
Buffer 2	20-50%
Buffer 3	100%
Buffer 4	0-20%
Buffer 5	20-50%

Applications: OPTIZYME™ NotI digests dsDNA with the recognition sequence indicated below.
Recognition Sequence: 5'...G C^G G C C G C...3'
Supplied With: 10X OPTIZYME Buffer 3
Conditions for 100% Activity: 1X OPTIZYME Buffer 3: 50mM Tris-HCl (pH 7.5 at 37°C), 10mM MgCl2 , 100mM NaCl and 0.1mg/ml BSA; Incubate at 37°C.
Concentration: 10u/μl
Storage Buffer Components: 20mM Tris-HCl (pH 7.8 at 25°C), 100mM NaCl, 0.1mM EDTA, 1mM DTT, 0.02% Triton X-100, 0.2mg/ml BSA and 50% (v/v) glycerol.
Unit Definition: One unit is defined as the amount of NotI required to digest 1 μg pTZ19RJL2 DNA-BseII fragments in 1 hour at 37°C in 50 μl of recommended reaction buffer.
Tested for: Overdigestion Assay, Ligation/Re-Cutting Assay, Blue/White Colony Certified, Genome Qualified
Recommended storage: -20°C

OPTIZYME™ NsiI Source: Moraxella phenylpyruvica RFL1103	
packaging	Mfr. No
1000 units Poly Tube	BP8058-1

Buffer 1	0-20%
Buffer 2	50-100%
Buffer 3	20-50%
Buffer 4	50-100%
Buffer 5	100%

Applications: OPTIZYME™ NsiI digests dsDNA with the recognition sequence indicated below.
Recognition Sequence: 5'...A T G C A^T...3'
Supplied With: 10X OPTIZYME Buffer 5
Conditions for 100% Activity: 1X OPTIZYME Buffer 5: 10mM Tris-HCl (pH 8.5 at 37°C), 10mM MgCl2, 100mM KCl and 0.1mg/ml BSA; Incubate at 37°C
Concentration: 10u/μl
Storage Buffer Components: 10mM Tris-HCl (pH 7.4 at 25°C), 200mM KCl, 1mM DTT, 1mM EDTA, 0.2mg/ml BSA and 50% (v/v) glycerol
Unit Definition: One unit is defined as the amount of NsiI required to digest 1 μg lambda DNA in 1 hour at 37°C in 50 μl of recommended reaction buffer.
Tested for: Overdigestion Assay, Ligation/Re-Cutting Assay, Blue/White Colony Certified
Recommended storage: -20°C

OPTIZYME™ PstI Source: Pseudomonas anquilliseptica RFL1	
packaging	Mfr. No
200 units Poly Tube	BP8073-1

Buffer 1	No Reaction%
Buffer 2	No Reaction%
Buffer 3	No Reaction%
Buffer 4	No Reaction%
Buffer 5	No Reaction%
Buffer Aarl, Scal, PstI	100%

Applications: OPTIZYME™ PstI digests dsDNA with the recognition sequence indicated below.
Recognition Sequence: 5'...C C^C W G G G...3'
Supplied With: 10X OPTIZYME Buffer Aarl, Scal, PstI
Conditions for 100% Activity: 1X OPTIZYME Buffer Aarl, Scal, PstI: 10mM Bis-Tris Propane-HCl (pH 6.5 at 37°C), 10mM MgCl2 , 100mM KCl and 0.1mg/ml BSA; Incubate at 55°C.
Concentration: 10u/μl
Storage Buffer Components: 10mM Tris-HCl (pH 7.4 at 25°C), 100mM KCl, 1mM DTT, 1mM EDTA, 0.2mg/ml BSA and 50% (v/v) glycerol
Unit Definition: One unit is defined as the amount of PstI required to digest 1 μg of lambda DNA-XagI fragments in 1 hour at 55°C in 50 μl of recommended reaction buffer.
Tested for: Overdigestion Assay, Ligation/Re-Cutting Assay, Blue/White Colony Certified, Genome Qualified
Recommended storage: -20°C

OPTIZYME™ PfoI Source: Pseudomonas fluorescens biovar 126	
packaging	Mfr. No
200 units Poly Tube	BP8077-1

Buffer 1	0-20%
Buffer 2	20-50%
Buffer 3	50-100%
Buffer 4	100%
Buffer 5	0-20%

Applications: OPTIZYME™ PfoI digests dsDNA with the recognition sequence indicated below.
Recognition Sequence: 5'...T^C C N G G A...3'
Supplied With: 10X OPTIZYME Buffer 4
Conditions for 100% Activity: 1X OPTIZYME Buffer 4: 33mM Tris-acetate (pH 7.9 at 37°C), 10mM Mg-acetate, 66mM K-acetate and 0.1mg/ml BSA; Incubate at 37°C
Concentration: 10u/μl
Storage Buffer Components: 10mM Tris-HCl (pH 7.4 at 25°C), 100mM KCl, 1mM DTT, 1mM EDTA, 0.2mg/ml BSA and 50% (v/v) glycerol
Unit Definition: One unit is defined as the amount of PfoI required to digest 1 μg of lambda DNA dam-, dcm- in 1 hour at 37°C in 50 μl of recommended reaction buffer.
Tested for: Overdigestion Assay, Ligation/Re-Cutting Assay, Genome Qualified
Recommended storage: -20°C

OPTIZYME™ PstI Source: Providencia stuarti	
packaging	Mfr. No
6000 units Poly Tube	BP8001-1

Buffer 1	50-100%
Buffer 2	50-100%
Buffer 3	100%
Buffer 4	50-100%
Buffer 5	100%

Applications: OPTIZYME™ PstI digests dsDNA with the recognition sequence indicated below.
Recognition Sequence: 5'...C T G C A^G...3'
Supplied With: 10X OPTIZYME Buffer 3
Conditions for 100% Activity: 1X OPTIZYME Buffer 3: 50mM Tris-HCl (pH 7.5 at 37°C), 10mM MgCl2 , 100mM NaCl and 0.1mg/ml BSA; Incubate at 37°C.
Concentration: 10u/μl
Storage Buffer Components: 10mM Tris-HCl (pH 7.4 at 25°C), 200mM NaCl, 1mM DTT, 0.1mM EDTA, 0.15% Triton X-100, 0.2mg/ml BSA and 50% (v/v) glycerol.
Unit Definition: One unit is defined as the amount of PstI required to digest 1 μg of lambda DNA in 1 hour at 37°C in 50 μl of reaction buffer.
Tested for: Overdigestion Assay, Ligation/Re-Cutting Assay, Blue/White Colony Certified, Genome Qualified
Recommended storage: -20°C

OPTIZYME™ PvuI Source: E. coli	
packaging	Mfr. No
300 units Poly Tube	BP8050-1

Buffer 1	0-20%
Buffer 2	20-50%
Buffer 3	50-100%
Buffer 4	50-100%
Buffer 5	100%

Applications: OPTIZYME™ PvuI digests dsDNA with the recognition sequence indicated below.
Recognition Sequence: 5'...C G A T^C G...3'
Supplied With: 10X OPTIZYME Buffer 5
Conditions for 100% Activity: 1X OPTIZYME Buffer 5: 10mM Tris-HCl (pH 8.5 at 37°C), 10mM MgCl2, 100mM KCl and 0.1mg/ml BSA; Incubate at 37°C
Concentration: 10u/μl
Storage Buffer Components: 10mM Tris-HCl (pH 7.5 at 25°C), 300mM KCl, 0.1mM EDTA, 1mM DTT, 0.2mg/ml BSA and 50% (v/v) glycerol
Unit Definition: One unit is defined as the amount of PvuI required to digest 1 μg lambda DNA-CpoI fragments in 1 hour at 37°C in 50 μl of recommended reaction buffer.
Tested for: Overdigestion Assay, Ligation/Re-Cutting Assay, Genome Qualified
Recommended storage: -20°C

OPTIZYME™ PvuII Source: Proteus vulgaris	
packaging	Mfr. No
2500 units Poly Tube	BP8022-1
5000 units Poly Tube	BP8022-5

Buffer 1	50-100%
Buffer 2	100%
Buffer 3	20-50%
Buffer 4	20-50%
Buffer 5	50-100%
Tested for	RNase activity, protease activity, and specific performance tests.

Applications: OPTIZYME™ PvuII digests dsDNA with the recognition sequence indicated below.
Recognition Sequence: 5'...C A G^C T G...3'
Supplied With: 10X OPTIZYME Buffer 2
Conditions for 100% Activity: 1X OPTIZYME Buffer 2: 10mM Tris-HCl (pH 7.5 at 37°C), 10mM MgCl2, 50mM NaCl and 0.1mg/ml BSA; Incubate at 37°C.
Concentration: 10U/μl
Storage Buffer Components: 10mM Tris-HCl (pH 7.4 at 25°C), 100mM KCl, 1mM DTT, 1mM EDTA, 0.2mg/ml BSA and 50% (v/v) glycerol
Unit Definition: One unit is defined as the amount of PvuII required to digest 1 μg of lambda DNA in 1 hour at 37°C in 50 μl of recommended reaction buffer.
Tested for: Overdigestion Assay, Ligation/Re-Cutting Assay
Recommended storage: -20°C

OPTIZYME™ RsaI Source: Rhodopseudomonas sphaeroides	
packaging	Mfr. No
1000 units Poly Tube	BP8000-1

Buffer 1	50-100%
Buffer 2	20-50%
Buffer 3	0-20%
Buffer 4	100%
Buffer 5	0-20%

Applications: OPTIZYME™ RsaI digests dsDNA with the recognition sequence indicated below.
Recognition Sequence: 5'...G T^A C...3'
Supplied With: 10X OPTIZYME Buffer 4
Conditions for 100% Activity: 1X OPTIZYME Buffer 4: 33mM Tris-acetate (pH 7.9 at 37°C), 10mM Mg-acetate, 66mM K-acetate and 0.1mg/ml BSA; Incubate at 37°C
Concentration: 10U/μl
Storage Buffer Components: 10mM Tris-HCl (pH 7.4 at 25°C), 100mM KCl, 1mM DTT, 1mM EDTA, 0.2mg/ml BSA and 50% (v/v) glycerol
Unit Definition: One unit is defined as the amount of RsaI required to digest 1 μg of lambda DNA in 1 hour at 37°C in 50 μl of recommended reaction buffer.
Tested for: Overdigestion Assay, Ligation/Re-Cutting Assay
Recommended storage: -20°C



OPTIZYME™ SacI
Source: Streptomyces achromogenes

packaging	Mfr. No
2000 units Poly Tube	BP8016-1

Buffer 1	50-100%
Buffer 2	20-50%
Buffer 3	0-20%
Buffer 4	50-100%
Buffer 5	0-20%

Applications: OPTIZYME™ SacI digests dsDNA with the recognition sequence indicated below.
Recognition Sequence: 5’...G A G C T^C...3’
Supplied With: 10X OPTIZYME Ecl136II, SacI Buffer
Conditions for 100% Activity: 1X OPTIZYME Buffer Ecl136II, SacI: 10mM Bis-Tris Propane-HCl (pH 6.5 at 37°C), 10mM MgCl2 and 0.1mg/ml BSA; Incubate at 37°C.
Concentration: 10u/μl
Storage Buffer Components: 10mM Tris-HCl (pH 7.4 at 25°C), 100mM KCl, 1mM DTT, 1 mM EDTA, 0.2mg/ml BSA and 50% (v/v) glycerol
Unit Definition: One unit is defined as the amount of SacI required to digest 1 μg of lambda DNA-HindIII fragments in 1 hour at 37°C in 50 μl of recommended reaction buffer.
Tested for: Overdigestion Assay, Ligation/Re-Cutting Assay, Blue/White Colony Certified, Genome Qualified
Recommended storage: -20°C

OPTIZYME™ SacII
Source: E. coli

packaging	Mfr. No
1200 units Poly Tube	BP8023-1

Buffer 1	100%
Buffer 2	50-100%
Buffer 3	0-20%
Buffer 4	50-100%
Buffer 5	0-20%

Applications: OPTIZYME™ SacII digests dsDNA with the recognition sequence indicated below.
Recognition Sequence: 5’...C C G C^G G...3’
Supplied With: 10X OPTIZYME Buffer 1
Conditions for 100% Activity: 1X OPTIZYME Buffer 1: 10mM Tris-HCl (pH 7.5 at 37°C), 10mM MgCl2 and 0.1mg/ml BSA; Incubate at 37°C.
Concentration: 10u/μl
Storage Buffer Components: 10mM potassium phosphate (pH 7.4 at 25°C), 100mM NaCl, 1mM EDTA, 7mM 2-mercaptoethanol, 0.2mg/ml BSA and 50% (v/v) glycerol
Unit Definition: One unit is defined as the amount of SacII required to digest 1 μg of lambda DNA in 1 hour at 37°C in 50 μl of recommended reaction buffer.
Tested for: Overdigestion Assay, Ligation/Re-Cutting Assay, Blue/White Colony Certified, Genome Qualified
Recommended storage: -20°C

OPTIZYME™ SalI
Source: Streptomyces albus G

packaging	Mfr. No
2000 units Poly Tube	BP8013-1

Buffer 1	0-20%
Buffer 2	0-20%
Buffer 3	100%
Buffer 4	0-20%
Buffer 5	20-50%

Applications: OPTIZYME™ SalI digests dsDNA with the recognition sequence indicated below.
Recognition Sequence: 5’...G^T C G A C...3’
Supplied With: 10X OPTIZYME Buffer 3
Conditions for 100% Activity: 1X OPTIZYME Buffer 3: 50mM Tris-HCl (pH 7.5 at 37°C), 10mM MgCl2, 100mM NaCl and 0.1mg/ml BSA; Incubate at 37°C
Concentration: 10u/μl
Storage Buffer Components: 10mM Tris-HCl (pH 7.4 at 25°C),100 mM KCl, 1mM DTT, 0.1mM EDTA, 0.2mg/ml BSA and 50% (v/v) glycerol
Unit Definition: One unit is defined as the amount of SalI required to digest 1 μg DNA-Eco8II fragments in 1 hour at 37°C in 50 μl of recommended reaction buffer.
Tested for: Overdigestion Assay, Ligation/Re-Cutting Assay, Blue/White Colony Certified, Genome Qualified
Recommended storage: -20°C

OPTIZYME™ Sau3AI
Source: Bacillus species RFL143

packaging	Mfr. No
300 units Poly Tube	BP8030-1

Buffer 1	20-50%
Buffer 2	20-50%
Buffer 3	0-20%
Buffer 4	50-100%
Buffer 5	0-20%

Applications: OPTIZYME™ Sau3AI digests dsDNA with the recognition sequence indicated below.
Recognition Sequence: 5’...^G A T C ...3’
Supplied With: 10X OPTIZYME Buffer Sau3AI
Conditions for 100% Activity: 1X OPTIZYME Buffer Sau3AI: 33mM Tris-acetate (pH 7.9 at 37°C), 10mM Mg-acetate, 66mM K-acetate, 0.02% Triton X-100 and 0.1mg/ml BSA; Incubate at 37°C
Concentration: 10u/μl
Storage Buffer Components: 10mM Tris-HCl (pH 7.4 at 25°C), 10mM KCl, 1mM DTT, 1mM EDTA, 0.2mg/ml BSA and 50% (v/v) glycerol
Unit Definition: One unit is defined as the amount of Sau3AI required to digest 1 μg of lambda DNA in 1 hour at 37°C in 50 μl of recommended reaction buffer.
Tested for: Overdigestion Assay, Ligation/Re-Cutting Assay
Recommended storage: -20°C

OPTIZYME™ Scal
Source: Streptomyces caespitosus

packaging	Mfr. No
1000 units Poly Tube	BP8037-1

Buffer 1	0-20%
Buffer 2	0-20%
Buffer 3	0-20%
Buffer 4	0-20%
Buffer 5	0-20%

Applications: OPTIZYME™ Scal digests dsDNA with the recognition sequence indicated below.
Recognition Sequence: 5’...A G T^A C T...3’
Supplied With: 10X OPTIZYME Aarl, Scal, Pasi Buffer
Conditions for 100% Activity: 1X OPTIZYME Buffer Aarl, Scal, Pasi: 10mM Bis-Tris Propane-HCl (pH 6.5 at 37°C), 10mM MgCl2 , 100mM KCl and 0.1mg/ml BSA; Incubate at 37°C.
Concentration: 10u/μl
Storage Buffer Components: 10mM Tris-HCl (pH 7.5 at 25°C), 50mM KCl, 1mM DTT, 0.1mM EDTA, 0.5mg/ml BSA and 50% (v/v) glycerol
Unit Definition: One unit is defined as the amount of Scal required to digest 1 μg of lambda DNA in 1 hour at 37°C in 50 μl of recommended reaction buffer.
Tested for: Overdigestion Assay, Ligation/Re-Cutting Assay
Recommended storage: -20°C

OPTIZYME™ SfaAI
Source: E. coli

packaging	Mfr. No
1000 units Poly Tube	BP8076-1

Buffer 1	50-100%
Buffer 2	0-20%
Buffer 3	0-20%
Buffer 4	100%
Buffer 5	0-20%

Applications: OPTIZYME™ SfaAI digests dsDNA with the recognition sequence indicated below.
Recognition Sequence: 5’...G C G A T^C G C...3’
Supplied With: 10X OPTIZYME Buffer 4
Conditions for 100% Activity: 1X OPTIZYME Buffer 4: 33mM Tris-acetate (pH 7.9 at 37°C), 10mM Mg-acetate, 66mM K-acetate and 0.1mg/ml BSA; Incubate at 37°C
Concentration: 10u/μl
Storage Buffer Components: 10mM Tris-HCl (pH 7.4 at 25°C), 100mM KCl, 1mM DTT, 1mM EDTA, 0.2mg/ml BSA and 50% (v/v) glycerol
Unit Definition: One unit is defined as the amount of SfaAI required to digest 1 μg of control DNA in 1 hour at 37°C in 50 μl of recommended reaction buffer. The control DNA is linearized pJET1 DNA with inserted SfaAI recognition site.
Tested for: Overdigestion Assay, Ligation/Re-Cutting Assay, Blue/White Colony Certified, Genome Qualified
Recommended storage: -20°C

OPTIZYME™ SmaI
Source: Serratia marcescens

packaging	Mfr. No
1200 units Poly Tube	BP8011-1
5000 units Poly Tube	BP8011-5

Buffer 1	50-100%
Buffer 2	0-20%
Buffer 3	0-20%
Buffer 4	100%
Buffer 5	0-20%

Applications: OPTIZYME™ SmaI digests dsDNA with the recognition sequence indicated below.
Recognition Sequence: 5’...C C C^G G G...3’
Supplied With: 10X OPTIZYME Buffer 4
Conditions for 100% Activity: 10mM Tris-HCl (pH 7.4 at 25°C), 100mM KCl, 1mM DTT, 1mM EDTA, 0.2mg/ml BSA and 50% (v/v) glycerol; Incubate at 30°C.
Concentration: 10u/μl
Storage Buffer Components: 10mM Tris-HCl (pH 7.4 at 25°C), 100mM KCl, 1mM DTT, 1mM EDTA, 0.2mg/ml BSA and 50% (v/v) glycerol
Unit Definition: One unit is defined as the amount of SmaI required to digest 1 μg of lambda DNA-Eco8II fragments in 1 hour at 30°C in 50μl of recommended reaction buffer.
Tested for: Overdigestion Assay, Ligation/Re-Cutting Assay, Blue/White Colony Certified, Genome Qualified
Recommended storage: -20°C

OPTIZYME™ SpeI
Source: Bacillus coagulans VS 29-022

packaging	Mfr. No
400 units Poly Tube	BP8018-1

Buffer 1	50-100%
Buffer 2	50-100%
Buffer 3	0-20%
Buffer 4	100%
Buffer 5	20-50%

Applications: OPTIZYME™ SpeI digests dsDNA with the recognition sequence indicated below.
Recognition Sequence: 5’...A^C T A G T...3’
Supplied With: 10X OPTIZYME Buffer 4
Conditions for 100% Activity: 1X OPTIZYME Buffer 4: 33mM Tris-acetate (pH 7.9 at 37°C), 10mM Mg-acetate, 66mM K-acetate and 0.1mg/ml BSA; Incubate at 37°C
Concentration: 10u/μl
Storage Buffer Components: 10mM Tris-HCl (pH 7.5 at 25°C), 300mM KCl, 1mM DTT, 0.1mM EDTA, 0.2mg/ml BSA and 50% (v/v) glycerol
Unit Definition: One unit is defined as the amount of SpeI required to digest 1 μg of control DNA in 1 hour at 37°C in 50 μl of recommended reaction buffer. The control DNA is pUC19 DNA with inserted SpeI recognition site-Psp1406I fragments.
Tested for: Overdigestion Assay, Ligation/Re-Cutting Assay, Blue/White Colony Certified, Genome Qualified
Recommended storage: -20°C



OPTIZYME™ SphI
Source: Pseudomonas aeruginosa

packaging	Mfr. No
500 units Poly Tube	BP8029-1

Buffer 1	100%
Buffer 2	50-100%
Buffer 3	0-20%
Buffer 4	50-100%
Buffer 5	0-20%

Applications: OPTIZYME™ SphI digests dsDNA with the recognition sequence indicated below.
Recognition Sequence: 5’...G C A T G^C...3’
Supplied With: 10X OPTIZYME Buffer 1
Conditions for 100% Activity: 1X OPTIZYME Buffer 1: 10mM Tris-HCl (pH 7.5 at 37°C), 10mM MgCl2 and 0.1mg/ml BSA; Incubate at 37°C.
Concentration: 10u/μl
Storage Buffer Components: 10mM potassium-phosphate (pH 7.0 at 25°C), 50mM KCl, 1mM DTT, 0.1mM EDTA, 0.15% Triton X-100, 0.5mg/ml BSA and 50% (v/v) glycerol
Unit Definition: One unit is defined as the amount of SphI required to digest 1 μg of lambda DNA in 1 hour at 37°C in 50 μl of recommended reaction buffer.
Tested for: Overdigestion Assay, Ligation/Re-Cutting Assay, Blue/White Colony Certified, Genome Qualified
Recommended storage: -20°C

OPTIZYME™ StuI
Source: E. coli

packaging	Mfr. No
1000 units Poly Tube	BP8027-1

Buffer 1	100%
Buffer 2	50-100%
Buffer 3	20-50%
Buffer 4	50-100%
Buffer 5	20-50%

Applications: OPTIZYME™ StuI digests dsDNA with the recognition sequence indicated below.
Recognition Sequence: 5’...A G G^C C T...3’
Supplied With: 10X OPTIZYME Buffer 1
Conditions for 100% Activity: 1X OPTIZYME Buffer 1: 10mM Tris-HCl (pH 7.5 at 37°C), 10mM MgCl2 and 0.1mg/ml BSA; Incubate at 37°C.
Concentration: 10u/μl
Storage Buffer Components: 10mM Tris-HCl (pH 7.4 at 25°C), 100mM KCl, 1mM DTT, 1mM EDTA, 0.2mg/ml BSA and 50% (v/v) glycerol
Unit Definition: One unit is defined as the amount of StuI required to digest 1 μg of lambda DNA-Eco81I in 1 hour at 37°C in 50 μl of recommended reaction buffer.
Tested for: Overdigestion Assay, Ligation/Re-Cutting Assay, Blue/White Colony Certified, Genome Qualified
Recommended storage: -20°C

OPTIZYME™ TaqI
Source: E. coli

packaging	Mfr. No
3000 units Poly Tube	BP8007-1

Buffer 1	0-20%
Buffer 2	20-50%
Buffer 3	20-50%
Buffer 4	20-50%
Buffer 5	20-50%
Buffer TaqI	100%

Applications: OPTIZYME™ TaqI digests dsDNA with the recognition sequence indicated below.
Recognition Sequence: 5’...T^C G A...3’
Supplied With: 10X OPTIZYME TaqI Buffer
Conditions for 100% Activity: 1X OPTIZYME Buffer TaqI: 10mM Tris-HCl (pH 8.0 at 37°C), 5mM MgCl2, 100mM NaCl and 0.1mg/ml BSA; Incubate at 65°C; To ensure higher efficiency digestion, perform the cleavage reaction under paraffin oil.
Concentration: 10u/μl
Storage Buffer Components: 10mM Tris-HCl (pH 7.5 at 25°C), 300mM KCl, 1mM DTT, 0.1mM EDTA, 0.5mg/ml BSA and 50% (v/v) glycerol
Unit Definition: One unit is defined as the amount of TaqI required to digest 1 μg of lambda DNA dam - in 1 hour at 65°C in 50 μl of recommended reaction buffer.
Tested for: Overdigestion Assay, Ligation/Re-Cutting Assay
Recommended storage: -20°C

OPTIZYME™ Tru9I
Source: Thermus ruber RFL1

packaging	Mfr. No
300 units Poly Tube	BP8064-1

Buffer 1	50-100%
Buffer 2	50-100%
Buffer 3	20-50%
Buffer 4	50-100%
Buffer 5	100%

Applications: OPTIZYME™ Tru9I digests dsDNA with the recognition sequence indicated below.
Recognition Sequence: 5’...T^T A A...3’
Supplied With: 10X OPTIZYME Buffer 5
Conditions for 100% Activity: 1X OPTIZYME Buffer 5: 10mM Tris-HCl (pH 8.5 at 37°C), 10mM MgCl2, 100mM KCl and 0.1mg/ml BSA; Incubate at 65°C; To ensure higher efficiency of digestion, perform the cleavage reaction under paraffin oil.
Concentration: 10u/μl
Storage Buffer Components: 10mM Tris-HCl (pH 7.4 at 25°C), 100mM KCl, 1mM DTT, 1mM EDTA, 0.2mg/ml BSA and 50% (v/v) glycerol
Unit Definition: One unit is defined as the amount of TruI required to digest 1 μg of lambda DNA in 1 hour at 65°C in 50 μl of recommended reaction buffer.
Tested for: Overdigestion Assay, Ligation/Re-Cutting Assay
Recommended storage: -20°C

OPTIZYME™ VspI
Source: Vibrio species

packaging	Mfr. No
1000 units Poly Tube	BP8055-1

Buffer 1	0-20%
Buffer 2	50-100%
Buffer 3	100%
Buffer 4	100%
Buffer 5	20-50%

Applications: OPTIZYME™ VspI digests dsDNA with the recognition sequence indicated below.
Recognition Sequence: 5’...A T^T A A T...3’
Supplied With: 10X OPTIZYME Buffer 3
Conditions for 100% Activity: 1X OPTIZYME Buffer 3: 50mM Tris-HCl (pH 7.5 at 37°C), 10mM MgCl2, 100mM NaCl and 0.1mg/ml BSA; Incubate at 37°C
Concentration: 10u/μl
Storage Buffer Components: 10mM Tris-HCl (pH 7.4 at 25°C), 100mM KCl, 1mM DTT, 1mM EDTA, 0.2mg/ml BSA and 50% (v/v) glycerol
Unit Definition: One unit is defined as the amount of VspI required to digest 1 μg of lambda DNA in 1 hour at 37°C in 50 μl of recommended reaction buffer.
Tested for: Overdigestion Assay, Ligation/Re-Cutting Assay, Genome Qualified
Recommended storage: -20°C

OPTIZYME™ XbaI
Source: Xanthomonas badrii

packaging	Mfr. No
2000 units Poly Tube	BP8008-1
10000 units Poly Tube	BP8008-5

Buffer 1	50-100%
Buffer 2	50-100%
Buffer 3	20-50%
Buffer 4	100%
Buffer 5	0-20%

Applications: OPTIZYME™ XbaI digests dsDNA with the recognition sequence indicated below.
Recognition Sequence: 5’...T^C T A G A...3’
Supplied With: 10X OPTIZYME Buffer 4
Conditions for 100% Activity: 1X OPTIZYME Buffer 4: 33mM Tris-acetate (pH 7.9 at 37°C), 10mM Mg-acetate, 66mM K-acetate and 0.1mg/ml BSA; Incubate at 37°C
Concentration: 10u/μl
Storage Buffer Components: 10mM Tris-HCl (pH 7.4 at 25°C), 100mM KCl, 7mM 2-mercaptoethanol, 1mM EDTA, 0.2mg/ml BSA and 50% (v/v) glycerol
Unit Definition: One unit is defined as the amount of XbaI required to digest 1 μg of lambda DNA dam --SmaI fragments in 1 hour at 37°C in 50 μl of recommended reaction buffer.
Tested for: Overdigestion Assay, Ligation/Re-Cutting Assay, Blue/White Colony Certified, Genome Qualified
Recommended storage: -20°C

OPTIZYME™ XhoI
Source: Xanthomonas holcicola

packaging	Mfr. No
3000 units Poly Tube	BP8010-1
10000 units Poly Tube	BP8010-5

Buffer 1	0-20%
Buffer 2	50-100%
Buffer 3	50-100%
Buffer 4	20-50%
Buffer 5	100%

Applications: OPTIZYME™ XhoI digests dsDNA with the recognition sequence indicated below.
Recognition Sequence: 5’....C^T C G A G...3’
Supplied With: 10X OPTIZYME Buffer 5
Conditions for 100% Activity: 1X OPTIZYME Buffer 5: 10mM Tris-HCl (pH 8.5 at 37°C), 10mM MgCl2 , 100mM KCl and 0.1mg/ml BSA; Incubate at 37°C.
Concentration: 10u/μl
Storage Buffer Components: 10mM Tris-HCl (pH 7.4 at 25°C), 100mM KCl, 1mM DTT, 1mM EDTA, 0.2mg/ml BSA and 50% (v/v) glycerol
Unit Definition: One unit is defined as the amount of XhoI required to digest 1 μg of lambda DNA-HindIII fragments in 1 hour at 37°C in 50 μl of recommended reaction buffer.
Tested for: Overdigestion Assay, Ligation/Re-Cutting Assay, Blue/White Colony Certified, Genome Qualified
Recommended storage: -20°C

OPTIZYME™ XmaI
Source: Xanthomonas maltophilia Jo 85-025

packaging	Mfr. No
200 units Poly Tube	BP8082-1

Buffer 1	20-50%
Buffer 2	50-100%
Buffer 3	50-100%
Buffer 4	100%
Buffer 5	50-100%

Applications: OPTIZYME™ XmaI digests dsDNA with the recognition sequence indicated below.
Recognition Sequence: 5’...C^C T A G G...3’
Supplied With: 10X OPTIZYME Buffer 4
Conditions for 100% Activity: 1X OPTIZYME Buffer 4: 33mM Tris-acetate (pH 7.9 at 37°C), 10mM Mg-acetate, 66mM K-acetate and 0.1mg/ml BSA; Incubate at 37°C
Concentration: 10u/μl
Storage Buffer Components: 10mM Tris-HCl (pH 7.4 at 25°C), 100mM KCl, 1mM DTT, 1mM EDTA, 0.2mg/ml BSA and 50% (v/v) glycerol
Unit Definition: One unit is defined as the amount of XmaI required to digest 1 μg of lambda DNA-SmaI fragments in 1 hour at 37°C in 50 μl of recommended reaction buffer.
Tested for: Overdigestion Assay, Ligation/Re-Cutting Assay, Blue/White Colony Certified, Genome Qualified
Recommended storage: -20°C

OPTIZYME™ XmnI Source: E. coli	
packaging	Mfr. No
500 units Poly Tube	BP8052-1
Buffer 1	20-50%
Buffer 2	50-100%
Buffer 3	0-20%
Buffer 4	100%
Buffer 5	0-20%

Applications: OPTIZYME™ XmnI digests dsDNA with the recognition sequence indicated below.
Recognition Sequence: 5’...G A A N N^N N T T C...3’
Supplied With: 10X OPTIZYME Buffer 4
Conditions for 100% Activity: 1X OPTIZYME Buffer 4: 33mM Tris-acetate (pH 7.9 at 37°C), 10mM Mg-acetate, 66mM K-acetate and 0.1mg/ml BSA; Incubate at 37°C
Concentration: 10u/μl
Storage Buffer Components: 10mM Tris-HCl (pH 7.5 at 25°C), 100mM KCl, 1mM DTT, 5mM MgCl2, 0.2mg/ml BSA and 50% (v/v) glycerol
Unit Definition: One unit is defined as the amount of XmnI required to digest 1 μg of lambda DNA in 1 hour at 37°C in 50 μl of recommended reaction buffer.
Tested for: Overdigestion Assay, Ligation/Re-Cutting Assay, Genome Qualified
Recommended storage: -20°C

OPTIZYME™ 10X Buffer 1	
packaging	Mfr. No
5X1mℓ Poly Tube	BP8084-1

Applications: The OPTIZYME™ Five Buffer System ensures the optimum reaction conditions for each restriction enzyme. This system consists of 10X Buffer 1, 10X Buffer 2, 10X Buffer 3, 10X Buffer 4 and 10X Buffer 5. The recommended buffer is supplied with each enzyme. To ensure consistent enzyme performance, OPTIZYME™ restriction enzyme buffers contain BSA, which enhances the stability of many enzymes and binds contaminants that may be present in DNA preparations. Multiple freeze-thaw cycles of the buffers will not cause BSA precipitation. OPTIZYME™ restriction enzymes exhibit 100% of their certified activity in the recommended buffer.
Components: 100mM Tris-HCl (pH 7.5 at 37°C), 100mM MgCl2, 1 mg/ml BSA
Tested for: Functional Test, Endo- and Exonucleases Assay, Nicking Activity Assay
Recommended storage: -20°C

OPTIZYME™ 10X Buffer 2	
packaging	Mfr. No
5X1mℓ Poly Tube	BP8085-1

Applications: The OPTIZYME™ Five Buffer System ensures the optimum reaction conditions for each restriction enzyme. This system consists of 10X Buffer 1, 10X Buffer 2, 10X Buffer 3, 10X Buffer 4 and 10X Buffer 5. The recommended buffer is supplied with each enzyme. To ensure consistent enzyme performance, OPTIZYME™ restriction enzyme buffers contain BSA, which enhances the stability of many enzymes and binds contaminants that may be present in DNA preparations. Multiple freeze-thaw cycles of the buffers will not cause BSA precipitation. OPTIZYME™ restriction enzymes exhibit 100% of their certified activity in the recommended buffer.
Components: 100mM Tris-HCl (pH 7.5 at 37°C), 100mM MgCl2, 500mM NaCl, 1mg/ml BSA
Tested for: Functional Test, Endo- and Exonucleases Assay, Nicking Activity Assay
Recommended storage: -20°C

OPTIZYME™ 10X Buffer 3	
packaging	Mfr. No
5X1mℓ Poly Tube	BP8086-1

Applications: The OPTIZYME™ Five Buffer System ensures the optimum reaction conditions for each restriction enzyme. This system consists of 10X Buffer 1, 10X Buffer 2, 10X Buffer 3, 10X Buffer 4 and 10X Buffer 5. The recommended buffer is supplied with each enzyme. To ensure consistent enzyme performance, OPTIZYME™ restriction enzyme buffers contain BSA, which enhances the stability of many enzymes and binds contaminants that may be present in DNA preparations. Multiple freeze-thaw cycles of the buffers will not cause BSA precipitation. OPTIZYME™ restriction enzymes exhibit 100% of their certified activity in the recommended buffer.
Components: 500mM Tris-HCl (pH 7.5 at 37°C), 100mM MgCl2, 1M NaCl, 1mg/ml BSA
Tested for: Functional Test, Endo- and Exonucleases Assay, Nicking Activity Assay
Recommended storage: -20°C

OPTIZYME™ 10X Buffer 4	
packaging	Mfr. No
5X1mℓ Poly Tube	BP8087-1

Applications: The OPTIZYME™ Five Buffer System ensures the optimum reaction conditions for each restriction enzyme. This system consists of 10X Buffer 1, 10X Buffer 2, 10X Buffer 3, 10X Buffer 4 and 10X Buffer 5. The recommended buffer is supplied with each enzyme. To ensure consistent enzyme performance, OPTIZYME™ restriction enzyme buffers contain BSA, which enhances the stability of many enzymes and binds contaminants that may be present in DNA preparations. Multiple freeze-thaw cycles of the buffers will not cause BSA precipitation. OPTIZYME™ restriction enzymes exhibit 100% of their certified activity in the recommended buffer.
Components: 330mM Tris-acetate (pH 7.9 at 37°C), 100mM magnesium acetate, 660mM potassium acetate, 1mg/ml BSA
Tested for: Functional Test, Endo- and Exonucleases Assay, Nicking Activity Assay
Recommended storage: -20°C

OPTIZYME™ 10X Buffer 5	
packaging	Mfr. No
5X1mℓ Poly Tube	BP8088-1

Applications: The OPTIZYME™ Five Buffer System ensures the optimum reaction conditions for each restriction enzyme. This system consists of 10X Buffer 1, 10X Buffer 2, 10X Buffer 3, 10X Buffer 4 and 10X Buffer 5. The recommended buffer is supplied with each enzyme. To ensure consistent enzyme performance, OPTIZYME™ restriction enzyme buffers contain BSA, which enhances the stability of many enzymes and binds contaminants that may be present in DNA preparations. Multiple freeze-thaw cycles of the buffers will not cause BSA precipitation. OPTIZYME™ restriction enzymes exhibit 100% of their certified activity in the recommended buffer.
Components: 100mM Tris HCl (pH 8.5 at 37°C), 100mM MgCl2, 1M KCl , 1mg/ml BSA
Tested for: Functional Test, Endo- and Exonucleases Assay, Nicking Activity Assay
Recommended storage: -20°C

OPTIZYME™ 10X Buffer AarI, Scal, PstI	
packaging	Mfr. No
1 mℓ Poly Tube	BP8096-1

Applications: The OPTIZYME™ Five Buffer System ensures the optimum reaction conditions for each restriction enzyme. This system consists of 10X Buffer 1, 10X Buffer 2, 10X Buffer 3, 10X Buffer 4 and 10X Buffer 5. The recommended buffer is supplied with each enzyme. To ensure consistent enzyme performance, OPTIZYME™ restriction enzyme buffers contain BSA, which enhances the stability of many enzymes and binds contaminants that may be present in DNA preparations. Multiple freeze-thaw cycles of the buffers will not cause BSA precipitation. OPTIZYME™ restriction enzymes exhibit 100% of their certified activity in the recommended buffer.
Components: 100mM Bis-Tris Propane-HCl (pH 6.5 at 37°C), 100mM MgCl2, 1M KCl, 1mg/ml BSA
Tested for: Functional Test, Endo- and Exonucleases Assay, Nicking Activity Assay
Recommended storage: -20°C

OPTIZYME™ 10X Buffer BamHI	
packaging	Mfr. No
5X1mℓ Poly Tube	BP8089-1

Applications: The OPTIZYME™ Five Buffer System ensures the optimum reaction conditions for each restriction enzyme. This system consists of 10X Buffer 1, 10X Buffer 2, 10X Buffer 3, 10X Buffer 4 and 10X Buffer 5. The recommended buffer is supplied with each enzyme. To ensure consistent enzyme performance, OPTIZYME™ restriction enzyme buffers contain BSA, which enhances the stability of many enzymes and binds contaminants that may be present in DNA preparations. Multiple freeze-thaw cycles of the buffers will not cause BSA precipitation. OPTIZYME™ restriction enzymes exhibit 100% of their certified activity in the recommended buffer.
Components: 100mM Tris-HCl (pH 8.0 at 37°C), 50mM MgCl2, 1M KCl , 0.2 % Triton X-100, 1mg/ml BSA
Tested for: Functional Test, Endo- and Exonucleases Assay, Nicking Activity Assay
Recommended storage: -20°C

OPTIZYME™ 10X Buffer Cfr9I	
packaging	Mfr. No
1 mℓ Poly Tube	BP8093-1

Applications: The OPTIZYME™ Five Buffer System ensures the optimum reaction conditions for each restriction enzyme. This system consists of 10X Buffer 1, 10X Buffer 2, 10X Buffer 3, 10X Buffer 4 and 10X Buffer 5. The recommended buffer is supplied with each enzyme. To ensure consistent enzyme performance, OPTIZYME™ restriction enzyme buffers contain BSA, which enhances the stability of many enzymes and binds contaminants that may be present in DNA preparations. Multiple freeze-thaw cycles of the buffers will not cause BSA precipitation. OPTIZYME™ restriction enzymes exhibit 100% of their certified activity in the recommended buffer.
Components: 100mM Tris-HCl (pH 7.2 at 37°C), 50mM MgCl2, 2M sodium glutamate, 1mg/ml BSA
Tested for: Functional Test, Endo- and Exonucleases Assay, Nicking Activity Assay
Recommended storage: -20°C

OPTIZYME™ 10X Buffer Ecl136II, SacI	
packaging	Mfr. No
1 mℓ Poly Tube	BP8095-1

Applications: The OPTIZYME™ Five Buffer System ensures the optimum reaction conditions for each restriction enzyme. This system consists of 10X Buffer 1, 10X Buffer 2, 10X Buffer 3, 10X Buffer 4 and 10X Buffer 5. The recommended buffer is supplied with each enzyme. To ensure consistent enzyme performance, OPTIZYME™ restriction enzyme buffers contain BSA, which enhances the stability of many enzymes and binds contaminants that may be present in DNA preparations. Multiple freeze-thaw cycles of the buffers will not cause BSA precipitation. OPTIZYME™ restriction enzymes exhibit 100% of their certified activity in the recommended buffer.
Components: 100mM Bis-Tris Propane-HCl (pH 6.5 at 37°C), 100mM MgCl2, 1mg/ml BSA
Tested for: Functional Test, Endo- and Exonucleases Assay, Nicking Activity Assay
Recommended storage: -20°C

OPTIZYME™ 10X Buffer EcoRI	
packaging	Mfr. No
5X1mℓ Poly Tube	BP8090-1

Applications: The OPTIZYME™ Five Buffer System ensures the optimum reaction conditions for each restriction enzyme. This system consists of 10X Buffer 1, 10X Buffer 2, 10X Buffer 3, 10X Buffer 4 and 10X Buffer 5. The recommended buffer is supplied with each enzyme. To ensure consistent enzyme performance, OPTIZYME™ restriction enzyme buffers contain BSA, which enhances the stability of many enzymes and binds contaminants that may be present in DNA preparations. Multiple freeze-thaw cycles of the buffers will not cause BSA precipitation. OPTIZYME™ restriction enzymes exhibit 100% of their certified activity in the recommended buffer.
Components: 500mM Tris-HCl (pH 7.5 at 37°C), 100mM MgCl2, 1M NaCl, 0.2 % Triton X-100, 1mg/ml BSA
Tested for: Functional Test, Endo- and Exonucleases Assay, Nicking Activity Assay
Recommended storage: -20°C

OPTIZYME™ 10X Buffer KpnI	
packaging	Mfr. No
1 mℓ Poly Tube	BP8094-1

Applications: The OPTIZYME™ Five Buffer System ensures the optimum reaction conditions for each restriction enzyme. This system consists of 10X Buffer 1, 10X Buffer 2, 10X Buffer 3, 10X Buffer 4 and 10X Buffer 5. The recommended buffer is supplied with each enzyme. To ensure consistent enzyme performance, OPTIZYME™ restriction enzyme buffers contain BSA, which enhances the stability of many enzymes and binds contaminants that may be present in DNA preparations. Multiple freeze-thaw cycles of the buffers will not cause BSA precipitation. OPTIZYME™ restriction enzymes exhibit 100% of their certified activity in the recommended buffer.
Components: 100mM Tris-HCl (pH 7.5 at 37°C), 100mM MgCl2, 0.2% Triton X-100, 1mg/ml BSA
Tested for: Functional Test, Endo- and Exonucleases Assay, Nicking Activity Assay
Recommended storage: -20°C



OPTIZYME™ 10X Buffer Sau3AI

packaging		Mfr. No
1 mℓ	Poly Tube	BP8091-1

Applications: The OPTIZYME™ Five Buffer System ensures the optimum reaction conditions for each restriction enzyme. This system consists of 10X Buffer 1, 10X Buffer 2, 10X Buffer 3, 10X Buffer 4 and 10X Buffer 5. The recommended buffer is supplied with each enzyme. To ensure consistent enzyme performance, OPTIZYME™ restriction enzyme buffers contain BSA, which enhances the stability of many enzymes and binds contaminants that may be present in DNA preparations. Multiple freeze-thaw cycles of the buffers will not cause BSA precipitation. OPTIZYME™ restriction enzymes exhibit 100% of their certified activity in the recommended buffer.

Components: 330mM Tris-acetate (pH 7.9 at 37°C), 100mM Mg-acetate, 660mM K-acetate, 0.2% Triton X-100 and 1mg/mℓ BSA

Tested for: Functional Test, Endo- and Exonucleases Assay, Nicking Activity Assay

Recommended storage: -20°C

OPTIZYME™ 10X Buffer TaqI

packaging		Mfr. No
1 mℓ	Poly Tube	BP8092-1

Applications: The OPTIZYME™ Five Buffer System ensures the optimum reaction conditions for each restriction enzyme. This system consists of 10X Buffer 1, 10X Buffer 2, 10X Buffer 3, 10X Buffer 4 and 10X Buffer 5. The recommended buffer is supplied with each enzyme. To ensure consistent enzyme performance, OPTIZYME™ restriction enzyme buffers contain BSA, which enhances the stability of many enzymes and binds contaminants that may be present in DNA preparations. Multiple freeze-thaw cycles of the buffers will not cause BSA precipitation. OPTIZYME™ restriction enzymes exhibit 100% of their certified activity in the recommended buffer.

Components: 100mM Tris-HCl (pH 8.0 at 37°C), 50mM MgCl2, 1M NaCl, 1mg/mℓ BSA

Tested for: Functional Test, Endo- and Exonucleases Assay, Nicking Activity Assay

Recommended storage: -20°C

Denhardt’s Reagent
50X Powder

packaging		Mfr. No
150 mg	AmberGlass	BP520-5

DNase	Not detected
Optical Absorbance of a 5X Solution (at 278nm)	0.6-0.8
pH of 1X solution (at 25°C)	5.0-7.0
Protease	Not detected
RNase	Not detected

Applications: Denhardt’s Reagent is used as a blocking agent in DNA hybridizations.

A 50X concentrate is obtained by reconstituting with 5mℓ of sterile water.

Components: 50mg each of Ficoll*, Polyvinylpyrrolidone, and Bovine Albumin. [26873-85-8 (Ficoll)] ; [9003-39-8 (Polyvinylpyrrolidone)] ; [9048-46-8 (Bovine Albumin)]

Recommended Storage: 2°C to 8°C

Denhardt’s Reagent
50X Solution

packaging		Mfr. No
5 mℓ	AmberGlass	BP515-5

DNase	Not detected
Optical Absorbance of a 5X Solution (at 278nm)	0.6-0.8
pH of 1X solution (at 25°C)	5.0-7.0
Protease	Not detected
RNase	Not detected

Applications: Denhardt’s Reagent is used as a blocking agent in DNA hybridizations.

Components: <1% each of Ficoll*, Polyvinylpyrrolidone, and Bovine Albumin; and Water. Each 5mℓ of 50X Solution contains 50mg each of Ficoll*, Polyvinylpyrrolidone, and Bovine Albumin. [26873-85-8 (Ficoll*)] ; [9003-39-8 (Polyvinylpyrrolidone)] ; [9048-46-8 (Bovine Albumin)]

Recommended Storage: -20°C

Dextran
White to Off-white Powder

packaging		Mfr. No
100 g	PolyBottle	BP1580-100

CAS: 9004-54-0

Heavy Metals (Pb)	<=0.0005%
Loss on drying	<=5%
pH of 1% Solution in 0.9% NaCl	6.5 ±0.3
Sulfated ash	<=1%

Applications: Dextran may be used with Polyethylene Glycol under certain conditions in the fractionation of proteins, cells, virus particles, and other macromolecules.

Dextran at neutral pH may be autoclaved.

Average M.W. of 500.000

Recommended Storage: RT

Dextran Sulfate Sodium Salt
White to Off-white Powder

packaging		Mfr. No
100 g	PolyBottle	BP1585-100
500 g	PolyBottle	BP1585-500

CAS: 9011-18-1

Moisture	<=7.5%
Sulfur Content	17-20%

Applications: Dextran Sulfate accelerates the rate of hybridization of DNA or RNA bound to membranes and may be used in Southern and Northern blotting procedures.

Average M.W. of 500.000

Recommended Storage: RT

Formamide

packaging		Mfr. No
100 mℓ	AmberGlass	BP227-100
500 mℓ	AmberGlass,EcoSafPak*	BP227-500

CH ₃ NO	H360D	
CAS: 75-12-7	P201, P308+P313	
MW: 45.04		
EINECS: 200-842-0		

Assay	>=99.0%
Conductivity	<=350µmhos/cm
Copper	<=0.0001%
Freezing Point	1.0°-3.0°C
Iron	<=0.0005%
Lead	<=0.0005%
Optical Absorbance at 280nm	<=0.1
Zinc	<=0.0005%

H2NCH2O; CH3NO; F.W. 45.04

Applications: Molecular Biology-grade Formamide is commonly used in nucleic acid hybridization and sequencing.

It requires pretreatment with a mixed-bed resin.

Recommended Storage: 4°C

EcoSafPak* is an environmentally friendly packaging system made of 100% recyclable material by an SFI certified supplier.

Vacuum distilled and packaged under Nitrogen

Hybridization Cocktail
50% Formamide

packaging		Mfr. No
50 mℓ	PolyBottle	BP1575-50

DNase	Not detected
RNase	Not detected
Suitability for Nucleic Acid Hybridization	To pass test

Applications: Hybridization Cocktail is formulated for use in Southern and Northern transfers, dot blots, and in situ DNA and RNA hybridizations.

Components: <85.0% Formamide, 2.2% Sodium Citrate Dihydrate, 4.4% Sodium Chloride, and <6.0% Dextran, Hydrogen Sulfate, Sodium Salt [75-12-7 (Formamide)] ; [6132-04-3 (Sodium Citrate Dihydrate)] ; [7647-14-5 (Sodium Chloride)] ; [9011-18-1 (Dextran, Hydrogen Sulfate, Sodium Salt)]

Recommended Storage: -20°C

Saline-Sodium Citrate
20X solution

packaging		Mfr. No
1 ℓ	PolyBottle	BP1325-1
4 ℓ	PolyPac*	BP1325-4
20 ℓ	PolyPac*	BP1325-20

CAS: [7647-14-5 (Sodium Chloride)], [6132-04-3 (Sodium Citrate Dihydrate)]

DNase	Not detected
NaCl Molarity	3 ±0.10M
pH (20X solution) (at 23°C)	7.0±0.1
Protease	Not detected
RNase	Not detected

Filtered through a 2-5 micron filter.

Applications: Saline-Sodium Citrate (SSC) is used as a buffer in Southern Blot transfer protocols.

Recommended storage: RT

Saline-Sodium Phosphate-EDTA
20X solution

packaging		Mfr. No
1 ℓ	PolyBottle	BP1328-1
4 ℓ	PolyPac*	BP1328-4
20 ℓ	PolyPac*	BP1328-20

0.2M NaH₂PO₄, 3.0M NaCl and 0.02M EDTA

CAS: [7647-14-5 (Sodium Chloride)], [7558-80-7 (Sodium Phosphate Monobasic, anhydrous)], [60-00-4 (Ethylenediamine Tetraacetic Acid)]

DNase	Not detected
Ethylenediaminetetraacetic acid	0.6%
Molarity of NaCl	3 ±0.10M
pH (20X solution) (at 23°C)	7.4±0.1
Protease	Not detected
RNase	Not detected
Sodium Biophosphate	2.8%
Sodium (Na)	2.8% phosphate
sodium chloride	17.5 %

Filtered through a 2-5 micron filter.

Applications: Saline-Sodium Phosphate-EDTA is used in bacterial screening and hybridization procedures.

0.2M NaH₂PO₄, 3.0M NaCl and 0.02M EDTA

Recommended storage: RT

8-Hydroxyquinoline
White to Light-buff Needles or Powder

packaging		Mfr. No
100 g	AmberGlass	BP436-100

C ₉ H ₇ NO	H312, H302, H341	
CAS: 148-24-3	P281, P260, P301+P312,	
MW: 145.16	P304+P340, P302+P352,	
EINECS: 205-711-1	P280, P305+P351+P338	

H319, H315, H335, H332,

Assay	>=99.0%
Insoluble in Alcohol	<=0.05%
Melting Point	72.5°-74.0°C
Residue after ignition	<=0.05%
Suitability for Magnesium Determination	To pass test
Sulfate (SO ₄)	To pass test (about 0.02%)

Applications: 8-Hydroxyquinoline is used to stabilize phenol.

Recommended Storage: RT

Ammonium Acetate			
White Crystals			
packaging		Mfr. No	
500 g	AmberGlass	BP326-500	
1 kg	AmberGlass	BP326-1	
C ₂ H ₇ NO ₂			
CAS: 631-61-8			
MW: 77.08			
Assay		>=97.0%	
Chloride (Cl)		<=5ppm	
Heavy Metals (Pb)		<=5ppm	
Insoluble matter		<=0.005%	
Iron		<=5ppm	
Nitrate		<=0.001%	
Optical Absorbance of a 1M Solution at 254nm		<=0.02	
Optical Absorbance of a 1M Solution at 280nm		<=0.01	
Optical Absorbance of a 1M Solution at 350nm		<=0.01	
pH of 5% Solution (at 25°C)		6.7-7.3	
Residue after ignition		<=0.01%	
Sulfate (SO ₄)		<=0.001%	

Applications: Ammonium Acetate can be used to precipitate DNA from enzymatic reactions.
Recommended Storage: 4°C

Cesium Chloride Crystalline Powder			
packaging		Mfr. No	
500 g	PolyBottle	BP1595-500	
1 kg	PolyBottle	BP1595-1	
CICs		EINECS: 231-600-2	
CAS: 7647-17-8			
MW: 168.36			
Assay		>=99.5%	
DNase		Not detected	
Lithium metal in liquid paraffin		<=0.0002%	
Potassium		<=0.03%	
Protease		Not detected	
Refractive Index (50% solution at 25°C)		1.369 ±0.005	
RNase		Not detected	
Sodium		<=0.02%	
UV A260nm of a 50% Solution (w/v)		<=0.08	

Applications: This economical grade of Cesium Chloride can be used where low levels of metal ion impurities do not interfere. It is recommended for isolating phage DNA, plasmid DNA, or total RNA.
Recommended Storage: RT

Cesium Chloride		
White Crystalline Powder		
packaging		Mfr. No
100 g	PolyBottle	BP210-100
500 g	PolyBottle	BP210-500
1 kg	PolyBottle	BP210-1
CICs		EINECS: 231-600-2
CAS: 7647-17-8		
MW: 168.36		
Assay		>=99.999%
DNase		Not detected
Lithium metal in liquid paraffin		<=1ppm
Potassium		<=2ppm
Protease		Not detected
Refractive Index of a 50% Solution		1.369 ±0.002
RNase		Not detected
Sodium		<=4ppm
UV A260nm of a 50% Solution (w/v)		<=0.02

Applications: This Molecular Biology-grade cesium chloride is ideal for UV analysis of nucleic acids during density-gradient centrifugation and for the isolation of plasmid DNA, phage DNA, virus particles, and RNA.
Recommended Storage: RT

Ethanol, Absolute (200 proof)		DNase-, RNase- and Protease-Free	
Molecular Biology Grade			
Packaging		Mfr. No	
100 mL	Amber Glass	BP2818-100	
500 mL	Amber Glass	BP2818-500	
4 L	Amber Glass	BP2818-4	
C ₂ H ₆ O		H225	
CAS: 64-17-5		P210, P240	
MW: 46.07			
Acetone, Isopropyl alcohol		To pass test	
Assay		>=99.5%	
Benzene		<=2 ppm	
DNase		Not detected	
Infrared spectrum		Conforms to ref.	
Methanol		<=0.1%	
Protease		Not detected	
Residue after evaporation		<=0.001%	
RNase		Not detected	
Solubility in water		Pass test	
Substances darkened			
by H ₂ SO ₄		Pass test	
Substances			
reducing KMnO ₄		Pass test	
Titratable acid		<=0.0005mEq /g	
Titratable base		<=0.0002mEq /g	
Water Content		<=0.2%	

Applications: For the purification and precipitation of biomolecules such as nucleic acids and proteins. It is also used in histology to prepare staining and destaining reagents and for dehydrating tissues prior to embedding.
Recommended storage: RT

Ethylenediamine Tetraacetic Acid 0.5 M Solution, pH 8.0 (for DNA Work) Clear, Colorless Liquid			
Packaging		Mfr. No	
100 mL	AmberGlass	BP2482-100	
500 mL	AmberGlass	BP2482-500	
1 L	AmberGlass	BP2482-1	
20 L	PolyPac*	BP2482-20	
<chem>HOOCCH2CH2NCH2CH2N(CH2COOH)2</chem> C ₁₀ H ₁₆ N ₂ O ₈ CAS: 60-00-4 MW: 292.25		H319 P280, P305+P351+P338	
			
Molarity		0.475-0.525 M	
pH		8.0±0.1	

Applications: EDTA is used extensively in molecular biology to minimize metal ion impurities in reaction buffers. Ideal for DNA work.
Recommended Storage: RT

Guanidine Hydrochloride			
Colorless-to-white Crystals			
Packaging		Mfr. No	
500 g	PolyBottle	BP178-500	
1 kg	PolyBottle	BP178-1	
CH ₅ N ₃ .HCl		H315, H302, H319	
CAS: 50-01-1		P301+P312, P302+P352,	
MW: 95.53		P280, P305+P351+P338	
EINECS: 200-002-3			
Absorbance of a 6M Solution at 260nm		<=0.03	
Absorbance of a 6M Solution at 280nm		<=0.01	
Arsenic		<=5ppm	
Assay		>=99.5%	
Iron		<=5ppm	
Lead		<=5ppm	

H₂NC(NH₂)₂:NH.HCl;CH₅N₃.HCl;F.W.95.53
Applications: Guanidine Hydrochloride is a strong protein-denaturing agent used in the isolation of RNA from cell extracts. RNase activity is inhibited in the presence of this reagent.
Recommended Storage: RT

Guanidine Thiocyanate White Crystalline Powder			
packaging		Mfr. No	
250 g	PolyBottle	BP221-250	
1 kg	PolyBottle	BP221-1	
CH ₅ N ₃ .CHNS CAS: 593-84-0 MW: 118.16 H312, H332, H302, H412, H314, EUH032		P280, P301+P330+P331, P305+P351+P338, P310, P302+P352, P273, P233	
		 	
A280nm (70% in H ₂ O at 25°C)		<=0.8	
Assay		>=99.0%	
Melting Point		117°-122°C	
Solubility (70% in H ₂ O at 25°C)		Clear and haze-free	
Applications: A more powerful denaturing agent than Guanidine Hydrochloride, Guanidine Thiocyanate is suitable for the isolation of intact RNA from cell extracts.			
Recommended Storage: RT			

Isopropanol		Molecular Biology-Grade	
packaging		Mfr. No	
500 ml	Amber Glass	BP2618-500	
1 ℓ	Amber Glass	BP2618-1	
2.5 ℓ	Amber Glass	BP2618-212	
4 ℓ	Amber Glass	BP2618-4	
C ₃ H ₈ O		H225, H336, H319	
CAS: 67-63-0		P261, P280,	
MW: 60.1		P305+P351+P338, P210,	
EINECS: 200-661-7		P240	
			
			
Acetone			<=0.002%
Assay			>=99.9%
Color (APHA)			<=5 APHA
DNase			Not detected
Fluorescence Background (as Quinine Sulfate)			<=1 ppb
Optical Abs. at 205nm			<=1
Optical Abs. at 220nm			<=0.20
Optical Abs. at 230nm			<=0.10
Optical Abs. at 254nm			<=0.015
Propionaldehyde			<=0.002%
Protease			Not detected
Ref. index at 25°C			Inclusive between 1.3740-1.3760
Residue after evaporation			<=1 ppm
RNase			Not detected
Solubility in water			To pass test
Substances			
reducing KMnO ₄			To pass test
Titrat. acid or base			<=0.0001 mEq/g.
Water Content			<=0.05%

Applications: Used in fundamental applications, such as purification and precipitation of nucleic acids and proteins and preservation of biological specimens.
Recommended storage: RT

Phenol/Chloroform/Isoamyl Alcohol			
125:24:1 Mixture, pH 4.3, Liquid			
Packaging		Mfr. No	
100 mL	Coated Amber Glass/Poison Pack	BP17541-100	
400 mL	Coated Amber Glass/Poison Pack	BP17541-400	
CAS: 108-95-2		P280, P301+P330+P331,	
EINECS: 203-632-7		P302+P350, P304+P340,	
H301, H314, H341, H318,		P305+P351+P338, P310,	
H331, H311, H351, H373		P260	
			
			
			
Assay		>99.0%	
Optical Absorbance at 330nm		<0.2	
Optical Absorbance at 405nm		<0.1	
Optical Absorbance at 510nm		<0.1	
pH		4.3±0.2	

Applications: This mixture is used to purify total RNA from mixtures containing RNA, DNA, and protein.
Recommended Storage: 4°C
UN 2821; DOT Class 6.1:Poison

Phenol/Chloroform/Isoamyl Alcohol 25:24:1 Mixture, pH 5.2, Liquid			
Packaging		Mfr. No	
100 mL	Coated Amber Glass/Poison Pack	BP17531-100	
400 mL	Coated Amber Glass/Poison Pack	BP17531-400	
CAS: 108-95-2		P280, P301+P330+P331,	
EINECS: 203-632-7		P302+P350, P304+P340,	
H301, H314, H341, H318, H331, H311, H351, H373		P305+P351+P338, P310, P260	
			
			
			
Assay		>99.0%	
Optical Absorbance at 330nm		<0.2	
Optical Absorbance at 405nm		<0.1	
Optical Absorbance at 510nm		<0.1	
pH (with buffer)		5.2±0.3	

Applications: This mixture is used for extracting protein from nucleic acid preparations.
Recommended Storage: 4°C
UN 2821; DOT Class 6.1:Poison



Phenol/Chloroform/Isoamyl Alcohol
25:24:1 Mixture, pH 6.7/8.0, Liquid

packaging	Mfr. No
100 mℓ CoatedAmberGlass/PoisonPack	BP1752I-100
400 mℓ CoatedAmberGlass/PoisonPack	BP1752I-400
CAS: 136112-00-0 H341, H373, H314, H351, H301, H311, H331 P301+P330+P331, P280, P305+P351+P338, P310, P281, P304+P340, P260	
  	
Assay	>99.0%
Optical Absorbance at 330nm	<0.2
Optical Absorbance at 405nm	<0.1
Optical Absorbance at 510nm	<0.1
pH	6.7±0.2
Applications: This mixture is used for extracting protein from crude nucleic acid preparations. Recommended Storage: 4°C UN 2821; DOT Class 6.1:Poison	

Phenol, Saturated
pH 4.3, Liquid

packaging	Mfr. No
100 mℓ CoatedAmberGlass/PoisonPack	BP1751I-100
400 mℓ CoatedAmberGlass/PoisonPack	BP1751I-400
C ₆ H ₅ OH MW: 94.04 H301, H311, H331, H314, H373, H341 P280, P302+P350,	
P301+P330+P331, P305+P351+P338, P304+P340, P301+P310, P260	
  	
Assay	>99.0%
Optical Absorbance at 330nm	<0.2
Optical Absorbance at 405nm	<0.1
Optical Absorbance at 510nm	<0.1
pH (with buffer)	4.3±0.2
Applications: Phenol saturated with Tris buffer is used for purifying RNA. Capped with a layer of buffer. Recommended Storage: 4°C UN 2821; DOT Class 6.1:Poison	

Phenol, Saturated
pH 6.6/7.9, Liquid

packaging	Mfr. No
100 mℓ CoatedAmberGlass/PoisonPack	BP1750I-100
400 mℓ CoatedAmberGlass/PoisonPack	BP1750I-400
H314, H331, H373, H341, H311 P280, P302+P350, P301+P330+P331, P305+P351+P338, P304+P340, P301+P310, P260	
  	
Assay	>99.0%
DNase	Not detected
Optical Absorbance at 330nm	<0.2
Optical Absorbance at 405nm	<0.1
Optical Absorbance at 510nm	<0.1
pH	6.6±0.2
pH (with buffer)	7.9±0.2
RNase	Not detected
Applications: Phenol saturated with Tris buffer is used in DNA extraction procedures. Capped with a layer of buffer. Recommended Storage: 4°C UN 2821; DOT Class 6.1:Poison	

Proteinase K
From Tritirachium album

packaging	Mfr. No
50 mg PolyBottle	BP1700-50
100 mg PolyBottle	BP1700-100
500 mg PolyBottle	BP1700-500
CAS: 39450-01-6 EINECS: 254-457-8 H315, H319, H335, H334	
P261, P280, P305+P351+P338	
 	
DNase	Not detected
RNase	Not detected
Specific Activity	>=30U/mg

Applications: Proteinase K is widely used in the isolation of DNA and RNA from cell and tissue preparations. It is a wide-spectrum protease, inactivating many enzymes, including endogenous nucleases.
Recommended Storage: -20°C
Not on TSCA inventory: for R and D use only; not for manufacturing or commercial purposes.

Sodium Acetate Anhydrous
White Crystals or Granular Powder

packaging	Mfr. No
500 g PolyBottle	BP333-500
1 kg PolyBottle	BP333-1
C ₂ H ₃ NaO ₂ CAS: 127-09-3 MW: 82.03	
EINECS: 204-823-8	
Assay	>=99.0%
Calcium (Ca)	<=0.005%
Chloride (Cl)	<=0.002%
Heavy Metals (Pb)	<=0.001%
Insoluble matter	<=0.01%
Iron	<=0.001%
Loss on Drying (at 120°C)	<=1.0%
Magnesium (Mg)	<=0.002%
pH of 5% Solution (at 25°C)	7.0-9.2
Phosphate (PO ₄)	<=0.001%
Sulfate (SO ₄)	<=0.003%

Applications: Sodium Acetate is added to DNA reactions to facilitate ethanol precipitation of the DNA.
Recommended Storage: RT

Sodium Acetate Trihydrate
Colorless-to-white Crystals

packaging	Mfr. No
500 g PolyBottle	BP334-500
C ₂ H ₃ NaO ₂ ·3H ₂ O CAS: 6131-90-4 MW: 136.08	
Assay	
99.0-101.0%	
Calcium (Ca)	<=0.005%
Chloride (Cl)	<=0.001%
Heavy Metals (Pb)	<=5ppm
Insoluble matter	<=0.005%
Iron	<=5ppm
Magnesium (Mg)	<=0.002%
pH of 5% Solution (at 25°C)	7.5-9.2
Phosphate (PO ₄)	<=5ppm
Potassium	<=0.005%
Substances Reducing Permanganate	To pass test
Sulfate (SO ₄)	<=0.002%

Applications: The trihydrate form of Sodium Acetate can be used for many of the same applications as the anhydrous form.
Recommended Storage: RT

Sodium Chloride-Tris-EDTA (STE) Buffer
10X, pH 8.0

packaging	Mfr. No
100 mℓ PolyBottle	BP2479-100
500 mℓ PolyBottle	BP2479-500
1 ℓ PolyBottle	BP2479-1
DNase	
Not detected	
pH of 1X solution (at 25°C)	7.9-8.1
Protease	Not detected
RNase	Not detected

Applications: STE buffer is suitable for biomolecular lab procedures. 10X Solution contains 100mM Tris-HCl, 10mM EDTA, and 1M NaCl. Autoclaved.
[77-86-1 (Tris Base)] ; [60-00-4 (EDTA)] ; [7647-01-0 (HCl)] ; [7647-14-5 (NaCl)]
Recommended Storage: RT
Filtered through a 5-micron filter.

Sodium Chloride-Tris-EDTA Buffer
1X, pH 8.0

packaging	Mfr. No
500 mℓ PolyBottle	BP2478-500
1 ℓ PolyBottle	BP2478-1
DNase	
Not detected	
pH (at 25°C)	7.9-8.1
Protease	Not detected
RNase	Not detected

Applications: STE buffer is suitable for biomolecular lab procedures. 1X solution contains 10mM Tris-HCl, 1mM EDTA and 100mM NaCl. Autoclaved.
[77-86-1 (Tris Base)] ; [60-00-4 (EDTA)] ; [7647-01-0 (HCl)] ; [7647-14-5 (NaCl)]
Recommended Storage: RT
Filtered through a 5-micron filter.

Spermidine Trihydrochloride
White Crystalline Powder

packaging	Mfr. No
1 g AmberGlass	BP2540-1
5 g AmberGlass	BP2540-5
25 g AmberGlass	BP2540-25
C ₇ H ₁₉ N ₃ ·3HCl CAS: 334-50-9 MW: 254.63 EINECS: 206-379-0	
H319, H335, H315 P261, P302+P352, P280, P305+P351+P338	
	
Assay	>=99%
FTIR	Conforms to standard
Melting Point	253°-263°C

Applications: Spermidine Trihydrochloride forms stable compounds with nucleic acids. It is used in biochemical research.
Recommended Storage: Store below 20°C, desiccate.

SurePrep* FFPE RNA Isolation Kit
For the rapid and efficient extraction and purification of RNA (including microRNA) from FFPE samples

packaging	Mfr. No
50 preps	BP2816-50

Average Yield .. Variable due to age of paraffin blocks; 2-3 µg of total RNA per 1mg of fresh FFPE hamster kidney
Binding capacity per column Up to 50µg RNA
Maximum amount of starting material. 5 slices of 20mm thick paraffin slices (25 mg of unsectioned block)
Maximum loading volume per spin column 650 µℓ
Size of RNA purified All sizes, including small RNA and microRNA (<200n.t.)
Time to complete 10 purifications 1 to 4 hours

Applications: Rapid method for the isolation and purification of total RNA from formaldehyde-fixed paraffin-embedded (FFPE) tissue samples.
Description: High quality and integrity of the isolated RNA - The purified total RNA is of the highest quality and integrity, and can be used in any sensitive downstream applications; Isolate a diversity of RNA species - All RNA species can be isolated, from large mRNA and ribosomal RNA down to microRNA (miRNA) and small interfering RNA (siRNA); High yields - The FFPE RNA Purification Kit allows for the purification of high yields of total RNA; No phenol:chloroform extractions - Total RNA is isolated from FFPE tissue samples without the use of harmful chemicals such as phenol or chloroform; Rapid procedure - Isolate total RNA from FFPE tissue sections using a rapid spin column format in as little as 1 hour
Components: Enzyme incubation buffer (6mℓ); Digestion buffer (20mℓ); Wash solution (22mℓ); Binding solution (20mℓ); Proteinase K (12mg); DNase I (1 vial); Elution buffer (6mℓ); Spin columns (50); Collection tubes (50); Elution tubes, 1.7mℓ (50).
Recommended storage: RT

SurePrep® Plant/Fungi Total RNA Purification Kit
For the rapid purification of total RNA (including microRNA) from plants and fungi

packaging	Mfr. No
50 preps	BP2817-50

Average yield: 50mg Grape leaves	35 µg
Average yield: 50mg Peach leaves	30 µg
Average yield: 50mg Plum leaves	32 µg
Average yield: 50mg Tobacco leaves	60 µg
Average yield: 50mg Tomato leaves	60 µg
Column binding capacity	50 µg
Maximum amount of starting material: Fungi (wet weight)	50 µg
Maximum amount of starting material: Plant cells	1 X 10 ⁶ cells
Maximum amount of starting material: Plant tissues	50 mg
Maximum column loading volume	650 µL
Size of RNA purified	All sizes, including small RNA and microRNA (<200n.t.)
Time to complete 10 purifications	30 min

Applications: Rapid method for the isolation and purification of total RNA from a wide range of plant and filamentous fungi species.
Description: No liquid nitrogen required for homogenization - Liquid nitrogen is not required for homogenization of samples, making RNA purification rapid and convenient; Isolate RNA from a wide range of samples - Total RNA can be isolated from a wide range of plant and filamentous fungi samples; Isolate a diversity of RNA species - All sizes of RNA are isolated from large mRNA down to microRNA, without the use of phenol or chloroform; Isolate total RNA, including viral RNA - RNA samples can be used for the detection of viral pathogens, as viral RNA is isolated with the total RNA; High yield of RNA - High yields of purified RNA can be isolated with this kit; No phenol:chloroform extractions - Total RNA is isolated without the use of harmful chemicals such as phenol or chloroform. RNA is of the highest quality and can be used in a number of downstream applications.
Components: Wash solution (22mL); Lysis solution (40mL); Elution buffer (6mL); Spin columns (50); Collection tubes (50); Elution tubes, 1.7mL (50).
Recommended storage: RT

SurePrep® Soil DNA Isolation Kit
For the rapid preparation of inhibitor-free DNA from common types of soil

packaging	Mfr. No
50 preps	BP2815-50

Column binding capacity	50 µL
Maximum column loading volume	650 µL
Maximum soil input	250 mg
Time to complete 10 purifications	30 min
Type of soil processed	All common soil types, including clay, loam and sandy soil

Applications: Isolation and purification of total DNA from microorganisms in soil and sediment samples, such as bacteria, fungi, and algae.
Description: Process all types of soil - All types of soil samples can be processed with this kit, including common soil samples such as clay, loam and sandy soil; Remove all humic acid from DNA samples - The kit removes all traces of humic acid using the provided Bead Tubes and a combination of chemical and physical homogenization and lysis. All the brown colour is removed from the samples; Rapid detection of microorganisms in soil samples - Isolate total DNA from all microorganisms found in soil, including bacteria, fungi and algae; Fast and easy processing - Rapid spin column format allows for the isolation of DNA in under 30 minutes; .Isolate high quality total DNA - The purified DNA is free from all inhibitors including humic acid, and can be used directly in downstream applications including PCR.
Components: Lysis additive (6mL); Binding solution (6mL); Lysis solution (45mL); Wash solution I (30mL); Wash solution II (18mL); Elution buffer (6mL); Spin columns (50); Collection tubes (50); Beads (30g); Screw cap tubes (50); Elution tubes, 1.7 mL (50).
Recommended storage: RT

SurePrep® Water RNA/DNA Purification Kit
For the rapid isolation of DNA/RNA from microorganisms in non-turbid water samples

packaging	Mfr. No
25 preps	BP2814-25

Maximum filter column loading volume	20 mL
Maximum spin column loading volume	650 µL
Maximum water input	100 mL
Elution volume	100 µL
safety	10 mL
Time to complete 10 purifications	45 min

Applications: Isolation and purification of total RNA and DNA simultaneously from all microorganisms found in environmental water samples, including bacteria, fungi, and algae.
Description: Rapid isolation of DNA/RNA from microorganisms in non-turbid water samples -Only three buffers are required to complete isolation of total DNA and RNA from all microorganisms found in water, including bacteria, fungi and algae; Complete column purification - The RNA and DNA are both column purified simultaneously using the same column; Fast and easy processing - Rapid spin column format allows for the isolation of both RNA and DNA in under 45 minutes; Isolate highly concentrated RNA and DNA - The purified RNA and DNA is eluted in low volumes, allowing for the direct use of the nucleic acid in downstream applications; Isolate a diversity of RNA species - All sizes of RNA are isolated, from large mRNA down to microRNA, without the use of phenol or chloroform
Components: Lysis solution (15mL); Enzyme incubation buffer (4mL); Nucleic acid elution buffer (4mL); Nucleic acid wash solution (11mL); Spin columns (25); Collection tubes (25); Elution tubes, 1.7mL (25); Bead tubes (25); Filter column, 0.45µm (25).
Recommended storage: RT

N-Methylpyrrolidinone

packaging	Mfr. No
4 L AmberGlass,EcoSafPak®	BP1172-4
C ₅ H ₉ NO	P201, P308+P313, P261,
CAS: 872-50-4	P302+P352, P280,
MW: 99.13	P305+P351+P338
H319, H315, H335, H360D	

Assay (GC)	>=99.8%
Color (APHA)	<=20
IR	Conforms to standard
Optical Absorbance at 265nm	<=1.0
Optical Absorbance at 270nm	<=0.7
Optical Absorbance at 300nm	<=0.3
Optical Absorbance at 400nm	<=0.02
Titrateable Base	<=0.001mEq/g
Water	<=0.05%

Recommended storage: RT
Also available in recyclable FisherPak® and NOWPak® containers. EcoSafPak® is an environmentally friendly packaging system made of 100% recyclable material by an SFI certified supplier.

Fisher BioReagents® SurePrep® Small RNA Purification Kit

Optimized for enrichment of small RNA by removing large RNA for studying regulation of gene expression

- ▶ Isolate all small RNA molecules (<200 nucleotides)
- ▶ Minimal contamination from large RNA molecules and genomic DNA
- ▶ Fast and easy processing using rapid spin-column format
- ▶ No phenol or chloroform extractions
- ▶ High quality small RNA such as miRNA and siRNA can be used in various downstream applications



BP2801-25

APPLICATIONS: Purified small RNA is suitable for use in RT-PCR, Northern blotting, and microarray analysis

INCLUDES: Lysis solution (2 x 20mL), Wash solution (22mL), RNA elution buffer (6mL), large RNA removal columns (25), small RNA enrichment columns (25), Collection tubes (50), Elution tubes (1.7mL) (50) and Product insert

STORAGE CONDITIONS: All solutions should be kept tightly sealed and stored at room temperature. Reagents should remain stable for at least two years in unopened containers.

Kit Specifications	
Column binding capacity	50µg
Maximum column loading volume	600µL
Size of RNA Purified	<200 nt
Maximum amount of starting material:	
Animal cells	3 × 10 ⁶ cells
Animal tissues	25mg
Bacteria	5 × 10 ⁶ cells
Plant tissues	50mg
Time to complete ten purifications	45 minutes

No. of Preps	MFR No.
25	BP2801-25

Fisher BioReagents® SurePrep® TrueTotal® RNA Purification Kit

Rapidly isolate and purify total RNA from a variety of biological sources

- ▶ Fast and easy processing using rapid spin-column format
- ▶ Isolate all sizes of RNA including large RNA and small RNA species
- ▶ No phenol or chloroform extractions
- ▶ Isolate and detect RNA from as little as a single animal cell
- ▶ Purify total RNA from a variety of biological sources such as cultured animal cells, tissue samples, blood, bacteria, yeast, fungi and plants



BP2800-50

APPLICATIONS: Purified intact RNA is suitable for Quantitative RT-PCR, Northern blotting, RNase protection, and expression array assays

INCLUDES: Lysis solution (40mL), Wash solution (22mL), Elution buffer (6mL), Micro spin columns (50), Collection tubes (50), Elution tubes (1.7mL) (50) and product insert

STORAGE CONDITIONS: All solutions should be kept tightly sealed and stored at room temperature. Reagents should remain stable for at least two years in unopened containers.

Fisher BioReagents® SurePrep® Leukocyte RNA Purification Kit

Rapid method for the isolation and purification of total leukocyte (white blood cell) RNA from mammalian blood samples

- ▶ Fast and easy processing using rapid spin-column format
- ▶ Differential red blood cell lysis allows for the removal of a majority of globin mRNA
- ▶ Isolate total leukocyte RNA, including all small RNA species
- ▶ No phenol or chloroform extractions
- ▶ High quality leukocyte RNA can be used in a number of downstream applications



BP2807-50

APPLICATIONS: Purified RNA is suitable for use in real time PCR, reverse transcription PCR, Northern blotting, RNase protection and primer extension, and expression array assays

INCLUDES: RBC Lysis Solution (180mL); Binding Solution (25mL); Wash Solution (22mL); RNA Elution Buffer (6mL); Spin Columns (50); Collection Tubes (50); Elution tubes (1.7mL) (50); Product Insert (1)

STORAGE CONDITIONS: All solutions should be kept tightly sealed and stored at room temperature. The RBC Lysis Solution should be stored at 4°C upon arrival. These reagents should remain stable for at least two years in their unopened containers.

Kit Specifications	
Column binding capacity	50µg
Maximum Column Loading Volume	600µL
Size of RNA Purified	All sizes, including small RNA (<200 nt)
Maximum Blood Input	2mL or 3 x 10 ⁶ Leukocytes
Minimum Blood Input	10µL
Time to complete ten purifications	45 minutes
Average Yield: 500µL human blood	1.5µg

No. of Preps	MFR No.
50	BP2807-50

Kit Specifications	
Column binding capacity	50µg
Maximum column loading volume	600µL
Size of RNA Purified	All sizes, including small RNA (<200 nt)
Maximum amount of starting material:	
Animal cells	3 × 10 ⁶ cells
Animal tissues	25mg
Blood	100µL
Bacteria	1 × 10 ⁶ cells
Yeast	1 × 10 ⁶ cells
Fungi	50mg
Plant tissues	50mg
Time to complete ten purifications	20 minutes

Average yields¹	
HeLa cells (1 × 10 ⁶ cells)	15µg
E. coli (1 × 10 ⁶ cells)	50µg

¹ Average yields will vary depending upon a number of factors including species, growth conditions, and developmental stage.

No. of Preps	MFR No.
50	BP2800-50

Fisher BioReagents* SurePrep* RNA Cleanup and Concentration Kit

Provides a rapid method to purify and concentrate RNA isolated using other methods such as phenol:chloroform extraction or after enzymatic reactions

- Fast and easy processing using rapid spin-column format
- Suitable for all sizes of RNA, from large mRNA to microRNA (miRNA)
- Able to clean up RNA isolated using different methods, including phenol/chloroform extractions
- Efficient RNA cleanup from various enzymatic reactions such as labeling, DNase treatment and in vitro transcription



BP2809-50

STORAGE CONDITIONS: All solutions should be kept tightly sealed and stored at room temperature. Reagents should remain stable for at least two years in unopened containers.

Kit Specifications	
Column binding capacity	50µg
Size of RNA Purified	All sizes, including small RNA (<200 nt)
Maximum Amount of Starting Material	50µg of RNA
Minimum Elution Volume	20µL
Time to complete ten purifications	20 minutes
Average Recovery	≥90%

No. of Preps	MFR No.
50	BP2809-50

APPLICATIONS: Purified RNA is suitable for Quantitative RT-PCR, Northern blotting, RNase protection, primer extension, and expression array assays

INCLUDES: Binding solution (2 x 20mL), Wash solution (22mL), RNA elution buffer (6mL), Collection tubes (50), Spin columns (50), Elution tubes (1.7mL) (50), Product insert

Fisher BioReagents* SurePrep* RNA/Protein Purification Kit

Rapid method for the isolation and purification of total RNA and proteins sequentially from a single sample using the same column in less than 20 minutes

- Fast and easy processing using rapid spin-column format
- All columns for RNA purification and protein purification provided
- Sequentially isolate high quality RNA and proteins from a single lysate – no need to split the lysate
- Isolate all sizes of RNA, from large mRNA and rRNA down to microRNA (miRNA)
- No phenol or chloroform extractions
- Purify total RNA and proteins from a variety of biological sources such as cultured animal cells, tissue samples, blood, bacteria, yeast, fungi or plants



BP2806-50

Kit Specifications	
Column binding capacity	50µg for RNA 200µg for protein
Maximum Column Loading Volume	600µL
Size of RNA Purified	All sizes, including small RNA (<200 nt)
Maximum Amount of Starting Material:	
Animal Cells	3 x 10 ⁶ cells
Animal Tissues	25mg
Blood	100µL
Bacteria	1 x 10 ⁸ cells
Yeast	1 x 10 ⁸ cells
Fungi	50mg
Plant Tissues	50mg
Time to complete ten purifications	30 minutes

Average Yield*	
HeLa Cells (1 x 10 ⁶ cells)	15µg RNA
HeLa Cells (1 x 10 ⁶ cells)	150µg protein

* Average yields will vary depending upon a number of factors including species, growth conditions used and developmental stage.

No. of Preps	MFR No.
50	BP2806-50

APPLICATIONS: Purified total RNA and protein from same biological sample for studies of gene expression including gene silencing experiments, mRNA knockdowns, or experiments correlating RNA and protein expression levels. RNA qualified for use in real time PCR, reverse transcription PCR, Northern blotting, RNase protection and expression array assays. Proteins suitable for use in SDS-PAGE analysis and Western blots.

INCLUDES: Lysis solution (40mL); Nucleic acid wash solution (22mL); Nucleic acid elution buffer (15mL); Protein column regeneration buffer (30mL); Protein column activation and wash buffer (60mL); Protein pH binding buffer (4mL); Protein elution buffer (8mL); Enzyme incubation buffer (6mL); Protein neutralizer (2mL); Protein loading dye (2mL); Spin columns (50); Collection tubes (150); Product insert (1)

STORAGE CONDITIONS: The Protein Loading Dye should be stored at -20°C upon arrival. All other solutions should be kept tightly sealed and stored at room temperature. These reagents should remain stable for at least one year in their unopened containers.

Fisher BioReagents* SurePrep* Nuclear or Cytoplasmic RNA Purification Kit

Rapid method for the isolation and purification of both nuclear and cytoplasmic RNA from cultured animal cells and animal tissue samples

- Fast and easy processing using rapid spin-column format
- High quality intact RNA can be used in various downstream applications
- Enrich for cytoplasmic or nuclear RNA depending on downstream application
- No genomic DNA contamination in the cytoplasmic fraction
- No phenol or chloroform extractions
- Isolate all sizes of RNA from each fraction, including all small RNA molecules (>200 nt)
- Kit provides sufficient reagents to perform either 50 cytoplasmic RNA preparations or 25 cytoplasmic and nuclear preparations



BP2805-50

APPLICATIONS: Purification of Nuclear and Cytoplasmic RNA suitable for real time PCR, reverse transcription PCR, Northern blotting, RNase protection and primer extension, and expression array assays.

INCLUDES: Lysis solution (20mL); Wash solution (22mL); Binding solution (25mL); RNA elution solution (6mL); Spin columns (50); Collection tubes (50); Elution tubes (1.7mL) (50); Product insert (1)

STORAGE CONDITIONS: All solutions should be kept tightly sealed and stored at room temperature. These reagents should remain stable for at least two years in their unopened containers.

Kit Specifications	
Column binding capacity	50µg
Maximum Column Loading Volume	600µL
Size of RNA Purified	All sizes, including small RNA (<200 nt)
Maximum Amount of Starting Material:	
Animal Cells	3 x 10 ⁶ cells
Animal Tissues	15mg
Time to complete ten purifications	45 minutes

Average Yields*	
HeLa Cells - cytoplasmic (1 x 10 ⁶ cells)	15µg
HeLa Cells - nuclear (1 x 10 ⁶ cells)	≤3.5µg

No. of Preps	MFR No.
50 (Cyto); 25 (Nuc & Cyto)	BP2805-50

Fisher BioReagents* MOPS 10X Buffer for Agarose RNA Electrophoresis

Convenient, easy to use 10X concentrated solution

A low ionic strength buffer used for electrophoretic separation of formaldehyde-denatured RNA.

- Easy to dilute solution, pH 7.0
- One-year shelf life
- 0.2-micron filtered
- RNase and DNase free
- Available in 500mL and 1L sizes, both in poly translucent bottles

Quantity	MFR No.
500mL	BP2900-500
1L	BP2900-1

Fisher BioReagents* SurePrep* Urine Bacterial RNA Purification Kit

Designed for the rapid preparation of bacterial RNA from human or animal urine samples

- Fast and easy processing using rapid spin-column format
- Isolate total RNA, from large mRNA down to microRNA (miRNA)
- No phenol or chloroform extractions
- Isolate RNA from both Gram negative and Gram positive bacteria, including *E. coli*, *Proteus spp.*, *Klebsiella spp.*, *Enterobacter spp.*, *Serratia spp.*, *Pseudomonas spp.*, *Clostridial spp.* and *Leptospiriosis spp.*, as well as *Chlamydia trachomatis* and *Neisseria gonorrhoeae*



BP2804-50


APPLICATIONS: Purified bacterial RNA suitable for real time PCR, reverse transcription PCR, Northern blotting, RNase protection, primer extension, and expression array assays.

INCLUDES: Lysis solution (40mL); Wash solution (22mL); RNA elution buffer (6mL); Spin columns (50); Collection tubes (50); Elution tubes (1.7mL) (50); Product insert (1)

STORAGE CONDITIONS: All solutions should be kept tightly sealed and stored at room temperature. These reagents should remain stable for at least two years in their unopened containers.

Kit Specifications	
Column binding capacity	50µg
Volume of Urine Processed	10-50mL
Size of RNA Purified	All sizes, including small RNA (<200 nt)
Time to complete ten purifications	30 minutes
Average Yield	~0.5µg RNA per 1 x 10 ⁷ cells (Varies due to cell density of sample)

No. of Preps	MFR No.
50	BP2804-50



Lanes 2 and 3 demonstrate the electrophoresis of 1µg of purified RNA from *E. coli*, resolved on a 1.5% denaturing formaldehyde-agarose gel, using the MOPS 10X Buffer (BP2900), diluted to 1X as a running buffer. Lanes 1 and 4 contain a 1kb RiboLadder RNA standard (BP2811-50). Electrophoresis was performed at 100V for 70 min., and the gel was stained with ethidium bromide (BPT302-10).

Fisher BioReagents SurePrep® Urine Exfoliated Cell RNA Purification Kit

Rapid method for the isolation and purification of total RNA from exfoliated cells that have been shed into the urine from the urinary tract

- Fast and easy processing using rapid spin-column format
- Isolate total RNA, from large mRNA down to microRNA (miRNA)
- RNA can be isolated and detected from as little as 100 exfoliated cells
- Isolate high quality total RNA from urine
- No phenol or chloroform extractions



BP2803-50

APPLICATIONS: Purified RNA from urine exfoliated cells suitable for real time PCR, reverse transcription PCR, Northern blotting, RNase protection and primer extension, and expression array assays

INCLUDES: Lysis solution (40mL); Wash solution (22mL); RNA elution solution (6mL); Spin columns (50); Collection tubes (50); Elution tubes (1.7mL) (50); Product insert (1)

STORAGE CONDITIONS: All solutions should be kept tightly sealed and stored at room temperature. These reagents should remain stable for at least two years in their unopened containers.

Kit Specifications	
Column binding capacity	50µg
Volume of Urine Processed	1-50mL
Maximum Input of Exfoliated Cells	1 x 10 ⁶ cells
Size of RNA Purified	All sizes, including small RNA (~200 nt)
Time to complete ten purifications	20 minutes
Average Yield	~1µg RNA per 1 x 10 ⁵ cells (varies due to cell density of sample)

No. of Preps	MFR No.
50	BP2803-50

Fisher BioReagents® SurePrep® RNA/DNA/Protein Purification Kit

Rapid method for the isolation and purification of total RNA, genomic DNA and proteins sequentially from a single sample of cultured animal cells, tissue samples, blood, bacteria, yeast, fungi or plants

The total RNA, genomic DNA and proteins are all column-purified using the same column. This kit is ideal for researchers who are interested in studying the genome, proteome and transcriptome of a single biological sample, such as studies of gene expression including gene silencing experiments or mRNA knockdowns, studies involving biomarker discovery, studies in epigenetics, and for characterization of cultured cell lines.



BP2802-50

- Fast and easy processing using rapid spin-column format
- All columns for RNA, DNA and protein purification provided
- Sequential isolation of nucleic acids and proteins from a single lysate, no need to split the lysate
- Isolate all sizes of RNA, including large mRNA and rRNA to microRNA (miRNA)
- Isolate high quality total RNA
- Isolate high quality genomic DNA with a molecular weight ≥ 30kb
- High yields of isolated proteins
- No phenol or chloroform extractions

APPLICATIONS: Purified RNA, DNA and protein suitable for assays such as PCR, real time PCR, reverse transcription PCR, Northern blotting, RNase protection, expression array assays, DNA sequencing, Southern blotting, SNP analysis, SDS-PAGE and Western Blots

INCLUDES: Lysis solution (40mL); RNA wash solution (20mL); RNA elution solution (6mL); gDNA wash solution (15mL); gDNA elution buffer (15mL); Protein column regeneration buffer (30mL); Protein column activation and wash buffer (60mL); Protein pH binding buffer (4mL); Protein elution buffer (8mL); Protein neutralizer (2mL); Protein loading dye (2mL); Spin columns (50); Collection tubes (200); Product inserts (1)

STORAGE CONDITIONS: The Protein Loading Dye should be stored at -20°C upon arrival. All other solutions should be kept tightly sealed and stored at room temperature. These reagents should remain stable for at least one year in their unopened containers

Kit Specifications	
Column binding capacity	50µg for RNA 20µg for DNA 200µg for protein
Maximum column loading volume	600µL
Size of RNA Purified	All sizes, including small RNA (~200nt)
Size of DNA Purified	≥30kb

Maximum amount of starting material:	
Animal cells	3 x 10 ⁶ cells
Animal tissues	25mg
Blood	100µL
Bacteria	1 x 10 ⁸ cells
Yeast	1 x 10 ⁸ cells
Fungi	50mg
Plant tissues	50mg
Time to complete ten purifications	30 minutes

Average Yields*	
HeLa Cells (1 x 10 ⁶ cells)	15µg RNA
HeLa Cells (1 x 10 ⁶ cells)	8µg DNA
HeLa Cells (1 x 10 ⁶ cells)	150µg protein

* Average yields will vary depending upon a number of factors including species, growth conditions used and developmental stage.

No. of Preps	MFR No.
50	BP2802-50

4-Dimethylaminopyridine

Peptide Synthesis

packaging	Mfr. No
100 g AmberGlass	BP596-100
C ₇ H ₁₀ N ₂ CAS: 1122-58-3 MW: 122.17 EINECS: 214-353-5 H310, H314, H301	P280, P361, P302+P350, P301+P330+P331, P305+P351+P338, P310



Assay	≥99%
IR	Conforms to standard
Melting Point	108°-111°C

Applications: 4-Dimethylaminopyridine is a catalyst for several acylation reactions, especially in oligonucleotide synthesis.

Recommended Storage: RT
UN 2928; DOT Class 6.1:Poison

5-Bromo-2'-Deoxyuridine

Off-white Powder

packaging	Mfr. No
250 mg AmberGlass	BP2508-250
500 mg AmberGlass	BP2508-500
1 g AmberGlass	BP2508-1
5 g AmberGlass	BP2508-5

C ₉ H ₁₁ BrN ₂ O ₅ CAS: 59-14-3 MW: 307.09 H335, H319, H341, H315	P281, P261, P302+P352, P280, P305+P351+P338
--	--



FTIR	Conforms to standard
Melting Point	183-193°C
UV/VIS: lambda max (H ₂ O)	279, (209)nm ± 3nm

Applications: 5-Bromo-2'-Deoxyuridine is a thymidine analog used as a mutagen in genetic research.
Recommended Storage: -20°C

8-Bromoadenine

Light Tan Powder

packaging	Mfr. No
100 mg AmberGlass	BP2507-100
500 mg AmberGlass	BP2507-500
1 g AmberGlass	BP2507-1

C ₅ C ₄ N ₆ Br CAS: 6974-78-3 MW: 213.97	EINECS: 230-225-1
---	-------------------

Applications: 8-Bromoadenine is used in several biochemical research applications.
Recommended Storage: Store below 0° C, desiccate
Not on TSCA inventory: for R and D use only; not for manufacturing or commercial purposes.

Acetic Acid, Glacial
Aldehyde-free

Sequencing

packaging	Mfr. No
500 mL Glass	BP1185-500
C ₂ H ₄ O ₂ CAS: 64-19-7 MW: 60.05 H226, H314	P280, P305+P351+P338, P302+P352, P304+P340, P210, P309+P311



Acetic Anhydride ((CH ₃ CO) ₂ O)	<=0.01%
Assay	≥99.7%
Chloride (Cl)	<=0.4ppm
Color (APHA)	<=10
Copper	<=0.1ppm
Dilution Test	To pass test
Heavy Metals (Pb)	<=0.5ppm
Impurities Reducing Dichromate	Not detected
Impurities Reducing Permanganate	Not detected
Iron	<=0.2ppm
Nickel	<=0.1ppm
Residue after evaporation	<=8ppm
Sulfate (SO ₄)	<=0.4ppm
Titrateable Base	<=0.0004mEq/g

Recommended storage: RT

Acetonitrile DNA Synthesis

packaging	Mfr. No
450 mL AmberGlass	BP1170-450
4 L AmberGlass	BP1170-4

C ₂ H ₃ N CAS: 75-05-8 MW: 41.04 EINECS: 200-835-2 H225, H319, H332, H302, H312	P280, P302+P352, P301+P312, P304+P340, P305+P351+P338, P210, P240
---	--



Assay (by GLC)	≥99.9%
Color (APHA)	<=10
Density (at 25°C)	0.775-0.780g/mL
Fluorescence Background (as Quinine Sulfate)	To pass test
IR	Conforms to standard
LC Gradient Suitability	To pass test
Optical Absorbance at 190nm	<=1.00
Optical Absorbance at 195nm	<=0.15
Optical Absorbance at 200nm	<=0.07
Optical Absorbance at 205nm	<=0.05
Optical Absorbance at 210nm	<=0.04
Optical Absorbance at 220nm	<=0.02
Optical Absorbance at 254nm	<=0.01
Refractive Index (at 25°C)	1.3405-1.3425
Residue after evaporation	<=1.0ppm
Titrateable Acid	<=0.008mEq/g
Titrateable Base	<=0.0006mEq/g
Water	<=10ppm

Also available in recyclable FisherPak® and NOWPak® containers.
Recommended storage: RT



Chloroform
Approx. 0.75% Ethanol as Preservative

Molecular Biology

packaging		Mfr. No
1 ℓ	AmberGlass,EcoSafPak*	BP1145-1
CHCl ₃	H302, H315, H351, H373	
CAS: 67-66-3	P281, P260, P301+P312,	
MW: 119.38	P302+P352	
EINECS: 200-663-8		
Acetone and Aldehyde		Not detected
Acid and Chloride		Not detected
Assay (GC)		>=99%
Color (APHA)		<=10
Density (at 25°C)		1.471-1.476g/ml
Fluorescence Background (as Quinine Sulfate)		To pass test
Free Chlorine		Not detected
Lead		<=0.05ppm
Optical Absorbance at 244nm		<=1.00
Optical Absorbance at 245nm		<=0.80
Optical Absorbance at 250nm		<=0.28
Optical Absorbance at 254nm		<=0.13
Peroxides		Not detected
Preservative (Ethanol)		0.5-1.0v/v%
Reactive Impurities with H ₂ SO ₄		Not detected
Refractive Index (at 25°C)		1.4420-1.4450
Residue after evaporation		<=2ppm
Water		<=0.02%

EcoSafPak* is an environmentally friendly packaging system made of 100% recyclable material by an SFI certified supplier.

Applications: Used in Phenol/Chloroform extractions to remove proteins from DNA or RNA samples.

Recommended storage: RT

Ethyl Acetate

Sequencing

packaging		Mfr. No
1 ℓ	AmberGlass,EcoSafPak*	BP1125-1
C ₄ H ₈ O ₂	H225, H319, H336, EUH066	
CAS: 141-78-6	P261, P280,	
MW: 88.11	P305+P351+P338, P210,	
EINECS: 205-500-4	P240	
Assay (GC)		>=99.9%
Color (APHA)		<=10
Density (at 25°C)		0.893 to 0.895g/ml
Fluorescence Background (as Quinine Sulfate)		To pass test
Optical Absorbance at 255nm		<=1.00
Optical Absorbance at 260nm		<=0.15
Optical Absorbance at 270nm		<=0.025
Peroxides		Not detected
Reactive Impurities with H ₂ SO ₄		Not detected
Refractive Index (at 25°C)		1.3680-1.3710
Residue after evaporation		<=10ppm
Titrateable Acid		<=0.0009mEq/g
Water		<=0.02%

Applications: Solvent that is commonly used in DNA or peptide synthesis.

Recommended storage: RT

EcoSafPak* is an environmentally friendly packaging system made of 100% recyclable material by an SFI certified supplier.

Formic Acid
Aldehyde-free

Sequencing

packaging		Mfr. No
500 mℓ	PolyBottle	BP1215-500
CH ₂ O ₂	P305+P351+P338,	
CAS: 64-18-6	P301+P330+P331,	
MW: 46.02	P302+P352, P210	
H226, H314		
P280, P260,		
Acetic acid glacial min.		<=0.4%
Aldehyde (HCHO)		Not detected
Ammonium (NH ₄)		<=0.005%
Assay (HCOOH w/w in H ₂ O)		>=88%
Chloride (Cl)		<=0.001%
Color (APHA)		<=15
Dilution Test		To pass test
Heavy Metals (Pb)		<=5ppm
Iron		<=5ppm
Residue after evaporation		<=0.002%
Sulfate (SO ₄)		<=0.002%
Sulfite		To pass test

Recommended storage: RT

Methylene Chloride
Pentene Preservative

packaging		Mfr. No
4 ℓ	AmberGlass	BP1186-4
CH ₂ Cl ₂	H351	
CAS: 75-09-2	P280	
MW: 84.93		
EINECS: 200-838-9		
Assay		>=99.9%
Boiling Range (1 drop to dryness)		<=1.0° incl. 39.8°C at 760mmHg
Color (APHA)		<=10
Density (at 25°C)		1.315-1.321g/ml
Fluorescence Background (as Quinine Sulfate)		<=1ppb
Free Halogens		To pass test
Optical Absorbance at 233nm		<=1.00
Optical Absorbance at 240nm		<=0.12
Optical Absorbance at 254nm		<=0.010
Refractive Index (at 25°C)		1.4200-1.4230
Residue after evaporation		<=1ppm
Titrateable Acid (as HCl)		<=0.0003mEq/g
Water		<=25ppm

Recommended storage: RT

Also available in recyclable FisherPak* and NOWPak* containers.

N,N-Dimethylformamide Synthesis

packaging		Mfr. No
500 mℓ	AmberGlass,EcoSafPak*	BP1160-500
4 ℓ	AmberGlass	BP1160-4
C ₃ H ₇ NO	P301+P310, P302+P352,	
CAS: 68-12-2	P304+P340, P280,	
MW: 73.09	P305+P351+P338, P210	
H319, H312, H332, H360D,		
H226		
Amines (as Dimethylamine)		<=5ppm
Assay (GC)		>=99.5%
Boiling Point		153.0° ±0.1°C
Boiling Range (1 drop to dryness)		2.0°C
Color (APHA)		<=15
Density (at 25°C)		0.942-0.946g/ml
IR		Conforms to standard
Optical Absorbance at 270nm		<=1.00
Optical Absorbance at 275nm		<=0.30
Optical Absorbance at 295nm		<=0.10
Optical Absorbance at 310nm		<=0.05
Optical Absorbance at 340-400nm		<=0.01
Residue after evaporation		<=0.005%
Titrateable Acid		<=0.0005mEq/g
Titrateable Base		<=0.003mEq/g
Water		<=0.03%

Recommended storage: RT

Also available in recyclable FisherPak*and NOWPak* containers. EcoSafPak* is an environmentally friendly packaging system made of 100% recyclable material by an SFI certified supplier.

n-Butyl Acetate

Sequencing

packaging		Mfr. No
500 mℓ	AmberGlass	BP1135-500
C ₆ H ₁₂ O ₂	H226, H336	
CAS: 123-86-4	P304+P340, P261, P210,	
MW: 116.16	P240	
EINECS: 204-658-1		
Assay (GC)		>=99.5%
Boiling Range (1 drop to dryness)		<=1.0°C
Color (APHA)		<=10
Density (at 25°C)		0.860-0.900 g/ml
Fluorescence Background (as Quinine Sulfate)		To pass test
Foreign Esters		Not detected
Optical Absorbance at 255nm		<=1.00
Optical Absorbance at 260nm		<=0.30
Optical Absorbance at 280nm		<=0.050
Reactive Impurities with H ₂ SO ₄		Not detected
Refractive Index (at 25°C)		1.3890-1.3950
Residue after evaporation		<=10 ppm
Titrateable Acid		<=0.003 meq/g
Water		<=0.05%

Applications: Solvent is commonly used in DNA or peptide synthesis.

Recommended storage: RT

n-Propanol
Peroxide-free

Sequencing

packaging		Mfr. No
500 mℓ	AmberGlass,EcoSafPak*	BP1130-500
C ₃ H ₈ O	P280, P305+P351+P338,	
CAS: 71-23-8	P310, P261, P301+P312,	
MW: 60.1	P210, P240,	
EINECS: 200-746-9	P304+P340	
H225, H336, H318		
Acidity (as Acetic Acid)		<=0.02%
Aldehyde (HCHO)		Not detected
Assay (GC)		>=99.5%
Boiling Range		96°-98°C
Color (APHA)		<=20
Density (at 25°C)		0.8016 ±0.005g/ml
Peroxides		Not detected
Residue after evaporation		<=0.01%

Also available in recyclable FisherPak* and NOWPak* containers. EcoSafPak* is an environmentally friendly packaging system made of 100% recyclable material by an SFI certified supplier.

Applications: Used for extracting ethidium bromide from DNA samples.

Recommended storage: RT

Pyridine

Sequencing

packaging		Mfr. No
500 mℓ	AmberGlass,EcoSafPak*	BP1155-500
C ₅ H ₅ N	H225, H312, H332, H302	
CAS: 110-86-1	P280, P302+P352,	
MW: 79.1	P304+P340, P210,	
EINECS: 203-809-9	P240	
Ammonia		<=0.002%
Assay (GC)		>=99.5%
Boiling Point		115.3° ±0.1°C
Boiling Range (1 drop to dryness)		<=2.0°C
Chloride (Cl)		<=0.001%
Copper		To pass test (about 5ppm)
Impurities Oxidized by Permanganate		Not detected
Residue after evaporation		<=0.002%
Solubility in H ₂ O		To pass test
Sulfate (SO ₄)		<=0.001%
Water		<=0.10%

Also available in recyclable FisherPak* and NOWPak* containers. EcoSafPak* is an environmentally friendly packaging system made of 100% recyclable material by an SFI certified supplier.

Applications: Commonly used in automated synthesis protocols.

Recommended storage: RT

PCR Nucleotide Mix		
10mM solution of each dATP, dCTP, dGTP, dTTP		
packaging		Mfr. No
200 µℓ	PolyMicroTube	BP2565-2
1 mL	PolyMicroTube	BP2565-1
Functional analysis		Pass Test

Applications: This premixed solution is used for PCR. Supplied as 10mM solution in purified water (pH 7.5). Functionally tested in PCR.
High purity:> 98% triphosphate.
[1927-31-7 (dATP)] ; [18423-43-3 (dTTP)] ; [93919-41-6 (dGTP)] ; [102783-51-7 (dCTP)]
Recommended Storage: -20°C
Not on TSCA inventory: for R and D use only; not for manufacturing or commercial purposes.

Mineral Oil	
Clear, Viscous Liquid	
packaging	Mfr. No
1 ℓ PolyBottle, EcoSafPak*	BP2629-1
CAS: 8012-95-1	
EINECS: 232-384-2	
Appearance	Clear, colorless, viscous liquid
DNase	None detected
Identification	Pass test
Limit of Polynuclear Compounds	Pass test
Neutrality	Pass test
Protease	None detected
Readily Carbonizable Substances	Pass test
RNase	None detected
Specific Gravity (at 25°C)	0.818-0.880
USP Protocol required	Pass test
Viscosity (at 40°C)	Not more than 33.5cs

Applications: Mineral Oil can be used in PCR applications.
Recommended Storage: RT

Taq DNA Polymerase		
Packaging		Mfr. No
100 units	Poly Tube	FB6000-10
5 x 100 units	Poly Tube	FB6000-15
500 units	Poly Tube	FB6000-25
5 x 500 units	Poly Tube	FB6000-30
2500 units	Poly Tube	FB6000-35
100 units	Poly Tube	FB6000-45
5 x 100 units	Poly Tube	FB6000-50
500 units	Poly Tube	FB6000-60
5 x 500 units	Poly Tube	FB6000-65
25 x 100 units	Poly Tube	FB6000-20
5 x 2500 units	Poly Tube	FB6000-40
25 x 100 units	Poly Tube	FB6000-55
2500 units	Poly Tube	FB6000-70
5 x 2500 units	Poly Tube	FB6000-75
CAS: 9012-90-2		

Tested for human genomic DNA using primers for the p53 gene.

Applications: Taq DNA Polymerase is a recombinant enzyme that is licensed for use in PCR. Provides superior performance in routine PCR and RT-PCR. Capable of amplifying large DNA target regions (up to 10kb). Well-suited for the following PCR applications:STR, RAPD, AFLP, and SNP analyses.
Description: Taq DNA Polymerase is a thermostable enzyme of approximately 94kDa size which was isolated originally from Thermus aquaticus. Replicates DNA optimally at 72°C and has a half-life of 40 minutes at 95°C allowing up to 45-50 thermal cycles without appreciable loss of enzyme activity.
The enzyme catalyzes the polymerization of nucleotides into double-stranded DNA in the 5'->3' direction in the presence of magnesium and it also has a 5'->3' exonuclease activity.
Source: Escherichia coli carrying a recombinant plasmid that encodes the T. aquaticus DNA polymerase gene.
Storage Buffer: 50mM Tris-HCl (pH 7.5), 0.1mM EDTA, 5mM DTT, added stabilizers, and 50% glycerol.
10X Buffer A: Supplied. Contains 500mM KCl, 15mM MgCl2, and 100mM Tris-HCl (pH 9.0 at 25°C).
10X Buffer B: Supplied. Contains 500mM KCl and 100mM Tris-HCl (pH 9.0 at 25°C). MgCl2 supplied separately.
Concentration: 5000 units/mL
Unit Definition: One unit catalyzes the incorporation of 10 nmole of total nucleotide into acid-insoluble product in 30 minutes at 70°C utilizing M13mp18 DNA as template.
Recommended Storage: -20°C

4-Chloro-1-naphthol Crystalline Powder		
packaging	Mfr. No	
50 g PolyBottle	BP1550-50	
C ₁₀ H ₇ ClO	P261, P302+P352, P280, P305+P351+P338	
CAS: 604-44-4		
MW: 178.62		
H335, H319, H315		
Assay (GC)	>=98%	
IR	Conforms to standard	
Melting Point	120°C ±3°C	

Applications: 4-Chloro-1-naphthol is a chromogenic substrate for Horseradish Peroxidase used in immunostaining, ELISA, and EIA.
Recommended Storage: -20°C



5-Bromo-4-chloro-3-indolyl Phosphate, p-Toluidine Salt White to Off-white Powder		
packaging	Mfr. No	
100 mg AmberGlass	BP1610-100	
500 mg AmberGlass	BP1610-500	
C ₁₅ H ₁₅ BrClN ₂ O ₆ P	EINECS: 229-506-1	
CAS: 6578-06-9		
MW: 433.61		
Assay	>=99%	
Lambda Max. (2mg in 100mℓ of 0.1N NaOH)	289nm ±2nm	
Protease	Not detected	
TLC	To pass test	

Applications: BCIP is used in conjunction with Nitro Blue Tetrazolium (NBT) in the colorimetric detection of Alkaline Phosphatase. It is suitable for use in immunohistochemistry, immunoblot staining, and ELISA applications.
Recommended Storage: <0°C

Bovine Serum Albumin (Fraction V)	
Cold-ethanol Precipitated	
packaging	Mfr. No
100g PolyBottle	BP1605-100
CAS: 9048-46-8	
Heavy Metals (Pb)	<=10ppm
pH (1% Solution in 0.9% NaCl)	7.0 ±0.3
Protease	Not Detected
Protein	>96%
Purity (albumin)	>98%
Solubility (0.4g/10mℓ H ₂ O)	Clear and haze-free
Sulfated ash	<=2%
Water Content	<5%

Applications: Cold-ethanol Precipitated Bovine Serum Albumin is used as a stabilizer for enzymes or enzymatic reactions, or as a blocker of nonspecific sites.
Recommended Storage: 4°C

Bovine Serum Albumin (Fraction V)	
Heat-shock Treated	
packaging	Mfr. No
100g PolyBottle	BP1600-100
CAS: 9048-46-8	
Heavy Metals (Pb)	<=10ppm
Loss on Drying (at 105°C)	<=5%
pH (1% Solution in 0.9% NaCl)	7.0 ±0.3
Protein	>96%
Purity (albumin)	>98%
Solubility (0.4g/10mℓ H ₂ O)	Clear and haze-free
Sulfated ash	<=2%
Water Content	<5%

Applications: Heat-shock Treated Bovine Serum Albumin is suitable for immunological studies.
Recommended Storage: 4°C

Bovine Serum Albumin (BSA)		Microbial Grade
packaging		Mfr. No
100 g Poly Bottle	BP9700-100	
CAS: 9048-46-8		
Appearance	Beige powder	
Heavy metals (as Pb)	<=20 ppm	
Loss on Drying (at 105°C)	<=6	
pH (1% in 0.15M NaCl)	Inclusive between 6.8-7.2	
Protein	>=95%	
Purity (albumin)	>=95%	
Sulfated ash	<=2%	

Applications: Immunological studies, Blocking non-specific sites
Description: Purified by a modified cold ethanol procedure.
Recommended storage: 4°C

Bovine Serum Albumin (BSA)		DNase- and Protease-Free Powder
Packaging		Mfr. No
100 g	Poly Bottle	BP8805-100
CAS: 9048-46-8		
Appearance	Pale Yellow Flakes	
Bioburden	<100 CFU/g	
DNase	Not detected	
pH (1% in 0.15M NaCl)	6.4-7.4	
Protease	Not detected	
Purity (albumin)	>=97%	
RNase	Not detected	

Applications: Immunological studies, Blocking non-specific sites
Description: Purified by proprietary heat shock process.
Recommended storage: 4°C

Bovine Serum Albumin (BSA)		Fatty Acid-Free Powder
packaging	Mfr. No	
100 g Poly Bottle	BP9704-100	
CAS: 9048-46-8		
Appearance	Beige to slightly yellow powder	
Free fatty acids	<=0.02%	
Loss on Drying (at 105°C)	<=5	
pH (1% in 0.15M NaCl)	Inclusive between 6.5-7.5	
Protein	>=96%	
Purity (albumin)	>=98%	
Sulfated ash	<=3%	

Applications: Immunological studies, Blocking non-specific sites
Description: Purified by a heat shock process with additional treatments; Fatty Acids: <=0.02%.
Recommended storage: 4°C

Bovine Serum Albumin (BSA)		Low Endotoxin Powder
packaging		Mfr. No
100 g Poly Bottle		BP9705-100
CAS: 9048-46-8		
Appearance	Beige to slightly yellow powder	
Endotoxin	<=1.0 EU/mg	
Heavy metals (as Pb)	<=20 ppm	
Loss on Drying (at 105°C)	<=5	
pH (1% in 0.15M NaCl)	Inclusive between 6.5-7.5	
Purity (albumin)	95%	
Sulfated ash	<=2%	

Applications: Immunological studies, Blocking non-specific sites
Description: Purified by a heat shock process with additional treatments;
Endotoxin: <=0.1 EU/mg.
Recommended storage: 4°C

Bovine Serum Albumin (BSA)		Protease-Free Powder
packaging		Mfr. No
100 g Poly Bottle		BP9703-100
CAS: 9048-46-8		
Appearance	Off-white crystalline powder	
Heavy metals (as Pb)	=<10 ppm	
Loss on Drying (at 105°C)	<=5	
pH (1% in 0.15M NaCl)	Inclusive between 6.8-7.2	
Protein	>=96%	
Purity (albumin)	>=98%	
Sulfated ash	<=2%	

Applications: Immunological studies, Blocking non-specific sites
Description: Purified by a heat shock process; Protease-Free.
Recommended storage: 4°C

Sulfosalicylic Acid Dihydrate		
Fine White Crystals		
packaging		Mfr. No
500 g	AmberGlass	BP177-500
C ₇ H ₆ O ₆ S.2H ₂ O	P301+P330+P331, P280,	
CAS: 5965-83-3	P305+P351+P338, P310,	
MW: 254.22	P303+P361+P353	
H314, H302		
Assay	99.0-101.0%	
Chloride (Cl)	<=0.001%	
Heavy Metals (Pb)	<=0.002%	
Insoluble matter	<=0.02%	
Iron	<=0.001%	
Residue after ignition	<=0.1%	
Salicylic Acid (HOC ₆ H ₄ COOH)	<=0.04%	
Sulfate (SO ₄)	To pass test (about 0.02%)	

Applications: Sulfosalicylic Acid is used for fixing proteins in agarose and polyacrylamide gels.
Recommended Storage: RT
UN 2585; DOT Class 8:Corrosive

p-Nitrophenyl Phosphate Disodium Salt, Hexahydrate	
packaging	Mfr. No
500 mg AmberGlass	BP2534-500
1 g AmberGlass	BP2534-1
5 g AmberGlass	BP2534-5
10 g AmberGlass	BP2534-10
50 g AmberGlass	BP2534-50
C ₆ H ₄ NNa ₂ O ₆ P.6H ₂ O	
EINECS: 224-246-5	
CAS: 4264-83-9	
MW: 371.13	
Assay	>99%
description	Pale Yellow Crystalline Powder
Free p-Nitrophenol	0.005%
Water Content (Hexahydrate)	27-31%

Applications: p-Nitrophenyl Phosphate is suitable for use as a substrate for alkaline and acid phosphatase.
Recommended Storage: Store below 0°C

2,2,2-Trifluoroethanol		Peptide Synthesis
Packaging		Mfr. No
100 mL	AmberGlass	BP622-100
C ₂ H ₃ F ₃ O	H312, H332, H361f	
CAS: 75-89-8	P280, P301+P351+P338,	
MW: 100.04	P310, P302+P352,	
EINECS: 200-913-6	P301+P312, P304+P340,	
H226, H315, H318, H302,	P210	
Assay	>=99%	
Boiling Range	74°-78°C	
Water	<=0.1%	

Applications: Trifluoroethanol is often used in peptide synthesis procedures.
Recommended Storage: RT
UN 1992; DOT Class 3:Flammable Liquid

Cacodylic Acid Crystalline	
packaging	Mfr. No
100 g PolyBottle/PoisonPack	BP324-100
C ₂ H ₇ AsO ₂	
CAS: 75-60-5	H331, H301, H410
MW: 138.01	P301+P310, P304+P340, P273
EINECS: 200-883-4	
Assay (dry basis)	>=99%
IR	Conforms to standard
Melting Point	195°-196°C
Solubility (10g/100mL H ₂ O)	Clear and colorless

Applications: Cacodylic Acid is used in buffers for DNA sequencing and recombinant DNA procedures.
Recommended Storage: RT
UN 1572; DOT Class 6.1:Poison

Dansyl Chloride Crystalline (Yellow)	
packaging	Mfr. No
1 g AmberGlass	BP427-1
C ₁₂ H ₁₂ ClNO ₂ S	
CAS: 605-65-2	P301+P330+P331, P280, P305+P351+P338, P310, P402+P404, P301+P312
MW: 269.75	
H302, H314, EUH029	
Assay	>=99.0%

Applications: Dansyl Chloride is a fluorochrome useful in the detection of N-terminal amino acids in proteins and peptides.
It is also used in the preparation of fluorescent derivatives of amino acids.
Recommended Storage: 0°C
UN 1759; DOT Class 8:Corrosive

N,N'-Diisopropylcarbodiimide		Peptide Synthesis
packaging		Mfr. No
100 mL	AmberGlass	BP590-100
C ₇ H ₁₄ N ₂	H330, H334, H318	   
CAS: 693-13-0	P305+P351+P338, P310,	
MW: 126.2	P260, P280, P302+P352,	
H226, H315, H317, H335,	P210	
Assay	≥99%	
Boiling Point	147°-149°C	
Color (APHA)	≤100	
IR	Conforms to standard	
Refractive Index (at 20°C)	1.4330 ±0.0015	

Applications: N,N'-Diisopropylcarbodiimide is used as an alternative to dicyclohexylcarbodiimide in peptide synthesis reactions.
Recommended Storage: RT
UN 1993; DOT Class 3:Flammable Liquid (n.o.s.)

N,N-Diisopropylethylamine		Peptide Synthesis
packaging		Mfr. No
500 mL AmberGlass		BP592-500
C ₈ H ₁₉ N	P301+P330+P331, P305+P351+P338, P280, P308+P313, P273, P210	
CAS: 7087-68-5		
MW: 129.24		
EINECS: 230-392-0		
H225, H302, H314, H412		
Assay		>=99%
Color (APHA)		<=20
IR		Conforms to standard
Refractive Index (at 25°C)		1.4133 ±0.002
Water (Karl Fischer)		<=0.2%

Applications: N,N-Diisopropylethylamine is commonly used in peptide synthesis reactions.
Recommended Storage: RT
UN 1993; DOT Class 3:Flammable Liquid (n.o.s.)

N-Methylmorpholine		Peptide Synthesis
packaging		Mfr. No
500 mL	AmberGlass	BP610-500
C ₄ H ₁₁ NO	P301+P330+P331, P305+P351+P338, P310, P302+P352, P304+P340, P210, P240, P280	
CAS: 109-02-4		
MW: 101.15		
EINECS: 203-640-0		
H225, H314, H312, H332, H302		
Assay (GC)		>=99.0%
Boiling Range (ACS)		114°-116°C

Applications: N-Methylmorpholine is used in peptide synthesis protocols.
Recommended Storage: RT
UN 2535; DOT Class 3:Flammable Liquid

p-Nitrophenol Yellow Crystals or Powder		Peptide Synthesis
packaging		Mfr. No
1 kg	AmberGlass	BP612-1
C ₆ H ₅ NO ₃	H373, H302, H312, H332	
CAS: 100-02-7	P280, P304+P340,	
MW: 139.11	P302+P352, P260	
EINECS: 202-811-7		
Assay (GC)		>=99%
Melting Point		112°-114°C

Applications: Nitrophenol is commonly used in peptide synthesis procedures.
Recommended Storage: RT
UN 1663; DOT Class 6.1:Poison

Sodium Cacodylate Trihydrate White Crystals or Crystalline Powder	
packaging	Mfr. No
50 g PolyBottle/PoisonPack	BP325-50
C ₂ H ₆ AsNaO ₂ .3H ₂ O	
CAS: 6131-99-3	P301+P310, P201, P308+P313, P304+P340, P273
MW: 214.03	
H331, H410, H301, H350	
Assay	>=97%
Loss on Drying (at 105°C)	25%±2%
pH of 5% Solution (at 25°C)	8.0 to 9.0
Solubility in H ₂ O (5g/100mL)	Clear and colorless

Applications: Sodium Cacodylate is a component of the DMS reaction buffer in Maxam-Gilbert sequencing and is also used as a buffer in terminal transferase reactions.
Recommended Storage: RT
UN 1688; DOT Class 6.1:Poison



Triethylamine		Peptide Synthesis	
packaging		Mfr. No	
500 mℓ	AmberGlass	BP616-500	
C ₆ H ₁₅ N		  	
CAS: 121-44-8	P280, P302+P350, P305+P351+P338,		
MW: 101.19	P304+P340, P308+P313,		
EINECS: 204-469-4	P261, P210		
H225, H311, H302, H332, H314, H335			
Assay			>=99%
Refractive Index (at 2°C)			1.400±0.003

Applications: Triethylamine is used as a solvent in automated synthesis protocols.
Recommended Storage: RT

Trifluoroacetic Acid		Peptide Synthesis Grade
packaging	Mfr. No	
500 mℓ AmberGlass	BP618-500	
C ₂ HF ₃ O ₂ CAS: 76-05-1 MW: 114.02 EINECS: 200-929-3 H314, H332, H412	P301+P330+P331, P280, P305+P351+P338, P310, P273, P304+P340	
		
		
Assay	>=99.5%	
Chloride (Cl)	<=5ppm	
Fluoride	<=14ppm	
Sulfate (SO ₄)	<=2ppm	
UV Absorbance at 275nm	<=0.06	
UV Absorbance at 300nm	<=0.03	
Water	<=0.02%	

Applications: Trifluoroacetic acid is used in peptide sythesis and sequencing protocols.

4-Methylumbelliferyl Phosphate Disodium Salt		Off-white Powder
packaging	Mfr. No	
25 mg AmberGlass	BP2666-250	
C ₁₆ H ₇ Na ₂ O ₆ P.3H ₂ O CAS: 22919-26-2 MW: 354.16	EINECS: 245-325-0	
FTIR	Conforms to Standard	
Loss on drying	16%	
Melting Point	193°-203°C	
Purity	>=97%	

Applications: 4-Methylumbelliferyl Phosphate is a biochemical substrate for membrane-bound phosphodiesterases.
Recommended Storage: Store below 0°C

α-Cyclodextrin		White Powder
packaging	Mfr. No	
5 g AmberGlass	BP2501-5	
25 g AmberGlass	BP2501-25	
C ₃₆ H ₆₀ O ₃₀ CAS: 10016-20-3 MW: 972.9	EINECS: 233-007-4	
Loss on drying	<=13%	
Optical Rotation (c = 1, H ₂ O)	+148° ±5°	

Applications: α-Cyclodextrin is used in chromatographic separations and purification methods. Useful for selective precipitation of enantiomeric, positional or structural isomers.
Recommended Storage: RT

β-Mercaptoethylamine Hydrochloride		
White Powder		
Packaging		Mfr. No
5 g	AmberGlass	BP2664-5
25 g	AmberGlass	BP2664-25
100 g	AmberGlass	BP2664-100
C ₂ H ₇ NS.HCl		
CAS: 156-57-0		
MW: 113.61		
EINECS: 205-858-1		
H315, H319, H302, H335		
Assay		98%
FTIR		Conforms to Standard
Melting Point		58°-73°C

Applications: β-Mercaptoethylamine HCl is used as an antioxidant in several biochemical procedures.
Recommended Storage: 4°C, desiccate

Forskolin From Coleus forskohlii		
Off-white Powder		
Packaging		Mfr. No
1 mg	Amberglass	BP2520-1
5 mg	AmberGlass	BP2520-5
10 mg	AmberGlass	BP2520-10
25 mg	AmberGlass	BP2520-25
50 mg	AmberGlass	BP2520-50
<chem>C22H34O7</chem> CAS: 66575-29-9 MW: 410.5 EINECS: 266-410-9		H312 P280, P302+P352
		

Applications: Forskolin is a molecular biology reagent. Forskolin functions as an antihypertensive and vasodilator. Adenylcyclase activator.
Recommended Storage: RT, desiccate
Not on TSCA inventory: for R and D use only; not for manufacturing or commercial purposes.

Polyvinylpyrrolidone		
packaging	Mfr. No	
100 g PolyBottle	BP431-100	
500 g PolyBottle	BP431-500	
(C ₆ H ₉ NO) _n CAS: 9003-39-8		
Arsenic	<=0.0002%	
Ash	<=0.05%	
Assay (dry basis)	>=95%	
Heavy Metals (Pb)	<=0.001%	
pH (5% in H ₂ O at Room Temperature)	3.0-7.0	
Viscosity (at 23°C)	1.2 to 1.3 Centistokes	

Applications: Polyvinyl Pyrrolidone is a component of Denhardt’s Reagent, commonly used in nucleic acid hybridization procedures.
Average Molecular Wt.:40.000
Recommended Storage: RT

Protease	
From <i>S. aureus</i> , strain V8	
packaging	Mfr. No
5 mg AmberGlass	BP2674-5
CAS: 9001-92-7	
Guaranteed analysis	
Activity	From <i>S. aureus</i> , strain V8
Applications: Protease is an enzyme that will cleave peptide bonds on the carboxyl-terminal side of either aspartate or glutamate residues.	
Unit Definition: One unit will change A280nm 0.001 per minute at 37°C, pH 7.8, using casein as the substrate.	
Activity:500 units/mg	
Recommended Storage: -20°C	

Pyridoxal 5-Phosphate Monohydrate		Off-white Crystalline Powder
packaging	Mfr. No	
1 g AmberGlass	BP2676-1	
5 g AmberGlass	BP2676-5	
25 g AmberGlass	BP2676-25	
C ₈ H ₁₀ NO ₆ P CAS: 54-47-7 MW: 247.14	EINECS: 200-208-3	
FTIR	Conforms to standard	
Purity	98%	

Applications: Pyridoxal 5-Phosphate is useful for modifying lysyl and valyl residues in proteins.
Recommended Storage: Store below 0°C

Sorbitol		
packaging	Mfr. No	
500 g AmberGlass	BP439-500	
C ₆ H ₁₄ O ₆ CAS: 50-70-4 MW: 182.17	EINECS: 200-061-5	
Arsenic	<=3ppm	
Assay	91.0-100.5%	
Chloride (Cl)	<=0.005%	
Heavy Metals (as Pb) (USP method)	<=0.001%	
Reducing Sugars	To pass test	
Sulfate (SO ₄)	<=0.01%	
Total Sugar	To pass test	

Applications: Sorbitol is used in the preparation of isoelectric focusing gels.
Recommended Storage: RT

Yeast DNA PFGE Markers		
packaging	Mfr. No	
5 blocks PolyMicroTube	BP2659-40	

Applications: For sizing DNA molecules from 225kb-1.9Mb. Consists of the 16 chromosomes of Saccharomyces cerevisiae YNN295. Supplied in 100μℓ agarose (0.5%) blocks in 1% lauroyl sarcosine, 0.5M EDTA (pH 9.0); each block contains 3μg DNA from 2 x 105 cells.
Recommended Storage: -20°C
Not on TSCA inventory: for R and D use only; not for manufacturing or commercial purposes.

Pyridoxal Hydrochloride		Off-white Powder
packaging	Mfr. No	
1 g AmberGlass	BP2675-1	
5 g AmberGlass	BP2675-5	
25 g AmberGlass	BP2675-25	
100 g AmberGlass	BP2675-100	
C ₈ H ₉ NO ₃ .HCl CAS: 65-22-5 MW: 203.62		
FTIR	Conforms to standard	
Melting Point	170°-180°C	
pH (0.1% aqueous solution)	2.7-3.5	
Solubility Test (0.1% aqueous solution, 15 min.)	Clear, colorless solution	
UV/VIS lambda max (water)	253-5, 290-1, 316-8nm	

Applications: Pyridoxal Hydrochloride is useful for labeling amino acids and for their detection in the picomole range.
Recommended Storage: 0°C, desiccate, protect from light.

Acrylamide Solution		Electrophoresis
40%		
packaging	Mfr. No	
1 ℓ AmberGlass	BP1402-1	
C ₃ H ₅ NO CAS: 79-06-1 MW: 71.08 H315, H319, H372, H311,	H350, H302, H340, H317, H361f P280, P312, P302+P350, P201, P308+P313, P260	
Acrylic Acid	<=0.003%	
Conductivity (at 25°C)	<=5μmhos/cm	
DNase	Not detected	
Electrophoresis	To pass test	
pH (at 25°C)	6.5 ±0.5	
Protease	Not detected	
RNase	Not detected	


Applications: 40% Acrylamide Solution is made from 3X recrystallized acrylamide and is suitable for electrophoretic separation of proteins and nucleic acids.
BP1402 is made from Acrylamide that meets BP170 specifications.
Recommended Storage: 4°C
UN 2810; DOT Class 6.1:Poison

Fisher BioReagents* EZ-Run* Protein Marker

Mixture of seven purified proteins optimized to produce uniform, compact bands when analyzed by SDS-PAGE and stained with Coomassie* Blue

- Supplied in gel loading buffer for direct application to SDS-PAGE gels
- Ideal for molecular weight determination of polypeptides analyzed in various gel systems

STORAGE CONDITIONS: Store at -20°C.



BP3600-500

Purified proteins contained in Marker		
Protein	Source	Molecular Weight (kDa)
β-galactosidase	<i>E. coli</i>	116.0
Bovine serum albumin	Bovine plasma	66.2
Ovalbumin	Chicken egg white	45.0
Lactate dehydrogenase	Porcine muscle	35.0
Restriction endonuclease <i>Bsp981</i>	<i>E. coli</i>	25.0
β-lactoglobulin	Bovine milk	18.4
Lysozyme	Chicken egg white	14.4

Molecular Weight Range	No. of Bands	Quantity	MFR No.
14.4 to 116.0 kDa	7	500μL (100 loadings) 2 × 500μL (200 loadings)	BP3600-500 BP3600-1

Fisher BioReagents* EZ-Run* Prestained Protein Marker

Ideal for monitoring protein separation during electrophoresis without staining and for verifying electroblotting of proteins onto membranes

- Comprised of six purified proteins covalently coupled to a blue chromophore
- Optimized to produce well-defined blue bands after SDS-PAGE
- Supplied in gel loading buffer


APPLICATIONS: The Prestained Protein Marker is designed for monitoring protein separation during electrophoresis without staining. It is also used to verify electroblotting of proteins onto membranes.

STORAGE CONDITIONS: Store at -20°C.

NOTE: The covalently coupled chromophore affects protein mobility; this prestained marker should only be used for approximating the molecular weight of separated proteins.

Purified proteins contained in Marker		
Protein	Source	Molecular Weight (kDa)
β-galactosidase	<i>E. coli</i>	118.0
Bovine serum albumin	Bovine plasma	87.0
Ovalbumin	Chicken egg white	49.0
Lactate dehydrogenase	Bovine erythrocytes	35.0
β-lactoglobulin	Bovine milk	26.0
Lysozyme	Chicken egg white	20.0

Molecular Weight Range	No. of Bands	Quantity	MFR No.
20 to 118 kDa	6	500μL (100 loadings) 2 × 500μL (200 loadings)	BP3601-500 BP3601-1




BP3601-1

Fisher BioReagents* EZ-Run* Rec Protein Ladder

Mixture of 14 highly purified recombinant proteins optimized for precise sizing of proteins by SDS-PAGE and Western blotting

- Ladder is ideal for precise molecular weight determination of polypeptides in denaturing PAGE and Western blots
- Each protein contains an integral Strep-tag* II sequence¹ and may be detected on Western blots using Strep-Tactin*²-AP conjugate² or an antibody against Strep-tag II
- 50kDa protein band has greater intensity than other bands to serve as a reference
- Supplied in gel loading buffer

STORAGE CONDITIONS: Store at -20°C.



BP3602 Series

Molecular Weight Range	No. of Bands	Quantity	MFR No.
10 to 200 kDa	14	500μL (100 loadings) 2 × 500μL (200 loadings)	BP3602-500 BP3602-1

¹ Strep-tag technology for protein purification and detection is covered by U.S. Patent No. 5,506,121; U.K. Patent No. 2,272,698; and French Patent No. 93 13 006.
² Strep-Tactin is covered by U.S. Patent No. 5,103,483.

Fisher BioReagents* EZ-Run* Prestained Rec Protein Ladder


Comprised of 10 recombinant proteins covalently coupled to a blue chromophore plus 10 and 72kDa reference bands tagged with green and orange dyes, respectively

- Ladder is ideal for monitoring protein separation during SDS-PAGE, verifying Western transfer efficiency, and approximating the molecular weight of blotted proteins
- Supplied in gel loading buffer

STORAGE CONDITIONS: Store at -20°C.

NOTE: The covalently coupled chromophore affects protein mobility; this prestained marker should only be used for approximating the molecular weight of separated proteins.

This product is manufactured under the license for Strep-tag* technology covered by U.S. Patent No. 5,506,121; U.K. Patent No. 2,272,698; and French Patent No. 93 13 006.



BP3603 Series

Molecular Weight Range	No. of Bands	Quantity	MFR No.
10 to 170 kDa	10	500μL (100 loadings) 2 × 500μL (200 loadings)	BP3603-500 BP3603-1

Fisher BioReagents* EZ-Run* Protein Gel Staining Solution

Ready-to-use solution for highly sensitive staining of proteins separated in polyacrylamide gels

EZ-Run Protein Gel Staining Solution contains Coomassie* Brilliant Blue G-250. Due to the property of the dye to form colloidal particles, proteins are preferentially stained without significant background staining of the gel matrix.

- Detection is linear over two orders of magnitude with as little as 5ng protein detected
- Sensitivity: 5ng BSA in 10% SDS-PAGE, 1.0mm thick gel, well width 3mm
- Does not contain toxic methanol or acetic acid
- Provides rapid destaining with simple water washes

- Can be reused twice without significant loss of sensitivity
- Tested for staining sensitivity in minigels using purified proteins and Protein Molecular Weight Standards

STORAGE CONDITIONS: Store at 4°C.

Quantity	MFR No.
1L	BP3620-1
4L	BP3620-4

Fisher BioReagents* EZ-Run* Protein Gel Solutions

Resolution equal to or better than gradient precast gels!

A ready-to-use premixed solution of acrylamide, bisacrylamide, buffer, and SDS that enables ultra-fine resolution of protein bands by denaturing PAGE. This novel formulation provides the resolution of a gradient gel with less preparative work than a standard discontinuous gel system.

- Ready-to-pour EZ-Run Protein Gel Solution for 1-D or 2-D PAGE—just add APS and TEMED
- Wide separation range—from 3 to 200kDa on same minigel
- Continuous gel system—no stacking gel necessary
- Proprietary (patent-pending) gel chemistry


APPLICATIONS: Compatible with post-electrophoresis applications such as Western blotting, MALDI analysis, and protein sequencing. Suitable for all commonly used gel staining methods including Coomassie* blue, silver staining and fluorescent dyes.

STORAGE CONDITIONS: Store at room temperature.

NOTE: EZ-Run Protein Gel Solution comes with an optimized proprietary running buffer; other traditional SDS-PAGE buffers should not be used with the EZ-Run gel solution.


Product Specifications	
Appearance	Clear, colorless solution
Conductivity (1:20, water)	
10% Gel Solution	460-520uS
12.5% Gel Solution	460-510uS
15% Gel Solution	440-500uS
Electrophoresis	To pass test
pH (at 25°C)	6.8-7.2
Polymerization Time	30 minutes
Protease	Not detected

Description	Quantity	MFR No.
10% EZ-Run Protein Gel Solution with buffer	30mL	BP7710-30
10% EZ-Run Protein Gel Solution with buffer	100mL	BP7710-100
10% EZ-Run Protein Gel Solution with buffer	500mL	BP7710-500
12.5% EZ-Run Protein Gel Solution with buffer	30mL	BP7712-30
12.5% EZ-Run Protein Gel Solution with buffer	100mL	BP7712-100
12.5% EZ-Run Protein Gel Solution with buffer	500mL	BP7712-500
15% EZ-Run Protein Gel Solution with buffer	100mL	BP7715-100
15% EZ-Run Protein Gel Solution with buffer	500mL	BP7715-500
EZ-Run Running Buffer 20X	500mL	BP7700-500



Resolution equal to or better than gradient precast gels!

EZ-Run Protein Gel Solution provides superior separation of closely spaced, small proteins (<20kDa) compared to a commercial gradient precast gel.



Separate wide range of protein sizes (3–200kDa) on the same minigel

The EZ-Run continuous gel system enables separation of small peptides and high MW proteins on the same minigel. For example, a commercial 12% precast discontinuous gel is not capable of resolving the 10 and 15kDa proteins compared to the 12% EZ-Run gel.



EZ-Run gel matrix compatible with common gel staining methods such as fluorescent dyes

Serial dilution of BSA (66kDa) and Ovalbumin (45kDa) are loaded in lanes 2 to 5 of an EZ-Run gel and detected with SYPRO* Ruby fluorescent protein stain. Protein standard in lane 1 is BP3602 EZ-Run Rec Protein Ladder.



Acrylamide		Electrophoresis
White, Needle-like Crystals		
packaging	Mfr. No	
100 g	PolyBottle	BP170-100
500 g	PolyBottle	BP170-500
5 kg	PolyPail	BP170-5
C ₃ H ₅ NO		P301+P310, P280,
CAS: 79-06-1		P302+P350, P201,
MW: 71.08		P308+P313, P260,
H315, H319, H301, H317,		P304+P340, P305+P351+P338
H332, H350, H340, H372,		
H361f, H312		
Acrylic Acid		<=0.002%
Appearance		To pass test
Arsenic		<=1 ppm
Assay		>=99.9%
Conductivity of a 50% Solution		<=2.50µmhos/cm
DNase		Not detected
Electrophoresis		To pass test
Insoluble in Methanol		<=0.005%
Insoluble in H ₂ O		<=0.005%
Iron		<=1 ppm
Lead		<=1 ppm
Mercury		<=1 ppm
pH of 10% Solution in 0.1M NaCl		6.5±0.5
Protease		Not detected
RNase		Not detected

Applications: This Acrylamide is 3X recrystallized to obtain a high-purity crystalline material suitable for electrophoretic separation of proteins and nucleic acids.
Recommended Storage: RT
UN 2074; DOT Class 6.1:Poison

Acrylamide: Bis-Acrylamide 19:1		Electrophoresis
Powder		
packaging	Mfr. No	
100 g	PolyBottle	BP1364-100
Conductivity (30% aqueous solution)		<=10µmhos/cm
DNase		Not detected
Electrophoresis		To pass test
Protease		Not detected
RNase		Not detected
Solubility (30% aqueous solution)		Clean and clear

Applications: This dry powder mix is premixed for 5% cross-linking and is suitable for DNA sequencing and separation of low-molecular-weight proteins. BP1364 is made from Acrylamide and Bis-Acrylamide that meet BP170 and BP171 specifications, including being DNase-, RNase-, and Protease-free. [79-06-1 (Acrylamide)] ; [110-26-9 (Methylene Diacrylamide)]
Recommended Storage: RT
UN 2074; DOT Class 6.1:Poison

Acrylamide: Bis-Acrylamide 19:1		Electrophoresis
40% Solution		
packaging	Mfr. No	
1 ℓ	PolyBottle	BP1406-1
H315, H319, H372, H311, H350, H302, H340, H317, H361f		
P280, P312, P302+P350, P201, P308+P313,		P260
Conductivity (at 25°C)		<=10µmhos/cm
DNase		Not detected
Electrophoresis		To pass test
Protease		Not detected
RNase		Not detected
Applications: This Acrylamide Solution is premixed for 5% cross-linking and is suitable for DNA sequencing and separation of low-molecular-weight proteins. BP1406 is made from Acrylamide and Bis-Acrylamide that meet BP170 and BP171 specifications, including being DNase-, RNase-, and Protease-free. [79-06-1 (Acrylamide)] ; [110-26-9 (Methylene Diacrylamide)]		
Recommended Storage: 4°C		
UN 2810; DOT Class 6.1:Poison		

Acrylamide: Bis-Acrylamide 29:1		Electrophoresis
Powder		
packaging	Mfr. No	
100 g	PolyBottle	BP1366-100
Conductivity (30% aqueous solution)		<=10µmhos/cm
DNase		Not detected
Electrophoresis		To pass test
Protease		Not detected
RNase		Not detected
Solubility (30% aqueous solution)		Clean and clear

Applications: This dry powder mix is premixed for 3.3% cross-linking and is suitable for DNA sequencing and separation of proteins. [79-06-1 (Acrylamide)] ; [110-26-9 (Methylene Diacrylamide)]
Recommended Storage: RT
UN 2074; DOT Class 6.1:Poison

Acrylamide: Bis-Acrylamide 29:1		Electrophoresis
40% Solution		
packaging	Mfr. No	
1 ℓ	PolyBottle	BP1408-1
H315, H319, H372, H311, H350, H302, H340, H317, H361f		
P280, P312, P302+P350, P201, P308+P313,		P260
Conductivity (at 25°C)		<=10µmhos/cm
DNase		Not detected
Electrophoresis		To pass test
Protease		Not detected
RNase		Not detected
Applications: This Acrylamide Solution is premixed for 3.3% cross-linking and is suitable for DNA sequencing and separation of proteins. BP1408 is made from Acrylamide and Bis-Acrylamide that meet BP170 and BP171 specifications, including being DNase-, RNase-, and Protease-free. [79-06-1 (Acrylamide)] ; [110-26-9 (Methylene Diacrylamide)]		
Recommended Storage: 4°C		
UN 2810; DOT Class 6.1:Poison		

Acrylamide: Bis-Acrylamide 37.5:1		Electrophoresis
Powder		
packaging	Mfr. No	
100 g	PolyBottle	BP1368-100
Conductivity (30% aqueous solution)		<=10µmhos/cm
DNase		Not detected
Electrophoresis		To pass test
Protease		Not detected
RNase		Not detected
Solubility (30% aqueous solution)		Clean and clear

Applications: This dry powder mix is premixed for 2.6% cross-linking and is suitable for separation of high-molecular-weight proteins. BP1368 is made from Acrylamide and Bis-Acrylamide that meet BP170 and BP171 specifications, including being DNase-, RNase-, and Protease-free. [79-06-1 (Acrylamide)] ; [110-26-9 (Methylene Diacrylamide)]
Recommended Storage: RT
UN 2074; DOT Class 6.1:Poison

Acrylamide: Bis-Acrylamide 37.5:1		Electrophoresis
40% Solution		
packaging	Mfr. No	
1 ℓ	PolyBottle	BP1410-1
H315, H319, H372, H311, H350, H302, H340, H317, H361f		
P280, P312, P302+P350, P201, P308+P313,		P260
Conductivity (at 25°C)		10µmhos/cm
DNase		Not detected
Electrophoresis		To pass test
Protease		Not detected
RNase		Not detected

Applications: This Acrylamide Solution is premixed for 2.6% cross-linking and is suitable for separation of high-molecular-weight proteins. BP1410 is made from Acrylamide and Bis-Acrylamide that meet BP170 and BP171 specifications, including being DNase-, RNase-, and Protease-free. [79-06-1 (Acrylamide)] ; [110-26-9 (Methylene Diacrylamide)]
Recommended Storage: 4°C
UN 2810; DOT Class 6.1:Poison

Agarose		IEF Grade
Isoelectric Focusing of Proteins		
packaging	Mfr. No	
10 g	PolyBottle	BP163-10
25 g	PolyBottle	BP163-25
C ₁₂ H ₁₈ O ₉		
CAS: 9012-36-6		
MW: 306.12		
EEO (-Mr)		Not detectable
Gel Strength		>500g/cm2
Gelation Temperature		<45°C
Moisture Content		<10%
Sulfate Content		<0.2%

Applications: Highly purified agarose for separating proteins by isoelectric focusing.
Recommended Storage: RT

Ammonium Persulfate		Electrophoresis
Colorless-to-white Crystals		
packaging	Mfr. No	
25 g	AmberGlass	BP179-25
100 g	AmberGlass	BP179-100
H ₈ N ₂ O ₈ S ₂		H317, H319, H302
CAS: 7727-54-0		P342+P311, P280,
MW: 228.19		P301+P312, P302+P352,
EINECS: 231-786-5		P305+P351+P338, P210
H272, H334, H335, H315,		
Assay		>=98.0%
Chloride and Chlorate (as Cl)		<=0.001%
Heavy Metals (Pb)		<=0.005%
Insoluble matter		<=0.005%
Iron		<=0.001%
Manganese		<=0.5ppm
Residue after ignition		<=0.05%
Titrateable Free Acid		<=0.04mEq/g

Applications: Electrophoresis-grade Ammonium Persulfate is an initiator of acrylamide polymerization for the preparation of polyacrylamide electrophoresis gels.
Recommended Storage: RT
UN 1444; DOT Class 5.1:Oxidizer

Ammonium Sulfate		
packaging	Mfr. No	
1 kg	Glass	BP212R-1
2.5 kg	PolyBottle	BP212-212
H ₈ N ₂ O ₄ S		
CAS: 7783-20-2		
MW: 132.13		
Arsenic		<=0.5ppm
Assay		>=99.5%
Heavy Metals (Pb)		<=5ppm
Iron		<=5ppm
Nitrate		<=0.001%
Optical Absorbance at 260nm		<=0.03
Optical Absorbance at 280nm		<=0.025
pH of 5% Solution (at 25°C)		5 to 6
Protease		Not detected
Residue on Ignition (sulfated)		<=0.005%

Applications: This enzyme-grade Ammonium Sulfate is suitable for use in the purification of proteins.
Recommended Storage: RT

Bis-Acrylamide Solution		Electrophoresis
2% w/v		
packaging	Mfr. No	
250 mL	AmberGlass	BP1404-250
C ₇ H ₁₀ N ₂ O ₂		EINECS: 203-750-9
CAS: 110-26-9		
MW: 154.17		
Electrophoresis		To pass test

Applications: Bis-Acrylamide Solution is a cross-linker commonly used in the preparation of polyacrylamide electrophoresis gels. BP1404 is made from Bis-Acrylamide that meets BP171 specifications.
Recommended Storage: 4°C

Bis-Acrylamide		Electrophoresis	
packaging		Mfr. No	
25 g	PolyBottle	BP171-25	
100 g	PolyBottle	BP171-100	
C ₇ H ₁₀ N ₂ O ₂ CAS: 110-26-9 MW: 154.17 EINECS: 203-750-9 H319, H315, H335, H302,	H312, H332 P280, P302+P352, P305+P351+P338, P261, P270, P308+P313		
Absorbance of a 1% Solution (at 290nm)			<=0.2
Assay			>=95%
DNase			Not detected
Electrophoresis			To pass test
pH of 1% Solution in 0.1M NaCl			>=5.0
Protease			Not detected
RNase			Not detected

Applications: Bis-Acrylamide is a cross-linker commonly used in the preparation of polyacrylamide electrophoresis gels.
Recommended Storage: 4°C

Brilliant Blue G-250	
packaging	Mfr. No
25 g AmberGlass	BP100-25
50 g AmberGlass	BP100-50
C ₂₇ H ₄₈ N ₃ NaO ₇ S ₂ CAS: 6104-58-1 MW: 854.02	
E1% _{1cm}	>=600g-1cm-1
Lambda Max.	610nm ±2nm
Loss on Drying (at 105°C)	<=5%
Solubility	To pass test

Applications: General protein stain.
Recommended storage: RT
C.I. 42655

Brilliant Blue R-250	
packaging	Mfr. No
25 g mberGlass	BP101-25
50 g AmberGlass	BP101-50
C ₄₅ H ₄₄ N ₃ NaO ₇ S ₂ CAS: 6104-59-2 MW: 825.96	
E1% _{1cm}	>=700g-1cm-1
Lambda Max.	592nm ±2nm
Loss on Drying (at 105°C)	<=5%
Solubility	To pass test

Applications: Staining proteins in polyacrylamide gels.
Recommended storage: RT
C.I. 42660

Glycine		White Crystals or Crystalline Powder
packaging	Mfr. No	
500 g PolyBottle	BP381-500	
1 kg PolyBottle	BP381-1	
5 kg PolyPail	BP381-5	
C ₂ H ₅ NO ₂ CAS: 56-40-6 MW: 75.07	EINECS: 200-272-2	
Ammonium (NH ₄)	<=0.02%	
Arsenic	<=1ppm	
Assay	98.5-101.0%	
Chloride (Cl)	<=0.003%	
Color of 10% Solution	<=10	
Conductivity of 1.5% Solution	<=25µmhos/cm	
Heavy Metals (Pb)	<=10ppm	
Iron	<=7ppm	
Loss on drying	<=0.20%	
Optical Absorbance of 1M Solution (at 280nm)	<=0.1 AU	
Other Amino Acids	Chromatographically not detectable	
Residue on Ignition (sulfated)	<=0.10%	

Applications: Glycine is used in the preparation of Tris-Glycine electrophoresis buffers and culture media.
Recommended Storage: RT

Lithium Dodecyl Sulfate		Molecular Biology
Powder or Flakes		
packaging	Mfr. No	
5 g	PolyBottle	BP1309-5
25 g	PolyBottle	BP1309-25
C ₁₂ H ₂₅ LiO ₄ S	P261, P301+P312,	 
CAS: 2044-56-6	P304+P340, P302+P352,	
MW: 272.33	P305+P351+P338, P233,	
H228, H335, H302, H315,	P210	
H332, H319		
Absorbance of a 3% Solution at 230nm	<=0.3	
Absorbance of a 3% Solution at 280nm	<=0.05	
Absorbance of a 3% Solution at 405nm	<=0.01	
Assay	>=99.0%	
Heavy Metals (Pb)	<=5ppm	
pH (at 25°C, 3% aqueous solution)	7.0±1.0	
Phosphate (PO ₄)	<=1ppm	

Recommended storage: RT

MES		
Fine White Crystals		
packaging	Mfr. No	
100 g PolyBottle	BP300-100	
C ₆ H ₁₃ NO ₄ S	H315, H335, H319	
CAS: 4432-31-9	P280, P305+P351+P338,	
MW: 195.24	P261, P302+P352, P264,	
EINECS: 224-632-3	P271	
Assay	>=98%	
DNase	Not detected	
Moisture	<=10%	
pKa (at 25°C)	6.09±0.2	
RNase	Not detected	

Applications: This biological buffer has a usable pH range of 5.5 to 6.7.
Recommended Storage: RT

Ponceau S	
packaging	Mfr. No
10 g AmberGlass	BP103-10
C ₂₂ H ₁₂ N ₄ Na ₄ O ₁₃ S ₄ CAS: 6226-79-5 MW: 760.54 H319, H335, H315	P261, P302+P352, P280, P305+P351+P338
E1% _{1cm}	>=350 g-1cm-1 min
Lambda Max.	514nm±2nm
Solubility	To pass test

Applications: Staining serum and plasma proteins on nitrocellulose and cellulose acetate membranes

Potassium Persulfate		
Fine, Colorless-to-white Crystals		
packaging		Mfr. No
100 g	PolyBottle	BP180-100
K₂O₈S₂ CAS: 7727-21-1 MW: 270.3 EINECS: 231-781-8 H272, H334, H302, H319, H335, H315, H317		P304+P341, P342+P311, P261, P301+P312, P302+P352, P280, P305+P351+P338, P210
		
		
		
Assay		>=99.0%
Chlorine Compounds (as Cl)		<=0.001%
Heavy Metals (Pb)		<=0.001%
Insoluble matter		<=0.005%
Iron		<=5ppm
Manganese		<=2ppm

Applications: Potassium Persulfate is a strong oxidizing agent and initiator of acrylamide polymerization. It can also be used as an initiator in the immobilization of cells and enzymes in an acrylamide matrix.
Recommended Storage: RT
UN 1492; DOT Class 5.1:Oxidizer

Protein Gel-Loading Dye for SDS-PAGE		Electrophoresis
Dark-purple Liquid - 2X		
packaging	Mfr. No	
1 mL PolyTube	BP637-1	
5 mL PolyBottle	BP637-5	
DNase	Not detected	
Optical Absorbance (dilution 1250 with deionized water) at 525nm	0.15-0.35	
Optical Absorbance (dilution 1250 with deionized water) at 588-594nm	0.32-0.47	
Optical Absorbance (dilution 1250 with deionized water) at 635-641nm	0.25-0.40	
Protease	Not detected	
RNase	Not detected	

Applications: This dye is added to protein samples prior to electrophoresis on SDS-PAGE gels.
Components: Sodium Dodecyl Sulfate (<5%), Glycerol (<12%), Tris (<1.0%), Proprietary Component (<2%), and Water (80%).
[151-21-3 (Sodium Dodecyl Sulfate)] ; [56-81-5 (Glycerol)] ; [126-72-7 (Tris)]
Recommended Storage: RT

Sodium Dodecyl Sulfate		Electrophoresis
White Powder		
packaging	Mfr. No	
100 g PolyBottle	BP166-100	
500 g PolyBottle	BP166-500	
5 kg PolyPail	BP166-5	
C ₁₂ H ₂₅ NaO ₄ S	H315, H319	 
CAS: 151-21-3	P280, P302+P352, P210,	
MW: 288.38	P270, P301+P312,	
EINECS: 205-788-1	P305+P351+P338	
H228, H311, H302, H335,		
Absorbance of a 0.1% Solution at 230nm	<=0.4	
Absorbance of a 0.1% Solution at 280nm	0.1	
Assay (as Fatty Alcohol Sulfate)	>=99%	
Assay (as C ₁₂ -OH)	98%	
DNase	Not detected	
Electrophoresis	To pass test	
Lead	<=5ppm	
Phosphate (PO ₄)	<=0.0001%	
RNase	Not detected	

Applications: Sodium Dodecyl Sulfate (SDS) is the most commonly used detergent in protein purification and electrophoresis. It is used as a protein denaturing agent during both extraction procedures and SDS polyacrylamide gel electrophoresis.
Recommended Storage: RT

Sodium Dodecyl Sulfate		Electrophoresis
10% solution		
packaging	Mfr. No	
200 mL PolyBottle	BP2436-200	
1 L PolyBottle	BP2436-1	
C ₁₂ H ₂₅ NaO ₄ S	H319, H335, H315	
CAS: 151-21-3	P280, P261, P302+P352,	
MW: 288.38	P305+P351+P338	
EINECS: 205-788-1		
Assay (as Fatty Alcohol Sulfate)		>=99%
Chloride (Cl)		<=0.1%
Concentration		10.0±0.5%
Copper		<=0.0005%
DNase		Not detected
Lead		<=0.0005%
Optical Absorbance of a 0.1% Solution at 230nm		<=0.4
Optical Absorbance of a 0.1% Solution at 280nm		<=0.1
Phosphate (PO ₄)		<=0.0005%
Protease		Not detected
RNase		Not detected

Applications: Sodium Dodecyl Sulfate (SDS) is the most commonly used detergent in protein purification and electrophoresis. It is used as a protein denaturing agent during both extraction procedures and SDS polyacrylamide gel electrophoresis.
Recommended Storage: RT
Filtered through a 0.2-micron filter.



Sodium Dodecyl Sulfate Electrophoresis/Molecular Biology
20% Solution

packaging		Mfr. No
200 mℓ	PolyBottle	BP1311-200
1 ℓ	PolyBottle	BP1311-1
C ₁₂ H ₂₅ NaO ₄ S	H319, H335, H315	
CAS: 151-21-3	P280, P261, P302+P352,	
MW: 288.38	P305+P351+P338	
EINECS: 205-788-1		
Assay (as Fatty Alcohol Sulfate)		>=99%
Chloride (Cl)		<=0.1%
Concentration		20 ±0.5%
Copper		<=0.0005%
DNase		Not detected
Lead		<=0.0005%
Optical Absorbance of a 0.1% Solution at 230nm		<=0.4
Optical Absorbance of a 0.1% Solution at 280nm		<=0.1
Phosphate (PO ₄)		<=0.0005%
Protease		Not detected
RNase		Not detected

Filtered through a 0.2-micron filter.
Recommended storage: RT

Tris-Glycine-SDS Electrophoresis
Dry Powder Mix of Tris/Glycine/SDS

packaging		Mfr. No
92 g	PolyBottle	BP1398-92
H335, H319, H315		
P261, P264, P302+P352, P280,		
P305+P351+P338		
Electrophoresis	To pass test	
pH of 1X Buffer (at 25°C)	8.3 ±0.2	
Solubility (buffer in final volume of 1 liter in H2O)	Clear solution	

Components in 1X buffer: 0.025M Tris, 0.192M glycine, and 0.1% SDS.
Components: Glycine (82.3%), Tris (17.1%), and SDS (0.6%) [56-40-6 (Glycine)] ; [77-86-1 (Tris)] ; [151-21-3 (SDS)]
Recommended Storage: RT
Makes 1 liter of 5X buffer. Final concentrations of components in 1X buffer: 0.025M Tris, 0.192M glycine, and 0.1% SDS.

Silver Nitrate Certified ACS
Crystalline

packaging		Mfr. No
100 g	AmberGlass	BP2546-100
AgNO ₃		
CAS: 7761-88-8		
MW: 169.87		
Assay	>=99.7%	
Chloride (Cl)	<=5ppm	
Clarity of Solution	To pass test	
Copper	<=2ppm	
Free Acid	To pass test	
Iron	<=2ppm	
Lead	<=0.001%	
Substances not Precipitated by HCl	<=0.01%	
Sulfate (SO ₄)	<=0.002%	

Applications: Silver Nitrate is used in biochemistry to stain agarose and acrylamide gels.
Recommended Storage: RT
UN 1493; DOT Class 5.1:Oxidizer

Sodium Persulfate, >98% Electrophoresis
White Crystalline Powder

Packaging		Mfr. No
1 kg	PolyBottle	BP2637-1
Na ₂ O ₈ S ₂	P261, P280, P308+P313,	  
CAS: 7775-27-1	P302+P352, P210,	
MW: 238.09	P305+P351+P338	
H272, H315, H302, H317,		
H334, H319, H335		
Assay	>=98.0%	

Applications: Sodium Persulfate is used as a promoter for polymerization reactions.
Recommended Storage: RT
UN 1505; DOT Class 5.1:Oxidizer

Sucrose Gel-Loading Dye for DNA Gels Electrophoresis
5X, Contains 40% Sucrose

packaging		Mfr. No
1 mℓ	PolyTube	BP655-1
5 mℓ	PolyBottle	BP655-5
DNase	Not detected	
Electrophoresis	To pass test	
Optical Absorbance (dilution 1:500 with deionized water) at 525nm	0.30-0.45	
Optical Absorbance (dilution 1:500 with deionized water) at 588-594nm	0.50-0.65	
Optical Absorbance (dilution 1:500 with deionized water) at 635-641nm	0.47-0.62	
RNase	Not detected	

Applications: This dye includes sucrose instead of glycerol. Add the dye to DNA samples prior to electrophoresis on agarose gels.
Components: Sucrose (40%), Proprietary Component (<3%), and Water. [57-50-1 (Sucrose)]
Recommended Storage: RT

TEMED Electrophoresis

packaging		Mfr. No
20 g	AmberGlass	BP150-20
100 g	AmberGlass	BP150-100
C ₆ H ₁₆ N ₂	EINECS: 203-744-6	
CAS: 110-18-9		
MW: 116.21		
Assay (by titration)	>=97%	
DNase	Not detected	
Electrophoresis	To pass test	
Protease	Not detected	
RNase	Not detected	

Applications: TEMED is a commonly used catalyst for the polymerization of acrylamide/bis-acrylamide.
Recommended Storage: RT
UN 2372; DOT Class 3:Flammable Liquid

Tris-Glycine-SDS Electrophoresis
10X Powder

packaging		Mfr. No
1	FoilPack	BP1342-1
H315, H319 P280, P305+P351+P338		
Electrophoresis	To pass test	
pH (1X solution) (at 25°C)	8.4±0.1	
Protease	Not detected	

Applications: Tris-Glycine-SDS is commonly used as a buffer in SDS-PAGE. Each pack contains preweighed powder to make 1ℓ of a 10X solution (0.25M Tris Base, 1.92M Glycine, and 1.0% w/v SDS).
[77-86-1 (Tris)] ; [56-40-6 (Glycine)] ; [151-21-3 (SDS)]
Recommended Storage: RT

Tris-Glycine Electrophoresis
10X Powder

packaging		Mfr. No
1	FoilPack	BP1307-1
CAS: [77-86-1 (Tris)], [56-40-6 (Glycine)]		
H315, H319		
P305+P351+P338, P280		
DNase		Not detected
Electrophoresis		To pass test
pH (1X solution) (at 25°C)		8.4±0.1
Protease		Not detected
RNase		Not detected

Recommended storage: RT

Tris-Glycine Electrophoresis
10X Solution

packaging		Mfr. No
1 ℓ	PolyBottle	BP1306-1
4 ℓ	PolyPac*	BP1306-4
CAS: [77-86-1 (Tris)], [56-40-6 (Glycine)]		
DNase	Not detected	
Electrophoresis	To pass test	
pH (1X solution) (at 25°C)	8.3-8.5	
Protease	Not detected	
RNase	Not detected	

Filtered through a 0.2-micron filter.
Recommended storage: RT

Tris-Glycine Electrophoresis
1X Solution

packaging		Mfr. No
4 ℓ	PolyPac*	BP2439-4
20 ℓ	PolyPac*	BP2439-20
DNase	Not detected	
Electrophoresis	To pass test	
pH (at 25°C)	8.4±0.1	
Protease	Not detected	
RNase	Not detected	

Applications: Tris-Glycine is commonly used as a buffer in native polyacrylamide gel electrophoresis.
0.025M Tris Base and 0.192M Glycine [77-86-1 (Tris)] ; [56-40-6 (Glycine)]
Recommended Storage: RT
Filtered through a 0.2-micron filter.

Tris-Glycine-SDS Electrophoresis
10X Solution

packaging		Mfr. No
1 ℓ	PolyBottle	BP1341-1
4 ℓ	PolyPac*	BP1341-4
DNase	Not detected	
Electrophoresis	To pass test	
pH (1X solution) (at 25°C)	8.4±0.1	
Protease	Not detected	
RNase	Not detected	

Applications: Tris-Glycine-SDS is commonly used as a buffer in SDS-PAGE.
Components: 0.25M Tris Base, 1.92M Glycine, and 1.0% (w/v) SDS. [77-86-1 (Tris)] ; [56-40-6 (Glycine)] ; [151-21-3 (SDS)]
Recommended Storage: RT
Filtered through a 0.2-micron filter.

Tris-Glycine-SDS Electrophoresis
1X Solution

packaging		Mfr. No
4 ℓ	PolyPac*	BP2440-4
20 ℓ	PolyPac*	BP2440-20
DNase	Not detected	
Electrophoresis	To pass test	
pH of 1X solution (at 25°C)	8.3±0.1	
Protease	Not detected	
RNase	Not detected	

Applications: Tris-Glycine-SDS is commonly used as a buffer in SDS-PAGE. 0.025M Tris Base, 0.192M Glycine, and 0.1% (w/v) SDS. [56-40-6 (Glycine)] ; [77-86-1 (Tris)] ; [151-21-3 (SDS)]
Recommended Storage: RT
Filtered through a 0.2-micron filter.

Tween* 20		
packaging	Mfr. No	
100 mℓ PolyBottle	BP337-100	
500 mℓ PolyBottle,EcoSafPak*	BP337-500	
CAS: 9005-64-5		
Arsenic	<=0.0003%	
Arsenic (As)	<=0.0003%	
Copper	<=0.0005%	
Copper (Cu)	<=0.0005%	
Hydroxyl Number	96-108	
Hydroxyl number	96-108	
Lead	<=0.0005%	
Lead (Pb)	<=0.0005%	
Protease	Not detected	
Residue after ignition	<=0.25%	
Water	<=3%	

Applications: This nonionic detergent is used in the isolation of water-soluble membrane proteins.
Recommended Storage: RT
EcoSafPak* is an environmentally friendly packaging system made of 100% recyclable material by an SFI certified supplier.

Tween* 80		
packaging		Mfr. No
500 mℓ	AmberGlass,EcoSafPak*	BP338-500
CAS: 9005-65-6		
Arsenic		<=0.0001%
Copper		<=0.0005%
Hydroxyl Number		65-80
Lead		<=0.0005%
Residue after ignition		<=0.25%
Water		<=3.0%

Applications: This nonionic detergent is used in the isolation of water-soluble membrane proteins.
Recommended Storage: RT
EcoSafPak* is an environmentally friendly packaging system made of 100% recyclable material by an SFI certified supplier.

Urea		Molecular Biology
Colorless-to-white Crystals or Crystalline Powder		
packaging	Mfr. No	
500 g PolyBottle	BP169-500	
2.5 kg PolyBottle	BP169-212	
10 kg PolyPail	BP169-10	
CH ₄ N ₂ O	H315	
CAS: 57-13-6	P302+P352	
MW: 60.06		
EINECS: 200-315-5		
Absorbance of a 5M Solution at 280nm	<=0.05	
Assay	>=99.0%	
Biuret Test	Negative	
CNO	<=5ppm	
Copper	<=1ppm	
DNase	Not detected	
Heavy Metals (Pb)	<=5ppm	
Iron	<=2ppm	
Melting Point	132°-135°C	
Protease	Not detected	
RNase	Not detected	
Solubility of an 8M Solution	To pass test	

Applications: Urea is used as a denaturing agent for DNA and proteins.
Recommended Storage: RT

2-Mercaptoethanol		Electrophoresis
packaging	Mfr. No	
100 g AmberGlass/PoisonPack	BP176-100	
C ₂ H ₆ OS	P280, P312, P302+P350, P261, P273, P301+P312, P305+P351+P338	
CAS: 60-24-2		
MW: 78.13		
H302, H411, H315, H335, H311, H319		
Assay	>=98%	
DNase	Not detected	
Electrophoresis	To pass test	
RNase	Not detected	
Total Disulfide (S-S) Absorbance of a 0.02M Solution (at 285nm)	<=0.05	

Applications: 2-Mercaptoethanol reduces disulfide linkages in proteins and peptides and is used in many enzymatic reactions as a protein stabilizer.
Recommended Storage: RT
UN 2966; DOT Class 6.1:Poison

4-(2-Aminoethyl)benzenesulfonyl Fluoride Hydrochloride		
White Powder		
packaging	Mfr. No	
10 mg AmberGlass	BP2644-10	
50 mg AmberGlass	BP2644-50	
100 mg AmberGlass	BP2644-100	
500 mg AmberGlass	BP2644-500	
1 g AmberGlass	BP2644-1	
C ₈ H ₁₀ FO ₂ S.CIH	P280, P301+P330+P331, P310, P305+P351+P338, P303+P361+P353, P210	
CAS: 30827-99-7		
MW: 239.7		
H314		
FTIR	Conforms to reference	
Melting Point	175°-185°C	
Purity (by HPLC)	>=97%	

Applications: 4-(2-Aminoethyl)benzenesulfonyl Fluoride Hydrochloride is a covalent binding serine protease inhibitor. Similar in structure to the commonly used inhibitor PMSF, AEBSF offers better solubility in water, lower toxicity, greater stability and higher inhibitory activity. AEBSF belongs to the family of irreversible sulfonyl fluoride inhibitors that block trypsin and chymotrypsin-type enzymes.
Recommended Storage: RT, protect from light.
Not on TSCA inventory: for R and D use only; not for manufacturing or commercial purposes.

AEBSF Hydrochloride		
White to Off-white Crystals		
packaging	Mfr. No	
50 mg PolyMicroTube	BP635-50	
100 mg PolyMicroTube	BP635-100	
500 mg PolyMicroTube	BP635-500	
C ₈ H ₁₀ FO ₂ S.CIH	P280, P301+P330+P331, P310, P305+P351+P338, P303+P361+P353, P210	
CAS: 30827-99-7		
MW: 239.7		
H314		
Assay (TLC, HPLC)	>=98.0%	
Carbon	40.09% ±0.5%	
Hydrogen	4.62% ±0.04%	
Melting Point	182.0°-190.0°C	
Nitrogen	5.84% ±0.2%	
Solubility	TIUTrypsin Inhibitor Units determined using Pefachrome TRY.	
Solubility 20g in 100mℓ Water	To pass test	
Solubility 7.5g in 100mℓ Absolute Alcohol	To pass test	
Solubility Inhibitor Activity	>250TIU/g	
Sulfur	13.38% ±0.1%	

Applications: This specific, irreversible inhibitor of serine proteases is a stable, nontoxic alternative to PMSF.
Recommended Storage: RT
Not on TSCA inventory: for R and D use only; not for manufacturing or commercial purposes.

Benzamidine Hydrochloride Hydrate		
White Crystalline Powder		
packaging	Mfr. No	
25 g PolyBottle	BP435-25	
1 kg AmberGlass	BP435-1	
C ₇ H ₈ N ₂ .HCl	H335, H315, H319	
CAS: 1670-14-0	P261, P302+P352, P280, P305+P351+P338	
MW: 156.61		
EINECS: 216-795-4		
Assay	>=99.0%	
Heavy Metals (Pb)	<=0.005%	
Lambda Max. (in Methanol)	230nm ±5nm	
Solubility (1% in H ₂ O at Room Temp.)	Clear and haze-free	
Water (internal)	<=15%	

Applications: Benzamidine Hydrochloride is commonly used as a protease inhibitor.
Recommended Storage: 4°C

Brij* 35		
packaging	Mfr. No	
500 g AmberGlass	BP345-500	
C ₅₈ H ₁₁₈ O ₂₄	P280, P305+P351+P338, P330, P264, P301+P312	
CAS: 9002-92-0		
MW: 1199.56		
H302, H315, H319, H335		
Hydroxyl Number	40-60	
Water	<=3.0%	

Applications: This nonionic detergent is used for solubilizing proteins and isolating membrane complexes.
Avg. Molecular Wt.:~600
Recommended Storage: RT

CHAPS		
White Crystalline Powder		
packaging	Mfr. No	
1 g AmberGlass	BP571-1	
5 g AmberGlass	BP571-5	
C ₃₂ H ₅₈ N ₂ O ₇ S	P261, P302+P352, P280, P305+P351+P338	
CAS: 75621-03-3		
MW: 614.9		
H335, H319, H315		
Assay (TLC)	>=98%	
Conductivity of a 10% Aqueous Solution (at 24°C)	<50µmhos/cm	
FTIR	Conforms to standard	
Optical Absorbance (1%) at 260nm	<=0.1	
Optical Absorbance (1%) at 280nm	<=0.1	
pH of 10% Solution	5.0-7.0	
Solubility (0.5g/50mℓ)	Clear and colorless liquid	

Applications: CHAPS is a Zwitterionic detergent used for solubilizing membrane proteins.
Recommended Storage: RT
Not on TSCA inventory: for R and D use only; not for manufacturing or commercial purposes.

CHAPSO		
White Powder		
packaging	Mfr. No	
500 mg AmberGlass	BP575-500	
C ₃₂ H ₅₈ N ₂ O ₈ S	H319, H335	
CAS: 82473-24-3	P280, P305+P351+P338	
MW: 630.88		
H302, H312, H332, H315,		
Assay (N ₂)	>=99%	
Conductivity of a 0.5M Aqueous Solution (at 24°C)	<=100µmhos/cm	
Optical Absorbance (1%) at 260nm	<=0.05	
pH of 1% Aqueous Solution	5.0-8.0	
Residue after ignition	<=0.1%	
Solubility (10% in Water)	Clear, haze-free	
TLC	Comparable to reference	
Water (Karl Fischer)	<=7.5%	

Applications: CHAPSO is a Zwitterionic detergent used for protein solubilization.
Recommended Storage: RT
Not on TSCA inventory: for R and D use only; not for manufacturing or commercial purposes.

Dithioerythritol		Electrophoresis
White Crystals or Powder		
packaging	Mfr. No	
5 g PolyBottle	BP173-5	
C ₄ H ₁₀ O ₂ S ₂	P261, P301+P312, P302+P352, P280, P305+P351+P338	
CAS: 6892-68-8		
MW: 154.24		
EINECS: 229-998-8		
H302, H335, H319, H315		
Assay	>=99%	
Corrected Absorbance of Disulfide (S-S) in 0.02M Solution (at 285nm)	<=0.13	
Electrophoresis	To pass test	
IR	Conforms to standard	
Melting Point	83° ±2°C	

Applications: Dithioerythritol is used to reduce disulfide bonds in proteins and peptides and is added in low concentrations as an enzyme stabilizer during protein purification and enzymatic reactions.
Recommended Storage: -20°C

Dithiothreitol		Electrophoresis
White Crystals or Powder		
packaging	Mfr. No	
5 g AmberGlass	BP172-5	
25 g AmberGlass	BP172-25	
C ₄ H ₁₀ O ₂ S ₂	P261, P301+P312, P302+P352, P280, P305+P351+P338	
CAS: 3483-12-3		
MW: 154.24		
H302, H335, H319, H315		
Assay	>=99%	
Corrected Absorbance of Disulfide (S-S) in 0.02M Solution (at 280nm)	<=0.04	
DNase	Not detected	
Electrophoresis	To pass test	
IR	Conforms to standard	
Melting Point	41° ±2°C	
Protease	Not detected	
RNase	Not detected	

Applications: Dithiothreitol is used to reduce disulfide bonds in proteins and peptides and is added in low concentrations as an enzyme stabilizer during protein purification and enzymatic reactions.
Recommended Storage: -20°C

Leupeptin Hemisulfate
White Powder

packaging	Mfr. No
1 mg AmberGlass	BP2662-1
5 mg AmberGlass	BP2662-5
25 mg AmberGlass	BP2662-25
100 mg AmberGlass	BP2662-100
C ₂₀ H ₃₈ N ₆ O ₄ ·0.5H ₂ O ₄ S CAS: 103476-89-7 MW: 475.6 H302, H312, H332	P280, P301+P312, P305+P351+P338, P302+P352



Applications: Leupeptin Hemisulfate is a biochemical used as an inhibitor for trypsin, plasmin, papain, and cathepsin B.
Recommended Storage: -20°C, desiccate, protect from light.
Not on TSCA inventory: for R and D use only; not for manufacturing or commercial purposes.

Thapsigargin
Clear or Whitish Oil

packaging	Mfr. No
1 mg GlassAmpule	BP629-1
C ₃₄ H ₅₀ O ₁₂ CAS: 67526-95-8 MW: 650.76 H311, H331, H301	P301+P310, P280, P312, P302+P350, P304+P340



Applications: This compound is useful as an inhibitor of microsomal calcium-ATPase and as a potent cell-permeable, intracellular calcium releaser that is independent of IP3.
Isolated from Thapsia garganica *ℓ*.
Recommended Storage: -20°C
Not on TSCA inventory: for R and D use only; not for manufacturing or commercial purposes.
UN 2811; DOT Class 6.1:Poison

Pepstatin A
White Powder

packaging	Mfr. No
5 mg AmberGlass	BP2671-5
10 mg AmberGlass	BP2671-10
25 mg AmberGlass	BP2671-25
100 mg AmberGlass	BP2671-100
250 mg AmberGlass	BP2671-250
C ₃₄ H ₆₃ N ₅ O ₉ CAS: 26305-03-3 MW: 685.9	EINECS: 247-600-0
FTIR	Conforms to Standard
Melting Point	230°-240°C

Applications: Pepstatin A is an inhibitor for pepsin, renin, cathepsin D, and other acid proteases.
Recommended Storage: RT
Not on TSCA inventory: for R and D use only; not for manufacturing or commercial purposes.

Triton* X-100

Electrophoresis

packaging	Mfr. No
100 mℓ PolyBottle	BP151-100
500 mℓ PolyBottle,EcoSafPak*	BP151-500
C ₃₄ H ₆₂ O ₁₁ CAS: 9002-93-1 MW: 646.85 H412, H318, H302	P280, P305+P351+P338, P310, P273, P301+P312



Applications: This nonionic detergent is used to solubilize membranes.
Recommended Storage: RT
EcoSafPak* is an environmentally friendly packaging system made of 100% recyclable material by an SFI certified supplier.

Sarkosyl

packaging	Mfr. No
500 g PolyBottle	BP234-500
C ₁₅ H ₂₈ NNaO ₃ CAS: 137-16-6 MW: 293.39	EINECS: 205-281-5
Assay (dry basis)	>=95%
pH of 10% Solution	7.5-8.5

Applications: This anionic detergent is used for the solubilization of membrane proteins.
Recommended Storage: RT

Actinomycin-D
Red Crystals

packaging	Mfr. No
5 mg PolyMicroTube	BP606-5
10 mg PolyMicroTube	BP606-10
C ₆₂ H ₈₆ N ₁₂ O ₁₆ CAS: 50-76-0 MW: 1255.44 EINECS: 200-063-6	H301 P301+P310



Applications: This antibiotic is a widely used transcription inhibitor. Isolated from Streptomyces parvuius.
Recommended Storage: 4°C
Not on TSCA inventory: for R and D use only; not for manufacturing or commercial purposes.

Amikacin From Kanamycin A
White Crystalline Powder

packaging	Mfr. No
250 mg AmberGlass	BP2643-250
1 g AmberGlass	BP2643-1
C ₂₂ H ₄₃ N ₅ O ₁₃ CAS: 37517-28-5 MW: 585.61 Melting point: White	Crystalline Powder EINECS: 253-538-5
Loss on drying	<=9%
pH (1% Aqueous Solution)	9.5-11.5
Potency	900µg/mg

Melting Point: 203°-204°C
Recommended Storage: 0° to 5°C, desiccate.
Not on TSCA inventory: for R and D use only; not for manufacturing or commercial purposes.

Amphotericin B
Yellow Powder

packaging	Mfr. No
250 mg AmberGlass	BP928-250
C ₄₇ H ₇₃ NO ₁₇ CAS: 1397-89-3 MW: 924.09 EINECS: 215-742-2	H315, H319, H335 P280, P305+P351+P338



Applications: Amphotericin B is an antibiotic used for controlling microbial contamination in tissue culture applications.
Recommended Storage: RT
Not on TSCA inventory: for R and D use only; not for manufacturing or commercial purposes.

Amphotericin B 250µg/mℓ
Yellow Solution

packaging	Mfr. No
20 mℓ PolyBottle	BP2645-20
50 mℓ PolyBottle	BP2645-50
C ₄₇ H ₇₃ NO ₁₇ CAS: 1397-89-3 MW: 924.1	EINECS: 215-742-2
Concentration	250 µg/mℓ
pH at 25°C	Inclusive between 6.5-7.5
Sterility	Sterile filtered

Applications: Amphotericin B is effective against yeast and fungi. An antibiotic solution used in cell biology.
Recommended Storage: -20°C
Not on TSCA inventory: for R and D use only; not for manufacturing or commercial purposes.

Ampicillin Sodium Salt
Crystalline Powder

packaging	Mfr. No
5 g PolyBottle	BP1760-5
25 g PolyBottle	BP1760-25
C ₁₆ H ₁₈ N ₃ NaO ₄ S CAS: 69-52-3 MW: 371.39 H317, H334	P261, P304+P341, P342+P311, P280
Dimethylaniline	<0.002%
Moisture	<2.0%
pH of 1% Aqueous Solution	8.0-10.0
Potency	845-988µg/mg



Recommended Storage: 4°C
Not on TSCA inventory: for R and D use only; not for manufacturing or commercial purposes.

Ampicillin Trihydrate
Off-white Powder

packaging	Mfr. No
25 g PolyBottle	BP902-25
C ₁₆ H ₁₉ N ₃ O ₄ S·3H ₂ O CAS: 7177-48-2 MW: 403.45 H317, H334	P261, P304+P341, P342+P311, P280
Assay	>=90%
pH of 1% Solution (at 25°C)	3.5-6.0
Potency (Anhydrous)	Inclusive between 950-1050µg/mg



Recommended Storage: 4°C
Not on TSCA inventory: for R and D use only; not for manufacturing or commercial purposes.

Antibiotic A23187 Free Acid
Off-white Crystalline Powder

packaging	Mfr. No
5 mg PolyMicroTube	BP2502-5
C ₂₉ H ₃₇ N ₃ O ₆ CAS: 52665-69-7 MW: 523.63 EINECS: 258-084-1 H302, H312, H332, H315,	H319, H335 P302+P352, P280, P261, P301+P312, P304+P340, P305+P351+P338
Di-n-butyl ether	>=99.0%
Purity (HPLC)	>=99.0%



Applications: Antibiotic A23187 is a biochemical tool used to study the role of divalent cations in various biological systems. This antibiotic demonstrates weak in vitro activity against gram-positive bacteria and fungi. A widely used ionophore. It has the ability to increase the intracellular concentrations of cations and is widely used to increase the intracellular levels of calcium.
Recommended Storage: -20°C
Not on TSCA inventory: for R and D use only; not for manufacturing or commercial purposes.

Kanamycin Sulfate
White Powder

packaging

Mfr. No

5 gPolyBottle

BP906-5

C₁₈H₃₆N₄O₁₁·H₂SO₄
CAS: 25389-94-0
MW: 582.59
EINECS: 246-933-9

H334, H335, H360D
P261, P304+P341,
P342+P311, P201,
P308+P313

!

Chromatographic Purity

To pass USP test

Ignition residue

<1.0%

Kanamycin B

<4.0%

pH of a 1% solution

Inclusive between 6.5-8.5

Potency

>750µg/mg

Recommended Storage: RT

Not on TSCA inventory: for R and D use only; not for manufacturing or commercial purposes.

Miconazole
Free Base

packaging

Mfr. No

1 gAmberGlass

BP2668-1

5 gAmberGlass

BP2668-5

25 gAmberGlass

BP2668-25

C₁₈H₁₄Cl₄N₂O
CAS: 22916-47-8
MW: 416.13

EINECS: 245-324-5

Chromatography purity

To pass test

Ignition residue

<0.2%

Loss on drying

<0.5%

Purity

Inclusive between 98-102%

Residual solvents

To pass test

Applications: Miconazole is used as an antifungal agent.

Recommended Storage: RT

Not on TSCA inventory: for R and D use only; not for manufacturing or commercial purposes.

Ketoconazole
White Powder

packaging

Mfr. No

50 mgAmberGlass

BP2734-50

100 mgAmberGlass

BP2734-100

500 mgAmberGlass

BP2734-500

1 gAmberGlass

BP2734-1

C₂₆H₂₈Cl₂N₄O₄
CAS: 65277-42-1
MW: 531.44
EINECS: 265-667-4

H360F, H301, H373, H410
P264, P270, P201, P281,
P308+P313, P273

FTIR

Conforms to standard

Melting Point

148°-152°C

Purity

>=99%

Applications: Ketoconazole is a cytochrome P450 inhibitor and demonstrates anti-psoriatic activity. An antifungal agent.

Recommended Storage: RT

UN 2811; DOT Class 6.1:Poison

Metolazone
White Crystalline Powder

packaging

Mfr. No

25 mgAmberGlass

BP2667-25

50 mgAmberGlass

BP2667-50

100 mgAmberGlass

BP2667-100

C₁₆H₁₆ClN₃O₃S
CAS: 17560-51-9
MW: 365.84

EINECS: 241-539-3

FTIR

Conforms to Reference

Melting Point

245°-260°C

Applications: Metolazone is an anti-hypertensive.

Recommended Storage: RT

Not on TSCA inventory: for R and D use only; not for manufacturing or commercial purposes.

Nalidixic Acid

packaging

Mfr. No

25 gPolyBottle

BP908-25

C₁₂H₁₂N₂O₃
CAS: 389-08-2
MW: 232.24
EINECS: 206-864-7

H317, H334, H351
P261, P304+P341,
P342+P311, P280,
P281

Assay

>=98%

Heavy Metals (Pb)

<=0.002%

Loss on Drying (at 105°C)

<=0.5%

Melting Range

225°-231°C

Applications: Nalidixic Acid functions as a DNA polymerase inhibitor by interfering with DNA gyrase activity.

Recommended Storage: RT

Not on TSCA inventory: for R and D use only; not for manufacturing or commercial purposes.

Neomycin Sulfate
White Powder

packaging

Mfr. No

5 gAmberGlass

BP2669-5

25 gAmberGlass

BP2669-25

C₂₃H₄₆N₆O₁₃·3H₂O₄S
CAS: 1405-10-3
MW: 908.87
EINECS: 215-773-1

H334, H317, H361d
P261, P304+P341,
P342+P311, P280,
P281

FTIR

Conforms to Standard

Loss on drying

<=8%

pH (3.3% aqueous solution)

5-7.5

Potency

>=630µg neomycin/mg

Purity

>=98% of Neomycin B

Applications: Neomycin Sulfate is an antibiotic and is used in tissue culture studies. It can act as an inhibitor of protein biosynthesis.

Recommended Storage: RT, desiccate, protect from light.

Nystatin
White or Almost White Powder

packaging

Mfr. No

5 gAmber Glass

BP2949-5

C₄₇H₇₅NO₁₇
CAS: 1400-61-9
MW: 926.1

EINECS: 215-749-0
Fungicidin Mycostatin

Loss on drying

<=5.0%

pH at 25°C

Inclusive between 6.0-8.0

Potency

>4400 U/mg

Source

Streptomyces noursei

Suspendibility

To pass test

Applications: Increases the permeability of the cell membrane of sensitive fungi by binding to sterols. Antimicrobial spectrum: Yeasts and molds.

Recommended storage: 0°C

Paromomycin Sulfate
Dry, White or Off-White Crystalline Powder

packaging

Mfr. No

1 gAmber Glass

BP2953-1

C₂₃H₄₅N₅O₁₄·H₂SO₄
CAS: 1263-89-4
MW: 713.71
EINECS: 215-031-7
Aminosidine Sulfate

H335, H319, H315
P261, P302+P352, P280,
P305+P351+P338

!

Loss on drying

<5%

Potency

>675 U/mg

Source

Streptomyces krestomuceticus

Applications: Used to study bacterial protein synthesis at the level of 16S ribosomal RNA and 30S ribosome assembly. Inhibits initiation and elongation during protein synthesis. Antimicrobial spectrum: Gram-negative and Gram-positive bacteria.

Recommended storage: RT

Penicillin G Sodium Salt
Dry, Fine White Powder

packaging

Mfr. No

5 gAmber Glass

BP2955-5

C₁₆H₁₇N₂NaO₄S
CAS: 69-57-8
MW: 356.38
Benzylpenicillin sodium salt

H317, H334
P261, P304+P341,
P342+P311, P280

Loss on drying

<1.5%

pH

Inclusive between 5.0-7.5

Potency

Inclusive between 1500-1750 U/mg

Source

Penicillum fungi

Applications: Inhibits bacterial cell wall synthesis; Antimicrobial spectrum: Gram-positive bacteria.

Recommended storage: RT

Penicillin/Streptomycin Mixture
Clear Solution

packaging

Mfr. No

50 mℓPlastic Bottle

BP2959-50

CAS: Penicillin (69-57-8); Streptomycin(3810-74-0)
H334, H317, H360FD
P261, P304+P341, P342+P311, P280, P201,
P308+P313

Penicillin concentration

10000 U/mℓ

pH

Inclusive between 5.0-7.5

sodium chloride

0.9%

Sterile-filtered

To pass test

Streptomycin concentration

10 mg/mℓ

Applications: Interferes with synthesis of bacterial cell wall. Inhibits protein synthesis; binds to 30S ribosomal subunit.

Description: Sterile-Filtered; Penicillin: 10.000U/ml; Streptomycin: 10mg/mℓ

Recommended storage: <0°C

Penicillin/Streptomycin/Glutamine Mixture
Clear Solution

packaging

Mfr. No

50 mℓPlastic Bottle

BP2960-50

CAS: Penicillin (69-57-8); Streptomycin(3810-74-0); Glutamine (56-85-9)
H317, H334, H360FD
P261, P304+P341, P342+P311, P280, P201,
P308+P313

Glutamine concentration

200 mM

Penicillin concentration

10000 U/mℓ

pH

Inclusive between 6.0-7.0

sodium chloride

0.9%

Sterile-filtered

To pass test

Streptomycin concentration

10 mg/mℓ

Applications: Interferes with synthesis of bacterial cell wall. Inhibits protein synthesis. Interferes with bacterial cell division.

Description: Sterile-Filtered; Penicillin: 10.000U/ml; Streptomycin: 10mg/ml; Glutamine: 200mM

Recommended storage: <0°C

Penicillin/Streptomycin/Neomycin Mixture
Clear Solution

packaging

Mfr. No

50 mℓPlastic Bottle

BP2961-50

CAS: Penicillin (69-57-8); Streptomycin(3810-74-0); Neomycin (1405-10-3)
H317, H334, H360FD
P261, P304+P341, P342+P311, P280, P201,
P308+P313

Neomycin concentration

10 mg/mℓ

Penicillin concentration

5000 U/mℓ

pH

Inclusive between 6.3-6.9

sodium chloride

0.9%

Sterile-filtered

To pass test

Streptomycin concentration

5 mg/mℓ

Applications: Interferes with synthesis of bacterial cell wall. Inhibits protein synthesis. Interferes with bacterial cell division.

Description: Sterile-Filtered; Penicillin: 5.000U/ml; Streptomycin: 5mg/ml; Neomycin: 10mg/mℓ

Recommended storage: <0°C

Penicillin-G Potassium Salt
White to Off-white Powder

packaging		Mfr. No
100 g	PolyBottle	BP914-100
C ₁₆ H ₁₇ KN ₂ O ₄ S CAS: 113-98-4 MW: 372.47 EINECS: 204-038-0 H334, H317	P261, P304+P341, P342+P311, P280, P302+P352, P308+P313	
Assay	>=85%	
Loss on Drying (at 105°C)	<=1.0%	
pH of 0.6% Solution	5.5-7.5	

Applications: Penicillin inhibits microbial contamination in tissue culture applications.
Not for drug use.
Recommended Storage: RT

Phleomycin, 2% solution in water

Clear Blue Solution

packaging		Mfr. No
100 mg	Plastic Bottle	BP2962-100
C ₅₅ H ₈₅ CuO ₂₁ N ₂₀ S ₂ .HCl CAS: 11006-33-0 MW: 1526.5 Phleomycin complex	H317 P280	
Concentration	20 mg/ml	
Source	Streptomyces verticillus	
Sterile-filtered	To pass test	
UV abs (200 ug/ml at 244 nm)	Inclusive between 0.75-1.05	

Applications: Disrupts the integrity of DNA, and blocks S-phase entry in the cell cycle.
Description: Concentration: 20mg/ml
Recommended storage: <0°C

Polymyxin B Sulfate

Dry, White or Off-White Powder

packaging		Mfr. No
1 g	Amber Glass	BP2954-1
C ₅₅ H ₉₆ N ₁₆ O ₁₃ .2H ₂ SO ₄ CAS: 1405-20-5 MW: 1385.61 EINECS: 215-774-7	H302 P301+P312	
Assay	>6000 U/mg	
Loss on drying	<7%	
pH	Inclusive between 5.0-7.5	
Source	Bacillus polymyxa	

Applications: Used as a immobilized (substrate bound) agent for removal of endotoxins. Mode of Action: Binds to and interferes with the permeability of the cytoplasmic membrane; Antimicrobial spectrum: Gram-negative bacteria.
Recommended storage: 4°C

Puromycin Dihydrochloride

Dry, White or Off-White Powder

packaging		Mfr. No
100 mg	Clear Glass	BP2956-100
C ₂₇ H ₂₆ N ₇ O ₅ .2HCl.xH ₂ O CAS: 58-58-2 MW: 544.44	H302 P301+P312	
Purity (HPLC)	>98%	
Source	Streptomyces alboniger	
Water (Karl Fischer)	<12%	

Applications: For selection of cells transfected with puromycin N-acetyl transferase gene (pac).
Recommended storage: <0°C

Rapamycin

Dry, White or Off-White Powder

packaging		Mfr. No
1 mg	Clear Glass	BP2963-1
C ₆₁ H ₇₆ NO ₁₃ CAS: 53123-88-9 MW: 914.2	(-)-Rapamycin; Sirolimus	
Melting point	183-184	
Purity (HPLC)	>98%	

Applications: Possesses immunosuppressive properties. Blocks a protein that is involved in cell division. A type of serine/threonine kinase inhibitor.
Recommended storage: <0°C

Rifampicin

Bright Orange Powder

packaging		Mfr. No
250 mg	AmberGlass	BP2679-250
1 g	AmberGlass	BP2679-1
5 g	AmberGlass	BP2679-5
25 g	AmberGlass	BP2679-25
C ₄₃ H ₅₈ N ₄ O ₁₂ CAS: 13292-46-1 MW: 822.94 EINECS: 236-312-0	H360Fd, H302 P264, P270, P201, P281, P308+P313	 

Heavy metals	<20 ppm
Loss on Drying (at 105°C)	<1.0%
Purity	Inclusive between 97-102%
Rifampicin quinine	<1.5%
Sulfated ash	<0.1%
UV/VIS (lambda max) methanol	237-239nm

Applications: Rifampicin inhibits DNA-dependent bacterial RNA Polymerase. Mammalian RNA polymerase is not affected.
Recommended Storage: 4°C
Not on TSCA inventory: for R and D use only; not for manufacturing or commercial purposes.

Spectinomycin Pentahydrate Dihydrochloride
Dry, White or Almost White Powder

packaging		Mfr. No
1 g	Amber Glass	BP2957-1
C ₁₄ H ₂₄ N ₇ O ₇ .2HCl.2H ₂ O CAS: 22189-32-8 MW: 495.35	H335, H315, H319 P261, P302+P352, P280, P305+P351+P338	
Spectinomycin Dihydrochloride Hydrate		
Assay	>603 U/mg	
pH at 25°C	Inclusive between 3.8-5.6	
Residue after ignition	<1%	
Source	Streptomyces spectabilis	
Specific rotation	Inclusive between 15-21	
Water (Karl Fischer)	Inclusive between 16-20%	

Applications: For selection of transformed cells containing the spectinomycin resistance gene and inhibition of bacterial protein synthesis at the level of peptidyl tRNA translocation.
Recommended storage: 4°C



Staurosporine
From Streptomyces sp.

packaging		Mfr. No
100 µg	PolyMicroTube	BP2541-100
500 µg	PolyMicroTube	BP2541-500
1 mg	PolyMicroTube	BP2541-1
C ₂₈ H ₂₆ N ₄ O ₃ CAS: 62996-74-1 MW: 466.5	H351 P281	
Assay	>=98%	

Applications: Staurosporine is a potent protein kinase inhibitor useful for studies on protein phosphorylation in the regulation of cellular functions.
Recommended Storage: 4°C, protect from light.
Not on TSCA inventory: for R and D use only; not for manufacturing or commercial purposes.

Streptomycin Sulfate
White to Light Yellow Powder

packaging		Mfr. No
50 g	PolyBottle	BP910-50
C ₄₂ H ₇₈ N ₁₄ O ₂₄ .3H ₂ SO ₄ CAS: 3810-74-0 MW: 1457.4 EINECS: 223-286-0	H317, H302, H334, H360FD P261, P280, P201, P308+P313, P285	 
Loss on Drying (at 105°C)	<=5%	
pH of a 20% solution	Inclusive between 4.5-7.0	
Potency	Inclusive between 650-850µg/mg	

Applications: Streptomycin Sulfate inhibits protein synthesis. It prevents microbial contamination in tissue culture applications and confers antibiotic resistance/sensitivity in molecular biology procedures.
Not for drug use.
Recommended Storage: 4°C

Tetracycline hydrochloride
Yellow Powder

packaging		Mfr. No
100 g	PolyBottle	BP912-100
C ₂₂ H ₂₄ N ₂ O ₈ .HCl CAS: 64-75-5 MW: 480.9	H319, H361d P281, P280, P305+P351+P338	 
4-epianhydrotetracycline	<=2.0%	
Heavy metals	<=0.005%	
Loss on Drying (at 105°C)	<=2.0%	
pH of a 1% solution	Inclusive between 1.8-2.8	
Potency (Assay)	>900µg/mg	
Specific rotation	Inclusive between -255 to -240 °	

Applications: Tetracycline blocks protein synthesis and is used to select tetracycline-resistant recombinant plasmids.
Not for drug use.
Recommended Storage: RT

Vancomycin Hydrochloride

Dry, White to Light Orange Powder

packaging		Mfr. No
1 g	Amber Glass	BP2958-1
C ₆₆ H ₇₅ Cl ₂ N ₆ O ₂₄ .HCl CAS: 1404-93-9 MW: 1485.73	H319, H335, H315, H317 P280, P261, P302+P352, P305+P351+P338	
Vancocine Hydrochloride; Vancocin Hydrochloride		
Potency	>900 U/mg	
Source	Streptomyces orientalis	
Water (Karl Fischer)	<10%	

Applications: Inhibits formation of peptidoglycan polymers of bacterial cell wall.
Recommended storage: <0°C

(-)-Nicotine Tartrate
White Solid

packaging		Mfr. No
250 mg	AmberGlass	BP2533-250
1 g	AmberGlass	BP2533-1
C ₁₀ H ₁₄ N ₂ .2C ₄ H ₆ O ₆ .2H ₂ O CAS: 65-31-6 MW: 498.44		
Melting Point	85°-95°C	

Applications: (-)-Nicotine Tartrate is used in biochemical research applications.
Recommended Storage: RT
Not on TSCA inventory: for R and D use only; not for manufacturing or commercial purposes.

2'-Deoxyguanosine-5'-triphosphate trisodium salt dihydrate

packaging		Mfr. No
10 mg	AmberGlass	BP419-10
C ₁₀ H ₁₃ N ₅ O ₁₃ P ₃ Na ₃ .2H ₂ O CAS: 2564-35-4 MW: 609.03	EINECS: 219-887-2	
Molar Extinction Coefficient (253nm, phosphate buffer, pH 7.0)	>=13.100	

Applications: The trisodium salt of dGTP is suitable for use in many molecular biology applications.
Recommended Storage: -20°C
Not on TSCA inventory: for R and D use only; not for manufacturing or commercial purposes.

Adenosine 5'-Triphosphate Disodium Salt Trihydrate

packaging		Mfr. No
25 g	PolyBottle	BP413-25
C ₁₀ H ₁₄ N ₅ Na ₂ O ₁₃ P ₃ .3H ₂ O CAS: 51963-61-2 MW: 605.19		
Lead	<=0.002%	
Solubility (5% in H ₂ O, 25°C)	Clear and haze-free	

Applications: The disodium salt of ATP is suitable for use in many molecular biology applications.
Recommended Storage: Store below 0°C

Aphidicolin
White to Off-white Solid

packaging	Mfr. No
1 mg PolyMicroTube	BP615-1
$C_{20}H_{34}O_4$ CAS: 38966-21-1 MW: 338.49	
Assay	>=98%
Carbon	70.97% ±1%
Hydrogen	10.12% ±0.2%
Melting Point	226°-235°C
Solubility (10mg in 1mL Methanol)	To pass test

Applications: This compound is a reversible inhibitor of eukaryotic nuclear DNA replication.
Isolated from *Nigrospora orycae*.
Recommended Storage: 4°C
Not on TSCA inventory: for R and D use only; not for manufacturing or commercial purposes.

Calcium Ionophore A23187
Free Acid

packaging	Mfr. No
1 mg PolyMicroTube	BP595-1
5 mg PolyMicroTube	BP595-5
10 mg PolyMicroTube	BP595-10
$C_{29}H_{37}N_3O_6$ EINECS: 258-084-1 CAS: 52665-69-7 MW: 523.63	
Assay (TLC, Ethylacetate:toluene (1:1))	>=98%
Carbon	66.52% ±1%
Hydrogen	-50.0° ±3° 7.12% ±0.2%
Melting Point	186°-192°C
Nitrogen	8.02% ±0.2%
Optical Rotation α^{25}_D	-50.0° ±3°
Solubility (in DMSO, Ethylacetate, Methanol, and Chloroform)	To pass test

Applications: This widely used ionophore increases the permeability of biological membranes for bivalent cations.
Isolated from *Streptomyces chartreusensis*.
Recommended Storage: 4°C
Not on TSCA inventory: for R and D use only; not for manufacturing or commercial purposes.

Citric Acid Anhydrous Crystalline

packaging		Mfr. No.
500 g	PolyBottle	BP339-500
C ₆ H ₈ O ₇	H319	
CAS: 77-92-9	P264, P305+P351+P338,	
MW: 192.13	P337+P313, P280	
EINECS: 201-069-1		
Assay		>=99.5%
Chloride (Cl)		<=0.001%
Insoluble matter		<=0.005%
Iron		<=3ppm
Lead		<=2ppm
Oxalate		To pass test (about 0.05%)
Phosphate (PO ₄)		<=0.001%
Residue after ignition		<=0.02%
Substances Carbonizable by H ₂ SO ₄ (Tartrates, etc.)		To pass test
Sulfur compounds (as SO ₄)		<=0.002%

Applications: Anhydrous citric acid is suitable for use in tissue culture systems requiring additives.
Recommended Storage: RT

Coenzyme A Trilithium Salt Dihydrate

White Crystalline Powder

packaging		Mfr. No.
25 mg	AmberGlass	BP2510-25
50 mg	AmberGlass	BP2510-50
100 mg	AmberGlass	BP2510-100
250 mg	AmberGlass	BP2510-250
1 g	AmberGlass	BP2510-1
$C_{21}H_{33}N_7O_{16}P_3SLi_3 \cdot 2H_2O$ CAS: 18439-24-2 MW: 821.4		
Assay		>=96%
OD250/OD260 (pH 2)		0.80-0.89

Applications: Coenzyme A is used in biochemical research applications.
Recommended Storage: -20°C
Not on TSCA inventory: for R and D use only; not for manufacturing or commercial purposes.

Copper (II) Sulfate Pentahydrate
Crystalline (Blue)

packaging		Mfr. No.
500 g	AmberGlass	BP346-500
CuO ₄ S ₂ .5H ₂ O	P273, P301+P312,	 
CAS: 7758-99-8	P302+P352, P280,	
MW: 249.68	P305+P351+P338	
H315, H410, H302, H319		
Assay		98.0-102.0%
Calcium (Ca)		<=0.005%
Chloride (Cl)		<=0.001%
Insoluble matter		<=0.005%
Iron		<=0.003%
Nickel		<=0.005%
Nitrogen Compounds (as N)		<=0.002%
Potassium		<=0.01%
Sodium		<=0.02%

Applications: Copper Sulfate is suitable for use in tissue culture systems requiring additives.
Recommended Storage: RT
 DOT Class 9:Miscellaneous

Cyanocobalamin
Crystalline (Dark Red)

packaging	Mfr. No.
1 g AmberGlass	BP862-1
$C_{63}H_{88}CoN_{14}O_{14}P$ CAS: 68-19-9 MW: 1355.38	
Assay	>=98%

Applications: Cyanocobalamin is used as a supplement in cell culture and tissue culture media.
Not for human use.
Recommended Storage: -4°C

Cyclopiazonic Acid
White to Off-white Crystals

packaging	Mfr. No.
5 mg PolyMicroTube	BP630-5
$C_{20}H_{20}N_2O_3$ CAS: 18172-33-3 MW: 336.39 H301	P301+P310, P302+P352, P280, P312 
Assay (TLC-Methylene Chloride:Methanol (9:1))	>=98.0%
Carbon	71.41% ±1%
Hydrogen	5.99% ±0.2%
Maximum Absorbance (lambda Max.)	221 and 281nm
Nitrogen	8.33% ±0.2%
Solubility (10mg/ml in Methanol and Methylene Chloride)	To pass test

Applications: Inhibitor of Ca²⁺-ATPases.
Recommended Storage: -20°C
Not on TSCA inventory: for R and D use only; not for manufacturing or commercial purposes.
UN 2811; DOT Class 6.1:Poison

D-(-)-Arabinose
White Powder

packaging	Mfr. No
25 g AmberGlass	BP2504-25
100 g AmberGlass	BP2504-100
1 kg AmberGlass	BP2504-1
$C_5H_{10}O_5$ CAS: 28697-53-2 MW: 150.13	
Assay	>=99%
FTIR	Conforms to standard
Melting Point	153°-164°C
Optical Rotation α_D^{25}	-97° to -105°

Applications: D-(-)-Arabinose is a component used for culture media.
Recommended Storage: RT
Not on TSCA inventory: for R and D use only; not for manufacturing or commercial purposes.

D-Glycogen

Beef Liver, White to Straw-yellow Crystals

packaging	Mfr. No
5 g AmberGlass	BP676-5
(C ₆ H ₁₀ O ₂) _n	
CAS: 9005-79-2	
EINECS: 9005-79-2	
Specific Rotation (c=0.25, H ₂ O)	+175° ±10°
TLC	Single spot

Applications: Glycogen is used as a carrier in precipitating DNA or RNA in place of sonicated DNA or tRNA.
Recommended Storage: 4°C

D-(-)-Ribose
White to Off-white Powder

packaging	Mfr. No
25 g AmberGlass	BP696-25
$C_5H_{10}O_5$	EINECS: 200-059-4
CAS: 50-69-1	
MW: 150.13	
Assay	>=98%
Specific Rotation α^{20}_D	-20° $\pm$ 1°
Melting Point	88° $\pm$ 2°C
TLC	Single spot

Applications: This sugar is suitable for use in tissue culture systems requiring additives.

Recommended Storage: RT

D-(+)-Maltose Monohydrate
White to Off-white Powder

packaging	Mfr. No
500 g AmberGlass	BP684-500
C ₁₂ H ₂₂ O ₁₁ ·H ₂ O	
CAS: 6363-53-7	
MW: 360.29	
Assay	>=97%
Specific Rotation α ²⁰ _D	+130° ±2°
Melting Point	119°-139°C
TLC	Single spot

Applications: Maltose is suitable for use in cell culture systems requiring sugar additives.

Recommended Storage: RT

D-(+)-Trehalose, Dihydrate
White Powder

packaging		Mfr. No
10 g	PolyBottle	BP2687-10
25 g	PolyBag	BP2687-25
100 g	PolyBottle	BP2687-100
1 kg	PolyBottle	BP2687-1
$C_{12}H_{22}O_{11} \cdot 2H_2O$ CAS: 6138-23-4 MW: 378.32		H315 P302+P352 
Assay		>=98%
FTIR		Conforms to standard
Loss on drying		8.5-10.5%

Applications: D-(+)-Trehalose is used in several biochemical techniques.
Recommended Storage: RT
Not on TSCA inventory: for R and D use only; not for manufacturing or commercial purposes.

D-(+)-Xylose
White to Off-white Powder

packaging	Mfr. No
250 g AmberGlass	BP708-250
<chem>C5H10O3</chem>	
CAS: 58-86-6	
MW: 150.13	
Assay	>=98%
Specific Rotation α^{20}_D	+20° $\pm$ 1°
TLC	Single spot

Applications: D-(+)-Xylose is often used as a cell culture additive.
Recommended Storage: RT

Dacarbazine
Pale Yellow Powder

packaging	Mfr. No
100 mg AmberGlass	BP2513-100
250 mg AmberGlass	BP2513-250
1 g AmberGlass	BP2513-1
5 g AmberGlass	BP2513-5

C₆H₁₀N₆O
CAS: 4342-03-4
MW: 182.19

EINECS: 224-396-1

Applications: Dacarbazine is an alkylating agent used in cancer research. Antineoplastic.
Recommended Storage: 4° C
Not on TSCA inventory: for R and D use only; not for manufacturing or commercial purposes.

D-Biotin

packaging	Mfr. No
1 g AmberGlass	BP232-1

C₁₀H₁₆N₂O₃S
CAS: 58-85-5
MW: 244.31

Applications: D-Biotin, which forms a stable complex with Avidin, is useful in the study of Avidin-Biotin interactions and in affinity chromatography. Not for human use.
Recommended Storage: -20°C

Dextrose Anhydrous
Crystalline Granules

Molecular Biology

packaging	Mfr. No
500 g AmberGlass	BP350-500
1 kg AmberGlass	BP350-1

C₆H₁₂O₆
CAS: 50-99-7
MW: 180.16

EINECS: 200-075-1

Applications: Dextrose is often used as a cell culture additive.
Recommended Storage: RT

DL-Thioctic Acid

packaging	Mfr. No
1 g AmberGlass	BP2682-1
5 g AmberGlass	BP2682-5
25 g AmberGlass	BP2682-25

C₈H₁₄O₂S₂
CAS: 1077-28-7
MW: 206.32
EINECS: 214-071-2
H302, H335, H319, H315

P261, P301+P312,
P302+P352, P280,
P305+P351+P338

Applications: DL-Thioctic Acid is a growth factor for many bacteria and protozoa; prosthetic group, coenzyme, or substrate in plants, microorganisms, and animal tissue.
Recommended Storage: RT

D-Mannitol
White to Off-white Powder

packaging	Mfr. No
500 g PolyBottle	BP686-500

C₆H₁₄O₆
CAS: 69-65-8
MW: 182.17

Assay >=98%
Specific Rotation α²⁰_D (c=8, HCl) +24° ±1°
Melting Point 165°-168°C
TLC Single spot

Applications: Mannitol is suitable for use in cell culture systems requiring sugar additives.
Recommended Storage: RT

D-Sucrose

Molecular Biology

packaging	Mfr. No
1 kg PolyBottle	BP220-1
2.5 kg PolyBottle	BP220-212
10 kg PolyPail	BP220-10

C₁₂H₂₂O₁₁
CAS: 57-50-1
MW: 342.29

EINECS: 200-334-9

Assay >=99.9%
Color of a 50% Solution (APHA) <=10
Specific Rotation α²⁰_D +66° to +67°
DNase Not detected
Insoluble matter <=0.005%
Invert Sugar <=0.1%
Iron <=5ppm
Lead <=5ppm
Loss on Drying (at 105°C) <=0.03%
Melting Point 160°-186°C
Protease Not detected
Residue after ignition <=0.01%
RNase Not detected
Titratable Acid <=0.0008mEq/g

Applications: D-Sucrose is suitable for the preparation of density gradients used in purification of proteins and nucleic acids by ultracentrifugation.
Recommended Storage: RT

Folic Acid
Yellow to Orange Powder

packaging	Mfr. No
5 g AmberGlass	BP2519-5
10 g AmberGlass	BP2519-10
25 g AmberGlass	BP2519-25
100 g AmberGlass	BP2519-100

C₁₉H₁₉N₇O₆
CAS: 59-30-3
MW: 441.4

Assay >=95%
FTIR Conforms to standard
Loss on drying <=8.5%
Residue on Ignition (sulfated) <=0.3%
UV/VIS lambda max (0.1N NaOH) 256 ± 3nm

Applications: Folic Acid is used in cell culture.
Recommended Storage: RT, protect from light

Upsilon -L-Glutamyl-p-Nitroanilide, Hydrochloride

packaging	Mfr. No
500 mg AmberGlass	BP2656-500
1 g AmberGlass	BP2656-1
5 g AmberGlass	BP2656-5
10 g AmberGlass	BP2656-10
25 g AmberGlass	BP2656-25

C₁₁H₁₃N₃O₅.ClH
CAS: 67953-08-6
MW: 303.7

Lambda max (ethanol) 317nm ±5nm
Melting Point 210°-222°C
Optical Rotation α²⁵_D +30° ±5°

Applications: Upsilon -L-Glutamyl-p-Nitroanilide, Hydrochloride is a water soluble form of GPNA suitable for use as substrate for Upsilon -glutamyl transpeptidase.
Recommended Storage: 0° to 5°C

Glutathione Reduced
White Crystalline Powder

packaging	Mfr. No
1 g AmberGlass	BP2521-1
5 g AmberGlass	BP2521-5
10 g AmberGlass	BP2521-10
25 g AmberGlass	BP2521-25
50 g AmberGlass	BP2521-50
100 g AmberGlass	BP2521-100

C₁₀H₁₇N₃O₆S
CAS: 70-18-8
MW: 307.32

EINECS: 200-725-4

Arsenic <=1.0%
Assay >=98.%
Heavy metals <=10 ppm
Purity >=98.0
Residue after ignition <=0.01%
Specific Rotation (at 20°C) -16° ±5°
Sulfate <=0.048
Trichloroethylene <=0.020%

Applications: Glutathione is a useful tripeptide involved in many aspects of metabolism, including transport of Upsilon -glutamyl amino acids and reductive cleavage of disulfide bonds.
Recommended Storage: 4°C

Glycyl-L-Glutamic Acid
White Fluffy Powder

packaging	Mfr. No
100 mg AmberGlass	BP2522-100
500 mg AmberGlass	BP2522-500
1 g AmberGlass	BP2522-1

C₇H₁₂N₂O₅
CAS: 7412-78-4
MW: 204.18

EINECS: 231-019-4

Applications: Glycyl-L-Glutamic Acid is used in biochemical research applications.
Recommended Storage: 0°C
Not on TSCA inventory: for R and D use only; not for manufacturing or commercial purposes.

Hematoxylin Trihydrate
Reddish-tan Powder

packaging	Mfr. No
5 g AmberGlass	BP2523-5
25 g AmberGlass	BP2523-25
100 g AmberGlass	BP2523-100

C₁₆H₁₄O₆
CAS: 517-28-2
MW: 302.27

EINECS: 208-237-3

Loss on drying Report
UV/VIS lambda max (MeOH) 290-292nm

Applications: Hematoxylin Trihydrate is used mainly as a stain in microscopy.
Recommended Storage: RT

Isopropyl-β-D-thiogalactopyranoside
Dioxane-free

packaging	Mfr. No
1 g AmberGlass	BP1755-1
10 g AmberGlass	BP1755-10

C₉H₁₈O₅S
CAS: 367-93-1
MW: 238.31

EINECS: 206-703-0

Appearance White powder or crystals
Assay (TLC) >99.0%
Specific Rotation α²⁰_D -31.5°±3.0°
Melting Point 110° to 114°C
NMR Analysis Dioxane-free

Applications: Because IPTG induces β-Galactosidase activity in many bacteria, it is used to maximize the expression of cloned genes under control of the lac operon.
Recommended Storage: 4°C

Isopropyl-β-D-thiogalactopyranoside
Dioxane-free

packaging	Mfr. No
100 g AmberGlass	BP1755-100

CAS: 367-93-1
EINECS: 206-703-0

Appearance White powder or crystals
Assay (TLC) >99.0%
Specific Rotation α²⁰_D -31.5°±3.0°
Melting Point 110° to 114°C
NMR Analysis Dioxane-free

Applications: Because IPTG induces β-Galactosidase activity in many bacteria, it is used to maximize the expression of cloned genes under control of the lac operon.
Recommended Storage: 4°C

Isopropyl-β-D-thiogalactopyranoside
White Powder, Dioxane Crystallized

packaging	Mfr. No
1 g AmberGlass	BP1620-1
10 g AmberGlass	BP1620-10
C ₉ H ₁₈ O ₅ S CAS: 367-93-1 MW: 238.3 EINECS: 206-703-0 H351, H319, H335,	EUH019, EUH066 P280, P305+P351+P338, P302+P350, P233
Optical Absorbance of a 5% Aqueous Solution at 300nm	<=0.13
Optical Absorbance of a 5% Aqueous Solution at 400nm	<=0.06
Solubility (5% in H ₂ O)	To pass test
Specific rotation	-27.0° ±1.0°
TLC	Single spot

Applications: Because IPTG induces β-Galactosidase activity in many bacteria, it is used to maximize the expression of cloned genes under the control of the lac operon.
Recommended Storage: -20°C

Kinetin (6-Furfurylaminopurine)
Off-white to Tan Powder

packaging	Mfr. No
1 g AmberGlass	BP942-1
C ₁₀ H ₉ N ₅ O CAS: 525-79-1 MW: 215.21 EINECS: 208-382-2 H315, H302, H335, H319	P261, P301+P312, P302+P352, P280, P305+P351+P338
Melting Point	265° ±3°C
TLC	Single spot
UV Extinction in 0.1N NaOH (=272nm)	17.000-18.000

Applications: Kinetin is a cytokine often used in plant cell culture research.
Recommended Storage: RT

L-Alanine
White Crystals or Crystalline Powder

packaging	Mfr. No
100 g AmberGlass	BP369-100
C ₃ H ₇ NO ₂ CAS: 56-41-7 MW: 89.09	EINECS: 200-273-8
Ammonium (NH ₄)	<=0.02%
Arsenic	<=1ppm
Assay	98.5-101.0%
Specific Rotation α ²⁰ _D	+14.3° to +15.2°
Heavy Metals (Pb)	<=10ppm
Loss on Drying (at 105°C)	<=0.20%
Other Amino Acids	Chromatographically not detectable
Residue on Ignition (sulfated)	<=0.10%
State of Solution	>=98.0% transmittance

Applications: This amino acid is suitable for use in tissue culture systems requiring additives.
Recommended Storage: RT

L-Arginine, Free Base
White Powder

packaging	Mfr. No
100 g PolyBottle	BP2505-100
500 g PolyBottle	BP2505-500
1 kg PolyBottle	BP2505-1
C ₆ H ₁₄ N ₄ O ₂ CAS: 74-79-3 MW: 174.2 EINECS: 200-811-1	H319 P280, P305+P351+P338
Assay	>=99%
FTIR	Conforms to standard
Moisture	<=0.3%
Optical Rotation (c=1.6, 6N HCl)	+26.3° to +27.7°

Applications: L-Arginine is an essential amino acid used in biochemical research applications.
Recommended Storage: RT

L-Arginine
White Crystals or Crystalline Powder

packaging	Mfr. No
100 g AmberGlass	BP370-100
C ₆ H ₁₄ N ₄ O ₂ CAS: 74-79-3 MW: 174.2 EINECS: 200-811-1	H319 P280, P305+P351+P338
Ammonium (NH ₄)	<=0.02%
Arsenic	<=1ppm
Assay	98.5-101.0%
Specific Rotation α ²⁰ _D	+26.9° to +27.9°
Heavy Metals (Pb)	<=10ppm
Loss on Drying (at 105°C)	<=0.5%
Other Amino Acids	Chromatographically not detectable
Residue on Ignition (sulfated)	<=0.10%
State of Solution	>=98.0% transmittance

Applications: This amino acid is suitable for use in tissue culture systems requiring additives.
Recommended Storage: RT

L-Arginine Hydrochloride
White Crystals or Crystalline Powder

packaging	Mfr. No
100 g AmberGlass	BP372-100
C ₆ H ₁₄ N ₄ O ₂ ·HCl CAS: 1119-34-2 MW: 210.67	EINECS: 214-275-1
Ammonium (NH ₄)	<=0.02%
Arsenic	<=1ppm
Assay	98.5-101%
Chloride (Cl)	16.58-17.00%
Specific Rotation α ²⁰ _D	+22° to +23° (+26.7° to +27.7° calc. as free base)
Heavy Metals (Pb)	<=10ppm
Loss on Drying (at 105°C)	<=0.20%
Other Amino Acids	<1%
Residue on Ignition (sulfated)	<=0.10%
State of Solution	>=98.0% transmittance

Applications: The hydrochloride salt of L-Arginine is suitable for use in tissue culture systems requiring additives and has many of the same applications as L-Arginine.
Recommended Storage: RT

L-Ascorbic Acid
White Crystalline Powder

packaging	Mfr. No
500 g AmberGlass	BP351-500
C ₆ H ₈ O ₆ CAS: 50-81-7 MW: 176.13	EINECS: 200-066-2
Assay	>=99.0%
Heavy Metals (Pb)	<=0.002%
Iron	<=0.001%
Residue after ignition	<=0.1%
Specific rotation	+20.5° to +21.5°

Applications: Ascorbic Acid is suitable for use in tissue culture systems requiring additives.
Not for human use.
Recommended Storage: RT

L-Asparagine Monohydrate
White Crystals or Crystalline Powder

packaging	Mfr. No
100 g AmberGlass	BP373-100
C ₅ H ₈ N ₂ O ₃ ·H ₂ O CAS: 5794-13-8 MW: 150.14	
Ammonium (NH ₄)	<=0.1%
Arsenic	<=1ppm
Assay	98.5-101.0%
Specific Rotation α ²⁰ _D	+33.5° to +36.5°
Heavy Metals (Pb)	<=10ppm
Loss on Drying (at 130°C)	11.5-12.5%
Other Amino Acids	<=1.0%
Residue on Ignition (sulfated)	<=0.10%
State of Solution	>=98.0% transmittance

Applications: This amino acid is suitable for use in tissue culture systems requiring additives.
Recommended Storage: RT

L-Aspartic Acid
White Crystals or Crystalline Powder

packaging	Mfr. No
100 g AmberGlass	BP374-100
C ₄ H ₇ NO ₄ CAS: 56-84-8 MW: 133.1	EINECS: 200-291-6
Ammonium (NH ₄)	<=0.02%
Arsenic	<=1ppm
Assay	98.5-100.5%
Specific Rotation α ²⁰ _D	+24.8° to +25.8°
Heavy Metals (Pb)	<=10ppm
Loss on drying	<=0.20%
Other Amino Acids	Chromatographically not detectable
Residue on Ignition (sulfated)	<=0.10%
State of Solution	>=98.0% transmittance

Applications: This amino acid is suitable for use in tissue culture systems requiring additives.
Recommended Storage: RT

L-Cysteine Hydrochloride Monohydrate
White Crystals

packaging	Mfr. No
100 g AmberGlass	BP376-100
C ₃ H ₇ NO ₂ ·S·HCl·H ₂ O CAS: 7048-04-6 MW: 175.64	H315, H319, H335 P280, P305+P351+P338
Ammonium (NH ₄)	<=0.02%
Arsenic (as As ₂ O ₃)	<=1ppm
Assay	98.5-101%
Chloride (Cl)	19.89-20.29%
Specific Rotation α ²⁰ _D (c=8, HCl)	+5.5° to +7.0° (+8.0° to +10.1° calc. as cysteine)
Heavy Metals (Pb)	<=10ppm
Loss on drying	8.5-12%
Other Amino Acids	Chromatographically not detectable
Residue on Ignition (sulfated)	<=0.10%
State of Solution	>=98.0% transmittance

Applications: This amino acid is suitable for use in tissue culture systems requiring additives.
Recommended Storage: RT

L-Cystine
White Crystals or Crystalline Powder

packaging	Mfr. No
100 g AmberGlass	BP377-100
C ₆ H ₁₂ N ₂ O ₄ S ₂ CAS: 56-89-3 MW: 240.3	P261, P301+P312, P302+P352, P280, P305+P351+P338
EINECS: 200-296-3 H315, H319, H302, H335	
Ammonium (NH ₄)	<=0.02%
Arsenic	<=1ppm
Assay	98.5-101.0%
Specific Rotation α ²⁰ _D	-215° to -225°
Heavy Metals (Pb)	<=10ppm
Loss on drying	<=0.30%
Other Amino Acids	Chromatographically not detectable
Residue on Ignition (sulfated)	<=0.10%
State of Solution	>=98.0% transmittance

Applications: This amino acid is suitable for use in tissue culture systems requiring additives.
Recommended Storage: RT

L-Glutamic Acid
White Crystals or Crystalline Powder

packaging	Mfr. No
100 g AmberGlass	BP378-100
C ₅ H ₉ NO ₄ CAS: 56-86-0 MW: 147.13	EINECS: 200-293-7
Ammonium (NH ₄)	<=0.02%
Arsenic	<=1ppm
Assay	99.0-100.5%
Specific Rotation α ²⁰ _D	+31.5° to +32.5°
Heavy Metals (Pb)	<=10ppm
Loss on Drying (at 105°C)	<=0.1%
Other Amino Acids	Chromatographically not detectable
Residue on Ignition (sulfated)	<=0.10%
State of Solution	>=98.0% transmittance

Applications: This amino acid is suitable for use in tissue culture systems requiring additives.
Recommended Storage: RT

L-Glutamine White Crystals or Crystalline Powder		
packaging	Mfr. No	
100 g AmberGlass	BP379-100	
C ₅ H ₁₀ N ₂ O ₃ CAS: 56-85-9 MW: 146.15	EINECS: 200-292-1	
Ammonium (NH ₄)	<=0.1%	
Arsenic	<=1ppm	
Assay	98.5-101.0%	
Specific Rotation α ²⁰ _D	+6.3° to +7.3°	
Heavy Metals (Pb)	<=10ppm	
Loss on drying	<=0.20%	
Other Amino Acids	<=1.0%	
Residue on Ignition (sulfated)	<=0.10%	
State of Solution	>=95.0% transmittance	

Applications: This amino acid is suitable for use in tissue culture systems requiring additives.
Recommended Storage: RT

L-Histidine White Crystals or Crystalline Powder		
packaging	Mfr. No	
100 g AmberGlass	BP382-100	
C ₆ H ₉ N ₃ O ₂ CAS: 71-00-1 MW: 155.16	EINECS: 200-745-3	
Ammonium (NH ₄)	<=0.02%	
Arsenic	<=1ppm	
Assay	98.5-101.0%	
Specific Rotation α ²⁰ _D	+12.5° ±0.5°	
Heavy Metals (Pb)	<=10ppm	
Loss on drying	<=0.20%	
Other Amino Acids	Chromatographically not detectable	
Residue on Ignition (sulfated)	<=0.10%	
State of Solution	>=98.0% transmittance	

Applications: This amino acid is suitable for use in tissue culture systems requiring additives.
Recommended Storage: RT

L-Histidine Hydrochloride Monohydrate White Crystals or Crystalline Powder		
packaging	Mfr. No	
100 g AmberGlass	BP383-100	
C ₆ H ₉ N ₃ O ₂ .HCl.H ₂ O CAS: 5934-29-2 MW: 209.64		
Ammonium (NH ₄)	<=0.02%	
Arsenic	<=1ppm	
Assay	98.5-101.0%	
Chloride (Cl)	16.66-17.08%	
Specific Rotation α ²⁰ _D	+8.9° to +9.5° (+12.0° to +12.8° calc. as free base)	
Heavy Metals (Pb)	<=10ppm	
Loss on drying	<=0.20%	
Other Amino Acids	Chromatographically not detectable	
Residue on Ignition (sulfated)	<=0.10%	
State of Solution	>=98.0% transmittance	

Applications: L-Histidine Hydrochloride is suitable for use in tissue culture systems or as an amino acid standard.
Recommended Storage: RT

L-Isoleucine White Crystals or Crystalline Powder		
packaging	Mfr. No	
100 g AmberGlass	BP384-100	
C ₆ H ₁₃ NO ₂ CAS: 73-32-5 MW: 131.17		
Ammonium (NH ₄)	<=0.02%	
Arsenic	<=1ppm	
Assay	98.5-101.0%	
Specific Rotation α ²⁰ _D	+39.5° to +41.5°	
Heavy Metals (Pb)	<=10ppm	
Loss on drying	<=0.20%	
Other Amino Acids	Chromatographically not detectable	
Residue on Ignition (sulfated)	<=0.10%	
State of Solution	>=98.0% transmittance	

Applications: This amino acid is suitable for use in tissue culture systems requiring additives.
Recommended Storage: RT

L-Leucine White Crystals or Crystalline Powder		
packaging	Mfr. No	
100 g AmberGlass	BP385-100	
C ₆ H ₁₃ NO ₂ CAS: 61-90-5 MW: 131.17		
Ammonium (NH ₄)	<=0.02%	
Arsenic	<=1ppm	
Assay	98.5-101.0%	
Specific Rotation α ²⁰ _D	+14.9° to +16.1°	
Heavy Metals (Pb)	<=10ppm	
Loss on Drying (at 105°C)	<=0.2%	
Other Amino Acids	Chromatographically not detectable	
Residue on Ignition (sulfated)	<=0.10%	
State of Solution	>=98.0% transmittance	

Applications: This amino acid is suitable for use in tissue culture systems requiring additives.
Recommended Storage: RT

L-Lysine Hydrochloride White Powder		
packaging	Mfr. No	
100 g AmberGlass	BP386-100	
C ₆ H ₁₄ N ₂ O ₂ .HCl CAS: 657-27-2 MW: 182.65	EINECS: 211-519-9	
Ammonium (NH ₄)	<=0.02%	
Arsenic	<=1ppm	
Assay	98.5-100.5%	
Chloride (Cl)	19.12-19.51%	
Specific Rotation α ²⁰ _D	+20.7° to +21.5° (+25.9° to 26.9° calc. as free base)	
Heavy Metals (Pb)	<=10ppm	
Loss on Drying (at 105°C)	<=0.40%	
Other Amino Acids	Chromatographically not detectable	
Residue on Ignition (sulfated)	<=0.10%	
State of Solution	>=98.0% transmittance	

Applications: This amino acid is suitable for use in tissue culture systems requiring additives.
Recommended Storage: RT

L-Methionine White Crystals or Crystalline Powder		
packaging	Mfr. No	
100 g AmberGlass	BP388-100	
C ₅ H ₁₁ NO ₂ S CAS: 63-68-3 MW: 149.21	EINECS: 200-562-9	
Ammonium (NH ₄)	<=0.02%	
Arsenic	<=1 ppm	
Assay	98.5-100.5%	
Specific Rotation α ²⁰ _D	+23.0° to +24.5°	
Heavy Metals (Pb)	<=10 ppm	
Loss on Drying (at 105°C)	<=0.20%	
Other Amino Acids	Chromatographically not detectable	
Residue on Ignition (sulfated)	<=0.10%	
State of Solution	>=98.0% transmittance	

Applications: This amino acid is suitable for use in tissue culture systems requiring additives.
Recommended Storage: RT

L-Ornithine Hydrochloride White Crystals or Crystalline Powder		
packaging	Mfr. No	
100 g AmberGlass	BP389-100	
C ₅ H ₁₂ N ₂ O ₂ .HCl CAS: 3184-13-2 MW: 168.62	EINECS: 221-678-6	
Ammonium (NH ₄)	<=0.02%	
Arsenic	<=1ppm	
Assay	98.5-101.0%	
Specific Rotation α ²⁰ _D	+23.0° to +25.0° (+29.3° to +31.9° calc. as free base)	
Heavy Metals (Pb)	<=10ppm	
Loss on Drying (at 105°C)	<=0.20%	
Other Amino Acids	Chromatographically not detectable	
Residue on Ignition (sulfated)	<=0.10%	
State of Solution	>=98.0% transmittance	

Applications: This amino acid is suitable for use in tissue culture systems requiring additives.
Recommended Storage: RT

L-Phenylalanine White Crystals or Crystalline Powder		
packaging	Mfr. No	
100 g AmberGlass	BP391-100	
C ₉ H ₁₁ NO ₂ CAS: 63-91-2 MW: 165.19	EINECS: 200-568-1	
Ammonium (NH ₄)	<=0.02%	
Arsenic	<=1ppm	
Assay	98.5-100.5%	
Specific Rotation α ²⁰ _D	-33.5° to -35.0°	
Heavy Metals (Pb)	<=10ppm	
Loss on Drying (at 105°C)	<=0.20%	
Other Amino Acids	Chromatographically not detectable	
Residue on Ignition (sulfated)	<=0.10%	
State of Solution	>=98.0% transmittance	

Applications: This amino acid is suitable for use in tissue culture systems requiring additives.
Recommended Storage: RT

L-Proline White Crystals or Crystalline Powder		
packaging	Mfr. No	
100 g AmberGlass	BP392-100	
C ₅ H ₉ NO ₂ CAS: 147-85-3 MW: 115.13	EINECS: 205-702-2	
Ammonium (NH ₄)	<=0.02%	
Arsenic	<=1ppm	
Assay	98.5-101.0%	
Specific Rotation α ²⁰ _D	-84.5° to -86.0°	
Heavy Metals (Pb)	<=10ppm	
Loss on Drying (at 105°C)	<=0.30%	
Other Amino Acids	Chromatographically not detectable	
Residue on Ignition (sulfated)	<=0.10%	
State of Solution	>=98.0% transmittance	

Applications: This amino acid is suitable for use in tissue culture systems requiring additives.
Recommended Storage: RT

L-Serine White Crystals or Crystalline Powder		
packaging	Mfr. No	
100 g AmberGlass	BP393-100	
C ₃ H ₇ NO ₃ CAS: 56-45-1 MW: 105.09	EINECS: 200-274-3	
Ammonium (NH ₄)	<=0.02%	
Arsenic	<=1ppm	
Assay	98.5-101.0%	
Specific Rotation α ²⁰ _D	+14.4° to +15.5°	
Heavy Metals (Pb)	<=10ppm	
Loss on Drying (at 105°C)	<=0.20%	
Other Amino Acids	Chromatographically not detectable	
Residue on Ignition (sulfated)	<=0.10%	
State of Solution	>=98.0% transmittance	

Applications: This amino acid is suitable for use in tissue culture systems requiring additives.
Recommended Storage: RT

L-Threonine White Crystals or Crystalline Powder		
packaging	Mfr. No	
100 g AmberGlass	BP394-100	
C ₄ H ₉ NO ₃ CAS: 72-19-5 MW: 119.12		
Ammonium (NH ₄)	<=0.02%	
Arsenic	<=1ppm	
Assay	98.5-100.5%	
Specific Rotation α ²⁰ _D	-27.6° to -29.0°	
Heavy Metals (Pb)	<=10ppm	
Loss on Drying (at 105°C)	<=0.20%	
Other Amino Acids	Chromatographically not detectable	
Residue on Ignition (sulfated)	<=0.10%	
State of Solution	>=98.0% transmittance	

Applications: This amino acid is suitable for use in tissue culture systems requiring additives.
Recommended Storage: RT

L-Tryptophan
White Crystals or Crystalline Powder

packaging	Mfr. No
100 g AmberGlass	BP395-100
C ₁₁ H ₁₂ N ₂ O ₂ CAS: 73-22-3 MW: 204.23	
Ammonium (NH ₄)	<=0.02%
Arsenic	<=1ppm
Assay	98.5-100.5%
Specific Rotation α _D ²⁰	-30° to -33°
Heavy Metals (Pb)	<=10ppm
Loss on Drying (at 105°C)	<=0.20%
Other Amino Acids	Chromatographically not detectable
Residue on Ignition (sulfated)	<=0.10%
State of Solution	>=95.0% transmittance

Applications: This amino acid is suitable for use in tissue culture systems requiring additives.
Recommended Storage: RT

L-Tyrosine
White Crystals or Crystalline Powder

packaging	Mfr. No
100 g AmberGlass	BP396-100
C ₉ H ₁₁ NO ₃ CAS: 60-18-4 MW: 181.19 H335, H319, H315	
Ammonium (NH ₄)	<=0.02%
Arsenic	<=1ppm
Assay	98.5-100.5%
Specific Rotation α _D ²⁰	-11.3° to -12.1°
Heavy Metals (Pb)	<=10ppm
Loss on drying	<=0.20%
Other Amino Acids	Chromatographically not detectable
Residue on Ignition (sulfated)	<=0.10%
State of Solution	>=95.0% transmittance

Applications: This amino acid is suitable for use in tissue culture systems requiring additives.
Recommended Storage: RT

L-Valine
White Crystals or Crystalline Powder

packaging	Mfr. No
100 g AmberGlass	BP397-100
C ₆ H ₁₁ NO ₂ CAS: 72-18-4 MW: 117.15	
Ammonium (NH ₄)	<=0.02%
Arsenic	<=1ppm
Assay	98.5-100.5%
Specific Rotation α _D ²⁰	+27.6° to +29.0°
Heavy Metals (Pb)	<=10ppm
Loss on Drying (at 105°C)	<=0.20%
Other Amino Acids	Chromatographically not detectable
Residue on Ignition (sulfated)	<=0.10%
State of Solution	>=98.0% transmittance

Applications: This amino acid is suitable for use in tissue culture systems requiring additives.
Recommended Storage: RT

Lymphocyte Separation Medium

packaging	Mfr. No
500 mℓ PolyBottles	BP2663-500

Applications: LSM is formulated for isolation of mononuclear cells from defribinated or heparinized whole human blood. One-step centrifugation permits separation of mononuclear lymphocytes from erythrocytes, polynuclear lymphocytes and most platelets. LSM is a sterile-filtered solution of a sucrose polymer and diatrizoate salts at a specific gravity of 1.077-1.080 g/mℓ at 20°C.
Recommended Storage: 0° to 5°C

Methotrexate
Yellow Powder

packaging	Mfr. No
10 mg AmberGlass	BP2665-10
100 mg AmberGlass	BP2665-100
C ₂₀ H ₂₂ N ₈ O ₅ ·2H ₂ O CAS: 59-05-2 MW: 490.47 H315, H319, H335 P281, P305+P351+P338, P308+P313	
Chromatographic Purity	Pass test
Moisture	<12%
Organic Volatile Impurities	Pass test
Purity	>98%
Residue after ignition	<0.1%
Specific rotation	+19° to +24°

Applications: Methotrexate is a folic acid antagonist.
Recommended Storage: 0°C
UN 2811; DOT Class 6.1:Poison

Methyl α-D-Mannopyranoside
White Crystalline Powder

packaging	Mfr. No
10 g AmberGlass	BP2530-10
25 g AmberGlass	BP2530-25
100 g AmberGlass	BP2530-100
500 g AmberGlass	BP2530-500
C ₇ H ₁₄ O ₆ CAS: 617-04-9 MW: 194.18	
Assay	>=99%
FTIR	Conforms to standard
Loss on drying	<0.5%
Melting Point	187-197°C
Optical Rotation α _D ²⁵	+76° to +85°

Applications: Methyl α-D-Mannopyranoside is used in biochemical research applications.
Recommended Storage: RT
Not on TSCA inventory: for R and D use only; not for manufacturing or commercial purposes.

Nitro Blue Tetrazolium

packaging	Mfr. No
1 g AmberGlass	BP108-1
C ₄₀ H ₃₀ Cl ₂ N ₁₀ O ₆ CAS: 298-83-9 MW: 817.65 EINECS: 206-067-4	
E1% _{1cm}	>=650 g-1cm-1 min
Lambda Max.	260nm±2nm
Solubility	To pass test

Applications: Used in conjunction with BCIP as a substrate for acid and alkaline phosphatases.
Recommended storage: 4 C

Nitrilotriacetic Acid
White Powder

packaging	Mfr. No
100 g PolyBottle	BP2670-100
500 g PolyBottle	BP2670-500
1 kg PolyBottle	BP2670-1
C ₆ H ₉ NO ₆ CAS: 139-13-9 MW: 191.14 EINECS: 205-355-7 H319, H302, H315, H351	
Purity	>=98%

Applications: Nitrilotriacetic Acid is a chelating and sequestering agent. A builder in synthetic detergents.
Recommended Storage: RT

o-Phenylenediamine
Tan to Brown Solid

packaging	Mfr. No
50 g AmberGlass	BP2537-50
250 g AmberGlass	BP2537-250
1 kg AmberGlass	BP2537-1
C ₆ H ₈ N ₂ CAS: 95-54-5 MW: 108.14	
EINECS: 202-430-6	
FTIR	Conforms to standard
Melting Point	Report

Recommended Storage: RT
UN 1673; DOT Class 6.1:Poison
Free Base

Phosphate Buffered Saline
10X Solution

packaging	Mfr. No
500 mℓ PolyBottle	BP399-500
1 ℓ PolyBottle	BP399-1
4 ℓ PolyPac*	BP399-4
20 ℓ PolyPac*	BP399-20

Chloride Concentration of a 1X Solution 0.140 ±0.004 moles/ℓ
Conductivity of a 1X Solution 14.000 to 17.800μmhos
DNase Not detected
pH of 10X solution (at 25°C) 6.7-6.9
pH of 1X solution (at 25°C) 7.3-7.5
Protease Not detected
RNase Not detected

Applications: A concentrated, standard phosphate buffer solution.
Components: 1.37M Sodium Chloride, 0.027M Potassium Chloride, and 0.119M Phosphate Buffer [7647-14-5 (Sodium Chloride)] ; [7447-40-7 (Potassium Chloride)] ; [7558-79-4 (Sodium Phosphate Dibasic)] ; [7778-77-0 (Potassium Phosphate Monobasic)]
Recommended Storage: RT

Phosphate Buffered Saline
1X Solution

packaging	Mfr. No
4 ℓ PolyPac*	BP2438-4
20 ℓ PolyPac*	BP2438-20

Chloride Concentration of a 1X Solution 0.140 ±0.004 moles/ℓ
Conductivity of a 1X Solution 14.000 to 17.800μmhos/cm
DNase Not detected
pH of 1X solution (at 25°C) 7.3-7.5
Protease Not detected
RNase Not detected

Applications: A standard phosphate buffer used in many biomolecular procedures.
Components:0.137M Sodium Chloride, 0.0027M Potassium Chloride, and 0.0119M Phosphates.
[7647-14-5 (Sodium Chloride)] ; [7447-40-7 (Potassium Chloride)] ; [7558-79-4 (Sodium Phosphate Dibasic)] ; [7778-77-0 (Potassium Phosphate Monobasic)]
Recommended Storage: RT
Filtered through a 0.2-micron filter.

Phorbol 12-Myristate 13-Acetate
Clear or Whitish Film

packaging	Mfr. No
1 mg GlassAmpule	BP685-1
5 mg GlassAmpule	BP685-5
10 mg GlassAmpule	BP685-10
C ₃₆ H ₅₆ O ₈ CAS: 16561-29-8 MW: 616.84 H315, H312 P302+P352, P280	
Assay (HPLC)	>=99.0%
NMR	Identical to standard reference
Solubility (in DMSO, Methanol, Ethanol, Acetone, Ether, and DMF)	To pass test

Applications: PMA is an activator of Protein Kinase C; it is a skin irritant and mouse skin tumor promoter. Film adheres to inside of vial.
Recommended Storage: -20°C
Not on TSCA inventory: for R and D use only; not for manufacturing or commercial purposes.

Pleiotrophin
Human, Recombinant

packaging	Mfr. No
50 µg PolyMicroTube	BP2536-50
Assay	>=97%
description	White Lyophilized Powder

Applications: Pleiotrophin is a biochemical that enhances neurite outgrowth of cerebral cortical neurons.
Sf 21 Expressed.
Recommended Storage: -20°C
Not on TSCA inventory: for R and D use only; not for manufacturing or commercial purposes.

Polyethylene Glycol 8000

packaging	Mfr. No
100 g AmberGlass	BP233-100
1 kg PolyBottle	BP233-1
H(OCH ₂ CH ₂)nOH CAS: 25322-68-3 MW: 7000-8000	
description	White powder
Identification	Pass test

Applications: PEG 8000 is suitable as a medium for the fusion of mammalian cells. It is also used to precipitate bacteriophage from lysed cell supernatants.
Recommended Storage: RT

Pyridoxine Hydrochloride
White Powder

packaging	Mfr. No
10 g AmberGlass	BP2677-10
50 g AmberGlass	BP2677-50
100 g AmberGlass	BP2677-100
C ₆ H ₁₁ NO ₃ .HCl	P261, P302+P352, P280, P305+P351+P338
CAS: 58-56-0	
MW: 205.64	
H335, H315, H319	
Chloride Content	16.9-17.6 on dried basis
FTIR	Conforms to standard
Heavy Metals (Pb)	<=0.003%
Loss on drying	<=0.5%
Melting Point	202°-214°C
pH (10% aqueous)	2.0-4.0
Purity	>=98%
Residue on Ignition (sulfated)	<=0.1%

Applications: Pyridoxine HCl is used as a vitamin (enzyme cofactor).
Recommended Storage: RT, protect from light.

Riboflavin
Fine, Orange-yellow Crystals

Electrophoresis

packaging	Mfr. No	
50 g AmberGlass	BP167-50	
C ₁₇ H ₂₆ N ₄ O ₆ CAS: 83-88-5 MW: 376.36		
EINECS: 201-507-1		
Assay (by UV-VIS)		>=95%
Electrophoresis		To pass test

Applications: Riboflavin is a photoinitiator of acrylamide/bis-acrylamide polymerization.
Not for human use.
Recommended Storage: RT

Sodium Oxalate
White Crystals or Powder

packaging	Mfr. No
500 g PolyBottle	BP353-500
C ₂ Na ₂ O ₄ CAS: 62-76-0 MW: 134	
Ammonium (NH ₄)	<=0.002%
Assay	>=99.0%
Chloride (Cl)	<=0.002%
Heavy Metals (Pb)	<=0.002%
Insoluble matter	<=0.005%
Iron	<=0.001%
Loss on Drying (at 105°C)	<=0.05%
Neutrality	To pass test
Potassium	<=0.005%
Substances Darkened by Hot H ₂ SO ₄	To pass test
Sulfate (SO ₄)	<=0.002%

Applications: This sodium salt is used in buffers for molecular biology and cell culture applications.
Recommended Storage: RT
UN 2811; DOT Class 6.1:Poison

Sodium Pyruvate
White Powder

packaging	Mfr. No	
100 g AmberGlass	BP356-100	
C ₃ H ₃ NaO ₃	H319, H315	
CAS: 113-24-6	P302+P352, P280,	
MW: 110.04	P305+P351+P338	
EINECS: 204-024-4		
Assay (by titration)		>=99%
Free Acid (as Pyruvic Acid)		<=1.0%
Solubility (6g/100mL H ₂ O)	Clear, haze-free	

Applications: This sodium salt is used in buffers for molecular biology and cell culture applications.
Recommended Storage: RT

Sodium Sulfate Anhydrous
White Granules

packaging	Mfr. No
500 g AmberGlass	BP354-500
Na ₂ O ₄ S CAS: 7757-82-6 MW: 142.04	

Assay	>=99.0%
Calcium (Ca)	<=0.01%
Chloride (Cl)	<=0.001%
Heavy Metals (Pb)	<=5ppm
Insoluble matter	<=0.01%
Iron	<=0.001%
Loss on ignition	<=0.5%
Magnesium (Mg)	<=0.005%
Nitrogen Compounds (as N)	<=5ppm
pH of 5% Solution (at 25°C)	5.2-9.2
Phosphate (PO ₄)	<=0.001%
Potassium	<=0.01%

Applications: This sodium salt is used in buffers for molecular biology and cell culture applications.
Recommended Storage: RT

Sodium Sulfite Anhydrous
White Crystals or Crystalline Powder

packaging	Mfr. No
500 g AmberGlass	BP355-500
Na ₂ O ₃ S CAS: 7757-83-7 MW: 126.04	
Assay	>=98.0%
Chloride (Cl)	<=0.02%
Free Acid	To pass test
Heavy Metals (Pb)	<=0.001%
Insoluble matter	<=0.005%
Iron	<=0.001%
Phosphate (PO ₄)	<=5ppm
Titrateable Free Base	<=0.03mEq/g

Applications: This sodium salt is used in buffers for molecular biology and cell culture applications.
Recommended Storage: RT

Sodium Tartrate Dihydrate
Colorless Crystals

packaging	Mfr. No
500 g AmberGlass	BP352-500
C ₄ H ₄ Na ₂ O ₆ .2H ₂ O CAS: 6106-24-7 MW: 230.08	
Ammonium (NH ₄)	<=0.003%
Assay	99.0-101.0%
Calcium (Ca)	<=0.01%
Chloride (Cl)	<=0.005%
Heavy Metals (Pb)	<=0.001%
Insoluble matter	<=0.005%
Iron	<=0.001%
Loss on Drying (at 150°C)	15.61-15.71%
pH of 5% Solution (at 25°C)	7.0 to 9.0
Phosphate (PO ₄)	<=5ppm
Sulfate (SO ₄)	<=0.005%

Applications: This sodium salt is used in buffers for molecular biology and cell culture applications.
Recommended Storage: RT

Succinic Acid
White Crystals

packaging	Mfr. No	
500 g AmberGlass	BP336-500	
C ₄ H ₆ O ₄	H335, H319, H315	
CAS: 110-15-6	P261, P302+P352, P280,	
MW: 118.09	P305+P351+P338	
EINECS: 203-740-4		
Arsenic		<=0.0003%
Ash		<=0.025%
Assay	>=99.5%	
Heavy Metals (Pb)	<=0.001%	
Melting Range	185°-190°C	

Applications: Succinic Acid is used in the preparation of cell culture and tissue culture media.
Recommended Storage: RT

Thiamine Hydrochloride
White Powder

packaging	Mfr. No	
100 g AmberGlass	BP892-100	
C ₁₂ H ₁₇ ClN ₄ OS.HCl	H335, H319, H315	
CAS: 67-03-8	P261, P302+P352, P280,	
MW: 337.26	P305+P351+P338	
EINECS: 200-641-8		
Assay		>=98%
Nitrate		None detected
Residue after evaporation	<=0.2%	
Water	<=6.0%	

Applications: This vitamin is suitable for use in systems requiring additives.
Not for human use.
Recommended Storage: RT

Thimerosal
White Powder

packaging	Mfr. No
10 g AmberGlass	BP2542-10
25 g AmberGlass	BP2542-25
100 g AmberGlass	BP2542-100
500 g AmberGlass	BP2542-500
C ₉ H ₉ HgNaO ₂ S	H310, H330, H373, H410, H300
CAS: 54-64-8	P260, P280, P301+P310,
MW: 404.8	P302+P350, P304+P340,
EINECS: 200-210-4	P273
Assay (anhydrous)	>=97%
Loss on drying	<=0.5%
Melting Point	225°-236°C
pH of 1% Aqueous Solution	6.0-7.4

Recommended Storage: RT, protect from light, desiccate.

Transferrin From Guinea Pig
Golden-Tan Frozen Liquid

packaging	Mfr. No
1 mg AmberGlass	BP2544-1
5 mg AmberGlass	BP2544-5
10 mg AmberGlass	BP2544-10
Assay	>=98%
Sterility	To pass test

Applications: Transferrin is an iron binding protein that displays bacteriostatic and fungistatic characteristics.
Purified by affinity chromatography; sterile filtered; concentrated.
Recommended Storage: -20°C
Not on TSCA inventory: for R and D use only; not for manufacturing or commercial purposes.

Vimentin
From Bovine Lens

packaging	Mfr. No
100 µg PolyMicroTube	BP2545-100
250 µg PolyMicroTube	BP2545-250
CAS: 156289-80-4	
MW: 57	

Recommended Storage: -20°C
Not on TSCA inventory: for R and D use only; not for manufacturing or commercial purposes.

5-Bromo-4-chloro-3-indolyl- β-D-galactopyranoside

packaging	Mfr. No
100 mg AmberGlass	BP1615-100
1 g AmberGlass	BP1615-1
C ₁₄ H ₁₅ BrClNO ₆	
CAS: 7240-90-6	
MW: 408.63	

IR	Conforms to standard
Molar Extinction (292±2nm)	>=4000ℓ g-1cm-1
Solubility (5% in DMF)	Clear and colorless
Specific Rotation (c=0.1, 50% DMF)	-60° ±4°
TLC	Single spot

Applications: X-Gal is a chromogenic substrate for the enzyme β-Galactosidase used in the detection of recombinant bacteriophage. It is also used for immunoblotting and immunocytochemical assays.
Recommended Storage: -20°C
Not on TSCA inventory: for R and D use only; not for manufacturing or commercial purposes.

Toluidine Blue O

packaging	Mfr. No
10 g AmberGlass	BP107-10
C ₁₅ H ₁₆ ClN ₃ S	H302
CAS: 92-31-9	P301+P312
MW: 305.82	
E1% _{1cm}	>=700ℓ g-1cm-1
Lambda Absorption Max.	620-638nm
Loss on Drying (at 105°C)	<=10%
Solubility	To pass test

Applications: Nuclear stain that may be used as an alternative to Methylene Blue
Recommended storage: RT
C.I. 52040

Ethylene Glycol

packaging	Mfr. No
1 ℓ AmberGlass,EcoSafPak*	BP230-1
4 ℓ AmberGlass,EcoSafPak*	BP230-4
C ₂ H ₆ O ₂	H302
CAS: 107-21-1	P280, P301+P310
MW: 62.06	
EINECS: 203-473-3	
Assay	>=99%
Boiling Range	196°-199°C
Color (APHA)	<=10
Iron	<=5ppm
Titratable Acid	0.002mEq/g
Water	<=0.5%

Applications: Due to its property as an antifreeze, Ethylene Glycol may be used for the storage of enzymes at low temperatures.
Recommended Storage: RT
EcoSafPak* is an environmentally friendly packaging system made of 100% recyclable material by an SFI certified supplier.

Acridine Orange Hemizinc Salt

packaging	Mfr. No
10 g AmberGlass	BP116-10
C ₁₇ H ₁₆ N ₃ .HCl.0.5ZnCl ₂	H341, H302, H312, H332
CAS: 10127-02-3	P280
MW: 369.94	
EINECS: 233-353-6	
E1% _{1cm}	>=1100ℓ g-1cm-1
Lambda Max.	490nm ±2nm
Loss on Drying (at 105°C)	<=5%
Solubility	To pass test

Recommended storage: RT
C.I. 46005

Eosin Y
Disodium Salt

packaging	Mfr. No
25 g AmberGlass	BP2419-25
100 g AmberGlass	BP2419-100
C ₂₀ H ₆ Br ₄ Na ₂ O ₅	H319
CAS: 17372-87-1	P280, P305+P351+P338
MW: 691.85	
EINECS: 241-409-6	
Assay	>=88%
Certification as Biological Stain	To pass test

Applications: Eosin Y is used in microbiological differential media or as a biological stain.
Recommended Storage: RT
C.I. 45380

Fast Green FCF

packaging	Mfr. No
10 g AmberGlass	BP123-10
C ₃₇ H ₃₄ N ₂ Na ₂ O ₁₀ S ₃	H341
CAS: 2353-45-9	P201, P202, P281,
MW: 808.84	P308+P313
EINECS: 219-091-5	
Assay	>=90.0%
Lambda Max.	622nm ±2nm

Recommended storage: RT
C.I. 42053

Giemsa Stain
Crystalline

packaging	Mfr. No
5 g AmberGlass	BP2422-5
10 g AmberGlass	BP2422-10
C ₁₄ H ₁₄ ClN ₃ S	EINECS: 257-438-2
CAS: 51811-82-6	
MW: 291.80	
Certification as Biological Stain	To pass test

Applications: Giemsa is a biological stain used for thin blood films to differentiate leukocytes, for staining malarial parasites in thin and thick blood films, and for staining bone marrow to show cell morphology.
Recommended Storage: RT

Glutaraldehyde Solution
25%

packaging	Mfr. No
1 ℓ AmberGlass,EcoSafPak*	BP2548-1
4 ℓ AmberGlass,EcoSafPak*	BP2548-4
C ₅ H ₈ O ₂	H334, H400
CAS: 111-30-8	P301+P330+P331,
MW: 100.12	P305+P351+P338, P310,
EINECS: 203-856-5	P280, P260,
H302, H331, H317, H314,	P273

Assay	24.0-26.5%
Color (APHA)	<=20

Recommended Storage: RT
UN 2810; DOT Class 6.1:Poison
EcoSafPak* is an environmentally friendly packaging system made of 100% recyclable material by an SFI certified supplier.

Glutaraldehyde Solution
50%

packaging	Mfr. No
1 ℓ AmberGlass,EcoSafPak*	BP2547-1
C ₅ H ₈ O ₂	H301, H331, H400
CAS: 111-30-8	P280, P305+P351+P338,
MW: 100.12	P270, P260, P273,
EINECS: 203-856-5	P308+P313
H317, H335, H314, H334,	

Assay	49.0-52.0%
Clarity of Solution	To pass test
Color (APHA)	<=125
pH (at 25°C)	3.7-4.5
Specific Gravity (at 20°C)	1.125-1.135

Recommended Storage: 4°C
UN 2927; DOT Class 6.1:Poison
EcoSafPak* is an environmentally friendly packaging system made of 100% recyclable material by an SFI certified supplier.

Hematoxylin
Crystalline

packaging	Mfr. No
25 g AmberGlass	BP2424-25
100 g AmberGlass	BP2424-100
C ₁₆ H ₁₄ O ₆	EINECS: 208-237-3
CAS: 517-28-2	
MW: 302.27	
Certification as Biological Stain	To pass test

Applications: Hematoxylin is used as a stain in microscopy.
Recommended Storage: RT
C.I. 75290

Methylene Blue

packaging	Mfr. No
100 g AmberGlass	BP117-100
C ₁₆ H ₁₈ ClN ₃ S	H302
CAS: 61-73-4	P301+P312
MW: 319.85	
E1% _{1cm}	>=2000g-1cm-1
Lambda Max.	662nm±2nm
Loss on Drying (at 105°C)	8-18%
Solubility	To pass test

Recommended storage: RT
C.I. 52015

2M Calcium Chloride Solution
Microbial Cell Culture

packaging	Mfr. No
5x10mℓ Poly Bottle	BP9742-10X5
CaCl ₂	MW: 110.99
CAS: 10035-04-08 (Calcium	H319
Chloride Dihydrate),	P280+P305+P351+P338
7732-18-5 (Water)	
Appearance	Clear solution
Cell culture tested-MTT	To pass test
Sterility	To pass test

Applications: Preparation of competent E. coli strains
Components: Calcium Chloride: 294g/ℓ
Recommended storage: RT

2M Magnesium Chloride Solution
Microbial Cell Culture

packaging	Mfr. No
5x10mℓ Poly Bottle	BP9741-10X5
MgCL ₂	MW: 95.21
CAS: 7791-18-6	H319
(Magnesium Chloride	P280+P305+P351+P338
Hexahydrate), 7732-18-5 (Water)	
Appearance	Clear solution
Cell culture tested-MTT	To pass test
Sterility	To pass test

Applications: Preparation of competent E. coli strains
Components: Magnesium Chloride: 406.6g/ℓ
Recommended storage: RT

2xYT Agar Molecular Genetics
Powder

packaging	Mfr. No
100 g PolyBottle	BP2466-100
500 g PolyBottle	BP2466-500
2 kg PolyBottle	BP2466-2
CAS: 9002-18-0	
Growth promotion ability E. coli	Found
Growth promotion ability S. aureus	Found

Applications: 2xYT Agar is used for growth of hosts for replication of M13 vectors and plating.
Components in g/ℓ:Tryptone: 16g; Yeast Extract: 10g; NaCl: 5g; Agar:15g
Recommended Storage: RT

2X YT Broth Molecular Genetics Powder

packaging	Mfr. No
100 g PolyBottle	BP2467-100
500 g PolyBottle	BP2467-500
2 kg PolyBottle	BP2467-2
10 kg PolyPail	BP2467-10

Growth promotion ability E. coli	Found
Growth promotion ability S. aureus	Found

Applications: 2X YT Broth is used for growth of hosts for replication of M13 vectors.
Components in g/ℓ:Tryptone: 16g; Yeast Extract: 10g; NaCl: 5g
Recommended Storage: RT

2 XYT Broth Capsules

packaging	Mfr. No
32x15.5 g Poly Bottle	BP9736-500
CAS: [Peptone: 91079-40-2; Sodium Chloride: 7647-14-5; Yeast Extract: 8013-01-2]	

Appearance	Capsules
Color (prepared media)	Light to medium amber
E.coli test	To pass test

Applications: Standard bacteriological media for the propagation and maintenance of E. coli strains in molecular biology procedures; dust-free granules allow for easy weighing.
Description: Contains 32 x 15.5g capsules; Use two capsules to make 1ℓ broth
Components: (per liter) Peptone: 16g; Yeast Extract: 10g; Sodium Chloride: 5g
Recommended storage: RT

Agar (80-100 Mesh)

packaging	Mfr. No
100 g PolyBottle	BP2641-100
500 g PolyBottle	BP2641-500
1 kg PolyBottle	BP2641-1
CAS: 9002-18-0	

Acid Insoluble Ash	<=0.5%
Arsenic	<=3ppm
Foreign Insoluble Matter	<=1.0%
Foreign Organic Matter	<=1.0%
Foreign Starch	To pass test
Gelatin	To pass test
Heavy Metals (Pb)	<=40ppm
Identity Test A	To pass test
Identity Test B	To pass test
Lead	<=10ppm
Microbial Limits (Salmonella)	Negative
Organic Volatile Impurities	To pass test
Total Ash (Dry Weight Basis)	<=6.5%
Water	<=20%
Water Absorption	To pass test

Agar Molecular Genetics Granulated

packaging	Mfr. No
500 g PolyBottle	BP1423-500
2 kg PolyPail	BP1423-2
CAS: 9002-18-0	

Gelling Point	32°-39°C
Growth promotion ability E. coli	Found
Growth promotion ability S. aureus	Found
Melting Point	83°-89°C

Applications: Agar, Granulated, is used as a solidifying agent in preparing microbiological culture media.
Recommended Storage: RT

LB Broth, Lennox Large Granules

packaging	Mfr. No
2 kg Poly Pail	BP9722-2
500 g Poly Bottle	BP9722-500
CAS: [Peptone: 91079-40-2; Sodium Chloride: 7647-14-5; Yeast Extract: 8013-01-2]	

Color (APHA)	Light beige
E.coli test	To pass test
particle size	1 to 2 mm

Applications: Standard bacteriological media for the propagation and maintenance of E. coli strains in molecular biology procedures; dust-free granules allow for easy weighing.
Components: (per liter) Peptone: 10g; Yeast Extract: 5g; Sodium Chloride: 5g
Recommended storage: RT

LB Broth, Lennox Pre-Buffered Capsules

packaging	Mfr. No
500 x 1 g Capsules	BP9733-500
CAS: [Peptone: 91079-40-2; Sodium Chloride: 7647-14-5; Yeast Extract: 8013-01-2; Tris/Tris HCl: 77-86-1; 1185-53-1] H315, H319 P280+P302+P352, P305+P351+P338	



Appearance	Capsules
Color (prepared media)	Yellow to gold
E.coli test	To pass test

Applications: Standard bacteriological media for the propagation and maintenance of E. coli strains in molecular biology procedures.
Components: (per liter) Peptone: 10g; Yeast Extract: 5g; Sodium Chloride: 5g; Tris/Tris-HCl: 1.5g
Recommended storage: RT

LB Broth, Lennox Capsules

packaging	Mfr. No
25x20 g Capsules	BP9730-500
CAS: [Peptone: 91079-40-2; Sodium Chloride: 7647-14-5; Yeast Extract: 8013-01-2; Tris/Tris HCl: 77-86-1; 1185-53-1] H315, H319 P280+P302+P352, P305+P351+P338	



Appearance	Capsules
Color (prepared media)	Yellow to gold
E.coli test	To pass test

Applications: Standard bacteriological media for the propagation and maintenance of E. coli strains in molecular biology procedures.
Components: (per liter) Peptone: 10g; Yeast Extract: 5g; Sodium Chloride: 5g; Tris/Tris-HCl: 1.5g
Recommended storage: RT

LB Broth, Lennox Molecular Genetics Powder

packaging	Mfr. No
500 g PolyBottle	BP1427-500
2 kg PolyPail	BP1427-2

Growth promotion ability E. coli	Found
Growth promotion ability S. aureus	Found

Applications: LB (Luria-Bertani) Broth, Lennox, is used for the culturing of E. coli in molecular biology procedures.
Components: (per liter) Tryptone 10g; Yeast Extract 5g; Sodium Chloride 5g
pH:7.0±0.2 at 25°C
Recommended Storage: RT

LB Broth, Miller Molecular Genetics Granular

packaging	Mfr. No
500 g PolyBottle	BP1426-500
2 kg PolyPail	BP1426-2

Growth promotion ability E. coli	Found
Growth promotion ability S. aureus	Found

Applications: LB (Luria-Bertani) Broth, Miller, is used for the culturing of E. coli in molecular biology procedures.
Components: (per liter) Tryptone 10g; Yeast Extract 5g; Sodium Chloride 10g
pH:7.0±0.2 at 25°C
Recommended Storage: RT

LB Agar, Miller Large Granules

packaging	Mfr. No
2 kg Poly Pail	BP9724-2
500 g Poly Bottle	BP9724-500
CAS: [Peptone: 91079-40-2; Sodium Chloride: 7647-14-5; Yeast Extract: 8013-01-2; Agar: 9002-18-0]	

Color (APHA)	Light beige
E.coli test	To pass test
particle size	1 to 2 mm

Applications: Standard bacteriological media for the propagation and maintenance of E. coli strains in molecular biology procedures; dust-free granules allow for easy weighing.
Components: (per liter) Peptone: 10g; Yeast Extract: 5g; Sodium Chloride: 10g; Agar: 12g
Recommended storage: RT

LB Agar, Miller Pre-Buffered Capsules

packaging	Mfr. No
500x1 g Capsules	BP9735-500
CAS: [Peptone: 91079-40-2; Sodium Chloride: 7647-14-5; Yeast Extract: 8013-01-2; Agar: 9002-18-0; Tris/Tris HCl: 77-86-1; 1185-53-1]	

Appearance	Capsules
Color (prepared media)	Yellow to gold
E.coli test	To pass test

Applications: Standard bacteriological media for the propagation and maintenance of E. coli strains in molecular biology procedures.
Components: (per liter) Peptone: 10g; Yeast Extract: 5g; Sodium Chloride: 10g; Agar: 15g; Tris/Tris-HCl: 1.5g
Recommended storage: RT

LB Agar, Miller Pre-Buffered Capsules

packaging	Mfr. No
25x20 g Capsules	BP9734-500
CAS: [Peptone: 91079-40-2; Sodium Chloride: 7647-14-5; Yeast Extract: 8013-01-2; Agar: 9002-18-0; Tris/Tris HCl: 77-86-1; 1185-53-1]	

Appearance	Capsules
Color (prepared media)	Yellow to gold
E.coli test	To pass test

Applications: Standard bacteriological media for the propagation and maintenance of E. coli strains in molecular biology procedures.
Components: (per liter) Peptone: 10g; Yeast Extract: 5g; Sodium Chloride: 10g; Agar: 15g; Tris/Tris-HCl: 1.5g
Recommended storage: RT

LB Broth, Miller Large Granules

packaging	Mfr. No
2 kg Poly Pail	BP9723-2
500 g Poly Bottle	BP9723-500
CAS: [Peptone: 91079-40-2; Sodium Chloride: 7647-14-5; Yeast Extract: 8013-01-2]	

Color (APHA)	Light beige
E.coli test	To pass test
particle size	1 to 2 mm

Applications: Standard bacteriological media for the propagation and maintenance of E. coli strains in molecular biology procedures; dust-free granules allow for easy weighing.
Components: (per liter) Peptone: 10g; Yeast Extract: 5g; Sodium Chloride: 10g
Recommended storage: RT



LB Broth, Miller		Pre-Buffered Capsules
packaging	Mfr. No	
500x1 g Capsules	BP9732-500	
CAS: [Peptone: 91079-40-2; Sodium Chloride: 7647-14-5; Yeast Extract: 8013-01-2; Tris/Tris HCl: 77-86-1; 1185-53-1] H315, H319 P280+P302+P352, P305+P351+P338		
Appearance		Capsules
Color (prepared media)		Yellow to gold
E.coli test		To pass test
Applications: Standard bacteriological media for the propagation and maintenance of E. coli strains in molecular biology procedures. Components: (per liter) Peptone: 10g; Yeast Extract: 5g; Sodium Chloride: 10g; Tris/Tris-HCl: 1.5g Recommended storage: RT		

NZY Broth Molecular Genetics Powder		
Packaging		Mfr. No
100 g	PolyBottle	BP2465-100
500 g	PolyBottle	BP2465-500
2 kg	PolyBottle	BP2465-2
10 kg	PolyPail	BP2465-10
Growth promotion ability E. coli		Found
Growth promotion ability S. aureus		Found
Applications: NZY Broth is used for the maintenance and propagation of bacteriophage lambda.		
Components: (per liter) Tryptone 10g; Yeast Extract 5g; NaCl 5g; MgSO ₄ 2g		
Recommended Storage: RT		

LB Broth, Miller		Pre-Buffered Capsules
packaging	Mfr. No	
25x20 g Capsules	BP9731-500	
CAS: [Peptone: 91079-40-2; Sodium Chloride: 7647-14-5; Yeast Extract: 8013-01-2; Tris/Tris HCl: 77-86-1; 1185-53-1] H315, H319 P280+P302+P352, P305+P351+P338		
Appearance		Capsules
Color (prepared media)		Yellow to gold
E.coli test		To pass test
Applications: Standard bacteriological media for the propagation and maintenance of E. coli strains in molecular biology procedures. Components: (per liter) Peptone: 10g; Yeast Extract: 5g; Sodium Chloride: 10g; Tris/Tris-HCl: 1.5g Recommended storage: RT		

Peptone Molecular Genetics Powder		
Packaging		Mfr. No
100 g	PolyBottle	BP1420-100
500 g	PolyBottle	BP1420-500
2 kg	PolyBottle	BP1420-2
Growth promotion ability E. coli Found		
Growth promotion ability S. aureus Found		
Applications: Peptone contains nitrogen in a form readily available to bacteria.		
Recommended Storage: RT		

NZCYM Broth Molecular Genetics Powder	
packaging	Mfr. No
100 g PolyBottle	BP2464-100
500 g PolyBottle	BP2464-500
2 kg PolyBottle	BP2464-2
10 kg PolyPail	BP2464-10
Growth promotion ability E. coli Found	
Growth promotion ability S. aureus Found	
Applications: NZCYM Broth is used for the maintenance and propagation of bacteriophage lambda. Components: (per liter) Tryptone 10g; Yeast Extract 5g; Casamino Acids 1g; NaCl 5g; MgSO ₄ 1g Recommended Storage: RT	

Peptone		Large Granules
Packaging		Mfr. No
500 g	Poly Bottle	BP9725-500
2 kg	Poly Pail	BP9725-2
CAS: 73049-73-7		
Ash		<=15%
Color (APHA)		Light beige
Moisture		<=5%
particle size		1 to 2 mm
pH of a 2% solution		Inclusive between 6.5-7.5
Total nitrogen content		>=10%
Applications: Standard bacteriological media for the propagation and maintenance of E. coli strains in molecular biology procedures; dust-free granules allow for easy weighing.		
Description: Mixture of peptides and free amino acids from pancreatic hydrolysis of animal tissue; pH (at 25°C, 2% Solution): 6.5 to 7.5		
Recommended storage: RT		

SOB Broth		Capsules
packaging	Mfr. No	
36x14 g Capsules	BP9737-500	
CAS: [Peptone: 91079-40-2; Sodium Chloride: 7647-14-5; Yeast Extract: 8013-01-2; Magnesium Sulfate: 7487-88-9; Potassium Chloride: 7447-40-7]		
Appearance	Capsules	
Color (prepared media)	Light to medium amber	
E.coli test	To pass test	
Applications: Standard bacteriological media for the propagation and maintenance of E. coli strains in molecular biology procedures.		
Components: (per liter) Peptone: 20g; Yeast Extract: 5g; Sodium Chloride: 0.5g; Magnesium Sulfate: 2.4g; Potassium Chloride: 0.186g		
Recommended storage: RT		

SOC Medium	
packaging	Mfr. No
5x10mℓ Poly Bottle	BP9740-10X5
CAS: [Peptone: 91079-40-2; Sodium Chloride: 7647-14-5; Yeast Extract: 8013-01-2; Magnesium Sulfate: 7487-88-9; Potassium Chloride: 7447-40-7; Magnesium Chloride: 7786-30-3; Glucose: 50-99-7]	
Color (APHA) Yellow to light brown	
pH test To pass test	
Sterility To pass test	
Applications: Standard bacteriological media for the propagation and maintenance of E. coli strains in molecular biology procedures. Components: Peptone: 20g/ℓ, Yeast Extract: 5g/ℓ, Sodium Chloride: 0.58g/ℓ, Potassium Chloride: 0.19g/ℓ, Magnesium Chloride: 0.95g/ℓ, Magnesium Sulfate: 1.2g/ℓ, Glucose: 3.6g/ℓ Recommended storage: 4°C	

Terrific Broth Molecular Genetics Powder	
packaging	Mfr. No
100 g PolyBottle	BP2468-100
500 g PolyBottle	BP2468-500
2 kg PolyBottle	BP2468-2
10 kg PolyPail	BP2468-10
Growth promotion ability E. coli Found	
Growth promotion ability S. aureus Found	
Applications: Terrific broth is a buffered enriched medium used to obtain high yields of plasmid DNA from transformed E. coli. Components: (per liter) Tryptone 12g; Yeast Extract 24g; K2HPO4 12.5g; KH2PO4 2.3g Recommended Storage: RT	

Terrific Broth, Modified		Capsules
packaging	Mfr. No	
25x24 g Poly Bottle	BP9729-600	
CAS: [Peptone: 91079-40-2; Yeast Extract: 8013-01-2; Potassium Phosphate, Monobasic: 7778-77-0; Potassium Phosphate, Dibasic: 7758-11-4]		
Appearance	Capsules	
Color (prepared media)	Light to medium amber	
E.coli test	To pass test	
Applications: Standard bacteriological media for the propagation and maintenance of E. coli strains in molecular biology procedures.		
Description: Contains 25 x 24g capsules; Use 2 capsules to make 1ℓ broth		
Components: (per liter) Peptone: 12g; Yeast Extract: 24g; Potassium Phosphate, Monobasic: 2.2g; Potassium Phosphate, Dibasic: 9.4g.		
Recommended storage: RT		

Terrific Broth, Modified		Large Granules
packaging		Mfr. No
2 kg	Poly Pail	BP9728-2
500 g	Poly Bottle	BP9728-500
CAS: [Peptone: 91079-40-2; Yeast Extract: 8013-01-2; Potassium Phosphate, Monobasic: 7778-77-0; Potassium Phosphate, Dibasic: 7758-11-4]		
Color (APHA)		Beige to yellow
E.coli test		To pass test
particle size		1 to 2 mm
Applications: Standard bacteriological media for the propagation and maintenance of E. coli strains in molecular biology procedures.		
Components: (per liter) Peptone: 12g; Yeast Extract: 24g; Potassium Phosphate, Monobasic: 2.2g; Potassium Phosphate, Dibasic: 9.4g.		
Recommended storage: RT		

Tryptone		Large Granules
packaging		Mfr. No
2 kg	Poly Pail	BP9726-2
500 g	Poly Bottle	BP9726-500
CAS: 91079-40-2		
Ash		<=15%
Color (APHA)		Light beige
Moisture		<=5%
particle size		1 to 2 mm
pH of a 2% solution		Inclusive between 6.5-7.5
Total nitrogen content		>=10%
Applications: Standard bacteriological media for the propagation and maintenance of E. coli strains in molecular biology procedures; dust-free granules allow for easy weighing.		
Description: An enzymatic digest of casein, a milk protein; pH (at 25°C, 2% Solution): 6.5 to 7.3		
Recommended storage: RT		

Yeast Extract Molecular Genetics Powder

packaging	Mfr. No
100 g PolyBottle	BP1422-100
500 g PolyBottle	BP1422-500
2 kg PolyBottle	BP1422-2

CAS: 8013-01-2
EINECS: 232-387-9

Growth promotion ability E. coli

Found

Growth promotion ability S. aureus

Found

Applications: Yeast Extract is the water-soluble portion of autolyzed yeast used in bacteriological culture media.

Recommended Storage: RT

n-Butanol Clear Colorless Liquid

packaging	Mfr. No
25 mℓ AmberGlass	BP505-25
500 mℓ AmberGlass,EcoSafPak*	BP505-500

C₄H₁₀O
CAS: 71-36-3
MW: 74.12
EINECS: 200-751-6
H226, H315, H302, H336,

Acidity <=0.005%
Aldehyde (HCHO) None
Assay >=99.8%
Chloride (Cl) None
Color (APHA) <=5
iso-Butanol <=0.1%
Water <=0.1%

Applications: n-Butanol is used as an analytical reagent and solvent.

Recommended Storage: RT

UN 1120; DOT Class 3:Flammable Liquid

EcoSafPak* is an environmentally friendly packaging system made of 100% recyclable material by an SFI certified supplier.

Yeast Extract Large Granules

packaging	Mfr. No
2 kg Poly Pail	BP9727-2
500 g Poly Bottle	BP9727-500

CAS: 8013-01-2
EINECS: 232-387-9

Color (APHA) Light beige to brown
Moisture <=6%
particle size 1 to 2 mm
pH of a 2% solution Inclusive between 5.5-7.5
Total nitrogen content >=8%

Applications: Standard bacteriological media for the propagation and maintenance of E. coli strains in molecular biology procedures; dust-free granules allow for easy weighing.

Description: Dried yeast autolysate with a high content of amino nitrogen and water soluble B complex vitamins

Recommended storage: RT

ACES Fine White Crystals

packaging	Mfr. No
100 g AmberGlass	BP303-100

C₄H₁₀N₂O₄S
CAS: 7365-82-4
MW: 182.19
EINECS: 230-908-4

Assay >=96%
IR Conforms to standard
PKa (at 25°C) 6.78 ±0.2

Applications: ACES biological buffer has a usable pH range of 6.1 to 7.5.

Recommended Storage: RT

YPD Broth Molecular Genetics Powder

packaging	Mfr. No
100 g PolyBottle	BP2469-100
500 g PolyBottle	BP2469-500
2 kg PolyBottle	BP2469-2
10 kg PolyPail	BP2469-10

Growth promotion ability E. coli Found
Growth promotion ability S. aureus Found

Applications: YPD Broth is used for the maintenance and propagation of yeast.

Recommended Storage: RT

Amido Black 10B

packaging	Mfr. No
10 g AmberGlass	BP124-10

C₂₂H₁₄N₆Na₂O₉S₂
CAS: 1064-48-8
MW: 616.48
EINECS: 213-903-1

E¹%_{1cm} >=450ℓ g-1cm-1
Lambda Max. 618nm±2nm
Loss on Drying (at 105°C) <=10%
Solubility To pass test

Recommended storage: RT

C.I. 20470

Ammonium Bicarbonate Fine, White Crystalline Powder

packaging	Mfr. No
500 g PolyBottle	BP2413-500

CH₃NO₃
CAS: 1066-33-7
MW: 79.06
EINECS: 213-911-5

Assay (as NH₃) 21.30-21.73
Chloride (Cl) <=5ppm
Heavy Metals (Pb) <=5ppm
Iron <=5ppm
Residue after ignition <=0.005%
Sulfate (SO₄) <=0.002%

Applications: Ammonium Bicarbonate is used for many biochemical applications.

Recommended Storage: RT

Ammonium Carbonate Certified ACS Lumps

packaging	Mfr. No
500 g Amber Glass	BP2414-500

CH₈N₂O₃
CAS: 8000-73-5
MW: 96.09
Hartshorn salt

Assay (as NH₃) >=30%
Chloride (Cl) <=5ppm
Heavy Metals (Pb) <=5ppm
Insoluble matter <=0.005%
Iron (Fe) <=3 ppm
Non-volatile matter <=0.005%
Sulfur compounds (as SO₄) <=0.002%

Applications: Ammonium Carbonate is used in many biochemical procedures.

Description: This product is a mixture of Ammonium bicarbonate and Ammonium carbonate.

Recommended storage: RT

Ammonium Phosphate, Monobasic White Crystalline Powder

packaging	Mfr. No
500 g PolyBottle	BP2427-500

H₂NO₄P
CAS: 7722-76-1
MW: 115.03
EINECS: 231-764-5

Assay >=98.0%
Calcium (Ca) <=0.001%
Chloride (Cl) <=5ppm
Heavy Metals (Pb) <=5ppm
Insoluble matter <=0.005%
Iron <=0.001%
Magnesium (Mg) <=0.0005%
Nitrate <=0.001%
pH of 5% Solution (at 25°C) 3.8-4.4
Potassium <=0.005%
Sodium <=0.005%
Sulfate (SO₄) <=0.01%

Applications: Ammonium Phosphate, Monobasic is used in several biochemical research applications.

Recommended Storage: RT

Ammonium Phosphate, Dibasic Colorless Crystals

packaging	Mfr. No
500 g PolyBottle	BP361-500

H₂N₂O₄P
CAS: 7783-28-0
MW: 132.06
H335, H315, H319

Assay 98.0-101.0%
Calcium (Ca) <=0.001%
Chloride (Cl) <=5ppm
Heavy Metals (Pb) <=0.001%
Insoluble matter <=0.005%
Iron <=0.001%
Magnesium (Mg) <=0.0005%
Nitrate <=0.003%
pH of 5% Solution (at 25°C) 7.7 to 8.1
Potassium <=0.005%
Sodium <=0.005%
Sulfate (SO₄) <=0.01%

Applications: Ammonium Phosphate Dibasic is suitable for use in many biological applications.

Recommended Storage: RT

BES White Crystalline Powder

packaging	Mfr. No
25 g PolyBottle	BP501-25
100 g PolyBottle	BP501-100
500 g PolyBottle	BP501-500

C₆H₁₃NO₃S
CAS: 10191-18-1
MW: 213.25
EINECS: 233-465-5

Assay (as free acid) >=99%
DNase Not detected
Loss on Drying (at 105°C) <=1%
PKa (at 25°C) 7.1 ±0.1
Protease Not detected
RNase Not detected

Applications: BES is a biological buffer with a usable pH range of 6.4 to 7.8.

Recommended Storage: RT

Bicine White Powder

packaging	Mfr. No
100 g PolyBottle	BP2646-100
250 g PolyBottle	BP2646-250
1 kg PolyBottle	BP2646-1

C₆H₁₃NO₄
CAS: 150-25-4
MW: 163.17
EINECS: 205-755-1

FTIR Conforms to standard
OD260nm (10% aqueous solution) <0.100
OD280nm (10% aqueous solution) <0.100
pH of 1% Aqueous Solution 4.5-5.5

Applications: Bicine is one of the zwitterionic amino acids known as "Good" buffers. Used as a biological buffer and chelating agent.

Recommended Storage: RT



Bis-Tris

Fine White Crystals

packaging

Mfr. No

100 g PolyBottle

BP301-100

C₈H₁₉NO₃
CAS: 6976-37-0
MW: 209.24
EINECS: 230-237-7
H290, H315, H319, H335

P234, P390, P261,
P302+P352, P280,
P305+P351+P338

!

Assay

>=99%

Melting Point

104°C ±2°C

PKa (at 25°C)

6.5 ±0.2

Applications: This biological buffer has a usable pH range of 5.8 to 7.2.

Recommended Storage: RT

CellPURE™ Bis-Tris Propane

Cell Culture-Tested

White crystalline powder

packaging

Mfr. No

100 g Poly Bottle

BP2943-100

C₁₁H₂₆N₂O₆
CAS: 64431-96-5
MW: 282.33

Absorbance 33% w/w at 290nm

<=0.05

Cell culture test

To pass test

DNase

Not detected

Endotoxin

<=1.0 EU/mg

Heavy metals (as Pb)

<=5 ppm

Infrared scan

Conforms to ref.

Purity (by titration)

>=99.0%

Nickase

Not detected

NMR

<=0.5%

Protease

Not detected

RNase

Not detected

Solubility

Clear and colorless

Standard plate count

<=100CFU/gm

Water (Karl Fischer)

<=0.2%

Applications: Cell cultivation, isolation of cells, enzyme assays, and other biochemical applications

Description: Cell culture tested; Analyzed for the absnce of nucleases and proteases; Tested for endotoxin and bioburden levels; pKa values mostly independent of temperature and concentration; High water solubility; Minimal permeability to biological membranes; Manufactured under strict quality control guidelines to ensure performance and reliability.

Recommended storage: RT

JustPURE™ BIS-TRIS Propane

Cellular and Molecular Biology - White Powder

packaging

Mfr. No

25 g Poly Bottle

BP2929-25

C₁₁H₂₆N₂O₆
CAS: 64431-96-5
MW: 282.33

Aluminum (Al)

<=0.0005%

Ammonia

<=0.05%

Calcium (Ca)

<=0.0005%

Chloride

<=0.05%

Copper (Cu)

<=0.0005%

Iron (Fe)

<=0.0005%

Lead (Pb)

<=0.001%

Magnesium (Mg)

<=0.0005%

Purity (by titration)

>=99.0%

Optical Absorbance at 260nm

<=0.1

Optical Absorbance at 280nm

<=0.05

pH at 20°C

Inclusive between 10.0-12.0

Phosphorus

<=0.005%

Residue after ignition

<=0.1%

Sodium (Na)

<=0.005%

Solubility

Clear and colorless

Sulfate (SO₄)

<=0.005%

Applications: For high efficiency transfection of mammalian cells, gel electrophoresis of RNA, protein isolation applications, cell cultures and enzyme assays and bioanalytical methods, such as IEF, 2-D electrophoresis and SDS-PAGE.

Description: Ultra-pure zwitterionic buffers; Optimized for research in cellular and molecular biology; Minimal permeability to biological membranes; Low interference with biological reactions; pKa values mostly independent of temperature and concentration; High water solubility; Minimal absorption in spectral range 240 to 700nm.

Recommended storage: RT

Calcium Chloride Dihydrate

White Crystals to Powder

packaging

Mfr. No

100 g PolyBottle

BP510-100

250 g PolyBottle

BP510-250

500 g PolyBottle

BP510-500

CaCl₂·2H₂O
CAS: 10035-04-8
MW: 147.02

H319
P280, P305+P351+P338

!

Assay

>=99-103.3%

DNase

Not detected

pH (5% solution in H₂O, at 25°C)

4.5-8.5

Protease

Not detected

RNase

Not detected

Applications: Calcium Chloride is suitable for use in the preparation of buffer solutions.

Recommended Storage: RT

Calcium Phosphate Dibasic Anhydrous

White Powder

packaging

Mfr. No

500 g AmberGlass

BP441-500

CaHO₄P
CAS: 7757-93-9
MW: 136.06

Arsenic

<=8ppm

Assay (as Ca)

30.0-31.7%

Barium

To pass test

Chloride (Cl)

<=0.25%

Fluoride

<=0.005%

Heavy Metals (Pb)

<=0.003%

Loss on ignition

6.6-8.5%

Sulfate (SO₄)

<=0.5%

Applications: Calcium Phosphate can be used for column chromatography of proteins.

Recommended Storage: RT

CAPS

Fine White Crystals

packaging

Mfr. No

100 g AmberGlass

BP321-100

500 g AmberGlass

BP321-500

C₈H₁₉NO₃S
CAS: 1135-40-6
MW: 221.31
EINECS: 214-492-1
H335, H319, H315

P261, P301+P312,
P302+P352, P280,
P305+P351+P338

!

Assay

>=98%

IR

Conforms to standard

PKa (at 25°C)

10.4 ±0.2

Applications: This biological buffer has a usable pH range of 9.7 to 11.1.

Recommended Storage: RT

CHES

Fine White Crystals

packaging

Mfr. No

100 g PolyBottle

BP318-100

C₈H₁₇NO₃S
CAS: 103-47-9
MW: 207.29
EINECS: 203-115-6

H335, H315, H319
P261, P302+P352, P280,
P305+P351+P338

!

Assay

>=98%

DNase

Not detected

PKa (at 25°C)

9.3 ±0.2

RNase

Not detected

Applications: This biological buffer has a usable pH range of 8.6 to 10.0.

Recommended Storage: RT

Cytoscint* Scintillation Cocktail

Molecular Biology

packaging

Mfr. No

4 ℓ AmberGlass,EcoSafPak*

BP458-4

3H Background

<=50cpm

3H Efficiency

>=50%

Aqueous PBS Buffer Capacity

>=4.5%

Distilled Water Capacity

>=8.5%

GC, IR, TLC

Conforms to standard

Specific Gravity

0.88 ±0.02

Applications: Universal LSC Cocktail.

Components: 65-95% Pseudocumene, 0-35% Ethoxylated Nonylphenol, <0.1% 2.5-Diphenyloxazole, <1% 1.4-bis[2-(2-methylphenyl)ethenyl]-benzene [95-63-6 (Pseudocumene)] ; [9016-45-9 (Ethoxylated Nonylphenol)] ; [92-71-7 (2.5-Diphenyloxazole)] ; [13280-61-0 (1.4-Bis[2-(2-methylphenyl)ethenyl]-benzene)]

UN 1993; DOT Class 3:Flammable Liquid

EcoSafPak* is an environmentally friendly packaging system made of 100% recyclable material by an SFI certified supplier.

Deoxycholic Acid Sodium Salt

packaging

Mfr. No

100 g PolyBottle

BP349-100

C₂₄H₃₉NaO₄
CAS: 302-95-4
MW: 414.58

EINECS: 206-132-7

Arsenic

<=0.001%

Assay (dry basis)

>=99%

Lead

<=0.001%

Moisture

<=5%

Sodium Choleate

<=2%

Applications: This anionic detergent is useful for solubilizing and isolating membrane-bound proteins in an active state.

Recommended Storage: RT

Dimethyl Sulfoxide

packaging

Mfr. No

100 mL AmberGlass

BP231-100

1 ℓ AmberGlass,EcoSafPak*

BP231-1

4 ℓ AmberGlass,EcoSafPak*

BP231-4

C₂H₆OS
CAS: 67-68-5
MW: 78.13

EINECS: 200-664-3

A280nm

<=0.10

Assay

>=99.7%

Color (APHA)

<=10

Density (at 25°C)

1.095 ±0.005g/mL

Freezing Point

>=18.0°C

Residue after evaporation

<=0.001%

Water

<=0.2%

Applications: DMSO is a cryopreservation agent for cells, and may also be used in bacterial transformations.

Note: Sold for nonmedical use only. Not to be used or sold as a drug for humans or animals.

Recommended Storage: RT

EcoSafPak* is an environmentally friendly packaging system made of 100% recyclable material by an SFI certified supplier.

480

481



Electrode Storage Solution			
Packaging			Mfr. No.
1 ℓ	PolyBottle		BP2418-1
pH (at 25°C)			3.9 ±0.1
Potassium chloride			6.2 ±0.3g/ℓ
Residue after evaporation			1.65 ±0.08%

Applications: Electrode Storage Solution is a safe and convenient method to maintain sensitivity by removing buildup from electrodes.
For use in storage of pH electrodes.
Recommended Storage: RT

JustPURE™ EPPS		Cellular and Molecular Biology	
White Powder			
packaging		Mfr. No	
100 g	Poly Bottle	BP2933-100	
C ₉ H ₂₀ N ₂ O ₄ S		EINECS: 240-198-8	
CAS: 16052-06-5			
MW: 252.33			
Ammonia		<=0.1%	
Bromide		<=0.1%	
Calcium (Ca)		<=0.002%	
Copper (Cu)		<=0.0005%	
Iron (Fe)		<=0.0005%	
Lead (Pb)		<=0.001%	
Magnesium (Mg)		<=0.0005%	
Purity (by titration)		>=99.5%	
Optical Absorbance at 260nm		<=0.1	
Optical Absorbance at 280nm		<=0.1	
pH at 20°C		Inclusive between 5.0-7.0	
Phosphorus		<=0.005%	
Residue after ignition		<=0.1%	
Sodium (Na)		<=0.005%	
Solubility		Clear and colorless	

Applications: For high efficiency transfection of mammalian cells, gel electrophoresis of RNA, protein isolation applications, cell cultures and enzyme assays and bioanalytical methods, such as IEF, 2-D electrophoresis and SDS-PAGE.
Description: Ultra-pure zwitterionic buffers; Optimized for research in cellular and molecular biology; Minimal permeability to biological membranes; Low interference with biological reactions; pKa values mostly independent of temperature and concentration; High water solubility; Minimal absorption in spectral range 240 to 700nm.
Recommended storage: RT

Ethylenediamine Tetraacetate, DEPC-Treated Solution			
0.5M, pH 8.0			
packaging		Mfr. No	
100 mℓ	PolyBottle	BP2483-100	
500 mℓ	PolyBottle	BP2483-500	
1 ℓ	PolyBottle	BP2483-1	
CAS: 60-00-4			
H319			
P280, P305+P351+P338			
			
Molarity			0.475-0.525M
pH			8.0 ± 0.1
Protease			Not detected
RNase			Not detected

Applications: EDTA, DEPC-Treated solution is ideal for RNA work.
[60-00-4 (Ethylenediamine Tetraacetic Acid, <14.0%)] ; [1310-73-2 (Sodium Hydroxide, <4.0%)] ; [1609-47-8 (Diethyl Pyrocarbonate, 1.0%)]
Recommended Storage: RT

Ethylenediamine Tetraacetic Acid			
White Crystalline Powder			
packaging		Mfr. No	
500 g	AmberGlass	BP118-500	
C ₁₀ H ₁₆ N ₂ O ₈		P264, P280,	
CAS: 60-00-4		P305+P351+P338, P337+P313	
MW: 292.23			
H319, H332			
Assay		>=99.5%	
Insoluble in Dilute NH ₄ OH		Clear and haze-free	
Iron		<=0.005%	
Magnesium (Mg)		<=0.0005%	
Nitrilotriacetic Acid (N(CH ₂ COOH) ₃)		<=0.1%	
pH of 5% Solution (at 25°C)		2.5-3.5	
Residue after ignition		<=0.2%	

Recommended storage: RT

Ethylenediamine Tetraacetic Acid, Dipotassium Salt Dihydrate			
White Crystalline Powder			
packaging		Mfr. No	
500 g	AmberGlass	BP119-500	
C ₁₀ H ₁₄ K ₂ N ₂ O ₈ ·2H ₂ O		P261, P302+P352, P280,	
CAS: 25102-12-9		P305+P351+P338	
MW: 404.46			
H335, H319, H315			
Assay		>=99.0%	
Chloride (Cl)		<=0.02%	
Heavy Metals (Pb)		<=0.001%	
Iron		<=0.001%	
pH of 5% Solution (at 25°C)		4-5	

Recommended storage: RT

Ethylenediamine Tetraacetic Acid, Disodium Salt Dihydrate			
Crystalline Powder - Electrophoresis			
packaging		Mfr. No	
500 g	PolyBottle	BP120-500	
1 kg	PolyBottle	BP120-1	
C ₁₀ H ₁₄ N ₂ Na ₂ O ₈ ·2H ₂ O		H332	
CAS: 6381-92-6		P304+P312	
MW: 372.23			
Arsenic		<=2ppm	
Assay		>=99%	
Calcium (Ca)		<=20ppm	
Electrophoresis		To pass test	
Insoluble matter		<=0.005%	
Iron		<=10ppm	
Lead		<=5ppm	
Mercury		<=2ppm	
Nitrilotriacetic Acid (N(CH ₂ COOH) ₃)		<=0.1%	
pH of 5% Solution in Water (at 25°C)		4.0-6.0	

Recommended storage: RT

Ethylenediamine Tetraacetic Acid, Tetrasodium Salt Dihydrate			
White Powder			
packaging		Mfr. No	
500 g	AmberGlass	BP121-500	
C ₁₀ H ₁₂ N ₂ Na ₄ O ₈ ·2H ₂ O		P261, P280,	
CAS: 10378-23-1		P301+P330+P331, P310,	
MW: 416.2		P305+P351+P338	
H302, H332, H318			
Assay		>=98.0%	
Assay (dry basis)		>=99%	
Insoluble matter		<=0.005%	
Nitrilotriacetic Acid (N(CH ₂ COOH) ₃)		<=0.2%	
Nitrolotriacetic acid		<=0.2%	
pH of 1% w/w Solution		10.7-11.7	
pH of a 1% w/w solution		Inclusive between 10.7-11.7	

Recommended storage: RT

Ficoll* 400		Molecular Biology
White to Off-white Powder		
packaging		Mfr. No
5 g	PolyBottle	BP525-5
10 g	PolyBottle	BP525-10
25 g	PolyBottle	BP525-25
100 g	PolyBottle	BP525-100
500 g	PolyBottle	BP525-500
CAS: 26873-85-8		
MW: 400.000		
Specific Rotation α^{20}_D (c=4, H ₂ O)		52°±2°
DNase		Not detected
Protease		Not detected
RNase		Not detected

Applications: Ficoll 400 is used as a component of Denhardt’s Reagent and in gel-loading buffers for DNA gels.
It is also used in forming density gradients to separate cells and subcellular components.
Recommended Storage: RT

Formamide Super Pure			
packaging		Mfr. No	
100 mℓ	AmberGlass	BP228-100	
CH ₃ NO		H360D	
CAS: 75-12-7		P201, P308+P313	
MW: 45.04			
EINECS: 200-842-0			
Assay		>=99.5%	
Conductivity		<=100µmhos/cm	
Copper		<=0.0001%	
Freezing Point		1.0°-3.0°C	
Iron		<=0.0005%	
Lead		<=0.0005%	
Optical Absorbance at 280nm		<0.05	
Zinc		<=0.0005%	

H2NCH:O: CH3NO; F.W. 45.04
Vacuum distilled and packaged under Nitrogen.
Applications: Super Pure Formamide has been pretreated with a mixed-bed resin so that it is ready for use in nucleic acid hybridization and sequencing.
Recommended Storage: 4°C

Glycerol		Molecular Biology	
packaging		Mfr. No	
1 ℓ	PolyBottle,EcoSafPak*	BP229-1	
4 ℓ	PolyBottle,EcoSafPak*	BP229-4	
C ₃ H ₈ O ₃		EINECS: 200-289-5	
CAS: 56-81-5			
MW: 92.09			
Assay		>=99.5%	
Calcium (Ca)		<=5ppm	
Color (APHA)		<=10	
Density (at 20°C)		1.260 ±0.005g/mℓ	
DNase		Not detected	
Heavy Metals (Pb)		<=5ppm	
Iron		<=5ppm	
Magnesium (Mg)		<=1ppm	
Protease		Not detected	
RNase		Not detected	
Zinc		<=2ppm	

Applications: Glycerol is used for low-temperature storage of enzymes and bacterial cultures, purification of proteins, and as a component in electrophoresis loading buffers.
Recommended Storage: RT
EcoSafPak* is an environmentally friendly packaging system made of 100% recyclable material by an SFI certified supplier.

JustPURE™ Gly-Gly		Cellular and Molecular Biology	
White Powder			
packaging		Mfr. No	
50 g	Poly Bottle	BP2930-50	
C ₄ H ₈ N ₂ O ₃		EINECS: 209-127-8	
CAS: 556-50-3		Diglycine, Glycyl-glycine	
MW: 132.12			
Ammonia		<=0.02%	
Arsenic (As)		<=0.00001%	
Bismuth (Bi)		<=0.0005%	
Cadmium (Cd)		<=0.0005%	
Calcium (Ca)		<=0.001%	
Chloride		<=0.005%	
Chromium (Cr)		<=0.0005%	
Cobalt (Co)		<=0.0005%	
Copper (Cu)		<=0.0005%	
Iron (Fe)		<=0.0005%	
Lead (Pb)		<=0.0005%	
Lithium (Li)		<=0.0005%	
Loss on Drying (at 110°C)		<=0.05	
Magnesium (Mg)		<=0.0005%	
Manganese (Mn)		<=0.0005%	
Purity (by titration)		>=99.5%	
Optical Absorbance at 260nm		<=0.075	
Optical Absorbance at 280nm		<=0.072	
pH at 20°C		Inclusive between 4.5-6.0	
Sodium (Na)		<=0.005%	
Solubility		Clear and colorless	
Strontium (Sr)		<=0.0005%	
Sulfate (SO ₄)		<=0.005%	

Applications: For high efficiency transfection of mammalian cells, gel electrophoresis of RNA, protein isolation applications, cell cultures and enzyme assays and bioanalytical methods, such as IEF, 2-D electrophoresis and SDS-PAGE.
Description: Ultra-pure zwitterionic buffers; Optimized for research in cellular and molecular biology; Minimal permeability to biological membranes; Low interference with biological reactions; pKa values mostly independent of temperature and concentration; High water solubility; Minimal absorption in spectral range 240 to 700nm.
Recommended storage: RT

HEPES Buffer		Molecular Biology
1M Solution, pH 7.3		
packaging		Mfr. No
100 mℓ	PolyBottle	BP299-100
500 mℓ	PolyBottle	BP299-500
1 ℓ	PolyBottle	BP299-1
C ₈ H ₁₈ N ₂ O ₄ S		EINECS: 230-907-9
CAS: 7365-45-9		
MW: 238.31		
Absorbance of a 0.1M Solution at 260nm		<=0.008
Absorbance of a 0.1M Solution at 280nm		<=0.004
DNase		Not detected
Molarity		1.0±0.05
pH (at 37°C)		7.3±0.1
Protease		Not detected
RNase		Not detected

Applications: HEPES is commonly used as a buffering agent.
Recommended Storage: RT

HEPES Sodium Salt
White Powder

Packaging		Mfr. No.
500 g	PolyBottle	BP410-500
1 kg	PolyBottle	BP410-1
2.5 kg	PolyBottle	BP410-212
C ₈ H ₁₇ N ₂ NaO ₄ S		EINECS: 278-169-7
CAS: 75277-39-3		
MW: 260.27		
Absorbance of a 0.1M Solution in Water (1cm) at 260nm		<=0.04
Absorbance of a 0.1M Solution in Water (1cm) at 280nm		<=0.02
Assay (dry weight basis)		>=99%
Endotoxin Assay (LAL)		<=0.06 EU/ml
Heavy Metals (Pb)		<=0.001%
Loss on Drying (at 110°C)		<=3.0%
pH of 1% Solution (at 25°C)		10.2 ±0.8
pKa (at 25°C)		7.55 ±0.1
Solubility (0.1M in water)		Clear, complete & colorless

Applications: HEPES Sodium Salt is used as a buffer for in vitro cell culture.
Recommended Storage: RT

HEPES
Fine White Crystals

packaging		Mfr. No
100 g	PolyBottle	BP310-100
500 g	PolyBottle	BP310-500
1 kg	PolyBottle	BP310-1
5 kg	PolyPail	BP310-5
C ₈ H ₁₈ N ₂ O ₄ S		EINECS: 230-907-9
CAS: 7365-45-9		
MW: 238.3		
Absorbance of a 0.1M Solution at 280nm		<=0.01
Assay		>=99%
DNase		Not detected
pKa (at 25°C)		7.5±0.2
Protease		Not detected
RNase		Not detected
Solubility (1.196g/100ml H2O)		Clear and colorless
TLC		Single spot

Applications: HEPES is commonly used as a buffering agent.
Recommended Storage: RT

Hyamine Hydroxide*
Clear Solution

Packaging		Mfr. No.
500 mL	AmberGlass	BP2658-500
C₂₈H₄₅NO₃ CAS: 26248-39-5 MW: 443.66		EINECS: 247-537-9

Recommended Storage: RT
Not on TSCA inventory: for R and D use only; not for manufacturing or commercial purposes.
UN 2924; DOT Class 3:Flammable Liquid

Hydrogen Peroxide
30% in Water

Packaging		Mfr. No
500 g	PolyBottle	BP2633-500
H ₂ O ₂	H318, H302	
CAS: 7722-84-1	P280, P305+P351+P338,	
MW: 34	P301+P312	
EINECS: 231-765-0		

Ammonium (NH ₄)	<=5µg/g
Assay (as H ₂ O ₂ , %wt)	29.0-32.0%
Chloride (Cl)	<=3µg/g
Color (APHA)	<=10
Heavy Metals (Pb)	<=1µg/g
Iron	<=0.5µg/g
Nitrate	<=2µg/g
Phosphate (PO ₄)	<=2µg/g
Residue after evaporation	<=0.002%
Sulfate (SO ₄)	<=5µg/g
Titrateable Acid	<=0.0006mEq/g

Applications: Hydrogen Peroxide is used as an antiseptic and for decomposition of inorganic and organic samples.
Recommended Storage: RT
UN 2014; DOT Class 5.1:Oxidizer

Imidazole
Molecular Biology

Packaging		Mfr. No.
50 g	polyBottle	BP305-50
C ₃ H ₄ N ₂ CAS: 288-32-4 MW: 68.08 H302, H314, H360D		P301+P330+P331, P280, P305+P351+P338, P310, P281, P301+P312

Absorbance of a 0.1M Solution at 260nm	<=0.2
Assay	>=99%
DNase	Not detected
Melting Point	90°±2°C
PKa (at 25°C)	6.99±0.2
Protease	Not detected
RNase	Not detected

HNCH:CHN:CH; C3H4N2; F.W. 68.08
Applications: Imidazole is an excellent biological buffer with a usable pH range of 6.2 to 7.8.
Recommended Storage: RT
UN 3263; DOT Class 8:Corrosive

Isoamyl Alcohol

Packaging		Mfr. No
500 mL	AmberGlass,EcoSafPak*	BP1150-500
C₅H₁₂O CAS: 123-51-3 MW: 88.15 EINECS: 204-633-5 H226, H335, H332, EUH066 P261, P304+P340, P210		 

Acidity (as Acetic Acid)	<=0.01%
Assay (GC)	>=99%
Carbonyl (as C=O)	<=0.1%
Color (APHA)	20
Peroxides	Not detected
Residue after evaporation	<=0.003%
Water	<=0.5%

EcoSafPak is an environmentally friendly packaging system made of 100% recyclable material by an SFI certified supplier.*
Applications: A mixture of Chloroform and Isoamyl Alcohol (24:1 v/v) is used to remove proteins during nucleic acid extraction procedures. It may also be used to extract Ethidium Bromide from nucleic acids.
Recommended storage: RT

Magnesium Acetate Tetrahydrate

packaging	Mfr. No
500 g AmberGlass	BP215-500
C₄H₆O₄Mg.4H₂O CAS: 16674-78-5 MW: 214.46	
Assay	>=99.0%
Calcium (Ca)	<=0.01%
Chloride (Cl)	<=0.001%
Heavy Metals (as Pb) (ACS method)	<=5ppm
Iron	<=5ppm
Nitrogen Compounds	<=0.001%
pH of 5% Solution (at 25°C)	6.5-8.5
Sulfate (SO ₄)	<=0.005%
Water-insoluble Matter	<=0.005%

Applications: Magnesium Acetate is suitable for use in the preparation of buffer solutions and is required as a stabilizer for all stages of bacteriophage purification.
Recommended Storage: RT

Magnesium Chloride Hexahydrate

Packaging		Mfr. No.
500 g	PolyBottle	BP214-500
Cl ₂ Mg.6H ₂ O		
CAS: 7791-18-6		
MW: 203.31		
H335, H319 P261, P280, P305+P351+P338		
Ammonium (NH ₄)		<=0.002%
Assay		>=99.0%
Calcium (Ca)		<=0.01%
Iron		<=0.0005%
Lead		<=5ppm
Nitrate		<=0.001%
Phosphate (PO ₄)		<=5ppm
Sulfate (SO ₄)		<=0.002%
Water-insoluble Matter		<=0.005%

Applications: Magnesium Chloride is suitable for use in the preparation of buffer solutions.
Recommended Storage: RT

Magnesium Sulfate Heptahydrate

Packaging		Mfr. No
1 kg	AmberGlass	BP213-1
MgO ₄ S.7H ₂ O		
CAS: 10034-99-8		
MW: 246.48		

Ammonium (NH ₄)	<=0.002%
Assay	98-102.0%
Calcium (Ca)	<=0.02%
Chloride (Cl)	<=5ppm
Heavy Metals (Pb)	<=5ppm
Insoluble matter	<=0.005%
Iron	<=5ppm
Manganese	<=5ppm
Nitrate	<=0.002%
pH of 0.5% solution (at 25°C)	5.0-8.2
Potassium	<=0.005%
Sodium	<=0.005%
Strontium	<=0.005%

Applications: Magnesium Sulfate is suitable for use in the preparation of buffer solutions.
Recommended Storage: RT

Manganese Chloride 1.0M Solution
Clear, Light-pink Solution

Packaging			Mfr. No
10 mL	PolyTube		BP541-1
100 mL	PolyBottle		BP541-100
Cl ₂ , Mn		H302, H373	
CAS: 7773-01-5		P260, P301+P312	
MW: 125.84			
Concentration			1.0M±0.01
DNase			Not detected
Heavy Metals (Pb)			<=5.0ppm
Protease			Not detected
RNase			Not detected

Applications: Manganese Chloride is used for dC and dG tailing of DNA.
0.2-micron filtered solution in 18-megohm water.
Recommended Storage: RT

JustPURE™ MES Hydrate Cellular and Molecular Biology White Powder

packaging		Mfr. No
250 g	Poly Bottle	BP2920-250
C ₆ H ₁₃ NO ₄ .xH ₂ O	H315, H319, H335	
CAS: 4432-31-9	P280, P305+P351+P338,	
MW: 195.23	P261, P302+P352, P264,	
EINECS: 224-632-3	P271	
Arsenic (As)	<=0.0001%	
Bismuth (Bi)	<=0.0005%	
Cadmium (Cd)	<=0.0005%	
Calcium (Ca)	<=0.002%	
Chloride	<=0.005%	
Chromium (Cr)	<=0.0005%	
Cobalt (Co)	<=0.0005%	
Copper (Cu)	<=0.0005%	
Iron (Fe)	<=0.0005%	
Lead (Pb)	<=0.0005%	
Lithium (Li)	<=0.0005%	
Magnesium (Mg)	<=0.0005%	
Purity (by titration)	>=99.5%	
Optical Absorbance at 260nm	<=0.025	
Optical Absorbance at 280nm	<=0.020	
pH at 20°C	Inclusive between 2.5-4.0	
Sodium (Na)	<=0.005%	
Strontium (Sr)	<=0.005%	

Applications: For high efficiency transfection of mammalian cells, gel electrophoresis of RNA, protein isolation applications, cell cultures and enzyme assays and bioanalytical methods, such as IEF, 2-D electrophoresis and SDS-PAGE.

Description: Ultra-pure zwitterionic buffers; Optimized for research in cellular and molecular biology; Minimal permeability to biological membranes; Low interference with biological reactions; pKa values mostly independent of temperature and concentration; High water solubility; Minimal absorption in spectral range 240 to 700nm.

Recommended storage: RT

Molecular Sieves, Grade 513 Type 4A; 4-8 Mesh

packaging		Mfr. No
500 g	PolyBottle	BP2634-500
CAS: 70955-01-0		
Mesh Size	To pass test	

Applications: Used for dehydration of gases.

Effective pore size:4Å...

Base:Alumina-Silicate

Cation:Sodium [1344-28-1 (Aluminum Oxide)] ; [14808-60-7 (Quartz)]

Recommended Storage: RT

Molecular sieves, Grade 514 (Type 4A; 8-12 Mesh) Type 4A; 8-12 Mesh

packaging		Mfr. No
500 g	PolyBottle	BP2631-500
CAS: 70955-01-0		
Mesh Size	To pass test	

Applications: Used for dehydration of liquids.

Effective pore size:4Å...

Base:Alumina-Silicate

Cation:Sodium [1344-28-1 (Aluminum Oxide)] ; [14808-60-7 (Quartz)]

Recommended Storage: RT

MOPS Fine White Crystals Molecular Biology

packaging		Mfr. No
100 g	PolyBottle	BP308-100
500 g	PolyBottle	BP308-500
C ₇ H ₁₅ NO ₄ S	H315, H335, H319	
CAS: 1132-61-2	P261, P302+P352, P280,	
MW: 209.26	P305+P351+P338	
EINECS: 214-478-5		
Assay	>=97%	
DNase	Not detected	
pKa (at 25°C)	7.2±0.2	
Protease	Not detected	
RNase	Not detected	

Applications: This biological buffer has a usable pH range of 6.5 to 7.9.

Recommended Storage: RT

CellPURE™ MOPS Sodium Salt Cell Culture-Tested White Crystalline powder

packaging		Mfr. No
25 g	Poly Bottle	BP2946-25
C ₇ H ₁₄ NNaO ₄ S		
CAS: 71119-22-7		
MW: 231.25		
Absorbance 33% w/w at 290nm	<=0.05	
Cell culture test	To pass test	
DNase	Not detected	
Endotoxin	<=0.1 EU/mg	
Heavy metals (as Pb)	<=5 ppm	
Infrared scan	Conforms to ref.	
Purity (by titration)	>=99.5%	
Nickase	Not detected	
Protease	Not detected	
RNase	Not detected	
Solubility	Clear and colorless	
Standard plate count	<=100CFU/gm	

Applications: Cell cultivation, isolation of cells, enzyme assays, and other biochemical applications

Description: Cell culture tested; Analyzed for the absence of nucleases and proteases; Tested for endotoxin and bioburden levels; pKa values mostly independent of temperature and concentration; High water solubility; Minimal permeability to biological membranes; Manufactured under strict quality control guidelines to ensure performance and reliability.

Recommended storage: RT

JustPURE™ MOPS Sodium Salt Cellular and Molecular Biology White Powder

packaging		Mfr. No
100 g	Poly Bottle	BP2926-100
C ₇ H ₁₄ NNaO ₄ S		
CAS: 71119-22-7		
MW: 231.25		
Ammonia	<=0.05%	
Calcium (Ca)	<=0.001%	
Chloride	<=0.05%	
Copper (Cu)	<=0.0005%	
Iron (Fe)	<=0.0005%	
Lead (Pb)	<=0.001%	
Magnesium (Mg)	<=0.0005%	
Purity (by titration)	>=99.5%	
Optical Absorbance at 260nm	<=0.1	
Optical Absorbance at 280nm	<=0.1	
pH at 20°C	Inclusive between 10.0-12.0	
Phosphorus	<=0.005%	
Solubility	Clear and colorless	
Sulfate (SO ₄)	<=0.05%	
Water (Karl Fischer)	<=5%	

Applications: For high efficiency transfection of mammalian cells, gel electrophoresis of RNA, protein isolation applications, cell cultures and enzyme assays and bioanalytical methods, such as IEF, 2-D electrophoresis and SDS-PAGE.

Description: Ultra-pure zwitterionic buffers; Optimized for research in cellular and molecular biology; Minimal permeability to biological membranes; Low interference with biological reactions; pKa values mostly independent of temperature and concentration; High water solubility; Minimal absorption in spectral range 240 to 700nm.

Recommended storage: RT

Heptane Sequencing Peroxide-free

packaging		Mfr. No
500 mL	AmberGlass,EcoSafPak*	BP1115-500
C ₇ H ₁₆	H400, H410	   
CAS: 142-82-5	P261, P301+P310, P331,	
MW: 100.2	P273, P302+P352,	
EINECS: 205-563-8	P210	
H225, H315, H336, H304,		
Aldehyde (HCHO)		Not detected
Assay		99.0%
Color (APHA)		<=10
Fluorescence Background (as Quinine Sulfate)		To pass test
Optical Absorbance at 200nm		<=0.75
Optical Absorbance at 215nm		<=0.20
Optical Absorbance at 254nm		<=0.014
Peroxides		Not detected
Refractive Index (at 25°C)		1.3845-1.3865
Residue after evaporation		<=5ppm
Water		<=0.02%

EcoSafPak is an environmentally friendly packaging system made of 100% recyclable material by an SFI certified supplier.*

Applications: Solvent for other organics and is used as a carrier in chromatography.

Recommended storage: RT

n-Octyl-β-D-Glucopyranoside White Crystalline Powder

packaging		Mfr. No
500 mg	PolyMicroTube	BP585-500
1 g	PolyMicroTube	BP585-1
5 g	AmberPolyBottle	BP585-5
25 g	AmberPolyBottle	BP585-25
C ₁₄ H ₂₈ O ₆	EINECS: 249-887-8	
CAS: 29836-26-8		
MW: 292.36		
Assay (GC)	>=98.0%	
Specific Rotation α _D ²⁰ (c=1, anhydrous MeOH)	-32.0° to -28.0°	

Applications: n-Octyl-β-D-Glucopyranoside is a highly purified, nonionic detergent.

Recommended Storage: -20°C

Not on TSCA inventory: for R and D use only; not for manufacturing or commercial purposes.

Phosphate Buffered Saline Tablets

packaging		Mfr. No
100 tablets g	Poly Bottle	BP2944-100
CAS: [Sodium chloride: 7647-14-5; Sodium phosphate dibasic: 7558-79-4; Dihydrogen potassium phosphate: 7778-77-0; Potassium chloride: 7447-40-7]		

Form Round tablet

pH at 25°C Inclusive between 7.2-7.6

Total weight Inclusive between 1814-2005mg/tab

Applications: Specifically developed for use in immunoassay procedures. The need to weigh and mix individual components is eliminated.

Description: One tablet dissolved in 200mL water yields 0.01M phosphate buffer, 0.0027M KCl, and 0.137M NaCl, pH 7.4 at 25°C

Recommended storage: RT

Phosphate Buffered Saline with Tween 20* White Powder

packaging		Mfr. No
10PK	Foil Pouches	BP2938-10
CAS: [Sodium chloride: 7647-14-5; Sodium phosphate dibasic: 7558-79-4; Tween 20: 9005-64-5; Di-hydrogen potassium phosphate: 7778-77-0; Potassium chloride: 7447-40-7]		
PBST		
pH at 25°C	Inclusive between 7.2-7.6	
Total weight	Inclusive between 9.55-10.55grams	

Applications: Used as a wash buffer and diluent for ELISA and Western blot.

Description: Dissolve the pre-weighed contents of one pouch in 1L of water.

Formulation per liter: 0.01M phosphate buffered saline, 0.138 M NaCl, 0.0027 M KCl, and 0.05% Tween 20, pH 7.4 at 25°C

Recommended storage: RT

Sodium Carbonate White Crystalline Powder		
packaging	Mfr. No	
1 kg PolyBottle	BP357-1	
CNa ₂ O ₃	H319	
CAS: 497-19-8	P280, P305+P351+P338	
MW: 105.99		
EINECS: 207-838-8		
Assay	>=99.5%	
Calcium (Ca)	<=0.03%	
Chloride (Cl)	<=0.001%	
Heavy Metals (Pb)	0.0005%	
Insoluble matter	<=0.01%	
Iron	0.0005%	
Loss on Heating (at 285°C)	<=1.0%	
Magnesium (Mg)	<=0.005%	
Phosphate (PO ₄)	<=0.001%	
Potassium	0.005%	
Silica	<=0.005%	
Sulfur Compounds (as SO ₄)	<=0.003%	

Applications: This sodium salt is used in buffers for molecular biology and cell culture applications.
Recommended Storage: RT

Sodium Chloride		
packaging	Mfr. No	
1 kg PolyBottle	BP358-1	
2.5 kg PolyBottle	BP358-212	
10 kg PolyPail	BP358-10	
ClNa	EINECS: 231-598-3	
CAS: 7647-14-5		
MW: 58.44		
Calcium (Ca)	<=0.005%	
Heavy Metals (Pb)	0.0005%	
Iodide	<=0.002%	
Iron	<=0.0002%	
Magnesium (Mg)	<=0.005%	
Phosphate (PO ₄)	0.0005%	
Sulfate (SO ₄)	<=0.004%	

Applications: Sodium Chloride is a main component of many electrophoresis and molecular biology buffers and media.
Recommended Storage: RT

Sodium Citrate Dihydrate Small, Colorless Granules		
packaging	Mfr. No	
500 g PolyBottle	BP327-500	
1 kg PolyBottle	BP327-1	
C ₆ H ₅ Na ₃ O ₇ ·2H ₂ O		
CAS: 6132-04-3		
MW: 294.09		
Ammonia	<=0.003%	
Assay	99.0-101.0%	
Calcium (Ca)	<=0.005%	
Chloride (Cl)	<=0.003%	
Free Acid (as Citric Acid)	<=0.15%	
Free Alkali	None	
Heavy Metals (Pb)	<=5ppm	
Insoluble matter	<=0.005%	
Iron	<=0.001%	
pH of 5% Solution (at 25°C)	7.0-9.0	
Sulfate (SO ₄)	<=0.005%	

Applications: Sodium Citrate is used primarily in the preparation of SSC buffer used in Southern transfers.
Recommended Storage: RT

Sodium Hydroxide White Pellets		
packaging	Mfr. No	
500 g PolyBottle	BP359-500	
2.5 kg PolyBottle	BP359-212	
HNaO	H314	
CAS: 1310-73-2	P280, P301+P330+P331,	
MW: 40	P305+P351+P338, P302+P352	
EINECS: 215-185-5		
Assay	>=97%	
Calcium (Ca)	<=0.005%	
Chloride (Cl)	<=0.005%	
Copper	0.001%	
Heavy Metals (as Ag)	0.002%	
Iron	<=0.001%	
Magnesium (Mg)	0.002%	
Mercury	<=0.1ppm	
Nickel	<=0.001%	
Nitrogen Compounds (as N)	<=0.001%	
Phosphate (PO ₄)	0.001%	
Potassium	<=0.02%	
Sodium carbonate	<=1.0%	
Sulfate (SO ₄)	<=0.003%	

Applications: Sodium Hydroxide is used to adjust the pH of buffers for molecular biology and cell culture applications.
Recommended Storage: RT
UN 1823; DOT Class 8:Corrosive

Sodium Iodide White Crystals		Molecular Biology
packaging	Mfr. No	
100 g AmberGlass	BP323-100	
INa	EINECS: 231-679-3	
CAS: 7681-82-5		
MW: 149.89		
Assay	>=99%	
Chloride and Bromide (as Cl)	<=0.01%	
Copper	<=3ppm	
DNase	Not detected	
Heavy Metals (Pb)	<=5ppm	
Insoluble matter	<=0.005%	
Iodate	<=0.003%	
Iron	<=3ppm	
Phosphate (PO ₄)	<=0.001%	
RNase	Not detected	
Sulfate (SO ₄)	<=0.005%	

Applications: This sodium salt is used in buffers for molecular biology and cell culture applications.
Recommended Storage: RT

Sodium Nitrate Colorless Crystals		
packaging	Mfr. No	
500 g PolyBottle	BP360-500	
NNaO ₃	P261, P301+P312,	
CAS: 7631-99-4	P302+P352, P280,	
MW: 84.99	P305+P351+P338, P210	
EINECS: 231-554-3		
H272, H315, H335, H302, H319		
Assay	>=99%	
Calcium (Ca)	<=0.005%	
Chloride (Cl)	<=0.001%	
Heavy Metals (Pb)	<=5ppm	
Insoluble matter	<=0.005%	
Iodate	<=5ppm	
Iron	<=3ppm	
Magnesium (Mg)	<=0.002%	
Nitrite	<=0.001%	
pH of 5% Solution (at 25°C)	5.5-8.3	
Phosphate (PO ₄)	5ppm	
Sulfate (SO ₄)	<=0.003%	

Applications: This sodium salt is used in buffers for molecular biology and cell culture applications.
UN 1498; DOT Class 5.1:Oxidizer

Sodium Phosphate Dibasic Anhydrous White Granules or Powder		
packaging	Mfr. No	
500 g PolyBottle	BP332-500	
1 kg PolyBottle	BP332-1	
HNa ₂ O ₄ P	EINECS: 231-448-7	
CAS: 7558-79-4		
MW: 141.96		
Assay	>=99.0%	
Chloride (Cl)	<=0.002%	
Heavy Metals (Pb)	<=0.001%	
Insoluble matter	<=0.010%	
Iron	<=0.001%	
Loss on Drying (at 105°C)	<=0.20%	
pH of 5% Solution (at 25°C)	8.7-9.3	
Sulfate (SO ₄)	<=0.005%	

Applications: This buffer component is commonly used in tissue culture and molecular biology applications.
Recommended Storage: RT

Sodium Phosphate Dibasic Heptahydrate Colorless-to-white Crystals		
packaging	Mfr. No	
500 g AmberGlass	BP331-500	
1 kg AmberGlass	BP331-1	
HNa ₂ O ₄ P·7H ₂ O		
CAS: 7782-85-6		
MW: 268.07		
Assay	98.0-102.0%	
Chloride (Cl)	<=0.001%	
Heavy Metals (Pb)	<=0.001%	
Insoluble matter	<=0.005%	
Iron	<=0.001%	
pH of 5% Solution (at 25°C)	8.7-9.3	
Sulfate (SO ₄)	<=0.005%	

Applications: This buffer component is commonly used in tissue culture and molecular biology applications.
Recommended Storage: RT

Sodium Phosphate Monobasic Anhydrous Colorless-to-white Crystals		
packaging	Mfr. No	
500 g AmberGlass	BP329-500	
1 kg AmberGlass	BP329-1	
H ₂ NaPO ₄	EINECS: 231-449-2	
CAS: 7558-80-7		
MW: 119.98		
Assay	>=99%	
Heavy Metals (Pb)	<=0.001%	
Insoluble matter	<=0.03%	
Moisture	0.5%	
pH of 1M Solution (at 25°C)	4.0-6.0	

Applications: Sodium Phosphate is commonly used in tissue culture and molecular biology applications.
Recommended Storage: RT

Sodium Phosphate Monobasic Monohydrate Colorless-to-white Crystals		
packaging	Mfr. No	
500 g AmberGlass	BP330-500	
1 kg AmberGlass	BP330-1	
H ₂ NaO ₄ P·H ₂ O		
CAS: 10049-21-5		
MW: 137.99		
Assay	98.0-102.0%	
Calcium (Ca)	<=0.005%	
Chloride (Cl)	<=5ppm	
Heavy Metals (Pb)	<=0.001%	
Insoluble matter	<=0.01%	
Iron	<=0.001%	
pH of 5% Solution (at 25°C)	4.1-4.5	
Potassium	0.01%	
Sulfate (SO ₄)	<=0.003%	

Applications: Sodium Phosphate is commonly used in tissue culture and molecular biology applications.
Recommended Storage: RT

Sodium Tetraborate Decahydrate White Crystals		Molecular Biology
packaging	Mfr. No	
500 g AmberGlass	BP175-500	
B ₄ Na ₂ O ₇ ·10H ₂ O	H360FD	
CAS: 1303-96-4	P201, P308+P313	
MW: 381.36		
Assay	99.5-105.0%	
Calcium (Ca)	<=0.005%	
Chloride (Cl)	<=0.001%	
DNase	Not detected	
Heavy Metals (Pb)	<=0.001%	
Insoluble matter	<=0.005%	
Iron	<=5ppm	
pH of 0.01M Solution (at 25°C)	9.15-9.20	
Phosphate (PO ₄)	<=0.001%	
RNase	Not detected	
Sulfate (SO ₄)	<=0.005%	

Applications: Sodium Borate can be used in the preparation of TBE electrophoresis buffer.
Recommended Storage: RT

STET Buffer, 1X Solution pH 8.0			
packaging			Mfr. No.
100 mℓ	PolyBottle		BP2480-100
1 ℓ	PolyBottle		BP2480-1
DNase			Not detected
pH (at 25°C)			7.9-8.1
Protease			Not detected
RNase			Not detected

Applications: STET buffer is used in biomolecular procedures. Saline/Tris/EDTA/Triton* X-100 Solution 1X solution contains 10mM Tris-HCl, 1mM EDTA, 100mM NaCl, and 5% Triton* X-100 (v/v). Autoclaved.
[77-86-1 (Tris Base)] ; [60-00-4 (EDTA)] ; [7647-01-0 (HCl)] ; [7647-14-5 (NaCl)] ; [9002-93-1 (Triton* X-100)]
Recommended Storage: RT
Filtered through a 5-micron filter.

Sudan Black	
packaging	Mfr. No
10 g AmberGlass	BP109-10
C ₂₉ H ₂₄ N ₆ CAS: 4197-25-5 MW: 456.54	EINECS: 224-087-1
E ¹ % _{1cm}	>=600ℓ g-1cm-1
Lambda Max. in Acetone	600nm ±2nm
Loss on Drying (at 105°C)	<=5%
Solubility	To pass test
Applications: Used to stain chromosomes and lipoproteins, as well as fat in bacteria.	
Recommended storage: RT	
C.I. 26150	

TAPS (Free Acid) White Powder		
packaging		Mfr. No
25 g	PolyBottle	BP551-25
100 g	PolyBottle	BP551-100
500 g	PolyBottle	BP551-500
C ₇ H ₁₇ NO ₆ S CAS: 29915-38-6 MW: 243.27		EINECS: 249-954-1
Assay (purity)		>=99%
DNase		Not detected
Heavy Metals (Pb)		<=5ppm
Moisture		1%
Protease		Not detected
RNase		Not detected
Applications: This biological buffer has a usable pH range of 7.7 to 9.1.		
Recommended Storage: RT		

TES White Crystals				
Packaging			Mfr. No	
100 g	AmberGlass		BP309-100	
C ₆ H ₁₅ NO ₆ S		H335, H315, H319		
CAS: 7365-44-8		P261, P302+P352, P280,		
MW: 229.25		P305+P351+P338		
EINECS: 230-906-3				
Assay				>=98%
Melting Point				219° ±5°C
PKa (at 25°C)				7.4±0.2

Applications: This biological buffer has a usable pH range of 6.8 to 8.2.
Recommended Storage: RT

Trichloroacetic Acid Flakes or Crystals			
Packaging		Mfr. No	
250 g	AmberGlass	BP555-250	
500 g	AmberGlass	BP555-500	
1 kg	PolyBottle	BP555-1	
C ₂ HCl ₃ O ₂		P301+P330+P331, P280,	
CAS: 76-03-9		P305+P351+P338, P310,	
MW: 163.39		P261, P273	
EINECS: 200-927-2			
H314, H335, H410			
Assay		>=99.0%	
Chloride (Cl)		<=0.002%	
Heavy Metals (Pb)		<=0.002%	
Iron		<=0.001%	
Nitrate		<=0.002%	
Phosphate (PO ₄)		<=5ppm	
Residue after ignition		0.03%	
Substances Darkened by Sulfuric Acid		To pass test	
Sulfate (SO ₄)		<=0.02%	
Applications: TCA is used to precipitate proteins and nucleic acids from solution.			
Recommended Storage: RT			
UN 1839; DOT Class 8:Corrosive			

Tricine White Crystals	
packaging	
100 g	PolyBottle
C ₆ H ₁₃ NO ₃	H315, H319, H335
CAS: 5704-04-1	P261, P302+P352, P280,
MW: 179.17	P305+P351+P338
EINECS: 227-193-6	
Assay	98%
Assay	>=98%
Melting Point	184°±2°C
PKa (at 25°C)	8.10±0.2
Applications: This biological buffer has a usable pH range of 7.4 to 8.8.	
Recommended Storage: RT	

Tris Base		
White Crystals or Crystalline Powder		
packaging		Mfr. No
1 kg	PolyBottle	BP154-1
C ₄ H ₁₁ NO ₃	H335, H315, H319	
CAS: 77-86-1	P280, P305+P351+P338,	
MW: 121.14	P261, P302+P352	
EINECS: 201-064-4		
Absorbance of a 2M Solution at 280nm		<=0.05
Arsenic.....		<=5ppm
Assay (corrected for H2O)		>=99.9%
Copper		<=1ppm
DNase		Not detected
Electrophoresis		To pass test
Iron		<=2ppm
Lead		<=1ppm
Magnesium (Mg)		<=10ppm
pH of 0.05M Solution (at 25°C)		10.3-10.5
Protease		Not detected
RNase		Not detected
Water		<=0.5%

Applications: This ultrapure, triple-crystallized buffer component is alcohol-free, so it is compatible with sensitive molecular biology, tissue culture, and electrophoresis methods.
Recommended Storage: RT

Tris Base			Molecular Biology
White Crystals or Crystalline Powder			
packaging			Mfr. No
500 g	PolyBottle		BP152-500
1 kg	PolyBottle		BP152-1
5 kg	PolyPail		BP152-5
10 kg	PolyPail		BP152-10
25 kg	FiberDrum		BP152-25
C ₄ H ₁₁ NO ₃		H335, H315, H319	
CAS: 77-86-1		P280, P305+P351+P338,	
MW: 121.14		P261, P302+P352	
EINECS: 201-064-4			
Absorbance of a 1M Solution at 280nm			<=0.05
Arsenic			<=5ppm
Assay (corrected for H2O)			99.8%
Copper			<=1ppm
DNase			Not detected
Electrophoresis			To pass test
Iron			<=2ppm
Lead			1ppm
Magnesium (Mg)			<=10ppm
pH of 0.05M Solution (at 25°C)			10.3-10.5
Protease			Not detected
RNase			Not detected
Water			<=2%

Applications: Tris Base is a buffer component in molecular biology, tissue culture, and electrophoresis procedures.
Recommended Storage: RT

Tris Buffered Saline 10X Solution		Molecular Biology
Packaging		Mfr. No.
100 mL	PolyBottle	BP2471-100
500 mL	PolyBottle	BP2471-500
1 L	PolyBottle	BP2471-1
H315, H319		
P305+P351+P338, P280, P264, P302+P352		
		
DNase	Not detected	
pH (1X solution) (at 25°C)	7.3-7.5	
Protease	Not detected	
RNase	Not detected	

Applications: Tris Buffered Saline (TBS) is used in biochemical and molecular biology applications. Filtered and autoclaved.
Components: 10X solution contains 1.37M Sodium Chloride, 0.027M Potassium Chloride, and 0.25M Tris/Tris-HCl.
pH:pH 7.4 [77-86-1 (Tris)] ; [1185-53-1 (Tris Hydrochloride)] ; [7447-40-7 (Potassium Chloride)] ; [7647-14-5 (Sodium Chloride)]
Recommended Storage: RT

Tris Buffered Saline 1X Solution, pH 7.4		Molecular Biology
Packaging		Mfr. No.
100 mL	PolyBottle	BP2472-100
1 L	PolyBottle	BP2472-1
H315, H319		
P305+P351+P338, P280, P264, P302+P352		
DNase	Not detected	
pH (at 25°C)	7.3-7.5	
Protease	Not detected	
RNase	Not detected	

Applications: Tris Buffered Saline (TBS) is used in biochemical and molecular biology applications. 1X solution contains 0.137M Sodium Chloride, 0.0027M Potassium Chloride and 0.025M Tris/Tris-HCl. Filtered and autoclaved.
[77-86-1 (Tris)] ; [1185-53-1 (Tris Hydrochloride)] ; [7447-40-7 (Potassium Chloride)] ; [7647-14-5 (Sodium Chloride)]
Recommended Storage: RT

Tris Hydrochloride Small White Flakes		Molecular Biology
packaging		Mfr. No
500 g	PolyBottle	BP153-500
1 kg	PolyBottle	BP153-1
C ₄ H ₁₁ NO ₃ .HCl CAS: 1185-53-1 MW: 157.6 EINECS: 214-684-5	H335, H319, H315 P261, P302+P352, P280, P305+P351+P338	
Absorbance of a 1M Solution at 280nm	<=0.05	
Arsenic	<=10ppm	
Assay	>=99%	
Barium	<=1ppm	
Calcium (Ca)	<=10ppm	
Copper	<=1ppm	
DNase	Not detected	
Electrophoresis	To pass test	
Iron	<=2ppm	
Lead	<=5ppm	
Magnesium (Mg)	10ppm	
Manganese	<=1ppm	
Melting Point	52°±2°C	
Protease	Not detected	
RNase	Not detected	

Applications: Tris Hydrochloride is a buffering medium for electrophoresis, molecular biology, and cell culture applications. It may be combined with Tris Base to simplify Tris buffer preparation.
Recommended Storage: RT

Tris-EDTA (TE), 10X Solution pH 7.4		Molecular Biology
packaging		Mfr. No
100 mℓ	PolyBottle	BP2477-100
500 mℓ	PolyBottle	BP2477-500
1 ℓ	PolyBottle	BP2477-1
CAS: 38641-82-6		
DNase		Not detected
pH (1X solution) (at 25°C)		7.3-7.5
Protease		Not detected
RNase		Not detected

Applications: Tris-EDTA (TE) is used for suspending nucleic acid samples. 10X solution contains 100mM Tris and 10mM EDTA. [60-00-4 (EDTA)] ; [77-86-1 (Tris)]
Recommended Storage: RT
Filtered through a 5-micron filter and autoclaved.

Tris-EDTA (TE), 10X Solution pH 7.6		Molecular Biology
packaging		Mfr. No
100 mℓ	PolyBottle	BP2475-100
500 mℓ	PolyBottle	BP2475-500
1 ℓ	PolyBottle	BP2475-1
CAS: 38641-82-6		
DNase		Not detected
pH (1X solution) (at 25°C)		.5-7.7
Protease		Not detected
RNase		Not detected

Applications: Tris-EDTA (TE) is used for suspending nucleic acid samples. 10X solution contains 100mM Tris and 10mM EDTA. [60-00-4 (EDTA)] ; [77-86-1 (Tris)]
Recommended Storage: RT
Filtered through a 5-micron filter and autoclaved.

Tris-EDTA 1X Solution, pH 7.4		Molecular Biology
packaging		Mfr. No
100 mℓ	PolyBottle	BP2476-100
500 mℓ	mℓPolyBottle	BP2476-500
CAS: 38641-82-6		
DNase		Not detected
pH (1X solution) (at 25°C)		7.4±0.1
Protease		Not detected
RNase		Not detected

Applications: Tris-EDTA (TE) is routinely used for suspending nucleic acid samples. 10mM Tris and 1mM EDTA [60-00-4 (EDTA)] ; [77-86-1 (Tris)] ; [7647-01-0 (Hydrogen Chloride)]
Recommended Storage: RT
Filtered through a 5-micron filter.

Tris-EDTA 1X Solution, pH 7.6		Molecular Biology
packaging		Mfr. No
100 mℓ	PolyBottle	BP2474-100
500 mℓ	PolyBottle	BP2474-500
1 ℓ	PolyBottle	BP2474-1
CAS: 38641-82-6		
DNase		Not detected
pH (1X solution) (at 25°C)		7.6±0.1
Protease		Not detected
RNase		Not detected

Applications: Tris-EDTA (TE) is routinely used for suspending nucleic acid samples. 10mM Tris and 1mM EDTA [60-00-4 (EDTA)] ; [77-86-1 (Tris)] ; [7647-01-0 (Hydrogen Chloride)]
Recommended Storage: RT
Filtered through a 5-micron filter.

Tris-EDTA 1X Solution, pH 8.0		Molecular Biology
packaging		Mfr. No
100 mℓ	PolyBottle	BP2473-100
500 mℓ	PolyBottle	BP2473-500
1 ℓ	PolyBottle	BP2473-1
CAS: 38641-82-6		
DNase		Not detected
pH (1X solution) (at 25°C)		8.0±0.1
Protease		Not detected
RNase		Not detected

Applications: Tris-EDTA (TE) is routinely used for suspending nucleic acid samples. 10mM Tris and 1mM EDTA [60-00-4 (EDTA)] ; [77-86-1 (Tris)] ; [7647-01-0 (Hydrogen Chloride)]
Recommended Storage: RT
Filtered through a 5-micron filter.

Water Sterile-Filtered		Molecular Biology
packaging		Mfr. No
100 mℓ	Poly Bottle	BP2819-100
1 ℓ	Poly Bottle	BP2819-1
4 ℓ	PolyPac*	BP2819-4
10 ℓ	PolyPac*	BP2819-10
20 ℓ	PolyPac*	BP2819-20
H ₂ O	EINECS: 231-791-2	
CAS: 7732-18-5		
MW: 18.02		
Cadmium (Cd)		<=10 ppb
Calcium (Ca)		<=10 ppb
Chromium (Cr)		<=10 ppb
Cobalt (Co)		<=10 ppb
Copper (Cu)		<=10 ppb
DNase		Not detected
Iron (Fe)		<=10 ppb
Lead (Pb)		<=10 ppb
Magnesium (Mg)		<=10 ppb
Manganese (Mn)		<=10 ppb
pH at 25°C		Inclusive between 5.4-7.0
Protease		Not detected
RNase		Not detected
Solubility in water		<2uS/cm
Vanadium (V)		<=10 ppb

Applications: Ideal for many fundamental procedures such as PCR, electrophoresis, DNA sequencing, and buffers for enzymatic analyses.
Description: 0.3µm filtered to ensure high purity
Recommended storage: RT

Water, Sterile		Biotech Grade
packaging		Mfr. No
20 ℓ	PolyPac*	BP2485-20
4 ℓ	PolyPac*	BP2485-4
H ₂ O	EINECS: 231-791-2	
CAS: 7732-18-5		
MW: 18.02		
Heavy Metals (Pb)		<=0.01ppm
pH (at 25°C)		5.4-7.0
Specific Conductance (at 25°C)		<=2°-10-6 ohm-1 cm-1
Substances Reducing Permanganate		To pass test

Applications: Biotech-grade water can be used for several biological test procedures.

Water, Sterile For DNA Work		DNA Grade
packaging		Mfr. No
1 ℓ	PolyBottle	BP2470-1
H ₂ O	EINECS: 231-791-2	
CAS: 7732-18-5		
MW: 18.02		
DNase		Not detected
Protease		Not detected

Applications: DNA-grade Water is DNase- and Protease-free and is suitable for all DNA work.
Filtered through a 0.2-micron filter and autoclaved.

Water Sterile-Filtered		Microbial Cell Culture-Grade
packaging		Mfr. No
100 mℓ	Poly Bottle	BP2820-100
500 mℓ	Poly Bottle	BP2820-500
1 ℓ	Poly Bottle	BP2820-1
H ₂ O	EINECS: 231-791-2	
CAS: 7732-18-5		
MW: 18.02		
Arsenic (As)		<=10 ppb
Cadmium (Cd)		<=10 ppb
Calcium (Ca)		<=20 ppb
Chromium (Cr)		<=10 ppb
Cobalt (Co)		<=10 ppb
Copper (Cu)		<=10 ppb
Endotoxin		<0.025 EU/mℓ
Iron (Fe)		<=10 ppb
Lead (Pb)		<=10 ppb
Magnesium (Mg)		<=10 ppb
Manganese (Mn)		<=10 ppb
Oxidizable substances		To pass test
pH at 25°C		Inclusive between 5.4-7.0
Silver (Ag)		<=10 ppb
Sodium (Na)		<=20 ppb
Solubility in water		<2uS/cm
Tin (Sn)		<=10 ppb
Vanadium (V)		<=10 ppb

Applications: For the preparation of microbiological growth media (functionally tested for promoting microbial growth).
Description: 0.3µm filtered and autoclaved
Recommended storage: RT

Water, Sterile Nuclease Free		
packaging		Mfr. No
50 mℓ	PolyBottle	BP2484-50
100 mℓ	PolyBottle	BP2484-100
H ₂ O	EINECS: 231-791-2	
CAS: 7732-18-5		
MW: 18.02		
DNase		Not detected
Protease		Not detected
RNase		Not detected

Applications: Diethylpyrocarbonate (DEPC)-treated Water is guaranteed nuclease-free (tested for both DNase and RNase) and is suitable for all RNA work.
Filtered through a 0.2-micron filter.

Water, Sterile For RNA Work, DEPC-treated and Nuclease-free		Molecular Biology
packaging		Mfr. No
1 ℓ	PolyBottle	BP561-1
H ₂ O	EINECS: 231-791-2	
CAS: 7732-18-5		
MW: 18.02		
DNase		Not detected
Protease		Not detected
RNase		.Not detected
Specific Conductance (at 25°C)		.<=18°-10-6ohm-1cm-1

Applications: Diethylpyrocarbonate (DEPC)-treated Water is guaranteed nuclease-free (tested for both DNase and RNase) and is suitable for all RNA work.
DEPC-treated and autoclaved.
Filtered through a 0.2-micron filter.

Hazard Class and category

Asp. Tox. 1	Aspiration hazard; Category 1
Acute Tox. 1	Acute toxicity (oral); Category 1
Acute Tox. 2	Acute toxicity (oral); Category 2
Acute Tox. 3	Acute toxicity (oral); Category 3
Acute Tox. 4	Acute toxicity (oral); Category 4
Acute Tox. 1	Acute toxicity (dermal); Category 1
Acute Tox. 2	Acute toxicity (dermal); Category 2
Acute Tox. 3	Acute toxicity (dermal); Category 3
Acute Tox. 4	Acute toxicity (dermal); Category 4
Acute Tox. 1	Acute toxicity (inhalation); Category 1
Acute Tox. 2	Acute toxicity (inhalation); Category 2
Acute Tox. 3	Acute toxicity (inhalation); Category 3
Acute Tox. 4	Acute toxicity (inhalation); Category 4
Skin Corr. 1A	Skin corrosion/irritation; Category 1A
Skin Corr. 1B	Skin corrosion/irritation; Category 1B
Skin Corr. 1C	Skin corrosion/irritation; Category 1C
Skin Irrit. 2	Skin corrosion/irritation; Category 2
Eye Dam. 1	Serious eye damage/eye irritation; Category 1
Eye Irrit. 2	Serious eye damage/eye irritation; Category 2
Resp. Sens. 1	Respiratory sensitisation; Category 1
Skin Sens. 1	Skin sensitisation; Category 1
Muta. 1A	Germ cell mutagenicity; Category 1A
Muta. 1B	Germ cell mutagenicity; Category 1B
Muta. 2	Germ cell mutagenicity; Category 2
Carc. 1A	Carinogenicity; Category 1A
Carc. 1B	Carcinogenicity; Category 1B
Carc. 2	Carcinogenicity; Category 2
Repr. 1A	Reproductive toxicity; Category 1A
Repr. 1B	Reproductive toxicity; Category 1B
Repr. 2	Reproductive toxicity; Category 2
STOT SE 1	Specific target organ toxicity - single exposure; Category 1
STOT SE 2	Specific target organ toxicity - single exposure; Category 2
STOT SE 3	Specific target organ toxicity - single exposure; Category 3
STOT RE 1	Specific target organ toxicity - repeated exposure; Category 1

STOT RE 2	Specific target organ toxicity - repeated exposure; Category 2
Aquatic Acute 1	Hazardous to the aquatic environment - acute hazard; Category 1
Aquatic Chronic 1	Hazardous to the aquatic environment - chronic hazard; Category 1
Aquatic Chronic 2	Hazardous to the aquatic environment - chronic hazard; Category 2
Aquatic Chronic 3	Hazardous to the aquatic environment - chronic hazard; Category 3
Aquatic Chronic 4	Hazardous to the aquatic environment - chronic hazard; Category 4
Ozone	Hazardous to the ozone layer; Category 1
Flam. Liq. 1	Flammable liquid; Category 1
Flam. Liq. 2	Flammable liquid; Category 2
Flam. Liq. 3	Flammable liquid; Category 3
Flam. Sol. 1	Flammable solid; Category 1
Flam. Sol. 2	Flammable solid; Category 2
Org. Perox. C	Organic peroxide; Type C
Org. Perox. D	Organic peroxide; Type D
Ox. Liq. 1	Oxidising liquid; Category 1
Ox. Liq. 2	Oxidising liquid; Category 2
Ox. Liq. 3	Oxidising liquid; Category 3
Ox. Sol. 1	Oxidising solid; Category 1
Ox. Sol. 2	Oxidising solid; Category 2
Ox. Sol. 3	Oxidising solid; Category 3
Pyr. Sol. 1	Pyrophoric solid; Category 1
Self-heat. 1	Self-heating substance/mixture; Category 1
Self-react. C	Self-reactive substance; Type C
Met. Corr. 1	Substance/mixture corrosive to metal; Category 1
Water-react. 1	Substance/mixture which, in contact with water, emits flammable gases; Category 1
Water-react. 2	Substance/mixture which, in contact with water, emits flammable gases; Category 2
Water-react. 3	Substance/mixture which, in contact with water, emits flammable gases; Category 3

Hazard and precautionary statements

EUH Phrases	
EUH001	Explosive when dry
EUH006	Explosive with or without contact with air
EUH014	Reacts violently with water
EUH018	In use, may form flammable/explosive vapor-air mixture
EUH019	May form explosive peroxides
EUH029	Contact with water liberates toxic gas
EUH031	Contact with acids liberates toxic gas
EUH032	Contact with acids liberates very toxic gas
EUH044	Risk of explosion if heated under confinement
EUH059	Hazardous to the ozone layer
EUH066	Repeated exposure may cause skin dryness or cracking
EUH070	Toxic by eye contact
EUH071	Corrosive to the respiratory tract
EUH201	Contains lead. Should not be used on surfaces liable to be chewed or sucked by children
EUH201A	Warning! Contains lead
EUH202	Cyanoacrylate. Danger. Bonds skin and eyes in seconds. Keep out of the reach of children
EUH203	Contains chromium (VI). May produce an allergic reaction
EUH204	Contains isocyanates. May produce an allergic reaction
EUH205	Contains epoxy constituents. May produce an allergic reaction
EUH206	Warning! Do not use together with other products. May release dangerous gases (chlorine)
EUH207	Warning! Contains cadmium. Dangerous fumes are formed during use. See information supplied by the manufacturer. Comply with the safety instructions
EUH208	May produce an allergic reaction
EUH209	Can become highly flammable in use
EUH210	Safety data sheet available on request
EUH401	To avoid risks to human health and the environment, comply with the instructions for use

H Phrases	
H200	Unstable explosive
H201	Explosive; mass explosion hazard
H202	Explosive; severe projection hazard
H203	Explosive; fire, blast or projection hazard
H204	Fire or projection hazard
H205	May mass explode in fire
H220	Extremely flammable gas
H221	Flammable gas
H222	Extremely flammable aerosol
H223	Flammable aerosol
H224	Extremely flammable liquid and vapor
H225	Highly flammable liquid and vapor
H226	Flammable liquid and vapor
H227	Combustible liquid
H228	Flammable solid
H240	Heating may cause an explosion
H241	Heating may cause a fire or explosion
H242	Heating may cause a fire
H250	Catches fire spontaneously if exposed to air
H251	Self-heating: may catch fire
H252	Self-heating in large quantities; may catch fire
H260	In contact with water releases flammable gases which may ignite spontaneously
H261	In contact with water releases flammable gases
H270	May cause or intensify fire; oxidizer
H271	May cause fire or explosion; strong oxidizer
H272	May intensify fire; oxidizer
H280	Contains gas under pressure; may explode if heated
H281	Contains refrigerated gas; may cause cryogenic burns or injury
H290	May be corrosive to metals
H300	Fatal if swallowed
H301	Toxic if swallowed
H302	Harmful if swallowed
H303	May be harmful if swallowed
H304	May be fatal if swallowed and enters airways
H305	May be harmful if swallowed and enters airways
H310	Fatal in contact with skin
H311	Toxic in contact with skin
H312	Harmful in contact with skin
H313	May be harmful in contact with skin
H314	Causes severe skin burns and eye damage
H315	Causes skin irritation
H316	Causes mild skin irritation
H317	May cause an allergic skin reaction
H318	Causes serious eye damage
H319	Causes serious eye irritation
H320	Causes eye irritation
H330	Fatal if inhaled
H331	Toxic if inhaled

H332	Harmful if inhaled
H333	May be harmful if inhaled
H334	May cause allergy or asthma symptoms or breathing difficulties if inhaled
H335	May cause respiratory irritation
H336	May cause drowsiness or dizziness
H340	May cause genetic defects
H341	Suspected of causing genetic defects
H350	May cause cancer
H350i	May cause cancer by inhalation
H351	Suspected of causing cancer
H360	May damage fertility or the unborn child
H360D	May damage the unborn child
H360Di	May damage the unborn child. Suspected of damaging fertility
H360F	May damage fertility
H360Fd	May damage fertility. Suspected of damaging the unborn child
H361	Suspected of damaging fertility or the unborn child
H361d	Suspected of damaging the unborn child
H361f	Suspected of damaging fertility
H361fd	Suspected of damaging fertility. Suspected of damaging the unborn child
H362	May cause harm to breast-fed children
H370	Causes damage to organs
H371	May cause damage to organs
H372	Causes damage to organs through prolonged or repeated exposure
H373	May cause damage to organs through prolonged or repeated exposure
H400	Very toxic to aquatic life
H401	Toxic to aquatic life
H402	Harmful to aquatic life
H410	Very toxic to aquatic life with long lasting effects
H411	Toxic to aquatic life with long lasting effects
H412	Harmful to aquatic life with long lasting effects
H413	May cause long lasting harmful effects to aquatic life
H420	Harms public health and the environment by destroying ozone in the upper atmosphere

P Phrases	
P101	If medical advice is needed, have product container or label at hand
P102	Keep out of reach of children
P103	Read label before use
P201	Obtain special instructions before use
P202	Do not handle until all safety precautions have been read and understood
P210	Keep away from heat/sparks/open flames/hot surfaces. - No smoking
P211	Do not spray on an open flame or other ignition source
P220	Keep/Store away from clothing/ combustible materials
P221	Take any precaution to avoid mixing with combustibles
P222	Do not allow contact with air
P223	Keep away from any possible contact with water, because of violent reaction and possible flash fire
P230	Keep wetted with .?
P231	Handle under inert gas
P231 + P232	Handle under inert gas. Protect from moisture
P232	Protect from moisture
P233	Keep container tightly closed
P234	Keep only in original container
P235	Keep cool
P235+ P410	Keep cool. Protect from sunlight
P240	Ground/Bond container and receiving equipment
P241	Use explosion-proof electrical/ ventilating/ lighting/ equipment
P242	Use only non-sparking tools
P243	Take precautionary measures against static discharge
P244	Keep reduction valves free from grease and oil
P250	Do not subject to grinding/shock/friction
P251	Pressurized container: Do not pierce or burn, even after use
P260	Do not breathe dust/fume/gas/mist/vapors/spray
P261	Avoid breathing dust/ fume/ gas/ mist/ vapors/ spray
P262	Do not get in eyes, on skin, or on clothing
P263	Avoid contact during pregnancy/ while nursing
P264	Wash face, hands and any exposed skin thoroughly after handling
P270	Do not eat, drink or smoke when using this product
P271	Use only outdoors or in a well-ventilated area
P272	Contaminated work clothing should not be allowed out of the workplace
P273	Avoid release to the environment
P280	Wear protective gloves/ protective clothing/ eye protection/ face protection
P281	Use personal protective equipment as required
P282	Wear cold insulating gloves/ face shield/ eye protection
P283	Wear fire/ flame resistant/ retardant clothing
P284	Wear respiratory protection
P285	In case of inadequate ventilation wear respiratory protection

Hazard and precautionary statements

P Phrases	
P301	IF SWALLOWED
P301 + P310	IF SWALLOWED: Immediately call a POISON CENTER or doctor/physician
P301 + P312	IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell
P301 + P330	IF SWALLOWED: Rinse mouth.
+ P331	Do NOT induce vomiting
P302	IF ON SKIN
P302 + P334	IF ON SKIN: Immerse in cool water/ wrap in wet bandages
P302 + P350	IF ON SKIN: Gently wash with plenty of soap and water
P302 + P352	IF ON SKIN: Wash with plenty of soap and water
P303	IF ON SKIN (or hair)
P303 + P361	IF ON SKIN (or hair): Remove/ Take off immediately all contaminated clothing. Rinse skin with water/ shower
+ P353	contaminated clothing. Rinse skin with water/ shower
P304	IF INHALED
P304 + P312	IF INHALED: Call a POISON CENTER or doctor/physician if you feel unwell
P304 + P340	IF INHALED: Remove to fresh air and keep at rest in a position comfortable for breathing
P304 + P341	IF INHALED: If breathing is difficult, remove to fresh air and keep at rest in a position comfortable for breathing
P305	IF IN EYES
P305 + P351	IF IN EYES: Rinse cautiously with water for several minutes.
+ P338	Remove contact lenses, if present and easy to do. Continue rinsing
P306	IF ON CLOTHING
P306 + P360	IF ON CLOTHING: Rinse immediately contaminated clothing and skin with plenty of water before removing clothes
P307	IF exposed
P307 + P311	IF exposed: Call a POISON CENTER or doctor/physician
P308	IF exposed or concerned
P308 + P313	IF exposed or concerned: Get medical advice/ attention
P309	IF exposed or if you feel unwell
P309 + P311	IF exposed or if you feel unwell: Call a POISON CENTER or doctor/physician
P310	Immediately call a POISON CENTER or doctor/physician
P311	Call a POISON CENTER or doctor/physician
P312	Call a POISON CENTER or doctor/physician if you feel unwell
P313	Get medical advice/ attention
P314	Get medical advice/ attention if you feel unwell
P315	Get immediate medical advice/ attention
P320	Specific treatment is urgent (see ... on this label)
P321	Specific treatment (see ... on this label)
P322	Specific measures (see ... on this label)
P330	Rinse mouth
P331	Do NOT induce vomiting
P332	If skin irritation occurs
P332 + P313	If skin irritation occurs: Get medical advice/ attention
P333	If skin irritation or rash occurs
P333 + P313	If skin irritation or rash occurs: Get medical advice/ attention
P334	Immerse in cool water/ wrap in wet bandages
P335	Brush off loose particles from skin
P335 + P334	Brush off loose particles from skin. Immerse in cool water/ wrap in wet bandages
P336	Thaw frosted parts with lukewarm water. Do no rub affected area
P337	If eye irritation persists
P337 + P313	If eye irritation persists: Get medical advice/ attention
P338	Remove contact lenses, if present and easy to do. Continue rinsing
P340	Remove to fresh air and keep at rest in a position comfortable for breathing

P341	If breathing is difficult, remove to fresh air and keep at rest in a position comfortable for breathing
P342	If experiencing respiratory symptoms
P342 + P311	If experiencing respiratory symptoms: Call a POISON CENTER or doctor/physician
P350	Gently wash with plenty of soap and water
P351	Rinse cautiously with water for several minutes
P352	Wash with plenty of soap and water
P353	Rinse skin with water/ shower
P360	Rinse immediately contaminated clothing and skin with plenty of water before removing clothes
P361	Remove/Take off immediately all contaminated clothing
P362	Take off contaminated clothing and wash before reuse
P363	Wash contaminated clothing before reuse
P370	In case of fire
P370 + P376	In case of fire: Stop leak if safe to do so
P370 + P378	In case of fire: Use ... for extinction
P370 + P380	In case of fire: Evacuate area
P370 + P380	In case of fire: Evacuate area. Fight fire remotely due to the risk of explosion
+ P375	In case of major fire and large quantities
P371	In case of major fire and large quantities: Evacuate area.
P371 + P380	In case of major fire and large quantities: Evacuate area.
+ P375	Fight fire remotely due to the risk of explosion
P372	Explosion risk in case of fire
P373	DO NOT fight fire when fire reaches explosives
P374	Fight fire with normal precautions from a reasonable distance
P375	Fight fire remotely due to the risk of explosion
P376	Stop leak if safe to do so
P377	Leaking gas fire: do not extinguish, unless leak can be stopped safely
P378	Use ... for extinction
P380	Evacuate area
P381	Eliminate all ignition sources if safe to do so
P390	Absorb spillage to prevent material damage
P391	Collect spillage
P401	Store ...
P402	Store in a dry place
P402 + P404	Store in a dry place. Store in a closed container
P403	Store in a well-ventilated place
P403 + P233	Store in a well-ventilated place. Keep container tightly closed
P403 + P235	Store in a well-ventilated place. Keep cool
P404	Store in a closed container
P405	Store locked up
P406	Store in corrosive resistant/ ... container with a resistant inner liner
P407	Maintain air gap between stacks/ pallets
P410	Protect from sunlight
P410 + P403	Protect from sunlight. Store in a well-ventilated place
P410 + P412	Protect from sunlight. Do not expose to temperatures exceeding 50 °C/ 122 °F
P411	Store at temperatures not exceeding ... °C/ ... °F
P411 + P235	Store at temperatures not exceeding ... °C/ ... °F. Keep cool
P412	Do not expose to temperatures exceeding 50 °C/ 122 °F
P413	Store bulk masses greater than ... kg/ ... lbs at temperatures not exceeding ... °C/ ... °F
P420	Store away from other materials
P422	Store contents under ...
P501	Dispose of contents/container to ...
P502	Refer to manufacturer/ supplier for information on recovery/ recycling

Name	Page
OPTIZYME™ Aarl	400
OPTIZYME™ Aatl	401
ACES	478
Acetic Acid, Glacial	431
Acetonitrile DNA Synthesis	431
Acridine Orange Hemizinc Salt	472
Acrylamide	444
Acrylamide Solution	441
Acrylamide: Bis-Acrylamide 19:1	444
Acrylamide: Bis-Acrylamide 29:1	444
Acrylamide: Bis-Acrylamide 37.5:1	445
Actinomycin-D	452
Adenosine 5'-Triphosphate Disodium Salt Trihydrate	459
AEBSF Hydrochloride	450
Agar	474
2xYT Agar Molecular Genetics	473
Agar Molecular Genetics	474
LB Agar, Miller	475
Agarose	386
Agarose	445
Agarose	386
Agarose	385
L-Alanine	464
OPTIZYME™ Alkaline Phosphatase	396
OPTIZYME™ Alol	401
OPTIZYME™ Alul	401
OPTIZYME™ Alw44l	401
Amido Black 10B	478
Amikacin From Kanamycin A	453
4-(2-Aminoethyl)benzenesulfonyl Fluoride Hydrochloride	450
Ammonium Acetate	422
Ammonium Bicarbonate	479
Ammonium Carbonate	479
Ammonium Persulfate	445
Ammonium Phosphate, Dibasic	479
Ammonium Phosphate, Monobasic	479
Ammonium Sulfate	445
Amphotericin B	453
Amphotericin B 250µg/mL	453
Ampicillin Sodium Salt	453
Ampicillin Trihydrate	453
Antibiotic A23187 Free Acid	453
OPTIZYME™ Apal	402
Aphidicolin	460
D-(-)-Arabinose	461
L-Arginine	464
L-Arginine Hydrochloride	464
L-Arginine, Free Base	464
L-Ascorbic Acid	465
L-Asparagine Monohydrate	465
L-Aspartic Acid	465
OPTIZYME™ Aval	402
OPTIZYME™ Avall	402
Bacitracin	454
OPTIZYME™ Ball	402
OPTIZYME™ BamHI	403
100 Base-Pair Ladder	391
250 Base-Pair Ladder	395
OPTIZYME™ BclI	403
Benzamidine Hydrochloride Hydrate	451
BES	479
OPTIZYME™ BgII	403
OPTIZYME™ BgIII	403
Bicine	479
D-Biotin	462
Bis-Acrylamide	446
Bis-Acrylamide Solution	445
Bis-Tris	480
Blasticidin S Hydrochloride	
From Streptomyces griseochromogenes	454
Boric Acid	386
Bovine Serum Albumin (BSA)	437
Bovine Serum Albumin (BSA)	438
Bovine Serum Albumin (Fraction V)	437
OPTIZYME™ Bpil	404
Brij® 35	451
Brilliant Blue G-250	446
Brilliant Blue R-250	446
8-Bromoadenine	431
5-Bromo-4-chloro-3-indolyl- β-D-galactopyranoside	472
5-Bromo-4-chloro-3-indolyl Phosphate, p-Toluidine Salt	437
5-Bromo-2'-Deoxyuridine	431
Bromophenol Blue	386
Bromophenol Blue Sodium Salt	386
2X YT Broth Molecular Genetics	474
LB Broth, Lennox	474
LB Broth, Lennox Molecular Genetics	475
LB Broth, Miller	476
LB Broth, Miller	475
LB Broth, Miller Molecular Genetics	475
OPTIZYME™ Bsh1236I	404
OPTIZYME™ BshTI	404
OPTIZYME™ BssHII	404
OPTIZYME™ BstEI	405
n-Butanol	478
n-Butyl Acetate	433

Name	Page
Cacodylic Acid	438
Calcium Chloride Dihydrate	480
2M Calcium Chloride Solution	473
Calcium Ionomphore A23187	460
Calcium Phosphate Dibasic Anhydrous	481
CAPS	481
Carbenicillin Disodium Salt	454
Cefotaxime, Sodium Salt	454
CellPURE™ Bis-Tris Propane	480
CellPURE™ MOPS Sodium Salt	486
CellPURE™ Phosphate Buffered Saline	488
Cesium Chloride	422
OPTIZYME™ Cfr9I	405
CHAPS	451
CHAPSO	451
CHES	481
Chloramphenicol Crystalline Powder	454
4-Chloro-1-naphthol	436
Chloroform	432
Citric Acid Anhydrous	460
OPTIZYME™ Clal	405
Coenzyme A Trilithium Salt Dihydrate	460
Collagenase From Clostridium histolyticum	395
Control Primer no. 1,	
for use with Fisher BioReagents exACTGene PCR Kit	435
Control Primer no. 2,	
for use with Fisher BioReagents exACTGene PCR Kit	435
Control Primer no. 3,	
for use with Fisher BioReagents exACTGene PCR Kit	435
Copper (II) Sulfate Pentahydrate	460
OPTIZYME™ Csp6I	405
Cyanocobalamin	460
α-Cyclodextrin	440
Cyclopiazonic Acid	461
L-Cysteine Hydrochloride Monohydrate	465
L-Cystine	465
Cytoscint® Scintillation Cocktail	481
Cytosine-β-D-Arabinofuranoside Hydrochloride	396
Dacarbazine	462
Dansyl Chloride	439
OPTIZYME™ Ddel	406
Denhardt's Reagent	420
2'-Deoxyadenosine 5'-Triphosphate (dATP)	434
Deoxycholic Acid Sodium Salt	481
2'-Deoxycytidine 5'-Triphosphate (dCTP)	434
2'-Deoxyguanosine 5'-Triphosphate (dGTP)	434
2'-Deoxyguanosine-5'-triphosphate trisodium salt dihydrate	459
2'-Deoxythymidine 5'-Triphosphate (dTTP)	434
Dextran	420
Dextran Sulfate Sodium Salt	420
Dextrose Anhydrous	462
N,N'-Diisopropylcarbodiimide	439
N,N-Diisopropylethylamine	439
Dimethyl Sulfoxide	481
4-Dimethylaminopyridine	431
N,N-Dimethylformamide Synthesis	433
Dithioerythritol	451
Dithiothreitol	451
OPTIZYME™ T4 DNA Ligase	399
DNA Marker Fisher BioReagents	
exACTGene 100 bp PCR DNA Ladder 100 Lanes	392
DNA Marker Fisher BioReagents	
exACTGene Cloning DNA Ladder 100 Lanes	393
DNA Marker Fisher BioReagents*	
exACTGene® 100 bp DNA Ladder 100 Lanes	393
DNA Marker Fisher BioReagents*	
exACTGene® 1kb DNA Ladder 100 Lanes	394
DNA Marker Fisher BioReagents*	
exACTGene® 1kb DNA Ladder 100 Lanes	393
DNA Marker Fisher BioReagents*	
exACTGene® 1kb Plus DNA Ladder 100 Lanes	394
DNA Marker Fisher BioReagents*	
exACTGene® 24kb Max DNA Ladder 100 Lanes	394
DNA Marker Fisher BioReagents*	
exACTGene® 50 bp Mini Ladder 100 Lanes	392
DNA Marker Fisher BioReagents*	
exACTGene® Low Range DNA Ladder 100 Lanes	392
DNA Marker Fisher BioReagents*	
exACTGene® Low Range Plus DNA Ladder 100 Lanes	392
DNA Marker Fisher BioReagents*	
exACTGene® Mid Range Plus DNA Ladder 100 Lanes	393
OPTIZYME™ T4 DNA Polymerase	399
OPTIZYME DNA Polymerase I	396
Ø X-174 RF DNA-Hae III Digest	395
DNA-Hind III Digest	395
DNA-Hind III/ϕX-174 RF DNA-Hae III Digest	395
OPTIZYME™ recombinant DNase I	398
rDNase I (RNase-free)	398
OPTIZYME™ DNase I, RNase-Free	397
dNTP Mix exACTGene® 10mM solution: dATP, dCTP, dGTP and dTTP ultrapure grade (>99% triphosphate)	435
dNTP Set 4 tubes/pack	434
dNTP Set exACTGene® 100mM solution: dATP, dCTP, dGTP and dTTP ultrapure grade (>99% triphosphate)	435
Doxorubicin Hydrochloride	454
Doxycycline Hydrochloride	455

Name	Page
OPTIZYME™ DpnI	406
OPTIZYME™ Dral	406
dTPP exACTGene® 100mM solution ultrapure grade	
(>99% triphosphate)	435
OPTIZYME™ Ecl136II	406
OPTIZYME™ Eco57I	407
OPTIZYME™ EcoRI	407
OPTIZYME™ EcoRV	407
Electrode Storage Solution	482
Eosin Y	472
Erythromycin	455
OPTIZYME™ Esp3I	408
Ethanol, Absolute (200 proof)	422
Ethidium Bromide	387
Ethidium Bromide	387
Ethyl Acetate	432
Ethylene Glycol	472
Ethylenediamine Tetraacetate, DEPC-Treated Solution	482
Ethylenediamine Tetraacetic Acid	482
Ethylenediamine Tetraacetic Acid 0.5 M Solution, pH 8.0 (for DNA Work)	422
Ethylenediamine Tetraacetic Acid, Dipotassium Salt Dihydrate	482
Ethylenediamine Tetraacetic Acid, Disodium Salt Dihydrate	482
Ethylenediamine Tetraacetic Acid, Tetrasodium Salt Dihydrate	483
PCR Buffer A, 10X, 1mL, 15mM MgCl2,	
for use with Fisher BioReagents exACTGene PCR Kit	435
PCR Buffer B, 10X, 1mL, contains no MgCl2,	
for use with Fisher BioReagents exACTGene PCR Kit	435
exACTGene Complete PCR Kit	435
exACTGene PCR Core Reagent Set B, 80 rxns, Includes: Taq, Buffer B, MgCl2, and dNTPs	435
exACTGene PCR core reagent Taq DNA polymerase	435
MgCl2, 25mM, for use with Fisher BioReagents exACTGene PCR Kit	435
dATP exACTGene® 100mM solution ultrapure grade (>99% triphosphate)	435
dCTP exACTGene® 100mM solution ultrapure grade (>99% triphosphate)	435
dGTP exACTGene® 100mM solution ultrapure grade (>99% triphosphate)	435
OPTIZYME Exonuclease III	397
EZ Run protein gel staining solution	443
EZ-run buffered protein gel solution	443
EZ-Run prestain protein marker	442
EZ-Run prestained recombinant protein ladder	442
10% EZ-run protein gel solution	443
12.5% EZ-run protein gel solution	443
15% EZ-run protein gel solution	443
EZ-Run recombinant protein ladder	442
Fast Green FCF	472
Ficoll® 400	483
Folic Acid	463
Formaldehyde	387
Formamide	421
Formamide Super Pure	483
Formic Acid	432
Forskolin From Coleus forskohlii	440
G418 Sulfate	455
Gentamycin Sulfate	455
Giemsa Stain	472
L-Glutamic Acid	465
L-Glutamine	466
Y-L-Glutamyl-p-Nitroanilide, Hydrochloride	463
Glutaraldehyde Solution	473
Glutathione Reduced	463
Glycerol	483
Glycerol Gel-Loading Dye, 5X	387
Glycine	446
D-Glycogen	461
Glycyl-L-Glutamic Acid	463
Glyoxal	387
Guanidine Hydrochloride	422
Guanidine Thiocyanate	423
OPTIZYME™ HaeII	408
Hematoxylin	473
Hematoxylin Trihydrate	463
HEPES	484
HEPES Buffer	484
HEPES Sodium Salt	484
Heptane Sequencing	487
OPTIZYME™ HincII	408
OPTIZYME™ HindIII	408
OPTIZYME™ HinfI	409
L-Histidine	466
L-Histidine Hydrochloride Monohydrate	466
OPTIZYME™ Hpal	409
OPTIZYME™ HpaII	409
OPTIZYME™ Hpy8I	409
Hyamine Hydroxide®	484
Hybridization Cocktail	421
Hydrogen Peroxide	484
8-Hydroxyquinoline	421
Hygromycin B	455
Imidazole	484
Ionomycin, Calcium Salt	455
Isoamyl Alcohol	485

Name	Page
L-Isoleucine	466
Isopropanol	423
Isopropyl-β-D-thiogalactopyranoside	463
Isopropyl-β-D-thiogalactopyranoside	464
JustPURE™ BIS-TRIS Propane	480
JustPURE™ EPPS	482
JustPURE™ Gly-Gly	483
JustPURE™ MES Hydrate	486
JustPURE™ MOPS Sodium Salt	487
Kanamycin Sulfate	456
Ketoconazole	456
Kilobase DNA Marker	395
Kinetin (6-Furfurylaminopurine)	464
OPTIZYME™ Klenow Fragment	397
OPTIZYME™ KpnI	410
50 Base-Pair Ladder	395
L-Leucine	466
Leupeptin Hemisulfate	452
OPTIZYME™ Lglul	410
Lithium Dodecyl Sulfate	446
Lymphocyte Separation Medium	468
L-Lysine Hydrochloride	466
Lysozyme, Egg White	396
Magnesium Acetate Tetrahydrate	485
Magnesium Chloride Hexahydrate	485
2M Magnesium Chloride Solution	473
Magnesium Sulfate Heptahydrate	485
D-(+)-Maltose Monohydrate	461
Manganese Chloride 1.0M Solution	485
D-Mannitol	462
2-Mercaptoethanol	450
β-Mercaptoethylamine Hydrochloride	440
MES	446
Methanol	388
L-Methionine	467
Methotrexate	468
Methyl α-D-Mannopyranoside	468
Methylene Blue	473
Methylene Chloride	432
N-Methylmorpholine	439
N-Methylpyrrolidinone	426
4-Methylumbelliferyl Phosphate Disodium Salt	440
Metolazone	456
Miconazole	456
Mineral Oil	436
Mitomycin C	456
OPTIZYME™ MluI	410
OPTIZYME™ M-MLV Reverse Transcriptase	397
Molecular Sieves, Grade 513	486
Molecular sieves, Grade 514 (Type 4A 8-12 Mesh)	486
MOPS	429
MOPS	486
OPTIZYME™ MspI	410
OPTIZYME™ NaeI	411
Nalidixic Acid	456
OPTIZYME™ NcoI	411
OPTIZYME™ NdeI	411
Neomycin Sulfate	457
OPTIZYME™ NheI	411
(-)-Nicotine Tartrate	459
Nitroliotriacetic Acid	469
Nitro Blue Tetrazolium	469
p-Nitrophenol	439
p-Nitrophenyl Phosphate Disodium Salt, Hexahydrate	438
Nonquantitative DNA Ladders	394
OPTIZYME™ NotI	412
OPTIZYME™ NsiI	412
Nystatin	457
NZCYM Broth Molecular Genetics	476
NZY Broth Molecular Genetics	476
n-Octyl-β-D-Glucopyranoside	487
OPTIZYME™ 10X Buffer 1	418
OPTIZYME™ 10X Buffer 2	418
OPTIZYME™ 10X Buffer 3	418
OPTIZYME™ 10X Buffer 4	418
OPTIZYME™ 10X Buffer 5	418
OPTIZYME™ 10X Buffer AarI, Scal, PstI	419
OPTIZYME™ 10X Buffer BamHI	419
OPTIZYME™ 10X Buffer Cfr9I	419
OPTIZYME™ 10X Buffer Ecl136II, SacI	419
OPTIZYME™ 10X Buffer EcoRI	419
OPTIZYME™ 10X Buffer KpnI	419
OPTIZYME™ 10X Buffer Sau3AI	420
OPTIZYME™ 10X Buffer TaqI	420
L-Ornithine Hydrochloride	467
Paromomycin Sulfate	457
OPTIZYME™ PstI	412
PCR Nucleotide Mix	436
Penicillin-G Potassium Salt	458
Penicillin G Sodium Salt	457
Penicillin/Streptomycin Mixture	457
Penicillin/Streptomycin/Glutamine Mixture	457
Penicillin/Streptomycin/Neomycin Mixture	457
Pepstatin A	452
Peptone	476
Peptone Molecular Genetics	476

Name	Page
Periodic Acid	488
OPTIZYME™ PfoI	412
Phenol	489
Phenol, Saturated	424
Phenol/Chloroform/Isoamyl Alcohol	424
Phenol/Chloroform/Isoamyl Alcohol	423
L-Phenylalanine	467
o-Phenylenediamine	469
Phleomycin, 2% solution in water	458
Phorbol 12-Myristate 13-Acetate	469
Phosphate Buffered 10X Powder Concentrate	488
Phosphate Buffered 1X Powder Concentrate	488
Phosphate Buffered Saline	469
Phosphate Buffered Saline Tablets	487
Phosphate Buffered Saline with Tween 20*	487
PIPES, Sesquisodium Salt	489
Pleiotrophin	470
Polyethylene Glycol 8000	470
Polymyxin B Sulfate	458
OPTIZYME™ T4 Polynucleotide Kinase	399
Polyvinylpyrrolidone	441
Ponceau S	447
Potassium Persulfate	447
Potassium Phosphate Dibasic	489
Potassium Phosphate Monobasic	489
L-Proline	467
n-Propanol	433
Protease	441
Protein Gel-Loading Dye for SDS-PAGE	447
Proteinase K	424
OPTIZYME™ PstI	413
Puromycin Dihydrochloride	458
OPTIZYME™ PvuI	413
OPTIZYME™ PvuII	413
Pyridine	433
Pyridoxal Hydrochloride	441
Pyridoxal 5-Phosphate Monohydrate	441
Pyridoxine Hydrochloride	470
Rapamycin	458
Riboflavin	470
Ribonuclease H	400
Ribonuclease Inhibitor	398
D-(-)-Ribose	461
Rifampicin	458
OPTIZYME™ SP6 RNA Polymerase	398
OPTIZYME™ T3 RNA Polymerase	399
OPTIZYME™ T7 RNA Polymerase	400
OPTIZYME™ RsaI	413
OPTIZYME™ SacI	414
OPTIZYME™ SacII	414
OPTIZYME™ SalI	414
Saline-Sodium Citrate	421
Saline-Sodium Phosphate-EDTA	421
Sarkosyl	452
OPTIZYME™ Sau3AI	414
OPTIZYME™ Scal	415
Sequencing Gel-Loading Dye, 3X	388
L-Serine	467
OPTIZYME™ SfaI	415
Silver Nitrate	448
OPTIZYME™ SmaI	415
SOB Broth	477
SOC Medium	477
Sodium Acetate Anhydrous	425
Sodium Acetate Trihydrate	425
Sodium Azide	489
Sodium Bicarbonate	489
Sodium Cacodylate Trihydrate	439
Sodium Carbonate	490
Sodium Chloride	490
Sodium Chloride-Tris-EDTA (STE) Buffer	425
Sodium Chloride-Tris-EDTA Buffer	425
Sodium Citrate Dihydrate	490
Sodium Dodecyl Sulfate	447
Sodium Dodecyl Sulfate	448
Sodium Dodecyl Sulfate	447
Sodium Hydroxide	490
Sodium Iodide	490
Sodium Nitrate	491
Sodium Oxalate	470
Sodium Persulfate, >98%	448
Sodium Phosphate Dibasic Anhydrous	491
Sodium Phosphate Dibasic Heptahydrate	491
Sodium Phosphate Monobasic Anhydrous	491
Sodium Phosphate Monobasic Monohydrate	491
Sodium Pyruvate	470
Sodium Sulfate Anhydrous	470
Sodium Sulfite Anhydrous	471
Sodium Tartrate Dihydrate	471
Sodium Tetraborate Decahydrate	491
Sorbitol	441
Spectinomycin Pentahydrate Dihydrochloride	458
OPTIZYME™ Spel	415
Spermidine Trihydrochloride	425
OPTIZYME™ SphI	416
Staurosporine	459

Name	Page
STET Buffer, 1X Solution	492
Streptomycin Sulfate	459
OPTIZYME™ Stul	416
Succinic Acid	471
D-Sucrose	462
Sucrose Gel-Loading Dye for DNA Gels	448
Sudan Black	492
Sulfosalicylic Acid Dihydrate	438
Sureprep Cleanup and Concentration Kit	428
Sureprep leukocyte RNA purification kit	427
Sureprep nuclear or cytoplasmic RNA purification kit	429
Sureprep RNA/DNA/protein purification kit	430
Sureprep RNA/protein purification kit	428
Sureprep Small RNA Purification Kit	427
Sureprep TrueTotal RNA Purification Kit	427
SurePrep Urine Bacterial RNA Purification Kit	429
Sureprep urine exfoliated cell RNA purification kit	430
SurePrep* FFPE RNA Isolation Kit	425
SurePrep* Plant/Fungi Total RNA Purification Kit	426
SurePrep* Soil DNA Isolation Kit	426
SurePrep* Water RNA/DNA Purification Kit	426
TAPS (Free Acid)	492
Taq DNA Polymerase	436
OPTIZYME™ TaqI	416
TBE Buffer Mix	388
TEMED	448
Template, Control, Lambda DNA, used to troubleshoot PCR parameters	435
OPTIZYME™ Terminal Deoxynucleotidyl Transferase	400
Terrific Broth Molecular Genetics	477
Terrific Broth, Modified	477
TES	492
Tetracycline hydrochloride	459
Tetrahydrofuran	434
Thapsigargin	452
Thiamine Hydrochloride	471
Thimerosal	471
DL-Thioctic Acid	462
L-Threonine	467
Toluidine Blue O	472
Transferrin From Guinea Pig	471
D-(+)-Trehalose, Dihydrate	461
Trichloroacetic Acid	492
Tricine	492
Triethylamine	440
Trifluoroacetic Acid	440
2,2,2-Trifluoroethanol	438
Tris Base	493
Tris Base	493
Tris Buffer	389
Tris Buffered Saline	493
Tris Hydrochloride	494
Tris Hydrochloride Solution pH 7.0	389
Tris Hydrochloride Solution pH 7.5	389
Tris Hydrochloride Solution pH 8.0	390
Tris-Acetate-EDTA	388
Tris-Acetate-EDTA	390
Tris-Acetate-EDTA	391
Tris-Borate-EDTA	388
Tris-Borate-EDTA	391
Tris-EDTA	494
Tris-EDTA	389
Tris-EDTA	494
Tris-EDTA (TE), 10X Solution	391
Tris-Glycine	449
Tris-Glycine-SDS	448
Tris-Glycine-SDS	449
Triton® X-100	452
OPTIZYME™ Tru9I	416
Tryptone	477
L-Tryptophan	468
Tween® 20	450
Tween® 80	450
L-Tyrosine	468
Urea	450
L-Valine	468
Vancomycin Hydrochloride	459
Vimentin	472
OPTIZYME™ VspI	417
Water	495
Water, Sterile	495
Water, Sterile	495
OPTIZYME™ XbaI	417
2 XYT Broth	474
OPTIZYME™ XhoI	417
OPTIZYME™ XmaJI	417
OPTIZYME™ XmnI	418
Xylene Cyanole FF	391
D-(+)-Xylose	461
Yeast DNA PFGE Markers	441
Yeast Extract	478
Yeast Extract Molecular Genetics	478
YPD Broth Molecular Genetics	478

Mfr. No	Page
A/0080/08	13
A/0080/17	13
A/0160/48	13
A/0200/50	13
A/0200/53	13
A/0360/25	15
A/0360/25	234
A/0360/27	15
A/0360/PB08	15
A/0360/PB15	15
A/0360/PB15	234
A/0360/PB17	15
A/0360/PB17	234
A/0400/25	15
A/0400/25	234
A/0400/26	15
A/0400/26	234
A/0400/PB08	15
A/0400/PB08	234
A/0400/PB15	15
A/0400/PB17	234
A/0400/PC15	15
A/0400/PC15	234
A/0400/PC17	15
A/0400/PC17	234
A/0406/PB08	13
A/0406/PB08	309
A/0406/PB15	13
A/0406/PB15	309
A/0407/PB08	13
A/0407/PB08	312
A/0411/PB07	14
A/0411/PB07	25
A/0411/PB08	14
A/0411/PB08	25
A/0411/PB15	14
A/0411/PB15	25
A/0411/PB17	14
A/0411/PB17	25
A/0420/PB17	15
A/0430/PB08	15
A/0440/PB08	16
A/0440/PB17	16
A/0480/PB08	16
A/0480/PB17	16
A/0520/17	18
A/0520/25	18
A/0520/PB17	18
A/0560/08	17
A/0560/15	17
A/0560/17	17
A/0560/21	17
A/0560/25	17
A/0560/27	17
A/0560/FP21	17
A/0560/PB17	17
A/0560/PC17	17
A/0560/PC21	17
A/0560/PC24	17
A/0560/PC25	17
A/0600/15	17
A/0600/15	234
A/0600/17	17
A/0600/17	234
A/0600/21	17
A/0600/21	234
A/0600/25	17
A/0600/25	234
A/0600/27	17
A/0600/27	234

Mfr. No	Page
A/0600/PB08	17
A/0600/PB17	17
A/0600/PC17	17
A/0600/PC17	234
A/0600/PC21	17
A/0600/PC21	234
A/0601/17	17
A/0601/21	17
A/0601/25	17
A/0602/15	17
A/0602/17	17
A/0603/15	16
A/0603/17	16
A/0604/17	16
A/0606/10RSS	16
A/0606/10RSS	309
A/0606/15	16
A/0606/15	309
A/0606/17	16
A/0606/17	309
A/0606/25RSS	16
A/0606/25RSS	309
A/0606/30RSS	16
A/0606/30RSS	309
A/0620/21	21
A/0620/25	21
A/0620/27	21
A/0620/PB08	21
A/0620/PB17	21
A/0621/21	21
A/0621/25	21
A/0621/27	21
A/0621/PB15	21
A/0621/PB17	21
A/0622/PB15	21
A/0622/PB17	21
A/0623/15	20
A/0623/17	20
A/0624/17	20
A/0626/10RSS	20
A/0626/10RSS	310
A/0626/15	20
A/0626/15	310
A/0626/17	20
A/0626/17	310
A/0626/21	20
A/0626/21	310
A/0626/25RSS	20
A/0626/25RSS	310
A/0626/30RSS	20
A/0626/30RSS	310
A/0626/PB17	20
A/0626/PB17	310
A/0627/10RSS	19
A/0627/10RSS	234
A/0627/10RSS	309
A/0627/15	19
A/0627/15	234
A/0627/17	19
A/0627/17	234
A/0627/17	309
A/0627/21	19
A/0627/21	234
A/0627/21	309
A/0627/25RSS	19
A/0627/25RSS	234
A/0627/25RSS	309
A/0627/30RSS	19
A/0627/30RSS	234
A/0627/30RSS	309

Mfr. No	Page
A/0627/PB17	19
A/0627/PB17	234
A/0627/PB17	309
A/0628/15	21
A/0628/17	21
A/0628/PB17	21
A/0630/PB15	20
A/0630/PB15	312
A/0630/PB17	20
A/0630/PB17	312
A/0632/PB15	19
A/0632/PB15	309
A/0632/PB17	19
A/0632/PB17	309
A/0635/PB15	20
A/0635/PB15	312
A/0635/PB17	20
A/0635/PB17	312
A/0636/15	19
A/0636/15	309
A/0636/17	19
A/0636/17	309
A/0640/PB08	22
A/0650/PB15	18
A/0650/PB15	309
A/0650/PB17	18
A/0650/PB17	309
A/0800/PB08	22
A/0802/PB05	22
A/0880/PB08	23
A/0881/99	23
A/0920/PB05	23
A/0920/PB08	23
A/0960/53	23
A/0976/50	26
A/0977/50	26
A/0977/60	26
A/0977/70	26
A/1000/53	26
A/1020/99	26
A/1040/53	116
A/1080/48	26
A/1080/53	26
A/1080/68	26
A/1160/48	26
A/1162/48	26
A/1278/46	26
A/1280/48	26
A/1280/50	26
A/1460/45	27
A/1460/45	158
A/1460/45	158
A/1520/46	27
A/1520/46	158
A/1552/15	27
A/1565/05	27
A/1566/47	27
A/1570/50	28
A/1570/53	28
A/1600/50	28
A/1605/50	28
A/1605/53	28
A/1610/50	28
A/1610/53	28
A/1690/50	28
A/1810/99	29
A/1840/60	29

Mfr. No	Page
A/1880/53	29
A/1960/53	29
A/1960/53	234
A/2000/53	29
A/2080/60	29
A/2160/45	184
A/2200/50	30
A/2200/53	30
A/2200/60	30
A/2200/65	30
A/2202/53	29
A/2202/60	29
A/2280/53	30
A/2280/60	30
A/2320/53	30
A/2360/60	30
A/2370/53	30
A/2372/53	30
A/2374/53	30
A/2380/53	30
A/2400/53	30
A/2400/60	30
A/2400/63	30
A/2520/61	31
A/2560/53	31
A/2560/60	31
A/2560/65	31
A/2600/53	31
A/2840/48	32
A/2930/43	32
A/2960/08	33
A/3040/46	32
A/3042/46	32
A/3150/43	33
A/3150/44	33
A/3222/21	35
A/3222/PB17	35
A/3240/PB15	36
A/3240/PB17	36
A/3280/21	35
A/3280/25	35
A/3280/PB15	35
A/3280/PB17	35
A/3285/PB17	35
A/3290/PB08	35
A/3290/PB15	35
A/3295/PB05	34
A/3295/PB05	310
A/3320/PB15	34
A/3320/PB17	34
A/3320/PB21	34
A/3360/PB15	34
A/3360/PB17	34
A/3366/15	34
A/3366/21	34
A/3367/15	34
A/3367/17	34
A/3400/53	37
A/3400/60	37
A/3400/65	37
A/3440/50	37
A/3440/53	37
A/3440/60	37
A/3440/65	37
A/3445/53	37
A/3445/60	37
A/3446/50	36
A/3446/50	310
A/3450/50	36
A/3450/50	312
A/3520/53	37

Mfr. No	Page	Mfr. No	Page	Mfr. No	Page	Mfr. No	Page	Mfr. No	Page	Mfr. No	Page	Mfr. No	Page	Mfr. No	Page	Mfr. No	Page
A/5360/53	42	A/7360/53	48	A463-1	152	A511-P500	24	B/3110/46	58	B/4791/46	310	BP1130-500	433	BP1408-1	444	BP175-500	491
A/5360/60	42	A/7486/50	48	A463-250	24	A511-P500	228	B/3200/48	58	B/4800/08	68	BP1135-500	433	BP1410-1	445	BP17501-100	424
A/5400/53	42	A/7490/PB05	48	A463-250	152	A512-P500	24	B/3620/44	58	B/4800/15	68	BP114-25	386	BP1420-100	476	BP17501-400	424
A/5680/48	43	A/7600/PB08	49	A463-500	24	A512-P500	33	B/3640/03	361	B/4800/17	68	BP1140-1	434	BP1420-2	476	BP17511-100	424
A/5680/50	43	A/7680/48	49	A463-500	152	A513-P500	24	B/3647/45	361	B/4800/21	68	BP1145-1	432	BP1420-500	476	BP17511-400	424
A/5680/70	43	A/7760/46	49	A465-1	14	A513-P500	152	B/3750/53	59	B/4800/25	68	BP115-25	386	BP1422-100	478	BP17521-100	424
A/5720/48	43	A/7770/50	49	A465-1	24	A955-1	18	B/3750/53	234	B/4800/MC15	68	BP1150-500	485	BP1422-2	478	BP17521-400	424
A/5720/50	43	A/7770/53	49	A465-250	14	A955-1	308	B/3750/60	59	B/4800/PB17	68	BP1155-500	433	BP1422-500	478	BP17531-100	423
A/5720/53	43	A/7770/60	49	A465-250	24	A955-212	18	B/3750/60	234	B/4850/08	68	BP116-10	472	BP1423-2	474	BP17531-400	423
A/5760/48	43	A/7840/53	50	A465-500	14	A955-212	308	B/3750/65	59	B/4850/15	68	BP1160-4	433	BP1423-500	474	BP17541-100	423
A/5760/50	43	A/8200/50	50	A465-500	24	A955-4	18	B/3750/65	234	B/4850/17	68	BP1160-500	433	BP1426-2	475	BP17541-400	423
A/5760/53	43	A/8202/50	50	A466-1	24	A955-4	308	B/3760/53	59	B/4850/25	68	BP117-100	473	BP1426-500	475	BP1755-1	463
A/5760/60	43	A/8240/53	50	A466-1	147	A955-500	18	B/3760/53	234	B/4850/PB17	68	BP1170-4	431	BP1427-2	475	BP1755-10	463
A/5760/70	43	A/8280/46	50	A466-250	24	A955-500	308	B/3760/60	59	B/4860/PB08	67	BP1170-450	431	BP1427-500	475	BP1755-100	463
A/5800/50	44	A/8300/46	50	A466-250	147	B/0052/50	52	B/3760/60	234	B/4900/25	69	BP1172-4	426	BP150-100	448	BP1756-100	389
A/5800/53	44	A/8440/53	51	A466-500	24	B/0100/48	52	B/3760/65	59	B/4900/PB08	69	BP118-500	482	BP150-20	448	BP1756-500	389
A/5880/53	44	A/8880/48	51	A466-500	147	B/0100/53	52	B/3760/65	234	B/4900/PB17	69	BP1185-500	431	BP151-100	452	BP1757-100	389
A/5880/60	44	A/8880/53	51	A467-1	24	B/0102/48	52	B/3760/70	59	B/4950/08	70	BP1186-4	432	BP151-500	452	BP1757-500	389
A/5880/65	44	A/8880/60	51	A467-1	215	B/0120/48	52	B/3760/70	234	B/4950/17	70	BP119-500	482	BP152-1	493	BP1758-100	390
A/5920/53	44	A/8882/48	51	A467-250	24	B/0150/46	53	B/3800/53	59	B/4950/21	70	BP120-1	482	BP152-10	493	BP1758-500	390
A/5920/60	44	A/8882/53	51	A467-250	215	B/0350/53	53	B/3800/60	59	B/4950/24	70	BP120-500	482	BP152-25	493	BP1759-100	389
A/5920/62	44	A/8961/53	51	A467-500	24	B/0400/48	53	B/3800/62	59	B/4950/25	70	BP121-500	483	BP152-5	493	BP1759-500	389
A/6000/53	44	A/9080/45	51	A467-500	215	B/0400/53	53	B/3800/70	59	B/4950/27	70	BP1215-500	432	BP152-500	493	BP176-100	450
A/6040/48	44	A/P470/46	44	A468-1	24	B/0500/53	54	B/3820/53	59	B/4951/15	70	BP123-10	472	BP153-1	494	BP1760-25	453
A/6040/53	44	A113-10X1AMP	13	A468-1	336	B/0500/60	54	B/3850/53	60	B/4951/17	70	BP124-10	478	BP153-500	494	BP1760-5	453
A/6040/60	44	A113-10X1AMP	308	A468-250	24	B/0500/65	54	B/4000/PB05	24	B/4951/24	70	BP1302-10	387	BP154-1	493	BP1761-1	389
A/6120/60	45	A113-1AMP	13	A468-250	336	B/0500/70	54	B/4000/PB08	60	B/4951/27	70	BP1306-1	449	BP1550-50	436	BP1761-500	389
A/6120/63	45	A113-1AMP	308	A468-500	24	B/0550/53	53	B/4050/PB05	60	B/4960/17	70	BP1306-4	449	BP1575-50	421	BP177-500	438
A/6160/53	45	A113-50	13	A468-500	336	B/0550/60	53	B/4050/PB08	60	B/4960/21	70	BP1307-1	449	BP1580-100	420	BP178-1	422
A/6160/60	45	A113-50	308	A469-1	24	B/0550/63	53	B/4055/PB05	60	B/5000/17	69	BP1309-25	446	BP1585-100	420	BP178-500	422
A/6200/44	45	A114-50	36	A469-1	228	B/0550/65	53	B/4300/07	60	B/5000/21	69	BP1309-5	446	BP1585-500	420	BP179-100	445
A/6200/46	45	A114-50	308	A469-250	24	B/0550/70	53	B/4310/PB05	61	B/5100/25	69	BP1311-1	448	BP1595-1	422	BP179-25	445
A/6220/46	45	A115-50	40	A469-250	228	B/0620/44	54	B/4320/43	61	B/5100/27	69	BP1311-200	448	BP1595-500	422	BP180-100	447
A/6280/53	45	A115-50	308	A469-500	24	B/0620/44	158	B/4320/43	158	B/5100/PB17	69	BP1325-1	421	BP160-100	385	BP210-1	422
A/6400/50	45	A116-05AMP	308	A469-500	228	B/0650/53	54	B/4320/44	61	B/5150/08	68	BP1325-20	421	BP160-500	385	BP210-100	422
A/6440/53	46	A116-05AMP	358	A470-1	24	B/0650/60	54	B/4320/44	158	B/5150/25	68	BP1325-4	421	BP1600-100	437	BP210-500	422
A/6440/60	46	A116-10X1AMP	308	A470-1	33	B/0700/53	54	B/4320/46	61	B/5150/PB17	68	BP1328-1	421	BP1605-100	437	BP212-212	445
A/6440/62	46	A116-10X1AMP	358	A470-250	24	B/0700/60	54	B/4320/46	158	B/5200/15	68	BP1328-20	421	BP161-100	386	BP212R-1	445
A/6440/65	46	A116-1AMP	308	A470-250	33	B/0800/53	55	B/4320L/08	61	B/5200/25	68	BP1328-4	421	BP1610-100	437	BP213-1	485
A/6480/53	46	A116-1AMP	358	A470-500	24	B/0800/70	55	B/4330/44	61	B/5200/PB17	68	BP1330-1	390	BP1610-500	437	BP214-500	485
A/6480/60	46	A116-2AMP	308	A470-500	33	B/0850/50	54	B/4330/44	158	B/5210/15	68	BP1331-1	388	BP1615-1	472	BP215-500	485
A/6480/62	46	A116-2AMP	358	A471-500	24	B/0850/53	54	B/4330/44	158	B/5210/PB17	68	BP1332-1	390	BP1615-100	472	BP220-1	462
A/6480/65	46	A116-50	308	A471-500	146	B/0850/60	54	B/4330/46	61	B/5250/08	69	BP1332-20	390	BP162-100	385	BP220-10	462
A/6480/70	46	A116-50	358	A507-P1	14	B/0850/70	54	B/4330/46	158	B/5250/17	69	BP1332-4	390	BP1620-1	464	BP220-212	462
A/6560/PB08	45	A117-05AMP	132	A507-P1	24	B/1080/53	55	B/4330/46	158	B/5250/25	69	BP1332-500	390	BP1620-10	464	BP221-1	423
A/6560/PB17	45	A117-05AMP	308	A507-P212	14	B/1120/53	55	B/4330L/08	61	B/5251/08	69	BP1333-1	391	BP163-10	445	BP221-250	423
A/6590/53	341	A117-10X1AMP	132	A507-P212	24	B/1120/61	55	B/4350/PB07	61	B/5251/17	69	BP1333-20	391	BP163-25	445	BP226-100	489
A/6640/50	46	A117-10X1AMP	308	A507-P500	14	B/1300/60	55	B/4350/PB08	61	B/5260/PB07	70	BP1333-4	391	BP164-100	385	BP226-500	489
A/6640/53	46	A117-1AMP	132	A507-P500	24	B/1400/PB08	55	B/4400/PB05	61	B/5300/08	71	BP1334-1	388	BP164-25	385	BP227-100	421
A/6640/60	46	A117-1AMP	308	A508-P1	24	B/1450/PB07	55	B/4450/PB07	61	B/5300/17	71	BP1335-1	390	BP164-500	385	BP227-500	421
A/6680/50	46	A117-2AMP	132	A508-P1	147	B/1460/48	55	B/4540/PB05	61	B/5500/PB08	71	BP1335-20	390	BP165-25	385	BP228-100	483
A/6680/53	46	A117-2AMP	308	A508-P212	24	B/1800/PB07	56	B/4550/PB05	62	B/5550/PB08	67	BP1335-4	390	BP166-100	447	BP229-1	483
A/6680/60	46	A117-50	132	A508-P212	147	B/1900/53	56	B/4630/44	62	B/P560/53	58	BP1335-500	390	BP166-5	447	BP229-4	483
A/6720/53	47	A117-50	308	A508-P500	24	B/1900/60	56	B/4630/44	158	B/P560/60	58	BP1338-1	391	BP166-500	447	BP230-1	472
A/6720/70	47	A456-1	199	A508-P500	147	B/1900/70	56	B/4630/44	158	B/P560/65	58	BP1338-4	391	BP167-50	470	BP230-4	472
A/6800/08	47	A456-1	308	A509-P1	24	B/1950/48	56	B/4630/46	62	B/P620/44	62	BP1339-1	389	BP168-1	386	BP231-1	481
A/6880/08	47	A456-212	199	A509-P1	216	B/1950/53	56	B/4630/46	158	B/P620/46	62	BP1341-1	449	BP168-500	386	BP231-100	481
A/6880/17	47	A456-212	308	A509-P212	24	B/1950/60	56	B/4630/46	158	BP100-25	446	BP1341-4	449	BP169-10	450	BP231-4	481
A/6920/08	47	A456-4	199	A509-P212	216	B/2050/50	56	B/4630L/05	62	BP100-50	446	BP1342-1	449	BP169-212	450	BP232-1	462
A/6920/17	47	A456-4	308	A509-P500	24	B/2150/48	56	B/4630L/08	62	BP101-25	446	BP1356-100	385	BP169-500	450	BP233-1	470
A/6960/08	47	A456-500	199	A509-P500	216	B/2300/PB07	56	B/4631L/08	62	BP101-50	446	BP1356-500	385	BP170-100	444	BP233-100	470
A/6960/PB17	47	A456-500	308	A510-P1	24	B/2300/PB08	56	B/4660/43	62	BP102-1	387	BP1360-100	386	BP170-5	444	BP234-500	452
A/7000/08	47	A461-1	167	A510-P1	336	B/2350/PB08	56	B/4690/44	62	BP102-5	387	BP1364-100	444	BP170-500	444	BP2410-100	385
A/7000/PB17	47	A461-1	308	A510-P212	24	B/2550/08	57	B/4690L/08	62	BP103-10	447	BP1366-100	444	BP1700-100	424	BP2413-500	479
A/7190/44	27	A461-212	167	A510-P212	336	B/2550/08	234	B/4760/77	64	BP107-10	472	BP1368-100	445	BP1700-50	424	BP2414-500	479
A/7240/PB08	48	A461-212	308	A510-P500	24	B/2550/PB17	57	B/4760/77	65	BP108-1							



Mfr. No	Page	Mfr. No	Page	Mfr. No	Page	Mfr. No	Page	Mfr. No	Page	Mfr. No	Page	Mfr. No	Page	Mfr. No	Page	Mfr. No	Page	Mfr. No	Page
BP2483-500	482	BP2534-10	438	BP2631-500	486	BP2677-100	470	BP299-500	484	BP3602-1	442	BP531-25	387	BP696-25	461	BP8060-1	406	BP928-250	453
BP2484-100	495	BP2534-5	438	BP2633-500	484	BP2677-50	470	BP300-100	446	BP3602-25	442	BP531-500	387	BP708-250	461	BP8060-5	406	BP942-1	464
BP2484-50	495	BP2534-50	438	BP2634-500	486	BP2679-1	458	BP301-100	480	BP3602-500	442	BP535-1	396	BP7700-500	443	BP8064-1	416	BP9700-100	437
BP2485-20	495	BP2534-500	438	BP2637-1	448	BP2679-25	458	BP303-100	478	BP3603-1	442	BP535-10	396	BP7710-100	443	BP8066-1	407	BP9703-100	438
BP2485-4	495	BP2536-50	470	BP2641-1	474	BP2679-250	458	BP304-100	489	BP3603-500	442	BP535-100	396	BP7710-30	443	BP8067-1	410	BP9704-100	437
BP2501-25	440	BP2537-1	469	BP2641-100	474	BP2679-5	458	BP304-500	489	BP361-500	479	BP535-25	396	BP7710-500	443	BP8068-1	405	BP9705-100	438
BP2501-5	440	BP2537-250	469	BP2641-500	474	BP2682-1	462	BP305-50	484	BP362-1	489	BP535-5	396	BP7712-100	443	BP8069-1	400	BP9722-2	474
BP2502-5	453	BP2537-50	469	BP2643-1	453	BP2682-25	462	BP308-100	486	BP362-500	489	BP541-1	485	BP7712-30	443	BP8070-1	408	BP9722-500	474
BP2504-1	461	BP2540-1	425	BP2643-250	453	BP2682-5	462	BP308-500	486	BP3620-1	443	BP541-100	485	BP7712-500	443	BP8071-1	404	BP9723-2	475
BP2504-100	461	BP2540-25	425	BP2644-1	450	BP2687-1	461	BP309-100	492	BP3620-4	443	BP551-100	492	BP7715-100	443	BP8072-1	404	BP9723-500	475
BP2504-25	461	BP2540-5	425	BP2644-10	450	BP2687-10	461	BP310-1	484	BP363-1	489	BP551-25	492	BP7715-500	443	BP8073-1	412	BP9724-2	475
BP2505-1	464	BP2541-1	459	BP2644-100	450	BP2687-100	461	BP310-100	484	BP363-500	489	BP551-500	492	BP8000-1	413	BP8075-1	401	BP9724-500	475
BP2505-100	464	BP2541-100	459	BP2644-50	450	BP2687-25	461	BP310-5	484	BP369-100	464	BP555-1	492	BP8001-1	413	BP8076-1	415	BP9725-2	476
BP2505-500	464	BP2541-500	459	BP2644-500	450	BP2734-1	456	BP310-500	484	BP370-100	464	BP555-250	492	BP8002-1	408	BP8077-1	412	BP9725-500	476
BP2507-1	431	BP2542-10	471	BP2645-20	453	BP2734-100	456	BP315-100	492	BP372-100	464	BP555-500	492	BP8002-5	408	BP8078-1	404	BP9726-2	477
BP2507-100	431	BP2542-100	471	BP2645-50	453	BP2734-50	456	BP318-100	481	BP373-100	465	BP561-1	495	BP8003-5	407	BP8079-1	409	BP9726-500	477
BP2507-500	431	BP2542-25	471	BP2646-1	479	BP2734-500	456	BP321-100	481	BP374-100	465	BP565-10	391	BP8004-1	412	BP8080-1	406	BP9727-2	478
BP2508-1	431	BP2542-500	471	BP2646-100	479	BP2800-50	427	BP321-500	481	BP376-100	465	BP571-1	451	BP8004-5	412	BP8081-1	405	BP9727-500	478
BP2508-250	431	BP2544-1	471	BP2646-250	479	BP2801-25	427	BP3215-1	400	BP377-100	465	BP571-5	451	BP8005-5	403	BP8082-1	417	BP9728-2	477
BP2508-5	431	BP2544-10	471	BP2647-100	454	BP2802-50	430	BP3222-1	398	BP378-100	465	BP575-500	451	BP8006-5	408	BP8083-1	410	BP9728-500	477
BP2508-500	431	BP2544-5	471	BP2647-25	454	BP2803-50	430	BP3222-5	398	BP379-100	466	BP581-100	488	BP8007-1	416	BP8084-1	418	BP9729-600	477
BP2510-1	460	BP2545-100	472	BP2647-50	454	BP2804-50	429	BP3226-1	398	BP381-1	446	BP581-25	488	BP8008-1	417	BP8085-1	418	BP9730-500	474
BP2510-100	460	BP2545-250	472	BP2648-1	454	BP2805-50	429	BP3226-2	398	BP381-5	446	BP585-1	487	BP8008-5	417	BP8086-1	418	BP9731-500	476
BP2510-25	460	BP2546-100	448	BP2648-250	454	BP2806-50	428	BP3227-1	398	BP381-500	446	BP585-25	487	BP8009-1	406	BP8087-1	418	BP9732-500	476
BP2510-250	460	BP2547-1	473	BP2648-5	454	BP2807-50	427	BP323-100	490	BP382-100	466	BP585-5	487	BP8010-1	417	BP8088-1	418	BP9733-500	474
BP2510-50	460	BP2548-1	473	BP2649-1	395	BP2809-50	428	BP324-100	438	BP383-100	466	BP585-500	487	BP8010-5	417	BP8089-1	419	BP9734-500	475
BP2512-10	396	BP2548-4	473	BP2653-1	455	BP2814-25	426	BP325-50	439	BP384-100	466	BP590-100	439	BP8011-1	415	BP8090-1	419	BP9735-500	475
BP2512-100	396	BP2550-100	395	BP2653-5	455	BP2815-50	426	BP326-1	422	BP385-100	466	BP592-500	439	BP8011-5	415	BP8091-1	420	BP9736-500	474
BP2512-5	396	BP2551-50	391	BP2656-1	463	BP2816-50	425	BP326-500	422	BP386-100	466	BP595-1	460	BP8012-1	407	BP8092-1	420	BP9737-500	477
BP2512-500	396	BP2552-100	395	BP2656-10	463	BP2817-50	426	BP327-1	490	BP388-100	467	BP595-10	460	BP8013-1	414	BP8093-1	419	BP9740-10X5	477
BP2513-1	462	BP2553-100	395	BP2656-25	463	BP2818-100	422	BP327-500	490	BP389-100	467	BP595-5	460	BP8014-1	403	BP8094-1	419	BP9741-10X5	473
BP2513-100	462	BP2554-25	395	BP2656-5	463	BP2818-4	422	BP328-1	489	BP391-100	467	BP596-100	431	BP8014-5	403	BP8095-1	419	BP9742-10X5	473
BP2513-250	462	BP2555-25	395	BP2656-500	463	BP2818-500	422	BP328-500	489	BP392-100	467	BP606-10	452	BP8015-1	401	BP8096-1	419	C/0080/48	72
BP2513-5	462	BP2556-200	395	BP2658-500	484	BP2819-1	495	BP329-1	491	BP393-100	467	BP606-5	452	BP8016-1	414	BP8097-1	396	C/0200/53	72
BP2516-1	454	BP2558-100	395	BP2659-40	441	BP2819-10	495	BP329-500	491	BP394-100	467	BP610-500	439	BP8017-1	411	BP8097-5	396	C/0280/48	73
BP2516-10	454	BP2558-100	395	BP2662-1	452	BP2819-100	495	BP330-1	491	BP395-100	468	BP6112	435	BP8017-5	411	BP8098-1	399	C/0480/48	73
BP2516-5	454	BP2560-4	434	BP2662-100	452	BP2819-20	495	BP330-500	491	BP396-100	468	BP6113	435	BP8018-1	415	BP8098-5	399	C/0480/53	73
BP2516-50	454	BP2561-4	434	BP2662-25	452	BP2819-4	495	BP331-1	491	BP397-100	468	BP6114	435	BP8019-1	411	BP8099-1	399	C/0520/48	73
BP2519-10	463	BP2562-4	434	BP2662-5	452	BP2820-1	495	BP331-500	491	BP399-1	469	BP6115	435	BP8020-1	411	BP8099-15	399	C/0640/46	82
BP2519-100	463	BP2563-4	434	BP2663-500	468	BP2820-100	495	BP332-1	491	BP399-20	469	BP6116	435	BP8021-1	410	BP8099-5	399	C/0640/48	82
BP2519-25	463	BP2564-1	434	BP2664-100	440	BP2820-500	495	BP332-500	491	BP399-4	469	BP6117	435	BP8022-1	413	BP8100-1	398	C/0680/46	81
BP2519-5	463	BP2564-4	434	BP2664-25	440	BP2900-1	429	BP333-1	425	BP399-500	469	BP6118	435	BP8022-5	413	BP8100-5	398	C/0680/48	81
BP2520-1	440	BP2565-1	436	BP2664-5	440	BP2900-500	429	BP333-500	425	BP410-1	484	BP612-1	439	BP8023-1	414	BP8101-1	399	C/0680/50	81
BP2520-10	440	BP2565-2	436	BP2665-10	468	BP2920-250	486	BP334-500	425	BP410-212	484	BP615-1	460	BP8024-1	405	BP8101-5	399	C/0680/60	81
BP2520-25	440	bp2565-2K	435	BP2665-100	468	BP2926-100	487	BP336-500	471	BP410-500	484	BP616-500	440	BP8024-5	405	BP8102-1	400	C/0682/48	81
BP2520-5	440	BP2570-100	392	BP2666-250	440	BP2929-25	480	BP337-100	450	BP413-25	459	BP618-500	440	BP8025-1	402	BP8102-5	400	C/0800/48	73
BP2520-50	440	BP2571-100	392	BP2667-100	456	BP2930-50	483	BP337-500	450	BP419-10	459	BP622-100	438	BP8026-1	406	BP8103-1	400	C/0840/46	74
BP2521-1	463	BP2572-100	392	BP2667-25	456	BP2933-100	482	BP338-500	450	BP427-1	439	BP629-1	452	BP8027-1	416	BP8103-5	400	C/0840/48	74
BP2521-10	463	BP2573-100	393	BP2667-50	456	BP2938-10	487	BP339-500	460	BP431-100	441	BP630-5	461	BP8029-1	416	BP8104-1	397	C/0840/50	74
BP2521-100	463	BP2574-100	392	BP2668-1	456	BP2940-4	488	BP345-500	451	BP431-500	441	BP633-5	386	BP8030-1	414	BP8104-5	397	C/0885/79	74
BP2521-25	463	BP2575-100	393	BP2668-25	456	BP2943-100	480	BP346-500	460	BP435-1	451	BP635-100	450	BP8032-1	409	BP8105-1	399	C/0920/50	74
BP2521-5	463	BP2577-100	393	BP2668-5	456	BP2944-100	487	BP349-100	481	BP435-25	451	BP635-50	450	BP8034-1	408	BP8105-5	399	C/0920/53	74
BP2521-50	463	BP2578-100	393	BP2669-25	457	BP2946-25	486	BP350-1	462	BP436-100	421	BP635-500	450	BP8035-1	402	BP8106-1	397	C/0960/50	74
BP2522-1	463	BP2578-100	394	BP2669-5	457	BP2949-5	457	BP350-500	462	BP439-500	441	BP637-1	447	BP8036-1	404	BP8106-5	397	C/1000/53	74
BP2522-100	463	BP2579-100	394	BP2670-1	469	BP2950-1	454	BP351-500	465	BP441-500	481	BP637-5	447	BP8037-1	415	BP8107-1	397	C/1040/60	75
BP2522-500	463	BP2580-100	394	BP2670-100	469	BP2950-5	454	BP352-500	471	BP458-4	481	BP639-1	388	BP8038-1	405	BP8108-1	397	C/1040/60	234
BP2523-100	463	BP2581-200	394	BP2670-500	469	BP2951-1	454	BP353-500	470	BP501-100	479	BP645-1	387	BP8039-1	402	BP8109-1	396	C/1040/63	75
BP2523-25	463	BP2582-200	394	BP2671-10	452	BP2952-1MU	455												



Mfr. No	Page	Mfr. No	Page	Mfr. No	Page	Mfr. No	Page	Mfr. No	Page	Mfr. No	Page	Mfr. No	Page	Mfr. No	Page	Mfr. No	Page
C/1400/53	75	C/4010/60	82	C/4966/PB17	310	C/7050/53	95	C/8921/17	100	D/0700/50	235	D/1856/15	310	D/2506/15	310	D/4120/25	116
C/1400/60	75	C/4020/60	82	C/4970/PB17	85	C/7120/48	95	C/8921/21	100	D/0700/53	126	D/1856/17	106	D/2506/17	109	D/4120/PB08	116
C/1400/62	75	C/4040/53	82	C/4990/08	87	C/7120/50	95	C/8921/25	100	D/0700/53	235	D/1856/17	310	D/2506/17	310	D/4120/PB17	116
C/1500/50	75	C/4040/60	82	C/5040/PB08	87	C/7120/53	95	C/8921/27	100	D/0700/60	126	D/1856/21	106	D/2506/PB17	109	D/4121/21	115
C/1500/53	75	C/4040/65	82	C/5360/48	87	C/7160/53	95	C/8921/PB17	100	D/0700/60	235	D/1856/21	310	D/2506/PB17	310	D/4121/24	115
C/1500/60	75	C/4040/70	82	C/5480/48	88	C/7160/60	95	C/8933/15	100	D/0700/65	126	D/1856/25RSS	106	D/2950/08	109	D/4121/25	115
C/1500/65	75	C/4080/50	82	C/5480/53	88	C/7170/50	95	C/8933/17	100	D/0700/68	126	D/1856/25RSS	310	D/2950/17	109	D/4121/PB08	115
C/1500/70	75	C/4120/53	83	C/5560/53	89	C/7210/53	95	C/8936/10RSS	99	D/0700/68	235	D/1856/30RSS	106	D/2950/25	109	D/4121/PB15	115
C/1610/53	75	C/4120/60	83	C/5680/53	89	C/7240/46	98	C/8936/10RSS	310	D/0700/70	126	D/1856/30RSS	310	D/2951/08	108	D/4121/PB17	115
C/1610/53	234	C/4280/53	83	C/5720/53	90	C/7320/PB08	98	C/8936/15	99	D/0700/70	235	D/1856/PB17	106	D/2951/17	108	D/4125/10RSS	115
C/1620/48	76	C/4280/53	234	C/5840/53	89	C/7360/25	98	C/8936/15	310	D/0750/50	127	D/1856/PB17	310	D/3000/46	111	D/4125/10RSS	310
C/1640/60	76	C/4280/68	83	C/5840/62	89	C/7360/PB08	98	C/8936/17	99	D/0755/46	103	D/1857/10RSS	106	D/3100/45	111	D/4125/25RSS	115
C/1640/65	76	C/4280/68	234	C/5880/50	89	C/7360/PB17	98	C/8936/17	310	D/0760/46	103	D/1857/10RSS	310	D/3200/04	111	D/4125/25RSS	310
C/1642/50	76	C/4320/50	83	C/5880/53	89	C/7440/44	99	C/8936/25RSS	99	D/0765/PB07	103	D/1857/15	106	D/3220/46	111	D/4125/30RSS	115
C/1642/53	76	C/4320/50	234	C/5920/50	89	C/7450L/05	98	C/8936/25RSS	310	D/0800/48	103	D/1857/15	310	D/3280/PB08	111	D/4125/30RSS	310
C/1642/60	76	C/4320/53	83	C/5920/53	89	C/7460/44	98	C/8936/30RSS	99	D/1250/PB08	103	D/1857/17	106	D/3400/08	112	D/4125/PB07	115
C/1720/53	76	C/4320/53	234	C/6002/46	90	C/7460/44	158	C/8936/30RSS	310	D/1338/PB07	103	D/1857/17	310	D/3410/08	112	D/4125/PB07	310
C/1720/60	76	C/4322/47	83	C/6098/08	90	C/7460L/08	99	C/8936/PB17	99	D/1338/PB07	310	D/1857/25RSS	106	D/3410/17	112	D/4125/PB08	115
C/1820/53	76	C/4345/20	83	C/6100/48	90	C/7480/43	99	C/8936/PB17	310	D/1340/PB08	104	D/1857/25RSS	310	D/3411/08	112	D/4125/PB08	310
C/1820/53	234	C/4360/50	84	C/6160/53	91	C/7480/44	99	C/9000/PB08	100	D/1341/05	104	D/1857/30RSS	106	D/3411/17	112	D/4125/PB17	115
C/1820/65	76	C/4360/53	84	C/6160/60	91	C/7551/45	99	C/9000/PB17	100	D/1341/05	310	D/1857/30RSS	310	D/3411/PB17	112	D/4125/PB17	310
C/1820/65	234	C/4400/53	84	C/6160/62	91	C/7720/48	99	C/9040/21	101	D/1345/53	104	D/1857/PB17	106	D/3450/08	112	D/4127/PB17	115
C/1880/53	77	C/4400/70	84	C/6160/65	91	C/7760/50	96	C/9040/25	101	D/1350/17	104	D/1857/PB17	310	D/3550/46	112	D/4128/PB17	115
C/1880/60	77	C/4600/48	84	C/6160/70	91	C/7760/53	96	C/9040/PB08	101	D/1350/25	104	D/1858/15	107	D/3550/48	112	D/4140/05	116
C/1880/63	77	C/4640/43	84	C/6180/53	90	C/7840/50	96	C/9040/PB17	101	D/1350/27	104	D/1858/17	107	D/3600/46	112	D/4500/08	116
C/1882/53	76	C/4680/15	84	C/6200/53	90	C/7879/53	96	C/9050/PB08	100	D/1500/08	104	D/1858/PB17	107	D/3600/48	112	D/4500/21	116
C/1882/60	76	C/4680/17	84	C/6200/53	234	C/7879/63	96	C/9050/PB17	100	D/1500/17	104	D/1859/15	106	D/3650/44	112	D/4500/25	116
C/1960/53	77	C/4680/25	84	C/6200/53	234	C/7880/50	96	C/9080/07	101	D/1500/27	104	D/1859/15	310	D/3802/PB05	113	D/4500/27	116
C/2100/60	77	C/4680/PB17	84	C/6200/60	90	C/7880/53	96	C/9080/08	101	D/1550/PB05	104	D/1859/17	106	D/3802/PB08	113	D/4500/PB17	116
C/2120/60	77	C/4681/15	84	C/6200/60	234	C/7880/70	96	C/9150/46	101	D/1550/PB15	104	D/1859/17	310	D/3820/PB05	113	D/4550/08	116
C/2160/50	77	C/4681/17	84	C/6200/60	234	C/7920/48	96	C/9152/48	101	D/1600/15	105	D/1860/PC21	108	D/3820/PB07	113	D/4550/15	116
C/2240/53	77	C/4681/PB17	84	C/6200/62	90	C/7920/50	96	C/9160/46	101	D/1600/17	105	D/1861/PB17	106	D/3830/PB08	113	D/4550/25	116
C/2280/63	78	C/4750/05	85	C/6200/62	234	C/7920/53	96	C/P330/53	90	D/1600/PB17	105	D/2050/08	108	D/3840/08	114	D/4550/PB15	116
C/2320/53	78	C/4756/15	85	C/6200/62	234	C/7920/65	96	C/P330/60	90	D/1610/08	104	D/2050/17	108	D/3840/15	114	D/4550/PB17	116
C/2360/50	78	C/4756/15	310	C/6200/65	90	C/7960/48	96	C/P540/46	94	D/1750/08	105	D/2100/08	108	D/3840/17	114	D/4551/PB15	116
C/2360/53	78	C/4756/17	85	C/6200/65	234	C/7960/60	96	C/P541/46	94	D/1750/17	105	D/2100/17	108	D/3840/21	114	D/4551/PB17	116
C/2360/60	78	C/4756/17	310	C/6200/65	234	C/8080/50	97	D/0008/44	101	D/1750/21	105	D/2110/08	108	D/3840/25	114	D/4556/PB15	116
C/2400/53	77	C/4758/17	85	C/6200/70	90	C/8080/53	97	D/0020/05	102	D/1750/25	105	D/2110/17	108	D/3840/27	114	D/4556/PB15	311
C/2440/53	77	C/4759/05	85	C/6200/70	234	C/8200/50	97	D/0020/08	102	D/1751/15	105	D/2250/48	108	D/3840/PB17	114	D/4556/PB17	116
C/2440/60	77	C/4880/PB08	85	C/6200/70	234	C/8200/53	97	D/0030/46	102	D/1751/17	105	D/2350/08	109	D/3841/08	114	D/4556/PB17	311
C/2520/53	78	C/4920/08	86	C/6230/53	90	C/8320/48	97	D/0035/44	102	D/1751/27	105	D/2400/17	110	D/3841/15	114	D/4600/48	117
C/2560/50	78	C/4920/15	86	C/6230/53	312	C/8320/53	97	D/0050/17	101	D/1756/15	105	D/2400/25	110	D/3841/17	114	D/4600/48	158
C/2640/48	78	C/4920/17	86	C/6240/53	90	C/8360/50	97	D/0052/17	102	D/1756/15	310	D/2400/27	110	D/3841/21	114	D/4700/48	117
C/2880/PB08	79	C/4920/21	86	C/6240/53	234	C/8480/53	98	D/0100/50	102	D/1756/17	105	D/2400/MS21	110	D/3841/24	114	D/4700/50	117
C/2880/PB17	79	C/4920/25	86	C/6280/05	91	C/8480/60	98	D/0122/48	102	D/1756/17	310	D/2400/PB08	110	D/3841/25	114	D/4750/46	117
C/2920/PB08	79	C/4920/27	86	C/6360/48	92	C/8480/63	98	D/0124/48	102	D/1820/44	105	D/2400/PB17	110	D/3841/27	114	D/4850/45	117
C/2920/PB17	79	C/4920/PB17	86	C/6400/50	92	C/8480/65	98	D/0130/48	102	D/1820/44	158	D/2450/15	110	D/3841/PB17	114	D/4850/45	158
C/3090/53	79	C/4920/PC21	86	C/6459/48	92	C/8520/50	98	D/0150/08	102	D/1850/08	107	D/2450/15	235	D/3842/15	114	D/4865/48	117
C/3091/53	79	C/4960/15	86	C/6500/48	93	C/8520/53	98	D/0250/05	103	D/1850/15	107	D/2450/17	110	D/3842/17	114	D/4865/65	117
C/3092/53	79	C/4960/17	86	C/6500/53	93	C/8520/60	98	D/0370/25	103	D/1850/17	107	D/2450/17	235	D/3842/PB17	114	D/4865/70	117
C/3093/53	79	C/4960/25	86	C/6520/48	93	C/8520/63	98	D/0370/PB08	103	D/1850/21	107	D/2450/25	110	D/3846/10RSS	113	D/5000/15	111
C/3100/53	79	C/4960/27	86	C/6520/53	93	C/8520/65	98	D/0370/PB17	103	D/1850/25	107	D/2450/25	235	D/3846/10RSS	310	D/5000/17	111
C/3160/44	79	C/4960/PB08	86	C/6520/60	93	C/8560/50	98	D/0450/50	126	D/1850/27	107	D/2450/27	110	D/3846/15	113	D/5000/21	111
C/3240/53	79	C/4960/PB15	86	C/6560/48	92	C/8560/53	98	D/0450/53	126	D/1850/MC15	107	D/2450/27	235	D/3846/15	310	D/5000/FP21	111
C/3280/08	80	C/4960/PB17	86	C/6560/53	92	C/8560/60	98	D/0450/68	126	D/1850/PB17	107	D/2450/MS21	110	D/3846/17	113	D/5000/PB17	111
C/3280/17	80	C/4961/17	86	C/6560/60	92	C/8560/63	98	D/0452/48	126	D/1850/PC21	107	D/2450/MS21	235	D/3846/17	310	D/5002/15	111
C/3320/48	80	C/4962/15	86	C/6560/65	92	C/8560/65	98	D/0452/50	126	D/1852/15	107	D/2450/PB08	110	D/3846/25RSS	113	D/5002/17	111
C/3320/53	80	C/4962/17	86	C/6600/48	92	C/8560/70	98	D/0500/48	126	D/1852/17	107	D/2450/PB08	235	D/3846/25RSS	310	D/5002/21	111
C/3360/44	80	C/4963/15	86	C/6640/48	93	C/8600/50	97	D/0500/70	126	D/1852/21	107	D/2450/PB08	235	D/3846/30RSS	113	D/5002/25	111
C/3400/08	80	C/4963/17	86	C/6640/53	93	C/8600/60	97	D/0600/53	127	D/1852/24	107	D/2450/PB15	110	D/3846/30RSS	310	D/5010/17	111
C/3440/05	80	C/4966/10RSS	85	C/6680/48	93	C/8680/48	96	D/0600/70	127	D/1852/25	107	D/2450/PB17	110	D/3846/PB17	113	D/5051/44	117
C/3680/48	81	C/4966/10RSS	310	C/6680/50	93	C/8800/48	97	D/0650/48	126	D/1852/27	107	D/2455/17	109	D/3846/PB17	310	D/5150/46	117
C/3720/48	81	C/4966/15	85	C/6840/50	94	C/8840/50	97	D/0650/50	126	D/1852/PB15	107	D/2500/08	110	D/3847/PB17	113	D/5290/05	118
C/3780/48	80	C/4966/15</															



Mfr. No	Page
E/0650/PC17	121
E/0650DF/08	121
E/0650DF/08	235
E/0650DF/15	121
E/0650DF/15	235
E/0650DF/17	121
E/0650DF/17	235
E/0650DF/21	121
E/0650DF/21	235
E/0650DF/24	121
E/0650DF/24	235
E/0650DF/25	121
E/0650DF/25	235
E/0650DF/28	235
E/0650DF/C17	121
E/0650DF/C17	235
E/0650DF/P17	121
E/0655/17	120
E/0655DF/17	120
E/0660/08	120
E/0665/15	120
E/0665/15	311
E/0665/17	120
E/0665/17	311
E/0665DF/10RS	120
E/0665DF/10RS	311
E/0665DF/15	120
E/0665DF/15	311
E/0665DF/17	120
E/0665DF/17	311
E/0665DF/25RS	120
E/0665DF/25RS	311
E/0665DF/30RS	120
E/0665DF/30RS	311
E/0700/08	123
E/0700/17	123
E/0700/25	123
E/0700/PB17	123
E/0701/08	123
E/0701/17	123
E/0750/15	123
E/0750/17	123
E/0750/25	123
E/0751/08	123
E/0800/PB17	124
E/0850/08	125
E/0850/15	125
E/0850/17	125
E/0850/21	125
E/0850/24	125
E/0850/25	125
E/0850/27	125
E/0850/FP21	125
E/0850/MC15	125
E/0850/PB17	125
E/0850/PC17	125
E/0900/15	124
E/0900/17	124
E/0900/21	124
E/0900/25	124
E/0900/27	124
E/0900/FP21	124
E/0900/PB17	124
E/0900/PC17	124
E/0903/15	124
E/0903/17	124
E/0906/10RSS	124
E/0906/10RSS	311
E/0906/15	124
E/0906/15	311
E/0906/17	124

Mfr. No	Page
E/0906/17	311
E/0906/25RSS	124
E/0906/25RSS	311
E/0906/30RSS	124
E/0906/30RSS	311
E/0906/PB17	124
E/0906/PB17	311
E/0907/17	124
E/0907/17	311
E/0908/17	124
E/0950/08	125
E/0950/17	125
E/1000/08	125
E/1100/08	125
E/1250/04	125
E/1300/17	109
E/1400/17	128
E/1450/08	128
E/1450/15	128
E/1450/17	128
E/1450/21	128
E/1450/25	128
E/1450/27	128
E/1450/PB17	128
E/1455/15	128
E/1455/17	128
E/1455/25	128
E/1455/27	128
E/1455/PB15	128
E/1455/PB17	128
E/1460/08	128
E/P140/48	126
E/P140/53	126
E/P140/65	126
E/P700/44	123
E/P800/03	123
E/P800/10	123
F/1170/50	129
F/1173/50	129
F/1200/46	130
F/1200/48	130
F/1250/44	130
F/1300/48	130
F/1300/48	159
F/1300/53	130
F/1300/53	159
F/1325/21	351
F/1350/08	130
F/1400/08	130
F/1451/25	131
F/1451/PB17	131
F/1501/25	130
F/1501/PB08	130
F/1501/PB15	130
F/1501/PB17	130
F/1510/21	131
F/1515/25	131
F/1520/21	131
F/1550/PB08	132
F/1550/PB17	132
F/1551/25	131
F/1551/27	131
F/1551/PB08	131
F/1551/PB17	131
F/1750/25	133
F/1750/PB15	133
F/1750/PB17	133
F/1800/27	133
F/1800/PB08	133
F/1800/PB17	133
F/1820/PB08	132

Mfr. No	Page
F/1820/PB17	132
F/1850/27	132
F/1850/PB08	132
F/1850/PB17	132
F/1900/25	132
F/1900/PB08	132
F/1900/PB15	132
F/1900/PB17	132
F/1950/48	133
F/1950/50	133
F/1952/48	133
F/1952/53	133
F/1952/60	133
F/2070/53	133
F/2080/PB17	134
F/2250/53	134
F/2250/60	134
F/P025/46	129
F/P565/PB05	131
F/P565/PB08	131
F/P770/46	133
FB6000-10	436
FB6000-15	436
FB6000-20	436
FB6000-25	436
FB6000-30	436
FB6000-35	436
FB6000-40	436
FB6000-45	436
FB6000-50	436
FB6000-55	436
FB6000-60	436
FB6000-65	436
FB6000-70	436
FB6000-75	436
FB6100	435
FB6111	435
FB6210	435
FB6225	435
FB6245	435
FB6260	435
FE/0100/17	309
FE/0100/17	370
FE/0200/17	22
FE/0200/17	309
FE/1100/17	22
FE/1100/17	309
FE/1200/17	309
FE/1200/17	370
FE/1300/17	202
FE/1300/17	309
G/0050/48	134
G/0150/50	134
G/0150/53	134
G/0150/53	235
G/0200/43	135
G/0275/46	135
G/0300/50	135
G/0300/53	135
G/0400/60	135
G/0400/60	235
G/0400/65	135
G/0400/65	235
G/0450/53	135
G/0450/60	135
G/0450/61	135
G/0450/65	135
G/0500/53	135
G/0500/60	135
G/0500/60	235

Mfr. No	Page
G/0500/61	135
G/0500/61	235
G/0500/65	135
G/0500/65	235
G/0500/70	135
G/0500/70	235
G/0505/55	135
G/0513/50	136
G/0514/50	136
G/0515/48	135
G/0518/PB08	136
G/0600/08	136
G/0600/17	136
G/0600/25	136
G/0650/08	136
G/0650/15	136
G/0650/17	136
G/0650/21	136
G/0650/25	136
G/0650/27	136
G/0800/48	137
G/0800/50	137
G/0800/60	137
G/0803/48	137
G/0850/90	137
G/0900/60	137
G/0910/07	137
G/0950/48	137
G/1050/53	138
G/1100/48	138
G/P450/05	136
G/P450/08	136
G/P460/53	136
G/P771/48	138
H/0010/46	138
H/0050/PB15	138
H/0095/45	138
H/0100/08	140
H/0100/15	140
H/0100/17	140
H/0100/27	140
H/0106/10RSS	138
H/0106/10RSS	311
H/0106/15	138
H/0106/15	311
H/0106/17	138
H/0106/17	311
H/0106/25RSS	138
H/0106/25RSS	311
H/0106/30RSS	138
H/0106/30RSS	311
H/0106/PB17	138
H/0106/PB17	311
H/0107/15	139
H/0107/15	312
H/0107/17	139
H/0107/17	312
H/0110/17	139
H/0155/08	139
H/0155/15	139
H/0155/17	139
H/0155/21	139
H/0155/25	139
H/0155/27	139
H/0155/FP21	139
H/0155/PB17	139
H/0160/08	139
H/0160/15	139
H/0160/17	139
H/0160/21	139
H/0160/27	139

Mfr. No	Page
H/0160/FP21	139
H/0166/07	140
H/0168/46	140
H/0168/48	140
H/0168/65	140
H/0169/46	140
H/0250/05	140
H/0262/43	140
H/0262/44	140
H/0270/05	140
H/0285/04	141
H/0300/53	141
H/0300/60	141
H/0300/65	141
H/0300/70	141
H/0320/53	141
H/0320/60	141
H/0320/65	141
H/0350/08	143
H/0350/15	143
H/0350/17	143
H/0350/21	143
H/0350/25	143
H/0350/27	143
H/0350/FP21	143
H/0350/PB17	143
H/0355/15	143
H/0355/17	143
H/0355/21	143
H/0355/25	143
H/0355/27	143
H/0400/08	143
H/0400/17	143
H/0401/15	167
H/0401/21	167
H/0401/24	167
H/0401/25	167
H/0401/27	167
H/0401/PB17	167
H/0402/PB15	143
H/0402/PB17	143
H/0403/15	142
H/0403/17	142
H/0405/10RSS	166
H/0405/10RSS	311
H/0405/21	166
H/0405/21	311
H/0405/25	166
H/0405/25	311
H/0405/25RSS	166
H/0405/25RSS	311
H/0405/30RSS	166
H/0405/30RSS	311
H/0405/PB15	166
H/0405/PB15	311
H/0405/PB17	166
H/0405/PB17	311
H/0405/26	149
H/0406/10RSS	149
H/0406/10RSS	311
H/0406/15	142
H/0406/15	311
H/0406/17	142
H/0406/17	311
H/0406/21	149
H/0406/21	311
H/0406/25RSS	142
H/0406/25RSS	311
H/0406/30RSS	142
H/0406/30RSS	311
H/0406/PB17	142
H/0406/PB17	311

Mfr. No	Page
H/0408/15	166
H/0408/17	166
H/0409/PB15	142
H/0409/PB15	312
H/0409/PB17	142
H/0409/PB17	312
H/0410/PB15	166
H/0410/PB15	312
H/0410/PB17	166
H/0410/PB17	312
H/0411/21	166
H/0411/27	166
H/0411/PB15	166
H/0411/PB17	166
H/0420/21	142
H/0420/25	142
H/0420/PB08	142
H/0420/PB17	142
H/0421/25	141
H/0421/27	141
H/0421/PB08	141
H/0421/PB17	141
H/0438/46	144
H/0438/48	144
H/0439/46	144
H/0440/07	144
H/0450/25	144
H/0450/27	144
H/0450/PB08	144
H/0450/PB17	144
H/0460/05	144
H/0467/48	144
H/0467/63	144
H/0480/17	144
H/0500/PB08	145
H/0550/PB07	145
H/0550/PB17	145
H/0560/PB05	145
H/0560/PB08	145
H/0564/21	339
H/0600/50	145
H/0650/53	146
H/0650/70	146
H/0700/48	146
H/0700/50	146
H/0750/05	146
H/0850/07	146
H/0900/27	147
H/0900/PB08	147
H/0900/PB17	147
H/0950/25	147
H/0950/PB08	147
H/0950/PB17	147
H/1000/PB17	149
H/1020/PB15	150
H/1020/PB17	150
H/1050/26	149
H/1050/PB17	149
H/1100/26	149
H/1100/PB15	149
H/1100/PB17	149
H/1111/25	148
H/1111/PB17	148
H/1150/21	149
H/1150/25	149
H/1150/26	149
H/1150/PB15	149
H/1150/PB17	149
H/1180/PB08	148
H/1180/PB15	148
H/1180/PB17	148

Mfr. No	Page
H/1185/PB17	148
H/1190/PB08	150
H/1190/PB17	150
H/1196/PB07	25
H/1196/PB07	148
H/1196/PB08	25
H/1196/PB08	148
H/1196/PB15	25
H/1196/PB15	148
H/1196/PB17	25
H/1196/PB17	148
H/1200/25	149
H/1200/PB08	149
H/1200/PB15	149
H/1200/PB17	149
H/1200/PC15	149
H/1200/PC17	149
H/1300/15	151
H/1300/17	151
H/1350/15	153
H/1350/17	153
H/1400/15	153
H/1400/17	153
H/1420/15	152
H/1420/17	152
H/1425/17	152
H/1430/90	153
H/1550/15	155
H/1550/17	155
H/1600/15	155
H/1600/17	155
H/1650/17	155
H/1750/15	155
H/1750/17	155
H/1800/15	154
H/1800/17	154
H/1800/25	154
H/1811/17	154
H/1811/21	154
H/1820/15	154
H/1860/15	153
H/1862/15	153
H/1900/48	156

Mfr. No	Page	Mfr. No	Page	Mfr. No	Page	Mfr. No	Page	Mfr. No	Page	Mfr. No	Page	Mfr. No	Page	Mfr. No	Page	Mfr. No	Page
J/2880/17	66	J/4200/PB08	133	J/4400/PB08	323	J/5655/PB08	311	J/7050/PB17	324	J/7610/21	237	J/7800/21	324	J/8030/05	162	J/8120/08	363
J/2880/21	63	J/4200/PB08	311	J/4400/PB15	161	J/5820/PB15	229	J/7300/PB08	272	J/7610/21	296	J/7800/24	296	J/8030/05	314	J/8150/05	314
J/2880/21	66	J/4250/17	150	J/4400/PB15	235	J/5820/PB15	323	J/7300/PB15	272	J/7610/21	324	J/7800/24	324	J/8030/08	162	J/8150/05	366
J/2880C/90	67	J/4270/17	150	J/4400/PB15	323	J/5820/PB17	229	J/7310/15	278	J/7610/PB15	237	J/7800/25	296	J/8030/08	314	J/8150/08	314
J/2885/08	63	J/4270/25	150	J/4400/PB17	161	J/5820/PB17	323	J/7310/17	278	J/7610/PB15	296	J/7800/25	324	J/8035/05	177	J/8150/08	366
J/2885/08	67	J/4310/15	151	J/4400/PB17	235	J/5900/PB08	240	J/7320/PB15	277	J/7610/PB15	324	J/7810/17	296	J/8035/05	314	J/8166/05	331
J/2885/15	63	J/4310/15	323	J/4400/PB17	323	J/6100/15	247	J/7320/PB15	324	J/7620/15	237	J/7810/25	296	J/8035/08	177	J/8166/08	314
J/2885/15	67	J/4310/17	151	J/4410/15	161	J/6100/15	323	J/7320/PB17	277	J/7620/15	296	J/7900/15	300	J/8035/08	314	J/8166/08	331
J/2885/17	63	J/4310/17	323	J/4410/15	235	J/6100/17	247	J/7320/PB17	324	J/7620/15	324	J/7900/15	325	J/8037/05	182	J/8200/08	314
J/2885/17	67	J/4310/21	151	J/4410/15	323	J/6100/17	323	J/7330/15	236	J/7620/17	237	J/7950/15	237	J/8037/05	314	J/8200/08	377
J/2885/21	63	J/4315/15	151	J/4410/17	161	J/6100C/05	247	J/7330/15	277	J/7620/17	296	J/7950/15	306	J/8037/08	182	J/8205/05	27
J/2885/21	67	J/4315/15	235	J/4410/17	235	J/6100C/05	326	J/7330/15	324	J/7620/17	324	J/7950/15	325	J/8037/08	314	J/8205/05	315
J/2885/24	63	J/4315/15	323	J/4410/17	323	J/6100C/90	247	J/7330/17	236	J/7620/21	237	J/7950/17	237	J/8040/05	186	J/8206/08	27
J/2885/24	67	J/4315/17	151	J/4450/15	161	J/6100C/90	326	J/7330/17	277	J/7620/21	296	J/7950/17	306	J/8040/05	314	J/8206/08	315
J/2885C/05	64	J/4315/17	235	J/4450/17	161	J/6212/15	236	J/7330/17	324	J/7620/21	324	J/7950/17	325	J/8040/08	186	J/8210/05	49
J/2885C/05	64	J/4315/17	323	J/4500/PB15	161	J/6212/15	251	J/7330/24	236	J/7620/24	237	J/7950/21	237	J/8040/08	314	J/8210/05	315
J/2885C/05	67	J/4315/24	151	J/4500/PB17	161	J/6212/15	323	J/7330/24	277	J/7620/24	296	J/7950/21	306	J/8045/05	191	J/8215/08	49
J/2885C/90	64	J/4315/24	235	J/4534/05	48	J/6212/17	236	J/7330/24	324	J/7620/24	324	J/7950/21	325	J/8045/05	314	J/8215/08	315
J/2885C/90	64	J/4320/15	151	J/4542/05	36	J/6212/17	251	J/7330/25	277	J/7620C/05	296	J/7950/24	237	J/8045/08	191	J/8220/05	52
J/2885C/90	67	J/4320/15	235	J/4544/05	60	J/6212/17	323	J/7330/25	324	J/7620C/05	326	J/7950/24	306	J/8045/08	314	J/8220/05	315
J/2892/15	249	J/4320/15	323	J/4546/05	83	J/6212C/05	252	J/7330C/05	278	J/7620C/90	296	J/7950/24	325	J/8047/05	194	J/8221/08	52
J/2980/05	78	J/4320/17	151	J/4548/05	130	J/6212C/05	326	J/7330C/05	326	J/7620C/90	326	J/7950C/05	306	J/8047/05	314	J/8221/08	315
J/2983/05	249	J/4320/17	235	J/4548/08	130	J/6212C/90	252	J/7330C/90	278	J/7625/17	296	J/7950C/05	326	J/8047/08	194	J/8230/05	59
J/2984/05	249	J/4320/17	323	J/4550/05	183	J/6212C/90	326	J/7330C/90	326	J/7625/17	325	J/7950C/90	306	J/8047/08	314	J/8230/05	315
J/2986/05	249	J/4320/21	151	J/4552/05	48	J/6220/15	252	J/7340/15	278	J/7625/PB15	296	J/7950C/90	326	J/8050/05	209	J/8231/08	59
J/3000/15	81	J/4320/21	235	J/4554/05	80	J/6220/17	252	J/7340/15	324	J/7625/PB15	325	J/7951C/05	306	J/8050/05	314	J/8231/08	315
J/3000/15	234	J/4320/21	323	J/4556/05	215	J/6610/15	257	J/7340/17	278	J/7630/15	296	J/7951C/05	326	J/8050/08	209	J/8235/05	72
J/3000/15	323	J/4320/24	151	J/4558/05	242	J/6610/15	324	J/7340/17	324	J/7630/15	325	J/7951C/90	306	J/8050/08	314	J/8235/05	315
J/3000/17	81	J/4320/24	235	J/4560/05	245	J/6610/17	257	J/7350/15	278	J/7630/17	296	J/7951C/90	326	J/8055/05	212	J/8236/08	71
J/3000/17	234	J/4320/24	323	J/4562/05	279	J/6610/17	324	J/7350/15	324	J/7630/17	325	J/8000/05	27	J/8055/05	314	J/8236/08	315
J/3000/17	323	J/4320C/05	151	J/4564/05	334	J/6612/15	257	J/7350/17	278	J/7630/21	296	J/8000/05	314	J/8055/08	212	J/8240/05	73
J/3000/24	81	J/4320C/05	326	J/4600/PB08	174	J/6612/15	324	J/7350/17	324	J/7630/21	325	J/8000/08	27	J/8055/08	314	J/8240/05	315
J/3010C/90	81	J/4320C/90	151	J/4600/PB15	174	J/6612/17	257	J/7370/07	278	J/7630/24	296	J/8000/08	314	J/8060/05	245	J/8241/08	73
J/3010C/90	326	J/4320C/90	326	J/4600/PB17	174	J/6612/17	324	J/7370/07	324	J/7630/24	325	J/8002/05	49	J/8060/05	314	J/8241/08	315
J/3500/PC15	87	J/4330/15	150	J/4620/PB08	170	J/6620/15	256	J/7370/15	278	J/7630/27	296	J/8002/05	314	J/8060/08	245	J/8245/05	88
J/3700/15	127	J/4330/15	323	J/4620/PB08	174	J/6620/15	324	J/7370/15	324	J/7630/27	325	J/8002/08	49	J/8060/08	314	J/8245/05	315
J/3700/15	235	J/4330/17	150	J/4620/PB15	170	J/6620/17	256	J/7370/17	278	J/7630C/05	296	J/8002/08	314	J/8061/05	273	J/8246/08	87
J/3700/15	323	J/4330/17	323	J/4620/PB15	174	J/6620/17	324	J/7370/17	324	J/7630C/05	326	J/8003/05	50	J/8061/05	314	J/8246/08	315
J/3700/17	127	J/4330/25	150	J/4620/PB17	170	J/6625/15	236	J/7380/PB08	277	J/7630C/90	296	J/8003/05	314	J/8061/08	273	J/8250/05	91
J/3700/17	235	J/4330/25	323	J/4620/PB17	174	J/6625/15	256	J/7400/PB17	279	J/7630C/90	326	J/8003/08	50	J/8061/08	314	J/8250/05	315
J/3700/17	323	J/4330C/05	151	J/5000/21	182	J/6625/15	323	J/7400/PC15	279	J/7640/15	297	J/8003/08	314	J/8062/05	275	J/8251/08	91
J/3700/24	127	J/4330C/05	326	J/5100/08	198	J/6630/17	236	J/7410/PB17	279	J/7640/15	325	J/8004/05	52	J/8062/05	314	J/8251/08	315
J/3700/24	235	J/4330C/90	151	J/5300/PB05	208	J/6630/17	257	J/7440C/05	281	J/7640/17	297	J/8004/05	314	J/8062/08	275	J/8255/05	94
J/3700/24	323	J/4330C/90	326	J/5500/PB08	211	J/6630/17	324	J/7440C/05	326	J/7640/17	325	J/8004/08	52	J/8062/08	314	J/8255/05	315
J/3700C/05	127	J/4340/15	150	J/5500/PB15	211	J/6630/PB15	236	J/7490/17	279	J/7640/21	297	J/8004/08	314	J/8064/05	276	J/8256/08	94
J/3700C/05	326	J/4340/15	323	J/5545/PB17	218	J/6630/PB15	257	J/7500/17	279	J/7640/21	325	J/8005/05	71	J/8064/05	314	J/8256/08	315
J/3700C/90	127	J/4340/17	150	J/5545/PB17	323	J/6630/PB15	324	J/7510/PB15	282	J/7640/24	297	J/8005/05	314	J/8064/08	276	J/8270/05	162
J/3700C/90	326	J/4340/17	323	J/5550/24	218	J/6630C/05	257	J/7510/PB15	324	J/7640/24	325	J/8005/08	71	J/8064/08	314	J/8270/05	315
J/3710/17	127	J/4350/15	150	J/5550/24	236	J/6630C/05	326	J/7510/PB17	282	J/7650/PB08	293	J/8005/08	314	J/8065/05	279	J/8271/08	162
J/3710/17	312	J/4350/15	235	J/5550/24	323	J/6630C/90	257	J/7510/PB17	324	J/7650/PB08	312	J/8006/05	59	J/8065/05	314	J/8271/08	315
J/3720/17	127	J/4350/15	323	J/5550/PB15	218	J/6630C/90	326	J/7520/17	285	J/7660/15	237	J/8006/05	314	J/8065/08	279	J/8275/05	178
J/3720/17	323	J/4350/17	150	J/5550/PB15	236	J/6640C/05	257	J/7520/17	324	J/7660/15	297	J/8006/08	59	J/8065/08	314	J/8275/05	315
J/3720/24	127	J/4350/17	235	J/5550/PB15	323	J/6640C/05	326	J/7520/PB15	285	J/7660/15	325	J/8006/08	314	J/8067/05	314	J/8276/08	178
J/3720/PB15	127	J/4350/17	323	J/5550/PB17	218	J/6640C/90	257	J/7520/PB15	324	J/7660/17	237	J/8010/05	73	J/8067/05	349	J/8276/08	315
J/3720/PB15	323	J/4350/21	150	J/5550/PB17	236	J/6640C/90	326	J/7580/17	284	J/7660/17	297	J/8010/05	314	J/8067/08	314	J/8280/05	183
J/3720C/05	127	J/4350/21	235	J/5550C/05	323	J/6690/15	236	J/7580/17	324	J/7660/17	325	J/8010/08	73	J/8067/08	349	J/8280/05	315
J/3720C/05	326	J/4350/21	323	J/5550C/05	218	J/6690/15	258	J/7580/PB15	284	J/7660/21	237	J/8010/08	314	J/8070/05	314	J/8281/08	183
J/3720C/90	127	J/4350/24	150	J/5550C/05	326	J/6690/15	324	J/7580/PB15	324	J/7660/21	297	J/8015/05	87	J/8070/05	373	J/8281/08	315
J/3720C/90	326	J/4350/24	235	J/5550C/90	218	J/6690/17	236	J/7580C/05	285	J/7660/21	325	J/8015/05	314	J/8070/08	314	J/8285/05	186
J/3800/15	129	J/4350/24	323	J/5550C/90	326	J/6690/17	258	J/7580C/05	326	J/7660/24	237	J/8015/08	87	J/8070/08	373	J/8285/05	315
J/3800/17	129	J/4350/25	150	J/5600/17	218	J/6690/17	324	J/7580C/90	285	J/7660/24	297	J/8015/08	314	J/8100/05	314	J/8286/08	186
J/3900/15	129	J/4350C/05	150	J/5600/17	323	J/6700/17	258	J/7580C/90	326	J/7660/24	325	J/8020/05	91	J/8100/05	342	J/8286/08	315
J/3900/17	129	J/4350C/05	326	J/5600/24	218	J/6700/17	324	J/7590C/05	286	J/766							

Mfr. No	Page	Mfr. No	Page	Mfr. No	Page	Mfr. No	Page	Mfr. No	Page	Mfr. No	Page	Mfr. No	Page	Mfr. No	Page	Mfr. No	Page
J/8430/24	238	J/9013/08	296	K/2560/04	173	L/1700/53	181	LS120-500	21	M/1400/48	190	M/3950/FP21	202	M/4101/25	203	M/5085/PB17	206
J/8430/24	325	J/9014/08	267	K/2710/99	170	L/1720/53	181	LS120-500	308	M/1450/48	191	M/3950/PB17	202	M/4101/PB17	203	M/5120/05	208
J/8430/24	339	J/9015/08	36	K/2710/99	173	L/1730/48	180	LS121-1	22	M/1450/53	191	M/3950/PC17	202	M/4150/08	203	M/5120/08	208
J/8430/25	325	J/9016/15	145	K/2720/99	170	L/1814/46	181	LS121-1	308	M/1660/53	192	M/3950/PC25	202	M/4150/17	203	M/5120/25	208
J/8430/25	339	J/9018/17	145	K/2720/99	174	L/1815/45	182	LS121-212	22	M/1665/50	192	M/4000/15	202	M/4200/08	203	M/5125/08	208
J/8430C/05	326	J/9019/17	145	K/2730/08	170	L/1815L/PB05	182	LS121-212	308	M/1665/53	192	M/4000/17	202	M/4350/17	204	M/5125/25	208
J/8430C/05	340	J/9021/17	298	K/2730/08	173	L/1819/48	182	LS121-4	22	M/1700/50	191	M/4000/21	202	M/4350/25	204	M/5150/45	158
J/8430C/90	326	K/0045/63	169	K/2740/99	170	L/1822/46	182	LS121-4	308	M/1800/50	192	M/4000/25	202	M/4350/30	204	M/5150/45	207
J/8430C/90	340	K/0050/60	169	K/2740/99	173	L/1840/46	182	LS121-500	22	M/1800/53	192	M/4000/27	202	M/4400/17	204	M/5150/46	158
J/8440/17	325	K/0050/63	169	K/2760/45	170	L/1860/07	182	LS121-500	308	M/1850/53	192	M/4000/FP21	202	M/4400/25	204	M/5150/46	207
J/8440/17	339	K/0050/70	169	K/2760/45	173	L/1900/17	182	M/0050/50	187	M/1850/60	192	M/4000/PB08	202	M/4400/PB17	204	M/5150L/08	207
J/8440/PB15	325	K/0090/60	174	K/2770/48	170	L/1950/08	182	M/0050/53	187	M/1850/70	192	M/4000/PB15	202	M/4400/PC17	204	M/5170/44	158
J/8440/PB15	339	K/0100/53	174	K/2770/48	173	L/1950/17	182	M/0200/50	187	M/2052/53	192	M/4000/PB17	202	M/4450/15	204	M/5170/44	207
J/8450/15	325	K/0100/60	174	L/0060/43	175	L/2000/45	183	M/0200/53	187	M/2250/53	193	M/4000/PC17	202	M/4450/17	204	M/5170/48	158
J/8450/15	339	K/0120/80	175	L/0100/PB08	176	L/2050/50	183	M/0300/53	187	M/2250/60	193	M/4000/PC25	202	M/4450/25	204	M/5170/48	207
J/8450/17	325	K/0120/81	175	L/0100/PB17	176	L/2080/48	184	M/0350/53	187	M/2300/50	192	M/4020/17	201	M/4450/27	204	M/5170L/08	207
J/8450/17	339	K/0125/81	175	L/0150/PB08	176	L/2080/53	184	M/0400/53	187	M/2300/53	192	M/4050/08	201	M/4450/30	204	M/5171L/08	207
J/8450/24	325	K/0130/80	175	L/0150/PB08	235	L/2100/50	184	M/0450/53	188	M/2350/60	193	M/4050/15	201	M/4450/MC15	204	M/5200/PB08	208
J/8450/24	339	K/0150/80	175	L/0150/PB17	176	L/2105/48	184	M/0450/53	235	M/2400/53	193	M/4050/17	201	M/4450/PB17	204	M/5200/PB17	208
J/8450/25	325	K/0150/81	175	L/0150/PB17	235	L/2200/48	184	M/0500/53	188	M/2405/53	193	M/4050/25	201	M/4460/PB08	204	M/5345/45	207
J/8450/25	339	K/0160/80	175	L/0160/80	176	L/2200/53	184	M/0500/53	235	M/2405/53	236	M/4050/PB17	201	M/4460/PB17	204	M/5350/48	207
J/8460/15	238	K/0160/81	175	L/0200/60	176	L/2200/65	184	M/0500/70	188	M/2420/45	193	M/4053/15	201	M/4475/46	204	M/5370/45	158
J/8460/15	238	K/0170/81	175	L/0200/65	176	L/2200/70	184	M/0500/70	235	M/2450/60	193	M/4053/17	201	M/4480/PB08	204	M/5370/45	208
J/8460/15	325	K/0200/15	175	L/0200/70	176	L/2201/50	184	M/0515/48	187	M/2450/63	193	M/4056/10RSS	200	M/4494/08	71	M/5370/46	158
J/8460/15	340	K/0200/21	175	L/0250/53	176	L/2201/65	184	M/0550/53	188	M/2451/60	193	M/4056/10RSS	311	M/4494/17	71	M/5370/46	208
J/8460/17	238	K/0210/PB17	175	L/0250/60	176	L/2300/48	185	M/0550/60	188	M/2451/63	193	M/4056/15	200	M/4494/21	71	M/5370L/08	208
J/8460/17	238	K/0220/15	175	L/0350/53	176	L/2300/50	185	M/0550/65	188	M/2510/45	193	M/4056/15	311	M/4494/25	71	M/5400/50	203
J/8460/17	325	K/0310/21	175	L/0350/60	176	L/2300/60	185	M/0600/53	188	M/2550/48	194	M/4056/17	200	M/4494/27	71	M/5650/48	209
J/8460/17	340	K/1900/15	170	L/0451/46	177	L/2300/68	185	M/0600/53	236	M/2650/50	196	M/4056/17	311	M/4494/PB17	71	M/5715/PB08	210
J/8460/24	238	K/1900/15	171	L/0451/48	177	L/2301/53	184	M/0600/60	188	M/2652/47	196	M/4056/21	200	M/4496/10RSS	71	M/5750/PB08	210
J/8460/24	238	K/1950/15	170	L/0455/46	177	L/2320/48	185	M/0600/60	236	M/2652/50	196	M/4056/21	311	M/4496/10RSS	310	M/5760/48	210
J/8460/24	325	K/1950/15	171	L/0455/48	177	L/2330/48	185	M/0600/61	188	M/2800/48	196	M/4056/25RSS	200	M/4496/17	71	M/P200/05	194
J/8460/24	340	K/1950/17	171	L/0460/07	176	L/2350/48	185	M/0600/61	236	M/2800/50	196	M/4056/25RSS	311	M/4496/17	310	M/P205/46	205
J/8460/25	238	K/2000/15	170	L/0460/15	176	L/2350/53	185	M/0600/65	188	M/2800/53	196	M/4056/30RSS	200	M/4496/25RSS	71	M/P206/46	205
J/8460/25	325	K/2000/15	171	L/0500/46	177	L/2350/65	185	M/0600/65	236	M/2850/48	196	M/4056/30RSS	311	M/4496/25RSS	310	MS/1010/50	209
J/8460/25	340	K/2000/17	171	L/0500/48	177	L/2350/68	185	M/0602/17	188	M/2850/50	196	M/4056/PB17	200	M/4496/30RSS	71	MS/1010/53	209
J/8460C/05	326	K/2100/15	170	L/0510/48	177	L/2400/53	185	M/0650/50	188	M/2850/53	196	M/4056/PB17	311	M/4496/30RSS	310	MS/1020/53	209
J/8460C/05	340	K/2100/15	171	L/0560/53	177	L/2450/48	185	M/0650/53	188	M/2950/46	197	M/4057/PB17	200	M/4650/17	204	MS/1020/62	209
J/8460C/90	326	K/2100/17	170	L/0565/48	177	L/2500/46	158	M/0750/53	189	M/2950/48	197	M/4057/PB17	312	M/4650/25	204	MS/1030/53	209
J/8460C/90	340	K/2100/17	171	L/0565/99	177	L/2500/46	185	M/0752/53	188	M/3012/48	197	M/4058/10RSS	200	M/4650/27	204	MS/1030/62	209
J/8470/17	325	K/2110/15	170	L/0605/60	179	L/2500L/08	185	M/0752/60	188	M/3050/48	197	M/4058/10RSS	309	M/4720/15	205	MS/1070/53	209
J/8470/17	340	K/2110/15	172	L/0610/53	178	L/2500L/17	185	M/0800/50	189	M/3100/48	197	M/4058/15	200	M/4850/48	205	N/0050/25	210
J/8470/24	325	K/2110/17	170	L/0700/53	179	L/2550/15	186	M/0800/53	189	M/3200/48	198	M/4058/15	309	M/4870/45	158	N/0050/PB17	210
J/8470/24	340	K/2110/17	172	L/0750/53	178	L/2550/17	186	M/0800/53	236	M/3200/50	198	M/4058/17	200	M/4870/45	205	N/0150/53	210
J/8470/PB15	325	K/2150/15	170	L/0800/53	178	L/2550/25	186	M/0850/53	189	M/3220/50	198	M/4058/17	309	M/4870/46	158	N/0211/45	210
J/8470/PB15	340	K/2150/15	172	L/0800/60	178	L/2700/48	186	M/0850/53	236	M/3300/46	159	M/4058/21	200	M/4870/46	205	N/0216/46	210
J/8490C/05	326	K/2150/17	172	L/0950/50	179	L/2700/60	186	M/0850/60	189	M/3300/46	198	M/4058/21	309	M/4870L/08	205	N/0300/50	211
J/8490C/05	377	K/2200/15	170	L/0950/53	179	LCMSKIT	18	M/0850/60	236	M/3450/48	196	M/4058/25RSS	200	M/4930/PB05	205	N/0350/48	210
J/8490C/90	326	K/2200/15	172	L/0950/60	179	LCMSKIT	308	M/0950/50	189	M/3500/53	197	M/4058/25RSS	309	M/4930/PB17	205	N/0400/53	211
J/8490C/90	377	K/2200/17	172	L/1000/50	179	LS118-1	308	M/0950/53	189	M/3550/48	197	M/4058/30RSS	200	M/4950/08	205	N/0500/46	158
J/8500/17	340	K/2250R/08	170	L/1000/53	179	LS118-1	369	M/0950/70	189	M/3650/50	195	M/4058/30RSS	309	M/4950/17	205	N/0500/46	211
J/8750/08	325	K/2250R/08	171	L/1000/60	179	LS118-212	308	M/0955/53	189	M/3650/60	195	M/4058/PB17	200	M/5000/46	158	N/0910/48	211
J/8750/08	344	K/2250R/15	170	L/1050/53	180	LS118-212	369	M/0955/53	236	M/3750/48	195	M/4058/PB17	309	M/5000/46	206	N/1250/45	211
J/8750/15	325	K/2250R/15	171	L/1050/70	180	LS118-4	308	M/1000/60	190	M/3750/50	195	M/4059/PB17	201	M/5000/48	158	N/1250/48	211
J/8750/15	344	K/2300R/08	170	L/1100/48	179	LS118-4	369	M/1000/60	236	M/3750/53	195	M/4059/PB17	312	M/5000/48	206	N/1252/44	211
J/8750/17	325	K/2300R/08	171	L/1100/50	179	LS118-500	308	M/1000/62	190	M/3750/60	195	M/4060/PB15	200	M/5000L/08	206	N/1252/45	211
J/8750/17	344	K/2300R/15	170	L/1150/50	179	LS118-500	369	M/1000/62	236	M/3751/48	195	M/4060/PB15	309	M/5020L/08	206	N/1270/44	158
J/8750/24	325	K/2300R/15	171	L/1150/53	179	LS119-1	308	M/1000/65	190	M/3754/90	195	M/4060/PB17	200	M/5050/08	206	N/1270/44	212
J/8750/24	344	K/2350/08	170	L/1200/50	180	LS119-1	370	M/1000/65	236	M/3755/48	195	M/4060/PB17	309	M/5078/08	206	N/1300/50	212
J/8780/17	325	K/2350/08	174	L/1350/50	181	LS119-212	308	M/1050/53	189	M/3851/53	198	M/4062/15	199	M/5078/15	206	N/1370/48	213
J/8780/17	353	K/2500/08	170	L/1350/53	181	LS119-212	370	M/1050/60	189	M/3874/48	198	M/4062/15	308	M/5078/17	206	N/1370/50	213
J/8800/PB08	312	K/2500/08	172	L/1400/48	180	LS119-4	308	M/1100/53	189	M/3876/46	198	M/4062/17	199	M/5078/21	206	N/1460/50	213
J/8800/PB08	359	K/2510/08	170	L/1445/53	181	LS119-4	370	M/1100/60	189	M/3900/17	202	M/4062/17	30				

Mfr. No	Page	Mfr. No	Page	Mfr. No	Page	Mfr. No	Page	Mfr. No	Page	Mfr. No	Page	Mfr. No	Page	Mfr. No	Page	Mfr. No	Page
O/0450/PB17	222	P/1038/48	226	P/1840/25	232	P/3200/PB05	243	P/4240/65	236	P/5080/62	254	P/6120/50	236	P/7480/27	265	P/7940/PB08	267
Q/0450/PB21	222	P/1039/46	226	P/1840/27	232	P/3298/44	243	P/4240/65	249	P/5120/53	254	P/6120/50	259	P/7481/15	265	P/7940/PB17	267
O/0500/25	222	P/1040/PB08	227	P/1880/17	233	P/3360/53	244	P/4240/70	249	P/5120/60	254	P/6120/53	236	P/7481/17	265	P/7950/PB15	267
O/0500/27	222	P/1040/PB17	227	P/1880/25	233	P/3360/60	244	P/4280/53	249	P/5240/53	254	P/6120/53	259	P/7481/21	265	P/7960/08	266
O/0500/PB08	222	P/1041/PB08	226	P/1880/27	233	P/3500/50	244	P/4280/60	249	P/5240/60	254	P/6120/60	236	P/7481/25	265	P/7960/15	266
O/0500/PB15	222	P/1041/PB17	226	P/2000/17	231	P/3520/PB05	244	P/4280/65	249	P/5240/65	254	P/6120/60	259	P/7481/27	265	P/7960/25	266
O/0500/PB17	222	P/1080/PB08	227	P/2040/17	232	P/3520/PB07	244	P/4280/68	249	P/5240/68	254	P/6120/65	236	P/7481/PB15	265	P/7960/PB08	266
O/0500/PC17	222	P/1100/07	227	P/2082/17	230	P/3600/87	244	P/4280/70	249	P/5245/53	254	P/6120/65	259	P/7486/15	264	P/7960/PB15	266
Q/0505/27	222	P/1120/46	227	P/2084/17	232	P/3600/88	244	P/4320/43	253	P/5245/60	254	P/6160/53	260	P/7486/15	311	P/7960/PB17	266
O/0505/PB17	222	P/1120/48	227	P/2084/PB08	232	P/3600/89	244	P/4360/50	250	P/5245/70	254	P/6162/50	259	P/7486/17	264	P/7965/PB05	266
O/0510/53	222	P/1120/53	227	P/2092/15	230	P/3610/88	244	P/4360/53	250	P/5248/53	255	P/6280/53	260	P/7486/17	311	P/7965/PB08	266
Q/0515/PB08	221	P/1160/48	227	P/2092/17	230	P/3640/87	244	P/4360/60	250	P/5248/60	255	P/6320/53	260	P/7490/08	169	P/7980/48	267
O/0515/PB08	312	P/1160/50	227	P/2094/17	231	P/3640/88	244	P/4400/50	249	P/5250/53	255	P/6320/60	260	P/7490/15	169	P/8000/43	267
Q/0580/53	222	P/1200/PB08	229	P/2100/PB17	233	P/3640/89	244	P/4400/53	249	P/5250/60	255	P/6360/50	260	P/7490/17	169	P/8000/44	267
O/0600/53	223	P/1200/PB17	229	P/2190/44	239	P/3652/43	244	P/4400/60	249	P/5250/70	255	P/6360/53	260	P/7490/21	169	P/8040/43	267
Q/0600/60	223	P/1240/25	229	P/1240/43	158	P/3675/08	245	P/4400/70	249	P/5320/48	255	P/6360/70	260	P/7490/25	169	P/8080/48	268
O/0600/62	223	P/1240/PB08	229	P/2200/43	239	P/3675/25	245	P/4440/53	250	P/5320/50	255	P/6400/53	260	P/7490/27	169	P/8120/48	267
O/0600/65	223	P/1240/PB17	229	P/2200/44	158	P/3676/08	245	P/4442/53	236	P/5320/53	255	P/6440/53	260	P/7490/FP21	169	P/8122/48	267
Q/0650/53	222	P/1280/PB08	228	P/2200/44	239	P/3676/25	245	P/4442/53	250	P/5320/60	255	P/6480/48	260	P/7490/PB08	169	P/8140/PB05	268
O/0650/60	222	P/1280/PB15	228	P/2200/48	239	P/3677/08	245	P/4442/63	236	P/5320/70	255	P/6520/53	261	P/7490/PB17	169	P/P560/45	261
Q/0650/65	222	P/1280/PB17	228	P/2220/PB05	239	P/3678/53	245	P/4442/63	250	P/5420/53	255	P/6520/65	261	P/7490/PC25	169	P/170-500	24
O/0650/70	222	P/1290/PB05	228	P/2315/PB17	239	P/3678/70	245	P/4520/53	250	P/5440/50	255	P/6560/50	236	P/7500/15	168	P/170-500	154
Q/0850/53	223	P/1290/PB05	311	P/2316/PB17	239	P/3679/53	245	P/4560/50	250	P/5440/53	255	P/6560/50	261	P/7500/17	168	PLBA12X	53
P/0120/43	223	P/1292/PB08	24	P/2320/53	240	P/3679/70	245	P/4600/50	250	P/5440/60	255	P/6560/53	236	P/7500/21	168	PLBA12X	316
P/0280/53	223	P/1292/PB08	228	P/2320/60	240	P/3680/53	245	P/4600/70	250	P/5480/60	255	P/6560/53	261	P/7500/24	168	PLCD12X	72
P/0295/48	223	P/1292/PB15	24	P/2360/53	236	P/2360/53	244	P/4640/60	251	P/5560/60	256	P/6560/65	236	P/7500/25	168	PLCD12X	317
P/0320/17	224	P/1292/PB15	228	P/2360/53	240	P/3681/70	244	P/4640/63	251	P/5600/53	256	P/6560/65	261	P/7500/27	168	PLCR13X	88
P/0320/25	224	P/1322/48	229	P/2360/60	236	P/3694/15	245	P/4680/50	251	P/5600/60	256	P/6560/70	236	P/7500/PB08	168	PLCR13X	317
P/0360/08	224	P/1360/17	233	P/2360/60	240	P/3700/48	245	P/4680/53	251	P/5600/61	256	P/6560/70	261	P/7500/PB17	168	PLRB23X	270
P/0360/17	224	P/1360/25	233	P/2360/61	236	P/3720/45	246	P/4680/60	251	P/5600/65	256	P/6600/50	262	P/7500/PC17	168	PLRB23X	320
P/0360/25	224	P/1400/08	230	P/2360/61	240	P/3720/48	246	P/4680/65	251	P/5600/68	256	P/6600/53	262	P/7500/PC25	168	PLRB23Y	270
P/0360/27	224	P/1400/17	230	P/2360/70	236	P/3760/53	247	P/4690/53	251	P/5640/53	236	P/6640/50	261	P/7501/15	168	PLRB23Y	320
P/0400/17	224	P/1400/25	230	P/2360/70	240	P/3760/60	247	P/4720/50	251	P/5640/53	256	P/6640/53	261	P/7501/17	168	PLSB53X	50
P/0400/25	224	P/1400/MS21	230	P/2390/43	158	P/3760/62	247	P/4720/53	251	P/5640/60	236	P/6720/53	262	P/7503/15	168	PLSB53X	316
P/0405/08	223	P/1440/08	231	P/2390/43	240	P/3760/65	247	P/4720/60	251	P/5640/60	256	P/6800/53	262	P/7503/17	168	PLSR13X	320
P/0440/53	224	P/1440/15	231	P/2390/44	158	P/3761/53	246	P/4720/70	251	P/5640/65	236	P/6800/60	262	P/7507/10RSS	167	PLSR13X	332
P/0480/53	224	P/1440/17	231	P/2390/44	240	P/3761/60	246	P/4760/53	252	P/5640/65	256	P/6800/65	262	P/7507/10RSS	311	PLTB23X	320
P/0600/90	224	P/1440/21	231	P/2392/82	240	P/3761/62	246	P/4760/60	252	P/5640/68	236	P/6800/70	262	P/7507/15	167	PLTB23X	343
P/0650/90	224	P/1440/25	231	P/2400/48	240	P/3880/50	247	P/4760/63	252	P/5640/68	256	P/6840/53	262	P/7507/15	311	PLTB23Y	320
P/0680/90	224	P/1440/27	231	P/2400/53	240	P/3880/53	247	P/4760/65	252	P/5640/70	236	P/6880/53	262	P/7507/17	167	PLTB23Y	343
P/0720/90	224	P/1440/PB17	231	P/2400/68	240	P/3920/50	247	P/4760/70	252	P/5640/70	256	P/6880/60	262	P/7507/17	311	PLTE32X	320
P/0840/53	224	P/1480/08	232	P/2400L/08	241	P/3920/53	247	P/4800/50	252	P/5645/17	256	P/6880/62	262	P/7507/25RSS	167	PLTE32X	342
P/0920/48	224	P/1480/15	232	P/2400L/17	241	P/3920/60	247	P/4800/53	252	P/5760/48	258	P/6960/53	263	P/7507/25RSS	311	PLTE53X	320
P/1006/10RSS	225	P/1480/17	232	P/2401L/08	241	P/3960/50	248	P/4800/60	252	P/5760/53	258	P/6960/60	263	P/7507/30RSS	167	PLTE53X	342
P/1006/10RSS	311	P/1480/21	232	P/2420/44	240	P/3960/53	248	P/4800/63	252	P/5760/70	258	P/6960/65	263	P/7507/30RSS	311	PLTL23X	321
P/1006/15	225	P/1480/25	232	P/2420/46	240	P/3960/60	248	P/4800/70	252	P/5800/48	257	P/7000/53	262	P/7507/PB17	167	PLTL23X	347
P/1006/15	311	P/1480/27	232	P/2420/53	240	P/3960/70	248	P/4806/50	252	P/5800/50	257	P/7000/60	262	P/7507/PB17	311	PLTL23Y	321
P/1006/17	225	P/1520/08	232	P/2420L/08	240	P/4000/50	248	P/4806/50	311	P/5800/53	257	P/7000/65	262	P/7508/17	167	PLTL23Y	347
P/1006/17	311	P/1520/17	232	P/2450/44	158	P/4000/53	248	P/4806/60	252	P/5800/60	257	P/7080/53	263	P/7508/17	309	Q/0120/48	268
P/1006/25RSS	225	P/1520/21	232	P/2450/44	159	P/4000/60	248	P/4806/60	311	P/5840/48	258	P/7140/53	263	P/7509/17	168	Q/0136/46	268
P/1006/25RSS	311	P/1520/25	232	P/2450/44	241	P/4020/48	247	P/4806/70	252	P/5840/50	258	P/7170/43	263	P/7509/17	312	Q/0150/50	268
P/1006/30RSS	225	P/1560/08	233	P/2460/44	158	P/4080/53	248	P/4806/70	311	P/5840/53	258	P/7200/50	263	P/7520/PB07	265	Q/0150/53	268
P/1006/30RSS	311	P/1560/15	233	P/2460/44	241	P/4080/60	248	P/4810/50	252	P/5840/60	258	P/7240/50	264	P/7555/17	168	Q/0152/53	268
P/1020/08	226	P/1560/17	233	P/2461/48	241	P/4080/62	248	P/4810/50	312	P/5840/65	258	P/7240/53	264	P/7555/21	168	R/0020/46	268
P/1020/17	226	P/1560/25	233	P/2520/48	241	P/4080/65	248	P/4840/50	254	P/5840/68	258	P/7280/50	263	P/7555/25	168	R/0030/46	269
P/1020/21	226	P/1600/17	233	P/2526/08	239	P/4080/70	248	P/4840/53	254	P/5880/48	258	P/7280/53	263	P/7560/25	265	R/0040/79	269
P/1020/25	226	P/1720/15	230	P/2530/PB05	239	P/4120/50	236	P/4880/50	254	P/5880/50	258	P/7300/PB08	264	P/7560/PB08	265	R/0040/99	269
P/1020/PB17	226	P/1720/17	230	P/2560/PB07	241	P/4120/50	248	P/4880/53	254	P/5880/53	258	P/7425/48	264	P/7560/PB17	265	R/0150/48	269
P/1021/15	226	P/1720/MS21	230	P/2600/50	241	P/4120/53	248	P/4880/60	254	P/5880/60	258	P/7440/08	236	P/7564/PB08	265	R/0150/50	269
P/1021/17	226	P/1760/15	230	P/2805/44	242	P/4120/53	248	P/4880/70	254	P/5880/63	258	P/7440/08	264	P/7600/08	266	R/0200/48	269
P/1021/PB17	226	P/1760/17	230	P/2840/50	242	P/4120/60	236	P/4920/53	253	P/5880/70	258	P/7440/17	236	P/7600/17	266	R/0200/50	269
P/1024/15	225	P/1760/21	230	P/2920/45	242	P/4120/60	248	P/4920/70	253	P/5920/53	259	P/7440/17	264	P/7640/08	266	R/0250/48	269
P/1024/17	225	P/1760/25	230	P/2920/46	242	P/4120/62	236	P/4960/53	253	P/5920/60	259	P/7440/25	236	P/7640/17	266	R/0300/44	269



Mfr. No	Page	Mfr. No	Page	Mfr. No	Page	Mfr. No	Page	Mfr. No	Page	Mfr. No	Page	Mfr. No	Page	Mfr. No	Page	Mfr. No	Page
S/1240/46	277	S/2720/62	284	S/3481/50	287	S/4440/53	292	S/5041/21	297	S/6400/70	303	S/8000/50	331	S/9231/PB08	25	SPCALMIX4500	322
S/1240/48	277	S/2720/65	284	S/3520/60	288	S/4440/60	292	S/5042/15	297	S/6480/53	303	S/8000/53	331	S/9231/PB08	337	SPCALMIX4500	327
S/1240/50	277	S/2760/53	284	S/3520/63	288	S/4440/65	292	S/5081/53	297	S/6490/48	303	S/8000/60	331	S/9231/PB15	25	SPCALMIX7100	322
S/1240/53	277	S/2760/60	284	S/3560/53	288	S/4450/53	292	S/5100/48	297	S/6490/68	303	S/8040/60	331	S/9231/PB15	337	SPCALMIX7100	327
S/1280/46	277	S/2760/65	284	S/3560/60	288	S/4450/60	292	S/5120/48	298	S/6520/60	304	S/8130/48	331	S/9231/PB17	25	SPCALMIX7500	322
S/1280/48	277	S/2840/62	284	S/3560/63	288	S/4450/65	292	S/5120/50	298	S/6560/53	304	S/8130/53	331	S/9231/PB17	337	SPCALMIX7500	327
S/1280/50	277	S/2840/70	284	S/3561/53	288	S/4450/68	292	S/5120/53	298	S/6560/60	304	S/8144/48	331	S/9240/17	338	SPCALMIX8100	322
S/1280/53	277	S/2880/53	284	S/3640/46	288	S/4450/70	292	S/5130/48	298	S/6600/53	304	S/8160/48	332	S/9240/25	338	SPCALMIX8100	327
S/1280/68	277	S/2880/60	284	S/3640/48	288	S/4460/53	292	S/5130/53	298	S/6600/60	304	S/8160/53	332	S/9240/26	338	SPCALMIX8500	322
S/1320/45	278	S/2880/62	284	S/3720/53	289	S/4460/53	312	S/5130/68	298	S/6600/63	304	S/8280/50	332	S/9240/PB15	338	SPCALMIX8500	327
S/1320/46	278	S/2880/65	284	S/3720/60	289	S/4480/53	292	S/5161/08	298	S/6600/65	304	S/8280/53	332	S/9240/PB17	338	SPCLAG22Y	276
S/1360/46	278	S/2880/70	284	S/3720/65	289	S/4480/60	292	S/5200/53	290	S/6600/70	304	S/8320/50	332	S/9240/PC15	338	SPCLAG22Y	320
S/1360/48	278	S/2920/53	284	S/3720/70	289	S/4480/65	292	S/5240/53	237	S/6640/53	304	S/8320/53	332	S/9240/PC17	338	SPCLAL22Y	28
S/1400/46	278	S/2920/60	284	S/3760/53	237	S/4520/53	237	S/5240/53	298	S/6640/60	304	S/8350/PB17	338	S/9360/PB17	338	SPCLAL22Y	316
S/1400/48	278	S/2920/62	284	S/3760/53	289	S/4520/53	292	S/5240/60	237	S/6640/63	304	S/8400/53	333	S/9440/PB08	339	SPCLAS22Y	51
S/1400/53	278	S/2920/65	284	S/3760/60	237	S/4520/60	237	S/5240/60	298	S/6640/65	304	S/8402/53	333	S/9500/PB17	339	SPCLAS22Y	31
S/1450/53	279	S/3000/60	285	S/3760/60	289	S/4520/60	292	S/5240/62	237	S/6640/68	304	S/8480/07	333	S/9560/PB08	340	SPCLAU11Y	137
S/1560/50	279	S/3000/63	285	S/3760/65	237	S/4520/63	237	S/5240/62	298	S/6650/53	304	S/8480/15	333	S/P500/53	287	SPCLAU11Y	318
S/1560/53	279	S/3040/53	285	S/3760/65	289	S/4520/63	292	S/5250/53	298	S/6650/60	304	S/8480/17	333	S/P500/65	287	SPCLBA22Y	53
S/1560/60	279	S/3040/60	285	S/3760/70	237	S/4520/65	237	S/5250/60	298	S/6650/65	304	S/8480/25	333	S/P520/53	288	SPCLBA22Y	316
S/1560/65	279	S/3045/50	285	S/3760/70	289	S/4520/65	292	S/5280/53	299	S/6660/60	304	S/8520/53	333	S/P520/63	288	SPCLBE22Y	57
S/1600/53	279	S/3100/63	286	S/3770/53	289	S/4520/70	237	S/5360/60	299	S/6840/53	305	S/8522/50	333	S/P525/53	292	SPCLBE22Y	317
S/1700/53	279	S/3105/63	286	S/3770/53	312	S/4520/70	292	S/5360/63	299	S/6840/60	305	S/8522/53	333	S/P525/63	292	SPCLBI21AY	58
S/1920/48	280	S/3105/70	286	S/3780/44	158	S/4690/53	293	S/5440/50	299	S/6840/63	305	S/8522/68	333	S/P530/48	289	SPCLBI21AY	317
S/2000/53	281	S/3120/53	286	S/3780/44	289	S/4720/53	293	S/5440/53	299	S/6850/53	305	S/8560/53	334	S/P530/53	289	SPCLBLKH2O	321
S/2000/60	281	S/3120/60	286	S/3800/53	289	S/4760/53	293	S/5480/48	299	S/6850/60	305	S/8560/60	334	S/P530/70	289	SPCLBLKH2O	369
S/2000/62	281	S/3120/63	286	S/3800/63	289	S/4800/60	294	S/5480/50	299	S/6920/53	342	S/8560/63	334	S/P770/53	333	SPCLBLKHCL	321
S/2000/65	281	S/3120/65	286	S/3880/50	290	S/4800/63	294	S/5520/53	299	S/6922/53	341	S/8560/65	334	S/P770/63	333	SPCLBLKHCL	369
S/2040/53	281	S/3120/70	286	S/3880/53	290	S/4800/65	294	S/5520/60	299	S/6960/48	305	S/8600/53	237	SPASBR92X	62	SPCLBLKHNO3	321
S/2040/60	281	S/3160/53	237	S/3920/50	290	S/4800/68	294	S/5520/63	299	S/6980/60	305	S/8600/53	334	SPASBR92Y	62	SPCLBLKHNO3	369
S/2040/62	281	S/3160/53	286	S/3920/53	290	S/4800/70	294	S/5560/53	299	S/7000/60	305	S/8600/60	237	SPASBR0392Y	60	SPCLCA22Y	74
S/2040/65	281	S/3160/60	237	S/3920/60	290	S/4840/60	294	S/5560/63	299	S/7000/63	305	S/8600/60	334	SPASCL91X	83	SPCLCA22Y	317
S/2046/50	281	S/3160/60	286	S/3920/63	290	S/4840/63	294	S/5600/53	300	S/7010/60	237	S/8600/63	237	SPASCL91Y	83	SPCLCAL1	322
S/2046/50	311	S/3160/63	237	S/3920/70	290	S/4845/60	294	S/5600/60	300	S/7010/60	305	S/8600/63	334	SPASCL92X	83	SPCLCAL1	327
S/2046/60	281	S/3160/63	286	S/4040/53	290	S/4845/70	294	S/5600/62	300	S/7010/63	237	S/8600/65	238	SPASCL92Y	83	SPCLCAL2	322
S/2046/60	311	S/3160/65	237	S/4040/70	290	S/4880/50	294	S/5640/53	300	S/7010/63	305	S/8600/65	334	SPASCL0292X	83	SPCLCAL2	327
S/2052/50	281	S/3160/65	286	S/4080/53	290	S/4880/53	294	S/5640/62	300	S/7010/70	237	S/8600/68	238	SPASCL0292Y	83	SPCLCAL3	322
S/2052/50	312	S/3160/68	237	S/4082/53	290	S/4880/60	294	S/5680/48	300	S/7010/70	305	S/8600/68	334	SPASCL0392X	83	SPCLCAL3	327
S/2080/53	281	S/3160/68	286	S/4086/50	290	S/4880/63	294	S/5720/48	300	S/7040/50	305	S/8600/70	238	SPASCL0392Y	83	SPCLCD22Y	72
S/2080/60	281	S/3160/70	237	S/4086/50	312	S/4880/65	294	S/5760/50	301	S/7040/53	305	S/8600/70	334	SPASF91AX	130	SPCLCD22Y	317
S/2080/65	281	S/3160/70	286	S/4100/53	291	S/4880/68	294	S/5760/53	301	S/7040/60	305	S/8620/46	334	SPASF91AY	130	SPCLCO22Y	92
S/2120/50	280	S/3161/53	237	S/4110/48	291	S/4920/53	237	S/5760/60	301	S/7040/65	305	S/8625/46	334	SPASF91X	130	SPCLCO22Y	317
S/2120/53	280	S/3161/53	285	S/4120/53	291	S/4920/53	294	S/5760/65	301	S/7120/53	306	S/8630/46	334	SPASF91Y	130	SPCLCR22Y	88
S/2120/60	280	S/3161/60	237	S/4120/60	291	S/4920/60	237	S/5760/70	301	S/7160/63	306	S/8640/53	335	SPASF92X	130	SPCLCR22Y	317
S/2120/70	280	S/3161/60	285	S/4160/53	291	S/4920/60	294	S/5761/53	300	S/7160/70	306	S/8640/70	335	SPASF92Y	130	SPCLCU22Y	95
S/2160/53	281	S/3161/63	237	S/4160/60	291	S/4920/63	237	S/5800/48	301	S/7200/53	306	S/8641/53	334	SPASNH3N92Y	36	SPCLCU22Y	317
S/2200/60	282	S/3161/63	285	S/4200/60	291	S/4920/63	294	S/5800/53	301	S/7200/60	306	S/8641/68	334	SPASN0292X	215	SPCLFE22Y	163
S/2272/50	282	S/3161/65	237	S/4200/63	291	S/4920/65	237	S/5840/50	301	S/7200/61	306	S/8660/50	335	SPASN0292Y	215	SPCLFE22Y	318
S/2330/48	282	S/3161/65	285	S/4200/65	291	S/4920/65	294	S/5840/53	301	S/7200/63	306	S/8680/53	335	SPASN02N92X	215	SPCLGE91AY	134
S/2360/48	282	S/3161/68	237	S/4200/70	291	S/4920/68	237	S/5880/53	301	S/7240/53	306	S/8720/48	335	SPASN02N92Y	215	SPCLGE91AY	318
S/2360/53	282	S/3161/68	285	S/4240/53	237	S/4920/68	294	S/5960/50	301	S/7240/60	306	S/8850/48	335	SPASN0392X	215	SPCLHG21AY	195
S/2360/68	282	S/3161/70	237	S/4240/53	291	S/4920/70	237	S/5960/53	301	S/7240/63	306	S/8850/53	335	SPASN0392Y	215	SPCLHG21AY	319
S/2380/48	282	S/3161/70	285	S/4240/60	237	S/4920/70	294	S/5965/53	301	S/7240/65	306	S/8920/60	336	SPASN03N92Y	215	SPCLIN21AY	160
S/2380/53	282	S/3165/53	285	S/4240/60	291	S/4930/05	295	S/5966/50	301	S/7240/70	306	S/8920/63	336	SPASPO492X	242	SPCLIN21AY	318
S/2440/53	236	S/3165/53	311	S/4240/63	237	S/4930/15	295	S/5966/50	312	S/7250/53	306	S/8970/60	336	SPASPO492Y	242	SPCLK22Y	246
S/2440/53	283	S/3200/50	286	S/4240/63	291	S/4930/17	295	S/5980/48	302	S/7250/60	306	S/8990/53	336	SPASPO4P92X	242	SPCLK22Y	319
S/2440/60	236	S/3240/53	286	S/4240/65	237	S/4930/21	295	S/5980/50	302	S/7250/63	306	S/9120/PB17	338	SPASPO4P92Y	242	SPCLMG22Y	187
S/2440/60	283	S/3280/53	287	S/4240/65	291	S/4930/25	295	S/5980/65	302	S/7320/48	307	S/9160/25	338	SPASSO492X	334	SPCLMG22Y	318
S/2440/70	236	S/3280/60	287	S/4240/70	291	S/4930/27	295	S/6040/53	302	S/7390/48	313	S/9160/26	338	SPASSO492Y	334	SPCLMN22Y	191
S/2440/70	283	S/3280/65	287	S/4246/50	291	S/4940/17	295	S/6041/53	302	S/7400/53	237	S/9160/PB15	338	SPASSO4S92X	334	SPCLMN22Y	319
S/2441/53	282	S/3320/53	237	S/4246/50	312	S/4940/17	312	S/6080/53	302	S/7400/53	313	S/9160/PB17	338	SPASSO4S92Y	334	SPCLMO92Y	209
S/2480/48	283	S/3320/53	287	S/4280/53	291	S/4945/31	295	S/6100/53	302	S/7400/63	237	S/9160/PB21	338	SPCALMIX10100	322	SPCLMO92Y	319
S/2560/45	283	S/3320/60	237	S/4360/53	293	S/4945/PB15	295	S/6120/60	303	S/7400/63	313	S/9200/PB08	338	SPCALMIX10100	328	SPCLMS1	322
S/2560/46	283	S/3320/60	287	S/4360/60	293	S/4945/PB21	295	S/6140/5									

Mfr. No	Page	Mfr. No	Page	Mfr. No	Page	Mfr. No	Page	Mfr. No	Page	Mfr. No	Page	Mfr. No	Page	Mfr. No	Page	Mfr. No	Page
SPMIXSTD3500	322	SPPLB93Y	317	SPPLC022Y	91	SPPLFE23X	318	SPPLLA23X	318	SPPLNB93Y	215	SPPLS93Y	336	SPPLSR23Y	320	SPPLW92Y	363
SPMIXSTD3500	329	SPPLBA22T	52	SPPLC022Y	317	SPPLFE23Y	163	SPPLNB93Y	319	SPPLSB72T	49	SPPLSR23Y	332	SPPLW93X	321	T/0481/07	312
SPMIXSTD3A100	322	SPPLBA22T	316	SPPLC023X	92	SPPLFE23Y	318	SPPLND22X	211	SPPLSB72T	316	SPPLTA92X	320	SPPLW93X	363	T/0481/07	343
SPMIXSTD3A100	329	SPPLBA22X	52	SPPLC023X	317	SPPLGA22X	134	SPPLLI12X	183	SPPLND22X	319	SPPLTA92X	341	SPPLW93Y	321	T/0486/07	312
SPMIXSTD3A500	322	SPPLBA22X	316	SPPLC023Y	92	SPPLGA22X	318	SPPLLI12X	318	SPPLND22Y	211	SPPLTB72X	316	SPPLW93Y	363	T/0486/07	343
SPMIXSTD3A500	329	SPPLBA22Y	52	SPPLC023Y	317	SPPLGA22Y	134	SPPLLI13X	183	SPPLND22Y	319	SPPLTB72Y	49	SPPLY22T	321	T/0490/48	344
SPMIXSTD4100	322	SPPLBA22Y	316	SPPLCR12X	88	SPPLGA22Y	318	SPPLLI13X	318	SPPLNI22T	212	SPPLTB72Y	316	SPPLY22T	372	T/0505/44	344
SPMIXSTD4100	329	SPPLBA23X	53	SPPLCR12X	317	SPPLGD22X	134	SPPLLI22X	183	SPPLNI22T	319	SPPLSB73X	49	SPPLY22X	321	T/0550/17	345
SPMIXSTD4500	322	SPPLBA23X	316	SPPLCR22T	88	SPPLGD22X	318	SPPLLI22X	318	SPPLNI22X	212	SPPLSB73X	316	SPPLY22X	372	T/0550/25	345
SPMIXSTD4500	329	SPPLBA23Y	53	SPPLCR22T	317	SPPLGD22Y	134	SPPLLI22Y	183	SPPLNI22X	319	SPPLSB73Y	49	SPPLY22Y	321	T/0550/PB17	345
SPMIXSTD4A100	322	SPPLBA23Y	316	SPPLCR22X	88	SPPLGD22Y	318	SPPLLI22Y	318	SPPLNI22Y	212	SPPLSB73Y	316	SPPLY22Y	372	T/0600/17	344
SPMIXSTD4A100	329	SPPLBE22T	57	SPPLCR22X	317	SPPLGD23X	134	SPPLLI23X	183	SPPLNI22Y	319	SPPLSC22T	272	SPPLY23X	321	T/0600/21	344
SPMIXSTD4A500	322	SPPLBE22T	316	SPPLCR22Y	88	SPPLGD23X	318	SPPLLI23X	318	SPPLNI23X	212	SPPLSC22T	320	SPPLY23X	372	T/0600/25	344
SPMIXSTD4A500	329	SPPLBE22X	57	SPPLCR22Y	317	SPPLGD23Y	134	SPPLLI23Y	183	SPPLNI23X	319	SPPLSC22X	272	SPPLY23Y	321	T/0600/PB17	344
SPMIXSTD5100	322	SPPLBE22X	316	SPPLCR23X	88	SPPLGD23Y	318	SPPLLI23Y	318	SPPLNI23Y	212	SPPLSC22X	320	SPPLY23Y	372	T/0601/17	344
SPMIXSTD5100	329	SPPLBE22Y	57	SPPLCR23X	317	SPPLGE92X	134	SPPLLU22X	186	SPPLNI23Y	319	SPPLSC22Y	272	SPPLY23Y	321	T/0601/25	344
SPMIXSTD5500	322	SPPLBE22Y	316	SPPLCR23Y	88	SPPLGE92X	318	SPPLLU22X	318	SPPLP92T	243	SPPLSC22Y	320	SPPLYB22X	372	T/0601/PB17	344
SPMIXSTD5500	329	SPPLBE23X	57	SPPLCR23Y	317	SPPLGE92Y	134	SPPLLU22Y	186	SPPLP92T	319	SPPLSC23X	272	SPPLYB22Y	321	T/0605/PB17	344
SPMIXSTD5A100	322	SPPLBE23X	317	SPPLCR92X	88	SPPLGE92Y	318	SPPLLU22Y	318	SPPLP92Y	243	SPPLSC23X	320	SPPLYB22Y	372	T/0611/46	312
SPMIXSTD5A100	329	SPPLBE23Y	57	SPPLCR92X	317	SPPLHF12X	138	SPPLMG12X	186	SPPLP92Y	319	SPPLSC23Y	272	SPPLYB22Y	321	T/0611/46	345
SPMIXSTD5A500	322	SPPLBE23Y	317	SPPLCR93X	88	SPPLHF12X	318	SPPLMG12X	318	SPPLP93X	243	SPPLSC23Y	320	SPPLYB22Y	372	T/0620/48	345
SPMIXSTD5A500	329	SPPLBI42X	58	SPPLCR93X	317	SPPLHF12Y	138	SPPLMG13X	187	SPPLP93X	319	SPPLSE22T	273	SPPLYB22X	321	T/0621/46	312
SPPLAG22T	276	SPPLBI42X	317	SPPLCS22X	81	SPPLHF12Y	318	SPPLMG13X	318	SPPLP93Y	243	SPPLSE22T	320	SPPLYB22X	372	T/0621/46	345
SPPLAG22T	320	SPPLBI42Y	58	SPPLCS22X	317	SPPLHF13X	318	SPPLMG22T	186	SPPLP93Y	319	SPPLSE22X	273	SPPLYB22Y	321	T/0700/08	346
SPPLAG22X	276	SPPLBI42Y	317	SPPLCS22Y	81	SPPLHG21AX	194	SPPLMG22T	318	SPPLPB22T	178	SPPLSE22X	320	SPPLYB22Y	372	T/0700/17	346
SPPLAG22X	320	SPPLBLKHCL	321	SPPLCS22Y	317	SPPLHG21AX	319	SPPLMG22Y	186	SPPLPB22T	318	SPPLSE22Y	273	SPPLYB22Y	321	T/0700/21	346
SPPLAG22Y	276	SPPLBLKHCL	369	SPPLCS23X	81	SPPLHG21AY	194	SPPLMG22Y	318	SPPLPB22X	178	SPPLSE22Y	320	SPPLYB22Y	372	T/0700/25	346
SPPLAG22Y	320	SPPLBLKHNO3	321	SPPLCS23X	317	SPPLHG21AY	319	SPPLMG23X	186	SPPLPB22X	318	SPPLSE23X	273	SPPLYB22Y	321	T/0700/27	346
SPPLAG23X	276	SPPLBLKHNO3	369	SPPLCS23Y	81	SPPLHG21X	194	SPPLMG23X	318	SPPLPB22Y	178	SPPLSE23X	320	SPPLYB22Y	372	T/0700/PB17	346
SPPLAG23X	320	SPPLC92X	78	SPPLCS23Y	317	SPPLHG21X	319	SPPLMG23Y	186	SPPLPB22Y	318	SPPLSE23Y	273	SPPLYB22Y	321	T/0700/PC21	346
SPPLAG23Y	276	SPPLC92X	317	SPPLCU12X	95	SPPLHG42T	194	SPPLMG23Y	318	SPPLPB23X	178	SPPLSE23Y	320	SPPLYB22Y	372	T/0701/15	346
SPPLAG23Y	320	SPPLC92Y	78	SPPLCU12X	317	SPPLHG42T	319	SPPLMN22T	191	SPPLPB23X	318	SPPLS92T	275	SPPLYB22Y	321	T/0701/17	346
SPPLAL12X	28	SPPLC92Y	317	SPPLCU13X	95	SPPLHG42X	194	SPPLMN22T	319	SPPLPB23Y	178	SPPLS92T	320	SPPLYB22Y	372	T/0701/21	346
SPPLAL12X	316	SPPLCA12X	74	SPPLCU13X	317	SPPLHG42X	319	SPPLMN22X	191	SPPLPB23Y	318	SPPLS92X	275	SPPLYB22Y	321	T/0701/25	346
SPPLAL13X	28	SPPLCA12X	317	SPPLCU22T	95	SPPLHG42Y	194	SPPLMN22X	319	SPPLPD32X	223	SPPLS92X	320	SPPLYB22Y	372	T/0701/PB15	346
SPPLAL13X	316	SPPLCA13X	74	SPPLCU22T	317	SPPLHG42Y	319	SPPLMN22Y	191	SPPLPD32X	319	SPPLS92Y	275	SPPLYB22Y	321	T/0701/PB17	346
SPPLAL22T	28	SPPLCA13X	317	SPPLCU22X	95	SPPLHG43X	194	SPPLMN22Y	319	SPPLPD32Y	223	SPPLS92Y	320	SPPLYB22Y	372	T/0702/PB15	346
SPPLAL22T	316	SPPLCA22T	73	SPPLCU22X	317	SPPLHG43X	319	SPPLMN23X	191	SPPLPD32Y	319	SPPLS93X	275	SPPLYB22Y	321	T/0702/PB17	346
SPPLAL22X	28	SPPLCA22T	317	SPPLCU22Y	95	SPPLHG43Y	194	SPPLMN23X	319	SPPLPR22X	264	SPPLS93X	320	SPPLYB22Y	372	T/0706/10RSS	312
SPPLAL22X	316	SPPLCA22Y	73	SPPLCU22Y	317	SPPLHG43Y	319	SPPLMN23Y	191	SPPLPR22X	319	SPPLS93Y	275	SPPLYB22Y	321	T/0706/10RSS	345
SPPLAL22Y	28	SPPLCA22Y	317	SPPLCU23X	95	SPPLHO22X	144	SPPLMN23Y	319	SPPLPR22Y	264	SPPLS93Y	320	SPPLYB22Y	372	T/0706/25RSS	312
SPPLAL22Y	316	SPPLCA23T	74	SPPLCU23X	317	SPPLHO22X	318	SPPLMO92T	209	SPPLPR22Y	319	SPPLS9A3X	276	SPPLYB22Y	321	T/0706/25RSS	345
SPPLAL23X	28	SPPLCA23T	317	SPPLCU23Y	95	SPPLHO22Y	144	SPPLMO92T	319	SPPLPT32X	244	SPPLS9A3X	320	SPPLYB22Y	372	T/0706/30RSS	312
SPPLAL23X	316	SPPLCA23X	74	SPPLCU23Y	317	SPPLHO22Y	318	SPPLMO92X	209	SPPLPT32X	319	SPPLSM22X	271	SPPLYB22Y	321	T/0706/30RSS	345
SPPLAL23Y	28	SPPLCA23X	317	SPPLDY22X	119	SPPLIN22X	160	SPPLMO92X	319	SPPLPT32Y	244	SPPLSM22X	320	SPPLYB22Y	372	T/0706/PB15	312
SPPLAL23Y	316	SPPLCA23Y	74	SPPLDY22X	317	SPPLIN22X	318	SPPLMO92Y	209	SPPLPT32Y	319	SPPLSM22Y	271	SPPLYB22Y	321	T/0706/PB15	345
SPPLAS12X	51	SPPLCA23Y	317	SPPLDY22Y	119	SPPLIN22Y	160	SPPLMO92Y	319	SPPLRB22X	270	SPPLSM22Y	320	SPPLYB22Y	372	T/0050/60	340
SPPLAS12X	316	SPPLCD22T	72	SPPLDY22Y	317	SPPLIN22Y	318	SPPLMO93X	209	SPPLRB22X	319	SPPLSM22Y	321	SPPLYB22Y	321	T/0050/62	340
SPPLAS22T	50	SPPLCD22T	317	SPPLER22X	119	SPPLIR32X	162	SPPLMO93X	319	SPPLRB22Y	270	SPPLSN22X	321	SPPLYB22Y	372	T/0050/68	340
SPPLAS22T	316	SPPLCD22X	72	SPPLER22X	317	SPPLIR32X	318	SPPLMO93Y	209	SPPLRB22Y	319	SPPLSN22X	321	SPPLYB22Y	321	T/0100/53	340
SPPLAS22X	50	SPPLCD22X	317	SPPLER22Y	119	SPPLIR32Y	162	SPPLMO93Y	319	SPPLRE92X	269	SPPLSN23X	350	SPPLYB22Y	372	T/0100/53	340
SPPLAS22X	316	SPPLCD22Y	72	SPPLER22Y	317	SPPLIR32Y	318	SPPLNA12X	280	SPPLRE92X	319	SPPLSN52T	321	SPPLYB22Y	321	T/0200/50	341
SPPLAS22Y	50	SPPLCD22Y	317	SPPLER23X	119	SPPLK12X	246	SPPLNA12X	320	SPPLRE92Y	269	SPPLSN52T	350	SPPLYB22Y	372	T/0200/53	341
SPPLAS22Y	316	SPPLCD23X	72	SPPLER23X	318	SPPLK12X	319	SPPLNA13X	280	SPPLRE92Y	319	SPPLSN52X	321	SPPLYB22Y	321	T/0200/60	341
SPPLAS23X	51	SPPLCD23X	317	SPPLER23Y	119	SPPLK13X	246	SPPLNA13X	320	SPPLRH32X	269	SPPLSN52X	350	SPPLYB22Y	372	T/0210/53	341
SPPLAS23X	316	SPPLCD23Y	72	SPPLER23Y	318	SPPLK13X	319	SPPLNA22T	280	SPPLRH32X	319	SPPLSN52Y	321	SPPLYB22Y	321	T/0210/65	341
SPPLAS23Y	51	SPPLCD23Y	317	SPPLER23Y	119	SPPLK22T	246	SPPLNA22T	320	SPPLRH32Y	269	SPPLSN52Y	350	SPPLYB22Y	372	T/0250/53	238
SPPLAS23Y	316	SPPLCE22X	80	SPPLER22Y	119	SPPLK22T	319	SPPLNA22X	280	SPPLRH32Y	319	SPPLSN53X	321	SPPLYB22Y	321	T/0250/53	340
SPPLAU32X	137	SPPLCE22X	317	SPPLER22Y	317	SPPLK22X	246	SPPLNA22X	320	SPPLRU32X	270	SPPLSN53X	350	SPPLYB22Y	372	T/0250/60	238
SPPLAU32X	318	SPPLCE22Y	80	SPPLER23X	119	SPPLK22X	319	SPPLNA22Y	280	SPPLRU32X	320	SPPLSN53Y	321	SPPLYB22Y	321	T/0250/60	340
SPPLAU32Y	137	SPPLCE22Y	317	SPPLER23X	318	SPPLK22Y	246	SPPLNA22Y	320	SPPLRU32Y	270	SPPLSN53Y	350	SPPLYB22Y	372	T/0300/53	342
SPPLAU32Y	318	SPPLCE23X	80	SPPLFE12X	162	SPPLK22Y	319	SPPLNA23X	280	SPPLRU32Y	320	SPPLSR12X	320	SPPLYB22Y	321	T/0420/08	343
SPPLB92T	59	SPPLCE23X	317	SPPLFE12X	318	SPPLK23X	246	SPPLNA23X	320	SPPLS92T	320	SPPLSR12X	332	SPPLYB22Y	372	T/0470/45	159
SPPLB92T	317	SPPLCE23Y	80	SPPLFE13X	318	SPPLK23X	319	SPPLNA23Y	280	SPPLS92T	336	SPPLSR22T	320	SPPLYB22Y	321	T/0470/45	343
SPPLB92X	59	SPPLCE23Y	317	SPPLFE22T	162	SPPLK23Y	246	SPPLNA23Y	320	SPPLS92X	320	SPPLSR22T	332	SPPLYB22Y	372	T/0470/47	159
SPPLB92X	317	SPPLCO12X	91	SPPLFE22T	318	SPPLK23Y	319	SPPLNB92X	215	SPPLS92X	336	SPPLSR22X	320	SPPLYB22Y	321	T/0470/47	343
SPPLB92Y	59	SPPLCO12X	317	SPPLFE22X	162	SPPLLA22X	176	SPPLNB92X	319	SPPLS92Y	320	SPPLSR22X	332	SPPLYB22Y	372	T/0476/46	312
SPPLB92Y	317	SPPLCO22T	91	SPPLFE22X	318	SPPLLA22X	318	SPPLNB92Y	215	SPPLS92Y	336	SPPLSR22Y	320	SPPLYB22Y	321	T/0476/46	343
SPPLB93X	60	SPPLCO22T	317	SPPLFE22Y	162	SPPLLA22Y	176	SPPLNB92Y	319	SPPLS93X	320	SPPLSR22Y	332	SP			

Mfr. No	Page	Mfr. No	Page	Mfr. No	Page	Mfr. No	Page	Mfr. No	Page
T/1470L/08	349	T/2323/17	354	T/3606/25RSS	312	U/P610/50	365	X/0200/21	372
T/1480/44	158	T/2500/48	355	T/3606/25RSS	360	U/P610/60	365	X/0200/25	372
T/1480/44	349	T/2500/53	355	T/3606/30RSS	312	UHPLCKIT	19	X/0200/27	372
T/1480/46	158	T/2650/53	355	T/3606/30RSS	360	UHPLCKIT	309	X/0200/PB17	372
T/1480/46	349	T/2660/46	355	T/3608/15	360	V/0011/PB05	362	X/0250/15	371
T/1480L/08	349	T/2750/79	355	T/3608/17	360	V/0100/48	366	X/0250/17	371
T/1500/43	349	T/2950/50	356	T/3622/46	361	V/0100/53	366	X/0250/21	371
T/1560/48	350	T/2950/53	356	T/3650/53	361	V/0100/70	366	X/0250/25	371
T/1560/53	350	T/2950/60	356	T/3652/48	361	V/0105/48	367	X/0250/27	371
T/1600/48	350	T/3000/50	356	T/3700/45	361	V/0200/53	367	X/0250/PB17	371
T/1600/53	350	T/3000/53	356	T/3700/46	361	W/0100/21	369	X/0255/17	371
T/1640/50	350	T/3000/60	356	T/3700/48	361	W/0100/24	369	X/0255/21	371
T/1645/48	351	T/3050/08	357	T/3710/48	362	W/0100/25	369	X/0255/27	371
T/1650/48	350	T/3050/25	357	T/3710/50	362	W/0100/27	369	X/0300/08	371
T/1652/48	351	T/3050/PB17	357	T/3710/60	362	W/0100/FP21	369	X/0300/17	371
T/1652/50	351	T/3100/08	356	T/3720/48	362	W/0106/10RSS	309	X/0300/27	371
T/1652/53	351	T/3100/17	356	T/3740/48	362	W/0106/10RSS	368	X/0350/08	371
T/1652/60	351	T/3100/21	356	T/3740/53	362	W/0106/15	309	X/0400/08	371
T/1654/48	351	T/3100/25	356	T/3751/08	363	W/0106/15	368	Y/0200/48	372
T/1654/50	351	T/3100/27	356	T/3751/17	363	W/0106/17	309	Z/0160/53	374
T/1654/53	351	T/3100/PB17	356	T/3753/44	362	W/0106/17	368	Z/0250/53	374
T/1654/60	351	T/3101/08	356	T/3757/08	363	W/0106/25RSS	309	Z/0250/61	374
T/1661/48	351	T/3101/PB17	356	T/3760/48	363	W/0106/25RSS	368	Z/0300/53	373
T/1661/53	351	T/3150/08	357	T/3770/46	363	W/0106/30RSS	309	Z/0300/60	373
T/1665/PB05	351	T/3150/17	357	T/3902/48	364	W/0106/30RSS	368	Z/0300/70	373
T/1670/48	351	T/3150/25	357	T/3940/48	364	W/0106/PB17	309	Z/0400/53	374
T/1800/PB15	353	T/3160/08	357	T/4050/15	364	W/0106/PB17	368	Z/0450/53	374
T/1840/48	352	T/3160/17	357	T/4200/17	364	W/0107/17	312	Z/0450/60	374
T/1900/53	353	T/3200/08	358	T/4200/25	364	W/0107/17	368	Z/0450/61	374
T/1950/PB08	352	T/3200/17	358	T/4200/PB15	364	W/0108/17	312	Z/0500/53	374
T/1958/PB08	352	T/3200/25	358	T/4200/PB17	364	W/0108/17	368	Z/0500/60	374
T/2000/PB08	353	T/3200/PB08	358	T/4203/60	364	W/0110/PB15	309	Z/0550/53	374
T/2060/46	158	T/3200/PB17	358	T/4205/60	364	W/0110/PB15	368	Z/0650/53	375
T/2060/46	353	T/3201/15	358	T/4206/60	364	W/0110/PB17	309	Z/0700/50	374
T/2150/17	355	T/3201/25	358	T/4210/48	364	W/0110/PB17	368	Z/0700/53	374
T/2150/25	355	T/3201/PB17	358	T/P010/15	361	W/0112/15	308	Z/0700/60	374
T/2150/27	355	T/3203/07	312	T/P050/15	361	W/0112/15	367	Z/0750/48	375
T/2150/PB15	355	T/3203/07	357	T/P050/24	361	W/0112/17	308	Z/0800/53	375
T/2150/PB17	355	T/3255/PB05	359	T/P190/04	342	W/0112/17	367	Z/0850/53	375
T/2200/08	355	T/3255/PB08	359	T/P190/17	342	W/0120/PB15	309	Z/0850/60	375
T/2200/17	355	T/3256/PB05	359	T/P500/50	356	W/0120/PB15	367	Z/0850/68	375
T/2200/21	355	T/3257/PB05	359	T/P500/60	356	W/0120/PB17	309	Z/0870/53	375
T/2200/25	355	T/3258/04	312	T/P505/48	357	W/0120/PB17	367	Z/0870/60	375
T/2200/27	355	T/3258/04	358	T/P630/48	362	W/0150/53	370	Z/1100/48	375
T/2200/FP21	355	T/3258/PB05	312	T/P630/53	362	W/0200/50	370	Z/1145/50	376
T/2200/PB17	355	T/3258/PB05	358	T/P630/60	362	W/0255/48	370	Z/1145/53	376
T/2250/17	355	T/3270/08	358	T/P631/48	362	W/0300/50	370	Z/1150/50	375
T/2250/25	355	T/3270/17	358	T/P631/53	362	W/0320/53	370	Z/1150/53	375
T/2250/PB17	355	T/3296/PB05	312	T/P631/60	362	W/0350/48	370	Z/1300/53	238
T/2300/15	354	T/3296/PB05	359	TAP-105-K	341	W/0350/53	370	Z/1300/53	376
T/2300/17	354	T/3300/48	359	U/0005/05	364	W6-1	308	Z/1300/60	238
T/2300/21	354	T/3350/48	359	U/0020/PB05	67	W6-1	367	Z/1300/60	376
T/2300/25	354	T/3510/10	359	U/0020/PB05	158	W6-212	308	Z/1350/53	376
T/2300/27	354	T/3520/50	87	U/0020/PB08	67	W6-212	367	Z/1350/60	376
T/2300/PB15	354	T/3600/08	360	U/0020/PB08	158	W6-4	308	Z/1350/70	376
T/2300/PB17	354	T/3600/17	360	U/0020/PB17	67	W6-4	367	Z/1400/53	376
T/2301/17	354	T/3600/21	360	U/0020/PB17	158	W6-500	308	Z/1550/53	376
T/2302/PB15	354	T/3600/25	360	U/0025/PB05	67	W6-500	367	Z/1550/60	376
T/2302/PB17	354	T/3600/27	360	U/0025/PB05	158	W9-1	24	Z/1550/61	376
T/2306/10RSS	312	T/3600/PB17	360	U/0025/PB08	67	W9-1	368	Z/1550/65	376
T/2306/10RSS	353	T/3601/15	360	U/0025/PB08	158	W9-2	24	Z/1600/53	376
T/2306/15	312	T/3601/17	360	U/0030/46	365	W9-2	368	Z/1600/65	376
T/2306/15	353	T/3601/25	360	U/0450/53	365	W9-500	24	Z/1600/70	376
T/2306/17	312	T/3606/10RSS	312	U/0450/60	365	W9-500	368	Z/1700/53	377
T/2306/17	353	T/3606/10RSS	360	U/0450/63	365	X/0100/17	372	Z/1850/50	377
T/2306/25RSS	312	T/3606/15	312	U/0500/53	365	X/0100/25	372	Z/1900/53	377
T/2306/25RSS	353	T/3606/15	360	U/0500/63	365	X/0100/27	372		
T/2306/30RSS	312	T/3606/17	312	U/0500/65	365	X/0100/MC15	372		
T/2306/30RSS	353	T/3606/17	360	U/0500/70	365	X/0100/PB17	372		
T/2310/PB17	353	T/3606/21	312	U/0650/78	366	X/0200/08	372		
T/2323/15	354	T/3606/21	360	U/0655/46	366	X/0200/17	372		



Ask your local account manager or Chemicals Specialist for a copy of our latest literature and technical resources.... or download them from www.acros.com. Register with your local Fisher Chemical distributor to receive new resources as they become available.

Optima LC/MS product range

Fisher Chemical OPTIMA LC/MS solvents have set the standard of excellence for consistent, reproducible performance in the mobile phase of LC/MS. These same high purity solvents are pre-blended with OPTIMA LC/MS modifiers such as Formic acid (FA) or Trifluoroacetic acid (TFA) to provide ready-to-use aqueous and organic mobile phase blends for LC/MS and LC/UV applications. Alternatively, solvents and modifiers can be purchased separately to be blended to your own specification. This brochure gives an overview of the product range and product specifications.



White paper – Optimizing Mobile Phase Solvent Purity for LC/MS

This paper presents comparative chromatographic data from four different instruments for LC/MS solvents offered by seven vendors. Data collected examine mass baseline, signal intensity, contamination by PEG and phthalates (plasticizers), and metal ion content which impacts ease of MS interpretation.



UHPLC product range

Our UHPLC gradient grade solvents have an extraordinarily low solvent filtration of 0.1 µm, ensuring low baseline noise and reducing pump blockages. Perfect for those working at high pressure, where high sensitivity and a fast run rate are key to successful analysis. This brochure gives an overview of the product range and product specifications.



Elemental Analysis - If you are looking for acids, look no further.

Whether you are analyzing environmental samples, etching glass, conducting routine testing or performing the most demanding laboratory procedures, Fisher Chemical can supply virtually any type of acid you need, in the grades, sizes and packaging that meet your requirements. This brochure gives an overview of the complete line of acid grades, detailing product specifications to ensure you can select the correct grade for your most challenging applications.



Specialized Chemical Services

This brochure gives an overview of our capabilities to support your chemistry when there is a need for products that are not listed amongst our catalogue products

- Semi-bulk and bulk requirements
- Tailored solvents and solvent blends or customized solutions
- Testing services or customized packaging and labeling





Fisher Scientific

Part of Thermo Fisher Scientific

AUSTRALIA
Thermo Fisher Scientific Australia
Pty Ltd
5 Caribbean Drive
Scoresby
Victoria 3179, Australia
Toll free customer service number:
1800 735 292
auinfo@thermofisher.com
www.thermofisher.au

CHINA
Thermo Fisher Scientific China
Toll free number:
800-810-5118, 400-650-5118
cnchemical@thermofisher.com
www.thermofisher.cn

Beijing office
Thermo Fisher Scientific
Units 702-715, 7/F,
Tower West, Yonghe Plaza,
No. 28, Andingmen East Street,
Beijing 100007, China
Tel: +86-10-8419-3588
Fax: +86-10-8419-3580
Shanghai office
Thermo Fisher Scientific
Building 7, No. 27 Xin Jinqiao Road
Shanghai 201206, China
Tel: +86-21-68654588

Guangzhou office
Thermo Fisher Scientific
2405-2406, Times Property Center,

INDIA
Thermo Fisher Scientific
403-404, B-Wing, Delphi
Hiranandani Business Park
Powai, Mumbai 400076, India
Tel: 18002097001 (Toll Free)
Fax: 022 6680 3001
qfc.customer@thermofisher.com

JAPAN
Thermo Fisher Scientific K.K.
4-24-11 Higashi-Ueno, Taito-ku
110-0015 Tokyo
Tel: 81 3 5826 1616
Fax: 81 3 5826 1650
sales.bid.jp@thermofisher.com
www.thermoscientific.jp/bid/

KOREA
Fisher Scientific Korea, Ltd.
Kookmin 1st Bldg. 7F, 1009-5
Daechi-dong, Gangnam-gu
Seoul, 135-851, Korea
Tel. (02) 3420.8700
Fax. (02) 3420.8710
sales.fsk@thermofisher.com
www.fishersci.co.kr

MALAYSIA
Fisher Scientific Malaysia Sdn Bhd
No. 3 Jalan Sepadu 25/123
Taman Perindustrian Axis Seksyen 25
40400 Shah Alam
Selangor Darul Ehsan Malaysia
Technical Service Hotline:
1-300-88-7868
Tel: (603) 51228888
Fax: (603) 51218899

NEW ZEALAND
Thermo Fisher Scientific
New Zealand Ltd
244 Bush Road
Albany, North Shore
Auckland, New Zealand
Toll free customer service number:
0800 933 966
nzinfo@thermofisher.com
www.thermofisher.co.nz

SINGAPORE
Fisher Scientific Singapore
Fisher Scientific Pte Ltd.
8 Pandan Crescent
LL4, #05-05 UE Tech Park
Singapore 128464
Tel : (65) 68736006
Fax: (65) 68735005
www.fishersci.com.sg



APEX CHEMICALS CO., LTD. (Head Office)

2034/87, Italthai Tower 19th Floor, New Petchburi Road, Bangkapi, Huaykwang, Bangkok 10310, Thailand

Tel. 02-038-9999 Fax. 02-038-0099



info@apexchemicals.co.th



ApexChemicals



@ApexChemicals



www.apexchemicals.co.th

